

OMNIS Coulometer Module



2.1019.0xx0

Product manual

8.1019.8002EN / 2025-05-16



Metrohm AG
Ionenstrasse
CH-9100 Herisau
Switzerland
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

OMNIS Coulometer Module

2.1019.0xx0

Product manual

8.1019.8002EN /
2025-05-16

Technical Communication
Metrohm AG
CH-9100 Herisau

This documentation is protected by copyright. All rights reserved.

This documentation is an original document.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Disclaimer

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g., conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

Table of contents

1	Overview	1
1.1	OMNIS Coulometer Module – Product description	1
1.2	OMNIS Coulometer Module – Product versions	1
1.3	About the documentation	2
1.4	Further information	3
1.5	Displaying the accessories	3
2	Safety	5
2.1	OMNIS Coulometer – Intended use	5
2.2	Responsibility of the operator	6
2.3	Requirements for operating personnel	6
2.4	Safety instructions	7
2.4.1	Danger from electrical potential	7
2.4.2	Danger from biological and chemical hazardous substances	7
2.4.3	Danger from highly flammable substances	8
2.4.4	Danger from leaking liquids	8
2.4.5	Danger during transport of the product	8
2.5	Design of warning messages	9
2.6	Meaning of warning signs	9
3	Functional description	11
3.1	OMNIS Coulometer Module – Overview	11
3.1.1	Coulometric Karl Fischer titration cell – Versions	12
3.1.2	Coulometric Karl Fischer titration cell – Overview	13
3.2	OMNIS Coulometer Module – Function	14
3.2.1	Magnetic stirrer – Functional description	14
3.2.2	Coulometric Karl Fischer titration cell – Functional description	14
3.3	OMNIS Module – Indicators	15
3.4	System – Signals	16
3.5	OMNIS Coulometer Module – Interfaces	17
4	Delivery and transport	19
4.1	Delivery	19
4.2	Packaging	19

12 Technical specifications	52
12.1 Ambient conditions	52
12.2 OMNIS Coulometer Module – Energy supply	52
12.3 OMNIS magnetic stirrer – Dimensions	52
12.4 OMNIS Coulometer Module – Dimensions	53
12.5 Magnetic stirrer – Dimensions	53
12.6 OMNIS Coulometer Module – Housing	53
12.7 Magnetic stirrer – Housing	54
12.8 OMNIS Coulometer Module – Connectors specifications	54
12.9 Display specifications	55
12.10 Current generator specifications	55
12.11 Measurement specifications	56
12.12 Magnetic stirrer – Specifications	58

1 Overview

1.1 OMNIS Coulometer Module – Product description

The OMNIS Coulometer Module is a module that is controlled by an OMNIS Coulometer or an OMNIS Titrator. The module works only together with an instrument that has an electrical and a network connection. Coulometric water content determination titrations according to Karl Fischer can be carried out with the aid of corresponding accessories.

You can request information on special applications in the "Application Bulletins" and "Application Notes", available free of charge via the regional Metrohm representative responsible. Various monographs on the subjects of titration techniques and electrodes are also available.


1.2 OMNIS Coulometer Module – Product versions

The product is available in the following versions:

Table 1 Product versions

Article number	Designation	Version feature
2.1019.0010	OMNIS Coulometer Module	without magnetic stirrer
2.1019.0110	OMNIS Coulometer Module	with integrated magnetic stirrer

The article number and serial number for identification of the product can be found on the type label:



1	(01) = Article number in accordance with GS1 standard	2	(21) = Serial number
3	(240) = Metrohm article number	4	Certification
5	Technical specifications	6	Certification
7	QR code		

1.3 About the documentation

Possible depictions in the documentation:

Depiction	Meaning
(5-12)	Cross-reference to figure legend (Figure number - Element in the figure)
1	Instruction step
Method	Parameters, menu items, tabs, and dialogs
File ► New	Menu path
[Continue]	Button or key
	Supplementary information to the descriptive text
	Note In graphics, orange arrows or frames indicate the reference to the descriptive text. The relevant elements may also be colored orange.



Movement

In graphics, blue arrows indicate the movement direction. The elements to be moved may also be colored blue.

1.4 Further information


Additional information on the OMNIS Coulometer Module is available on the following pages:

- Metrohm website <https://www.metrohm.com> – Overview of product family, document as PDF, details of the accessories, and information on applications.
- Software help of the OMNIS Software <https://guide.metrohm.com> – Thematically filtered individual content, videos, information on control software.

1.5 Displaying the accessories

Up-to-date information on the scope of delivery and on optional accessories can be found on the Metrohm website.

1 Searching for a product on the website


- Go to <https://www.metrohm.com>.
- Click on .
- Enter the article number of the product into the search field and press **[Enter]**.
 - Article number: See [OMNIS Coulometer Module – Product versions, chapter 1.2, page 1](#)
- In the result list, click on the desired product.

Detailed information regarding the product is displayed.


2 Displaying the accessories

- Scroll down (accessories subject to availability):
 - Included parts
 - Optional parts

3 Downloading the accessories list (included and optional parts)

- Click on  to download the accessories list as a PDF.



 Metrohm recommends keeping the downloaded PDF for reference purposes.

2 Safety

2.1 OMNIS Coulometer – Intended use

Metrohm products are used for the analysis and handling of chemicals.

Usage therefore requires the user to have basic knowledge and experience in handling chemicals. Knowledge regarding the application of fire prevention measures prescribed for laboratories is also mandatory.

Adherence to this technical documentation and compliance with the maintenance specifications make up an important part of intended use.

Any utilization in excess of, or deviating from, the intended use is regarded as misuse.

Specifications regarding the operating values and limit values of individual products are contained in the "Technical specifications" section, if relevant.

Exceeding and/or not observing the mentioned limit values during operation puts people and components at risk. The manufacturer assumes no liability for damage due to non-observance of these limit values.

The EU declaration of conformity loses its validity as soon as modifications are carried out on the products and/or the components.

The OMNIS Coulometer is the central instrument of an OMNIS titration system for coulometric Karl Fischer titrations. Coulometric water content determination titrations according to Karl Fischer can be carried out with the aid of corresponding accessories.

The OMNIS Coulometer Module is a module that is controlled by an OMNIS Coulometer. The module works only together with an instrument that has an electrical and a network connection.

2.4 Safety instructions

2.4.1 Danger from electrical potential

Contact with electrical potential can cause serious injuries or death. To avoid danger from electrical potential, observe the following:

- Operate the product only if it is in perfect condition. The housing must also be intact.
- Only use the product with the covers fitted. If covers are damaged or missing, disconnect the product from the energy supply and contact the regional Metrohm service representative.
- Protect live components (e.g. power supply unit, power cord, connection sockets) against moisture.
- Always have maintenance work and repairs on electrical components carried out by a regional Metrohm service representative.
- Disconnect the product from the energy supply immediately if at least one of the following cases occurs:
 - The housing is damaged or open.
 - Live parts are damaged.
 - Moisture penetrates.

2.4.2 Danger from biological and chemical hazardous substances

Contact with biological hazardous substances may cause poisoning from toxins or infections from microorganisms. Contact with aggressive chemical substances may cause poisoning or chemical burns. To avoid danger from biological or chemical hazardous substances, observe the following:

- Label the product according to regulations if it is used for substances that have a potential for chemical hazards and are generally subject to the Hazardous Substances Ordinance.
- Wear personal protective equipment (e.g. protective glasses, gloves).
- Use exhaust equipment when working with vaporizing hazardous substances.
- Dispose of hazardous substances in accordance with regulations.
- Clean and disinfect contaminated surfaces.
- Only use detergents that do not cause any unwanted side reactions with the materials to be cleaned.
- Dispose of chemically contaminated materials (e.g. cleaning material) in accordance with regulations.
- Proceed as follows in case of a return shipment to Metrohm AG or a regional Metrohm representative:
 - Decontaminate the product or product component.
 - Remove the labeling for hazardous substances.
 - Create a declaration of decontamination and enclose it with the product.



2.4.3 Danger from highly flammable substances

Using highly flammable substances or gases may cause fires or explosions. To avoid danger from highly flammable substances, observe the following:

- Avoid ignition sources.
- Use protective grounding.
- Use exhaust equipment.

2.4.4 Danger from leaking liquids

Leaking liquids may cause injuries and may damage the product. To avoid danger from leaking liquids, observe the following:

- Check the product and its accessories for leakages and loose connections.
- Replace leaking components and connecting elements without delay.
- Tighten loose connecting elements.
- Do not loosen tubing connections under pressure.
- Do not remove aspiration tubing under pressure.
- Carefully pull the tubing ends out of the vessels.
- Carefully allow liquids from the tubing to drain into suitable vessels.
- Insert the tubing tips completely into the vessels.
- Remove and dispose of escaping liquids in accordance with regulations.
- If you suspect that liquid has penetrated the instrument, disconnect the instrument from the energy supply. Then have the instrument checked by a regional Metrohm service representative.

2.4.5 Danger during transport of the product

Chemical or biological substances may be spilled during the transport of the product. Parts of the product may fall down or may be damaged. There is a risk of injury from chemical or biological substances and pieces of broken glass. To ensure safe transport, observe the following:

- Remove loose parts (e.g. sample racks, sample vessels, bottles) before transport.
- Remove liquids.
- Lift and transport the product with both hands on the base plate.
- Lift and transport heavy products only according to instructions.

2.5 Design of warning messages

The present documentation uses warning messages as follows.

Structure

1. Severity of the danger (signal word)
2. Type and source of danger
3. Consequences of disregarding the danger
4. Measures for averting the danger

Hazard levels

Signal color and signal word designate the hazard level.

DANGER

Indicates an immediate danger. It will result in serious injuries or death if not avoided.

WARNING

Indicates a potential danger. Failure to avoid the danger may result in death or serious injury.

CAUTION

Indicates a potential danger. If not avoided, it may result in light or minor injuries.

NOTICE

Indicates a potentially damaging situation. If not avoided, the product or something in the surrounding area could be damaged.












2.6 Meaning of warning signs

Warning signs on the product or in the documentation indicate potential dangers or draw attention to certain behaviors in order to avoid accidents or damage.

Depending on the application purpose, the operating company attaches additional warning signs to the product. The corresponding instructions of the operator must be followed.

Table 2 Warning signs according to ISO 7010 (examples)



Warning signs / meaning		Warning signs / meaning	
	General warning sign		Warning of hot surface
	Warning of sharp object (cut/puncture)		Warning of hand injuries (crushing)
	Warning of electrical voltage		Warning of corrosive substances
	Warning of optical radiation		Warning of a laser beam
	Warning of flammable materials		Warning of biological hazard
	Warning of toxic materials		

3 Functional description

3.1 OMNIS Coulometer Module – Overview




Figure 1 OMNIS Coulometer Module (with magnetic stirrer) – Overview

1	Lid	2	Internal measuring interface
3	Support rod (6.2016.050)	4	Titration cell holder (6.02047.020)
5	Clamping ring (6.2013.010)	6	Magnetic stirrer
7	Control bar of the magnetic stirrer	8	Status display



3.1.1 Coulometric Karl Fischer titration cell – Versions




Figure 2 3 versions of the coulometric Karl Fischer titration cell

- | | |
|---|--|
| <p>1 Karl Fischer titration cell / 80–250 mL / coulometric (6.1464.320)</p> | <p>2 Karl Fischer titration cell made of amber glass / 80–250 mL / coulometric (6.1464.323)</p> |
| <p>3 Karl Fischer titration cell with 2 side openings / 80–250 mL / coulometric (6.1465.320)</p> | |



3.1.2 Coulometric Karl Fischer titration cell – Overview




Figure 3 Coulometric Karl Fischer titration cell (equipped) – Overview

1 Adsorber tube (6.1403.030)

With ground-joint sleeve (6.2713.020)

In combination with an OMNIS Sample Robot Oven, if necessary with tubing olive (6.1808.310) for a tubing for discharging gases.

2 Indicator electrode cable (6.02104.040)

With blue coding for polarizable metal electrodes

- Coulometric Karl Fischer titration cell (6.1464.323) made of amber glass,
- Coulometric Karl Fischer titration cell (6.1465.320) with 2 side openings.

The coulometric Karl Fischer titration cell is attached to a titration cell holder on the support rod. The amber glass version is recommended for light-sensitive materials.

3.3 OMNIS Module – Indicators




Figure 4 OMNIS Module – Indicators







1 Status display
Multi-colored

The status of the instrument is displayed with the status display (4-1) using different colors. (see "System – Signals", chapter 3.4, page 16)



3.4 System – Signals

System components with status indicators show their operating status with colors and/or flashing patterns. The meaning of the colors and flashing patterns is explained in the following table.

Visual signal		Meaning
	LED lights up yellow.	System start or initialization
	LED flashes yellow (slowly).	Ready for connection setup or locking
	LED flashes yellow (fast).	Connection setup started or locking underway
	LED lights up green.	Ready for operation
	LED flashes green (slowly).	In operation
	LED flashes red (fast).	Malfunction or error

Some system components only use part of the explained flashing patterns.



3.5 OMNIS Coulometer Module – Interfaces




Figure 5 OMNIS Coulometer Module – Connector

1 MDL connector

MDL = Metrohm Device Link

Connection socket for connecting cables
between OMNIS instruments.




Figure 6 OMNIS Coulometer or OMNIS Coulometer Module – Internal measuring interface

<p>1 Indicator Indicator for the GENERATOR outlet</p>	<p>2 GENERATOR Outlet for generator electrode (gray coding)</p>
<p>3 INPUT 1 Temperature sensor (red coding) or Measuring input for polarizable metal electrode (blue coding) or Measuring input for potentiometric sensor (green coding)</p>	<p>4 INPUT 2 Temperature sensor (red marking) or Measuring input for potentiometric sensor (green coding)</p>
<p>5 Indicator Indicator for the internal measuring interface</p>	

Measuring inputs INPUT 1 and INPUT 2

The measuring inputs **INPUT 1** and **INPUT 2** are marked with colored circle segments. The markings indicate that only certain types of electrode cables may be plugged into the connection socket:

Table 3 Meaning of the colors

Red	The connector supports temperature sensors.
Blue	The connector supports polarizable metal electrodes.
Green	The connector supports potentiometric sensors.
Gray	The connector supports a generator electrode.

4 Delivery and transport


4.1 Delivery

Inspect the delivery immediately upon receipt:


- Check the delivery against the delivery note to ensure completeness.
- Check the product for damage.
- If the delivery is incomplete or damaged, contact your regional Metrohm representative.

4.2 Packaging

The product and accessories are supplied in protective special packaging. Keep this packaging to ensure safe transportation of the product. If a transport locking device is present, keep this as well for future reuse.

1 Mounting the support rod


- Screw the support rod onto the stand attachment.

2 Mounting the clamping ring

- Push the clamping ring over the support rod with the groove facing upward.



3 Mounting the titration cell holder



- Press the green locking lever on the titration cell holder.
- Push the titration cell holder over the support rod.
- To fix in place, release the green locking lever at the desired height.

5.4 Replacing the adsorber material





Depending on the OMNIS product, different adsorber cartridges or adsorber tubes are available.

Table 4 Available adsorber cartridges or adsorber tubes

Adsorber cartridge / adsorber tube	Article number	Figure
------------------------------------	----------------	--------


Adsorber cartridge for OMNIS Solvent Module	6.01807.000	
---	-------------	--



Adsorber cartridge / adsorber tube	Article number	Figure
Adsorber tube for coulometric Karl Fischer titration cell	6.1403.030	
Adsorber tube for volumetric Karl Fischer titration cell	6.01406.010	
Adsorber tube for a cylinder unit OMNIS	6.1619.020	
Adsorber tube for waste bottle for OMNIS Dosing Module	6.1609.000	



1 Removing the lid from the housing



- Adsorber cartridge: Unlatch and remove the lid including the seal from the housing.
- Adsorber tube: Remove the lid by pulling it out of the housing.

2 Removing the molecular sieve (if present)


- Remove the molecular sieve and regenerate for at least 24 hours in the drying oven at 300 °C. Place in a desiccator to cool down and then seal airtight in a glass bottle, see also [FAQ on Karl Fischer titration](#).

3 Filling the molecular sieve

- Adsorber cartridge: Place a base-covering cotton plug loosely into the housing at the bottom. Do not pack the wad of cotton too tightly as sufficient gas flow must be possible. Use a molecular sieve to fill the housing to approx. 1 cm under the housing edge.
- Adsorber tube: Place a small cotton plug on the molecular sieve. Do not pack the wad of cotton too tightly as sufficient gas flow must be possible.

4 Sealing the housing with the lid

- i Make sure that the sealing surface between the housing and the lid is clean and dry, and that there no residual filling material whatsoever is present.



- Adsorber cartridge: Hook the lid including the seal into the housing side and close it by clicking it into place.
- Adsorber tube: Seal the housing with the lid.

i At moderate humidity, replace the molecular sieve approx. every 6 weeks.

An increase in drift is an indication that the molecular sieve is saturated and that the humidity is therefore entering the Karl Fischer titration cell.

Hint:

After replacing the molecular sieve, write the date on the adsorber housing.

5.5 Equipping the coulometric Karl Fischer titration cell

CAUTION

Risk of cuts from sharp edges

Cut injuries due to damaged glass parts and pieces of broken glass.

- Handle glass parts (e.g. electrodes, bottles) with care.
- Only use undamaged glass parts.
- Dispose of damaged glass parts immediately.

Preparing the coulometric Karl Fischer titration cell

Prerequisite:

- 7 Fill the Karl Fischer titration cell. (see "Filling the coulometric Karl Fischer titration cell", chapter 5.6, page 27).
- 8 Depending on the application, insert the required adapter into the opening for application versions.

See also

Coulometric Karl Fischer titration cell – Overview (chapter 3.1.2, page 13)

5.6 Filling the coulometric Karl Fischer titration cell



WARNING

Contact with chemicals

Chemicals may cause chemical burns.

- Wear personal protective equipment (e.g. protective glasses, gloves).
- Use exhaust equipment when working with vaporizing hazardous substances.

Use of a generator electrode with diaphragm

Prerequisite:

- The Karl Fischer titration cell is completely equipped with a generator electrode with diaphragm. Pre-dry the components in the oven at 50 °C.

- 1 Remove the adsorber tube from the generator electrode.
- 2 Fill approximately 5 mL of catholyte into the generator electrode.
- 3 Insert the adsorber tube into the generator electrode.
- 4 Remove the right ground-joint stopper from the Karl Fischer titration cell.
- 5 Using a funnel, fill as much anolyte into the Karl Fischer titration cell until the 100 mL marking on the titration cell is reached.
The level of the anolyte should be roughly 1 to 2 mm above the level of the catholyte.
- 6 Close the ground-joint opening on the right with the ground-joint stopper (with ground-joint sleeve attached).


Use of a generator electrode without diaphragm**Prerequisite:**

- The Karl Fischer titration cell is completely equipped with a generator electrode without diaphragm. Pre-dry the components in the oven at 50 °C.

- 1** Remove the right ground-joint stopper from the Karl Fischer titration cell.
- 2** Use a funnel to fill approximately 100 mL of reagent into the Karl Fischer titration cell.
- 3** Close the ground-joint opening on the right with the ground-joint stopper (with ground-joint sleeve attached).

5.7 Mounting the coulometric Karl Fischer titration cell**Prerequisite:**

- The support rod is mounted with clamping ring and titration cell holder (see *Mounting the magnetic stirrer accessories, page 20*).


1 Placing the titration cell

- Slide the metal clip of the coulometric Karl Fischer titration cell over the titration cell holder.

2 Checking the position of the titration cell

- The titration cell is positioned in the center of the magnetic stirrer.


3 Setting the clamping ring



- Push the clamping ring under the titration cell holder.
- Rotate the clamping ring in such a way that the wedge on the titration cell holder fits into the groove of the clamping ring.
- Fix the clamping ring with the knurled screw in position at the desired height.

The clamping ring is used as the lower stop for the titration cell holder. The stop facilitates the correct positioning of the titration cell on the magnetic stirrer.

4 Guiding out the cable



Guide the cables out going under the bar.

5 Close the lid.



6 Start-up

6.1 Initial start-up by Metrohm

As a basic rule, the initial start-up of the system is carried out by the regional Metrohm service representative.

7 Coulometric titration

The **coulometric Karl Fischer titration** is a variation of the classic water content determination method according to Karl Fischer.

7.1 OMNIS Coulometer Module – Principle of coulometry according to Karl Fischer

The **coulometric Karl Fischer titration** is a variation of the classic water content determination method according to Karl Fischer. The conventional method works with a methanolic solution of iodine, sulfur dioxide and a base as buffer substance. If an aqueous sample is titrated, then several reactions take place that can be summarized in the following sum equation:



According to the equation above, the I_2 reacts quantitatively with H_2O . This chemical equation serves as a basis for the water content determination.

With the **coulometric Karl Fischer titration**, the necessary iodine is directly and electrochemically generated in the electrolyte containing iodine. Between the amount of electric charge and the amount of generated iodine, there is a strictly quantitative relationship, which is used for high-precision dosing of the iodine. Because the coulometric Karl Fischer method is an **absolute determination**, no titer needs to be determined. It must only be ensured that the reaction generating the iodine runs with a 100% current efficiency. All of the reagents available today ensure this.

The endpoint indication is effected voltametrically by modulating an alternating current of constant strength to a double Pt electrode. This results in a voltage difference between the Pt wires. This is drastically reduced as soon as even the slightest amounts of free iodine are present. This circumstance is used for detecting the endpoint of the titration.

7.2 OMNIS Coulometer Module – Working with water standards

Certified water standards

Commercially available, certified water standards with water contents of 1.00 ± 0.003 mg/g and/or 0.10 ± 0.005 mg/g should be used for validating the instrument as a whole, integrated system.

i The 1.0 mg/g water standard is easier to handle and is therefore preferred.

Table 5 Recommended weighing ranges

1.0 mg/g water standard	0.2–2.0 g
0.1 mg/g water standard	0.5–5.0 g

7.3 OMNIS Coulometer Module – Sample addition

This chapter contains a few notes concerning sample addition. An exhaustive discussion of this topic is not possible here. Further notes can be found in the literature by the reagent manufacturers and in the following **Metrohm Application Bulletins**:

Bulletin no.	Title
No. 137	Coulometric water content determination according to Karl Fischer
No. 142	Karl Fischer water content determination in non-explosive gases
No. 145	Determination of low water contents in plastics using the KF oven method
No. 209	Coulometric water content determinations according to the Karl Fischer method in insulating oils, hydrocarbons, and their products

Values for the sample sizes

The sample weight should be small in order to be able to titrate as many samples as possible in the same electrolyte solution and to keep the titration time short. However, ensure that the sample contains at least 50 µg of H₂O. The following table helps you determine the appropriate sample size.

Table 6 Recommended sample sizes

Water content of the sample	Sample size	Resulting water content
10,000 ppm = 1%	10–100 mg	100–1,000 µg
1,000 ppm = 0.1%	100 mg–1 g	100–1,000 µg
100 ppm = 0.01%	1 g	100 µg
10 ppm = 0.001%	5 g	50 µg

Working with liquid samples

Liquid samples are added with a syringe. The samples can be injected two different ways:

- Use a syringe with a long needle that you immerse in the reagent during the injection.
- Use a syringe with a short needle and aspirate the last drop back into the needle.

The best way for you to determine the injected sample amount is to weigh the sample by difference.

Glass syringes should be used for the **determination of traces and validations**. Metrohm recommends obtaining these from a specialized syringe manufacturer.

Highly volatile samples and samples of low viscosity should be cooled before sampling. Doing so avoids losses while working. The syringe must, however, not be cooled directly, as condensation could be formed. For the same reason, no air may be aspirated into a syringe into which a cooled sample has been aspirated beforehand.

Samples of high viscosity can be thinned by heating. The syringe must be heated as well. The same target can be reached by diluting with suitable solvents. In this case, the water content of the solvent has to be determined and subtracted as a blank value.

If samples contain only **traces of water**, then the syringe has to be predried well. If possible, the syringe should be rinsed with the sample solution by filling in and discarding the solution several times.

Working with solid samples

Solid samples such as powders, pastes, fats, and oils are extracted or dissolved in a suitable solvent if possible. The resulting solution is injected, during which a blank value correction for the solvent must be carried out.

If no suitable solvent can be found for a solid sample, or if the sample reacts with the Karl Fischer reagent, then a Karl Fischer oven should be used.

8.2.1 Switching the magnetic stirrer on and off

1 Switching the magnetic stirrer on


Press the  key.

The magnetic stirrer stirs with the most recently used stirring rate.

2 Switching the magnetic stirrer off

Press the  key once again.

The magnetic stirrer stops.

 If the magnetic stirrer is running at a high stirring rate, reduce the stirring rate before switching it off.

As an alternative, switch the magnetic stirrer on and off in the OMNIS Software under *Manual control*.

8.2.2 Setting the magnetic stirrer

The stirring rate can be adjusted in 15 steps.

Prerequisite:

The magnetic stirrer is switched on.

1 Increasing the stirring rate in steps

Press the  key.


Each pressing of the key increases the stirring rate by 1 step. The current stirring rate appears in the OMNIS Software under **Manual control**.

2 Reducing the stirring rate

Press the  key.

Each pressing of the key reduces the stirring rate by 1 step. The current stirring rate appears in the OMNIS Software under **Manual control**.

Alternatively, the stirring rate can also be set in the OMNIS Software under *Manual control*.

 The stirring direction can be exclusively set in the OMNIS Software under **Manual control**.




Figure 8 Reagent replacement with OMNIS Dosing Module – Example

<p>1 Antidiffusion tip (6.1543.200) without antidiffusion valve</p>	<p>2 Stoppers for reagent replacement (6.1446.060) with ground-joint sleeve (6.2713.000)</p>
<p>3 Adapter for reagent replacement (6.2730.030) with nipple and O-ring</p>	<p>4 FEP tubings (6.1805.100)</p>
<p>5 OMNIS Liquid Adapter (6.01600.010) on a bottle with KF reagent</p>	<p>6 Adsorber tube for a cylinder unit OMNIS (6.1619.020)</p>
<p>7 Waste bottle (6.1608.030)</p>	<p>8 Cylinder unit OMNIS 50 mL (6.01503.250)</p>

For a reagent replacement with an OMNIS Dosing Module, a tubing setup such as that shown in the figure is required. Proceed as follows:

Preparing the reagent replacement

Prerequisite:

- The ground-joint stopper has been removed from the ground-joint opening on the right.

8.3.2 Reagent replacement with OMNIS Solvent Module

When replacing reagents with the OMNIS Solvent Module, the used reagent is aspirated through the PTFE tubing from the Karl Fischer titration cell and pumped into the waste bottle.

i If necessary, clean the titration cell with multiple rinsing cycles.

The fresh reagent is pumped into the titration cell via the OMNIS Solvent Module.




Figure 9 Reagent replacement with OMNIS Solvent Module

- | | |
|--|---|
| 1 Siphon Breaker (6.01600.200)
on bottle with KF reagent | 2 PTFE tubing (6.1805.200) |
| 3 Waste bottle. (6.1608.030) | 4 Adapter for reagent replacement
(6.1446.220)
with ground-joint sleeve (6.2713.000) |

For a reagent replacement with the OMNIS Solvent Module, a tubing setup such as that shown in the figure is required. Proceed as follows:



Preparing reagent replacement

Prerequisite:

- The ground-joint stopper has been removed from the ground-joint opening on the right.

Required accessories:

- *(see figure 9, page 43)*

1 Inserting the adapter

- Place the ground-joint sleeve over the adapter.
- Insert the adapter together with the ground-joint sleeve into the right-hand ground-joint opening of the titration cell with the aspiration tips facing downwards.

2 Connecting the titration cell with the waste bottle

- Screw the first PTFE tubing onto the closed aspiration tip.
- Screw the other end of the PTFE tubing onto the waste bottle to pump the used reagent into the waste bottle.

3 Connecting the titration cell with the Siphon Breaker

- Screw the second PTFE tubing onto the open aspiration tip.
- Screw the other end of the PTFE tubing onto the Siphon Breaker on the fresh reagent to pump the fresh reagent into the titration cell.

8.3.3 Manual reagent replacement




Figure 10 OMNIS Solvent Module – Indicators and controls

- | | |
|--|---|
| <p>1 Add key
Deliver liquid (Solvent) into the titration cell</p> | <p>2 Aspirate key
Aspirate waste from the titration cell</p> |
| <p>3 Status display
Multi-colored</p> | |

Prerequisite:

- OMNIS Solvent Module is connected.
- The solvent bottle, waste bottle and Karl Fischer titration cell are fully mounted and connected with the corresponding tubings.

1 Manually emptying the KF titration cell

Press the  key on the OMNIS Solvent Module:

The OMNIS Solvent Module starts aspirating waste out of the Karl Fischer titration cell into the waste bottle.

Different versions are possible:

- Long pressing (> 1 s): Aspirating takes place until the key is released. The pumping duration is saved.
- Short pressing (\leq 1 s): The aspiration takes place during the saved pumping duration. Press the key again to stop the procedure prematurely.

9 Maintenance

9.1 Maintenance

Regularly perform maintenance work on the product to prevent malfunctions and to ensure a long service life.

- Metrohm recommends having the products maintained by specialist personnel of Metrohm AG as part of an annual service. Shorter maintenance intervals may be necessary if you frequently work with caustic and corrosive chemicals.
- Only perform maintenance work that is described in this instruction. Contact your regional Metrohm service representative for further maintenance work and repairs. The regional Metrohm service representative offers every form of technical advice for maintenance and service of all Metrohm products.
- Only use spare parts that meet the technical requirements of the manufacturer. Original spare parts always meet these requirements.

9.2 Cleaning the product surface

Regularly clean the product to prevent malfunctions and to ensure a long service life.

- Remove spilled chemicals immediately.
- Protect plug connections against contamination.



WARNING

Chemical hazardous substances

Contact with aggressive chemical substances may cause poisoning or chemical burns.

- Wear personal protective equipment (e.g. protective glasses, gloves).
- Use exhaust equipment when working with vaporizing hazardous substances.
- Clean contaminated surfaces.
- Only use detergents that do not cause any unwanted side reactions with the materials to be cleaned.
- Dispose of chemically contaminated materials (e.g. cleaning material) in accordance with regulations.



WARNING

Health hazards from electrical potential.

Severe injuries with possibly fatal consequences.

- Operate the product only if it is in perfect condition. The housing must also be intact.
- Only use the product with the covers fitted.
- Protect live components (e.g. power supply unit, power cord, connection sockets) against moisture.
- Always have maintenance work and repairs on electrical components carried out by a regional Metrohm service representative.

Prerequisite:

- The product is switched off and disconnected from the energy supply.

Required accessories:

- Cleaning cloth (soft, lint-free)
- Water or ethanol

1 Clean the surface with a damp cloth. Remove persistent contamination with ethanol.

2 Wipe the surface with a dry cloth.

3 Clean the connectors with a dry cloth.

10 Troubleshooting

Messages on malfunctions and errors are displayed in the control software or in the embedded software (e.g. on the display of an instrument) and contain the following information:

- Descriptions of causes of malfunctions (e.g. jammed drive)
- Descriptions of problems with the control (e.g. missing or invalid parameter)
- Information on how to solve the problem

System components with status display elements also indicate malfunctions and errors with a red flashing LED.

Troubleshooting on the product is often only possible with the control software or the embedded software (e.g. initializing, moving to a defined position).

See also

System – Signals (chapter 3.4, page 16)

10.1 Karl Fischer titration


Problem	Cause	Remedy
The drift is very high during conditioning.	The titration cell is leaking.	<ul style="list-style-type: none"> ▪ Check the seals and the septum. Replace if necessary. ▪ Replace the molecular sieve.
The drift becomes greater after each titration.	The sample releases water very slowly.	<ul style="list-style-type: none"> ▪ Adjust the method. ▪ Add solubility promoter. ▪ Increase the temperature (possibly using a KF oven). ▪ See technical literature.
	A side reaction is taking place.	<ul style="list-style-type: none"> ▪ Use special reagents. ▪ Adjust the method (increase/decrease the temperature, external extraction). ▪ See technical literature.
	The pH value is no longer in the optimum range.	Add buffer (see technical literature).



Problem	Cause	Remedy
The titration does not finish.	The titration cell is leaking.	<ul style="list-style-type: none"> Check the seals and the septum. Replace if necessary. Replace the molecular sieve.
	The stop criterion is unsuitable.	Adjust the control parameters (see manual/help of the software used): <ul style="list-style-type: none"> Increase the stop drift. Select a short delay time.
	See also: The drift becomes greater after each titration.	
The sample is overtitrated.	The amount of methanol in the working medium is too low.	<ul style="list-style-type: none"> Replace the working medium. Reduce the amount of solubility promoter, if working with solvent mixtures, see technical literature.
	The electrode may be covered.	Wipe off the electrode with ethanol or a suitable solvent.
The solution becomes darker after each titration.		Replace the working medium.
	The electrode may be covered.	Wipe off the electrode with ethanol or a suitable solvent.
	The electrode has a short circuit.	<ul style="list-style-type: none"> Check the Pt wires. Activate the electrode check.
The endpoint is reached too quickly.	The dosing rate outside the control range is too high.	Select the user-defined titration rate and reduce the dosing rate (see manual/help of the software used).



11 Disposal



Properly dispose of chemicals and of the product to reduce negative effects on the environment and public health. Local authorities, waste disposal companies or dealers provide more detailed information on disposal. Observe the WEEE EU directive (WEEE = Waste Electrical and Electronic Equipment) for the proper disposal of waste electronic equipment within the European Union.



12 Technical specifications

12.1 Ambient conditions

Nominal function range	+5 to +45 °C	at max. 80% relative humidity, non-condensing
Storage	+5 to +45 °C	at max. 80% relative humidity, non-condensing
Altitude / Pressure range	Max. 3,000 m.a.s.l. / min. 700 mbar	
Overvoltage category	II	
Pollution degree	2	

12.2 OMNIS Coulometer Module – Energy supply

Nominal voltage	24 VDC
------------------------	--------

12.3 OMNIS magnetic stirrer – Dimensions

Measurements

<i>Width</i>	142 mm
<i>Height</i>	70 mm
<i>Depth</i>	116 mm

Weight

0.7 kg

12.4 OMNIS Coulometer Module – Dimensions

Measurements

<i>Width</i>	142 mm
<i>Height</i>	358 mm
<i>Depth</i>	
Without magnetic stirrer	284 mm
With magnetic stirrer	400 mm

Weight

<i>Type</i>	
Without magnetic stirrer	4.0 kg
With magnetic stirrer	4.7 kg

12.5 Magnetic stirrer – Dimensions

Measurements

<i>Width</i>	142 mm
<i>Height</i>	70 mm
<i>Depth</i>	116 mm

Weight	700 g
---------------	-------

12.6 OMNIS Coulometer Module – Housing

Materials

<i>Lid</i>	PET	Poly(ethylene terephthalate)
<i>Back panel</i>	AW-5754 H12/H22	Aluminum, coated
<i>Base</i>	1.4301	Stainless steel
<i>Enclosure</i>	PBT	Poly(butylene terephthalate)

Magnetic stirrer – Housing



Front foils PET Poly(ethylene terephthalate), mat

IP degree of protection IP 40

12.7 Magnetic stirrer – Housing

Materials

Cover PBT Poly(butylene terephthalate)

Base Chromium steel sheet

Enclosure PBT Poly(butylene terephthalate)

Front foils PET Poly(ethylene terephthalate)

IP degree of protection IP 40

12.8 OMNIS Coulometer Module – Connectors specifications

Energy supply via MDL
Socket Round plug 6-pin, size 1, 0° (ODU MINI-SNAP)

MDL Metrohm Device Link 1 connector

Internal measuring interface

INPUT 1

Socket round plug 7-pin, size 0, 45°

Potentiometric pH, ISE, Redox Measuring input for potentiometric electrodes



Temperature	Temp.	Measuring input for temperature sensors of the Pt1000 or NTC type for automatic temperature compensation
Polarizer	Pol.	Measuring input for polarizable electrodes
<i>INPUT 2</i>		
Socket		round plug 7-pin, size 0, 45°
Potentiometric	pH, ISE, Redox	Measuring input for potentiometric electrodes
Temperature	Temp.	Measuring input for temperature sensors of the Pt1000 or NTC type for automatic temperature compensation
<i>GENERATOR</i>		
	Socket	Round plug 2-pin
	Current generator output	For reagent generation

12.9 Display specifications

Status display	LED	multi-colored
----------------	-----	---------------

12.10 Current generator specifications

Low voltage generator (for bromine 1492)

<i>Current range</i>	0.5– 60.0 mA
<i>Voltage range</i>	0.0– 29.0 V

High voltage generator (for KFC water and BRC bromine index)

<i>Current range</i>	50.0– 400.0 mA
<i>Low voltage range</i>	0.0– 29.0 V
<i>High voltage range</i>	0.0– 39.0 V



Iodine production for Karl Fischer water determination

<i>Determination range</i>	0.01– 200.0 mg H ₂ O	Recommended water quantity
<i>Resolution</i>	0.1 µg H ₂ O	
<i>Titration rate</i>	max. 2.24 mg H ₂ O/min	
<i>Reproducibility</i>	±3 µg H ₂ O	at 10 µg– 1000 µg H ₂ O
		Sample: Standard from the reagent manufacturer
	≤0.3%	>1000 µg H ₂ O

12.11 Measurement specifications

Potentiometric

<i>Measuring range</i>	-2,400 to +2,400 mV -13 to +20 pH	
<i>Resolution</i>	1.56 µV 0.001 pH	
<i>Measuring accuracy</i>	±0.5 mV ±0.003 pH	in the measuring range -2,000 mV to +2,000 mV
<i>Input resistance</i>	≥ 1*10 ¹² Ω	
<i>Offset current</i>	≤ ±1*10 ⁻¹² A	

Temperature

<i>Pt1000</i>		
Measuring range	-150 to +250 °C	
Resolution	approx. 0.002 °C	
Measuring accuracy	±0.4 °C	in the measuring range -20.0 to +150.0 °C
<i>NTC 30 kOhm</i>		
Measuring range	-5 to +250 °C	
Measuring resolution	approx. 0.002 °C	





Measuring accuracy	± 0.6 °C	in the measuring range +10.0 °C to +40.0 °C
--------------------	--------------	--

Polarizer

I_{pol DC}

Polarization current	-200.0 to +200.0 μ A	adjustable in 0.5 μ A steps
Measuring range	-2,400 to +2,400 mV	
Measuring resolution	0.1 mV	

I_{pol AC}

Polarization current	5 μ A, 10 μ A, 20 μ A, 30 μ A	Effective values
Measuring range	0—+1,700 mV	Effective value
Measuring resolution	0.1 mV	Effective value
Frequency	10 Hz	

Upol DC

<i>Polarization voltage</i>	-2,000 mV to +2,000 mV	adjustable in 5 mV steps
<i>Measuring range</i>	-200.0 μ A to +200.0 μ A	
<i>Measuring resolution</i>	0.01 μ A	

Load of the measuring input I_{pol}

<i>R_{L max.} ±10μA</i>	240 k Ω
<i>R_{L max.} ±50μA</i>	48 k Ω
<i>R_{L max.} ±100μA</i>	24 k Ω

Load of the measuring input Upol

<i>R_{L min.} ±300 mV</i>	1.5 k Ω
<i>R_{L min.} ±600 mV</i>	3 k Ω
<i>R_{L min.} ±1,000 mV</i>	5 k Ω

Measuring accuracy



applies for all measuring ranges without sensor error, under reference conditions, measuring interval 100 ms

Reference conditions

<i>Relative humidity</i>	≤ 60%
<i>Ambient temperature</i>	+25 °C (±3 °C)
<i>Instrument status</i>	

min. 30 minutes in operation

12.12 Magnetic stirrer – Specifications

Adjustment range for rotational speed	+1 to +15	Rotation in counter-clockwise direction (seen from above)
	-1 to -15	Rotation in clockwise direction (seen from above)
Rotational speed change per step	120 rpm	
Maximum rotational speed	1,800 rpm	
Lengths of stirring bar	8, 12, 16, 25, 30 mm	