Metrohm IC Driver for Empower™ Instructions for Use

Handbuch
8.102.8050DE / 2013-12-06

Metrohm AG
CH-9100 Herisau
Switzerland
Phone +41 71 353 85 85
Fax +41 71 353 89 01
info@metrohm.com
www.metrohm.com

Metrohm IC Driver for Empower™ Instructions for Use

Handbuch

Teachware
Metrohm AG
CH-9100 Herisau
teachware@metrohm.com

This documentation is protected by copyright. All rights reserved.

Although all the information given in this documentation has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please send us your comments using the address given above.

Change Control

Version	Date	Summary of Changes
1.0	July 2013	First Edition
1.1	December 2013	Addition of new IC devices 930 / 940

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IC – MSM

IC	Sample Processor	Dosino 59290 (2 ml)	Main Time Program				
Pump	Injector	MSM	Peristaltic Pump	Degasser	MCS	Column Oven	Conductivity Detector
<input checked="" type="checkbox"/> Active							
<input checked="" type="checkbox"/> Automatic stepping to next position during equilibration							
Interval	<input type="text" value="10.0"/> min						

IC – Peristaltic Pump

IC	Sample Processor	Dosino 59290 (2 ml)	Main Time Program				
Pump	Injector	MSM	Peristaltic Pump	Degasser	MCS	Column Oven	Conductivity Detector
<input checked="" type="checkbox"/> Active							
<input checked="" type="checkbox"/> On							
Rate	<input type="text" value="1"/>						

IC – Degasser

IC	Sample Processor	Dosino 59290 (2 ml)	Main Time Program				
Pump	Injector	MSM	Peristaltic Pump	Degasser	MCS	Column Oven	Conductivity Detector
<input checked="" type="checkbox"/> Active							
<input checked="" type="checkbox"/> On							

IC – MCS

IC	Sample Processor	Dosino 59290 (2 ml)	Main Time Program				
Pump	Injector	MSM	Peristaltic Pump	Degasser	MCS	Column Oven	Conductivity Detector
<input checked="" type="checkbox"/> Active							
<input checked="" type="checkbox"/> On							

IC – Column Oven

IC	Sample Processor	Dosino 59290 (2 ml)	Main Time Program				
Pump	Injector	MSM	Peristaltic Pump	Degasser	MCS	Column Oven	Conductivity Detector
<input checked="" type="checkbox"/> Active							
<input checked="" type="checkbox"/> On							
<input type="checkbox"/> Wait for stable temperature							
Temperature <input type="text" value="30.0"/> °C							

IC – Detector

IC	Sample Processor	Dosino 59290 (2 ml)	Main Time Program				
Pump	Injector	MSM	Peristaltic Pump	Degasser	MCS	Column Oven	Conductivity Detector
<input checked="" type="checkbox"/> Active							
Temperature Coefficient <input type="text" value="2.3"/> %/°C							
Warning Limit <input type="text" value="9999"/> µS/cm							
Thermostat <input type="text" value="40"/> °C							

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

Sample Processor – Peristaltic Pump

IC	Sample Processor	Dosino 59290 (2 ml)	Main Time Program
Tower	Peristaltic Pump		
<input checked="" type="checkbox"/> Active			
<input type="checkbox"/> On			
Rate <input type="text" value="3"/>			

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



IC Sample Processor Dosino 59290 (2 ml) Main Time Program					
Time	Device	Module	Command	Parameter	
0.00	Compact IC pro - Ani...	Injector	Switch	SwitchTo = Fill	
---	Professional Sample ...	Tower	Move rack	Move = 0, Numbe...	
---	Professional Sample ...	Tower	Lift	LiftPositionSelecti...	
0.00	Professional Sample ...	Peristaltic P...	On/Off	On = True, Rate ...	
2.00	Professional Sample ...	Peristaltic P...	On/Off	On = False, Rate ...	
2.00	Compact IC pro - Ani...	Injector	Switch	SwitchTo = Inject	
2.00	Compact IC pro - Ani...	Conductivity...	Start Data Acquisition		

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#	Name	Peak Label	Retention Time (min)	RT Window (min)	Peak Match	Channel	Y Value	X Value	Fit	Weighting	Internal Std	RT Reference	Rel RT Reference	Rel Resol Reference
1	Fluoride		4.263	0.213	Closest		Area	Amount	Quadratic	None				
2	Chloride		5.138	0.257	Closest		Area	Amount	Quadratic	None				
3	Nitrite		5.849	0.292	Closest		Area	Amount	Quadratic	None				
4	Bromide		7.600	0.360	Closest		Area	Amount	Quadratic	None				
5	Nitrate		8.600	0.408	Closest		Area	Amount	Quadratic	None				
6	Sulfate		10.381	0.519	Closest		Area	Amount	Quadratic	None				
7	Phosphate		11.399	0.570	Closest		Area	Amount	Quadratic	None				

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

#	Time (min)	Type	Value
1	0.000	Inhibit Integration	
2	0.000	Valley to Valley	
3	11.700	Inhibit Integration	

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


Status	Manual Control
Instrument Status:	IDLE
Pump Flow:	0.000 mL/min
Pump Pressure:	0.0 bar
Injector Position:	INJECT
Degasser State:	READY
MCS State:	READY
Last Stepping:	252.1 min
Column Heater Temperature:	25.1 °C

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The screenshot shows the 'Run Samples' window in the Anion software. The main window contains a table with columns: Vial, Inj Vol (uL), # of Injs, Label, SampleName, Level, Function, Method Set / Report Method, Label Reference, Processing, Run Time (Minutes), Data Start (Minutes), Next Inj. Delay (Minutes), SampleWeight, and Dilution. A dropdown menu is open over the 'Function' column, showing options like 'Inject Samples', 'Inject Controls', 'Inject RF Internal Standards', etc. Below the table, there are status indicators for 'Sample Set Time Remaining' and 'Total Samples Time Remaining'. A 'Sample Set' dropdown is also visible. At the bottom, a 'Status Manual Control' panel displays various instrument parameters such as 'Instrument Status: IDLE', 'Pump Flow: 0.000 mL/min', 'Pump Pressure: 0.0 bar', 'Injector Position: INJECT', 'Degasser State: READY', 'MCS State: READY', 'Last Stepping: 254.6 min', 'Column Heater Temperature: 25.1 °C', 'Column Heater Temp. OK: INSTABLE', 'Detector Conductivity: ---- µS/cm', and 'Rack Name: 6.2041.430'. A chromatogram plot is also visible, showing a flat line at 0.00 AU over a 1.00 minute period.

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
Run Sample Set

Name for this sample set : Anion smpleset

Sample set method name : Anion smpleset

Settings for this Sample Set

Wait For User

Run Mode : Run and Process


Suitability Mode : Continue on Fault

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Abort Now!

Abort after Vial is completed.

Abort after Injection is completed.

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

The screenshot shows a software configuration window for a method set. On the left, a tree view shows a folder named "Method Set" containing sub-items "Data Channels" and "Derived Channels". The main area contains several dropdown menus and a table.

Instrument Method: Anion equilibrate instmet

Default Processing Method: Anion procmct


Default Report Method: [Empty]

Channel Name	Processing Method	Report Method

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nsf_022013 - Background Processing and Reporting

Processing

Process

Use acquisition method set (i.e. from the sample set used to acquire data)

Use specified method set

Use specified processing method An1 procmet

Clear Calibration How: Calibrate and Quantitat

Reporting

Print \\srv-print\Farblaser-ICLAB-1

Use acquisition method set (i.e. from the sample set used to acquire data)

Use specified method set

Use specified report method

Exporting

Export

Use acquisition method set (i.e. from the sample set used to acquire data)

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



E	Name	Retention Time (min)	Area (μV*sec)	% Area	Height (μV)	Int Type	Amount	Units	Peak Type
1	Fluoride	2.118	17	50.85	2	BB	0.759	mg/L	Found
2	Chloride	2.675	11	32.27	1	BB	0.755	mg/L	Found
3	Nitrate	3.905	6	16.88	1	BB	0.785	mg/L	Found

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The screenshot shows the software interface with a data table and a chromatogram. The table has the following data:

SampleName	Vial	Injection	Sample Type	Processed Channel Descr.	Date Acquired	Date Processed	Processing Method	Result Id
Standard2	24	1	Standard	Conductivity Detector	06.03.2013 09:25:34 CET	06.03.2013 13:21:06 CET	An1 procmet	2505
Standard1	25	1	Standard	Conductivity Detector	06.03.2013 09:33:02 CET	06.03.2013 13:21:06 CET	An1 procmet	2506
Standard3	26	1	Standard	Conductivity Detector	06.03.2013 09:40:28 CET	06.03.2013 13:21:06 CET	An1 procmet	2507
Check 0.8mg/L	27	1	Unknown	Conductivity Detector	06.03.2013 09:47:54 CET	06.03.2013 13:21:06 CET	An1 procmet	2508

The chromatogram shows a red line representing the signal from the Conductivity Detector. The x-axis is labeled 'Minutes' and ranges from 0.00 to 10.00. The y-axis is labeled 'µS/cm' and ranges from 0.00 to 3.00. There are three distinct peaks at approximately 2.5, 3.5, and 4.5 minutes.

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
The 'Open Report Method' dialog box contains the following text and options:


Please select the Report Method that you would like to use to preview the data that you have selected.

- Use the Report Method in the acquisition Method sna 4uLmetset.
- Use the Report Method named Default.
- Use a Report Method that was generated to be appropriate for the selected data.
- Use the following Report Method: [Dropdown menu]
- Use the currently open Report Method named Unbenannt.

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
	SampleName	Sample Type	Vial	Inj #	Run Time (Minutes)	Injection Volume (ul)	Acquisition Method Set	Sample Weight	Processed Channel Descr.	Dilution
1	Check 0.8mg/L	Unknown	27	1	5.00	20.00	Anion metset	1.00000	Conductivity Detector	1.00000

Peak Results




	Name	RT	Area	Height	Amount	Units
1	Fluoride	2.118	17	2	0.759	mg/L
2	Chloride	2.675	11	1	0.755	mg/L
3	Nitrate	3.905	6	1	0.785	mg/L

Metrohm IC Post Run Report

Sample: Anionen Vial: 1 Inj: 1

IC Properties

Device Type 0940.1500
 Program Version 5.940.0100
 Serial Number 2106

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
The screenshot shows a software interface for configuring a method set. On the left is a tree view with the following structure:

- Method Set
 - Data Channels
 - Conductivity Detector
 - Derived Channels
 - cation

At the top right, there are three dropdown menus:

- Instrument Method: Conductivity Cation instmet
- Default Processing Method: (empty)
- Default Report Method: (empty)

Below these is a table with the following data:

	Channel Name	Processing Method	Report Method
1	cation	Conductivity Cation procmeth	