


Mira Cal



Tutorial
8.105.8020EN



Metrohm AG

CH-9100 Herisau

Switzerland

Phone +41 71 353 85 85

Fax +41 71 353 89 01

info@metrohm.com

www.metrohm.com

Mira Cal

Tutorial

Technische Dokumentation
Metrohm AG
CH-9100 Herisau
techdoc@metrohm.com

This documentation is protected by copyright. All rights reserved.

This document has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Table of contents

1	Introduction	1
1.1	Structure of the tutorial	1
1.2	Program description	1
1.3	Spectral libraries	2
1.4	Symbols and conventions	3
2	Preparation	4
2.1	Install the software	4
2.2	Start the software	5
2.3	Establish connection	5
2.4	Firmware Update	7
2.5	Update S.T. Japan Library	8
3	Software description	10
3.1	Overview	10
3.2	Menu items	11
3.3	Spectra view	12
3.4	Sample list	14
3.4.1	Save spectra	14
3.4.2	Load spectra	15
3.4.3	Manipulate spectra	16
3.4.4	Generate reports	16
4	Data acquisition	18
4.1	Data acquisition with Mira	18
4.2	Data acquisition with Mira Cal	19
5	Spectral libraries	22
5.1	Create and manage custom libraries	22
5.1.1	Create and load local libraries	22
5.1.2	Modify libraries	24
5.1.3	Upload library to spectrometer	26
5.2	Activate S.T. Japan libraries	27
5.3	Enable or disable libraries for search	28
6	Troubleshooting	30
6.1	Support Request	30



7 System requirements	31
Index	32

Table of figures

Figure 1	Menu item „Connect“	6
Figure 2	Dialog box „Device Connection“	6
Figure 3	Menu item „Disconnect“	7
Figure 4	Menu item „Update Firmware“	8
Figure 5	Dialog box „Update Device Firmware“	8
Figure 6	Menu item „Upload Encrypted Library“	9
Figure 7	Overview graphic	10
Figure 8	Menu item „Acquisition“	11
Figure 9	Menu item „Library Management“	11
Figure 10	Context menu	12
Figure 11	Zoom	12
Figure 12	Selection „Detail Pane“	13
Figure 13	Window „Zoom“	13
Figure 14	Selection „Overlay“	13
Figure 15	Selection „Legend“	14
Figure 16	Context menu	15
Figure 17	Context menu	15
Figure 18	Context menu	16
Figure 19	Context menu	18
Figure 20	Dialog box „Delete Data“	19
Figure 21	Arming Laser	20
Figure 22	Starting the measurement	20
Figure 23	Disarm Laser	21
Figure 24	Button „New“	23
Figure 25	Button „Open“	24
Figure 26	Context menu	25
Figure 27	Dialog box „Group directory“	25
Figure 28	Context menu	26
Figure 29	Context menu	26
Figure 30	Menu item „Upload Access File“	27
Figure 31	Dialog box „Device Library“	29
Figure 32	Context menu	29

1 Introduction

1.1 Structure of the tutorial

The present tutorial guides you through your first steps using the **Mira Cal** software.

The tutorial concentrates on four parts:

- Prepare the software
- A software description
- Perform measurements
- Create and manage spectra libraries

1.2 Program description

With the Mira Cal software you can conveniently store, manage and exchange data acquired on a Mira spectrometer. The software allows not only to activate/deactivate purchased libraries, but also the data transfer between the Mira spectrometer and the computer and also between different Mira spectrometers.

Using the software you can build your own library and save it on the instrument. This “open library” principle allows you to use the Mira spectrometer in the most flexible way possible.

Available Spectrometer

The following spectrometer can be used with Mira Cal:

- Mira M-1 Basic (User Manual: 8.923.8001)
- Mira M-1 Advanced (User Manual: 8.923.8001)
- Mira M-2 Basic (User Manual: 8.923.8002)
- Mira M-2 Advanced (User Manual: 8.923.8002)

Communication

Communication between the Mira spectrometer and the host computer is established by a USB Mini-B cable.



1.3 Spectral libraries

In addition to the option of creating custom libraries, the following libraries are additionally available and can be purchased:




Table 1 Spectral libraries

Article number	Designation	Number of spectra
6.6071.601	Complete Raman spectral library.	> 8690
6.6071.602	Raman spectra of active substances and auxiliary materials that are relevant to the pharmaceutical industry and medical research.	> 1170
6.6071.603	Raman spectra of solvents.	> 460
6.6071.604	Raman spectra of polymers, polymer additives, plastics, plasticizers and packaging materials.	> 920
6.6071.605	Raman spectra of aliphatic and aromatic aldehydes and ketones.	> 1070
6.6071.606	Raman spectra of alcoholic and phenolic compounds.	> 890
6.6071.607	Raman spectra of esters, lactones and anhydrides.	> 2930
6.6071.608	Raman spectra of hydrocarbons and halogenated hydrocarbons.	> 560
6.6071.609	Raman spectra of chemical substances that are used in the semiconductor industry.	> 370
6.6071.610	Raman spectra of selected hazardous substances that are listed in the "EPA Cameo Database for Chemical Emergencies and Responders" and the "USCG CHRIS Hazardous Chemicals Database".	> 1360
6.6071.611	Raman spectra of selected hazardous substances that are listed in the "EPA Cameo Database for Chemical Emergencies and Responders", "USCG CHRIS Hazardous Chemicals Database", and "NIOSH Guide to Chemical Hazards Databases", as well as chemicals that are regulated by the "Toxic Substances Control Act"	> 3030
6.6071.612	Raman spectra of substances that are relevant to forensic analysis.	> 740
6.6071.613	Raman spectra of pesticides, insecticides, herbicides, fungicides, algicides and similar agricultural chemicals.	> 460
6.6071.614	Raman spectra of selected dyes, colorants, pigments and indicators.	> 300
6.6071.615	Raman spectra of sulfur and phosphorus compounds.	> 970
6.6071.616	Raman spectra of substances with a high production volume, as listed in the "HPV Challenge Program Chemical List".	> 690
6.6071.620	Raman spectra of food additives, including FDA-controlled substances. Additionally, spectra of indirect food additives and substances that come into contact with foodstuffs, such as packaging materials and associated processing chemicals.	> 1070

Article number	Designation	Number of spectra
6.6071.621	Raman spectra of biochemicals, including vitamins, resins, starches, glycerins, fatty acids, sugars, carbohydrates, proteins, and peptides.	> 1900
6.6071.622	Raman spectra of flavors, scents, and other substances that are used for manufacturing cosmetics.	> 1030

1.4 Symbols and conventions

The following icons and formatting are used in this documentation:

	Instruction step Carry out these steps in the sequence shown.
Method	Dialog text, parameter in the software
File ▶ New	Menu or menu item
[Continue]	Button or key
	WARNING This symbol draws attention to a possible life-threatening hazard or risk of injury.
	Note This symbol highlights additional information and tips.



2 Preparation

2.1 Install the software

Check System requirements

- 1 Make sure, the used system fulfils the system requirements (*see Chapter 7, page 31*).

Check Windows User rights

- 1 Make sure, to have administrator rights before installing Mira Cal.
- 2 Make sure, that User rights settings are set to „Always notify“ or „Notify me only when programs try to make changes to my computer“.



NOTE

This step depends on the installed operating system and gis described here:

<https://support.microsoft.com/en-us/kb/975787>

Install the software

- 1 Double-click on „Mira_Cal.exe“.
- 2 Click on **[Next]** and follow the instructions on the screen.

Install spectrometer driver.

1



NOTE

As soon as the spectrometer connects for the first time, the driver for the spectrometer is installed automatically. Therefore an internet connection is needed.

2

If the installation was not successful:

- **A:** Open Mira Cal (see Chapter 2.2, page 5).
- **B:** Go to **File ▶ Advanced Updates ▶ Install ▶ Device Drivers**.

2.2 Start the software

1



NOTE

As soon as the software is installed, Mira Cal icon appears on the desktop of the computer.

Double-click on , to start the program.

2.3 Establish connection

Physically connect the spectrometer to the host

1

Connect the spectrometer to the PC using the USB Mini-B cable.

2

Turn on the instrument and enter the PIN code (default: 1234) (see instrument user manual "Spectrometer variants", Page 1).



Connect the spectrometer with Mira Cal




Figure 1 Menu item „Connect“

In Mira Cal choose Menu item **File ▶ Connect**.

2 Click on **[Refresh]**.

The spectrometer is listed in „Available COM Devices“.



NOTE

If the spectrometer is not listed, repeat step 2 until it is listed.







Figure 2 Dialog box „Device Connection“

Select the desired spectrometer and click on **[Select]**.

**NOTE**

A successful connection is indicated by the icon   in Mira Cal and the icon  on the screen of the spectrometer.

Disconnect the spectrometer

1




Figure 3 Menu item „Disconnect“

Click on **File ► Disconnect**, to disconnect the spectrometer.

2 Disconnect the USB Mini-B cable from the computer.

2.4 Firmware Update**NOTE**

After a Firmware Update, the instrument must be calibrated again. Information regarding calibration can be found in the user manual of the instrument.

**NOTE**

When you installed a new version of Mira Cal, make sure to update the instruments firmware as well.

Load a library to the spectrometer

Make sure, a spectrometer is connected.




Figure 6 Menu item „Upload Encrypted Library“

Click on **File ▶ Advanced ▶ Updates ▶ Upload Encrypted Library**.

- 2 Select the .zip file with the new library and click on **[Open]**.

The library will be loaded to the spectrometer. Subsequently the spectrometer will reboot automatically.



NOTE

The reboot can last several minutes. Do not disconnect the spectrometer within the reboot.

- 3 Activate library (see Chapter 5.2, page 27).

Sample list	<p>The sample list contains information about the collected samples (for example date and time).</p> <p>Spectra are divided in groups and can be opened by clicking on the triangle ▸ left of the group name.</p> <p>Possible matches with a library, can be opened by clicking on the triangle ▸ left of the spectrum name.</p>
--------------------	--

3.2 Menu items

Mira Cal has two different program parts which can be opened by clicking on the corresponding menu item symbol in the vertical bar on the left.

The menu item symbol for the opened program part is shown in white with black icons. The symbols for the other program part in green with white icons.

The menus, tool bars and content of the main window depend on the menu item currently opened.

Menu item „Acquisition“




Figure 8 Menu item „Acquisition“

- Spectra view (with **Overlay** and **DetailPane**).
- Acquire samples with remote controlling via Mira Cal.
- Match measured spectra to device library.
- Save and load acquired spectra to and from a local drive.
- Generate reports.

Menu item „Library Management“




Figure 9 Menu item „Library Management“

- Download libraries from Mira spectrometer to PC and import measured spectra.
- Create custom libraries.
- Add/remove substances from libraries.



- Upload changes to the connected Mira spectrometer.
- Manage libraries of connected Mira spectrometer.

3.3 Spectra view

Spectra view is a tool for visual inspection of acquired spectra. Spectra View can be used to load and display spectra in different formats (*.rmn, *.spc, *.txt).

Copy displayed spectrum to clipboard

The currently displayed spectrum can be copied to the clipboard.

1 Right-click into the Spectra View.




Figure 10 Context menu

Click on **Copy Image** in the dialog box.

Zoom and DetailPane

Zoom

1 Click and hold in Spectra view with left mouse key.




Figure 11 Zoom

Click and hold left mouse key, draw a rectangle and release left mouse key.

- 3 Double-click in Spectra View to get back to original size.

DetailPane


- 1 

Figure 12 Selection „Detail Pane“

Mark check box „DetailPane“.





Figure 13 Window „Zoom“

The opening window zooms into displayed spectrum.



NOTE

The same zoom scales can be used in different areas of a spectrum.

- 3 Click on  to get back to original size.

Overlay spectra


- 1 

Figure 14 Selection „Overlay“

Mark Check box „Overlay“ that spectra can be overlaid.

Spectra can be selected to display the names.

**NOTE**

A multiple selection of several groups is possible.




Figure 16 Context menu

Click on **Save As...** in the opening context menu.

- 3 In the dialog box select the file format and file location and click on **[Save]**.

3.4.2 Load spectra

Spectra in the above file formats can be loaded for visual inspections and transfer to libraries.

Load spectra from the local drive

- 1 Right-click on the Sample list pane.




Figure 17 Context menu

Choose **Load spectra ► Local data** from the context menu.

- 3 From the appearing dialog box choose the spectra to be loaded and click on **[Open]**.

4 Delete Data

Delete data saved on device?



Figure 20 Dialog box „Delete Data“

Click on **[Delete]** to delete Spectra from spectrometer or click on **[Keep]** to keep data on the spectrometer.

The spectra can be found in a new group called “Import from Instrument”.

4.2 Data acquisition with Mira Cal

Acquire data

- 1 Connect computer to an instrument (*see Chapter 2.3, page 5*).
- 2 Select menu item „Acquisition“.
- 3 Position the sample for measurement.



NOTE

Refer to the Mira user manual (*see "Available Spectrometer", page 1*) for detailed description.



4



WARNING

Risk of injury by laser radiation

The Mira Advanced spectrometer are class 3B laser which can emit laser light during data acquisition.

The laser radiation can cause serious eye injuries if the safety measures are not followed.

- The provided protective glasses must be worn when using an attachment lens.
- The safety regulations must be observed (refer to Mira user manual (see "Available Spectrometer", page 1)).

When the laser is armed, pressing **[Auto Acquire]** will fire it. Arming the laser is not required before each measurement.



Figure 21 Arming Laser

1 Button „Arm Laser“

When using a Mira Advanced, click on **[Arm Laser]** (21-1) to arm the laser.

5




Figure 22 Starting the measurement

2 Button „Auto Acquire“

Click on **[Auto Acquire]** (22-2) to start measurement.

The sample is now listed in the group „Recent Acquisitions“.




Figure 23 Disarm Laser

3 Button „Laser Armed“

Click on **[Laser Armed]** (23-3) to disarm the laser.



5 Spectral libraries



NOTE

A spectrometer must be connected in order to perform the following operations.

5.1 Create and manage custom libraries



NOTE

Custom libraries are built from measurements of user samples. It is of utmost importance that these samples are representative of the statistical population. Failing to select appropriate samples will lead to poor performance of the library matching. Consider measuring many samples and average the spectra.

5.1.1 Create and load local libraries

Create a new library

- 1 Select menu item „Library Management“.




Figure 24 Button „New“

In the „Local Library“ section click on **[New]**.

3 Rename the new library.

Load a library

1 Select menu item „Library Management“.




Figure 26 Context menu

Click on **Copy to...** in the context menu.




Figure 27 Dialog box „Group directory“

Select the library you want to add the spectrum to and click on **[Copy]**.

Remove samples from a library

- 1 Select menu item „Library Management“.
- 2 Select the desired library with the sample you want to delete.
- 3 In the Sample list right-click on the sample you want to remove.



Figure 28 Context menu

Click on **Remove** in the context menu.

Save library local

- 1 Select menu item „Library Management“.
- 2 In „Local Library“ right-click on the library you want to save.



Figure 29 Context menu

Click on **Save As...** in the opening context menu.

- 4 Choose the file location and click on **[Save]**.



NOTE

It is recommended to always save libraries in the .rmn format.

5.1.3 Upload library to spectrometer

- 1 Select menu item „Library Management“.
- 2 Right-click on the library you want to upload to the spectrometer.
- 3 Click on **Upload to Device** in the context menu.

5.2 Activate S.T. Japan libraries

Buy library

To purchase libraries, contact your local Metrohm representative. You will be provided with an activation file.

Activate library

The spectrometer must be connected.




Figure 30 Menu item „Upload Access File“

Click on **File** ► **Advanced** ► **Updates** ► **Upload Access File**.

2 Choose the file and click on **[Open]**.



5.3 Enable or disable libraries for search



NOTE

- The fewer libraries are activated for the spectral comparison, the shorter the search.
- If many libraries are activated, the sample may have a high match rate for several library substances and thus be difficult to identify. It is advisable to incorporate previous knowledge of the sample into the identification.
For this reason, the libraries are arranged according to subjects. We recommend activating only those substance classes that fit the sample.
- It is recommended not to use libraries with big and small data at the same time.
- The settings described in this chapter only determine whether a library is included in the search or not.
Newly purchased libraries first have to be activated in the Mira Cal.

Libraries can be enabled /disabled either directly on the spectrometer or using Mira Cal.

Enable or disable libraries on the spectrometer

- 1 Refer to the Mira user manual (*see "Available Spectrometer", page 1*) for a detailed description how to enable or disable libraries on the spectrometer.

Enable or Disable libraries using Mira Cal

- 1 Select menu item „Library Management“.
- 2 Click on **[List Libraries]**.

Device Library		List Libraries
Enabled	GroupName	
<input checked="" type="checkbox"/>	Raman_Demo	
<input checked="" type="checkbox"/>	Inorganics_Full_Range	
<input type="checkbox"/>	Inorganics_Partial_Range	
<input type="checkbox"/>	Minerals_Full_Range	
<input type="checkbox"/>	Minerals_Partial_Range	

Figure 31 Dialog box „Device Library“

The dialog box „Device Library“ opens.



NOTE

All libraries are listed. The checkbox in the „Enabled“ column indicates, whether libraries are enabled for search or not.

3 Right-click on a library.

4

Device Library		List Libraries
Enabled	GroupName	
<input type="checkbox"/>	Raman_Demo	
<input type="checkbox"/>	Inorganics Full Range	
<input type="checkbox"/>	In	
<input type="checkbox"/>	M	
<input type="checkbox"/>	Minerals_Partial_Range	

Enable/Disable

Download Spectra

Delete

Figure 32 Context menu

Click on **Enable** in the context menu, to enable the library for search or click on **Disable** to disable the library for search.

7 System requirements

<i>Processor</i>	Pentium 4 equivalent or higher
<i>Operating system 32 or 64 bit version of:</i>	<ul style="list-style-type: none">▪ Windows 8.1▪ Windows 8▪ Windows 7
<i>Screen resolution</i>	Resolution 1024 x 768 or higher, EGA equivalent or higher
<i>RAM</i>	All operating systems: 1 GB minimum; 4 GB recommended
<i>hard disk capacity</i>	Minimum 100 GB; 500 GB recommended
<i>USB connectors</i>	USB 2.0, two free ports recommended
<i>Mouse</i>	Required, USB, PS/2 style, or wireless
<i>Keyboard</i>	Compatible keyboard using USB / PS/2, or similar connection
<i>Printer compatibility</i>	Must be connected to one of the following: <ul style="list-style-type: none">▪ USB printer▪ Network printer within immediate work area▪ Alternative print connection, including parallel, or other reliable methods (A test page should be printed to verify correct operation).
<i>System backup</i>	Network or CD backup required for data archiving



Index

A

Acquisition 11

D

Data

 Acquire 18

 Transfer 18

Data acquisition 19

DetailPane 12

Driver 4

I

Installation 4

L

Legend 12

Library

 Activate 27, 28

 Buy 27

 Create 22

 Deactivate 28

 Load 22

 Modify 24

 Upload 26

Library Management 11

O

Overview 10

S

S.T. Japan 8, 27

Search 28

Spectra

 Copy 12

 Load 15

 Manipulate 16

 Overlay 12

 Save 14

Spectrometer

 Connect 5

 Disconnect 5

 Update 7

Structure 10

System Requirements 4

U

User rights 4

V

View 10

Z

Zoom 12