

Optrode M2



6.1125.000

Manual

8.109.8066EN / v1 / 2025-06-06



Metrohm AG
CH-9100 Herisau
Switzerland
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

Optrode M2

6.1125.000

Manual

Technical Communication
Metrohm AG
CH-9100 Herisau

This documentation is protected by copyright. All rights reserved.

This documentation is an original document.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Disclaimer

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g., conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

Table of contents

1	Introduction	1
1.1	Description	1
1.2	Symbols and conventions	1
2	Safety	3
2.1	General notes on safety	3
2.2	Intended use	3
2.3	Danger caused by breakage of glass	3
3	Overview	4
4	Installation	5
4.1	Example of connecting to a device with a USB connector	5
4.2	Example of connecting to a device with a USB Mini (OTG) connector	6
4.3	Energy supply with a USB power supply unit	6
4.4	Installing the Optrode M2	7
5	Operation and maintenance	9
5.1	Setting the wavelength	9
5.2	Sensor care	11
6	Troubleshooting	12
6.1	Possible problems and information on how to deal with them	12
7	Recycling and disposal	14
8	Technical specifications	15
8.1	Measuring mode	15
8.2	Measuring interval	15
8.3	Resolution	15
8.4	Ambient conditions	15
8.5	Energy supply	15
8.6	Reference conditions	16
8.7	Light sources (LEDs)	16



8.8	Dimensions	16
8.9	Material	16
	Index	17

Table of figures

Figure 1	Optrode M2 – Overview	4
Figure 2	Energy supply on a device with USB connector	5
Figure 3	Connect Optrode M2 to the indicator measuring input	5
Figure 4	Energy supply on a device with USB Mini (OTG) connector	6
Figure 5	Connect Optrode M2 to the indicator measuring input	6
Figure 6	Install the Optrode M2	7
Figure 7	Diagrams showing propeller stirrer, sensor and titration tip during a titration. a) clockwise stirring direction, b) counterclockwise stirring direction.	8
Figure 8	Wavelength display of the Optrode M2	9
Figure 9	Switching the LED	10

1 Introduction

1.1 Description

The **Optrode M2** is used as a photometric sensor for various titration methods that require the use of indicators. Its enclosed glass shaft enables it to be safely used in nonaqueous or corrosive media. It does not require any time-consuming care or conditioning.

Equipped with 8 LED light sources with different wavelength ranges, it can be used with a large variety of indicators. The sensor does not have a mechanical switch.


The LED light source can be selected either via the control software *tiamo* (version 2.5 and higher) or without the need for contact using the provided stirring bar.

The energy supply to the Optrode M2 is ensured via USB cable. Metrohm titrators are equipped with USB connectors for this purpose. In addition, a USB power adapter is available to power the Optrode M2 independently.

The Optrode M2 provides an analog measuring signal and can be connected reliably to any Metrohm titrator at the **Ind.** connector. This allows the Optrode M2 to be used for titration like other sensors.

1.2 Symbols and conventions

The following symbols and formatting may appear in this documentation:

(5-12)	Cross-reference to figure legend The first number refers to the figure number, the second to the instrument part in the figure.
1	Instruction step Perform the steps one after the other.
Method	Dialog text, parameter in the software
File ► New	Menu or menu item
[Continue]	Button or key
	WARNING This symbol draws attention to a possible life-threatening hazard or risk of injury.



WARNING

This symbol draws attention to a possible hazard due to electrical current.



WARNING

This symbol draws attention to a possible hazard due to heat or hot instrument parts.



WARNING

This symbol draws attention to a possible biological hazard.



WARNING

Warning of optical radiation



CAUTION

This symbol draws attention to possible damage to instruments or instrument parts.



NOTICE

This symbol highlights additional information and tips.

2 Safety

2.1 General notes on safety



WARNING

This sensor may be operated only in accordance with the information contained in this documentation.

2.2 Intended use

The Optrode M2 is designed for use in indoor environments (analytical laboratories). It is intended solely for use as an optical sensor for photometric titrations (together with a suitable titrator).

This sensor can be used in chemicals and flammable solvents. Usage of the Optrode M2 therefore requires the user to have basic knowledge and experience in handling toxic and caustic substances. Knowledge of the application of the fire prevention measures prescribed for laboratories or production operation is also mandatory.

2.3 Danger caused by breakage of glass

The Optrode M2 has a main body made of glass (*see chapter 8, page 15*).

If this is damaged, then the function may become impaired and there is a danger of cut injuries.



WARNING

Risk of cuts from sharp edges

Cut injuries due to damaged glass parts and pieces of broken glass.

- Handle glass parts (e.g. electrodes, bottles) with care.
- Only use undamaged glass parts.
- Dispose of damaged glass parts immediately.
- In the event of uncertainties, contact the regional Metrohm service representative.

3 Overview




Figure 1 Optrode M2 – Overview

1 USB plug For energy supply of the sensor.	2 Optical circuit board With wavelength display.
3 Glass shaft	4 Measuring gap Light path = 10 mm.
5 Photodiode For measuring the light intensity.	6 Storage vessel (6.2008.040)
7 Ground-joint sleeve	8 Sensor F plug For connecting to a titrator (connector Ind.).

4 Installation

Also refer to the manual for the respective instrument to connect it.

The Optrode M2's USB cable is used for the energy supply to the electronics and to the sensor's light source. The actual sensor cable (measuring signal) has a conventional F plug for connecting to the indicator measuring input (**Ind.**) of a titrator on the.

4.1 Example of connecting to a device with a USB connector




Figure 2 Energy supply on a device with USB connector




Figure 3 Connect Optrode M2 to the indicator measuring input



4.2 Example of connecting to a device with a USB Mini (OTG) connector




Figure 4 Energy supply on a device with USB Mini (OTG) connector

In order to connect the Optrode M2's USB cable, an adapter USB MINI (OTG) - USB A (6.2151.100) is required (see figure).




Figure 5 Connect Optrode M2 to the indicator measuring input

4.3 Energy supply with a USB power supply unit

A 6.02117.000 USB power supply unit is available if the power supply to the Optrode M2 cannot be provided via a titrator's USB port. This USB power supply unit can be connected to a power outlet directly.

The cable for energy supply to the Optrode M2 can be connected to the USB port on the USB power supply unit.

4.4 Installing the Optrode M2



NOTE

Ensure that the sensor is not exposed to direct sunlight and avoid disruptive light reflections.

Install the Optrode M2 as shown in the figure below.




Figure 6 Install the Optrode M2

The sensor must be firmly installed on the titration head. Particularly when automation is involved, it is important to ensure that the cable has sufficient play.

During the titration, it is important that the solution is mixed well. Select a stirring rate high enough to form a small vortex. If the stirring rate is too high, then air bubbles will be aspirated. This results in incorrect measured values. If the stirring rate is too low, then the solution at the sensor will not yet be mixed correctly.

In order for the measurement to be taken in a well-mixed solution after the addition of the titrant, position the antidiffusion tip where turbulence is high. Furthermore, the distance between the addition of the titrant and the sensor should be as large as possible. Take into account the stirring direction (counterclockwise or clockwise) when positioning the electrode and antidiffusion tip.

The sensor's optical window is to be positioned facing the flow direction. This removes small air bubbles from the measuring gap.




Figure 7 Diagrams showing propeller stirrer, sensor and titration tip during a titration. a) clockwise stirring direction, b) counter-clockwise stirring direction.

1	Propeller stirrer
3	Antidiffusion tip

2	Sensor
----------	---------------

5 Operation and maintenance

5.1 Setting the wavelength

The Optrode M2 has eight LEDs (LED = light-emitting diode) on its optical circuit board that serve as light sources. Each LED emits light in a different wavelength range. The LEDs are labeled with their primary wavelength on the optical circuit board. It is easy to recognize which LED is active and which wavelength is set by the fact that the label is lit in the corresponding color.



Figure 8 Wavelength display of the Optrode M2

Table 1 Wavelength ranges

LED	Color	Usable wavelength range / nm
470	blue	460–480
502	blue-green	485–520
520	green	505–535
574	yellow-green	560–585
590	yellow-orange	575–605
610	orange	595–625
640	light red	620–655
660	red	650–670

The magnetic switch in the sensor head is used to switch between LEDs or it is done automatically via **tiamo** (version 2.5 and higher).

Switching between LEDs with magnet

The wavelength is changed manually as follows:

- 1 Hold a magnet, e.g. a stirring bar, against the magnet symbol.




Figure 9 Switching the LED

- 2 Briefly move the magnet up and down again.
Moving the magnet switches to the next LED.
- 3 Move the magnet again until the desired wavelength is set.

Switching between LEDs via tiamo

The Optrode M2 is recognized automatically by the control software **tiamo** (version 2.5 and higher). However, this only works if the Optrode M2 is connected directly to the titrator or control instrument via USB.

The wavelength is changed as follows in **tiamo**:

- 1 Double-click on the command in the method run.
- 2 Open the **General/Hardware** tab.
- 3 Under **Sensor**, select the option **Optrode type 2**.
- 4 Under **Wavelength**, select the desired wavelength.

(A wavelength of 610 nm is suggested by default, no matter what wavelength the Optrode M2 is currently set to.)

The selected wavelength can later be used in the formula editor as variable (**.WVL**).

**NOTE**

- The LED only switches to the selected wavelength once the determination starts.
- The wavelength cannot be changed during a running determination (neither with a magnet nor via **tiamo**).

**CAUTION**

The light intensity of the LED must have stabilized sufficiently before use. Each time you switch on and each time you change the wavelength, wait at least ten minutes before starting the first determination.

**NOTE**

The wavelength remains set even if the Optrode M2 is switched off or its energy supply is interrupted (USB connection).


5.2 Sensor care

The Optrode M2 does not require any special care. Store the sensor dry in the storage vessel provided.

If it becomes dirty, the Optrode M2's glass shaft can be cleaned with warm water and a small amount of detergent or ethanol. Use a soft cloth for this.

Problem	Cause	Remedy
	<i>The Optrode M2 is defective.</i>	Replace the Optrode M2.
No signal. The measured value is constant at approximately 0 mV.	<i>The solution is too dark.</i>	<ul style="list-style-type: none"> ▪ Use less indicator. ▪ Dilute the solution.
	<i>The Optrode M2 is dirty.</i>	Clean the sensor, particularly the smooth surfaces in the measuring gap (<i>see chapter 5.2, page 11</i>).
Signal is too high. Constant measured value at approximately 900 mV.	<i>The solution is too bright.</i>	Use more indicator.
	<i>Too much ambient light present (e.g. direct sunlight).</i>	Protect the measurement setup against light.
Optrode M2 is not recognized automatically in <i>tiamo</i> (version 2.5 and higher).	<i>The Optrode M2 firmware is outdated.</i>	Have the firmware updated by a regional Metrohm service representative.

7 Recycling and disposal



Properly dispose of chemicals and of the product to reduce negative effects on the environment and public health. Local authorities, waste disposal companies or dealers provide more detailed information on disposal. Observe the WEEE EU directive (WEEE = Waste Electrical and Electronic Equipment) for the proper disposal of waste electronic equipment within the European Union.

8 Technical specifications

8.1 Measuring mode

Measured quantity Absorption
 (For titration only)

8.2 Measuring interval

Measuring cycle 3 ms

8.3 Resolution

Voltage 0.1 mV

8.4 Ambient conditions

Nominal function range 0–40 °C
 at max. 80% relative humidity, non-condensing

Altitude max. 3,000 m. above sea level

pH range 0–14

8.5 Energy supply

Voltage 4.75–5.25 VDC

Power consumption max. 85 mA



8.6 Reference conditions

Ambient temperature +25 °C (± 3 °C)

Device status At operating temperature, LED used operating for at least 10 minutes

8.7 Light sources (LEDs)

Typical wavelengths

- 470 nm
- 502 nm
- 520 nm
- 574 nm
- 590 nm
- 610 nm
- 640 nm
- 660 nm

8.8 Dimensions

Shaft diameter 12 mm

Sleeve diameter 14.2 mm

Height 177 mm

Installation length 135 mm

Cable length 1.2 m fixed cable (plug F and USB)

8.9 Material

Glass tube Borosilicate glass

Sleeve PBT

Index

A

Antidiffusion tip 7

C

Care 11

Cleaning 11

Connect 4

D

Dimensions 16

E

Energy supply 4, 6, 11, 15

Error search 12

I

Install 7

L

LED 9

Light intensity 11

Light reflections 7

Light sources 9

M

Magnet 10

Magnet symbol 10

Magnetic switch 10

P

Photodiode 4

Power outlet 6

Power supply unit 6

Propeller stirrer 8

R

Resolution 15

S

Sensor cable 5

Sensor F plug 4

Stirring direction 7

Stirring rate 7

Storage 11

Sunlight 7

Switch LED 10

Symbol 10

T

tiamo

 Switch LED 10

Troubleshooting 12

Turbulence 7

U

USB

 Adapter 6

 Cable 5

 Plug 4

 Power supply unit 6

W

wavelength 10

Wavelength 4, 9, 10, 16

Wavelength ranges 9