

# IC Equipment Set



Flow Cell for Bioscan – 6.5331.0X0

Manual

8.110.8003EN / 2019-09-25





Metrohm AG  
CH-9100 Herisau  
Switzerland  
Phone +41 71 353 85 85  
Fax +41 71 353 89 01  
info@metrohm.com  
www.metrohm.com

# **IC Equipment Set**

## **Flow Cell for Bioscan – 6.5331.0X0**

### **Manual**

Technical Communication  
Metrohm AG  
CH-9100 Herisau  
techcom@metrohm.com

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

# Table of contents

<b>1</b>	<b>Introduction</b>	<b>1</b>
1.1	Description .....	1
1.2	About the Documentation .....	1
1.2.1	Symbols and conventions .....	1
<b>2</b>	<b>Installation</b>	<b>3</b>
<b>3</b>	<b>Technical Data</b>	<b>6</b>
3.1	General .....	6
3.2	Gold flow cell .....	6
3.3	Glassy Carbon flow cell .....	7
3.4	Platinum flow cell .....	7
3.5	Silver flow cell .....	8
<b>4</b>	<b>Accessories</b>	<b>9</b>
	<b>Index</b>	<b>10</b>



## Table of figures

Figure 1	Parts and connections of flow cell .....	3
----------	--	---

# 1 Introduction

## 1.1 Description

The flow cell is used for amperometric detection with the Bioscan measuring instruments. The flow cell's construction is based on the Wall-Jet-Principle. According to the Wall-Jet-Principle the inlet is positioned exactly opposite to the working electrode, which allows for optimum measuring performance.

The flow cell is available in four versions:

- Gold flow cell (6.1254.010)
- Glassy carbon flow cell (6.1254.050)
- Platinum flow cell (6.1254.060)
- Silver flow cell (6.1254.070)

For a more detailed description of the specification and intended use of the individual cells, please see (*see chapter 3, page 6*).

## 1.2 About the Documentation




### CAUTION






Please study this documentation carefully before you start to use the instrument. The documentation contains information and warnings that must be observed by the user in order to guarantee the safe use of the instrument.

### 1.2.1 Symbols and conventions

The following symbols and styles are used in this documentation:

(5-12)	<p><b>Cross-reference to figure legend</b></p> <p>The first number refers to the figure number, the second to the instrument part in the figure.</p>
1	<p><b>Instruction step</b></p> <p>Carry out these steps in the sequence shown.</p>
	<p><b>Warning</b></p> <p>This symbol draws attention to a possible life hazard or risk of injury.</p>



	<p><b>Warning</b></p> <p>This symbol draws attention to a possible hazard due to electrical current.</p>
	<p><b>Warning</b></p> <p>This symbol draws attention to a possible hazard due to heat or hot instrument parts.</p>
	<p><b>Warning</b></p> <p>This symbol draws attention to a possible biological hazard.</p>
	<p><b>Caution</b></p> <p>This symbol draws attention to a possible damage of instruments or instrument parts.</p>
	<p><b>Note</b></p> <p>This symbol marks additional information and tips.</p>

## 2 Installation



### CAUTION

Never switch on the flow cell when it is..

- ..not being rinsed through with a conducting eluent at the time; or
- ..not completely connected; or
- ..moist on the outside so that a short circuit could occur between the connections of working electrode and auxiliary electrode.




Figure 1 Parts and connections of flow cell

**1 Reference electrode connection**  
For black electrode cable (6.2156.000).

**2 Inlet for eluent**  
Connection of PEEK capillary by 6.2744.130 KELLF pressure screw only. PEEK pressure screws may damage the threads of the cell.

**3 Auxiliary electrode connection**  
For blue electrode cable (6.2156.000).

**4 Working electrode connection**  
For red electrode cable (6.2156.000).

**5 Outlet for eluent**  
Connection of PEEK capillary by 6.2744.130 KELLF pressure screw only. PEEK pressure screws may damage the threads of the cell.

**6 Mounting screws (4x)**

**7 Distance piece 50  $\mu\text{m}$  (6.1254.020)**  
Defines the volume of the cell.

**8 Check marks**  
For correct alignment of the top and bottom part of the flow-through cell.





## **5 Mount the flow cell**

- Attach the flow cell to the measuring cell holder of the Bioscan.
- Rotate flow cell in its holder so that the outlet is located as high as possible. This allows possibly occurring air bubbles to escape from the cell.



## 3 Technical Data

### 3.1 General

<i>Construction</i>	Flow cell with working, reference and auxiliary electrode.								
<i>Cell volume</i>	The cell volume depends on the thickness of the selected distance piece:								
	<table> <thead> <tr> <th>Distance piece</th> <th>Cell volume</th> </tr> </thead> <tbody> <tr> <td>25 µM</td> <td>0.15 µL</td> </tr> <tr> <td>50 µM</td> <td>0.29 µL</td> </tr> <tr> <td>120 µM</td> <td>0.71 µL</td> </tr> </tbody> </table>	Distance piece	Cell volume	25 µM	0.15 µL	50 µM	0.29 µL	120 µM	0.71 µL
Distance piece	Cell volume								
25 µM	0.15 µL								
50 µM	0.29 µL								
120 µM	0.71 µL								
<i>Working temperature</i>	The flow cell should not be operated permanently at temperatures higher than 45 °C.								
<i>Reference electrode</i>									
<i>Type</i>	Solid phase reference electrode								
<i>Conversion to Ag/AgCl</i>	$E_{\text{Hy-REF}} = E_{\text{Ag/AgCl}} - 328 + 29.9 \text{ pH (mV)}$								

### 3.2 Gold flow cell

<i>Arbeitselektrode</i>	
<i>Material</i>	Gold
<i>Durchmesser</i>	3 mm
<i>Working electrode</i>	
<i>Material</i>	Gold
<i>Diameter</i>	3 mm
<i>Applications</i>	<p>Sugar and amino acids</p> <ul style="list-style-type: none"> <li>▪ mono-, di-, oligo- and polysaccharides</li> <li>▪ sugar alcohols</li> <li>▪ sugar amines</li> <li>▪ sugar acids</li> <li>▪ amino acids</li> <li>▪ antibiotics</li> </ul>
<i>Operating modes</i>	DC and PAD

*Operating range*

<i>Acidic medium</i>	-0.35 V...+1.1 V
<i>Alkaline medium</i>	-1.25 V...+0.75 V

**3.3 Glassy Carbon flow cell***Working electrode*

<i>Material</i>	Glassy-Carbon
<i>Diameter</i>	3 mm

*Applications*

- Aromatics and amines
- catecholamines, aromatic amines
  - inorganic anions (nitrite, sulfite ...)
  - phenols
  - vitamins
  - few amino acids

*Operating modes* DC only (not PAD)

*Operating range*

<i>Acidic medium</i>	-0.8 V... +1.3 V
<i>Alkaline medium</i>	-1.5 V... +0.6 V

**3.4 Platinum flow cell***Working electrode*

<i>Material</i>	Platinum
<i>Diameter</i>	3 mm

*Applications*

- Special applications
- alcohols
  - glycols
  - hydrogen peroxide
  - hydrazine
  - arsenite, hypochlorite

*Operating modes* DC and PAD

*Operating range*

<i>Acidic medium</i>	-0.2 V... +1.3 V
<i>Alkaline medium</i>	-0.9 V... +0.65 V



## 3.5 Silver flow cell

### *Working electrode*

*Material* Silver

*Diameter* 2 mm

*Applications* Environmental applications

- halides
- cyanide, sulfide
- thiosulfate
- pharmaceuticals

*Operating modes* DC and PAD

### *Operating range*

*Acidic medium* -0.55 V... +0.4 V

*Alkaline medium* -1.2 V... +0.1 V

## 4 Accessories

Up-to-date information on the scope of delivery and optional accessories for your product can be found on the Internet. You can download this information using the article number as follows:

### Downloading the accessories list

- 1** Enter <https://www.metrohm.com/> into your Internet browser.
- 2** Enter the article number (e.g. **6.5331.0X0**) into the search field.  
The search result is displayed.
- 3** Click on the product.  
Detailed information regarding the product is shown on various tabs.
- 4** On the **Included parts** tab, click on **Download the PDF**.  
The PDF file with the accessories data is created.



### NOTICE

Once you have received your new product, we recommend downloading the accessories list from the Internet, printing it out and keeping it together with the manual for reference purposes.



# Index

**A**

auxiliary electrode ..... 3

**C**

Cell volume ..... 6

**D**

Distance piece ..... 3

Distand piece ..... 6

**G**Glassy Carbon flow cell  
Technical data ..... 7

Gold flow cell

Technical data ..... 6

**I**

Installation ..... 3

**P**Platinum flow cell  
Technical data ..... 7**R**

reference electrode ..... 3

**S**Silver flow cell  
Technical data ..... 8**T**

Technical data ..... 6

Temperature ..... 6

**W**

working electrode ..... 3