

# 708 Sampling Unit

Series 01 ...



**Instructions for use 8.708.1003**

91.12 Ti/mm

# Instructions for Use of 708 Sampling Unit

## Table of contents

<b>1. Overview</b> .....	2
<b>2. Working with the 684 KF Coulometer and 698 Autosampler</b> .....	4
2.1 Cable connections .....	4
2.2 Tubing connections when 703 Ti Stand is used .....	4
2.3 Tubing connections when 649 Magnetic Stirrer is used .....	5
2.4 Operation .....	6
2.5 Calibration of the sample loop .....	7
<b>3. Working with the 701 KF Titrino and 698 Autosampler</b> .....	8
3.1 Cable connections .....	8
3.2 Tubing connections .....	8
3.3 Operation .....	9
3.4 Calibration of the sample loop .....	10
<b>4. Working with the 670 Titroprocessor and 677 Drive Unit</b> .....	11
4.1 Cable connections .....	11
4.2 Tubing connections .....	11
4.3 Example of PAGE 3 of the 670 Titroprocessor .....	13
4.4 Operation .....	14
4.5 Calibration of the sample loop .....	14
<b>5. Working with a titrator and the 664 Control Unit</b> .....	15
5.1 Cable connections .....	15
5.2 Tubing connections .....	15
5.3 Program of the 664 Control Unit and setting the waiting and valve opening times ..	17
5.4 Operation .....	18
5.5 Calibration of the sample loop .....	18
<b>6. Troubleshooting, rectification of malfunctions</b> .....	19
<b>7. Appendix</b> .....	20
7.1 Technical specifications .....	20
7.2 Pin assignment .....	21
7.3 Warranty .....	23
7.4 Scope of delivery and ordering designations .....	24
<b>Index</b> .....	27



# 1. Overview

## Front

1


### Valve block

For the following tubing connections:


"waste"	Waste tubing.
"smpl"	Sample feed.
"loop"	Sample loop, metering system.
"transf.sol."	Transfer solution. Transfers the sample from the sample loop to the titration vessel.
"to cell"	Tubing leading to titration vessel.

### Flow path

#### "fill" position



#### "inject" position



2

### Tubing cassettes

in working position. If the pump is switched off for a longer duration, the tubing cassettes should be raised. This measure prolongs lifetime of the pump tubings.


2a) Setting the contact pressure, i.e. control of the flow rate.

Lever up: highest pressure = highest flow rate.

Lever down: lowest contact pressure = lowest flow rate.

**Important:** An excessive contact pressure shortens the lifetime of the pump tubing. This must be especially noted in KF titrations.

Inserting the pump tubing in the tubing cassette:



3

### On/off switching of pump

Manual operation of the pump. Must be "on" when the pump is switched on and off by remote control.

4

### Switching the valve

Manual switching of the valve, for the flow path see ①.

Warning: Do not press key if a remote control cable is attached to the socket "Valve"! (The valve switching could be damaged.)

5

**On/off remote control of pump**

Must be "on" if the pump is switched on and off by remote control.

**Rear**

6

**Socket for on/off switching of pump**

Switching via remote control. Switches "Remote" and "Pump" must be "on".

7

**Socket for valve switching**

8

**Socket for connection to Titrino**

Ensures that a start is made only when the titration vessel is conditioned.

9

**Rating plate**

with manufacturing, serial and instrument number.

10

**Display of the set mains voltage**

Before switching the apparatus on for the first time, ensure that the set mains voltage matches the voltage of your supply system. If not, disconnect mains cable and alter the voltage setting.

11

**Mains connection**

With supply systems in which severe HF disturbances are superimposed on the mains voltage, the Sampling Unit should be operated via an additional mains filter, e.g. Metrohm model 615.

When the apparatus is opened or parts removed from it, certain components may be live if it is connected to the supply system. It is thus essential to disconnect the mains cable before altering certain settings or replacing parts.


12

**Power switch**

## 2. Working with the 684 Coulometer and 698 Autosampler

### 2.1 Cable connections


The instruments are interconnected as follows:



Configure the 3.540.2441 RS232 Interface of the 684 KF Coulometer so that calculations are performed without confirmation of the sample weight: DIP switch 8 of interface to ON.

### 2.2 Tubing connections when 703 Ti Stand is used

The following setup is recommended (see cover picture):



#### Tubing connections

- Connection 698 Autosampler - "smpl" of the valve block: This connection should be as short as possible (loss of sample!). With the recommended setup, a length of ca. 40 cm is sufficient. Shorten 6.1822.010 tubing to ca. 40 cm. Attach the 6.2744.000 Compression Fitting from the accessories of the 698 Autosampler to the free end of the tubing. Connect tubing between 698 Autosampler and "smpl".
- Attach 6.1822.010 waste tubing to "waste" of the valve block and route to waste bottle.


- Fit 6.1808.020 Tubing Connectors to both ends of 6.1826.000 pump tubing. Hold the tubing with a cloth to facilitate insertion of the tubing connectors. Keep a pump tubing on hand as a replacement. You should always have spare pump tubing in reserve.
- Remove both tubing cassettes from the holder of the 708 Sampling Unit. You need only one tubing cassette.  
Lay pump tubing in the guide of the cassette so that the stopper is on the left (see page 2).  
Replace cassette in the holder ensuring the pump tubing is free from kinks.
- Connect right end of the pump tubing using 6.1805.300 tubing to "trans.sol" of the valve block.
- Connect left end of pump tubing to the 6.1824.010 Twin Burette Tip.
- Connect twin burette tip to "to cell" of the valve block. Shorten this end of the burette tip somewhat so that the solvent flow can be observed.  
Fix the twin burette tip in the titration vessel upper half with the 6.1446.060 Stopper and 6.2730.030 Screw Fitting.
- Attach the desired 6.1825.XXX Sample Loop to the two "loop" connections of the valve block. For work with the 698 Autosampler, sample loops with a volume  $\leq 250 \mu\text{L}$  can be used.
- Attach 6.2021.020 Electrode Holder to the stand rod of the 703 Ti Stand and use 6.1228.000 Holder for storage of the aspiration tube. During aspiration, the twin burette tip can be stored in the 6.1228.000 Holder.

### Function check

Switch the pump on. The flow rate should be regulated by means of the contact pressure (see page 2) to ca. 5 - 7 mL/min ("fast drop rate"). In case of problems, see page 19.

## 2.3 Tubing connections when 649 Magnetic Stirrer is used

The following setup is recommended:



### Tubing connections

- Connection 698 Autosampler – "smpl" of the valve block: This connection should be as short as possible (loss of sample!). With the recommended setup, a length of ca. 40 cm is sufficient. Shorten 6.1822.010 tubing to ca. 40 cm. Attach the 6.2744.000 Compression Fitting from the accessories of the 698 Autosampler to the free end of the tubing.  
Connect tubing between 698 Autosampler and "smpl".
- Attach 6.1822.010 waste tubing to "waste" of the valve block and route to waste bottle.
- Fit 6.1808.020 Tubing Connectors to both ends of 6.1826.000 pump tubing. Hold the tubing with a cloth to facilitate insertion of the tubing connectors. Keep a pump tubing on hand as a replacement.
- Remove both tubing cassettes from the holder of the 708 Sampling Unit. You need only one tubing cassette.  
Lay pump tubing in the guide of the cassette so that the stopper is on the left (see page 2).  
Replace cassette in the holder ensuring the pump tubing is free from kinks.

- Connect right end of the pump tubing using 6.1805.300 tubing to "trans.sol" of the valve block.
- Connect left end of pump tubing with the aid of the 6.1805.310 tubing to the longer piece of the 6.1439.010 Dispensing/aspiration Tube of the coulometry cell.
- Take the second 6.1805.310 tubing, pull a length of ca. 11 cm out of the nipple at one end: Cut off brown seal, then cut tubing diagonally to a length of ca. 1 cm. Slide the seal over the tubing with the white part facing the tubing end. Hold tubing with flat-nose pliers and move seal ca. 11 cm up the tubing. Finally, cut tubing off level using a razor blade.  
Connect this end to the shorter part of the dispensing/aspiration tube. The tubing runs through the tube. Attach other end of 6.1805.310 tubing to "to cell" of the valve block.
- Attach the desired 6.1825.XXX Sample Loop to the two "loop" connections of the valve block. For work with the 698 Autosampler, sample loops with a volume  $\leq 250 \mu\text{L}$  can be used.

### Function check

Switch the pump on. The flow rate should be regulated by means of the contact pressure (see page 2) to ca. 5-7 mL/min ("fast drop rate"). In case of problems, see page 19.

## 2.4 Operation

1. Prepare 684 KF Coulometer: Add reagent to coulometry cell, condition cell (press <GO>). Switch on send and enter sample weight (<sample>) in mg, see also page 7.
2. Switch on the pump of the 708 Sampling Unit. Check flow rate ("fast drop rate"). If need be, move pump tubing a little to the right. In case of problems, see page 19.
3. Fill the sample vials of the 698 Autosampler with ca. 0.9 mL sample and press in the cover fully with the opening downwards. Load sample turntable of 698 Autosampler with vials. The sample vials should be dried for 1 h at 105°C and stored in a desiccator.  
**Warning:** The samples must not contain any solids! (danger of blockage)
4. "Program" 698 Autosampler:
 

<b>&lt;STOP&gt;</b>	Press stop key
VIAL	Display
XX <ENTER>	Enter first sample vial
THRU	Display
XX <ENTER>	Enter last sample vial
5. Starting 698 Autosampler:
 

<b>&lt;START&gt;</b>	Press start key
INIT V	Display
XX <ENTER>	Enter position of first sample vial
RINSE	Rinse cycle between samples (0 = no)
0 <ENTER>	Enter 0
LAST V	Position of the last sample vial
XX <ENTER>	This starts the first determination

A start is made as soon as the drift measured by the 684 KF Coulometer is less than 100  $\mu\text{g}/\text{min}$ . If you wish to start at lower drift values, you will have to wait until these are reached. A good drift value when the 708 Sampling Unit is running is ca. 12  $\mu\text{g}/\text{min}$ .

### Switching off

If the 708 Sampling Unit is switched off, lift up tubing cassette. This measure prolongs the life of the pump tubing!

If the 708 Sampling Unit is not used for some time, the wet system should be rinsed with methanol: Before switching off, siphon in methanol "as a sample" during ca. 3 min.

## 2.5 Calibration of the sample loop

The exact volume of the sample loop is determined with, e.g. a solution of water in methanol (standard methanol).

- Using a microlitre syringe or a Dosimat, inject a known volume  $V_1$  of your solution into the coulometry cell. Perform, e.g. 5 determinations and calculate the mean value of the titrated water  $\rightarrow t_1$ .
- Now add the same solution to the sample vials of the 698 Autosampler and start so that the solution is injected via the sample loop of the 708 Sampling Unit. Again perform 5 determinations and determine the mean value of the titrated water  $\rightarrow t_2$ . Same conditions and same titration parameters as above!
- Calculate the volume of the sample loop  $V_{loop}$  using the following formula:

$$V_{loop} = \frac{V_1 * t_2}{t_1}$$

- Write the volume on the rough surface of the small plate of the sample loop.
- Conversion of the volume into the sample weight:  
Weigh an exactly known volume  $V_0$  of your sample  $\rightarrow m_0$ . The weight is calculated using


$$m_{sample} = \frac{m_0 * V_{loop}}{V_0}$$

### 3. Working with the 701 Titrino and 698 Autosampler


It is also possible to use the 702 SM Titrino in place of the 701 KF Titrino for work in the SET mode with conditioning.

#### 3.1 Cable connections

The instruments are interconnected as follows:



#### 3.2 Tubing connections



#### Tubing connections

- Connection 698 Autosampler - "smpl" of the valve block: This connection should be as short as possible (loss of sample!). With the recommended setup, a length of ca. 40 cm is sufficient. Shorten 6.1822.010 tubing to ca. 40 cm. Attach the 6.2744.000 Compression Fitting from the accessories of the 698 Autosampler to the free end of the tubing. Connect tubing between 698 Autosampler and "smpl".
- Attach 6.1822.010 waste tubing to "waste" of the valve block and route to waste bottle.
- Fit 6.1808.020 Tubing Connectors to both ends of 6.1826.000 pump tubing. Hold the tubing with a cloth to facilitate insertion of the tubing connectors. Keep a pump tubing on hand as a replacement. You should always have spare pump tubing in reserve.

- Remove both tubing cassettes from the holder of the 708 Sampling Unit. You need only one cassette.  
Lay pump tubing in the guide of the cassette so that the stopper is on the left (see page 2).  
Replace cassette in the holder ensuring the pump tubing is free from kinks.
- Connect right end of the pump tubing using 6.1805.300 tubing to "trans.sol" of the valve block.
- Connect left end of pump tubing to the 6.1824.010 Twin Burette Tip.
- Connect twin burette tip to "to cell" of the valve block. Shorten this end of the burette tip somewhat so that the solvent flow can be observed.  
Fix the twin burette tip in the titration vessel upper half with the 6.2730.050 Stopper.
- Attach the desired 6.1825.XXX Sample Loop to the two "loop" connections of the valve block.  
For work with the 698 Autosampler, sample loops with a volume  $\leq 250 \mu\text{L}$  can be used.

### Function check

Switch the pump on. The flow rate should be regulated by means of the contact pressure (see page 2) to ca. 5 - 7 mL/min ("fast drop rate"). In case of problems, see page 19.

## 3.3 Operation

1. Prepare 701 KF Titrino: Add solvent and condition titration cell (press <START>). Select the desired report and enter sample size.  
After the start, the titration method of the 701 KF Titrino should not contain any inquiries for identification or sample size.
2. Switch on the pump of the 708 Sampling Unit. Check flow rate ("fast drop rate"). If need be, move pump tubing a little to the right. In case of problems, see page 19.
3. Fill the sample vials of the 698 Autosampler with ca. 0.9 mL sample and press in the cover fully with the opening downwards. Load sample turntable of 698 Autosampler with glass vials. The sample vials should be dried for 1 h at 105°C and stored in a desiccator.  
**Warning:** The vials must not contain any solids! (danger of blockage)
4. "Program" 698 Autosampler:
 

<b>&lt;STOP&gt;</b>	Press stop key
VIAL	Display
XX <ENTER>	Enter first sample vial
THRU	Display
XX <ENTER>	Enter last sample vial
5. Starting 698 Autosampler:
 

<b>&lt;START&gt;</b>	Press start key
INIT V	Display
XX <ENTER>	Enter position of first sample vial
RINSE	Rinse cycle between samples (0 = no)
0 <ENTER>	Enter 0
LAST V	Position of the last sample vial
XX <ENTER>	This starts the first determination

A start is made as soon as the 701 KF Titrino is conditioned.

**Switching off**

If the 708 Sampling Unit is switched off, lift up tubing cassette. This measure prolongs the life of the pump tubing!

If the 708 Sampling Unit is not used for some time, the wet system should be rinsed with methanol: Before switching off, siphon in methanol "as a sample" during ca. 3 min.

**3.4 Calibration of the sample loop**

The exact volume of the sample loop is determined with, e.g. a solution of water in methanol (standard methanol).

- Using a microlitre syringe or a Dosimat, inject a known volume  $V_1$  of your solution into the titration cell. Perform, e.g. 5 determinations and calculate the mean value of the titrant consumption  $\rightarrow t_1$ .
- Now add the same solution to the sample vials of the 698 Autosampler and start so that the solution is now injected via the sample loop of the 708 Sampling Unit. Again perform 5 determinations and determine the mean value of the titrant consumption  $\rightarrow t_2$ . Same conditions and same titration parameters as above!
- Calculate the volume of the sample loop  $V_{loop}$  using the following formula:

$$V_{loop} = \frac{V_1 * t_2}{t_1}$$


- Write the volume on the rough surface of the small plate of the sample loop.
- Conversion of the volume into the sample weight:  
Weigh an exactly known volume  $V_0$  of your sample  $\rightarrow m_0$ . The weight is calculated using

$$m_{sample} = \frac{m_0 * V_{loop}}{V_0}$$

## 4. Working with the 670 Titroprocessor and 677 Drive Unit


### 4.1 Cable connections

The instruments are interconnected as follows:




### 4.2 Tubing connections

- Fit 6.1808.020 Tubing Connectors to both 6.1826.000 pump tubings. Hold the tubing with a cloth to facilitate insertion of the tubing connectors. Keep a set of pump tubings on hand as a replacement.
- Remove both tubing cassettes from the holder of the 708 Sampling Unit.
- Setting up 1<sup>st</sup> pump circuit:



- . Lay pump tubing in the guide of a cassette so that the stopper is on the left (see page 2). Replace cassette in the holder ensuring the pump tubing is free from kinks.
- . Connect right end of the pump tubing using 6.1805.300 tubing to "trans.sol" of the valve block.
- . The left end of the pump tubing leads to the transfer solution.
- . Cut 6.1803.020 tubing using a sharp cutting tool (e.g. razor blade) to the correct length. Attach 6.1811.020 Nipple to other end of tubing: Slide nipple with thread facing end of tubing over the tubing. Cut tubing diagonally to a length of ca. 1 cm. Slide seal with the white part facing the end of the tubing over the tubing. Hold tubing with flat-nose pliers and move seal. Finally, cut off tubing level with seal using a razor blade.
- . Connect nipple to left end of pump tubing and route to transfer solution.
- . Fit 6.1811.020 Nipples to both ends of a ca. 60 cm length of 6.1803.020 tubing. Make connection between "to cell" of valve block and burette tip in the titration vessel and position burette tip in the titration vessel using 6.1446.030 Stopper.
- Setting up 2<sup>nd</sup> pump circuit:



- . Insert second pump tubing in the tubing cassette and place the cassette in the holder of the 708 Sampling Unit.
- . Use 6.1805.300 tubing to connect the right end of the pump tubing to "smpl" of the valve block.
- . The left end of the pump tubing leads to the sample.
- . Attach a 6.1811.020 Nipple to another piece of 6.1803.020 tubing. Connect this to the left end of the pump tubing and route to sample.
- . If you supply the samples using a sample changer, also fit a 6.1811.020 Nipple to the other end of the 6.1803.020 tubing. Attach a 6.1543.060 Burette Tip to this and insert it in the sample head of the sample changer.
- . Use another piece of 6.1803.020 tubing for the waste tubing and fit a 6.1811.020 Nipple to one end. Connect this tubing to "waste" of the valve block and lead it off to the waste bottle.
- Attach the desired 6.1825.XXX Sample Loop to both "loop" connections of the valve block.

**Function check**

Switch on all instruments.

Set the two switches "Remote" and "Pump" on the 708 Sampling Unit to "on".

Set the "Remote" switch on the two 683 Pump Units to "off".

Switch the 670 Titroprocessor to "manual control" (MC). Test the control lines:

- 1: Solvent addition
- 2: Aspiration
- 4: Switch on pump of 708 Sampling Unit
- 5: Switch valve of 708 Sampling Unit to "inject"
- 6: Switch valve of 708 Sampling Unit to "fill"

Activate stirrer.

The flow rate should be regulated by means of the contact pressure (see page 2) to ca. 5 - 7 mL/min (ca. 2 drops per second). In case of problems, see page 19.

**4.3 Example of PAGE 3 of the 670 Titroprocessor**

The following method includes commands for control of the 708 Sampling Unit and a sample changer. The latter are marked with "SC". If the sample is supplied without a sample changer, simply omit these commands.

1	OMOVE1;		SC	18	MEAS 1;	
2	CTR 12;	5.0 s	Aspirate	18a		
3	OCTR12;	100 ms		18b		
4	WPOS 1;		SC	18c		
5	CTR 16;	100 ms	708 "Fill"	19	DYNT 1;	Titration
6	OCTR16;	100 ms		19a		
7	CTR 14;	100 ms	708 Pump on	19b		
8	NOP ;	XX s 1)	Fill sample loop	19b		
9	CTR 11;	YY s 2)	Add solvent	20	TSTOP ;	
10	OCTR11;	100 ms		20a		
11	CTR 15;	100 ms	708 "inject"	20b		
12	OCTR15;	100 ms		20c		
13	STIR 1;		Stirrer on	21	OCTIR1;	Stirrer off
14	RPOS 1;		SC	22	CTR 12;	Aspirate
15	SHIFT1;		SC	23	OCTR12;CTR 11;	Rinse
16	NOP ;	XX s 1)	Rinse sample loop	24	OCTR11;CTR 12;	Aspirate
17	OCTR14;	100 ms	708 Pump off	25	OCTR12;CTR 11;	Rinse
26	OCTR11;CTR 12;	5.0 s	Aspirate			
27	OCTR12;CTR 11;	3.0 s	Rinse			
28	OCTR11;					
29	END ;					

**Remarks**

- 1): The waiting times depend on the size of the sample loop, the flow rate and the viscosity of the solutions. With a flow rate of ca. 5 - 7 mL/min and dilute aqueous solutions, the following waiting times result:

Sample loop	Waiting time
6.1825.060 5 ml	120 s
6.1825.050 1 ml	35 s
6.1825.040 0.5 ml	20 s

- 2): When the sample loop is being filled, transfer solution (solvent) is conveyed to the titration vessel. The waiting time for the addition of more solvent depends on this volume and the size of the titration vessel.

## 4.4 Operation

Switch the two switches "Pump" and "Remote" on the 708 Sampling Unit to "on".

Adjust the flow rate by means of the contact pressure (see page 2) to 5 – 7 mL/min (ca. 2 drops per second). In case of problems, see page 19.

Start method on 670 Titroprocessor.

**Warning:** The samples must not contain any solids! (danger of blockage)

### Switching off

If the 708 Sampling Unit is switched off, lift up tubing cassette. This measure prolongs the life of the pump tubing!

If the 708 Sampling Unit is not used for some time, the wet system should be rinsed with solvent: Before switching off, siphon in solvent "as a sample" for ca. 3 min.

## 4.5 Calibration of the sample loop

The exact volume of the sample loop is determined with a titratable solution (e.g. NaOH).

- Pipette a known volume  $V_1$  of your solution into the titration vessel. Perform, e.g. 5 determinations and calculate the mean value of the titrant consumption  $\rightarrow t_1$ .
- Now add the same solution via the sample loop of the 708 Sampling Unit. Again perform 5 determinations and determine the mean value of the titrant consumption  $\rightarrow t_2$ . Same conditions and same titration parameters as above!
- Calculate the volume of the sample loop  $V_{loop}$  using the following formula:

$$V_{loop} = \frac{V_1 * t_2}{t_1}$$


- Write the volume on the rough surface of the small plate of the sample loop.
- Conversion of the volume into the sample weight:  
Weigh an exactly known volume  $V_0$  of your sample  $\rightarrow m_0$ . The weight is calculated using

$$m_{sample} = \frac{m_0 * V_{loop}}{V_0}$$

## 5. Working with a titrator and the 664 Control Unit


### 5.1 Cable connections

The instruments are interconnected as follows:




### 5.2 Tubing connections

- Fit 6.1808.020 Tubing Connectors to both 6.1826.000 pump tubings. Hold the tubing with a cloth to facilitate insertion of the tubing connectors. Keep a set of pump tubings on hand as a replacement.
- Remove both tubing cassettes from the holder of the 708 Sampling Unit.
- Setting up 1<sup>st</sup> pump circuit:



- Lay a pump tubing in the guide of a cassette so that the stopper is on the left (see page 2). Replace cassette in the holder ensuring the pump tubing is free from kinks.
- Connect right end of the pump tubing using 6.1805.300 tubing to "trans.sol" of the valve block.

- The left end of the pump tubing leads to the transfer solution.  
Cut 6.1803.020 tubing using a sharp cutting tool (e.g. razor blade) to the correct length. Attach 6.1811.020 Nipple to one end of tubing: Slide nipple with thread facing end of tubing over the tubing. Cut tubing diagonally to a length of ca. 1 cm. Slide seal with the white part facing the end of the tubing over the tubing. Hold tubing with flat-nose pliers and move seal. Finally, cut off tubing level with seal using a razor blade.  
Connect nipple to left end of pump tubing and route to transfer solution.
- Fit 6.1811.020 Nipples to both ends of a ca. 60 cm length of 6.1803.020 tubing. Make connection between "to cell" of valve block and burette tip in the titration vessel and position burette tip in the titration vessel using 6.1446.030 Stopper.
- Setting up 2<sup>nd</sup> pump circuit:



- Insert second pump tubing in the tubing cassette and place the cassette in the holder of the 708 Sampling Unit.
- Use 6.1805.300 tubing to connect the right end of the pump tubing to "smpl" of the valve block.
- The left end of the pump tubing leads to the sample.  
Attach a 6.1811.020 Nipple to another piece of 6.1803.020 tubing. Connect this to the left end of the pump tubing and route to sample.  
If you supply samples using a sample changer, also fit a 6.1811.020 Nipple to the other end of the 6.1803.020 tubing. Attach a 6.1543.060 Burette Tip to this and insert it in the sample head of the sample changer.
- Use another piece of 6.1803.020 tubing fitted at one end with a 6.1811.020 Nipple for the waste tubing. Connect this tubing to "waste" of the valve block and lead off to the waste bottle.
- Attach the desired 6.1825.XXX Sample Loop to the two "loop" connections of the valve block.

### Function check

Set the "Remote" switch on the two 683 Pump Units to "on".

Switch on all instruments.

Set the two switches "Remote" and "Pump" on the 708 Sampling Unit to "on".

Manually activate the different switches on the 664 Control Unit:

V1: 683 Rinsing Pump is activated.

V2: Pump of 708 Sampling Unit is switched on.

V3: 683 Aspiration Pump is activated.

Switch on stirrer.



- 2): When the sample loop is being filled, transfer solution (solvent) is added to the titration vessel. The waiting time  $\Delta t_3$  for the addition of more solvent depends on this volume and the size of the titration vessel.

## 5.4 Operation

Switch the two switches "Pump" and "Remote" on the 708 Sampling Unit to "on".

Adjust the flow rate by means of the contact pressure (see page 2) to 5 - 7 mL/min (ca. 2 drops per second). In case of problems, see page 19.

Start method on 664 Control Unit.

**Warning:** The samples must not contain any solids! (danger of blockage)

### Switching off

If the 708 Sampling Unit is switched off, lift up tubing cassette. This measure prolongs the life of the pump tubing!

If the 708 Sampling Unit is not used for some time, the wet system should be rinsed with solvent: Before switching off, siphon in solvent "as a sample" for ca. 3 min.

## 5.5 Calibration of the sample loop

The exact volume of the sample loop is determined with a titratable solution (e.g. NaOH).


- Pipette a known volume  $V_1$  of your solution into the titration vessel. Perform, e.g. 5 determinations and calculate the mean value of the titrant consumption  $\rightarrow t_1$ .
- Now add the same solution via the sample loop of the 708 Sampling Unit. Again perform 5 determinations and determine the mean value of the titrant consumption  $\rightarrow t_2$ . Same conditions and same titration parameters as above!
- Calculate the volume of the sample loop  $V_{loop}$  using the following formula:

$$V_{loop} = \frac{V_1 \cdot t_2}{t_1}$$

- Write the volume on the rough surface of the small plate of the sample loop.
- Conversion of the volume into the sample weight:  
Weigh an exactly known volume  $V_0$  of your sample  $\rightarrow m_0$ . The weight is calculated using


$$m_{sample} = \frac{m_0 \cdot V_{loop}}{V_0}$$

## 6. Troubleshooting, rectification of malfunctions




### Insufficient or non-existent flow rate

- Increase contact pressure of pump tubing, see page 2.
- Check all connections for tightness.
- Use new pump tubing.
- Undo all connections of the pump circuit. First undo the first connection in the flow direction after siphoning on of the liquid. If liquid still flows here, check the next connection, etc. If liquid transport stops at the valve block:
  - . Valve block obstructed: Unscrew valve block. Clean with solvent, compressed air and a needle, then screw it on again (tighten evenly and without much pressure).
- If liquid transport stops "at intervals":
  - Tubing ends squeezed: Cut off seal.
    - . Cut tubing diagonally to a length of ca. 1 cm. Slide seal with white part facing end of tubing over tubing. Cut off tubing level with seal using a razor blade.




### Leaking of tubing nipples

- Cut off seal. Cut tubing diagonally to a length of ca. 1 cm. Slide seal with white part facing end of the tubing over the tubing. Cut tubing off level with seal using a razor blade.
- Use new 6.1811.020 Nipple.




### Tubing leaking or high drift in KF titrations

- Use new tubing or new 6.1826.000 pump tubing.



### Leaks at valve block

- Tighten 3 Phillips screws of valve block more. Tighten screws evenly.



### Air bubbles in pump circuit

- Immerse tubing which siphons the liquid on completely.
- Check connections for tightness.
- When connected to 698 Autosampler and on use of a "light" solvent (e.g. methanol), use 6.1822.010 waste tubing.


## 7. Appendix

### 7.1 Technical specifications

<b>Pump type</b>	2-channel peristaltic pump
<b>Pump capacity</b> with 6.1826.000 pump tubing, internal diameter 2.54 mm	5 ...7 ml/min (depending on contact pressure)
<b>Pressure</b>	max. 4 bar
<b>Sample requirements</b> Nature Sample loops, volume Viscosity (with standard tubing)	Clear liquids 50 µl...5000 µl max. ca. 200 cp
<b>Materials</b> Tubing connections Pump tubing Valve block Housing	PTFE FPM (Viton) PTFE, PCTFE Light-metal injection moulding, multicoat stove enamel
<b>Remote control</b> via control lines	Pump: on/off, Valve: fill/inject
<b>Safety requirements</b>	Designed and tested in accordance to IEC publication 348, safety class I. For indoor use only. These instructions for use contain information and warnings which must be observed by the user to ensure safe operation and to retain the apparatus in safe condition.
<b>Ambient temperature</b> Nominal operating range Storage transport	5 ... 40 °C – 20 ... 70 °C
<b>Mains connection</b> Mains voltage Mains frequency Greatest power consumption Fusing	100...125 V and 200...250 V, adjustable 50...60 Hz ca. 50 VA (when valve switched) Slow-blow fusible cutout (1 A), 5 x 20 mm
<b>Dimensions</b> Width Height Depth	250 mm 127 mm ca. 350 mm
<b>Weight</b>	ca. 5.7 kg

## 7.2 Pin assignment

	external	Function
<p><b>"Pump" socket</b></p>		<p>ON : 24 V AC OFF: 0 V AC</p> <p>ON: +5 V DC OFF: 0 V DC</p> <p>} not used</p>
<p><b>"Valve" socket Inputs</b></p> <p><b>Outputs</b></p>		<p>Inject = low (active). Additional condition for "Inject", see socket "Titrino", page 22.</p> <p>Fill = low (active)</p> <p>} External start: Pulses at a frequency of 100... 200 Hz, as long as line 2 of the "Titrino" socket is active (if the socket is activated at all) and sets the 708 to "Inject".</p>

	external	Function
<p><b>"Titrimo" socket Inputs</b></p>  <p><b>Outputs</b></p> <p>LSTTL</p>		<p>Socket is activated (+ 5 V)</p> <p>Gives external start (on "Valve" socket, pins 4 and 5).</p> <p>If this line is "low" (active) and Inject (on "Valve" socket, pin 1) is "low" too,</p> <p>then the 708 is set to "Inject" and a start pulse is outputted (start pulse of ca. 150 ms duration).</p>
<p>No liability whatsoever will be accepted for damage caused by improper interconnection of instruments.</p>		

## 7.3 Warranty

The warranty regarding our products is limited to rectification free of charge in our workshops of defects that can be proved to be due to material, design or manufacturing faults which appear within 12 months from the day of delivery. Transport costs are chargeable to the orderer.

For day and night operation, the warranty is valid for 6 months.

Glass breakage in the case of electrodes or other glass parts is not covered by the warranty. Checks which are not a result of material or manufacturing faults are also charged during the warranty period. For parts of outside manufacture insofar as these constitute an appreciable part of our instrument, the warranty stipulations of the manufacturer in question apply.

With regard to the guarantee of accuracy, the technical specifications in the Instructions for Use are authoritative.

Concerning defects in material, construction or design as well as the absence of guaranteed features, the orderer has no rights or claims except those mentioned above.

If damage of the packaging is evident on receipt of a consignment or if the goods show signs of transport damage after unpacking, the carrier must be informed immediately and a written damage report demanded. Lack of an official damage report releases METROHM from any liability to pay compensation.

If any instruments and parts have to be returned, the original packaging should be used if at all possible. This applies above all to instruments, electrodes, burette cylinders and PTFE pistons. Before embedment in wood shavings or similar material, the parts must be packed in a dustproof package (for instruments, use of a plastic bag is imperative). If open assemblies are enclosed in the scope of delivery that are sensitive to electromagnetic voltages (e.g. data interfaces etc.) these must be returned in the associated original protective packaging (e.g. conductive protective bag). (Exception: assemblies with built-in voltage source belong in a non-conductive protective packaging). For damage which arises as a result of non-compliance with these instructions, no warranty responsibility whatsoever will be accepted by METROHM.

## 7.4 Scope of delivery and ordering designations

<b>708 Sampling Unit</b>	<b>2.708.0010</b>
including the following accessories:	
1 PTFE tubing with M6 thread, length 15 cm	6.1805.300
4 Tubing connectors	6.1808.020
1 Nipple with M6 thread	6.1811.020
2 Capillary tubings with M6 thread, length 1 m	6.1822.010
1 Sample loop, volume ca. 100 µl	6.1825.020
1 Sample loop, volume ca. 250 µl	6.1825.030
2 Pump tubings, FPM (Viton)	6.1826.000
1 Mains cable with cable socket, type CEE(22),V	
Cable plug to customer's specifications:	
type SEV 12 (Switzerland...)	6.2122.020
type CEE(7),VII (Germany...)	6.2122.040
type NEMA/ASA (USA...)	6.2122.070
1 Instructions for Use for 708 Sampling Unit	8.708.1003

For instrument combinations ready for use the following additional instruments and accessories are needed:

### 684 KF Coulometer, 698 Autosampler and 703 Ti Stand

1 KF Coulometer	2.684.0140
with printer and corresponding cable	
1 Autosampler	2.698.0010
1 Ti Stand	2.703.0010
(Thread adaptor for Riedel de-Haën and Fluka reagents)	6.1618.000)
1 Cable 698 Autosampler – 708 Sampling Unit	6.2128.040
1 Cable 698 Autosampler – 684 KF-Coulometer	6.2128.050
1 Twin burette tip	6.1824.010
1 Stopper for twin burette tip	6.1446.060
1 Stopper with nipple and O ring	6.2730.030
1 Electrode holder	6.2021.020
1 Holding tube	6.1228.000

### 684 KF Coulometer, 698 Autosampler und 649 Magnetic Swing-out Stirrer

1 KF Coulometer	2.684.0140
with printer and corresponding cable	
1 Autosampler	2.698.0010
1 Magnetic Swing-out Stirrer	2.649.0040
1 Cable 698 Autosampler – 708 Sampling Unit	6.2128.040
1 Cable 698 Autosampler – 684 KF Coulometer	6.2128.050
1 Connector for glass tube with M6 thread	6.1811.000
2 PTFE tubings with M6 thread, length 40 cm	6.1805.310

**701 KF Titrino and 698 Autosampler**

1 Autosampler	2.698.0010
1 KF Titrino	2.701.0010
with printer and corresponding cable and exchange unit	
1 Ti Stand	2.703.0010
1 Cable 698 Autosampler – 708 Sampling Unit	6.2128.060
1 Cable 701 KF Titrino – 708 Sampling Unit	6.2128.070
1 Twin burette tip	6.1824.010
1 Stopper for twin burette tip	6.2730.050

Instead of the 703 Ti Stand, the 649 Magentic Swing-out Stirrer can be used (ordering number 2.649.0040).

**670 Titroprocessor, 677 Drive Unit and 2 x 683 Pump Unit**

Instruments, exchange unit, electrode	
1 Cable 670 Titroprocessor – 677 Drive Unit and 708 Sampling Unit	3.980.3580
1 Titration vessel	6.1415.220
1 Titration vessel lid	6.1414.010
1 Rotor nozzle	6.2740.000
1 Guide sleeve for rotor nozzle	6.1823.000
1 SGJ sleeve for guide sleeve	6.1236.040
1 Burette tip (sample transport into the titration vessel)	6.1543.060
1 Aspiration tip	6.1543.120
2 Movable ball stoppers	6.1446.030
1 Aspiration tube	6.1805.200
1 Screw connector for aspiration tube on 683 Pump Unit	6.1820.030
evtl. 1 sample loop, see "General accessories", page 26	6.1825.XXX
1 PTFE tubing for various connections, length 5 m	6.1803.020
4 Nipples with M6 thread	6.1811.020
1 PTFE tubing with M6 thread, length 15 cm	6.1805.300

If the samples are transported with a sample changer, the following additional accessories are needed:

1 Burette tip (for siphoning on the sample)	6.1543.060
1 Nipple with M6 thread	6.1811.020
evtl. PTFE tubing for various connections, length 5 m	6.1803.020

**Titration, 664 Control Unit and 2 x 683 Pump Unit**

Instruments, exchange unit, electrode and cable Titration – 664 Control Unit	
Special EPROM for 664 Control Unit	6.9960.002
1 Cable 664 Control Unit – 708 Sampling Unit (Valve)	3.980.3590
3 Cables 664 Control Unit – Pumps (708 and 683)	6.2121.010
1 Titration vessel	6.1415.220
1 Titration vessel lid	6.1414.010
1 Rotor nozzle	6.2740.000
1 Guide sleeve for rotor nozzle	6.1823.000
1 SGJ sleeve for guide sleeve	6.1236.040
1 Burette tip (sample transport into the titration vessel)	6.1543.060
1 Aspiration tip	6.1543.120
2 Movable ball stoppers	6.1446.030
1 Aspiration tube	6.1805.200
1 Screw connector for aspiration tubing on 683 Pump Unit	6.1820.030
evtl. 1 sample loop, see "General accessories", page 26	6.1825.XXX
1 PTFE tubing for various connections, length 5 m	6.1803.020
4 Nipples with M6 thread	6.1811.020
1 PTFE tubing with M6 thread, length 15 cm	6.1805.300

If the samples are transported with a sample changer, the following additional accessories are needed:

1 Burette tip (for siphoning on the sample)	6.1543.060
1 Nipple with M6 thread	6.1811.020
evtl. PTFE tubing for various connections, length 5 m	6.1803.020

### **General accessories for 708 Sampling Unit**

Sample loop, volume ca. 50 $\mu$ l	6.1825.010
Sample loop, volume ca. 100 $\mu$ l	6.1825.020
Sample loop, volume ca. 250 $\mu$ l	6.1825.030
Sample loop, volume ca. 500 $\mu$ l	6.1825.040
Sample loop, volume ca. 1 ml	6.1825.050
Sample loop, volume ca. 5 ml	6.1825.060
Pump tubing	6.1826.000
Tubing connectors for pump tubing	6.1808.020
PTFE tubing for various connections, 5 m	6.1803.020
Nipple with M6 thread for 6.1803.020 Tubing	6.1811.020
PTFE tubing with M6 thread, length 15 cm	6.1805.300
PTFE tubing with M6 thread, length 40 cm	6.1805.310
Capillary tubing with M6 thread, length 1 m	6.1822.010
Twin burette tip	6.1824.010

# Index

## A

Accessories	24ff
Air bubbles	19
Aspirating	13, 17
Autosampler	4, 8
- operation	6, 9
- cable connections	4, 8
- tubing connections	4, 5, 8

## C

Cables	24ff
Cable connections	
- 698 Autosampler	4, 8
- 678 EP/KF-Processor	15
- 684 KF Coulometer	4
- 701 KF Titrino	8
- 702 SM Titrino	15
- 670 Titroprocessor	11
- 672, 682, 686 Titroprocessor	15
Calibration of sample loop	7
Cassette	2
Connections	
- 698 Autosampler	4, 8
- 664 Control Unit	15ff
- 678 EP/KF Processor	15
- 684 KF Coulometer	4ff
- 701 KF Titrino	8
- 702 SM Titrino	15
- 670 Titroprocessor	11ff
- 672, 682, 686 Titroprocessor	15
Contact pressure	2
Control Unit	15ff
- accessories	25

## E

EP/KF Processor	15ff
Equipping the titration vessel	
- with 684 KF Coulometer	4, 5
- with 701 KF-Titrino	8
- with Titroprocessor or SM Titrino	11, 15

## F

Fill	2
Flow rate	2
Function control	
- with 698 Autosampler	5, 6, 9
- with 664 Control Unit	16
- with 684 KF Coulometer	5, 6
- with 701 KF Titrino	9
- with 670 Titroprocessor	13

## I

Inject	2
Interconnections	
- 698 Autosampler	4, 8
- 664 Control Unit	15ff
- 684 KF Coulometer	4ff
- 701 KF Titrino	8ff
- 670 Titroprocessor	11ff

## K

KF Coulometer	4ff
- accessories	23
- cable connections	4
- calibration of sample loop	7
- function control	5, 6
- tubing connections	4, 5
KF Titrino	8ff
- accessories	24
- cable connections	8
- calibration of sample loop	10
- function control	9
- tubing connections	8

## L

Leaking	19
Loop	26
- calibration	7, 10, 14, 18

## M

Mains connection	3
Manufacturing number	3

## O

Operation	2
- with 698 Autosampler	6, 9
- with 684 Control Unit	18
- with 684 KF Coulometer	6
- with 701 KF Titrino	9
- with 670 Titroprocessor	14
Ordering designations	24ff
Overview	2

## P

Pin assignment	21
Power connection	3
Problems	19
Pump tubing	26

## R

Rectification of malfunctions	19
Remote control	3, 21
Rinsing of the titration vessel	13, 17
Rinsing of the tubings	
- with 698 Autosampler	6, 9
- with 664 Control Unit	18
- with 684 KF Coulometer	6
- with 701 KF Titrino	9
- with 670 Titroprocessor	14

## S

Sample changer	11, 15
Sample loop	26
- calibration	7, 10, 14, 18
Scope of delivery	23
Serial number	3
SM Titrino	15ff
Sockets	3

Switching off	
- with 664 Control Unit	18
- with 684 KF-Coulometer	7
- with 701 KF-Titrino	9
- with 670 Titroprocessor	14
Switching the valve	2, 3

## T

Technical specifications	20
Titroprocessor	
- 670	11ff
- 672, 682, 686	15ff
Transfer solution	2, 11, 15
Transport path	2
Troubleshooting	19
Tubing cassette	2
Tubing connections	
- 698 Autosampler	4, 5, 8
- 678 EP/KF Processor	15
- 684 KF Coulometer	4, 5
- 701 KF Titrino	8
- 702 SM Titrino	15
- 670 Titroprocessor	11
- 672, 682, 686 Titroprocessor	15
Tubings	26

## V

Valves of 664 Control Unit	17
Valve block	2
Volume of the sample loop	7

## W

Waiting times	
- with 664 Control Unit	17
- with 670 Titroprocessor	13
Warranty	23