



CH-9101 Herisau/Switzerland
Tel. +41 71 353 85 85
Fax +41 71 353 89 01
E-Mail sales@metrohm.ch
<http://www.metrohm.ch>

711 Liquino

Remote Control Reference

8.711.1013

98.11 dm

1 The Metrohm Remote Control Language	1
1.1 General rules	1
1.2 Call up of objects	2
1.3 Triggers	3
1.4 Status messages	4
1.5 Detailed status conditions	4
1.6 Error messages	5
2 The 711 Liquino Remote Control tree	9
2.1 &Mode – Operating modes	9
2.1.1 XDOS – Extended Dosing	10
2.1.2 PIP – Pipetting	12
2.1.3 CONT – Preparation of standard solutions	13
2.1.4 SMPL - Sampling mode	14
2.1.5 SEQ - Sequence Mode	16
2.1.6 GLP - Dosing test mode	18
2.2 &Configuration	19
2.3 &Info - Instrument data	21
2.4 &Setup - Instrument setup	22
2.5 &UserMeth - user defined methods	24
2.6 &Assembly - Instrument components	25
2.7 &MemoryCard - File operations	27
2.8 &CurrData - Identifications for a determination	28
2.9 &Diagnose - Instrument diagnosis	28
3 Properties of the RS 232 Interface	29
3.1 Data Transfer Protocol	29
3.2 Handshake	29
3.3 Pin Assignment	32
3.4 What can you do if data transfer fails?	34
4 The Metrohm Remote Interface	35
4.1 Pin assignment of the Liquino Remote socket	35
4.2 Remote functions of the Liquino Remote socket	37
4.3 Liquino Remote lines during a determination	41
5 Technical specifications	43
6 Index	45

1 The Metrohm Remote Control Language

1.1 General rules

Metrohm instruments have an extensive remote control facility that allows full control of the instrument via the RS 232 interface, i.e. the instrument can receive data from an external controller or send data to an external controller. C_R and L_F are used as terminators for the data transfer. The instrument sends $2xC_R$ and L_F as termination of a data block, to differentiate between a data line which has C_R and L_F as terminators. The controller terminates its commands with C_R and L_F . If more than one command per line is sent by the controller, ";" is used as a separator between the individual commands.

The data are grouped logically and easy to understand. Thus e.g., for the selection of the dialog language, the following must be sent

&Config.Aux.Language "english"

whereby it is sufficient to only transmit the boldface characters, thus:

&C.A.L "english"

The quantities of the commands above are:

Config	configuration data
Aux	auxiliaries, various data
Language	setting the dialog language

The data are hierarchically structured (tree form). The quantities that occur in this tree are called **objects** in the following. The dialog language is an object which can be called up with the

&Config.Aux.Language

command.

If one is in the desired location in the tree, the value of the object can be queried.

&Config.Aux.Language \$Q Q means Query


The query command \$Q initiates the issuing of the value on the instrument and the value emission is triggered. Entries which start with \$, trigger something. They are thus called **triggers**.

Values of objects can not only be queried, they can also be modified. Values are always entered in quotes, for example:

&Config.Aux.Language "english"

1.2 Call up of objects

An excerpt from the object tree is represented below:



Rules	Example
The root of the tree is designated by &.	
The branches (levels) of a tree are marked with a dot (.) when calling up an object.	
When calling up an object, it is sufficient to give only as many letters as necessary to uniquely assign the object. If the call is not unequivocal, the first object in the series will be recognized.	Calling up the dialog language &Config.Aux.Language or &C.A.L
Upper- or lowercase letters may be used.	&C.A.L or &c.a.l
An object can be assigned a value. Values are signified at the beginning and end by quotes ("). They may contain up to 24 ASCII characters.	Entering the dialog language: &C.A.L"english"
Numerical values can contain up to 6 digits, a negative sign, and a decimal point. Numbers with more than 6 characters are not accepted; more than 4 decimal places are rounded off. For numbers < 1, it is necessary to enter leading zeros.	correct entry of numbers: "0.1" incorrect entry of numbers "1,5" or "+3" or ".1"
The current object remains until a new object is called.	entry of another dialog language: "deutsch"
New objects can be addressed relative to the old object: A preceding dot leads forwards to the next level in the tree.	From the root to node 'Aux': &C.A Forward from node 'Aux' to 'Prog': .P
More than one preceding dot leads one level backwards in the tree. n node backwards require n+1 preceding dots.	Jump from node 'Prog' to node 'Aux' and select a new object 'Language' at this level: ..L
If you must jump back to the root, enter a preceding &.	Change from node 'Language' via the root to node 'Mode': &M

1.3 Triggers

Triggers initiate an action on the instrument, for example, starting a process or sending data. Triggers are marked by the introductory symbol \$.

The following triggers are possible:

\$G	Go	Starts processes, for ex. starting the mode run or setting the RS 232 interface parameters
\$S	Stop	Stops processes
\$H	Hold	Holds processes
\$C	Continue	Continues processes after Hold
\$Q	Query	Queries all information from the current node in the tree forward up to and including the values
\$Q.P	Path	Queries the path from the root of the tree up to the current node
\$Q.H	Highest Index	Queries the number of son nodes of the current node
\$Q.N"i"	Name	Queries the name of the son node with index i, i = 1 - n
\$D	Detail-Info	Queries the detailed status information
\$U	qUit	Aborts the data flow of the instrument, for example, after \$Q

The triggers \$G and \$S are linked to particular objects, see the summary table.

All other triggers can be used at any time and at all locations on the object tree.

Examples:

Querying the value of the baud rate: **&Config.RSSet.Baud \$Q**
 Querying all values of the node "RSSet": **&Config.RSSet \$Q**
 Querying the path of the node "RSSet": **&Config.RSSet \$Q.P**
 Start mode: **&Mode \$G**
 Querying the detailed status: **\$D**

1.4 Status messages

In order to have an efficient control by an external control device, it must also be possible to query status conditions; they provide information on the status of the instrument. The trigger \$D initiates output of the status. Status messages consist of the global status, the detailed status and eventual error messages, e.g. \$S.Mode.SET;E26. The global status informs on the activity of the process, while the detailed status conditions show the exact activity within the process.

The following **global status conditions** are possible:

\$G	Go:	The instrument is executing the last command.
\$H	Hold:	The instrument has been held (\$H, key <meas/hold> or by an error which effects the hold status)
\$C	Continue:	The instrument has been restarted actively after hold
\$R	Ready:	The instrument has executed the last command and is ready
\$S	Stop:	A process has been aborted in an "unnatural manner". e.g. stopped or aborted because there was an error.

1.5 Detailed status conditions

Status conditions of the global \$G:

(Example: Liquino, XDOS mode)

\$G. Mde. XDOS . Inac:	instrument at the beginning or at the end of a mode.
. Run:	instrument in the XDOS mode, running.

\$G. Assembly. Fill:	Buret in filling process
\$G. Assembly. Prep:	Preparing a dosing unit
\$G. Assembly. Empty:	Emptying a dosing unit

Manual modes:

\$G.Manual.DOSING	Manual dosing
\$R.Manual.DIS	Manual distribution

Status conditions of the global \$H:

The status message of the action which has been held appears.

Status conditions of the global \$C:

The status conditions of the global \$C are identical with the ones of the global status \$G. They appear when the process has been restarted actively from the status "Hold" (\$C, key <CONTINUE/HOLD> or automatically after elimination of an error).

Status conditions of the global \$R:

\$R. Mde. XXXX. Inac: instrument in the XXXX mode, inactive.

Status conditions of the global \$S:

\$S. Mde. XXXX. Inac: The instrument gives the status from which it has been stopped. The detailed status information is therefore identical to those of the global status \$G.

1.6 Error messages

Error messages are added to the status messages and separated from them by the sign ";". These list of error messages applies to all instruments provided with the Metrohm remote control language.

E8	Card read/write error. Exit: Send new command.
E9	Wrong card, i.e. memory card from an other instrument or card has been removed/inserted during the inquiry. Exit: Replace card and send new command.
E10	The card has lost data. Exit: Replace card and end new command.
E18	Card battery low (it is between 2.37...2.64 V). Exit: Send new command. Change battery.
E20	Check exchange unit. Exit: Mount Exchange Unit (properly) or &m \$\$.
E21	Check electrode, short circuit. Replace electrode if necessary. Exit: Rectify fault or &m \$\$.
E22	Check electrode, break. See above. Exit: Rectify fault or &m \$\$.
E23	Division by zero. Exit: The error message disappears on next startup or on recalculation.
E24	Check drive unit. Exit: Connect drive unit (correctly) or &m \$\$.
E26	Manual stop. Exit: The error message disappears on next startup.
E27	Stop V reached in SET, STAT, DOS or DOC. Exit: The error message disappears on next startup.
E28	Wrong object call up Exit: Send correct path for object. Start path at root.
E29	Wrong value or no value allowed. Exit: Send correct value or call up new object.
E30	Wrong trigger, this trigger is not allowed or carrying-out of action not possible. Exit: Send correct trigger (exception: \$D) or call up new object.
E31	Command is not possible in active status. Repeat command in inactive status. Exit: Send new command.
E32	Command is not possible during titration. Repeat command during the conditioning phase or in inactive status. Exit: Send new command.
E33	Value has been corrected automatically. Exit: Send new command.
E34	instrument at the end of the titration and sample data is edited; the instrument at rest or editing during filling. Exit: &m \$\$.

RS receive errors:

- E36** Parity
Exit: <QUIT> and ensure settings of appropriate parameters at both devices are the same.
- E37** Framing error
Exit: <QUIT> and ensure settings of appropriate parameters at both devices are the same.
- E38** Overrun error. At least 1 character could not be read.
Exit: <QUIT>
- E39** The internal working-off buffer of the instrument is full (>82 characters).
Exit: <QUIT>

RS send errors:

- E42** CTS=OFF No proper handshake for more than 1 s.
Exit: <QUIT> Is the receiver switched on and ready to receive?
- E43** The transmission of the instrument has been interrupted with XOFF for at least 6 s.
Exit: Send XON or <QUIT>
- E45** The receive buffer of the instrument contains an incomplete command (L_F missing). Sending from the instrument is therefore blocked.
Exit: Send L_F or <QUIT>.
- E120** Overrange of the primary measured value (pH, U, I_{pol} , U_{pol} or T with MEAS T). The secondary measured value (temperature) may be instable as well.
Exit: Correct error or $\&m\ \$\$$.
- E121** Measuring point list overflow (more than 500 measuring points).
Exit: The error message disappears on next startup.
- E122** EP overflow.
Exit: The error message disappears on next startup or on recalculation.
- E123** Missing EP for calculation.
Exit: The error message disappears on next startup or on recalculation.
- E124** Number of EP does not correspond with the set windows.
Exit: The error message disappears on next startup or on recalculation.
- E125** Missing fix EP for calculation, has not been defined.
Exit: The error message disappears on next startup or on recalculation.
- E126** Fix-EP outside of measuring point list.
Exit: The error message disappears on next startup or on recalculation.
- E128** No new mean.
Exit: The error message disappears on next startup or on recalculation.
- E129** No new common variable, old value remains.
Exit: The error message disappears on next startup or on recalculation.
- E130** Wrong sample. For SET, KFT or DOC with preset titration direction the first measured value lies behind the endpoint.
Exit: The error message disappears on next startup.
- E131** No EP set for SET, STAT.
Exit: The error message disappears on next startup.
- E132** Silo empty and it has been started with open silo or empty silo has been opened.
Exit: Send a silo entry.
- E133** Silo full.
Exit: Send new command.

E134	No method. A method, which is required from the silo memory or in TIP, does not exist. Exit: The error message disappears on next startup.
E135	Check temp.sensor in MEAS T or with activated temperature monitoring. Temp. sensor may not be connected. Exit: Correct error or &m \$\$.
E136	Same buffer in CAL. Measured value of the second buffer differs less than 6 mV from the measured value of the first buffer. Exit: Rectify error or &m \$\$.
E137	XXX Bytes are missing so that the method, the silo line could not be stored or not enough RAM for running TIP. Exit: Send new command.
E149	Rate missing for calculation. Exit: The error message disappears on next start or on recalculation.
E150	Number of rates does not agree with the set windows or the mean rate C80 could not be calculated. Exit: The error message disappears on next start or on recalculation.
E151	Fix volume missing for calculation. Exit: The error message disappears on next startup or on recalculation.
E155	No new silo result (C24 or C25). Exit: The error message disappears on next start or on recalculation.
E157	No sequence defined in TIP. Exit: The error message disappears on next start.
E158	A second TIP has been called up in TIP. Exit: The error message disappears on next start.
E160	No new temporary variable. Exit: The error message disappears on next start.
E161	Measurement range of the secondary measured value (temperature) exceeded. The primary measured value (pH, U, Ipol, Upol) can also be unstable. Exit: Rectify error or &m \$\$.
E162	No ramp end set in DOC. Exit: The error message disappears on next start or &m \$\$.
E166	Save lines is "OFF" although a submethod of TIP includes an assignment to C24 or C25. Exit: The error message disappears on next start. Attention: The data of this sample will not be stored.
E167	Rate too high in DOS. No dispensing possible with the Exchange Unit currently mounted. Exit: The error message disappears on next start or &m \$\$.
E170	Fix time missing for calculation. Exit: The error message disappears on next start or on recalculation.
E171	Rate too low in DOS. No dispensing possible with the Exchange Unit currently mounted. Exit: The error message disappears on next start or &m \$\$.
E172	In TIP or DOS a QuickMeas was started, without defining a measuring quantity. Exit: The error message disappears on next start or &Mode.QuickMeas \$\$.
E173	The warning interval of the internal buret D0 called. Exit: Execute prep &a.b.p \$\$G or start next method.
E174	The warning interval of the external buret D1 called. Exit: Execute prep &a.b.p \$\$G or start next method.
E175	The warning interval of the external buret D2 called. Exit: Execute prep &a.b.p \$\$G or start next method.
E176	The function &Assembly.Buret.Prepare or &Assembly.Buret.Empty was interrupted manually. Exit: The error message disappears on next start.

- E177** Accessing to the memory card, the card was not (properly) inserted.
Exit: The error message disappears on next start.
- E178** The date of changing the battery of the card is expired.
Exit: The error message disappears on next start.
- E180** Memory card write-protected.
Exit: Send new command.
- E181** Memory card not formatted.
Exit: Send new command.
- E182** Memory card not accessible.
Exit: Send new command.
- E183** A directory with the same name exists already on the memory card.
Exit: Send new command.
- E184** Measured value below lower limit.
Exit: The error message disappears when again within the limits or &m \$\$.
- E185** Measured value above upper limit.
Exit: The error message disappears when again within the limits or &m \$\$.
- E186** Temperature below lower limit.
Exit: The error message disappears when again within the limits or &m \$\$.
- E187** Temperature above upper limit.
Exit: The error message disappears when again within the limits or &m \$\$.
- E188** Rate below lower limit.
Exit: The error message disappears when again within the limits or &m \$\$.
- E189** Rate above upper limit.
Exit: The error message disappears when again within the limits or &m \$\$.

2 The 711 Liquino remote control tree

2.1 &Mode – Operating modes

The remote control tree can be divided into the following main branches:

&	Root
. <u>M</u> ode	Method parameters
. <u>C</u> onfig	Instrument configuration
. <u>I</u> nfo	Current data
. <u>S</u> etup	Instrument setup
. <u>U</u> serMeth	User-defined methods
. <u>A</u> ssembly	Instrument components
. <u>C</u> urrData	Determination data
. <u>D</u> iagnosis	Instrument diagnosis

The following table lists all objects of the remote control tree. For the unambiguous designation of the objects, the underlined characters suffice. The meaning of the individual object is described here in brief, for more detailed information please refer to 'Instructions for use'. The default values of the objects are printed in boldface.

Object	Meaning	Entry range/Selection
&<u>M</u>ode	Mode	
. <u>D</u> IS	Start, stop or hold of the current method. Interrupted methods may be continued.	\$G; \$S; \$H; \$C
. <u>V</u> olume	Manual Distribution	
. <u>R</u> ate	Start or stop of the manual distribution.	
. <u>S</u> elect	Volume to be distributed	0.0.. 1.000 ..99999.9; infinite
. <u>M</u> ode	Distribution rate	0.001..166.66; max
. <u>N</u> ame	Selection of operating mode	XDOS ; PIP; CONT; SMPL; SEQ; GLP
. <u>P</u> arameter	Name of current method	read only
*	Parameters depend on the selected mode, see above.	

The structure of the **&Mode.Parameter** sub-branch depends on the selected operation mode (**&Mode.Select.Mode**). Changing the operation mode (**XDOS**, **PIP**, **CONT**, **SEQ**, **GLP**) or loading a method modifies the **&Mode.Parameter** sub-branch. A method (that can be saved or loaded at will) consists of the entire **&Mode.Parameter** sub-branch of the remote control tree.

2.1.1 XDOS – Extended Dosing

Object	Meaning	Entry range/Selection
&Mode	Mode	
:	Start, stop or hold of the current method. Interrupted methods may be continued.	\$G; \$S; \$H; \$C
L . Parameter		
. DosPara	Dosing parameters	
. Type	Selection of given parameters	volume&rate ; volume&time; time&rate
. Volume	Volume to be distributed	0.0.. 1.0 ..99999.9 mL; infinite
. Time	Dosing duration	00:00:01..99:59:59 (hh:mm:ss)
. Rate	Feeding rate	0.001..166.66 mL/min ; max
. TempMeas	Temperature measuring	on; off
. TDelta	Measuring interval	00:00:01..99:59:59 (hh:mm:ss)
. FillRate	Filling rate	0.001..166.66 mL/min ; max
. AutoFill	Automatic filling of empty burets	on ; off
. Burets	Buret setup	
. Setup	Use of one or two burets	single ; tandem
. Buret1	Buret 1	
. Drive	Drive address of first buret	DOS1 ; DOS2; DOS3; DOS4
. Duid	Descriptive name of dosing unit	8 ASCII characters
. PortUse	Type of tubing	default ; custom
. PortSpec	Port Specifications	
. FeedPort	Feeding port (outlet)	1..4
. FillPort	Filling port (inlet)	1..2..4
. WastePort	Waste port (outlet)	1..4
. Buret2	Buret 2	
. Drive	Drive address of second buret	DOS1; DOS2 ; DOS3; DOS4
. Duid	Descriptive name of dosing unit	8 ASCII characters
. PortUse	Type of tubing	default ; custom
. PortSpec	Port Specifications	
. FeedPort	Feeding port (outlet)	1..4
. FillPort	Filling port (inlet)	1..2..4
. WastePort	Waste port (outlet)	1..4
. Monitoring	Monitoring functions	
. Temp	Temperature monitoring affects the remote output lines according to the selection of the ... AlarmPins node. (upper limit reached=output 1 active, lower limit reached=output 2 active)	
. Status	Temperature monitoring	on; off
. LowLim	Lower temperature limit	-70.0 ..200.0 °C
. UpLim	Upper temperature limit	-70.0.. 200.0 °C
. Action	Action if limit is violated	none ; stop; hold; wait
. AlarmPins	Condition for alarm pin activation	none ; under; over; all
. Volume	The volume monitoring effects halting of a running method, when the given volume limit is reached.	
. Status	Volume monitoring	on; off
. CritVol	Volume limit	0.0.. 10000.0 ..99999.9

<ul style="list-style-type: none"> └─ . Prese1 └─ . IReq └─ . Prompts <ul style="list-style-type: none"> └─ . Id1 └─ . Id2 └─ . Id3 └─ . Pause 	<p>Preselections</p> <hr/> <p>Requests at method start off; id1; id1&2; all</p> <p>Request prompt for ID1 14 ASCII characters</p> <p>Request prompt for ID2 14 ASCII characters</p> <p>Request prompt for ID3 14 ASCII characters</p> <p>Start delay 00:00:00..99:59:59 (hh:mm:ss)</p>
<ul style="list-style-type: none"> └─ . Reports <ul style="list-style-type: none"> └─ . TAutoScale └─ . SLow └─ . SHigh └─ . 1 <ul style="list-style-type: none"> └─ . Select └─ . 8 <ul style="list-style-type: none"> └─ . Select 	<p>Report definitions</p> <hr/> <p>Autoscaling of temperature curve on; off</p> <p>Lower limit of temperature scale -70..200</p> <p>Upper limit of temperature scale -70..200</p> <p>Selection of report type off; result; parameter; curve; V-list; T-list; comb. list; ff</p>

2.1.2 PIP – Pipetting

Object	Meaning	Entry range/Selection
&Mode	Mode	
:	Start, stop or hold of the current method. Interrupted methods may be continued.	\$G; \$S; \$H; \$C
L . P parameter		
L . P ipPara	Pipetting parameters	
L . P ipVol	Sample volume to be pipetted	0.050.. 1.000 ...10.000 mL
L . D ilVol	Solvent volume for dilution	0.000 ..99.999 mL
L . A irGap	Air gap length	5.. 10 ..100 mm
L . B uret	Buret setup	
L . D rive	Drive address	DOS1 ; DOS2; DOS3; DOS4
L . D Uid	Descriptive name of dosing unit	8 ASCII characters
L . P ortUse	Type of tubing	default ; custom
L . P ortSpec	Port Specifications	
L . F eedPort	Feeding port (outlet)	1 ..4
L . F illPort	Filling port (inlet)	1.. 2 ..4
L . W astePort	Waste port (outlet)	1.. 4
L . R ates	Pipetting rates	
L . F eedRate	Feeding rate	0.001..166.66 mL/min; max
L . F illRate	Filling rate	0.001..166.66 mL/min; max
L . P resel	Preselections	
L . V Req	Volume request	off ; on
L . R eports	Report definitions	
L . 1		
L . S elect	Selection of report type	off ; result; parameter; ff
L . 3		
L . S elect		

2.1.3 CONT – Preparation of standard solutions

<i>Object</i>	<i>Meaning</i>	<i>Entry range/Selection</i>
&Mode	Mode	
:	Start, stop or hold of the current method. Interrupted methods may be continued.	\$G; \$S; \$H; \$C
L . Parameter		
. ContPara	Parameters for CONT-Mode	
. Type	Type of solution preparation	subst.conc ; mass conc.; mass fraction; volume fraction
. Unit	Unit of concentration	Mol/L ; ppm; g/L; %
. Content	Target concentration	0.001..1.0..9999.9
. TargetVol	Target volume	0.01..10.0..9999
. Standard	Concentration of the source solution	0.001..1.0..9999.9
. Factor	Correction factor	0.000001..1.0..999999.0
. MolarMass	Molar mass of the dissolved substance	0.0001..1.0..999999
. Density	Density of the solvent	0.1..1.0..99.9999 g/mL
. Buret	Buret setup	
. Drive	Drive address	DOS1 ; DOS2; DOS3; DOS4
. Duid	Descriptive name of dosing unit	read only
. PortUse	Type of tubing	default ; custom
. PortSpec	Port Specifications	
. FeedPort	Feeding port (outlet)	1..4
. FillPort	Filling port (inlet)	1..2..4
. WastePort	Waste port (outlet)	1..4
. Rates	Distribution rates	
. FeedRate	Feeding rate	0.001..166.66 mL/min; max
. FillRate	Filling rate	0.001..166.66 mL/min; max
. Presel	Preselections	
. IReq	Requests at method start	off ; id1; id1&2; all
. SReq	Sample weight request	off; on
. CntReq	Target content request	off; on
. Prompts		
. Id1	Request prompt for ID1	14 ASCII characters
. Id2	Request prompt for ID2	14 ASCII characters
. Id3	Request prompt for ID3	14 ASCII characters
. Reports	Report definitions	
. 1		
. Select	Selection of report type	off ; result;parameter; ff
. 3		
. Select		

2.1.4 SMPL - Sampling mode

Object	Meaning	Entry range/Selection
&Mode	Mode	
:	Start, stop or hold of the current method. Interrupted methods may be continued.	\$G; \$S; \$H; \$C
L		
. Parameter		
. Smp1Para	Sampling parameters	
. NumSamples	Number of samples to be processed	1..99999; infinite
. SMPLmode	Manual or automated operation	manual ; auto
. IntervalT	Interval time between each sample	00:00:00..99:59:59 (hh:mm:ss)
. Sampling	Sampling operation	
. SmpLVol	Sample volume	0.001.. 10.0 ..50.0 mL
. SolvVol	Solvent volume for sample dilution	0.001.. 10.0 ..50.0 mL
. SmpLRinseVol	Rinsing volume of sample solution	0.001.. 10.0 ..50.0 mL
. SmpLRinseRep	Repetitions of sample rinsing	1.. 3 ..10
. SmpLAspRate	Aspiration rate of sample solution	0.001..166.66 mL/min; max
. SmpLDosRate	Dosing rate of sample solution	0.001..166.66 mL/min; max
. SolvFillRate	Filling rate of solvent aspiration	0.001..166.66 mL/min; max
. AirgapLen	Air gap length	1.. 10 ..100 mm
. CellRinsing	Measuring cell rinsing	
. PumpExt	Use of external pumps for rinsing the cell	on; off
. SolvRinseRep	Repetitions of solvent rinsing	1.. 3 ..10
. External	External pumps settings	
. RinseTime	Duration of cell rinsing (solvent distribution)	00:00:00..99:59:59 (hh:mm:ss)
. AspTime	Duration of cell emptying (solvent aspiration)	00:00:00..99:59:59 (hh:mm:ss)
. Internal	Settings for rinsing without external pumps	
. WasteAspVol	Solvent aspiration volume (emptying)	0.1.. 10 ..500 mL
. WasteAspRate	Aspiration rate	0.001..166.66 mL/min max
. SolvRinseVol	Solvent rinsing volume (distribution)	0.1.. 50 ..250 mL
. SolvDosRate	Distribution rate	0.001..166.66 mL/min ; max
. Burets	Buret setup	
. Smp1DUid	Descriptive name of dosing unit (DOS1)	8 ASCII characters
. BeforeSmp1	Dosing of auxiliar solution before sample distribution	
. 1		
. Drive	Drive address	off ; DOS2; DOS3; DOS4
. DUid	Identification of dosing unit	8 ASCII characters
. DosVol	Volume of auxiliar solution	0.0001.. 1.0 ..99999.9 mL
. FeedRate	Feeding rate	0.001..166.66 mL/min; max
. FillRate	Filling rate	0.001..166.66 mL/min; max
. WaitT	Waiting time after dosing	00:00:00 ..99:59:59 (hh:mm:ss)

L . 3			
. Drive
.
. FillRate
. WaitT
L . AfterSmpl	Dosing of aux. solution after sample distribution		
. 1			
. Drive	Drive address		off ; DOS2; DOS3; DOS4
. Duid	Descriptive name of dosing unit		8 ASCII characters
. DosVol	Volume of auxiliar solution		0.0001..1.0..99999.9 mL
. FeedRate	Feeding rate		0.001..166.66 mL/min; max
. FillRate	Filling rate		0.001..166.66 mL/min; max
. WaitT	Waiting time after dosing		00:00:00..99:59:59 (hh:mm:ss)
. . .			
. 3			
. Drive
.
. FillRate
. WaitT
L . AutoCtrl	Control of external instruments		
. SmplCtrl	Use of a 730 Sample Changer		off ; on
. OnSmplCtrl	Sample Changer controls		
. StartSmplChng	Interface for Sample changer remote control		remote ; RS
. RSstartStr	Text string for starting the Sample Changer via RS232		Default: &M;\$G
. AwaitAck	Interface for Sample changer acknowledgement		remote; RS; off
. RSawaitStr	Acknowledge string to receive		Default: * (= any characters)
. ExtDevice	Use of an external device		on; off
. OnExtDev	Control of a measuring or titration instrument		
. StartExtDev	Interface for starting the external device		remote ; RS; off
. RSstartStr	Command string for starting via RS232		Default: &M;\$G
. AwaitAck	Interface for external device acknowledgement		remote; RS ; off
. RSawaitStr	Acknowledge string to receive		Default: * (= any characters)
L . Presel	Preselections		
. IReq	Requests at method start		off ; id1; id1&2; all
. Prompts			
. Id1	Request prompt for ID1		14 ASCII characters
. Id2	Request prompt for ID2		14 ASCII characters
. Id3	Request prompt for ID3		14 ASCII characters
L . Reports	Report definitions		
. 1			
. Select	Selection of report type		off ; result; parameter; ff
. 3			
. Select			

2.1.5 SEQ - Sequence Mode

Object	Meaning	Entry range/Selection
&Mode :	Mode Start, stop or hold of the current method. Interrupted methods may be continued.	\$G; \$S; \$H; \$C
L . P arameter		
. S equence	Sequence parameters	
. A utoStart	Number of sequences to be executed	1..99999, infinite
. 1	Sequence node	
. S elect	Command selection	NO OPERATION; METHOD; PAUSE; TIMER; SCAN; CONTROL; ENDSEQ
. 75	<i>Up to 75 nodes can make up a sequence. Each sequence node carries one of the following commands or child-nodes respectively. The content of the 'Select'-node determines the command and its parameters.</i>	
. S elect		
L . M ethod	Method to be executed	
L . N ame	Method name	8 ASCII characters
L . P ause	Sequence delay	
L . T ime	Waiting time before continuation of the sequence	hh:mm:ss
L . T imer	Timer function	
. D ate	Date for continuation of the sequence	YYYY-MM-DD
. T ime	Time for continuation of the sequence	hh:mm:ss
L . S can	Interface scan	
. A ddress	Interface address	Remote; RS
. P attern	Bit pattern for Remote scan	5 bits (1, 0 or *)
. S tring	Text string for RS232 scan	..15 ASCII characters Default: * (= any characters)
L . C trl	Output data or signals	
. A ddress	Interface address	Remote; RS
. P attern	Bit pattern for setting the Remote lines	4 bits (1, 0 or *)
. S tring	String to be sent via RS232 interface	15 ASCII characters

Example of a single command in a sequence:

```

: :
| | . 5 — Step 5 of the sequence
: : | . Select "METHOD" — Command selection
| | . Method
| | | . Name "Chloride" — Method name
    
```

<ul style="list-style-type: none"> ⋮ └─ .Prese1 └─ .IReq └─ .Prompts <ul style="list-style-type: none"> └─ .Id1 └─ .Id2 └─ .Id3 └─ .Reports <ul style="list-style-type: none"> └─ .1 ⋮ └─ .Select └─ .2 └─ .Select 	<p>Preselections</p> <hr/> <p>Requests at method start off; id1; id1&2; all</p> <p>Request prompt for ID1 8 ASCII characters</p> <p>Request prompt for ID2 8 ASCII characters</p> <p>Request prompt for ID3 8 ASCII characters</p> <p>Report definitions</p> <hr/> <p>Selection of report type off; parameter; ff</p>
---	--

2.1.6 GLP - Dosing test mode

Object	Meaning	Entry range/Selection
&Mode	Mode	
<ul style="list-style-type: none"> ⋮ * 	Start, stop or hold of the current method. Interrupted methods may be continued.	\$G; \$S; \$H; \$C
L . Parameter		
<ul style="list-style-type: none"> └ . GlpPara 	Dosing test parameters	
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . Drive 	Drive address	DOS1 ; DOS2; DOS3; DOS4
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . IDid 	Descriptive name of dosing unit	8 ASCII characters
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . PortUse 	Type of tubing	default ; custom
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . PortSpec 	Port Specifications	
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . FeedPort 	Feeding port (outlet)	1..4
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . FillPort 	Filling port (inlet)	1..2..4
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . WastePort 	Waste port (outlet)	1..4
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . LiqTemp 	Temperature of the distributed liquid	-70.0.. 20.0 ..200.0
<ul style="list-style-type: none"> └ . Options 	Dosing options	
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . Liquid 	Type of liquid	water ; other
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . OtherLiq 	Settings for non-aqueous liquids	
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . Name 	Name of liquid	13 ASCII characters
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . Density 	Density of the liquid	0.01.. 1.0 ..20.0 g/L
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . FeedRate 	Feeding rate	0.001..166.66 mL/min; max
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . FillRate 	Filling rate	0.001..166.66 mL/min; max
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . WeightFrom 	Sample weight input by	balance ; keyboard
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . BalID 	Balance identification	8 ASCII characters
<ul style="list-style-type: none"> └ . Presel 	Preselections	
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . IReq 	Requests at method start	off; id1; id1&2; all
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . Prompts 		
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . Id1 	Request prompt for ID1	read only
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . Id2 	Request prompt for ID2	read only
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . Id3 	Request prompt for ID3	read only
<ul style="list-style-type: none"> └ . Reports 	Report definitions	
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . 1 		
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . Select 	Selection of report type	off ; result; parameter; curve; detail; ff
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . 7 		
<ul style="list-style-type: none"> └ <ul style="list-style-type: none"> └ . Select 		

2.2 &Configuration

Object	Meaning	Entry range/Selection
&Config		
└─ P eriphUnit	Settings of peripheral units	
└─└─ C harSet	External printer definition	IBM ; Epson; Seiko; Citizen; HP DeskJ
└─└─ B alance	Selection of balance	Sartorius ; Mettler; Mettler AT; AND; Precisa
└─ A ux	Basic settings	
└─└─ L anguage	Dialog language	english ; deutsch
└─└─ S et	Date and time setting	\$G
└─└─└─ D ate	Date	YYYY-MM-DD
└─└─└─ T ime	Time	hh:mm:ss
└─└─ D evName	Device label	8 ASCII characters
└─└─ P rog	Program version	711.011(read only)
└─└─ C ontrast	Display contrast	0..3..7
└─└─ B eeper	Beeper	off; on
└─└─ 1 stTimePrep	Automatic PREP warning	off; on
└─ M anual	Settings for manual operation of the dosing drives	
└─└─ 1	Drive address (DOS1)	
└─└─└─ G eneral		
└─└─└─└─ D osUnit	Identification of the mounted Dosing unit	default ; any 12 ASCII characters
└─└─└─└─ F eedTo	Feeding port of the dosing unit	1..4
└─└─└─└─ F illFrom	Filling port of the Dosing unit	1..2..4
└─└─└─└─ F illRate	Filling rate	0.001..166.666 mL/min; max
└─└─└─└─ A utoFill	Automatic filling after distribution	on; off
└─└─└─ D osing	Settings for manual dosing	
└─└─└─└─ T ype	Increasing or constand feeding rate	const; increasing
└─└─└─└─ F eedRate	Nominal feeding rate	0.001..166.666 mL/min; max
└─└─└─└─ S lope	Slope for increasing feeding rate	1..4..10
└─└─└─		
└─└─└─ 4	Drive address (DOS4)	
└─└─└─└─ G eneral	..	
└─└─└─└─		
└─└─└─└─ S lope	..	
└─ D osUnits	Configuration of the Dosing units	
└─└─ D USelect	Selected Dosing unit	12 ASCII characters
└─└─ 1	Number of the Dosing unit	
└─└─└─ I D	Identification of the Dosing unit	default ; any 12 ASCII characters
└─└─└─ T ubes	Tube definitions	
└─└─└─└─ L en1	Tube length (port 1)	0.400..9999 mm
└─└─└─└─ D ia1	Tube diameter (port 1)	0.1..2.0..9.0 mm

<ul style="list-style-type: none"> ⋮ ├─ Len2 ├─ Di a2 ├─ Len3 ├─ Di a3 ├─ Len4 ├─ Di a4 ⋮ ├─ 10 ├─ ID ├─ Tubes ⋮ ├─ RSset ├─ Baud ├─ DataBit ├─ StopBit ├─ Parity ├─ Handsh ├─ Reports ├─ Header1 ├─ Header2 ├─ Header3 ├─ GLPErrLim ├─ 2ml ├─ ErrMax ├─ Slope ├─ Inter ├─ 5ml ├─ ErrMax ├─ Slope ├─ Inter ├─ 10ml ├─ ErrMax ├─ Slope ├─ Inter ├─ 20ml ├─ ErrMax ├─ Slope ├─ Inter ├─ 50ml ├─ ErrMax ├─ Slope ├─ Inter 	<ul style="list-style-type: none"> Tube length (port 2) 0..250..9999 mm Tube diameter (port 2) 0.1..2.0..9.0 mm Tube length (port 3) 0..9999 mm Tube diameter (port 3) 0.1..0.0..9.0 mm Tube length (port 4) 0..9999 mm Tube diameter (port 4) 0.1..0..9.0 mm ⋮ Number of the Dosing unit Identification of the Dosing unit default; any 12 ASCII characters Tube definitions ⋮ ⋮ Settings of the RS232 serial interface \$G <hr/> Baud rate 57600; 38400; 19200; 9600; 4800; 2400; 1200; 600; 300 Data bits 7; 8 Stop bits 1; 2 Parity even; odd; none Handshake Hardware; SWchar; SWline; none ⋮ Configuration of printed reports <hr/> First line of the report header 21 ASCII characters Second line of the report header 21 ASCII characters Third line of the report header 21 ASCII characters ⋮ Configuration of the GLP limits <hr/> Limits for 2 mL-Dosing units Maximum volume error 1..6..999 μL Maximum slope deviation 0.001..0.003..0.1 Maximum intercept (Y-axis) 1..6..99 μL Limits for 5 mL-Dosing units Maximum volume error 1..15..999 μL Maximum slope deviation 0.001..0.003..0.1 Maximum intercept (Y-axis) 1..15..99 μL Limits for 10 mL-Dosing units Maximum volume error 1..30..999 μL Maximum slope deviation 0.001..0.003..0.1 Maximum intercept (Y-axis) 1..20..99 μL Limits for 20 mL-Dosing units Maximum volume error 1..60..999 μL Maximum slope deviation 0.001..0.003..0.1 Maximum intercept (Y-axis) 1..30..99 μL Limits for 50 mL-Dosing units Maximum volume error 1..150..999 μL Maximum slope deviation 0.001..0.003..0.1 Maximum intercept (Y-axis) 1..50..99 μL
---	---

2.3 &Info - Instrument data

<i>Object</i>	<i>Meaning</i>	<i>Entry range/Selection</i>
&Info		
. Report	Report output	\$G
. Select	Report selection	result; parameter; curve: V-list; T-list; comb. list; config; user methods; current dir.; memory card; ff
. ActualInfo	Actual data	
. Buret	Buret data	
. 1	DOS1	
. Type	Type of dosing device	read only
. State	Drive status	read only
. CylVol	Cylinder volume	read only
. ActVol	Actual volume	read only
. ActPort	Actual port	read only
. LastPistDir	Piston direction	read only
. PistPos	Actual piston position	read only
. 4	DOS4	
. Type	..	
. PistPos	..	
. SMPL	SMPL mode data	
. Current	Current sample counter	read only
. SEQ	SEQ mode data	
. Current	Current sample counter (auto start)	read only
. Inputs	Remote interface data	
. Status	Status of input lines	read only
. Outputs	Remote interface data	
. Status	Status of the output lines	read only
. Display	Display content	
. L1	Display line 1	24 ASCII characters
. L2	Display line 2	24 ASCII characters

2.4 &Setup - Instrument setup

Object	Meaning	Entry range/Selection																																																																				
&Setup																																																																						
K eycode	Sending of key codes	on; off																																																																				
:																																																																						
:																																																																						
	<table border="1"> <thead> <tr> <th>Code</th> <th>Key</th> <th>Code</th> <th>Key</th> </tr> </thead> <tbody> <tr><td>1</td><td><HOLD / CONTINUE></td><td>16</td><td><7 / DOSING></td></tr> <tr><td>2</td><td><STOP></td><td>17</td><td><4 / PREP></td></tr> <tr><td>3</td><td><START></td><td>18</td><td><1 / CUR DATA></td></tr> <tr><td>4</td><td><CONFIG></td><td>19</td><td><0 / MEM CARD></td></tr> <tr><td>5</td><td><PARAM></td><td>20</td><td><HOME></td></tr> <tr><td>6</td><td><USER METHOD></td><td>21</td><td><END></td></tr> <tr><td>7</td><td></td><td>22</td><td><CLEAR></td></tr> <tr><td>8</td><td><9 / TEMP></td><td>23</td><td><ENTER / YES></td></tr> <tr><td>9</td><td><6 / EXCH></td><td>24</td><td><↑></td></tr> <tr><td>10</td><td><3 / EMPTY></td><td>25</td><td><↓></td></tr> <tr><td>11</td><td><* / MODE></td><td>26</td><td><SELECT / DOS / NR></td></tr> <tr><td>12</td><td><8 / DISPENSE></td><td>27</td><td><QUIT / NO></td></tr> <tr><td>13</td><td><5 / FILL></td><td>28</td><td><MANUAL / MAN OFF></td></tr> <tr><td>14</td><td><2 ></td><td>29</td><td><ALPHA></td></tr> <tr><td>15</td><td><./ PRINT></td><td>30</td><td><INSERT></td></tr> <tr><td></td><td></td><td>31</td><td><DELETE></td></tr> </tbody> </table>	Code	Key	Code	Key	1	<HOLD / CONTINUE>	16	<7 / DOSING>	2	<STOP>	17	<4 / PREP>	3	<START>	18	<1 / CUR DATA>	4	<CONFIG>	19	<0 / MEM CARD>	5	<PARAM>	20	<HOME>	6	<USER METHOD>	21	<END>	7		22	<CLEAR>	8	<9 / TEMP>	23	<ENTER / YES>	9	<6 / EXCH>	24	<↑>	10	<3 / EMPTY>	25	<↓>	11	<* / MODE>	26	<SELECT / DOS / NR>	12	<8 / DISPENSE>	27	<QUIT / NO>	13	<5 / FILL>	28	<MANUAL / MAN OFF>	14	<2 >	29	<ALPHA>	15	<./ PRINT>	30	<INSERT>			31	<DELETE>	
Code	Key	Code	Key																																																																			
1	<HOLD / CONTINUE>	16	<7 / DOSING>																																																																			
2	<STOP>	17	<4 / PREP>																																																																			
3	<START>	18	<1 / CUR DATA>																																																																			
4	<CONFIG>	19	<0 / MEM CARD>																																																																			
5	<PARAM>	20	<HOME>																																																																			
6	<USER METHOD>	21	<END>																																																																			
7		22	<CLEAR>																																																																			
8	<9 / TEMP>	23	<ENTER / YES>																																																																			
9	<6 / EXCH>	24	<↑>																																																																			
10	<3 / EMPTY>	25	<↓>																																																																			
11	<* / MODE>	26	<SELECT / DOS / NR>																																																																			
12	<8 / DISPENSE>	27	<QUIT / NO>																																																																			
13	<5 / FILL>	28	<MANUAL / MAN OFF>																																																																			
14	<2 >	29	<ALPHA>																																																																			
15	<./ PRINT>	30	<INSERT>																																																																			
		31	<DELETE>																																																																			
:	The key code is sent together with a log time which is measured from the power-up of the instrument. Format: Log time (in sec.) - key code																																																																					
T race	Message output on changed values	on; off																																																																				
T ree	Sending format of path info																																																																					
S hort	Short format of path info	on; off																																																																				
C hangedOnly	Paths of modified values only	on; off																																																																				
L ock	Locking of key functions																																																																					
K eyboard	Lock all keys	on; off																																																																				
C onfig	Lock <CONFIG> key	on; off																																																																				
P arameter	Lock <PARAM> key	on; off																																																																				
U serMeth	User method functions																																																																					
R ecall	Lock method loading	on; off																																																																				
S tore	Lock method storing	on; off																																																																				
D elete	Lock method deletion	on; off																																																																				
D isplay	Lock display functions	on; off																																																																				
M ode	Waiting interval settings																																																																					
S tartWait	Waiting time after method start	on; off																																																																				
F inWait	Waiting time after run	on; off																																																																				
S endMeas	Automatic sending of measured values																																																																					
S tatus	Connect/disconnect sending	on; off																																																																				
I nterval	Time interval of automatic sending	1..10..16200 s																																																																				
:																																																																						

2.5 &UserMeth - user defined methods

Object	Meaning	Entry range/Selection
&UserMeth		
File operations		
.FreeMem	Memory available	read only
.Recall	Loading of methods	\$G
.Name	Method name	8 ASCII characters
.Store	Saving of methods	\$G
.Name	Method name	8 ASCII characters
.Delete	Deleting of methods	\$G
.Name	Method name	8 ASCII characters
.DeIA11	Deleting of all methods	\$G
.List	List of stored methods	
.1		
.Name	Method name	read only
.Mode	Mode	read only
.Bytes	Method size in bytes	read only
.Checksum	Checksum of method	read only
.x	up to x methods, depending on memory available	
.Name
.Checksum

Remark on file operations

First set the file names (node: .Name), then send the trigger \$G for the corresponding function:.

Example: Loading a method

```
&U.R.N"Chloride"
&U.R;$G
```

2.6 &Assembly - Instrument components

Object	Meaning	Entry range/Selection
&Assembly		
. Drive	Drive address	1..4
. DosingUnit	Identification of Dosing unit	12 ASCII characters
. Rates	Rates of the dosing drive	
. Forward	Forward rates	
. Select	Type of rate control	digital ; waitticks
. Digital	Digital rate	0.001..166.666 mL/min; max
. Waitticks	Rate in wait ticks	225 ..65535
. Reverse	Reverse rates	
. Select	Type of rate control	digital ; waitticks
. Digital	Digital rate	0.001..166.666 mL/min; max
. Waitticks	Rate in wait ticks	225 ..65535
. CompiledSteps	Compiled buret functions	
. Dos	Dosing	\$G; \$S; \$H
. Volume	Volume to be distributed	0.0 ..99999.9 mL
. FillFrom	Filling port (inlet)	1..2..4
. FeedTo	Dosing port (outlet)	1..4
. Fill	Buret filling	\$G; \$S; \$H
. FillFrom	Filling port (inlet)	1..2..4
. ExchPos	Stopcock to exchange position	\$G; \$S; \$H
. FillFrom	Filling port (inlet)	1..2..4
. Prep	PREP function	\$G; \$S; \$H
. Waste	Waste port	1..4
. FillFrom	Filling port (inlet)	1..2..4
. Empty	Buret emptying	\$G; \$S; \$H
. ToPort	Empty port (outlet)	1..4
. BasicSteps	Basic buret functions	
	<i>Use these functions with great care. Improper use of the piston functions may lead to hardware damage of the dosing unit or the dosing drive.</i>	
. Cock	Stopcock positioning	\$G
. Pos	Target port	1..4
. Dir	Direction of stopcock rotation	desc ; asc
. GoPos	Piston positioning	\$G; \$S; \$H
. Pos	Target position	0 ..10000 pulses
. MoveRel	Relative piston positioning	\$G; \$S; \$H
. Delta	Number of pulses to move the piston	-50.0.. 0 ..50.0
. ToEnd	Stopcock to end position	\$G; \$S; \$H
. ZeroAdjust	Zero adjustment	\$G; \$S; \$H

⋮			
├	.Temp	Temperature read-out	read only
├	.Outputs		
├	.SetLines	Setting the output lines (Remote socket)	\$G
├	├ .L1	Line 1 = output 0	off ; active; inactive
├	├ .L2	Line 2 = output 1	off ; active; inactive
├	├ .L3	Line 3 = output 2	off ; active; inactive
├	├ .L4	Line 4 = output 3	off ; active; inactive
├	├ .L5	Line 5 = output 4	off ; active; inactive
├	├ .L6	Line 6 = output 5	off ; active; inactive
├	├ .L7	Line 7 = output 6	off ; active; inactive
├	├ .L8	Line 8 = output 7	off ; active; inactive
├	.ResetLines	Reset all output lines	\$G

How to use the remote output lines

During the run of a single mode (XDOS, SMPL, PIP, CONT, GLP) the remote output lines will be set automatically to indicate the instrument status.

In a sequence run (SEQ mode) the remote output lines may be set via the CONTROL command or by remote access via the RS232 interface.

One can set the output lines separately. First define the line status of any line (nodes **&Assembly.Outputs.SetLines.Lx**), then set the status of all lines by sending the trigger **'&Assembly.Outputs.SetLines;\$G'**.

The selections of the output lines are:

<i>off</i>	lock the current line status
<i>active</i>	set line active (1)
<i>inactive</i>	reset line to inactive (0)

Locked lines (selection: *off*) are not affected by the trigger **'...SetLines;\$G'**.

The trigger **'&Assembly.Outputs.ResetLines;\$G'** resets all output lines (even locked lines) to inactive.

2.7 &MemoryCard - File operations

&MemoryCard

├─ . Recall	Loading of methods	\$G
├─┬─ . Name	Method name	8 ASCII characters
├─ . Store	Saving of methods	\$G
├─┬─ . Name	Method name	8 ASCII characters
├─ . Delete	Deleting of methods	\$G
├─┬─ . Name	Method name	8 ASCII characters
├─ . ChangeDir	Change directory	\$G
├─┬─ . Name	Directory name	10 ASCII characters
├─┬─ . Checksum	Checksum of directory	\$G
├─┬─┬─ . Value	Checksum value	read only
├─ . CreateDir	Create directory	\$G
├─┬─ . Name	Directory name	10 ASCII characters
├─ . DelDir	Delete directory	\$G
├─┬─ . Name	Directory name	10 ASCII characters
├─ . Backup	Backup of internal memory	\$G
├─┬─ . Name	Directory name	10 ASCII characters
├─ . Reload	Reload backup from card	\$G
├─┬─ . Name	Directory name	10 ASCII characters
├─ . Format	Formatting the card	\$G
├─┬─ . CardLabel		
├─┬─┬─ . Name	Card label	8 ASCII characters
├─ . FreeMem	Memory available	read only
├─ . BatteryChange	Setting the date for next battery change	\$G
├─┬─ . Date	Date	YYYY:MM:DD
├─ . List	Lists of directories	
├─┬─ . Card	Content of the memory card	
├─┬─┬─ . 1	First directory	
├─┬─┬─┬─ . Name	Directory name	read only
├─┬─┬─┬─ . Bytes	Directory size in bytes	read only
├─┬─┬─┬─		
├─┬─┬─┬─ . x	Last directory	
├─┬─┬─┬─┬─ . Name	Directory name	read only
├─┬─┬─┬─┬─ . Bytes	Directory size in bytes	read only
├─┬─ . ActDir	Content of the current directory	
├─┬─┬─ . 1		
├─┬─┬─┬─ . Name	Directory name	read only
├─┬─┬─┬─ . Mode	Mode	read only
├─┬─┬─┬─ . Bytes	File size in bytes	read only
├─┬─┬─┬─ . Checksum	Checksum of method	read only
├─┬─┬─┬─		
├─┬─┬─┬─ . x		
├─┬─┬─┬─		

2.8 &CurrData - Identifications for a determination

&CurrData

• <u>I</u>d1	Identification 1	8 ASCII characters
• <u>I</u>d2	Identification 2	8 ASCII characters
• <u>I</u>d3	Identification 3	8 ASCII characters
• <u>S</u>mp1	Sample size (CONT mode only)	0.00001..1.000..1000.000
• <u>C</u>nt	Content (CONT mode only)	0.001..1.000..9999.9

2.9 &Diagnose - Instrument diagnosis

&Diagnose

• <u>I</u>nit	RAM initialization	\$G
• <u>S</u>elect	Selection of memory range to initialize	all , param, config, user meth, setup, assembly
• <u>R</u>amTest	RAM integrity test	\$G
• <u>L</u>cdTest	Display test	\$G, \$\$
• <u>C</u>ontrastTest	Display contrast test	\$G, \$\$
• <u>K</u>eyTest	Keyboard test	\$G, \$\$
• <u>I</u>oTest	Remote socket test	\$G
• <u>R</u>sTest	RS232 interface test	\$G
• <u>E</u>busTest	Dosino interface test	\$G
• <u>B</u>eeperTest	Beeper test	\$G, \$\$
• <u>C</u>ardTest	Memory card test	\$G
• <u>C</u>ardIfTest	Memory card interface test	\$G
• <u>D</u>osino	Dosino test	\$G, \$\$
• <u>T</u>empAdjust	Temperatur measuring test	\$G
• <u>S</u>imulateKey	Keycode simulation	1..31, see page 14
• <u>P</u>owerOn	'Power on' simulation	\$G
• <u>I</u>ns<trno< b=""></trno<>	Instrument number	not accessible via RS232 interface

3 Properties of the RS 232 Interface

3.3 Data Transfer Protocol

Metrohm instruments are configured as DTE (Data Terminal Equipment).

The RS 232 interface has the following technical specifications:

- Data interface according to the RS 232C standard, adjustable transfer parameters, see pages 8 and 154.
- Max. line length: 512 characters
- Control characters:
 - C_R (ASCII DEC 13)
 - L_F (ASCII DEC 10)
 - XON (ASCII DEC 17)
 - XOFF (ASCII DEC 19)
- Cable length: max. 15 m

Start	7 or 8 Data Bit	Parity Bit	1 or 2 Stop Bit
-------	-----------------	------------	-----------------


Only a shielded data cable (for example, METROHM D.104.0201) may be used to couple Metrohm instruments with foreign devices. The cable shield must be properly grounded on both instruments (pay attention to current loops; always ground in a star-head formation). Only plugs with sufficient shielding may be used (for example, METROHM K.210.0381 with K.210.9045).

3.4 Handshake


Software-Handshake, SWchar

Handshake inputs on the instrument (CTS) are not checked.
 Handshake outputs (DTR, RTS) are set by the instrument.
 The instrument sends XOFF when its input buffer contains 384 characters. After this it can receive 128 extra characters (including L_F).

Metrohm instruments as Receiver :



Metrohm instrument as Sender :




max. characters: 2 characters at 300...9600 baud
 16 characters at ≥ 19200 baud


Software-Handshake, SWiine

Handshake input ports on the instrument (CTS) are not checked.
 Handshake output ports (DTR, RTS) are set by the instrument.
 The instrument has an input buffer which can accept up to 512 characters.

Metrohm instrument as Receiver :




Metrohm instrument as Sender:




The instrument's transmission can be stopped by the external device with XOFF. After XOFF is received the Metrohm instrument completes sending the line already started. If data output is disabled for more than 6 s by XOFF, E43 appears in the display.

Hardware-Handshake, HWs

Metrohm instrument as Receiver :




Metrohm instrument as Sender:



The data flow can be interrupted by deactivating the CTS line.

3.5 Pin Assignment

RS232C Interface



Protective earthing
Direct connection from cable plug to the protective ground of the instrument.


Polarity allocation of the signals

- Data lines (TxD, RxD)
 - voltage negative (<-3 V): signal state "ON"
 - voltage positive (> +3 V): signal state "ZERO"
- control or message lines (CTS, RTS, DTR)
 - voltage negative (<-3 V): OFF state
 - voltage positive (> +3 V): ON state

In the transitional range from +3 V to -3 V the signal state is undefined.

Driver 14C88 according to EIA RS 232C specification
Receiver 14C89 " "

Contact arrangement at plug (female) for RS 232C socket (male)



View of soldered side of plug

Ordering numbers:
K.210.0381 and K.210.9045

No liability whatsoever will be accepted for damage or injury caused by improper interconnection of instruments.


3.6 What can you do if data transfer fails?

Problem	Questions for remedial action
No characters can be received on a connected printer.	<ul style="list-style-type: none"> - Are the instruments switched on and cables plugged in correctly? - Is the printer set to "on-line"? - Are baud rate, data bit and parity the same on both instruments? - Is the handshake set properly? <p>If everything seems to be OK, try to print a report with the key sequence <PRINT> <ENTER>.</p>
No data transmission and the display of the instrument shows an error message.	<ul style="list-style-type: none"> - error 42: Transmission error. Is the printer set to "on-line"? Is the connection cable properly wired? - error 43: Data output of the instrument disabled for longer than 6 s by XOFF. - error 36-39: Receive error. Are the RS settings the same on both devices?
The received characters are garbled.	<ul style="list-style-type: none"> - Are the RS settings the same on both devices? - Has the correct printer been selected? - Data transfer has been interrupted on the hardware side during the printout of a curve. Re-establish connections and switch printer off/on.
Wrong line spacing.	The printer does not emulate completely the preset mode. Usually these problems arise with the IBM mode. Set the printer to a different mode (e.g. Epson).
Printout of a curve is not ok. Other reports are printed ok.	<p>Handshake is necessary for the printout of curves.</p> <ul style="list-style-type: none"> - Is your cable correctly wired? (The DTR of the printer has to be connected to the CTS of the instrument.) - Set "HWs" for the handshake of the instrument. Configure the printer such that its DTR is set (possibly with DIP switches).

4 The Metrohm Remote Interface

4.1 Pin assignment of the Liquino Remote socket

	external	Function
<p>Inputs</p>	pin 21 (Input 0)	Start
	pin 9 (Input 1)	Stop
	pin 22 (Input 2)	Hold/Cont.
	pin 10 (Input 3)	ext. ready
	pin 23 (Input 4)	Sample changer ready
	pin 11 (Input 5)	not used in single modes
	pin 24 (Input 6)	not used in single modes
	pin 12 (Input 7)	not used in single modes
<p>Remote functions for remote control in single modes</p>		
<p>Scannable remote lines in SEQ mode</p>		
<p>Voltage</p>	pin 15	$I \leq 200 \text{ mA}$ 0 V: active 5 V: inactive
	pin 14	
	pin 25	



Ordering numbers for plug:
K.210.9004 (shell) and K.210.002

No liability whatsoever will be accepted for damage caused by improper interconnection of instruments.

4.2 Remote functions of the Liquino Remote socket

Pin 1	Output 6	external Start (in SMPL) **
Pin 2	Output 7	Start Sample Changer (in SMPL) **
Pin 3	Output 4	activate 1 (external pump in SMPL) **
Pin 4	Output 2	lower limit reached (in SMPL)
Pin 5	Output 0	ready
Pin 6	—	not used
Pin 7	—	not used
Pin 8	—	not used
Pin 9	Input 1	Stop
Pin 10	Input 3	external ready (in SMPL) *
Pin 11	Input 5	scannable remote line *
Pin 12	Input 7	scannable remote line *
Pin 13	—	not used
Pin 14	0 Volt	ground
Pin 15	+5 Volt	source voltage
Pin 16	Output 5	activate 2 (external pump in SMPL) **
Pin 17	Output 3	End of determination EOD
Pin 18	Output 1	upper limit reached (in SMPL)
Pin 19	—	not used
Pin 20	—	not used
Pin 21	Input 0	Start
Pin 22	Input 2	Hold / Cont.
Pin 23	Input 4	Sample Changer ready (in SMPL) *
Pin 24	Input 6	scannable remote line
Pin 25	0 Volt	ground

* Input 3...7 are scannable in SEQ mode (SCAN command)


** Output 4...7 may be freely set in SEQ mode (CONTROL command)

Use of the Remote socket


Particular instrument functions can be triggered by activating the individual pins (connection pins) of a remote socket. The instrument itself activates certain pins according to the instrument's condition and transmits signal pulses when particular events occur.

The functions which are used for remote control in an automated system are allocated to logical remote lines. They are subdivided into input lines (signal inputs to trigger functions, e. g. START or STOP) and output lines (signal outputs to display the instrument's condition or events, e. g. 'ready' or 'End of Determination' =EOD).

The pin assignment of a 25-pin remote plug can be seen from the information given below.



Remote plug (with 25 pins)



Electrical switching of the inputs and outputs of the remote socket

Standard allocation of the remote lines in Metrohm instruments:

Remote lines	Pin	Function
Input 0	21	Start
Input 1	9	Stop
Input 2	22	Enter
Input 3	10	Clear
Input 4	23	Sample Ready
Input 5	11	not used
Input 6	24	not used
Input 7	12	not used
Output 0	5	Ready
Output 1	18	Conditioning ok
Output 2	4	Determination busy
Output 3	17	EOD (End of Determination)
Output 4	3	L1
Output 5	16	Error
Output 6	1	Activate L3
Output 7	2	Pulse for recorder

Output lines 8...13 are mostly not used as standard.

Remark: Standard pin assignments are only activated in an instrument's single mode (not in sequence functions).


The standard Remote lines

Assignment of the remote socket according to pins:

Pin	I/O lead	Function
1	Output 6	L3 activate
2	Output 7	Pulse for recorder
3	Output 4	L1
4	Output 2	Determination busy
5	Output 0	Ready
6	Output 8	not used
7	Output 9	not used
8	Output 10	not used
9	Input 1	Stop
10	Input 3	Clear
11	Input 5	not used
12	Input 7	not used
13	Output 11	not used
14	0 Volt	
15	+5 Volt	
16	Output 5	Error
17	Output 3	EOD (End of Determination)
18	Output 1	Conditioning ok
19	Output 12	not used
20	Output 13	not used
21	Input 0	Start
22	Input 2	Enter
23	Input 4	Sample Ready
24	Input 6	not used
25	0 Volt	

The pin assignment or the allocation of functions of the remote lines may differ slightly from one Metrohm instrument to another. Consult the 'Instructions for use' of your instrument.

Wiring of the Metrohm Standard Remote cable e.g. 6.2141.020 (25-pin)




Output 0 (Pin 5)	(Pin 21) Input 0
Output 1 (Pin 18)	(Pin 9) Input 1
Output 2 (Pin 4)	(Pin 22) Input 2
Output 3 (Pin 17)	(Pin 10) Input 3
Output 4 (Pin 3)	(Pin 23) Input 4
Output 5 (Pin 16)	(Pin 11) Input 5
Output 6 (Pin 1)	(Pin 24) Input 6
Output 7 (Pin 2)	(Pin 12) Input 7
Output 8 (Pin 6)	(Pin 6)
Output 9 (Pin 7)	(Pin 7)
Output 10 (Pin 8)	(Pin 8)
Output 11 (Pin 13)	(Pin 13)
Output 12 (Pin 19)	(Pin 19)
Output 13 (Pin 20)	(Pin 20)
Input 0 (Pin 21)	(Pin 5) Output 0
Input 1 (Pin 9)	(Pin 18) Output 1
Input 2 (Pin 22)	(Pin 4) Output 2
Input 3 (Pin 10)	(Pin 17) Output 3
Input 4 (Pin 23)	(Pin 3) Output 4
Input 5 (Pin 11)	(Pin 16) Output 5
Input 6 (Pin 24)	(Pin 1) Output 6
Input 7 (Pin 12)	(Pin 2) Output 7
(Pin 6)	(Pin 6) Output 8
(Pin 7)	(Pin 7) Output 9
(Pin 8)	(Pin 8) Output 10
(Pin 13)	(Pin 13) Output 11
(Pin 19)	(Pin 19) Output 12
(Pin 20)	(Pin 20) Output 13
0 V (Pin 14)	0 V (Pin 14)
0 V (Pin 25)	0 V (Pin 25)
+5 Volt (Pin 15)	+5 Volt (Pin 15)




4.3 Liquino Remote lines during a determination

XDOS Mode




The output lines 1 and 2 are used with the temperature monitoring function. Activation of Output 1 or Output 2 indicates a violation of the temperature limits.

PIP Mode




The remote lines Output 7 and Input 6 are used for remote communication of the pipetting set. They can also be used in automated systems.

CONT Mode



SMPL Mode




Assignments of the Output lines:

- Output 4 (Activate 1) : use for an aspirating pump (emptying cell)
- Output 5 (Activate 2) : use for an dispensing pump (rinsing cell)
- Output 6 (External start) : use for the start of a titration or measuring instrument
- Output 7 (Sample changer) : may be used for remote control of a sample changer

- Input 3 (External ready) : use for the detection of the EOD pulse of the titration or measuring instrument
- Input 4 (Sample changer ready) : the ready signal of a 730 Sample changer has to be programmed in the Sample changer method (use Output 4).

SEQ Mode



5 Technical specifications

Short description	Compact, multi functional dosing and control instrument for analytical laboratory use and synthesis. Connection of 4 dosing drives and 1 temperature sensor possible.	
Dimensions	<i>Height:</i>	208 mm
	<i>Width:</i>	96 mm
	<i>Depth:</i>	208 mm
Weight	2,36 kg (without accessories and keyboard)	
Housing	Metal case, multiple stove-enamel Upper part PUR (Polyurethane) non inflammable Fire classification UL94HB(IEC707), lacquered	
	<i>Keyboard case:</i>	Crastine (PBTB), aluminized inside
	<i>Keayboard foil:</i>	Polyester, chemical resistant
Power supply	<i>Voltage selector</i>	100/115/230/240V
	Settings with other supply voltages:	
		110 V --> 115 V
		117 V --> 115 V
		120 V --> 115 V
		220 V --> 230 V
	<i>Tolerance</i>	± 10%
	<i>Frequency</i>	50 - 60 Hz
	<i>Power input</i>	approx. 28 VA
	<i>Fuses</i>	T315 mA for 100 and 115 V T120 mA for 230 and 240 V
Ambient temperature	<i>Nominal operating range</i>	5...40°C (at 20...80 % atmospheric humidity)
	<i>Storage, transport</i>	-20...60°C at 60°C < 50 % atmospheric humidity at 50°C < 85 % atmospheric humidity at 40°C < 95 % atmospheric humidity
Working memory (RAM)	Battery backed, 2 x 32 KB x 8	
Real time clock	Battery backed	
Supporting battery	Lithium 3.6V, 0.85Ah, Order No. Q.005.0001 <i>Not to be exchanged by customer!</i>	
Memory card	<i>Type</i>	Memory Card, JEIDA 4.x / PCMCIA 2.x (68 Pins)
	<i>Memory capacity:</i>	128 KByte to 2 Mbyte (SRAM) Order No. 6.2245.010 (128 Kbyte)
External display	LCD, 2 x 24 characters, character height 5 mm	
Keyboard	Tactile touch keyboard with action point, splashwater-proof	
RS232-Interface	Socket D-subminiature 9-pin Connection for PC, printer or balance programmable for serial data communication with external devices.	
Option	Cable RS232C - IBM PC AT, D-subminiature 9-pin --> D-subminiature 25-pin, Order No. 6.2125.010	

6 Index

- \$C**3; 4
\$D3
\$G3; 4
\$H3; 4
\$Q1; 3
\$Q.H3
\$Q.N"i"3
\$Q.P3
\$R4
\$S3; 4
\$U3
- &**2; 9
&Assembly25
&Configuration19
&CurrData28
&Diagnose28
&Info21
&MemoryCard27
&Mode9
&Setup22
&UserMeth24
- 25-pin remote plug37
- abort.....3
 action3
 Activate37; 42
 active26
 ASCII characters2
 aspirating pump.....42
 automated system37
 Automatic messages23
 automatic sending of
 measured values22
 auxiliar solution14
- backwards2
 basic buret functions25
 basic settings.....19
 bit pattern.....16
 branches2; 9
 buret data21
 buret emptying.....25
 buret filling25
 buret functions25
 buret setup..... 10; 12; 13; 14
 busy37
- cable29; 32
 call up of objects2
 card label27
 change directory27
 command selection16
 compiled buret functions25
 conditioning37
 configuration19
 configuration of printed
 reports20
 configuration of the Dosing
 units19
CONT13
- CONT** Mode41
 Continue3
 control characters29
 control of a measuring
 instrument15
 control of external
 instruments15
 C_R 1; 29
 CTS29; 31; 32
 current directory27
 curve output23
- data bit29
 data block1
 data line1
 data transfer failure34
 data transfer protocol29
 decimal point2
 detail-Info3
 diagnostics28
 digits2
 DIP switches34
 directory27
 dispensing pump42
 display content21
 dosing25
 dosing options18
 dosing parameters10
 dosing test mode18
 dosing test parameters18
 drive address19; 25
 DTE29
 DTR29; 32
- E4331
 earthing32
 emptying4
 End of Determination37; 39
 EOD37
 EOD pulse42
 Epson34
 Error37
 error 36-3934
 error 4234
 error 4334
 error messages6
 events37
 exchange position25
 extended dosing10
 external pumps14
 External ready42
 External start42
- file operations24; 27
 filling4
 forwards2
- general rules1
GLP18
GLP limits20
 GND32
 Go3
 handshake29
- hardware-handshake31
 hierarchically structured1
 Highest index3
 Hold3
 hold processes3
 HWs31
 I/O leads39
- IBM34
 identifications28
 inactive4; 26
 input37; 39
 input buffer29
 instrument components25
 instrument condition37
 instrument data21
 instrument diagnosis28
 instrument functions37
 instrument setup22
 interface scan16
- key codes22
- L1 37
 levels2
 L_F 1; 29
 limits20
 locking of key functions22
 lower case2
- main branches9
 manual distribution4; 9
 manual dosing4
 measuring cell rinsing14
 measuring instrument42
 memory27
 method to be executed16
 monitoring functions10
- name3
 name of the son node3
 negative sign2
 node2
 number of son nodes3
 numerical values2
- objects1
 on-line34
 operating modes9
 output37; 39
 output data16
 output lines26
- parameters for **CONT**-Mode13
 parity29
 path3
 peripheral units19
 pin37; 39
 pin assignment32; 35
PIP12
PIP Mode41
 pipetting12

pipetting parameters.....	12	rinsing	14	status messages	4
pipetting set	41	root.....	2; 9	Stop	3
plug	33	RS 232 interface	29	stop bit.....	29
polarity	32	RS232 serial interface.....	20	stop processes.....	3
port specifications.....	12; 13; 18	RS232C Interface.....	32	stopcock.....	25
power on simulation.....	28	RTS	29; 32	stopped	4
preceding dot.....	2	RxD	32	stored methods	24
PREP function.....	25			SWchar	29
preparation of standard		sample changer.....	42	SWline.....	30
solutions	13	Sample Ready	37		
preparing	4	sample size.....	28	technical specifications	43
preselections		sampling mode.....	14	temperature	26
.....	11; 12; 13; 15; 17; 18	sampling operation.....	14	temperature limits	41
printer.....	34	sampling parameters.....	14	temperature monitoring	
properties of the RS 232		send errors.....	7	10; 41
interface.....	29	sending format.....	22	temperature read-out.....	26
protective earthing	32	SEQ	16	tests	28
pump.....	42	SEQ Mode.....	42	timer function.....	16
		SEQ mode data.....	21	titration instrument	15
Query	3	sequence.....	16	transitional range.....	32
query command.....	1	sequence delay	16	transmission	31
qUit	3	sequence mode.....	16	tree.....	1
quotes	1	sequence node.....	16	trigger	1; 3
		sequence parameters.....	16	trigger functions.....	37
RAM initialization.....	28	setting the output lines.....	26	triggers	3
rates of the dosing drive	25	settings for manual operation...	19	TxD	32
receive errors	7	signal outputs	37		
relative address.....	2	signal pulses.....	37	upper case.....	2
Remote cable.....	40	signals	16	user defined methods.....	24
remote control language.....	1	SMPL	14		
remote control tree.....	9	SMPL Mode	42	value	2
Remote interface.....	35	SMPL mode data	21	voltage.....	32
Remote interface data.....	21	socket	33	volume monitoring.....	10
remote lines	16; 37	software-handshake	29		
Remote socket	26; 37	standard allocation	37	XDOS	10
report definitions		standard Remote cable	40	XDOS Mode.....	41
.....	11; 12; 13; 15; 17; 18	start of a titration	42	XOFF	29
report header	20	start processes	3	XON	29
report output	21	status conditions.....	4		
requests	11	status information	3		