719 S Titrino

Version de programme 5.719.0020

Mode d'emploi

Table des matières

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3 Operation via RS232 Interface

3.1 General rules

The Titrino has an extensive remote control facility that allows full control of the Titrino via the RS 232 interface, i.e. the Titrino can receive data from an external controller or send data to an external controller. C_R and L_F are used as terminators for the data transfer. The Titrino sends $2xC_R$ and L_F as termination of a data block, to differentiate between a data line which has C_R and L_F as terminators. The controller terminates its commands with C_R and L_F . If more than one command per line is sent by the controller, “;” is used as a separator between the individual commands.

The data are grouped logically and easy to understand. Thus e.g., for the selection of the dialog language, the following must be sent

&Config.Aux.Language "english"

whereby it is sufficient to only transmit the boldface characters, thus:

&C.A.L "english"

The quantities of the commands above are:

Config	configuration data
Aux	auxiliaries, various data
Language	setting the dialog language

The data are hierarchically structured (tree form). The quantities that occur in this tree are called **objects** in the following. The dialog language is an object which can be called up with the

&Config.Aux.Language

command.

If one is in the desired location in the tree, the value of the object can be queried.

&Config.Aux.Language \$Q Q means Query


The query command \$Q initiates the issuing of the value on the instrument and the value emission is triggered. Entries which start with \$, trigger something. They are thus called **triggers**.

Values of objects can not only be queried, they can also be modified. Values are always entered in quotes, for example:

&Config.Aux.Language "english"

3.1.1 Call up of objects

An excerpt from the object tree is represented below:

Rules	Example
The root of the tree is designated by &.	
The branches (levels) of a tree are marked with a dot (.) when calling up an object.	
When calling up an object, it is sufficient to give only as many letters as necessary to uniquely assign the object. If the call is not unequivocal, the first object in the series will be recognized.	Calling up the dialog language &Config.Aux.Language or &C.A.L
Upper- or lowercase letters may be used.	&C.A.L or &c.a.l
An object can be assigned a value. Values are signified at the beginning and end by quotes ("). They may contain up to 24 ASCII characters. Numerical values can contain up to 6 digits, a negative sign, and a decimal point. Numbers with more than 6 characters are not accepted; more than 4 decimal places are rounded off. For numbers <1, it is necessary to enter leading zeros.	Entering the dialog language: &C.A.L"english" correct entry of numbers: "0.1" incorrect entry of numbers "1,5" or "+3" or ".1"
The current object remains until a new object is called.	entry of another dialog language: "deutsch"
New objects can be addressed relative to the old object: A preceding dot leads forwards to the next level in the tree.	From the root to node 'Aux': &C.A Forward from node 'Aux' to 'Prog': .P
More than one preceding dot leads one level backwards in the tree. n node backwards require n+1 preceding dots.	Jump from node 'Prog' to node 'Aux' and select a new object 'Language' at this level: ..L
If you must jump back to the root, enter a preceding &.	Change from node 'Language' via the root to node 'Mode': &M

3.1.2 Triggers

Triggers initiate an action on the Titrino, for example, starting a process or sending data. Triggers are marked by the introductory symbol \$.

The following triggers are possible:

\$G	Go	Starts processes, for ex. starting the mode run or setting the RS 232 interface parameters
\$S	Stop	Stops processes
\$H	Hold	Holds processes
\$C	Continue	Continues processes after Hold
\$Q	Query	Queries all information from the current node in the tree forward up to and including the values
\$Q.P	Path	Queries the path from the root of the tree up to the current node
\$Q.H	Highest Index	Queries the number of son nodes of the current node
\$Q.N"i"	Name	Queries the name of the son node with index i, $i = 1 - n$
\$D	Detail-Info	Queries the detailed status information
\$U	qUit	Aborts the data flow of the instrument, for example, after \$Q

The triggers \$G and \$S are linked to particular objects, see the summary table page 63ff.

All other triggers can be used at any time and at all locations on the object tree.

Examples:

Querying the value of the baud rate: **&Config.RSSet.Baud \$Q**

Querying all values of the node "RSSet": **&Config.RSSet \$Q**

Querying the path of the node "RSSet": **&Config.RSSet \$Q.P**

Start mode: **&Mode \$G**

Querying the detailed status: **\$D**

3.1.3 Status messages

In order to have an efficient control by an external control device, it must also be possible to query status conditions; they provide information on the status of the Titrino. The trigger \$D initiates output of the status. Status messages consist of the global status, the detailed status and eventual error messages, e.g. \$S.Mode.SET;E26. The global status informs on the activity of the process, while the detailed status conditions show the exact activity within the process.

The following **global status conditions** are possible:

\$G	Go:	The Titrino is executing the last command.
\$H	Hold:	The Titrino has been held (\$H, key <meas/hold> or by an error which effects the hold status)
\$C	Continue:	The Titrino has been restarted actively after hold
\$R	Ready:	The Titrino has executed the last command and is ready
\$S	Stop:	A process has been aborted in an "unnatural manner". e.g. stopped or aborted because there was an error.

Detailed status conditions

Status conditions of the global \$G:

\$G	.Mode.SET	.Inac:	Instrument at the beginning or at the end of a titration.
		.Req.Id1:	Instrument in the SET mode, requesting Id1 after start.
		.Id2:	Instrument in the SET mode, requesting Id2 after start.
		.Id3:	Instrument in the SET mode, requesting Id3 after start.
		.Smp1:	Instrument in the SET mode, requesting sample size after start.
		.Unit:	Instrument in the SET mode, requesting unit of sample size after start.
		.Start:	Instrument in the SET mode, processing the start conditions.
		.SET1:	Instrument in the SET mode, titrating to the first endpoint.
		.SET2:	Instrument in the SET mode, titrating to the second endpoint.
		.Cond.Ok:	Instrument in the SET, conditioning, endpoint reached (after the first startup from the standby mode).
		.Cond.Prog:	Instrument in the SET mode, conditioning, endpoint not reached (Conditioning progressing).
\$G	.Mode.MEAS	.Inac:	Instrument at the beginning or at the end of a titration.
		.Req.Id1:	Instrument in the MEAS mode, requesting Id1 after start.
		.Id2:	Instrument in the MEAS mode, requesting Id2 after start.
		.Id3:	Instrument in the MEAS mode, requesting Id3 after start.
		.Smp1:	Instrument in the MEAS mode, requesting sample size after start.
		.Unit:	Instrument in the MEAS mode, requesting unit of sample size after start.
		.Meas:	Instrument in the MEAS mode, measuring.
\$G	.Mode.CAL	.Inac:	Instrument at the beginning or at the end of a calibration
		.Req.Temp:	Instrument in the CAL mode, requesting calibration temperature.
		.Meas.Temp:	Instrument in the CAL mode, measuring calibration temperature.
		.Req.Buf1:	Instrument in the CAL mode, requesting pH of buffer 1.
		.Meas.Buf1:	Instrument in the CAL mode, measures buffer 1.
		.Req.Buf2:	Instrument in the CAL mode, requesting pH of buffer 2.
		.Meas.Buf2:	Instrument in the CAL mode, measures buffer 2.
		etc.	
\$G	.Assembly.Bur	.Fill:	Buret in filling process
		.ModeDis:	Buret in DIS mode

Status conditions of the global \$H:

The status message of the action which has been held appears.
If the process is held because a monitored limit has been violated, its status message is \$H.Mode.XXX.XXX.

Status conditions of the global \$C:

The status conditions of the global \$C are identical with the ones of the global status \$G. They appear when the process has been restarted actively from the status "Hold" (\$C, key <meas/hold> or automatically after elimination of an error).

Status conditions of the global \$R:

\$R .Mode.XXXX.QuickMeas: Quick manual measurement from the initial status in mode XXXX.

\$R .Mode.SET .Inac: Instrument in the SET mode, inactive.
.Cond.Ok: Instrument in the SET mode, conditioning, endpoint reached.
.Cond.Prog: Instrument in the SET mode, conditioning, endpoint not reached.
\$R .Mode.MEAS .Inac: Instrument in the MEAS mode, inactive.
\$R .Mode.CAL .Inac: Instrument in the CAL mode, inactive.
\$R .Assembly.Bur.ModeDis: Buret in the DIS mode, inactive.

Status conditions of the global \$\$:

\$\$.Mode.XXXX.QuickMeas: Quick manual measurement from the initial status in mode XXXX.

The instrument gives the status from which it has been stopped. The detailed status information is therefore identical to for the global status \$G. Violation of monitored limits with action "end" give the status message **\$\$.Mode.XXX.Inac;EYYY**.

3.1.4 Error messages

Error messages are added to the status messages and separated from them by the sign ";".

E20	Check exchange unit. Exit: Mount Exchange Unit (properly) or &m \$\$.
E21	Check electrode, short circuit. Exit: Rectify fault or &m \$\$.
E22	Check electrode, break. Exit: Rectify fault or &m \$\$.
E23	Division by zero. Exit: The error message disappears on next startup or on recalculation.
E26	Manual stop. Exit: The error message disappears on next startup.
E27	Stop V reached in SET. Exit: The error message disappears on next startup.
E28	Wrong object call up Exit: Send correct path for object. Start path at root.
E29	Wrong value or no value allowed. Exit: Send correct value or call up new object.
E30	Wrong trigger, this trigger is not allowed or carrying-out of action not possible. Exit: Send correct trigger (exception: \$D) or call up new object.
E31	Command is not possible in active status. Repeat command in inactive status. Exit: Send new command.
E32	Command is not possible during titration. Repeat command during the conditioning phase or in inactive status. Exit: Send new command.
E33	Value has been corrected automatically. Exit: Send new command.
E34	Instrument at the end of the titration and sample data is edited; the instrument at rest or editing during filling. Exit: &m \$\$.

RS receive errors:

- E36** Parity
Exit: <QUIT> and ensure settings of appropriate parameters at both devices are the same.
- E37** Stop Bit
Exit: <QUIT> and ensure settings of appropriate parameters at both devices are the same.
- E38** Overrun error. At least 1 character could not be read.
Exit: <QUIT>
- E39** The internal working-off buffer of the Titrino is full (>82 characters).
Exit: <QUIT>

RS send errors:

- E40** DSR=OFF No proper handshake for more than 1 s.
Exit: <QUIT> Is the receiver switched on and ready to receive?
- E41** DCD=ON No proper handshake for more than 1 s.
Exit: <QUIT> Is the receiver switched on and ready to receive?
- E42** CTS=OFF No proper handshake for more than 1 s.
Exit: <QUIT> Is the receiver switched on and ready to receive?
- E43** The transmission of the Titrino has been interrupted with XOFF for at least 6 s.
Exit: Send XON or <QUIT>
- E44** The RS interface parameters are no longer the same for both devices. Reset.
- E45** The receive buffer of the Titrino contains an incomplete command (L_f missing). Sending from the Titrino is therefore blocked.
Exit: Send L_f or <QUIT>.

E120	Overrange of the primary measured value (pH, U, Ipol, Upol or T with MEAS T). The secondary measured value (temperature) may be instable as well. Exit: Correct error or &m \$\$.
E123	Missing EP for calculation. Exit: The error message disappears on next startup or on recalculation.
E128	No new mean. Exit: The error message disappears on next startup or on recalculation.
E129	No new common variable, old value remains. Exit: The error message disappears on next startup or on recalculation.
E130	Wrong sample. For SET with preset titration direction the first measured value lies behind the endpoint. Exit: The error message disappears on next startup.
E131	No EP set for SET. Exit: The error message disappears on next startup.
E132	Silo empty and it has been started with open silo or empty silo has been opened. Exit: Send a silo entry.
E133	Silo full. Exit: Send new command.
E134	No method. A method, which is required from the silo memory, does not exist. Exit: The error message disappears on next startup.
E135	Check temp.sensor in MEAS T. Exit: Correct error or &m \$\$.
E136	Same buffer in CAL. Measured value of the second buffer differs less than 6 mV from the measured value of the first buffer. Exit: Correct error or &m \$\$.
E137	XXX Bytes are missing so that the method, the silo line could not be stored. Exit: Send new command.

3.2 Remote control commands

3.2.1 Overview

The internal object tree can be divided into the following branches:

&	Root
Mode	Method parameters
UserMeth	Administration of the internal user-memory for methods
Config	Instrument configuration
SmplData	Sample specific data
Info	Current Data
Assembly	Component data
Setup	Setting the operating mode
Diagnose	Diagnostics program

&Mode

Object	Description	Input range	Reference
& Root			
Mode	Mode	\$G, \$\$, \$H, \$C	3.2.2.1.
.QuickMeas	Rapid meas. in basic mode	\$G, \$\$	3.2.2.2.
.Select	Mode selection	SET , MEAS,CAL	3.2.2.3.
.SETQuantity	Measured quantity for SET	pH , U, Ipol, Upol	ditto
.MEASQuantity	Measured quantity for MEAS	pH , U, Ipol, Upol, T	ditto
.Name	Name of current method	read only/read+write	3.2.2.4.
.Parameter*	Parameter of current mode, page 65ff		
.Def	Definitions for data output		
.Formulas	Calculation formulas		
.1	for result 1		
.Formula	Calculation formula	special	3.2.2.5.
.TextRS	Text for result output	up to 8 ASCII char	ditto
.Decimal	Number of decimal places	0... 2 ...5	ditto
.Unit	Unit for result output	up to 6 ASCII char	ditto
:	up to 9 results		
.ComVar	Assignment of common variables		
.C30	for C30	RSX,EPX,CXX,MNX	3.2.2.6.
up to C39			
.Report	Reports at the end of determination		
.Assign	Assignment	depends on mode	3.2.2.7.
.Mean	Assignment for mean calculation		
.1	MN1		
.Assign	Input of variable	RSX, EPX, CXX	3.2.2.8.
up to 9			
.CFmla	Calculation constants		
.1	Calculation constant C01		
.Value	Input of value	0 ...±999 999	3.2.2.9.
up to C19			

*Parameter	Tree part "Parameters for SET"		
.SET1	Control parameters for EP1		
.EP	Endpoint 1	depends on meas.quant.	3.2.2.10.
.UnitEp	Unit of endpoint	read only	ditto
.Dyn	Dynamics	depends on meas.quant.	3.2.2.11.
.UnitDyn	Unit of dynamics	read only	ditto
.MaxRate	Maximum dosing rate	0.01... 10 ...150, max.	ditto
.MinRate	Minimum dosing rate	0.01... 25.0 ...9999	ditto
.Stop	Titration stop		
.Type	Type of stop criterion	drift , time	3.2.2.12.
.Drift	Stop drift	1... 20 ...999	ditto
.Time	Switch-off delay time	0... 10 ...999, inf	ditto
.StopT	Stop time	0...999 999, OFF	ditto
.SET2	Control parameters for EP2, as for EP1		
.TitrPara	Titration parameters		
.Direction	Titration direction	+, -, auto	3.2.2.13.
.XPause	Pause 1	0 ...999 999	3.2.2.14.
.Start V	Start volume		
.Type	Type of start volume	abs., rel., OFF	3.2.2.15.
.V	Volume for absolute start volume	0 ...999.99	ditto
.Factor	Factor for relative start volume	0 ...±999 999	ditto
.Rate	Dispensing rate for start volume	0.01...150.0, max.	ditto
.Pause	Pause 2	0 ...999 999	3.2.2.16.
.ExtrT	Extraction time	0 ...999 999	3.2.2.17.
.MeasInput	Measuring Input	1 , 2, diff.	3.2.2.18.
.Ipol	Polarization current	0... 1 ...±127	ditto
.Upol	Polarization voltage	0... 400 ...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON, OFF	ditto
.Temp	Titration temperature	-170.0... 25.0 ...500.0	3.2.2.19.
.StopCond	Stop conditions		
.VStop	Stop volume		
.Type	Type of stop volume	abs. , rel., OFF	3.2.2.20.
.V	Volume for absolute stop volume	0... 99.99 ...9999.99	ditto
.Factor	Factor for relative stop volume	0...± 999 999	ditto
.FillRate	Filling rate	0.01...150.0, max.	3.2.2.21.
.Statistics	Statistics		
.Status	Status of statistics calculation	ON, OFF	3.2.2.22.
.MeanN	No. of individual determinations	2 ...20	ditto
.ResTab	Result table		
.Select		original , delete n, delete all	ditto
.DelN	Deletion of individual results	1 ...20	ditto
.Presel	Preselections		
.Cond	Conditioning	ON, OFF	3.2.2.23.
.DriftDisp	Display of drift during cond.	ON , OFF	ditto
.DCor	Drift Correction		
.Type	Drift type	auto, man., OFF	3.2.2.24.
.Value	Value for manual drift correction	0.0 ...99.9	ditto
.IReq	Request of Id's after start	id1, id1&2, all, OFF	3.2.2.25.
.SReq	Request of smpl size after start	value, unit, all, OFF	ditto
.ActPulse	Output of a pulse	first, all, cond., OFF	3.2.2.26.

*Parameter	Tree part "Parameters for MEAS"		
Measuring parameters			
.Measuring	Drift for meas.value acquisition	depends on meas.quant.	3.2.2.27.
.SignalDrift	Unit of measured value drift	read only	ditto
.UnitSigDrift	Equilibrium time	0...9999, OFF	ditto
.EquTime	Measuring input	1, 2, diff.	3.2.2.28.
.MeasInput	Polarization current	0...1...±127	ditto
.Ipol	Polarization voltage	0...400...±1270	ditto
.Upol	Test for polarized electrodes	ON, OFF	ditto
.PolElectrTest	Titration temperature	-170.0...25.0...500.0	3.2.2.29.
.Temp			
Statistics			
.Statistics	Status of statistics calculation	ON, OFF	3.2.2.22.
.Status	No. of individual determinations	2...20	ditto
.MeanN	Result table		
.ResTab		original ,delete n,delete all	ditto
.Select	Deletion of individual results	1...20	ditto
.DeIN			
Preselections			
.Presel	Request of Id's after start	id1, id1&2, all, OFF	3.2.2.24.
.IReq	Request of sample size after start	value, unit, all, OFF	ditto
.SReq	Output of a pulse	ON, OFF	3.2.2.25.
.ActPulse			

*Parameter	Tree part "Parameters for CAL"		
Calibration parameters			
.Calibration			
.MeasInput	Measuring input	1, 2, diff.	3.2.2.30.
.CalTemp	Calibration temperature	-20.0... 25.0 ...120.0	3.2.2.31.
.Buffer			
.1			
.Value	pH value of buffer 1	0... 7.00 ...±20.00	3.2.2.32.
.2			
.Value	pH value of buffer 2	0... 4.00 ...±20.00, OFF	ditto
.	up to 9 buffers		
.SignalDrift	Drift for meas.value acquisition	0.5... 4.00 ...999, OFF	3.2.2.33.
.EquTime	Equilibrium time	0... 110 ...9999, OFF	ditto
.ElectrodeId	Electrode identification	8 ASCII char.	3.2.2.34.
.SmplChanger	Calibration on a Titrimo	ON, OFF	3.2.2.35.
.ActPulse	Output of a pulse	first, all, OFF	3.2.2.36.
Statistics			
.Statistics			
.Status	Status of statistics calculation	ON, OFF	3.2.2.22.
.MeanN	No. of individual determinations	2 ...20	ditto
.ResTab	Result table		
.Select		original , delete n, delete all	ditto
.DeIN	Deletion of individual results	1 ...20	ditto

&UserMeth

Object	Description	Input range	Reference
& Root			
<pre> : : : - UserMeth : : </pre>			
Method memory			
.FreeMemory	Memory available	read only	3.2.2.37.
.Recall	Load method	\$G	3.2.2.38.
.Name	Method name	8 ASCII characters	ditto
.Store	Save method	\$G	ditto
.Name	Method name	8 ASCII characters	ditto
.Delete	Delete method	\$G	ditto
.Name	Method name	8 ASCII characters	ditto
.DelAll	Delete all methods	\$G	ditto
.List	List of methods		
.1	Method 1		
.Name	Method name	read only	3.2.2.39.
.Mode	Mode	read only	ditto
.Quantity	Measured quantity	read only	ditto
.Bytes	Method size in bytes	read only	ditto
.Checksum	Checksum of method	read only	ditto
.2	for each method		

&Config

Object	Description	Input range	Reference
& Root			
Config	Instrument configuration		
.PeriphUnit	Selection of peripheral units		
.CharSet	External printer	Epson, Seiko, Citizen IBM , HP	3.2.2.40.
.Balance	Selection of balance	Sartorius , Mettler, Mettler AT AND, Precisa	3.2.2.41.
.Plot	Selection of plot at analog output	U , dU/dt, V, dV/dt U(rel), T	ditto
.Aux	Miscellaneous		
.Language	Dialog language	english , deutsch, francais, espanol, italiano, portugese, svenska	3.2.2.42.
.Set	Setting of date and time	\$G	3.2.2.43.
.Date	Date	XXXX-XX-XX	
.Time	Time	XX:XX	
.RunNo	Run number	0 ...9999	3.2.2.44.
.AutoStart	Automatic start	1...9999, OFF	3.2.2.45.
.StartDelay	Start delay time	0 ...999 999	3.2.2.46.
.DevName	Device label	8 ASCII char.	3.2.2.47.
.Prog	Program version	read only	3.2.2.48.
.RSSet	Settings RS232	\$G	3.2.2.49.
.Baud	Baud rate	300,600,1200,2400,4800, 9600	
.DataBit	Number of data bits	7, 8	ditto
.StopBit	Number of stop bits	1 , 2	ditto
.Parity	Parity	even, odd, none	ditto
.Handsh	Handshake	HWs , HWf, SWchar, SWline, none	ditto
.ComVar	Values of common variables		
.C30	C30	0 ... ±999 999	3.2.2.50.
up to C39	0 ... ±999 999		

&SmplData

Object	Description	Input range	Reference
& Root			
SmplData	Sample data		
.Status	Status of silo memory	ON, OFF	3.2.2.51.
.OFFSilo	Current sample data		
.Id1	Sample identification 1	up to 8 ASCII char	3.2.2.52.
.Id2	Sample identification 2	up to 8 ASCII char	ditto
.Id3	Sample identification 3	up to 8 ASCII char	ditto
.ValSmpl	Sample size	±X.XXXXX	ditto
.UnitSmpl	Unit of sample size	up to 5 ASCII char	ditto
.ONSilo	Current sample data		
.Counter	Counter of silo memory		
.MaxLines	Maximum lines	read only	3.2.2.53.
.FirstLine	First line	read only	ditto
.LastLine	Last line	read only	ditto
.EditLine	Editing silo lines		
.1	1 st silo line		
.Method	Method name	up to 8 ASCII char	3.2.2.54.
.Id1	Sample identification 1	up to 8 ASCII char	ditto
.Id2	Sample identification 2	up to 8 ASCII char	ditto
.Id3	Sample identification 3	up to 8 ASCII char	ditto
.ValSmpl	Sample size	±X.XXXXX	ditto
.UnitSmpl	Unit of sample size	up to 5 ASCII char	ditto
up to 99 lines			
.DelLine	Delete silo line	\$G	3.2.2.55.
.LineNum	Line number	1...99, OFF	ditto
.DelAll	Delete silo line	\$G	3.2.2.56.
.CycleLines	Cycle lines	ON, OFF	3.2.2.57.

&Info

Object	Description	Input range	Reference
& Root			
- Info - Report - Select - Current data - Transmission of formatted reports - Report type \$G - configuration, parameters, smpl data, statistics, silo, calib C-fmla, def, user method, full, short, calc, all, ff			
CalibrationData	pH calibration data	\$G	3.2.2.58.
- Inp1 - pHas - Asymmetry pH - Slope - Slope of electrode - Temp - Calibration temperature - Date - Date of calibration - ElectrodeId - Id of the calibrated electrode - Inp2 - For measuring input 2, as for input 1 - Diff - For differential input, as for input 1	- For measuring input 1 0...7.00...±20.00 0...1.000...±9.999 -170.0...25.0...500.0 - read only - read only	- ditto - ditto - ditto - ditto - ditto	3.2.2.59.
Checksums	Checksums	\$G	3.2.2.60.
- ActualMethod - Checksum of current method	read only	ditto	
DetermData	Determination data	\$G	3.2.2.61.
- Write - Read/write for several nodes - ExV - Volume of Exchange/Dosing unit	- ON, OFF - read only/read + write	ditto	
TitrationResults	Titration results		
- RS - Calculated results 1st result - Value - Value - EP - Endpoint 1st result - Value - Value - Meas - Measured value - Var - Variables C4X - C40 - Start measured value - C41 - Titration end volume - C42 - Titration time - C43 - Volume drift in SET - C44 - Titration temperature - C45 - Start volume - C46 - Asymmetry pH - C47 - Slope of electrode - DTime - Time for drift correction	- up to 9 results - read only - read only - read only - read only/read + write - read only/read + write - read only/read + write - read only/read + write - read only/read + write - read only - read only - read only/read + write	3.2.2.62. - ditto - ditto - ditto	
StatisticsVal	Statistics values		
- ActN - Number of results in chart 1 1st mean - Mean - Mean - Std - Absolute standard deviation - RelStd - Relative standard deviation	- read only - read only - read only	- ditto - ditto - ditto	3.2.2.63.

└ up to 9 mean values			
.ActualInfo	Current data		
└ .Inputs	I/O Inputs		
└ .Status	Line status	read only	3.2.2.64.
└ .Change	Change of line status	read only	ditto
└ .Clear	Clear change	\$G	ditto
└ .Outputs	as for I/O Inputs		ditto
└ .Status	Line status	read only	ditto
└ .Change	Change of line status	read only	ditto
└ .Clear	Clear change	\$G	ditto
└ .Assembly	From Assembly		
└ .CyclNo	Cycle number	read only	3.2.2.65.
└ .Counter	Assembly counter	read only	3.2.2.66.
└ .V	Volume counter	read only	ditto
└ .Clear	Clears counter	\$G	ditto
└ .PistonPos	Measured value	read only	3.2.2.67.
└ .Meas	Measured value	read only	ditto
└ .Titrator	From Titrator		
└ .CyclNo	Cycle number	read only	3.2.2.68.
└ .V	Volume	read only	ditto
└ .Meas	Measured indicator voltage	read only	ditto
└ .dVdt	Volume drift dV/dt	read only	ditto
└ .dMeasdt	Measured value drift	read only	ditto
└ .dMeasdV	1st deviation of titration curve	read only	ditto
└ .EP	EP entry		
└ .Index	Index of entry	read only	3.2.2. 69.
└ .X	X coordinate	read only	ditto
└ .Y	Y coordinate	read only	ditto
└ .Display	Display		
└ .L1	Text line 1	up to 24 ASCII char	3.2.2.70.
└ .L2	Text line 2	up to 24 ASCII char	ditto
└ .Assembly	Assembly		
└ CycleTime	Cycle time	read only	3.2.2.71.
└ ExV	Volume of Exchange unit	read only	ditto

&Assembly

Object	Description	Input range	Reference
& Root			
- Assembly - Bur - Rates - Forward - .Select - .Digital - Reverse - .Select - .Digital - Fill - ModeDis - .Select - .V - .Time - .VStop - .AutoFill - Meas - .Status - .MeasInput - .Ipol - .Upol - Outputs - .AutoEOD - .SetLines - .LO - up to L 3 - .ResetLines	Assembly control Buret Rates Forward rate Type of rate control Digital rate as for forward rate Type of rate control Digital rate Fill Dispensing Type of dispensing control Volume to be dispensed Time to dispense Limit volume Filling after each increment Measuring Measuring ON/OFF Selection of measuring input Polarization current Polarization voltage I/O outputs Automatic output of EOD Set I/O lines Signal on LO Reset I/O lines	digital, analog 0...150, max. digital, analog 0...150, max. \$G,\$H,\$C \$G,\$S,\$H,\$C volume , time 0.0001... 0.1 ...9999 0.25... 1 ...86 400 0.0001...9999, OFF ON, OFF ON, OFF 1 , 2, diff., Ipol, Upol, Temp 0... 1 ...±127 0... 400 ...±1270 ON , OFF \$G active,inactive,pulse, OFF \$G	3.2.2.72. ditto ditto ditto 3.2.2.73. 3.2.2.74. ditto ditto ditto ditto ditto 3.2.2.75. ditto ditto ditto 3.2.2.76. ditto ditto ditto

&Setup

Object	Description	Input range	Reference
& Root			
Setup	Settings for the operating mode		
.Keycode	Send key code	ON, OFF	3.2.2.77.
.Tree	Sending format of path info		
.Short	Short format of path	ON, OFF	3.2.2.78.
.ChangedOnly	Paths of modified nodes only	ON, OFF	ditto
.Trace	Message on changed values	ON, OFF	3.2.2.79.
.Lock	Lock key functions		
.Keyboard	Lock all keyboard keys	ON, OFF	3.2.2.80.
.Config	Lock <CONFIG> key	ON, OFF	ditto
.Parameter	Lock <PARAM> key	ON, OFF	ditto
.SmplData	Lock <SMPL DATA> key	ON, OFF	ditto
.UserMeth	Lock functions		
.Recall	Lock "loading"	ON, OFF	ditto
.Store	Lock "saving"	ON, OFF	ditto
.Delete	Lock "deletion"	ON, OFF	ditto
.Display	Lock display function	ON, OFF	ditto
.Mode	Setting waiting intervals		
.StartWait	Waiting time after start	ON, OFF	3.2.2.81.
.FinWait	Waiting time after run	ON, OFF	ditto
.SendMeas	Automatic sending of measured values		
.SendStatus	Connect/disconnect sending	ON, OFF	3.2.2.82.
.Interval	Time interval	0.08...4...16200,	ditto
.Select	Selection	Assembly, Titration	3.2.2.83.
.Assembly	From assembly		
.CyclNo	Cycle number	ON, OFF	3.2.2.84.
.V	Volume	ON, OFF	ditto
.Meas	Measured indicator voltage	ON, OFF	ditto
.Titration	From Titration		
.CyclNo	Cycle number	ON, OFF	3.2.2.85.
.V	Volume	ON, OFF	ditto
.Meas	Measured indicator voltage	ON, OFF	ditto
.dVdt	Volume drift dV/dt	ON, OFF	ditto
.dMeasdt	Measured value drift	ON, OFF	ditto
.dMeasdV	1st deviation of titration curve	ON, OFF	ditto
.AutoInfo	Automatic message for changes		3.2.2.86.
.Status	Switch AutoInfo on/off	ON, OFF	ditto
.P	When mains is switched on	ON, OFF	ditto
.T	Titration infos		
.R	When "ready"	ON, OFF	ditto
.G	When method started	ON, OFF	ditto
.GC	When start is initiated	ON, OFF	ditto
.S	When stopped	ON, OFF	ditto
.B	Begin of method	ON, OFF	ditto
.F	End of process	ON, OFF	ditto
.E	Error	ON, OFF	ditto
.H	When "hold"	ON, OFF	ditto

-.C	Continue after "hold"	ON, OFF	ditto
-.O	Conditioning OK	ON, OFF	ditto
-.N	Conditioning not OK	ON, OFF	ditto
-.Re	Request after start	ON, OFF	ditto
-.Si	Silo empty	ON, OFF	ditto
-.EP	Entry in EP list	ON, OFF	ditto
-.RC	Recalculation of results done	ON, OFF	ditto
-.I	Changing an I/O input	ON, OFF	ditto
-.O	Changing an I/O output	ON, OFF	ditto
- .PowerOn	RESET (power on)	\$G	3.2.2.87.
- .Initialise	Set default values	\$G	3.2.2.88.
- .Select	Selection of branch	ActMeth , Config, Silo, Calib Assembly, Setup, All	ditto
- .RamInit	Initialization of working mem.	\$G	3.2.2.89.
- .InstrNo	Device Identification	\$G	3.2.2.90.
- .Value	Input of device identification	8 ASCII characters	ditto

&Diagnose

Object	Description	Input range	Reference
& Root			
. - Diagnose - .Report	Diagnose		
	Output of adjustment parameters	\$G	3.2.2.91.

3.2.2 Description of the remote control commands

3.2.2.1. Mode \$G, \$S, \$H, \$C

Start and stop (\$G, \$S) or hold of the current method (3.2.2.3) with \$H and continue with \$C.

\$G also serves to continue after inquiries of identifications and sample size after the start (see 3.2.2.25) as well as after inquiries of calibration temperature and pH values of buffers (see 3.2.2.31 and 3.2.2.32).

3.2.2.2. Mode.QuickMeas \$G, \$S

Start and stop of a measurement in the basic mode with the parameters (measured quantity, measuring input) of the current method. Corresponds to the <meas/hold> key.

With an ongoing measurement, the current mode can be started. This stops the measurement automatically.

3.2.2.3. Mode.Select SET, MEAS, CAL Mode.SETQuantity pH, U, Ipol, Upol Mode.MEASQuantity pH, U, Ipol, Upol, T

Selection of the standard mode. Mode and the measured quantity belong to the complete selection.

If a method is selected from the method memory, the nodes &Mode.Select and &Mode.XXXQuantity are overwritten with mode and measured quantity of the corresponding user method.

3.2.2.4. Mode.Name read only

Name of the current method in the working memory. \$Q sends 8 ASCII characters. Standard methods carry the name *****. The node can be set read + write, see 3.2.2.61.

3.2.2.5. Mode.Def.Formulas.1.Formula EPX, CXX, RSX, +, -, *, /, (,) Mode.Def.Formulas.1.TextRS up to 8 ASCII characters Mode.Def.Formulas.1.Decimal 0...2...5 Mode.Def.Formulas.1.Unit up to 6 ASCII characters Mode.Def.Formulas.2.Formula etc. up to .9

Entry of formulas. Rules for formula entry, see page 39f.

Example: "(EP2-EP1)*C01/C00"

In addition to the formula, a text for result output, the number of decimal places and a unit for the result output can be selected. "No unit" is selected with the blank string.

In place of "RSX", a result name may be entered (.TextRS). This name is outputted in the report full, short. It is used for the result and the corresponding mean value.

3.2.2.6. Mode.Def.ComVar.C30 RSX, MNX, EPX, CXX Mode.Def.ComVar.C31 etc., up to .C39

Assignment of common variables.

The values of the common variables are to be found in &Config.ComVar. They can be viewed and entered there, see 3.2.2.50.

3.2.2.7. Mode.Def.Report.Assign

SET, MEAS, CAL: full, short, calc, param, calib, ff

Definition of the report sequence, which is outputted automatically at the end of the determination. Entries of more than one block have to be separated with " , " .

3.2.2.8. Mode.Def.Mean.1.Assign

RS1, RSX, EPX, CXX

Mode.Def.Mean.2.Assign

etc., up to .9

Assignment of the statistics calculations. Valid assignments are a requirement for statistics calculations. In addition, the statistics calculation must be switched on, see 3.2.2.22. Rules for statistics calculations see page 41.

3.2.2.9. Mode.CFmla

Mode.CFmla.1.Value

0...±999 999

Mode.CFmla.2.Value

etc., up to .19

Calculation constants specific to a method. Stored in the method memory of the Titrimo. Operands specific to the sample (3.2.2.52 and 3.2.2.53) and values of common variables (3.2.2.50) on the other hand are not stored with the methods.

3.2.2.10. Mode.Parameter.SET1.EP

pH: 0...±20.00, OFF

U, Ipol: 0...±2000, OFF

Upol: 0...±200.0, OFF

Mode.Parameter.SET1.UnitEp read only

Parameters for SET: Setting the 1st endpoint as pH value, in mV (with U and Ipol) resp. μ A (with Upol). The corresponding unit can be read with .UnitEP. If the value is on "OFF", no further nodes will appear from SET1.

3.2.2.11. Mode.Parameter.SET1.Dyn

pH: 0.01...20.00, OFF

U, Ipol: 1...2000, OFF

Upol: 0.1...200.0, OFF

Mode.Parameter.SET1.UnitDyn read only

Mode.Parameter.SET1.MaxRate 0.01...10...150, max.

Mode.Parameter.SET1.MinRate 0.01...25.0...9999.9

Parameters for SET: Control parameters, see page 25.

.Dyn: Dynamics, control range in pH, mV (with U and Ipol) or μ A (with Upol). The corresponding unit can be read with .UnitDyn.

.MaxRate: Maximum allowed titration rate in mL/min. Max. means maximum possible rate with the Exchange Unit in current use.

.MinRate: Minimum titration rate in μ L/min.

3.2.2.12. Mode.Parameter.SET1.Stop.Type drift, time
Mode.Parameter.SET1.Stop.Drift 1...20...999
Mode.Parameter.SET1.Stop.Time 0...10...999, inf
Mode.Parameter.SET1.Stop.StopT 0...99 999, **OFF**

Parameters for SET: Type and size of the stop criterion of the titration.

.Type: Type of stop criterion after stop drift or switch-off delay time.

.Drift: Stop drift in $\mu\text{L}/\text{min}$. Applies when "drift" has been selected.

.Time: Switch-off delay time in s. Applies when "time" has been selected. "inf" means infinite.

.StopT: Stop time in s. Applies when "time" has been selected and the value of .Time is set to "inf".

3.2.2.13. Mode.Parameter.TitrPara.Direction +, -, **auto**

Parameters for SET: Titration direction.

"auto" means the titration direction is determined automatically by the instrument. If 2 EP's have been set in a SET titration, the titration direction is given by the two EP's. The entry of the titration direction is then invalid.

3.2.2.14. Mode.Parameter.TitrPara.XPause 0...999 999

Parameter for SET: Pause time 1 in s. Is waited off before the dispensing of the start volume.

3.2.2.15. Mode.Parameter.TitrPara.StartV.Type abs., rel., **OFF**
Mode.Parameter.TitrPara.StartV.V 0...999.99
Mode.Parameter.TitrPara.StartV.Factor 0... $\pm$ 999 999
Mode.Parameter.TitrPara.StartV.Rate 0.01...150, **max.**

Parameters for SET: Start volume.

If an absolute start volume (abs.) has been selected, the volume in mL is valid.

A relative start volume (rel.) is dispensed as a function of the sample size:

Start volume in mL = $\text{smpI size} \cdot \text{factor}$

The factor is valid.

The dispensing rate in mL/min applies to both cases. Max. means maximum possible dispensing rate with the Exchange Unit in current use.

3.2.2.16. Mode.Parameter.TitrPara.Pause 0...999 999

Parameter for SET: Pause time 2 in s. Is waited off after the dispensing of the start volume.

3.2.2.17. Mode.Parameter.TitrPara.ExtrT 0...999 999

Parameter for SET: Extraction time in s.

3.2.2.18. Mode.Parameter.TitrPara.MeasInput 1, 2, diff.
Mode.Parameter.TitrPara.Ipol -127...1...+127
Mode.Parameter.TitrPara.Upol -1270...400...+1270
Mode.Parameter.TitrPara.PolElectrTest ON, **OFF**

Parameters for SET:

Selection of the measuring input; valid with measured quantities pH and U.
 "diff." means differential amplifier, see page 124.

With Ipol, the inquiries for the polarization current in μA (Ipol) and .PolElectrTest are valid.

With Upol, the inquiry for the polarization voltage in mV (Upol) is valid. Entry in steps of 10 mV.

Besides .PolElectrTest is valid.

If the test for polarized electrodes is switched on, it is performed on change-over from the inactive state to an active state (titration or conditioning).

3.2.2.19. Mode.Parameter.TitrPara.Temp -170.0...25.0...500.0

Parameter for SET: Titration temperature in $^{\circ}\text{C}$. If a Pt100 or Pt1000 is connected, the temperature is measured continuously and the parameter .Temp is updated.

The temperature is used for the temperature correction in pH measurements.

3.2.2.20. Mode.Parameter.StopCond.VStop.Type abs., rel., OFF Mode.Parameter.StopCond.VStop.V 0...99.99...9999.99 Mode.Parameter.StopCond.VStop.Factor 0...±999 999

Parameter for SET: Stop volume.

If an absolute stop volume (abs.) has been selected, the volume in mL is valid.

A relative stop volume (rel.) is dispensed as a function of the sample size:

Stop volume in mL = smpl size * factor

The factor is valid.

OFF means that the criterion is not monitored.

3.2.2.21. Mode.Parameter.StopCond.FillRate 0.01...150, max.

Parameter for SET: Filling rate in the titration in mL/min. Max. means maximum possible filling rate with the Exchange Unit in current use.

3.2.2.22. Mode.Parameter.Statistics.Status ON, OFF Mode.Parameter.Statistics.MeanN 2...20 Mode.Parameter.Statistics.ResTab.Selected original, delete n, delete all Mode.Parameter.Statistics.ResTab.DeIN 1...20

Entries for the statistics calculations.

.Status: On/off switching. Requirement for statistics calculations is a valid assignment, see 3.2.2.8.

.MeanN: Number of individual results for statistics calculations.

.ResTab.Select: Selection of the table for the statistics calculations.

original: Original table. The original table is (again) set up, i.e. any individual results which have been deleted are reincorporated in the statistics calculations.

delete n: Single result lines are removed from the statistics calculation. All results of the corresponding line in the statistics table are deleted. Specification of the line number in .ResTab.DeIN.

delete all: Clear entire statistics table. The results can not be reactivated.

.ResTab.DelN: Specification of the line number to be deleted.

3.2.2.23.	Mode.Parameter.Presel.Cond	ON, OFF
	Mode.Parameter.Presel.DriftDisp	ON, OFF

Parameters for SET:

.Cond: Conditioning ON/OFF

.DriftDisp: Drift display during conditioning ON/OFF.

3.2.2.24. Mode.Parameter.Presel.DCor.Type auto, man, **OFF**
Mode.Parameter.Presel.DCor.Value **0.0**, 99.9

Parameters for SET:

.Type: Mode of drift correction

.Value: drift value in $\mu\text{L}/\text{min}$.

If "auto" is set, the actual drift value at the start of titration is used for drift correction.

Formula for drift correction: $EP - (\text{drift value} * \text{titration time})$

3.2.2.25. Mode.Parameter.Presel.IReq	id1, id1&2, all, OFF
Mode.Parameter.Presel.SReq	value, unit, all, OFF

Parameters for SET, MEAS: Automatic inquiry after the start of the determination. From such an inquiry, the determination continues if the requested entry/entries is/are made, e.g. &SmpIData.OFFSilo.Id1 (see 3.2.2.73) or with &M\$G, see 3.2.2.1.

\$H\$ is not possible in requests.

3.2.2.26. Mode.Parameter.Presel.ActPuls ON, OFF
Parameter for SET: first, all, cond., OFF

Output of a pulse on the I/O line "Activate", see page 131.

3.2.2.27. Mode.Parameter.Measuring.SignalDrift

pH, U, Ipol, T:	0.5...999, OFF
Upol:	0.05...99.9, OFF

Mode.Parameter.Measuring.UnitSigDrift read only

Mode.Parameter.Measuring.EquTime 0...9999, **OFF**

Parameters for MEAS: Criteria for the measured value acquisition. Measured value drift in mV/min (with pH, U, I_{pol}), μ A/min (with U_{pol}), resp. °C/min (with T). Equilibration time in s. OFF means that the corresponding criterion is switched off. If both criteria are OFF, the measurement continues indefinitely. If the equilibration time has never been edited, it is automatically calculated by the instrument to match the drift, see page 34. After it has been edited once, it remains in force with the set value.

3.2.2.28. Mode.Parameter.Measuring.MeasInput	1, 2, diff.
Mode.Parameter.Measuring.Ipol	±127...1...+127
Mode.Parameter.Measuring.Upol	±1270...400...+1270

Mode.Parameter.Measuring.PolElectrTest ON, OFF

Parameters for MEAS:

Selection of the measuring input; valid with measured quantities pH and U. "diff." means differential amplifier, see page 124.

With Ipol, the inquiries for the polarization current in μA (Ipol) and .PolElectrTest are valid.

With Upol the inquiry for the polarization voltage in mV (Upol) is valid. Entry in steps on 10 mV.

Besides .PolElectrTest is valid.

If the test for polarized electrodes is switched on, it is performed on change-over from the inactive state to the measurement.

3.2.2.29. Mode.Parameter.Measuring.Temp -170.0...25.0...500.0

Parameters for MEAS: Measurement temperature in $^{\circ}\text{C}$. If a Pt100 or Pt1000 is connected, the temperature is measured.

The temperature is used for the temperature correction in pH measurements.

3.2.2.30. Mode.Parameter.Calibration.MeasInput 1, 2, diff.

Parameters for CAL: Selection of the measuring input. "diff." means differential amplifier, see page 124.

3.2.2.31. Mode.Parameter.Calibration.CalTemp -20.0...25.0...120.0

Parameters for CAL: Calibration temperature in $^{\circ}\text{C}$. If a Pt 100 or Pt1000 is connected, the temperature is measured.

3.2.2.32. Mode.Parameter.Calibration.Buffer.1.Value 0...7.00...±20.00

Mode.Parameter.Calibration.Buffer.2.Value 0...4.00...±20.00, OFF etc. up to 9 buffers

Parameters for CAL: pH of buffers. The first buffer which is set to "OFF" determines the number of buffers in the calibration.

3.2.2.33. Mode.Parameter.Calibration.SignalDrift 0.5...2...999, OFF

Mode.Parameter.Calibration.EquTime 0...110...9999, OFF

Parameters for CAL: Criteria for measured value acquisition. Measured value drift in mV/min, equilibration time in s. OFF means that the corresponding criterion is switched off. If both criteria are on OFF, the measured value is acquired immediately.

If the equilibration time has never been edited, it is automatically calculated by the instrument to match the drift, see page 34. After it has been edited once, it remains in force with the set value.

3.2.2.34. Mode.Parameter.Calibration.ElectrodeId up to 8 ASCII char

Parameters for CAL: Electrode identification. It is classified under calibration data, see 3.2.2.59.

3.2.2.35. Mode.Parameter.Calibration.SmplChanger ON, OFF

Parameters for CAL: Calibration at Titrino.

With "ON", there are no hold points in the calibration sequence for entries, the first buffer is measured directly.

3.2.2.36. Mode.Parameter.Calibration.ActPulse first, all, **OFF**

Parameters for CAL: Output of a pulse on the I/O line "Activate", see page 131.

3.2.2.37. UserMeth.FreeMem read only

Memory space, available for user methods or silo lines. \$Q sends the number of free bytes, e.g. "4928".

3.2.2.38. UserMeth.Recall	\$G
UserMeth.Recall.Name	up to 8 ASCII characters
UserMeth.Store	\$G
UserMeth.Store.Name	up to 8 ASCII characters
UserMeth.Delete	\$G
UserMeth.Delete.Name	up to 8 ASCII characters
UserMeth.DeIAI	\$G

Management of the internal method memory: Load, store and delete methods. An action is performed if "\$G" is sent to the corresponding node just after entering the name.

Do not use blank characters before and after method name!

.DeIAI: Deletes all methods in the user memory.

3.2.2.39. UserMeth.List.1.Name	read only
UserMeth.List.1.Mode	read only
UserMeth.List.1.Quantity	read only
UserMeth.List.1.Bytes	read only
UserMeth.List.1.Checksum	read only
for each method	

List of the methods in the user method memory with the following characteristics:

.Name: Name of the method

.Mode: Mode

.Quantity: Measured quantity

.Bytes: Number of bytes of the user memory used by the method

.Checksum: Checksum of the method, see 3.2.2.59.

3.2.2.40. Config.PeriphUnit.CharSet Epson, Seiko, Citizen, HP, **IBM**

Selection of the character set and the graphics control characters of the Titrino.

IBM means the IBM character set following character set table 437 and IBM graphics control characters. Select 'IBM' for work with the computer.

3.2.2.41. Config.PeriphUnit.Balance	Sartorius, Mettler, Mettler AT,
	AND, Precisa
Config.PeriphUnit.Plot	U, dU/dt, V, dV/dt, U(rel), T

Selection of the balance type and the signal at the analog output.

3.2.2.42. Config.Aux.Language **english**, deutsch, francais, espanol, italiano, portugese, svenska

Selection of the dialog language.

3.2.2.43. Config.Aux.Set \$G
Config.Aux.Set.Date YYYY-MM-DD
Config.Aux.Set.Time HH:MM

Date and time.

Input format of the date: Year-month-day, two-digit, enter leading zeros.

Input format for the time: Hours:minutes, two-digit, enter leading zeros.

Date and time have to be set with &Config.Aux.Set \$G just after entry of the value.

3.2.2.44. Config.Aux.RunNo 0...9999

Current sample number.

Set to 0 on power on and initialization. After 9999, counting starts again at 0.

3.2.2.45. Config.Aux.AutoStart 1...9999, **OFF**

Number of automatic, internal starts.

3.2.2.46. Config.Aux.StartDelay 0...999 999

Start delay time in s. During this time, the data of the preceding determination are retained.

3.2.2.47. Config.Aux.DevName up to 8 ASCII characters

Name of the instrument for connections with several units. It is advisable to use only the letters A...Z (ASCII No. 65...90), a...z (ASCII No. 97...122) and the numbers 0...9 (ASCII No. 48...57) when the function Setup.AutoInfo (3.2.2.86) is used at the same time.

If a name has been entered, it will be printed out in the result report (full, short).

3.2.2.48. Config.Aux.Prog read only

Output of the program version.

The Titrino sends "719.0021" on requests with \$Q.

3.2.2.49. Config.RSSet \$G
Config.RSSet.Baud 300, 600, 1200, 2400, 4800, **9600**
Config.RSSet.DataBit 7, **8**
Config.RSSet.StopBit 1, **2**
Config.RSSet.Parity even, odd, **none**
Config.RSSet.Handsh **HWs**, HWf, SWchar, SWline, none

\$G sets all RS settings. The changes are performed only if the instrument is inactive. After the setting of the interface parameters, wait at least 2 s to allow the components to equilibrate.

Settings of the values for the data transmission via the RS interface: baud rate, data bit, stop bit, parity and type of handshake, see also page 95ff.

The setting of the values must be initiated with \$G immediately after entry of the values.

3.2.2.50. Config.ComVar.C30

with up to .C39, etc.

0... $\pm 999\,999$

Values of the common variables from C30 up to C39. Insert the common variables directly or describe the determination results directly from the method, see 3.2.2.6.

3.2.2.51. SmpIData.Status

ON, OFF

On/off switching of silo memory. When the silo memory is switched on, the sample data are fetched from the lowest valid silo line.

3.2.2.52. SmpIData.OFFSilo.Id1

up to 8 ASCII characters

SmpIData.OFFSilo.Id2

up to 8 ASCII characters

SmpIData.OFFSilo.Id3

up to 8 ASCII characters

SmpIData.OFFSilo.ValSmpI

6-digits, sign and decimal point

SmpIData.OFFSilo.UnitSmpI

up to 5 ASCII characters

Current sample data.

The identifications Id1...Id3 can be used in formulas as sample-specific calculation constants C21...C23.

If "no unit" is desired for the unit of the sample size, the blank string must be entered.

3.2.2.53. SmpIData.ONSilo.Counter.MaxLines

read only

SmpIData.ONSilo.Counter.FirstLine

read only

SmpIData.ONSilo.Counter.LastLine

read only

Information on silo memory.

.MaxLines: Maximum possible number of silo lines.

.FirstLine: Lowest valid silo line.

.LastLine: Last occupied silo line.

3.2.2.54. SmpIData.ONSilo.EditLine.1.Method

up to 8 ASCII characters

SmpIData.ONSilo.EditLine.1.Id1

up to 8 ASCII characters

SmpIData.ONSilo.EditLine.1.Id2

up to 8 ASCII characters

SmpIData.ONSilo.EditLine.1.Id3

up to 8 ASCII characters

SmpIData.ONSilo.EditLine.1.ValSmpI

6-digits, sign and dec.point

SmpIData.ONSilo.EditLine.1.UnitSmpI

up to 5 ASCII characters

etc., up to .99

Contents of a silo line.

.Method: Method used to process the sample, from the method memory or from the card.

.Id: The identifications Id1...Id3 can also be used as sample-specific calculation constants C21...C23 in formulas.

.UnitSmpl: If "no unit" is desired for the sample size, the blank string must be entered.

3.2.2.55. SmpData.ONSil.DelLine \$G
SmpData.ONSil.DelLine.LineNum 1...99, **OFF**

Deletion of a silo line. The line # is deleted with &SmpData.ONSil.DelLine \$G. If a formerly deleted line is edited again, it becomes valid (function "undelete").

3.2.2.56. SmpData.ONSil.DelAll \$G
Deletes the entire silo memory. Must be triggered with \$G.

3.2.2.57. SmpData.ONSil.CycleLines ON, **OFF**
Silo data cycling.
With "ON", executed lines are copied to the next free silo lines, see page 53.
Exercise caution if you edit the silo memory during the determinations!

3.2.2.58. Info.Report \$G
Info.Report.Select configuration, parameters, smpl data,
statistics, silo, calib, C-fmla, def,
user method, **full**, short, calc, all, ff

\$G sends the selected report to the COM which is set in
&Config.PeriphUnit.RepToComport:
configuration: Configuration report. Is not accessible during a running determination.

parameters: Parameter report of the current method. During a running determination only "live"-parameters are accessible.

smpl data: Current sample data.

statistics: Statistics table with the individual results.

silo: Contents of the silo memory.

calib: Calibration data of the measuring input in the current method.

C-fmla: Contents of the <C-fmla> key.

def: Contents of the <def> key.

user method: Contents of the method memory.

full: Full result report of the last completed determination.

short: Short result report of the last completed determination.

calc: Calculation report of the current method.

all: All reports.

ff: Form feed on printer.

Reports which are sent from the Titrino are marked with space (ASCII 32) and ' at the beginning. Then an individual identifier for each report follows. Reports which are triggered by RS232 (\$G) have the same introducer but without preceding space, i.e. they start with '.

3.2.2.59. Info.CalibrationData \$G
Info.CalibrationData.Inp1.pHas -20.00...**7.00**...+20.00
Info.CalibrationData.Inp1.Slope -9.999...**1.000**...+9.999
Info.CalibrationData.Inp1.Temp -170.0...**25.0**...+500.0

Info.CalibrationData.Inp1.Date read only
Info.CalibrationData.Inp1.ElectrodeId read only
 identical for .Inp2 and .Diff

pH calibration data for measuring input 1. After the calibration, the data are entered automatically together with the date of the calibration and the electrode identification, see 3.2.2.34.

Calibration data can be entered. They are accepted with &Info.CalibrationData \$G. If calibration data are entered, the calibration date is deleted.

3.2.2.60. Info.Checksums \$G

Info.Checksums.ActualMethod read only

The checksums can be used to identify the content of a file unequivocally, e.g. files with identical content

have identical results of the checksums. An empty file has checksum "0". The calculation of the checksums is triggered with \$G.

.ActualMethod: Result of the checksum of the current method in the working memory. Identical methods with different method names have the same results of the checksum.

3.2.2.61. Info.DetermData \$G

Info.DetermData.Write ON, OFF

Info.DetermData.ExV read only/read + write

Recalculation of the measuring data is triggered with \$G.

.Write: With "ON", the following nodes can be overwritten:

&Info.TitrResults.Var.C4X (X = 0...5) and &Mode.Name.

.ExV: Volume of the exchange unit, with which the determination was executed

3.2.2.62. Info.TitrResults.RS.1.Value read only

etc., up to .9

Info.TitrResults.EP.1.V read only

Info.TitrResults.EP.1.Meas read only

etc., up to .2

Info.TitrResults.Var.C40 read only/read + write

etc., up to .C47

.RS: Values of the calculated results.

.EP: Endpoints with SET:

Volume coordinate in mL, e.g. "1.2340"

Measured value coordinate in pH "5.12", mV (with U and Ipol) "-241" or μ A (with Upol) "43.7".

.Var: Various variables. You may overwrite the variables C40...C45, see 3.2.2.95.

C40: Initial measured value in pH "5.12", mV (with U and Ipol)

"41", μ A (with Upol) "43.7" or °C (with T) "25.0". In MEAS final measured value.

C41: End volume with SET in ml, "12.5360".

C42: Time from start of titration to end in s, "62".

C43: Volume drift on start of a SET titration from the conditioning in ul/min, "3.5".

C44: (Last measured) temperature in °C. Used for the temperature compensation in pH measurements.

C45: Start volume with SET in ml, "2.800".

C46: Asymmetry pH of CAL, "6.89".

C47: Relative electrode slope of CAL, "0.9950".

3.2.2.63.	Info.StatisticsVal.ActN	read only
	Info.Statistics.1.Mean	read only
	Info.Statistics.1.Std	read only
	Info.Statistics.1.RelStd	read only
	etc. up to .9	

The current values of the statistics calculation.

\$Q sends, e.g.

ActN: Current value of the individual results "3"

Data for MN1:

Mean: Mean value (decimal places as in result) "3.421"

Std: Standard deviation (1 decimal place more than in result) "0.0231"

RelStd: Relative standard deviation (in %, 2 decimal places) "0.14"

3.2.2.64.	Info.ActualInfo.Inputs.Status	read only
	Info.ActualInfo.Inputs.Change	read only
	Info.ActualInfo.Inputs.Clear	\$G
	Info.ActualInfo.Outputs.Status	read only
	Info.ActualInfo.Outputs.Change	read only
	Info.ActualInfo.Outputs.Clear	\$G

Status sends the current status of the I/O lines, Change sends the information regarding whether a change in status of a line has taken place since the last clearing, Clear clears the change information. For the output, there is a conversion from binary to decimal, e.g.

	0 0 0 0 1 0 1 0
Line No.	7 6 5 4 3 2 1 0

Output: $2^1 + 2^3 = "10"$

1 means ON or change; 0 means OFF or no change.

The lines are assigned as follows (see also pages 129ff):

Inputs:		Outputs:	
0	Start (pin 21)	0	Ready (pin 5)
1	Stop (pin 9)	1	Cond. ok (pin 18)
2	Enter (pin 22)	2	Titration (pin 4)
3	Clear (pin 10)	3	EOD (pin 17)
4	Smpl Ready (pin 23)	4	Freely selectable (pin 3)
5	pin 11	5	Error (pin 16)
6	pin 24	6	Activate (pin 1)
7	pin 12	7	Pulse for recorder (pin 2)

3.2.2.65.	Info.ActualInfo.Assembly.CyclNo	read only
------------------	--	-----------

\$Q sends the current cycle number of the voltage measurement cycle, e.g. "127". From the cycle number and the cycle time (see 3.2.2.71), a time frame can be set up.

The cycle number is set to 0 on switching on the instrument, on every start and for QuickMeas. It is incremented as long as the instrument remains switched on.

3.2.2.66. Info.ActualInfo.Assembly.Counter.V read only
Info.ActualInfo.Assembly.Counter.Clear \$G
 \$Q sends the volume. With the function &Info.Assembly.Counter.Clear \$G, the volume counter is set to zero.

3.2.2.67. Info.ActualInfo.Assembly.Meas read only
Info.ActualInfo.Assembly.PistonPos read only
 \$Q sends the current measured value from the assembly.
 Piston position: 0 = initial position (filled)
 10 000 = end position (empty)

3.2.2.68. Info.ActualInfo.Titrator.CyclNo read only
Info.ActualInfo.Titrator.V read only
Info.ActualInfo.Titrator.Meas read only
Info.ActualInfo.Titrator.dVdt read only
Info.ActualInfo.Titrator.dMeasdt read only
Info.ActualInfo.Titrator.dMeasdV read only
 \$Q sends the current values in the following formats:

	SET	MEAS	CAL
CyclNo	127	127	127
V(ml)	1.2345	-	-
Meas:			
pH	3.6(mV)	3.345	3.345
U, Ipol (mV)	-345.6	-345.6	-
Upol (uA)	-12.5	-12.5	-
T (°C)	-	25.0	-
dVdt (ul/s)	2.5142	-	-
dMeasdt			
pH,U,Ipol mV/s	0.7957	0.7957	0.7957
Upol (uA/s)	0.7957	0.7957	-
T (°C/s)	-	0.7957	-
dMeasdV (mV/ul)	10.6326	-	-

NV: Not Valid. If in the signal drift is OFF in modes MEAS and CAL, the signal drift is NV.

OV will be sent for "overrange".

A time frame can be set up from the cycle number and the cycle time (see 3.2.2.71). The cycle number is set to 0 at the start of a method and it is incremented until the end of the method.

3.2.2.69. Info.ActualInfo.EP.Index read only
Info.ActualInfo.EP.X read only
Info.ActualInfo.EP.Y read only
 .EP.X"1.234" Volume coordinate of the EP
 .EP.Y"5.34" Measured value coordinate of the EP

3.2.2.70. Info.ActualInfo.Display.L1 up to 24 ASCII characters

719 S Titrimo

Info.ActualInfo.Display.L2 up to 24 ASCII characters
 Lines of the display. The display can be written to from the computer. Proceed as follows:
 Lock the display, see 3.2.2.80.
 \$Q sends the contents of the corresponding display line.

3.2.2.71. Info.Assembly.CycleTime read only
Info.Assembly.ExV read only
 Inquiries regarding basic variables of the assembly: Cycle time in s, volume of the active Exchange Unit in mL.

3.2.2.72. Assembly.Bur.Rates.Forward.Selected digital, **analog**
Assembly.Bur.Rates.Forward.Digital 0...150, **max.**
Assembly.Bur.Rates.Reverse.Selected digital, **analog**
Assembly.Bur.Rates.Reverse.Digital 0...150, **max.**

Expel and aspirating rate.
 Digital or analog control. With digital control, the inputted value applies (in mL/min). "max." means maximum possible rate with the Exchange Unit in current use.
 Analog means rate control with the analog potentiometer on Titrino.

3.2.2.73. Assembly.Bur.Fill \$G, \$H, \$C
 \$G starts the 'FILL' mode of the buret function.

3.2.2.74. Assembly.Bur.ModeDis \$G, \$S, \$H, \$C
Assembly.Bur.ModeDis.Selected **volume**, time
Assembly.Bur.ModeDis.V 0.0001...0.1...9999
Assembly.Bur.ModeDis.Time 0.25...1...86400
Assembly.Bur.ModeDis.VStop 0.0001...9999, **OFF**
Assembly.Bur.ModeDis.AutoFill **ON**, **OFF**

Dispensing mode with parameters. The dispensing mode can only be started and stopped via the RS Control. During a running dosification, no method can be started at the Titrino.

.Selected: Dispensing of volume increments or during a preset time.
 .Volume, .Time: Size of the volume increments or entry of time.
 .VStop: Limit volume for the dispensing.
 .AutoFill: ON means automatic filling after every dispensing.

3.2.2.75. Assembly.Meas.Status **ON**, **OFF**
Assembly.Meas.MeasInput 1, 2, Diff., Ipol, Upol, Temp
Assembly.Meas.Ipol $\pm 127 \dots 1 \dots +127$
Assembly.Meas.Upol $\pm 1270 \dots 400 \dots +1270$

Measurement in assembly. The measuring function can only be started via RS Control. When the measuring function is switched on, no method can be started at the Titrino.

.Input: Selection of the potentiometric measuring input 1, 2, diff., polarized electrodes or temperature.

.Ipol: Polarization current in μA .
 .Upol: Polarization potential in mV, entry in steps of 10 mV.

3.2.2.76. Assembly.Outputs.AutoEOD **ON, OFF**
Assembly.Outputs.SetLines **\$G**
Assembly.Outputs.SetLines.L0 active, inactive, pulse, **OFF**
 up to .L 3
Assembly.Outputs.ResetLines **\$G**

Setting the I/O output lines.

.AutoEOD: The automatic output of the EOD (End of Determination) at the end of the determination can be switched off. Thus, for example, in conjunction with a Titrino several determinations can be performed in the same beaker. Before AutoEOD is switched on, line 3 must be set to "OFF".

.SetLines: With \$G, all lines are set.

.SetLines.LX: Set the line LX. "active" means setting of a static signal, "inactive" means resetting of the signal, "pulse" means output of a pulse of app. 150 ms, "OFF" means the line is not operated, see also page 130.

Warnings:

- L2 is the EOD line. If you have "AutoEOD" set to "ON", an active line 2 is set to "inactive" by the EOD pulse.
- L3 is the line of the activate pulse. An active line 3 is set to "inactive" by the activate pulse.

.ResetLines: Lines are set to the inactive status (= high).

3.2.2.77. Setup.Keycode **ON, OFF**

ON means the key code of a key pressed on the Titrino is outputted. The key code comprises 2 ASCII characters; table of the keys with their code, see page 108. A keystroke of key 11 is sent as follows:

#11

The beginning of the message is marked by a space (ASCII 32).

3.2.2.78. Setup.Tree.Short **ON, OFF**
Setup.Tree.ChangedOnly **ON, OFF**

Definition of the type of answer to \$Q.

.Short: With "ON", each path is sent with only the necessary amount of characters in order to be unequivocal (printed in bold in this manual). A combination of .Short and .ChangedOnly is not possible.

.ChangedOnly: Sends only the changed values, i.e. values which have been edited. All paths are sent absolute, i.e. from the root.

3.2.2.79. Setup.Trace **ON, OFF**

The Titrino automatically reports when a value has been confirmed with <enter> at the Titrino. Message, e.g.:

&SmplData.OFFSilo.Id1"Trace"

The beginning of the message is marked by a space (ASCII 32).

3.2.2.80.	Setup.Lock.Keyboard	ON, OFF
	Setup.Lock.Config	ON, OFF
	Setup.Lock.Parameter	ON, OFF
	Setup.Lock.SmplData	ON, OFF
	Setup.Lock.UserMeth.Recall	ON, OFF
	Setup.Lock.UserMeth.Store	ON, OFF
	Setup.Lock.UserMeth.Delete	ON, OFF
	Setup.Lock.Display	ON, OFF

ON means disable the corresponding function:

- .Keyboard: Disable all keys of the Titrimo
- .Config: Disable the <configuration> key
- .Parameter: Disable the <parameter> key
- .SmplData: Disable the <smpl data> key
- .UserMeth.Recall: Disable "recall" in <user meth> key
- .UserMeth.Store: Disable "store" in <user meth> key
- .UserMeth.Delete: Disable "delete" in <user meth> key
- .Display: Disable the display, i.e. it will not be written to by the device program of the Titrimo and can be operated from the computer.

3.2.2.81.	Setup.Mode.StartWait	ON, OFF
	Setup.Mode.FinWait	ON, OFF

Holding points in the method sequence. If they are "ON", the sequence stops until "OFF" is sent. Switching the instrument on sets both nodes to OFF:

- .StartWait: Holding point right after starting a method or submethod in TIP (holding point after AutoInfo !"T.GC").
- .FinWait: Holding point at the end a method or submethod in TIP (holding point after AutoInfo !"T.F").

3.2.2.82.	Setup.SendMeas.SendStatus	ON, OFF
	Setup.SendMeas.Interval	0.08...4...16200, MPList

.SendStatus: ON means the automatic transmission of measured values (see 3.2.2.84 and 3.2.2.85) in the inputted interval is active.

.Interval: Time interval (in s) for the automatic transmission of associated measured values defined under points 3.2.2.84 and 3.2.2.85. The inputted value is rounded off to a multiple of 0.08. The smallest possible time interval depends on the number of measured values which have to be sent, on the baud rate, on the load on the interface and on the type of device connection.

The automatic transmission is switched on/off with 'SendStatus'.

3.2.2.83.	Setup.SendMeas.Select	Assembly, Titrator
------------------	------------------------------	---------------------------

Selection of the unit of which the measured values should be sent (3.2.2.84 or 3.2.2.85).

3.2.2.84.	Setup.SendMeas.Assembly.CyclNo	ON, OFF
	Setup.SendMeas.Assembly.V	ON, OFF
	Setup.SendMeas.Assembly.Meas	ON, OFF

Selection of the values from Assembly for the output in the set time interval (see 3.2.2.82):

- .CyclNo: Cycle number of the potential measurement. Together with the cycle time (3.2.2.81), a time frame can be set up.
The cycle number is set to 0 on switching on the instrument and it is always incremented as long as the instrument remains switched on.
- .V: Volume
- .Meas: Measured value associated to the cycle number.
The unit "assembly" must be preset (see 3.2.2.83).

3.2.2.85.	Setup.SendMeas.Titrator.CyclNo	ON, OFF
	Setup.SendMeas.Titrator.V	ON, OFF
	Setup.SendMeas.Titrator.Meas	ON, OFF
	Setup.SendMeas.Titrator.dVdt	ON, OFF
	Setup.SendMeas.Titrator.dMeasdt	ON, OFF
	Setup.SendMeas.Titrator.dMeasdV	ON, OFF

Selection of the values from the titrator which are sent in the set time interval (see 3.2.2.82, formats see 3.2.2.68):

- .CyclNo: Cycle number. Together with the cycle time (3.2.2.71), a time frame can be set up. The other data belong to the corresponding cycle number. The cycle number is set to 0 at the start of a method and it is incremented until the end of the method.
- .V: Volume.
- .Meas: Measuring value
- .dVdt: associated volume drift.
- .dMeasdt: associated measured value drift.
- .dMeasdV: associated 1st derivative of the titration curve.
The unit "titrator" must be preset (see 3.2.2.83).

3.2.2.86.	Setup.AutoInfo.Status	ON, OFF
	Setup.AutoInfo.P	ON, OFF
	Setup.AutoInfo.T.R	ON, OFF
	Setup.AutoInfo.T.G	ON, OFF
	Setup.AutoInfo.T.GC	ON, OFF
	Setup.AutoInfo.T.S	ON, OFF
	Setup.AutoInfo.T.B	ON, OFF
	Setup.AutoInfo.T.F	ON, OFF
	Setup.AutoInfo.T.E	ON, OFF
	Setup.AutoInfo.T.H	ON, OFF
	Setup.AutoInfo.T.C	ON, OFF
	Setup.AutoInfo.T.O	ON, OFF
	Setup.AutoInfo.T.N	ON, OFF
	Setup.AutoInfo.T.Re	ON, OFF
	Setup.AutoInfo.T.Si	ON, OFF
	Setup.AutoInfo.T.M	ON, OFF
	Setup.AutoInfo.T.EP	ON, OFF
	Setup.AutoInfo.T.RC	ON, OFF
	Setup.AutoInfo.I	ON, OFF
	Setup.AutoInfo.O	ON, OFF

ON means that the Titrino reports automatically the moment the corresponding change occurs.

.Status: Global switch for all set AutoInfo.

.P PowerOn: Simulation of power on (3.2.2.87). Not from mains.

Messages from node .T, Titrator:

.T.R Ready: Status 'Ready' has been reached.

.T.G Go: Instrument has been started.

.T.GC GoCommand: Instrument has received a go command.

.T.S Stop: Status 'Stop' has been reached.

.T.B Begin of sequence.

.T.F Final: End of determination, the final steps will be carried out.

.T.E Error. Message together with error number, see page 60ff.

.T.H Hold: Status 'Hold' has been reached.

.T.C Continue: Continue after hold.

.T.O Conditioning OK: EP reached (in SET with conditioning).

.T.N Conditioning Not OK: EP not reached (in SET with conditioning).

.T.Re Request: In the inquiry of an identification or the sample size after start of titration.

.T.Si SiloEmpty: Silo empty, i.e. the last line has been removed from the silo memory.

.T.EP EPList: Entry into EP list (with SET)

.T.RC Results have been recalculated.

Messages for changings in the I/O lines. If the changings are made simultaneously, there is 1 message. Pulses receive 2 messages: one message each for line active and inactive.

.I Input: Change of an input line.

.O Output: Change of an output line (except 7, pin 2, for recorder pulses).

If a change occurs that requires a message, the Titrino sends space (ASCII 32) and ! as an introducer. This is followed by the name of the device (see 3.2.2.47). Special ASCII characters in the device name are ignored. If no device name has been entered, only ! is sent. Finally the Titrino sends the information which node has triggered the message.

Example: !John".T.Si": The message was triggered from instrument "John", node .T.Si

3.2.2.87. Setup.PowerOn \$G

Simulation of 'power on'. The device has the same status as after power on: The cylinder is filled, error messages deleted and the current sample number set to 0. The method last used is ready for operation.

3.2.2.88. Setup.Initialise \$G

Setup.Initialise.Select **ActMeth, Silo, Calib, Config, Assembly, Setup, All**

Setting of default values for the following areas:

ActMeth: Current method. Parameters, calculations, and assignments for the data output, operands C01...C19.

Silo: The silo memory is deleted. Same function as delete entire silo.

Calib: pH calibration data for all measuring inputs.

Config: All values under &Config.

Assembly: All values under &Assembly.

Setup: All values under &Setup.
All: Values of the entire tree (except silo and method memory).
The action must be triggered with &Setup.Initalize \$G.

3.2.2.89. Setup.RamInit \$G

Initializes instrument, see page 116. All parameters are set to their default value and error messages are cleared. The user and silo memories will be deleted. The user memory contains the default user methods from Metrohm.

3.2.2.90. Setup.InstrNo \$G **Setup.InstrNo.Value** **serial number**, 8 ASCII characters

Instrument identification for report output.
Set the value with &Setup.InstrNo \$G .

3.2.2.91. Diagnose.Report \$G

Output of the report containing the adjustment parameters. The Titrino has to be in its inactive basic state.

3.3 Properties of the RS 232 Interface

Data Transfer Protocol

The Titrino is configured as DTE (Data Terminal Equipment).

The RS 232 interface has the following technical specifications:

- Data interface according to the RS 232C standard, adjustable transfer parameters, see page 10.
- Max. line length: 512 characters
- Control characters: C_R (ASCII DEC 13)
L_F (ASCII DEC 10)
XON (ASCII DEC 17)
XOFF (ASCII DEC 19)
- Cable length: max. approx. 15 m

Start	7 or 8 Data Bit	Parity Bit	1 or 2 Stop Bit
-------	-----------------	------------	-----------------

Only a shielded data cable (for example, METROHM D.104.0201) may be used to couple the Titrino with foreign devices. The cable shield must be properly grounded on both instruments (pay attention to current loops; always ground in a star-head formation). Only plugs with sufficient shielding may be used (for example, METROHM K.210.0381 with K.210.9045).

3.3.1 Handshake

Software-Handshake, SWchar


Handshake inputs on the Titrino (CTS, DSR, DCD) are not checked.
Handshake outputs (DTR, RTS) are set by the Titrino.

As soon as a L_F is recognized, the Titrino sends XOFF. It can then receive 6 extra characters and store them.


However, the Titrino also sends XOFF if its input buffer contains 60 characters. After this, it can receive maximum 22 extra characters (incl. L_F).

If the transmission is interrupted for the time of 4 characters after the Titrino has sent XOFF, the string received earlier is processed even if no L_F has been sent.

Titrimo as Receiver :

Titrimo as Sender :


Software-Handshake, SWline


Handshake input ports on the Titrimo (CTS, DSR, DCD) are not checked.
Handshake output ports (DTR, RTS) are set by the Titrimo.

The Titrimo is equipped with an input buffer that can accommodate a string of up to 80 characters + $C_R L_F$. As soon as an L_F is recognized, the Titrimo sends XOFF. After this, it can receive maximum 6 extra characters and store them. The string sent previously is now processed by the Titrimo. Afterwards, the Titrimo sends XON and is again ready to receive.

Titrimo as Receiver :


Titrimo as Sender:


Titrimo transmission can be stopped by external instruments with XOFF. After XOFF is received the Titrimo completes sending the line already started. If data output is disabled for more than 3 s by XOFF, E43 appears in the display.

Hardware-Handshake, HWs

Titrimo as Receiver :

Titrimo as Sender:


The data flow can be interrupted by deactivating the CTS line.


Hardware-Handshake, HWf

All handshake inputs are checked at the Titrimo, handshake outputs are set.

Titrimo as Receiver :


Titrimo as Sender:

The data flow can be interrupted by deactivating the CTS line.

3.3.2 Pin Assignment

RS232C Interface

Protective earthing


Direct connection from cable plug to the protective ground of the instrument.

Polarity allocation of the signals

- Data lines (TxD, RxD)
 - voltage negative (< -3 V): signal state "ON"
 - voltage positive (> +3 V): signal state "ZERO"
- control or message lines (CTS, DSR, DCD, RTS, DTR)
 - voltage negative (< -3 V): OFF state
 - voltage positive (> +3 V): ON state

In the transitional range from +3 V to -3 V the signal state is undefined.

Driver 14C88 according to EIA RS 232C specification
Receiver 14C89 " "

Contact arrangement at plug (female) for RS 232C socket (male)

View of soldered side of plug

Ordering numbers:
K.210.9004 and K.210.0001

No liability whatsoever will be accepted for damage or injury caused by improper interconnection of instruments.

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Certificate of Conformity and System Validation

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

Name of commodity:	719 S Titrino
System software:	Stored in ROMs
Name of manufacturer:	Metrohm Ltd., Herisau, Switzerland

This Metrohm instrument has been built and has undergone final type testing according to the standards:

Electromagnetic compatibility: Emission
EN50081-1, EN50081-2, EN55022 (class B)

Electromagnetic compatibility: Immunity
EN50082-1, IEC61000-6-2, Namur, IEC61000-4-2, IEC61000-4-3,
IEC61000-4-5, IEC61000-4-6, IEC61000-4-11

Safety specifications
IEC61010-1, EN61010-1

It has also been certified by the Swiss Electrotechnical Association (SEV), which is member of the International Certification Body (CB/IEC).

The technical specifications are documented in the instruction manual.

The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance. The features of the system software are documented in the instruction manual.

Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.

Herisau, October 30, 2001

Dr. J. Frank

Ch. Buchmann

Development Manager

Production and
Quality Assurance Manager

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