


# Rod Stirrer

# 722

Series 01 ...



- 1 Rod Stirrer**  
Connect stirrer cable to the shorter cable of the controller.
- 2 6.2137.000 Controller**  
Route longer cable of the controller to the power supply.  
The stirring speed is controlled with the scaled rotary wheel.  
If the rotary stop latches at one end of the scale, turn to the other end until it snaps in to reset the scale.
- 3 6.1909.010 Stirrer Propeller**  
Lead rod stirrer through the opening in the holder and attach stirrer propeller from below.
- 4 6.1910.000 Switch Ring**  
For use without 6.2021.030 Electrode Holder.  
Place switch ring on controller. Position "ON" at back means "stirrer on". In all other positions, the stirrer is switched off.

## Technical Data

<b>Stirrer speed</b> measured in 30 ml H <sub>2</sub> O, with 1 electrode	440...2800 rpm
<b>Power supply</b>	7...12 V DC/ 0.2 A
<b>Power consumption</b>	ca. 3 VA
<b>Material</b>	PP
<b>Dimensions</b>	
Length with stirrer propeller	ca. 25 cm
Largest diameter	27 mm
<b>Weight</b>	
incl. stirrer propeller and controller	ca. 300 g

## Standard equipment

<b>722 Rod Stirrer</b>	<b>2.722.0010</b>
including the following accessories	
1 Stirrer propeller	6.1909.010
1 Controller	6.2137.000
1 Instructions for use	8.722.1003
<b>722 Rod Stirrer</b>	<b>2.722.0011</b>
Accessories as for 2.722.0010, but in addition:	
1 Switch ring	6.1910.000
1 Mains adapter 110 V/9V DC, USA...	6.2119.020
<b>722 Rod Stirrer</b>	<b>2.722.0014</b>
Accessories as for 2.722.0010, but in addition:	
1 Switch ring	6.1910.000
1 Mains adapter 220 V/9V DC, Europe...	6.2117.020

Ionenanalytik • Analyse des ions • Ion analysis • Análisis iónico

**722 Propeller Rod Stirrer**



Metrohm Ltd.  
CH-9101 Herisau  
Switzerland  
Phone +41 71 53 85 85  
Fax +41 71 53 89 01

## EC Declaration of Conformity

The METROHM AG company, Herisau, Switzerland hereby certifies, that the instrument:

### **722 Propeller Rod Stirrer**

meets the requirements of EC Directives 89/336/EWG and 73/23/EWG.

---

#### Source of the specifications:

EN 50082-1 Electromagnetic compatibility, basic specification  
Interference Immunity

#### Description of the instrument:

Rod stirrer for 727 Titration Stand and for general application on titration heads with standard ground joint; with separate controller and stirring propeller

---

Herisau, December 6, 1995

Dr. J. Frank

Development Manager

Ch. Buchmann

Production and  
Quality Assurance Manager