

728 Magnetic Stirrer



Manual

8.728.8001EN / 2019-11-29



Metrohm AG

CH-9100 Herisau

Switzerland

Phone +41 71 353 85 85

Fax +41 71 353 89 01

info@metrohm.com

www.metrohm.com

728 Magnetic Stirrer

Manual

Technical Communication
Metrohm AG
CH-9100 Herisau
techcom@metrohm.com

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Table of contents

1	Introduction	1
1.1	Instrument description	1
1.1.1	Model versions	1
1.1.2	Connectors	1
1.1.3	Intended use	1
1.2	About the documentation	2
1.2.1	Symbols and conventions	2
1.3	Safety instructions	3
1.3.1	General notes on safety	3
1.3.2	Electrical safety	3
1.3.3	Tubing and capillary connections	4
1.3.4	Flammable solvents and chemicals	4
1.3.5	Recycling and disposal	4
2	Overview of the instrument	6
3	Installation	8
3.1	Setting up the instrument	8
3.1.1	Packaging	8
3.1.2	Checks	8
3.1.3	Location	8
3.2	Installing the magnetic stirrer	8
3.3	Connecting the magnetic stirrer	8
4	Operation and maintenance	9
4.1	General notes	9
4.1.1	Care	9
4.1.2	Maintenance by Metrohm Service	9
5	Technical specifications	10
5.1	Magnetic stirrer	10
5.2	Power supply	10
5.3	Interfaces and connectors	10
5.4	Ambient temperature	10
5.5	Reference conditions	10
5.6	Dimensions/material	11
6	Accessories	12
	Index	13



Table of figures

Figure 1	Front	6
Figure 2	Rear with stand	7

1 Introduction

1.1 Instrument description

The 728 Magnetic Stirrer with base plate, support rod and electrode holder supplements Titrino and Dosino units.

1.1.1 Model versions

The 728 Magnetic Stirrer is available in the following six versions:

2.728.0010	728 Magnetic Stirrer	
2.728.0021	728 Magnetic Stirrer	With 115 V US power supply unit
2.728.0024	728 Magnetic Stirrer	With 230 V EURO power supply unit
2.728.0031	728 Magnetic Stirrer	With 115 V US power supply unit and stand
2.728.0034	728 Magnetic Stirrer	With 230 V EURO power supply unit and stand
2.728.0040	728 Magnetic Stirrer	Complete for Titrino/Dosino

Each version includes different accessories according to its use (*see chapter 6, page 12*).

1.1.2 Connectors

The 728 Magnetic Stirrer needs 5 to 12 V DC voltage.

The 728 Magnetic Stirrer is supplied with power via a direct connection to a Titrino, Dosino or via a power supply unit.

1.1.3 Intended use

The 728 Magnetic Stirrer is designed for usage in analytical laboratories.

This instrument is suitable for stirring chemicals and flammable samples. Usage of the 728 Magnetic Stirrer therefore requires the user to have basic knowledge and experience in handling toxic and caustic substances. Knowledge with respect to the application of the fire prevention measures prescribed for laboratories is also mandatory.



1.2 About the documentation



CAUTION

Please read through this documentation carefully before putting the instrument into operation. The documentation contains information and warnings which the user must follow in order to ensure safe operation of the instrument.

1.2.1 Symbols and conventions

The following symbols and formatting may appear in this documentation:

(5-12)

Cross-reference to figure legend

The first number refers to the figure number, the second to the instrument part in the figure.

1

Instruction step

Carry out these steps in the sequence shown.

Method

Dialog text, parameter in the software

File ▶ New

Menu or menu item

[Next]

Button or **key**



WARNING

This symbol draws attention to a possible life-threatening hazard or risk of injury.



WARNING

This symbol draws attention to a possible hazard due to electrical current.



WARNING

This symbol draws attention to a possible hazard due to heat or hot instrument parts.



WARNING

This symbol draws attention to a possible biological hazard.



CAUTION

This symbol draws attention to possible damage to instruments or instrument parts.

**NOTE**

This symbol highlights additional information and tips.

1.3 Safety instructions

1.3.1 General notes on safety

**WARNING**

Operate this instrument only according to the information contained in this documentation.

This instrument left the factory in a flawless state in terms of technical safety. To maintain this state and ensure non-hazardous operation of the instrument, the following instructions must be observed carefully.

1.3.2 Electrical safety

The electrical safety when working with the instrument is ensured as part of the international standard IEC 61010.

**WARNING**

Only personnel qualified by Metrohm are authorized to carry out service work on electronic components.

**WARNING**

Never open the housing of the instrument. The instrument could be damaged by this. There is also a risk of serious injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.

Supply voltage

**WARNING**

An incorrect supply voltage can damage the instrument.

Only operate this instrument with a supply voltage specified for it (see rear panel of the instrument).



Protection against electrostatic charges



WARNING

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Do not fail to pull the power cord out of the power socket before you set up or disconnect electrical plug connections at the rear of the instrument.

1.3.3 Tubing and capillary connections



CAUTION

Leaks in tubing and capillary connections are a safety risk. Tighten all connections well by hand. Avoid applying excessive force to tubing connections. Damaged tubing ends lead to leakage. Appropriate tools can be used to loosen connections.

Check the connections regularly for leakage. If the instrument is used mainly in unattended operation, then weekly inspections are mandatory.

1.3.4 Flammable solvents and chemicals



WARNING

All relevant safety measures are to be observed when working with flammable solvents and chemicals.

- Set up the instrument in a well-ventilated location (e.g. fume cupboard).
- Keep all sources of flame far from the workplace.
- Clean up spilled liquids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.

1.3.5 Recycling and disposal



This product is covered by European Directive 2012/19/EU, WEEE – Waste Electrical and Electronic Equipment.

The correct disposal of your old instrument will help to prevent negative effects on the environment and public health.



More details about the disposal of your old instrument can be obtained from your local authorities, from waste disposal companies or from your local dealer.




Figure 2 Rear with stand

1 Connector plug	2 Stand clamping screw
3 Stand	4 Cable cut-out



3 Installation

3.1 Setting up the instrument

3.1.1 Packaging

The instrument is supplied in protective packaging together with the separately packed accessories. Keep this packaging, as only this ensures safe transportation of the instrument.

3.1.2 Checks

Immediately after receipt, check whether the shipment has arrived complete and without damage by comparing it with the delivery note.

3.1.3 Location

The instrument has been developed for operation indoors and may not be used in explosive environments.

Place the instrument in a location of the laboratory which is suitable for operation and free of vibrations and which provides protection against corrosive atmosphere and contamination by chemicals.

The instrument should be protected against excessive temperature fluctuations and direct sunlight.

3.2 Installing the magnetic stirrer

The 728 Magnetic Stirrer is fastened to a support rod with a diameter of 10 mm.

3.3 Connecting the magnetic stirrer

Power supply

The 728 Magnetic Stirrer needs 5 to 12 V DC voltage.

It can be operated directly at a Metrohm instrument with a stirrer connector.

4 Operation and maintenance

4.1 General notes

4.1.1 Care

The 728 Magnetic Stirrer requires appropriate care. Excess contamination of the instrument may result in functional disruptions and a reduction in the lifetime of the otherwise sturdy mechanics and electronics.

Spilled chemicals and solvents should be removed immediately. Above all, the plug connections on the rear of the instrument (in particular the connection socket) should be protected from contamination.



CAUTION

Although this is extensively prevented by design measures, the connection cable should be unplugged immediately if aggressive media have found their way into the interior of the instrument to prevent serious damage to the instrument electronics. In such cases, Metrohm Service must be informed.

4.1.2 Maintenance by Metrohm Service

Maintenance of the 728 Magnetic Stirrer is best carried out as part of annual service, which is performed by specialist personnel from Metrohm. If working frequently with caustic and corrosive chemicals, a shorter maintenance interval could be necessary.

Metrohm Service offers every form of technical advice for maintenance and service of all Metrohm instruments.



5 Technical specifications

5.1 Magnetic stirrer

Direction of rotation Counterclockwise

Rotational speed 200 - 1,900 rpm

5.2 Power supply

Power supply 5 - 12 V DC

Power consumption ≤ 2.5 VA

5.3 Interfaces and connectors

The 728 Magnetic Stirrer does not have any control interfaces. It is operated and controlled via the power supply.

5.4 Ambient temperature

Nominal function range +5 - +45 °C

Transport and storage -40 - +70 °C

5.5 Reference conditions

Ambient temperature +25 °C (± 3 °C)

Relative humidity $\leq 60\%$

5.6 Dimensions/material

<i>Width</i>	90 mm
<i>Height</i>	76 mm
<i>Depth</i>	150 mm
<i>Weight</i>	660 g (without accessories)
<i>Material</i>	
<i>Housing, upper part</i>	PBT
<i>Housing, bottom part</i>	Steel sheet, stove-enameled



Index

C

Connect
Magnetic stirrer 8

D

Dimensions 11

E

Electrostatic charge 4

I

Instrument description 1

M

Magnetic stirrer
Direction of rotation 10
Rotational speed 10
Material 11
Model versions 1

P

Power supply
Power consumption 10
Supply voltage 10

S

Safety instructions 3
Service 3
Supply voltage 3