

732 IC Detector

733 IC Separation Center

5.732.0012 Program



8.732.1033 Instructions for Use

Table of contents

1	Introduction	1
1.1	Instrument description	1
1.2	Parts and controls	3
1.2.1	732 IC Detector	3
1.2.2	733 IC Separation Center	5
1.3	Information on the Instructions for Use	8
1.3.1	Organization	8
1.3.2	Notation and pictograms	9
1.4	Safety notes	10
1.4.1	Electrical safety	10
1.4.2	General precautionary rules	10
2	Installation	11
2.1	Flow chart	11
2.2	Setting up the instruments	12
2.2.1	Packaging	12
2.2.2	Check	12
2.2.3	Location	12
2.2.4	Arrangement of the instruments	12
2.3	Connection of 733 IC Separation Center	13
2.3.1	733.0010/733.0X30 IC Separation Center	13
2.3.2	733.0X20 IC Separation Center	14
2.3.3	Connection of syringe and suction tubing	16
2.3.4	Connection of the drain tube	16
2.3.5	Connection of the 6.5324.000 Bottle rack (option)	16
2.4	Mains connection	17
2.4.1	Setting the mains voltage	17
2.4.2	Fuses	18
2.4.3	Mains cable and mains connection	18
2.4.4	On/off switching of the instruments	18
2.5	Capillary connections	19
2.5.1	Capillaries	19
2.5.2	Steel connectors	19
2.5.3	PEEK connectors	20
2.6	Connection of 709 IC Pump	21
2.6.1	Electrical connection	21
2.6.2	Pulsation dampener	21
2.6.3	Filter unit PEEK	22
2.6.4	Filter unit Manufit	23
2.6.5	Connection to injection valve with PEEK capillaries	24
2.6.6	Connection to injection valve with steel capillaries	25
2.6.7	Passivation of the IC system	27
2.7	Precolumns	28
2.7.1	General information on precolumns	28
2.7.2	Precolumns with twin cartridge holder	28
2.7.3	Precolumns with cartridge head	30
2.7.4	IC anion precolumn SUPERSEP	31

2.8	Separating columns and suppressor module	32
2.8.1	General information on separating columns	32
2.8.2	Selection of the sample loop	32
2.8.3	General information on suppressor module	33
2.8.4	One-channel system without suppressor module.....	35
2.8.5	Two-channel system without suppressor module	35
2.8.6	One-channel system with suppressor module	37
2.8.7	Leak testing and conditioning	40
2.9	Connection of external devices	41
2.9.1	Connection of a recorder	41
2.9.2	Connection of «IC Metrodata for Win95»	41
2.9.3	Connection of the 750 Autosampler	41
2.9.4	Connection of the 766 IC Sample Processor	41
2.9.5	Connection of the 791 VA Detector	41
2.9.6	Connection of a printer	42
2.9.7	Connection of a PC	44
2.9.8	Connection of devices to the remote interface	44
3	Operating tutorial	45
3.1	Requirements	45
3.2	Preparations	46
3.3	Putting into operation	48
3.4	Calibration	59
3.5	Sample determination	62
3.6	Storing as a method	64
4	Operation	67
4.1	Operating sequences	67
4.1.1	General flow chart	67
4.1.2	Flow chart for basic settings	68
4.1.3	Flow chart for putting into operation	69
4.1.4	Flow chart for injection	70
4.2	Fundamentals of the operation	71
4.2.1	Display	71
4.2.2	Overview of key functions	72
4.2.3	Instrument dialog	77
4.2.4	Data entry	79
4.2.5	Text entry	80
4.3	Displays in the standby mode	81
4.3.1	Measured value and current time	81
4.3.2	Status messages	82
4.4	Basic settings	84
4.4.1	Setup	84
4.4.2	Configuration, <CONFIG> key	88
4.5	Measurement parameters	97
4.5.1	<PARAM> key	97
4.5.2	<FULL SCALE> key	102

4.6	Triggering of functions	103
4.6.1	<FILL> keys	103
4.6.2	<INJECT> keys	104
4.6.3	<ZERO> key	105
4.6.4	<ZERO OFF> key	106
4.6.5	<MARK> key	106
4.6.6	<PUMP R/S> key	106
4.7	Programming	107
4.7.1	<PROGRAM> key	107
4.7.2	<PROG R/S> key	114
4.7.3	<EVENT> key	114
4.7.4	<METHOD> key	117
4.8	Data output	119
4.8.1	<PRINT> key	119
4.8.2	<PLOT> key	122
4.8.3	<REPORT> key	124
4.9	Examples of methods	126
4.9.1	Cation determination with Metrosep Cation 1-2	126
4.9.2	Anion determination with Metrosep Anion Dual 2	129

5 Notes – Maintenance – Faults **133**

5.1	Practical notes on ion chromatography	133
5.1.1	Separating columns	133
5.1.2	Pumps	134
5.1.3	Eluents	135
5.1.4	Suppressor module	136
5.1.5	Connections	136
5.2	Maintenance and servicing	136
5.2.1	General information	136
5.2.2	Passivation	137
5.2.3	Recycling	137
5.2.4	Shutdown	137
5.2.5	Changing separating columns	138
5.2.6	Regeneration of the suppressor	140
5.2.7	Cleaning the suppressor	141
5.2.8	Replacing the suppressor	143
5.3	Faults and malfunctions	145
5.3.1	Error messages	145
5.3.2	Malfunctions and their rectification	148
5.4	Diagnostic tests	149
5.4.1	General	149
5.4.2	Preparing instruments	150
5.4.3	Checking working memory (RAM)	151
5.4.4	Check keypad	151
5.4.5	Check display	152
5.4.6	Check RS232 interface	153
5.4.7	Check remote interface	154
5.4.8	Internal hardware test	155
5.4.9	Initialize data memory	157
5.5	Validation / GLP	159

6	Interfaces	161
6.1	RS232 interfaces	161
6.1.1	General rules for remote control	161
6.1.2	Call-up of objects	162
6.1.3	Triggers	163
6.1.4	Status messages	164
6.1.5	Error messages	165
6.1.6	Remote control commands	166
6.1.7	Data transmission protocol	177
6.1.8	Handshake	178
6.1.9	Pin assignment	181
6.1.10	RS232 error rectification	182
6.2	Remote interfaces	183
6.2.1	"Remote" interface	183
6.2.2	"733 IC Separation Center" interface	185
6.3	Analog output	187
6.4	External power supply for 733 IC Separation Center	187
6.5	Valve interfaces	188
7	Appendix	189
7.1	Technical data	189
7.1.1	732 IC Detector	189
7.1.2	733 IC Separation Center	192
7.2	Standard equipment	193
7.2.1	732 IC Detector	193
7.2.2	733 IC Separation Center	194
7.3	Optional accessories	196
7.3.1	Accessories for 733 IC Separation Center	196
7.3.2	Separating columns and precolumns	198
7.3.3	Additional devices and cables	202
7.4	Warranty and conformity	204
7.4.1	Warranty	204
7.4.2	EU Declaration of conformity	205
7.4.3	Certificate of conformity and system validation	207
7.5	Index	209

List of figures

<u>Fig. 1:</u>	Block diagram of the ion chromatography system	2
<u>Fig. 2:</u>	Front of the 732 IC Detector	3
<u>Fig. 3:</u>	Rear of the 732 IC Detector	4
<u>Fig. 4:</u>	Front of the 733 IC Separation Center.....	5
<u>Fig. 5:</u>	Rear of the 733 IC Separation Center.....	6
<u>Fig. 6:</u>	Connection 732 – 2.733.0010/2.733.0X30	13
<u>Fig. 7:</u>	Connection 732 – 2.733.0X20	15
<u>Fig. 8:</u>	Setting the mains voltage	18
<u>Fig. 9:</u>	Connectors for capillaries.....	20
<u>Fig. 10:</u>	Connection of 709 IC Pump	21
<u>Fig. 11:</u>	6.2821.100 Filter unit PEEK	22
<u>Fig. 12:</u>	6.2821.000 Filter unit Manufit.....	23
<u>Fig. 13:</u>	Connection to injection valve with PEEK capillaries.....	25
<u>Fig. 14:</u>	Connection to injection valve with steel capillaries	26
<u>Fig. 15:</u>	Installing precolumn cartridges	29
<u>Fig. 16:</u>	Interior of the 733.0010 IC Separation Center	34
<u>Fig. 17:</u>	Interior of the 733.0X30 IC Separation Center.....	36
<u>Fig. 18:</u>	Connections at suppressor module	38
<u>Fig. 19:</u>	Ion chromatogram of the calibration	61
<u>Fig. 20:</u>	Ion chromatogram of the drinking water sample	63
<u>Fig. 21:</u>	Schematic representation of the instrument dialog.....	78
<u>Fig. 22:</u>	Ion chromatogram for cation standard with Metrosep Cation 1-2	128
<u>Fig. 23:</u>	Ion chromatogram for anion standard with Metrosep Anion Dual 2.....	131
<u>Fig. 24:</u>	Assembling the suppressor.....	142

List of numbered parts and controls

1	Display	3	54	Housing for filter unit	22
2	Main function keys	3	55	Connector	22
3	Numeric keys	3	56	Inlet capillary	23,26
4	Auxiliary function keys	3	57	Manufit pressure screw	23,29
5	Program status display	3	58	Counterpart end	23
6	Auto-zero display	3	59	PTFE gasket	23,29
7	Thermostat display	3	60	4 Steel meshes	23,29
8	Overload display	3	61	Steel mesh holding end	23
9	Mains switch	4,18	62	Manufit housing	23,29
10	Serial number	4	63	Outlet capillary	23,26
11	Analog output 0... 1 V	4	64	Filter unit PEEK	25
12	Analog output 0... 10 mV	4	65	PEEK capillary	25,34,36
13	Connection for detector block	4	66	Pulsation dampener	25,34,36
14	Connection for 733 IC Separation Center	4	67	Column connection capillary	25,29,34,36
15	Connection for 709 IC Pump	4	68	Injection valve	25,34,36
16	RS232 interface	4	69	Filter unit Manufit	26
17	Remote interface	4	70	Coupling	26
18	Mains connection plug	4,18	71	Steel capillary	26
19	Fuse holder	4,18	72	Outlet capillary	29
20	Door to interior	5	73	2 Steel meshes	29
21	Connection for syringe	5	74	Precolumn cartridge	29
22	Feedthrough for aspirating tubing	5	75	Inlet capillary	29
23	„FILL“ for valve A	5	76	IC separating column	29,34,36
24	„INJECT“ for valve A	5	77	Manufit pressure screw	29
25	„FILL“ for valve B	5	78	Steel spacer	29
26	„INJECT“ for valve B	5	79	Steel connector for ferrule	29
27	Connection/Feedthrough	5	80	Manufit housing	29
28	Feedthrough	5	81	Detector block	34,36
29	Opening for detector cable B	6	82	Inlet capillary for detector block	34,36
30	Opening for outlet capillary B	6	83	Mounting rail	34,36
31	Rear panel opening	6	84	Column holder	34,36
32	Opening for outlet capillary A	6	85	Capillary for syringe	34,36
33	Opening for detector cable A	6	86	Sample loop	34,36
34	Knurled screw	6	87	Inlet capillary for injector	34,36
35	Detachable rear panel	6	88	PTFE aspirating tubing	34,36
36	Connection for 732 IC Detector A	6	89	Suppressor inlet capillary for eluent	36,38
37	Knurled screw	6	90	Suppressor inlet capillary for H ₂ SO ₄	36,38
38	Terminal block for valve A	6	91	Suppressor outlet capillary for H ₂ SO ₄	36,38
39	Detachable rear panel	6	92	Suppressor outlet capillary for H ₂ O	36,38
40	Opening for inlet capillary A	6	93	Suppressor inlet capillary for H ₂ O	36,38
41	Rear panel opening	6	94	Suppressor outlet capillary for eluent	36,38
42	Connection for drain tube	6	95	Suppressor module	36
43	Rear panel opening	6,36	96	Coupling	36
44	Opening for inlet capillary B	6	97	Screw nut	142
45	Connection for 732 IC Detector B	6	98	Connection piece	142
46	Terminal block for valve B	6	99	Suppressor rotor	142
47	Model plate	6	100	Suppressor holder	142
48	Connection for external supply	6			
49	Ferrule	20,29			
50	Pressure screw	20,29			
51	Capillary	20			
52	Compression fitting	20			
53	Connector with filter	22			

1 Introduction

1.1 Instrument description

The **732 IC Detector** is a conductivity detector especially designed for ion chromatography with an extensive operating range and high sensitivity for the recording of chromatograms with and without chemical suppression. The associated thermostatable detector block is normally installed in the 733 IC Separation Center, but can also be used as a separate detector. The two following versions are available:

- **2.732.0010** IC Detector with standard detector block
- **2.732.0110** IC Detector with metal-free detector block

The 732 IC Detector is operated using the keypad with operator guidance via the two-line LCD. In addition to setting of the measurement parameters, time programs can be generated which can be used to initiate a large number of instrument functions for each of the maximum 20 program steps. Further, the same functions can be executed at a specific time with 4 programmable "events".

The 732 IC Detector is equipped with various interfaces for communication purposes. Recorders, integrators or the «IC Metrodata for Win95» chromatography data system can be connected to the analog output (1 V or 10 mV) for the plotting and evaluation of chromatograms. The two RS232 interfaces are used for the connection of a 709 IC Pump, a printer or a PC for remote control of the IC system. Finally, programmable signals at a "remote" interface can be employed to control any external devices which, in turn, can start functions at the IC system.

The **733 IC Separation Center** is a thermally and electronically isolated wet part which accommodates injectors, columns, detectors, suppressor module and pulsation dampener and is controlled by the 732 IC Detector. The following versions are available:

- **2.733.0010** IC wet part with 1 injector for a one-channel system with electronic suppression
- **2.733.0020** IC wet part with 2 injectors for a two-channel system with electronic suppression
- **2.733.0120** IC wet part with 2 injectors for a two-channel system with electronic suppression, metal-free
- **2.733.0030** IC wet part with 1 injector and 1 Metrohm Suppressor Module MSM for a one-channel system with chemical suppression
- **2.733.0130** IC wet part with 1 injector and 1 Metrohm Suppressor Module MSM for a one-channel system with chemical suppression, metal-free

The 732 IC Detector and 733 IC Separation Center are the main components of a modular **ion chromatography system** that can be expanded to meet the wishes of the individual user (see Fig. 1). The minimum configuration of the one-channel system also includes a 709 IC Pump, a separating column and a recorder. The two-channel system requires at least a second 732 IC Detector and a second 709 IC Pump. Printers, integrators, data recording devices, PC and autosamplers can be attached to both systems. Further, practically all HPLC peripherals and parts available on the market such as precolumns, additional separating columns, additional detectors and other injection systems can be seamlessly integrated in the system.

However, the individual IC units can also be freely combined with common HPLC instruments. This offers the possibility of expanding your system to a standalone ion chromatograph.




Fig. 1: Block diagram of the ion chromatography system

C	Separating column	M	Suppressor module (one-channel system only)	709	IC Pump
D	Detector	PC	PC	732	IC Detector
E	Eluent	PL	Recorder	733	IC Separation Center
I	Injector	PR	Printer	750	Autosampler
IF	Interface	S	Sample	752	Pump Unit
IN	Integrator	W	Waste	754	Dialysis Unit
				766	IC Sample Processor

1.2 Parts and controls

1.2.1 732 IC Detector




Fig. 2: Front of the 732 IC Detector

<p>1 Display (LCD) comprising 2 lines each of 24 characters</p>	<p>5 Program status display (LED) LED dark: Program inactive LED lit up: Program ready LED flashes: Program running</p>
<p>2 Main function keys Auto-zero, Parameters, Full Scale, Program start/stop</p>	<p>6 Auto-zero display (LED) LED dark: Auto-zero switched off LED lit up: Auto-zero switched on</p>
<p>3 Numeric keys Edit mode: Numeric keys Basic mode: Function keys</p>	<p>7 Thermostat display (LED) LED dark: Heating switched off LED lit up: Heating switched on</p>
<p>4 Auxiliary function keys Select, Clear, Quit, Enter</p>	<p>8 Overload display (LED) LED lit up: Meas. signal > 150% of full-scale range LED flashes: Meas. signal > 180% of full-scale range</p>




Fig. 3: Rear of the 732 IC Detector

<p>9 Mains switch switch to switch instrument on and off: I = ON 0 = OFF</p>	<p>15 Connection for 709 IC Pump RS232 interface, can also be used for connection of a printer</p>
<p>10 Serial number</p>	<p>16 RS232 interface connection of a printer, PC, etc.</p>
<p>11 Analog output 0¼ 1 V red socket: live black socket: common</p>	<p>17 Remote interface remote I/O lines for connection of external devices</p>
<p>12 Analog output 0¼ 10 mV red socket: live black socket: common</p>	<p>18 Mains connection plug mains connection, see <i>section 2.4</i></p>
<p>13 Connection for detector block</p>	<p>19 Fuse holder changing the fuses, see <i>section 2.4</i></p>
<p>14 Connection for 733 IC Separation Center</p>	

1.2.2 733 IC Separation Center




Fig. 4: Front of the 733 IC Separation Center

20	Door to interior	25	"FILL" key for valve B or "STEP" for suppressor module
21	Connection for 6.2816.020 Syringe	26	"INJECT" key for valve B
22	Feedthrough for aspirating tubing	27	733.0010: Feedthrough for capillary 733.0X20: Connection for 6.2816.020 Syringe 733.0X30: Feedthrough for suppressor inlet cap.
23	"FILL" key for valve A	28	733.0010: Feedthrough 733.0X20: Feedthrough for aspirating tubing 733.0X30: Feedthrough for suppressor inlet capillary
24	"INJECT" key for valve A		




Fig. 5: Rear of the 733 IC Separation Center

29 Opening for detector cable B opening for connecting cable detector block B – 732	39 Detachable rear panel access to bottom part of inner compartment
30 Opening for outlet capillary B discharge of the eluent of column B to waste	40 Opening for inlet capillary A supply of the eluent for column A
31 Rear panel opening (closed with plastic stopper) for additional supply and discharge lines to and from the inner compartment	41 Rear panel opening (closed with plastic stopper) for additional supply and discharge lines to and from the inner compartment
32 Opening for outlet capillary A discharge of the eluent of column A to waste	42 Connection for drain tube for discharge of spilled liquid from the inner compartment
33 Opening for detector cable A opening for connecting cable detector block A – 732	43 Rear panel opening (closed with plastic stopper) for additional supply and discharge lines to and from the inner compartment
34 Knurled screw for fastening the rear panel 35	44 Opening for inlet capillary B supply of the eluent for column B
35 Detachable rear panel access to top part of the inner compartment	45 Connection for 732 IC Detector B
36 Connection for 732 IC Detector A	46 Terminal block for valve B Ground, Fill, Inject: inputs for control of the valve Pos.Fill: output signal on switching of the valve to position "FILL" Integr.Start: output signal on switching of the valve to position "INJECT"
37 Knurled screw for fastening rear panel 39	47 Model plate with serial number
38 Terminal block for valve A Ground, Fill, Inject: inputs for control of the valve Pos.Fill: output signal on switching of the valve to position "FILL" Integr.Start: output signal on switching of the valve to position "INJECT"	48 Connection for external supply connection of power supply unit (5 V, 0.5 A / 24 V, 2 A) in operation without 732 IC Detector

1.3 Information on the Instructions for Use



Please read through these Instructions for Use carefully before you put the 732 IC Detector and 733 IC Separation Center into operation. The Instructions for Use contain information and warnings which must be heeded the user to assure safe operation of the instruments.

1.3.1 Organization

These **8.732.1033 Instructions for Use** for the 732 IC Detector and 733 IC Separation Center provide a comprehensive overview of the installation, startup procedure, operation, fault rectification and technical specifications of these instruments. The Instructions for Use are organized as follows:





- Section 1 Introduction**
General description of instruments, parts and controls and safety notes
- Section 2 Installation**
Installation of 732 IC Detector / 733 IC Separation Center, attachment of accessories and external devices
- Section 3 Operating tutorial**
Introduction to the operation using an example
- Section 4 Operation**
Detailed description of the operation and explanation of functions of all keys
- Section 5 Notes – Maintenance – Faults**
Notes on ion chromatography, maintenance, fault rectification, diagnostic tests, validation
- Section 6 Interfaces (green pages)**
Description of RS232 interfaces, remote interfaces, valve interface and analog output
- Section 7 Appendix**
Technical data, standard equipment, options, warranty, declarations of conformity, index

To find the required information on the instruments, you will find it an advantage to use either the **Table of contents** or the **Index** at the back. The **8.732.1043 Quick Reference Guide** is suitable for use as a reference work for daily use as it explains the most important parameters and key functions.

As a supplement to the Instructions for Use, the **Metrohm Monograph 8.732.2003 "Ion chromatography"** is also supplied. This provides an introduction to the theoretical fundamentals and general information on separating columns and sample pretreatment. You will find detailed information on the separating columns available from Metrohm and on special IC applications in the relevant **"Application Bulletins"**, which are available on request free of charge from your Metrohm agency.

1.3.2 Notation and pictograms

The following notations and pictograms (symbols) are used in these Instructions for Use:

<PARAM>	Key
"Range"	Parameter or entry value
35	Part or control of 732/733
<u>22</u>	Part or control of 709
>PARAM/detector range: 1.00 mS/cm	Display Text in display 1 of the 732 IC Detector
	Hazard This symbol draws attention to a possible danger to life or of injury if the associated directions are not followed correctly.
	Warning This symbol draws attention to possible damage to instruments or instrument parts if the associated directions are not followed correctly.
	Caution This symbol marks important information. First read the associated directions before you continue.
	Comment This symbol marks additional information and tips.

1.4 Safety notes

1.4.1 Electrical safety

While electrical safety in the handling of the 732 IC Detector and 733 Separation Center is assured in the context of the specifications IEC 1010-1 (protection class 1, degree of protection IP40), the following points should be noted:

- **Mains connection**



*Setting of the **mains voltage**, checking the **mains fuse** and the **mains connection** must be effected in accordance with the instructions in section 2.4.*

- **Opening the 732 IC Detector**



*If the 732 IC Detector is connected to the power supply, the instrument must not be opened nor must parts be removed from it, otherwise there is a danger of coming into contact with components which are live. Hence, always disconnect the instrument from all voltage sources before you open it and ensure that the **mains cable is disconnected from mains connection 18** !*

- **Opening the 733 IC Separation Center**



*Disconnect **connecting cable to the 732 IC Detector from connector 14** before you remove the middle housing panel with connectors.*

- **Protection against static charges**



Electronic components are sensitive to static charging and can be destroyed by discharges. Before you touch any of the components inside the 732 IC Detector or 733 IC Separation Center, you should earth yourself and any tools you are using by touching an earthed object (e.g. housing of the instrument or a radiator) to eliminate any static charges which exist.

1.4.2 General precautionary rules

- **Handling of solvents**




Check all lines of the IC system periodically for possible leaks. Follow the relevant instructions regarding the handling of flammable and/or toxic solvents and their disposal.

2 Installation

2.1 Flow chart

The following flow chart provides an overview of all installation work. You will find more detailed information in the relevant sections.



2.2 Setting up the instruments

2.2.1 Packaging

The 732 IC Detector and 733 IC Separation Center are supplied together with the separately packed accessories in special packagings containing shock-absorbing foam linings designed to provide excellent protection. The actual instruments are packed in an evacuated polyethylene bag to prevent the ingress of dust. Please store all these special packagings as only they assure transport of the instruments free from damage.

2.2.2 Check

After receipt, immediately check whether the shipment is complete and has arrived without damage (compare with delivery note and list of accessories in *section 7.2*). In the case of transport damage, see instructions in *section 7.4.1 "Warranty"*.

2.2.3 Location

Position the instruments in the laboratory at a location convenient for operation, free from vibrations and protected against a corrosive atmosphere and contamination by chemicals. The same applies to all other components of the IC system.



To avoid disturbing temperature influences on the insulated column compartment, the entire system including pump and eluent reservoir must be protected against direct sunlight.

2.2.4 Arrangement of the instruments

In one-channel operation, the 709 IC Pump, 733 IC Separation Center and 732 IC Detector are best stacked on top of one another in this order.

In two-channel operation (2.733.0X20 IC Separation Center), the optimum arrangement (1, 2 or 3 towers) depends on the laboratory space available. However, the 709 IC Pumps should be set up at the very bottom and the 732 IC Detectors at the very top.



To ensure that the arrangement of pumps and detectors for the two channels A and B is clearly apparent in two-channel operation, it is advantageous to mark the instruments. The 6.2248.000 Magnetic plate is enclosed with the 732 IC Detector for this purpose. It can be cut to the desired size, labeled (e.g. with "A" or "B") and affixed to the appropriate instrument.

2.3 Connection of 733 IC Separation Center

2.3.1 2.733.0010/2.733.0X30 IC Separation Center

The instrument versions 2.733.0010 and 2.733.0030 of the IC Separation Center are operated with a 732 IC Detector whose standard equipment also includes the **1.732.0100 Detector block**. For the metal-free 2.733.0130 instrument version the **metal-free 1.732.0110 Detector block** must be used. It is best to proceed as follows when connecting the two instruments and the detector block:

1 Install detector block

- Unscrew the four knurled screws **34** from the top rear panel **35** of the 733 IC Separation Center and remove rear panel (see *Fig. 5*).
- Position detector block from the back in the space provided in the 733 IC Separation Center on the right and push fully to the front (see *Fig. 16*).
- Insert the cable permanently attached to the detector block in opening **33** and the outlet capillary in opening **32** "Waste A" of the rear panel **35**.
- Replace rear panel **35** and screw to the 733 IC Separation Center using the four knurled screws **34**.

2 Connect detector block

- Plug the gray connecting cable permanently attached to the detector block into connection **13** "Detector Block" of the 732 IC Detector and fasten to the instrument by tightening the screws in the cable connector (see *Fig. 6*).




Fig. 6: Connection 732 – 2.733.0010/2.733.0X30

3 Connect waste container

- Lead the outlet capillary of the detector block to a sufficiently large waste container and fix in place.

4 Connect 732 to 733

- Plug one end of the 6.2125.090 Connecting cable into connection **14** "733 IC Separation Center" of the 732 IC Detector and fasten to the instrument by tightening the screws in the cable connector (see *Fig. 6*).
- Plug the other end of the 6.2125.090 Connecting cable into connection **36** "732 IC Detector" of the 733 IC Separation Center and fasten to the instrument by tightening the screws in the cable connector (see *Fig. 6*).

2.3.2 2.733.0X20 IC Separation Center

For operation of the 2.733.0020 or the metal-free 2.733.0120 instrument version of the 733 Separation Center, two 1.732.0100 or metal-free 1.732.0110 IC Detectors are needed. It is best to proceed as follows when connecting the instruments and the two detector blocks:

1 Install detector blocks

- Unscrew the four knurled screws **34** from the top rear panel **35** of the 733 IC Separation Center and remove rear panel (see *Fig. 5*).
- Position first detector block A from the back in the space provided in the 733 IC Separation Center on the right and push fully to the front (see *Fig. 16*).
- Position second detector block B from the back in the space provided in the 733 IC Separation Center on the left and push fully to the front (see *Fig. 16*).
- Insert the cable permanently attached to detector block A in opening **33** and the outlet capillary in opening **32** "Waste A" of the rear panel **35**.
- Insert the cable permanently attached to detector block B in opening **29** and the outlet capillary in opening **30** "Waste B" of the rear panel **35**.
- Replace rear panel **35** and fasten to 733 IC Separation Center with the four knurled screws **34**.

2 Connect detector blocks

- Plug the gray connecting cable permanently attached to detector block A into connection **13** "Detector block" of the first 732 IC Detector and fasten to the instrument by tightening the screws in the cable connector (see *Fig. 7*).
- Plug the gray connecting cable permanently attached to detector block B into connection **13** "Detector Block" of the second 732 IC Detector Block and fasten to the instrument by tightening the screws in the cable connector.

3 Connect waste container

- Lead the outlet capillaries of the two detector blocks to a sufficiently large waste container and fix in place.

4 Connect 732 to 733

- Plug one end of the 6.2125.090 Connecting cable into connection **14** "733 IC Separation Center" of the **first 732 IC Detector** and fasten to the instrument by tightening the screws in the cable connector (see *Fig. 7*).
- Plug other end of the 6.2125.090 Connecting cable into connection **36** "732 IC Detector" of the 733 IC Separation Center and fasten to the instrument by tightening the screws in the cable connector.
- Plug one end of the 6.2125.090 Connecting cable into connection **14** "733 IC Separation Center" of the **second 732 IC Detector** and fasten to the instrument by tightening the screws in the cable connector.
- Plug the other end of the 6.2125.090 Connecting cable into connection **45** "732 IC Detector" of the 733 IC Separation Center and fasten to the instrument by tightening the screws in the cable connector.




Fig. 7: Connection 732 – 2.733.0X20

2.3.3 Connection of syringe and aspirating tubing

For manual filling of the sample loops mounted on the injection valves, the 6.2816.020 Syringe and the PTFE aspirating tubing already screwed to the valve are needed. These accessories are mounted or adjusted as follows:

1 Connect syringe

- Push 6.2816.020 Syringe (without needle) as far as it will go into connection socket **21** (for valve A) or **27** (for valve B) (see *Fig. 4*).

2 Adjust aspirating tubing

- Loosen the rotary nipple screwed onto the interior side of connection **22** or **28**.
- Pull PTFE aspirating tubing **88** (see *Fig. 16* and *Fig. 17*) by hand out of connection **22** or **28** as far as desired.
- Retighten rotary nipple on the interior side of connection **22** or **28** to fix the aspirating tubing in place.

2.3.4 Connection of the drain tube

The 733 IC Separation Center has a connection at the rear to which a drain tube for discharged liquids can be attached. Proceed as follows:

1 Connect drain tube

- Mount 6.1816.00 Silicone tubing on connection nipple **42** (see *Fig. 5*).

2 Lead drain tube to collecting vessel

- Lead the other end of the drain tube to a suitable collecting vessel and fix in place.

2.3.5 Connection of the 6.5324.000 Bottle rack (option)

The optional available 6.5324.000 Bottle rack for supply vessels can be placed on top of the IC system tower. The accessories include the supply vessels for eluent (2 L), regeneration solution (1 L) and rinsing solution (1 L). For the connection of the supply capillaries leading to the 709 IC Pump and the suppressor module, see the instructions given on the enclosed leaflet.

2.4 Mains connection



Follow the instructions below for connecting to the power supply. If the instrument is operated with a mains voltage set wrongly and/or wrong mains fuse, there is a danger of fire!

2.4.1 Setting the mains voltage

Before switching on the 732 IC Detector for the first time, check that the mains voltage set on the instrument (see *Fig. 8*) matches the local mains voltage. If this is not the case, you must reset the mains voltage on the instrument as follows:

1 Disconnect mains cable

Disconnect mains cable from mains connection plug **18** of the 732 IC Detector.

2 Remove fuse holder

Using a screwdriver, loosen fuse holder **19** below the mains connection plug **18** and take out completely.

3 Check and change fuse if necessary

Carefully take the fuse installed for the desired mains voltage out of fuse holder **19** and check its specifications (the position of the fuse in the fuse holder is marked by the white arrow imprinted next to the mains voltage range):

100¼120 V 0.63 A (slow-blow) Metrohm No. U.600.0014

220¼240 V 0.315 A (slow-blow) Metrohm No. U.600.0011

4 Insert fuse

Change fuse if necessary and reinsert in fuse holder **19**.

5 Install fuse holder

Depending on the desired mains voltage, insert fuse holder **19** in the 732 IC Detector so that the corresponding mains voltage range can be read normally and the adjacent white arrow points to the white bar imprinted below the fuse holder (see *Fig. 8*).




Fig. 8: Setting the mains voltage

2.4.2 Fuses

One of the two fuses 0.63 A/slow-blow for 100...120 V or 0.315 A/slow-blow for 220...240 V is installed in fuse holder **19** of the 732 IC Detector as standard.



Ensure that the instrument is never put into operation with fuses of another type, otherwise there is danger of fire!

For checking or changing fuses, process as described in *section 2.4.1*.

2.4.3 Mains cable and mains connection

Mains cable

The instrument is supplied with one of three mains cables

- 6.2122.020 with plug SEV 12 (Switzerland, ...)
- 6.2122.040 with plug CEE(7), VII (Germany, ...)
- 6.2133.070 with plug NEMA 5-15 (USA, ...)

which are three-cored and fitted with a plug with an earthing pin. If a different plug has to be fitted, the yellow/green lead (IEC standard) must be connected to protective earth (protection class 1).



Any break in the earthing inside or outside the instrument can make it a hazard!

Mains connection

Plug the mains cable into mains connection plug **18** of the 732 IC Detector (see *Fig. 8*).

2.4.4 On/off switching of the instruments

The 732 IC Detector is switched on and off using mains switch **9**. When the instrument is switched on, display **1** lights up.

2.5 Capillary connections

2.5.1 Capillaries

Some of the connections under high pressure between the feed pump and the detector block must be set up by the user. For **metal-free systems**, the **6.1831.010 PEEK capillary** (i.d. = 0.25 mm, e.d. = 1/16", length = 3 m) must be used.

For **non metal-free systems**, the 6.1831.010 PEEK capillary can also be used in the pressure range of 0...25 MPa (0...250 bar); on the other hand, in the pressure range of 25...50 MPa (250...500 bar), which is permissible only together with the non metal-free 2.709.0010 version of the 709 IC Pump, the **6.2620.020 Steel capillary** (i.d. = 0.25 mm, e.d. = 1/16", length = 3 m) must be used.

With PEEK capillaries the connection is made preferably with the 6.2744.010 PEEK connectors, with steel capillaries, the connection is made preferably with 6.2620.000 and 6.2620.010 steel connectors (see section 2.5.2 and section 2.5.3).



Capillaries fitted with new connectors must have a perfectly flat cut surface. For PEEK capillaries it is best to use the **6.2621.080 Capillary tubing cutter**, for steel capillaries the **6.2621.040 Capillary tubing cutter** (both available as an option).

2.5.2 Steel connectors

For the connection of steel capillaries, the steel connectors **6.2620.010 Ferrule** and **6.2620.000 Pressure screw** available as an option can be used. Proceed as follows:

1 Mount connectors

Slide a pressure screw **50** (6.2620.000) and a ferrule **49** (6.2620.010) over the end of the capillary **51** to be fastened as shown in Fig. 9.

2 Insert capillary in connection

Push capillary end into the corresponding connection as far as it will go (to avoid dead volume).

3 Tighten compression fitting

Tighten pressure screw **50** with the open-end spanner 1/4" (6.2621.050) supplied.



For the connection of capillaries to the injection valves, use only the special steel connections contained in a plastic bag affixed to the valve (or as an alternative the 6.2744.010 PEEK compression fittings). If other steel connectors are used (e.g. 6.2620.000 and 6.2620.010), the valve connection may be damaged!




Fig. 9: Connectors for capillaries

49	Ferrule (6.2620.010)
50	Pressure screw (6.2620.000)
51	Capillary 6.2620.020 Steel capillary or 6.1831.010 PEEK capillary
52	Compression fitting (6.2744.010)

2.5.3 PEEK connectors

For the connection of 6.1822.010 PEEK capillaries (i.D. = 0.3 mm) or 6.1822.010 PTFE microcapillaries, the **6.2744.010 PEEK compression fittings** are used. Proceed as follows:

1 Mount compression fitting

Slide a compression fitting **52** (6.2744.010) over the end of the capillary **51** to be fastened as shown in *Fig. 9*.

2 Insert capillary in connection

Push capillary end in the corresponding connection as far as it will go (to avoid dead volume).

3 Tighten compression fitting

Tighten compression fitting **52** by hand (never use tools).

2.6 Connection of 709 IC Pump

2.6.1 Electrical connection

For operation of the 732 IC Detector and 733 IC Separation Center you can use any commercial HPLC pump. However, as the attainable sensitivity depends to a large extent on the quality of the pump, Metrohm advises use of the **709 IC Pump**, which has been specially developed for the demands of ion chromatography and has minimal pump pulsation and an outstanding flow constancy.

Startup and operation of the 709 IC Pump are described in the *709 Instructions for Use*. The eluent, which must be degassed and filtered (cf. *section 5.1.3*), is selected on the basis of the separating column installed in the 733 IC Separation Center and the current separation problem (see 8.732.2003 IC Monograph).

The **connection of the 709 IC Pump** at the connection **15** of the 732 IC Detector is shown in *Fig. 10*. For this you can use the 6.2125.060 Cable available as an option or another RS cable specified as a "null modem" cable. To ensure proper functioning of the communication between the 732 IC Detector and 709 IC Pump, the sliding switch **36** on the IC pump must be set to "RS 232" and the external control switched on with key **8** <EXT.> (see *709 Instructions for Use*).




Fig. 10: Connection of 709 IC Pump

2.6.2 Pulsation dampener

To protect the column material against pressure drops caused by the injector, the use of a pulsation dampener connected between the pump and the injection valve of the 733 IC Separation Center is recommended. The optional **6.2620.150 Pulsation dampener MF** (see *section 7.3.1*) is very well suited to this purpose.

The metal-free 6.2620.150 Pulsation Dampener is supplied fully assembled and has two connections for capillaries, for which either the connectors supplied or two 6.2744.010 PEEK compression fittings can be used. The flow direction is arbitrary. The pulsation dampener is positioned in the interior of the 733 IC Separation Center on the base below the injection valve (see *Fig. 16 and Fig. 17*).

2.6.3 Filter unit PEEK

The **6.2821.100 Filter unit PEEK** (see Fig. 11) supplied with the 709 IC Pump is installed between the 709 IC Pump and the injection valve at the 733 IC Separation Center. This filter unit serves to avoid contamination of the piston seals of the 709 IC Pump by abrasive particles and can be used in the pressure range 0...25 MPa (0...250 bar).

The two filter units PEEK supplied with the 2.733.0X30 IC Separation Center (with suppressor) are installed between the 752 Pump Unit and the inlet capillaries for regeneration and rinsing solution. These filter units serve to protect the suppressor module from foreign particles and bacterial growth.

The 6.2821.100 Filter unit PEEK consists of the housing **54** and the two connectors **53** (with filter) and **55** (without filter) screwed into the housing **54**. For the connection of capillaries **51** PEEK compression fittings **52** (6.2744.010) must be used. New connectors **53** with filter are available as an option with the ordering number 6.2821.110 (set of 10).



For the connection of the filter unit, please note the flow direction arrow printed on the housing.




Fig. 11: 6.2821.100 Filter unit PEEK

51 Capillary
6.1831.010 PEEK capillary

54 Housing for filter unit
Part of 6.2824.100 Filter unit

52 Compression fitting (6.2744.010)

55 Connector without filter
Part of 6.2824.100 Filter unit

53 Connector with filter (6.2824.110)
Part of 6.2824.100 Filter unit

2.6.4 Filter unit Manufit

The optional **6.2821.000 Filter unit Manufit** (see *section 7.3.1*) is installed between the 709 IC Pump and the injection valve at the 733 IC Separation Center. The filter unit serves to avoid contamination of the piston seals by abrasive particles and can be used in the pressure range 0...50 MPa (0...500 bar) for non metal-free systems together with steel capillaries. It is installed as follows (see *Fig. 12*):

1 Prepare Manufit housing

- Insert outlet capillary **63** with steel mesh holding end **61** into Manufit housing **62**.
- Insert the 4 steel meshes **60** provided into the steel mesh holding end **61**.
- Press the PTFE gasket **59** into the steel mesh holding end **61**.

2 Prepare Manufit pressure screw


- Insert inlet capillary **56** with counterpart end **58** into Manufit pressure screw **57**.

3 Assembly

- Fit the two capillary end pieces **58** and **61** together.
- Screw Manufit pressure screw **57** and Manufit housing **62** firmly together.



To replace contaminated steel meshes, proceed in the reverse order.



56	Inlet capillary
57	Manufit pressure screw
58	Counterpart end
59	PTFE Gasket (6.2821.010)
60	Steel meshes (6.2821.020)
61	Steel mesh holding end
62	Manufit housing
63	Outlet capillary

Fig. 12: 6.2821.000 Filter unit Manufit

2.6.5 Connection to injection valve with PEEK capillaries

For metal-free systems and for non metal-free systems in the pressure range 0...25 MPa (0...250 bar) we recommend to use 6.1831.010 PEEK capillaries, a 6.2620.150 Pulsation dampener (see *section 2.6.2*) and a 6.2821.100 Filter unit PEEK (see *section 2.6.3*) to connect the 709 IC Pump and the injection valve of the 733 IC Separation Center. Proceed as follows:



For the connection of capillaries to the injection valve, use only the 6.2744.010 PEEK compression fittings. If other steel connectors are used (e.g. 6.2620.000 and 6.2620.010), the valve connection may be damaged!

1 Connection to 709 IC Pump

- Cut connection capillary **22** (6.1831.010 PEEK capillary) to the required length and equip with the necessary connectors.
- Attach one end of connection capillary **22** to connection **23** of the 709 IC Pump (see *709 Instructions for Use*).
- Attach the other end of connection capillary **22** to the connector **53** (with filter) of the filter unit **64** (see *Fig. 13*).
- Attach the PEEK capillary **65** cut to the required length and equipped with connectors to the connector **55** of the filter unit **64**.

2 Installation of the capillary in the IC Separation Center

- Unscrew the four knurled screws **37** of the bottom rear panel **39** of the 733 IC Separation Center and remove rear panel (see *Fig. 5*).
- Insert PEEK capillary **65** from the back into the inner compartment of the IC Separation Center.
- Install rear panel **39** so that the capillary is positioned in opening **40** „Inlet A“ or **44** „Inlet B“ and screw on with the four knurled screws **37**.

3 Connection to injection valve

- Connect PEEK capillary **65** to pulsation dampener **66** (see *section 2.6.2*). Using another PEEK capillary **65**, connect pulsation dampener **67** to connection "5" of injection valve **68**.

4 Mount column connection capillary

- Connect column connection capillary **67** (ca. 20 cm of 6.1831.010 PEEK capillary) to connection "4" of injection valve **68**.




Fig. 13: Connection to injection valve with PEEK capillaries

22	Connection capillary 6.1831.010 PEEK capillary	66	Pulsation dampener (6.2620.150)
64	Filter unit PEEK (6.2824.100)	67	Column connection capillary 6.1831.010 PEEK capillary
65	PEEK capillary (6.1831.010)	68	Injection valve

2.6.6 Connection to injection valve with steel capillaries

For the pressure range 25...50 mPa (250...500 bar) we recommend to use 6.2620.020 Steel capillaries, a 6.2620.150 Pulsation dampener (see section 2.6.2) and a 6.2821.000 Filter unit Manufit (see section 2.6.4) to connect the 709 IC Pump and the injection valve of the 733 IC Separation Center. Proceed as follows:



For the connection of capillaries to the injection valve, use only the special steel connectors contained in a plastic bag affixed to the valve. If other steel connectors are used (e.g. 6.2620.000 and 6.2620.010), the valve connection may be damaged!

1 Connection to 709 IC Pump

- Attach inlet capillary **56** of the Manufit filter unit **69** to connection **23** of the 709 IC Pump (see *709 Instructions for Use*).
- Attach outlet capillary **63** of the filter unit Manufit **69** using coupling **70** to a steel capillary **71** cut to the required length (see *Fig. 14*).

2 Installation of the capillary in the IC Separation Center

- Unscrew the four knurled screws **37** of the bottom rear panel **39** of the 733 IC Separation Center and remove rear panel (see *Fig. 5*).
- Insert steel capillary **71** from the back into the inner compartment of the IC Separation Center.
- Install rear panel **39** so that the capillary is positioned in opening **40** „Inlet A“ or **44** „Inlet B“ and screw on with the four knurled screws **37**.

3 Connection to injection valve

- Connect steel capillary **71** to pulsation dampener **66** (see *section 2.6.2*). Using another steel capillary **71**, connect pulsation dampener **66** to connection "5" of injection valve **68**.

4 Mount column connection capillary

- Connect column connection capillary **67** (ca. 20 cm of 6.1831.010 PEEK capillary) to connection "4" of injection valve **68**.




Fig. 14: Connection to injection valve with steel capillaries

56 Inlet capillary of Filter unit Manufit 69	68 Injection valve
63 Outlet capillary of Filter unit Manufit 69	69 Filter unit Manufit (6.2821.000)
66 Pulsation dampener (6.2620.150)	70 Coupling (6.2620.060)
67 Column connection capillary 6.2620.020 Steel capillary	71 Steel capillary (6.2620.020)

2.6.7 Passivation of the IC system

With non metal-free IC systems, the entire IC system (without precolumn, separating column and suppressor module) must be passivated with nitric acid before being put into operation for the first time. Proceed as follows:

1 Connect detector block to injection valve

- Connect column connection capillary **67** using a 6.2620.060 Coupling directly to inlet capillary **82** of the detector block (see *Fig. 16*).

2 Set injection valve to "INJECT"

- Switch on 732 IC Detector using mains switch **9**.
- Press key **24** or **26** <INJECT> on 733 IC Separation Center. The green LED in the key should light up to show that the injection valve is in the position "INJECT".

3 Rinse with HNO₃

- Immerse aspirating capillary of the 709 IC Pump in $c(\text{HNO}_3) = 0.2 \text{ mol/L}$.
- Set a flow rate of 1 mL/min on the 709 IC Pump.
- Switch on 709 IC Pump and rinse IC system for ca. 10 min. During this time, check all capillaries and their connections between the 709 IC Pump and the detector block for escaping liquid. Should liquid escape at any position, the appropriate compression fitting must be tightened more or changed.
- Switch off 709 IC Pump.

4 Rinse with dist. H₂O

- Immerse aspirating capillary of the 709 IC Pump in distilled or demineralised water.
- Switch on 709 IC Pump and rinse the IC system for ca. 10 min.
- Switch off 709 IC Pump.

5 Rinse with eluent

- Immerse aspirating capillary of the 709 IC Pump in the eluent which will be needed for the separating column used later.
- Switch on 709 IC Pump and rinse IC system until the absolute conductivity shown on the 732 Detector is stable.
- Switch off 709 IC Pump.

6 Remove coupling

- Remove 6.2620.060 Coupling between column outlet capillary **71** and inlet capillary **82** (see *Fig. 16*). The IC system is now ready for the installation of precolumns, separating columns and suppressor module.

2.7 Precolumns

2.7.1 General information on precolumns

The use of easily exchangeable precolumns protects the separating columns and prolongs their lifetime. The precolumns available from Metrohm (see *section 7.3.2*) are either real precolumns or precolumn cartridges, which are used together with the 6.2821.050 Twin cartridge holder or a 6.2821.040 Cartridge head.



New IC precolumns are normally filled with solution and sealed at both ends. Before the precolumn is installed in the system, it must be ensured that this solution is freely miscible with the eluent used (check manufacturer's specifications).

2.7.2 Precolumns with twin cartridge holder

The precolumn cartridges are installed in the 6.2821.050 Twin cartridge holder as follows (see *Fig. 15*):

1 Install cartridge

- Insert inlet capillary **75** with end piece for precolumn cartridge in Manufit housing **62**.
- Insert outlet capillary **72** with end piece for precolumn cartridge in Manufit pressure screw **57**.
- Remove end caps from the precolumn cartridge **74** (the steel mesh **73** and gaskets **59** are already inserted in the cartridge).
- Fit the two capillary end pieces onto the precolumn cartridge **74** (comply with the flow direction if specified on the column).
- Screw Manufit pressure screw **57** and Manufit housing **62** firmly together.

2 Connect precolumn

- Fit connectors to inlet capillary **75** of the assembled precolumn (see *section 2.5*).
- Connect inlet capillary **75** either using 6.2620.060 Coupling to the column connection capillary **67** mounted on the injection valve (see *section 2.6.4*) or directly to connection "4" of injection valve A or B.
- Shorten outlet capillary **72** of the precolumn to ca. 5 cm and mount the connectors (see *section 2.5*).

3 Rinse the precolumn

- Place a beaker beneath outlet capillary **72**.
- Switch on 709 IC Pump, rinse precolumn for ca. 10 min with eluent and then switch off 709 IC Pump.




Fig. 15: Installing precolumn cartridges

49	Ferrule (6.2620.010)	73	2 Steel meshes (6.2821.020)
50	Pressure screw (6.2620.000)	74	Precolumn cartridge
57	Manufit pressure screw	75	Inlet capillary
59	PTFE gasket (6.2821.010)	76	IC separating column
60	4 Steel meshes (6.2821.020)	77	Manufit pressure screw
62	Manufit housing	78	Steel spacer (6.2821.080)
67	Column connection capillary of injector	79	Steel connector for ferrule (of IC separating column)
72	Outlet capillary	80	Manufit housing

2.7.3 Precolumns with cartridge head

The precolumn cartridges are installed in the 6.2821.040 Cartridge head as follows (see *Fig. 15*):

1 Prepare separating column

- Remove end caps from separating column **76**.
- Unscrew fastening screw from column inlet.
- Take steel connector **79** for ferrule out of fastening screw.

2 Install cartridge

- Remove end caps for precolumn cartridge **74** (the steel mesh **73** and gaskets **59** are already installed in the cartridge).
- Mount steel spacer **78** on separating column **76**.
- Mount precolumn cartridge **74** on the steel spacer (comply with flow direction if specified on the precolumn).
- Screw Manufit pressure screw **77** firmly to separating column **76**.
- Mount steel connector **79** for ferrule on the inlet side of the precolumn cartridge **74**.
- Screw on Manufit housing **80** with Manufit pressure screw **77**.

3 Connect precolumn

- Fit connectors to the column connection capillary **67** mounted on the injection valve (see *section 2.5*).
- Screw column connection capillary **67** to Manufit housing **80**.

2.7.4 IC anion precolumn SUPERSEP

The 6.1009.010 IC Anion Precolumn SUPERSEP has two connections for steel or PEEK capillaries and is installed as follows:

1 Connect precolumn

- Remove end caps from the precolumn.
- Fit connectors to column connection capillary **67** mounted on the injection valve (see *section 2.5*).
- Screw precolumn to column connection capillary **67**.
- Cut a piece from the 6.2620.020 Steel capillary or the 6.1831.010 PEEK capillary as short as possible and fit with connectors (see *section 2.5*).
- Fasten the prepared capillary to the other end of the precolumn.

2 Rinse precolumn

- Place a beaker beneath the outlet capillary of the precolumn.
- Switch on 709 IC Pump and rinse precolumn with eluent for ca. 10 min.
- Switch off IC pump.

2.8 Separating columns and suppressor module

2.8.1 General information on separating columns



New IC separating columns are normally filled with solution and sealed at both ends. Before the column is installed in the system, it must be ensured that this solution is freely miscible with the eluent used (check manufacturer's specifications).

The IC separating columns and precolumns currently available from Metrohm are listed in *section 7.3.2*. A test chromatogram and an information leaflet is provided with each column. You will find additional information concerning these columns in the 8.732.2003 Metrohm Monograph «Ion chromatography» and in special "Application Bulletins", which are available on request free of charge from your local Metrohm agency.



When you install the column, always ensure that this is inserted correctly in accordance with the flow direction shown (arrow must point upwards).

2.8.2 Selection of the sample loop

Selection of the sample loop depends on the separating column used. Normally, the following sample loops are used:

Anion columns	100 µL
Cation columns	10 µL
Columns for suppressor technique	20 µL

Depending on the instrument version, the following sample loops are installed in the 733 IC Separation Center:

Version	Valve	Sample loop	Volume
2.733.0010	A	6.2620.120 (steel)	100 µL
2.733.0020	A	6.2620.120 (steel)	100 µL
	B	6.2620.100 (steel)	10 µL
2.733.0030	A	6.1825.210 (PEEK)	20 µL
2.733.0120	A	6.1825.220 (PEEK)	100 µL
	B	6.1825.230 (PEEK)	10 µL
2.733.0130	A	6.1825.210 (PEEK)	20 µL

If desired, the built-in sample loop can be replaced by one of the sample loops available as an option (see *section 7.3.1*).

2.8.3 General information on suppressor module

The **Metrohm Suppressor Module MSM** for chemical suppression installed in the 2.733.0X30 IC Separation Center comprises a total of 3 suppressor units which are in turn used for suppression, regenerated with sulfuric acid and rinsed with water. To record every new chromatogram under comparable conditions, work is normally carried out with freshly regenerated suppressor. Switching is either automatic together with the valve switching or manual.



The suppressor units must never be regenerated with H_2SO_4 in the same flow direction used for the eluent. You should thus always install the inlet and outlet capillaries as described in section 2.8.6 according to the scheme shown in Fig. 18.

For operation of the suppressor module, a **two-channel peristaltic pump** is needed which conveys the regeneration solution (normally **20 mmol/L H_2SO_4**) and the rinsing solution (normally **dist. H_2O**) to the suppressor units. We advise working with a **flow rate of 0.5 mL/min**.



*The **752 Pump Unit** is available from Metrohm as an option. Two lengths of pump tubing (6.1826.050) are enclosed with this pump (flow rate 0.5 mL/min). Startup and operation of the 752 Pump Unit is described in the Instructions for Use enclosed with the pump.*

The three inlets and outlets numbered 1...3 on the suppressor module each have 2 permanently mounted PTFE capillaries, which must be connected as described in section 2.8.6 (see Fig. 16 and Fig. 17).

To avoid contamination of the suppressor module by foreign particles or bacterial growth, it is advantageous to install an **in-line filter** between the 752 Pump Unit and the suppressor module. The two **6.2821.100 Filter units PEEK** (see section 2.6.3) supplied with the 2.733.0X30 IC Separation Center are eminently suitable for this purpose.



The suppressor module must never be switched in the dry state as there is a danger of blocking. Before every switching operation of the suppressor module, the three suppressor units must have been rinsed for at least ½ h with eluent, regeneration and rinsing solution.




Fig. 16: Interior of the 2.733.0010 IC Separation Center

65 Capillary to 709 IC Pump	83 Mounting rail
66 Pulsation dampener (6.2620.150)	84 Column holder (6.2027.0X0)
67 Column connection capillary	85 Capillary for syringe PEEK capillary, fixed mounting
68 Injection valve	86 Sample loop (6.2620.120) 100 μ L, fixed mounting, steel
76 IC separating column	87 Inlet capillary for injector Steel or PEEK capillary
81 Detector block	88 PTFE aspirating tubing fixed mounting
82 Inlet capillary for detector block Steel capillary, fixed mounting	

2.8.4 One-channel system without suppressor module

With the one-channel system without suppressor module, the IC separating column is installed in the 2.733.0010 IC Separation Center as follows (see *Fig. 16*):

1 Connect column to injector

- Remove end caps from column **76**.
- *without precolumn:*
Screw inlet end of separating column **76** (note flow direction) to column connection capillary **67** mounted on the injector.
- *with precolumn in twin cartridge holder:*
Screw inlet end of separating column **76** (note flow direction) to outlet capillary **72** (see *Fig. 15*).
- *with precolumn in cartridge head:*
Install separating column **76** (note flow direction) in the cartridge head (see *Fig. 15*) as described in *section 2.7.3*.

2 Rinse column

- Place a beaker beneath the column outlet.
- Switch on 709 IC Pump and rinse the column for ca. 10 min with eluent.
- Switch off IC Pump.

3 Connect column to detector block

- Screw outlet end of separating column **76** to the inlet capillary **82** permanently mounted on the detector block.

4 Fix column

- Insert one or two column holders **84** (6.2027.030, 6.2027.040 or 6.2027.050) in the mounting rails **83** and fasten separating column in the column holder.

2.8.5 Two-channel system without suppressor module

With the two-channel system without suppressor module (2.733.0X20 IC Separation Center), the first IC separating column is connected on the left side to injection valve A and detector block A as with the one-channel system (see *section 2.8.4* and *Fig. 16*). The second column is connected on the right side to injection valve B and detector block B in an analogous manner.




Fig. 17: Interior of the 2.733.0X30 IC Separation Center

43 Rear panel opening	83 Mounting rail
65 Capillary to 709 IC Pump	84 Column holder (6.2027.0X0)
66 Pulsation dampener (6.2620.150)	85 Capillary for syringe PEEK capillary, fixed mounting
67 Column connection capillary PEEK capillary	86 Sample loop (6.1825.210) 20 μ L, fixed mounting, PEEK
68 Injection valve	87 Inlet capillary for injector Steel or PEEK capillary
76 IC separating column	88 PTFE aspirating tubing fixed mounting
81 Detector block	89 Suppressor inlet capillary for eluent
82 Inlet capillary for detector block (fixed mounting)	90 Suppressor inlet capillary for H_2SO_4

91	Suppressor outlet capillary for H ₂ SO ₄	94	Suppressor outlet capillary for eluent
92	Suppressor outlet capillary for H ₂ O	95	Suppressor module
93	Suppressor inlet capillary for H ₂ O	96	Coupling (6.2620.060; steel) or Coupling (6.2744.040; PEEK)

2.8.6 One-channel system with suppressor module

In the case of the one-channel system with suppressor module, first the IC separating column is installed in the 2.733.0X30 IC Separation Center (see Fig. 17) and then the suppressor module is connected to the 752 Pump Unit needed for operation (see Fig. 18). Proceed as follows:

1 Connect column to injector

- Remove end caps from column **76**.
- *without precolumn:*
Screw inlet end of separating column **76** (note flow direction) to column connection capillary **67** mounted on the injector.
- *with precolumn in twin cartridge holder:*
Screw inlet end of separating column **76** (note flow direction) to outlet capillary **72** (see Fig. 15).
- *with precolumn in cartridge head:*
Install separating column **76** (note flow direction) in the cartridge head as described in section 2.7.3 (see Fig. 15).

2 Rinse column

- Place a beaker beneath the column outlet.
- Switch on 709 IC Pump and rinse column with eluent for ca. 10 min.
- Switch off IC pump.

3 Connect column to suppressor module


- Cut inlet capillary **89** marked with "Eluent" of suppressor module **95** to the desired length using a sharp cutting tool (e.g. razor blade).
- Screw inlet capillary **89** to outlet end of separating column **76** using a 6.2744.010 Compression fitting.

4 Fix column

- Insert one or two column holders **84** (6.2027.030, 6.2027.040 or 6.2027.050) in mounting rails **83** and fasten separating column in the column holder.

5 Connect suppressor module to detector block

- Cut outlet capillary **94** marked with "Detector" of suppressor module **95** to the desired length using a sharp cutting tool (e.g. razor blade).



89 Suppressor inlet capillary for eluent

90 Suppressor inlet capillary for H₂SO₄

91 Suppressor outlet capillary for H₂SO₄

92 Suppressor outlet capillary for H₂O

93 Suppressor inlet capillary for H₂O

94 Suppressor outlet capillary for eluent

Fig. 18: Connections at suppressor module

- Screw outlet capillary **94** to coupling **96** using a 6.2744.010 Compression fitting.
- Screw inlet capillary **82** of detector block **81** to other end of coupling **96**.

6 Fix connection suppressor – detector block

- Insert one of the column holders **84** (6.2027.030, 6.2027.040 or 6.2027.050) in mounting rail **83** and fasten coupling **96** in the column holder.

7 Prepare 752 Pump Unit

- Take two tubing cartridges out of the holder of the 752 Pump Unit.
- Lay a length of 6.1826.050 Pump tubing in each of the two tubing cartridges and reinsert these in the holder without kinking the pump tubing.
- Regulate contact pressure on pump tubing in accordance with the instructions printed on the pump.



Pump tubing is consumable material with a lifetime which depends on the contact pressure. For this reason, raise the tubing cartridges completely by loosening the clip on the right if the pump is switched off for some considerable time (the set contact pressure remains unchanged).

8 Suppressor connection 2: H₂SO₄

- Loosen rotary nipple screwed onto the interior side of connection **27**. Pull inlet capillary **90** marked with "H₂SO₄" (see Fig. 17 and Fig. 18) by hand out of the opening of connection **27** as far as required. Retighten nipple on the interior side of connection **27** to fix inlet capillary **90**.

- Attach inlet capillary **90** using a 6.2744.010 Compression fitting to the connector **55** of the filter unit PEEK (see *section 2.6.3*).
- Attach a piece of the 6.1803.020 PTFE tubing (from 752 accessories) cut to the required length using a 6.2744.010 Compression fitting to the connector **53** with filter at the other end of the Filter unit PEEK.
- Attach the other end of the PTFE tubing using a 6.2744.010 Compression fitting to the 6.2744.030 Coupling (from 752 accessories) and mount it on the outlet end of the first length of pump tubing.
- Mount a 6.2744.030 Coupling to the inlet end of the first length of pump tubing. Attach a piece of the 6.1803.020 PTFE tubing (aspirating tubing, from 752 accessories) cut to the required length using a 6.2744.010 Compression fitting to the other end of this coupling.
- Immerse the other end of the aspirating tubing in a vessel containing regeneration solution (normally 20 mmol/L H₂SO₄) and fix in place.
- Pull outlet capillary **91** marked with "Waste" of the suppressor module through the rear panel opening **43**, lead it to a sufficiently large waste container and fix it in place.

9 Suppressor connection 3: H₂O

- Loosen rotary nipple screwed onto the interior side of connection **28**. Pull inlet capillary **93** marked with "H₂O" (see *Fig. 17* and *Fig. 18*) by hand out of the opening of connection **28** as far as required. Retighten nipple on the interior side of connection **28** to fix inlet capillary **93**.
- Attach inlet capillary **93** using a 6.2744.010 Compression fitting to the connector **55** of the filter unit PEEK (see *section 2.6.3*).
- Attach a piece of the 6.1803.020 PTFE tubing (from 752 accessories) cut to the required length using a 6.2744.010 Compression fitting to the connector **53** with filter at the other end of the Filter unit PEEK.
- Attach the other end of the PTFE tubing using a 6.2744.010 Compression fitting to the 6.2744.030 Coupling (from 752 accessories) and mount it on the outlet end of the second length of pump tubing.
- Mount a 6.2744.030 Coupling to the inlet end of the second length of pump tubing. Attach a piece of the 6.1803.020 PTFE tubing (aspirating tubing, from 752 accessories) cut to the required length using a 6.2744.010 Compression fitting to the other end of this coupling.
- Immerse the other end of the aspirating tubing in a vessel containing rinsing solution (normally dist. H₂O) and fix in place.
- Pull outlet capillary **92** marked with "Waste" of the suppressor module through the rear panel opening **43**, lead it to a sufficiently large waste container and fix it in place.

2.8.7 Leak testing and conditioning

Before sample solutions can be injected in the IC system, the entire system must be tested for leaks and then conditioned with eluent until the baseline is stable. Proceed as follows:

1 Switch on 709 IC Pump

- Immerse aspirating capillary of the 709 IC Pump in eluent.
- Set flow rate recommended on the 709 IC Pump for the separating column used (normally 0.5...2 mL/min).
- Set maximum shutoff pressure on the 709 IC Pump (normally ca. 3 MPa above the pressure observed with the column used).
- Switch on 709 IC Pump.

2 Check for leaks

- Check all capillaries and their connections between the 709 IC Pump and the detector block for escaping liquid. If eluent escapes anywhere, the appropriate compression fitting must be tightened further or changed.

3 Switch on 732 IC Detector

- Switch on 732 IC Detector with mains switch **9**.
- Set **operating temperature**: Under the <CONFIG> key, enter the desired value for the parameter "thermostat" (default value = 35 °C).
- Enter **cell constant**: Under the <CONFIG> key, enter the value printed on the detector block for the parameter "cell constant".
- Set **measuring range**: Under the <PARAM> key, set the parameter "range" so that the displayed absolute conductivity value of the eluent lies within the selected range (default value = 1 mS/cm).
- Set **full-scale range**: Under the <PARAM> or <FULL SCALE> key, enter the desired value for the parameter "full scale" (default value = range/1). To start with, it is advisable not to set too narrow a full-scale range as the conductivity of the eluent in the conditioning phase can change greatly until attainment of a stable temperature.

4 Condition IC system

- Rinse the IC system with eluent until the desired stability of the baseline is reached (normally 30...60 min; if the eluent is changed, the establishment of the ion exchanger equilibrium on the separating column can take longer).

5 Condition suppressor (if present)

- Load method "Prep-MSM" and start this method. The suppressor module is switched to the next position every 20 min and in this way conditioned.

2.9 Connection of external devices

2.9.1 Connection of a recorder

For the connection of a recorder the 732 IC Detector has the two analog outputs **11** (0...1 V) and **12** (0...10 mV) available (see *Fig. 3*). A diagram of the circuit of the two analogue outputs can be found in *section 6.3*. The polarity of the output signals at the analog output sockets can be switched at any time under the <PARAM> key (see *section 4.5.1*).

For the connection of recorders with banana connection sockets, the 6.2115.010 Cable available as an option can be used.

2.9.2 Connection of «IC Metrodata for Win95»

With «IC Metrodata for Win95» (ordering number 2.714.0310), Metrohm offers a chromatography data system for the automatic evaluation of chromatograms using a PC. It comprises a PC board and the associated integration software. You will find all relevant information about the installation of hardware and software in the *714 Instructions for Use*.

2.9.3 Connection of the 750 Autosampler

The 750 Autosampler available from Metrohm as an option is an automatic sampler for ion chromatography. The apparatus accommodates max. 128 samples each of 730 μL , which are automatically transferred to the sample loops attached to the injection valves of the 733 IC Separation Center. The electrical connection of the 750 Autosampler and the tubing connections for the sample feed is described in the *750 Instructions for Use*.

2.9.4 Connection of the 766 IC Sample Processor

The 766 IC Sample Processor available from Metrohm as an option is an automatic sampler for ion chromatography. The apparatus accommodates max. 127 samples each of 11 mL, which are automatically transferred to the sample loops attached to the injection valves of the 733 IC Separation Center. The electrical connection of the 766 IC Sample Processor and the tubing connections for the sample feed is described in the *766 Instructions for Use*.

2.9.5 Connection of the 791 VA Detector

The 791 VA Detector available from Metrohm as an option enables the electrochemical (amperometric) detection with carbon or metal electrodes. The accessories include a flow-through cell which can be fixed to the conductivity detector in the 733 IC Separation Center. You will find detailed information about installation and operation of the 791 VA Detector in the *791 Instructions for Use*.

2.9.6 Connection of a printer

External printers are normally connected to the 732 IC Detector via the RS232 interfaces **16** (see *Fig. 3*). If no 709 IC Pump is connected, a printer can also be attached to connection **15** "709 IC Pump". You will find more detailed information on the RS232 interfaces in *section 6.1*.








Before a printer is connected to the RS232 interfaces **16** or **15**, the 732 IC Detector always be **switched off** with mains switch **9** !

Printers with the following printer drivers can be connected to the two RS232 interfaces of the 732 IC Detector (see also *section 4.4.2*):

- IBM IBM Proprinter and printers with IBM emulation
- Epson EPSON printers and printers with EPSON emulation
- Seiko Seiko printer DPU-411/414
- Citizen Citizen printer IDP562 RS
- HP HP printer and printers with PCL3 emulation

The following table provides information on a few selected printers.

Printer	Cable	Settings on 732 IC Detector	Settings on printer
IBM Proprinter	6.2125.050	>CONFIG/printer send to: IBM >CONFIG/RS settings baud rate: 9600 data bit: 8 stop bit: 1 parity: none handshake: HWS	see printer manual
Epson with 6-pin round connector	6.2125.040	>CONFIG/printer send to: Epson >CONFIG/RS settings baud rate: 9600 data bit: 8 stop bit: 1 parity: none handshake: HWS	Settings of the DIP switches: 
Epson with additional serial interface #8148	6.2125.050	>CONFIG/printer send to: Epson >CONFIG/RS settings baud rate: 9600 data bit: 8 stop bit: 1 parity: none handshake: HWS	Settings of the DIP switches on the <u>interface</u> : 
Epson LX-300	6.2125.050	>CONFIG/printer send to: Epson >CONFIG/RS settings baud rate: 9600 data bit: 8 stop bit: 1 parity: none handshake: HWS	see printer manual
Epson and Canon printers with parallel interface	6.2125.020 + 2.145.0300 Serial/ parallel converter	>CONFIG/printer send to: Epson >CONFIG/RS settings baud rate: 9600 data bit: 8 stop bit: 1 parity: none handshake: HWS	see printer manual

Printer	Cable	Settings on 732 IC Detector	Settings on printer																																				
Seiko DPU-411	6.2125.020	<pre>>CONFIG/printer send to: Seiko >CONFIG/RS settings baud rate: 9600 data bit: 8 stop bit: 1 parity none handshake: HWS</pre>	<p>Settings of the DIP switches:</p> <p>DIP01 DIP02</p>  <p>The switchable 7-bit ASCII character set of the printer will be automatically set to the national character sets in accordance with the set dialog language.</p>																																				
Seiko DPU-414	6.2125.130	<pre>>CONFIG/printer send to: Seiko >CONFIG/RS settings baud rate: 9600 data bit: 8 stop bit: 1 parity none handshake: HWS</pre>	<p>Settings of the DIP switches:</p> <table border="1"> <thead> <tr> <th></th> <th>Dip SW-1</th> <th>Dip SW-2</th> <th>Dip SW-3</th> </tr> </thead> <tbody> <tr><td>1</td><td>OFF</td><td>ON</td><td>ON</td></tr> <tr><td>2</td><td>ON</td><td>OFF</td><td>ON</td></tr> <tr><td>3</td><td>ON</td><td>ON</td><td>ON</td></tr> <tr><td>4</td><td>OFF</td><td>ON</td><td>ON</td></tr> <tr><td>5</td><td>ON</td><td>ON</td><td>OFF</td></tr> <tr><td>6</td><td>OFF</td><td>ON</td><td>ON</td></tr> <tr><td>7</td><td>ON</td><td>OFF</td><td>ON</td></tr> <tr><td>8</td><td>ON</td><td>OFF</td><td>ON</td></tr> </tbody> </table> <p>The switchable 7-bit ASCII character set of the printer will be automatically set to the national character sets in accordance with the set dialog language.</p>		Dip SW-1	Dip SW-2	Dip SW-3	1	OFF	ON	ON	2	ON	OFF	ON	3	ON	ON	ON	4	OFF	ON	ON	5	ON	ON	OFF	6	OFF	ON	ON	7	ON	OFF	ON	8	ON	OFF	ON
	Dip SW-1	Dip SW-2	Dip SW-3																																				
1	OFF	ON	ON																																				
2	ON	OFF	ON																																				
3	ON	ON	ON																																				
4	OFF	ON	ON																																				
5	ON	ON	OFF																																				
6	OFF	ON	ON																																				
7	ON	OFF	ON																																				
8	ON	OFF	ON																																				
Citizen IDP562-RS	6.2125.050	<pre>>CONFIG/printer send to: Citizen >CONFIG/RS settings baud rate: 9600 data bit: 8 stop bit: 1 parity none handshake: HWS</pre>	<p>Settings of the DIP switches:</p>  <p>The switchable 7-bit ASCII character set of the printer can be changed to the national character sets only by setting the DIP switches 4 and 5 in the printer:</p> <table border="1"> <thead> <tr> <th>4</th> <th>5</th> <th>Character set</th> </tr> </thead> <tbody> <tr><td>open</td><td>open</td><td>USA</td></tr> <tr><td>closed</td><td>closed</td><td>Great Britain</td></tr> <tr><td>closed</td><td>open</td><td>France</td></tr> <tr><td>open</td><td>closed</td><td>Germany</td></tr> </tbody> </table> <p>Spanish does not have its own character set (it is best to select French).</p>	4	5	Character set	open	open	USA	closed	closed	Great Britain	closed	open	France	open	closed	Germany																					
4	5	Character set																																					
open	open	USA																																					
closed	closed	Great Britain																																					
closed	open	France																																					
open	closed	Germany																																					
HP Deskjet with serial interface	6.2125.050 or adapter cable 25-pin neg./9-pin pos. (e.g. HP C2933A)	<pre>>CONFIG/printer send to: HP >CONFIG/RS settings baud rate: 9600 data bit: 8 stop bit: 1 parity none handshake: HWS</pre>	<p>Settings of the DIP switches:</p> <p>A B</p> 																																				
HP Laserjet with serial interface	Adapter cable 25-pin neg./9-pin pos. (e.g. HP C2933A)	<pre>>CONFIG/printer send to: HP >CONFIG/RS settings baud rate: 9600 data bit: 8 stop bit: 1 parity none handshake: HWS</pre>	see printer manual																																				
HP Deskjet/Laserjet with parallel interface	6.2125.020 + 2.145.0300 Serial/parallel converter	<pre>>CONFIG/printer send to: HP >CONFIG/RS settings baud rate: 9600 data bit: 8 stop bit: 1 parity none handshake: HWS</pre>	see printer manual																																				



If you connect other printers, ensure that these emulate a printer mode supported by the 732 IC Detector. Most printers with a serial interface are connected using the 6.2125.050 Cable. Printers with a parallel interface need a serial/parallel converter (e.g. 2.145.0300) and the 6.2125.020 Cable.

2.9.7 Connection of a PC

IBM-compatible PCs are connected to the 732 IC Detector via RS232 interface **16** (see Fig. 3). You will find more detailed information on the RS232 interface in section 6.1, which also describes remote control of the 732 IC Detector via the RS interface.



*Before a PC is connected to RS232 interface **16**, the 732 IC Detector must always be switched off using mains switch **9** !*

The following table provides information on the connection of PCs. It lists the required cables and details on the configuration of the 732 IC Detector 732 and PC.

PC	Cable	Settings on 732 Detector	Settings on PC
PC with 25-pin RS232 connector	6.2125.060	>CONFIG/printer send to: IBM >CONFIG/RS settings. settings same as on PC	setting of the RS parameters dependent on control program
PC with 9-pin RS232 connector	6.2125.060 + 6.2125.010	>CONFIG/printer send to: IBM >CONFIG/RS settings settings same as on PC	setting of the RS parameters dependent on control program

2.9.8 Connection of devices to the remote interface

Any external devices can be connected to the 25-pin remote interface **17** (see Fig. 3). The 732 IC Detector can be remote controlled via the 8 input lines, the 8 output lines can be used to control external devices.



*Before an external device is connected to the remote interface **17**, the 732 IC Detector must always be switched off using mains switch **9** !*

The pin assignment of the remote interface, its functions and the electrical requirements and the conditions are described in section 6.2, the assignment of the remote input lines can be found in section 4.4.1.

3 Operating tutorial



This section introduces you to the operation of the 732 IC Detector and 733 IC Separation Center by means of a brief operating tutorial which describes the basic operating steps needed for the recording of an ion chromatogram and the development of a method.

The determination of the anionic content of a drinking water sample with the PRP-X100 IC Anion Column using the single column technique is used as an illustrative example. Please note that the steps and parameter settings described apply only to this example. If you perform a different determination, use a different separating column or other peripheral units, the procedures described in the tutorial must be modified appropriately.

For further explanations of the operation, please refer to section 4, which describes the functions of the individual keys and the programming in detail.

3.1 Requirements

For the determination of anions in drinking water described in this tutorial, the following instruments, accessories and solutions are required:

- **2.732.0X10 IC Detector**
- **2.733.0XX0 IC Separation Center**
without connection of the suppressor module
- **100 mL Sample loop (6.2620.120 steel or 6.1825.220 PEEK)**
already integrated in the 2.733.0010 and 2.733.0X20 IC Separation Center
- **2.709.0X10 IC Pump**
In the operating tutorial, the pump is connected to the 732 IC Detector and remote controlled by this.
- **6.2125.060 Cable** (connecting cable 732 – 709)
RS232 connecting cable for remote control of the 709 IC Pump by the 732 IC Detector.
- **6.2620.150 Pulsation dampener**
Use of the pulsation dampener available as an option is optional, but is advisable to protect the separating column.

- **6.2821.100 Filter unit PEEK**
Filter between 709 IC Pump and injection valve to avoid contamination.
- **6.1005.000 PRP-X100 IC Anion Column**
- **Eluent**
2 mmol/L phthalic acid / 8% acetone / pH 5.0 in dist. H₂O
flow: 2 mL/min
- **Standard**
Standard solution with 5 ppm Cl⁻ and 10 ppm each of NO₃⁻ and SO₄²⁻ (in dist. H₂O)
- **2.714.0310 IC Metrodata for Win95**
The «IC Metrodata for Win95» PC integration system, which comprises a PC board and the associated PC integration software, is used to record and evaluate chromatograms on a PC. The installation, operation and method development are not described in this brief tutorial. You will find all relevant information in the *Instructions for Use of IC Metrodata for Win95*.

Instead of «IC Metrodata for Win95», you can also use other data recording systems, integrators, recorders or printers, but their use is also not described in this tutorial.

3.2 Preparations

Before you start this brief tutorial, the entire IC system must be correctly installed as described in *section 2*. In what follows, the most important points for the installation are described once again (for details, see the sections mentioned).

1 Install 732 IC Detector and 733 IC Separation Center

- | | |
|--|------------------------|
| ⇒ Setting up instruments | <i>section 2.2</i> |
| ⇒ Installing and connecting detector block | <i>section 2.3.1/2</i> |
| ⇒ Mounting syringe and aspirating tubing | <i>section 2.3.3</i> |
| ⇒ Mounting drain tube | <i>section 2.3.4</i> |
| ⇒ Power supply connection | <i>section 2.4</i> |

2 Prepare eluent

- | | |
|--|----------------------|
| ⇒ Preparing eluent:
2 mmol/L phthalic acid / 8% acetone / pH 5.0 in dist. H ₂ O
(pH value of solution adjusted with NaOH) | |
| ⇒ Microfiltering and degassing eluent | <i>section 5.1.3</i> |
| ⇒ Stirring eluent in supply vessel | <i>section 5.1.3</i> |

3	Install 709 IC Pump	<i>*(see 709 Instructions for Use)</i>
	⇒ Setting up pump	<i>section 2.1*</i>
	⇒ Mounting pump head	<i>section 2.2*</i>
	⇒ Tubing connections	<i>section 2.4*</i>
	⇒ Power supply connection	<i>section 2.6*</i>
	⇒ Switching on pump	<i>section 2.6.4*</i>
	⇒ Degassing pump	<i>section 2.7*</i>

4	Connect 709 IC Pump	
	⇒ Electrical connection at 732 IC Detector (requires 6.2125.060 RS Cable)	<i>section 2.6.1</i>
	⇒ Installing pulsation dampener (option)	<i>section 2.6.2</i>
	⇒ Installing filter unit PEEK	<i>section 2.6.3</i>
	⇒ Making connection to injection valve	<i>section 2.6.5</i>
	⇒ Passivating the IC system	<i>section 2.6.7</i>

5	Connect separating column	
	⇒ Connecting PRP-X100 IC Anion Column	<i>section 2.8</i>
	⇒ Testing for leaks	<i>section 2.8.7</i>
	⇒ Conditioning	<i>section 2.8.7</i>

6	Connect external devices	
	⇒ Connecting recorder (if used)	<i>section 2.9.1</i>
	⇒ Connecting «IC Metrodata for Win95» or other recording or evaluation system (if used)	<i>section 2.9.2</i>
	⇒ Connecting printer (if used)	<i>section 2.9.6</i>

3.3 Putting into operation

After the complete IC system has been installed as described in *section 3.2*, it can be put into operation. In what follows, all operating steps up to the first calibration with the standard solution are described in sequence.

Please note that all displays refer to the condition in which the instrument was put into operation for the first time (initial condition). If you do not run through this tutorial until a later date, differences in regard to the dialog language and the parameter values can appear. How you return to the initial condition in such a case is described in *section 5.4.9*.

In this description of putting the instruments into operation, it is assumed that all instruments are shut down and must first be switched on. If this is not the case (e.g. if you start the operating tutorial immediately following the conditioning), you can skip steps **1 – 3**.

1 Switch on external devices

- ⇒ Switch on recorder (if used).
- ⇒ Switch on «IC Metrodata for Win95» and PC (if used).
- ⇒ Switch on other recording and evaluation system (if used).

2 Switch on 709 IC Pump

- ⇒ Switch on 709 IC Pump with mains switch **29** (see *section 2.6.4 of the 709 Instructions for Use*).


3 Switch on remote control for 709 IC Pump

(the pump must be connected to the 732 IC Detector with a **6.2125.060 RS Cable** for this)

- ⇒ Set sliding switch **36** at the rear of the 709 IC Pump to "RS 232" (see *Fig. 2 of the 709 Instructions for Use*).
- ⇒ Press key **8** <EXT.> on the 709 IC Pump to switch on the external control (see *section 3.6 of the 709 Instructions for Use*).

When the external control is activated, LED **7** above the <EXT.> key lights up. All keys of the 709 IC Pump with the exception of the <SELECT> key are blocked, operation is now possible only via the 732 IC Detector.

4 Switch on 732 IC Detector




⇒ Switch on 732 IC Detector with mains switch **9** at the rear of the instrument.

151.2 mS/cm	0.0 min
method	DEFAULT

After the instrument has been switched on, display **1** lights up and shows the status messages for the basic instrument condition.

The 732 IC Detector is now in the standby mode of conductivity measurement. The displayed values have the following meaning:



Instead of the loaded method, other status messages can be selected which are then continuously displayed in the standby mode in the bottom line of display **1**. Proceed as follows:

5 Select status message



⇒ Press the <SELECT> key.

151.2 mS/cm	0.4 min
1995-09-11	14:15:27

The bottom line of the display shows the current **date** and the current **time**. Please refer to point **6** to see how to change the date and time.



⇒ Press the <SELECT> key again.

151.2 mS/cm	0.6 min
full scale	1.00 mS/cm

The bottom line of the display now shows the **full-scale range**.



⇒ Press the <SELECT> key again.

151.2 mS/cm	0.9 min
abs. cond.	151.2 mS/cm

The bottom line of the display now shows the current measurement of the **absolute conductivity**, which is identical to the value already shown in the top line on the left.



⇒ Press the <SELECT> key again.

```
151.2 mS/cm   1.2 min
pump ready
```

The bottom line of the display shows the status message for the **709 IC Pump**. If the remote operation is switched on correctly (see point 4), the status message "pump ready" now appears.

If the IC Pump is operated independently, the message "pump not responding" appears instead.



⇒ Press the <SELECT> key again.

```
151.2 mS/cm   0.9 min
method        DEFAULT
```

The bottom line of the display shows the status message for the loaded **method** once again. The "DEFAULT" method in which all parameters are set to the initial values is loaded as standard. You have now returned to the status messages displayed at the start.

The next point of the tutorial is concerned with the fundamentals of data entry using the settings of the instrument configuration required for the example.

6 Set configuration



⇒ Press the <CONFIG> key to open the main menu for the instrument configuration.

```
CONFIG
>CONFIG/detector
```

The name of the selected key appears in the top line of the display, the title of the submenu "detector", which contains various inquiries for the conductivity detector, is shown in the bottom line.

The main menu of the <CONFIG> key contains several such submenus, which are selected by repeated pressing of this key. Each submenu has a title marked by ">".



⇒ Now press the <ENTER> key.

```
>CONFIG/detector
thermostat:    35 °C
```

This moves you from the title to the individual inquiries of the submenu "detector"; the title continues to be shown in the top line. The **operating temperature** of the conductivity cell appears in the display as the first inquiry.

Please note the character ":" with this display. It always appears when values can not be entered using the numeric keys, but must be selected from preset values using the <SELECT> key. However, for our example the preset value of 35 °C does not have to be changed.



⇒ Confirm the set operating temperature by pressing the <ENTER> key. The displayed value will be accepted and a switch made to the next inquiry.

```
>CONFIG/detector
'zero' unit:      mS/cm
```

The next inquiry of the submenu "detector" appears, namely the **unit for display of the auto-zero value**.

This parameter also contains the character ":". The <SELECT> key can be used to select the other values "%fs" (% full scale) or "mV". For our example, however, the preset value "µS/cm" need not be changed.



⇒ Confirm the set unit by pressing the <ENTER> key.

```
>CONFIG/detector
cell constant  16.7 /cm
```

The **cell constant** of the conductivity cell in the detector block appears as the next inquiry of the submenu "detector". Each detector block has a characteristic cell constant which is determined in the factory and printed on the block. To ensure that the absolute conductivity is correctly displayed, this value must be entered on first time startup.

Please note that this parameter does not contain the character ":". This means that here values must be entered using the numeric keys. Hence press the corresponding numeric keys in the correct sequence to enter the value printed on the detector block. During the entry you can always return to the default value by pressing the <CLEAR> key and restart the entry.

```
>CONFIG/detector
cell constant  17.1 /cm
```

A cell constant of 17.1 /cm is entered here as an example.



⇒ Confirm the cell constant just entered by pressing the <ENTER> key.

```
CONFIG
>CONFIG/printer
```

As the previous inquiry was the last of the submenu "detector", a switch is now made automatically to the title of the next submenu "printer", which contains various inquiries concerning the printout on an external printer. As no printer is attached in our example, the next submenu can be selected directly.



⇒ To do this, press the <CONFIG> key.

```
CONFIG
>CONFIG/print meas. value
```

The title of the submenu "print meas.value" appears in the bottom display line. The submenu contains various inquiries for the measured value printout on an external printer, but again we are not interested in this at the present time.



⇒ Press the <CONFIG> key again.

```
CONFIG
>CONFIG/auxiliaries
```

The title of the submenu "auxiliaries" appears in the bottom line of the display. This submenu contains, among other things, the inquiries for entry of the date, time and dialog language.



⇒ Press the <ENTER> key to move to the inquiries.

```
>CONFIG/auxiliaries
run number          0
```

The display shows the selection of the run number as the first inquiry.

As we have no further interest in this here, we can proceed directly to the next inquiry.



⇒ Press the <ENTER> key.

```
>CONFIG/auxiliaries
number of cycles    1
```

This display shows the selection of the number of cycles for loop programs.

As this is yet another parameter which does not interest us at present, we can proceed to the next inquiry.



⇒ Press the <ENTER> key again.

```
>CONFIG/auxiliaries
>CONFIG/aux/event
```

The display shows the submenu for EVENT settings, which has not to be opened.



⇒ Press the <CONFIG> key.

```
>CONFIG/auxiliaries
date                1995- 09- 11
```

The display now shows the inquiry regarding the current **date** with the numeric data for year, month and day.

If the displayed date is the same as the current date, you need only press the <ENTER> key to confirm it.

However, if you wish to change this date, enter the new numeric values in the order year – month – day using the numeric keys, e.g. "1995-10-05" for October 5, 1995.



⇒ Confirm the inputted date by pressing the <ENTER> key.

```
>CONFIG/auxiliaries
time                16: 43: 27
```

The current **time** with the numeric data for hours, minutes and seconds appears in the display as the next inquiry.

If the displayed time matches the current time, you need only press the <ENTER> key to confirm this.

If you wish to change the displayed time, enter the new numeric values in the order hours – minutes – seconds using the numeric keys, e.g. "08:32:00". The new time does not become active until you confirm it by pressing the <ENTER> key.



⇒ Confirm the new time by pressing the <ENTER> key.

```
>CONFIG/auxiliaries
dialog:      english
```

The display now shows the selection of the **dialog language**, which is set to "english", as the next inquiry. The character ":" again appears with this display. This means that values must be selected from preset values using the <SELECT> key.



⇒ Press the <SELECT> key to select the next language setting.

```
>CONFIG/auxiliaries
dialog:      deutsch
```

"deutsch" is now selected as the dialog language.



⇒ You could confirm this language setting by pressing the <ENTER> key to immediately switch the dialog language to "deutsch" (or any other language you have selected), but instead:

```
>CONFIG/auxiliaries
device label
```

Use the <SELECT> key to move through the selection back to "english" and confirm with the <ENTER> key. The display shows the device label as the next inquiry. However, at present we are not interested in the further inquiries of this submenu.



⇒ Press the <QUIT> key.

```
CONFIG
>CONFIG/auxiliaries
```

You now quit the inquiries and return to the title of the submenu "auxiliaries".



⇒ Press the <QUIT> key again to return to the standby mode of the instrument.

```
151.2 mS/cm    4.3 min
method         DEFAULT
```

The 732 IC Detector is now again in the standby mode and the display shows the status messages.



Fundamentals of the instrument dialog

Main menu Each key of the 732 IC Detector opens a main menu whose thematically arranged submenus can be selected by repeated pressing of this key. The name of the key always appears in the top line of the display.

Submenu Each submenu has a title which is marked by ">" and appears in the bottom line of the display. <ENTER> is used to move from the title to the individual inquiries where the most important settings of the instrument can be changed. <QUIT> effects a return to the standby mode of the instrument.

Inquiries With inquiries without ":" the valves must be entered using the numeric keys. The inputted value is accepted with <ENTER> and the next inquiry appears.

In the case of inquiries with ":" the admissible values must be selected using the <SELECT> key. The set value is accepted with <ENTER> and the next inquiry appears.

Depending on the parameter, <CLEAR> is used to reset the displayed value to the smallest possible value or the default value. The <CLEAR> key is also used to abort wrong entries.

<QUIT> is used to quit the inquiries and return to the title of the submenu.

7 Enter parameters for the conductivity detector



⇒ Press the <PARAM> key to open the main menu for the parameter settings.

```
PARAM
>PARAM/detector
```

The title of the submenu "detector" appears that contains various inquiries for the conductivity detector.



⇒ Press the <ENTER> key to move to the inquiries.

```
>PARAM/detector
range:      1.00 mS/cm
```

The first inquiry shows the **measuring range** for which 7 stages from 0...100 $\mu\text{S}/\text{cm}$ to 0...10 mS/cm are available. Select the measuring range to ensure that the conductivity value of the eluent used easily lies within the selected range.



⇒ Press the <SELECT> key until the desired measuring range appears in the display.

```
>PARAM/detector
range:      200 mS/cm
```

The eluent used for the PRP-X100 Anion Column has a conductivity of ca. 150 $\mu\text{S}/\text{cm}$. A range of 200 $\mu\text{S}/\text{cm}$ is thus selected as the measuring range.



⇒ Confirm the selected measuring range by pressing the <ENTER> key.

```
>PARAM/detector
full scale:      200 mS/cm
```

The next inquiry which appears concerns the **full-scale range**. The full-scale range (operating range) is used to set the desired sensitivity for the display and analog output of the measuring signal during the recording of the chromatogram. This value is set using the <FULL SCALE> key, but not until later (see section 3.4, point 2).



⇒ Confirm the preset value corresponding to the measuring range by pressing the <ENTER> key.

```
>PARAM/detector
temp. coeff.:   2.5 %/°C
```

The next inquiry concerns the **temperature coefficient** for the automatic conversion of the conductivity from the operating temperature of the measuring cell to the reference temperature of 20° C. The preset value of "2.5 %/°C" applies to anions, the value applicable to cations "1.5 %/°C" can be selected by pressing the <SELECT> key.



⇒ Confirm the preset value for anions by pressing the <ENTER> key.

8 Enter parameters for analog output

```
PARAM
>PARAM/anal og output
```

After the last action of point 7, the title of the submenu "analog output" appears.



⇒ Press the <ENTER> key to move to the inquiries.

```
>PARAM/anal og output
pol ari ty:      +
```

The **polarity** of the analog output signal can be switched with the <SELECT> key between "+" and "-". To ensure that a positive signal is always outputted for peaks, normally "+" is used for the determination of anions and "-" for that of cations.



⇒ Confirm the correct polarity "+" for our example with the <ENTER> key.

```
>PARAM/anal og output
offset:          0 %fs
```

The second inquiry concerns the **offset** of the analog output signal. This offset of the zero point of the conductivity can be set to 10 % or 50 % of the full-scale range using the <SELECT> key. Such an offset of the zero point is primarily advisable when the integrator or the evaluation system do not accept any negative voltage values. As this is not the case for «IC Metrodata for Win95», the value can be left at "0 %fs".



⇒ Confirm the preset value with the <ENTER> key.

```
>PARAM/anal og output
damp ing:                off
```



The electronic **damping** of the analog output signal must normally not be switched on.

⇒ Confirm the switched-off status of the damping with the <ENTER> key.

```
PARAM
>PARAM/pl ot
```

The title of the submenu "p1ot" appears. This submenu contains various inquiries regarding the graphics printout on an external printer. As we shall not deal with the settings for such a printer in this brief tutorial, we can switch directly to the next inquiry.

9 Enter parameters for 709 IC Pump



If you do not operate the 709 IC Pump by remote control via the 732 IC Detector, the parameters below must be entered directly on the pump (see 709 Instructions for Use).



⇒ Press the <PARAM> key.

```
PARAM
>PARAM/709 IC Pump
```

The submenu "709 IC Pump" contains various inquiries concerning the 709 IC Pump.



⇒ Press the <ENTER> key to move to these inquiries.

```
>PARAM/709 IC Pump
fl ow                0.5 mL/mi n
```

The first inquiry concerns the **flow** of the 709 IC Pump, which can be set to a value between 0.01 and 5.00 mL/min using the numeric keys.

```
>PARAM/709 IC Pump
fl ow                2.0 mL/mi n
```

For the PRP-X100 IC Anion Column used in the example, the recommended flow rate is 2 mL/min. Enter this value using the numeric keys.



⇒ Confirm the inputted value using the <ENTER> key.

```
>PARAM/709 IC Pump
Pmax                10.0 MPa
```

The parameter "Pmax" denotes the **maximum shutoff pressure** for the 709 IC Pump and is used to protect the separating column against excessive pressure. If the pressure exceeds this value, the pump is automatically shut off. Use the numeric keys to enter the desired value here. This limit value should be ca. 3 MPa above the normal operating pressure of the separating column used or correspond to the maximum admissible operating pressure of the column (1 MPa = 10 bar).

>PARAM/709 IC Pump
Pmax 10.0 MPa

The maximum admissible pressure for the PRP-X100 IC Anion Column is 34 MPa, the normal operating pressure is 2 mL/min at ca. 7 MPa. A value of 10 MPa is thus entered as the maximum shutoff pressure.



⇒ Confirm the set maximum shutoff pressure with the <ENTER> key.

>PARAM/709 IC Pump
Pmin 0.0 MPa

The parameter "Pmin" denotes the **minimum shutoff pressure** for the 709 IC Pump. If the pressure drops below this value for a considerable length of time (e.g. owing to leaks or if the eluent supply is interrupted), the pump is automatically switched off.

Use the numeric keys to enter here the desired value. This limit value should be sufficiently far below the normal operating pressure of the separating column used (1 MPa = 10 bar).

>PARAM/709 IC Pump
Pmin 1.0 MPa

A value of 1 MPa can be entered as the minimum shutoff pressure for the PRP-X100 IC Anion Column.



⇒ Confirm the set minimum shutoff pressure using the <ENTER> key.

>PARAM/709 IC Pump
flow corr. 1.00

The parameter "flow corr." denotes the **correction factor for the flow** of the 709 IC Pump. This correction factor is used to match the displayed flow to the actual flow. If you are interested in the exact display of the flow, this correction factor must be determined by measurement of the actual flow and then entered here (see *section 3.3 of the 709 Instructions for Use*).



⇒ Confirm the inputted flow correction with the <ENTER> key. As this was the last inquiry of the last submenu of the <PARAM> key, an automatic switch is made to the standby mode.

151.2 mS/cm 8.5 min
method DEFAULT

The 732 IC Detector is again in the standby mode and the display shows the status messages.

10 Start 709 IC Pump



If you do not operate the 709 IC Pump by remote control via the 732 IC Detector, the pump must be switched on directly using the <PUMP RS> key (see 709 Instructions for Use).



⇒ Press the <PUMP R/S> key. The pump drive of the 709 IC Pump is started.

151.2 mS/cm	8.9 min
pump running	6.9 MPa

After the start of the 709 IC Pump, an automatic switch is made to the status message for the pump. If this is running properly, the message "pump running" and the current measured pressure appear.

11 Condition IC system

⇒ Rinse IC system with eluent until the desired stability of the baseline is achieved.

To assess the baseline stability, it is advantageous if the measured conductivity is continuously recorded with a recorder or printer or (e.g. with «IC Metrodata for Win95») shown on the PC.

It normally takes 30...60 min until the IC system is ready for analysis. If the eluent is changed, the establishment of the ion exchanger equilibrium on the separating column may take longer.

151.2 mS/cm	53.2 min
pump running	6.9 MPa

Example: With the PRP-X100 IC Anion Column a stable conductivity value of ca. 150 μ S/cm is reached after around 45 min and scarcely changes after this time.

3.4 Calibration

After putting into operation and conditioning the IC system as described in *section 3.3*, the first calibration can now be performed. This requires a standard solution which contains the analyte substances in about the same concentration as expected in the sample.

For our example dealing with the determination of drinking water with the PRP-X100 IC Anion Column, a 100 µL sample loop is used which is filled with the following standard solution:

5 ppm Cl⁻, 10 ppm NO₃⁻, 10 ppm SO₄²⁻
(as Na⁺ or K⁺ salts in dist. water)

1 Select status message for full-scale range



⇒ Repeatedly press the <SELECT> key until the status message for the full-scale range appears.

151.2 mS/cm	54.9 min
full scale	200 mS/cm

The **full-scale range** currently in force appears in the bottom line of the display.

2 Set full-scale range




⇒ Press the <FULL SCALE> key.

151.2 mS/cm	55.2 min
full scale:	<u>200 mS/cm</u>

The numeric value of the full-scale range start to flash. This indicates that the setting here must be made with the <←> and <→> keys and not as is normal with the <SELECT> key.

The **full-scale range** (operating range) is used to set the desired sensitivity for the display and the analog output of the measuring signal during the recording of the chromatogram. The range should be so selected that the largest peaks which need to be evaluated are clearly within the full-scale range.



⇒ Change the full-scale range by pressing the <←> key (value decremented) or the <→> key (value incremented) until the desired value appears in the display.

151.2 mS/cm	55.8 min
full scale:	<u>4.00 mS/cm</u>

For the determination of anions in drinking water, a full-scale range of 4 µS/cm is selected.



⇒ Confirm the desired full-scale range with the <ENTER> key.

151.2 mS/cm	56.1 min
full scale:	4.00 mS/cm

The 732 IC Detector is again in the standby mode and the display shows the status message for the full-scale range, but now the numeric value no longer flashes.



As the selected full-scale range is clearly less than the measured absolute conductivity, red overload display **8** on the 732 IC Detector now lights up.

3 Trigger auto-zero



⇒ Press the <ZERO> key. This starts the automatic electronic background compensation. During zeroing, the green LED in the key starts to flash.

+0.000 mS/cm	56.4 min
full scale: 4.00 mS/cm	

As soon as zeroing is at an end, the green LED in the <ZERO> key is permanently on. The top line of the display shows the **auto-zero value** in place of the absolute conductivity. In contrast to the absolute conductivity, the auto-zero value is always outputted with the sign in front of the number.

With the **auto-zero function** the current conductivity value becomes the new zero point of the selected full-scale range. Each time the <ZERO> key is pressed, the measuring signal is always reset to zero. The auto-zero function can be switched off with the <ZERO OFF> key.

4 Set injection valve A to "FILL"



⇒ Press the <FILL> key for injection valve A at the 733 IC Separation Center. This switches the injection valve to the "FILL" position. This position is indicated by lighting up of the green LED in the <FILL> key.

5 Fill sample loop

- ⇒ Immerse the aspirating tubing **88** attached to connection **22** in the standard solution.
- ⇒ Using the syringe fixed to connection **21** siphon in ca. 1 mL standard solution.

6 Set injection valve A to "INJECT"



⇒ Press the <INJECT> key for injection valve A on the 733 IC Separation Center and the chromatogram is started. This switches the injection valve to the "INJECT" position. The "INJECT" position is indicated by lighting up of the green LED in the <INJECT> key.

⇒ At the same time start the recorder (if used).


An integrator attached to the 733 IC Separation Center or a PC with evaluation software (e.g. «IC Metrodata for Win95») is started automatically when injection valve A is switched to the "INJECT" position.

+0.000 mS/cm 0.0 min
 full scale: 4.00 mS/cm

As soon as injection valve A is in the "INJECT" position, the current time is reset to "0.0 min". This marks the start time of the chromatogram.

The ion chromatogram of the standard is plotted in the next 20 min or so.

Fig. 19 shows an example of a calibration with 5 ppm Cl⁻, 10 ppm NO₃⁻ and 10 ppm SO₄²⁻.



FULL REPORT

Ret Time (Min)	Component Name	Concentr. ppm	Area (uV*Sec)	Height (uV)
1.758	Chloride	5.000000	2791676.00	416181.718
3.164	Nitrate	10.000000	2405831.75	202204.046
10.575	Sulphate	10.000000	4652037.50	146991.609

Fig. 19: Ion chromatogram of the calibration

3.5 Sample determination

Following the calibration of the IC system as described in *section 3.4* the first sample solution can be injected.

1 Filter drinking water sample

⇒ Filter the drinking water sample with a 0.45 µm microfilter.

2 Trigger auto-zero



⇒ Press the <ZERO> key. This starts the automatic electronic background compensation. During zeroing, the green LED in the key starts to flash. As soon as zeroing is at an end, the green LED in the <ZERO> key is permanently on.

3 Set injection valve A to "FILL"



⇒ Press the <FILL> key for injection valve A on the 733 IC Separation Center. This switches the injection valve to the "FILL" position. This position is indicated by lighting up of the green LED in the <FILL> key.

4 Fill sample loop

- ⇒ Immerse the suction tubing **88** attached to connection **22** in the vessel containing the drinking water sample.
- ⇒ Use the syringe fastened to connection **21** to siphon in ca. 1 mL sample.

5 Set injection valve A to "INJECT"



⇒ Press the <INJECT> key for injection valve A on the 733 IC Separation Center. This switches the injection valve to the "INJECT" position and the chromatogram is started. The "INJECT" position is indicated by lighting up of the green LED in the <INJECT> key.

⇒ At the same time start the recorder (if used).


An integrator attached to the 733 IC Separation Center or a PC with evaluation software (e.g. «IC Metrodata for Win95») is started automatically when injection valve A is switched to the "INJECT" position.

+0.000 mS/cm	0.0 min
full scale:	4.00 mS/cm

As soon as injection valve A is in the "INJECT" position, the current time is reset to "0.0 min". This marks the start time of the chromatogram.

The ion chromatogram of the drinking water sample is recorded during the next 20 min or so.

Fig. 20 shows an example of a drinking water sample.



FULL REPORT

Ret Time (Min)	Component Name	Concentr. ppm	Area (uV*Sec)	Height (uV)
1.749	Chloride	5.487972	3064128.00	381971.937
3.160	Nitrate	9.647280	2320973.25	189443.562
10.323	Sulphate	7.142281	3322615.75	112246.203

Fig. 20: Ion chromatogram of the drinking water sample


In comparison with the standard, in the drinking water sample the system peak which appears after the sulfate peak is very large. The next sample can be started at the earliest 18 min following injection of the previous sample.

3.6 Storing as a method

The parameter settings used for the drinking water determination can be stored in the 732 IC Detector as a method and hence later recalled at any time. Proceed as follows to store the inputted settings under the name "Water":




⇒ Press the <METHOD> key to open the main menu for method management.



The title of the submenu "METHOD/recall" appears. This is used to load a method already stored in the working memory. Proceed further to the next submenu.




⇒ Press the <METHOD> key again.



The title of the submenu "METHOD/store" appears. This is used to store a method loaded in the working memory.




⇒ Press the <ENTER> key.



The inquiry for entry of the method name which can have a length of max. 8 characters appears. The name "Water" selected for this example is entered as follows:




⇒ Press the <→> key to start the text entry.



The first 8 letters of the alphabet "ABCDEFGH" appear with the first letter A flashing. The letter at the flashing position can now be selected with the <←> and <→> keys.




⇒ Press the <→> key repeatedly until the flashing letter "w" appears.



The first letter of the name is now the desired "w".



⇒ Press the <ENTER> key to confirm the first letter.



The second letter of the name starts to flash and can again be selected with the <←> and <→> keys.



⇒ Press the <→> key repeatedly until the letter "a" appears at the flashing position.

>METHOD/store
name: Wabcdefg



The second letter of the name is now the desired "a".
⇒ Press the <ENTER> key to confirm the second letter.

>METHOD/store
name: Waabcdef



The third letter of the name starts to flash and can be selected with the <←> and <→> keys.
⇒ Press the <→> key repeatedly until the letter "t" appears at the flashing position.

>METHOD/store
name: Watuvwxyz



The third letter of the name is now the desired "t".
⇒ Press the <ENTER> key to confirm the third letter.

>METHOD/store
name: Watuvw



The fourth letter of the name starts to flash can be selected with the <←> and <→> keys.
⇒ Press the <←> key repeatedly until the letter "e" appears at the flashing position.

>METHOD/store
name: Watefghi



The fourth letter of the name is now the desired "e".
⇒ Press the <ENTER> key to confirm the fourth letter.

>METHOD/store
name: Wateefgh



The fifth letter of the name starts to flash and can be selected with the <←> and <→> keys.
⇒ Press the <→> key repeatedly until the letter "r" appears at the flashing position.

>METHOD/store
name: Waterrstu



The fifth letter of the name is now the desired "r".
⇒ Press the <ENTER> key to confirm the fifth letter.

>METHOD/store
name: Waterrst



The desired name "water" has now been entered in full.
⇒ Press the <QUIT> key to end the text entry.

```
>METHOD/store
name:      Water
```



```
+0.000 mS/cm   28.4 min
full scale:   4.00 mS/cm
```

The display shows the complete method name "Water" for confirmation.

⇒ Press the <ENTER> key to confirm the method name. The parameters are now stored in the instrument under the name "Water".

The 732 IC Detector is now again in the standby mode and the display shows the status messages.

A method report of the method "Water" just stored can be outputted with the <REPORT> key (see *section 4.8.3*) and has the following appearance:

```
732 IC Detector      01104   732.0012
date 1995-09-18   time 11:03:47

METHOD
Method name      Water
date: 1995-09-18 10:39:10

PARAM
>PARAM/detector
range:           200 µS/cm
full scale:      4.00 µS/cm
temp.coeff.:     2.5 %/°C
>PARAM/analog output
polarity:        +
offset:          0 %fs
damping:         off
>PARAM/plot
autom.start:     off
time interval    1.0 s
time scale       10.0 mm/min
time scale label: rel
stop time        off min
left:            0.000 µS/cm
right:           10 µS/cm
>PARAM/709 IC Pump
flow             2.00 mL/min
Pmax             10.0 MPa
Pmin             1.0 MPa
flow corr.       1.00

PROGRAM
-----
```

4 Operation




This section provides a detailed description of the operation of the 732 IC Detector and 733 IC Separation Center using the keypad and dialog display. After an overview of the various operating sequences (section 4.1), the fundamentals of the operation are explained (section 4.2). This is followed by a detailed description of the display (section 4.3) and key functions (section 4.4 – 4.8). A few selected illustrative methods complete this section 4 (section 4.9).

4.1 Operating sequences


4.1.1 General flow chart

The following flow diagram shows the general sequence of the determination of an ion chromatogram. Additional detailed charts and in-depth information can be found in the sections mentioned.




4.1.2 Flow chart for basic settings

The following flow diagram shows the basic SETUP and CONFIG settings which must be selected for the 732 IC Detector as a function of the attached devices. You will find detailed information in *section 4.4*.




4.1.3 Flow chart for putting into operation

The following flow diagram shows how the IC system can be put back into operation after a separating column has been removed and all devices have been switched off (shutdown, see *section 5.2.4*). A requirement is that the installation (*section 2.1*) is complete and the basic settings (*section 4.1.2*) have been made. You will find further information in the sections mentioned.



4.1.4 Flow chart for injection

The following flow diagram shows how an injection is started on the IC system with and without a program. If the points <ZERO>, <FILL> and "Siphoning in sample" are not programmed or must be performed by the Autosampler, they must also be manually triggered. You will find further information in the sections mentioned.




4.2 Fundamentals of the operation

4.2.1 Display

LCD **1** on the 732 IC Detector comprises two lines each of 24 characters. What appears on the display depends on whether the instrument is in the **standby mode** or the **edit mode**.


Standby mode

After the 732 IC Detector has been switched on, the instrument is always automatically in the standby mode of conductivity measurement. The first display line shows the measured value and the current time, the second line displays status messages selectable with the <SELECT> key.




Edit mode

Pressing the appropriate keys switches from the standby mode to the edit mode. The first line of the display then shows the name of the pressed key, the second line the title of the first submenu:




Pressing the <ENTER> key switches from the menu title to the individual inquiries. The first line shows the menu title, the second line is used as an entry line for parameters.




4.2.2 Overview of key functions

Key field of 732 IC Detector



The key field of the 732 IC Detector contains 4 colored main function keys, 12 gray numeric keys and 4 gray auxiliary function keys. In the standby mode, the numeric keys can also be used to trigger the specified functions by means of their numbers.





Key field of 733 Separation Center









The key field of the 733 IC Separation Center contains 4 keys for operation of the injection valves and the suppressor module.







The following table provides a brief overview of the different functions of the individual keys of the 732 IC Detector and 733 Separation Center in the standby mode and the edit mode. You will find more detailed information on the key functions in *sections 4.4 ¼ 4.8*.

Main function keys of 732 IC Detector





Key	Standby mode	Edit mode
	Trigger auto-zero <ul style="list-style-type: none"> Automatic zeroing of the current conductivity (the green LED in the key lights up). 	Trigger auto-zero <ul style="list-style-type: none"> Automatic zeroing of the current conductivity (the green LED in the key lights up).
	Open parameter menu <ul style="list-style-type: none"> Opens the main menu for setting the parameters for the conductivity detector, analog output, graphics plot and 709 IC Pump. All settings which are set in the parameter menu can be stored together with the program as a method. 	Select submenu or parameters <ul style="list-style-type: none"> Selection of the next submenu in the main menu. Selection of the next parameter in the submenu.
	Set full-scale range <ul style="list-style-type: none"> Direct selection of the full-scale range (value flashes). Selection of the value with the < ← > and < → > keys. 	Return to standby mode <ul style="list-style-type: none"> Acceptance of the displayed value for the full-scale range and return to the standby mode.
	Start/stop program <ul style="list-style-type: none"> <i>LED in key dark:</i> Program inactive; No function (program start not possible until program status = "active"). <i>LED in key lit up:</i> Program active; Program start with program type "cycle". <i>LED in key flashes:</i> Program running; Program stop with all program types. 	No function

Numeric keys of 732 IC Detector





Key	Standby mode	Edit mode
	Open configuration menu <ul style="list-style-type: none"> Opens the main menu for configuration of the 732 IC Detector. The settings in the configuration menu are retained until they are changed or the working memory (RAM) is initialised. 	Numeric entry ('7') or select submenu or parameter <ul style="list-style-type: none"> Selection of the next submenu in the main menu. Selection of the next parameter in the submenu.
	Start/stop pump <ul style="list-style-type: none"> On/off switching of the pump drive of the 709 IC Pump (if pump remote controlled by 732 IC Detector). 	Numeric entry ('8')
	Open method menu <ul style="list-style-type: none"> Opens the main menu for the recall, storage and deletion of user-defined methods. 	Numeric entry ('9') or select submenu or parameter <ul style="list-style-type: none"> Selection of the next submenu in the main menu. Selection of the next parameter in the submenu.
	Switch off auto-zero <ul style="list-style-type: none"> Switches off the auto-zero function (green LED in the <ZERO> key goes out). 	Numeric entry ('4') or switch off auto-zero <ul style="list-style-type: none"> Switches off the auto-zero function (green LED in the <ZERO> key goes out).
	Open event menu <ul style="list-style-type: none"> Opens the main menu for the programming, editing and deletion of events. 	Numeric entry ('5') or select submenu or parameter <ul style="list-style-type: none"> Selection of the next submenu in the main menu. Selection of the next parameter in the submenu.
	Open program menu <ul style="list-style-type: none"> Opens the main menu for the development, editing and deletion of time programs. Definition of program type. Set programs active/inactive. 	Numeric entry ('6') or select submenu or parameter <ul style="list-style-type: none"> Selection of the next submenu in the main menu. Selection of the next parameter in the submenu.

Key	Standby mode	Edit mode
	Start graphics plot <ul style="list-style-type: none"> Output of measured value curve on an external printer (abort with <QUIT>). 	Numeric entry ('1')
	Start report output <ul style="list-style-type: none"> Selection of report and output on an external printer (abort with <QUIT>). 	Numeric entry ('2') or return to standby mode
	Start measured value output <ul style="list-style-type: none"> Output of single measured values or start of the continuous measured value output on an external printer (abort with <QUIT>). 	Numeric entry ('3')
	Trigger marking <ul style="list-style-type: none"> Sets a marking signal of ca. 10% of the full-scale range at the analog output. 	Numeric entry ('0')
	No function	Numeric entry ('.') or narrowing of the full-scale range (see <FULL SCALE>) or text entry <ul style="list-style-type: none"> Activation of the text entry mode. Shift of the moving character string to the left (see <i>section. 4.2.5</i>).
	No function	Change in the sign ('-/+') or widening of the full-scale range (see <FULL SCALE>) or text entry <ul style="list-style-type: none"> Activation of the text entry mode. Shift of the moving character string to the right (see <i>section 4.2.5</i>).

Auxiliary function keys of 732 IC Detector

Key	Standby mode	Edit mode
	Select status message <ul style="list-style-type: none"> Selection of the status messages for full-scale range, absolute conductivity, 709 IC Pump, program, method, date and time. 	Select preset parameter values <ul style="list-style-type: none"> Selection from preset parameter values for parameters which are marked by a colon " : ".
	No function	Delete parameter values <ul style="list-style-type: none"> Deletion of the displayed parameter values and resetting to default values or "0".
	Acknowledge error messages <ul style="list-style-type: none"> Acknowledgment of the error messages displayed on the first line. or abort printing operations <ul style="list-style-type: none"> Abort printing operations started with <PLOT>, <REPORT> or <PRINT>. 	Abort entry <ul style="list-style-type: none"> Abort the parameter entry (parameter reset to original value). Exit from rolling inquiries and selection of the next higher menu level or the standby mode.
	No function	Confirmation of parameters <ul style="list-style-type: none"> Confirmation of existing parameter values or those just entered. Selection of the next menu line.

Operating keys of 733 IC Separation Center

Key	Injection valve	Suppressor module
	Switching to "FILL" <ul style="list-style-type: none"> Switching of injection valve A to "FILL" position. 	No function
	Switching to "INJECT" <ul style="list-style-type: none"> Switching of injection valve A to "INJECT" position. 	No function
	Switching to "FILL" <ul style="list-style-type: none"> Switching of injection valve B to "FILL" position. 	Rotate suppressor module <ul style="list-style-type: none"> Switch suppressor module to next position.
	Switching to "INJECT" <ul style="list-style-type: none"> Switching of injection valve B to "INJECT" position. 	No function

4.2.3 Instrument dialog

The instrument dialog of the 732 IC Detector is organized in the form of so-called rolling inquiries which are arranged in menu levels in a hierarchical manner and are subject to the following rules:

Main menu

The main function keys as well as most numeric keys of the 732 IC Detector open a main menu whose thematically arranged submenus are selected by repeated pressing of this key. The name of the key always appears in the top line of the display.

Submenu

Each submenu has a title marked by ">" which appears in the bottom line of the display. <ENTER> is used to move from the title to the individual inquiries, which are used to change the most important settings of the instrument. With <QUIT> a switch is made back to the standby mode.

Inquiries

In inquiries without " : ", the values must be entered using the numeric keys. <ENTER> is used to accept the inputted value and the next inquiry appears.

In the case of inquiries with " : ", the admissible values must be selected with the <SELECT> key. <ENTER> is used to accept the set value and the next inquiry appears.

Depending on the parameter, <CLEAR> is used to reset the displayed value to the smallest possible value or the default value. The <CLEAR> key also serves to abort wrong entries.

<QUIT> is used to quit the inquiries and return to the main menu.

You will find a schematic representation of the instrument dialog in *Fig. 21*.




Fig. 21: Schematic representation of the instrument dialog

4.2.4 Data entry

Numeric entry




With inquiries in which the parameter title has no colon ":", the parameter values can be entered directly with the numeric keys.

<Select> selection



With inquiries in which the parameter title has a colon ":", the <SELECT> key can be used to display a preset selection of data. This <SELECT> selection has a cyclic structure like a selector drum.

Example:



Confirmation



Confirmation of the inputted or selected parameter value by pressing the <ENTER> key.

Correction



Deletion of the parameter value just entered or selected and resetting to the default value or "0".



Abort of the parameter entry with <QUIT> (e.g. as a wrong value has been entered). The parameter value is reset to the original value, a beep sounds and a direct return to the submenu is effected.

4.2.5 Text entry

With certain inquiries which require the entry of a text, any ASCII characters can be entered to write the texts.

Opening the text editor

The "←" or "→" keys open the text editor.



"←" is used to delete an existing character string and to position the text cursor at the left edge of the entry field.



With "→" an existing character string is retained, the text cursor is positioned on the last character of the existing text.

After the text editor has been opened, a moving string is displayed which comprises an alphabetical list of all characters which can be entered. The character which flashes is that which is currently selected and can be changed (text cursor).

Character selection



The "←" and "→" keys move the string of the selectable characters (uppercase and lowercase letters, numbers and special characters, arranged alphabetically) in the selected direction past the text cursor.



Pressing these keys once causes a shift in the string by one position to the left or right. By pressing and holding the keys, the string can be moved rapidly.

Numbers can also be entered directly using the numeric keys.

Confirmation of the character selection



The <ENTER> key appends the character currently at the text cursor position to the existing text line. When the entire width of the text field is full, the text entry mode is quit and the text line accepted with <ENTER>.

Deleting characters



The <CLEAR> key is used to delete the rearmost character of the existing text line. As a result, the text cursor automatically moves one character to the left.

Closing text entry



The text entry mode is exited with <QUIT>. The displayed text line can then be accepted with <ENTER> or discarded by pressing <QUIT> again.

You will find an example illustrating text entry in *section 3.6*.

4.3 Displays in the standby mode

4.3.1 Measured value and current time

In the standby mode of the conductivity measurement, the top line of the display always shows the current measured value (left) and the current time (right).

151.2 mS/cm method	13.5 min DEFAULT
-----------------------	---------------------

Measured value

The following possibilities exist for the display of the measured value:

Measured value without sign

151.2 $\mu\text{S/cm}$ Absolute conductivity in $\mu\text{S/cm}$ or mS/cm.

This value is also shown automatically if the auto-zero value lies outside the full-scale range.

Measured value with sign

+0.004 $\mu\text{S/cm}$ Auto-zero value in $\mu\text{S/cm}$ (default setting)

+0.0 %fs Auto-zero value in %fs (% of the full-scale range, can be set under ">CONFIG/detector/'zero' unit")

+0.4 mV Auto-zero value in mV (can be set under ">CONFIG/detector/'zero' unit")

No measured value

----- Measured value is more than 10% outside the set measuring range.

151.2 mS/cm method	13.5 min DEFAULT
-----------------------	---------------------

Current time

The following possibilities exist for the display of the current time:

Without "INJECT A" or program start

Current time after instrument switched on

With "INJECT A"

Current time after the last switching of injection valve A to "INJECT" position (does not apply to an "INJECT A" command within a program)

With program start

Current time from last program start

4.3.2 Status messages



Pressing the <SELECT> key in the standby mode selects the status messages in turn. These are shown in the bottom line of the display. Depending on the instrument status, the following messages can be displayed:

+0.004 mS/cm	13.5 min
method	XXXXXXXX

Status messages for a method

- method XXXXXXXX
An unchanged method is loaded in the working memory (see *section 4.7.4*).
- method XXXXXXXX modif.
The method loaded in the working memory has been changed (see *section 4.7.4*).
- method DEFAULT
The default method "DEFAULT" is loaded in the working memory (see *section 4.7.4*).

+0.004 mS/cm	13.5 min
1995-10-12	11:24:31

Display of date and time

The current date and the current time are displayed. Date and time can be changed under <CONFIG> (see *section 4.4.2*).

151.2 mS/cm	13.5 min
full scale	1.00 mS/cm

Display of the full-scale range

The full-scale range currently set is displayed (operating range, see *section 4.5.1*) and can be changed with the <FULL SCALE> key (see *section 4.5.2*).

+0.004 mS/cm	13.5 min
abs. cond.	151.2 mS/cm

Display of the absolute conductivity

The absolute conductivity currently measured is displayed. This is primarily useful when the auto-zero value is displayed as the measured value.

+0.004 mS/cm	13.5 min
pump ready	

Status messages for 709 IC Pump

These status messages appear only with "SETUP/peripherals/operation with 709 = on".

- pump ready
The 709 IC Pump is switched on and ready for start with the <PUMP R/S> key (see *section 4.6.6*).
- pump running XX.X MPa
The 709 IC Pump is running, the current pressure P_{actual} is shown on the right. The pump can be stopped using the <PUMP R/S> key (see *section 4.6.6*).

stopped: YY-MM-DD HH:MM

The 709 IC Pump was stopped at the time indicated owing to violation of the shutoff limit "Pmin" or "Pmax". To restart a pump stopped in this manner, the pump keys <EXT.> and <R/S> must be pressed in succession and then <EXT.> pressed again (see also section 4.5.1).

pump not responding

This message appears when

- the 709 IC Pump is not connected to the 732 IC Detector,
- the 709 IC Pump is not switched on,
- the remote control of the 709 IC Pump is not switched on.

+0.004 mS/cm	13.5 min
prog. type	cycle no. XXX

Status messages for program type

These status messages appear only with "PROGRAM/parameters/status = active".

prog. type cycle no. XXX

Display of the program type "cycle" and the cycle number (see section 4.7.1).

prog. type inject

Display of the program type "inject" (see section 4.7.1).

prog. type remote

Display of the program type "remote" (see section 4.7.1).

+0.004 mS/cm	13.5 min
next step	ready

Status messages for program step

These status messages appear only with "PROGRAM/parameters/status = active".

next step ready

The program is ready for a new start (see section 4.7.1).

next step XXX.X min

Display of the time when the next program step will be executed.

next step ---

The program has been set to "inactive", a new start is not possible.

next step ∞

The program continues to run indefinitely (no "end" command available).

4.4 Basic settings

4.4.1 Setup

Switch on +



Under "Setup", various basic settings of the 732 IC Detector are combined which normally seldom need to be changed. They are accessible only when the <CONFIG> key is pressed at the same time as the instrument is switched on. The "SETUP" part which then appears comprises the following main menu:

```
SETUP
>SETUP/output
```

Setting of the remote output lines

```
SETUP
>SETUP/i nput assi gn
```

Assignment of the remote input lines

```
SETUP
>SETUP/graphi cs
```

Setting of general graphics parameters

```
SETUP
>SETUP/peri pheral s
```

Setting of the peripherals 733 and 709

Repeated pressing of the <CONFIG> key selects the submenus in turn. The individual inquiries of a submenu are accessed using the <ENTER> key and exited with the <QUIT> key. The following listings show all dialog items which appear under "SETUP". The values shown in the displays are the default values, the possible entry values or ranges are shown below the display.

```
>SETUP
>SETUP/output
```

Setting of the remote output lines

For further details on the remote interface, see *section 6.2*.

```
>SETUP/output
remote      00000000
```

Set the remote output lines 1...8

Sets the basic setting for the remote output lines when the instrument is switched on.

Each of the remote output lines 1...8 can be set manually from left to right in turn:

- 0 Line off, inactive (open)
- 1 Line on, active (0 V)

Each time the instrument is switched on, the remote output lines are reset in accordance with the parameters entered here.

0, 1
<CLEAR> sets all values to 0

<p>>SETUP >SETUP/i nput assi gn</p>	<p>Assignment of the remote input lines</p> <p>The 732 IC Detector has 4 programmable remote input lines which can be used to trigger a total of 15 different instrument functions (for details, see <i>section 6.2.1</i>). The default assignment of the individual functions to the status of the 4 input lines 2...5 (defined by the corresponding decimal code 1...15) can be changed if need be. Several functions can thus also be assigned the same status.</p>
<p>>SETUP/i nput assi gn PROG R/S 1</p> <p>1...15</p>	<p>Start/stop program</p> <p>An active remote program is started, any programs running are stopped.</p>
<p>>SETUP/i nput assi gn PUMP R/S 2</p> <p>1...15</p>	<p>Start/stop 709 IC Pump</p> <p>Equivalent to pressing the <PUMP R/S> key.</p>
<p>>SETUP/i nput assi gn FILL A 4</p> <p>1...15</p>	<p>Trigger "FILL" at valve A</p> <p>Equivalent to pressing the <FILL A> key on the 733 IC Separation Center.</p>
<p>>SETUP/i nput assi gn INJECT A 8</p> <p>1...15</p>	<p>Trigger "INJECT" at valve A</p> <p>Equivalent to pressing the <INJECT A> key on the 733 IC Separation Center.</p>
<p>>SETUP/i nput assi gn FILL B/STEP 3</p> <p>1...15</p>	<p>Trigger "FILL" at valve B or "STEP" on the suppressor module</p> <p>Equivalent to pressing the <FILL B> key or the <STEP> key on the 733 IC Separation Center.</p>
<p>>SETUP/i nput assi gn INJECT B 12</p> <p>1...15</p>	<p>Trigger "INJECT" at valve B</p> <p>Equivalent to pressing the <INJECT B> key on the 733 IC Separation Center.</p>
<p>>SETUP/i nput assi gn REPORT 5</p> <p>1...15</p>	<p>Report output</p> <p>Equivalent to pressing the <REPORT> key.</p>
<p>>SETUP/i nput assi gn ZERO 6</p> <p>1...15</p>	<p>Switch on auto-zero</p> <p>Equivalent to pressing the <ZERO> key.</p>
<p>>SETUP/i nput assi gn MARK 7</p> <p>1...15</p>	<p>Switch on marking</p> <p>Equivalent to pressing the <MARK> key.</p>

>SETUP/i nput assi gn ZERO OFF 9	Switch off auto-zero Equivalent to pressing the <ZERO OFF> key.
1...15	
>SETUP/i nput assi gn PLOT 10	Plot output Equivalent to pressing the <PLOT> key.
1...15	
>SETUP/i nput assi gn pol ari ty 11	Switch the polarity The polarity of the signal at the analog output is switched.
1...15	
>SETUP/i nput assi gn SELECT 13	Select Equivalent to pressing the <SELECT> key.
1...15	
>SETUP/i nput assi gn QUIT 14	Quit Equivalent to pressing the <QUIT> key.
1...15	
>SETUP/i nput assi gn ENTER 15	Enter Equivalent to pressing the <ENTER> key.
1...15	

>SETUP >SETUP/graphi cs	General graphics parameters General parameters for the printout of graphics plots on an external printer.
----------------------------	---

>SETUP/graphi cs gri d: off	Grid lines for graphics plot on Dotted grid lines are plotted on the graphics printout. off No grid lines are plotted on the graphics printout.
on, off	

>SETUP/graphi cs frame: off	Frame for graphics plot on A frame is plotted on the graphics printout. off Only the y and x axes are plotted on the graphics printout, but no frame.
on, off	

>SETUP/graphi cs wi dt h 0.8	Relative width of the graphics printout The width of the graphics printout must be matched to the attached printer.
0.4...1.0	

>SETUP
>SETUP/peripherals

Setting the peripherals

>SETUP/peripherals
operation with 733: on
on, off

Operation with 733 IC Separation Center

- on The 732 IC Detector is operated with the 733 IC Separation Center (normal case).
- off The 732 IC Detector is operated without the 733 IC Separation Center. All parameters which affect the 733 IC Separation Center are suppressed in the dialog.

>SETUP/peripherals
operation with 709: on
on, off

Operation with 709 IC Pump

- on The 732 IC Detector is operated with the 709 IC Pump (normal case).
- off The 732 IC Detector is operated without the 709 IC Pump. All parameters which affect the 709 IC Pump are suppressed in the dialog. In this case, the RS interface "709 IC Pump" can be used for the connection of an external printer.

4.4.2 Configuration, <CONFIG> key



The <CONFIG> key is used for the entry of basic settings of general applicability which are also retained when the instrument is switched off. The key opens the following main menu:

```
CONFIG
>CONFIG/detector
```

Settings for the conductivity detector

```
CONFIG
>CONFIG/printer
```

Settings for external printer

```
CONFIG
>CONFIG/print meas. value
```

Settings for measured value printout

```
CONFIG
>CONFIG/auxiliaries
```

Various general instrument settings

```
CONFIG
>CONFIG/RS settings
```

Settings for interface "RS 232"

```
CONFIG
>CONFIG/RS settings 709
```

Settings for interface "709 IC Pump" (RS232)

```
CONFIG
>CONFIG/733 IC Sep. Cent.
```

Settings for 733 IC Separation Center

Repeated pressing of the <CONFIG> key selects the submenus in turn. The individual inquiries of a submenu are accessed using the <ENTER> key and exited with the <QUIT> key. The following listings show all dialog items which appear under <CONFIG>. The values shown in the display are the default values, the possible entry values or ranges are shown below the display.

```
>CONFIG
>CONFIG/detector
```

Settings for the conductivity detector

```
>CONFIG/detector
thermostat:      35 °C
```

25, 30, 35, 40, 45 °C, off

Thermostatting of the conductivity detector

Setting of the operating temperature of the conductivity cell.

At constant ambient temperature the heating built into the 732 IC Detector regulates the selected measuring cell temperature to an accuracy of ± 0.01 °C and thus establishes the precondition for highly sensitive determinations. It normally takes 30..60 min after the instrument has been switched on until this temperature stability is attained.

When the heating is switched on, the green display lamp **7** "THERMOSTAT" is on. Thermostatting of the conductivity detector is switched off by pressing the <CLEAR> key. The display shows the value "off °C", at the same time the green display lamp **7** goes out.



Thermostatting functions only if the ambient temperature is at least 5 °C lower than the operating temperature.

>CONFIG/detector
'zero' unit: mS/cm

 μS/cm, %fs, mV

Unit for display of the auto-zero value

The following units can be selected for the auto-zero value displayed in the instrument standby mode (see section 4.6.3):

- μS/cm Display in conductivity unit μS/cm.
- %fs Display in % of the full-scale range which has been set under the <PARAM> key (see section 4.5.1) or <FULL SCALE> key (see section 4.5.2).
- mV Display in mV.

>CONFIG/detector
cell constant 16.7/cm

 13.0...21.0 /cm

Cell constant of the conductivity cell

Setting of the cell constant of the measuring cell for correct display of the absolute conductivity.

In ion chromatography, interest is normally centered only on relative changes in the conductivity and not on the absolute value. With the cell constant of 16.7/cm set in the factory, the error in the display of the absolute conductivity is ca. ±10%.

If you wish to have a more accurate display of the absolute conductivity, the cell constant must be determined by means of a calibration solution. To do this, pump a solution of known conductivity through the IC system, observe the display of the absolute conductivity and change the cell constant until the displayed value matches the actual value.

>CONFIG
>CONFIG/printer

Settings for external printer

>CONFIG/printer
id. 1

 18 ASCII characters

Identification 1 for 1st line of print header

Freely selectable ASCII character string for the first line of the print header (text entry, section 4.2.5).

>CONFIG/printer
id. 2

18 ASCII characters

Identification 2 for 2nd line of print header

Freely selectable ASCII character string for the second line of the print header (text entry, see *section 4.2.5*).

>CONFIG/printer
print header: once

once, always, off

Printout of the print header

The print header comprises title line (with device name, serial number and program number), identification 1 and 2 as well as date and time if selected and is printed out before a measured value, report or curve plot.

- once Printout of the print header only once after the instrument has been switched on before each measured value, report or curve plot.
- always Printout of the print header before every measured value, report or curve plot.
- off No printout of the print header.

>CONFIG/printer
date&time: on

on, off

Printout of date and time in the print header

- on Date and time are printed out in the print header.
- off Date and time are not printed out in the print header.

>CONFIG/printer
send to: IBM

IBM, Epson, Seiko,
Citizen, HP

Selection of the character set/printer type

- IBM Character set 437 for IBM PC or printer type IBM Proprinter
- Epson EPSON printers
- Seiko Seiko printer DPU-411/414
- Citizen Citizen printer IDP562 RS
- HP HP printers (Deskjet..., Laserjet... etc.)

You will find further details on the printer connection in *section 2.9.5*.

>CONFIG
>CONFIG/print meas. value

Settings for measured value printout

Output of measured values via RS232 interface.

>CONFIG/print meas. value
print crit.: immed.

immed., time, off

Criterion for measured value printout

- immed. Measured value printout each time the <PRINT> is pressed.
- time Automatic measured value printout at selectable time intervals.
- off No measured value printout.

<pre>>CONFIG/print meas. value time interval 1.0 s</pre>	<p>Time interval for measured value printout</p> <p><i>This inquiry appears only with "print crit. = time".</i></p> <p>Time interval in seconds between printout of the individual measured values.</p>
0.4...99999 s	

<pre>>CONFIG/print meas. value stop time off min</pre>	<p>Stop time for measured value printout</p> <p><i>This inquiry appears only with "print crit. = time".</i></p> <p>Time to stop of the measured value printout.</p>
1...999 min, off	<p>off No time restriction (press <CLEAR>).</p>

<pre>>CONFIG/print meas. value date&time: off</pre>	<p>Printout of date and time</p> <p><i>This inquiry appears only with "print crit. = immed." or "print crit. = time".</i></p>
on, off	<p>on Date and time are printed out for every measured value.</p> <p>off Date and time are not printed out.</p>

>CONFIG
 >CONFIG/auxiliaries

General instrument settings

<pre>>CONFIG/auxiliaries run number 0</pre>	<p>Run number</p> <p>The run number is incremented by +1 for each new determination as follows:</p> <p style="margin-left: 40px;"><i>without program</i> Increase with every "INJECT A".</p> <p style="margin-left: 40px;"><i>with program</i> Increase with every new start and with every "return" in cycle programs.</p> <p>The parameter value has the following meaning:</p>
0...999, off	<p>0...999 Start point for numbering.</p> <p>off No numbering (press <CLEAR>).</p>

<pre>>CONFIG/auxiliaries number of cycles 1</pre>	<p>Number of cycles</p> <p>This parameter is identical to the parameter of the same name under the <PROGRAM> key (see <i>section 4.7.1</i>) and determines how many times cycle programs should be run. As access to the <PROGRAM> key is blocked while a program is running, the number of cycles can subsequently be changed here.</p>
1...999	

>CONFIG/auxiliaries >CONFIG/aux/event	<h3>Settings for events</h3>								
<i>Submenu</i>									
>CONFIG/aux/event enabled in program: off on, off	<h3>Perform events in a program</h3> <p>on Events are performed during a running program.</p> <p>off Events are <u>not</u> performed during a running program.</p>								
>CONFIG/auxiliaries date YY-MM-DD YYYY: 1995...2094; MM: 01...12; DD: 01...31	<h3>Date</h3> <p>Current date with numeric data for year (YYYY), month (MM) and day (DD).</p>								
>CONFIG/auxiliaries time HH:MM:SS HH: 00...23 MM: 00...59 SS: 00...59	<h3>Time</h3> <p>Current time with numeric data for hours (HH), minutes (MM) and seconds (SS). If a new time is entered, it becomes active when the <ENTER> key is pressed.</p>								
>CONFIG/auxiliaries dialog: english english, deutsch, français, español	<h3>Dialog language</h3> <table border="0"> <tr> <td>english</td> <td>English</td> </tr> <tr> <td>deutsch</td> <td>German</td> </tr> <tr> <td>français</td> <td>French</td> </tr> <tr> <td>español</td> <td>Spanish</td> </tr> </table>	english	English	deutsch	German	français	French	español	Spanish
english	English								
deutsch	German								
français	French								
español	Spanish								
>CONFIG/auxiliaries device label 8 ASCII characters	<h3>Device label</h3> <p>Freely selectable ASCII character string for the device label (text entry, see <i>section 4.2.5</i>).</p>								
>CONFIG/auxiliaries program 732.0012 <i>display only</i>	<h3>Number of the program version</h3> <p>Display only (no entry possibility). Please specify this number in inquiries to Metrohm.</p>								
>CONFIG/auxiliaries >CONFIG/aux/beep	<h3>Settings for the beeper</h3>								
<i>submenu</i>									
>CONFIG/aux/beep status: on on, only error, off	<h3>Status of the beeper</h3> <p>on Single beep on entry errors and on program return, three beeps with error messages.</p> <p>only error Beep (three times) only with error messages.</p> <p>off No beep.</p>								


```
>CONFIG
>CONFIG/RS settings 709
```

Settings for RS interface "709 IC Pump"

This RS232 interface is normally used for connection of the 709 IC Pump. However, it can also be used to attach an external printer if the instrument is appropriately configured in the Setup program ("operation with 709 = off", see *section 4.4.1*). For detailed information on the RS232 interface, see *section 6.1*.

```
>CONFIG/RS settings 709
connected: printer
```

```
printer, off
```

Connection to RS interface

This inquiry appears only with "operation with 709 = off".

printer Connection of a printer (the RS232 parameters are displayed).

off No device connected (no RS232 parameters appear).

```
>CONFIG/RS settings 709
baud rate: 9600
```

```
9600, 4800, 2400,
1200, 600, 300
```

Data transmission rate (baud rate)

This inquiry appears only with "operation with 709 = on" or "operation with 709 = off" and "connected = printer".

Data transmission rate in bit/s

```
>CONFIG/RS settings 709
data bit: 8
```

```
7, 8
```

Data bits

This inquiry appears only with "operation with 709 = off" and "connected = printer".

```
>CONFIG/RS settings 709
stop bit: 1
```

```
1, 2
```

Stop bits

This inquiry appears only with "operation with 709 = off" and "connected = printer".

```
>CONFIG/RS settings 709
parity: none
```

```
none, odd, even
```

Parity

This inquiry appears only with "operation with 709 = off" and "connected = printer".

none Parity is not checked.

odd Odd parity.

even Even parity.

```
>CONFIG/RS settings 709
handshake: Hws
```

```
Hws, Hwf,
SWchar, SWline, none
```

Handshake

This inquiry appears only with "operation with 709 = off" and "connected = printer".

Hws Reduced hardware handshake.

Hwf Full hardware handshake.

SWchar Software handshake with character stop.

SWline Software handshake with line stop.

none No handshake.

For detailed information on the handshake, see *section 6.1.8*.

>CONFIG >CONFIG/733 IC Sep. Cent.	<h3>Settings for 733 IC Separation Center</h3> <p><i>This submenu appears only with "operation with 733 = on".</i></p>
--------------------------------------	--

>CONFIG/733 IC Sep. Cent. valve A	<h3>Configuration 733 IC Separation Center</h3> <p>Depending on the instrument configuration, the following messages appear:</p> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="padding: 2px;">valve A</td> <td style="padding: 2px;">733.0010 (1 injector).</td> </tr> <tr> <td style="padding: 2px;">valve A + valve B</td> <td style="padding: 2px;">733.0020 (2 injectors).</td> </tr> <tr> <td style="padding: 2px;">valve A + suppressor</td> <td style="padding: 2px;">733.0030 (1 injector + 1 suppressor module).</td> </tr> <tr> <td style="padding: 2px;">733 not accessible</td> <td style="padding: 2px;">733 not connected or cable connection faulty.</td> </tr> </table>	valve A	733.0010 (1 injector).	valve A + valve B	733.0020 (2 injectors).	valve A + suppressor	733.0030 (1 injector + 1 suppressor module).	733 not accessible	733 not connected or cable connection faulty.
valve A	733.0010 (1 injector).								
valve A + valve B	733.0020 (2 injectors).								
valve A + suppressor	733.0030 (1 injector + 1 suppressor module).								
733 not accessible	733 not connected or cable connection faulty.								

display only

>CONFIG/733 IC Sep. Cent. control: no restriction	<h3>Valve control of the 733 IC Separation Center</h3> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="padding: 2px;">no restriction, 732 only</td> <td style="padding: 2px;">no restriction</td> <td style="padding: 2px;">Suppressor module and valves can be operated with the <FILL> and <INJECT> or <STEP> keys at 733 IC Separation Center and also via the 732 IC Detector.</td> </tr> <tr> <td style="padding: 2px;"></td> <td style="padding: 2px;">732 only</td> <td style="padding: 2px;">Suppressor module and valves can be operated only via 732 IC Detector (the keys on the 733 are locked).</td> </tr> </table>	no restriction, 732 only	no restriction	Suppressor module and valves can be operated with the <FILL> and <INJECT> or <STEP> keys at 733 IC Separation Center and also via the 732 IC Detector.		732 only	Suppressor module and valves can be operated only via 732 IC Detector (the keys on the 733 are locked).
no restriction, 732 only	no restriction	Suppressor module and valves can be operated with the <FILL> and <INJECT> or <STEP> keys at 733 IC Separation Center and also via the 732 IC Detector.					
	732 only	Suppressor module and valves can be operated only via 732 IC Detector (the keys on the 733 are locked).					

>CONFIG/733 IC Sep. Cent. >CONFIG/733/valve A	<h3>Settings for injection valve A</h3>
--	---

submenu

>CONFIG/733/valve A status	<h3>Display of the valve setting</h3> <p>Depending on the position of the valve, one of the following messages appears:</p> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="padding: 2px;">fill</td> <td style="padding: 2px;">Filling of sample loop.</td> </tr> <tr> <td style="padding: 2px;">inject</td> <td style="padding: 2px;">Injection of loop contents.</td> </tr> <tr> <td style="padding: 2px;">undefined</td> <td style="padding: 2px;">Undefined position.</td> </tr> </table>	fill	Filling of sample loop.	inject	Injection of loop contents.	undefined	Undefined position.
fill	Filling of sample loop.						
inject	Injection of loop contents.						
undefined	Undefined position.						

display only

>CONFIG/733/valve A trigger:	<h3>Switching of the valve</h3> <p>After confirmation with <ENTER>, the is moved to the selected position:</p> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="padding: 2px;">fill</td> <td style="padding: 2px;">Filling of sample loop.</td> </tr> <tr> <td style="padding: 2px;">inject</td> <td style="padding: 2px;">Injection of loop contents.</td> </tr> </table>	fill	Filling of sample loop.	inject	Injection of loop contents.
fill	Filling of sample loop.				
inject	Injection of loop contents.				

inject, fill

>CONFIG/733 IC Sep. Cent.
>CONFIG/733/valve B

submenu

Settings for injection valve B

This submenu appears only with the 733.0020 IC Separation Center with 2 injection valves.

>CONFIG/733/valve B
status fill

display only

Display of the valve setting

Depending on the position of the valve, one of the following messages appears:

fill	Filling of sample loop.
inject	Injection of loop contents.
undefined	Undefined position.

>CONFIG/733/valve B
trigger: fill

inject, fill

Switching of the valve

After confirmation with <ENTER>, the valve is moved to the selected position:

fill	Filling of sample loop.
inject	Injection of loop contents.

>CONFIG/733 IC Sep. Cent.
>CONFIG/733/suppressor

submenu

Settings for the suppressor module

This submenu appears only with the 733.0030 IC Separation Center with suppressor module.

>CONFIG/733/suppressor
auto step: fill

fill, inject, off

Automatic triggering of "STEP"

Switching of the suppressor module to the next position can be triggered automatically as follows:

fill	Triggering of "STEP" after each "FILL A".
inject	Triggering of "STEP" after each "INJECT A".
off	No automatic triggering.

>CONFIG/733/suppressor
status in position

display only

Display of the suppressor position

Depending on the position of the suppressor module, one of the following messages appears:

in position	Suppressor module is in one of the three possible positions.
undefined	Undefined position.

>CONFIG/733/suppressor
trigger: ---

---, step


Switching of the suppressor module

After confirmation of the entry with <ENTER>, the following action results:

---	No action.
step	Switching of the suppressor module to the next position.

4.5 Measurement parameters

4.5.1 <PARAM> key



The <PARAM> key combines the most important parameter settings for the 732 IC Detector. The parameter values entered here are available for the measurement immediately following their confirmation. They can be stored and recalled together with the current existing program as a method. The key opens the following main menu:

PARAM >PARAM/detector	Settings for the conductivity detector
PARAM >PARAM/anal og output	Settings for analog output
PARAM >PARAM/pl ot	Settings for graphics plot
PARAM >PARAM/709 IC Pump	Settings for 709 IC Pump

Repeated pressing of the <PARAM> key selects the submenus in turn. The individual inquiries are accessed using the <ENTER> key and exited with the <QUIT> key. The following listings show all dialog items appearing under <PARAM>. The values shown in the displays are default values, possible entry values or ranges are shown below the display.

>PARAM
>PARAM/detector

Settings for the conductivity detector

>PARAM/detector
range: 1.00 mS/cm
100, 200, 500 µS/cm
1, 2, 5, 10 mS/cm

Measuring range

7 stages from 0...100 µS/cm to 0...10 mS/cm are available for the measuring range. Select the measuring range so that the conductivity value of the eluent used is always within the selected range.

If the measured conductivity value exceeds the upper range limit by more than 10%, "-----" appears in display **1** instead of a measured value. In such a case, set the next higher range.

The selected measuring range limits the possible settings of the full-scale range, which can be selected to be maximum 2000 times more sensitive than the measuring range.



If the measuring range is changed, the auto-zero function is automatically switched off.


```
>PARAM/detector
full scale: 1.00 mS/cm
```

0.05 $\mu\text{S}/\text{cm}$...10 mS/cm

Full-scale range (operating range)

The full-scale range (operating range) sets the desired sensitivity for the display and analog output of the measuring signal during the recording of a chromatogram. The following points should be noted:

- The possible values of the full-scale range which can be selected with the <SELECT> key depend on the preset measuring range "range". This value can be reduced in a total of 11 stages by a factor 1 (lowest sensitivity) up to maximum a factor of 2000 (highest sensitivity).
- The full-scale range is primarily used to increase the sensitivity for the recording of the chromatogram following the electronic background compensation of the eluent sensitivity (auto-zero function, see also description in *section 4.6.3*).
- For the selected full-scale range (FS), the linearity of the conductivity measurement is assured in the range $-1.5 \text{ FS} \dots +1.5 \text{ FS}$. If the measured value violates these limits, overload display **8** lights up, but the measured value is still displayed. If the measured value lies outside $\pm 1.8 \text{ FS}$, overload display **8** starts to flash, at the same time display **1** shows the absolute conductivity in place of the auto-zero value.
- The full-scale range determines the limits for the output of the analog signal at analog outputs **11** (0...1 V) and **12** (0...10 mV). The polarity of the output signal can be changed at any time with the parameter "polarity" to "+" (0...+FS) or "-" (0...-FS).



>PARAM/detector
temp. coeff. : 2.5 %/°C
1.5, 2.5 %/°C

Temperature coefficient

The temperature coefficient is used to take into account the temperature dependence of the conductivity measurement. Two values can be selected:

- 1.5 %/°C Temperature coefficient for cations
- 2.5 %/°C Temperature coefficient for anions

When the thermostating of the measuring cell is switched on (thermostat display **7** lit up), the conductivity k_T measured at the set operating temperature (25...45 °C) is automatically converted to the conductivity k_{20} at the reference temperature 20 °C.

>PARAM
>PARAM/anal og output

Settings for the analog output

The output signal at the two analog outputs **11** (0...1 V) and **12** (0...10 mV) is determined by the parameter "full scale" (see above).

>PARAM/anal og output
pol arity: +
+,-

Polarity of the output signal

Selection of whether a positive or negative full-scale range should be outputted

- + Positive polarity (0...+FS)
- Negative polarity (0...-FS)

>PARAM/anal og output
off set: 0 %fs
0, 10, 50 %fs

Zero point offset

Offset of the zero point of the conductivity.

- 0 %fs No zero point offset.
- 10 %fs Offset by 10% of the full-scale range; at the analog output the signal -0.1 FS...+0.9 FS (+) or 0.1 FS...-0.9 FS (-) is outputted.
- 50 %fs Offset by 50% of the full-scale range; at the analog output the signal -0.5 FS...+0.5 FS (+) or 0.5 FS...-0.5 FS (-) is outputted.

>PARAM/anal og output
damp ing: off
on, off

Damping

- on Damping to reduce disturbing pulsations or a high noise level. Standard solutions and samples must both be recorded with damping as this can influence the peak heights.
- off Damping switched off (normal condition with pulsation-free and constant flow pumps).

```
>PARAM
>PARAM/pl ot
```

Settings for graphics plot

These settings are effective with the graphics plot on an external printer triggered using the <PLOT> key (see section 4.8.2).

```
>PARAM/pl ot
auto start:
```

```
off
on, off
```

Automatic start

- on Automatic start of the graphics plot on every "INJECT A".
- off No automatic start of the graphics plot.

```
>PARAM/pl ot
time interval
```

```
1.0 s
0.4...99999 s
```

Time interval

Time interval after which a new measured value is plotted.

```
>PARAM/pl ot
time scale
```

```
10.0 mm/min
0.1, 0.2, 0.5, 1, 2, 5,
10, 20, 60, 120 mm/min
```

Time scale (paper speed)

The entry value corresponds to the paper feed of the printer in mm/min (the effective paper feed depends on the printer used).

If values other than those specified are entered, these will be automatically rounded to the next preset value.

```
>PARAM/pl ot
time scale label:
```

```
rel
rel, abs
```

Labeling of the time axis

The time axis is labeled at every main division.

- abs Labeling with absolute (current) time (e.g. "08:34").
- rel Labeling with relative time beginning with the start of the curve printout (e.g. "2m40s" for 2 min 40 s).

```
>PARAM/pl ot
stop time
```

```
off min
1...999 min, off
```

Stop time for graphics plot

Time until graphics plot is stopped.

- off No time limit (press <CLEAR>).

```
>PARAM/pl ot
left:
```

```
0.000 mS/cm
-2000...2000 µS/cm
-2000...2000 mS/cm
-2...2 S/cm
```

Left boundary of conductivity axis

Left limit value for graphics plot of the conductivity on an external printer.

The unit ($\mu\text{S/cm}$, mS/cm , S/cm) can be selected with the <SELECT> key.

```
>PARAM/pl ot
right:
```

```
10.00 mS/cm
-2000...2000 µS/cm
-2000...2000 mS/cm
-2...2 S/cm
```

Right boundary of the conductivity axis

Right limit value for graphics plot of the conductivity on an external printer.

The unit ($\mu\text{S/cm}$, mS/cm , S/cm) can be selected with the <SELECT> key.

>PARAM >PARAM/709 IC Pump	<h3>Settings for 709 IC Pump</h3> <p><i>This submenu appears only with "SETUP/peripherals/operation with 709 = on".</i></p> <p>These settings are effective only if the 709 IC Pump is connected to the 732 IC Detector and the remote control is activated with the <EXT.> key (see <i>section 2.6.1</i>).</p>
------------------------------	---

>PARAM/709 IC Pump flow 0.50 mL/min 0.01...5.00 mL/min	<h3>Flow rate</h3> <p>Flow rate of the 709 IC Pump.</p>
---	---

>PARAM/709 IC Pump Pmax 10.0 MPa 0.1...50 MPa	<h3>Maximum shutoff pressure</h3> <p>Maximum shutoff pressure for the automatic shutoff of the pump drive of the 709 IC Pump to protect the separating column against excessive pressure.</p> <p>This limit value should be ca. 3 MPa (= 30 bar) above the normal operating pressure of the separating column used or correspond to the maximum admissible operating pressure of the separating column.</p> <p>If the 709 Pump exceeds the specified limit value, the pump drive is switched off within a pump cycle and blocked, at the same time the bottom line of the display of the 732 IC Detector shows the message "stopped: YY-MM-DD HH:MM" with details of the time when the pump was stopped. The LED "P_{max}" on the 709 IC Pump flashes, the shutoff pressure can be requested under "P_{actual}".</p> <p>A 709 IC Pump stopped in this manner is put back into operation as follows:</p> <ul style="list-style-type: none"> • Press <EXT.> key to switch off remote operation. • Press <R/S> key to remove block. The LED "P_{max}" no longer flashes. • Press <EXT.> key to switch remote operation back on.
--	---

>PARAM/709 IC Pump Pmin 0.0 MPa 0.0...50 MPa	<h3>Minimum shutoff pressure</h3> <p>Minimum shutoff pressure for the automatic shutoff of the pump drive of the 709 IC Pump to protect the IC system when leaks occur or the eluent supply is interrupted.</p> <p>This limit value should be sufficiently far below the normal operating pressure of the separating column used. The entry of "0.0" switches the automatic shutoff off.</p> <p>If the pressure value of the 709 IC Pump is below the specified limit value for a considerable length of time, the pump drive of the pump is switched off and blocked, at the same time the bottom line of the display</p>
---	--

of the 732 IC Detector shows the message "stopped: YY-MM-DD HH:MM" with details of the time when the pump was stopped. The LED "P_{min}" on the 709 IC Pump flashes, the shutoff pressure can be requested under "P_{actual}".

A 709 IC Pump stopped in this manner is put back into operation as follows:

- Press <EXT.> key to switch off remote operation.
- Press <R/S> key to remove the block. The LED "P_{min}" no longer flashes.
- Press <EXT.> key to switch the remote operation back on.

>PARAM/709 IC Pump flow corr.	1.00
----------------------------------	------

0.90...1.10

Correction factor of the flow rate

Here a correction factor for the flow rate can be entered to match the displayed flow rate with the flow rate actually measured. You will find further information in *section 3.3* of the 709 Instructions for Use.

4.5.2 <FULL SCALE>key



The <FULL SCALE> key can be used for direct entry of the full-scale range to avoid the circuitous route via parameter menu (see *section 4.5.1*). You will find a detailed description of the parameter "full scale" in *section 4.5.1*.

XXXXX mS/cm	XXXX min
full scale:	1.00 mS/cm

0.05 µS/cm...10 mS/cm

Full-scale range (operating range)



After the <FULL SCALE> key has been pressed, the bottom line of the display shows the parameter "full scale" with a flashing numeric value. In contrast to the entry under the <PARAM> key, here the parameter value is not entered using the <SELECT> key, but selected with the <←> and <→> keys. Proceed as follows:

- Press the <←> key to lower the full-scale range stepwise to the lowest possible value.
- Press the <→> key to increase the full-scale range stepwise to the highest possible value.
- Confirm the desired full-scale range with the <ENTER> or <FULL SCALE> key.

4.6 Triggering of functions


4.6.1 <FILL> keys

The two <FILL> keys at the 733 IC Separation Center each have the following function, which depends on the instrument version:

	733.0010 (1 valve)	733.0X20 (2 valves)	733.0X30 (valve+ suppr.)
 <p>A</p>	Switching of injection valve A to "FILL" position	Switching of injection valve A to "FILL" position	Switching of injection valve A to "FILL" position
 <p>B</p>	No function	Switching of injection valve B to "FILL" position	Switching of suppressor module to next position ("STEP" function)

If the injection valve is in the "FILL" position, the green LED in the key lights up. In this position the injection valve is connected as shown in the diagram opposite:

- The sample loop is filled by siphoning in the sample solution using the syringe.
- The eluent flows directly to the separating column.





The "FILL" function can be triggered at any time with the <FILL> key – also in the edit mode or when a program is running. If this is not required, the key can be locked. For this, the parameter ">CONFIG/733 IC Sep.Cent./contro1" must be set to "732 on1y" under the <CONFIG> key (see section 4.4.2). Switching of the injection valve is then possible only using the 732 IC Detector via a programmed "FILL" command or with ">CONFIG/733 IC Sep.Cent./trigger = fill".

With the 733.0030 IC Separation Center, switching of the suppressor module to the next position can be triggered automatically by the switching of injection valve A to the "FILL" position. For this, the parameter ">CONFIG/733/suppressor/auto step" must be set to "fill" under the <CONFIG> key (see section 4.4.2).


4.6.2 <INJECT> keys

The two <INJECT> keys at the 733 IC Separation Center have the following function, depending on the instrument version:

	733.0010 (1 valve)	733.0X20 (2 valves)	733.0X30 (valve+ suppr.)
 A	Switching of injection valve A to "INJECT" position	Switching of injection valve A to "INJECT" position	Switching of injection valve A to "INJECT" position
 B	No function	Switching of injection valve B to "INJECT" position	No function

If the injection valve is in the "INJECT" position, the green LED in the key lights up. In this position, the injection valve is connected as shown in the diagram opposite:


- The eluent flows via the sample loop to the separating column. This causes injection of the sample previously added to the sample loop.
- Sample inlet tubing and syringe are directly connected.




The "INJECT" function can be triggered at any time with the <INJECT> key – also in the edit mode or when a program is running. If this is not required, the key can be locked. For this, the parameter ">CONFIG/733 IC Sep.Cent./control" must be set to "732 only" under the <CONFIG> key (see section 4.4.2). Switching of the injection valve is then possible only using the 732 IC Detector via a programmed "INJECT" command or with ">CONFIG/733 IC Sep.Cent./trigger = inject". Switching of injection valve A to the "INJECT" position can also be used for the automatic start of a program of the "inject" type.

With the 733.0X30 Separation Center, switching of the suppressor module to the next position can be triggered automatically by the switching of injection valve A to the "INJECT" position. For this, the parameter ">CONFIG/733/suppressor/auto step" must be set to "inject" under the <CONFIG> key (see section 4.4.2).

4.6.3 <ZERO> key



The <ZERO> key is used to trigger the auto-zero function. "Auto-zero" is the name used for the automatic electronic background compensation, i.e. the current measured conductivity value is set to zero and hence lies in the middle of the selected full-scale range. The diagram below shows the operating principle of the auto-zero function. Before the auto-zero function is activated for the first time or after it has been switched off, the measuring signal lies outside the full-scale range (FS); it is automatically set to 0 $\mu\text{S}/\text{cm}$ each time the <ZERO> key is pressed.



The auto-zero function can be triggered at any time using the <ZERO> key – also in the edit mode or when a program is running. Lighting up of green LED **6** in the <ZERO> key shows that the electronic background compensation is active.



The auto-zero function works properly only when the measured value remains relatively stable (e.g. not during the appearance of the injection peak at the start of the chromatogram).

Lighting up of the red overload display **8** shows that the measured value lies outside $\pm 150\%$ of the selected full-scale range. If the measured value is outside $\pm 180\%$ of the selected full-scale range, the overload display **8** starts to flash and the display shows the absolute conductivity (see also section 4.5.1).

4.6.4 <ZERO OFF> key



The <ZERO OFF> key is used to switch off the auto-zero function. The auto-zero function can be triggered using the <ZERO OFF> key in the standby mode or when a program is running, but not in the edit mode.

When the green LED **6** in the <ZERO> key goes out, this shows that the electronic background compensation is no longer active.


4.6.5 <MARK> key



The <MARK> key is used to trigger a marking signal. In addition to the measuring signal, a signal of ca. 10% of the full-scale range is outputted briefly at the analog output. This allows, e.g. the start of a chromatogram to be marked during recording on a recorder.

Marking can be triggered with the <MARK> key in the standby mode or when a program is running, but not in the edit mode.

4.6.6 <PUMP R/S> key



The <PUMP R/S> key is used for on/off switching of the pump drive of the 709 IC Pump. A requirement for this is that the 709 IC Pump is connected to the 732 IC Detector and the remote control is switched on at the pump with the <EXT.> key (see *section 2.6.1*).

With the <PUMP R/S> key the pump drive of the 709 IC Pump can be switched on or off in the standby mode or when a program is running, but not in the edit mode. When the pump drive is switched on, the red LED above the <R/S> key on the pump lights up.

The settings of the pump parameters on the 732 IC Detector 732 are described under the <CONFIG> key in *section 4.4.2*, the pump status messages in the standby mode in *section 4.3.2*.

4.7 Programming

4.7.1 <PROGRAM> key



The <PROGRAM> key can be used to develop and edit time programs which can comprise maximum 20 program steps each with up to 14 individual program points. This allows the chronological sequence of actions at the 732 IC Detector and 733 IC Separation Center to be automated virtually at will.

The program loaded in the working memory can be stored together with the parameter values currently under the <PARAM> key as a method and recalled.

The <PROGRAM> key opens the following main menu:

PROGRAM >PROGRAM/parameters	General program parameters
PROGRAM >PROGRAM/edit	Edit program
PROGRAM >PROGRAM/delete all	Delete all program steps

Repeated pressing of the <PROGRAM> key selects the submenus in turn. The individual inquiries of a submenu are accessed using the <ENTER> key and exited with the <QUIT> key. The following listings show all dialog items appearing under <PROGRAM>. The values shown in the display are default values, the possible entry values or ranges are shown below the display.

```
>PROGRAM
>PROGRAM/parameters
```

General program parameters

```
>PROGRAM/parameters
type:                cycle
```


cycle, remote, inject

Program type

The program type is used to define the type of program start and program sequence.

cycle **Loop program with start via <PROG R/S>**

A loop program can be started by using the <PROG R/S> key, by external start via a remote input line (see *section 6.2.1*), or via EVENT (see *section 4.7.3*). It is automatically run several times if it includes a "return". The number of repetitions is defined by the parameter "number of cycles". The basic structure of a loop program has the following appearance:



remote Single program with start via remote
 A remote program is started externally via a remote input line (see section 6.2.1) and executed once only.

inject Single program with start via "INJECT"
 An inject program is started each time injection valve A is switched to the "INJECT" position and executed once only.

```

>PROGRAM/parameters
number of cycles      1
                      1...999
  
```

Number of cycles

This inquiry appears only with "PROGRAM/parameters/type = cycle".

Number of executions of loop programs.



This parameter is identical to that of the same name in the submenu ">CONFIG/auxiliaries" of the <CONFIG> key (see section 4.4.2). As access to the <PROGRAM> key is blocked while a program is running, the number of cycles can be changed only there.

```

>PROGRAM/parameters
status:              di sabled
                    enabled, disabled, test
  
```

Program status

The three following program statuses are possible:

enabled Program enabled
 Irrespective of the program type, a program can be started only when it is enabled. This status is shown by the lighting up of the green LED **5** in the <PROG R/S> key. If a method with program steps is loaded, the program status is automatically set to "enabled".

disabled Program disabled
 A disabled program can not be started. In this status, the green LED **5** in the <PROG R/S> key remains off. When the instrument is switched on, the program status is always set to "disabled".

test

Program test

The test status is used to test a program. After the program status has been set to "test" and return to the standby mode, a switch is automatically made to the test mode, which runs as follows:

- The bottom line of display shows the message "test step" with details of the time "x.xxx min" when the first program step should be executed.
- Pressing the <ENTER> key executes the actions associated with the program step.
- The display shows the next program steps in turn which are triggered with <ENTER>.
- After execution of the last program step, the message "program end" appears and must be confirmed with <ENTER>. The instrument then returns to the standby mode and the program status is set to "disabled".

```
>PROGRAM
>PROGRAM/edi t
```

Edit program

This submenu is used for the entry of new program steps or to edit existing program steps and their program points.

```
>PROGRAM/edi t
time:          x. x mi n
              0.0...999.9 min
```

Edit program step

Under this inquiry new program steps can be entered or existing ones selected, changed and deleted. Each program step includes a definite time (two program steps with the same time can not be entered). A program can contain up to 20 program steps.

Enter new program step

A new program step can be entered when "x.x min" for the time appears in the display. The new time entered with the numeric keys must be confirmed with the <ENTER> key, the individual program points which should be executed at this time can then be entered.

Select existing program steps

Program steps already entered can be selected with the <SELECT> key. After the <ENTER> key has been pressed, the individual program points of the selected

program step can be edited.

Change time for existing program steps

For a program step selected with <SELECT>, a new time can be entered directly with the numeric keys and must be confirmed with <ENTER>. The individual program points of the program step can then be edited (if no change is necessary, press <QUIT>).

Delete existing program steps

A program step selected with <SELECT> can be deleted by pressing the two keys <CLEAR> and <ENTER> simultaneously. If you wish to delete all program steps at once, it is better to use the submenu "PROGRAM/delete all".

```
>PROGRAM/edi t. XXX. X mi n
XXXXXXXXXX      ---
```

Edit program points

After confirmation of the newly entered or selected time with <ENTER>, it is possible to edit the maximum 14 possible program points of a program step. The time of the program step then appears in the top line of the display, whereas in the bottom line the program points listed below associated with the program step can be edited in turn. For all program points, the default value " --- " means that no action results.

```
>PROGRAM/edi t. XXX. X mi n
fl ag:          ---
---, return, reset, end
```

Set program flag

After setting of a program flag, no additional program points can be entered.

return Return with loop programs

Each time the "return" flag is reached, a return is made to the start of the program until the inputted number of cycles has been reached.

reset Reset the parameters

With this flag all parameters of the <PARAM> key are reset to the values valid before the start of the program.

end Program end

This flag is used to end a program, additional program steps are no longer executed. Programs without "end" are not stopped until <PROG R/S> is pressed.

```
>PROGRAM/edi t. XXX. X mi n
FS change:     ---
---, on
```

Change full-scale range

on The full-scale range should be changed. As the possible settings of the full-scale range depend on the measuring range, both the measuring range and the full-scale range must be reentered.

>PROGRAM/edit. XXX. X min
range: 100 mS/cm

100, 200, 500 μ S/cm
1, 2, 5, 10 mS/cm

Reset measuring range

This inquiry appears only with "FS change = on".

The measuring range is set to the inputted value (for further details of the parameter "range", see section 4.5.1).



If the measuring range just entered differs from the measuring range valid at the start of the program, an auto-zero function that is switched on will be switched off. If the auto-zero function should be switched on again, this can not be done in the same program step.

>PROGRAM/edit. XXX. X min
full scale: 100 mS/cm

0.05 μ S/cm...10 mS/cm

Reset full-scale range

This inquiry appears only with "FS change = on".

The full-scale range is set to the inputted range (for further details of the parameter "full scale", see section 4.5.1).

>PROGRAM/edit. XXX. X min
zero: ---

---, on, off

Switch auto-zero on/off

This inquiry does not appear with "FS change = on".

- on Trigger the auto-zero function (equivalent to the <ZERO> key, see section 4.6.3).
- off Switch off the auto-zero function (equivalent to the <ZERO OFF> key, see section 4.6.4).

>PROGRAM/edit. XXX. X min
polarity: ---

---, +, -

Reset the polarity of the output signal

- + Output of the positive full-scale range (0...+FS; see section 4.5.1).
- Output of the negative full-scale range (0...-FS; see section 4.5.1).

>PROGRAM/edit. XXX. X min
mark: ---

---, on

Trigger marking signal

- on Output of a marking signal at the analog output (equivalent to the <MARK> key, see section 4.6.5).

>PROGRAM/edit. XXX. X min
valve A: ---

---, fill, inject

Switch injection valve A

This inquiry appears only with "SETUP/peripherals/operation with 733 = on".

- fill Switching of injection valve A to "FILL" position (equivalent to the <FILL> key, see section 4.6.1).
- inject Switching of injection valve A to "INJECT" position (equivalent to the <INJECT> key, see section 4.6.2).

```
>PROGRAM/edit. XXX. X min
valve B:      ---
```

```
---, fill, inject
```

Switch injection valve B

This inquiry appears only with "SETUP/peripherals/operation with 733 = on" and with the 733.0020 IC Separation Center with two injection valves.

- fill** Switching of injection valve B to "FILL" position (equivalent to the <FILL> key, see section 4.6.1).
- inject** Switching of injection valve B to "INJECT" position (equivalent to the <INJECT> key, see section 4.6.2).

```
>PROGRAM/edit. XXX. X min
suppressor:  ---
```

```
---, step
```

Switch suppressor module

This inquiry appears only with "SETUP/peripherals/operation with 733 = on" and with the 733.0X30 IC Separation Center with suppressor module.

- step** Switching of the suppressor module to the next position (equivalent to the key function <STEP>, see section 4.6.1).

```
>PROGRAM/edit. XXX. X min
remote       *****
```

```
*, 0, 1
```

```
<CLEAR> sets all values to *
```

Set remote output lines 1...8

Sets the remote output lines 1...8 (from left to right) to the following values:

- 0** Line off, inactive (open).
- 1** Line on, active (0 V).
- *** Leave line in current status (entry of "*" with numeric key <. >).

You will find further details of the remote output lines in section 6.2.

```
>PROGRAM/edit. XXX. X min
pump R/S:    ---
```

```
---, on, off
```

Start/stop pump drive

This inquiry appears only with "SETUP/peripherals/operation with 709 = on".

- on** Switch on pump drive of the 709 IC Pump (equivalent to the <PUMP R/S> key, see section 4.6.6).
- off** Switch off 709 IC Pump Drive (equivalent to the <PUMP R/S> key, see section 4.6.6).

```
>PROGRAM/edit. XXX. X min
flow:        --- mL/min
```

```
---, 0.01...5 mL/min
```

Reset flow rate

This inquiry appears only with "SETUP/peripherals/operation with 709 = on".

The flow rate of the 709 IC Pump is set to the inputted value (for further details of the "flow" parameter, see section 4.5.1).

```
>PROGRAM/edit. XXX.X min
Pmax:      --- MPa
          ---, 0.1...50 MPa
```

Reset maximum shutoff pressure

This inquiry appears only with "SETUP/peripherals/operation with 709 = on".

The maximum shutoff pressure for automatic shutoff of the 709 IC Pump is set to the inputted value (for further details of the parameter "Pmax", see *section 4.5.1*).

```
>PROGRAM/edit. XXX.X min
Pmin:      --- MPa
          ---, 0.0...50 MPa
```

Reset minimum shutoff pressure

This inquiry appears only with "SETUP/peripherals/operation with 709 = on".

The minimum shutoff pressure for the automatic shutoff of the 709 IC Pump is set to the inputted value (for further details of the parameter "Pmin", see *section 4.5.1*).

```
>PROGRAM
>PROGRAM/delete all
```

Delete all program steps

```
>PROGRAM/delete all
delete all?
<ENTER>, <QUIT>
```

Confirm deletion of all program steps

- <ENTER> All program steps of the program will be deleted.
- <QUIT> Abort of the deletion operation.

4.7.2 <PROG R/S> key



The <PROG R/S> key is used to start or stop the program loaded in the working memory. The status of the LED in the key indicates what it can trigger:

- LED off** **Program disabled**
No action is triggered with the <PROG R/S> key. A loop program can not be started until the parameter "PROGRAM/parameters/status = enabled" has been set.
- LED lit up** **Program enabled**
The <PROG R/S> key is used to start a loop program. The current time is automatically set to "0.0 min".
- LED flashes** **Program running**
The <PROG R/S> key is used to stop an ongoing loop, remote or inject program.



If a program is running, the <PARAM>, <FULL SCALE>, <METHOD> and <PROGRAM> keys are locked.

4.7.3 <EVENT> key



Under the <EVENT> key maximum 4 different events, which can be triggered (once or daily) at a particular time, can be entered. These events can each comprise up to 14 individual program points.

The <EVENT> key opens the following main menu:

EVENT >EVENT/new	Enter new event
EVENT >EVENT/edi t	Edit existing event
EVENT >EVENT/del ete	Delete individual event
EVENT >EVENT/del ete all	Delete all events

Repeated pressing of the <EVENT> key selects the submenus in turn. The individual inquiries of the submenus are accessed using the <ENTER> key and exited with the <QUIT> key. The following listings show all dialog items which appear under <EVENT>. The values shown in the displays are the default values, the possible entry values and ranges are shown below the display.

```
>EVENT
>EVENT/new
```

Enter new event

```
>EVENT/new
format:          date
                date, daily
```

Event format

date The event is triggered once at a particular time.

daily The event is triggered daily at a particular time.

```
>EVENT/new
event: YY- MM- DD  HH: MM: SS
      YY: 00...99  HH: 00...23
      MM: 01...12  MM: 00...59
      DD: 01...31  SS: 00...59
```

Date and time for a new event

This inquiry appears only with "format = date".

Entry of date and time when the event should be triggered (YY: 00...94 = 2000...2094, 95...99 = 1995...1999).

```
>EVENT/new
event:          HH: MM: SS
              HH: 00...23
              MM: 00...59
              SS: 00...59
```

Time for new event

This inquiry appears only with "format = daily".

Entry of the time when the event should be triggered daily.

```
EVENT YY- MM- DD  HH: MM: SS
XXXXXXXXXX      ---
```

Edit program points

After confirmation of the inputted date or time with <ENTER>, it is possible to edit the maximum 14 possible program points of a new event. Date and time of the event then appear in the top line of the display, whereas in the bottom line the program points associated with the event can be edited in turn. For all program points, the default value " - - - " means that no action follows.

Except for the missing parameter "flag", the program points for events are identical with the program points which can be entered under the <PROGRAM> key and are thus not described further here (see *section 4.7.1*).

In addition to this points, the following program point can be entered:

```
EVENT YY- MM- DD  HH: MM: SS
prog R/S:        ---
                ---, on, off
```

Start/stop program

on Start of the program loaded in the working memory (corresponds to the <PROG R/S> key, see *section 4.7.2*).

off Stop of the program loaded in the working memory (corresponds to the <PROG R/S> key, see *section 4.7.2*).

```
>EVENT
>EVENT/edi t
```

Edit existing event

```
>EVENT/edi t
event: YY-MM-DD HH:MM:SS
```

Select existing event

Selection of the event with the <SELECT> key.



Date and time of an existing event can not be changed later.

```
EVENT YY-MM-DD HH:MM:SS
XXXXXXXXX ---
```

Edit program points

After confirmation of the selected date or time with <ENTER>, it is possible to edit the maximum 14 possible program points of an existing event. Date and time of the event then appear in the top line of the display, whereas in the bottom line the program points associated with the event can be edited in turn. For all program points, the default value " - - - " means that no action follows.

Except for the missing parameter "flag", the program points for events are identical with the program points which can be entered under the <PROGRAM> key and are thus not described further here (see *section 4.7.1*).

In addition to this points, the program point "prog R/S" can be entered (see under ">EVENT/new").

```
>EVENT
>EVENT/delete
```

Delete individual event

```
>EVENT/delete
event: YY-MM-DD HH:MM:SS
```

Select event for deletion

Selection of the event with the <SELECT> key.

<ENTER>, <QUIT>

<ENTER> The event will be deleted

<QUIT> Abort of the deletion operation.

```
>EVENT
>EVENT/delete all
```

Delete all events

```
>EVENT/delete all
delete all?
```

Confirm deletion of all events

<ENTER> All events will be deleted.

<QUIT> Abort of the deletion operation.

<ENTER>, <QUIT>

4.7.4 <METHOD> key



The <METHOD> key is used for the management of maximum 9 user-defined methods, which can be stored and recalled under freely selectable names. A method comprises the measurement parameters defined under the <PARAM> key and the program defined under the <PROGRAM> key.

The <METHOD> key opens the following main menu:

METHOD >METHOD/recal l	Recall method
METHOD >METHOD/store	Store method
METHOD >METHOD/del et e	Delete method

Repeated pressing of the <EVENT> key selects the submenus in turn. The individual inquiries of a submenu are accessed using the <ENTER> key and exited with the <QUIT> key. The following listings show all dialog items which appear under <METHOD>. The values shown in the displays are the default values, the possible entry values and ranges are shown below the display.

```
>METHOD
>METHOD/recal l
```

Recall existing method

```
>METHOD/recal l
name:                DEFAULT
                    DEFAULT, Prep-MSM,
                    XXXXXXXXX
```

Selection of the method

Use of the <SELECT> key to select the existing method which should be recalled.

- DEFAULT Default method with empty program and parameters set to the default values.
- Prep-MSM Method for conditioning the suppressor. The suppressor module is switched to the next position every 20 min.
- xxxxxxx Name of the user method stored.

After confirmation of the name with the <ENTER> key, the method is loaded in the working memory. The parameters under the <PARAM> and <PROGRAM> keys will be overwritten.

```
>METHOD
>METHOD/store
```

Store method

```
>METHOD/store
name:          XXXXXXXX
```

8 ASCII characters

Enter or select method name

Entry of a new name (text entry, see *section 4.2.5*) under which the method should be stored or selection of an existing method name with the <SELECT> key.

After confirmation of the name with the <ENTER> key, the parameters loaded in the working memory under the <PARAM> and <PROGRAM> keys are stored as a method.

```
>METHOD/store
overwrite XXXXXXXX?
```

<ENTER>, <QUIT>

Confirmation of overwriting

This inquiry appears only if a method should be stored under a name which already exists.

<ENTER> The method is stored under the selected name.

<QUIT> Abort of the storage operation.

```
>METHOD
>METHOD/del ete
```

Delete method

```
>METHOD/del ete
name:          XXXXXXXX
```

Selection of the method

Selection of the method which should be deleted using the <SELECT> key.

```
>METHOD/del ete
del ete XXXXXXXX?
```

<ENTER>, <QUIT>

Confirmation of the deletion

<ENTER> The selected method will be deleted.

<QUIT> Abort of the deletion operation.

4.8 Data output

4.8.1 <PRINT> key



The <PRINT> key is used for the output of single measured values or to start the automatic measured value output on an external printer or PC connected to the RS232 interface. The output is triggered directly at a keystroke, no inquiries appear.

Preparations

- For the output on an **external printer**, printer type and RS232 settings must be correctly selected under <CONFIG> (see *section 4.4.2*).
- For the output on a **PC** the character set under <CONFIG> ("send to: IBM") and the RS232 settings must be correctly selected (see *section 4.4.2*). For receipt of the data, a suitable program must be started on the PC (e.g. Metrodata VESUV 2.1 or terminal program in Windows accessories).

Elements of the measured value output

Depending on the parameter settings under ">CONFIG/print meas. value", the measured value output can comprise the following elements:

- **Print header**
The print header comprises a title line (device name, serial number, program number), identifications 1 and 2 (if used) as well as the date and time if desired (see *section 4.4.2*) and is outputted once, always or never before a measured value output, depending on the setting of the parameter ">CONFIG/printer/print header".
- **Run number**
The run number (parameter ">CONFIG/auxiliaries/run number", see *section 4.4.2*) is incremented by +1 and outputted on every "INJECT A" (without program) or on every new start or "return" (with program).
- **Run time**
The current time after the program start, after the last "INJECT A" or after the instrument has been switched on is outputted before the measured value.
- **Measured value**
The measured value shown in the display is outputted:

Measured value without sign	Absolute conductivity in $\mu\text{S}/\text{cm}$ or mS/cm
Measured value with sign	Auto-zero value in $\mu\text{S}/\text{cm}$, %fs or mV
" - - - - - "	Measured value is >10% outside the measuring range

- **Identification "ZEROREF"**
Reference value for the auto-zero function (measured value which is compensated).
- **Identification "OVERLOAD"**
The measured value is 50...80% outside the selected full-scale range.
- **Identification "BACKGROUND"**
The measured value is more than 80% outside the selected full-scale range. The absolute conductivity (background) is outputted in place of the auto-zero value.

Output of individual measured values

- For the output of individual measured values, ">CONFIG/print meas.value = immed." must be selected (see *section 4.4.2*). Each time the <PRINT> key is pressed, a measured value will be outputted via the RS232 interface.
- Example:
Individual measured value output with print header output (with date and time), run time, measured value

```

>CONFIG/printer
  id.1
  id.2
  print header:      always
  date&time:         on
  send to:           Seiko
>CONFIG/print meas.value
  print crit.:       immed.
  date&time:         off

```

<pre> 732 IC Detector 01104 732.0012 date 1995-10-09 time 13:24:58 12.3 138.9 uS/cm ===== </pre>	Print header Date and time Run time, measured value
<pre> 732 IC Detector 01104 732.0012 date 1995-10-09 time 13:25:37 139.1 uS/cm ZEROREF 13.0 +0.0 uS/cm ===== </pre>	Print header Date and time Auto-zero reference value Run time, auto-zero value

Automatic measured value output

- For the automatic measured value output of measured values at selectable time intervals, ">CONFIG/print meas.value/print crit. = time" must be selected (see *section 4.4.2*). When the <PRINT> key is pressed, the measured value output defined by the parameters "time interval" and "stop time" is started via the RS232 interval. The measured value output can always be aborted with <QUIT>.
- Example:
Automatic measured value output with single print header output (with date and time), run time, measured value

```
>CONFIG/printer
id.1 column PRP-X100
id.2 conditioning
print header:      once
date&time:         on
send to:           Seiko
>CONFIG/print meas.value
print crit.:       time
time interval      600.0 s
stop time          60 min
date&time:         off
```

```
732 IC Detector      01104   732.0012
id.1 column PRP-X100
id.2 conditioning
date 1995-10-09   time 14:36:43

0.0  109.6 uS/cm
10.0 125.1 uS/cm
20.0 133.9 uS/cm
30.0 138.2 uS/cm
40.0 140.6 uS/cm
50.0 141.0 uS/cm
60.0 141.2 uS/cm
```

Print header
 Identification 1
 Identification 2
 Date and time

 Run time, 1st meas. value
 Run time, 2nd meas. value
 Run time, 3rd meas. value
 Run time, 4th meas. value
 Run time, 5th meas. value
 Run time, 6th meas. value
 Run time, 7th meas. value

4.8.2 <PLOT> key



The <PLOT> key is used to start the graphics plot on an external printer connected to the RS232 interface. The output is started directly at a keystroke, no inquiries appear.

Preparations

- Under "SETUP" the general graphics parameters (frame, grid, print width) must be set to the desired values (see *section 4.4.1*).
- Under <CONFIG> printer type and RS232 settings for the output on an external printer must be correctly selected (see *section 4.4.2*).
- Under <PARAM> the parameters effective for the graphics plot must be set (see *section 4.5.1*).

Elements of the measured value output

Depending on the parameter settings selected under "SETUP", <CONFIG> and <PARAM>, the graphics plot can comprise the following elements:

- **Print header**
The print header comprises title line (device name, serial number, program number), identifications 1 and 2 (if used) as well as the date and time if required (see *section 4.4.2*) and is outputted once, always or never before a graphics plot, depending on the setting of the parameter ">CONFIG/printer/print header".
- **Grid**
With "SETUP/grid = on" grid lines are drawn in.
- **Frame**
With "SETUP/frame = on" a frame is drawn in.
- **Labeling of the y axis**
Conductivity in $\mu\text{S}/\text{cm}$, mS/cm or S/cm
- **Labeling of the x axis**
Absolute time "HH:MM:SS" or relative time "XXXmXXs"
- **Conductivity measurement value**
Solid line


Example

- Graphics plot of a chromatogram with the PRP-X100 IC Anion column on a Citizen printer IDP562 RS

```
>SETUP/graphics
grid:                off
frame:               off
width                1.0
```

```
>CONFIG/printer
id.1 column PRP-X100
id.2 standard
print header:        always
date&time:           on
send to:             Citizen
```

```
>PARAM/printer
auto start:          on
time interval        1.0 s
time scale           10.0 mm/min
time scale label:    rel
stop time            10 min
left:                -0.200 µS/cm
right:               1.800 µS/cm
```



4.8.3 <REPORT> key



The <REPORT> key is used for the output of reports on an external printer or PC connected to the RS232 interface.

The <REPORT> key opens the following inquiry:

XXXXXX mS/cm report:	XXXX mi n all	Selection of the report
all, config, event, method, method list, param, program, pump		Selection of the report which should be outputted:
	all	All reports in the order "config, event, method, method list, pump"
	config	Configuration report (parameters of the <CONFIG> key)
	event	Event report (parameters of the <EVENT> key)
	method	Method report (comprises method name, parameter and program report)
	method list	Method list (details of stored methods)
	param	Parameter report (parameters of the <PARAM> key)
	program	Program report (parameters of the <PROGRAM> key)
	pump	Report of the pump parameters

Example of "config" report

Configuration for 732 IC Detector with attached Seiko printer and 733.0010 IC Separation Center with 1 injection valve

<pre> CONFIG >CONFIG/detector thermostat: 35 °C 'zero' unit: µS/cm cell constant 16.7 /cm >CONFIG/printer id.1 PRP-X100 id.2 drinking water print header: once date&time: on send to: Seiko >CONFIG/print meas.value print crit.: time time interval 300.0 s stop time 60 min date&time: off >CONFIG/auxiliaries run number 0 number of cycles 1 >CONFIG/aux/event enabled in program: off date 1995-10-10 time 11:38:20 dialog: english device label program 732.0012 </pre>	<pre> >CONFIG/aux/beeper status: on repeat time 60 s >CONFIG/RS settings baud rate: 9600 data bit: 8 stop bit: 1 parity: none handshake: HWS control via RS: on >CONFIG/RS settings 709 baud rate: 9600 >CONFIG/733 IC Sep.Cent. valve A control: no restriction >CONFIG/733/valve A status inject trigger: inject ----- </pre>
---	---

Example of "event" report

Two daily events (adjust pump flow) and a date event (switch off pump supply)

```
EVENT
daily:          07:00:00
               pump flow      2.00 mL/min
daily:          19:00:00
               pump flow      0.30 mL/min
date:   95-10-13 16:00:00
               pump           off
               -----
```

Example of "method" report

Method for the determination of drinking water with PRP-X100 Anion Column

```
METHOD
method name      PRP-X100
date: 1995-10-10 15:15:41

PARAM
>PARAM/detector
  range:          200 µS/cm
  full scale:     4.00 µS/cm
  temp.coeff.:    2.5 %/°C
>PARAM/analog output
  polarity:       +
  offset:         0 %fs
  damping:        off
>PARAM/plot
  auto start:     off
  time interval   1.0 s
  time scale      10.0 mm/min
  time scale label: rel
  stop time       off min
  left:           0.000 µS/cm
  right:          4.000 µS/cm
>PARAM/709 IC Pump
  flow            2.00 mL/min
  Pmax            10.0 MPa
  Pmin            1.0 MPa
  flow corr.      1.00

PROGRAM
  progr.typ       inject
  17.9 zero       on
  18.0 flag       end
  -----
```

Example of "method list" report

```
Method list
method name      'PRP-X100'
  progr.typ       inject
  date 1995-10-10 15:15:41
  number of steps      3

  free memory      8(9)
  -----
```

Example of "param" report

Measurement parameters for drinking water determination with PRP-X100 IC Anion Column

```
PARAM
>PARAM/detector
  range:          200 µS/cm
  full scale:     4.00 µS/cm
  temp.coeff.:    2.5 %/°C
>PARAM/analog output
  polarity:       +
  offset:         0 %fs
  damping:        off
>PARAM/plot
  auto start:     off
  time interval   1.0 s
  time scale      10.0 mm/min
  time scale label: rel
  stop time       off min
  left:           0.000 µS/cm
  right:          4.000 µS/cm
>PARAM/709 IC Pump
  flow            2.00 mL/min
  Pmax            10.0 MPa
  Pmin            1.0 MPa
  flow corr..     1.00
  -----
```

Example of "program" report

Program parameters for drinking water determination with the PRP-X100 IC Anion Column

```
PROGRAM
  progr.type:     inject
  17.9 zero:      on
  18.0 flag:      end
  -----
```

Example of "pump" report

Report of the parameters for the 709 IC Pump

```
PUMP
  status          running
  Pactual         6.9 MPa
  error           no error
  -----
```

4.9 Examples of methods

4.9.1 Cation determination with Metrosep Cation 1-2

This example shows a method for the determination of monovalent and divalent cations with the 6.1010.000 IC Cation Column Metrosep Cation 1-2, which can be used with or without an autosampler.

General conditions

- **Separating column:** 6.1010.000 Metrosep Cation 1-2
- **Sample loop:** 10 μ L
- **Eluent:** 4 mmol/L tartaric acid,
1 mmol/L dipicolinic acid
- **Flow:** 1 mL/min
- **Eluent conductivity:** ca. 670 μ S/cm
- **Full scale:** 5 μ S/cm

Configuration

```

CONFIG
>CONFIG/detector
  thermostat:          35 °C
  'zero' unit:          $\mu$ S/cm
  cell constant       17.1 /cm
>CONFIG/printer
  id.1
  id.2
  print header:       once
  date&time:          on
  send to:             IBM
>CONFIG/print meas.value
  print crit.:        off
>CONFIG/auxiliaries
  run number          0
  number of cycles    1
>CONFIG/aux/event
  enabled in program: off
  date                1995-10-16
  time                10:59:20
  dialog:             english
  device label
  program             732.0012
>CONFIG/aux/beep
  status:              on
  repeat time         60 s
>CONFIG/RS settings.
  baud rate:          9600
  data bit:           8
  stop bit:           1
  parity:             none
  handshake:          HWS
  control via RS:     on
>CONFIG/RS settings 709
  baud rate:          9600
>CONFIG/733 IC Sep.Cent.
  valve A
  control:            no restriction
>CONFIG/733/valve A
  status              inject
  trigger:            inject
  -----

```

Method

```

METHOD
method name      Cat.1-2
date: 1995-10-16 10:57:42


PARAM
>PARAM/detector
range:          1.00 mS/cm
full scale:     5.00 µS/cm
temp.coeff.:    1.5 %/°C
>PARAM/analog output
polarity:       -
offset:         0 %fs
damping:       off
>PARAM/plot
auto start:     off
time interval   1.0 s
time scale      10.0 mm/min
time scale label: rel
stop time       off min
left:           1.000 µS/cm
right:          -4.000 µS/cm
>PARAM/709 IC Pump
flow            1.00 mL/min
Pmax            10.0 MPa
Pmin            1.0 MPa
flow corr.     1.00

PROGRAM
progr.type:     inject
11.9 zero:      on
12.0 flag:      end
-----
    
```

Triggering of auto-zero function
Program end

The program of the type "inject" is started each time injection valve A is switched to the "INJECT" position. After 12.0 min, the IC system is ready for the next injection.

Chromatogram for cation standard



FULL REPORT

Ret Time (Min)	Component Name	Concentr. ppm	Area (uV*Sec)	Height (uV)
2.682	Li	1.000000	3100993.50	537565.125
3.252	Na	5.000000	4215109.00	627987.062
3.624	NH ₄	5.000000	4517200.50	508054.437
4.566	K	10.000000	4022451.00	313609.843
6.937	Ca	10.000000	6971803.50	396256.875
9.982	Mg	10.000000	15063622.0	501869.562

Fig. 22: Ion chromatogram for cation standard with Metrosep Cation 1-2

4.9.2 Anion determination with Metrosep Anion Dual 2

This example shows a method for the determination of monovalent and divalent anions with the 6.1006.100 IC Anion Column Metrosep Anion Dual 2 and the Metrohm suppressor module, which can be used with or without an autosampler.

General conditions

- **Separating column:** 6.1006.100 Metrosep Anion Dual 2
- **Sample loop:** 20 µL
- **Eluent:** 1.3 mmol/L sodium carbonate,
2.0 mmol/L sodium bicarbonate
- **Flow:** 0.8 mL/min
- **Eluent conductivity:** ca. 14 µS/cm
- **Full scale:** 50 mS/cm and 5 mS/cm
- **Suppressor module:** Regeneration solution 20 mmol/L H₂SO₄,
rinsing solution Milli-Q water

Configuration

```

CONFIG
>CONFIG/detector
  thermostat:          35 °C
  'zero' unit:         µS/cm
  cell constant        17.1 /cm
>CONFIG/printer
  id.1
  id.2
  print header:        once
  date&time:           on
  send to:              IBM
>CONFIG/print meas.value
  print crit.:         off
>CONFIG/auxiliaries
  run number           0
  number of cycles     1
>CONFIG/aux/event
  enabled in program:  off
  date                 1995-10-16
  time                 10:59:20
  dialog:              english
  device label
  program              732.0012
>CONFIG/aux/beep
  status:              on
  repeat time          60 s
>CONFIG/RS settings
  baud rate:           9600
  data bit:            8
  stop bit:            1
  parity:              none
  handshake:           HWS
  control via RS:      on
>CONFIG/RS settings 709
  baud rate:           9600
>CONFIG/733 IC Sep.Cent.
  valve A + suppressor
  control:             no restriction
>CONFIG/733/valve A
  status               inject
  trigger:             inject
>CONFIG/733/suppressor
  auto step:           fill
  status               in position
  trigger:             ---
  
```

Method

METHOD	
method name	An.Dual
date:	1995-10-16 14:30:39
PARAM	
>PARAM/detector	
range:	100 $\mu\text{S}/\text{cm}$
full scale:	50.00 $\mu\text{S}/\text{cm}$
temp.coeff.:	2.5 %/ $^{\circ}\text{C}$
>PARAM/analog output	
polarity:	+
offset:	0 %fs
damping:	off
>PARAM/Plot	
auto start:	off
time interval	1.0 s
time scale	10.0 mm/min
time scale label:	rel
stop time	off min
left:	-5.000 $\mu\text{S}/\text{cm}$
right:	45.00 $\mu\text{S}/\text{cm}$
>PARAM/709 IC Pump	
flow	0.80 mL/min
Pmax	7.0 MPa
Pmin	1.0 MPa
flow corr.	1.00
PROGRAM	
progr.type:	cycle
number of cycles	1
9.0 range:	100 $\mu\text{S}/\text{cm}$
full scale:	5.00 $\mu\text{S}/\text{cm}$
mark:	on
14.7 range:	100 $\mu\text{S}/\text{cm}$
full scale:	50.0 $\mu\text{S}/\text{cm}$
mark:	on
19.9 zero:	on
49.9 valve A:	fill
remote	0*****
pump R/S:	off
50.0 flag:	end

Switching of full-scale range to 5 $\mu\text{S}/\text{cm}$
 Triggering of a marking signal
 Switching of full-scale range to 50 $\mu\text{S}/\text{cm}$
 Triggering of a marking signal
 Triggering of the auto-zero function
 Switching of valve A to "FILL"
 Switching off of 752 Pump Unit
 Switching off of 709 IC Pump
 Program end

The program of the type "inject" is started each time injection valve A is switched to the "INJECT" position. For the determination of nitrate and phosphate, the full-scale range is briefly narrowed to 5 $\mu\text{S}/\text{cm}$ and this range marked by marking signals. After 20.0 min the IC system is ready for the next injection. If there is no new program start within 50 min after the last injection, valve A is automatically switched to the "FILL" position, the 709 IC Pump and 752 Pump Unit are stopped and the program ended.


Settings for 709 IC Pump

- Installation: as shown in *Fig. 10* (see *section 2.6*)
- <EXT.> key: Switched on (remote operation activated)

Settings for 752 Pump Unit

- Installation: as described in *section 2.8.6*

Chromatogram for anion standard



FULL REPORT

Ret Time (Min)	Component Name	Concentr. ppm	Area (uV*Sec)	Height (uV)
4.918	Chloride	19.661419	4187692.25	331573.781
11.505	Nitrate	2.010445	1906272.75	68276.7421
13.399	Phosphate	1.973833	1011560.75	30180.7343
16.617	Sulphate	19.895790	2717087.00	69003.4375

Fig. 23: Ion chromatogram for anion standard with Metrosep Anion Dual 2

5 Notes – Maintenance – Faults

5.1 Practical notes on ion chromatography

5.1.1 Separating columns

Separation efficiency

The attainable quality of analyses with the 732/733 IC system depends to a large extent on the separation efficiency of the column used. When purchasing an IC column you should ensure that the separation efficiency suffices for the analysis problems at hand. Ascertain the **characteristic data of the IC column** on the standard chromatogram enclosed with the column such as capacity factors, selectivity, plate number and resolution and check these data with your own measurements. If any difficulties arise, you should always first check the quality of the column by recording a **standard chromatogram**.

You will find additional general tips on handling IC separating columns in the **8.732.2003 Metrohm Monograph "Ion chromatography"**, as well as detailed information on the separating columns available from Metrohm (see *section 7.3.2*) in the leaflets supplied and in the special **Application Bulletins** available free of charge from your nearest Metrohm agency.

Protection

To protect the column against foreign particles which could have an adverse influence on the separation efficiency, we advise you to subject both the eluents and all samples to **microfiltration** (0.45 µm filter) and to siphon the eluent through the **6.2821.090 Aspirating Filter**.

To avoid contamination by abrasive particles arising from piston seals of the 709 IC Pump, it is advantageous to install an **in-line filter** between the pump and the 733 IC Separation Center. **The 6.2821.100 Filter unit PEEK** (see *section 2.6.3*) supplied with the 709 IC Pump for operation with PEEK capillaries in the pressure range 0...25 MPa and the optional **6.2821.000 Filter unit Manufit** (see *section 2.6.4*) for operation with steel capillaries in the pressure range 0...50 MPa are eminently suitable for this purpose.

The use of readily interchangeable **precolumns** serves to protect the actual separating columns and increase their service life appreciably. The precolumns available from Metrohm (see *section 7.3.2*) are either actual precolumns or so-called precolumn cartridges which are used together with the 6.2821.050 Twin cartridge holder or the 6.2828.010 Precolumn cartridge holder (see *section 2.7*).

Storage

Always store the separating columns closed when not in use and filled in accordance with the manufacturer's specifications.

Dead volume

Dead volume at the end of a column can be the cause of extreme peak broadening or splitting (appearance of double peaks). Filling the column with glass beads ($\varnothing \leq 100 \mu\text{m}$) frequently improves the separation efficiency.

Regeneration

If the separation properties of the column have deteriorated, it can be regenerated in accordance with the column manufacturer's specifications. With the separating columns available from Metrohm (see *section 7.3.2*), the instructions for regeneration can be found on the leaflet enclosed with every column.



*In the case of separating columns with carrier material based on silica, **only solutions with pH 2½7** may be used for regeneration, otherwise the columns could be damaged.*

5.1.2 Pumps

Pulsation dampener

For sensitive measurements with the 732/733 IC System, high pressure pumps as free from pulsation as possible and with a very constant capacity are required. If the pulsations that appear are too large, a pulsation dampener between the pump and the 733 IC Separation Center may possibly help. The optional **6.2620.150 Pulsation dampener MF** (see *section 7.3.1*) whose installation is described in *section 2.6.2* is eminently suitable for this purpose. A pulsation dampener also offers protection against pressure shocks on the column material caused by injection.

Maintenance

To protect the pump against foreign particles, we advise you to subject the eluent to **microfiltration** (0.45 μm filter) and siphon the eluent through the **6.2821.090 Aspirating Filter**.

In many cases, an unstable baseline (pulsation, flow fluctuations) can be traced to contaminated valves or faulty, leaky piston seals.

Contaminated valves are cleaned by rinsing with water, RBS solution or organic solvents. The rinsing effect is improved by brief treatment in an ultrasonic bath. When the cleaned valves are reinstalled, you must ensure that the flow direction is correct.

The **replacement of piston seals** has to be done in accordance with the pump manufacturer's directions. The corresponding maintenance work for the 709 IC Pump is described in *section 4.2* of the *709 Instructions for Use*.

Salt crystals between the piston and the seal are the cause of abrasive particles, which can enter the eluent. These lead to contaminated valves, pressure rise and in extreme cases to scratched pistons. It is thus essential to ensure that **no precipitates** can appear (see also *section 5.1.3*).

5.1.3 Eluents
Treatment

For the preparation of the eluents only chemicals of a purity degree of at least "p.a." should be used. For diluting please use only **high purity water**.

Fresh eluents should always be **microfiltered** (0.45 µm filter) and **de-gassed** (with N₂, He or vacuum). For high sensitive measurements, the eluent should be continuously **stirred** with a magnetic stirrer, particularly when the recycling procedure is employed or when alkaline eluents are used. For alkaline eluents and eluents with low buffering capacity one should preferably use a **CO₂ absorber** (e.g. the absorber supplied with the optional 6.5324.000 Bottle rack).

The supply vessel containing the eluent must be closed as tightly as possible to avoid excessive evaporation. This is primarily important with eluents containing organic solvents (e.g. acetone), the evaporation of which can lead to drifts in the long term. If work is performed in a very sensitive range, even if one drop of condensate falls back in the eluent this can cause a noticeable change in the background conductivity.

Influence of various parameters on anion columns

- *Concentration:* An increase in the concentration usually leads to shorter retention times and quicker separation, but also to a higher background conductivity.
- *pH:* pH alterations lead to shifts in the dissociation equilibrium and thus to changes in the retention times.
- *Organic modifiers:* Addition of an organic solvent (e.g. methanol, acetone, acetonitrile) to aqueous eluents generally accelerates lipophilic ions.

Eluent change

When the eluent is changed, it must be ensured that **no precipitates** can be formed. Solutions used in direct succession must therefore be miscible. If the system has to be rinsed with an organic solution, several solvents with increasing or decreasing lipophilic character may possibly have to be used (e.g. water ↔ acetone ↔ chloroform).

5.1.4 Suppressor module

To avoid contamination of the suppressor module by foreign particles or bacterial growth, it is advantageous to install an **in-line filter** between the 752 Pump Unit and the suppressor module. The two **6.2821.100 Filter units PEEK** (see *section 2.6.3 and section 2.8.3*) supplied with the 2.733.0X30 IC Separation Center are eminently suitable for this purpose.

5.1.5 Connections

All connections between injector, column and detector must be as short as possible, have a low dead volume and be absolutely tight. The PEEK capillary after the detector block must be free from constriction (the measuring cell is tested to 5 MPa = 50 bar back pressure).

5.2 Maintenance and servicing

5.2.1 General information

Care

The 732 IC Detector and 733 IC Separation Center require proper care and attention. Excessive contamination of the instruments could possibly lead to malfunctions and a shorter service life of the inherently rugged mechanical and electronic parts.

Spilled chemicals and solvents should be wiped up immediately. It is especially important to protect the plug connections at the rear of the instrument (particular the mains plug) against contamination.



Although constructional measures have been designed to virtually eliminate such a situation, should corrosive media penetrate the interior of the instruments the mains plug of the 732 IC Detector must be immediately disconnected to prevent extensive damage to the instrument electronics. Inform Metrohm service if your instrument(s) have been damaged in such a way.



The instrument must not be opened by untrained personnel. Please comply with the safety notes in section 1.4.1.

Maintenance by Metrohm service

Maintenance of the 732 IC Detector and 733 IC Separation Center is best done as part of an annual service performed by specialists from the Metrohm company. If work is frequently performed with caustic and

corrosive chemicals, it may be necessary to shorten the interval between servicing.

The Metrohm service department is always willing to offer expert advice on the maintenance and servicing of all Metrohm instruments.

5.2.2 Passivation

Passivation of the entire IC system (without column) by rinsing with 20...50 mL 0.2 mol/L HNO₃ is appropriate in the following cases:

- in the installation
(only with non metal-free IC systems, see *section 2.6.7*)
- if the separation system is changed
(only with non metal-free IC systems, see *section 5.2.5*)
- if exceptional changes are observed in the measurement properties of the cell
(with metal-free and non metal-free IC systems).

For passivation the separating column is removed from the 733 IC Separation Center. The two capillaries **67** and **82** (see *Fig. 16* and *Fig. 17*) are directly connected to each other using the 6.2620.060 Coupling enclosed with the accessories.

5.2.3 Recycling

To keep the eluent consumption between injections to a minimum when the system is at rest (e.g. overnight), the so-called recycling procedure can be used. In recycling the eluent exiting the outlet capillary of the detector block is led back directly to the eluent supply vessel. The IC system is thus quickly ready for new injections without a long conditioning period.



The recycling procedure must **not** be used

- in operation with the suppressor module,
- with alkaline eluents,
- with the 6.1010.000 IC Cation column METROSEP Cation 1-2.

5.2.4 Shutdown

If the 732/733 IC-System is shut down for a considerable length of time, the entire IC system (without column and suppressor) must be **rinsed free from salt** with methanol/water (1:4) to avoid crystallization of eluent salts with the corresponding subsequent damage.

For rinsing the separating column is removed; the two capillaries **67** and **82** (see *Fig. 16* and *Fig. 17*) are directly connected to each other with the 6.2620.060 Coupling enclosed with the accessories. Rinse with methanol/water (1:4) until the conductivity drops below 10 µS/cm.

5.2.5 Changing separating columns

Identical separation system

If you wish to replace an IC separating column by a column of the same type, proceed as follows (see *Fig. 16* and *Fig. 17*):

1 Remove old column

- Switch off pump drive of the 709 IC Pump.
- Unscrew column from inlet capillary **82** of the detector block or from suppressor inlet capillary **89**.
- Unscrew column from column connection capillary **67** or the precolumn.

2 Connect new column to injector

- Remove end caps from column **76**.
- Screw inlet end of separating column **76** (note flow direction) to column connection capillary **67** or to the precolumn (see *section 2.8.4/section 2.8.6*).

3 Rinse column

- Place beaker beneath the column outlet.
- Switch on 709 IC Pump and rinse column with eluent for ca. 10 min, then switch off pump.

4 Connect column to detector block

- Screw outlet end of separating column **76** to inlet capillary **82** or suppressor inlet capillary **89**.

Changing the separation system

If you wish to replace an IC separating column by a column of a different type, proceed as follows (see *Fig. 16* and *Fig. 17*):

1 Remove old column

- Switch off pump drive of 709 IC Pump.
- Unscrew column from inlet capillary **82** of the detector block or from the suppressor inlet capillary **89**.
- Unscrew column from column connection capillary **67** or the precolumn.

2 Connect detector block to injection valve

- Connect column connection capillary **67** using a 6.2620.060 Coupling directly to inlet capillary **82** of the detector block (see *Fig. 16*).

3 Rinse with HNO₃ (passivation)

- Rinse IC system with $c(\text{HNO}_3) = 0.2 \text{ mol/L}$ (flow rate 1 mL/min) for ca. 10 min.

4 Rinse with dist. H₂O

- Rinse IC system with dist. Or demin. Water (flow rate 1 mL/min) for ca. 10 min.

5 Rinse with eluent

- Rinse IC system with the eluent needed for the separating column you intend to install until the absolute conductivity displayed on the 732 IC Detector is stable.

6 Remove coupling

- Remove 6.2620.060 Coupling between column outlet capillary **67** and inlet capillary **82** (see *Fig. 16*).

7 Connect new column to injector

- Remove end caps from column **76**.
- Screw inlet end of separating column **76** (note flow direction) to column connection capillary **67** or to the precolumn (see *section 2.8.4/section 2.8.6*).

8 Rinse column

- Place beaker beneath the column outlet.
- Switch on 709 IC Pump and rinse column with eluent for ca. 10 min, then switch off pump.

9 Connect column to detector block

- Screw outlet end of separating column **76** to inlet capillary **82** or suppressor inlet capillary **89**.

5.2.6 Regeneration of the suppressor

Regenerating a suppressor operating at reduced capacity

If the suppressor units are exposed to certain heavy metals (e.g. iron) or organic contaminants for long periods of time, these can no longer be completely removed by the regeneration solution normally used (20 mmol/L H_2SO_4). This diminishes the capacity of the suppressor units, which, in milder cases, results in a reduced sensitivity to phosphates and, in severe cases, in a strong increase in the baseline. If such capacity problems occur at one or several positions, the suppressor units must be treated as follows:

1 Disconnect suppressor from IC system

- Disconnect suppressor from separating column and detector.

2 Regenerate suppressor

- Rinse each suppressor unit for about 15 min with one of the following solutions:

Contamination with heavy metals

1 mol/L H_2SO_4

Contamination with organic cationic complexing agents

0.1 mol/L H_2SO_4 / 0.1 mol/L oxalic acid / acetone 5%

Severe contamination with organic substances

0.2 mol/L H_2SO_4 / acetone \geq 20%



The 6.1826.050 pump tubing is made of PVC and must not be used for rinsing with solutions which contain acetone. In such cases, rinse with different pump tubing or a different pump.

3 Connect suppressor to IC system

- Reconnect suppressor to the IC system. If capacity problems persist, replace the suppressor rotor (see section 5.2.8).

Regenerating the suppressor in case of high counterpressure

If excessive counterpressure is observed in one or several suppressor units, treat the units as follows:

1 Disconnect suppressor from IC system

- Disconnect suppressor from separating column and detector.

2 Regenerate suppressor

- Connect the inlet capillary **90** marked "H₂SO₄" to the 709 IC Pump using the 6.2620.060 or 6.2744.040 coupling.
- Set the flow at the 709 IC Pump to 0.5 mL/min and rinse the suppressor unit with 1 mol/L H₂SO₄ for 5 to 10 min.
- As the pressure falls, slowly increase the flow at the 709 IC Pump to 2 mL/min. Do not exceed a maximum pressure of 2 MPa (20 bar).
- Switch suppressor to the next position using the <STEP> key.
- Connect the inlet capillary **93** marked "H₂O" to the 709 IC Pump using coupling 6.2620.060 or 6.2744.040.
- Set the flow at the 709 IC Pump to 0.5 mL/min and rinse the suppressor unit with distilled H₂O for 5 to 10 min.
- As the pressure falls, slowly increase the flow at the 709 IC Pump to 2 mL/min. Do not exceed a maximum pressure of 2 MPa (20 bar).

3 Connect suppressor to IC system

- Reconnect the suppressor to the IC system and the 752 Pump Unit. If the pressure problems persist, replace the suppressor rotor (see *section 5.2.8*).

5.2.7 Cleaning the suppressor

It may be necessary to clean the suppressor in the following cases:

- High counterpressure on the suppressor connection tubing
- Irremediable blockage of the suppressor (the suppressor can no longer deliver solutions)
- Irremediable obstruction of the suppressor (the suppressor can no longer be switched to next position)

To clean the connection piece and the rotor, proceed as follows (see *Fig. 24*):

1 Disconnect suppressor from IC system

- Disconnect input lead of suppressor from the separating column and output lead to the detector.

2 Dismantle suppressor

- Unscrew nut **97** from suppressor holder **100**.
- Pull out connection piece **98** and suppressor rotor **99** from suppressor holder **100** (the connection piece and the rotor normally stick together).
- Loosen connection piece **98** from suppressor rotor **99**.

3 Clean input and output leads

- Connect each of the 6 capillary tubings attached to connection piece **98** to the pump one after another, and pump through ultrapure water.
- Check whether solution emerges from connection piece **98**. If one of the input or output leads remains blocked, replace the connection piece (order number 6.2832.010).

4 Clean suppressor rotor

- Clean the sealing surface of suppressor rotor **99** using a lint-free cloth and ethanol.

5 Insert suppressor rotor

- Insert suppressor rotor **99** in suppressor holder **100** in such a way that the tubing connections at the rear of the rotor fit in the corresponding openings inside the rotor, and that one of the three holes in the rotor can be seen from below in one of the openings of the holder.
- If the rotor has been inserted correctly, its sealing surface will be about 4 mm inside the holder. If this is not the case, bring the rotor into the correct position from below with the aid of a sharp object (e.g. a screwdriver).

6 Clean connection piece

- Clean the sealing surface of connection piece **98** using a lint-free cloth and ethanol.




Fig. 24: Assembling the suppressor

97	Screw nut	99	Suppressor rotor (6.2832.000)
98	Connection piece (6.2832.010) with input and output leads	100	Suppressor holder

7 Insert connection piece

- Insert connection piece **98** in suppressor holder **100** in such a way that connection "1" is at the top, and that the three lugs on the connection piece fit in the corresponding openings of the holder.
- Screw nut **97** onto the thread of suppressor holder **100** manually (do not use tools).

8 Connect and condition the suppressor

- Reconnect the suppressor to the IC system.
- Before switching the suppressor to the next position for the first time, rinse all 3 suppressor units with solution for 5 min.

5.2.8 Replacing the suppressor

The suppressor in the suppressor block may have to be replaced in the following cases:

- Irremediable loss of suppressor capacity (reduced phosphate sensitivity and/or strong rise in baseline)
- Irremediable blockage of the suppressor (the suppressor can no longer deliver solutions)

Both the 6.2832.000 Suppressor rotor and the 6.2832.010 Connection piece with the input and output leads can be replaced. To replace these components proceed as follows (see *Fig. 24*):

1 Disconnect suppressor from IC system

- Disconnect all input and output leads of the suppressor from IC system and the pump of the 752 Pump Unit.

2 Dismantle suppressor

- Unscrew nut **97** from suppressor holder **100**.
- Pull out connection piece **98** and suppressor rotor **99** from suppressor holder **100** (the connection piece and the rotor normally stick together).
- Loosen connection piece **98** from suppressor rotor **99**.

3 Clean suppressor rotor

- Clean the sealing surface of new suppressor rotor **99** (6.2832.000) using a lint-free cloth and ethanol.

4 Insert suppressor rotor

- Insert new suppressor rotor **99** in suppressor holder **100** in such a way that the tubing connections at the rear of the rotor fit in the corresponding openings inside the rotor and that one of the three holes in the rotor can be seen from below in one of the openings of the holder.

- If the rotor has been inserted correctly, its sealing surface will be about 4 mm inside the holder. If this is not the case, bring the rotor into the correct position from below with the aid of a sharp object (e.g. a screwdriver).

5 Clean connection piece

- Clean the sealing surface of new connection piece **98** (6.2832.010) with the aid of a lint-free cloth and ethanol.

6 Insert connection piece

- Insert new connection piece **98** in suppressor holder **100** in such a way that connection "1" is at the top and that the three lugs on the connection piece fit in the corresponding openings of the holder.
- Screw nut **97** onto the thread of suppressor holder **100** manually (do not use tools).

7 Connect and condition the suppressor

- Reconnect the suppressor to the IC system.
- Before switching the suppressor to the next position for the first time, rinse all 3 suppressor units with solution for 5 min.

5.3 Faults and malfunctions

5.3.1 Error messages

If any type of malfunction occurs during operation of the 732 IC Detector, this is shown by an error message which appear in the top line of display **1** and (if switched on under <CONFIG>) by three beeps. To return to the instrument standby mode, the error messages must be confirmed with the <QUIT> key.

You will find further details of the error messages, their possible causes and the procedure for their rectification in the following table, which is arranged alphabetically. If the errors can not be rectified, please follow the procedure described in *section 5.4*.

Error message	Meaning/Causes	Rectification
adjustment error	Instrument adjustment invalid.	Inform Metrohm service.
all methods faulty	The stored methods can no longer be recalled.	Re-enter methods and store. If this error occurs frequently, inform Metrohm service.
connect 709 IC Pump	The active program contains a command for the 709 IC Pump, which can not be addressed at present.	Check connection to 709 IC Pump, switch on pump and remote operation with <EXT.>.
connect 733 Sep.Center	The active program contains a command for the 733 IC Separation Center, which can not be addressed at present.	Check connection to 733 IC Separation Center.
DAC coarse error #	Error in internal hardware test.	Inform Metrohm service.
DAC fine error #	Error in internal hardware test.	Inform Metrohm service.
EEPROM error 233	Error on storage of the EEPROM block "Program".	Re-enter program.
EEPROM error 234	Error on storage of the EEPROM block "Parameters".	Re-enter parameters.
EEPROM error 235	Error on storage of the EEPROM block "Method".	Re-enter method (program and parameters).
EEPROM error 236	Error on storage of the EEPROM block "Event".	Re-enter events.
EEPROM error 237	Error on storage of the EEPROM block "Config".	Re-enter configuration.
EEPROM error 238	Error on storage of the EEPROM block "Device number".	Re-enter device number.
EEPROM error 239	Error on storage of the EEPROM block "Adjustment".	Inform Metrohm service.
EEPROM error 240	Error on storage of the EEPROM blocks.	Re-enter all values.
Err1	Program checksum wrong.	Inform Metrohm service.
Err3	Instrument adjustment faulty.	<QUIT>. If new instrument adjustment necessary, inform Metrohm service.
Err4	Timer interrupt for multitasking missing.	Inform Metrohm service.

<i>Error message</i>	<i>Meaning/Causes</i>	<i>Rectification</i>
Err5	RS232 interface at 732 IC Detector faulty.	Inform Metrohm service.
Err6	Clock faulty.	Inform Metrohm service.
Err7	LCD write/read error.	Inform Metrohm service.
Err8	AD converter faulty.	Inform Metrohm service.
ErrC	RS232 interface at 709 IC Pump faulty.	Inform Metrohm service.
faulty method #	Method # (0...9) can no longer be recalled.	Re-enter method and store.
inject in inject program	The active inject program contains an "inject" command. As this itself triggers the program start, the program will be executed endlessly.	Delete "inject" command in the program.
I0 error 50	Error at input 1 or output 1.	Inform Metrohm service.
I0 error 51	Error at input 2 or output 2.	Inform Metrohm service.
I0 error 52	Error at input 3 or output 3.	Inform Metrohm service.
I0 error 53	Error at input 4 or output 4.	Inform Metrohm service.
I0 error 54	Error at input 5 or output 5.	Inform Metrohm service.
I0 error 55	Error at input 6 or output 6.	Inform Metrohm service.
I0 error 56	Error at input 7 or output 7.	Inform Metrohm service.
I0 error 57	Error at input 8 or output 8.	Inform Metrohm service.
I0 error 58	Frequency less than lower limit frequency of an RC element.	Inform Metrohm service.
I0 error 59	Frequency greater than upper limit frequency of an RC element.	Inform Metrohm service.
method memory full	The memory for the user-defined methods is full.	Delete unused or seldom needed methods.
missing 'end'	An "end" command is missing in the active program.	Insert program step with "Flag = end" in the active program. Without an "end" command, the program continues to run indefinitely until stopped by <PROG R/S>.
missing 'return'	In the active loop program there is no "return" command so that the program runs only once.	Enter "return" command in the program.
multiple 'end'	Active program contains several "end" commands.	Delete superfluous "end" commands in the program.
multiple 'return'	Active program contains several "return" commands.	Delete superfluous "return" commands in the program.
plot data overflow	Volume of data too great for graphics plot.	<QUIT> (graphics plot aborted).
program step after 'end'	The active program contains additional program steps after an "end" command which can no longer be executed.	Delete program steps after "end" or increase time for "end" command.
pump not responding	No connection to 709 IC Pump or baud rate set wrongly.	Connect 709 IC Pump to 732 IC Detector, switch on pump and remote control with <EXT.> or check baud rate of 709 and 732.
RAM error	Error in RAM test.	Inform Metrohm service.

Error message	Meaning/Causes	Rectification
range error #	Error in internal hardware test.	Inform Metrohm service.
RS error 36	RS receive error; wrong parity.	<QUIT>, set same parity for both devices.
RS error 37	RS receive error; wrong stop bit.	<QUIT>, set same stop bit for both devices.
RS error 38	RS-receive error; overflow (at least 1 character could not be read).	<QUIT>, set same baud rate for both devices, restart sender.
RS error 39	RS receive error; overflow of internal receive buffer (>82 characters).	<QUIT>.
RS error 40	RS transmission error; DSR = OFF. Handshake unsatisfactory for more than 1 s.	<QUIT>, check receiver (switched on and ready?).
RS error 41	RS transmission error; DCD = ON. Handshake unsatisfactory for more than 1 s.	<QUIT>, check receiver (switched on and ready?).
RS error 42	RS transmission error; CTS = OFF. Handshake unsatisfactory for more than 1 s.	<QUIT>, check receiver (switched on and ready?).
RS error 43	RS transmission error; transmission of the IC Detector interrupted with XOFF for at least 3 s.	<QUIT> or send XON.
RS error 44	RS transmission error; the RS232 parameters for the two devices are no longer the same.	<QUIT>, reset RS232 parameters for both devices.
RS error 45	RS transmission error; the receive buffer of the IC Detector contains an incomplete string (L _F missing), the transmission is thus blocked.	<QUIT> or send L _F .
RS error 60 ... 89	Error in RS232 test.	<QUIT>, check connection between the two RS232 interfaces.
RS error XX at 709	RS error 36...89 at the RS232 interface "709 IC Pump".	See individual RS errors.
RS error XX at RS 232	RS error 36...89 at the RS232 interface "RS 232".	See individual RS errors.
sensitivity error #	Error in internal hardware test.	Inform Metrohm service.
stopped:YY-MM-DD HH:MM	709 IC Pump stopped at specified time owing to violation of shutoff limits.	Press the 709 IC Pump keys <EXT.>, <R/S> in succession and then <EXT.> again.
store operation failed	Storage of the method failed.	Switch instrument off/on. Store method again. If error reappears, inform Metrohm service.
U cell error #	Error in internal hardware test.	Inform Metrohm service.

5.3.2 Malfunctions and their rectification

If difficulties appear with the IC system during analyses, their causes are best investigated in the order **separating column** ® **pump** ® **eluent** ® **732/733 IC System**. Several of the malfunctions which may appear are listed in the following table with details of possible causes and countermeasures.

<i>Malfunction</i>	<i>Cause</i>	<i>Rectification</i>
Baseline with high noise level, pulsation	<ul style="list-style-type: none"> Contaminated pump valves Faulty piston seals Quality of the pump does not suffice for the selected sensitivity 	<ul style="list-style-type: none"> Clean the valves (see <i>section 5.1.2</i>) Replace the piston seals (see <i>section 5.1.2</i>) Use pulsation dampener, use more powerful pump or lower the sensitivity
Drift of the baseline	<ul style="list-style-type: none"> Thermal equilibrium not yet reached Leak in system Evaporation of organic solvent in eluent 	<ul style="list-style-type: none"> Condition system with heating switched on Check connections and make leakproof Ensure better closure of eluent supply vessel
Considerable pressure drop	<ul style="list-style-type: none"> Leak in system 	<ul style="list-style-type: none"> Check connections and make leakproof
Considerable pressure rise	<ul style="list-style-type: none"> Contamination of the filter in the 6.2821.100 Filter unit PEEK Contamination of the filter in the 6.2821.000 Filter unit Manufit Contamination of the column inlet filter Change of column packing by injection of contaminated samples 	<ul style="list-style-type: none"> Replace the 6.2821.110 Filter (see <i>section 2.6.3</i>) Clean or replace 6.2821.020 Steel mesh(es) Clean or replace 6.2821.020 Steel mesh(es) Regenerate the column (see <i>section 5.1.1</i>) or replace column <p><i>Note:</i> Samples should always be microfiltered.</p>
Chromatograms with poor resolution, change in the retention times	<ul style="list-style-type: none"> Deterioration in separation efficiency of the IC column 	<ul style="list-style-type: none"> Regenerate the column (see <i>section 5.1.1</i>) or replace column
Extreme peak broadening, splitting (double peaks)	<ul style="list-style-type: none"> Dead volume at the column ends 	<ul style="list-style-type: none"> Fill column with glass beads ($\varnothing \leq 100 \mu\text{m}$) or replace column
No feed of regeneration or rinsing solution for suppressor	<ul style="list-style-type: none"> Contamination of the filter in the 6.2821.100 Filter unit PEEK Counterpressure too high in suppressor module 	<ul style="list-style-type: none"> Replace the 6.2821.110 Filter (see <i>section 2.6.3</i>) Clean or replace suppressor module (see <i>section 5.2.6</i> / <i>5.2.8</i>)

5.4 Diagnostic tests

5.4.1 General

The 732 IC Detector and 733 IC Separation Center are very precise and dependable instruments. Thanks to their rugged construction, it is highly unlikely that external mechanical or electrical influences will have any adverse effect on their functions.

Although a fault in the instruments can not be excluded with certainty, the possibility is greater that malfunctions are caused by wrong operation or handling, through improper connections and the operation with third-party devices.

Whatever the case, it is always advisable to localize the fault with the diagnostic tests, which can be performed quickly and simply. The customer need call Metrohm service only when the instrument really has a fault. Further, he can use the results of the specific diagnostic function to provide the service engineer with much more precise information.

In the case of inquiries, always quote the serial number on model plate **10** of the 732 IC Detector (see *Fig. 3*) or on model plate **47** of the 733 IC Separation Center (see *Fig. 5*), program version (see *section 4.4.2*) and if applicable the error message.

Procedure

The diagnostic menu listed in *section 5.4.2* shows all components for which detailed instructions (diagnostic steps) are available for checking the functionality.

In the case of a possible malfunction, we advise you to perform either the corresponding diagnostic step or all diagnostic steps as a routine step on the instrument.

The reactions of the instruments to the instructions must be compared with the descriptions in the diagnostic step. If the instruments do not show the expected reaction ("No" case), the appropriate diagnostic step must be repeated to exclude operating errors. However, it is highly probable that repeated wrong reactions indicate a malfunction.

Instruments needed:

Necessary only if RS232 or remote should be checked:

- **3.496.8510** Test plug (at "Remote" socket)
- **3.496.8480** Test plug (at "RS 232" socket)

5.4.2 Prepare instruments

- Power off.
- Disconnect all external connections at the 732 IC Detector (except mains cable).
- Power on and immediately press and hold <9> key until

```
i n i t i a l i s a t i o n  E E P R O M ?
ENTER>yes '9' >di agnosi s
```

- Press <9> key to open the following diagnostic menu:

di agnosi s >di ag/EEPROM i n i t i a l .	Initialize data memory	1)	<i>sectn 5.4.9</i>
di agnosi s >di ag/run mode	Parameters for sequence control of "internal test" and "external test"	2)	
di agnosi s >di ag/RAM test	Check working memory (RAM)	3)	<i>sectn 5.4.3</i>
di agnosi s >di ag/keyboard test	Check keypad	3)	<i>sectn 5.4.4</i>
di agnosi s >di ag/di splay test	Check display	3)	<i>sectn 5.4.5</i>
di agnosi s >di ag/RS test	Check RS232 interfaces	3)	<i>sectn 5.4.6</i>
di agnosi s >di ag/I O test	Check remote interfaces	3)	<i>sectn 5.4.7</i>
di agnosi s >di ag/i n t e r n a l t e s t	Internal hardware test	3)	<i>sectn 5.4.8</i>
di agnosi s >di ag/external test	External hardware test	2)	
di agnosi s >di ag/i n s t r . a d j u s t m e n t	Instrument adjustment	2)	
di agnosi s >di ag/i n s t r u m e n t n u m b e r	Reading and writing the instrument number	2)	
di agnosi s >di ag/power on reset	Power on reset	2)	

Pressing the <9> key again selects the submenus in turn. The individual inquiries of a submenu are accessed using the <ENTER> key, exit is with the <QUIT> key.

1) Initialization: This diagnostic step is intended for the initialization of the data memory.
 2) Care: These diagnostic steps are reserved for the service engineer and are not described.
 3) Test functions: These diagnostic steps are test functions which the customer can perform.

5.4.3 Checking working memory (RAM)

This diagnostic step executes a non-destructive test over the entire area of the RAM content (working memory).

- Prepare instrument for diagnostic test (see *section 5.4.2*).
- If need be, press <9> key repeatedly until

```
di agnosi s
>di ag/RAM test
```

- <ENTER>

If no faults are found, the display shows:

```
>di ag/RAM test
RAM test ok
```

- <ENTER>

```
di agnosi s
>di ag/keyboard test
```

5.4.4 Check keypad

With this diagnostic step all keys of the keypad are checked for functionality.

- Prepare instrument for diagnostic test (see *section 5.4.2*).
- If need be, press <9> key repeatedly until

```
di agnosi s
>di ag/keyboard test
```

- <ENTER>

```
>di ag/keyboard test
```

- Press all keys in succession and check the reaction on the display.

The display shows the corresponding matrix code for the pressed key. The matrix code can be seen in the key table (e.g. the following display should appear if the <CONFIG> key has been pressed).

```
di ag/keyboard test
code: 2
```

- The test is quit by pressing the <CLEAR> key twice.

```
di agnosi s
>di ag/di spl ay test
```

Key table

Code	Key	Code	Key
0	<9 / METHOD>	12	<CLEAR>
1	<8 / PUMP R/S>	13	<SELECT>
2	<7 / CONFIG>	14	<3 / PRINT>
3	<ZERO>	15	<2 / REPORT>
4	<PARAM>	16	<1 / PLOT>
5		17	
6	<6 / PROGRAM>	18	<ENTER>
7	<5 / EVENT>	19	<QUIT>
8	<4 / ZERO OFF>	20	<-/+ / →>
9	<PROG R/S>	21	<./ ←>
10	<FULL SCALE>	22	<0 / MARK>
11		23	

5.4.5 Check display

This diagnostic step is used to check all light emitting diodes and the display for functionality.

- Prepare instrument for diagnostic test (see section 5.4.2).
- If need be, press <9> key repeatedly until

```
di agnosi s
>di ag/di spl ay test
```

- <ENTER>

After the <ENTER> key has been pressed, the program automatically runs through a test sequence for a visual check on the light emitting diodes and the display.

⇒ The light emitting diodes for OVERLOAD, THERMOSTAT, ZERO and PROG R/S are switched on for a certain time.

P Power on test pattern (each pixel active) appears.

⇒ Both lines of the display are cleared.

⇒ Both lines of the display are written to with the characters "#", "H" and finally with "I".

⇒ Both lines are written to from right to left with the continuous moving display „0123456789ABCDEFGHIJKLMNOPQRSTUVWXYZ“.

P The light emitting diodes for OVERLOAD, THERMOSTAT, ZERO and PROG R/S flash in succession for a brief time.

- The test sequence can be held by pressing any key (with the exception of <9>) and restarted.
- The test is quit with the <9> key.

```
di agnosi s
>di ag/RS test
```

5.4.6 Check RS232 interfaces

With this diagnostic step all outputs and inputs of the two RS232 interfaces **16** "RS 232" and **15** "709 IC Pump" can be checked for functionality.

- Prepare instrument for diagnostic test (see *section 5.4.2*).
- If need be, press <9> key repeatedly until

```
di agnosi s
>di ag/RS test
```

- <ENTER>

```
>di ag/RS test
>di ag/RS test/RS232
```

- <ENTER>

```
>di ag/RS test/RS232
RS connector?
```

- Without switching off the instrument, plug the **3.496.8480** Test plug into socket "RS232".
- <ENTER>

The test runs automatically. If no fault is found, the following appears:

```
>di ag/RS test/RS232
RS test ok
```

- Remove test plug and <ENTER>.

```
>di ag/RS test
>di ag/RS test/709 Pump
```

- <ENTER>

```
>di ag/RS test/709 Pump
RS connector?
```

- Without switching off the instrument, plug the **3.496.8480** Test plug into socket "709 IC Pump".
- <ENTER>

The test runs automatically. If no fault is found, the following appears:

```
>di ag/RS test/709 Pump
RS test ok
```

- Remove test plug and <ENTER>.

```
di agnosi s
>di ag/10 test
```

5.4.7 Check remote interfaces

This diagnostic step can be used to check all outputs and inputs of the two remote interfaces **17** "Remote" and **14** "733 IC Separation Center" for functionality.

- Prepare instrument for diagnostic test (see *section 5.4.2*).
- If need be, press <9> key repeatedly until

```
di agnosi s
>di ag/I0 test
```

- <ENTER>

```
>di ag/I0 test
>di ag/I0 test/remote
```

- <ENTER>

```
>di ag/I0 test/remote
I0 connector?
```

- Without switching off the instrument, plug the **3.496.8510** Test plug into socket "Remote".
- <ENTER>

The test runs automatically. If no fault is found, the following appears:

```
>di ag/I0 test/remote
I0 test ok ..
```

- Remove test plug and <ENTER>.

```
>di ag/I0 test
>di ag/I0 test/733 Sep. C.
```

- <ENTER>

```
>di ag/I0 test/733 Sep. C.
I0 connector?
```

- Without switching off the instrument, plug **3.496.8510** Test plug into socket "733 IC Separation Center".
- <ENTER>

The test runs automatically. If no fault is found, the following appears:

```
>di ag/I0 test/733 Sep. C.
I0 test ok ..
```

- Remove test plug and <ENTER>.

```
di agnosi s
>di ag/i nternal test
```

5.4.8 Internal hardware test

This diagnostic step can be used to check the hardware for measured value acquisition in the 732 IC Detector for functionality. The test is divided into separate tests performed in succession. All these tests are carried out with the possibilities included in the IC Detector and are thus not dependent on working aids such as test plugs or test adapters.

- The instrument must be in operational readiness for the following tests (switched on for at least 15 min).
- Prepare instrument for diagnostic tests (see *section 5.4.2*).
- If need be, press <9> key repeatedly until

```
di agnosi s  
>di ag/i nternal test
```

- <ENTER>

```
>di ag/i nternal test  
>di ag/i nt/U cell
```

Check Ucell

With this diagnostic step the voltage generator (generation of the cell voltage) is checked for functionality.

- <ENTER>

The test runs automatically and lasts ca. 3 min. During this time various measurement results appear in the display, but these should be ignored. If the test discovers results which lie outside the specified tolerance zones, the number of deviations found appears following completion of the test. If all results of the test are within the specified tolerance zones, the following appears:

```
U cell test ok
```

- <ENTER>

```
>di ag/i nternal test  
>di ag/i nt/range
```

Check range

This diagnostic step is used to test the range amplifier for functionality.

- <ENTER>

The test runs automatically and lasts ca. ½ min. During this time various measurement results appear in the display, but these should be ignored. If the test discovers results which lie outside the specified tolerance zones, the number of deviations found appears following completion of the test. If all results of the test are within the specified tolerance zones, the following appears:

```
range test ok
```

- <ENTER>

```
>di ag/i nternal test
>di ag/i nt/sensi tivity
```

Check sense

With this diagnostic step the sense amplifier is checked for functionality.

- <ENTER>

The test runs automatically and lasts ca. 2 min. During this time various measurement results appear in the display, but these should be ignored. If the test discovers results which lie outside the specified tolerance zones, the number of deviations found appears following completion of the test. If all results of the test are within the specified tolerance zones, the following appears:

```
sensitivity test ok
```

- <ENTER>

```
>di ag/i nternal test
>di ag/i nt/DAC coarse
```

Check CompCoarse

With this diagnostic step the DA converter for the coarse adjustment is checked for functionality.

- <ENTER>

The test runs automatically and lasts ca. ½ min. During this time various measurement results appear in the display, but these should be ignored. If the test discovers results which lie outside the specified tolerance zones, the number of deviations found appears following completion of the test. If all results of the test are within the specified tolerance zones, the following appears:

```
DAC coarse test ok
```

- <ENTER>

```
>di ag/i nternal test
>di ag/i nt/DAC fine
```

Check CompFine

With this diagnostic step the DA converter for the fine adjustment is checked for functionality.

- <ENTER>

The test runs automatically and lasts ca. 3 min. During this time various measurement results appear in the display, but these should be ignored. If the test discovers results which lie outside the specified tolerance zones, the number of deviations found appears following completion of the test. If all results of the test are within the specified tolerance zones, the following appears:

```
err: 0
DAC fine test ok
```

- <ENTER>

```
>diagnosi s
>diag/external test
```

5.4.9 Initialize data memory

This diagnostic step can be used to write default values to the instrument parameters using the keypad and thus switch the instrument to the original condition. This measure is important with the following two points:



The setting of certain instrument parameters such as the locking of keys is possible only via RS232, i.e. with the aid of a PC. If such instrument parameters are set and no PC is available to cancel the settings, full use can not be made of the instrument.




In rare cases, it is possible that major interference signals such as line spikes and lightning can have an adverse influence on the contents of the data memory. If the contents of the data memory are undefined, this may lead to a system crash.

With the initialization step, default values can be written into parts of the data memory (EEPROM). Although the adjustment data are retained, the initialization should be performed only if necessary as the stored user data (etc.) are deleted.

- Prepare instrument for diagnostic test (see section 5.4.2).
- If need be, press <9> key repeatedly until

```
diagnosi s
>diag/EEPROM initial.
```

- Press <ENTER> key to open the following diagnostic menu:



Pressing the <SELECT> key selects the submenus in turn. The individual initialization alternatives are accessed using the <ENTER> key, exit is with the <QUIT> key.

The table shows which parts of the data memory are affected by the corresponding initialization alternatives.

- If need be, press the <SELECT> key repeatedly until the desired initialization alternative appears (such as “USER“)

```

>di ag/EEPROM i n i t i a l .
bl ock: USER
    
```

- <ENTER>

After the initialization has run automatically, the instrument quits the diagnostic menu and runs a power on reset.

¹⁾ Care: This initialization alternative is reserved for the service engineer!

5.5 Validation / GLP

The requirements of **GLP** (**G**ood **L**aboratory **P**ractice) include a periodic check of analytical measuring instruments with regard to their reproducibility and accuracy using **S**tandard **O**perating **P**rocedures, **SOP**.

The 732 IC Detector and 733 IC Separation Center as part of the complete ion chromatography system, whose most important components also include separating column, pump and evaluation system, must be incorporated in its comprehensive validation.

Testing of the electronic and mechanical function groups of Metrohm instruments can and should be performed as part of a regular service by trained personnel of the manufacturing company (see *section 5.2.1*). All Metrohm instruments are equipped with start-up-test routines which check for perfect functioning of the relevant assemblies when the instrument is switched on. If no error message is displayed, it may be assumed the instrument is operating without faults.

The Metrohm company also supplies its instruments with an integrated diagnostic program (see *section 5.4*) which, in the case of possible malfunctions or faulty behavior, allows the user to check the functioning of certain assemblies and localize the fault. Diagnostic programs can also be integrated in a validation procedure.

6 Interfaces

6.1 RS232 interfaces

The 732 IC Detector has two RS232 interfaces **15** and **16** (see *Fig. 3*) which have the same features (see *section 6.1.7 – 6.1.9*), but can be used differently.

RS232 interface **16** is reserved for the connection of a PC or a printer (see *section 2.9.4 – 2.9.5*) and allows remote operation of the 732 IC Detector 732 and 733 Separation Center 733 by an attached PC (see *section 6.1.1 – 6.1.6*).

RS232 interface **15** is normally used for the attachment of a 709 IC Pump, but a printer can also be connected in place of the pump (see *section 2.9.6*).

6.1.1 General rules for remote control

The 732 IC Detector is equipped with the comprehensive Metrohm remote control language, which allows full control over the instrument via an RS232 interface, i.e. the IC Detector can receive data from an external device or send data to an external device. The 732 IC Detector sends $2 \times C_R$ and L_F as the terminator of a requested **data block**. In contrast, C_R and L_F are used as the terminator of a **data line**. On receipt of data from an external device, this must always close its commands with C_R and L_F . If more than one command is sent on a line, ';' must be used as a delimiter between the individual commands.

The data are grouped logically and readily understandable. For example, to select the dialog language the command

```
&Config.Aux.Dialog"english"
```

must be sent with entry of the boldface characters sufficing, in other words

```
&C. A. D"english"
```

All quantities of the 732 IC Detectors are collected in **groups**. The entries for the configuration, for example, are located in the group

```
&Config
```

The 'Config' group contains sub-groups, e.g. for setting the RS232 interface parameters

```
&Config.RSset
```

or for various settings

```
&Config.Aux
```

The data have a hierarchical structure (tree structure). The quantities which appear in this tree are called **objects** in what follows. The dialog language is that object which is called up with the command

&Config.Aux.Dialog

Once you are at the desired location in the tree, you can request the value of the object:

&Config.Aux.Dialog \$Q Q for query


The inquiry '\$Q' initiates the output of the value on the instrument, in other words the value output is triggered. Entries which begin with the character '\$' always trigger something and are hence referred to as **triggers** in what follows.

However, values of objects can not only be requested, they can also be modified. Values are always inputted in inverted commas, e.g.

&Config.Aux.Dialog"english"

6.1.2 Call-up of objects

A section from the object tree is shown below:



The following **rules** apply to the call-up of objects:

Rules	Example
The root of the tree is designated by '&'.	
For the call-up of an object, the nodes (levels) of the tree are marked by a point (.).	
For the call-up of the objects, as many characters as necessary to allow unambiguous allocation of the object suffice. If the call-up is not unambiguous, the first object in the series is identified.	Call up of the dialog language: &Config.Aux.Dialog or &C. A. D
Uppercase and lowercase letters can be used.	&C. A. D or &c. a. d
A value can be assigned to an object. Values are marked at both their beginning and end by inverted commas ("). They can contain maximum 24 ASCII characters. In the case of parameters with specified text expressions (e.g. on, off), only the English expressions may be used. Numeric values can contain up to 6 digits, a negative sign and a decimal point. Numbers with more than 6 digits are not accepted; more than 4 decimal places are rounded off. With numbers <1, leading zeros must be entered.	Entry of the dialog language: &C. A. D"english" Correct numeric entries: "0.1" Incorrect numeric entries: "1,5" or "+3" or ".1"

<i>Rules</i>	<i>Example</i>
Until a new object is called up, the old object remains in force.	Entry of a different dialog language: " deutsch "
New objects can be addressed relative to the old object: A leading point leads one node forwards in the tree. More than one leading point leads one node backwards in the tree. n nodes backwards require n+1 leading points.	From the root to the node 'Aux': &C. A Forwards from the node 'Aux' to 'Prog': .P Jump from node 'Prog' to the node 'Aux' and selection of the new object 'Dialog' at this node: ..D
To return to the root, enter a leading '&'.	Switch from node 'Dialog' via the root to node 'Parameter': &P

6.1.3 Triggers

Triggers initiate an action at the 732 IC Detector, e.g. sequence start or data transmission. Triggers are marked by the introducer '\$'.

The following triggers are possible:

\$G	Go	Starts processes, e.g. initiation of the auto-zero function or program start
\$S	Stop	Stops processes
\$Q	Query	Used to request all information from the current node in the tree forwards up to and including the values
\$Q.P	Path	Used to request the path from the root of the tree up to the current node
\$Q.H	Highest Index	Used to request the number of daughter nodes of the current node
\$Q.N"i"	Name	Used to request the name of the daughter node with index i, i = 1...n
\$D	Detail-Info	Used to request detailed status information
\$U	qUit	Used to terminate the data flow of the instrument, e.g. after \$Q

The triggers '\$G' and '\$S' are linked to particular objects, see Overview table in *section 6.1.6*.

All other triggers can always be used at all locations in the data tree.

Examples:

Request of the baud rate value: &Config.RSset.Baud \$Q

Request of all values of the node 'RSset': &Config.RSset \$Q

Request of the path of the node 'RSset': &Config.RSset \$Q.P

Triggering the auto-zero function: &Zero \$G

Request of the detailed status: \$D

6.1.4 Status messages

To ensure appropriate control of an external control unit is possible, it must also be possible to request status conditions which provide information on the status of the 732 IC Detector. The output of a status message is triggered by the trigger '\$D'. Status messages comprise the global status '\$R' and the detailed status messages listed below, which can also appear in combination, e.g. '\$R.Cond.ProgWait'. If a fault appears, the error message is also appended to the status message, e.g. '\$R.Cond.ZeroOK;E42'.

\$R	Ready: The 732 IC Detector is ready to measure.
.Cond	Basic status of the conductivity measurement
.Zero	Auto-zero function active, value not yet valid
.ZeroOK	Auto-zero function active, value valid
.Overflow	The measured value is 50...80% outside the full-scale range
.Background	The measured value is more than 80% outside the full-scale range
.ProgWait	Program active (can be started)
.ProgGo	Program running
.ProgEnd	Program ended
.ProgTest	Program test mode active
.Input	Change in the signal of an input line
.Output	Change in the signal of an output line
.Diagnosis	Instrument in diagnostic status

6.1.5 Error messages

Error messages 'exxx' are appended to the status message and separated from this by a ';'.

<i>Error</i>	<i>Meaning</i>	<i>Exit/Corrective action</i>
E28	Wrong object call-up.	Correct path.
E29	Wrong value.	Enter correct value or new path.
E30	Wrong trigger.	Enter correct trigger or new path.
E36	RS232 receive error; parity.	<QUIT>, set same parity for both devices.
E37	RS232 receive error; stop bit.	<QUIT>, set same stop bit for both devices.
E38	RS232 receive error; overflow (at least 1 character could not be read).	<QUIT>, set same baud rate for both devices, restart sender.
E39	RS232 receive error; overflow of internal receive buffer (>82 characters).	<QUIT>.
E40	RS232 send error; DSR=OFF. Handshake not answered for longer than 1 s.	<QUIT>, check receiver (switched on and ready?).
E41	RS232 send error; DCD=ON. Handshake not answered for longer than 1 s.	<QUIT>, check receiver (switched on and ready?).
E42	RS232 send error; CTS=OFF. Handshake not answered for more than 1 s.	<QUIT>, check receiver (switched on and ready?).
E43	RS232 send error; transmission of the IC Detector was interrupted with XOFF for at least 3 s.	<QUIT> or send XON.
E44	RS232 send error ; the RS parameters are no longer the same for both devices.	<QUIT>, reset RS parameters for both devices.
E45	RS232 send error; the receive buffer of the IC Detector contains an incomplete string (L _F missing), the transmission is thus blocked.	<QUIT> or send L _F .
E50...E59	Error in IO test.	<QUIT>, inform Metrohm Service.
E60...E89	Error in RS232 test.	<QUIT>, check connection between the RS interfaces.
E200	Instrument adjustment invalid.	Inform Metrohm Service.
E202	Wrong entry via RS232 remote control.	Enter correct value, correct trigger or new path.
E240	Error in storage of the EEPROM blocks.	<QUIT>, re-enter all values.

6.1.6 Remote control commands

The remote control tree can be divided into the following main branches:

&	Root
C onfig	Instrument configuration
M ethods	Methods
P arameter	Parameters
P rogram	Program
S electEvent	Events
Z ero	On/off switching of auto-zero
M arker	Triggering marking signal
P lot	Start/stop graphics plot
P rint	Start/stop measured value printout
I nf	Instrument information
S etup	Settings of operating mode
A ssembly	Data of the assemblies
D iagnosis	Diagnostic tests

The following table lists all objects of the remote control tree. For the unambiguous designation of the objects, the boldface characters suffice. The meaning of the individual object is described here only in brief, for more detailed information please refer to *section 4*. The default values of the objects are printed in boldface.

<i>Object</i>	<i>Meaning</i>	<i>Entry range/Selection</i>
&C onfig	Instrument configuration	
D etector	Conductivity detector	
T hermostat	Operating temperature of conductivity cell	25, 30, 35 , 40, 45°C, off
Z eroUnit	Unit for display of auto-zero value	µS/cm , %fs, mV
C ellConst	Cell constant of conductivity cell	13.0... 16.7 ...21.0 /cm
P rinter	External printer	
I d1	1 st line of printout header	18 ASCII characters
I d2	2 nd line of printout header	18 ASCII characters
P rintHead	Output of printout header: once, always, never	once , always, off
D ateTime	Output of date and time in printout header	on , off
S endTo	Selection of printer driver	IBM , Epson, Seiko, Citizen, HP
P rintMeasVal	Measured value printout (function same as <PRINT> key)	\$G, \$S
P rintCrit	Criterion for measured value printout: immediately (single value), time (continuous output), off	immed. , time, off
T ime	Timed output of measured values	
I nterval	Time interval for output of measured values	0.4... 1.0 ...99999 s
S topTime	Stop time (off = infinite)	1...999 min, off
D ateTime	Output of date and time of measured value	on, off

Object	Meaning	Entry range/Selection
. Aux	General instrument settings	
. RunNo	Run number (automatically incremented by 1 after every determination)	0 ...999, off
. Cycles	Number of cycles in loop programs	1 ...999
. CEvent	Settings for events	
. CEenable	Perform events during running program	on, off
. Set	Setting of date and time	\$G, \$S
. Date	Date	YYYY-MM-DD
. Time	Time	HH:MM:SS
. Dialog	Dialog language	english , deutsch, francais, español
. DevLabel	Device label	8 ASCII characters
. Prog	Number of program version	read only
. Beeper	Beeper	
. Status	Status: on, error only, off	on , only error, off
. Repeat	Repeat time for error message	5... 60 ...999 s, off
. RSset	RS232 interface	
. Baud	Baud rate in bit/s	9600 , 4800, 2400, 1200, 600, 300
. DataBit	Data bits	7, 8
. StopBit	Stop bits	1 , 2
. Parity	Parity: none, odd, even	none , odd, even
. Handsh	Handshake: hardware simple, hardware full, software character, software line, none	HWs , HWf, SWchar, SWline, none
. RSset709	2nd RS232 interface "709 IC Pump"	
. Status	Connection in operation without 709 IC Pump	printer , off
. 709	Connected: 709 IC Pump	
. Baud	Baud rate in bit/s	9600 , 4800, 2400, 1200, 600, 300
. Printer	Connected: printer	
. Baud	Baud rate in bit/s	9600 , 4800, 2400, 1200, 600, 300
. DataBit	Data bits	7, 8
. StopBit	Stop bits	1 , 2
. Parity	Parity: none, odd, even	none , odd, even
. Handsh	Handshake: hardware simple, hardware full, software character, software line, none	HWs , HWf, SWchar, SWline, none
. SepCenter	733 IC Separation Center	
. Config	Display of the configuration	read only
. Control	Control: without restriction, via 732 only	no restriction , 732 only
. AValve	Injection valve A	
. StatusA	Display of valve setting: inject, fill	read only
. TriggerA	Switching the valve	inject, fill
. BValve	Injection valve B	
. StatusB	Display of the valve setting: inject, fill	read only
. TriggerB	Switching the valve	inject, fill
. Suppressor	Suppressor module	
. Autostep	Automatic triggering of "step"	fill , inject, off
. Status	Display of the suppressor setting: in position or undefined	read only
. Trigger	Switching the suppressor	---, step

Object	Meaning	Entry range/Selection
&Methods	Methods	
. FreeMemory	Free memory (1...9)	read only
. Recall	Load method	\$G
. Name	Method name	8 ASCII characters
. Store	Store method	\$G
. Name	Method name	8 ASCII characters
. Delete	Delete method	\$G
. Name	Method name	8 ASCII characters
. AllDelete	Delete all methods	\$G
&Parameter	Parameters	
. Detector	Conductivity detector	
. Range	Measurement range	100, 200, 500 $\mu\text{S/cm}$, 1 , 2, 5, 10 mS/cm
. FullScale	Full-scale range (operating range)	0.05 $\mu\text{S/cm}$... 1 ...10 mS/cm
. TempCoeff	Temperature coefficient	1.5, 2.5 $\%^\circ\text{C}$
. AnalogOut	Analog output	
. Polarity	Polarity	+ , -
. Offset	Zero point shift in % of full-scale range	0 , 10, 50 %fs
. Damping	Damping	on, off
. PlotPara	Graphics plot	
. Autostart	Automatic start of graphics plot	on, off
. Interval	Time interval for graphics plot	0.4... 1.0 ...99999 s
. TimeScale	Time scale (paper speed)	0.1, 0.2, 0.5, 1, 2, 5, 10 , 20, 60, 120 mm/min
. TLabel	Labeling of time axis: relative, absolute	rel , abs
. StopTime	Stop time for graphics plot	1...999 min, off
. Left	Left boundary for conductivity axis	-2000... 0 ...2000 $\mu\text{S/cm}$, -2000... 0 ...2000 mS/cm , -2... 0 ...2 S/cm
. Right	Right boundary for conductivity axis	-2000... 0 ... 10 ...2000 $\mu\text{S/cm}$, -2000... 0 ...2000 mS/cm , -2... 0 ...2 S/cm
. Pump709	709 IC Pump	
. Flow	Flow rate	0.01... 0.50 ...5.00 mL/min
. Pmax	Maximum shutoff pressure	0.1... 10.0 ...50 MPa
. Pmin	Minimum shutoff pressure	0.0 ...50 MPa
. Flowcorr	Correction factor for flow rate	0.90... 1.00 ...1.10
&Program	Program	\$G, \$S

Object	Meaning	Entry range/Selection
. Parameter	General program parameters	
. Type	Program type: cycle, remote, inject	cycle , remote, inject
. Cycle	Loop program	
. No	Number of cycles	1...999
. Status	Program status: inactive, active, test	disabled , enabled, test
. Edit	Edit program	
. Select	Selection of program steps (max. 20)	
. 1	Program step 1	
. Time	Time for program step 1	x.x , 0.0...999.9 min
. Action	Actions of program step 1	
. Flag	Program flag	--- , return, reset, end
. FSChange	Change full-scale range	--- , on
. Range	Reset measurement range	100 , 200, 500 μ S/cm, 1, 2, 5, 10 mS/cm
. FullScale	Reset full-scale range	0.05... 100 ...999 μ S/cm, 1...10 mS/cm
. Zero	On/off switching of auto-zero	--- , on, off
. Polarity	Change polarity	--- , +, -
. Mark	Trigger marking signal	--- , on
. ValveA	Switch injection valve A	--- , fill, inject
. ValveB	Switch injection valve B	--- , fill, inject
. Suppressor	Switch suppressor module	--- , step
. Remote	Set remote output lines 1...8: * (leave unchanged), 0 (off, inactive, open), 1 (on, active, 0 V)	* , 0, 1
. PumpRS	On/off switching of pump drive	--- , on, off
. Flow	Change flow rate	--- , 0.01...5.00 mL/min
. Pmax	Change maximum shutoff pressure	--- , 0.1...50 MPa
. Pmin	Change minimum shutoff pressure	--- , 0.0...50 MPa
. 20	Program step 20	
. Time	Time for program step 20	x.x , 0.0...999.9 min
. Action	Actions of program step 20	
:	see Program step 1	
. DeleteAll	Delete all program steps possible only in program status "disabled"	\$G
&SselectEvent	Events	
. 1	Event 1	
. Date	Date/time or only time (daily) for event 1	YY-MM-DD HH:MM:SS
. Action	Actions of event 1 1	
. FSChange	Change full-scale range	--- , on
. Range	Reset measurement range	100 , 200, 500 μ S/cm, 1, 2, 5, 10 mS/cm
. FullScale	Reset full-scale range	0.05... 100 ...999 μ S/cm, 1...10 mS/cm
. Zero	On/off switching of auto-zero	--- , on, off
. Polarity	Change polarity	--- , +, -

Object	Meaning	Entry range/Selection
<ul style="list-style-type: none"> └ . Mark └ . AValve └ . BValve └ . Suppressor └ . Remote └ . ProgRS └ . PumpRS └ . Flow └ . Pmax └ . Pmin 	Trigger marking signal Switch injection valve A Switch injection valve B Advance suppressor module Set remote output lines 1...8: * (leave unchanged), 0 (off, inactive, open), 1 (on, active, 0 V) Start/stop program On/off switching of pump drive Change flow rate Change maximum shutoff pressure Change minimum shutoff pressure	---, on ---, fill, inject ---, fill, inject ---, step *, 0, 1 ---, on, off ---, on, off ---, 0.01...5.00 mL/min ---, 0.1...50 MPa ---, 0.0...50 MPa
. 4 <ul style="list-style-type: none"> └ . Date └ . Action └ : 	Event 4 Date/time or only time (daily) for event 4 Actions of event 4 <i>see Event 1</i>	YY-MM-DD HH:MM:SS
&Zero <ul style="list-style-type: none"> └ . RefValue 	Start/stop of the auto-zero function Auto-zero ref. value (compensated conductivity)	\$G, \$S read only
&Marker	Triggering a marking signal	\$G, \$S
&Plot <ul style="list-style-type: none"> └ . State 	Start/stop of the graphics plot Status of the graphics plot: on, off	\$G, \$S read only
&Print <ul style="list-style-type: none"> └ . State 	Start/stop of the measured value printout Status of the measured value printout: on, off	\$G, \$S read only
&Info <ul style="list-style-type: none"> └ . Report <ul style="list-style-type: none"> └ . Select └ . ActualInfo <ul style="list-style-type: none"> └ . Inputs <ul style="list-style-type: none"> └ . State 	Instrument information Send formatted reports Selection of the report Current information Remote inputs Status of the input lines of the two remote inter- faces "733 IC Separation Center" and "Remote" in byte form (1=on, low, active; 0=off, high, inactive): "733 IC Sep. Center" "Remote" n 0 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 Pin 21 9 22 10 23 11 24 12 21 9 22 10 23 11 24 12 $Status = \sum_{n=0}^{14} 2^n - 2^{15}$ <u>Example:</u> 1 0 1 0 0 0 1 0 0 0 0 0 1 0 0 1 Status = $2^0 + 2^2 + 2^6 + 2^{12} - 2^{15} =$ = 1 + 4 + 64 + 4096 - 32768 = -28603	\$G all , config, event, method, method list, param, program, pump read only

Object	Meaning	Entry range/Selection																																																																			
. C hange	Change in the status of the input lines since the last "clear", in byte form (1 = change, 0 = no change, see above)	read only																																																																			
	. C lear	Reset "change" byte to 0	\$G																																																																		
. O utputs	Remote outputs																																																																				
. S tate	Status of the output lines of the two remote interfaces "733 IC Separation Center" and "Remote" in byte form: 1=on, low, active; 0=off, high, inactive	read only																																																																			
	<table style="margin-left: auto; margin-right: auto;"> <tr> <td></td> <td colspan="7" style="text-align: center;">"733 IC Sep. Center "</td> <td colspan="8" style="text-align: center;">"Remote"</td> </tr> <tr> <td>n</td> <td>0</td><td>1</td><td>2</td><td>3</td><td>4</td><td>5</td><td>6</td><td>7</td> <td>8</td><td>9</td><td>10</td><td>11</td><td>12</td><td>13</td><td>14</td><td>15</td> </tr> <tr> <td>Pin</td> <td>5</td><td>18</td><td>4</td><td>17</td><td>3</td><td>16</td><td>1</td><td>2</td> <td>5</td><td>18</td><td>4</td><td>17</td><td>3</td><td>16</td><td>1</td><td>2</td> </tr> <tr> <td>Output</td> <td></td><td></td><td></td><td></td><td></td><td></td><td></td><td>8</td> <td>7</td><td>6</td><td></td><td>1</td><td>2</td><td></td><td>3</td><td>4</td><td>5</td> </tr> </table> $Status = \sum_{n=0}^{14} 2^n - 2^{15}$ <p><u>Example:</u> 0 0 1 0 0 0 0 0 1 0 0 0 1 0 0 1 Status = $2^2 + 2^8 + 2^{12} - 2^{15} = 4 + 256 + 4096 - 32768 = -28412$</p>		"733 IC Sep. Center "							"Remote"								n	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	Pin	5	18	4	17	3	16	1	2	5	18	4	17	3	16	1	2	Output								8	7	6		1	2		3	4	5
	"733 IC Sep. Center "							"Remote"																																																													
n	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15																																																					
Pin	5	18	4	17	3	16	1	2	5	18	4	17	3	16	1	2																																																					
Output								8	7	6		1	2		3	4	5																																																				
. C hange	Change in the status of the output lines since the last "clear", in byte form (1 = change, 0 = no change, see above)	read only																																																																			
	. C lear	Reset "change" byte to 0	\$G																																																																		
. M easValue	Measured values																																																																				
. C onductivity	Current measured value in S/cm	read only																																																																			
	. B ackground	Absolute conductivity in S/cm	read only																																																																		
	. R unTime	Run time in min	read only																																																																		
. P ump709	709 IC Pump																																																																				
. S tate	Status: ready, running, stopped, unknown	read only																																																																			
	. P actual	Current pressure in MPa	read only																																																																		
	. E rror	Error status: no error, stopped: YY-MM-DD XX:MM (stopped with date)	read only																																																																		
. M ethod	Method in working memory																																																																				
. N ame	Method name	read only																																																																			
	. S tatus	Status: original, modified	read only																																																																		
	. I ndex	Method index (1...10)	read only																																																																		
. D isplay	Display																																																																				
. S elect	Status display: full scale, absolute, prog.type, prog.step, method, date	read only																																																																			
	. L 1	Display in LCD line 1	read only																																																																		
	. L 2	Display in LCD line 2	read only																																																																		
	. L ED	Status of the LED displays in byte form: 1=on, 0=off																																																																			
	<table style="margin-left: auto; margin-right: auto;"> <tr> <td>n</td> <td>0</td><td>1</td><td>2</td><td>3</td> </tr> <tr> <td>LED</td> <td>OVERLOAD</td><td>THERMOSTAT</td><td>PROG R/S</td><td>ZERO</td> </tr> </table> $Status = \sum_{n=0}^3 2^n$ <p><u>Example:</u> 0 1 1 1 Status = $2^1 + 2^2 + 2^3 = 2 + 4 + 8 = 14$</p>	n	0	1	2	3	LED	OVERLOAD	THERMOSTAT	PROG R/S	ZERO																																																										
n	0	1	2	3																																																																	
LED	OVERLOAD	THERMOSTAT	PROG R/S	ZERO																																																																	

Object	Meaning	Entry range/Selection
&Setup	Operating mode	
. IdReport	Send identification before report Elements of the message: «Space (dec 32), ', report identification» " 'co" config " 'ev" event " 'me" method " 'ml" method list " 'cp" parameter " 'pr" program " 'pu" pump " 'mp" measured value output with <PRINT>	on, off
. Keycode	Send key code of pressed keys Elements of the message: «Space (dec 32), # or ù, two-place code» 0 9 (METHOD) 12 CLEAR 1 8 (PUMP R/S) 13 SELECT 2 7 (CONFIG) 14 3 (PRINT) 3 ZERO 15 2 (REPORT) 4 PARAM 16 1 (PLOT) 6 6 (PROGRAM) 18 ENTER 7 5 (EVENT) 19 QUIT 8 4 (ZERO OFF) 20 -/+ (→) 9 PROG R/S 21 . (←) 10 FULL SCALE 22 0 (MARK)	on, off
. Tree	Definition of the response to \$Q	
. Short	Path names are sent with just the required number of characters (boldface characters)	on, off
. ChangedOnly	Only path names and their values which have been edited once are sent.	on, off
. Trace	Send path and value on changes Element of the change message: «Space (dec 32), path, "value"»	on, off
. Lock	Lock functions	
. Keyboard	Lock all keys	on, off
. Config	Locking of <CONFIG> key	on, off
. Event	Locking of <EVENT> key	on, off
. FullScale	Locking of <FULL SCALE> key	on, off
. Mark	Locking of <MARK> key	on, off
. Method	Locking of <METHOD> key	on, off
. Param	Locking of <PARAM> key	on, off
. Plot	Locking of <PLOT> key	on, off
. Print	Locking of <PRINT> key	on, off
. ProgRS	Locking of <PROG R/S> key	on, off
. Program	Locking of <PROGRAM> key	on, off
. PumpRS	Locking of <PUMP R/S> key	on, off
. Report	Locking of <REPORT> key	on, off
. Select	Locking of <SELECT> key	on, off
. Zero	Locking of <ZERO> key	on, off
. Zerooff	Locking of <ZERO OFF> key	on, off

Object	Meaning	Entry range/Selection
<ul style="list-style-type: none"> . Methods <ul style="list-style-type: none"> . Recall . Store . Delete . Display . Remote 	<p>Locking of functions of the <METHOD> key</p> <p>Locking of the "load" functions</p> <p>Locking of the "store" functions</p> <p>Locking of the "delete" functions</p> <p>Locking of the LCD (measured values are no longer displayed)</p> <p>Locking of all remote lines</p>	<p>on, off</p> <p>on, off</p> <p>on, off</p> <p>on, off</p> <p>on, off</p>
<ul style="list-style-type: none"> . AutoInfo 	<p>Automatic messages on change</p> <p><i>Note: While the settings of the AutoInfo branch are stored after saving with "&Setup.Save", they are not initialised when the instrument is switched on again. They do not become active until one of the settings has been re-set.</i></p>	
<ul style="list-style-type: none"> . Message <ul style="list-style-type: none"> . DateTime . Error . Ready . Stopped . Wait . PowerOn . Inputs . Outputs 	<p>Elements of the automatic message: «Space (dec 32), !, instrument designation, message from sub-nodes, date (optional), time (optional)» <u>Example:</u> « !732B".R;.0"95-10-25 16:30:47»</p> <p>Output of date and time of appearance</p> <p>Message "E" on error</p> <p>Message "R" on attainment of basic status and on program end</p> <p>Message "S" on stop of a sequence</p> <p>Message "W" on entry into wait status</p> <p>Message "P" when instrument switched on</p> <p>Message "I" on change in an input line</p> <p>Message "O" on change in an output line</p>	<p>on, off</p> <p>on, off</p> <p>on, off</p> <p>on, off</p> <p>on, off</p> <p>on, off</p> <p>on, off</p> <p>on, off</p>
<ul style="list-style-type: none"> . Save 	<p>Saving of all parameters</p> <p>(without saving the parameters changed by remote control will be lost when the instrument is switched off)</p>	\$G
<ul style="list-style-type: none"> . InstrNo <ul style="list-style-type: none"> . Value 	<p>Instrument number</p> <p>Serial number</p>	8 ASCII characters
<ul style="list-style-type: none"> . Interface <ul style="list-style-type: none"> . OutputAssign <ul style="list-style-type: none"> . RemOut . InputAssign <ul style="list-style-type: none"> . ProgRS . PumpRS . FillA . InjectA . FillB . InjectB . Report . Zero 	<p>Interfaces</p> <p>Setting of the remote output lines</p> <p>Basic settings for remote output lines 1...8 when instrument switched on: 1=on, low, active; 0=off, high, inactive</p> <p>Assignment of the remote input lines (details, see section 6.2.1)</p> <p>Start/stop program (= <PROG R/S>)</p> <p>Start/stop 709 IC Pump (= <PUMP R/S>)</p> <p>Switch injection valve A to "FILL"</p> <p>Switch injection valve A to "INJECT"</p> <p>Switch injection valve B to "FILL"</p> <p>Switch injection valve B to "INJECT"</p> <p>Trigger report (= <REPORT> key)</p> <p>Switch on auto-zero (= <ZERO> key)</p>	<p>0000000...11111111</p> <p>1...15</p> <p>1...2...15</p> <p>1...4...15</p> <p>1...8...15</p> <p>1...3...15</p> <p>1...12...15</p> <p>1...5...15</p> <p>1...6...15</p>

Object	Meaning	Entry range/Selection
. Mark	Trigger marking signal (= <MARK> key)	1... 7 ...15
. Zerooff	Switch off auto-zero (= <ZERO OFF> key)	1... 9 ...15
. Plot	Trigger graphics plot (= <PLOT> key)	1... 10 ...15
. Polarity	Switch polarity at analog output	1... 11 ...15
. Select	Select (= <SELECT> key)	1... 13 ...15
. Quit	Quit (= <QUIT> key)	1... 14 ...15
. Enter	Enter (= <ENTER> key)	1... 15 ...15
. Graphics	General graphics parameters	
. Grid	Grid lines on graphics output	on, off
. Frame	Frame on graphics output	on, off
. Recorder	Setting of printout	
. Width	Relative width of printout	0.4... 0.8 ...1.0
. Peripherals	Peripheral instruments	
. With733	Operation with 733 IC Separation Center	on , off
. With709	Operation with 709 IC Pump	on , off

&Assembly

Basic elements of the assembly

. Meas	Measurement	
. State	On/off switching of measurement	on , off
. Outputs	Remote lines	
. SmpIX	Activate (on) or deactivate (off) remote output lines	on , off
. AutoEOD	Automatic output of EOD on program end and "return" with loop programs	on , off
. SetLines	Set remote outputs 1...16:	\$G
	Interface "709 IC Pump": outputs 1...8	
	Interface "Remote": outputs 9...16	
	active setting of a static signal	
	inactive resetting of the static signal	
	pulse output a pulse (length ca. 150 ms)	
	off do not set line	
. L1	Signal of output line 1 (pin 5)	active, inactive, pulse, off
. L2	Signal of output line 2 (pin 18)	active, inactive, pulse, off
. L3	Signal of output line 3 (pin 4)	active, inactive, pulse, off
. L4	Signal of output line 4 (pin 17)	active, inactive, pulse, off
. L5	Signal of output line 5 (pin 3)	active, inactive, pulse, off
. L6	Signal of output line 6 (pin 16)	active, inactive, pulse, off
. L7	Signal of output line 7 (pin 1)	active, inactive, pulse, off
. L8	Signal of output line 8 (pin 2)	active, inactive, pulse, off
. L9	Signal of output line 9 (pin 5)	active, inactive, pulse, off
. L10	Signal of output line 10 (pin 18)	active, inactive, pulse, off
. L11	Signal of output line 11 (pin 4)	active, inactive, pulse, off
. L12	Signal of output line 12 (pin 17)	active, inactive, pulse, off
. L13	Signal of output line 13 (pin 3)	active, inactive, pulse, off
. L14	Signal of output line 14 (pin 16)	active, inactive, pulse, off
. L15	Signal of output line 15 (pin 1)	active, inactive, pulse, off
. L16	Signal of output line 16 (pin 2)	active, inactive, pulse, off
. ResetLines	Set remote outputs 1...16 to "off"	\$G

Object	Meaning	Entry range/Selection
&D Diagnosis	Diagnostic tests	
<ul style="list-style-type: none"> . EEPROMInit <ul style="list-style-type: none"> . BlockSelect 	<p>Initialize data memory</p> <p>Selection of blocks to be initialized:</p> <p>USER All blocks except instrument adjustment and instrument number</p> <p>CONFIG Block &Config</p> <p>INFO Block &Info</p> <p>SETUP Block &Setup</p> <p>ASSEMBLY Block &Assembly</p> <p>EVENT Block &SelectEvent</p> <p>METHOD Blocks &Parameter and &Program</p> <p>ALL All blocks incl. instrument adjustment (security code necessary)</p>	<p>\$G</p> <p>USER, CONFIG, INFO, SETUP, ASSEMBLY, EVENT, METHOD, ALL</p>
<ul style="list-style-type: none"> . RunMode <ul style="list-style-type: none"> . Report . Stop 	<p>Parameters for sequence control</p> <p>Output of a diagnostic report on printer</p> <p>Stop after error</p>	<p>on, off</p> <p>on, off</p>
. R AMTest	Check working memory	\$G
. K eyTest	Check keypad	\$G, \$\$
. D isplayTest	Check display	\$G, \$\$
. RS Test	Check RS232 interfaces (possible only with 3.496.8480 Test plug, see section 5.4.6)	
<ul style="list-style-type: none"> . RS232 . IC Pump 	<p>Test of "RS232" interface</p> <p>Test of "709 IC Pump" interface</p>	<p>\$G, \$\$</p> <p>\$G, \$\$</p>
. I O Test	Check remote interfaces (possible only with 3.496.8510 Test plug, see section 5.4.7)	
<ul style="list-style-type: none"> . Remote . SepCenter 	<p>Test of "Remote" interface</p> <p>Test "733 IC Separation Center" interface</p>	<p>\$G, \$\$</p> <p>\$G, \$\$</p>
. I ntTest	Internal hardware test	
<ul style="list-style-type: none"> . UcellTest . RangeTest . SenseTest . DACCoarse . DACFine 	<p>Check voltage generator</p> <p>Check range amplifier</p> <p>Check sense amplifier</p> <p>Check DA converter for coarse adjustment</p> <p>Check DA converter for fine adjustment</p>	<p>\$G, \$\$</p> <p>\$G, \$\$</p> <p>\$G, \$\$</p> <p>\$G, \$\$</p> <p>\$G, \$\$</p>
. E xtTest	External hardware test	
<ul style="list-style-type: none"> . LowpassTest . OffsetTest . PolarTest . TCoeffTest . FreqTest . ThermostatTest 	<p>Check low-pass function</p> <p>Check offset</p> <p>Check polarity</p> <p>Check temperature coefficient</p> <p>Check frequency</p> <p>Check thermostat</p>	<p>\$G, \$\$</p> <p>\$G, \$\$</p> <p>\$G, \$\$</p> <p>\$G, \$\$</p> <p>\$G, \$\$</p> <p>\$G, \$\$</p>

<i>Object</i>	<i>Meaning</i>	<i>Entry range/Selection</i>
. SimulateKey	Simulate keystroke	0...22
	0 9 (METHOD)	12 CLEAR
	1 8 (PUMP R/S)	13 SELECT
	2 7 (CONFIG)	14 3 (PRINT)
	3 ZERO	15 2 (REPORT)
	4 PARAM	16 1 (PLOT)
	6 6 (PROGRAM)	18 ENTER
	7 5 (EVENT)	19 QUIT
	8 4 (ZERO OFF)	20 -/+ (→)
	9 PROG R/S	21 . (←)
	10 FULL SCALE	22 0 (MARK)
. Adjust	Instrument adjustment (security code necessary)	\$G
. InstrNo	Serial number (accessible by remote control only under &Setup.InstrNo)	
. PowerOn	Simulation "Power on "	\$G


6.1.7 Data transmission protocol

The two RS232 interfaces of the 732 IC Detectors are configured as DTE (Data Terminal Equipment) with the following technical specifications:

- *Standard* Data interface in accordance with EIA standard RS 232C (DIN 66020, page 1), transmission parameters can be set under ">CONFIG/RS settings", see section 4.4.2.
- *Control characters*

C_R	DEC 13	HEX 0D
L_F	DEC 10	HEX 0A
XON	DEC 17	HEX 11
XOFF	DEC 19	HEX 13
- *Max. line length* 80 characters + $C_R L_F$
- *Mode* Full duplex (simultaneous transmission and receiving)

Restriction:
If data are first received by the interface, transmission is not started until the receipt is at an end.
- *Cable* For interconnection of the 732 IC Detector with third-party devices, only a shielded data cable may be used (e.g. Metrohm D.104.0201). The cable shield must be perfectly earthed at both devices (pay attention to current loops: always use star earthing). Use only connectors with sufficient shielding (e.g. Metrohm K.210.0001 with K.210.9004).
- *Cable length* max. ca. 15 m




6.1.8 Handshake

No handshake (none)


The 732 IC Detector neither checks handshake inputs (CTS, DSR, DCD) nor sets handshake outputs (DTR, RTS).

Reduced hardware handshake (HWs)

732 IC Detector as **receiver**:



IC Detector 732 as **sender**:




The data flow can be interrupted by deactivating the CTS line.


Full hardware handshake (HWf)

All handshake inputs are checked, all handshake outputs are set.

732 IC Detector as **receiver**:



732 IC Detector as **sender**:



The data flow can be interrupted by deactivation of the CTS line.

Software handshake with character stop (SWchar)


Handshake inputs (CTS, DSR, DCD) are not checked at the 732 IC Detector, handshake outputs (DTR, RTS) are set.

As soon as an L_F is recognized, the 732 IC Detector sends XOFF. After this time, it can still receive 6 characters and store them temporarily.


However, the 732 IC Detector also sends XOFF when its input buffer contains 60 characters. After this time, it can still receive maximum 22 characters (incl. L_F).

If the transmission is interrupted for the time of 4 characters after the 732 IC Detector has sent XOFF, the string previously received will be processed even if no L_F has been sent.

732 IC Detector as **receiver**:



732 IC Detector as **sender**:




Software handshake with line stop (SWline)


Handshake inputs (CTS, DSR, DCD) are not checked at the 732 IC Detector, handshake outputs (DTR, RTS) are set.

The 732 IC Detector has an input buffer which can accept a string of up to 80 characters + $C_R L_F$. As soon as an L_F is recognized, the 732 IC Detector sends XOFF. After this time, it can still receive max. 6 characters and store them temporarily. The string sent previously is now processed by the 732 IC Detector. It then sends XON and is again ready to receive.

732 IC Detector as **receiver**:



732 IC Detector as **sender**:



The transmission of the 732 IC Detector can be stopped by the external device with XOFF. After the receipt of XOFF, the IC Detector finishes transmission of the line already started. If the data output is blocked for more than 3 s by XOFF, "RS error 43" appears in the display.

6.1.9 Pin assignment

RS232C interface	external
<p>Transmitted Data (TxD) If no data transmission takes place, the line is maintained in the "ON" condition. Data are sent only when CTS and DSR are in the "ON" condition and DCD is in the "OFF" condition.</p> <p>Received Data (RxD) Data are received only when DCD is "ON".</p> <p>Request to Send (RTS) ON condition: 732 IC Detector is ready to send data.</p> <p>Clear to Send (CTS) ON condition: Remote station is ready to receive data.</p> <p>Data Set Ready (DSR) ON condition: The transmission line is connected.</p> <p>Signal Ground (GND)</p> <p>Data Carrier Detect (DCD) ON condition: The received signal level is within the tolerance range (remote station is ready to send data).</p> <p>Data Terminal Ready (DTR) ON condition: 732 IC Detector is ready to receive data.</p>	
<p>Protective earth Direct connection from cable connector to protective earth of the device.</p> <p>Polarity assignment of the signals</p> <ul style="list-style-type: none"> Data lines (TxD, RxD) <ul style="list-style-type: none"> Voltage negative (<math>< -3\text{ V}</math>): signal status "ONE" Voltage positive (> +3 V): signal status "ZERO" Control or message lines (CTS, DSR, DCD, RTS, DTR) <ul style="list-style-type: none"> Voltage negative (<math>< -3\text{ V}</math>): OFF status Voltage positive (> +3 V): ON status <p>In the transition region from +3 V to -3 V, the signal status is undefined.</p> <p>Driver 14C88 to EIA RS 232C specification</p> <p>Receiver 14C89 to EIA RS 232C specification</p>	<p>Contact arrangement at the sockets "RS 232" and "709 IC Pump" (male)</p> <p>Connection cables to external devices must have an appropriate 25-pin connector (female). Metrohm offers connector sockets (K.210.9004) and the associated housing (K.210.0001) as an option for the preparation of such cables.</p>
<p><i>No liability whatsoever will be accepted for damage arising from the improper connection of devices.</i></p>	

6.1.10 RS232 error rectification









<i>Problem</i>	<i>Questions for corrective action</i>
<p>No characters can be received on an attached printer.</p>	<p>⇒ Are the devices switched on and the connecting cables plugged in properly?</p> <p>⇒ Is the printer set to "on-line"?</p> <p>⇒ Do the baud rate, data bits and parity for the two devices have the same setting?</p> <p>⇒ Is the handshake set properly?</p> <p>If everything appears ok, try to print out a report by pressing the <REPORT> key.</p>
<p>No data transmission occurs and an error message appears in the display of the 732 IC Detector.</p>	<p>ⓘ RS error 36¼439: Receive errors. Do the RS232 data transmission parameters of the two devices have the same setting?</p> <p>⇒ RS error 40¼42: transmission errors. Is the cable used correctly wired and plugged in? Is the printer switched on and set to "on-line"?</p> <p>⇒ RS error 43: Data output of the 732 IC Detector blocked for longer than 3 s by XOFF.</p>
<p>The received characters are garbled.</p>	<p>⇒ Do the data bits and parity of the two devices have the same setting?</p> <p>⇒ Does the baud rate of the two devices have the same setting?</p> <p>⇒ Is the correct printer selected?</p> <p>⇒ The data transfer has been interrupted by the hardware during a printout. Set up connections again, switch printer off then on.</p>
<p>The graphics plot is not printed out properly. Other reports are ok.</p>	<p>A handshake is necessary for the graphics plot.</p> <p>⇒ Is your cable correctly wired? (The DTR of the printer must be wired with the CTS of the 732 IC Detectors.)</p> <p>⇒ Set the handshake on the 732 IC Detector to "HwS". The printer must be so configured that its DTR is set (usually with DIP switches).</p>

6.2 Remote interfaces

The 732 IC Detector has the two remote interfaces **14** "733 IC Separation Center" and **17** "Remote" (see *Fig. 3*) which each have 8 input and 8 output lines.

6.2.1 "Remote" interface

The remote interface **17** "Remote" is used for the connection of any type of external device (see *section 2.9.8*) and has the following pin assignment:

732 IC Detector	Pin	Function																																																																																																									
Inputs																																																																																																											
	21	Print Triggering the print command Equivalent to pressing the <PRINT> key (see <i>section 4.8.1</i>). 																																																																																																									
	9																																																																																																										
	22	2¹ Triggering of functions (remote control) Assignment of the input lines see <i>section 4.4.1</i> <table border="1" data-bbox="746 1064 1396 1608"> <thead> <tr> <th>Input</th> <th>Decimal</th> <th>Function</th> </tr> <tr> <th>23</th> <th>10</th> <th>22</th> <th>9</th> <th>code</th> <th>(default settings)</th> </tr> </thead> <tbody> <tr><td>H</td><td>H</td><td>H</td><td>H</td><td>0</td><td>Inactive</td></tr> <tr><td>H</td><td>H</td><td>H</td><td>L</td><td>1</td><td>PROG R/S</td></tr> <tr><td>H</td><td>H</td><td>L</td><td>H</td><td>2</td><td>PUMP R/S</td></tr> <tr><td>H</td><td>H</td><td>L</td><td>L</td><td>3</td><td>FILL B/STEP</td></tr> <tr><td>H</td><td>L</td><td>H</td><td>H</td><td>4</td><td>FILL A</td></tr> <tr><td>H</td><td>L</td><td>H</td><td>L</td><td>5</td><td>REPORT</td></tr> <tr><td>H</td><td>L</td><td>L</td><td>H</td><td>6</td><td>ZERO</td></tr> <tr><td>H</td><td>L</td><td>L</td><td>L</td><td>7</td><td>MARK</td></tr> <tr><td>L</td><td>H</td><td>H</td><td>H</td><td>8</td><td>INJECT A</td></tr> <tr><td>L</td><td>H</td><td>H</td><td>L</td><td>9</td><td>ZERO OFF</td></tr> <tr><td>L</td><td>H</td><td>L</td><td>H</td><td>10</td><td>PLOT</td></tr> <tr><td>L</td><td>H</td><td>L</td><td>L</td><td>11</td><td>Switch polarity</td></tr> <tr><td>L</td><td>L</td><td>H</td><td>H</td><td>12</td><td>INJECT B</td></tr> <tr><td>L</td><td>L</td><td>H</td><td>L</td><td>13</td><td>SELECT</td></tr> <tr><td>L</td><td>L</td><td>L</td><td>H</td><td>14</td><td>QUIT</td></tr> <tr><td>L</td><td>L</td><td>L</td><td>L</td><td>15</td><td>ENTER</td></tr> </tbody> </table>	Input	Decimal	Function	23	10	22	9	code	(default settings)	H	H	H	H	0	Inactive	H	H	H	L	1	PROG R/S	H	H	L	H	2	PUMP R/S	H	H	L	L	3	FILL B/STEP	H	L	H	H	4	FILL A	H	L	H	L	5	REPORT	H	L	L	H	6	ZERO	H	L	L	L	7	MARK	L	H	H	H	8	INJECT A	L	H	H	L	9	ZERO OFF	L	H	L	H	10	PLOT	L	H	L	L	11	Switch polarity	L	L	H	H	12	INJECT B	L	L	H	L	13	SELECT	L	L	L	H	14	QUIT	L	L	L	L	15	ENTER
Input	Decimal		Function																																																																																																								
23	10	22	9	code	(default settings)																																																																																																						
H	H	H	H	0	Inactive																																																																																																						
H	H	H	L	1	PROG R/S																																																																																																						
H	H	L	H	2	PUMP R/S																																																																																																						
H	H	L	L	3	FILL B/STEP																																																																																																						
H	L	H	H	4	FILL A																																																																																																						
H	L	H	L	5	REPORT																																																																																																						
H	L	L	H	6	ZERO																																																																																																						
H	L	L	L	7	MARK																																																																																																						
L	H	H	H	8	INJECT A																																																																																																						
L	H	H	L	9	ZERO OFF																																																																																																						
L	H	L	H	10	PLOT																																																																																																						
L	H	L	L	11	Switch polarity																																																																																																						
L	L	H	H	12	INJECT B																																																																																																						
L	L	H	L	13	SELECT																																																																																																						
L	L	L	H	14	QUIT																																																																																																						
L	L	L	L	15	ENTER																																																																																																						
	10																																																																																																										
	23	2³ 																																																																																																									
	11 24 12																																																																																																										
		not assigned																																																																																																									

732 IC Detector	Pin	Function
Outputs		
	5	Ready L = ready (program inactive or not started) H = not ready (program running) $V_{CEO} = 40\text{ V}$ $I_C = 20\text{ mA}$
	18	Remote 1 Programmable remote output line 1 (command "remote" in program or event) $V_{CEO} = 40\text{ V}$ $I_C = 20\text{ mA}$
	4	Remote 2 Programmable remote output line 2 (command "remote" in program or event) $V_{CEO} = 40\text{ V}$ $I_C = 20\text{ mA}$
	17	Advance pulse (EOD) The EOD signal is outputted in the following cases: <ul style="list-style-type: none"> without program: on every "INJECT A" with program: each time a program is ended and on every return to the loop start with loop programs
	3	Remote 3 Programmable remote output line 3 (command "remote" in program or event) $V_{CEO} = 40\text{ V}$ $I_C = 20\text{ mA}$
	16	Error L = instrument error (as long as it exists) H = instrument ok $V_{CEO} = 40\text{ V}$ $I_C = 20\text{ mA}$
	1	Remote 4 Programmable remote output line 4 (command "remote" in program or event) $V_{CEO} = 40\text{ V}$ $I_C = 20\text{ mA}$
	2	Remote 5 Programmable remote output line 5 (command "remote" in program or event) $V_{CEO} = 40\text{ V}$ $I_C = 20\text{ mA}$
Voltages		
	15	$I \leq 40\text{ mA}$ $R_i \cong 12\ \Omega$
	14	
	25	
Contact arrangement at socket "Remote" (female)		
		Connecting cables to external devices must have an appropriate 25-pin connector (male). Metrohm offers connector sockets (K.210.9060) and the associated housing (K.210.0002) as an option for the preparation of such cables.
<p><i>No liability whatsoever will be accepted for damage arising from the improper connection of devices.</i></p>		

6.2.2 "733 IC Separation Center" interface


The "733 IC Separation Center" is normally attached to remote interface **14** "733 IC Separation Center" (see *section 2.3*). If the 732 IC Detector is operated without the 733 IC Separation Center, this remote interface can also be used for the connection of external devices. The remote interface "733 IC Separation Center" has the following pin assignment:

732 IC Detector	Pin	Function																						
Inputs																								
	21	<p>Position "FILL" at injection valve A L = Valve A in position "FILL" H = Valve A not in position "FILL"</p>																						
	9	<p>Position "INJECT" at injection valve A L = Valve A in position "INJECT" H = Valve A not in position "INJECT"</p>																						
	22	<p>Position "FILL" at injection valve B or "in position" at suppressor L = Valve B in position "FILL" or suppressor set to "in position" H = Valve B not in position "FILL" or suppressor not set to "in position"</p>																						
	10	<p>Position "INJECT" at injection valve B L = Valve B in position "INJECT" H = Valve B not in position "INJECT"</p>																						
	23	<p>Configuration of the 733 IC Separation Center</p> <table border="1"> <thead> <tr> <th colspan="2">Input</th> <th rowspan="2">Decimal code</th> <th rowspan="2">Status</th> </tr> <tr> <th>11</th> <th>23</th> </tr> </thead> <tbody> <tr> <td>H</td> <td>H</td> <td>0</td> <td>undefined</td> </tr> <tr> <td>H</td> <td>L</td> <td>1</td> <td>valve A</td> </tr> <tr> <td>L</td> <td>H</td> <td>2</td> <td>valve A + suppressor</td> </tr> <tr> <td>L</td> <td>L</td> <td>3</td> <td>valve A + valve B</td> </tr> </tbody> </table>	Input		Decimal code	Status	11	23	H	H	0	undefined	H	L	1	valve A	L	H	2	valve A + suppressor	L	L	3	valve A + valve B
Input			Decimal code	Status																				
11	23																							
H	H	0	undefined																					
H	L	1	valve A																					
L	H	2	valve A + suppressor																					
L	L	3	valve A + valve B																					
	11																							
	24	<p>Connection 733 IC Separation Center L = 733 IC Separation Center connected H = 733 IC Separation Center not accessible</p>																						
	12	not assigned																						

732 IC Detector	Pin	Function
Outputs		
	5	Valve A ® "FILL" Valve A is switched to position "FILL".
	18	Valve A ® "INJECT" Valve A is switched to position "INJECT".
	4	Lock 733 keypad L = <FILL> and <INJECT> keys are locked H = <FILL> and <INJECT> are not locked $V_{CE0} = 40\text{ V}$ $I_C = 20\text{ mA}$
	17	Valve B ® "FILL" or suppressor ® "STEP" Valve B is switched to position "FILL" or suppressor module is switched to next position.
	3	Valve B ® "INJECT" Valve B is switched to position "INJECT".
	2	Remote 6 Programmable remote output line 6 (command "remote" in program or event) $V_{CE0} = 40\text{ V}$ $I_C = 20\text{ mA}$
	1	Remote 7 Programmable remote output line 7 (command "remote" in program or event) $V_{CE0} = 40\text{ V}$ $I_C = 20\text{ mA}$
	16	Remote 8 Programmable remote output line 8 (command "remote" in program or event) $V_{CE0} = 40\text{ V}$ $I_C = 20\text{ mA}$
Voltages		
		5 V: inactive High (H) $I \leq 40\text{ mA}$ 0 V: active Low (L) $R_i \cong 12\ \Omega$
Contact arrangement at socket "733 IC Separation Center" (female)		
		Connecting cables to external devices must have an appropriate 25-pin connector (male). Metrohm offers connector sockets (K.210.9060) and the associated housing (K.210.0002) as an option for the preparation of such cables.
<i>No liability whatsoever will be accepted for damage arising from the improper connection of devices.</i>		

6.3 Analog output

For the attachment of recorders or integration systems, the 732 IC Detector offers the two analog outputs **11** (0...1 V) and **12** (0...10 mV) (see *section 3*). The attachment of these instruments is described in detail in *section 2.9*. The polarity of the output signal at the analog output sockets can always be switched under the <PARAM> key (see *section 4.5.1*). The circuit of the two analog outputs has the following appearance:








6.4 External power supply for 733 IC Separation Center

If the 733 IC Separation Center is operated without the 732 IC Detector, an external power supply must be attached to connection **48** "External Power". Power supply and connecting cable must comply with the following conditions:

"External Power" socket	Function	
	1 +5 V DC / 0.5 A	
	2 not assigned	
	3 +24 V DC / 2 A	
	4 0 V (digital) *	* Pin 4 and pin 5 must be connected in the vicinity of the power source ("zero point of star")
	5 0 V (analog) *	
	6 Ground	
	7 Ground	

6.5 Valve interfaces

The two valve interfaces **38** and **46** of the 733 IC Separation Center can be used for the attachment of external devices (e.g. 750 Autosampler, «IC Metrodata for Win95») which can be employed to control the injection valves or the suppressor module (see *section 2.3*). The two valve interfaces for the injection valves A and B (or the suppressor module) are identical and have the following connections

733 IC Sep. Center	Connection	Function
Inputs		
	Ground	Ground
	Fill	Valve @ "FILL" Valve is switched to the position "FILL" (or suppressor module is switched further)
	Inject	Valve @ "INJECT" Valve is switched to the position "INJECT"
Outputs		
	COM RUN	Position "Fill" A pulse is outputted when the valve is switched to the position "FILL".
	COM RUN	Integrator Start A pulse is outputted when the valve is switched to the position "INJECT".
<i>No liability whatsoever will be accepted for damage arising from the improper connection of devices.</i>		

7 Appendix

7.1 Technical data

7.1.1 732 IC Detector

Conductivity measurement

<i>Measurement ranges</i>	100, 200, 500, 1000, 2000, 5000, 10 000 $\mu\text{S/cm}$
<i>Full scale</i>	0.05...10 000 $\mu\text{S/cm}$
<i>Temperature correction</i>	Automatic correction of the measured conductivity to reference temperature 20°C with the adjustable temperature coefficients (1.5 %/°C or 2.5 %/°C)
<i>Accuracy, absolute value</i>	< 4 % of full-scale value for all ranges under the following reference conditions: Adjustment with reference resistance to 1.000 V (range 1 mS/cm, cell constant 16.7 /cm, temp. coeff. 2.5 %/°C, frequency 5 kHz, room temperature) to accuracy < ± 2 mV, measurement accuracy < 0.5 % after 30 min warm-up time
<i>Linearity</i>	Up to ± 150 % of each full-scale range Deviations < 0.5 % of full-scale range
<i>Temperature dependence</i>	Typically 25 ppm/°C
<i>Dependence of cell const..</i>	< ± 1 % over entire adjustment range
<i>Drift (electronic)</i>	< 0.0013 % of selected measurement range /h/°C
<i>Noise (electronic)</i>	Typically < 0.0003 % of selected meas. range
<i>Damping</i>	2-stage damping (Bessel 4 th order) "damping = off" 0.25 s (10...90 %) "damping = on" 2.00 s (10...90 %)

Auto-zero

<i>Function</i>	Automatic zero setting (electronic background compensation) over entire selectable measurement range
<i>Initiation</i>	Manually or externally (RS232, remote)
<i>Maximum error</i>	± 0.6 % of full-scale range
<i>Linear range</i>	± 150 % of full-scale range

Marker

<i>Signal</i>	Ca. 10 % of full-scale range
<i>Initiation</i>	Manual or externally (RS232, remote)

Conductivity detector

Construction Thermostatted conductivity detector with 2 ring-shaped steel electrodes

Measurement principle Alternating current measurement with following cell supply:

<i>Meas. Range</i>	<i>Amplitude</i>	<i>Frequency</i>
100 $\mu\text{S/cm}$	1 V peak	714 Hz
200 $\mu\text{S/cm}$	1 V peak	5000 Hz
500 $\mu\text{S/cm}$	1 V peak	5000 Hz
1 mS/cm	1 V peak	5000 Hz
2 mS/cm	0.5 V peak	5000 Hz
5 mS/cm	0.2 V peak	5000 Hz
10 mS/cm	0.1 V peak	5000 Hz

Cell volume 1.5 μL

Cell constant 16.7 /cm; adjustable 13.0...21.0 /cm

Maximum back pressure for measuring cell 5.0 MPa (50 bar)

Thermostating Connectable dynamic control to adjustable operating temperature, indicated by LED

Operating temperature Adjustable in steps of 5°C from 25...45°C

Max. temperature deviation $\pm 1.5^\circ\text{C}$

Temperature stability $\leq 0.01^\circ\text{C}$ at constant ambient temperature

Parts and controls

Display LCD, 2 lines, each of 24 characters (character height 5 mm)

Indicators LED for auto-zero display
LED for display of program status
LED for overload display
LED for display of block heating

Keypad Chemically resistant membrane keypad made of polyester with function and numeric keys

Mains switch At rear of instrument

Mains connection

Voltage 115 V: 100...120 V $\pm 10\%$
230 V: 220...240 V $\pm 10\%$
Switching with mains voltage selector in fuse holder (see *section 2.4.1*)

Frequency 50...60 Hz

Power consumption 70 VA

Fuse 5 mm \varnothing , 20 mm length
100...120 V: 0.63 A (slow-blow)
220...240 V: 0.315 A (slow-blow)

Interfaces

<i>RS232 interfaces</i>	Specifications, see <i>section 6.1</i>
<i>Remote interfaces</i>	Specifications, see <i>section 6.2</i>
<i>Analog output</i>	Specifications, see <i>section 6.3</i>

Safety specifications

<i>Construction/testing</i>	According to IEC 1010 / EN 61010 / UL 3101-1, protection class 1, degree of protection IP40
<i>Safety directions</i>	The Instructions for Use include information and warnings which must be heeded by the user to assure safe operation of the instrument.

Electromagnetic compatibility (EMC)

<i>Emitted interference</i>	Standards met: EN 55011 (class B), EN 55022 (class B), EN 50081-1/2
<i>Immunity to interference</i>	Standards met: IEC801-2/IEC1000-4-2 (class 4), IEC801-3/ IEC1000-4-3 (class 3), IEC801-4/IEC1000-4-4 (class 4), IEC801-5/IEC1000-4-5 (class 2/3), IEC801-6/IEC1000-4-6 (class 3), EN50082-2, EN61000-3-2/3/IEC1000-3-2/3, EN50093/ IEC1000-4-11

Ambient temperature

<i>Nominal operating range</i>	+5...+45°C (at 20...80 % atmospheric humidity)
<i>Storage, transport</i>	-20...+70°C

Diagnostic tests

<i>Self-diagnostic test</i>	Automatic self-diagnosis when instrument switched on
<i>User diagnostic test</i>	Built-in diagnostic program (see <i>section 5.4</i>)

Housing

<i>Material of cover</i>	Polyurethane rigid foam (PUR) with fire protection for fire class UL94VO, FCH-free
<i>Material of base</i>	Steel, enameled
<i>Width</i>	255 mm
<i>Height</i>	130 mm
<i>Depth</i>	343 mm
<i>Weight</i>	8.2 kg (with detector block, without accessories)

7.1.2 733 IC Separation Center

Parts and controls

<i>Keypad</i>	Chemically resistant membrane keypad made of polyester with function keys
<i>Indicators</i>	LEDs for display of valve position

Power supply

<i>Internal</i>	Internal supply of 732 IC Detector via connecting cable
<i>External</i>	External supply via DIN connector: 5 V / 0.5 A 24 V / 2 A (transient, 200...300 ms)

Interfaces

<i>Valve interfaces</i>	Specifications, see <i>section 6.5</i>
-------------------------	--

Safety specifications

<i>Construction/testing</i>	According to IEC 1010 / EN 61010 / UL 3101-1, protection class 1, degree of protection IP40
<i>Safety notes</i>	The Instructions for Use include information and warnings which must be heeded by the user to assure safe operation of the instrument.

Electromagnetic compatibility (EMC)

<i>Emitted interference</i>	Standards met: EN 55011 (class B), EN 55022 (class B), EN 50081-1/2
<i>Immunity to interference</i>	Standards met: IEC801-2/IEC1000-4-2 (class 4), IEC801-3/ IEC1000-4-3 (class 3), IEC801-4/IEC1000-4-4 (class 4), IEC801-5/IEC1000-4-5 (class 2/3), IEC801-6/IEC1000-4-6 (class 3), EN50082-2

Ambient temperature

<i>Nominal operating range</i>	+5...+45°C (at 20...80 % atmospheric humidity)
<i>Storage, transport</i>	-20...+70°C

Housing

<i>Material</i>	Polyurethane rigid foam (PUR) with fire protection for fire class UL94VO, FCH-free
<i>Width</i>	255 mm
<i>Height</i>	385 mm
<i>Depth</i>	343 mm
<i>Weight</i>	6.4 kg (without accessories)

7.2 Standard equipment






Subject to changes !
All dimensions are given in mm.

7.2.1 732 IC Detector

The 732 IC Detector is available in two versions:

- **2.732.0010** IC Detector with standard detector block
- **2.732.0110** IC Detector with metal-free detector block

These instruments include the following parts:








Quant.		Order No.	Description												
2.732.0010	2.732.0110														
1	-	1.732.0100	Detector block with permanently attached connecting cable to 732 IC Detector												
-	1	1.732.0110	Detector block, metal-free with permanently attached connecting cable to 732 IC Detector												
1	1	6.2125.060	Connecting cable Connecting cable 732 IC Detector (RS232) – PC 												
1	1	6.2125.090	Connecting cable Connecting cable 732 IC Detector – 733 IC Separation Center 												
1	1	6.2248.000	Magnetic plate for identification of the 732 IC Detector 												
1	1	6.2122.0X0	Mains cable to customer's specifications: <table style="width: 100%; border: none;"> <tr> <td style="border: none;"><u>Cable socket</u></td> <td style="border: none;"><u>Cable connector</u></td> <td style="border: none;"></td> </tr> <tr> <td style="border: none;">Type IEC 320/C 13</td> <td style="border: none;">Type SEV 12 (CH...)</td> <td style="border: none;">6.2122.020</td> </tr> <tr> <td style="border: none;">Type IEC 320/C 13</td> <td style="border: none;">Type CEE (7), VII (D...)</td> <td style="border: none;">6.2122.040</td> </tr> <tr> <td style="border: none;">Type CEE (22), V</td> <td style="border: none;">Type NEMA 5-15 (USA...)</td> <td style="border: none;">6.2122.070</td> </tr> </table>	<u>Cable socket</u>	<u>Cable connector</u>		Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020	Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040	Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070
<u>Cable socket</u>	<u>Cable connector</u>														
Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020													
Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040													
Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070													
1	1	8.732.1033	Instructions for Use (English) for 732 IC Detector and 733 IC Separation Center												
1	1	8.732.1043	Quick reference guide (English) for 732 IC Detector and 733 IC Separation Center												









7.2.2 733 IC Separation Center

The 733 IC Separation Center is available in the following versions:

- **2.733.0010** Instrument with 1 injection valve
- **2.733.0020** Instrument with 2 injection valves
- **2.733.0120** Instrument with 2 injection valves, metal-free
- **2.733.0030** Instrument with 1 injection valve and 1 suppressor module
- **2.733.0130** Instrument with 1 inj. valve and 1 suppr. module, metal-free








These instruments contain the following parts:










Quant.					Order No.	Description
2.733.0010	2.733.0020	2.733.0120	2.733.0030	2.733.0130		
1	1	1	1	1	6.1816.000	Silicone tubing Drain tube for inner compartment, length = 0.5 m 
1	1	1	1	1	6.1831.010	PEEK capillary Length = 3 m 
1	2	2	1	1	6.2027.030	Column holder Diameter d = 8.5 mm 
1	2	2	1	1	6.2027.040	Column holder Diameter d = 11.3 mm 
1	2	2	1	1	6.2027.050	Column holder Diameter d = 15.0 mm 
1	1	-	1	-	6.2620.020	Steel capillary Length = 3 m 
2	4	-	2	-	6.2620.060	Coupling 1/16" - 1/16" for connection of 1/16" steel or PEEK capillaries; incl. two 6.2620.000 Compression screws and two 6.2620.010 Ferrules 


Quant.					Order No.	Description
2.733.0010	2.733.0020	2.733.0120	2.733.0030	2.733.0130		
1	1	-	1	-	6.2621.000	Adjustable spanner 
1	1	-	1	-	6.2621.050	Open-end spanner 1/4" 
1	1	-	1	-	6.2621.060	Open-end spanner 5/16" 
1	2	2	2	2	6.2744.010	PEEK compression fitting For the connection of 6.1831.010 PEEK capillaries or 6.1822.010 PTFE microcapillaries, set of 5 
1	2	2	1	1	6.2744.020	Coupling 1/16" – Luer Coupling for connection of a 6.1803.000 PTFE capillary to connection 21 or 27 of the 733 IC Separation Center when a 750 Autosampler or a 766 IC Sample Processor is used 
-	-	4	-	2	6.2744.040	PEEK coupling for the connection of 1/16" capillaries 
1	2	2	1	1	6.2816.020	Syringe made of PP, volume = 10 mL; for manual filling of the sample loop 
-	-	-	2	2	6.2821.100	Filter unit PEEK 2 mm To avoid contamination of the suppressor module. Spare part: 6.2821.110 Filter 

7.3 Optional accessories






7.3.1 Accessories for 733 IC Separation Center









Order No.	Description	
6.1825.XXX	Sample loop, made of PEEK For injection valve; incl. 2 PEEK compression fittings 6.2744.010 6.1825.230: Volume = 10 μ L 6.1825.210: Volume = 20 μ L 6.1825.220: Volume = 100 μ L	
6.2620.000	Pressure screw, made of steel For connection of 6.2620.020 Steel capillaries, set of 5.	
6.2620.010	Ferrule, made of steel For connection of 6.2620.020 Steel capillaries, set of 5.	
6.2620.040	Coupling $\frac{1}{4}$" - $\frac{1}{4}$" Connector for plastic separating columns with $\frac{1}{4}$ "-28 thread.	
6.2620.XXX	Sample loop, made of steel For injection valve. 6.2620.100: Volume = 10 μ L 6.2620.110: Volume = 20 μ L 6.2620.120: Volume = 100 μ L	
6.2620.150	Pulsation dampener MF Metal-free pulsation dampener to reduce pulsations and prolong the life of separating columns.	
6.2621.040	Capillary tubing cutter for steel capillaries	











Order No.	Description	
6.2621.080	Capillary tubing cutter for plastic capillaries For 6.1831.010 PEEK capillaries and 6.1822.010 PTFE microcapillaries incl. 5 additional blades	
6.2744.030	PEEK Coupling Connection between 6.2744.010 PEEK compression fitting and 6.1826.0X0 Pump tubing; set of 4	
6.2821.000	Filter unit Manufit To avoid contamination through abrasive particles from piston seals. For the pressure range 25...50 MPa with steel capillaries.	
6.2821.010	PTFE gasket Spare part for 6.2821.000 Filter unit Manufit and 6.1005.000, 6.1005.030, 6.1007.000 and 6.1010.000 IC Separating Columns; set of 4.	
6.2821.020	Steel mesh Spare part for 6.2821.000 Filter unit Manufit and 6.1005.000, 6.1005.030, 6.1007.000 and 6.1010.000 IC Separating Columns; set of 4.	
6.2821.040	Cartridge head To hold precolumn cartridges; mounted directly on the separating column.	
6.2821.050	Twin cartridge holder To hold precolumn cartridges; installed in the inlet capillary of the separating column.	
6.2821.080	Steel spacer Spare part for 6.2821.040 Cartridge head.	
6.2821.090	Aspirating filter PE 20 mm For 6.1834.000 Aspirating tubing (supplied with the 709 IC Pump); set of 5.	









Order No.	Description
6.2821.110	<p>Filter for Filter unit PEEK 2 mm Spare part for 6.2821.100 Filter unit PEEK. Set of 10.</p> 
6.5324.000	<p>Bottle rack For holding 3 supply bottles with eluent, regeneration solution and rinsing solution, incl. accessories (bottles, siphon, etc.)</p>

7.3.2 Separating columns and precolumns


Order No.	Description
6.1005.000	<p>IC anion column PRP-X100 (125 mm) For the determination of anions without chemical suppression. Column dimensions: 125 × 4.0 mm Precolumn: 6.1005.020</p> 
6.1005.010	<p>IC anion column PRP-X100 (250 mm) For the determination of anions without chemical suppression. Column dimensions: 250 × 4.0 mm Precolumn: 6.1005.020</p> 
6.1005.020	<p>IC precolumn cartridge PRP-X100 To prolong the service life of 6.1005.000 and 6.1005.010 IC anion columns PRP-X100. Column dimensions: 20 × 4.0 mm Installation using 6.2821.040 Cartridge head or 6.2821.050 Twin cartridge holder.</p> 
6.1005.030	<p>IC exclusion column PRP-X300 For the determination of organic acids without chemical suppression. Column dimensions: 250 × 4.0 mm Precolumn: 6.1005.040</p> 
6.1005.040	<p>IC precolumn cartridge PRP-X300 To prolong the service life of the 6.1005.030 IC exclusion column PRP-X300. Column dimensions: 20 × 4.0 mm Installation using 6.2821.040 Cartridge head or 6.2821.050 Twin cartridge holder.</p> 










Order No.	Description	
6.1005.050	IC precolumn cartridge PRP-1 To prolong the service life of the 6.1009.000 IC anion column SUPERSEP and the 6.1006.100 IC anion column METROSEP Anion Dual 2. Column dimensions: 20 × 4.0 mm Installation using 6.2821.040 Cartridge head or 6.2821.050 Twin cartridge holder.	
6.1005.100	IC anion column Star-Ion A300 For the determination of anions with chemical suppression. Column dimensions: 100 × 4.6 mm	
6.1005.200	IC anion column Organic Acids For the determination of organic acids. Column dimensions: 250 × 7.5 mm	
6.1006.020	IC column cartridge METROSEP Anion Dual 1 For the determination of anions with and without chemical suppression. Column dimensions: 150 × 3.0 mm Installation using 6.2828.000 Glass cartridge holder.	
6.1006.030	IC precolumn cartridge METROSEP Anion Dual 1 Set of 3 To prolong the service life of the 6.1006.020 IC column cartridge METROSEP Anion Dual 1. Column dimensions: 30 × 3.0 mm Installation using 6.2828.010 Precolumn cartridge holder.	
6.1006.040	IC precolumn cartridge METROSEP Anion Dual 1 Set of 3 To prolong the service life of the 6.1006.020 IC column cartridge METROSEP Anion Dual 1. Column dimensions: 70 × 3.0 mm Installation using 6.2828.020 Glass cartridge holder.	
6.1006.100	IC anion column METROSEP Anion Dual 2 For the determination of anions with and without chemical suppression. Column dimensions: 75 × 4.6 mm Precolumn: 6.1005.050 (installation with 6.2821.050 Twin cartridge holder)	
6.1006.200	Preconcentration cartridge METROSEP Anion For anion preconcentration. Column dimensions: 30 × 3.0 mm Installation using 6.2828.010 Precolumn cartridge holder.	

Order No.	Description	
6.1007.000	IC cation column Nucleosil 5SA For the determination of divalent cations without chemical suppression. Column dimensions: 125 × 4.0 mm Precolumn: 6.1007.010	
6.1007.010	IC precolumn cartridge Nucleosil 5SA To prolong the service life of the 6.1007.000 IC cation column Nucleosil 5SA. Column dimensions: 20 × 4.0 mm Installation using 6.2821.040 Cartridge head or 6.2821.050 Twin cartridge holder.	
6.1008.010	IC cation column Hyperrez Monovalent For the determination of monovalent cations without chemical suppression. Column dimensions: 50 × 4.6 mm	
6.1009.000	IC anion column SUPERSEP For the determination of anions without chemical suppression. Column dimensions: 100 × 4.6 mm Precolumns: 6.1009.010 IC anion column SUPERSEP or 6.1005.010 IC precolumn cartridge	
6.1009.010	IC anion precolumn SUPERSEP To prolong the service life of the 6.1009.000 IC anion column SUPERSEP.	
6.1010.000	IC cation column METROSEP Cation 1-2 For the determination of monovalent and divalent cations without chemical suppression. Column dimensions: 125 × 4.0 mm Precolumn: 6.1010.010	
6.1010.010	IC precolumn cartridge METROSEP Cation 1-2 To prolong the service life of the 6.1010.000 IC cation column METROSEP Cation 1-2. Column dimensions: 20 × 4.0 mm Installation using 6.2821.040 Cartridge head or 6.2821.050 Twin cartridge holder.	
6.1012.X00	Sample pretreatment cartridge IC-RP For non-polar solid phase extraction. Removes organic substances; for the enrichment of heavy metals. With Luer connection. 6.1012.000: set of 50 6.1012.100: set of 10	
6.1012.X10	Sample pretreatment cartridge IC-H Cation exchanger in H ⁺ form. Removes interfering cations, CO ₃ ²⁻ , HCO ₃ ⁻ or for alkaline samples. With Luer connection. 6.1012.010: set of 50 6.1012.110: set of 10	

Order No.	Description	
6.1012.X20	Sample pretreatment cartridge IC-Ag Cation exchanger in Ag ⁺ form. Removes halides. With Luer connection. 6.1012.020: set of 50 6.1012.120: set of 10	
6.1012.X30	Sample pretreatment cartridge IC-OH Cation exchanger in OH ⁻ form. For highly acidic samples. With Luer connection. 6.1012.030: set of 50 6.1012.130: set of 10	
6.1012.200	Sample pretreatment cartridge Chromafix C18 Removes organic substances (not suitable for fluoride determinations). With Luer connection. 6.1012.200: set of 50	
6.2828.000	Glass cartridge holder To hold the 6.1006.0020 Column cartridge METROSEP Anion Dual 1.	
6.2828.010	Precolumn cartridge holder To hold the 6.1006.0030 Precolumn cartridge METROSEP Anion Dual 1.	
6.2828.020	Glass cartridge holder To hold the 6.1006.0040 Column cartridge METROSEP Anion Dual 1.	
6.2832.000	Suppressor rotor Replacement cartridge for Metrohm suppressor module	
6.2832.010	Connection piece for suppressor rotor with input and output leads	

7.3.3 Additional devices and cables

Order No.	Description
2.145.0300	Serial/parallel converter For the connection of printers with a parallel interface to the RS232 interface of 732 IC Detector
2.709.0X10	709 IC Pump Serial, high-performance dual piston pump with minimal residual pulsation and outstanding flow constancy; incl. accessories. Versions: 2.709.0010 Standard version 2.709.0110 Metal-free version
2.714.0310	IC Metrodata for Win95 Chromatography data system comprising a PC board and evaluation software, for the automatic evaluation of chromatograms with a PC. Versions: 2.714.0310 For 1 chromatography system with 2 channels 2.714.0320 For 2 chromatography systems each with 2 channels
2.750.0010	750 Autosampler Sampler for the automation of sample injection. Capacity: 128 sample vessels each with an effective volume of ca. 700 µL; incl. accessories. Accessories: 6.2413.000 Glass Sample Vessels, set of 1000 6.2743.000 PP Sample Vessels, set of 1000 6.2743.010 Polyethylene Stoppers, transparent, set of 1000 6.2743.020 Polyethylene Stoppers, red, set of 1000 6.2743.030 Filter Stoppers, set of 100
2.752.0010	752 Pump Unit Dual-channel peristaltic pump for the regeneration and rinsing solutions for the suppressor module; incl. 2 lengths of 6.1826.050 Pump tubing.
2.754.0010	754 Dialysis Unit Instrument for automatic sample dialysis for sample pretreatment, incl. dialysis cell and dual-channel peristaltic pump.
2.766.0010	766 IC Sample Processor Sampler for the automation of sample injection. Capacity: 127 sample vessels each with an effective volume of ca. 11 mL; incl. accessories. Accessories: 6.2743.050 PP Sample vessels, set of 2000 6.2743.060 Polyethylene stoppers, transparent, set of 1000
2.791.0020	791 VA Detector Current measuring instrument for electrochemical detection; complete with detection cell and accessories (without working electrode). Working electrodes: 6.0807.000 Mini carbon paste electrode 6.0807.010 Mini glassy carbon electrode 6.1204.100 Ultra Trace graphite tip * 6.1204.120 Platinum tip * 6.1204.130 Silver tip * 6.1204.140 Gold tip * (* 6.2103.110 Contact pin needed)
6.2115.010	Cable Connecting cable 732 IC Detector (analog output) – recorder 

Order No.	Description
6.2115.060	<p>Cable Connecting cable 732 IC Detector (analog output) – integrator or data acquisition system (analog cable)</p> 
6.2115.070	<p>Cable Connecting cable between 733 IC Separation Center and integrator or data acquisition system (start cable)</p> 
6.2125.010	<p>Cable Connecting cable 732 IC Detector (RS232 interface) – PC Adapter cable 25-pin to 9-pin.</p> 
6.2125.020	<p>Cable Connecting cable 732 IC Detector (RS232 interface) – Seiko printer DPU-411 Can also be used as an RS232 extension cable.</p> 
6.2125.040	<p>Cable Connecting cable 732 IC Detector (RS232 interface) – Epson printer with 9-pin circular connector</p> 
6.2125.050	<p>Cable Connecting cable 732 IC Detector (RS232 interface) – printer (IBM Proprinter, Citizen IDP562-RS, Epson with serial interface, Epson LX-300, HP Deskjet/Laserjet with serial interface, etc.)</p> 
6.2125.060	<p>Cable Connecting cable 732 IC Detector (RS232 interface) – PC or 709 IC Pump</p> 
6.2125.110	<p>Cable Connecting cable 732 IC Detector (RS232 interface) – PC</p> 
6.2125.130	<p>Cable Connecting cable 732 IC Detector (RS232 interface) – Seiko printer DPU-414 Adapter cable 25-pin to 9-pin.</p> 

7.4 Warranty and conformity

7.4.1 Warranty

The warranty on our products is limited to defects that are traceable to material, construction or manufacturing error which occur within 12 months from the day of delivery. In this case, the defects will be rectified in our workshops free of charge. Transport costs are to be paid by the customer.

For day and night operation, the warranty is limited to 6 months.

Glass breakage in the case of electrodes or other parts is not covered by the warranty. Checks which are not a result of material or manufacturing faults are also charged during the warranty period. For parts of outside manufacture insofar as these constitute an appreciable part of our instrument, the warranty stipulations of the manufacturer in question apply.

With the regard to the guarantee of accuracy, the technical specifications in the instruction manual are authoritative.

Concerning defects in material, construction or design as well as the absence of guaranteed features, the orderer has no rights or claims except those mentioned above.

If damage of the packaging is evident on receipt of a consignment or if the goods show signs of transport damage after unpacking, the carrier must be informed immediately and a written damage report demanded. lack of an official damage report releases Metrohm from any liability to pay compensation.

If any instruments and parts have to be returned, the original packaging should be used if at all possible. This applies above all to instruments, electrodes, burette cylinders and PTFE pistons. Before embedment in wood shavings or similar material, the parts must be packed in a dust-proof package (for instruments, use of a plastic bag is imperative). If open assemblies are enclosed in the scope of delivery that are sensitive to electromagnetic voltages (e.g. data interfaces etc.) these must be returned in the associated original protective packaging (e.g. conductive protective bag). (Exception: assemblies with built-in voltage source belong in a non-conductive protective packaging). For damage which arises as a result of non-compliance with these instructions, no warranty responsibility whatsoever will be accepted by Metrohm.

7.4.2 EU Declaration of conformity



EU Declaration of Conformity

The METROHM AG company, Herisau, Switzerland hereby certifies, that the instrument:

732 IC Detector

meets the requirements of EC Directives 89/336/EWG and 73/23/EWG.

Source of the specifications:

- | | |
|--------------|--|
| EN 50081-1/2 | Electromagnetic compatibility, basic specification
Emitted Interference |
| EN 50082-2 | Electromagnetic compatibility, basic specification
Interference Immunity |
| EN 61010 | Safety requirements for electrical laboratory measurement
and control equipment |

Description of the instrument:

Thermostatted conductivity detector for ion chromatographic analyses with electronic or chemical suppression.

Herisau, September 30, 1995



Dr. J. Frank
Development Manager

Ch. Buchmann
Production and
Quality Assurance Manager



EU Declaration of Conformity

The METROHM AG company, Herisau, Switzerland hereby certifies, that the instrument:

733 IC Separation Center

meets the requirements of EC Directives 89/336/EWG and 73/23/EWG.

Source of the specifications:

EN 50081-1/2	Electromagnetic compatibility, basic specification Emitted Interference
EN 50082-2	Electromagnetic compatibility, basic specification Interference Immunity

Description of the instrument:

Thermally and electronically isolated wet part for ion chromatography.

Herisau, September 30, 1995



Dr. J. Frank

Ch. Buchmann

Development Manager

Production and
Quality Assurance Manager

7.4.3 Certificate of conformity and system validation

Certificate of Conformity and System Validation	
<p>This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.</p>	
<p>Name of commodity:</p> <p>System software:</p> <p>Name of manufacturer:</p> <p>Principal technical information:</p>	<p>732 IC Detector</p> <p>Stored in ROMs</p> <p>Metrohm Ltd., Herisau, Switzerland</p> <p>Voltages: 100...120, 220...240 V Frequency: 50...60 Hz</p>
<p>This Metrohm instrument has been built and has undergone final type testing according to the standards:</p> <p style="margin-left: 40px;">IEC801-2/IEC1000-4-2 (class 4), IEC801-3/ IEC1000-4-3 (class 3), IEC801-4/IEC1000-4-4 (class 4), IEC801-5/IEC1000-4-5 (class 2/3), IEC801-6/IEC1000-4-6 (class 3), EN50082-2, EN61000-3-2/3/IEC1000-3-2/3, EN50093/IEC1000-4-11, EN55011 (class B), EN55022 (class B), EN50081-1/2 — <i>Electromagnetic compatibility</i></p> <p style="margin-left: 40px;">IEC1010, EN61010, UL3101-1 — <i>Security specifications</i></p> <p>It has also been certified by the Swiss Electrotechnical Association (SEV), which is member of the International Certification Body (CB/IEC).</p> <p>The technical specifications are documented in the instruction manual.</p> <p>The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance. The features of the system software are documented in the instruction manual.</p>	
<p>Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.</p>	
<p>Herisau, September 30, 1995</p> <div style="display: flex; justify-content: space-around; margin-top: 20px;"> <div style="text-align: center;">  <p>Dr. J. Frank</p> <p>Development Manager</p> </div> <div style="text-align: center;">  <p>Ch. Buchmann</p> <p>Production and Quality Assurance Manager</p> </div> </div>	

Certificate of Conformity and System Validation

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

Name of commodity:	733 IC Separation Center
Name of manufacturer:	Metrohm Ltd., Herisau, Switzerland
Principal technical information:	Voltages: 5 V DC, 24 V DC

This Metrohm instrument has been built and has undergone final type testing according to the standards:

IEC801-2/IEC1000-4-2 (class 4), IEC801-3/ IEC1000-4-3 (class 3),
 IEC801-4/IEC1000-4-4 (class 4), IEC801-5/IEC1000-4-5 (class 2/3),
 IEC801-6/IEC1000-4-6 (class 3), EN50082-2, EN55011 (class B),
 EN55022 (class B), EN50081-1/2 — *Electromagnetic compatibility*

IEC1010, EN61010, UL3101-1 — *Security specifications*

It has also been certified by the Swiss Electrotechnical Association (SEV), which is member of the International Certification Body (CB/IEC).

The technical specifications are documented in the instruction manual.

Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.

Herisau, September 30, 1995



Dr. J. Frank

Ch. Buchmann

Development Manager

Production and
Quality Assurance Manager

7.5 Index

7

709see IC Pump 709
 732 see IC Detector 732
 733 see IC Separation Center 733
 733 not accessible 95
 750 see Autosampler 750
 752 see Pump Unit 752
 754 see Dialysis Unit 754
 766 see IC Sample Processor 766
 791 see VA Detector 791

A

Abort entry 76
 Abort of the
 deletion operation 113,116,118
 Abort of the storage operation 118
 Abort printing operations 76
 Abrasive particles 134
 abs. cond. 82
 Absolute conductivity 49,81,82,89,119
 Absolute time 100
 Accuracy 189
 Acknowledge error messages 76
 Additional devices and cables 202
 Adjustable spanner (6.2621.000) 195
 Advance pulse 184
 all 124
 Alternating current measurement 190
 Ambient temperature 89,191,192
 Analog output **11**
 Circuit 187
 Connection of a recorder 41
 Figure 4
 Settings 99
 Signal output 98
 Analog output **12**
 Circuit 187
 Connection of a recorder 41
 Figure 4
 Settings 99
 Signal output 98
 Anion determination 129
 Appendix 189
 Application Bulletins 8,133
 Arrangement of the instruments 12
 Aspirating Filter
 (6.2821.090) 133,134,197
 Aspirating tubing
 see PTFE aspirating tubing
 Assembling the suppressor 142
 Assembly 174
 Assignment of the
 remote input lines 85
 auto start 100
 auto step 96
 Automatic measured value output 121
 Automatic measured value printout 90
 Automatic messages on change 173
 Automatic start 100

Automatic triggering of "STEP" 96
 Autosampler 750
 Connection 41
 Ordering designation 202
 Auto-zero 60,62,111,189
 Auto-zero display **6**
 Figure 3
 Function 105,106
 Auto-zero
 function 97,98,105,106,111,170
 Auto-zero reference value 120
 Auto-zero value 81,89,119
 Auxiliary function keys **4**
 Figure 3
 Overview 72,76

B

Back pressure 190
BACKGROUND 120
 Background compensation
 see Auto-zero
 Baseline with high noise level 148
 Basic elements of the assembly 174
 Basic settings 68,84
 baud rate 93,94
 Beep 79
 beep 92
 Beeper 92
 Block diagram 2
 Blockage 141,143
 Bottle rack (6.5324.000)
 Connection 16
 Ordering designation 198

C

Cable (6.2115.010) 41,202
 Cable (6.2115.060) 203
 Cable (6.2115.070) 203
 Cable (6.2125.010) 44,203
 Cable (6.2125.020) 43,203
 Cable (6.2125.040) 42,203
 Cable (6.2125.050) 42,43,203
 Cable (6.2125.060) 21,44,45,193,203
 Cable (6.2125.090) 193
 Cable (6.2125.110) 203
 Cable (6.2125.130) 43,203
 Calibration 59
 Call-up of objects 162
 Capacity factor 133
 Capacity problems 140
 Capillaries 19
 Capillary **51**
 Figure 20,22
 Steel connectors 19
 Capillary **65**
 Connection 24
 Figure 25,34,36
 Capillary **85**

 Figure 34,36
 Capillary connections 19
 Capillary tubing cutter
 (6.2621.040) 19,196
 Capillary tubing cutter
 (6.2621.080) 19,197
 Care 136
 Cartridge head (6.2821.040)
 Installation 30
 Ordering designation 197
 Cation determination 126
 Caution 9
 Cell constant 40,51,89,190
 cell constant 89
 Cell supply 190
 Cell volume 190
 Certificate of conformity
 and system validation 207
 Change full-scale range 110
 Change in the sign 75
 Change time for
 program steps 109,110
 Changing of separating column 138
 Changing the separation system 138
 Character selection 80
 Character set 90
 Check 12
 Check CompCoarse 156
 Check CompFine 156
 Check display 152
 Check keypad 151
 Check range 155
 Check sense 156
 Check Ucell 155
 Checking working memory 151
 Chemical suppression 33
 Chromatography data system 41
 Citizen printer IDP562 RS 42,90,123
 Cleaning the suppressor 141
 <CLEAR> 76,80
 Closing text entry 80
 Column see Separating column
 Column connection capillary **67**
 Figure 25,26,29,34,36
 Installation 24,26,28,30,31,35,37
 Passivation 27,138
 Column holder **84**
 Figure 34,36
 Installation 37
 Ordering designation 194
 Comment 9
 Communication 1
 Compression fitting **52**
 Figure 20,22
 Installation 20
 Ordering designation 195
 Conditioning 40,58
 Conductivity 81,89
 Conductivity cell 89
 Conductivity detector 88,97
 Conductivity measurement 189
 Conductivity measurement value 122
 <CONFIG> 74,88

config.....	124	Connection of	
CONFIG/733 IC Sep. Cent.	95	733 IC Separation Center.....	13
CONFIG/733/suppressor.....	96	Connection of 750 Autosampler.....	41
CONFIG/733/valve A.....	95	Connection of	
CONFIG/733/valve B.....	96	766 IC Sample Processor	41
CONFIG/aux/beep.....	92	Connection of 791 VA Detector.....	41
CONFIG/aux/event.....	92	Connection of a PC.....	44
CONFIG/auxiliaries.....	91	Connection of a printer.....	42
CONFIG/detector.....	88	Connection of a recorder.....	41
CONFIG/print meas. value.....	90	Connection of devices to the remote	
CONFIG/printer.....	89	interface.....	44
CONFIG/RS settings.....	93	Connection of external devices.....	41
CONFIG/RS settings 709.....	94	Connection of syringe and aspirating	
Configuration.....	50,74,88,166	tubing.....	16
Configuration report.....	124	Connection of waste container.....	14,15
Confirm of deletion.....	113,116	Connection piece 98	
Confirmation.....	79	Cleaning.....	142,144
Confirmation of deletion.....	118	Figure.....	142,143
Confirmation of overwriting.....	118	Ordering designation.....	201
Confirmation of parameters.....	76	Replacement.....	143
Conformity.....	205	Connection to injection valve.....	24,25
connected.....	94	Connection to RS interface.....	94
Connection 13		Connections.....	136
Connection of detector block.....	13,14	Connections at	
Figure.....	4	suppressor module.....	38
Connection 14		Connector 53	
Connection of 733.....	14,15	Connection of	
Figure.....	4	inlet capillary 90	38,39
Connection 15		Figure.....	22
Connection of 709.....	21	Installation.....	22
Figure.....	4	Ordering designation.....	198
Connection 21		Connector 55	
Connection of syringe.....	16	Connection of inlet capillary 90 ..	38,39
Figure.....	5	Figure.....	22
Connection 22		Installation.....	22
Adjust aspirating tubing.....	16	Connectors for capillaries.....	19
Figure.....	5	Contaminated valves.....	134
Connection 23		Contamination with heavy metals.....	140
Attachment of		Contamination with organic	
connection capillary 22	24	cationic complexing agents.....	140
Connection of filter unit		Contamination with	
Manufit 58	25	organic substances.....	140
Connection 27		control.....	95
Connection of syringe.....	16	Control characters.....	177
Figure.....	5	Control via RS232 interface.....	93
Installation of suppressor		Correction.....	79
inlet capillary.....	38	Correction factor.....	57,102
Connection 28		Counterpart end 58	
Adjust aspirating tubing.....	16	Figure.....	23
Figure.....	5	Installation.....	23
Installation of suppressor		Coupling 70	
inlet capillary.....	39	Figure.....	26
Connection 36		Installation.....	25
Connection of 732.....	14,15	Coupling 96	
Figure.....	7	Figure.....	37
Connection 42		Installation.....	38
Connection of drain tube.....	16	Coupling (6.2620.040).....	196
Figure.....	7	Coupling (6.2620.060).....	27,137,138,194
Connection 45		Coupling (6.2744.020).....	195
Connection of 732.....	15	Coupling (6.2744.030).....	197
Figure.....	7	Coupling (6.2744.040).....	195
Connection 48		Criterion for measured value printout..	90
Figure.....	7	Current information.....	170
Pin assignment.....	187	Current time.....	81
Connection capillary 22		cycle.....	83,107
Figure.....	25	Cycle program.....	91,107,110
Installation.....	24		
Connection of			
«IC Metrodata for Win95».....	41		
Connection of 709 IC Pump.....	21		

D

daily.....	115
Damping.....	56,99,189
dampng.....	99
data bit.....	93,94
Data entry.....	79
Data output.....	119
Data receipt.....	93
Data selection.....	79
Data transmission protocol.....	177
Data transmission rate.....	93
Date.....	49,52,82,90,91,92,115,119,122
date.....	92,115
date&time.....	90,91
Dead volume.....	148
Declaration of conformity.....	205
DEFAULT.....	82
Default value.....	77
Degassing of eluents.....	135
Degree of protection.....	10
Delete all events.....	116
Delete all program steps.....	113
delete all?.....	113,116
Delete event.....	116
Delete method.....	118
Delete parameter values.....	76
Delete program steps.....	109,110
delete?.....	118
Deleting characters.....	80
Detector block 81	
Connection at 732.....	13
Connection of	
suppressor module.....	35
Figure.....	34,36
Ordering designation.....	193
Detector block (1.732.0100).....	13,193
Detector block (1.732.0110).....	13,193
Determination of	
anions in drinking water.....	45
Device label.....	92
device label.....	92
Device name.....	90,119,122
Diagnostic steps.....	149
Diagnostic tests.....	149,175
dialog.....	92
Dialog language.....	53,92
Dialysis Unit 754	
Connection at IC system.....	2
Ordering designation.....	202
Dimensions.....	191,192
diabl ed.....	108
Display 1	
Description.....	71
Error messages.....	145
Figure.....	3
Representation.....	9
Standby mode.....	81
Test.....	152
Display of date and time.....	82
Display of the absolute conductivity ...	82
Display of the auto-zero value.....	89
Display of the full-scale range.....	82
Display of the suppressor position.....	96
Display of the valve setting.....	95,96
Displays in the standby mode.....	81

Disposal 10
 Door **20**
 Figure 5
 Double peaks 134,148
 Drain tube 16
 Drift 135,148,189

E

Earthing 10
 Edit event 116
 Edit mode 71,73
 Edit program 109
 Edit program points 110,115,116
 Edit program step 109
 Electrical safety 10
 Electromagnetic compatibility 191,192
 Elements of the
 measured value output 119,122
 Eluent change 135
 Eluents 135
 EMC 191,192
 Emitted interference 191
enabled 108
enabled in program 92
end 110
 <ENTER> 76,79,80
 ENTER 86
 Enter new event 115
 Enter new program step 109
 EOD 184
 Epson 90
 EPSON printers 42,90
 Error messages
 Beep 92
 Display 145
 RS232 interface 165,182
 Table of messages 145
 Error signal 184
 EU Declaration of conformity 205
 Evaporation 135,148
 <EVENT> 74,114
 event 115,116,124
 Event format 115
 Event report 124
 EVENT/delete 116
 EVENT/delete all 116
 EVENT/edit 116
 EVENT/new 115
 Events 114,169
 Examples 125,126
 Exit from rolling inquiries 76
 External devices 41,44,185,188
 External power supply for 733 187

F

Faults 133,145,149
 Feedthrough **22**
 Adjust aspirating tubing 16
 Figure 5
 Feedthrough **28**
 Adjust aspirating tubing 16
 Figure 5

Installation of suppressor
 inlet capillary 39
 Ferrule **49**
 Figure 20,29
 Installation 19
 Ordering designation 196
 <FILL> 76,103
fill 96,112
FILL A 85
FILL B/STEP 85
 Filling of sample loop 96,103
 Filter unit Manufit **69**
 Figure 23,26
 Installation 23
 Notes 133
 Ordering designation 197
 Filter unit PEEK **64**
 Connection at 752 Pump Unit 38,39
 Figure 22,25
 Installation 22
 Notes 133,136
 Ordering designation 195
 Filtration 62
flag 110
flow 56,101,112
 Flow chart
 Basic settings 68
 General 67
 Injection 70
 Installation 11
 Putting into operation 69
flow corr. 102
 Flow rate 56,101,112
format 115
 Frame 86,122
frame 86
FS change 110
 Full duplex 177
 Full hardware handshake 178
 <FULL SCALE> 73,102
full scale 40,82,98,99,102,111
 Full-scale range
 40,49,55,59,82,98,102,110
 Fundamentals of the
 instrument dialog 54
 Fundamentals of the operation 71
 Fuse 17,18,190
 Fuse holder **19**
 Figure 4,18
 Replacement of fuses 17

G

General flow chart 67
 General graphics parameters 86
 General information on
 maintenance 136
 General information on
 precolumns 28
 General information on
 separating columns 32
 General information on
 suppressor module 33
 General instrument settings 91
 General precautionary rules 10
 General program parameters 107
 Glass beads 134

Glass cartridge holder
 (6.2828.000) 201
 GLP 159
 Graphics parameters 86,122
 Graphics plot 86,100,122,168
grid 86
 Grid lines 86,122

H

Handling of solvents 10
 Handshake 93,94,178
handshake 93,94
 Hazard 9
 Heating 88
 Housing 191,192
 Housing **54** for filter unit
 Figure 22
 Installation 22
 HP printers 42,43,90

I

IBM 90,119
 IBM Proprinter 42,90
 IC anion column
 Metrosep Anion Dual 2 129,199
 IC anion column Organic Acids 199
 IC anion column PRP-X100 123,198
 IC anion column Star-Ion A300 199
 IC anion column SUPERSEP 200
 IC anion precolumn SUPERSEP
 Installation 31
 Ordering designation 200
 IC cation column
 Hyperrez Monovalent 200
 IC cation column
 Metrosep Cation 1-2 126,137,200
 IC cation column Nucleosil 5SA 200
 IC column cartridge
 METROSEP Anion Dual 1 199
 IC Detector 732
 Certificate of Conformity
 and System Validation 207
 Connection of 709 21
 Connection of 733 13
 Connection of 750 41
 Connection of 766 41
 Connection of a PC 44
 Connection of a printer 42
 EU Declaration of Conformity 205
 Front 3
 Key field 72
 Opening 10
 Parts and controls 3
 Rear 4
 Standard equipment 193
 Technical data 189
 IC exclusion column PRP-X300 198
 IC Metrodata for Win95
 Connection 41
 Ordering designation 202
 IC precolumn cartridge
 METROSEP Anion Dual 1 199
 IC precolumn cartridge
 METROSEP Cation 1-2 200
 IC precolumn cartridge
 Nucleosil 5SA 200

IC precolumn cartridge PRP-1.....	199		
IC precolumn cartridge PRP-X100	198		
IC precolumn cartridge PRP-X300	198		
IC Pump 709			
Automatic shutoff.....	101		
Connection at 732	21		
Electrical connection	21		
Maintenance.....	134		
Ordering designation.....	202		
Practical notes.....	134		
RS232 settings	94		
Settings.....	101,131		
Start/stop pump drive.....	106,112		
Status messages.....	82		
Switch on remote control.....	48		
IC Sample Processor 766			
Connection	41		
Ordering designation.....	202		
IC separating column 76			
Changing	138		
Figure	29,34,36		
Install cartridge head.....	30		
Installation.....	35,37		
IC separating columns			
.....	see <i>Separating columns</i>		
IC Separation Center 733			
Certificate of Conformity			
and System Validation.....	208		
Configuration	95		
Connection of 732	13		
Connection of 750	41		
Connection of 766	41		
EU Declaration of conformity.....	206		
External power supply	187		
Front	5		
Interior.....	34,36		
Key field	72		
Opening.....	10		
Operating keys	76		
Optional accessories.....	196		
Parts and controls	5		
Rear	6		
Settings.....	95		
Standard equipment.....	194		
Technical data	192		
Valve interfaces	188		
i d. 1	89		
i d. 2	90		
Identical separation system.....	138		
Identification	89,122		
Immunity to interference	191		
i n p o s i t i o n	96		
Increase in the baseline.....	140		
Individual measured values.....	120		
Information on the			
Instructions for Use.....	8		
Initial condition.....	48		
Initialize data memory.....	157		
<INJECT>	76,104		
i n j e c t	83,96,108,112		
INJECT A	85		
INJECT B	85		
Inject program	108		
Injection	96,104		
Injection valve 68			
Connection of 709	24,26		
Figure	25,26,34,36		
Operation.....	76		
Settings.....	95,96		
Switching	103,104,111		
Inlet capillary 56			
Figure	23,26		
Installation	23		
Inlet capillary 75			
Figure	29		
Installation	28		
Inlet capillary 82			
Figure	34,36		
Installation	35		
Passivation	27,138		
Inlet capillary 87			
Figure	34,36		
Inquiries.....	54,77		
Installation	11,46		
Installing precolumn cartridges.....	29		
Installing the Manufit filter unit.....	23		
Installing the PEEK Filter unit.....	22		
Instructions for Use			
(8.732.1033)	8,193		
Instrument adjustment.....	176		
Instrument description.....	1		
Instrument dialog.....	54,77		
Instrument information.....	170		
Instrument number	173		
Instrument settings.....	91		
Integrator start	188		
Interfaces.....	161		
Interior of the 2.733.0010			
IC Separation Center	34		
Interior of the 2.733.0X30			
IC Separation Center	36		
Internal hardware test.....	155		
Introduction	1		
Ion chromatogram			
of the calibration.....	61		
Ion chromatogram of the			
drinking water sample	63		
Ion chromatography.....	133		
Ion chromatography system	2		
K			
Key 23 "FILL"			
Figure	5		
Function.....	103		
Key 24 "INJECT"			
Figure	5		
Function.....	104		
Key 25 "FILL"			
Figure	5		
Function.....	103		
Key 26 "INJECT"			
Figure	5		
Function.....	104		
Key code	152,172		
Key field.....	72		
Key name	71		
Key table	152		
Keypad	190,192		
Keypad test	151		
Keys			
<←>	80,102		
<→>	80,102		
<CLEAR>	76,80		
<CONFIG>	74,88		
<ENTER>	76,79,80		
<EVENT>	74,114		
<FILL>	76,103		
<FULL SCALE>	73,102		
<INJECT>	76,104		
<MARK>	75,106		
<METHOD>	74,117		
<PARAM>	73,97		
<PLOT>	75,122		
<PRINT>	75,119		
<PROG R/S>	73,114		
<PROGRAM>	74,107		
<PUMP R/S>	74,106		
<QUIT>	76,79,80		
<REPORT>	75,124		
<SELECT>	76,79		
<ZERO OFF>	74,106		
<ZERO>	73,105		
Knurled screw 34			
Figure	7		
Open rear panel	14		
Open in rear panel	13		
Knurled screw 37			
Figure	7		
Open rear panel	24,26		
L			
Labeling of the time axis	100		
Labeling of the x axis.....	122		
Labeling of the y axis	122		
Leak testing.....	40		
Leaks.....	10,101,148		
l e f t	100		
Left boundary	100		
Linearity.....	98,189		
Location	12		
Lock functions.....	172		
Loop program	91,107,110		
Loss of suppressor capacity.....	143		
M			
Magnetic plate (6.2248.000)			
Installation	12		
Ordering designation	193		
Main function keys 2			
Figure	3		
Overview.....	72,73		
Main menu	54,77		
Mains cable			
Connection.....	18		
Ordering designation	193		
Mains connection			
Procedure.....	17,18		
Safety notes.....	10		
Technical data.....	190		
Mains connection plug 18			
Connection of mains cable	18		
Figure	4,18		
Mains frequency.....	190		
Mains switch 9			
Figure	4,18		
On/off switching of			
the instruments.....	18		
Mains voltage			
Setting	17		
Technical data.....	190		
Maintenance.....	133,136		
Malfunctions.....	145,148		
Manufit filter unit.....	see <i>Filter unit Manufit</i>		
Manufit housing 62			

Figure 23,29
 Installation 23,28

Manufit housing **80**
 Figure 29
 Montieren 30

Manufit pressure screw **57**
 Figure 23,29
 Installation 23,28

Manufit pressure screw **77**
 Figure 29
 Installation 30

<MARK> 75,106

MARK 85,111

Marking signal 106,111,170,189

Maximum back pressure 190

Maximum shutoff pressure 56,101,113

Measured value 71,81,119

Measured value output 119,122

Measured value printout 90,166

Measurement parameters 97,117,168

Measuring cell temperature 88

Measuring range 40,54,97,110,111

Menu title 71

Method 50,82,107,117,168

<METHOD> 74,117

method 124

Method examples 126

Method list 124

method list 124

Method name 117,118

Method report 124

METHOD/delete 118

METHOD/recall 117

METHOD/store 118

Metrohm Monograph
 "Ion chromatography" 8,133

Metrohm service 136

Metrohm suppressor
 module *see Suppressor module*

Microfiltration 133,134,135

Minimum shutoff pressure 57,101,113

Model plate **47**
 Figure 7

modif. 82

Monograph
 "Ion chromatography" 8,133

Mounting rail **83**
 Figure 34,36
 Installation of column holder 35

MSM *see Suppressor module*

N

name 117,118

next step 83

Noise 99,189

Notation 9

Notes 133

Number of cycles 91,108,110

number of cycles 91,107,108

Number of the program version 92

Numeric entry 79

Numeric keys **3**
 Figure 3
 Numeric entry 80
 Overview 72,74

O

Object 161,166

Obstruction 141

offset 55,99

On/off switching of the instruments 18

One-channel system 2,35,37

Open configuration menu 74

Open event menu 74

Open method menu 74

Open parameter menu 73

Open program menu 74

Open-end spanner (6.2621.050) .. 19,195

Open-end spanner (6.2621.060) 195

Opening **29**
 Figure 7
 Insertion of detector cable 14

Opening **30**
 Figure 7
 Insertion of outlet capillary 14

Opening **32**
 Figure 7
 Insertion of outlet capillary 14

Opening **33**
 Figure 7
 Insertion of detector cable 14

Opening **40**
 Figure 7
 Insert capillary 24,26

Opening **44**
 Figure 7
 Insert capillary 24,26

Opening the text editor 80

Operating mode 172

Operating range 98,102

Operating sequences 67

Operating temperature .. 40,50,88,99,190

Operating tutorial 45

Operation 67

operation with 709 87

Operation with 709 IC Pump 87

operation with 733 87

Operation with
 733 IC Separation Center 87

Optional accessories 196

Organic modifiers 135

Outlet capillary **63**
 Connection 25
 Figure 23,26
 Installation 23

Outlet capillary **72**
 Figure 29
 Installation 28,35,37

Output of individual
 measured values 120

OVERLOAD 120

Overload display **8**
 Figure 3
 Function 98,105

Overview of key functions 72

overwrite? 118

P

Packaging 12

Paper feed 100

Paper speed 100

<PARAM> 73,97

param 124

PARAM/709 IC Pump 101

PARAM/analog output 99

PARAM/detector 97

PARAM/plot 100

Parameter name 71

Parameter report 124

Parameter settings 97

Parameter value 71

Parameters for 709 IC Pump 56

Parameters for analog output 55

Parameters for the
 conductivity detector 54

Parity 93,94

pari ty 93,94

Parts and controls 3,190

Passivation 27,137,138

Path 172

PC 119,161

PCL3 emulation 42

Peak broadening 134,148

PEEK capillary (6.1831.010) 19,20,194

PEEK connectors 20

Perform events in a program 92

Peripherals 87

pH 135

Pictograms 9

Pin assignment 181

Plate number 133

Plot *see Graphics plot*

<PLOT> 75,122

PL0T 86

Pmax 82,83,101,113

Pmi n 82,83,101,113

Polarity 55,86,98,99,111,187

pol ari ty 86,98,99,111

Power consumption 190

Practical notes 133

Precautionary rules 10

Precipitates 134,135

Precolumn cartridge **74**
 Figure 29
 Installation 28,30

Precolumn cartridge
 holder (6.2828.010) 201

Precolumns
 Installation 28
 Notes 133
 Ordering designations 198

Precolumns with cartridge head 30

Precolumns with
 twin cartridge holder 28

Preconcentration
 cartridge METROSEP Anion 199

Preparations for operating tutorial 46

Prep- MSM 40,117

Pressure drop 148

Pressure rise 148

Pressure screw **50**
 Figure 20,29
 Installation 19
 Ordering designation 196

<PRINT> 75,119

print crit. 90,121

Print header 89,119,122

print header	90	Example.....	131	Remote program.....	108
Printer	42,89,94,119,122,124,161,182	General information.....	33	repeat time.....	93
Printer drivers.....	42	Ordering designation.....	202	Repeat time for beep.....	93
Printer settings.....	42	Preparation.....	38	Replacing the suppressor.....	143
Printer type	90,122	Putting into operation	48,69	<REPORT>	75,124
Printout of date and time	90,91	Q		REPORT	85,124
Printout of the print header	90	Quick reference guide		report	124
<PROG R/S>	73,114	(8.732.1043)	8,193	Report identification	172
PROG R/S	85,115,116	<QUIT>	76,79,80	Report selection	124
prog. type	83	QUIT	86	reset.....	110
Program.....	107,117,169	R		Reset the parameters.....	110
<PROGRAM>	74,107	RAM test.....	151	Resolution	133,148
program.....	92,124	range.....	40,97,111	return.....	107,110,119
Program disabled.....	108,114	Ready signal	184	Return with loop programs.....	110
Program enabled.....	108,114	Rear panel 35		right.....	100
Program end.....	110	Figure	7	Right boundary.....	100
program end.....	108,109	Open rear panel	13,14	Rinse with dist. H ₂ O.....	27,139
Program flag.....	110	Rear panel 39		Rinse with eluent	27,139
Program number	90,119,122	Figure	7	Rinse with HNO ₃	27,138
Program parameters	107	Open rear panel	24,26	Rinsing of IC system	137
Program point.....	107,110,114,115,116	Rear panel opening 31		Rinsing solution.....	33,39
Program report	124	Figure	7	Rolling inquiries.....	77
Program running.....	114	Open rear panel	24,26	RS control	93
Program start.....	73	Rear panel opening 41		RS232 error rectification.....	182
Program status display 5		Figure	7	RS232 interface 15	
Figure	3	Rear panel opening 43		Connection of a printer.....	42
Function	108	Figure	7,36	Data transmission protocol	177
Program step.....	83,107,109,113	Insert suppressor		Description	161
Program stop.....	73	outlet capillary	38,39	Error rectification	182
Program test	108,109	Recall existing method	117	Figure	4
Program type.....	107	Recorder.....	41	Handshake	178
Program version	92	Rectification of malfunctions	148	Pin assignment.....	181
PROGRAM/del et e all	113	Recycling.....	135,137	Settings	94
PROGRAM/edi t	109	Reduced hardware handshake	178	Test.....	153
PROGRAM/parameters	107	Reference temperature.....	99,189	RS232 interface 16	
Programming.....	107	Reference value for the		Connection of a PC	44
Protection class.....	10,18	auto-zero function.....	120	Connection of a printer.....	42
Protective earth.....	18	Regeneration	134,140	Data transmission protocol	177
PTFE aspirating tubing 88		Regeneration solution.....	33,38,39,140	Description	161
Adjusting.....	16	Relative time	100	Error rectification	182
Figure	34,36	Relative width of the		Figure	4
PTFE gasket 59		graphics printout	86	Handshake	178
Figure	23,29	remote.....	83,108,112	Pin assignment.....	181
Installation.....	23,28,30	Remote control commands.....	166	Remote control commands.....	166
Ordering designation.....	197	Remote control for 709 IC Pump.....	48	Settings	93
PTFE microcapillary (6.1822.010)	20	Remote control language.....	161	Test.....	153
Pulsation dampener 66		Remote control tree.....	166	Rules for remote	
Connection	21	Remote control via RS232.....	161	control language	162
Figure	25,26,34,36	Remote input lines		Run number	91,119
Ordering designation.....	196	Assignment.....	85,173	run number.....	91,119
Practical notes	134	Functions.....	183	Run time	119
Pulsations	99,148	Status	170	S	
Pump.....	see IC Pump 709	Remote interface 14		Safety notes	10
remote.....	124	Description	185	Safety specifications	191,192
Pump drive	112	Pin assignment.....	185	Sample determination	62
pump not respondi ng.....	82,83	Test.....	154	Sample loop 86	
<PUMP R/S>	74,106	Remote interface 17		Figure	34,36
PUMP R/S	85,112	Connection of external devices	44	Filling.....	60,62,95,96,103
pump ready	82	Description	183	Ordering designation	196
Pump report.....	124	Figure	4	Selection.....	32
pump runni ng	82	Pin assignment.....	183	Sample pretreatment cartridge	
Pump tubing (6.1826.050)		Test.....	154	Chromafix C18	201
General information.....	33	Remote output lines		Sample pretreatment cartridge	
Installation.....	38	Basic settings	84,174	IC-Ag.....	201
Pump Unit 752		Functions.....	184	Sample pretreatment cartridge	
Connection of aspirating tubing .	38,39	Settings	112	IC-H.....	200

Sample pretreatment cartridge IC-OH 201	Sign change 75	Suppressor holder 100 Figure 142,143 Mounting 144
Sample pretreatment cartridge IC-RP 200	Silicone tubing (6.1816.000) 194	Suppressor inlet capillary 89 Figure 36,38 Installation 37
Screw nut 97 Figure 142 Mounting 144	Single-channel system 2,35,37	Suppressor inlet capillary 90 Figure 36,38 Installation 38
Seiko printer DPU-411/414 42,90	Software handshake with character stop 179	Suppressor inlet capillary 93 Figure 37,38 Installation 39
<SELECT> 76,79	Software handshake with line stop 180	Suppressor module 95 Cleaning 141 Conditioning 40 Connection 1 37 Connection 2 38 Connection 3 39 Connection of column 37 Display of position 96 Figure 37 General information 33 Installation 37 Protection 136 Regeneration 140 Replacement 143 Rotation 76 Settings 96 Switching 96,103,112
SELECT 86	SOP 159	Suppressor outlet capillary 91 Figure 37,38 Installation 38,39
Select event 116	Spilled chemicals 136	Suppressor outlet capillary 92 Figure 37,38 Installation 39
Select preset parameter values 76	Splitting 134,148	Suppressor outlet capillary 94 Figure 37,38 Installation 37
Select program steps 109	Standard chromatogram 133	Suppressor rotor 99 Cleaning 143 Figure 142 Inserting 143 Ordering designation 201 Replacement 143
Select status message 76	Standard equipment 193	Switch off auto-zero 74,106
Select program steps 109	Standard operating procedures 159	Switch off data receipt 93
Select status message 76	Standard solution 59	Switch on data receipt 93
Selection of data 79	Standby mode 71,73,81	Switching of the suppressor module 96,103,104,112
Selection of the character set/printer type 90	Start graphics plot 75	Switching of the valve 95,96
Selection of the method 117,118	Start measured value output 75	Switching to "FILL" 76
Selection of the report 124	Start of a program 114	Switching to "INJECT" 76
Selection of the sample loop 32	Start report output 75	Syringe (6.2816.020) Connection of syringe 16 Ordering designation 195
Selectivity 133	Start the graphics plot 122	
send to 90	Start/stop program 73,115	
Sensitivity 98	Start/stop pump drive 74,106,112	
Separating column Changing 138	Static charges 10	
Characteristic data 133	status 92,95,96,108	
Dead volume 134	Status messages 49,71,82,164	
Fixing 35,37	Status messages for 709 IC Pump 82	
General information 32	Status messages for method 82	
Installation 35,37	Status messages for program step 83	
Ordering designations 198	Status of the beeper 92	
Practical notes 133	Steel capillary 71 Connection 25 Figure 26	
Protection 133	Steel capillary (6.2620.020) 19,194	
Regeneration 134	Steel connector 79 Figure 29 Installation 30	
Rinsing 35,37	Steel connectors 19	
Separation efficiency 133	Steel mesh holding end 61 Figure 23 Installation 23	
Storage 133	Steel meshes 60 Figure 23,29 Installation 23 Ordering designation 197	
Separation efficiency 133	Steel meshes 73 Figure 29 Installation 28,30 Ordering designation 197	
Serial number 90,119,122,149	Steel spacer 78 Figure 29 Installation 30 Ordering designation 197	
Serial number 10 Figure 4	step 96,112	
Serial/parallel converter (2.145.0300) 43,44,202	Stirring of the eluent 135	
Service 136,159	stop bit 93,94	
Servicing 136	Stop of a program 73,114	
Setting up the instruments 12	stop time 91,100,121	
Settings 709 IC Pump 101	Stop time for graphics plot 100	
733 IC Separation Center 95	Stop time for measured value printout 91	
Analog output 99	stopped 82,83,101	
Beeper 92	Store method 64,118	
Conductivity detector 88,97	Submenu 54,71,77	
Events 92	Sunlight 12	
External printer 89	Suppression 33	
Graphics plot 100	suppressor 112	
Injection valve A 95		
Injection valve B 96		
Mains voltage 17		
Measured value printout 90		
Peripherals 87		
Remote output lines 84		
RS232 interface 93		
Suppressor module 96		
Setup 84		
SETUP/graphi cs 86		
SETUP/i nput assi gn 85		
SETUP/out put 84		
SETUP/per i pheral s 87		
Shutdown 137		
Shutoff limits 82,83		

T

Technical data 189
temp. coeff. 99
Temperature coefficient 55,99,189
Temperature correction 189
Temperature dependence 99,189
Temperature deviation 190
Temperature of conductivity cell 88
Temperature stability 88,190
Terminal block 38 Figure 7 Pin assignment 188

Terminal block 46	
Figure	7
Pin assignment	188
test	108,109
Test chromatogram	32
Test mode.....	108,109
Test plug (3.496.8480)	149,153
Test plug (3.496.8510)	149,154
Test status	108,109
test step	108,109
Text cursor	80
Text editor	80
Text entry	75,80
thermostat	40,88
Thermostat display 7	
Figure	3
Function	88,89,99
Thermostating.....	88,99,190
Time.....	49,52,82,90,91,92,115,119,122
time	92,109
Time axis	100
Time interval	100,121
Time interval for measured value printout	91
time interval	91,100
Time program.....	107
Time restriction	91
Time scale	100
time scale label	100
time scale	100
Transport	12
Transport damage.....	204
Treatment of eluents.....	135
Trigger	162,163
trigger	95,96
Trigger a marking signal.....	75,106
Trigger auto-zero	73
Triggering of functions.....	103,183
Triggering the print command.....	183
Tubing cartridge	38
Twin cartridge holder (6.2821.050)	
Figure	29
Install precolumn cartridge.....	28
Ordering designation.....	197
Two-channel system.....	2,35
type.....	107

U

Unit for display of the auto-zero value	89
--	----

V

VA Detector 791	
Connection	41
Ordering designation.....	202
Validation.....	159
Valve.....	<i>see Injection valve</i>
val ve A.....	95,111
val ve A + suppressor.....	95
val ve A + val ve B.....	95
val ve B.....	112
Valve control	95
Valve interfaces	188
Valve setting	95

W

Warning	9
Warranty	204
Waste container.....	14,15
width.....	86
Width of the graphics printout	86

Z

<ZERO>	73,105
ZERO.....	85,111
<ZERO OFF>.....	74,106
ZERO OFF	86
Zero point offset	99
'zero' unit	89
ZEROREF	120