CH-9101 Herisau/Schweiz
Telefon +41 71 353 85 85
Fax +41 71 353 89 01
EMail sales@metrohm.ch
Internet <http://www.metrohm.ch>

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8.748.1001

99.01 mst, ti

CH-9101 Herisau/Schweiz

Telefon +41 71 353 85 85

Fax +41 71 353 89 01

E-Mail sales@metrohm.ch

Internet <http://www.metrohm.ch>

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DH= Dipping Head

1 Einführung

1.1 Übersicht

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- 1** **Kabel-/Schlauchhalter 6.2053.010**
(3 Stück)

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3 Operation via RS232 Interface

3.1 General rules

The DH Sample Changer has a remote control facility that allows its control via TiNet. The parameters of the Sample Changer may be easily set by the method editor of TiNet.

The following description offers special possibilities of low level control via the transfer command of the calculation block of the TiNet method editor.

3.1.1 Tree of commands

The commands are hierarchically ordered with the following rules:

Rules	Example
The root of the tree is designated by &.	
The branches (levels) of a tree are marked with a dot (.) when calling up an object.	
Upper- or lowercase letters may be used.	&STIRER.ON or &Stirer.On
An object can be assigned a value. Values are marked at the beginning and end by quotes ("). They may contain ASCII characters or numerical values.	&Stirer.On"True" &Init.SpeedXY"3"
The current object remains valid until a new object is called.	
New objects can be addressed relative to the old object: A preceding dot leads forwards to the next level in the tree.	
More than one preceding dot leads one (or several) level(s) backwards in the tree. n node backwards require n+1 preceding dots.	&Stirer.On"True";..\$G Leads back from nod &Stirer.On to nod &Stirer
If you want to jump back to the root, enter a preceding &.	

3.1.2 Triggers

\$G **Go** Commands are executed with \$G.

Example:

&Stirer.On"True" Sets command.

..\$G Starts stirrer.

3.2 Remote control commands

3.2.1 Overview

Object	Description	Input range	Reference
& Root			
Stirer	Stirrer control	\$G	3.2.2.1.
├ .On	Switching stirrer	True, False	ditto
GoTo	Go to position	\$G	3.2.2.2.
├ .Vessel	Designation of vessel	special	ditto
GoAbove	Go to position	\$G	3.2.2.3.
├ .Vessel	Designation of vessel	special	ditto
GoIn	Go into vessel	\$G	3.2.2.4.
├ .Vessel	Designation of vessel	special	ditto
GoToHome	Goes to home position	\$G	3.2.2.5.
WashElectrodes	Wash electrodes	\$G	3.2.2.6.
├ .Times	Number of wash cycles	1...9999	ditto
SetCtrlLine	Setting I/O line	\$G	3.2.2.7.
├ .On	Setting line	True, False	ditto
├ .Line	Selection of line	1...3	ditto
WaitForInput	Waiting for input on I/O line	\$G	3.2.2.8.
├ .On	Line state	True , False	ditto
├ .Line	Selection of line	1...3	ditto
Init	Setup of Sample Changer	\$G	3.2.2.9.
├ .SpeedXY	Speed in XY direction	0... 5 ...7	ditto
├ .SpeedXYFast	Speed in XY direction	0... 2 ...7	ditto
├ .SpeedZ	Speed in Z direction	0... 3 ...7	ditto
├ .SpeedZFast	Speed in Z direction	0 ...7	ditto
├ .SpeedZSlow	Speed in Z direction	0... 7	ditto
StopAll	Emergency stop	\$G	3.2.2.10.
Assembly	Basic assembly functions		
├ .GoToXYZ	Go to position	\$G	3.2.2.11.
├ .X	X coordinate	0 ...630	ditto
├ .Y	Y coordinate	0 ...272	ditto
├ .Z	Z coordinate	0 ...120	ditto
├ .SpeedXY	Speed in XY direction	0... 3 ...7	ditto
├ .SpeedZ	Speed in Z direction	0... 3 ...7	ditto
├ .PumpOn	Switching on/off pump	True, False , \$G	3.2.2.12.

3.2.2 Description of the remote control commands

3.2.2.1 Stirer \$G
 Stirer.On True, **False**

&Stirer.On"True";..\$G Switching the stirrer on.
 &Stirer.On"False";..\$G Switching the stirrer off.

3.2.2.2 GoTo \$G
 GoTo.Vessel 1...XXX, Home, Soak, Wash

Makro command which moves to the predefined beaker using the speeds from branch "Init" (see 3.2.2.9):

- Moves to shift position (set in TiNet configuration), SpeedZSlow
- Moves to beaker position, SpeedXY
- Moves to Z = 0 (as defined in the rack definition file ..\bin\PW74xDef\748*r.def of TiNet) position.

3.2.2.3 GoAbove \$G
 GoAbove.Vessel 1...XXX, Home, Soak, Wash

Makro command which moves to the predefined beaker (without dipping) using the speeds from branch "Init" (see 3.2.2.9):

- Moves to shift position (set in TiNet configuration), SpeedZSlow
- Moves to beaker position, SpeedXY

3.2.2.4 GoIn \$G
 GoIn.Vessel 1...XXX, Home, Soak, Wash

Makro command which moves to the lowest possible position in the predefined beaker using the speeds from branch "Init" (see 3.2.2.10):

- Moves to shift position (set in TiNet configuration), SpeedZSlow
- Moves to beaker position, SpeedXY
- Moves to lowest possible Z position, Zmax. (defined in the corresponding rack definition file), SpeedZ

3.2.2.5 GoToHome \$G

Makro command which moves to home position, readjusting the Sample Changer and using the speeds from branch "Init" (see 3.2.2.9):

3.2.2.6 WashElectrodes \$G
 WashElectrodes.Times 1...9999

Makro command which moves to wash position position, washes the electrodes and uses the speeds from branch "Init" (see 3.2.2.9):

- Moves up to Z=0, SpeedZ
- Moves to wash position, SpeedXY
- Moves down to wash position, SpeedZ
- Switches pump on
- Moves down to lowest possible position Zmax, SpeedZSlow
- Switches pump off

- Stirrs during 3 s
 - Moves up to Z=0, SpeedZFast
- .Times: The number of washing cycles can be set. If the number is >1, the following additional steps will be carried out (n times); the repetitive steps beginning at the marked position ➤ in the sequence:
- Moves up to Z=0, SpeedZFast
 - Switches pump on
 - Moves down to Zmax., Speed Zslow
 - Switches pump off

3.2.2.7 SetCtrlLine \$G
 SetCtrlLine.On True, False
 SetCtrlLine.Line 1...3

A remote control line can be set. Select socket 1...3 and set line as follows:

&SetCtrlLine.On"True";..\$G Sets an active line (0 V), see also page 20.

&SetCtrlLine.On"False";..\$G Sets the line to inactive (5 V)

The beaker sensor is normally plugged into socket 1.

3.2.2.8 WaitForInput \$G
 WaitForInput.On True, False
 WaitForInput.Line 1...3

The Sample Changer waits for an active line (with WaitForInput.ON"True") or for an inactive line (with WaitForInput.ON"False"). No other commands are worked off during this waiting time.

3.2.2.9 Init \$G
 Init.SpeedXY 0...5...7
 Init.SpeedXYFast 0...2...7
 Init.SpeedZ 0...3...7
 Init.SpeedZFast 0...7
 Init.SpeedZSlow 0...7

Sets the global speeds for all makro commands of the Sample Changer. A small number means a high speed, i.e. 0 is the highest possible speed.

The speeds are only set if the init movement is carried out (command &Init \$G).

&Init \$G moves to home position and Z=0.

3.2.2.10 StopAll \$G
 Emergency stop. Stops all commands immediatelay. May be used e.g. to set remote control lines to inactive, or to stop stirrer or pump.

3.2.2.11	Assembly.GoToXYZ	\$G
	Assembly.GoToXYZ.X	0...630
	Assembly.GoToXYZ.Y	0...272
	Assembly.GoToXYZ.Z	0...120
	Assembly.GoToXYZ.SpeedXY	0...3...7
	Assembly.GoToXYZ.SpeedZ	0...3...7

Moves to the predefined position. The above speeds are only valid for Assembly commands.

The extreme positions which are still in the moving range of the Sample Changer are given in the corresponding rack definition file (`..\bin\PW74xDef\748*r.def` of your TiNet program).

You may enter higher values for X, Y, and Z but the Sample Changer has to be reinitialized after being moved to a position out of its moving range.

3.2.2.12	Assembly.PumpOn	True, False , \$G
	&Assembly.Pump"True";\$G	Switching the pump on.
	&Assembly.Pump"False";\$G	Switching the pump off.

3.3 Properties of the RS 232 Interface

The DH Sample Changer is configured as DTE (Data Terminal Equipment).

The RS 232 interface has the following technical specifications:

- Data interface according to the RS 232C standard
- RS232 parameters: 9600 baud
8 data bit
1 stop bit
no parity
SW handshake
- Control characters: C_R (ASCII DEC 13)
L_F (ASCII DEC 10)
XON (ASCII DEC 17)
XOFF (ASCII DEC 19)
- Cables 6.2125.060 25/25 pol
6.2125.110 25/9 pol
length: max. approx. 15 m

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Die Firma Metrohm AG, Herisau, Schweiz bescheinigt hiermit, dass das Gerät: