CH-9101 Herisau/Schweiz

Tel. +41 71 353 85 85

Fax +41 71 353 89 01

CompuServe 100031,3703

Internet <http://www.metrohm.com>

751 GPD Titrino

ab Serie 9

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8.751.1101

99.07 Ti

Inhalt

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3 Operation via RS232 Interface

3.1 General rules

The Titrino has an extensive remote control facility that allows full control of the Titrino via the RS 232 interface, i.e. the Titrino can receive data from an external controller or send data to an external controller. C_R and L_F are used as terminators for the data transfer. The Titrino sends $2xC_R$ and L_F as termination of a data block, to differentiate between a data line which has C_R and L_F as terminators. The controller terminates its commands with C_R and L_F . If more than one command per line is sent by the controller, “;” is used as a separator between the individual commands.

The data are grouped logically and easy to understand. Thus e.g., for the selection of the dialog language, the following must be sent

&Config.Aux.Language "english"

whereby it is sufficient to only transmit the boldface characters, thus:

&C.A.L "english"

The quantities of the commands above are:

Config	configuration data
Aux	auxiliaries, various data
Language	setting the dialog language

The data are hierarchically structured (tree form). The quantities that occur in this tree are called **objects** in the following. The dialog language is an object which can be called up with the

&Config.Aux.Language

command.

If one is in the desired location in the tree, the value of the object can be queried.

&Config.Aux.Language \$Q Q means Query


The query command \$Q initiates the issuing of the value on the instrument and the value emission is triggered. Entries which start with \$, trigger something. They are thus called **triggers**.

Values of objects can not only be queried, they can also be modified. Values are always entered in quotes, for example:

&Config.Aux.Language "english"

3.1.1 Call up of objects

An excerpt from the object tree is represented below:

Rules	Example
The root of the tree is designated by &.	
The branches (levels) of a tree are marked with a dot (.) when calling up an object.	
When calling up an object, it is sufficient to give only as many letters as necessary to uniquely assign the object. If the call is not unequivocal, the first object in the series will be recognized.	Calling up the dialog language &Config.Aux.Language or &C.A.L
Upper- or lowercase letters may be used.	&C.A.L or &c.a.l
An object can be assigned a value. Values are signified at the beginning and end by quotes ("). They may contain up to 24 ASCII characters. Numerical values can contain up to 6 digits, a negative sign, and a decimal point. Numbers with more than 6 characters are not accepted; more than 4 decimal places are rounded off. For numbers <1, it is necessary to enter leading zeros.	Entering the dialog language: &C.A.L"english" correct entry of numbers: "0.1" incorrect entry of numbers "1,5" or "+3" or ".1"
The current object remains until a new object is called.	entry of another dialog language: "deutsch"
New objects can be addressed relative to the old object: A preceding dot leads forwards to the next level in the tree.	From the root to node 'Aux': &C.A Forward from node 'Aux' to 'Prog': .P
More than one preceding dot leads one level backwards in the tree. n node backwards require n+1 preceding dots.	Jump from node 'Prog' to node 'Aux' and select a new object 'Language' at this level: ..L
If you must jump back to the root, enter a preceding &.	Change from node 'Language' via the root to node 'Mode': &M

3.1.2 Triggers

Triggers initiate an action on the Titrino, for example, starting a process or sending data. Triggers are marked by the introductory symbol \$.

The following triggers are possible:

\$G	Go	Starts processes, for ex. starting the mode run or setting the RS 232 interface parameters
\$S	Stop	Stops processes
\$H	Hold	Holds processes
\$C	Continue	Continues processes after Hold
\$Q	Query	Queries all information from the current node in the tree forward up to and including the values
\$Q.P	Path	Queries the path from the root of the tree up to the current node
\$Q.H	Highest Index	Queries the number of son nodes of the current node
\$Q.N"i"	Name	Queries the name of the son node with index i, $i = 1 - n$
\$D	Detail-Info	Queries the detailed status information
\$U	qUit	Aborts the data flow of the instrument, for example, after \$Q

The triggers \$G and \$S are linked to particular objects, see the summary table page 116ff.

All other triggers can be used at any time and at all locations on the object tree.

Examples:

Querying the value of the baud rate: **&Config.RSSet.Baud \$Q**

Querying all values of the node "RSSet": **&Config.RSSet \$Q**

Querying the path of the node "RSSet": **&Config.RSSet \$Q.P**

Start mode: **&Mode \$G**

Querying the detailed status: **\$D**

3.1.3 Status messages

In order to have an efficient control by an external control device, it must also be possible to query status conditions; they provide information on the status of the Titrino. The trigger \$D initiates output of the status. Status messages consist of the global status, the detailed status and eventual error messages, e.g. \$S.Mode.SET;E26. The global status informs on the activity of the process, while the detailed status conditions show the exact activity within the process.

The following **global status conditions** are possible:

\$G	Go:	The Titrino is executing the last command.
\$H	Hold:	The Titrino has been held (\$H, key <meas/hold> or by an error which effects the hold status)
\$C	Continue:	The Titrino has been restarted actively after hold
\$R	Ready:	The Titrino has executed the last command and is ready
\$S	Stop:	A process has been aborted in an "unnatural manner". e.g. stopped or aborted because there was an error.

Detailed status conditions

Status conditions of the global \$G:

\$G	.Mode.DET	.Inac:	Instrument at the beginning or at the end of a titration.
		.Req .Id1:	Instrument in the DET mode, requesting Id1 after titration start.
		.Id2:	Instrument in the DET mode, requesting Id2 after titration start.
		.Id3:	Instrument in the DET mode, requesting Id3 after titration start.
		.Smp1:	Instrument in the DET mode, requesting sample size after titration start.
		.Unit:	Instrument in the DET mode, requesting unit of sample size after titration start.
		.Start:	Instrument in the DET mode, processing the start conditions.
		.Titr:	Instrument in the DET mode, titrating.
\$G	.Mode.MET...		As DET.
\$G	.Mode.SET	.Inac:	Instrument at the beginning or at the end of a titration.
		.Req .Id1:	Instrument in the SET mode, requesting Id1 after start.
		.Id2:	Instrument in the SET mode, requesting Id2 after start.
		.Id3:	Instrument in the SET mode, requesting Id3 after start.
		.Smp1:	Instrument in the SET mode, requesting sample size after start.
		.Unit:	Instrument in the SET mode, requesting unit of sample size after start.
		.Start:	Instrument in the SET mode, processing the start conditions.
		.SET1:	Instrument in the SET mode, titrating to the first endpoint.
		.SET2:	Instrument in the SET mode, titrating to the second endpoint.
		.Cond.Ok:	Instrument in the SET, conditioning, endpoint reached (after the first startup from the standby mode).
		.Cond.Prog:	Instrument in the SET mode, conditioning, endpoint not reached (Conditioning progressing).
\$G	.Mode.KFT...		As SET.
\$G	.Mode.STAT	.Inac:	Instrument at the beginning or at the end of a determination.
		.Req .Id1:	Instrument in the STAT mode, requesting Id1 after start.
		.Id2:	Instrument in the STAT mode, requesting Id2 after start.
		.Id3:	Instrument in the STAT mode, requesting Id3 after start.
		.Smp1:	Instrument in the STAT mode, requesting sample size after start.
		.Unit:	Instrument in the STAT mode, requesting unit of sample size after start.
		.Start:	Instrument in the STAT mode, processing the start conditions.
		.Titr:	Instrument in the STAT mode, controlling.
		.Wait:	Instrument in the STAT mode, waiting because a monitored limit has been violated.
\$G	.Mode.DOS...		As STAT.

\$G .Mode.DOC...	As STAT.
\$G .Mode.MEAS .Inac:	Instrument at the beginning or at the end of a titration.
.Req .Id1:	Instrument in the MEAS mode, requesting Id1 after start.
.Id2:	Instrument in the MEAS mode, requesting Id2 after start.
.Id3:	Instrument in the MEAS mode, requesting Id3 after start.
.Smp1:	Instrument in the MEAS mode, requesting sample size after start.
.Unit:	Instrument in the MEAS mode, requesting unit of sample size after start.
.Meas:	Instrument in the MEAS mode, measuring.
\$G .Mode.CAL .Inac:	Instrument at the beginning or at the end of a calibration
.Req.Temp:	Instrument in the CAL mode, requesting calibration temperature.
.Meas.Temp:	Instrument in the CAL mode, measuring calibration temperature.
.Req.Buf1:	Instrument in the CAL mode, requesting pH of buffer 1.
.Meas.Buf1:	Instrument in the CAL mode, measures buffer 1.
.Req.Buf2:	Instrument in the CAL mode, requesting pH of buffer 2.
.Meas.Buf2:	Instrument in the CAL mode, measures buffer 2.
etc.	
\$G .Assembly.Bur .Fill:	Buret in filling process
.ModeDis:	Buret in DIS mode
\$G .Prep.X.Active:	Preparing buret X, X=0, 1, 2.
.Empty.X.Active:	Emptying buret X, X=1, 2.

In TIP, its global status as well as the step number (X) is available.

\$G .TIP.X .Inac:	Instrument at the beginning or at the end of a TIP.
.Req .Id1:	Instrument in the TIP mode, requesting Id1 after start.
.Id2:	Instrument in the TIP mode, requesting Id2 after start.
.Id3:	Instrument in the TIP mode, requesting Id3 after start.
.Smp1:	Instrument in the TIP mode, requesting sample size after start.
.Unit:	Instrument in the TIP mode, requesting unit of sample size after start.
.Pause:	Instrument in the TIP mode, in pause.
.Info:	Instrument in the TIP mode, in info.
.Mode...:	Instrument in the TIP mode, working off a submethod. The detailed status messages of the submethod appear, see above.

Status conditions of the global \$H:

The status message of the action which has been held appears.

If the process is held because a monitored limit has been violated, its status message is \$H.Mode.XXX.Titr.

Status conditions of the global \$C:

The status conditions of the global \$C are identical with the ones of the global status \$G. They appear when the process has been restarted actively from the status "Hold" (\$C, key <meas/hold> or automatically after elimination of an error).

Status conditions of the global \$R:

\$R .Mode.XXXX.QuickMeas: Quick manual measurement from the initial status in mode XXXX.

\$R .Mode.DET .Inac: Instrument in the DET mode, inactive.

\$R .Mode.MET .Inac: Instrument in the MET mode, inactive.

\$R .Mode.SET .Inac: Instrument in the SET mode, inactive.

.Cond.Ok: Instrument in the SET mode, conditioning, endpoint reached.

.Cond.Prog: Instrument in the SET mode, conditioning, endpoint not reached.

\$R .Mode.KFT... As SET.

\$R .Mode.STAT .Inac: Instrument in the STAT mode, inactive.

\$R .Mode.DOS .Inac: Instrument in the DOS mode, inactive.

\$R .Mode.DOC .Inac: Instrument in the DOC mode, inactive.

\$R .Mode.MEAS .Inac: Instrument in the MEAS mode, inactive.

\$R .Mode.CAL .Inac: Instrument in the CAL mode, inactive.

\$R .Assembly.Bur.ModeDis: Buret in the DIS mode, inactive.

\$R .TIP.Inac: Instrument in TIP, inactive.

Status conditions of the global \$\$:

\$\$.Mode.XXXX.QuickMeas: Quick manual measurement from the initial status in mode XXXX.

The instrument gives the status from which it has been stopped. The detailed status information is therefore identical to for the global status \$G. Violation of monitored limits with action "end" give the status message \$\$S.Mode.XXX.Inac;EYYY.

3.1.4 Error messages

Error messages are added to the status messages and separated from them by the sign ";".

E8	Card read/write error. Exit: Send new command.
E9	Wrong card, a card has been removed/inserted during the inquiry. Exit: Send new command.
E10	The card has lost data. Exit: Send new command.
E18	Card battery low (it is between 2.37...2.64 V). Exit: Send new command.
E20	Check exchange unit. Exit: Mount Exchange Unit (properly) or &m \$\$.
E21	Check electrode, short circuit. Exit: Rectify fault or &m \$\$.
E22	Check electrode, break. Exit: Rectify fault or &m \$\$.
E23	Division by zero. Exit: The error message disappears on next startup or on recalculation.
E24	Check drive unit. Exit: Connect drive unit (correctly) or &m \$\$.
E26	Manual stop. Exit: The error message disappears on next startup.
E27	Stop V reached in SET, STAT, DOS or DOC. Exit: The error message disappears on next startup.
E28	Wrong object call up Exit: Send correct path for object. Start path at root.
E29	Wrong value or no value allowed. Exit: Send correct value or call up new object.
E30	Wrong trigger, this trigger is not allowed or carrying-out of action not possible. Exit: Send correct trigger (exception: \$D) or call up new object.
E31	Command is not possible in active status. Repeat command in inactive status. Exit: Send new command.
E32	Command is not possible during titration. Repeat command during the conditioning phase or in inactive status. Exit: Send new command.
E33	Value has been corrected automatically. Exit: Send new command.
E34	Instrument at the end of the titration and sample data is edited; the instrument at rest or editing during filling. Exit: &m \$\$.

RS receive errors:

- E36** Parity
Exit: <QUIT> and ensure settings of appropriate parameters at both devices are the same.
- E37** Framing error
Exit: <QUIT> and ensure settings of appropriate parameters at both devices are the same.
- E38** Overrun error. At least 1 character could not be read.
Exit: <QUIT>
- E39** The internal working-off buffer of the Titrino is full (>82 characters).
Exit: <QUIT>

RS send errors:

- E42** CTS=OFF No proper handshake for more than 1 s.
Exit: <QUIT> Is the receiver switched on and ready to receive?
- E43** The transmission of the Titrino has been interrupted with XOFF for at least 6 s.
Exit: Send XON or <QUIT>
- E45** The receive buffer of the Titrino contains an incomplete command (L_f missing). Sending from the Titrino is therefore blocked.
Exit: Send L_f or <QUIT>.
- E120** Overrange of the primary measured value (pH, U, I_{pol} , U_{pol} or T with MEAS T). The secondary measured value (temperature) may be instable as well.
Exit: Correct error or &m \$\$.
- E121** Measuring point list overflow (more than 500 measuring points).
Exit: The error message disappears on next startup.
- E122** EP overflow.
Exit: The error message disappears on next startup or on recalculation.
- E123** Missing EP for calculation.
Exit: The error message disappears on next startup or on recalculation.
- E124** Number of EP does not correspond with the set windows.
Exit: The error message disappears on next startup or on recalculation.
- E125** Missing fix EP for calculation, has not been defined.
Exit: The error message disappears on next startup or on recalculation.
- E126** Fix-EP outside of measuring point list.
Exit: The error message disappears on next startup or on recalculation.
- E128** No new mean.
Exit: The error message disappears on next startup or on recalculation.
- E129** No new common variable, old value remains.
Exit: The error message disappears on next startup or on recalculation.

- E130** Wrong sample. For SET, KFT or DOC with preset titration direction the first measured value lies behind the endpoint.
Exit: The error message disappears on next startup.
- E131** No EP set for SET, STAT.
Exit: The error message disappears on next startup.
- E132** Silo empty and it has been started with open silo or empty silo has been opened.
Exit: Send a silo entry.
- E133** Silo full.
Exit: Send new command.
- E134** No method. A method, which is required from the silo memory or in TIP, does not exist.
Exit: The error message disappears on next startup.
- E135** Check temp.sensor in MEAS T or with activated temperature monitoring.
Exit: Correct error or &m \$\$.
- E136** Same buffer in CAL. Measured value of the second buffer differs less than 6 mV from the measured value of the first buffer.
Exit: Correct error or &m \$\$.
- E137** XXX Bytes are missing so that the method, the silo line could not be stored or not enough RAM for running TIP.
Exit: Send new command.
- E149** Rate missing for calculation.
Exit: The error message disappears on next start or on recalculation.
- E150** Number of rates does not agree with the set windows or the mean rate C80 could not be calculated.
Exit: The error message disappears on next start or on recalculation.
- E151** Fix volume missing for calculation.
Exit: The error message disappears on next startup or on recalculation.
- E155** No new silo result (C24 or C25).
Exit: The error message disappears on next start or on recalculation.
- E157** No sequence defined in TIP.
Exit: The error message disappears on next start.
- E158** A second TIP has been called up in TIP.
Exit: The error message disappears on next start.
- E160** No new temporary variable.
Exit: The error message disappears on next start.
- E161** Measurement range of the secondary measured value (temperature) exceeded. The primary measured value (pH, U, Ipol, Upol) can also be unstable.
Exit: Rectify error or &m \$\$.
- E162** No ramp end set in DOC.
Exit: The error message disappears on next start or &m \$\$.
- E166** Save lines is "OFF" although a submethod of TIP includes an assignment to C24 or C25.
Exit: The error message disappears on next start. Attention: The data of this sample will not be stored.
- E167** Rate too high in DOS. No dispensing possible with the Exchange Unit currently mounted.
Exit: The error message disappears on next start or &m \$\$.
- E170** Fix time missing for calculation.
Exit: The error message disappears on next start or on recalculation.

E171	Rate too low in DOS. No dispensing possible with the Exchange Unit currently mounted. Exit: The error message disappears on next start or &m \$\$.
E172	In TIP or DOS a QuickMeas was started, without defining a measuring quantity. Exit: The error message disappears on next start or &Mode.QuickMeas \$\$.
E173	The warning interval of the internal buret D0 called. Exit: Execute prep &a.b.p \$G or start next method.
E174	The warning interval of the external buret D1 called. Exit: Execute prep &a.b.p \$G or start next method.
E175	The warning interval of the external buret D2 called. Exit: Execute prep &a.b.p \$G or start next method.
E176	The function &Assembly.Buret.Prepare or &Assembly.Buret.Empty was interrupted manually. Exit: The error message disappears on next start.
E177	Accessing to the memory card, the card was not (properly) inserted. Exit: The error message disappears on next start.
E178	The date of changing the battery of the card is expired. Exit: The error message disappears on next start.
E180	Memory card write-protected. Exit: Send new command.
E181	Memory card not formatted. Exit: Send new command.
E182	Memory card not accessible. Exit: Send new command.
E183	A directory with the same name exists already on the memory card. Exit: Send new command.
E184	Measured value below lower limit. Exit: The error message disappears when again within the limits or &m \$\$.
E185	Measured value above upper limit. Exit: The error message disappears when again within the limits or &m \$\$.
E186	Temperature below lower limit. Exit: The error message disappears when again within the limits or &m \$\$.
E187	Temperature above upper limit. Exit: The error message disappears when again within the limits or &m \$\$.
E188	Rate below lower limit. Exit: The error message disappears when again within the limits or &m \$\$.
E189	Rate above upper limit. Exit: The error message disappears when again within the limits or &m \$\$.
E196	Result is out of limits. Exit: The error message disappears on next start or on recalculation.
E197	Sample size is out of limits. Exit: The error message disappears on next start or on introduction of new sample size.
E198	Validation interval is expired. Exit: The error message disappears on next start or clear counter with &Config.Monitoring.Validation.ClearCount \$G.
E199	Service date is reached. Exit: The error message disappears on next start or change date in &Config.Monitoring.Service.Date.

- E203** No Oven parameters: Oven not (correctly) connected.
Exit: The error message disappears on next start. If you don't wish oven parameters in your report, select `&Mode.Parameter.Presel.Oven "no"` in your method(s).
- E205** Calibration interval is expired.
Exit: The error message disappears on next calibration or if you delete the calibration.
- E212** Transmission error from Remote Box. Unknown characters.
Exit: Rectify error and switch Titrino off and on again.
- E213** Time-out error from PC keyboard (Remote Box)
Exit: Rectify error and switch Titrino off and on again.
- E214** Check Remote Box. Remote Box not (properly) connected but activated in `&Config.Periph.RemoteBox`.
Exit: Rectify error and switch Titrino off and on again.
- E270** Overload in dosing element: piston is blocked.
Exit: `<STOP>`, `<QUIT>`. Rectify fault at the wet part. Mount it again onto the dos driver and press `<CLEAR>`. The dosing element will be initialized.
- E282** Overload in dosing element: cock is blocked.
Exit: `<STOP>`, `<QUIT>`. Rectify fault at the wet part. Mount it again onto the dos driver and press `<CLEAR>`. The dosing element will be initialized.

3.2 Remote control commands

3.2.1 Overview

The internal object tree can be divided into the following branches:

&	Root
Mode	Method parameters
UserMeth	Administration of the internal user-memory for methods
MemoryCard	Administration of the memory card
Config	Instrument configuration
SmplData	Sample specific data
Hotkey	Keys with direct access
Info	Current Data
Assembly	Component data
Setup	Setting the operating mode
Diagnosis	Diagnostics program

&Mode

Object	Description	Input range	Reference
& Root			
├─ Mode	Mode	\$G, \$\$, \$H, \$C	3.2.2.1.
│ ├─ .QuickMeas	Rapid meas. in basic mode	\$G, \$\$	3.2.2.2.
│ ├─ .Select	Mode selection	DET , MET, SET, KFT, STAT, DOS, DOC, MEAS, CAL, TIP	3.2.2.3.
│ ├─ .DETQuantity	Measured quantity for DET	pH , U, lpol, Upol	ditto
│ ├─ .METQuantity	Measured quantity for MET	pH , U, lpol, Upol	ditto
│ ├─ .SETQuantity	Measured quantity for SET	pH , U, lpol, Upol	ditto
│ ├─ .MEASQuantity	Measured quantity for MEAS	pH , U, lpol, Upol, T	ditto
│ ├─ .KFTQuantity	Measured quantity for KFT	lpol , Upol	ditto
│ ├─ .STATQuantity	Measured quantity for STAT	pH , U, lpol, Upol	ditto
│ ├─ .DOCQuantity	Measured quantity for DOC	pH , U, lpol, Upol	ditto
│ ├─ .Name	Name of current method	read only/read + write	3.2.2.4.
│ └─ .Parameter*	Parameter of current mode, see below		
├─ Def	Definitions for data output		
│ ├─ .Formulas	Calculation formulas		
│ │ ├─ .1	for result 1		
│ │ │ ├─ .Formula	Calculation formula	special	3.2.2.5.
│ │ │ ├─ .TextRS	Text for result output	up to 8 ASCII char	ditto
│ │ │ ├─ .Decimal	Number of decimal places	0...2...5	ditto
│ │ │ ├─ .Unit	Unit for result output	up to 6 ASCII char	ditto
│ │ │ ├─ .Limits	Limits for result	ON, OFF	ditto
│ │ │ ├─ .LoLim	Lower limit	0...±999 999	ditto
│ │ │ ├─ .UpLim	Upper limit	0...±999 999	ditto
│ │ │ └─ .Output	Output on L13	active, pulse, OFF	ditto
│ │ └─ :	up to 9 results		
│ ├─ .SiloCalc	Silo calculations		
│ │ ├─ .Assign	Assignment		
│ │ │ ├─ .C24	Store as variable C24	RSX, EPX, CXX	3.2.2.6.
│ │ │ └─ .C25	Store as variable C25	RSX, EPX, CXX	
│ │ └─ .MatchId	Matching of Id's	id1, id1&2, all, OFF	
│ ├─ .ComVar	Assignment of common variables		
│ │ ├─ .C30	for C30	RSX, EPX, CXX, MNX	3.2.2.7.
│ │ └─ up to C39			
│ ├─ .Report	Reports at the end of determination		
│ │ ├─ .Assign1	Output to COM1	special	3.2.2.8.
│ │ └─ .Assign2	Output to COM 2	as COM1	
│ ├─ .Mean	Assignment for mean calculation		
│ │ ├─ .1	MN1		
│ │ │ └─ .Assign	Input of variable	RSX, EPX, CXX	3.2.2.9.
│ │ └─ :			
│ ├─ .TempVar	Assignment of temporary variables		
│ │ ├─ .C70	for C70	RSX, EPX, CXX	3.2.2.10.
│ │ └─ up to C79			
└─ .CFmla	Calculation constants		
│ ├─ .1	Calculation constant C01		
│ │ └─ .Value	Input of value	0...±999 999	3.2.2.11.
│ └─ up to C19			

*Parameter	Tree part "Parameters for DET"		
.TitrPara			
Titration parameters			
.MptDensity	Measuring point density	0... 4 ...9	3.2.2.12.
.MinIncr	Minimum increment	0... 10.0 ...999.9	ditto
.DosRate	Dispensing rate for increments	0.01...150.0, max.	3.2.2.13.
.SignalDrift	Drift for meas. value acquisition	depends on meas.quant.	3.2.2.14.
.UnitSigDrift	Unit of measured value drift	read only	ditto
.EquTime	Equilibrium time	0... 26 ...9999, OFF	ditto
.StartV	Start volume		
.Type	Type of start volume	abs., rel., OFF	3.2.2.15.
.V	Volume for absolute start volume	0 ...999.99	ditto
.Factor	Factor for relative start volume	0 ...±999 999	ditto
.Rate	Dispensing rate for start volume	0.01...150.0, max.	ditto
.Pause	Waiting time	0 ...999 999	3.2.2.16.
.DosUnit	Selection of the dosing unit	internal D0 , external D1, external D2	3.2.2.17.
.MeasInput	Measuring input	1 , 2, diff.	3.2.2.18.
.Ipol	Polarization current	0... 1 ...±127	ditto
.Upol	Polarization voltage	0... 400 ...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON, OFF	ditto
.Temp	Titration temperature	-170.0... 25.0 ...500.0	3.2.2.19.
.StopCond			
Stop conditions			
.VStop	Stop volume		
.Type	Type of stop volume	abs. , rel., OFF	3.2.2.20.
.V	Volume for absolute stop volume	0... 99.99 ...9999.99	ditto
.Factor	Factor for relative stop volume	0...± 999 999	ditto
.MeasStop	Stop measured value pH, U, I	depends on meas.quant.	3.2.2.21.
.UnitMStop	Unit of stop measured value	read only	ditto
.EPStop	Stop after a number of EP's	1... 9 , OFF	3.2.2.22.
.FillRate	Filling rate	0.01...150.0, max.	3.2.2.23.
.Statistics			
Statistics			
.Status	Status of statistics calculation	ON, OFF	3.2.2.24.
.MeanN	No. of individual determinations	2 ...20	ditto
.ResTab	Result table		
.Select		original , delete n, delete all	ditto
.DeIN	Deletion of individual results	1 ...20	ditto
.Evaluation			
Evaluation			
.EPC	EP criterion	0... 5 ...200	3.2.2.25.
.Recognition	EP recognition		
.Select	Type of EP recognition	all , greatest, last, window, OFF	ditto
.Window	Window		
.1	up to 9 windows		
.LowLim	Lower limit window 1	depends on meas.quant.	ditto
.UpLim	Upper limit window 1	depends on meas.quant.	ditto
.FixEP	Fix endpoints		
.1	up to 9 fix EP's		
.Value	Measured value for fix EP1	depends on meas.quant.	ditto
.pK	pK or HNP evaluation	ON, OFF	ditto
.Presel			
Preselections			
.IReq	Request of Id's after start	id1, id1&2, all, OFF	3.2.2.26.
.Sreq	Request of smpl size after start	value, unit, all, OFF	ditto
.LimSmplSize	Limits for sample size		3.2.2.27.
.Status	Status of limit control	ON, OFF	ditto
.LoLim	Lower limit	0.0 ...999 999	ditto
.UpLim	Upper limit	0.0... 999 999	ditto
.ActPulse	Output of a pulse	ON, OFF	3.2.2.28.

*Parameter	Tree part "Parameters for MET"		
.TitrPara			
Titration parameters			
.VStep	Volume increment	0... 0.10 ...999.9	3.2.2.12.
.DosRate	Dispensing rate for increments	0.01...150.0, max.	3.2.2.13.
.SignalDrift	Drift for meas. value acquisition	depends on meas.quant.	3.2.2.14.
.UnitSigDrift	Unit of measured value drift	read only	ditto
.EquTime	Equilibrium time	0... 26 ...9999, OFF	ditto
.StartV			
Start volume			
.Type	Type of start volume	abs., rel., OFF	3.2.2.15.
.V	Volume for absolute start volume	0...999.99	ditto
.Factor	Factor for relative start volume	0...±999 999	ditto
.Rate	Dispensing rate for start volume	0.01...150.0, max.	ditto
.Pause	Waiting time	0...999 999	3.2.2.1.
.DosUnit	Selection of the dosing unit	internal D0 , external D1, external D2	3.2.2.17.
.MeasInput	Measuring input	1, 2, diff.	3.2.2.18.
.Ipol	Polarization current	0...1...±127	ditto
.Upol	Polarization voltage	0... 400 ...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON, OFF	ditto
.Temp	Titration temperature	-170.0... 25.0 ...500.0	3.2.2.19.
.StopCond			
Stop conditions			
.VStop			
Stop volume			
.Type	Type of stop volume	abs. , rel., OFF	3.2.2.20.
.V	Volume for absolute stop volume	0... 99.99 ...9999.99	ditto
.Factor	Factor for relative stop volume	0...± 999 999	ditto
.MeasStop	Stop measured value pH, U, I	depends on meas.quant.	3.2.2.21.
.UnitMStop	Unit of stop measured value	read only	ditto
.EPStop	Stop after a number of EP's	1... 9 , OFF	3.2.2.22.
.FillRate	Filling rate	0.01...150.0, max.	3.2.2.23.
.Statistics			
Statistics			
.Status	Status of statistics calculation	ON, OFF	3.2.2.24.
.MeanN	No. of individual determinations	2 ...20	ditto
.ResTab			
Result table			
.Select		original , delete n, delete all	ditto
.DelN	Deletion of individual results	1 ...20	ditto
.Evaluation			
Evaluation			
.EPC	EP criterion	depends on meas.quant.	3.2.2.25.
.Recognition			
EP recognition			
.Select	Type of EP recognition	all , greatest, last, window, OFF	ditto
.Window			
Window			
.1	up to 9 windows		
.LowLim	Lower limit window 1	depends on meas.quant.	ditto
.UpLim	Upper limit window 1	depends on meas.quant.	ditto
.FixEP			
Fix endpoints			
.1	up to 9 fix EP's		
.Value	Measured value for fix EP1	depends on meas.quant.	ditto
.pK	pK or HNP evaluation	ON, OFF	ditto
.Presel			
Preselections			
.IReq	Request of Id's after start	id1, id1&2, all, OFF	3.2.2.26.
.SReq	Request of sample size after start	value, unit, all, OFF	ditto
.LimSmplSize			
Limits for sample size			
.Status	Status of limit control	ON, OFF	ditto
.LoLim	Lower limit	0.0 ...999 999	ditto
.UpLim	Upper limit	0.0... 999 999	ditto
.ActPulse	Output of a pulse	ON, OFF	3.2.2.28.

*Parameter	Tree part "Parameters for SET"		
.SET1	Control parameters for EP1		
.EP Endpoint 1	depends on meas.quant.	3.2.2.29.	
.UnitEp	Unit of endpoint	read only	ditto
.Dyn	Dynamics	depends on meas.quant.	3.2.2.30.
.UnitDyn	Unit of dynamics	read only	ditto
.MaxRate	Maximum dosing rate	0.01... 10 ...150, max.	ditto
.MinRate	Minimum dosing rate	0.01... 25.0 ...9999	ditto
.Stop	Titration stop		
.Type	Type of stop criterion	drift , time	3.2.2.31.
.Drift	Stop drift	1... 20 ...999	ditto
.Time	Switch-off delay time	0... 10 ...999, inf	ditto
.StopT	Stop time	0...999 999, OFF	ditto
.SET2	Control parameters for EP2, as for EP1		
.TitrPara	Titration parameters		
.Direction	Titration direction	+, -, auto	3.2.2.32.
.XPause	Waiting time before start volume	0 ...999 999	3.2.2.33.
.Start V	Start volume		
.Type	Type of start volume	abs., rel., OFF	3.2.2.15.
.V	Volume for absolute start volume	0 ...999.99	ditto
.Factor	Factor for relative start volume	0 ...±999 999	ditto
.Rate	Dispensing rate for start volume	0.01...150.0, max.	ditto
.Pause	Waiting time after start volume	0 ...999 999	3.2.2.16.
.ExtrT	Extraction time	0 ...999 999	3.2.2.34.
.DosUnit	Selection of the dosing unit	internal D0 , external D1, external D2	3.2.2.17. 3.2.2.18.
.MeasInput	Measuring input	1 , 2 , diff.	3.2.2.18.
.Ipol	Polarization current	0... 1 ...±127	ditto
.Upol	Polarization voltage	0... 400 ...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON , OFF	ditto
.Temp	Titration temperature	-170.0... 25.0 ...500.0	3.2.2.19.
.TDelta	Time interv. for meas.acquisition	1... 2 ...999 999	3.2.2.35.
.StopCond	Stop conditions		
.VStop	Stop volume		
.Type	Type of stop volume	abs. , rel., OFF	3.2.2.20.
.V	Volume for absolute stop volume	0... 99.99 ...9999.99	ditto
.Factor	Factor for relative stop volume	0...± 999 999	ditto
.FillRate	Filling rate	0.01...150.0, max.	3.2.2.23.
.Statistics	Statistics		
.Status	Status of statistics calculation	ON , OFF	3.2.2.24.
.MeanN	No. of individual determinations	2 ...20	ditto
.ResTab	Result table		
.Select		original , delete n, delete all	ditto
.DelN	Deletion of individual results	1 ...20	ditto
.Presel	Preselections		
.Cond	Conditioning	ON , OFF	3.2.2.36.
.DriftDisp	Display of drift during cond.	ON , OFF	ditto
.DCor	Drift correction		
.Type	Type of drift acquisition	auto, man., OFF	ditto
.Value	Drift value for manual drift corr.	0.0 ...99.9	ditto
.IReq	Request of Id's after start	id1, id1&2, all, OFF	3.2.2.26.
.SReq	Request of smpl size after start	value, unit, all, OFF	ditto
.LimSmplSize	Limits for sample size		3.2.2.27
.Status	Status of limit control	ON , OFF	ditto
.LoLim	Lower limit	0.0 ...999 999	ditto
.UpLim	Upper limit	0.0... 999 999	ditto
.ActPulse	Output of a pulse	first, all, cond., OFF	3.2.2.28.

*Parameter	Tree part "Parameters for KFT"		
CtrlPara			
.EP	Endpoint	depends on meas.quant.	3.2.2.37.
.UnitEp	Unit of endpoint		read only ditto
.Dyn	Dynamics	depends on meas.quant.	3.2.2.38.
.UnitDyn	Unit of dynamics		read only ditto
.MaxRate	Maximum dosing rate		0.01...150, max. ditto
.MinIncr	Minimum increment		0.1...9.9, min. ditto
Stop			
.Type	Type of stop criterion		drift , time 3.2.2.30.
.Drift	Stop drift		1...20...999 ditto
.Time	Switch-off delay time		0...10...999, inf ditto
.StopT	Stop time		0...999 999, OFF ditto
Titration parameters			
Titration parameters			
.Direction	Titration direction		+, -, auto 3.2.2.31.
.XPause	Waiting time before start volume		0...999 999 3.2.2.32.
StartV			
.Type	Type of start volume		abs., rel., OFF 3.2.2.15.
.V	Volume for absolute start volume		0...999.99 ditto
.Factor	Factor for relative start volume		0...±999 999 ditto
.Rate	Dispensing rate for start volume		0.01...150.0, max. ditto
.Pause	Waiting time after start volume		0...999 999 3.2.2.16.
.ExtrT	Extraction time		0...999 999 3.2.2.33.
.DosUnit	Selection of the dosing unit		internal D0 , external D1, external D2 3.2.2.17.
.MeasInput		without meaning	
.Ipol	Polarization current		0...50...±127 3.2.2.18.
.Upol	Polarization voltage		0...400...±1270 ditto
.PolElectrTest	Test for polarized electrodes		ON, OFF ditto
.Temp	Titration temperature		-170.0...25.0...500.0 3.2.2.19.
.TDelta	Time interv. for meas.acquisition		1...2...999 999 3.2.2.34.
Stop conditions			
Stop conditions			
VStop			
.Type	Type of stop volume		abs. , rel., OFF 3.2.2.20.
.V	Volume for absolute stop volume		0...99.99...9999.99 ditto
.Factor	Factor for relative stop volume		0...±999 999 ditto
.FillRate	Filling rate		0.01...150.0, max. 3.2.2.23.
Statistics			
Statistics			
.Status	Status of statistics calculation		ON, OFF 3.2.2.24.
.MeanN	No. of individual determinations		2...20 ditto
ResTab			
.Select			original , delete n, delete all ditto
.DelN	Deletion of individual results		1...20 ditto
Preselections			
Preselections			
Cond			
.Cond	Conditioning		ON , OFF 3.2.2.35.
.DriftDisp	Display of drift during cond.		ON , OFF ditto
DCor			
.Type	Type of drift acquisition		auto, man., OFF ditto
.Value	Drift value for manual drift corr.		0.0...99.9 ditto
.IReq	Request of Id's after start		id1, id1&2, all, OFF 3.2.2.26.
.SReq	Request of smpl size after start		value, unit, all, OFF ditto
LimSmplSize			
.Status	Status of limit control		ON, OFF 3.2.2.27
.LoLim	Lower limit		0.0...999 999 ditto
.UpLim	Upper limit		0.0...999 999 ditto
.Oven	KF Oven connected		COM1, COM2, no 3.2.2.39.
.ActPulse	Output of a pulse		first, all, cond., OFF 3.2.2.28.

*Parameter	Tree part "Parameters for STAT"		
.CtrlPara Control parameters			
.EP	Measuring point	depends on meas.quant.	3.2.2.36.
.UnitEp	Unit of endpoint	read only	ditto
.Dyn	Dynamics	depends on meas.quant.	3.2.2.37.
.UnitDyn	Unit of dynamics	read only	ditto
.MaxRate	Maximum dosing rate	0.01... 10 ...150, max.	ditto
.MinRate	Minimum dosing rate	0.01... 25.0 ...9999	3.2.2.40.
.TitrPara Titration parameters			
.StartV	Start volume		
.Type	Type of start volume	abs., rel., OFF	3.2.2.15.
.V	Volume for absolute start volume	0 ...999.99	ditto
.Factor	Factor for relative start volume	0 ...±999 999	ditto
.Rate	Dispensing rate for start volume	0.01...150.0, max.	ditto
.Pause	Waiting time after start volume	0 ...999 999	3.2.2.16.
.Tstart	Start time	0 ...999 999	3.2.2.41.
.StartMeas	Start measuring value	depends on meas.quant.	ditto
.RateStart	Start rate	0.01...150, OFF	ditto
.TDelta	Time interv.for meas.acquisition	1... 2 ...999 999	3.2.2.34.
.Direction	Titration direction	+, -, auto	3.2.2.31.
.DosUnit	Selection of the dosing unit	internal D0 , external D1, external D2	3.2.2.17.
.MeasInput	Measuring input	1 , 2, diff.	3.2.2.18.
.Ipol	Polarization current	0... 1 ...±127	ditto
.Upol	Polarization voltage	0... 400 ...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON, OFF	ditto
.Temp	Titration temperature	-170.0... 25.0 ...500.0	3.2.2.19.
.StopCond Stop conditions			
.TStop	Stop time		
.Type	Type of stop time	abs.,rel.,delta,delay, OFF	3.2.2.42.
.Time	Time for absolute stop time	0... 999 999	ditto
.Factor	Factor for relative stop time	0...± 999 999	ditto
.Delta	Time after EP is once reached	0... 999 999	ditto
.Delay	Time after last dosing step	0... 999 999	ditto
.VStop	Stop volume		
.Type	Type of stop volume	abs. , rel., OFF	3.2.2.20.
.V	Volume for absolute stop volume	0... 99.99 ...9999.99	ditto
.Factor	Factor for relative stop volume	0...± 999 999	ditto
.RateStop	Stop rate	0.01...150.0, OFF	3.2.2.43.
.FillRate	Filling rate	0.01...150.0, max.	3.2.2.23.
.Statistics Statistics			
.Status	Status of statistics calculation	ON, OFF	3.2.2.24.
.MeanN	No. of individual determinations	2 ...20	ditto
.ResTab	Result table		
.Select		original ,delete n,delete all	ditto
.DelN	Deletion of individual results	1 ...20	ditto

*Parameter	Tree part "Parameters for STAT", continuation		
Evaluation			
.TimeWin	Time windows for rate evaluation		
.1	up to 9 windows		
.LowLim	Lower limit window 1	0...999 999, OFF	3.2.2.44.
.UpLim	Upper limit window 1	0...999 999, OFF	ditto
.FixVol	Fix volumes		
.1	up to 9 fix volumes		
.Value	Value for fix volume 1	0...999 999, OFF	3.2.2.45.
.FixTime	Fix times		
.1	up to 9 fix times		
.Value	Value for fix time 1	0.01...1, OFF	3.2.2.46.
Monitoring			
.MeasVal	Monitoring of measured values		
.Status	Status	ON, OFF	3.2.2.47.
.LowLim	Lower limit	depends on meas.quant.	ditto
.UpLim	Upper limit	depends on meas.quant.	ditto
.Action	Action if out of limits	end, hold, wait, none	ditto
.Rate	Monitoring of rates		
.Status	Status	ON, OFF	3.2.2.48.
.LowLim	Lower limit	0.000...150	ditto
.UpLim	Upper limit	0.000...150	ditto
.Action	Action if out of limits	end, hold, wait, none	ditto
.Temp	Monitoring of temperatures		
.Status	Status	ON, OFF	3.2.2.49.
.LowLim	Lower limit	-170.0...500.0	ditto
.UpLim	Upper limit	-170.0...500.0	ditto
.Action	Action if out of limits	end, hold, wait, none	ditto
.L10Output	Assignment of output L10		
.AssignOutput	To type of limit violation	meas,temp,rate,all, none	3.2.2.50.
.Limit	Assign output to violated limit	upper, lower, any	ditto
.Output	Signal to be set on output L10	active, pulse	ditto
.L11Output	as for output L10		
.L12Output	as for output L10		
.L13Output	as for output L10		
Preselections			
.IReq	Request of Id's after start	id1, id1&2, all, OFF	3.2.2.26.
.SReq	Request of sample size after start	value, unit, all, OFF	ditto
.LimSmplSize	Limits for sample size		3.2.2.27
.Status	Status of limit control	ON, OFF	ditto
.LoLim	Lower limit	0.0...999 999	ditto
.UpLim	Upper limit	0.0...999 999	ditto
.RateDisp	Display of drift during cond.	ON, OFF	3.2.2.51.
.ActPulse	Output of a pulse	ON, OFF	3.2.2.28.

*Parameter	Tree part "Parameters for DOS"		
.DosPara	Dosing parameters		
.Type	Type of dosing	volume , time, rate	3.2.2.52.
.Volume	Volume dosing		
.Volume	Volume	0... 10 ...99 999.99	ditto
.DisType	Second dosing criterion	time, rate	ditto
.Rate	Rate	0.001...150, max.	ditto
.Time	Dosing time	1... 100 ...999 999	ditto
.Time	Time dosing		
.Time	Dosing time	1... 100 ...999 999	ditto
.DisType	Second dosing criterion	volume , rate	ditto
.Rate	Rate	0.001...150, max.	ditto
.Volume	Volume	0... 10 ...99 999.99	ditto
.Rate	Rate dosing		
.Rate	Rate	0.001...150, max.	ditto
.StopType	Second dosing criterion	volume , time	ditto
.Time	Dosing time	1... 100 ...999 999	ditto
.Volume	Volume	0... 10 ...99 999.99	ditto
.Pause	Waiting time after start volume	0 ...999 999	3.2.2.53.
.TDelta	Time interv.for meas.acquisition	1... 10 ...999 999	3.2.2.54.
.DosUnit	Selection of the dosing unit	internal D0 , external D1, external D2	3.2.2.17.
.Temp	Titration temperature	-170.0... 25.0 ...500.0	3.2.2.55.
.StopCond	Stop conditions		
.VStop	Stop volume		
.Type	Type of stop volume	abs., rel., OFF	3.2.2.20.
.V	Volume for absolute stop volume	0... 99.99 ...9999.99	ditto
.Factor	Factor for relative stop volume	0...± 999 999	ditto
.FillRate	Filling rate	0.01...150.0, max.	3.2.2.23.
.Statistics	Statistics		
.Status	Status of statistics calculation	ON, OFF	3.2.2.24.
.MeanN	No.of individual determinations	2 ...20	ditto
.ResTab	Result table		
.Select		original , delete n, delete all	ditto
.DeIN	Deletion of individual results	1 ...20	ditto
.Monitoring	Monitoring		
.MeasVal	Monitoring of measured values		
.MeasMode	Measured value	pH, U, Ipol, Upol, OFF	3.2.2.56.
.MeasInput	Measuring input	1 , 2, diff.	ditto
.PCurrent	Polarization current	0... 1 ...±127	ditto
.PVoltage	Polarization voltage	0... 400 ...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON, OFF	ditto
.pH	Limits for monitoring of pH values		
.LowLim	Lower limit	-20.00 ...20.00	ditto
.UpLim	Upper limit	-20.00 ...20.00	ditto
.U	Limits for monitoring of U values		
.LowLim	Lower limit	-2000 ...2000	ditto
.UpLim	Upper limit	-2000 ...2000	ditto
.Ipol	Limits for monitoring of Ipol values		
:	as for U		
.U	Limits for monitoring of Upol values		
:	as above		
.Action	Action if out of limits	end, hold, wait, none	ditto

*Parameter	Tree part "Parameters for DOS, monitoring", continuation		
<ul style="list-style-type: none"> - .Monitoring - .Temp - .Status - .LowLim - .UpLim - .Action - .L10Output - .AssignOutput - .Limit - .Output - .L11Output - .L12Output - .L13Output - .Presel - .IReq - .SReq - .LimSmplSize - .Status - .LoLim - .UpLim - .ActPulse 	<p>Monitoring</p> <p>Monitoring of temperatures</p> <p>Status</p> <p>Lower limit</p> <p>Upper limit</p> <p>Action if out of limits</p> <p>Assignment of output L10</p> <p>To type of limit violation</p> <p>Assign output to violated limit</p> <p>Signal to be set on output L10</p> <p>as for output L10</p> <p>as for output L10</p> <p>as for output L10</p> <p>Preselections</p> <p>Request of Id's after start</p> <p>Request of smpl size after start</p> <p>Limits for sample size</p> <p>Status of limit control</p> <p>Lower limit</p> <p>Upper limit</p> <p>Output of a pulse</p>	<p>ON, OFF</p> <p>-170.0...500.0</p> <p>-170.0...500.0</p> <p>end, hold, wait, none</p> <p>meas,temp,all,none</p> <p>upper, lower, any</p> <p>active, pulse</p> <p>id1, id1&2, all, OFF</p> <p>value, unit, all, OFF</p> <p>ON, OFF</p> <p>0.0...999 999</p> <p>0.0...999 999</p> <p>ON, OFF</p>	<p>3.2.2.47.</p> <p>ditto</p> <p>ditto</p> <p>ditto</p> <p>3.2.2.48.</p> <p>ditto</p> <p>ditto</p> <p>3.2.2.26.</p> <p>ditto</p> <p>3.2.2.27</p> <p>ditto</p> <p>ditto</p> <p>3.2.2.28.</p>

*Parameter	Tree part "Parameters for DOC"		
.DosPara			
Dosing parameters			
.BeginMeas	Begin of measuring	depends on meas.quant.	3.2.2.57.
.EndMeas	End of measuring	depends on meas.quant.	ditto
.UnitMeas	Unit of measured quantity	read only	ditto
.SweepTime	Sweep time	0... 300 ...999 999	3.2.2.58.
.Dyn	Dynamics	depends on meas.quant.	3.2.2.59.
.UnitDyn	Unit of dynamics	read only	ditto
.MaxRate	Maximum dosing rate	0.01... 150, max.	ditto
.MinRate	Minimum dosing rate	0.01... 5.0 ...9999	ditto
.Direction	Direction	+, -, auto	3.2.2.60.
.StartV	Start volume		
.Type	Type of start volume	abs., rel., OFF	3.2.2.61.
.V	Volume for absolute start volume	0...999.99	ditto
.Factor	Factor for relative start volume	0...±999 999	ditto
.Rate	Dispensing rate for start volume	0.01...150.0, max.	ditto
.Pause	Waiting time after start volume	0...999 999	ditto
.TDelta	Time interv. for meas.acquisition	1... 10 ...999 999	3.2.2.52.
.DosUnit	Selection of the dosing unit	internal D0 , external D1, external D2	3.2.2.17.
.MeasInput	Measuring input	1, 2, diff.	3.2.2.62.
.Ipol	Polarization current	0... 1 ...±127	ditto
.Upol	Polarization voltage	0... 400 ...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON, OFF	ditto
.Temp	Titration temperature	-170.0... 25.0 ...500.0	3.2.2.53.
.StopCond			
Stop conditions			
.VStop	Stop volume		
.Type	Type of stop volume	abs. , rel., OFF	3.2.2.20.
.V	Volume for absolute stop volume	0... 99.99 ...9999.99	ditto
.Factor	Factor for relative stop volume	0...± 999 999	ditto
.FillRate	Filling rate	0.01...150.0, max.	3.2.2.23.
.Statistics			
Statistics			
.Status	Status of statistics calculation	ON, OFF	3.2.2.24.
.MeanN	No. of individual determinations	2 ...20	ditto
.ResTab	Result table		
.Select		original , delete n, delete all	ditto
.DeIN	Deletion of individual results	1 ...20	ditto
.Monitoring			
Monitoring			
.MeasVal	Monitoring of measured values		
.Status	Status	ON, OFF	3.2.2.45.
.LowLim	Lower limit	depends on meas.quant.	ditto
.UpLim	Upper limit	depends on meas.quant.	ditto
.Action	Action if out of limits	end, hold, wait, none	ditto
.Temp	Monitoring of temperatures		
.Status	Status	ON, OFF	3.2.2.47.
.LowLim	Lower limit	-170.0 ...500.0	ditto
.UpLim	Upper limit	-170.0... 500.0	ditto
.Action	Action if out of limits	end, hold, wait, none	ditto

*Parameter	Tree part "Parameters for DOC, monitoring", continuation		
<ul style="list-style-type: none"> - .Monitoring <ul style="list-style-type: none"> - .L10Output - .AssignOutput - .Limit - .Output - .L11Output - .L12Output - .L13Output - .Presel <ul style="list-style-type: none"> - .IReq - .SReq - .LimSmplSize <ul style="list-style-type: none"> - .Status - .LoLim - .UpLim - .RateDisp - .ActPulse 	<p>Monitoring</p> <p>Assignment of output L10</p> <p>To type of limit violation</p> <p>Assign output to violated limit</p> <p>Signal to be set on output L10</p> <p>as for output L10</p> <p>as for output L10</p> <p>as for output L10</p> <p>Preselections</p> <p>Request of Id's after start</p> <p>Request of sample size after start</p> <p>Limits for sample size</p> <p>Status of limit control</p> <p>Lower limit</p> <p>Upper limit</p> <p>Display of drift during cond.</p> <p>Output of a pulse</p>	<p>meas,temp,all,none</p> <p>upper, lower, any</p> <p>active, pulse</p> <p>id1, id1&2, all, OFF</p> <p>value, unit, all, OFF</p> <p>ON, OFF</p> <p>0.0...999 999</p> <p>0.0...999 999</p> <p>ON, OFF</p> <p>ON, OFF</p>	<p>3.2.2.48.</p> <p>ditto</p> <p>ditto</p> <p>3.2.2.26.</p> <p>ditto</p> <p>3.2.2.27</p> <p>ditto</p> <p>ditto</p> <p>ditto</p> <p>3.2.2.49.</p> <p>3.2.2.28.</p>

*Parameter	Tree part "Parameters for MEAS"		
.Measuring	Measuring parameters		
.SignalDrift	Drift for meas.value acquisition	depends on meas.quant.	3.2.2.63.
.UnitSigDrift	Unit of measured value drift	read only	ditto
.EquTime	Equilibrium time	0...9999, OFF	ditto
.MeasInput	Measuring input	1, 2, diff.	3.2.2.64.
.Ipol	Polarization current	0...1...±127	ditto
.Upol	Polarization voltage	0...400...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON, OFF	ditto
.Temp	Titration temperature	-170.0...25.0...500.0	3.2.2.65.
.TDelta	Time interv.for meas.acquisition	1...2...999 999	3.2.2.66.
.Statistics	Statistics		
.Status	Status of statistics calculation	ON, OFF	3.2.2.24.
.MeanN	No. of individual determinations	2...20	ditto
.ResTab	Result table		
.Select		original , delete n, delete all	ditto
.DelN	Deletion of individual results	1...20	ditto
.Presel	Preselections		
.IReq	Request of Id's after start	id1, id1&2, all, OFF	3.2.2.26.
.SReq	Request of sample size after start	value, unit, all, OFF	ditto
.LimSmplSize	Limits for sample size		3.2.2.27
.Status	Status of limit control	ON, OFF	ditto
.LoLim	Lower limit	0.0...999 999	ditto
.UpLim	Upper limit	0.0...999 999	ditto
.ActPulse	Output of a pulse	ON, OFF	3.2.2.28.

*Parameter	Tree part "Parameters for CAL"		
.Calibration	Calibration parameters		
.MeasInput	Measuring input	1, 2, diff.	3.2.2.67.
.CalTemp	Calibration temperature	-20.0...25.0...120.0	3.2.2.68.
.Buffer			
.1			
.Value	pH value of buffer 1	0...7.00...±20.00	3.2.2.69.
.2			
.Value	pH value of buffer 2	0...4.00...±20.00, OFF	ditto
.	up to 9 buffers		
.SignalDrift	Drift for meas.value acquisition	depends on meas.quant.	3.2.2.70.
.EquTime	Equilibrium time	0...26...9999, OFF	ditto
.Electrodel	Electrode identification	8 ASCII char.	3.2.2.71.
.SmplChanger	Calibration on a Titrino	ON, OFF	3.2.2.72.
.ActPulse	Output of a pulse	first, all, OFF	3.2.2.73.
.Statistics	Statistics		
.Status	Status of statistics calculation	ON, OFF	3.2.2.24.
.MeanN	No. of individual determinations	2...20	ditto
.ResTab	Result table		
.Select		original , delete n, delete all	ditto
.DelN	Deletion of individual results	1...20	ditto

*Parameter	Tree part "Parameters for TIP"		
.Sequence	Sequence		
.1	Step 1		
.Select	Step selection	method,pause,L4 output,L6 output, info,prep,stirrer, OFF	3.2.2.74.
.Method	Method from mem.or card	special	3.2.2.75.
.Pause	Waiting time	0 ...999 999, INF	ditto
.L4Output	Line L4	active,inactive,pulse, OFF	ditto
.L6Output	Line L6	active,inactive,pulse, OFF	ditto
.Info	Display information	up to 16 ASCII char.	ditto
.Prep	Preparation of titrating buret	internal D0 , external D1, external D2	ditto
.Stirrer	Stirrer	ON , OFF	ditto
:	up to 30 steps		
.Statistics	Statistics		
.Status	Status of statistics calculation	ON , OFF	3.2.2.24.
.MeanN	No. of individual determinations	2 ...20	ditto
.ResTab	Result table		
.Select		original ,delete n,delete all	ditto
.DelN	Deletion of individual results	1 ...20	ditto
.Presel	Preselections		
.IReq	Request of Id's after start	id1, id1&2, all, OFF	3.2.2.26.
.SReq	Request of sample size after start	value, unit, all, OFF	ditto
.LimSmplSize	Limits for sample size		3.2.2.27
.Status	Status of limit control	ON , OFF	ditto
.LoLim	Lower limit	0.0 ...999 999	ditto
.UpLim	Upper limit	0.0 ... 999 999	ditto
.MeasMode	Measuring mode for man.meas.	pH,U,lpol,Upol,T, OFF	3.2.2.76.
.MeasInput	Measuring input	1 , 2 , diff.	ditto
.Ipol	Polarization current	0 ... 1 ...±127	ditto
.Upol	Polarization voltage	0 ... 400 ...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON , OFF	ditto
.Temp	Titration temperature	-170.0... 25.0 ...500.0	ditto

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&MemoryCard

Object	Description	Input range	Reference
& Root			
MemoryCard	Administration of the memory card		
Recall	Load method	\$G	3.2.2.80.
Name	Method name	8 ASCII characters	ditto
Store	Save method	\$G	ditto
Name	Method name	8 ASCII characters	ditto
Delete	Delete method	\$G	ditto
Name	Method name	8 ASCII characters	ditto
ChangeDir	Change directory	\$G	3.2.2.81.
Name	Directory name	10 ASCII characters	ditto
Checksum	Checksum of directory	\$G	ditto
Value	Value of checksum	read only	ditto
CreateDir	Create new directory	\$G	3.2.2.82.
Name	Directory name	10 ASCII characters	ditto
DelDir	Delete directory	\$G	3.2.2.83.
Name	Directory name	10 ASCII characters	ditto
Backup	Backup of internal memory	\$G	3.2.2.84.
Name	Directory name	10 ASCII characters	ditto
Reload	Reload backup from the card	\$G	3.2.2.85.
Name	Directory name	10 ASCII characters	ditto
Format	Format the card	\$G	3.2.2.86.
CardLabel	Directory name		
Name	Value of checksum	8 ASCII characters	ditto
FreeMemory	Memory available	read only	3.2.2.87.
BatteryChange	Date for battery change	\$G	3.2.2.88.
Date	Date	YYYY-MM-DD	ditto
List	Lists		
Card	Directories on the card	read only	
.1	Method 1		
Name	Directory name	read only	3.2.2.89.
Bytes	Directory size in bytes	read only	ditto
.2	for each directory		
ActDir	Methods in the current directory		
.1	Method 1		
Name	Method name	read only	3.2.2.90.
Mode	Mode	read only	ditto
Quantity	Measured quantity	read only	ditto
DosUnit	Dosing element	read only	ditto
Bytes	Method size in bytes	read only	ditto
Checksum	Checksum of method	read only	ditto
.2	for each method		

- .RSSet1		Settings RS232, 1	\$G	3.2.2.108.
-	.Baud	Baud rate	300,600,1200,2400,4800, 9600 ,19200,38400,57600, 115200	ditto
-	.DataBit	Number of data bits	7, 8	ditto
-	.StopBit	Number of stop bits	1 , 2	ditto
-	.Parity	Parity	even, odd, none	ditto
-	.Handsh	Handshake	HWs , SWchar, SWline, none	ditto
-	.RSSet2	as for RS1		
- .ComVar		Values of common variables		
-	.C30	C30	0... ±999 999	3.2.2.109.
-	up to C39	0... ±999 999		
- .DosPrep		Parameters for the preparation of burets		
-	.PowerOnPrep	Warning after power on	ON, OFF	3.2.2.110.
-	.Report	Preparation report	ON, OFF	3.2.2.111.
-	.Select	Selection of dosing unit	internal D0 , external D1, external D2	3.2.2.112.
-	.D0	Internal Titrino buret		
-	.WarnInterval	Warning interval for preparation	5...9999, OFF	3.2.2.113.
-	.V	Volume	0... 3.5 ...99999.99	3.2.2.114.
-	.Repeat	Number of cycles	1... 2 ...9	3.2.2.115.
-	.DosRate	Dosing rate	0.01...150, max.	3.2.2.116.
-	.FillRate	Filling rate	0.01...150, max.	ditto
-	.D1	Buret D1		
-	.WarnInterval	Warning interval for preparation	5...9999, OFF	3.2.2.104.
-	.Select	Type of dosing unit	Dosimat , Dosino	3.2.2.117.
-	.Dosimat	Parameters for Dosimats		
-	.V	Volume	0... 3.5 ...99999.99	3.2.2.107.
-	.Repeat	Number of cycles	1... 2 ...9	3.2.2.108.
-	.DosRate	Dosing rate	0.01...150, max.	3.2.2.109.
-	.FillRate	Filling rate	0.01...150, max.	ditto
-	.Dosino	Parameters for Dosinos		
-	.Outlet	Expelling outlet	tip , flask	3.2.2.118.
-	.DosTubing	Size of dosing tubes		
-	.Length	Length	1... 40.0 ...999.9	3.2.2.119.
-	.Diam	Diameter	0.1... 2.0 ...9.9	ditto
-	.AspirTubing	Size of aspirating tubes		
-	.Length	Length	1... 25.0 ...999.9	ditto
-	.Diam	Diameter	0.1... 2.0 ...9.9	ditto
-	.DosRate	Dosing rate	0.01...150, max.	3.2.2.120.
-	.FillRate	Filling rate	0.01...150, max.	ditto
-	.D2	Buret D2 as for D1		

&Info

Object	Description	Input range	Reference
& Root			
Info	Current data		
- .Report	Transmission of formatted reports	\$G	3.2.2.130.
- .Select	Report type	configuration, parameters, smpl data, statistics, silo, calib C-fmla, def, user method, full , short, mplist, curve, derive, comb, meas crv, temp crv, adj para, scalc full, scalc srt, prep, calc, act dir, mem card, all, ff	ditto
- .CalibrationData	pH calibration data	\$G	3.2.2.131.
- .Inp1	For measuring input 1		
- .pHas	Asymmetry pH	0... 7.00 ...±20.00	ditto
- .Slope	Slope of electrode	0... 1.000 ...±9.999	ditto
- .Temp	Calibration temperature	-170.0... 25.0 ...500.0	ditto
- .Date	Date of calibration	read only	ditto
- .ElectrodeId	Id of the calibrated electrode	read only	ditto
- .Inp2	For measuring input 2, as for input 1		
- .Diff	For differential input, as for input 1		
- .PrepData	Preparation data		
- .D0	For the internal Titrino buret		
- .Date	Preparation date	read only	3.2.2.132.
- .Time	Preparation time	read only	ditto
- .D1	For buret D1		
- .Type	Buret type	read only	ditto
- .Date	Preparation date	read only	ditto
- .Time	Preparation time	read only	ditto
- .D2	For buret D2, as for D1		
- .Checksums	Checksums	\$G	3.2.2.133.
- .MPList	Checksum of meas.point list	read only	ditto
- .ActualMethod	Checksum of current method	read only	ditto
- .DetermData	Determination data	\$G	3.2.2.134.
- .Write	Read/write for several nods	ON, OFF	
- .ExV	Volume of Exchange/Dosing unit	read only/read + write	ditto
- .MPList	Measuring point list		
- .1	Measuring point 1		
- .Attribute	Attribute	read only/read + write	ditto
- .X	X coordinate	read only/read + write	ditto
- .Y	Y coordinate	read only/read + write	ditto
- .Z1	Z1 coordinate	read only/read + write	ditto
- .Z2	Z2 coordinate	read only/read + write	ditto
	for each measuring point		
- .TitrResults	Titration results		
- .RS	Calculated results		
- .1	1 st result		
- .Value	Value	read only	3.2.2.135.
	up to 9 results		

"Info",

"Info", continuation			
.EP	Endpoint		
.1	1 st result		
.V	Value	read only	ditto
.Meas	Measured value	read only	
.Mark	Mark if more than 1 EP per window	read only	
up to 9 EP's			
.Var	Variables C4X		
.C40	Start measured value	read only/read + write	ditto
.C41	Titration end volume	read only/read + write	
.C42	Titration time	read only/read + write	
.C43	Volume drift in SET/KFT	read only/read + write	
.C44	Titration temperature	read only/read + write	
.C45	Start volume	read only/read + write	
.C46	Asymmetry pH	read only	
.C47	Slope of electrode	read only	
.C48	Volume at maximum voltage	read only/read + write	
.C49	Volume at minimum voltage	read only/read + write	
.DTime	Time for drift corr.or dosing time	read only/read + write	
.FixEP	Fix EP or fix volume		
.51	C51		
.Value	Value	read only	3.2.2.136.
up to 59			
.pK	pK/HNP or fix times		
.61	C61		
.Value	Value	read only	ditto
up to 69			
.TempVar	Temporary variables C7X		
.C70	up to C79	read only/read + write	ditto
.TimeWin	Rates in time windows		
.81	C81		
.Mean	Rate	read only	3.2.2.137.
.Dev	Standard deviation of the rate	read only	ditto
up to 89			
.MeanRateC80	Mean rate C80		
.Mean	Rate	read only	ditto
.Dev	Standard deviation of the rate	read only	ditto
Statistics values			
.StatisticsVal			
.ActN	Number of results in chart	read only	3.2.2.138.
.1	1 st mean		
.Mean	Mean	read only	ditto
.Std	Absolute standard deviation	read only	ditto
.RelStd	Relative standard deviation	read only	ditto
up to 9 mean values			
Values of silo calculations			
.SiloCalc			
.C24	Values of variable C24		
.Name	Name	read only	3.2.2.139.
.Value	Value	read only	ditto
.Unit	Unit	read only	ditto
.C25	as for C24		
.C26	Values of variable C26		
.ActN	Number of single values	read only	ditto
.Mean	Mean value	read only	ditto
.Std	Absolute standard deviation	read only	ditto
.RelStd	Relative standard deviation	read only	ditto
.C27	as for C26		

- .ActualInfo	"Info", continuation		
- .Inputs	Current data		
- .Status	I/O Inputs		
- .Change	Line status	read only	3.2.2.140.
- .Clear	Change of line status	read only	ditto
- .Outputs	Clear change	\$G	ditto
- .Assembly	as for I/O Inputs		
- .CyclNo	From Assembly		
- .Counter	Cycle number	read only	3.2.2.141.
- .V	Assembly counter	read only	3.2.2.142.
- .Clear	Volume counter	read only	ditto
- .Meas	Clears counter	\$G	ditto
- .Titrator	Measured value	read only	3.2.2.143.
- .CyclNo	From Titrator		
- .V	Cycle number	read only	3.2.2.144.
- .Meas	Volume	read only	ditto
- .dVdt	Measured indicator voltage	read only	ditto
- .dMeasdt	Volume drift dV/dt	read only	ditto
- .dMeasdV	Measured value drift	read only	ditto
- .ERC	1st deviation of titration curve	read only	ditto
- .T	ERC from DET	read only	ditto
- .MeasPt	Temp.as secondary meas.	read only	ditto
- .Index	Entry in measuring point list		
- .X	Index of entry	read only	3.2.2.145.
- .Y	X coordinate	read only	ditto
- .Z1	Y coordinate	read only	ditto
- .Z2	Z1 coordinate	read only	ditto
- .EP	Z2 coordinate	read only	ditto
- .Index	EP entry		
- .X	Index of entry	read only	ditto
- .Y	X coordinate	read only	ditto
- .Oven	Y coordinate	read only	ditto
- .HeatTime	Oven data		
- .SampleTemp	Heating time	read only	3.2.2.146.
- .LowTemp	Sample temperature	read only	ditto
- .HighTemp	Lowest temperature	read only	ditto
- .GasFlow	Highest temperature	read only	ditto
- .UnitFlow	Gas flow	read only	ditto
- .Display	Unit of gas flow	read only	ditto
- .L1	Display		
- up to line 8	Text line 1	up to 32 ASCII char	3.2.2.147.
- .DelAll	Delete display	\$G	ditto
- .Comport	Comport		
- .Number	COM where PC is connected	read only	3.2.2.148.
- .Assembly	Assembly		
- CycleTime	Cycle time	read only	3.2.2.149.
- ExV	Volume of Exchange/Dosing unit	read only	ditto

&Assembly

Object	Description	Input range	Reference
& Root			
-			
- Assembly	Assembly control		
- .Bur	Buret		
- .Select	Selection of buret	internal D0 , external D1 external D2	3.2.2.150.
- .Empty	Empties the buret	\$G,\$S,\$H,\$C	3.2.2.151.
- .Prep	Prepares the buret	\$G,\$S,\$H,\$C	3.2.2.152.
- .Rates	Rates		
- .Forward	Forward rate		
- .Select	Type of rate control	digital, analog	3.2.2.153.
- .Digital	Digital rate	0...150, max.	ditto
- .Reverse	as for forward rate		
- .Select	Type of rate control	digital, analog	ditto
- .Digital	Digital rate	0...150, max.	ditto
- .Fill	Fill	\$G,\$H,\$C	3.2.2.154.
- .ModeDis	Dispensing	\$G,\$S,\$H,\$C	3.2.2.155.
- .Select	Type of dispensing control	volume , time	ditto
- .V	Volume to be dispensed	0.0001... 0.1 ...9999	ditto
- .Time	Time to dispense	0.25... 1 ...86 400	ditto
- .VStop	Limit volume	0.0001...9999, OFF	ditto
- .AutoFill	Filling after each increment	ON, OFF	ditto
- .Meas	Measuring		
- .Status	Measuring ON/OFF	ON, OFF	3.2.2.156.
- .MeasInput	Selection of measuring input	1 , 2, diff., lpol, Upol, Temp	ditto
- .Ipol	Polarization current	0... 1 ...±127	ditto
- .Upol	Polarization voltage	0... 400 ...±1270	ditto
- .Outputs	I/O outputs		
- .AutoEOD	Automatic output of EOD	ON , OFF	3.2.2.157.
- .SetLines	Set I/O lines	\$G	ditto
- .LO	Signal on LO	active,inactive,pulse, OFF	ditto
- up to L13			
- .ResetLines	Reset I/O lines	\$G	ditto
- .Stirrer	Stirrer	ON, OFF	3.2.2.158.

&Setup

Object	Description	Input range	Reference
& Root			
.			
└ Setup			
	Settings for the operating mode		
└ .Comport	Output of automatic info	1,2,1&2	3.2.2.159.
└ .Keycode	Send key code	ON, OFF	3.2.2.160.
└ .Tree	Sending format of path info		
└ .Short	Short format of path	ON, OFF	3.2.2.161.
└ .ChangedOnly	Paths of modified nodes only	ON, OFF	ditto
└ .Trace	Message on changed values	ON, OFF	3.2.2.162.
└ .Lock	Lock key functions		
└ .Keyboard	Lock all keyboard keys	ON, OFF	3.2.2.163.
└ .Config	Lock <CONFIG> key	ON, OFF	ditto
└ .Parameter	Lock <PARAM> key	ON, OFF	ditto
└ .SmplData	Lock <SMPL DATA> key	ON, OFF	ditto
└ .UserMeth	Lock functions		
└ .Recall	Lock "loading"	ON, OFF	ditto
└ .Store	Lock "saving"	ON, OFF	ditto
└ .Delete	Lock "deletion"	ON, OFF	ditto
└ .Display	Lock display function	ON, OFF	ditto
└ .Mode	Setting waiting intervals		
└ .StartWait	Waiting time after start	ON, OFF	3.2.2.164.
└ .FinWait	Waiting time after run	ON, OFF	ditto
└ .SendMeas	Automatic sending of measured values		
└ .SendStatus	Connect/disconnect sending	ON, OFF	3.2.2.165.
└ .Interval	Time interval	0.08...4...16200, MPList	ditto
└ .Select	Selection	Assembly, Titration	3.2.2.166.
└ .Assembly	From assembly		
└ .CyclNo	Cycle number	ON, OFF	3.2.2.167.
└ .V	Volume	ON, OFF	ditto
└ .Meas	Measured indicator voltage	ON, OFF	ditto
└ .Titration	From Titration		
└ .CyclNo	Cycle number	ON, OFF	3.2.2.168.
└ .V	Volume	ON, OFF	ditto
└ .Meas	Measured indicator voltage	ON, OFF	ditto
└ .dVdt	Volume drift dV/dt	ON, OFF	ditto
└ .dMeasdt	Measured value drift	ON, OFF	ditto
└ .dMeasdV	1st deviation of titration curve	ON, OFF	ditto
└ .ERC	ERC from DET	ON, OFF	ditto
└ .T	Temp.as secondary meas.	ON, OFF	ditto

"Setup", continuation			
.AutolInfo	Automatic message for changes		3.2.2.169.
.Status	Switch AutolInfo on/off	ON, OFF	ditto
.P	When mains is switched on	ON, OFF	ditto
.T	Titration infos		
.R	When "ready"	ON, OFF	ditto
.G	When method started	ON, OFF	ditto
.GC	When start is initiated	ON, OFF	ditto
.S	When stopped	ON, OFF	ditto
.B	Begin of method	ON, OFF	ditto
.F	End of process	ON, OFF	ditto
.E	Error	ON, OFF	ditto
.H	When "hold"	ON, OFF	ditto
.C	Continue after "hold"	ON, OFF	ditto
.O	Conditioning OK	ON, OFF	ditto
.N	Conditioning not OK	ON, OFF	ditto
.Re	Request after start	ON, OFF	ditto
.Si	Silo empty	ON, OFF	ditto
.M	Entry in measuring point list	ON, OFF	ditto
.EP	Entry in EP list	ON, OFF	ditto
.RC	Recalculation of results done	ON, OFF	ditto
.C	Comport infos		
.B1	When COM1 sends a report	ON, OFF	ditto
.R1	When COM1 is ready again	ON, OFF	ditto
.B2	When COM2 sends a report	ON, OFF	ditto
.R2	When COM2 is ready again	ON, OFF	ditto
.I	Changing an I/O input	ON, OFF	ditto
.O	Changing an I/O output	ON, OFF	ditto
.Graphics	Changing the curve output		
.Grid	Grid on curve	ON, OFF	3.2.2.170.
.Frame	Frame on curve	ON, OFF	ditto
.Scale	Type of depending axis	Full, Auto	ditto
.Recorder	Length of axes		
.Right	Length of meas value axis	0.2... 0.5 ...1.00	ditto
.Feed	Length of paper drive axis	0.01... 0.05 ...1.00	ditto
.PowerOn	RESET (power on)	\$G	3.2.2.171.
.Initialise	Set default values	\$G	3.2.2.172.
.Select	Selection of branch	ActMeth, Config, Silo, Calib	
		Assembly, Setup, All	ditto
.RamInit	Initialization of working mem.	\$G	3.2.2.173.
.InstrNo	Device Identification	\$G	3.2.2.174.
.Value	Input of device identification	8 ASCII characters	ditto

&Diagnose

Object	Description	Input range	Reference
& Root . ├ Diagnose └ .Report	Diagnose Output of adjustment parameters	\$G	3.2.2.175.

3.2.2 Description of the remote control commands

3.2.2.1. Mode \$G, \$S, \$H, \$C
 Start and stop (\$G, \$S) or hold of the current method (3.2.2.3) with \$H and continue with \$C.

\$G also serves to continue after inquiries of identifications and sample size after the start (see 3.2.2.26) as well as after inquiries of calibration temperature and pH values of buffers (see 3.2.2.68 and 3.2.2.69).

3.2.2.2. Mode.QuickMeas \$G, \$S
 Start and stop of a measurement in the basic mode with the parameters (measured quantity, measuring input) of the current method. Corresponds to the <meas/hold> key. In TIP and DOS, the measured quantity is selected with &Mode.Parameter, see 3.2.2.76 and 3.2.2.56.
 With an ongoing measurement, the current mode can be started. This stops the measurement automatically.

3.2.2.3. Mode.Select DET, MET, SET, KFT, STAT, DOS, DOC, MEAS, CAL, TIP

Mode.DETQuantity	pH, U, Ipol, Upol
Mode.METQuantity	pH, U, Ipol, Upol
Mode.SETQuantity	pH, U, Ipol, Upol
Mode.MEASQuantity	pH, U, Ipol, Upol, T
Mode.KFTQuantity	Ipol, Upol
Mode.STATQuantity	pH, U, Ipol, Upol
Mode.DOCQuantity	pH, U, Ipol, Upol

Selection of the standard mode. Mode and the measured quantity belong to the complete selection.

If a method is selected from the method memory, the nodes &Mode.Select and &Mode.XXXQuantity are overwritten with mode and measured quantity of the corresponding user method.

3.2.2.4. Mode.Name read only
 Name of the current method in the working memory. \$Q sends 8 ASCII characters. Standard methods carry the name *****. The node can be set read + write, see 3.2.2.134.

3.2.2.5. Mode.Def.Formulas.1.Formula EPX, CXX, RSX, +, -, *, /, (,)
Mode.Def.Formulas.1.TextRS up to 8 ASCII characters
Mode.Def.Formulas.1.Decimal 0...2...5
Mode.Def.Formulas.1.Unit up to 6 ASCII characters
Mode.Def.Formulas.1.Limits ON, OFF
Mode.Def.Formulas.1.LoLim 0...±999 999
Mode.Def.Formulas.1.UpLim 0...±999 999
Mode.Def.Formulas.1.Output active, pulse, OFF
Mode.Def.Formulas.2.Formula
 etc. up to .9

Entry of formulas. Rules for formula entry, see page 74.
 Example: "(EP2-EP1)*C01/C00"

In addition to the formula, a text for result output, the number of decimal places and a unit for the result output can be selected. "No unit" is selected with the blank string.

In place of "RSX", a result name may be entered (.TextRS). This name is outputted in the report full, short, scalc full and scalc srt. It is used for the result and the corresponding mean value.

The limit control for results can also be activated. If a result is out of limit, a message appears in the result report, E196 is sent, and output line L13 can be set.

3.2.2.6. Mode.Def.SiloCalc.Assign.C24 RSX, EPX, CXX
Mode.Def.SiloCalc.Assign.C25 RSX, EPX, CXX
Mode.Def.SiloCalc.MatchId id1, id1&2, all, **OFF**

.Assign.C2X: Assignment to store results in the silo as C2X.

.MatchId: Indication which sample identification(s) have to match so that the results can be combined.

3.2.2.7. Mode.Def.ComVar.C30 RSX, MNX, EPX, CXX
Mode.Def.ComVar.C31
 etc., up to **.C39**

Assignment of common variables.

The values of the common variables are to be found in &Config.ComVar. They can be viewed and entered there, see 3.2.2.109.

3.2.2.8. Mode.Def.Report.Assign1
Mode.Def.Report.Assign2
 DET: full, short, mplist, curve, derive comb, scalc full, scalc srt, calc, param, calib, ff
 MET, SET, MEAS: full, short, mplist, curve, scalc full, scalc srt, calc, param, calib, ff
 KFT: full, short, mplist, curve, scalc full, scalc srt, calc, param, ff
 CAL: full, short, scalc full, scalc srt, calc, param, calib, ff
 STAT, DOC: full, short, mplist, curve, meas crv, temp crv, scalc full, scalc srt, calc, param, calib, ff
 DOS: full, short, mplist, curve, meas crv, temp crv, scalc full, scalc srt, calc, param, ff
 TIP: full, short, scalc full, scalc srt, calc, param, ff

Definition of the report sequence, which is outputted automatically at the end of the determination. Entries of more than one block have to be separated with ",".

.Assign1: Output to COM1 of the Titrimo. Identical for COM2.

3.2.2.9. Mode.Def.Mean.1.Assign RS1, RSX, EPX, CXX
Mode.Def.Mean.2.Assign
 etc., up to **.9**

Assignment of the statistics calculations. Valid assignments are a requirement for statistics calculations. In addition, the statistics calculation must be switched on, see 3.2.2.24. Rules for statistics calculations see page 77.

3.2.2.10. Mode.Def.TempVar.C70 RSX, EPX, CXX
etc. up to **.C79**

Assignment of temporary variables in a submethod for calculations in TIP.

3.2.2.11. Mode.CFmla
Mode.CFmla.1.Value 0...±999 999
Mode.CFmla.2.Value
etc., up to **.19**

Calculation constants specific to a method. Stored in the method memory of the Titrino. Operands specific to the sample (3.2.2.117 and 3.2.2.135) and values of common variables (3.2.2.109) on the other hand are not stored with the methods.

3.2.2.12. Mode.Parameter.TitrPara.MptDensity 0...4...9
Mode.Parameter.TitrPara.MinIncr 0...10.0...999.9
Mode.Parameter.TitrPara.VStep 0...0.10...9.999

.MptDensity: Parameter for DET: Measuring point density.

.MinIncr: Parameter for DET: Minimum increment in uL. If the minimum increment is set to 0, measured values are stored vs. time.

.VStep: Parameter for MET: Volume increment in mL. With "0", there is no dispensing and measured values vs. time are entered in the measuring point list.

3.2.2.13. Mode.Parameter.TitrPara.DosRate 0.01...150, **max.**
Parameters for DET and MET: Dispensing rate for the volume increments in mL/min. Max. means maximum possible dispensing rate with the Exchange Unit in current use.

3.2.2.14. Mode.Parameter.TitrPara.SignalDrift pH,U,Ipol:0.5...50...999, OFF
Upol: 0.05...50...99.9, OFF
Mode.Parameter.TitrPara.UnitSigDrift read only
Mode.Parameter.TitrPara.EquTime 0...26...9999, OFF

Parameters for DET and MET: Criteria for the measured value acquisition.

Measured value drift in mV/min (with pH, U, Ipol) or uA/min (with Upol), equilibration time in s. OFF means that the corresponding criterion is switched off. If both criteria are OFF, the measured values are acquired immediately after dispensing.

If the equilibration time has never been edited, it is automatically calculated by the instrument to match the drift, see page 20. After it has been edited once, it remains in force with the set value.

3.2.2.15. Mode.Parameter.TitrPara.StartV.Type abs., rel., **OFF**
Mode.Parameter.TitrPara.StartV.V 0...999.99
Mode.Parameter.TitrPara.StartV.Factor 0...±999 999
Mode.Parameter.TitrPara.StartV.Rate 0.01...150, **max.**

Parameters for DET, MET, SET, KFT, STAT: Start volume.

If an absolute start volume (abs.) has been selected, the volume in mL is valid.

A relative start volume (rel.) is dispensed as a function of the sample size:

Start volume in mL = $\text{smpl size} \cdot \text{factor}$

The factor is valid.

The dispensing rate in mL/min applies to both cases. Max. means maximum possible dispensing rate with the Exchange Unit in current use.

3.2.2.16. Mode.Parameter.TitrPara.Pause 0...999 999
Parameters for DET, MET, SET, KFT, STAT: Pause time in s. Is waited off after the dispensing of the start volume.

3.2.2.17. Mode.Parameter.TitrPara.DosUnit internal D0, external D1/D2
Mode.Parameter.DosPara.DosUnit internal D0, external D1/D2
.TitrPara.DosUnit: Parameter for DET, MET, SET, KFT, STAT. Selection of dosing unit.
.DosPara.DosUnit: Parameter for DOS, DOC. Selection of dosing unit.

3.2.2.18. Mode.Parameter.TitrPara.MeasInput 1, 2, diff.
Mode.Parameter.TitrPara.Ipol -127...1...+127
Mode.Parameter.TitrPara.Upol -1270...400...+1270
Mode.Parameter.TitrPara.PolElectrTest ON, OFF

Parameters for DET, MET, SET, KFT STAT:

Selection of the measuring input; valid with measured quantities pH and U.

"diff." means differential amplifier, see page 215.

With Ipol, the inquiries for the polarization current in uA (Ipol) and

.PolElectrTest are valid.

With Upol, the inquiry for the polarization voltage in mV (Upol) is valid. Entry in steps of 10 mV.

Besides .PolElectrTest is valid.

If the test for polarized electrodes is switched on, it is performed on change-over from the inactive state to an active state (titration or conditioning).

3.2.2.19. Mode.Parameter.TitrPara.Temp -170.0...25.0...500.0
Parameters for DET, MET, SET, KFT, STAT: Titration temperature in °C. If a Pt100 or Pt1000 is connected, the temperature is measured continuously and the parameter .Temp is updated.
The temperature is used for the temperature correction in pH measurements.

3.2.2.20. Mode.Parameter.StopCond.VStop.Type abs., rel., OFF
Mode.Parameter.StopCond.VStop.V 0...99.99...9999.99
Mode.Parameter.StopCond.VStop.Factor 0...±999 999

Parameters for DET, MET, SET, KFT, STAT, DOS, DOC: Stop volume. (With DOS the default value for .Type is "OFF" and the input range for V is 0...99999.99.)

If an absolute stop volume (abs.) has been selected, the volume in mL is valid.

A relative stop volume (rel.) is dispensed as a function of the sample size:

Mode.Parameter.Evaluation.Recognition.Selected all, greatest,
last, window, OFF

Mode.Parameter.Evaluation.Recognition.Window.1.LowLim
pH: 0...±20.00, OFF
U, Ipol: 0...±2000, OFF
Upol: 0...±200.0, OFF

Mode.Parameter.Evaluation.Recognition.Window.1.UpLim
Input range as LowLim

etc. up to 9 windows

Mode.Parameter.Evaluation.FixEP.1.Value pH: 0...±20.00, OFF
U; Ipol: 0...±2000, OFF
Upol: 0...±200.0, OFF

etc. up to 9 fix EP's

Mode.Parameter.Evaluation.pK ON, OFF

Parameters for DET and MET: Evaluation of the EP's, see page 21.

.EPC: EP criterion in pH, in mV (with U and Ipol) or in uA (with Upol).

.Recognition.Selected: EP recognition.

all: All endpoints found are recognized.

great: Only the largest EP is recognized.

last: Only the last EP is recognized.

window: Only EP's that lie within set windows are recognized.

OFF: The EP evaluation is switched off.

.Recognition.Window.1.LowLim: Lower limit for window in pH, mV (with U and Ipol) or uA (with Upol).

.Recognition.Window.1.UpLim: Upper limit for window in pH, mV (with U and Ipol) or uA (with Upol).

Windows are opened until the lower limit is set to OFF. For every expected EP, an individual window must be set, see page 23.

.FixEP.1.Val: Fix-EP's in pH, mV (for U, Ipol) resp. uA (for Upol). Fix EP's are evaluated until the setting OFF is found.

.pK: pK or HNP evaluation. Possible only in pH and U titrations.

3.2.2.26. Mode.Parameter.Presel.IReq id1, id1&2, all, OFF
Mode.Parameter.Presel.SReq value, unit, all, OFF

Parameters for DET, MET, SET, KFT, STAT, DOS, DOC, MEAS: Automatic inquiry after the start of the determination. From such an inquiry, the determination continues if the requested entry/entries is/are made, e.g.

&SmplData.OFFSilo.Id1 (see 3.2.2.122) or with &M \$G, see 3.2.2.1.

\$H is not possible in requests.

3.2.2.27. Mode.Parameter.Presel.LimSmplSize.Status ON, OFF
Mode.Parameter.Presel.LimSmplSize.LoLim 0.0...999 999
Mode.Parameter.Presel.LimSmplSize.UpLim 0.0...999 999

Limit control for the sample size.

3.2.2.28. Mode.Parameter.Presel.ActPuls ON, OFF
for SET, KFT: first, all, cond., OFF

Output of a pulse on the I/O line "Activate", see page 224.

3.2.2.29. Mode.Parameter.SET1.EP pH: 0...±20.00, OFF

751 GPD Titrimo

U, Ipol: 0...±2000, **OFF**
 Upol: 0...±200.0, **OFF**
Mode.Parameter.SET1.UnitEp read only
 Parameters for SET: Setting the 1st endpoint as pH value, in mV (with U and Ipol) resp. uA (with Upol). The corresponding unit can be read with .UnitEP. If the value is on "OFF", no further nodes will appear from SET1.

3.2.2.30. Mode.Parameter.SET1.Dyn pH: 0.01...20.00, **OFF**
 U, Ipol: 1...2000, **OFF**
 Upol: 0.1...200.0, **OFF**
Mode.Parameter.SET1.UnitDyn read only
Mode.Parameter.SET1.MaxRate 0.01...**10**...150, max.
Mode.Parameter.SET1.MinRate 0.01...**25.0**...999.9
 Parameters for SET: Control parameters, see page 31.
 .Dyn: Dynamics, control range in pH, mV (with U and Ipol) or uA (with Upol). The corresponding unit can be read with .UnitDyn.
 .MaxRate: Maximum allowed titration rate in mL/min. Max. means maximum possible rate with the Exchange Unit in current use.
 .MinRate: Minimum titration rate in ul/min.

3.2.2.31. Mode.Parameter.SET1.Stop.Type drift, time
Mode.Parameter.SET1.Stop.Drift 1...**20**...999
Mode.Parameter.SET1.Stop.Time 0...**10**...999, inf
Mode.Parameter.SET1.Stop.StopT 0...99 999, **OFF**
 Parameters for SET: Type and size of the stop criterion of the titration.
 .Type: Type of stop criterion after stop drift or switch-off delay time.
 .Drift: Stop drift in ul/min. Applies when "drift" has been selected.
 .Time: Switch-off delay time in s. Applies when "time" has been selected. "inf" means infinite.
 .StopT: Stop time in s. Applies when "time" has been selected and the value of .Time is set to "inf.".

3.2.2.32. Mode.Parameter.TitrPara.Direction +, -, **auto**
 Parameters for SET, KFT, STAT: Titration direction.
 "auto" means the titration direction is determined automatically by the instrument. If 2 EP's have been set in a SET titration, the titration direction is given by the two EP's. The entry of the titration direction is then invalid.

3.2.2.33. Mode.Parameter.TitrPara.XPause 0...999 999
 Parameter for SET, KFT: Pause time in s. Runs before dosing the start volume.

3.2.2.34. Mode.Parameter.TitrPara.ExtrT 0...999 999
 Parameter for SET, KFT: Extraction time in s.

3.2.2.35. Mode.Parameter.TitrPara.TDelta 1...2...999 999

Parameter for SET, KFT, STAT: Time interval in s for the entry of a measurement point in the list of measured points.

3.2.2.36. Mode.Parameter.Presel.Cond ON, **OFF**
Mode.Parameter.Presel.DriftDisp ON, **OFF**
Mode.Parameter.Presel.DCor.Type auto, man., **OFF**
Mode.Parameter.Presel.DCor.Value 0.0...99.9

Parameters for SET, KFT:

.Cond: Conditioning ON/OFF

.DriftDisp: Drift display during conditioning ON/OFF.

.DCor.Type: Type of drift take-over for the drift correction. auto: Take-over of the drift value at start.

.DCor.Value: Drift value for the manual drift correction.

3.2.2.37. Mode.Parameter.CtrlPara.EP pH: 0...±20.00, **OFF**
U, Ipol: 0...±2000, **OFF**
Upol: 0...±200.0, **OFF**
Mode.Parameter.CtrlPara.UnitEp read only

Parameters for KFT, STAT: Setting of the EP, resp. Control point as a pH value, in mV (with U and Ipol) or uA (with Upol). The corresponding unit can be read with .UnitEP.

With KFT the standard value for Ipol is 250 mV, for Upol 25 uA.

3.2.2.38. Mode.Parameter.CtrlPara.Dyn pH: 0.01...20.00, **OFF**
U, Ipol: 1...2000, **OFF**
Upol: 0.1...200.0, **OFF**
Mode.Parameter.CtrlPara.UnitDyn read only
Mode.Parameter.CtrlPara.MaxRate 0.01...10...150, max.
Mode.Parameter.CtrlPara.MinIncr 0.1...9.9, **min.**

Parameters for KFT, STAT: Control parameters.

.Dyn: Dynamics (control range) in pH, mV (with U and Ipol) or uA (with Upol). The corresponding unit can be read with .UnitDyn. Standard with KFT Ipol 100, with KFT Upol 10.

.MaxRate: Maximum allowed titration rate in ml/min. Max. means maximum possible rate with mounted Exchange Unit.

.MinIncr: Parameter for KFT. Minimum volume increment in ul.

3.2.2.39. Mode.Parameter.Presel.Oven COM1, COM2, **no**

Parameter for KFT.

If an Oven is connected, its results will be incorporated into the result report of the Titrino.

If there is no Oven connected via RS232, this parameter has to be on "no".

3.2.2.40. Mode.Parameter.CtrlPara.MinRate 0.01...25.0...9999

Parameter for STAT: Control parameters.

.MinRate: Minimum titration rate in ul/min.

3.2.2.41. Mode.Parameter.TitrPara.TStart 0...999 999
Mode.Parameter.TitrPara.StartMeas pH: 0.00...±20.00, **OFF**
 U, Ipol: 0...±2000, **OFF**
 Upol: 0.0...±200.0, **OFF**
Mode.Parameter.TitrPara.RateStart 0.01...150, **OFF**

Parameters for STAT: Start of the measured point acquisition after the start condition has been met.

.TStart: Start time in s.

.StartMeas: Start measured value.

.RateStart: Start rate in mL/min. The current rate must be less than the start rate. This criterion does not become effective until 10 s after the start.

3.2.2.42. Mode.Parameter.StopCond.TStop.Type abs.,rel.,delta,delay,**OFF**
Mode.Parameter.StopCond.TStop.Time 0...999 999
Mode.Parameter.StopCond.TStop.Factor 0...±999 99
Mode.Parameter.StopCond.TStop.Delta 0...999 999
Mode.Parameter.StopCond.TStop.Delay 0...999 999

Parameters for STAT: Stop time.

If an absolute stop time (abs.) has been selected, the time is in s (Time).

A relative stop time (rel.) is calculated as a function of the sample size:

Stop time in s = sample size * factor

The factor selected (Factor) applies.

delta: Time in s after control point is reached for the first time.

delay: Time in s after last dosing step.

OFF means that the criterion is not monitored.

3.2.2.43. Mode.Parameter.StopCond.RateStop 0.01...150, **OFF**

Parameters for STAT: Stop rate in ml/min. The current rate must be less than the start rate. This criterion does not become effective until 10 s after the start.

3.2.2.44. Mode.Parameter.Evaluation.TimeWin.1.LowLim 0...999 999, **OFF**
Mode.Parameter.Evaluation.TimeWin.1.UpLim 0...999 999, **OFF**

Parameters for STAT: Evaluation of rates in time windows. Entry of the times in s. The evaluated rates are available as variables C8X (X = 1...9).

3.2.2.45. Mode.Parameter.Evaluation.FixVol.1.Value 0...999 999, **OFF**

Parameters for STAT: Interpolation of the volume at preset times. Entry of the times in s. The interpolated volumes are available as variables C5X.

3.2.2.46. Mode.Parameter.Evaluation.FixTime.1.Value 0.01...1.00, **OFF**

Parameters for STAT: Interpolation of the time at preset fraction of the final volume. Entry as part of V(tot). The interpolated volumes are available as variable C6X.

3.2.2.47. Mode.Parameter.Monitoring.MeasVal.Status ON, **OFF**

Mode.Parameter.Monitoring.MeasVal.LowLim pH: 0...±20.00
Mode.Parameter.Monitoring.MeasVal.UpLim U, Ipol: 0...±2000
 Upol: 0...±200.0

Mode.Parameter.Monitoring.MeasVal.Action end, hold, wait, **none**

Parameters for STAT, DOC: Monitoring of limit values. Measured values (pH, U or I) are entered in the list of measured points only if the monitoring is active.

On limit value violation, an action can be selected:

end: Abort of the determination.

hold: Interrupt dosing until the sequence is continued by manual switching (& m \$C).

wait: Interrupt dosing until the limit value is no longer violated, then continue automatically.

3.2.2.48. Mode.Parameter.Monitoring.Rate.Status ON, **OFF**
Mode.Parameter.Monitoring.Rate.LowLim **0.000...150**
Mode.Parameter.Monitoring.Rate.UpLim **0.000...150**
Mode.Parameter.Monitoring.Rate.Action end, hold, wait, **none**

Parameters for STAT: Monitoring of limit values of the rate. Entry in ml/min.

On limit value violation, an action can be selected:

end: Abort of the determination.

hold: Interrupt dosing until the sequence is continued by manual switching (& m \$C).

wait: Interrupt dosing until the limit value is no longer violated, then continue automatically. On violation of the lower limit, the sequence can no longer be continued as the rate would become even less than that when dosing was interrupted.

3.2.2.49. Mode.Parameter.Monitoring.Temp.Status ON, **OFF**
Mode.Parameter.Monitoring.Temp.LowLim **-170.0...500.0**
Mode.Parameter.Monitoring.Temp.UpLim **-170.0...500.0**
Mode.Parameter.Monitoring.Temp.Action end, hold, wait, **none**

Parameters for STAT, DOS, DOC: Monitoring of the limit values of the temperature in °C. Temperature values are entered in the list of measured points only if the monitoring is active.

On limit value violation, an action can be selected:

end: Abort of the determination.

hold: Interrupt dosing until the sequence is continued by manual switching (& m \$C).

wait: Interrupt dosing until the limit value is no longer violated, then continue automatically.

3.2.2.50. Mode.Parameter.Monitoring.L10Output.AssignOutput
 meas, temp, rate, all, **none**
Mode.Parameter.Monitoring.L10Output.Limit upper, lower, **any**
Mode.Parameter.Monitoring.L10Output.Output active, **pulse**

Parameters for STAT, DOS, DOC: Output of a signal on L10 output (pin 8) of the remote socket on limit value violation.

The output can be assigned a definite or all monitoring value(s) as well as the violated limit and the type of signal to be outputted determined.

3.2.2.51. Mode.Parameter.Presel.RateDisp ON, OFF
 Parameters for STAT, DOC: Display of the current rate during the determination.

3.2.2.52. Mode.Parameter.DosPara.Type volume, time, rate
Mode.Parameter.DosPara.Volume.Volume 0.000...10...99999.99
Mode.Parameter.DosPara.Volume.DisType time, rate
Mode.Parameter.DosPara.Volume.Rate 0.001...150, max.
Mode.Parameter.DosPara.Volume.Time 1...100...999 999
Mode.Parameter.DosPara.Time.Time 1...100...999 999
Mode.Parameter.DosPara.Time.DisType volume, rate
Mode.Parameter.DosPara.Time.Rate 0.001...150, max.
Mode.Parameter.DosPara.Time.Volume 0.000...10...99999.99
Mode.Parameter.DosPara.Rate.Rate 0.001...150, max.
Mode.Parameter.DosPara.Rate.StopType volume, time
Mode.Parameter.DosPara.Rate.Time 1...100...999 999
Mode.Parameter.DosPara.Rate.Volume 0.000...10...99999.99

Parameters for DOS: In the equation Rate = Volume/Time, two parameters are preset, the third is calculated. Depending on the entry of the type, the following information is required: .Volume, .Time or .Rate.
 Entries: .Volume in ml, .Time in s, .Rate in ml/min.

3.2.2.53. Mode.Parameter.DosPara.Pause 0...999 999
 Parameters for DOS: Pause time in s.

3.2.2.54. Mode.Parameter.DosPara.TDelta 1...10...999 999
 Parameters for DOS, DOC: Time interval for the entry of a measured point in the list of measured points.

3.2.2.55. Mode.Parameter.DosPara.Temp -170.0...25.0...500.0
 Parameters for DOS, DOC: Temperature in °C. The temperature is measured continuously and updated if a T sensor is attached.
 The temperature is used for temperature compensation in pH measurements.

3.2.2.56. Mode.Parameter.Monitoring.MeasVal.MeasModepH,U,Ipol,Upol OFF
Mode.Parameter.Monitoring.MeasVal.MeasInput 1, 2, diff.
Mode.Parameter.Monitoring.MeasVal.PCurrent 0...1...±127
Mode.Parameter.Monitoring.MeasVal.PVoltage 0...400...±1270
Mode.Parameter.Monitoring.MeasVal.PolElectrTest ON, OFF
Mode.Parameter.Monitoring.MeasVal.pH.LowLim 0...±20.00
Mode.Parameter.Monitoring.MeasVal.pH.UpLim 0...±20.00
Mode.Parameter.Monitoring.MeasVal.U.LowLim 0...±2000
Mode.Parameter.Monitoring.MeasVal.U.UpLim 0...±2000
Mode.Parameter.Monitoring.MeasVal.Ipol.LowLim 0...±2000
Mode.Parameter.Monitoring.MeasVal.Ipol.UpLim 0...±2000

Mode.Parameter.Monitoring.MeasVal.Upol.LowLim 0...±200.0
Mode.Parameter.Monitoring.MeasVal.Upol.UpLim 0...±200.0
Mode.Parameter.Monitoring.MeasVal.Action end, hold, wait, none

Parameters for DOS: Selection of the measured quantity for measured value monitoring. Measured values (pH, U and I) are entered in the list of measured points only when the monitoring is active.

Selection of the measurement input (MeasInput) applies with measured quantities pH and U. diff. means differential amplifier, see page 215.

With Ipol, the requests for the polarization current in uA (Ipol) and .PolElectrTest apply.

With Upol, the request for the polarization voltage in mV (Upol) applies. Entry in steps of 10 mV. .PolElectrTest also applies. If the test for polarized electrodes is switched on (ON), it is performed on the change from the inactive basic mode to an active mode.

On limit value violation, an action can be selected:

stop: Abort of the determination.

hold: Interrupt dosing until the sequence is continued by manual switching (&m \$C).

wait: Interrupt dosing until the limit value is no longer violated, then continue automatically.

3.2.2.57. Mode.Parameter.DosPara.BeginMeas pH: 0...±20.00, **init**
 U, Ipol: 0...±2000, **init**
 Upol: 0...±200.0, **init**

Mode.Parameter.DosPara.EndMeas

as with .Begin, but OFF instead of init

Mode.Parameter.DosPara.UnitMeas

read only

Parameters for DOS: Setting of ramp start and ramp end. init means ramp start with default measured value.

3.2.2.58. Mode.Parameter.DosPara.SweepTime 0...**300**...999 999
 Sweep time for the ramp in s.

3.2.2.59. Mode.Parameter.DosPara.Dyn pH: 0.01...**0.25**...20.00, OFF
 U: 1...**15**...2000, OFF
 Ipol: 1...**6**...2000, OFF
 Upol: 0.1...**10**...200.0, OFF

Mode.Parameter.DosPara.UnitDyn

read only

Mode.Parameter.DosPara.MaxRate

0.01...150, **max.**

Mode.Parameter.DosPara.MinRate

0.01...**5.0**...9999

Parameters for DOC: Control parameters, see page 68.

.Dyn: Dynamics in pH, mV (with U and Ipol) or uA (with Upol). The corresponding unit can be read with .UnitDyn.

.MaxRate: Maximum allowed titration rate in ml/min. Max. means maximum possible rate with the mounted Exchange Unit.

.MinRate: Minimum titration rate in ul/min.

3.2.2.60. Mode.Parameter.DosPara.Direction +, -, **auto**

Parameters for DOC: Titration direction.

auto means the titration direction is determined automatically by the instrument. The direction is used only when the ramp start is the same as init.

3.2.2.61.	Mode.Parameter.DosPara.StartV.Type	abs., rel., OFF
	Mode.Parameter.DosPara.StartV.V	0...999.99
	Mode.Parameter.DosPara.StartV.Factor	0...±999.999
	Mode.Parameter.DosPara.StartV.Rate	0.01...150, max.
	Mode.Parameter.DosPara.Pause	0...999.999

Parameters for DOC: Start volume and pause time in s. Runs after dispensing of the start volume.

If an absolute start volume (abs.) has been selected, the volume is in ml (Volume).

A relative start volume (rel.) is dispensed as a function of the sample size: Start volume in ml = sample size * factor (the factor selected (Factor))

The dosing rate in ml/min holds for both cases. Max. means the maximum possible dosing rate with the Exchange Unit currently mounted.

3.2.2.62.	Mode.Parameter.DosPara.MeasInput	1, 2, diff.
	Mode.Parameter.DosPara.Ipol	0...1...±127
	Mode.Parameter.DosPara.Upol	0...400...±1270
	Mode.Parameter.DosPara.PolElectrTest	ON, OFF

Parameters for DOC:

Selection of the measurement input (MeasInput) applies to the measured quantities pH and U. diff. means differential amplifier, see page 215.

With Ipol, the requests for the polarization current in uA (Ipol) and .PolElectrTest apply.

With Upol, the request for the polarization voltage in mV (Upol) applies. Entry in steps of 10 mV. .PolElectrTest also applies.

If the test for polarized electrodes is switched on (ON), it will be performed on the change from the inactive basic mode to an active mode.

3.2.2.63.	Mode.Parameter.Measuring.SignalDrift	
	pH, U, Ipol, T:	0.5...999, OFF
	Upol:	0.05...99.9, OFF
	Mode.Parameter.Measuring.UnitSigDrift	read only
	Mode.Parameter.Measuring.EquTime	0...9999, OFF

Parameters for MEAS: Criteria for the measured value acquisition. Measured value drift in mV/min (with pH, U, Ipol, T), uA/min (with Upol), resp. °C/min (with T). Equilibration time in s. OFF means that the corresponding criterion is switched off. If both criteria are OFF, the measurement continues indefinitely. If the equilibration time has never been edited, it is automatically calculated by the instrument to match the drift, see page 69. After it has been edited once, it remains in force with the set value.

3.2.2.64.	Mode.Parameter.Measuring.MeasInput	1, 2, diff.
	Mode.Parameter.Measuring.Ipol	±127...1...+127
	Mode.Parameter.Measuring.Upol	±1270...400...+1270
	Mode.Parameter.Measuring.PolElectrTest	ON, OFF

Parameters for MEAS:

Selection of the measuring input; valid with measured quantities pH and U.
"diff." means differential amplifier, see page 215.

With Ipol, the inquiries for the polarization current in uA (Ipol) and .PolElectrTest are valid.

With Upol the inquiry for the polarization voltage in mV (Upol) is valid. Entry in steps on 10 mV.

Besides .PolElectrTest is valid.

If the test for polarized electrodes is switched on, it is performed on change-over from the inactive state to the measurement.

3.2.2.65. Mode.Parameter.Measuring.Temp -170.0...**25.0**...500.0

Parameters for MEAS: Measurement temperature in °C. If a Pt100 or Pt1000 is connected, the temperature is measured.

The temperature is used for the temperature correction in pH measurements.

3.2.2.66. Mode.Parameter.Measuring.TDelta 1...**2**...999 999

Parameters for MEAS: Time interval in s for the acquisition of a measured value into the measuring point list.

3.2.2.67. Mode.Parameter.Calibration.MeasInput 1, 2, diff.

Parameters for CAL: Selection of the measuring input. "diff." means differential amplifier, see page 215.

3.2.2.68. Mode.Parameter.Calibration.CalTemp -20.0...**25.0**...120.0

Parameters for CAL: Calibration temperature in °C. If a Pt 100 or Pt1000 is connected, the temperature is measured.

3.2.2.69. Mode.Parameter.Calibration.Buffer.1.Value 0...**7.00**...±20.00

Mode.Parameter.Calibration.Buffer.2.Value 0...**4.00**...±20.00, OFF
etc. up to 9 buffers

Parameters for CAL: pH of buffers. The first buffer which is set to "OFF" determines the number of buffers in the calibration.

3.2.2.70. Mode.Parameter.Calibration.SignalDrift 0.5...**2**...999, OFF

Mode.Parameter.Calibration.EquTime 0...**110**...9999, OFF

Parameters for CAL: Criteria for measured value acquisition. Measured value drift in mV/min, equilibration time in s. OFF means that the corresponding criterion is switched off. If both criterions are on OFF, the measured value is acquired immediately.

If the equilibration time has never been edited, it is automatically calculated by the instrument to match the drift, see page 20. After it has been edited once, it remains in force with the set value.

3.2.2.71. Mode.Parameter.Calibration.Electrodeld up to 8 ASCII char

Parameters for CAL: Electrode identification. It is classified under calibration data, see 3.2.2.131.

3.2.2.72. Mode.Parameter.Calibration.SmplChanger ON, OFF
 Parameters for CAL: Calibration at Titrino.
 With "ON", there are no hold points in the calibration sequence for entries, the first buffer is measured directly.

3.2.2.73. Mode.Parameter.Calibration.ActPulse first, all, OFF
 Parameters for CAL: Output of a pulse on the I/O line "Activate", see page 224.

3.2.2.74. Mode.Parameter.Sequence.X.Select method, pause, L4 output, L6 output, info, prep, stirrer, OFF
 Parameters for TIP: Selection of an element for step X (X = 1...30). For the parameters of the elements see 3.2.2.75.

3.2.2.75. Mode.Parameter.Sequence.X.Method Method name
Mode.Parameter.Sequence.X.Pause 0...999 999, INF
Mode.Parameter.Sequence.X.L4Output active, inactive, pulse, OFF
Mode.Parameter.Sequence.X.L6Output as for L4
Mode.Parameter.Sequence.X.Info up to 16 ASCII characters
Mode.Parameter.Sequence.X.Prep internal D0, external D1, external D2
Mode.Parameter.Sequence.X.Stirrer ON, OFF

Parameters for TIP: Parameters of the elements of TIP.

.Method: Method name of a method available in the user memory or on the card. Up to 8 ASCII characters.

.Pause: Pause time in s. INF means infinite. Continue the sequence with &m \$G.

.L4 Output: Warning: A pulse triggered by the limit value monitoring at L4 (pin 3) in a submethod sets an output set to active in TIP to inactive.

.L6 Output: Warning: An activate pulse at L6 output (pin 1) in a submethod sets an output set to active in TIP to inactive.

.Info: Entry of a message which is written into the display. The sequence remains in the display with the corresponding message. Continue with &m \$G.

.Prep: Preparation of titrating burette.

.Stirrer: Switching stirrer on/off.

3.2.2.76. Mode.Parameter.Presel.MeasMode pH, U, Ipol, Upol, OFF
Mode.Parameter.Presel.MeasInput 1, 2, diff.
Mode.Parameter.Presel.Ipol 0...1...±127
Mode.Parameter.Presel.Upol 0...400...±1270
Mode.Parameter.Presel.PolElectrTest ON, OFF
Mode.Parameter.Presel.Temp -170...25.0...500.0

Parameters for TIP: Selection of the measured quantity for manual measurements in the inactive state, see 3.2.2.2. Selection of the measuring input

(MeasInput) applies to measured quantities pH and U. "diff." means differential amplifier, see page 215. With Ipol the requests for the polarization current in uA (Ipol) and .PolElectrTest apply. With Upol the request for the polarization voltage in mV (Upol) applies. Entry in steps of 10 mV. .PolElectrTest also applies. If the test for polarized electrodes is switched on (ON), it will be performed on the change from the inactive state to an active state. The temperature applies to pH measurements.

3.2.2.77. UserMeth.FreeMem read only
Memory space, available for user methods or silo lines. \$Q sends the number of free bytes, e.g.
"4928".

3.2.2.78. UserMeth.Recall \$G
UserMeth.Recall.Name up to 8 ASCII characters
UserMeth.Store \$G
UserMeth.Store.Name up to 8 ASCII characters
UserMeth.Delete \$G
UserMeth.Delete.Name up to 8 ASCII characters
UserMeth.DeIAI \$G

Management of the internal method memory: Load, store and delete methods. An action is performed if "\$G" is sent to the corresponding node just after entering the name.

Do not use blank characters before and after method name!

.DeIAI: Deletes all methods in the user memory.

3.2.2.79. UserMeth.List.1.Name read only
UserMeth.List.1.Mode read only
UserMeth.List.1.Quantity read only
UserMeth.List.1.DosUnit read only
UserMeth.List.1.Bytes read only
UserMeth.List.1.Checksum read only
for each method

List of the methods in the user method memory with the following characteristics:

.Name: Name of the method
.Mode: Mode
.Quantity: Measured quantity
.DosUnit: Buret of the method
.Bytes: Number of bytes of the user memory used by the method
.Checksum: Checksum of the method, see 3.2.2.133.

3.2.2.80. MemoryCard.Recall \$G
MemoryCard.Recall.Name up to 8 ASCII characters
MemoryCard.Store \$G
MemoryCard.Store.Name up to 8 ASCII characters
MemoryCard.Delete \$G
MemoryCard.Delete.Name up to 8 ASCII characters

Administration of the method memory of the current directory of the memory card: load, store and delete methods. The action is carried out, if "\$G" is transmitted to the corresponding node.

Do not use blank characters before and after method name!

3.2.2.81. MemoryCard.ChangeDir \$G
MemoryCard.ChangeDir.Name up to 10 ASCII characters
MemoryCard.ChangeDir.Checksum \$G
MemoryCard.ChangeDir.Checksum.Value read only

Changing the current directory. The action is carried out if "\$G" is transmitted to the corresponding node.

Do not use blank characters before and after name!

The checksum characterizes the content of the directory.

3.2.2.82. MemoryCard.CreateDir \$G
MemoryCard.CreateDir.Name up to 10 ASCII characters

Delete directory. The action is carried out if "\$G" is transmitted to the corresponding node.

Do not use blank characters before and after name!

3.2.2.83. MemoryCard.DeIDir \$G
MemoryCard.DeIDir.Name up to 10 ASCII characters

Delete directory. The action is carried out if "\$G" is transmitted to the corresponding node.

Do not use blank characters before and after name!

3.2.2.84. MemoryCard.Backup \$G
MemoryCard.Backup.Name up to 10 ASCII characters

Backup of the internal memory onto the card. The action is carried out if "\$G" is transmitted to the corresponding node.

Do not use blank characters before and after name!

3.2.2.85. MemoryCard.Reload \$G
MemoryCard.Reload.Name up to 10 ASCII characters

Reload a backup from the memory card into the internal memory. The action is carried out if "\$G" is transmitted to the corresponding node.

Do not use blank characters before and after name!

3.2.2.86. MemoryCard.Format \$G
MemoryCard.CardLabel.Name up to 8 ASCII characters

Format the memory card. The action is carried out if "\$G" is transmitted to the corresponding node.

Do not use blank characters before and after name!

3.2.2.87. MemoryCard.FreeMem read only

Free memory on the card. \$Q sends number of free bytes, e.g. "4928".

3.2.2.88. MemoryCard.BatteryChange \$G
MemoryCard.BatteryChange.Date YYYY-MM-DD
 Date for changing battery. The date will be set with "\$G".

3.2.2.89. MemoryCard.List.Card.1.Name read only
MemoryCard.List.Card.1.Bytes read only
 for each directory

List of all directories on the memory card with the following information:

.Name: Name of the directory
 .Bytes: Number of bytes used by the directory

3.2.2.90. MemoryCard.List.ActDir.1.Name read only
MemoryCard.List.ActDir.1.Mode read only
MemoryCard.List.ActDir.1.Quantity read only
MemoryCard.List.ActDir.1.DosUnit read only
MemoryCard.List.ActDir.1.Bytes read only
MemoryCard.List.ActDir.1.Checksum read only
 for each method

List of all methods of the current card directory with the following information:

.Name: Name of the method
 .Mode: Mode
 .Quantity: Measured quantity
 .DosUnit: Buret of the method
 .Bytes: Number of bytes used by the method
 .Checksum: Checksum of the method, see 3.2.2.133.

3.2.2.91. Config.Monitoring.Validation.Status ON, OFF
Config.Monitoring.Validation.Interval 1...365...9999
Config.Monitoring.Validation.Counter 0...9999
Config.Monitoring.Validation.ClearCount \$G

Monitoring of validation.

.Interval: Time interval in days for validation.
 .Counter: Time counter in days since last validation.
 .ClearCount: Clears the above counter.

3.2.2.92. Config.Monitoring.Calibration.Status ON, OFF
Config.Monitoring.Calibration.MeasInput 1, 2, diff
Config.Monitoring.Calibration.Interval 1...7...9999
Config.Monitoring.Calibration.Counter 0...9999

Monitoring of pH calibration.

.MeasInput: Measuring input.
 .Interval: Time interval in days for calibration from the last calibration date.
 .Counter: Time counter in days since last calibration. The counter is reset to zero if a new calibration is carried out or if the calibration for the corresponding measuring input is entered manually.

- 3.2.2.93.** **Config.Monitoring.Service.Status** ON, OFF
Config.Monitoring.Service.Date XXXX-XX-XX
 Monitoring of service interval.
- 3.2.2.94.** **Config.Monitoring.DiagRep** ON, OFF
 Printing of system test report after each switching on of the Titrino.
- 3.2.2.95.** **Config.PeriphUnit.CharSet1** Epson, Seiko, Citizen, HP, IBM
Config.PeriphUnit.CharSet2
 Selection of the character set and the graphics control characters for COM1 resp. COM2 of the Titrino.
 IBM means the IBM character set following character set table 437 and IBM graphics control characters. Select 'IBM' for work with the computer.
- 3.2.2.96.** **Config.PeriphUnit.RepToComport** 1, 2, 1&2
 Selection of COM of the Titrino where manually triggered reports should be outputted.
- 3.2.2.97.** **Config.PeriphUnit.Balance** Sartorius, Mettler, Mettler AT, AND, Precisa
 Selection of the balance type.
- 3.2.2.98.** **Config.PeriphUnit.Stirrer** ON, OFF
 Automatic stirrer control. With "ON" the stirrer will be switched on after starting a method. At the end of the method it is switched off again. Valid for DET, MET, SET, KFT, STAT, DOC.
- 3.2.2.99.** **Config.PeriphUnit.RemoteBox.Status** ON, OFF
Config.PeriphUnit.RemoteBox.Keyboard US, deutsch, francais, español, schweiz.
Config.PeriphUnit.RemoteBox.Barcode input, method, id1, id2, id3, smpl size
 Connections via Remote Box.
 .Status: Select if a Remote Box is connected.
 .Keyboard: Type of keyboard which is connected to the Remote Box.
 .Barcode: Select target in Titrino where you wish to have the string from the barcode reader. "input" means that the string comes into the field where the cursor is currently placed.
- 3.2.2.100.** **Config.Aux.Language** english, deutsch, francais, espanol, italiano, portugese, svenska
 Selection of the dialog language.
- 3.2.2.101.** **Config.Aux.Set** \$G
Config.Aux.Set.Date YYYY-MM-DD

- Config.Aux.Set.Time** HH:MM
 Date and time.
 Input format of the date: Year-month-day, two-digit, enter leading zeros.
 Input format for the time: Hours:minutes, two-digit, enter leading zeros.
 Date and time have to be set with &Config.Aux.Set \$G just after entry of the value.
- 3.2.2.102. Config.Aux.RunNo** 0...9999
 Current sample number.
 Set to 0 on power on and initialization. After 9999, counting starts again at 0.
- 3.2.2.103. Config.Aux.AutoStart** 1...9999, **OFF**
 Number of automatic, internal starts.
- 3.2.2.104. Config.Aux.StartDelay** 0...999 999
 Start delay time in s. During this time, the data of the preceding determination are retained.
- 3.2.2.105. Config.Aux.ResDisplay** **bold**, standard
 Character set for the result display at the end of the determination.
- 3.2.2.106. Config.Aux.DevName** up to 8 ASCII characters
 Name of the instrument for connections with several units. It is advisable to use only the letters A...Z (ASCII No. 65...90), a...z (ASCII No. 97...122) and the numbers 0...9 (ASCII No. 48...57) when the function Setup.AutoInfo (3.2.2.169) is used at the same time.
 If a name has been entered, it will be printed out in the result report (full, short).
- 3.2.2.107. Config.Aux.Prog** read only
 Output of the program version.
 The Titrimo sends "751.0020" on requests with \$Q.
- 3.2.2.108. Config.RSSet1** \$G
Config.RSSet1.Baud 300, 600, 1200, 2400, 4800, **9600**,
19200, 38400, 57600, 115200
Config.RSSet1.DataBit 7, **8**
Config.RSSet1.StopBit **1**, 2
Config.RSSet1.Parity even, odd, **none**
Config.RSSet1.Handsh **HWs**, SWchar, SWline, none
 \$G sets all RS settings. The changes are performed only if the instrument is inactive. After the setting of the interface parameters, wait at least 2 s to allow the components to equilibrate.
 Settings of the values for the data transmission via the RS interface: baud rate, data bit, stop bit, parity and type of handshake, see also page 180 ff.

Baud rates >9600 need a PC which is equipped accordingly (e.g. with 16550 component).

The setting of the values must be initiated with \$G immediately after entry of the values.

3.2.2.109. Config.ComVar.C30 0... ±999 999
with up to **.C39**, etc.

Values of the common variables from C30 up to C39. Insert the common variables directly or describe the determination results directly from the method, see 3.2.2.7

3.2.2.110. Config.DosPrep.PowerOnPrep ON, OFF
Warning after power on, that the preparation of the buret has to be executed.

3.2.2.111. Config.DosPrep.Report ON, OFF
Report output after preparation of burets.

3.2.2.112. Config.DosPrep.Select internal D0, external D1, external D2
Selection of buret.

3.2.2.113. Config.DosPrep.DX.WarnInterval 5...9999, OFF
Time interval in min for automatic warning that a preparation has to be executed.

3.2.2.114. Config.DosPrep.DX(.Dosimat).V 0...3.5...99999.99
Volume in ml, which will be expelled at the preparation. Start of sequence see 3.2.2.152.

3.2.2.115. Config.DosPrep.DX(.Dosimat).Repeat 1...2...9
Number of cycles for expelling the volume at the preparation. Start of sequence see 3.2.2.152.

3.2.2.116. Config.DosPrep.DX(.Dosimat).DosRate 0.01...150, max.
Config.DosPrep.DX(.Dosimat).FillRate 0.01...150, max.
Dosing and filling rate in ml/in for the preparation. Start of sequence see 3.2.2.152.

3.2.2.117. Config.DosPrep.DX.Select Dosimat, Dosino
Selection of buret type. Start of sequence see 3.2.2.152.

3.2.2.118. Config.DosPrep.DX.Dosino.Outlet tip, flask
Place where the liquid of the preparation is expelled. Start of sequence see 3.2.2.152.

tip: To the tip.
flask: Back to the flask.

3.2.2.119. **Config.DosPrep.DX.Dosino.DosTubing.Length** 1...**40**...999.9
Config.DosPrep.DX.Dosino.DosTubing.Diameter 0.1...**2**...9.9
Config.DosPrep.DX.Dosino.AspirTubing.Length 1...**25**...999.9
Config.DosPrep.DX.Dosino.AspirTubing.Diameter 0.1...**2**...9.9

.DosTubing: Dosing tube. Length in cm, diameter in mm.

.AspirTubing: Aspiration tube. Length in cm, diameter in mm.

Start of sequence see 3.2.2.152.

3.2.2.120. **Config.DosPrep.DX.Dosino.DosRate** 0.01...150, **max.**
Config.DosPrep.DX.Dosino.FillRate 0.01...150, **max.**

Dosing and filling rate in ml/min for preparation. Start of sequence see 3.2.2.152.

3.2.2.121. **SmplData.Status** ON, **OFF**

On/off switching of silo memory. When the silo memory is switched on, the sample data are fetched from the lowest valid silo line.

3.2.2.122. **SmplData.OFFSilo.Id1** up to 8 ASCII characters
SmplData.OFFSilo.Id2 up to 8 ASCII characters
SmplData.OFFSilo.Id3 up to 8 ASCII characters
SmplData.OFFSilo.ValSmpl 6-digits, sign and decimal point
SmplData.OFFSilo.UnitSmpl up to 5 ASCII characters

Current sample data.

The identifications Id1...Id3 can be used in formulas as sample-specific calculation constants C21...C23.

If "no unit" is desired for the unit of the sample size, the blank string must be entered.

3.2.2.123. **SmplData.ONSil.Counter.MaxLines** read only
SmplData.ONSil.Counter.FirstLine read only
SmplData.ONSil.Counter.LastLine read only

Information on silo memory.

.MaxLines: Maximum possible number of silo lines.

.FirstLine: Lowest valid silo line.

.LastLine: Last occupied silo line.

3.2.2.124. **SmplData.ONSil.EditLine.1.Method** up to 8 ASCII characters
SmplData.ONSil.EditLine.1.Id1 up to 8 ASCII characters
SmplData.ONSil.EditLine.1.Id2 up to 8 ASCII characters
SmplData.ONSil.EditLine.1.Id3 up to 8 ASCII characters
SmplData.ONSil.EditLine.1.ValSmpl 6-digits, sign and dec.point
SmplData.ONSil.EditLine.1.UnitSmpl up to 5 ASCII characters
SmplData.ONSil.EditLine.1.C24 read only
SmplData.ONSil.EditLine.1.C25 read only

SmplData.ONSilO.EditLine.1.Mark read only
 etc., up to **.255**
 Contents of a silo line.
 .Method: Method used to process the sample, from the method memory or from the card.
 .Id: The identifications Id1...Id3 can also be used as sample-specific calculation constants C21...C23 in formulas.
 .UnitSmpl: If "no unit" is desired for the sample size, the blank string must be entered.
 .C24, .C25: Results which have been assigned to C24 and C25.
 .Mark: Mark of the silo line: "*" = deleted line, "+" = line which is worked off, "-" = line which is worked off and not valid for silo calculations (deleted), "/" = last worked-off line, where recalculation can still be done. Silo lines which have been worked off are "read only".

3.2.2.125. SmplData.ONSilO.DeLine \$G
SmplData.ONSilO.DeLine.LineNum 1...255, **OFF**
 Deletion of a silo line. The line # is deleted with &SmplData.ONSilO.DeLine \$G. If a formerly deleted line is edited again, it becomes valid (function "undelete").

3.2.2.126. SmplData.ONSilO.DeAll \$G
 Deletes the entire silo memory. Must be triggered with \$G.

3.2.2.127. SmplData.ONSilO.CycleLines **ON, OFF**
 Silo data cycling.
 With "ON", executed lines are copied to the next free silo lines, see page 96.
 Exercise caution if you edit the silo memory during the determinations!

3.2.2.128. SmplData.ONSilO.SaveLines **ON, OFF**
 Silo lines are not deleted when they are worked off. Assigned results are stored as C24 and C25. "Save lines" can only be set to "ON" if the silo is completely empty. Delete the silo, see 3.2.2.126.

3.2.2.129. HotKey.User.Name up to 10 ASCII characters
HotKey.User.Delete \$G
HotKey.User.Delete.Name up to 10 ASCII characters
HotKey.User.DeAll \$G
HotKey.User.List.1.Name read only
 Management of user names.
 .Name: Input of user names.
 .Delete.Name: Deletes selected user name with &HotKey.User.Delete \$G.
 .List: List of all user names.

3.2.2.130. Info.Report \$G
Info.Report.Select configuration, parameters, smpl data,

statistics, silo, calib, C-fmla, def, user method, **full**, short, mplist, curve, deriv, comb, meas crv, temp crv, adj para, scalc full, scalc srt, prep, calc, act dir, mem card, all, ff

\$G sends the selected report to the COM which is set in &Config.PeriphUnit.RepToComport:

configuration: Configuration report. Is not accessible during a running determination.

parameters: Parameter report of the current method. During a running determination only "live"-parameters are accessible.

smpl data: Current sample data.

statistics: Statistics table with the individual results.

silo: Contents of the silo memory.

calib: Calibration data of the measuring input in the current method.

C-fmla: Contents of the <C-fmla> key.

def: Contents of the <def> key.

user method: Contents of the method memory.

full: Full result report of the last completed determination.

short: Short result report of the last completed determination.

mplist: Measuring point list of the running determination.

curve: Titration curve of the last determination.

derive: 1st derive of titration curve of the last determination (with DET).

comb: 1st derive combined with the titration curve of the last determination (with DET).

meas crv: Curve measured value vs. time (with STAT, DOS, DOC and activated meas monitoring).

temp crv: Curve temperature vs. time (with STAT, DOS, DOC and activated temperature monitoring).

adj para: Adjustment parameters.

scalc full: Full report of the silo calculations.

scalc srt: Short report of the silo calculations.

prep: Preparation report.

calc: Calculation report of the current method.

act dir: Methods of the current directory of the memory card.

mem card: All methods of the memory card.

all: All reports.

ff: Form feed on printer.

Reports which are sent from the Titrino are marked with space (ASCII 32) and ' at the beginning. Then an individual identifier for each report follows. Reports which are triggered by RS232 (\$G) have the same introducer but without preceding space, i.e. they start with '.

3.2.2.131. Info.CalibrationData \$G

Info.CalibrationData.Inp1.pHas	-20.00... 7.00 ... +20.00	
Info.CalibrationData.Inp1.Slope	-9.999... 1.000 ... +9.999	
Info.CalibrationData.Inp1.Temp	-170.0... 25.0 ... +500.0	
Info.CalibrationData.Inp1.Date		read only
Info.CalibrationData.Inp1.ElectrodeId		read only

identical for .Inp2 and .Diff

pH calibration data for measuring input 1. After the calibration, the data are entered automatically together with the date of the calibration and the electrode identification, see 3.2.2.71.

Calibration data can be entered. They are accepted with &Info.CalibrationData \$G. If calibration data are entered, the calibration date is deleted.

3.2.2.132. Info.PrepData.D0.Date read only
 Info.PrepData.D0.Time read only
 Info.PrepData.D1.Type read only
 Info.PrepData.D1.Date read only
 Info.PrepData.D1.Time read only
 identical to .D2

Data from the last correctly executed preparation.

3.2.2.133. Info.Checksums \$G
 Info.Checksums.MPList read only
 Info.Checksums.ActualMethod read only

The checksums can be used to identify the content of a file unequivocally, e.g. files with identical content

have identical results of the checksums. An empty file has checksum "0". The calculation of the checksums is triggered with \$G.

.MPList: Result of the checksum of the current measuring point list.

.ActualMethod: Result of the checksum of the current method in the working memory. Identical methods with different method names have the same results of the checksum.

3.2.2.134. Info.DetermData \$G
 Info.DetermData.Write ON, OFF
 Info.DetermData.MPList.1.Attribute read only/read + write
 Info.DetermData.ExV read only/read + write
 Info.DetermData.MPList.1.X read only/read + write
 Info.DetermData.MPList.1.Y read only/read + write
 Info.DetermData.MPList.1.Z1 read only/read + write
 Info.DetermData.MPList.1.Z2 read only/read + write
 for every measuring point

Determination data in hexadecimal format. A measuring point list is available in mode DET, MET, STAT, DOS, DOC, SET, KFT, and MEAS.

Recalculation of the measuring data is triggered with \$G.

.Write: With "ON", the following nodes can be overwritten:
 &Info.DetermData.MP.List, &Info.TitrResults.Var.C4X (X = 0...5),
 &Info.TitrResults.TempVar.C7X (X = 0...9), and &Mode.Name.

.ExV: Volume of the exchange unit, with which the determination was executed

.MPList.1.Attribute: Attribute

.MPList.X: X coordinate, time

.MPList.Y: Y coordinate, volume

.MPList.Z1: Z1 coordinate, measuring value

.MPList.Z2: Z2 coordinate, temperature

3.2.2.135. Info.TitrResults.RS.1.Value read only
 etc., up to .9
 Info.TitrResults.EP.1.V read only

Info.TitrResults.EP.1.Meas read only
 etc., up to **.2**
Info.TitrResults.Var.C40 read only/read+ write
 etc., up to **.C47**
Info.TitrResults.Var.DTime read only/read+ write
.RS: Values of the calculated results.
.EP: Endpoints with DET, MET, SET, KFT:
 Volume coordinate in mL, e.g. "1.2340"
 Measured value coordinate in pH "5.12", mV (with U and Ipol) "-
 241" or uA (with Upol) "43.7".
.Var: Various variables. You may overwrite the variables C40...C45, see
 3.2.2.134.
 C40: Initial measured value in pH "5.12", mV (with U and Ipol)
 "41", uA (with Upol) "43.7" or °C (with T) "25.0". In MEAS final
 measured value.
 C41: End volume with SET, STAT, DOS, DOC in ml, "12.5360".
 C42: Time from start of titration to end in s, "62".
 C43: Volume drift on start of a SET/KFT titration from the condi-
 tioning in ul/min, "3.5".
 C44: (Last measured) temperature in °C. Used for the temperature
 compensation in pH measurements.
 C45: Start volume with SET, STAT, DOC in ml, "2.800".
 C46: Asymmetry pH of CAL, "6.89".
 C47: Relative electrode slope of CAL, "0.9950".
 C48: Voltage at maximum voltage (not valid in CAL and TIP)
 C49: Voltage at minimum voltage (not valid in CAL and TIP)
 DTime: Dosing time in s in DOS/DOC or time for the drift correc-
 tion in SET/KFT with conditioning, resp.

3.2.2.136. Info.TitrResults.FixEP.51.Value read only
 etc. up to **.59**
Info.TitrResults.pK.61.Value read only
 etc. up to **.69**
Info.TitrResults.TempVar.C70 read only/read+ write
 etc. up to **.C79**
.FixEP: Fix EP with DET, MET resp. Fix V with STAT. C5X corresponds to the
 fix volumes X, X = 1...9.
.pK: With DET, MET resp. time with given part of the end volume in s in STAT.
 C6X corresponds to X = 1...9.
.TempVar: Temporary variables in TIP corresponding to the assignments in
 the submethods.

3.2.2.137. Info.TitrResults.TimeWin.81.Mean read only
Info.TitrResults.TimeWin.81.Dev read only
 etc. up to **.89**
Info.TitrResults.MeanRateC80.Mean read only
Info.TitrResults.MeanRateC80.Dev read only

For STAT: Mean rates and standard deviations in mL/min which have been
 calculated with linear regression. C80 is the rate calculated over all points of
 the measuring point list, C8X are the rates in the time windows X, X = 1...9.

3.2.2.138.	Info.StatisticsVal.ActN	read only
	Info.Statistics.1.Mean	read only
	Info.Statistics.1.Std	read only
	Info.Statistics.1.RelStd	read only
	etc. up to .9	

The current values of the statistics calculation.

\$Q sends, e.g.

ActN: Current value of the individual results "3"

Data for MN1:

Mean: Mean value (decimal places as in result) "3.421"

Std: Standard deviation (1 decimal place more than in result) "0.0231"

RelStd: Relative standard deviation (in %, 2 decimal places) "0.14"

3.2.2.139.	Info.SiloCalc.C24.Name	read only
	Info.SiloCalc.C24.Value	read only
	Info.SiloCalc.C24.Unit	read only
	for .C25 as for .C24	
	Info.SiloCalc.C26.ActN	read only
	Info.SiloCalc.C26.Mean	read only
	Info.SiloCalc.C26.Std	read only
	Info.SiloCalc.C26.RelStd	read only
	for .C27 as for .C26	

The current values from the silo calculations. C26 is the mean value out of the C24 variables; C27 comes from C25.

\$Q sends:

C24.Name: Name of the assigned value "RS1"

C24.Value: Value "2.222"

C24.Unit: Unit of the assigned value "%"

C26.ActN: Number of single results "3"

C26.Mean: Mean (decimal places as for the result itself) "3.421"

C26.Std: Standard deviation (decimal places as for the result + 1) "0.0231"

C26.RelStd: Relative standard deviation (in %, 2 decimal places) "0.14"

3.2.2.140.	Info.ActualInfo.Inputs.Status	read only
	Info.ActualInfo.Inputs.Change	read only
	Info.ActualInfo.Inputs.Clear	\$G
	Info.ActualInfo.Outputs.Status	read only
	Info.ActualInfo.Outputs.Change	read only
	Info.ActualInfo.Outputs.Clear	\$G

Status sends the current status of the I/O lines, Change sends the information regarding whether a change in status of a line has taken place since the last clearing, Clear clears the change information. For the output, there is a conversion from binary to decimal, e.g.

	0 0 0 0 0 0 0 0 0 0 0 1 0 1 0
Line No.	13 12 11 10 9 8 7 6 5 4 3 2 1 0
Output:	$2^1 + 2^3 = "10"$

1 means ON or change; 0 means OFF or no change.

The lines are assigned as follows (see also pages 221ff):

Inputs:

0 Start (pin 21)

Outputs:

0 Ready (pin 5)

1	Stop (pin 9)	1	Cond. ok (pin 18)
2	Enter (pin 22)	2	Titration (pin 4)
3	Clear (pin 10)	3	EOD (pin 17)
4	Smpl Ready (pin 23)	4	Monitoring, line L4 (pin 3)
5	pin 11	5	Error (pin 16)
6	pin 24	6	Activate, line L6 (pin 1)
7	pin 12	7	Pulse for recorder (pin 2)
		8	not used (pin 6)
		9	not used (pin 7)
		10	Monitoring (pin 8)
		11	Monitoring (pin 13)
		12	Monitoring (pin 19)
		13	Monitoring (pin 20)

3.2.2.141. Info.ActualInfo.Assembly.CyclNo read only
 \$Q sends the current cycle number of the voltage measurement cycle, e.g. "127". From the cycle number and the cycle time (see 3.2.2.149), a time frame can be set up.

The cycle number is set to 0 on switching on the instrument, on every start and for QuickMeas. It is incremented as long as the instrument remains switched on.

3.2.2.142. Info.ActualInfo.Assembly.Counter.V read only
Info.ActualInfo.Assembly.Counter.Clear \$G
 \$Q sends the volume. With the function &Info.Assembly.Counter.Clear \$G, the volume counter is set to zero.

3.2.2.143. Info.ActualInfo.Assembly.Meas read only
 \$Q sends the current measured value from the assembly.

3.2.2.144. Info.ActualInfo.Titrator.CyclNo read only
Info.ActualInfo.Titrator.V read only
Info.ActualInfo.Titrator.Meas read only
Info.ActualInfo.Titrator.dVdt read only
Info.ActualInfo.Titrator.dMeasdt read only
Info.ActualInfo.Titrator.dMeasdV read only
Info.ActualInfo.Titrator.ERC read only
Info.ActualInfo.Titrator.T read only

\$Q sends the current values in the following formats:

	DET	MET	SET	KFT	STAT, DOC	DOS	MEAS	CAL
CyclNo	127	127	127	127	127	127	127	127
V(ml)	1.2345	1.2345	1.2345	1.2345	1.2345	1.2345	-	-
Meas:								
pH	3.345	3.345	3.6(mV)	-	3.6(mV)	3.6(mV)	3.345	3.345
U, Ipol (mV)	-345.6	-345.6	-345.6	-345.6	-345.6	-345.6	-345.6	-
Upol (uA)	-12.5	-12.5	-12.5	-12.5	-12.5	-12.5	-12.5	-
T (°C)	-	-	-	-	-	-	25.0	-

dVdt (ul/s)	-	-	2.5142	2.5142	2.5142	-	-	-
dMeasdt								
pH,U,IpH mV/s	0.7957	0.7957	0.7957	0.7957	0.7957	-	0.7957	0.7957
UpH (uA/s)	0.7957	0.7957	0.7957	0.7957	0.7957	-	0.7957	-
T (°C/s)	-	-	-	-	-	-	0.7957	-
dMeasdV (mV/ul)	-	-	10.6326	10.6326	10.6326	-	-	-
ERC	34	-	-	-	-	-	-	-
T (°C)	25.9	25.9	-	-	25.9	25.9	25.9	-

NV: Not Valid. If in the signal drift is OFF in modes MEAS and CAL, the signal drift is NV.

OV will be sent for "overrange".

A time frame can be set up from the cycle number and the cycle time (see 3.2.2.149). The cycle number is set to 0 at the start of a method and it is incremented until the end of the method.

3.2.2.145. Info.ActualInfo.MeasPt.Index read only
 Info.ActualInfo.MeasPt.X read only
 Info.ActualInfo.MeasPt.Y read only
 Info.ActualInfo.MeasPt.Z1 read only
 Info.ActualInfo.MeasPt.Z2 read only
 Info.ActualInfo.EP.Index read only
 Info.ActualInfo.EP.X read only
 Info.ActualInfo.EP.Y read only

\$Q sends the last entry into the measuring point list (.MeasPt) or the last entry into the list of EP's with DET, MET, SET, KFT.

.MeasPt.X"165" Volume (DET, MET), time (STAT, DOS, DOC, SET, KFT, MEAS), resp. of the MPList in s

.MeasPt.Y"3.654" Measured value (DET, MET), volume (STAT, DOS, DOC, SET, KFT), resp. of the MPList in mL

.MeasPt.Z1"6.34" Measured value (STAT, DOS, DOC, SET, KFT, MEAS) of the MPList, format depends on the measured quantity

.MeasPt.Z2"25.8" Temperature of the MPList in °C

.EP.X"1.234" Volume coordinate of the EP

.EP.Y"5.34" Measured value coordinate of the EP

3.2.2.146. Info.ActualInfo.Oven.HeatTime read only
 Info.ActualInfo.Oven.SampleTemp read only
 Info.ActualInfo.Oven.LowTemp read only
 Info.ActualInfo.Oven.HighTemp read only
 Info.ActualInfo.Oven.GasFlow read only
 Info.ActualInfo.Oven.UnitFlow read only

\$Q sends the current values from a connected KF Oven. If no Oven is connected, the values are empty.

.HeatTime: Heating time of sample in s.

.SampleTemp: Nominal sample temperature in °C.

.LowTemp: Lowest temperature during the sample heating time in °C.

.HighTemp: Highest temperature during the sample heating time in °C.

.GasFlow: Average gas flow during sample heating time.

.UnitFlow: Unit of gas flow.

3.2.2.147. **Info.ActualInfo.Display.L1** up to 32 ASCII characters
Info.ActualInfo.Display.L8 up to 32 ASCII characters
Info.ActualInfo.Display.DelAll \$G

Lines of the display. The display can be written to from the computer. Proceed as follows:

1. Lock the display, see 3.2.2.163.
2. Delete the whole display (.DelAll).
3. For writing onto the display, the standard character set will be used.
4. Unlock the display, see 3.2.2.163
5. Delete the whole display (.DelAll).
6. Send a value to nod &Config.Aux.ResDisplay (see 3.2.2.105) to refresh the display.

\$Q sends the contents of the corresponding display line.

3.2.2.148. **Info.ActualInfo.Comport.Number** read only
 \$Q sends the comport number of the Titrimo where the PC is connected.

3.2.2.149. **Info.Assembly.CycleTime** read only
Info.Assembly.ExV read only

Inquiries regarding basic variables of the assembly: Cycle time in s, volume of the active Exchange Unit in mL.

3.2.2.150. **Info.Assembly.Bur.Select** **internal D0**, external D1, external D2
 Selection of buret for assembly functions.

3.2.2.151. **Info.Assembly.Bur.Empty** \$G, \$S, \$H, \$C
 Starts the function "empty". Only possible with the burette type "Dosino". The parameters are under the function &Config.DosPrep, see 3.2.2.118 up to 3.2.2.120.

3.2.2.152. **Info.Assembly.Bur.Prep** \$G, \$S, \$H, \$C
 Starts the function "preparation". The parameters are under the function &Config.DosPrep, see 3.2.2.110 up to 3.2.2.120.

3.2.2.153. **Assembly.Bur.Rates.Forward.Selected** digital, **analog**
Assembly.Bur.Rates.Forward.Digital 0...150, **max.**
Assembly.Bur.Rates.Reverse.Selected digital, **analog**
Assembly.Bur.Rates.Reverse.Digital 0...150, **max.**

Expel and aspirating rate.

Digital or analog control. With digital control, the inputted value applies (in mL/min). "max." means maximum possible rate with the Exchange Unit in current use.

Analog means rate control with the analog potentiometer on Titrimo.

3.2.2.154. Assembly.Bur.Fill \$G, \$H, \$C
 \$G starts the 'FILL' mode of the burette function.

3.2.2.155. Assembly.Bur.ModeDis \$G, \$\$, \$H, \$C
Assembly.Bur.ModeDis.Selected volume, time
Assembly.Bur.ModeDis.V 0.0001...0.1...9999
Assembly.Bur.ModeDis.Time 0.25...1...86400
Assembly.Bur.ModeDis.VStop 0.0001...9999, OFF
Assembly.Bur.ModeDis.AutoFill ON, OFF

Dispensing mode with parameters. The dispensing mode can only be started and stopped via the RS Control. During a running dosification, no method can be started at the Titrino.

.Selected: Dispensing of volume increments or during a preset time.

.Volume, .Time: Size of the volume increments or entry of time.

.VStop: Limit volume for the dispensing.

.AutoFill: ON means automatic filling after every dispensing.

3.2.2.156. Assembly.Meas.Status ON, OFF
Assembly.Meas.MeasInput 1, 2, Diff., Ipol, Upol, Temp
Assembly.Meas.Ipol $\pm 127...1...+127$
Assembly.Meas.Upol $\pm 1270...400...+1270$

Measurement in assembly. The measuring function can only be started via RS Control. When the measuring function is switched on, no method can be started at the Titrino.

.Input: Selection of the potentiometric measuring input 1, 2, diff., polarized electrodes or temperature.

.Ipol: Polarization current in uA.

.Upol: Polarization potential in mV, entry in steps of 10 mV.

3.2.2.157. Assembly.Outputs.AutoEOD ON, OFF
Assembly.Outputs.SetLines \$G
Assembly.Outputs.SetLines.L0 active, inactive, pulse, OFF
 up to .L13
Assembly.Outputs.ResetLines \$G

Setting the I/O output lines.

.AutoEOD: The automatic output of the EOD (End of Determination) at the end of the determination can be switched off. Thus, for example, in conjunction with a Titrino several determinations can be performed in the same beaker. Before AutoEOD is switched on, line 3 must be set to "OFF".

.SetLines: With \$G, all lines are set.

.SetLines.LX: Set the line LX. "active" means setting of a static signal, "inactive" means resetting of the signal, "pulse" means output of a pulse of app. 150 ms, "OFF" means the line is not operated, see also page 223.

Warnings:

- If you have "AutoEOD" to "ON", an active line 3 is set to "inactive" by the EOD pulse.

- L6 is the line of the activate pulse. An active line 6 is set to "inactive" by the activate pulse.
- L5 is the error line. It is continuously controlled by the Titrimo program and can therefore not be set freely.

Line assignments in Titrimo program:

L0	Ready, inactive state
L1	Conditioning OK
L2	Titration in progress
L3	EOD (End Of Determination)
L4	Can be set in TIP
L5	Error
L6	Activate pulse + can be set in TIP
L7	Buret volume pulses
L8, 9	---
L10-13	Monitoring in modes STAT, DOS, DOC

.ResetLines: Lines are set to the inactive status (= high).

3.2.2.158. Assembly.Stirrer **ON, OFF**
Switching stirrer ON/OFF.

3.2.2.159. Setup.Comport **1, 2, 1&2**
Selects the Titrimo COM for the output of automatic info:
&Setup.Keycode
&Setup.Trace
&Setup.SendMeas
&Setup.AutoInfo

3.2.2.160. Setup.Keycode **ON, OFF**
ON means the key code of a key pressed on the Titrimo is outputted. The key code comprises 2 ASCII characters; table of the keys with their code, see page 195. A keystroke of key 11 is sent as follows:
#11
The beginning of the message is marked by a space (ASCII 32).

3.2.2.161. Setup.Tree.Short **ON, OFF**
Setup.Tree.ChangedOnly **ON, OFF**
Definition of the type of answer to \$Q.
.Short: With "ON", each path is sent with only the necessary amount of characters in order to be unequivocal (printed in bold in this manual). A combination of .Short and .ChangedOnly is not possible.
.ChangedOnly: Sends only the changed values, i.e. values which have been edited. All paths are sent absolute, i.e. from the root.

3.2.2.162. Setup.Trace **ON, OFF**
The Titrimo automatically reports when a value has been confirmed with <enter> at the Titrimo. Message, e.g.:
&SmpData.OFFSilo.Id1"Trace"

The beginning of the message is marked by a space (ASCII 32).

3.2.2.163. Setup.Lock.Keyboard	ON, OFF
Setup.Lock.Config	ON, OFF
Setup.Lock.Parameter	ON, OFF
Setup.Lock.SmplData	ON, OFF
Setup.Lock.UserMeth.Recall	ON, OFF
Setup.Lock.UserMeth.Store	ON, OFF
Setup.Lock.UserMeth.Delete	ON, OFF
Setup.Lock.Display	ON, OFF

ON means disable the corresponding function:

- .Keyboard: Disable all keys of the Titrimos
- .Config: Disable the <configuration> key
- .Parameter: Disable the <parameter> key
- .SmplData: Disable the <smpl data> key
- .UserMeth.Recall: Disable "recall" in <user meth> key
- .UserMeth.Store: Disable "store" in <user meth> key
- .UserMeth.Delete: Disable "delete" in <user meth> key
- .Display: Disable the display, i.e. it will not be written to by the device program of the Titrimo and can be operated from the computer.

3.2.2.164. Setup.Mode.StartWait	ON, OFF
Setup.Mode.FinWait	ON, OFF

Holding points in the method sequence. If they are "ON", the sequence stops until "OFF" is sent. Switching the instrument on sets both nodes to OFF:

- .StartWait: Holding point right after starting a method or submethod in TIP (holding point after AutoInfo !" .T.GC").
- .FinWait: Holding point at the end a method or submethod in TIP (holding point after AutoInfo !" .T.F").

3.2.2.165. Setup.SendMeas.SendStatus	ON, OFF
Setup.SendMeas.Interval	0.08...4...16200, MPList

- .SendStatus: ON means the automatic transmission of measured values (see 3.2.2.167 and 3.2.2.168) in the inputted interval is active.
- .Interval: Time interval (in s) for the automatic transmission of associated measured values defined under points 3.2.2.167 and 3.2.2.168. The inputted value is rounded off to a multiple of 0.08. The smallest possible time interval depends on the number of measured values which have to be sent, on the baud rate, on the load on the interface and on the type of device connection. With "MPList" the measured values are sent at the time of their entry into the measured point list.

The automatic transmission is switched on/off with 'SendStatus'.

3.2.2.166. Setup.SendMeas.Select	Assembly, Titrator
---	---------------------------

Selection of the unit of which the measured values should be sent (3.2.2.167 or 3.2.2.168).

3.2.2.167. Setup.SendMeas.Assembly.CyclNo	ON, OFF
Setup.SendMeas.Assembly.V	ON, OFF
Setup.SendMeas.Assembly.Meas	ON, OFF

Selection of the values from Assembly for the output in the set time interval (see 3.2.2.165):

- .CyclNo:** Cycle number of the potential measurement. Together with the cycle time (3.2.2.149), a time frame can be set up. The cycle number is set to 0 on switching on the instrument and it is always incremented as long as the instrument remains switched on.
- .V:** Volume
- .Meas:** Measured value associated to the cycle number. The unit "assembly" must be preset (see 3.2.2.166).

3.2.2.168. Setup.SendMeas.Titrator.CyclNo	ON, OFF
Setup.SendMeas.Titrator.V	ON, OFF
Setup.SendMeas.Titrator.Meas	ON, OFF
Setup.SendMeas.Titrator.dVdt	ON, OFF
Setup.SendMeas.Titrator.dMeasdt	ON, OFF
Setup.SendMeas.Titrator.dMeasdV	ON, OFF
Setup.SendMeas.Titrator.ERC	ON, OFF
Setup.SendMeas.Titrator.T	ON, OFF

Selection of the values from the titrator which are sent in the set time interval (see 3.2.2.165, formats see 3.2.2.144):

- .CyclNo:** Cycle number. Together with the cycle time (3.2.2.149), a time frame can be set up. The other data belong to the corresponding cycle number. The cycle number is set to 0 at the start of a method and it is incremented until the end of the method.
- .V:** Volume.
- .Meas:** Measuring value with STAT, DOS, DOC with activated measured value monitoring
- .dVdt:** associated volume drift.
- .dMeasdt:** associated measured value drift.
- .dMeasdV:** associated 1st derivative of the titration curve.
- .ERC:** ERC in DET.
- .T:** In STAT, DOS, DOC, with activated temperature monitoring. The unit "titrator" must be preset (see 3.2.2.166).

3.2.2.169. Setup.AutoInfo.Status	ON, OFF
Setup.AutoInfo.P	ON, OFF
Setup.AutoInfo.T.R	ON, OFF
Setup.AutoInfo.T.G	ON, OFF
Setup.AutoInfo.T.GC	ON, OFF
Setup.AutoInfo.T.S	ON, OFF
Setup.AutoInfo.T.B	ON, OFF
Setup.AutoInfo.T.F	ON, OFF
Setup.AutoInfo.T.E	ON, OFF
Setup.AutoInfo.T.H	ON, OFF
Setup.AutoInfo.T.C	ON, OFF
Setup.AutoInfo.T.O	ON, OFF
Setup.AutoInfo.T.N	ON, OFF

Setup.AutoInfo.T.Re	ON, OFF
Setup.AutoInfo.T.Si	ON, OFF
Setup.AutoInfo.T.M	ON, OFF
Setup.AutoInfo.T.EP	ON, OFF
Setup.AutoInfo.T.RC	ON, OFF
Setup.AutoInfo.C.B1	ON, OFF
Setup.AutoInfo.C.R1	ON, OFF
Setup.AutoInfo.C.B2	ON, OFF
Setup.AutoInfo.C.R2	ON, OFF
Setup.AutoInfo.I	ON, OFF
Setup.AutoInfo.O	ON, OFF

ON means that the Titrino reports automatically the moment the corresponding change occurs.

.Status: Global switch for all set AutoInfo.

.P PowerOn: Simulation of power on (3.2.2.171). Not from mains.

Messages from node .T, Titrator:

.T.R Ready: Status 'Ready' has been reached.

.T.G Go: Instrument has been started.

.T.GC GoCommand: Instrument (or submethod in TIP) has received a go command.

.T.S Stop: Status 'Stop' has been reached.

.T.B Begin of sequence (or submethod).

.T.F Final: End of determination (or submethod), the final steps will be carried out.

.T.E Error. Message together with error number, see page 111ff.

.T.H Hold: Status 'Hold' has been reached.

.T.C Continue: Continue after hold.

.T.O Conditioning OK: EP reached (in SET, KFT with conditioning).

.T.N Conditioning Not OK: EP not reached (in SET, KFT with conditioning).

.T.Re Request: In the inquiry of an identification or the sample size after start of titration.

.T.Si SiloEmpty: Silo empty, i.e. the last line has been removed from the silo memory.

.T.M MeasList: Entry in the measuring point list (with DET, MET, STAT, DOS, DOC, SET, KFT, MEAS).

.T.EP EPList: Entry into EP list (with DET, MET, SET, KFT)

.T.RC Results have been recalculated.

Messages from node .C, Comport:

.C.B1 COM1: A report is outputted on COM1. During this time, COM2 will be blocked. COM2 is generally blocked, if COM1 is busy.

.C.R1 COM1 is ready again. (Comes also when you <QUIT> an error.)

.C.B2, .R2 Identical for COM2.

Messages for changings in the I/O lines. If the changings are made simultaneously, there is 1 message. Pulses receive 2 messages: one message each for line active and inactive.

.I Input: Change of an input line.

.O Output: Change of an output line (except 7, pin 2, for recorder pulses).

If a change occurs that requires a message, the Titrino sends space (ASCII 32) and ! as an introducer. This is followed by the name of the device (see 3.2.2.106). Special ASCII characters in the device name are ignored. If no de-

vice name has been entered, only ! is sent. Finally the Titrino sends the information which node has triggered the message.

Example: !John".T.Si": The message was triggered from instrument "John", node .T.Si

3.2.2.170. Setup.Graphics.Grid **ON, OFF**
Setup.Graphics.Frame **ON, OFF**
Setup.Graphics.Scale **Full, Auto**
Setup.Graphics.Recorder.Right 0.2...**0.5**...1.00
Setup.Graphics.Recorder.Feed 0.01...**0.05**...1.00

Change in the appearance and the format of the curve for the output. The settings are valid for both Titrino COM ports.

.Grid: On/off switching of grid over curve.
.Frame: On/off switching of frame surrounding the curve. If grid and frame are switched off, the curve is printed faster as the print head does not have to move to the end of the paper.
.Scale: Type of scaling of the measured value axis: Full means that the scale runs from the smallest up to the greatest measured point. With auto, the smallest measured value is taken and the next smaller tick defines the beginning of the scale; the next greater tick to the greatest measured value is the end of the scale.
.Right: Relative specification of the width of the output medium (e.g. paper width) for the length of the measured value axis. 1 means the measured value axis is plotted over the entire width of the paper (largest possible width). In extreme cases, the writing of the right tick may lie outside.
.Feed: Length of the volume axis referred to the burette cylinder volume, V(B) per cm (0.1 means, e.g. 1 mL/cm with a 10 mL Exchange Unit). Depending on the printer, the measure in cm may not always be correct.

3.2.2.171. Setup.PowerOn **\$G**
Simulation of 'power on'. The device has the same status as after power on: The cylinder is filled, error messages deleted and the current sample number set to 0. The method last used is ready for operation.

3.2.2.172. Setup.Initialise **\$G**
Setup.Initialise.Select **ActMeth, Silo, Calib, Config, Assembly, Setup, All**

Setting of default values for the following areas:

ActMeth: Current method. Parameters, calculations, and assignments for the data output, operands C01...C19.
Silo: The silo memory is deleted. Same function as delete entire silo.
Calib: pH calibration data for all measuring inputs.
Config: All values under &Config.
Assembly: All values under &Assembly.
Setup: All values under &Setup.
All: Values of the entire tree (except silo and method memory).
The action must be triggered with &Setup.Initalise \$G.

3.2.2.173. Setup.RamInit \$G

Initializes instrument, see page 204. All parameters are set to their default value and error messages are cleared. The user and silo memories will be deleted. The user memory contains the default user methods from Metrohm.

3.2.2.174. Setup.InstrNo \$G
Setup.InstrNo.Value **serial number**, 8 ASCII characters

Instrument identification for report output.
Set the value with &Setup.InstrNo \$G .

3.2.2.175. Diagnose.Report \$G

Output of the report containing the adjustment parameters. The Titrino has to be in its inactive basic state.

3.3 Properties of the RS 232 Interface

Data Transfer Protocol

The Titrino is configured as DTE (Data Terminal Equipment).

The RS 232 interface has the following technical specifications:

- Data interface according to the RS 232C standard, adjustable transfer parameters, see page 10.
- Max. line length: 512 characters
- Control characters: C_R (ASCII DEC 13)
 L_F (ASCII DEC 10)
XON (ASCII DEC 17)
XOFF (ASCII DEC 19)
- Cable length: max. approx. 15 m

Start	7 or 8 Data Bit	Parity Bit	1 or 2 Stop Bit
-------	-----------------	------------	-----------------

Only a shielded data cable (for example, METROHM D.104.0201) may be used to couple the Titrino with foreign devices. The cable shield must be properly grounded on both instruments (pay attention to current loops; always ground in a star-head formation). Only plugs with sufficient shielding may be used (for example, METROHM K.210.0381 with K.210.9045).

3.3.1 Handshake


Software-Handshake, SWchar

Handshake inputs on the Titrino (CTS) are not checked.


Handshake outputs (DTR, RTS) are set by the Titrino.

The Titrino sends XOFF when its input buffer contains 384 characters. After this it can receive 128 extra characters (including L_F).

Titrino as Receiver :

Titrimo as Sender :


max. characters: 2 characters at 300...9600 baud
16 characters at ≥ 19200 baud


Software-Handshake, SWline

Handshake input ports on the Titrimo (CTS) are not checked.
Handshake output ports (DTR, RTS) are set by the Titrimo.
The Titrimo has an input buffer which can accept up to 512 characters.

Titrimo as Receiver :


Titrimo as Sender:


Titrimo transmission can be stopped by external instruments with XOFF. After XOFF is received the Titrimo completes sending the line already started. If data output is disabled for more than 6 s by XOFF, E43 appears in the display.

Hardware-Handshake, HWs

Titrimo as Receiver :


Titrimo as Sender:

The data flow can be interrupted by deactivating the CTS line.

3.3.2 Pin Assignment

RS232C Interface

Protective earthing

Direct connection from cable plug to the protective ground of the instrument.

Polarity allocation of the signals


- Data lines (TxD, RxD)
 - voltage negative (< -3 V): signal state "ON"
 - voltage positive (> +3 V): signal state "ZERO"
- control or message lines (CTS, RTS, DTR)
 - voltage negative (< -3 V): OFF state
 - voltage positive (> +3 V): ON state

In the transitional range from +3 V to -3 V the signal state is undefined.

Driver 14C88 according to EIA RS 232C specification

Receiver 14C89 " "

Contact arrangement at plug (female) for RS 232C socket (male)

View of soldered side of plug

Ordering numbers:
K.210.0381 and K.210.9045

No liability whatsoever will be accepted for damage or injury caused by improper interconnection of instruments.

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6.4.2 Certificate of Conformity and System Validation

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

Name of commodity:	751 GPD Titrimo
System software:	Stored in ROMs
Name of manufacturer:	Metrohm Ltd., Herisau, Switzerland

This Metrohm instrument has been built and has undergone final type testing according to the standards:

Electromagnetic compatibility: Emission

EN50081-1/92, EN55022/class B, EN55011/class B	Generic emission
---	------------------

Electromagnetic compatibility: Immunity

EN50082-2/95, EN50082-1/92	Immunity
IEC801-2/91 (level 4), NAMUR/93	Static discharge
IEC801-3, ENV50140/93+ENV50204/93 (level 3)	Radiated rf electromag.field immunity
IEC801-4, IEC1000-4-4/95 (level 4)	EI.fast transient requirements
IEC801-5, IEC1000-4-5/95 (level 2/3)	"Surges" immunity
IEC801-6, IEC1000-4-6/96, ENV50141/93 (level 3)	Immunity to conducted disturbances
NAMUR/93 Paragr. 3.2.2.	Voltage dips, short interruptions

Security specifications

IEC1010 class1, EN61010 class1, UL3101-1, EN60947:IP31

The technical specifications are documented in the instruction manual.

The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance.

The features of the system software are documented in the instruction manual.

Metrohm Ltd. is holder of the SQS certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.

Herisau, March 14. 1997

Dr. J. Frank
Development Manager

Ch. Buchmann
Production and
Quality Assurance Manager

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