


 **Metrohm**  
Ion analysis  
CH-9101 Herisau/Switzerland  
Phone ++41 71 353 85 85  
Fax ++41 71 353 89 01  
Internet www.metrohm.ch  
E-Mail info@metrohm.ch

# 763 PVC Thermomat

with «763 PVC Thermomat 1.0» PC program



8.763.1003 Instructions for Use

26.10.1999 / dö

# Table of contents

<b>1</b>	<b>Introduction</b> .....	<b>1</b>
1.1	Determination of thermal stability.....	1
1.2	Instrument description.....	2
1.3	Operating elements.....	3
1.4	Information about these Instructions for Use.....	6
1.4.1	Organization.....	6
1.4.2	Notation and pictograms.....	7
1.5	Safety notes.....	8
1.5.1	Electrical safety.....	8
1.5.2	General precautionary rules.....	8
<b>2</b>	<b>Installation</b> .....	<b>9</b>
2.1	Setting up the instrument.....	9
2.1.1	Packaging.....	9
2.1.2	Check.....	9
2.1.3	Location.....	9
2.2	Attaching the accessories.....	10
2.2.1	Accessories for external gas supply.....	10
2.2.2	Reaction and measuring vessels.....	10
2.2.3	Exhaust collection tube (option).....	11
2.3	Mains connection.....	13
2.3.1	Checking the mains voltage.....	13
2.3.2	Exchange of fuses.....	13
2.3.3	Mains cable and mains connection.....	14
2.3.4	Switching the instrument on/off.....	14
2.4	Connection to the PC.....	15
2.4.1	Connection PVC Thermomat – PC.....	15
2.4.2	Software installation.....	15
2.4.3	Basic settings.....	16
<b>3</b>	<b>Operating tutorial</b> .....	<b>19</b>
3.1	Determinations.....	19
3.1.1	Installation.....	19
3.1.2	Preparing for determinations.....	20
3.1.3	Carrying out determinations.....	22
3.2	Results.....	24
3.2.1	Determination overview.....	24
3.2.2	Determination and method data.....	26
3.2.3	Graphics.....	28

<b>4</b>	<b>Operation</b>	<b>31</b>
<b>4.1</b>	<b>Fundamentals of the operation</b>	<b>31</b>
4.1.1	Starting/closing the program	31
	Start the «763 PVC Thermomat» program	31
	Close the «763 PVC Thermomat» program	31
4.1.2	Glossary	32
4.1.3	Control window	33
	Open	33
	Arrangement	33
	Menus	33
	Icons	34
	Close	34
4.1.4	Results window	35
	Open	35
	Arrangement	35
	Menus	35
	Icons	36
	Close	36
4.1.5	Overview of file types	37
4.1.6	Context-sensitive menus	37
4.1.7	Mouse functions	37
4.1.8	Help	37
<b>4.2</b>	<b>Instrument and software settings</b>	<b>38</b>
4.2.1	Communication	38
4.2.2	User permissions	39
	Functions for groups	39
	Functions for users	40
4.2.3	Timer	41
4.2.4	Gas flow control	41
4.2.5	Recording temperature	42
4.2.6	Options	43
	Save settings	43
	Clear sample ID1 history	43
<b>4.3</b>	<b>Program information</b>	<b>44</b>
4.3.1	Unit information	44
4.3.2	Status overview	45
4.3.3	Events overview	46
	Event window	46
	Menus	47
	Icons	47
	Filter	47
	Parameter for event recording	47

<b>4.4</b>	<b>Calibration functions</b>	<b>48</b>
4.4.1	Determination of cell constants	48
	Manual entry of cell constants	48
	Automatic determination of cell constants	48
4.4.2	Determination of Delta T	50
	Preparations	50
	Carrying out	52
<b>4.5</b>	<b>Methods</b>	<b>54</b>
4.5.1	Method handling	54
	Create method	54
	Open method	54
	Method manager	54
	Save method	55
	Print method report	55
4.5.2	Method parameters	55
	Overview	55
	Parameters	56
	Evaluation	58
	Induction time	59
	Stability time	59
	Curves	60
	Formulas	61
	Standards	62
	Documentation	63
	Description	64
<b>4.6</b>	<b>Determinations</b>	<b>65</b>
4.6.1	Flow chart	65
4.6.2	Sample preparation	66
4.6.3	Prepare instruments and accessories	66
	Clean instruments and accessories	66
	Prepare measuring vessels	67
	Prepare reaction vessels	68
4.6.4	Carry out determinations	69
	Select method	69
	Switch on heating	69
	Enter sample identifiers	70
	Insert reaction vessels	70
	Start determination	70
	Live curve	71
	Live parameters	71
	Stop determination	72
<b>4.7</b>	<b>Results</b>	<b>73</b>
4.7.1	Determination overview	73
	Open determination overview	73
	Format determination overview	74
	Mark determinations	78
	Find determination	78
	Delete determinations	79
	Sort and filter determinations	80

<b>4.7.2</b>	<b>Determination and method data</b> .....	<b>84</b>
	Open view .....	84
	Determination .....	84
	Results .....	85
	Method .....	86
	Parameters .....	87
	Evaluation .....	88
	Formulas .....	89
	Standards .....	90
	Curves .....	90
	Documentation .....	91
	Warnings .....	92
<b>4.7.3</b>	<b>Graphics and reevaluation</b> .....	<b>93</b>
	Single graph .....	93
	Multiple graph .....	93
	Live graph .....	94
	Switch zoom on/off .....	95
	Switch autoscaling on/off .....	95
	Switch offset on/off .....	95
	Copy graph to clipboard .....	95
	Print graph .....	95
	Reevaluation .....	96
<b>4.7.4</b>	<b>Extrapolation</b> .....	<b>99</b>
	Extrapolation window .....	99
	Print extrapolation results .....	100
	Settings for extrapolation .....	101
	Parameters for extrapolation .....	101
<b>4.7.5</b>	<b>Recalculation</b> .....	<b>102</b>
	Recalculation window .....	102
	Method .....	102
	Formulas .....	103
	Standards .....	104
<b>4.7.6</b>	<b>Data export and printing</b> .....	<b>106</b>
	Copy to clipboard .....	106
	Export to other database .....	106
	Export to Microsoft Word .....	106
	Export to Microsoft Excel .....	106
	Export of measured values .....	107
	Export of determination and method data .....	107
	Printing of results in the control window .....	108
	Printing of data in the results window .....	109
<b>4.7.7</b>	<b>Options</b> .....	<b>110</b>
	Program settings .....	110
	Report settings .....	110
	Company settings .....	111
	Results .....	112
	Scaling .....	112
	Extrapolation .....	113
	Live graph .....	114
	Options for the administrator .....	114
	Save settings .....	115
<b>4.7.8</b>	<b>Windows</b> .....	<b>115</b>
	Tile windows vertically .....	115
	Tile windows horizontally .....	115
	Cascade windows .....	115
	Arrange icons .....	115
	Close all windows .....	115

<b>4.8</b>	<b>GLP</b> .....	<b>116</b>
4.8.1	General information about GLP and validation.....	116
4.8.2	GLP properties .....	117
4.8.3	GLP status .....	117
4.8.4	GLP tests.....	118
	GLP temperature test .....	118
	GLP conductivity test.....	123
	GLP gas flow test .....	126
4.8.5	GLP results.....	128
	Print GLP certificate.....	128
	View GLP results.....	128
	Print GLP results.....	130
<b>5</b>	<b>Maintenance – Faults</b> .....	<b>131</b>
<b>5.1</b>	<b>Maintenance and service</b> .....	<b>131</b>
5.1.1	Care.....	131
5.1.2	Maintenance by Metrohm service .....	131
<b>5.2</b>	<b>Faults and malfunctions</b> .....	<b>132</b>
5.2.1	Error messages in the PC program .....	132
5.2.2	Status and error messages on the instrument.....	132
<b>6</b>	<b>Appendix</b> .....	<b>135</b>
<b>6.1</b>	<b>Technical data</b> .....	<b>135</b>
6.1.1	General information .....	135
6.1.2	Temperature regulation and measurement.....	135
6.1.3	Conductivity measurement.....	136
6.1.4	Gas flow regulation .....	136
6.1.5	GLP test set .....	136
6.1.6	RS232 interface.....	136
6.1.7	Mains connection .....	137
6.1.8	Safety specifications.....	137
6.1.9	Electromagnetic compatibility (EMC) .....	137
6.1.10	Ambient temperature .....	138
6.1.11	Housing .....	138
<b>6.2</b>	<b>Standard equipment</b> .....	<b>139</b>
<b>6.3</b>	<b>Optional accessories</b> .....	<b>142</b>
<b>6.4</b>	<b>Warranty and conformity</b> .....	<b>143</b>
6.4.1	Warranty.....	143
6.4.2	Software license .....	143
6.4.3	EU Declaration of conformity .....	144
6.4.4	Certificate of conformity and system validation.....	145
<b>6.5</b>	<b>Index</b> .....	<b>146</b>

## List of figures

<u>Fig. 1:</u>	Schematic measuring arrangement .....	1
<u>Fig. 2:</u>	Front of the 763 PVC Thermomat .....	3
<u>Fig. 3:</u>	Rear of the 763 PVC Thermomat .....	4
<u>Fig. 4:</u>	Installation of reaction and measuring vessels .....	12
<u>Fig. 5:</u>	Accessories for Delta-T determination .....	51
<u>Fig. 6:</u>	Accessories for GLP temperature test.....	121
<u>Fig. 7:</u>	Indicator lamps .....	132


## List of numbered parts and controls

<b>1</b>	Reaction vessel.....	3,12,51	<b>22</b>	Connection piece .....	12
<b>2</b>	Reaction vessel cover.....	3,12,51	<b>23</b>	Opening "Out" .....	12
<b>3</b>	Connecting tubing .....	3,13	<b>24</b>	Opening "In" .....	12
<b>4</b>	Measuring vessel.....	3,13	<b>25</b>	Marking field.....	12
<b>5</b>	Measuring vessel cover .....	3,13	<b>26</b>	Connection plug.....	12
<b>6</b>	Connecting tubing .....	3,12,51	<b>27</b>	PTFE tube.....	12
<b>7</b>	Mains pilot lamp .....	3,132	<b>28</b>	Electrodes .....	12
<b>8</b>	Gas flow display .....	3,132	<b>29</b>	Protective ring .....	12
<b>9</b>	Temperature display .....	3,132	<b>30</b>	Connection nipple.....	12,51
<b>10</b>	Error display .....	3,132	<b>31</b>	Connection.....	12,51
<b>11</b>	Instrument display .....	3,132	<b>32</b>	Gas inlet tube (short) .....	12
<b>12</b>	Mains switch .....	5	<b>33</b>	Sealing ring .....	12,51
<b>13</b>	Mains connection plug .....	5	<b>34</b>	Connection.....	12,51
<b>14</b>	Fuse holder.....	5	<b>35</b>	Gas inlet tube (long).....	51
<b>15</b>	Model plate .....	5	<b>36</b>	Temperature sensor Pt100.....	51,121
<b>16</b>	PC connection .....	5	<b>37</b>	Distance piece .....	51,121
<b>17</b>	Pt100 connection.....	5	<b>38</b>	PTFE cylinder .....	121
<b>18</b>	Connection "N <sub>2</sub> in".....	5	<b>39</b>	Hexagonal screw.....	121
<b>19</b>	Gas supply connection.....	5	<b>40</b>	Allen key .....	121
<b>20</b>	Electrode connection.....	5	<b>41</b>	Aluminium cylinder.....	121
<b>21</b>	Mounting for exhaust collection tube .....	5			

# 1 Introduction

## 1.1 Determination of thermal stability

The determination of the thermal stability of polyvinyl chloride according to DIN 53 381, part 1, is based on the fact that PVC decomposes at higher temperatures with the release of HCl. The hydrochloric acid formed is taken up by a gas stream (normally N<sub>2</sub>) and transferred into the measuring vessel where it is absorbed in the measuring solution (dist. water). The decomposition process is monitored by measuring the conductivity of the aqueous HCl solution. Both the stability time and the induction time are determined. The **stability time** is defined as the time required to achieve a conductivity difference of 50 μS/cm. The **induction time** is the time that elapses before the conductivity starts to rise sharply. The method is used for testing PVC in all stages of processing as well as for testing stabilizers.



**Fig. 1:** Schematic measuring arrangement

## 1.2 Instrument description

The **763 PVC Thermomat** is a PC-controlled instrument for determining the thermal stability of PVC and similar halogenated plastics and is available in the following two versions:

- **2.763.0014 PVC Thermomat for 230 V, 50¼60 Hz**
- **2.763.0015 PVC Thermomat for 115 V, 50¼60 Hz**

The 763 PVC Thermomat is equipped with two **heating blocks** each with 4 measuring positions. Each block can be individually heated, i.e. 2 sets of 4 samples can be measured at two different temperatures or 8 samples can be measured at the same temperature. Measurements at the individual measuring positions can be started individually.


The complete **operation** of the 763 PVC Thermomat is carried out via a PC connected to the RS232 interface with the aid of the **«763 PVC Thermomat»** control and evaluation program. Up to 4 instruments can be connected to 1 PC allowing a maximum of 32 samples to be analyzed at the same time. The evaluation algorithm of the PC program determines the point of inflection of the Thermomat curve and therefore the **induction time** fully automatically. Apart from the induction time, the so-called **stability time**, i.e. the time taken until a certain alteration in the conductivity is reached, can also be determined. The results obtained can be subjected to further mathematical processes. In particular, induction times can be converted to the standard temperatures contained in the relevant standards.

Each Thermomat curve can also be evaluated **manually**. A PC-supported tangent method is available for this in which the tangents can be applied to the curve as required. This means that evaluations are possible even in extreme cases.

The results of the determinations are stored in a **database** together with all the data concerning the method and determination. In the results display part of the program determinations can be searched for, sorted, filtered, exported and printed out. Apart from the graphical display of single and multiple curves it is also possible to carry out recalculations with altered parameters and to extrapolate the results to a particular temperature.


**GLP** (Good Laboratory Practice) and instrument validation are continuously increasing in importance. In the 763 PVC Thermomat GLP tests are provided for the measurement of temperature, conductivity and gas flow. They determine whether and which tests are to be carried out. The time interval between the tests and the required accuracy can also be laid down as required. If the GLP function is selected then on each result report a remark will be made as to whether the GLP test requirements have been fulfilled. As an option Metrohm offers a GLP test set (6.5616.000), with which the most important tests can be carried out.

## 1.3 Operating elements



**Fig. 2:** Front of the 763 PVC Thermomat

<b>1</b>	<b>Reaction vessel (6.1429.040)</b>	<b>7</b>	<b>Mains pilot lamp</b> lit up when instrument switched on
<b>2</b>	<b>Reaction vessel cover (6.2753.100)</b>	<b>8</b>	<b>Gas flow display</b> lit up when gas flow switched on
<b>3</b>	<b>Connecting tubing (6.1816.010)</b> between reaction vessel and measuring vessel	<b>9</b>	<b>Temperature display</b> flashes when heating switched on lit up when temperature reached
<b>4</b>	<b>Measuring vessel (6.1428.100)</b>	<b>10</b>	<b>Error display (red)</b> lit up at error
<b>5</b>	<b>Measuring vessel cover (6.0913.130)</b> with built-in conductivity measuring cell	<b>11</b>	<b>Instrument display</b> indicates numbers (1...4) of units connected
<b>6</b>	<b>Connecting tubing (6.1805.080)</b> for gas supply into reaction vessel		



**Fig. 3: Rear of the 763 PVC Thermomat**

**12 Mains switch**

to switch instrument on and off:

I = ON    0 = OFF

---

**13 Mains connection plug**

mains connection, see *section 2.3*

---

**14 Fuse holder**

changing the fuses, see *section 2.3*

---

**15 Model plate**

indicating mains voltage and serial number

---

**16 PC connection**

RS232 interface

---

**17 Pt100 connection**

Connection for external temperature sensor

---

**18 Connection "N<sub>2</sub> in"****19 Gas supply connection**

for connection tubing **6**

---

**20 Electrode connection**

for measuring vessel cover **5**

---

**21 Mounting**

for 6.2757.000 Exhaust collection tube (option)

## 1.4 Information about these Instructions for Use



*Please read through these Instructions for Use carefully before you put the 763 PVC Thermomat into operation. The Instructions for Use contain information and warnings to which the user must pay attention in order to assure safe operation of the instrument.*

### 1.4.1 Organization

These **8.763.1003 Instructions for Use** for the 763 PVC Thermomat provide a comprehensive overview of the installation, startup procedure, operation, fault rectification and technical specifications of this instrument. The Instructions for Use are organized as follows:





- Section 1 Introduction**  
Measuring method, description of instrument, operating elements and safety notes
- Section 2 Installation**  
Installation of accessories, mains connection, connection to PC
- Section 3 Operating tutorial**  
Introduction to the operation using an example
- Section 4 Operation**  
Detailed description of the operation
- Section 5 Maintenance – Faults**  
Maintenance, fault display and fault rectification
- Section 6 Appendix**  
Technical data, standard equipment, options, warranty, declarations of conformity, index

To find the required information on the instruments, you will find it an advantage to use either the **Table of contents** or the **Index** at the back.

You will find additional information on the PVC Thermomat method in the relevant "**Application Bulletins**", which are available on request free of charge from your Metrohm agency, as well as in the "**Application Notes**" concerning the PVC Thermomat method, which can be downloaded from the Internet under [«www.metrohm.ch»](http://www.metrohm.ch).

**1.4.2 Notation and pictograms**

The following notations and pictograms (symbols) are used in these Instructions for Use:

Range	Menu item, parameter or entry value
763 PVC THERMOMAT RESULTS	<b>Program window</b>
<OK>	<b>Button</b>
[ Ctrl ]	<b>Key</b>
<b>35</b>	<b>Operating element of 763</b>
	<b>Hazard</b> This symbol draws attention to a possible danger to life or of injury if the associated directions are not followed correctly.
	<b>Warning</b> This symbol draws attention to possible damage to instruments or instrument parts if the associated directions are not followed correctly.
	<b>Caution</b> This symbol marks important information. First read the associated directions before you continue.
	<b>Comment</b> This symbol marks additional information and tips.

## 1.5 Safety notes

### 1.5.1 Electrical safety

While electrical safety in the handling of the 763 PVC Thermomat is assured in the context of the specifications IEC 1010-1 (protection class 1, degree of protection IP20), the following points should be noted:

- **Mains connection**



*Setting of the **mains voltage**, checking the **mains fuse** and the **mains connection** must be effected in accordance with the instructions in section 2.3.*

- **Opening the 763 PVC Thermomat**



*If the 763 PVC Thermomat is connected to the power supply, the instrument must not be opened nor must parts be removed from it, otherwise there is a danger of coming into contact with components which are live. Hence, always disconnect the instrument from all voltage sources before you open it and ensure that the **mains cable is disconnected from mains connection 13 !***

- **Protection against static charges**



*Electronic components are sensitive to static charging and can be destroyed by discharges. Before you touch any of the components inside the 763 PVC Thermomat, you should earth yourself and any tools you are using by touching an earthed object (e.g. housing of the instrument or a radiator) to eliminate any static charges which exist.*

### 1.5.2 General precautionary rules

- **Take care with hot reaction vessels**



*Avoid any contact with hot reaction vessels. Place them in the supplied vessel holder to cool them down.*

## 2 Installation

### 2.1 Setting up the instrument

#### 2.1.1 Packaging

The 763 PVC Thermomat is supplied together with the separately packed accessories in special packaging containing shock-absorbing foam linings designed to provide excellent protection. The instrument itself is packed in an evacuated polyethylene bag to prevent the ingress of dust. Please store all these special packaging as only they assure transport of the instrument free from damage.

#### 2.1.2 Check

After receipt, immediately check whether the shipment is complete and has arrived without damage (compare with delivery note and list of accessories in *section 6.2*). In the case of transport damage, see instructions in *section 6.4.1 "Warranty"*.

#### 2.1.3 Location

Position the instrument in the laboratory at a location convenient for operation, free from vibrations and protected against a corrosive atmosphere and contamination by chemicals. It is advisable to set up the instrument in a fume cupboard to guard against odor emission, but the PC can be positioned outside on the regular laboratory bench.



*In order to improve accessibility to the measuring places the instrument can also be placed on the 6.2059.000 Rotation ring, which is available as an option.*

## 2.2 Attaching the accessories

### 2.2.1 Accessories for external gas supply

Nitrogen (N<sub>2</sub>) is required for external gas supply to the 763 PVC Thermomat. The following accessories must be mounted on the rear panel of the PVC Thermomat:

---

#### 1 Fitting the connecting tubing

- Screw one end of connecting tubing 6.1805.030 (150 cm) onto connection **18 "N<sub>2</sub> in"** on the rear panel of the PVC Thermomat (see *Fig. 3*).
- Screw the 6.1808.020 connection piece onto the other end of connecting tubing 6.1805.030.

---

#### 2 Connect up gas supply

- Attach the gas supply from the cylinder containing nitrogen to the tubing nozzle of 6.1808.020 Connection piece.
- Set gas pressure at the nitrogen cylinder to ca. 1.5 bar using the reducing valve.

### 2.2.2 Reaction and measuring vessels

*Fig. 4* shows the details of how the accessories for measuring the thermal stability are to be attached and connected to each other. Proceed as follows:

---

#### 1 Fitting the measuring vessel cover

- Insert PTFE tube **27** from above into opening **24 "In"** of measuring vessel cover **5**.
- Screw connection piece **22** into opening **24** of measuring vessel cover **5**.

---

#### 2 Inserting the measuring vessel in the instrument

- Fill 60 mL dist. water into measuring vessel **4**.
- Place measuring vessel cover **5** on measuring vessel **4**.
- Insert measuring vessel **4** with attached measuring vessel cover **5** into the opening provided for it on the 763 PVC Thermomat and connect measuring vessel cover **5** with connection plug **26** to electrode connection **20** (see *Fig. 3*).

---

#### 3 Fitting the reaction vessel

- Insert gas inlet tube **32** from below into connection **34** of reaction vessel cover **2**.
- Push sealing ring **33** over the upper end of gas inlet tube **32**.
- Screw connection nipple **30** loosely into connection **34** and press gas inlet tube **32** against connection nipple **30** from below; then fix gas inlet tube **32** to reaction vessel cover **2** by firmly tightening connection nipple **30**.

- Place reaction vessel cover **2** on reaction vessel **1**.
- Connect the white connecting tubing **3** to connection **31** of reaction vessel cover **2** (see *Fig. 3*).

---

**4 Inserting the reaction vessel into the instrument**

- Place the sample in the reaction vessel (see section 4.6.2).
- When the required temperature has been reached insert reaction vessel **1** with attached reaction vessel cover **2** in the openings provided for it on the 763 PVC Thermomat.

---

**5 Making the tubing connections**

- Connect the white connecting tubing **3** attached to reaction vessel cover **2** to connection piece **22** of measuring vessel cover **5**.
- Screw one end of brown connecting tubing **6** to connection nipple **30** of reaction vessel cover **2**.
- Screw the other end of brown connecting tubing **6** to the corresponding connection **19** on the 763 PVC Thermomat (see *Fig. 3*).



*The 6.1428.020 transparent glass measuring vessel available as an option can be used instead of the 6.1428.100 polycarbonate measuring vessel; in contrast to the polycarbonate vessel this can be cleaned with acetone.*

**2.2.3 Exhaust collection tube (option)**

For selective exhaust lead off the **6.2757.000 Exhaust collection tube** available as an option can be fitted to the 763 PVC Thermomat. In addition to the exhaust collection tube 8 x **6.1816.010 Silicone tubing** (22 cm) must also be ordered. Attach the exhaust collection tube as described below:

---


**1 Positioning the exhaust collection tube**

- Place the exhaust collection tube with the two supports in the two mountings **21** (see *Fig. 3*) on the 763 PVC Thermomat so that the connection for leading off the exhaust is located to the rear.

---

**2 Connecting the measuring vessel**

- Screw a connection piece **22** into opening **23** of the measuring vessel cover **5**.



**Fig. 4:** Installation of reaction and measuring vessels

<b>1</b>	Reaction vessel (6.1429.040)	<b>26</b>	Connection plug
<b>2</b>	Reaction vessel cover (6.2753.100)	<b>27</b>	PTFE tube (6.1819.080) for air supply into measuring vessel
<b>3</b>	Connecting tubing (6.1816.010)	<b>28</b>	Electrodes
<b>4</b>	Measuring vessel (6.1428.100)	<b>29</b>	Protective ring (4.422.0520)
<b>5</b>	Measuring vessel cover (6.0913.130) with built-in conductivity measuring cell	<b>30</b>	Connection nipple (6.1808.090)
<b>6</b>	Connecting tubing (6.1805.080) for gas supply to reaction vessel	<b>31</b>	Connection for connecting tubing <b>3</b>
<b>22</b>	Connection piece (6.1808.050)	<b>32</b>	Gas inlet tube (6.2418.110)
<b>23</b>	Opening "Out" for measuring vessel exhaust	<b>33</b>	Sealing ring (6.1454.040)
<b>24</b>	Opening "In" for measuring vessel inlet	<b>34</b>	Connection for connection nipple <b>30</b>
<b>25</b>	Marking field for entering the cell constant		

- Connect one end of the 6.1816.010 Silicone tubing to connection piece **22**.
- Plug the other end of the 6.1816.010 Silicone tubing into the corresponding exhaust collection tube opening.
- Close the unused openings on the exhaust collection tube with the included E.400.0010 stoppers.

### 3 Connecting the exhaust collection tube

- Attach a piece of suitable tubing to the exhaust collection tube connection and connect it to an active aspirating device (e.g. water-jet pump).

## 2.3 Mains connection



*Follow the instructions below for connecting to the power supply. If the instrument is operated with a mains voltage set wrongly and/or wrong mains fuse, there is a danger of fire!*

### 2.3.1 Checking the mains voltage

Before switching on the 763 PVC Thermomat for the first time, check that the mains voltage indicated on the model plate **15** (see *Fig. 3*) matches the local mains voltage. If this is not the case, call the Metrohm service.

### 2.3.2 Exchange of fuses

Two fuses 4 A/slow-blow for 115 V or 2 A/slow-blow for 230 V are installed in fuse holder **14** of the 763 PVC Thermomat as standard.



*Ensure that the instrument is never put into operation with fuses of another type, otherwise there is danger of fire!*

For changing fuses, proceed as follows:

#### 1 Disconnect mains cable

Disconnect mains cable from mains connection plug **13** of the 763 PVC Thermomat.

#### 2 Remove fuse holder

Using a screwdriver, loosen fuse holder **14** above the mains connection plug **13** and take out completely.

**3 Change fuses**

Remove old fuses carefully from fuse holder **14** and replace them with two new fuses which are suitable for the set mains voltage, type TH (slow-blow, with high switching capacity):

<b>115 V</b>	<b>4 A (TH)</b>	Metrohm-No. U.600.0022
<b>230 V</b>	<b>2 A (TH)</b>	Metrohm-No. U.600.0107

**4 Install fuse holder**

Push fuse holder **14** back into the instrument until it clicks into position.

**2.3.3 Mains cable and mains connection****Mains cable**

The instrument is supplied with one of three mains cables

- 6.2122.020 with plug SEV 12 (Switzerland, ...)
- 6.2122.040 with plug CEE(7), VII (Germany, ...)
- 6.2133.070 with plug NEMA 5-15 (USA, ...)

which are three-cored and fitted with a plug with an earthing pin. If a different plug has to be fitted, the yellow/green lead (IEC standard) must be connected to protective earth (protection class 1).



*Any break in the earthing inside or outside the instrument can make it a hazard!*

**Mains connection**

Plug the mains cable into mains connection plug **13** of the 763 PVC Thermomat (see *Fig. 3*).

**2.3.4 Switching the instrument on/off**

The 763 PVC Thermomat is switched on and off using mains switch **12**. When the instrument is switched on, the mains pilot lamp **7** "ON" (see *Fig. 2*) lights up.

## 2.4 Connection to the PC

### 2.4.1 Connection PVC Thermomat – PC



*Always switch off 763 PVC Thermomat and PC before you connect the two instruments with the 6.2134.100 Cable.*

The PC program «763 PVC Thermomat» allows a maximum of 4 PVC Thermomats to be controlled. The following possibilities exist for connecting these PVC Thermomats to serial PC interfaces:

- Connection to built-in COM interfaces on the PC
- Connection to an additional built-in interface extension board
- Connection to a softswitch (order number: 2.145.0014 for 230 V, 2.145.0015 for 115 V) which is connected to one of the COM interfaces.

Connect the PC connection **16** at the 763 PVC Thermomat to one of the serial COM ports at the PC or to a softswitch using the 6.2134.100 Cable (9 pin/9 pin). If only a 25-pin COM interface is available on the PC then the 6.2125.110 Adapter cable or a commercially available adapter must be used.

### 2.4.2 Software installation

The PC program «**763 PVC Thermomat 1.0**» is required for the operation of the 763 PVC Thermomat; this is contained on the 6.6037.000 CD included in the accessories. This program runs under Windows NT (recommended, especially if several instruments are connected), Windows 95 and Windows 98 operating systems and is installed as follows:

#### 1 Install program

- For Windows NT: log in as operator with administrator rights.
- Insert 6.6037.000 Installation CD into CD drive.
- Select <Start> and **Run**. Browse for the **setup.exe** file on the installation CD and click on <OK>.
- Click on "763" and follow the instructions given in the setup program.

The software package will be installed in the desired directory. In addition to the program files, the following folders are installed:

- Database** Folder for database files (\*.mrd)
- Log** Folder for event log files (\*.mel) and temperature recording files (\*.txt)
- Template** Folder with empty database (**Default.mrd**) and copy of the result program (**Nachaus.prg**)

#### 2 Registration

- Please send us your 8.763.8007 Registration card as soon as possible.

### 2.4.3 Basic settings


When the program is started for the first time several basic settings must be made for the 763 PVC Thermomat. Proceed as follows:

#### 1 Switch on instruments

- Check if the PVC Thermomat has been connected correctly to the PC (see *section 2.4.1*).
- Switch on the PVC Thermomat using mains switch **12**.
- Switch on the PC.

#### 2 Start program and login

- Double-click the software icon to start the program. The program window with the opening picture is opened and the **Login** window appears on the screen:




- Do not enter any password here, just click on **<Change password>**. The following window appears:




- Leave the **Old password** field empty. Under **New password** enter the required password which will in future allow you to log in as **Administrator**. Confirm the new password by renewed entry in the **Confirm new password** field and click on **<OK>**.
- Enter the new **Password** again in the **Login** window. You are now logged in as **Administrator**.
- Confirm the message **There are no units configured yet** with **<OK>**. The **763 PVC THERMOMAT CONTROL** window appears.

### 3 Communication settings

- Click on **763 PVC THERMOMAT CONTROL / Options / Communication** to open the **Select number of units** window:



- Enter the number of PVC Thermomats (1...4) connected to the PC and click on <OK>. The **Unit configuration** window appears:



- Under **COM port** select for each instrument the number of the serial interface (1...8) to which this instrument is connected. If a softswitch is used then the option **Softswitch in use** must be switched on. Under **COM port** the number of the serial interface to which the softswitch is connected is selected; under **Softswitch port** the number of the interface on the softswitch to which the instrument is connected is selected.
- Click on <OK> to close the **Unit configuration** window.
- Confirm the message which then appears with <OK>.

### 4 Close the program and start it up again

- Close the program with **763 PVC THERMOMAT CONTROL / File / Exit**.
- Restart the «PVC Thermomat 1.0» program by double-clicking the program icon.
- In the **Login** window enter the **Password** for the administrator. The control window appears. In the status line the message **Unit#: Ready** should appear for all connected instruments.



## 3 Operating tutorial



*This section introduces you to the operation of the 763 PVC Thermomat by means of a brief operating tutorial which describes the basic operating steps needed for performing determinations and display of results.*

*For further explanations of the operation, please refer to section 4.*

### 3.1 Determinations



*The determination of induction and stability time of a PVC sample at 200 °C and a gas flow of 7 L/h is used as an illustrative example. Please note that the steps and parameter settings described apply only to this example. If you use a different sample and different parameters, the procedures described in the tutorial must be modified appropriately.*

#### 3.1.1 Installation

Before you start this brief tutorial, the 763 PVC Thermomat must be correctly installed as described in *section 2*. In what follows, the most important points for the installation are listed once again (for details, see the sections mentioned).

- |   |                      |
|---|----------------------|
| ⇒ Setting up instrument                     | <i>section 2.1</i>   |
| ⇒ Installing air supply accessories         | <i>section 2.2.1</i> |
| ⇒ Installing reaction and measuring vessels | <i>section 2.2.2</i> |
| ⇒ Mains connection                          | <i>section 2.3</i>   |
| ⇒ Connection to PC                          | <i>section 2.4</i>   |

### 3.1.2 Preparing for determinations

Before starting determinations, the method must be selected, the heating must be started and the reaction and measuring vessels must be prepared. Proceed as follows:

#### 1 Switch on 763 PVC Thermomat

⇒ Switch on 763 PVC Thermomat with mains switch **12** on the rear of the instrument. After the instrument has been switched on the mains pilot lamp **7** lights up.


#### 2 Switch on PC

⇒ Switch on PC and start «763 PVC Thermomat» program.

#### 3 Select method "Default"

⇒ In the opened control window, select the determination method **Default** for Block A (and – if more than 4 samples are to be determined – also for Block B) in the **Method** field.




⇒ Click on  , if you want to change the parameters of the default method (details see *section 4.5.2*).



*The determinations described in this tutorial have been carried out using the default settings (Temperature = 200 °C, Delta T = 0.6 °C, Gas flow = 7 L/h, etc.).*


#### 4 Start heating

⇒ To start the heating of Block A and – if desired – of Block B click on the following button:



As soon as the heating is switched on the color of the frame of this button switches to **red**. At the same time the temperature display **9 "TEMPERATURE"** LED on the 763 PVC Thermomat starts to blink for the selected block.

The actual temperature is shown digitally beside the button. A bar below shows the absolute variation of the actual temperature from the target temperature; the bar length corresponds to 50 °C.



**5 Prepare measuring vessels**


- ⇒ Clean the used measuring vessels and their covers (details see *section 4.6.3*).
- ⇒ Fill each cleaned measuring vessel **4** with **60 mL distilled water**.
- ⇒ Place the measuring vessel cover **5** fitted with a PTFE tube **27** on the measuring vessel **4** (see *Fig. 4*).
- ⇒ Place the measuring vessel **4** with its cover **5** in one of the openings provided for it on the 763 PVC Thermomat and connect the electrode plug **26** on the cover to the corresponding socket **20** on the 763 PVC Thermomat.

**6 Prepare reaction vessels**

- ⇒ Weigh out **0.5 g PVC** into the reaction vessels **1**.
- ⇒ Take the upper rim of the reaction vessel **1** in your hand (e.g. in the space between thumb and index finger) and rotate the glass through  $360^\circ$ . This provides the degreased glass with a thin grease film, without which it is very difficult to remove the cover from the vessel after the determination.
- ⇒ Insert an gas inlet tube **32** into the reaction vessel cover **2** and fasten it by screwing down the connection nipple **30** (see *Fig. 4*).
- ⇒ Place the reaction vessel cover **2** on the reaction vessel **1**. Rotate the cover so that the air inlet tube **32** is as close as possible to the vessel wall.
- ⇒ Fasten the white connecting tubing **3** on to the connection **31** of the reaction vessel cover **2**.
- ⇒ Place the prepared reaction vessel into the vessel holder.

**7 Enter sample identifiers**

- ⇒ Enter the sample identifiers **ID1** and **ID2** in the corresponding fields of the operating area in the control window.




The screenshot shows a control window titled "ID1 & ID2". It contains two input fields. The first field is a dropdown menu with "PVC" selected. The second field is a text box containing "KA356".

### 3.1.3 Carrying out determinations

When the preparations are complete the determinations can be started. Proceed as follows:

#### 1 Wait until the temperature has been reached

⇒ Wait until the heating block reaches the temperature defined in the method. This condition is indicated by temperature display **9** "TEMPERATURE" which stops blinking and remains lit up. At the same time or after a slight delay the frame color of the <Heat to...> button turns **green**.



#### 2 Insert reaction vessels


- ⇒ Place the prepared reaction vessels **1** in the heating block openings on the 763 PVC Thermomat. Close any openings which are not used with stoppers or empty reactions vessels to protect against contamination.
- ⇒ Connect the white connecting tubing **3** fastened to the reaction vessel cover **2** to the corresponding connection **22** on the measuring vessel cover **5**.
- ⇒ Screw the brown gas supply tubing **6** onto the connection nipple **30** of the reaction vessel cover **2** and the corresponding connection **19** on the 763 PVC Thermomat (see *Fig. 3*).

#### 3 Start determination


⇒ In the operating field of the control window for every determination you want to start you must click on the button



After the start of a determination the **live curve** of the started determination is shown in the operating window beside the channel number. The actual measured value appears above the curve.



As soon as a point of inflection is found during the measurement the determination is terminated. The window with the live curve becomes gray again and the message **finished** appears.



**4 Automatic report output**

At the end of the determination the report defined in the method is printed out automatically.

<b>Unit</b>	1	<b>Serial number</b>	01109	<b>Druckdatum</b>	09.09.1999 09:34:03
<b>Block</b>	A	<b>Cell constant</b>	1	<b>Bestimmungsdatum</b>	23.06.1999 13:25:52
<b>Channel</b>	1	<b>User</b>	Administrator		

---

<b>ID 1</b>	<b>ID 2</b>
PVC	KA356

---

<b>Induction time</b>	1.27 h		
<b>Stability time</b>	1.33 h	<b>Delta Kappa</b>	50 µS/cm
<b>Standard time</b>	6.45 h	(F: 2.254; T: 180°C)	

---

<b>Method name</b>	PVC 200°C	<b>Creation date</b>	23.06.1999 08:27:20
<b>Creator</b>	Administrator		

---

<b>Temperature</b>	200 °C	<b>Stop time</b>	0.00 h
<b>Delta T</b>	0.60 °C	<b>Stop at conductivity</b>	200 µS/cm
<b>Gas flow</b>	7 L/h	<b>Stopat endpoint</b>	<input checked="" type="checkbox"/>

---

<b>Start delay</b>	0 min	<b>Delta Kappa</b>	50 µS/cm
<b>Start mode</b>	<input checked="" type="radio"/> per channel <input type="radio"/> per 2 channels <input type="radio"/> per block	<b>Evaluation delay</b>	0.00 h
		<b>Evaluation suppression start</b>	0.00 h
		<b>Evaluation suppression end</b>	0.00 h

---

**Send to file**

GLP test passed.

## 3.2 Results




The demo database **Demo.mrd** is used as an example for the display and further processing of the results of recorded determinations. Please note that the steps and parameter settings described apply only to this database. If you use a different database, the procedures described in the tutorial must be modified appropriately.

### 3.2.1 Determination overview



The determination overview displays an overview table with selected data fields of all the determinations contained in the opened database.

#### 1 Open database "demo.mrd"

- ⇒ Click on  or **File / Results** in the control window in order to open the **763 PVC THERMOMAT RESULTS** window. This loads the **Repos.mrd** database and opens the determination overview as standard.
- ⇒ Click on **File / Open database** in the results window. Select the **demo.mrd** database in the **Open database** window and click on **<Open>**. A sub-window opens in the results window which shows the data from 5 determinations in an overview table.

Determination data					
	ID 1	ID 2	Determination date	Temperature	Induction time automatic
	PVC	KA356	24.06.1999 13:32:43	200°C	1.27 h
▶	PVC	KA356	24.06.1999 13:32:41	190°C	2.90 h
	PVC	KA356	22.06.1999 13:47:30	170°C	14.14 h
	PVC	KA356	21.06.1999 17:53:40	160°C	33.43 h
	PVC	KA356	21.06.1999 17:53:38	180°C	6.33 h


#### 2 Format determination overview


- ⇒ **Select fields:** Click on  or **Format / Select fields** in the results window. In the **Display fields** window move the fields until the fields required for the determination overview are visible in the **Display** field and then click on **<OK>**.
- ⇒ **Adjust field width:** Mark the column in the determination overview whose width is to be altered. Click on **Format / Adjust width** in the results window. Enter the column width in the **Column width** window or click on **<Best fit>** in order to automatically adapt the column width to the longest field contents. Alternatively the column width can also be set directly by using the **mouse**. This is done by moving the cursor between two columns in the title line of the table until it appears as . The column can now be pulled to the required width with the **left-hand mouse key** pressed down. A **double-click** on the

selected column automatically adapts the column width to the longest field contents (including title).

- ⇒ **Format fields:** Click on **Format / Fields** in the results window. Set down the format for the fields in the determination overview in the **Cells effects** window and click on **<OK>**.
- ⇒ **Select fonts:** Click on **Format / Fonts** in the results window. Select the font to be used in the determination overview in the **Font** window and click on **<OK>**.
- ⇒ **Save settings:** Click on **Options / Save settings now** to save the new format of the determination overview permanently.




### 3 Print determination overview

- ⇒ Click on  or **File / Print** in the results window to open the **Printing of** window.



- ⇒ Select **Overview table** if the determination overview is to be printed out in tabular form.
- ⇒ Select **Overview list** if the determination overview is to be printed out as a list with all the fields selected in the determination overview.
- ⇒ Click on **<Print>**.
- ⇒ In the **Print** window select the printer, printing area and number of copies and click on **<Print>**.


#### 4 Export determination overview


- ⇒ **Copy to clipboard:** Mark the lines or columns in the determination overview which are to be copied. Click on  or **E**dit / **C**opy. The marked texts are loaded into the clipboard from where they can be further processed with other programs.
- ⇒ **Export to other database:** Mark the determinations (lines) which are to be exported in the determination overview. Click on **D**etermination / **E**xport / **E**xport to database. In the **Export data to...** window select the database to which the determinations are to be exported and click on <Save>.
- ⇒ **Export to Word:** Mark the lines or columns in the determination overview which are to be exported. Click on  or **D**etermination / **E**xport / **E**xport to MS **W**ord. The marked texts are loaded into an RTF file which is automatically opened with Microsoft Word.
- ⇒ **Export to Excel:** Mark the lines or columns in the determination overview which are to be exported. Click on  or **D**etermination / **E**xport / **E**xport to MS **E**xcel. The marked texts are loaded into an XLS file which is automatically opened with Microsoft Excel.

### 3.2.2 Determination and method data


All the determination and method data of each determination in a determination overview can be shown, printed out and exported.

#### 1 Display determination and method data


- ⇒ In the opened determination overview mark the determination for which all parameters and results are to be displayed.
- ⇒ Click on  or **V**iew / **D**etermination and method data in the results window or double-click on the required line in the determination overview. The **Determination and method data** window opens in which the data are shown on several register cards.



## 2 Print determination and method data

- ⇒ Click on  or **File / Print** in the results window to open the **Printing of** window (see *section 3.2.1*).
- ⇒ Select the **All determination and method data** option.
- ⇒ Click on **<Print>**.
- ⇒ In the **Print** window select the printer, printing area and number of copies and click on **<Print>**.


## 3 Export determination and method data


- ⇒ Mark the determinations (lines) in the determination overview from which all data is to be exported.
- ⇒ Click on  or **Determination / Export / Determination and method data**.
- ⇒ In the **Save to export file** window select the required folder and the file name of the \*.txt export file under which the data are to be stored and click on **<Save>**.

### 3.2.3 Graphics




Curves for each determination in a determination overview can be shown, printed out and copied.


#### 1 Open single graph

- ⇒ Mark the determination in the determination overview whose single curve is to be displayed.
- ⇒ Click on  or **Determination / Graphics / Single graph** in the results window. A new window is opened in which the **single graph** for the determination selected in the determination overview is shown. The sample identifiers **ID 1** and **ID 2** are shown in the window title. If present then the determined **Induction time** (red vertical line and numbers) and **Stability time** (black vertical line and numbers) will be shown in the graphics window as well as the **measuring curve** (green).





#### 2 Open multiple graph

- ⇒ Mark the determinations in the determination overview whose curves are to be displayed in a multiple graph.
- ⇒ Click on  or **Determination / Graphics / Multiple graph** in the results window. A new window is opened in which **all curves** for the determinations marked in the determination overview are shown. The individual measuring curves are shown in color, the legends referring to the colors (**ID 1** and **ID 2**) are located in the graph header. If present the determined **Induction time** (  ) and **Stability time** (  ) will also be shown.





### 3 Open live graph

- ⇒ Mark the determination in the determination overview whose curve is to be displayed.
- ⇒ Click on  or **Determination / Graphics / Live graph with reevaluation** in the results window. A new window is opened in which the **measuring curve** of the marked determination is shown. If a different determination is marked in the determination overview then the curve is automatically updated. If present then the determined **Induction time** (red vertical line and numbers) and **Stability time** (black vertical line and numbers) will be shown in the graphics window as well as the **measuring curve** (green). With **<2nd derivative>** the second derivative of the measuring curve can be shown as a blue line. The live graph offers the possibility of the manual **reevaluation** of the measuring curves (see section 4.7.3).



---

#### 4 Print graphics

- ⇒ Open the required graph in the results window which is to be printed.
- ⇒ Click on  or **File / Print** in the results window. The curve shown in the active graph window will be printed out.
- ⇒ If you want to print several single or multiple curves at the same time first mark the required determinations in the determination overview. Click on  or **File / Print** (see section 3.2.1). In the **Printing of...** window select the option **Single graph** or **Multiple graph** and click on <Print>.

---

#### 5 Copy graphics

- ⇒ Open the required graph in the results window which is to be copied.
- ⇒ Use the right-hand mouse key to click on the graph window and select the option **Copy** in the context-sensitive menu. The selected curve is copied into the clipboard where it can be further processed with other programs.

# 4 Operation



*This section describes the most important points concerning the operation of the 763 PVC Thermomat. For further details please refer to the on-line help in the PC program which can provide you with the required information rapidly and conveniently from any place in the program.*

## 4.1 Fundamentals of the operation


### 4.1.1 Starting/closing the program

#### Start the «763 PVC Thermomat» program



#### Start the program

Double-click this icon or the **PvcTm.exe** file to start the «763 PVC Thermomat 1.0» program. The **Login** window appears:



Select the desired **User** name and enter the **Password**. Then click on <OK>.



*After software installation, the program can be started without entering a **Password**. For the definition of users, see section 4.2.2.*

#### Close the «763 PVC Thermomat» program

**763 PVC THERMOMAT CONTROL / File / Exit**

Exit the «763 PVC Thermomat» program.

The program is also quitted by clicking on  in the upper right part of the **763 PVC THERMOMAT CONTROL** window.

## 4.1.2 Glossary

### **Control window**

The main window **763 PVC THERMOMAT CONTROL** is known as the control window. It contains all the functions required for controlling the 763 PVC Thermomats connected to the PC.

### **Method**

A method includes all parameters for carrying out and evaluating a determination.

### **Determination**

A determination is understood to be the automatic determination of the induction time and/or stability time of a sample. In order to carry out a determination a suitable method must be selected for the samples.

### **Results**

The results of a determination are automatically stored in the **Repos.mrd** database and can be observed in the results window.

### **Reevaluation**

Reevaluation is the subsequent reprocessing of determinations, in particular the manual determination of the induction time with the aid of tangents.

### **Recalculation**

Recalculation means that results, equations and standards can be calculated again at a later date.

### **Extrapolation**

Extrapolation can be used to convert the results obtained at various temperatures to a required target temperature. In addition, this method can be used for converting the induction time to the standard time.

### 4.1.3 Control window


The **763 PVC THERMOMAT CONTROL** window is used for controlling the PVC Thermomats, managing methods, displaying the live curves and accessing various program functions.

#### Open

The **763 PVC THERMOMAT CONTROL** window is opened automatically when the program is started (see *section 4.1.1*).

#### Arrangement

The elements of the control window are the menu bar, the toolbar and the status bar, indicating prompts, the logged-in user and the instrument status. The operating area contains a diagram of the measuring unit of the 763 PVC Thermomat which can be used to start, display and stop determinations.



#### Menus














The control window contains the following main menus:

- File** Management of method files, opening the results window, printing, new log in.
- View** View of toolbar and status bar, unit information, live parameters, status overview, event log overview, unit selection.


<b>Tools</b>	Cell constants, GLP test, timer, gas flow control, start all channels, determination of Delta T, temperature recording, service diagnosis.
<b>Options</b>	General settings, instrument configuration, user permissions.
<b>Help</b>	Program-specific on-line help.

## Icons

The following icons are displayed in the control window:

	Create new method (see <i>section 4.5.1</i> ).
	Open existing method (see <i>section 4.5.1</i> ).
	Method manager (see <i>section 4.5.1</i> ).
	Print results (see <i>section 4.7.6</i> ).
	Open results window (see <i>section 4.1.4</i> ).
	Live parameters, which can be altered during a run (see <i>section 4.6.4</i> ).
	Status overview of the connected instruments (see <i>section 4.3.2</i> ).
	Unit information (see <i>section 4.3.1</i> ).
	Display GLP status (see <i>section 4.8.3</i> ).
	Select unit 1...4.
	Start all channels of Block A (see <i>section 4.6.4</i> ).
	Start all channels of Block B (see <i>section 4.6.4</i> ).
	Call up help.


## Close

The control window is closed by clicking on **763 PVC THERMOMAT CONTROL / File / Exit** or  in the upper right part of the window.

### 4.1.4 Results window


The **763 PVC THERMOMAT RESULTS** window is used to display, output and recalculate the results of determinations recorded by the 763 PVC Thermomat. The determination data are stored in the \*.mrd database files and can be shown in this window as an overview table and as curves. The results window can run even when the control window is closed.

#### Open

The **763 PVC THERMOMAT RESULTS** window is opened by clicking on  or **763 PVC THERMOMAT CONTROL / File / Results** in the control window.

#### Arrangement

The elements of the results window are the menu bar, the toolbar and the status bar. Within the results window sub-windows can be opened showing determination overviews, single, multiple and live graphs.



#### Menus


























The results window contains the following main menus:

- File**            Open database, print, close window.
- Edit**            Copy, mark, delete filter.
- View**            Select display: determination overview, method and determination data, GLP.


<b>Determination</b>	Seek, sort, filter, graphics, extrapolation, recalculation, export, delete.
<b>Format</b>	Format determination overview.
<b>Options</b>	General program settings.
<b>Window</b>	Display of windows.
<b>Help</b>	Program-specific on-line help.

## Icons

The following icons are shown in the results window:

	Print results, curves and overview lists (see <i>section 4.7.6</i> ).		Filter selection (see <i>section 4.7.1</i> ).
	Copy to clipboard (see <i>section 4.7.6</i> ).		Selection not in filter (see <i>section 4.7.1</i> ).
	Select fields for determination overview (see <i>section 4.7.1</i> ).		Single graph (see <i>section 4.7.3</i> ).
	Show determination overview (see <i>section 4.7.1</i> ).		Multiple graph (see <i>section 4.7.3</i> ).
	Show all method and determination data (see <i>section 4.7.2</i> ).		Single graph with possibility of reevaluation (see <i>section 4.7.3</i> ).
	Find (see <i>section 4.7.1</i> ).		Export to MS Word (see <i>section 4.7.6</i> ).
	Sort ascending (see <i>section 4.7.1</i> ).		Export to MS Excel (see <i>section 4.7.6</i> ).
	Sort descending (see <i>section 4.7.1</i> ).		Export measured values (see <i>section 4.7.6</i> ).
	Selection based filter (see <i>section 4.7.1</i> ).		Export determination and method data (see <i>section 4.7.6</i> ).
	Special filter/sorting (see <i>section 4.7.1</i> ).		Tile windows vertically (see <i>section 4.7.8</i> ).
	Apply filter/sorting (see <i>section 4.7.1</i> ).		Tile windows horizontally (see <i>section 4.7.8</i> ).
	Remove filter/sorting (see <i>section 4.7.1</i> ).		Close results window.
	Delete filter (see <i>section 4.7.1</i> ).		

## Close

The results window is closed by clicking on  or **763 PVC THERMO-MAT RESULTS / File / Return**.

### 4.1.5 Overview of file types

The following file types are produced by the «763 PVC Thermomat» program:

- \*.mrd      Database file**  
Contains measuring data and results. \*.mrd files are automatically stored in the **Database** folder.
- \*.mel      Event log file**  
This file contains a protocol of all events which have occurred with the connected PVC Thermomats. \*.mel files are automatically stored in the **Log** folder.
- \*.txt      Measurement table file**  
Contains measurement tables in ASCII format. \*.txt files are automatically stored in the **Log** folder.


### 4.1.6 Context-sensitive menus

A lot of menu functions of the program window can also be selected by clicking the required window or element with the **right-hand mouse key**. The menu options which are then opened depend on the selected active window or element

### 4.1.7 Mouse functions

The **mouse** can be used to carry out the normal program operating functions such as the selection of menu items and fields. It can additionally be used for magnifying a section of a curve (**zooming**). To zoom a portion of the curve place the mouse cursor to the upper left corner of the square to zoom, press the left mouse button and drag the cursor to the lower right corner of the rectangle. After releasing of the left mouse button the selected region will be zoomed full-screen.

### 4.1.8 Help

By clicking on  or **Help / Help** or by pressing the [ F1 ] key you can get on-line help on the current topic anywhere in the program.

- Green texts*                      can be clicked to jump to a different Help topic.
- Violet texts*                     identify dialog items, parameters or buttons in the corresponding window.
- Blue texts*                        Identify titles and important information.


## 4.2 Instrument and software settings

### 4.2.1 Communication


#### 763 PVC THERMOMAT CONTROL / Options / Communication

This menu item is used for the configuration of the communication between the PC and up to 4 connected 763 PVC Thermomats.

The **Select number of units** window appears first; the number of 763 PVC Thermomats connected to the PC must be entered here.



Then the **Unit configuration** window opens; the following parameters can be selected here for each instrument which is connected:



**COM port** Selects the COM port (1...8) to which the instrument or the softswitch is connected.

#### Softswitch

**Softswitch in use** Use a softswitch between PC and instrument.

**Softswitch port** Selects the port (1...4) on the softswitch to which the instrument is connected.




*Two softswitches are available as options from Metrohm: 2.145.0014 (for 230 V) or 2.145.0015 (for 115 V). 4 PVC Thermomats can be connected to a softswitch.*

### 4.2.2 User permissions

The «763 PVC Thermomat» program is equipped with a comprehensive **password protection** system which allows each menu item to be individually provided with **access rights**. The access rights are issued for various **groups**. The individual **user** can then be assigned to a group

#### 763 PVC THERMOMAT CONTROL / Options / User permissions


This menu item is used to issue access rights and is reserved for the **Administrator** only. The **Users and groups manager** window opens in which the already defined groups with their users are listed.



The following menu items are available by clicking on the entries with the right-hand mouse key.

#### Functions for groups

- Add user**                      Adds a new user to the selected group. The name of the new user must be entered in the **Enter new user name** window.
- Add group**                      Adds a new group to the previous groups. The name of the new group must be entered in the **Enter new group name** window.
- Delete group**                      Deletes the selected group.
- Properties**                      The **Group permissions** window opens in which access to all program functions can be switched on or off individually.



### Functions for users

**Delete user**                      Deletes the selected user.

It is recommended to make user lists and enter passwords as a first action after software installation.




*The **Administrator** is not included in the list of operators but possesses all rights of access and cannot be deleted.*

After configuration of user permissions the program prompts for the selection of the user name and entry of the password in the **Login** window every time the program starts. The user name stamps all methods, determinations and reports. It is possible at any time to change the user with the menu item **763 PVC THERMOMAT CONTROL / File / New login**.

### 4.2.3 Timer

763 PVC THERMOMAT CONTROL / Tools / Timer function...

The **Automatic heater start** window opens in which the automatic start of heating at a predetermined time can be entered independently for block A or block B.



**Start heater**

- Never** No automatic start of heating.
- Once** Automatic start of heating at the time defined in the **Start at...** field.
- Daily** Automatic start of heating on the day selected in the **<Which days>** field at the time defined in the **Start at...** field.

---

**<Which days>** This button opens the **Select days** window in which the required day can be selected.

---

**Start at...** If **Once** has been selected then the time and date for the automatic start of heating can be entered, if **Daily** has been selected the time.

### 4.2.4 Gas flow control


The flow of gas through the reaction vessels to the measuring vessels can be switched on and off manually and shown in a separate window.

763 PVC THERMOMAT CONTROL / Tools / Gas flow control / Gas flow on  
Switches on gas flow.

763 PVC THERMOMAT CONTROL / Tools / Gas flow control / Gas flow off  
Switches off gas flow.

**763 PVC THERMOMAT CONTROL / Tools / Gas flow control / Display gas flow...**

The **Gas flow display** window opens in which the actual gas flow is shown. This function is only available if no determination is running.



#### 4.2.5 Recording temperature

**763 PVC THERMOMAT CONTROL / Tools / Temp. logging / Block A on**

Switches on temperature recording for block A and external temperature sensor.

**763 PVC THERMOMAT CONTROL / Tools / Temp. logging / Block A off**

Switches off temperature recording for block A and external temperature sensor.

**763 PVC THERMOMAT CONTROL / Tools / Temp. logging / Block B on**

Switches on temperature recording for block B and external temperature sensor.

**763 PVC THERMOMAT CONTROL / Tools / Temp. logging / Block B off**

Switches off temperature recording for block B and external temperature sensor.

The temperature can be individually recorded at any time for each heating block. The measured values are written in a text file which contains the following values:

<b>t(s)</b>	Time in s
<b>Block [°C]</b>	Measured temperature in heating block
<b>Ext. Sensor [°C]</b>	Measured temperature with the external temperature sensor

The text file is stored under the name **U#X-YYMMDD-HHMMSS.txt** in the folder **Log**; # stands for the instrument number (1...4), X for the block (A or B) and **YYMMDD-HHMMSS** for the date and time of the start of data recording.

## 4.2.6 Options

### Save settings

The program and window settings can be saved with the two following menu items:

#### 763 PVC THERMOMAT CONTROL / Options / Save settings on exit

If this option is switched on then the current program and window settings will be automatically saved each time the program is shut down.

#### 763 PVC THERMOMAT CONTROL / Options / Save settings now

The current program and window settings will be saved immediately.

### Clear sample ID1 history

#### 763 PVC THERMOMAT CONTROL / Options / Clear sample ID history

This menu item is used to delete the list automatically formed from the entered sample identifiers **ID1**.


## 4.3 Program information

### 4.3.1 Unit information



763 PVC THERMOMAT CONTROL / View / Unit information

The Unit information and Service diagnosis window opens.



<b>Serial number</b>	Serial number of the selected 763 PVC Thermomat.
<b>Program number</b>	Number of the EEPROM program of the selected 763 PVC Thermomat.
<b>Last maintenance</b>	Date of the last maintenance with the signature of the service technician who carried out the maintenance work.
<b>Last calibration</b>	Date of the last adjustment with signature and code of the office which carried it out.
<b>Unit operating hours</b>	Number of hours that the instrument has been switched on.
<b>Pump operating hours</b>	Number of hours that the gas flow has been switched on.
<b>&lt;Print report&gt;</b>	Print out unit information report.



*The Service Diagnosis block is password-protected and is only accessible to trained service personnel.*

4.3.2 Status overview



763 PVC THERMOMAT CONTROL / View / Status overview

The **STATUS - 763 PVC THERMOMAT CONTROL** window opens and displays an overview of all available channels of the instruments connected to the PC.

S Status - 763 PVC Thermomat Control									
Unit	Block	Channel	Method	Sample ID1	Status	Stab. time	Induc. time	Set temp.	Current temp.
1	A	1	Default		ready			121.60°C	23.29°C
1	A	2	Default		ready			121.60°C	23.29°C
1	A	3	Default		ready			121.60°C	23.29°C
1	A	4	Default		ready			121.60°C	23.29°C
1	B	1	Default		ready			121.60°C	22.88°C
1	B	2	Default		ready			121.60°C	22.88°C
1	B	3	Default		ready			121.60°C	22.88°C
1	B	4	Default		ready			121.60°C	22.88°C

This window contains the following columns:

<b>Unit</b>	Instrument number (1...4).
<b>Block</b>	Instrument block (A, B).
<b>Channel</b>	Channel (1...4).
<b>Method</b>	Method name.
<b>Sample ID1</b>	Sample identifier 1.
<b>Status</b>	Status of channel:
<b>ready</b>	No active measurement. The channel is ready for the start of a determination.
<b>running</b>	Running measurement.
<b>finished</b>	Determination finished. The channel is ready for the start of a new determination.
<b>error</b>	Communication error between instrument and PC.
<b>Stab. time</b>	Determined stability time.
<b>Induc. time</b>	Determined induction time.
<b>Set temp.</b>	Set temperature defined in method.
<b>Current temp.</b>	Actual block temperature.


The width of the columns in the **STATUS - 763 PVC THERMOMAT CONTROL** window can be altered by pulling the right-hand field margin of the column title with the mouse. These settings together with the window size and window position can be permanently stored with **763 PVC THERMOMAT CONTROL / Options / Save settings now.**













### 4.3.3 Events overview

#### Event window




##### 763 PVC THERMOMAT CONTROL / View / Event log overview

The **PVCTM.MEL - MLOGVIEW** window opens in which an overview of all the events concerning the 763 PVC Thermomats connected to the PC is displayed.



Type	Date	Time	User	Unit	Event
	29.09.1999	16:19:06	Administrator	1	Heater A off
	29.09.1999	16:18:46	Administrator	1	Heater A on
	29.09.1999	16:18:16	Administrator	1	Connection established (serial no. 1111)
	29.09.1999	16:12:27	Administrator	1	Gas flow off
	29.09.1999	16:12:22	Administrator	1	Gas flow on
	29.09.1999	16:09:50	Administrator	1	Connection established (serial no. 1111)
	29.09.1999	16:09:05	Administrator	n.a.	Permission settings edited
	29.09.1999	16:08:52	Administrator	n.a.	User "Administrator" logged in
	29.09.1999	16:08:51	n.a.	n.a.	Program start
	29.09.1999	16:08:31	Administrator	n.a.	Program closed
	29.09.1999	15:31:46	Administrator	1	Connection established (serial no. 1111)
	29.09.1999	15:31:44	Administrator	n.a.	User "Administrator" logged in

This window contains the following columns:

<b>Type</b>	Type of event:  Information about an event which has taken place properly.  Warning about a special event.  Error message about the faulty event.
<b>Date</b>	Date of the event. The format depends on the settings defined under <b>Regional Settings / Date</b> in Windows.
<b>Time</b>	Time of the event. The format depends on the settings defined under <b>Regional Settings / Time</b> in Windows.
<b>User</b>	Name of the operator logged in when the event occurred.
<b>Unit</b>	Instrument number(1...4).
<b>Event</b>	Description of the event.






**Menus**

The events window contains the following main menus:

- File**           Open and save event log files.
- View**           View toolbar and status bar, update.
- Filter**          Filter the events.
- Help**           Program-specific on-line help.

**Icons**

The following icons are shown in the events window:

-  Open existing event log file (\*.mel).
-  Save event log file as text file (\*.mel).
-  Update event overview.
-  Selection-based filter.
-  Remove filter.

**Filter**

 **PVCTM.MEL - MLOGVIEW / Filter / Selection based filter**

The selection-based filter is a method for filtering records rapidly and conveniently. This is done by first marking an entry in a field in the determination overview. After activating the selection-based filter the table is filtered so that only records which contain the marked entry are shown.


 **PVCTM.MEL - MLOGVIEW / Filter / Remove filter**

The filter used on the event overview is removed and all events are shown again.

**Parameter for event recording**

763 PVC THERMOMAT CONTROL / Options / Event log

The Log settings window opens:



**Delete entries older than...**

Number of days before events are automatically deleted.  
Range: 1 ... 365

## 4.4 Calibration functions


### 4.4.1 Determination of cell constants

As in PVC Thermomat measurements only the alterations in the conductivity are measured and evaluated, the cell constants of the conductivity cells used in the PVC Thermomat, which are normally in the range 1.00...1.20, only need to be entered exactly if the absolute value for the conductivity is to be shown correctly. The exact entry is also required before carrying out the GLP test for conductivity. The cell constants can either be entered manually or determined automatically with the aid of a defined standard solution.

#### Manual entry of cell constants

763 PVC THERMOMAT CONTROL / Tools / Cell constants / Manual...


The **Enter/check cell constants** window appears in which the cell constant for the conductivity cell can be entered for each **Channel**.



#### Automatic determination of cell constants

763 PVC THERMOMAT CONTROL / Tools / Cell constants / Automatic...

The **Determine cell constants** window appears:



**Block A**

**Channel 1...4**

Each channel of Block A for which the cell constant of the conductivity cell is to be determined can be marked.

**Block B**

**Channel 1...4**

Each channel of Block B for which the cell constant of the conductivity cell is to be determined can be marked.

**Conductivity of standard solution**


Conductivity of the standard solution in  $\mu\text{S}/\text{cm}$  which is filled into all measuring vessels for which the cell constant of the conductivity cell is to be determined.

A solution of  $c(\text{KCL}) = 1 \text{ mmol}/\text{L}$ , prepared from the 6.2301.060 Conductivity standard (KCl 0.1 mol/L) available as an option, can be used as standard solution. The conductivity of this standard solution is:

Temperature	Conductivity
18 °C	127 $\mu\text{S}/\text{cm}$
19 °C	130 $\mu\text{S}/\text{cm}$
20 °C	133 $\mu\text{S}/\text{cm}$
21 °C	136 $\mu\text{S}/\text{cm}$
22 °C	138 $\mu\text{S}/\text{cm}$
23 °C	141 $\mu\text{S}/\text{cm}$
24 °C	144 $\mu\text{S}/\text{cm}$
25 °C	147 $\mu\text{S}/\text{cm}$

**<Start>**

Starts the automatic determination of the cell constants. When the determination has ended the **Enter/check cell constants** window appears in which the automatically determined cell constants have been entered:



Click on **<OK>** to accept these values or on **<Cancel>** to retain the previous values.

## 4.4.2 Determination of Delta T

The **temperature correction Delta T** describes the variation of the actual sample temperature from the temperature in the heating block and is a parameter of the method (see *section 4.5.2*). It can be automatically determined with the 6.1111.010 Calibrated external temperature sensor (part of the 6.5616.000 GLP test set available from Metrohm as an option). Since the direct temperature measurement in PVC samples is impossible, silicone oil must be used as sample instead of PVC.

### Preparations

*Fig. 5* shows in detail how to install the accessories for Delta T determination. Proceed as follows:

---

#### 1 Prepare reaction vessel cover

- Install gas inlet tube **35** (part of the 6.5616.000 GLP test set, length = 148 mm) on the reaction vessel cover **2** (see *section 2.2.2*).
- Fasten distance piece **37** to the gas inlet tube **35** so that it is located at a distance of approx. 2 cm from the upper end of the reaction vessel.
- Fasten distance piece **37** to the gas inlet tube **35** so that it is located at a distance of approx. 5 cm from the lower end of the reaction vessel.
- Insert temperature sensor **36** from above in opening **31** of the reaction vessel cover **2** and fasten it in the corresponding openings of the distance pieces **37**.

---

#### 2 Prepare reaction vessel

- Fill reaction vessel **1** with **5 g silicone oil** (e.g. FLUKA 85409).
- Place the reaction vessel cover **2** with temperature sensor **36** on the reaction vessel **1**.
- Slide the temperature sensor **36** right down (the sensor must touch the base of the vessel).

---


#### 3 Insert and connect reaction vessel

- Place the prepared reaction vessel **1** in channel 2 or 3 of the heating block.
- Attach one end of the brown connecting tubing **6** at connection nipple **30** on the reaction vessel cover **2**.
- Attach the other end of the brown connecting tubing **6** at the connection **19** on the 763 PVC Thermomat (see *Fig. 3*).

---


#### 4 Modify method parameters

- In the control window select the required **Method** for the heating block for which the temperature correction **Delta T** is to be determined.



<b>1</b>	<b>Reaction vessel (6.1429.040)</b>
<b>2</b>	<b>Reaction vessel cover (6.2753.100)</b>
<b>6</b>	<b>Connecting tubing (6.1805.080)</b> for gas supply to reaction vessel
<b>30</b>	<b>Connection nipple (6.1808.090)</b>
<b>31</b>	<b>Connection</b> Feedthrough for temperature sensor <b>36</b>
<b>33</b>	<b>Sealing ring (6.1454.040)</b>
<b>34</b>	<b>Connection</b> for connection nipple <b>30</b>
<b>35</b>	<b>Gas inlet tube (6.2418.000)</b> length = 148 mm
<b>36</b>	<b>Temperature sensor Pt100 (6.1111.010)</b>
<b>37</b>	<b>Distance piece (6.2042.040)</b>

**Fig. 5:** Accessories for Delta T determination

- Click on  to open the **Method parameters** window. Enter the required set temperature for the sample under **Temperature** and set **Delta T** to **0°C**. Set **Gas flow** to the desired value and click on **<Save>**.

### **5 Start gas flow and heating**


- Switch on gas flow manually.
- Start heating for the selected block manually.

**Carrying out**

763 PVC THERMOMAT CONTROL / Tools / Determine Delta T / Block A


763 PVC THERMOMAT CONTROL / Tools / Determine Delta T / Block B

These menu items are used to start the automatic determination of **Delta T** for Block A or Block B. The **External temperature sensor** window first appears in which the calibration data of the temperature sensor used can be entered (once data have been entered they are retained):



<b>Sensor name</b>	Name of temperature sensor.
<b>Point 1...3</b>	Calibration points 1...3 of the temperature sensor.
<b>T</b>	Temperature at the calibration point in °C.
<b>R</b>	Resistance at the calibration point in ohms.

When the external temperature sensor data have been entered the **Determination of Delta T** window appears.



If the temperature defined in the method has not been reached then the heating is started automatically. The message **Waiting until temperature is reached...** appears in the window.

When the temperature defined in the method has been reached the message **Measuring temperature. Please wait...** appears in the window. The three following values are shown below this message:

<b>Block temperature</b>	The temperature in the heating block measured by the internal temperature sensor.
<b>Temperature ext. sensor</b>	The temperature in the reaction vessel measured by the external temperature sensor.
<b>Delta T</b>	The temperature correction determined from the difference between <b>Block temperature – Temperature ext. sensor</b> .

As soon as both the measured temperatures are stable the message **Measured temperature is now stable!** appears. The determined temperature correction **Delta T** can now be included in the method with **<Accept>**. The message **Delta T is now set in the corresponding method** then appears.

## 4.5 Methods

### 4.5.1 Method handling

#### Create method



763 PVC THERMOMAT CONTROL / File / New method


Creates a new method. A window opens for the input of the method parameters (see *section 4.5.2*). When all parameters have been entered the method must be stored under a new name with <Save as...>.

#### Open method



763 PVC THERMOMAT CONTROL / File / Open method

Opens an existing method. The **Open method** window opens for method selection.



**Method name** Name of the method.

**Creation date** Date and time at method creation.

**Creator** Name of the logged-in user at method creation.


<Open> Open the method parameters window (see *section 4.5.2*).

#### Method manager



763 PVC THERMOMAT CONTROL / File / Method manager

The **Method manager** window opens for renaming, deleting and opening methods.



<b>Method name</b>	Name of the method.
<b>Creation date</b>	Date and time at method creation.
<b>Creator</b>	Name of the logged-in user at method creation.
<b>&lt;Rename&gt;</b>	Rename the selected method.
<b>&lt;Delete&gt;</b>	Delete the selected method.
<b>&lt;Open&gt;</b>	Open the method parameters window (see <i>section 4.5.2</i> ).

### Save method

In the **Method parameters** window (see *section 4.5.2*) the two following buttons are available for storing an opened method:

<b>&lt;Save&gt;</b>	Saves the method under the same name.
<b>&lt;Save as...&gt;</b>	Saves the method under a new name. The old method is retained. Methods can also be copied in this way.

### Print method report

The following button is available in the **Method parameters** window (see *section 4.5.2*) for printing out all the method parameters of the opened method:

<b>&lt;Print&gt;</b>	Prints out method report.
----------------------	---------------------------

## 4.5.2 Method parameters

### Overview



763 PVC THERMOMAT CONTROL / File / Open method

By clicking on this menu item or on  beside the method name in the operating field of the **763 PVC THERMOMAT CONTROL** window the window with the method parameters opens; it contains the following tabs:


<b>Parameters</b>	Measuring parameters.
<b>Evaluation</b>	Parameters for curve evaluation.
<b>Curves</b>	Parameters for scaling the axes.
<b>Formulas</b>	Equation input for further calculation of the induction and stability times.
<b>Standards</b>	Parameters for calculating the standard time.
<b>Documentation</b>	Parameters for report output.
<b>Description</b>	Freely-definable comments about method.

The **Method parameters** window contains the following buttons:

<Save>	Saves method under the same name.
<Save as>	Saves method under a new name; the old method is retained.
<Print>	Prints out method report.
<Cancel>	Closes method parameter window without saving the modified parameters.

## Parameters

Tab **Parameters** of the **Method parameters** window with measuring parameters.



**Temperature**

**Set temperature** of the sample.  
Range: 50...220 °C

**Delta T**

**Temperature correction:** Correction value for increasing the block temperature in order to reach the set temperature in the sample. This value is determined by measuring the variation of the actual sample temperature from the block temperature using a calibrated external temperature sensor (see *section 4.4.2*).  
Range: -9.9...9.9 °C

If no temperature sensor is available for the **Delta T** determination, a value of **+0.6 °C** can be entered for temperatures of **160 °C to 200 °C** (this value was determined using 5 g silicone oil and a gas flow of 7 L/h nitrogen).

---

**Gas flow** **Gas flow** through sample. It cannot be set individually for each block; the flow rate for both blocks depends always on the value entered in the method loaded for Block A.  
Range: **4...15 L/h**

---

**Start options**

**Start delay** Delay period until start of measurement.  
Range: **0...1000 min**

**Start mode per channel** Mode for start and stop of the channels:  
All 8 channels can be started and stopped individually.

**per 2 channels** Each 2 channels can be started and stopped together.

**per block** The 4 channels in one block can be started and stopped together.

---

**Stop criteria**

**Time** Time until automatic termination of the determination.  
Range: **0.1...100 h**

**Conductivity** Limit value for conductivity. If this value is reached then the determination is terminated automatically.  
Range: **1...400 mS/cm**

**Endpoint(s)** Automatic termination of the determination when each active channel has reached all endpoints.

---

**Shut off at end of determinations**

**Heater** Automatically switches off the heating at the end of the determination.


**Gas flow** Automatically switches off the gas flow at the end of the determination.



*If the gas flow is automatically switched off, the **reaction vessels must be removed immediately** from the PVC Thermomat after the end of the determination. If the reaction vessels are not removed, the instrument can be damaged.*

## Evaluation

Tab **Evaluation** of the **Method parameters** window with parameters for curve evaluation.



**Evaluate induction time** Determines the **induction time** (see below).

**Evaluate stability time** Determines the **stability time** (see below).

**Delta Kappa** Conductivity alteration for the determination of the **stability time**.  
Range: **1...400 mS/cm**

**Evaluation delay** (option for the induction time)  
Delay period until start of the curve evaluation.  
Range: **0...100 h**

**Evaluation suppression** (option for the induction time)


**Start** Time from start of the determination up to start of the evaluation suppression.  
Range: **0...100 h**

**End** Time from start of the determination up to end of the evaluation suppression.  
Range: **0...100 h**

### Induction time

The **induction time** is the time up to the inflection point of the *conductivity vs. time* curve recorded by the PVC Thermomat. The induction time is a characteristic value for the thermal stability of the sample under investigation.


For automatic evaluation of the induction time the second derivative of the measured curve is used. This derivative has a maximum at the inflection point and can be displayed in the **Live graph** window. In order to recognize an inflection point certain criteria concerning the height and width of the peak in the second derivative must be fulfilled. The inflection point can also be evaluated manually as point of intersection of the two tangents of the extended straight arms of the curve.



With certain samples it is possible for the conductivity to increase in small steps a long time before the induction time itself. In order that this increase is not evaluated as the endpoint the operator can delay the evaluation or suppress it in a particular range.

### Stability time


The **stability time** is the time required to achieve a given conductivity alteration of **Delta Kappa**.



The input of an evaluation delay or evaluation suppression has no influence on this evaluation mode; the calculation is always carried out from the time that the measurement was started.

### Curves

Tab **Curves** of the **Method parameters** window with parameters for scaling the curve axes.



#### Time axis(x)

- Auto scaling** Automatic fixing of the initial and final values of the time axis.
- Manual scaling** Manual fixing of the final value of the time axis.
- End value** Final value of the time axis.  
Range: **0.5...100 h**


---

#### Conductivity axis(y)

- Auto scaling** Automatic fixing of the initial and final values of the conductivity axis.
- Manual scaling** Manual fixing of the final value of the conductivity axis.
- End value** Final value of the conductivity axis.  
Range: **10...400 mS/cm**

### Formulas

Tab **Formulas** of the **Method parameters** window with parameters for entering equations for further calculations with the induction and stability times.



#### Formula 1 for induction time

If this option is switched on then the automatically determined **Induction time** will be calculated from the values entered in the equation by the operator and then stored as the **Formula 1** result under the name entered in the first field.

#### Formula 2 for stability time


If this option is switched on then the automatically determined **Stability time** will be calculated from the values entered in the equation by the operator and then stored as the **Formula 2** result under the name entered in the first field.

#### Formula 3 for stability time

If this option is switched on then the automatically determined **Stability time** will be calculated from the values entered in the equation by the operator and then stored as the **Formula 3** result under the name entered in the first field.

## Standards

Tab **Standards** of the **Method parameters** window with parameters for calculating the standard time.



**Calculate standard time** Switches the calculation of the standard time according to the following equation on/off:


$$\text{Standard time} = \text{Factor}^{(\text{Temperature} - \text{Target temperature})/10} \cdot \text{Induction time}$$

### Formula parameters

<b>Factor</b>	Parameter for calculating the standard time which gives the factor by which the time changes by 10 °C. This value can be determined by the <b>extrapolation</b> (see <i>section 4.7.8</i> ) of results obtained at different temperatures.
<b>Target temperature</b>	Target temperature for calculating the standard time.

### Documentation

Tab **Documentation** of the **Method parameters** window with parameters for the report output.



#### Include in report


<b>Header</b>	Output of the header with company settings (if defined), printout date, sample identifications, instrument information, determination date, user, cell constant and serial number.
<b>Method parameters</b>	Output of the method parameters.
<b>Curve(s)</b>	Output of curves.
<b>Method description</b>	Output of the method description given on the <b>Description</b> tab.
<b>Results</b>	Output of the selected results: <b>Induction time</b> , <b>Stability time</b> , <b>Standard time</b> , <b>Formula 1...3</b> .

#### Automatic report

<b>Send to printer immediately</b>	Automatic printout of the reports defined under <b>Include in report</b> following the end of the determination.
<b>Send to file</b>	Automatic output of determination and method data following the end of the determination as a text file *.txt for further processing in other programs. The contents of this report does not depend on the report elements defined under <b>Include in report</b> .
<b>&lt;Browse&gt;</b>	Selection of the <b>Directory</b> for storing this file.

## Description

Tab **Description** of the **Method parameters** window with freely-definable **comments about the method**.




**Enter description text**    Comments about the method.

## 4.6 Determinations

### 4.6.1 Flow chart

The following flow chart provides an overview of all steps for carrying out determinations of the oxidative stability.



## 4.6.2 Sample preparation

PVC or other plastics which decompose through the action of heat with cleavage of hydrogen halides are prepared in accordance with DIN 53'381 as follows:

### 1 Powdery samples

- Polymers, dry blends and coating powders need no special sample pretreatment.

### 2 Pastes

- Pastes are fused on a glass plate in a drying oven to a sheet with a thickness of about 0.5 mm. After cooling, the sheet is cut into pieces of edge length about 2 mm.

### 3 Solid samples

- Molding compounds in granular form, molded and finished parts, semi-finished products as well as coatings separated from the substrate are pulverized in accordance with DIN 53'733 and screened to separate the screening fraction between the analytical screens with wire mesh screen DIN 4188-1.4 (wire mesh of stainless steel with 1.4 mm inside diameter mesh size) and wire mesh screen DIN 4188-2 (wire mesh of stainless steel with 2 mm inside diameter mesh size) for the experiments.

### 4 Sheets

- Sheets are cut into pieces with an edge length of about 2 mm.



*According to DIN 53'581, the samples are to be added directly into the prewarmed reaction vessels inserted in the heating block. Nevertheless, comparative measurements have shown identical results for sample addition to cold or hot reaction vessels.*

## 4.6.3 Prepare instruments and accessories

### Clean instruments and accessories

The **cleanness** of instruments and accessories is an absolute necessity for achieving reliable, reproducible and correct analytical results. Even minute contamination can catalytically accelerate the thermal decomposition and lead to completely false results. Make sure that the following checks and cleaning procedures are carried out before the start of a determination:

### 1 Openings for reaction vessels

- Check whether the openings for the reaction vessels are clean and empty. Blow out any dust in the openings with nitrogen.

- When the instrument is not in use always close the openings with stoppers.

## 2 Measuring vessels and accessories

- After the contents have been discarded clean used **measuring vessels 4** several times with alcohol and distilled water (**do not use acetone!**).
- Rinse the **measuring vessel cover 5** complete with electrodes **28** and PTFE tube **27** several times with distilled water and acetone. If necessary, remove the protective ring **29** for better access to the electrodes **28** (see *Fig. 4*).
- Rinse the **connecting tubing 3** between reaction vessel and measuring vessel several times with distilled water and acetone.

## 3 Reaction vessels and accessories

- It is recommended to use **new reaction vessels 1** and **gas inlet tubes 32** for each determination. Blow out new reaction vessels with nitrogen.
- Used and not too strongly contaminated reaction vessels and gas inlet tubes can be cleaned by immersion in boiling RBS solution for 1 h or a similar laboratory flushing agent. They must then be thoroughly rinsed with distilled water and acetone.



*RBS attacks glass. Reaction vessels and gas inlet tubes therefore cannot be cleaned ad infinitum.*

- Remove the gas inlet tube **32** from the **reaction vessel cover 2** and rinse it several times with distilled water and acetone. Dry the cover at 80 °C in a drying oven.



*Reaction vessel covers which have not been dried properly can falsify the results of the following determinations.*

### Prepare measuring vessels

The cleaned measuring vessels must always be prepared as follows before the start of a determination (see *Fig. 4*):

#### 1 Fill measuring vessels

- Fill each cleaned **measuring vessel 4** with **60 mL distilled water**. For analysis times >24 h more than 60 mL distilled water must be added to take the evaporation losses of approx. 7 mL/d into account and to ensure that the electrodes remain adequately immersed.

## 2 Insert measuring vessels

- Place the **measuring vessel cover 5** fitted with a gas inlet tube on the measuring vessel **4** (see *Fig. 4*).
- Place the measuring vessel **4** with its cover in one of the openings provided for it on the 763 PVC Thermomat and connect the electrode plugs **26** on the cover to the corresponding sockets **20** on the 763 PVC Thermomat (see *Fig. 3*).

## Prepare reaction vessels

Before the start of a determination the reaction vessels are filled with the samples and prepared as follows (see *Fig. 4*):

### 1 Weighing in the samples

- Weigh out **0.5 g** of each of the prepared samples into the cold reaction vessels **1**.



*According to DIN 53'581, the samples are to be added directly into the prewarmed reaction vessels inserted in the heating block. Nevertheless, comparative measurements have shown identical results for sample addition to cold or hot reaction vessels.*

## 2 Mount accessories

- Take the upper rim of the reaction vessel **1** in your hand (e.g. in the space between thumb and index finger) and rotate the glass through 360°. This provides the degreased glass with a thin **grease film**, without which it is very difficult to remove the cover from the vessel after the determination.
- Insert an **gas inlet tube 32** (length = 98 mm) into the reaction vessel cover **2** and fasten it by screwing down the connection nipple **30**.
- Place the **reaction vessel cover 2** on the reaction vessel **1**. Rotate the cover so that the gas inlet tube **32** is as close as possible to the vessel wall.
- Fasten the **white connecting tubing 3** to the measuring vessel **4** on to the corresponding connection **31** on the reaction vessel cover **2**.
- Place the prepared reaction vessel **1** into the vessel holder.


### 4.6.4 Carry out determinations

#### Select method

The determination method is selected in the **Method** field in the control window.




Different methods with different temperatures can be selected for Block A and Block B. In contrast, the gas flow which is defined in the method loaded for Block A also applies to Block B.


By clicking on  the method parameters (see *section 4.5.2*) can be altered as required.

#### Switch on heating

The heating of the two blocks A and B to the temperature set in the method can be switched on individually for each block. This is done by clicking on the following button.




As soon as the heating is switched on the color of the frame of this button switches to red. At the same time the temperature display **9** "TEMPERATURE" on the 763 PVC Thermomat starts to blink for the selected block.




If this button is clicked on again during the heating-up phase then the heating will be switched off.

The actual temperature is shown digitally beside the button. A bar below shows the absolute variation of the actual temperature from the target temperature; the bar length corresponds to 50 °C.




When the selected temperature has been reached the temperature display **9** "TEMPERATURE" stops blinking and the LED remains lit up. At the same time the color of the frame of the button switches to **green**. The block is now ready for receiving the samples and starting the determination.



If this button is clicked on the heating is switched off.

### Enter sample identifiers

The sample identifiers **ID1** and **ID2** can be entered in the corresponding fields of the operating area in the control window before or after starting the determination.



The entry for **ID1** can be selected from the automatically formed list of identifiers which have already been entered. This list can be deleted again with **763 PVC THERMOMAT CONTROL / Options / Clear sample ID history**.

### Insert reaction vessels

When the heating block reaches the temperature defined in the method (frame color of the <Heat to...> button turns **green**) the reaction vessels can be inserted into the PVC Thermomat 763 as follows (see *Fig. 4*):

#### 1 Insert reaction vessels







- Place the prepared reaction vessels **1** in the heating block openings on the 763 PVC Thermomat. Close any openings which are not in use with stoppers or empty reaction vessels to prevent contamination.

#### 2 Connect reaction vessels

- Connect the white connecting tubing **3** fastened to the reaction vessel cover **2** to the corresponding connection piece **22** on the measuring vessel cover **5**.
- Screw the brown gas supply tubing **6** onto the reaction vessel cover **2** and the corresponding connection **19** on the 763 PVC Thermomat (see *Fig. 3*).


### Start determinations

Determinations can be started by clicking on the following buttons in the operating area of the **763 PVC THERMOMAT CONTROL** window:


	Starts a single determination for the corresponding channel.
	Starts a double determination for channels 1 and 2.
	Starts a double determination for channels 3 and 4.
	Starts a multiple determination for channels 1 to 4.
	Starts all Block A channels.
	Starts all Block B channels.

### Live curve


The live curve of a started determination is shown in the operating window beside the channel number. Either a status message or the actual measured value appears above the curve. A differentiation is made between the following conditions:




**No active measurement.**  
The channel is ready for the start of a determination.




**No active measurement.**  
Communication error between instrument and PC.



**Determination started.**  
The first value is being measured.



**Running measurement.**  
Measured value and time are shown above the live curve.




**Determination finished.**  
The channel is ready for the start of a new determination.

### Live parameters



763 PVC THERMOMAT CONTROL / View / Live parameters

The **Live parameters** window opens for altering method parameters during a running determination.



The following parameters for Block A or Block B can be altered live:

---

#### Stop criteria

<b>Time</b>	Time up to automatic termination of the determination. Range: <b>0.1...100 h</b>
<b>Conductivity</b>	Conductivity limit value. When this value is reached the determination is automatically terminated. Range: <b>1...400 mS/cm</b>
<b>Endpoint(s)</b>	Automatic termination of the determination when each active channel has reached all endpoints.

---

#### Shut off at end of determinations





<b>Heater</b>	Automatically switches off the heating at the end of the determination.
<b>Gas flow</b>	Automatically switches off the gas flow at the end of the determination.



*If the gas flow is automatically switched off, the **reaction vessels must be removed** from the PVC Thermomat within 1 h after the end of the determination. If the reaction vessels are not removed, a vacuum is generated in the course of the time by condensation of the further produced HCl in condensed water vapor from the measuring vessel and the measuring solution is siphoned back into the hot reaction vessel.*

#### Stop determinations

Determinations are normally terminated automatically when the stop criteria defined in the method have been achieved. Running determinations can also be stopped manually by clicking on the following buttons in the operating area of the **763 PVC THERMOMAT CONTROL** window:

	Stops a single determination for the corresponding channel.
	Stops a double determination for channels 1 and 2.
	Stops a double determination for channels 3 and 4.
	Stops a multiple determination for channels 1 to 4.

If a determination is stopped manually an inquiry **Do you want to save this determination?** appears. With **No** neither measuring data will be saved nor will results be automatically outputted.

## 4.7 Results

### 4.7.1 Determination overview

#### Open determination overview

The determination overview shows selected information in tabular form about all the determinations which are stored in the opened database. To open a determination overview the results window must be opened first and then the database opened.



#### 763 PVC THERMOMAT CONTROL / File / Results

This menu item is used to open the **763 PVC THERMOMAT RESULTS** window. It is used for the display, output and recalculation of **results** of the determinations recorded by the 763 PVC Thermomat. The determination data are stored in the \*.mrd database files and can be shown in this window as an overview table and as curves. The results window can run even when the control window is closed.

When the results window is opened the **Repos.mrd** database is automatically loaded; all the recorded determinations are loaded in this database as standard.



*The program part for display of the results window is organized as database and saved in the **Nachaus.prog** file, which normally has approx. 5 MB. If the size of this file exceeds 10 MB it should be newly optimized and compressed. To do this, close the program and select **Optimize / Program** from the PVC Thermomat program folder in the **Start / Program** menu.*



#### 763 PVC THERMOMAT RESULTS / File / Open database

Opens an existing database. The **Open database** window opens for the selection of the database file \*.mrd. When the database has been selected the corresponding **determination overview** is automatically opened. Only one database can be opened at any time.



*The **Repos.mrd** database has to be optimized and compressed from time to time, especially if data records are deleted frequently. To do this, close the program and select **Optimize / Database** from the PVC Thermomat program folder in the **Start / Program** menu.*



#### 763 PVC THERMOMAT RESULTS / View / Determination overview


Displays an overview table with selected data fields of all the determinations contained in the opened database.

## Format determination overview








### 763 PVC THERMOMAT RESULTS / Format / Select fields

The **Display fields** window opens for selection of the **fields** which are to be shown in the determination overview.




<b>Selection</b>	Selection of all the fields which can be used for the determination overview. The following fields are available:
<b>Block</b>	Number of the block used for the determination.
<b>Cell constant</b>	Cell constant of the conductivity cell used for the determination.
<b>Channel</b>	Number of the channel used for the determination.
<b>Delta Kappa</b>	Alteration in conductivity for the determination of the <b>Stability time</b> .
<b>Determination date</b>	Date and time at which the determination was started.
<b>Formula 1 for induction time</b>	Result for <b>Formula 1</b> .
<b>Formula 2 for stability time</b>	Result for <b>Formula 2</b> .
<b>Formula 3 for stability time</b>	Result for <b>Formula 3</b> .
<b>Gas flow</b>	Gas flow through sample.
<b>ID 1</b>	Sample identifier 1.
<b>ID 2</b>	Sample identifier 2.

<b>Induction time automatic</b>	Automatically determined <b>induction time</b> .
<b>Induction time manual</b>	<b>Induction time</b> manually determined with the re-evaluation.
<b>Method name</b>	Name of method.
<b>Send to file</b>	Folder in which the *.txt determination file is to be stored.
<b>Serial number</b>	Serial number of the instrument used for the determination.
<b>Stability time</b>	Automatically determined <b>stability time</b> .
<b>Standard time</b>	Automatically determined <b>standard time</b> .
<b>Temperature</b>	Set temperature for the method.
<b>Unit</b>	Number of the instrument used for the determination.
<b>User</b>	Name of the user logged in at the start of the determination.
<hr/>	
<b>Display</b>	Fields which have been selected for the determination overview.
<hr/>	
	Moves all fields in <b>Selection</b> to <b>Display</b> .
	Moves the field selected in <b>Selection</b> to <b>Display</b> . The field is inserted at the end of the list.
<hr/>	
	Moves the field selected in <b>Selection</b> to <b>Display</b> . The field is inserted above the field which is marked in <b>Display</b> .
<hr/>	
	Moves the field selected in <b>Display</b> to <b>Selection</b> .
	Moves all fields in <b>Display</b> to <b>Selection</b> .



**763 PVC THERMOMAT RESULTS / Format / Adjust width**

The **Column width** window opens in which the width of the columns selected for the determination overview can be fixed.




<b>Column Width</b>	Column width in characters.
<b>Standard width</b>	Resets the column width to the default value (18.5 characters).
<b>&lt;Best fit&gt;</b>	Adapts the column width to the longest field content (including title).

Alternatively the column width can be set directly by using the mouse. This is carried out by moving the cursor in the title line of the table between two columns until the cursor appears as . With the **left-hand mouse key** the column can then be drawn to the required size. A **double-click** on the selected column automatically adapts the column width to the longest field content (including title).



**763 PVC THERMOMAT RESULTS / Format / Fields...**

The **Cells effects** window opens in which the format of cells shown in the determination overview can be fixed.



<b>Gridlines shown</b>	
<b>Horizontal</b>	Shows horizontal grid lines.
<b>Vertical</b>	Shows vertical grid lines.
<b>Cell effect</b>	
<b>Flat</b>	Shows flat cells.
<b>Raised</b>	Shows raised cells.
<b>Sunken</b>	Shows depressed cells.
<b>Gridline color</b>	Color of horizontal and vertical grid lines.
<b>Background color</b>	Color of cell background.
<b>Sample</b>	Display of the selected format for the cells.

**A** 763 PVC THERMOMAT RESULTS / Format / Fonts...

The **Font** window opens in which the font used for the determination overview can be fixed.



<b>Font</b>	Selects the font from those available in the PC.
<b>Font style</b>	Selects the typestyle for the selected font.
<b>Size</b>	Selects the type size for the selected font in points.
<b>Effects</b>	
<b>Underline</b>	Selects that the characters are to be underlined.
<b>Color</b>	Selects the color for the selected font.
<b>Sample</b>	Display of the selected font.

## Mark determinations

### 763 PVC THERMOMAT RESULTS / Edit / Select record

Marks the record (line) in the determination overview in which the cursor is located.

Alternatively the record can be marked by using the **mouse**; this is done by clicking on the gray field to the left of the required line. If the left-hand mouse key is kept pressed down then several lines in sequence can be marked.

### 763 PVC THERMOMAT RESULTS / Edit / Select all records

Marks all records (lines) in the determination overview.

This function can be triggered by clicking on the first empty field in the title line of the determination overview.

## Find determination



### 763 PVC THERMOMAT RESULTS / Determination / Find

The **Find** window opens in which the search template for the search in the determination overview can be entered.



#### Find What

#### Search template.

The following wildcards are also allowed:

\* any character sequence (e.g. \*ab, \*ab\*, ab\*)

? individual character (e.g. a?a, ??a, ab?)

# individual number (e.g. ##20, ab###, #23##)

#### Search

All

#### Search mode.

Searches through the whole table.

Down

Searches downwards to the end of the table.

Up

Searches upwards to the start of the table.

#### Match

Whole field

#### Search area.

The search is only successful when the whole field agrees with the search template.

Start of field

The search is successful when the start of the field contents agrees with the search template.

Any part of field

The search is successful when part of the field contents agrees with the search template.

<b>Match Case</b>	If this option is switched on then the search will only be successful when the upper and lower case characters of the search template and field contents coincide.
<b>Search Fields As Formatted</b>	If this option is switched on then, in searches for the whole field content in fields with units (e.g. °C, h), the unit will also be taken into account. If only the numerical value is being looked for then this option must be switched off (standard setting).
<b>Search Only Current Field</b>	If this option is switched on then the search will be limited to the selected field. The field name is shown in the window title. This option is not available if one or more complete lines are marked.
<Find First>	Begins the search at the start of the determination overview.
<Find Next>	Continues the search from the location of the cursor in the determination overview.

### Delete determinations

763 PVC THERMOMAT RESULTS / Determination / Delete

Deletes the marked determinations from the database.



*This menu item is reserved for the **Administrator**.*

## Sort and filter determinations




### 763 PVC THERMOMAT RESULTS / Determination / Sort / Ascending


Sorts the records in the determination overview according to the field contents in the selected column (marked column or column in which the cursor is located) in increasing sequence. If several columns are marked then the records will be sorted beginning with the first marked column from the left. The sorting conditions are stored under **Special filter/sorting**.



### 763 PVC THERMOMAT RESULTS / Determination / Sort / Descending

Sorts the records in the determination overview according to the field contents in the selected column (marked column or column in which the cursor is located) in decreasing sequence. If several columns are marked then the records will be sorted beginning with the first marked column from the left. The sorting conditions are stored under **Special filter/sorting**.

The ascending or descending sort can be removed again by clicking on  or **Determination / Filter / Remove filter/sorting**. If the sort is not removed then it is retained until a new sorting process is carried out.


If sorting has been switched off it can be used again by clicking on  or **Determination / Filter / Apply filter/sorting**.


For complex sorting processes with several columns a **Special sorting** must be used.



### 763 PVC THERMOMAT RESULTS / Determination / Filter / Selection based filter

The selection-based filter is a method for filtering records rapidly and conveniently. This is done by first marking a part of an entry or the complete entry in a field in the determination overview. After activating the selection-based filter the table is filtered so that only records which contain the marked entry are shown. The filtering conditions are stored under **Special filter/sorting**.

The filter can be removed again by clicking on  or **Determination / Filter / Remove filter/sorting**. If the filter is not removed then it is retained until a new filtering process is carried out.


If filtering has been switched off it can be used again by clicking on  or **Determination / Filter / Apply filter/sorting**.

For complex filters a **Special filter** must be used.



763 PVC THERMOMAT RESULTS / Determination / Filter / Special filter/sorting

The **Special filter/sorting** window opens in which the complex conditions for filtering and sorting records can be defined.



The lower, important part of the window contains a table for entry of conditions with the following lines:

<b>Field</b>	Selects a field (column) in the determination overview.
<b>Sort</b>	Sorts the selected field: <b>Ascending</b> , <b>Descending</b> or <b>(not sorted)</b> .
<b>Criteria</b>	Criteria for filtering records.
<b>or</b>	Further criteria for filtering records.

In the lines **Criteria** and **or** and in all following lines the number of records to be searched for can be limited by the use of expressions. An **expression** is a combination of operators, functions and field names for laying down a condition for filtering records. Only the most important elements and rules for entering expressions and a few selected examples are given below. For more information please refer to the documentation about Microsoft Access.

**Operators, functions, field names**

<b>Mathematical operators</b>	+, -, =, *, /
<b>Comparison operators</b>	<, >, <>, <=, >=, <b>Between ... And ...</b>
<b>Logical operators</b>	<b>And, Or, Not</b>
<b>Wildcards</b>	*, ?, # (see <b>Find</b> )
<b>Current date</b>	<b>Date()</b>
<b>Current time</b>	<b>Now()</b> (date and time)
<b>Other fields</b>	["field name"]

---

**Linkages with different fields**

If conditions are fixed in different fields then these conditions can also be linked with **And** and **Or**. The following conditions apply:


Criteria in the **same line** are **And-conditions**, i.e. only those records are shown which fulfill all conditions.


Criteria in **different lines** are **Or-conditions**, i.e. all records which fulfill one of the conditions will be shown.


---


**Examples**


>10	Values larger than 10.
*PVC*	All data with "PVC" as part of their name.
##.##.1998	All data for 1998.
<b>Between Date() And (Date() - 7)</b>	all data of previous week.

Newly entered conditions are only saved when they are used with the **Special filter/sorting** window open by clicking on  or **Determination / Filter / Apply filter/sorting**.

By clicking on  or **Edit / Delete filter** all the conditions entered in the **Special filter/sorting** window are deleted again.

By clicking on  or **Edit / Undo** the last alterations made in the **Special filter/sorting** window are canceled.


The filter/sort can be removed again by clicking on  or **Determination / Filter / Remove filter/sorting**. If the filter/sort is not removed then it is retained until a new filtering/sorting process is carried out.

If filtering has been switched off it can be used again by clicking on  or **Determination / Filter / Apply filter/sorting**.



**763 PVC THERMOMAT RESULTS / Determination / Filter / Apply filter/sorting**


The conditions for filtering and sorting records stored under **Special filter/sorting** are applied to the determination overview.

The filter/sorting can be removed again by clicking on  or **Determination / Filter / Remove filter/sorting**. If the filter/sorting is not removed then it is retained until a new filtering/sorting process is carried out.



**763 PVC THERMOMAT RESULTS / Determination / Filter / Remove filter/sorting**

The filter/sorting applied to the determination overview is removed and all records are shown again.

If filter/sorting is switched off it can be switched on again by clicking on  or **Determination / Filter / Apply filter/sorting**.



**763 PVC THERMOMAT RESULTS / Edit / Undo**

The last alteration made in the **Special filter/sorting** window is canceled.



**763 PVC THERMOMAT RESULTS / Edit / Delete filter**

All the conditions entered in the **Special filter/sorting** window are deleted again.



**763 PVC THERMOMAT RESULTS / Determination / Filter / Filter selection**

The lines marked in the determination overview are incorporated in the filter, i.e. only those lines which have previously been marked will be shown.



**763 PVC THERMOMAT RESULTS / Determination / Filter / Selection not in filter**

The lines marked in the determination overview are removed from the filter, i.e. only those lines which were not previously marked will be shown.

## 4.7.2 Determination and method data

### Open view



#### 763 PVC THERMOMAT RESULTS / View / Determination and method data


The **Determination and method data** window is opened which contains all parameters and results for the selected determination on the following tabs:

<b>Determination</b>	Information about the determination.
<b>Results</b>	Overview of results of determination.
<b>Method</b>	Information about the method.
<b>Parameters</b>	Measuring parameters.
<b>Evaluation</b>	Parameters for curve evaluation.
<b>Formulas</b>	Formula input for further calculation of induction and stability times.
<b>Standards</b>	Parameters for calculating the standard time.
<b>Curves</b>	Parameters for scaling the axes.
<b>Documentation</b>	Parameters for report output.
<b>Warnings</b>	Warnings about the determination.

If the determination in the determination overview is changed by marking a different determination then the parameters in the **Determination and method** window will be updated automatically.

### Determination

**Determination** tab of the **Determination and method data** window with information about the determination.




The screenshot shows the 'Determination and method data' window with the 'Determination' tab selected. The window contains the following fields and controls:

- Standards** (selected): ID 1: PVC, ID 2: KA356
- Curves**: Thermomat 1 (4x2 grid of radio buttons, top-left selected)
- Documentation**: Determination date: 24.06.1999, User: Administrator
- Warnings**: Serial number: 01109, Cell constant: 1
- Formulas**: GLP test: GLP test passed.

An 'OK' button is located at the bottom right of the window.

<b>ID 1</b>	Sample identifier 1.
<b>ID 2</b>	Sample identifier 2.
<b>Determination date</b>	Date and time at which the determination was started.
<b>User</b>	Name of the user logged in at the start of the determination.
<b>Serial number</b>	Serial number of the instrument used for the determination.
<b>Cell constant</b>	Cell constant of the conductivity cell used for the determination.
<b>GLP test</b>	GLP status at the time the determination was carried out. The following results are possible:  <b>GLP test not switched on</b> GLP monitoring was not switched on for either temperature or conductivity or gas flow.  <b>Due GLP test</b> The next GLP test was due.  <b>GLP test not passed</b> The last GLP test carried out was failed.  <b>GLP test passed</b> The last GLP test carried out was passed.  Display of instrument, block and channel used for the determination



**Results**

**Results** tab of the **Determination and method data** window with results overview for the determination.

**Induction time automatic** Automatically determined **Induction time**.

**Induction time manual** **Induction time** determined manually by re-evaluation.

**Stability time** Determined **Stability time** for **Delta Kappa**.


**Delta Kappa** Alteration in conductivity for the determination of the **Stability time**.

**Formula 1 for induction time**  
Result for **Formula 1**.

**Formula 2 for induction time**  
Result for **Formula 2**.

**Formula 3 for induction time**  
Result for **Formula 3**.


**Standard time** Result for **standard time**.



Standards	Curves	Documentation	Warnings
Determination	Results	Method	Parameters
Induction time automatic	<input type="text" value="3.26"/>	h	
Induction time manual	<input type="text" value="0.00"/>	h	
Stability time	<input type="text" value="0.00"/>	h	
Delta Kappa	<input type="text" value="0"/>	μS/cm	
Formula 1 for induction time	<input type="text" value="0.00"/>	h	
Formula 2 for stability time	<input type="text" value="0.00"/>	h	
Formula 3 for stability time	<input type="text" value="0.00"/>	h	
Standard time	<input type="text" value="0.00"/>	h	

## Method

**Method** tab of the **Determination and method data** window with information about the method used for the determination.




Standards	Curves	Documentation	Warnings
Determination	Results	Method	Parameters
Method name	<input type="text" value="PVC 200°C"/>		
Creation date	<input type="text" value="16.06.1999"/>	Creator	<input type="text" value="Administrator"/>
Method description	<input type="text"/>		

<b>Method name</b>	Name of the method.
<b>Creation date</b>	Date on which the method was last saved.
<b>Creator</b>	Name of the logged-in used the last time that the method was saved.
<b>Method description</b>	<b>Description</b> of the method (see <i>section 4.5.2</i> ).

### Parameters

Parameters tab of the **Determination and method data** window with the display of the measuring parameters used in the determination.




<b>Temperature</b>	<b>Set temperature</b> of sample.
<b>Delta T</b>	<b>Temperature correction:</b> Correction value for increasing the block temperature in order to reach the set temperature in the sample.
<b>Gas flow</b>	<b>Gas flow</b> through sample.
<b>Start options</b>	
<b>Start delay</b>	Delay time until start of measurement.
<b>Start mode</b>	Mode for starting the channels:
<b>per channel</b>	All channels have been started individually.
<b>per 2 channels</b>	Each 2 channels have been started together.
<b>per block</b>	4 channels in one block have been started together.
<b>Stop criteria</b>	
<b>Time</b>	Time until automatic termination of determination.
<b>Conductivity</b>	Conductivity limit. When this value is reached the determination is automatically terminated.
<b>Endpoint(s)</b>	Automatic termination of the determination when each active channel has reached all endpoints.

**Shut off at end of determinations**

<b>Heater</b>	Heating was switched off automatically at end of determination.
<b>Gas flow</b>	Gas flow was switched off automatically at end of determination.

**Evaluation**


**Evaluation** tab of the **Determination and method data** window with the display of the parameters used for curve evaluation.



<b>Evaluate induction time</b>	Determination of <b>induction time</b> .
<b>Evaluate stability time</b>	Determination of <b>stability time</b> .
<b>Delta Kappa</b>	Alteration in conductivity for the determination of the <b>stability time</b> .
<b>Evaluation delay</b> (Option for induction time)	Delay time before the start of curve evaluation.
<b>Evaluation suppression</b> (Option for induction time)	
<b>Start</b>	Time from start of determination until start of evaluation suppression.
<b>End</b>	Time from start of determination until end of evaluation suppression.

### Formulas

Formulas tab of the **Determination and method data** window with the display of the formulas used for the calculation.



#### Formula 1 for induction time

The automatically determined **Induction time** was calculated from the values entered in the equation by the operator and then stored as the **Formula 1 result** under the name entered in the first field.

#### Formula 2 for induction time


The automatically determined **Stability time** was calculated from the values entered in the equation by the operator and then stored as the **Formula 2 result** under the name entered in the first field.

#### Formula 3 for induction time

The automatically determined **Stability time** was calculated from the values entered in the equation by the operator and then stored as the **Formula 3 result** under the name entered in the first field.

### Standards

**Standards** tab of the **Determination and method data** window with the parameters used for the calculation of the **standard time**.



#### Calculation of standard time


Shows whether the standard time was calculated or not.

#### Formula parameters

- Factor** Parameter for calculating the standard time which gives the factor by which the time changes per 10 °C.
- Target temperature** Target temperature for calculating the standard time.

### Curves

**Curves** tab of the **Determination and method data** window with the display of the parameters used for scaling the axes of the curve.



**Time axis (x)**


<b>Autoscaling</b>	Automatic fixing of the initial and final values of the time axis.
<b>Manual scaling</b>	Manual fixing of the final value of the time axis.
<b>End value</b>	Final value of the time axis.

**Conductivity axis (y)**

<b>Autoscaling</b>	Automatic fixing of the initial and final values of the conductivity axis.
<b>Manual scaling</b>	Manual fixing of the final value of the conductivity axis.
<b>End value</b>	Final value of the conductivity axis.

**Documentation**

**Documentation** tab of the **Determination and method data** window with display of the settings used for the report output.



**Include in report**

<b>Header</b>	Output of the header.
<b>Curve</b>	Output of curves.
<b>Method description</b>	Output of the method description entered on the <b>Description</b> tab.
<b>Results</b>	Output of the selected results: <b>Induction time</b> , <b>Stability time</b> , <b>Standard time</b> , <b>Formula 1...3</b> .
<b>Method parameters</b>	Output of method parameters.

**Automatic report output****Send to printer immediately**

Automatic printout of the reports defined under **Include in report** following the end of the determination.

**Send to file**


Automatic output of determination and method data following the end of the determination as a text file \*.txt.

**Directory**

Directory in which the report will be stored.

**Warnings**

**Warnings** tab of the **Determination and method data** window with the display of warnings about the determination.

**Description**

If the determination was not carried out properly then the following warnings are possible:

**Temperature not reached**

Set temperature was not reached at start of determination.

**Missing measuring points (start time)**

Measuring points are missing from the given time (in s) onwards.

**Missing measuring points (end time)**

Measuring points are again present from the given time (in s) onwards.

**Manual stop**

The determination was stopped manually.

**Abortion after communication error**

Automatic determination stop because connection between the 763 PVC Thermomat and PC was interrupted.

**Value**

Value of the warning (e.g. temperature, time).


### 4.7.3 Graphics and reevaluation

#### Single graph



763 PVC THERMOMAT RESULTS / Determination / Graphics / Single graph

Opens a new window in which the **single graph** for the determination selected in the determination overview is shown. The sample identifiers **ID 1** and **ID 2** are shown in the window title.



If present then the determined **induction time** (red vertical line and numbers) and **stability time** (black vertical line and numbers) will be shown in the graphics window as well as the **measuring curve** (green). A maximum of 16 single graph windows can be displayed.

By clicking the graph using the right mouse key the context-sensitive menu with the following menu items is opened:

- Zoom off**                      Switch off zooming.
- Copy**                              Copy graph to clipboard.
- Autoscaling**                      Switch on/off autoscaling.



*In contrast to the **Multiple graph** and **Live graph** windows, the window with the **single graph** is **not** automatically updated when the determination in the determination overview is changed.*


#### Multiple graph



763 PVC THERMOMAT RESULTS / Determination / Graphics / Multiple graph

Opens a new **Multiple graph** window in which **all curves** for the determinations marked in the determination overview are shown.

The individual measuring curves are shown in color, the legends referring to the colors (sample identifiers **ID 1** and **ID 2**) are located in the graph header. If present the determined **induction time** ( ) and **stability time** ( ) will also be shown.



By clicking the graph using the right mouse key the context-sensitive menu with the following menu items is opened:


<b>Offset</b>	Switch on/off curve offset.
<b>Zoom off</b>	Switch off zooming.
<b>Copy</b>	Copy graph to clipboard.

### Live graph



763 PVC THERMOMAT RESULTS / Determination / Graphics / Live graph with reevaluation

Opens a new **Live graph** window in which the **measuring curve** of the marked determination is shown. If a different determination is marked in the determination overview then the curve is automatically updated.



If present then the determined **induction time** (red vertical line and numbers) and **stability time** (black vertical line and numbers) will be shown in the graphics window as well as the **measuring curve**

(green). With <2<sup>nd</sup> derivative> the second derivative of the measuring curve can be shown as a blue line.

By clicking the graph using the right mouse key the context-sensitive menu with the following menu items is opened:

<b>Reevaluation</b>	Open reevaluation window (see below).
<b>Zoom off</b>	Switch off zooming.
<b>Copy</b>	Copy graph to clipboard.
<b>Autoscaling</b>	Switch on/off autoscaling.

### Switch zoom on/off

**Switch on zoom:** move the cursor to the left-hand upper corner of the area to be enlarged, press the left-hand mouse key and move the cursor to the lower right-hand corner of the required zoom area. When the left-hand mouse key is released the selected area will be enlarged to fill the whole window.

**Switch off zoom:** click on the graph window with the right-hand mouse key and select the menu item **Zoom off** in the context-sensitive menu.

### Switch autoscaling on/off

With the **Autoscaling** option in the context-sensitive menu of the graph window the manual scaling of the axes made under **Curves** can be canceled. If this option is switched on then the initial and final values of the axes will be fixed automatically so that all measuring points are visible.

### Switch offset on/off

The **Offset** option in the context-sensitive menu of the **Multiple graph** window can be used to switch the offset for the displayed measuring curves on and off. If offset is switched on the curves will be shown offset from each other. In this case the difference of the whole axis is shown in  $\mu\text{S/cm}$  instead of the normal conductivity axis marking.

### Copy graph to clipboard

The **Copy** option in the context-sensitive menu of the graph window can be used to copy the displayed graph into the clipboard.

### Print graph




**763 PVC THERMOMAT RESULTS / File / Print**

If one of the **Single graph**, **Multiple graph** or **Live graph** windows is open and selected (active) then this menu item can be used to print out the graph shown in the graph window directly.

For printing several single or multiple graphs at the same time, see *section 4.7.6*.

## Reevaluation

The **Reevaluation** option in the context-sensitive menu of the **Live graph** window opens the **Reevaluation** window for the manual reevaluation of measuring curves.



### Induction time

<b>Automatic</b>	Automatically determined <b>Induction time</b> .
<b>Manual</b>	Induction time determined manually in the reevaluation.

### Stability time

<b>Stability time</b>	Determined <b>Stability time</b> for <b>Delta Kappa</b> .
<b>Delta Kappa</b>	Alteration in conductivity for the determination of the <b>Stability time</b> .

**ID 1** Sample identifier 1.

**ID 2** Sample identifier 2.

**Determination date** Date and time at start of determination.

**Temperature** Sample temperature in the determination.

The alterations in the results achieved by means of the reevaluation can be stored, canceled or reset with the following buttons:

**<Save>** Saves the results altered in the reevaluation.

**<Undo>** Cancels all the alterations made in the reevaluation since the last time the results were saved.


**<Reset>** Deletes all the alterations produced with the reevaluation. The determination will then only contain the results which were automatically determined during data recording.

The following reevaluation possibilities are offered:




### Set induction time manually

If this option is switched on then the **induction time** can be set manually. This is done by pressing the [ **Ctrl** ] key and at the same time moving the red vertical line in the **Live graph** window with **pressed-down left-hand mouse key** to the required location. The new manually set induction time is displayed in the **Manual** field and in the **Live graph** window.



### Set stability time manually

If this option is switched on then the alteration in conductivity **Delta Kappa** for the calculation of the **stability time** can be altered manually. After alterations in the **Delta Kappa** field the newly calculated **Stability time** is shown both in the field with the same name and in the **Live graph** window.




### Set tangents manually

If this option is switched on then tangents can be set manually to determine the **induction time**. The induction time is then defined as being the intersection point of the tangents. The tangents are set as follows:

#### Set first tangent


Press the [ **Ctrl** ] key and use the left-hand mouse key to

click on the required starting point of the first tangent on the measuring curve and then draw the tangent to the required location with the pressed-down left-hand mouse key.



### Set second tangent

Press the [ **Shift** ] key and use the left-hand mouse key to click on the required starting point of the second tangent on the measuring curve and then draw the tangent to the required location with the pressed-down left-hand mouse key.



### Alter tangents

Each of the two tangents has two points indicated by small black squares; these can be used as follows to alter the tangents after they have been created:

- Move the mouse close to such a point until the cursor appears as a crosswire.
- Press the [ **Ctrl** ] key and move tangent to the required location with the left-hand mouse key pressed down.

The new induction time determined with the tangents is displayed both in the **Manual** field and in the **Live graph** window.



### Delete manually set tangents

Deletes the manually set tangents; the newly determined induction time is retained.

### 4.7.4 Extrapolation


#### Extrapolation window

##### 763 PVC THERMOMAT RESULTS / Determination / Extrapolation

Opens the **Extrapolation** window in which the extrapolation values calculated from the results of the determinations marked in the determination overview are displayed numerically and graphically. The extrapolation is based on an exponential equation with the help of the result quantities defined under **<Settings>** and **Target temperature**. Extrapolation can be used to convert results measured at different temperatures to the required target temperature. In addition this method can be used to determine the **Factor** for converting the **induction time** into the **standard time**.



*The time determined for a particular sample by extrapolation to low target temperatures (e.g. 25 °C) is a purely theoretical value. The relationship between this time and, for example, the working life of the sample must always be determined experimentally.*



#### Extrapolation graph

In the **Temperature vs. Time** graph the line  $T = \ln(t)/B - \ln(A)/B$ , calculated with the help of a linear regression, is shown. The measured values of the selected determination are shown in **black**, the regression line in **green** and the time calculated for the **Target temperature** in **red**.

With the **Copy** option in the context-sensitive menu, which is opened by clicking on the graph with the right-hand mouse key, the displayed extrapolation graph can be copied into the clipboard.

### Extrapolation results

The results of the extrapolation calculation are displayed beneath the graph:

<b>Formula</b>	Exponential equation for calculating the extrapolated time for the target temperature from the measuring values.
<b>A</b>	Calculated coefficient <b>A</b> of the exponential equation.
<b>B</b>	Calculated coefficient <b>B</b> of the exponential equation.
<b>r<sup>2</sup></b>	Calculated coefficient of correlation. If this value is below the <b>limit</b> (see below) set by the administrator then a warning appears when the extrapolation window is opened.
<b>Time</b>	Result of the extrapolation for the target temperature in h and years.
<b>Target temperature</b>	Target temperature for the extrapolation.
<b>Number of determinations</b>	Number of determinations used for the extrapolation.
<b>Standard factor</b>	Value determined from the regression line which gives the factor by which the time changes per 10 °C. This value can be entered in the equation for determination of the <b>Standard time</b> as a <b>Factor</b> (see <i>section 4.5.2</i> ).
<Settings>	Opens the <b>Program settings</b> window with the <b>Extrapolation</b> tab for the entry of parameters for extrapolation (see below).
<Update>	Updates the display of graph and results after an alteration in the selection of the determinations in the determination overview.
<Detail>	Shows only the measured values in the graph.

### Print extrapolation results




#### 763 PVC THERMOMAT RESULTS / File / Print

If the **Extrapolation** window is open then this menu item prints out the extrapolation graphs (normal and detail), the extrapolation results and an overview of the determinations used for the extrapolation with the **ID 1**, **ID 2**, **Temperature** and **Results** fields.

### Settings for extrapolation

763 PVC THERMOMAT RESULTS / Options / Settings / Extrapolation

Extrapolation tab of the Program settings window (see section 4.7.7) with settings for extrapolation.




<b>Induction time</b>	Extrapolation with the results for the <b>induction time</b> .
<b>Stability time</b>	Extrapolation with the results for the <b>stability time</b> .
<b>Formula 1..3 ...</b>	Extrapolation with the <b>Formula results 1...3</b> .
<b>Target temperature</b>	Temperature to which the selected result is to be converted with the help of the regression line obtained from the extrapolation.

### Parameters for extrapolation

763 PVC THERMOMAT RESULTS / Options / Administrator / Extrapolation parameters

The Extrapolation parameters window opens.



<b>Correlation limit</b>	Limit for the correlation coefficient $r^2$ determined in the extrapolation. If this value is below the limit a warning will be given.
--------------------------	--



*This menu item is reserved for the Administrator.*

## 4.7.5 Recalculation

### Recalculation window


#### 763 PVC THERMOMAT RESULTS / Determination / Recalculation

The **Recalculation** window opens in which the results of the determinations marked in the determination overview can be recalculated; it contains the following three tabs:

<b>Method</b>	Recalculation of the induction and stability times.
<b>Formulas</b>	Recalculation of formula 1...3.
<b>Standards</b>	Recalculation of standard time.

### Method

**Method** tab of the **Recalculation** window for the recalculation of induction and stability times.



#### Recalculate induction time

Switches recalculation of the **induction time** on/off. If recalculation is switched off the original result is retained.

**Calculation active** If this option is **switched on** then the parameters for the calculation of the induction time can be newly entered.

If this option is **switched off** then no result will be produced (the original result is deleted).

**Evaluation delay** Delay time until start of the curve evaluation.

#### Evaluation suppression start

Time from start of determination until start of evaluation suppression.

#### Evaluation suppression end

Time from start of determination until end of evaluation suppression.

**Recalculate stability time**

Switches recalculation of the **stability time** on/off. If recalculation is switched off the original result is retained.

**Calculation active**

If this option is **switched on** then **Delta Kappa** can be newly entered to calculate the stability time.


If this option is **switched off** then no result will be produced (the original result is deleted).

**Delta Kappa**

Alteration in conductivity for the determination of the stability time in  $\mu\text{S/cm}$ .

**Formulas**

**Formulas** tab of the **Recalculation** window for the recalculation of formula results.



**Recalculate formula 1**

Switches recalculation of **Formula 1** for the **induction time** on/off. If recalculation is switched off the original result is retained.

**Calculation active**


If this option is **switched on** then the parameters for formula 1 can be newly entered.

If this option is **switched off** then no result will be produced (the original result is deleted).

<b>Recalculate formula 2</b>	Switches recalculation of <b>Formula 2</b> for the <b>stability time</b> on/off. If recalculation is switched off the original result is retained.
<b>Calculation active</b>	If this option is <b>switched on</b> then the parameters for formula 2 can be newly entered. If this option is <b>switched off</b> then no result will be produced (the original result is deleted).
<b>Recalculate formula 3</b>	Switches recalculation of <b>Formula 3</b> for the <b>stability time</b> on/off. If recalculation is switched off the original result is retained.
<b>Calculation active</b>	If this option is <b>switched on</b> then the parameters for formula 3 can be newly entered. If this option is <b>switched off</b> then no result will be produced (the original result is deleted).

## Standards

**Standards** tab of the **Recalculation** window for the recalculation of the standard time.



**Recalculation**

Method Formulas Standards

**Recalculate standard time**

Calculation active

Formula parameters

Factor

Target temperature  °C

Formula: Standard time = Factor ^ ((T - Ttarget) / 10) x Induction time

Cancel OK

### Recalculate standard time

Switches the recalculation of the **standard time** from the **induction time** on/off. If recalculation is switched off the original result is retained.

**Calculation active** If this option is **switched on** then the parameters for calculating the standard time can be newly entered. If this option is **switched off** then no result will be produced (the original result is deleted).

**Formula parameters****Factor**

Parameter for calculating the standard time which gives the factor by which the time changes per 10 °C. This value can be determined by the extrapolation of results obtained at different temperatures.

**Target temperature**

Target temperature for calculating the standard time.

## 4.7.6 Data export and printing

### Copy to clipboard



763 PVC THERMOMAT RESULTS / Edit / Copy

Copies the marked lines or columns in the determination overview into the clipboard.

### Export to other database

763 PVC THERMOMAT RESULTS / Determination / Export / Export to database

Exports the marked determinations from the opened database to a different database. The **Export data to...** window opens in which either an existing \*.mrd database file can be selected or a new database can be created by entering a name in the **File name** field.

When a determination is exported all the current GLP tests (see *section 4.8.4*) which apply to the particular determination are also exported. The GLP results always remain stored in the original database.

### Export to Microsoft Word



763 PVC THERMOMAT RESULTS / Determination / Export / Export to MS Word

Copies the lines or columns marked in the determination overview to **Microsoft Word**. The copied text is stored in an RTF file which is automatically opened by Microsoft Word.

### Export to Microsoft Excel



763 PVC THERMOMAT RESULTS / Determination / Export / Export to MS Excel

Copies the lines or columns marked in the determination overview to **Microsoft Excel**. The copied text is stored in an XLS file which is automatically opened by Microsoft Excel.

### Export of measured values



#### 763 PVC THERMOMAT RESULTS / Determination / Export / List of measured values

Exports the measured values of the determinations marked in the determination overview. The **Save to export file** window opens for the selection of the folder and the entry of the **File name** for the export file. After clicking on **<Save>** the list of measurements for each of the marked determinations is stored in a **\*.txt** text file. This text file contains the following values for each determination:

<b>ID 1</b>	Sample identifier 1.
<b>ID 2</b>	Sample identifier 2.
<b>Determination date</b>	Date of determination in format DD.MM.YYYY hh:mm:ss.
<b>Output date</b>	Output date in format DD.MM.YYYY hh:mm:ss.
<b>Time [s]</b>	Time for every measurement in s.
<b>Measured value [mS/cm]</b>	Measured conductivity in $\mu\text{S/cm}$ .

### Export of determination and method data



#### 763 PVC THERMOMAT RESULTS / Determination / Export / Determination and method data


Exports all the determination and method data for the determinations marked in the determination overview. The **Save to export file** window opens for the selection of the directory and entry of the **File name** for the export file. After clicking on **<Save>** the determination and method data for each of the marked determinations is stored in a **\*.txt** text file. The first line of this text file contains the **parameter title**, separated by tabulator characters. The following lines contain the **parameter values**, also separated by tabulator characters.

## Printing of results in the control window



### 763 PVC THERMOMAT CONTROL / File / Print

The **Print results** window opens and provides an overview of all running and finished determinations of all channels of the connected units.



<b>Unit</b>	Instrument number (1...4).
<b>Block</b>	Instrument block (A, B).
<b>Channel</b>	Channel number (1...4).
<b>Method</b>	Method name.
<b>Sample ID1</b>	Sample identifier 1.
<b>Stab. time</b>	Determined stability time.
<b>Induc. time</b>	Determined induction time.

The width of the columns in the **Print results** window can be altered by pulling the right-hand field margin of the column title with the mouse. These settings are permanently stored when the window is closed.


The required determinations must be marked for the **printout of results**. After clicking on <Print> the results are printed out according to the settings defined in the method on the **Documentation** tab.

### Printing of data in the results window



763 PVC THERMOMAT RESULTS / File / Print

If the **Determination data** window is open and selected (active) then the **Printing of...** window opens.



<b>Overview table</b>	Prints out the determination overview in <b>tabular form</b> . Only those columns which are within the maximum printing width (15 cm) will be printed out.
<b>Overview list</b>	Print out all fields for all determinations in the form of a <b>list</b> .
<b>All determination and method data</b>	Prints out all determination and method data for all the marked determinations. A two-page report is printed out for each determination.
<b>Single graph</b>	Prints out the single graphs for all marked determinations.
<b>Multiple graph</b>	Prints out the multiple graphs with all marked determinations.
<b>Determination acc. to method</b>	Prints out the report defined in the method (see <i>section 4.5.2</i> ) for all marked determinations.
<b>Determination acc. to settings</b>	Prints out the report defined under <b>Settings</b> for all marked determinations.
<b>&lt;Settings&gt;</b>	Opens the <b>Program settings</b> window for fixing the report settings (see <i>section 4.7.7</i> ).
<b>&lt;Print&gt;</b>	Prints all selected options for the determinations selected in the determination overview.

## 4.7.7 Options

### Program settings


763 PVC THERMOMAT RESULTS / Options / Settings

The **Program settings** window opens; it contains the following tabs:

<b>Report settings</b>	Settings for report output.
<b>Company settings</b>	Settings for company logo and text.
<b>Results</b>	Settings for results output.
<b>Scaling</b>	Settings for scaling the curve axes.
<b>Extrapolation</b>	Settings for extrapolation.
<b>Live graph</b>	Settings for displaying the live graph.

### Report settings

Tab of the **Program settings** window with settings for the report output.



#### Definition of elements for report output

<b>Header</b>	Output of the header with company settings (if defined), date of printout, sample identifiers, instrument information, date of determination, user, cell constant and serial number.
<b>Results</b>	Output of the results selected on the <b>Results</b> tab.
<b>Curve</b>	Output of curves.
<b>Method parameters</b>	Output of method parameters.
<b>Method description</b>	Output of the method description entered under <b>Description</b> .


**Definition for multiple printing**

**with offset**

If this option is switched on then in a printout of multiple graphs the curves will be shown offset from each other. In this case the difference of the whole axis is shown in  $\mu\text{S}/\text{cm}$  instead of the normal conductivity axis marking.

**Company settings**

Tab of the **Program settings** window with settings for the output of the company logo and company header in the report header.



<Company logo...>

The **Find company logo** window appears in which the required graphics file (\*.wmf or \*.bmp) with the company logo can be selected. After clicking on <Open> the logo is shown in the adjacent window.

<Remove logo>

The included company logo is removed.

**Logo on the...**

**left side**

The company logo is printed out at the left-hand margin of the header.

**right side**


The company logo is printed out at the right-hand margin of the header.

**Company text**

The text entered here will be printed out at the opposite side of the header to the company logo.

### Results

Tab of the **Program settings** window with settings for the output of results.




Output of the following results on the reports

<b>Induction time</b>	Induction time.
<b>Stability time</b>	Stability time.
<b>Formula 1..3 ...</b>	Formula results.
<b>Standard time</b>	Standard time.

### Scaling


Tab of the **Program settings** window with settings for scaling the curve axes.



<b>Acc. to method</b>	Axis scaling according to the settings defined in the method under <b>Curves</b> .
<b>Autoscaling</b>	Automatic fixing of the initial and final values of the curve axes.
<b>Manual scaling</b>	Manual fixing of the initial and final values of the curve axes:
<b>X Start</b>	Starting value of time axis in h.
<b>X End</b>	Finishing value of time in h.
<b>Y Start</b>	Starting value of conductivity axis in $\mu\text{S/cm}$ .
<b>Y End</b>	Finishing value of conductivity axis in $\mu\text{S/cm}$ .

### Extrapolation


Tab of the **Program settings** window with settings for extrapolation.



<b>Induction time</b>	Extrapolation with the results for the <b>induction time</b> .
<b>Stability time</b>	Extrapolation with the results for the <b>stability time</b> .
<b>Formula 1..3 ...</b>	Extrapolation with the <b>Formula results 1...3</b> .
<b>Target temperature</b>	Temperature to which the selected result is to be converted with the help of the regression line obtained from the extrapolation.

## Live graph

Tab of the **Program settings** window with settings for the display of the **Live graph**.



### Retain zoom at determination change

The selected zoom range is retained after a change in determination in the determination overview.

### Live graph always on top

The **Live graph** window is always shown in the foreground.


## Options for the administrator

### 763 PVC THERMOMAT RESULTS / Options / Administrator / Database password

The **Open** window opens in which the required \*.mrd database file can be selected. After clicking on <Open> the **Password** window appears in which the **Password** for opening the database with **Microsoft Access 97** is given.

### 763 PVC THERMOMAT RESULTS / Options / Administrator / Extrapolation parameters

The **Extrapolation parameters** window opens.



### Correlation limit

Limit for the correlation coefficient  $r^2$ . If this value is below the limit a warning will be given.

### Save settings

The current window settings can be saved using the following two menu items:

#### 763 PVC THERMOMAT RESULTS / Options / Save settings on exit

If this option is switched on then the current window settings will be saved automatically each time the results window is closed.

#### 763 PVC THERMOMAT RESULTS / Options / Save settings now

The current window settings are saved immediately.

## 4.7.8 Windows

### Tile windows vertically



#### 763 PVC THERMOMAT RESULTS / Window / Tile vertically

All the subwindows opened in the results window (determination overview, curves, GLP window) are arranged one above the other.

### Tile windows horizontally



#### 763 PVC THERMOMAT RESULTS / Window / Tile horizontally

All the subwindows opened in the results window (determination overview, curves, GLP window) are arranged side by side.

### Cascade windows

#### 763 PVC THERMOMAT RESULTS / Window / Cascade

All the subwindows opened in the results window (determination overview, curves, GLP window) are arranged so that they slightly overlap horizontally and vertically.

### Arrange icons

#### 763 PVC THERMOMAT RESULTS / Window / Arrange icons

Arranges all the minimized subwindows contained in the results window (determination overview, curves, GLP window) so that all are visible side by side.

### Close all windows

#### 763 PVC THERMOMAT RESULTS / Window / Close all

Closes all the subwindows contained in the results window (determination overview, curves, GLP window).

## 4.8 GLP

### 4.8.1 General information about GLP and validation

The requirements of **GLP** (**Good Laboratory Practice**) include a periodic check of analytical measuring instruments with regard to their reproducibility and accuracy using **Standard Operating Procedures, SOP**. Under the title «**Application Bulletin No. 278 – Validation of Metrohm PVC Thermomats**» an example of such a standard operating procedure for the 763 PVC Thermomat is available from Metrohm.

The operating software for the 763 PVC Thermomat contains prepared **GLP tests** for temperature, conductivity and gas flow measurements. The operator can determine whether and which tests are to be carried out; in addition the intervals between the tests and the required accuracy can also be set down. If the GLP function is activated then each result report will contain a remark stating whether the GLP test requirements have been fulfilled. Metrohm offers the **6.5616.000 GLP test set** as an option for carrying out the most important tests.

Further information on the subjects of QA, GLP and validation can also be found in the brochure «**Quality management with Metrohm**», which is available from your local Metrohm agency.


Testing of the electronic and mechanical function groups of Metrohm instruments can and should be performed as part of a regular service by trained personnel of the manufacturing company (see *section 5.1.1*). All Metrohm instruments are equipped with start-up-test routines which check for perfect functioning of the relevant assemblies when the instrument is switched on. If no error message is displayed, it may be assumed the instrument is operating without faults.

The Metrohm company also supplies its instruments with an integrated diagnostic program which, in the case of possible malfunctions or faulty behavior, allows the service technician to check the functioning of certain assemblies and localize the fault.

### 4.8.2 GLP properties

763 PVC THERMOMAT CONTROL / Tools / GLP test / Properties

This menu item is used to open the **GLP properties** window for the selected instrument; this contains the tabs for temperature, conductivity and gas flow.



- Enable GLP monitoring**    Activates automatic GLP monitoring.

---

- GLP alert interval**    Number of days before next GLP test is due.

---

- Tolerance**    States the tolerance values for the particular GLP test. If the result is within the **Set value ± Tolerance** range then the test has been passed. By clicking on **<Set default>** the preset default value is entered in this field.


**Additional alert message text**  
 Additional text which is to be added to the request to carry out the next GLP test.

### 4.8.3 GLP status



763 PVC THERMOMAT CONTROL / Tools / GLP test / Status

Opens the **GLP status** window for the selected instrument.



Test type	Type of GLP test ( <b>Temperature</b> , <b>Conductivity</b> or <b>Gas flow</b> ).
Next test due date	Date on which the next GLP test is to be carried out. If GLP monitoring is not switched on then <b>Disabled</b> will appear in this column.
Due in	Number of days before next GLP test is due.
Last test result	States whether the previous test was passed ( <b>Ok</b> ) or not ( <b>Not Ok!</b> ).

#### 4.8.4 GLP tests

##### GLP temperature test

###### 763 PVC THERMOMAT CONTROL / Tools / GLP test / Temperature

This menu item is used to start the automatic GLP test for **temperature measurement**. This test requires the use of the 6.1111.010 calibrated temperature sensor and the 6.1253.000 GLP insert; these are contained in the 6.5616.000 GLP test set available from Metrohm as an option. The temperature sensor must be connected to the "**Pt100**" connection on the rear panel of the 763 PVC Thermomat.



*Before starting the GLP test, the temperature of the GLP insert must have been stabilized. If the GLP insert is placed into a cold heating block, this condition is fulfilled after the normal heating time of heating block and GLP insert. On the other hand, if the GLP insert is placed in an already heated heating block, wait at least 30 min before starting the GLP test.*


In the GLP test a separate window opens which guides you through the test step-by-step. The next step is reached with **<Next>**, the previous step with **<Back>** and the test is aborted with **<Cancel>**. Proceed as follows to carry out the GLP temperature test:

##### 1 Select instrument and method

- Select unit 1...4 in the control window for which the GLP temperature test has to be carried out.
- Select the determination method whose temperature has to be used for the GLP test in the control window (see section 4.6.4) and alter the temperature if necessary.


##### 2 Start GLP test

- Click on **763 PVC THERMOMAT CONTROL / Tools / GLP test / Temperature** and click on **<Next>** in the GLP test window.



<b>Sensor name</b>	Name of the temperature sensor.
<b>Point 1...3</b>	Calibration points 1...3 of the temp. sensor.
<b>T</b>	Temperature at calibration point in °C.
<b>R</b>	Resistance at calibration point in ohms.


- Enter the calibration data of the temperature sensor used (once data have been entered they are retained) and click on <Next>.



<b>Channels used for test</b>	
<b>Block A</b>	At least one channel must be selected for Block A.
<b>Block B</b>	At least one channel must be selected for Block B.


**Test duration per channel** Measuring time per channel in min.

- Select the channels, for which the temperature test has to be carried out. Enter the desired **Test duration per channel** and click on **<Next>**.



### 3 Insert GLP test set

- Insert the temperature sensor **36** into the opening of the PTFE cylinder **38** from above (see *Fig. 6*).
- Position aluminium cylinder **41** vertically and introduce temperature sensor **36** into the aluminium cylinder from above until it reaches the stop.
- Fasten temperature sensor **36** to aluminium cylinder **41**, by gently tightening hexagonal screw **39** with the help of Allen key **40**.
- Place temperature sensor **36** together with PTFE and aluminium cylinder in the first channel selected for Block A.
- Connect temperature sensor **36** to connection **17** "Pt100" on the rear panel of the 763 PVC Thermomat and click on **<Next>** to start the measurement.



- 
- 36** Temperature sensor Pt100 (6.1111.010)

---

  - 38** PTFE cylinder (4.422.1370, part of 6.1253.000)

---

  - 39** Hexagonal screw

---

  - 40** Allen key (6.2042.040)


---

  - 41** Aluminium cylinder (3.422.1400, part of 6.1253.000)


**Fig. 6:** Accessories for GLP temperature test

#### 4 Carry out GLP temperature test

- If the temperature defined in the method has not been reached then the heating is started automatically. The message **Waiting until temperature is reached...** appears in the **Measure status** window.



- When the temperature defined in the method has been reached the message **Measuring temperature. Please wait...** appears in the window.




**Ext. temperature** The temperature measured by the external temperature sensor.

**Measuring time** The elapsed time since the start of the temperature measurement.

- When the measurement has been carried out you are requested to insert the external temperature sensor together with the GLP insert in the next selected channel and to start the measurement. The change must be carried out quickly.
- Carry out the GLP temperature measurement for each selected channel in sequence.

## 5 Confirm/cancel test results

- When the last temperature measurement has been completed the **TEST RESULTS** for the measured channels are shown.



**Test passed:** The variation of the measured temperature is within the preset tolerance.



**Test failed:** The variation of the measured temperature is larger than the preset tolerance.

- Check the results and click on <Finish> to save the GLP test results or on <Cancel> to abort the GLP test without saving the results.

## GLP conductivity test

### 763 PVC THERMOMAT CONTROL / Tools / GLP test / Conductivity

This menu item is used to start the automatic GLP test for **conductivity measurement**. This test can either be carried out with a standard solution with known conductivity or with the 6.2109.030 Test resistor which is included in the 6.5616.000 GLP test set available from Metrohm as an option. The test resistor must be connected to the electrode connections **20** of the 763 PVC Thermomat.

In the GLP test a separate window opens which guides you through the test step-by-step. The next step is reached with **<Next>**, the previous step with **<Back>** and the test is aborted with **<Cancel>**.



The **cell constants** of the conductivity cells have to be determined before carrying out the GLP conductivity test with standard solutions for the first time (see section 4.4.1).


Proceed as follows to carry out the GLP conductivity test:

#### 1 Select instrument

- Select unit 1...4 in the control window for which the GLP conductivity test has to be carried out.
- If the GLP test is carried out with the test resistor carry on from step **2**, otherwise continue with step **3**.

#### 2 Carry out GLP test with test resistor

- Click on **763 PVC THERMOMAT CONTROL / Tools / GLP test / Conductivity** and click on **<Next>** in the GLP test window.




<b>Use test resistor</b>	Carries out the GLP test with 6.2109.030 Test resistor.
<b>Resistance</b>	Resistance in ohms (10 k $\Omega$ for 6.2109.030).

- Enable the **Use test resistor** option and enter the **Resistance** of the 6.2109.030 Test resistor in ohms. Afterwards, click on **<Next>**.
- For each channel in sequence, connect the 6.2109.030 Test resistor to the electrode connections **20** and click on **<Next>**.
- Continue with step **4**.

### 3 Carry out GLP test with standard solution

- Insert a measuring vessel filled with 60 mL distilled water together with a measuring vessel cover into all channels (see *section 4.6.3*).
- Click on **763 PVC THERMOMAT CONTROL / Tools / GLP test / Conductivity** and click on **<Next>** in the GLP test window.



GLP test for unit with serial no. 01111

**PREPARATION OF TEST**

Use test resistor   
 Use standard solution

Test with resistor

Resistance  Ohms

Test with standard solution

Enter approx. 60 mL of a standard solution into the vessel of each measuring channel.  
Make sure the electrode covers are inserted properly!

Conductivity of standard solution   $\mu\text{S}/\text{cm}$

<b>Use standard solution</b>	Carries out the GLP test with standard solutions.
------------------------------	---

#### Conductivity of standard solution

Conductivity of the standard solution filled into the measuring vessels in  $\mu\text{S}/\text{cm}$ .


A solution of  $c(\text{KCl}) = 1 \text{ mmol}/\text{L}$ , prepared from the 6.2301.060 Conductivity standard (KCl 0.1 mol/L) available as an option, can be used as standard solution. The conductivity of this standard solution is:

Temperature	Conductivity
18 °C	127 µS/cm
19 °C	130 µS/cm
20 °C	133 µS/cm
21 °C	136 µS/cm
22 °C	138 µS/cm
23 °C	141 µS/cm
24 °C	144 µS/cm
25 °C	147 µS/cm

- Click on <Next> in the GLP test window. The conductivity of each channel is measured in sequence.

**4 Confirm/cancel test results**

- When the last measurement has been completed the **TEST RESULTS** (set values and measured values) for the measured channels will be shown.



**Test passed:** The variation of the measured conductivity is within the preset tolerance.



**Test failed:** The variation of the measured conductivity is larger than the preset tolerance.

- Check the results and click on <Finish> to save the GLP test results or on <Cancel> to abort the GLP test without saving the results.

## GLP gas flow test

### 763 PVC THERMOMAT CONTROL / Tools / GLP test / Gas flow

This menu item is used to start the automatic GLP test for **gas flow measurement**. A gas flow meter is required for this test (not available from Metrohm) so that the gas flow can be measured at each channel.

In the GLP test a separate window opens which guides you through the test step-by-step. The next step is reached with <Next>, the previous step with <Back> and the test is aborted with <Cancel>.


Proceed as follows to carry out the GLP gas flow test:

#### 1 Select instrument

- Select unit 1...4 in the control window for which the GLP gas flow test has to be carried out.

#### 2 Start GLP test


- Click on 763 PVC THERMOMAT CONTROL / Tools / GLP test / Gas flow and click on <Next> in the GLP test window.




Gas flow rate

Gas flow in L/h.

- Enter the **Gas flow rate** and click on <Next>.



- Measure the gas flow for each channel in sequence and enter the measured values in the corresponding fields.
- Afterwards, click on <Next>.



**Test passed:** The variation of the measured gas flow is within the preset tolerance.



**Test failed:** The variation of the measured gas flow is larger than the preset tolerance.

- Check the results and click on <Finish> to save the GLP test results or on <Cancel> to abort the GLP test without saving the results.

### 4.8.5 GLP results

#### Print GLP certificate

763 PVC THERMOMAT CONTROL / Tools / GLP test / Print GLP results

This menu item is used to print out the **GLP certificate** for the selected instrument; this contains the printout date, the **Serial number** of the instrument and the results of the previous GLP tests carried out on this instrument. Each of the GLP tests for temperature, conductivity and gas flow contains the following information:

<b>Date</b>	Date of the last GLP test carried out.
<b>Tester</b>	Operator who carried out the GLP test.
<b>Test result</b>	Information about whether the test was passed or failed.



*In order to print out the test results **GLP monitoring** must be enabled (see section 4.8.2).*

#### View GLP results



763 PVC THERMOMAT RESULTS / View / GLP

The **GLP** selection window opens to provide three possible selections:

**<Temperature>** Opens the overview of the results contained in the current database of the **GLP test for temperature measurement**.

Date	Serial number	Channel								Test done / test passed		
		A1	A2	A3	A4	B1	B2	B3	B4			
▶ 04.10.1999	01111	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>

**<Conductivity>** Opens the overview of the results contained in the current database of the **GLP test for conductivity measurement**.

Date	Serial number	Channel								Test passed
		A1	A2	A3	A4	B1	B2	B3	B4	
▶ 04.10.1999	01111	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

<Gas flow>

Opens the overview of the results contained in the current database of the **GLP test for gas flow measurement**.

GLP Gas flow overview										
Date	Serial number	Channel				Test passed				
		A 1	A 2	A 3	A 4	B 1	B 2	B 3	B 4	
▶ 04.10.1999	01111	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

**Date** Date on which the GLP test was carried out.

**Serial number** Serial number of the instrument on which the GLP test was carried out.

**Channel** Whether the GLP test was carried out and whether it was passed or failed is shown for each of the 8 channels **A1...A4** and **B1...B4**.

The **GLP Temperature**, **GLP Conductivity** or **GLP Gas flow** window is opened by a double-click on a single GLP test within the **Date**, **Serial number** or **Channel** fields; the detailed results of the GLP test are listed in this window.




*If you require an overview of all the GLP tests carried out with the connected instruments then the **Repos.mrd** database must be open. If a different database is opened then only those GLP tests will be shown which were exported to this database at the export of the determinations (see section 4.7.6).*

## Print GLP results



### 763 PVC THERMOMAT RESULTS / File / Print

If a GLP window is open then the GLP results can be printed out with this menu item. If the database contains GLP tests carried out on different instruments then the **Select unit** window opens first in which the **Serial number** of the instrument for which the GLP results are to be printed out must be selected.



One of the following options for the printout can then be selected in the **GLP Printing selection** window.



<b>Last certificate</b>	Prints out the last GLP certificate which in addition to the <b>Serial number</b> of the instrument contains the <b>Date</b> , the <b>Tester</b> and the <b>Test result (Test passed or Test not passed)</b> .
<b>Certificate of date</b>	Prints out the GLP certificate for the given <b>Date</b> which in addition to the <b>Serial number</b> of the instrument contains the <b>Date</b> , the <b>Tester</b> and the <b>Test result (Test passed or Test not passed)</b> .
<b>All temperature GLP's</b>	Prints out the detailed results of all GLP tests for temperature measurement.
<b>All conductivity GLP's</b>	Prints out the detailed results of all GLP tests for conductivity measurement.
<b>All gas flow GLP's</b>	Prints out the detailed results of all GLP tests for gas flow measurement.

# 5 Maintenance – Faults

## 5.1 Maintenance and servicing

### 5.1.1 General information

#### Care

The 763 PVC Thermomat requires proper care and attention. Excessive contamination of the instrument could possibly lead to malfunctions and a shorter service life of the inherently rugged mechanical and electronic parts.

Spilled chemicals and solvents should be wiped up immediately. It is especially important to protect the plug connections at the rear of the instrument (particular the mains plug) against contamination.



*Although constructional measures have been designed to virtually eliminate such a situation, should corrosive media penetrate the interior of the instruments the mains plug of the 763 PVC Thermomat must be immediately disconnected to prevent extensive damage to the instrument electronics. Inform Metrohm service if your instrument(s) have been damaged in such a way.*



*The instrument must not be opened by untrained personnel. Please comply with the safety notes in section 1.5.1.*

### 5.1.2 Maintenance by Metrohm service


Maintenance of the 763 PVC Thermomat is best done as part of an annual service performed by specialists from the Metrohm company. If work is frequently performed with caustic and corrosive chemicals, it may be necessary to shorten the interval between servicing.

The Metrohm service department is always willing to offer expert advice on the maintenance and servicing of all Metrohm instruments.

## 5.2 Faults and malfunctions

### 5.2.1 Error messages in the PC program

If any type of malfunction occurs during operation of the 763 PVC Thermomat, this is shown by error messages in the PC program, which appear in an **error window**, e.g.:




Follow the instructions given in the **error window** and close this window with <OK>.

### 5.2.2 Status and error messages on the instrument

The 763 PVC Thermomat has 10 different LEDs to indicate the status of the instrument (see *Fig. 7*). The status and error messages which could occur on the 763 PVC Thermomat are listed in the following tables. The symbols have the following meanings:

- LED switched off
- LED switched on
- ☼ LED flashes




---

**7 Mains pilot lamp**  
lit up when instrument switched on

---

**8 Gas flow display**  
lit up when gas flow switched on

---

**9 Temperature display**  
flashes when heating switched on  
lit up when temperature reached

---

**10 Error display (red)**  
lit up at error

---

**11 Instrument display**  
indicates numbers (1...4) of units connected

**Fig. 7:** Indicator lamps

**Operating status**

<i>Display</i>	<i>Status/Error</i>	<i>Rectification</i>
ON ○	Instrument switched on.	
ON ●	Instrument switched off.	

**Self-test following switch-on**

After the instrument has been switched on an internal hardware self-test is carried out. Each of the 7 test steps only lasts about a second and is shown on instrument display **11** by a particular combination pattern of the 4 display lamps. If there is a fault then the corresponding LEDs will flash.

<i>Display</i>	<i>Status/Error</i>	<i>Rectification</i>
DEVICE ● ● ● ●	LED test. If all LEDs do not light up during this test then an LED is faulty.	Inform Metrohm service.
DEVICE ○ ●●●●	RAM test error.	Switch instrument off and on again. If the error reoccurs please contact Metrohm service.
DEVICE ●●○●●	PROM check sum test error.	Switch instrument off and on again. If the error reoccurs please contact Metrohm service.
DEVICE ●●●○●	BUSY test error.	Switch instrument off and on again. If the error reoccurs please contact Metrohm service.
DEVICE ●●●○	ADC test error.	Switch instrument off and on again. If the error reoccurs please contact Metrohm service.
DEVICE ○●●○	EEPROM test error.	Switch instrument off and on again. If the error reoccurs please contact Metrohm service.
DEVICE ●○●○	Incorrect mains voltage (Variation > 10%).	Check power supply.

**Gas flow display during operation**

<i>Display</i>	<i>Status/Error</i>	<i>Rectification</i>
FLOW ○	Gas flow switched off.	
FLOW ●●	Gas flow control active, set gas flow not yet reached.	
FLOW ●	Gas flow control active, set gas flow reached.	

### Temperature display during operation

The heating status is shown separately for Block A and Block B.

Display	Status/Error	Rectification
A B ○-TEMPERATURE-○	Heating switched off.	
A B ☀-TEMPERATURE-☀	Heating active, set temperature not yet reached.	
A B ●-TEMPERATURE-●	Heating active, set temperature reached.	

### Error displays during operation

Errors are shown separately for Block A and Block B.

Display	Status/Error	Rectification
○—ERROR—○	System working properly.	
☀—ERROR—☀	Operating error. Possible causes: <ul style="list-style-type: none"> <li>• Heating fault</li> <li>• Maximum temperature (220°C) exceeded by 10 °C.</li> </ul>	<ul style="list-style-type: none"> <li>• Switch instrument off and on again. If the error reoccurs please contact Metrohm service.</li> <li>• Switch off instrument and contact Metrohm service.</li> </ul>
●—ERROR—●	Self-test error (see above).	Switch instrument off and on again. If the error reoccurs please contact Metrohm service.

### Instrument displays during operation

The status of each connected instrument is shown separately.

Display	Status/Error	Rectification
DEVICE ○ ○ ○ ○	Instrument must be logged in to PC.	Set up instrument communication (see <i>section 2.4.3</i> ).
DEVICE ☀ ☀ ☀ ☀	Connection to PC interrupted.	Check connection between PC and PVC Thermomat.
DEVICE ● ○ ○ ○	Instrument 1 is logged in and connected to the PC.	
DEVICE ○ ● ○ ○	Instrument 2 is logged in and connected to the PC.	
DEVICE ○ ○ ● ○	Instrument 3 is logged in and connected to the PC.	
DEVICE ○ ○ ○ ●	Instrument 4 is logged in and connected to the PC.	

# 6 Appendix

## 6.1 Technical data

### 6.1.1 General information

<i>No. of instruments which can be connected</i>	1 ... 4
<i>No. of samples per instrument</i>	1 ... 8 (4 per heating block)
<i>Sample quantity</i>	a few grams / a few milliliters

### 6.1.2 Temperature regulation and measurement

<i>Temperature setting range</i>	50 ... 220 °C
<i>Temperature correction (Delta T)</i>	0 ... ± 9.9 °C (can be entered manually or determined automatically with the help of 6.1111.010 Temperature sensor)
<i>Temperature measuring range</i>	0 ... 250 °C
<i>Resolution</i>	0.1 °C
<i>Max. variation of heating block temperature from set value (50 ... 220 °C)</i>	± 0.3 °C
<i>Reproducibility of set temperature</i>	typ. < 0.2 °C
<i>Max. temperature difference between different measuring positions in one block</i>	typ. < 0.3 °C
<i>Temperature variation (when operating temperature has been reached, with inserted and identically filled reaction vessels with an air flow of 20 L/h)</i>	typ. < 0.1 °C
<i>Switch-off temperature</i>	260 ± 11 °C (if this temperature is exceeded the heating is switched off, it switches on again at 210 ± 11 °C)
<i>Heating-up time of instrument from 20 °C</i>	(up to a temperature stability of ± 0.1 °C)
<i>to 120 °C</i>	approx. 45 min
<i>to 220 °C</i>	approx. 60 min
<i>External instrument temperature</i>	< 50 °C (at operating temperature 220 °C)

### 6.1.3 Conductivity measurement

<i>Sensor</i>	6.0913.130 Conductivity cell built into measuring vessel cover
<i>Construction</i>	Conductivity measurement cell with 2 steel electrodes
<i>Measurement principle</i>	Alternating current measurement with 1 kHz frequency and approx. 1.7 V amplitude (peak to peak).
<i>Cell constant</i>	1.00 ... 1.20 (the exact value can be manually entered or automatically determined)
<i>Measurement range</i>	0 ... 400 $\mu\text{S}/\text{cm}$
<i>Resolution</i>	0.1 $\mu\text{S}/\text{cm}$
<i>Maximum error</i>	$\pm 1 \% \pm 0.5 \mu\text{S}/\text{cm}$


### 6.1.4 Gas flow regulation

<i>Gas flow supply</i>	External supply from N <sub>2</sub> cylinder
<i>Adjustable range</i>	7 ... 25 L/h
<i>Max. variation from set range</i>	$\pm 10 \%$


### 6.1.5 GLP test set

<i>External temperature sensor</i>	
<i>Type</i>	Pt100 jacketed element 1.9 mm dia., with 4-lead technology; DIN IEC 751 class B
<i>Calibration accuracy</i>	0 ... 90 °C: $\pm 0.02 \text{ }^\circ\text{C}$ 90 ... 220 °C: $\pm 0.03 \text{ }^\circ\text{C}$
<i>Max. variation of measured temperature from true value (50 ... 220 °C)</i>	$\pm 0.3 \text{ }^\circ\text{C}$
<i>Test resistor</i>	
<i>Resistance</i>	10 k $\Omega$

### 6.1.6 RS232 interface

<i>Connector</i>	Dsub 9-pin (male)
	
<i>Function</i>	TxD and RxD signal for connection with software handshake
<i>Default settings</i>	9600 baud, 8 bit, 1 stop bit, no parity, XON/XOFF

*Pin assignment*



**6.1.7 Mains connection**

<i>Voltage</i>	2.763.0015: 115 V: 100 ... 120 V ± 10 % 2.763.0014: 230 V: 220 ... 240 V ± 10 %
<i>Frequency</i>	50 ... 60 Hz
<i>Power consumption</i>	450 VA
<i>Fuse</i>	5 mm dia., 20 mm length 100 ... 120 V: 2.0 ATH (slow-blow) 220 ... 240 V: 4.0 ATH (slow-blow)

**6.1.8 Safety specifications**

<i>Construction/testing</i>	According to IEC 1010 / EN 61010 / UL 3101-1, protection class 1, degree of protection IP20
<i>Safety directions</i>	The Instructions for Use include information and warnings which must be heeded by the user to assure safe operation of the instrument.

**6.1.9 Electromagnetic compatibility (EMC)**

<i>Emitted interference</i>	Standards met: EN 55011 (class B) EN 55022 (class B) EN 50081-1 EN 50081-2 EN 61000-3-2 (class A) EN 61000-3-3
<i>Immunity to interference</i>	Standards met: EN/IEC 61000-4-2 (class 3) EN/IEC 61000-4-3 (class 3) EN/IEC 61000-4-4 (class 4) EN/IEC 61000-4-5 (class 2/3) EN/IEC 61000-4-6 (class 3) EN/IEC 61000-4-11 (class 3) ENV 50204

**6.1.10 Ambient temperature**

<i>Nominal operating range</i>	+5...+45°C (at 20...80 % atmospheric humidity)
<i>Storage</i>	-20...+70°C
<i>Transport</i>	-40...+70°C

**6.1.11 Housing**

<i>Material of cover</i>	Polyurethane rigid foam (PUR) with fire protection for fire class UL94VO, CFC-free
<i>Material of base</i>	Steel, enameled
<i>Width</i>	405 mm
<i>Height</i>	268 mm (without accessories) 353 mm (with accessories)
<i>Depth</i>	402 mm
<i>Weight</i>	19.4 kg (without accessories) 26.2 kg (with accessories)

## 6.2 Standard equipment











*Subject to changes !  
All dimensions are given in mm.*





The 763 PVC Thermomat is available in two versions:

- **2.763.0014** PVC Thermomat for mains voltage 230 V
- **2.763.0015** PVC Thermomat for mains voltage 115 V






These instruments include the following parts:

Quant.	Order No.	Description
8	6.0913.130	<b>Measuring vessel cover</b> with integrated conductivity measurement cell
8	6.1428.100	<b>Measuring vessel</b> made of polycarbonate
1	6.1429.040	<b>Reaction vessel</b> made of clear glass  Set of 117
2	6.1454.040	<b>Sealing ring for air inlet tube</b> made of nitrile rubber  Set of 6

Quant.	Order No.	Description	
1	6.1801.000	<b>PVC tubing</b> for connection of nitrogen L = 1 m, d1 = 4 mm, d2 = 7 mm	
1	6.1805.030	<b>FEP tubing</b> with two screw nipples with M6 thread, Length = 150 cm	
8	6.1805.080	<b>FEP tubing</b> with two screw nipples with M6 thread, Length = 25 cm	
1	6.1808.020	<b>Connection piece</b> with M6 thread and tubing stub	
16	6.1808.050	<b>Connection piece</b> with M8 thread and tubing stub	
8	6.1808.090	<b>Connection nipple</b> with M6 and M8 threads	
8	6.1816.010	<b>Connecting tubing</b> made of silicone Length = 22 cm	
8	6.1819.080	<b>PTFE tube</b> for air supply to the measuring vessel	

Quant.	Order No.	Description
1	6.2041.190	<b>Reaction vessel holder</b> for holding 8 6.1429.040 Reaction vessels 
1	6.2122.0X0	<b>Mains cable</b> to customer's specifications: <u>Cable socket</u> <u>Cable connector</u> Type IEC 320/C 13                  Type SEV 12 (CH...) ..... 6.2122.020 Type IEC 320/C 13                  Type CEE (7), VII (D...) ..... 6.2122.040 Type CEE (22), V                    Type NEMA 5-15 (USA...) ..... 6.2122.070
1	6.2134.100	<b>Connecting cable</b> Connecting cable 763 PVC Thermomat (RS232) – PC 
1	6.2418.110	<b>Gas inlet tube</b> for 6.1429.040 Reaction vessel Length = 98 mm Set of 117 
24	6.2753.100	<b>Reaction vessel cover</b> for 6.1429.040 Reaction vessel 
1	6.6037.000	<b>Software CD "Metrodata 763 PC Software 1.0"</b>
1	8.763.1003	<b>Instructions for Use (English)</b> for 763 PVC Thermomat
1	8.763.8007	<b>Registration card (German/English)</b> for PC program «Metrodata 763 PVC Thermomat 1.0»

## 6.3 Optional accessories

Order No.	Description
2.145.001X	<b>Softswitch</b> for connection of 4 PVC Thermomats to 1 COM interface 2.145.0014: for 230 V 2.145.0015: for 115 V
6.1428.020	<b>Measuring vessel</b> made of clear glass 
6.1816.010	<b>Silicone tubing</b> for connection of measuring vessels to the 6.2757.000 Exhaust collection tube Length = 22 cm 
6.2059.000	<b>Rotation ring</b> for rotatable storage of a 763 PVC Thermomat. 
6.2125.110	<b>Connecting cable</b> for connection of the 763 PVC Thermomat to a 25-pin COM interface at the PC. 
6.2301.060	<b>Conductivity standard 250 mL</b> $c(\text{KCl}) = 0.1 \text{ mol/L}$
6.2757.000	<b>Exhaust collection tube</b> for leading off waste air from the measuring vessels. With E.400.0010 Stoppers for closing the openings. 
6.5616.000	<b>GLP test set</b> for carrying out GLP tests, incl. the following accessories: <ul style="list-style-type: none"> <li>1 × 6.1111.010 Temperature sensor Pt100</li> <li>1 × 6.1253.000 GLP insert, incl. the following accessories:               <ul style="list-style-type: none"> <li>1 × 3.422.1400 Aluminium cylinder</li> <li>1 × 4.422.1370 PTFE cylinder</li> </ul> </li> <li>2 × 6.2042.040 Distance piece</li> <li>1 × 6.2109.030 Test resistor 10 kΩ</li> <li>1 × 6.2418.000 Gas inlet tube, length = 148 mm, set of 12</li> <li>1 × 6.2621.120 Hexagon key 1.5 mm</li> </ul>

## 6.4 Warranty and conformity

### 6.4.1 Warranty

The warranty on our products is limited to defects that are traceable to material, construction or manufacturing error which occur within 12 months from the day of delivery. In this case, the defects will be rectified in our workshops free of charge. Transport costs are to be paid by the customer. For day and night operation, the warranty is limited to 6 months.

Glass breakage in the case of electrodes or other parts is not covered by the warranty. Checks which are not a result of material or manufacturing faults are also charged during the warranty period. For parts of outside manufacture insofar as these constitute an appreciable part of our instrument, the warranty stipulations of the manufacturer in question apply.

With the regard to the guarantee of accuracy, the technical specifications in the instruction manual are authoritative.

Concerning defects in material, construction or design as well as the absence of guaranteed features, the orderer has no rights or claims except those mentioned above.

If damage of the packaging is evident on receipt of a consignment or if the goods show signs of transport damage after unpacking, the carrier must be informed immediately and a written damage report demanded. Lack of an official damage report releases Metrohm from any liability to pay compensation.

If any instruments and parts have to be returned, the original packaging should be used if at all possible. This applies above all to instruments, electrodes, burette cylinders and PTFE pistons. Before embedment in wood shavings or similar material, the parts must be packed in a dust-proof package (for instruments, use of a plastic bag is imperative). If open assemblies are enclosed in the scope of delivery that are sensitive to electromagnetic voltages (e.g. data interfaces etc.) these must be returned in the associated original protective packaging (e.g. conductive protective bag). (Exception: assemblies with built-in voltage source belong in a non-conductive protective packaging). For damage which arises as a result of non-compliance with these instructions, no warranty responsibility whatsoever will be accepted by Metrohm.

### 6.4.2 Software license

As the purchaser of the «763 PVC Thermomat 1.0» PC software you also receive the license to use this software on your PC. The PC software has been validated by Metrohm in accordance with the requirements of the ISO 9001 quality system concerning its functioning, analytical performance and accuracy of results. The software functions are documented in the 'Instructions for Use' provided.

Please send us your 8.763.8007 Registration card as soon as possible so that we can register you as an official purchaser. As a registered purchaser you will receive any updated program versions for a preferential price.

### 6.4.3 EU Declaration of conformity



## EU Declaration of Conformity

The METROHM AG company, Herisau, Switzerland hereby certifies, that the instrument:

### 763 PVC Thermomat

meets the requirements of EC Directives 89/336/EEC and 73/23/EEC.

#### Source of the specifications:

EN 50081-1	Electromagnetic compatibility, basic specification; Emitted Interference
EN 50082-1	Electromagnetic compatibility, basic specification; Interference Immunity
EN 61010	Safety requirements for electrical laboratory measurement and control equipment

#### Description of the instrument:

Instrument for automatic determination of the thermal stability of PVC samples.

Herisau, September 21, 1999





Dr. J. Frank

Ch. Buchmann

Development Manager

Production and  
Quality Assurance Manager

**6.4.4 Certificate of conformity and system validation**

<b>Certificate of Conformity and System Validation</b>	
<p>This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.</p>	
Name of commodity:	763 PVC Thermomat
System software:	Stored in ROMs
Name of manufacturer:	Metrohm Ltd., Herisau, Switzerland
Principal technical information:	Voltages: 100...120, 220...240 V Frequency: 50...60 Hz
<p>This Metrohm instrument has been built and has undergone final type testing according to the standards:</p> <p style="text-align: center;">                         EN/IEC 61000-4-2 (class 3), EN/IEC 61000-4-3 (class 3),                          EN/IEC 61000-4-4 (class 4), EN/IEC 61000-4-5 (class 2/3),                          EN/IEC 61000-4-6 (class 3), EN/IEC 61000-4-11, ENV50204  <span style="float: right;"><i>— Electromagnetic compatibility</i></span> </p> <p style="text-align: center;">                         IEC1010, EN61010, UL3101-1  <span style="float: right;"><i>— Security specifications</i></span> </p> <p>It has also been certified by the Swiss Electrotechnical Association (SEV), which is member of the International Certification Body (CB/IEC).</p> <p>The technical specifications are documented in the instruction manual.</p> <p>The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance. The features of the system software are documented in the instruction manual.</p>	
<p>Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.</p>	
<p>Herisau, September 21, 1999</p> <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p>Dr. J. Frank Development Manager</p> </div> <div style="text-align: center;">  <p>Ch. Buchmann Production and Quality Assurance Manager</p> </div> </div>	

## 6.5 Index

### 2

<2<sup>nd</sup> derivative> ..... 95

### A

A 100

Acc. to method ..... 113

Accessories

Cleaning ..... 66

Installation ..... 10

Standard equipment ..... 139

Additional alert message text ..... 117

Administrator ..... 39,40,114

All conductivity GLP's ..... 130

All determination and method data ..... 109

All gas flow GLP's ..... 130

All temperature GLP's ..... 130

Allen key **40**

Figure ..... 121

Tightening hexagonal screw **39** ..... 120

Aluminium cylinder **41**

Figure ..... 121

Fitting ..... 120

Ambient temperature ..... 138

And-condition ..... 82

Appendix ..... 135

Application Bulletins ..... 6

Application Notes ..... 6

Apply filter/sorting ..... 83

Arrangement

Control window ..... 33

Results window ..... 35

Automatic report ..... 63

Automatic report output ..... 92

Autoscaling ..... 60,91,93,95,113

### B

B ..... 100

Background color ..... 76

Basic settings ..... 16

<Best fit> ..... 76

Block ..... 45,74,108

Block A ..... 119

Block B ..... 119

Block temperature ..... 53

<Browse> ..... 63

### C

Cable (6.2125.110) ..... 15,142

Cable (6.2134.100) ..... 15,141

Calculate standard time ..... 62

Calculation active ..... 102,103,104

Calculation of standard time ..... 90

Calibration functions ..... 48

Care ..... 131

Caution ..... 7

Cell constant

Automatic determination ..... 48

Manual entry ..... 48

Technical data ..... 136

Cell constant ..... 74,85

Cell effect ..... 76

Cells effects ..... 76

Certificate of conformity and system

validation ..... 145

Certificate of date ..... 130

Channel ..... 45,48,74,108,129

Channels used for test ..... 119

Check ..... 9

Cleaning ..... 66

Clipboard ..... 26,30,95,106

Coefficient of correlation ..... 100

Color ..... 77

Column width ..... 75,76

COM port ..... 38

COM port settings ..... 17

Comment ..... 7

Communication ..... 38

Company logo ..... 111

<Company logo> ..... 111

Company settings ..... 111

Company text ..... 111

Comparison operators ..... 81

Compress

Database ..... 73

Program ..... 73

Conductivity

GLP test ..... 123

Maximum error ..... 136

Measurement ..... 136

Measurement range ..... 136

Resolution ..... 136

Standard solution ..... 124

<Conductivity> ..... 128

Conductivity axis ..... 60,91

Conductivity measurement cell ..... 136

Conductivity of standard solution ..... 49

Conductivity standard

(6.2301.060) ..... 124,142

Conformity ..... 144

Connecting tubing **3**

Cleaning ..... 67

Connection to measuring

vessel cover **5** ..... 22,68,70

Figure ..... 3,12

Fitting ..... 10,11,21

Ordering designation ..... 140,142

Connecting tubing **6**

Connection to reaction

vessel cover **2** ..... 22,70

Figure ..... 3,12,51

Fitting ..... 11,50

Ordering designation ..... 140

Connecting tubing (6.1805.030)

Ordering designation ..... 140

Connection

to the PC ..... 15

Connection **17**

Connecting temperature

sensor **36** ..... 120

Figure ..... 5

Connection **18**

Figure ..... 5

Fitting connecting tubing ..... 10

Connection **19**

Attach connecting tubing **6** ..... 11,50,70

Figure ..... 5

Connection **31**

Attach connecting tubing **3** ..... 10,11,68

Figure ..... 12,51

Insert temperature sensor **36** ..... 50

Connection **34**

Figure ..... 12,51

Inserting gas inlet tube **32** ..... 10

Connection nipple **30**

Attach connecting tubing **6** ..... 11,50

Figure ..... 12,51

Fitting ..... 10

Ordering designation ..... 140

Connection piece **22**

Attach connecting tubing **3** ..... 11,70

Figure ..... 12

Fitting ..... 10,11

Ordering designation ..... 140

Connection piece (6.1808.020) ..... 10,140

Connection plug **26**

Connection ..... 10,68

Figure ..... 12

Context-sensitive menus ..... 37

Control window

Arrangement ..... 33

Close ..... 34

Definition ..... 32

Icons ..... 34

Menus ..... 33

Open ..... 33

Overview ..... 33

Copy

Determination ..... 106

Extrapolation graph ..... 100

Graphics ..... 95

to clipboard ..... 26

Copy ..... 93,94,95

Correlation coefficient ..... 101,114

Correlation limit ..... 101,114

Creation date ..... 54,55,86

Creator ..... 54,55,86

Criteria ..... 81

Current temp. ..... 45

Curve ..... 91,110

Curves ..... 63

Parameter display ..... 90

Parameter entry ..... 60

Database

Compress ..... 73

Delete determinations ..... 79

File type ..... 37

Open ..... 24,73

Optimize ..... 73

Password for opening with Access 97114

Date ..... 46,128,129

Declaration of conformity ..... 144

Definition for multiple printing ..... 111

Definition of elements for report output 110

Degree of protection ..... 8

<Delete> ..... 55

Delete filter ..... 83

Delta Kappa ..... 58,59,74,85,88,96,97,103

Delta T ..... 53,56,87

Delta-T determination ..... 50

Description ..... 64,92

Database

Compress ..... 73

Delete determinations ..... 79

File type ..... 37

Open ..... 24,73

Optimize ..... 73

Password for opening with Access 97114

<Detail> ..... 100

Determination

- Carry out ..... 22,69
- Cell constant ..... 48
- Copy ..... 106
- Definition ..... 32
- Delete ..... 79
- Delta T ..... 50
- Export to Excel ..... 106
- Export to other database ..... 106
- Export to Word ..... 106
- Filter ..... 80
- Find ..... 78
- Flow chart ..... 65
- Information ..... 84
- Mark all ..... 78
- Marking ..... 78
- Operating tutorial ..... 19
- Preparation ..... 20
- Sort ..... 80
- Start ..... 22,70
- Stop ..... 72

**Determination acc. to method** ..... 109

**Determination acc. to settings** ..... 109

Determination and method data

- Curves ..... 90
- Determination ..... 84
- Display ..... 26
- Documentation ..... 91
- Evaluation ..... 88
- Export ..... 27,107
- Formulas ..... 89
- Method ..... 86
- Open view ..... 84
- Operating tutorial ..... 26
- Parameters ..... 87
- Print ..... 27
- Results ..... 85
- Standards ..... 90
- Warnings ..... 92

**Determination date** ..... 74,85,96,107

**Determination of Delta T** ..... 52

Determination overview

- Display ..... 73
- Export ..... 26
- Format ..... 24,74
- Open ..... 73
- Operating tutorial ..... 24
- Print ..... 25

**Diagnosis** ..... 44

Diagnostic program ..... 116

Dimensions ..... 138

**Display** ..... 75

**Display fields** ..... 74

Distance piece **37**

- Figure ..... 51
- Installation ..... 50

Documentation

- Parameter display ..... 91
- Parameter entry ..... 63

**Due in** ..... 118

**E**

Earthing ..... 8

**Effects** ..... 77

Electrical safety ..... 8

Electrode connection **20**

- Connect test resistor ..... 124
- Connecting measuring vessel cover **510,68**
- Figure ..... 5

Electrodes **28**

- Cleaning ..... 67
- Figure ..... 12

Electromagnetic compatibility ..... 137

EMC ..... 137

Emitted interference ..... 137

**Enable GLP monitoring** ..... 117

**End value** ..... 60,91

**Enter description text** ..... 64

**Enter/check cell constants** ..... 48,49

Error display **10**

- Figure ..... 3,132
- Status display ..... 134

Error messages ..... 46,132

Error window ..... 132

EU Declaration of conformity ..... 144

**Evaluate induction time** ..... 58,88

**Evaluate stability time** ..... 58,88

Evaluation

- Parameter display ..... 88
- Parameter entry ..... 58

**Evaluation delay** ..... 58,88,102

**Evaluation suppression** ..... 58,88

**Evaluation suppression end** ..... 58,102

**Evaluation suppression start** ..... 58,102

**Event** ..... 46

Event log file ..... 37

Event recording ..... 47

Event window ..... 46

Events overview ..... 46

Events window

- Icons ..... 47
- Menus ..... 47

Excel ..... 106

Exhaust collection tube (6.2757.000)

- Fitting ..... 11
- Ordering designation ..... 142

**Exit** ..... 31

Export

- Determination and method data 27,107
- Measured values ..... 107
- to Excel ..... 26,106
- to other database ..... 26,106
- to Word ..... 26,106

Expression ..... 81

**Ext. temperature** ..... 122

External temperature ..... 135

**External temperature sensor** ..... 52

Extrapolation ..... 101

- Definition ..... 32
- Graph ..... 99
- Parameters ..... 101,114
- Print results ..... 100
- Results ..... 100
- Settings ..... 101,113
- Window ..... 99

**F**

**Factor** ..... 62,90,100,105

Faults ..... 132

**Field** ..... 81

Field names ..... 81

Fields

- Adjust width ..... 24,75
- Format ..... 24,25
- Select ..... 24
- Selection ..... 74

File types ..... 37

Filter

- Apply ..... 83
- Delete ..... 83
- Remove ..... 47,83
- Selection based filter ..... 80
- Special filter ..... 81

**Filter selection** ..... 83

**Find** ..... 78

<Find First> ..... 79

<Find Next> ..... 79

**Find What** ..... 78

**Font** ..... 77

**Font style** ..... 77

Fonts ..... 24,25,77

Format cells ..... 76

**Formula** ..... 100

**Formula 1** ..... 101,112,113

Formula 1 for induction time ..... 61,74,85,89

**Formula 2** ..... 101,112,113

**Formula 2 for induction time** ..... 85,89

Formula 2 for stability time ..... 61,74

**Formula 3** ..... 101,112,113

**Formula 3 for induction time** ..... 85,89

Formula 3 for stability time ..... 61,74

**Formula parameters** ..... 62,90,105

Formulas

- Entry ..... 61
- Recalculation ..... 103
- View ..... 89

Front ..... 3

Functions ..... 81

Fundamentals of the operation ..... 31

Fuse holder **14**

- Changing fuses ..... 13
- Figure ..... 5

Fuses ..... 13,14,137

**G**

Gas flow

- Automatic stop ..... 57,72
- Display ..... 42
- Entry ..... 57
- GLP test ..... 126
- Manual switching off ..... 41
- Manual switching on ..... 41
- Regulation ..... 136
- Status display ..... 133

**Gas flow** ..... 57,72,74,87,88

<Gas flow> ..... 129

Gas flow control ..... 41

Gas flow display **8**

- Figure ..... 3,132
- Status display ..... 133

**Gas flow rate** ..... 126

Gas inlet tube **32**

- Cleaning ..... 67
- Figure ..... 12
- Fitting ..... 10,21
- Insert ..... 68
- Ordering designation ..... 141

Gas inlet tube **35**

- Figure ..... 51
- Fitting ..... 50

General information on maintenance ..... 131

General precautionary rules ..... 8

Glossary ..... 32

GLP ..... 116

Certificate .....	128
Conductivity test .....	123
Gas flow test .....	126
General information .....	116
Monitoring .....	117
Print results .....	130
Properties .....	117
Selection windows .....	128
Status .....	117
Temperature test .....	118
View detailed results .....	129
View results .....	128
<b>GLP alert interval</b> .....	117
<b>GLP Printing selection</b> .....	130
<b>GLP test</b> .....	85
GLP test set (6.5616.000)	
Conductivity test .....	123
Ordering designation .....	142
Technical data .....	136
Temperature test .....	118
Graphics	
Autoscaling .....	95
Copy to clipboard .....	30,95
Extrapolation .....	99
Live curve .....	29
Live graph .....	94
Multiple graph .....	28,93
Offset .....	95
Operating tutorial .....	28
Print .....	30,95
Single graph .....	28,93
Zoom .....	95
<b>Gridline color</b> .....	76
<b>Gridlines shown</b> .....	76
Group	
Add .....	39
Delete .....	39
Properties .....	39
<b>Group permissions</b> .....	39
<b>H</b>	
Hazard .....	7
<b>Header</b> .....	63,91,110
<b>Heater</b> .....	57,72,88
Heating	
Automatic start .....	41
Automatic stop .....	57,72
Start .....	20
Status display .....	134
Switch on .....	69
Heating-up time .....	135
Help .....	37
Hexagonal screw <b>39</b>	
Figure .....	121
Tightening .....	120
Housing .....	138
<b>I</b>	
Icons	
Control window .....	34
Events window .....	47
Results window .....	36
<b>ID 1</b> .....	21,28,43,45,70,74, 85,93,96,100,107
<b>ID 2</b> .....	21,28,70,74,85,93,96,100,107
Immunity to interference .....	137
<b>Include in report</b> .....	63,91
<b>Indicator lamps</b> .....	132
<b>Induc. time</b> .....	45,108
Induction time	
Definition .....	59
Graphical display .....	93,94
Recalculation .....	102
Set manually .....	97
<b>Induction time</b> .....	28,29,61,96, 97,101,112,113
<b>Induction time automatic</b> .....	75,85
<b>Induction time manual</b> .....	75,85
Information about the	
Instructions for Use .....	6
Installation	
Accessories for external	
gas supply .....	10
Exhaust collection tube .....	11
Measuring vessel .....	10
Operating tutorial .....	19
Procedure .....	9
Reaction vessel .....	10
Instructions for Use (8.763.1003) ...	6,141
Instrument description .....	2
Instrument display <b>11</b>	
Figure .....	3,132
Self-test .....	133
Status display .....	134
Instrument settings .....	38
Introduction .....	1
<b>L</b>	
Last calibration .....	44
<b>Last certificate</b> .....	130
Last maintenance .....	44
<b>Last test result</b> .....	118
License .....	143
Linkages .....	82
Live curve .....	22,71
Live graph	
Copy .....	30
Open .....	29
Print .....	30
Settings .....	114
Window .....	94
<b>Live graph always on top</b> .....	114
Live parameters .....	71
Location .....	9
Logical operators .....	81
<b>Login</b> .....	31,40
Login window .....	16
<b>Logo on the</b> .....	111
<b>M</b>	
Mains cable	
Connection .....	14
Ordering designation .....	141
Mains connection	
Procedure .....	13,14
Safety notes .....	8
Technical data .....	137
Mains connection plug <b>13</b>	
Connection of mains cable .....	14
Figure .....	5
Mains frequency .....	137
Mains pilot lamp <b>7</b>	
Figure .....	3,132
Function .....	14
Status display .....	133
Mains switch <b>12</b>	
Figure .....	5
On/off switching of the instrument ...	14
Mains voltage	
Check .....	13
Error display .....	133
Technical data .....	137
Maintenance .....	131
<b>Manual scaling</b> .....	60,91,113
Marking .....	78
Marking field <b>24</b>	
Figure .....	12
<b>Match</b> .....	78
<b>Match Case</b> .....	79
Mathematical operators .....	81
<b>Measure status</b> .....	121
<b>Measured value</b> .....	107
Measurement table file .....	37
<b>Measuring time</b> .....	122
Measuring vessel <b>4</b>	
Cleaning .....	67
Figure .....	3,12
Filling .....	67
Fitting .....	10
Insert in instrument .....	68
Ordering designation .....	139
Preparation .....	21,67
Measuring vessel (6.1428.020)	
Fitting .....	11
Ordering designation .....	142
Measuring vessel cover <b>5</b>	
Cleaning .....	67
Connecting exhaust	
collection tube .....	11
Figure .....	3,12
Fitting .....	10,21
Ordering designation .....	139
Place on measuring vessel <b>4</b> .....	68
Technical data .....	136
Menu	
Context-sensitive .....	37
Control window .....	33
Events window .....	47
Results window .....	35
Menu bar .....	33,35
Method .....	45
Create .....	54
Creation date .....	54,55
Creator .....	54,55
Curves .....	60
Definition .....	32
Delete .....	55
Description .....	64
Documentation .....	63
Evaluation .....	58
Formulas .....	61
Handling .....	54
Information .....	86
Manager .....	54
Open .....	54,55
Parameters .....	55,56
Print method report .....	55
Recalculation .....	102
Rename .....	55
Save .....	55
Selection .....	20,69
Standards .....	62
<b>Method</b> .....	20,69,108
<b>Method description</b> .....	63,86,91,110
<b>Method manager</b> .....	54
<b>Method name</b> .....	54,55,75,86
<b>Method parameters</b> .....	63,91,110
Metrohm service .....	131
Microsoft Excel .....	106
Microsoft Word .....	106
Model plate <b>15</b>	
Check mains voltage .....	13

Figure ..... 5  
**Mounting 21**  
 Figure ..... 5  
 Fitting exhaust collection tube ..... 11  
**Mouse** ..... 37  
**Multiple graph**  
 Copy ..... 30  
 Open ..... 28  
 Print ..... 30  
 Window ..... 93  
**Multiple graph** ..... 109

**N**

**New login** ..... 40  
**Next test due date** ..... 118  
**Notation** ..... 7  
**Number of determinations** ..... 100

**O**

**Offset** ..... 94,95  
**On-line help** ..... 37  
**<Open>** ..... 54,55  
**Open method** ..... 54  
**Opening 23**  
 Connecting exhaust collection tube ..... 11  
 Figure ..... 12  
**Opening 24**  
 Figure ..... 12  
 Inserting PTFE needle **27** ..... 10  
**Opening the 763 PVC Thermomat** ..... 8  
**Operating elements** ..... 3  
**Operating status** ..... 133  
**Operating tutorial** ..... 19  
**Operation** ..... 31  
**Operators** ..... 81  
**Optimize**  
 Database ..... 73  
 Program ..... 73  
**Optional accessories** ..... 142  
**Options** ..... 110  
**Options for the administrator** ..... 114  
**or** ..... 81  
**Or-condition** ..... 82  
**Output date** ..... 107  
**Overview list** ..... 25,109  
**Overview table** ..... 25,109

**P**

**Packaging** ..... 9  
**Parameters**  
 Entry ..... 56  
 Extrapolation ..... 101,114  
 View ..... 87  
**Password** ..... 31,114  
**Password protection** ..... 39  
**PC connection 16**  
 Connection to the PC ..... 15  
 Figure ..... 5  
 Technical data ..... 136  
**Pictograms** ..... 7  
**Point** ..... 119  
**Point 1...3** ..... 52  
**Power consumption** ..... 137  
**Print**  
 Extrapolation results ..... 100

GLP certificate ..... 128,130  
 GLP results ..... 130  
**Graphics** ..... 95  
**Results** ..... 108  
 Selection window ..... 109  
**<Print>** ..... 55,56,109  
**<Print report>** ..... 44  
**Print results** ..... 108  
**Printing of** ..... 109  
**Program**  
 Close ..... 31  
 Compress ..... 73  
 Installation ..... 15  
 Optimize ..... 73  
 Program icon ..... 16,17  
 Program number ..... 44  
 Settings ..... 38  
 Start ..... 31  
**Program information** ..... 44  
**Program number** ..... 44  
**Prompts** ..... 33  
**Protection class** ..... 8,14  
**Protective earth** ..... 14  
**Protective ring 29**  
 Figure ..... 12  
 Removing ..... 67  
**PTFE cylinder 38**  
 Figure ..... 121  
 Installation ..... 120  
**PTFE tube 27**  
 Cleaning ..... 67  
 Figure ..... 12  
 Fitting ..... 10  
 Ordering designation ..... 140  
**Pump operating hours** ..... 44  
**PVC tubing (6.1801.000)**  
 Ordering designation ..... 140

**R**

**R** ..... 52,119  
**r^2** ..... 100  
**Reaction vessel 1**  
 Cleaning ..... 67  
 Figure ..... 3,12,51  
 Fitting ..... 10,11,50  
 Insert ..... 11,22,50,70  
 Ordering designation ..... 139  
 Precautionary rules ..... 8  
 Preparation ..... 21,68  
**Reaction vessel cover 2**  
 Cleaning ..... 67  
 Figure ..... 3,12,51  
 Fitting ..... 21  
 Fitting gas inlet tube **32** ..... 10  
 Fitting gas inlet tube **35** ..... 50  
 Ordering designation ..... 141  
 Place on reaction vessel **1** ..... 68  
**Reaction vessel holder (6.2041.190)** ..... 141  
**Rear** ..... 4  
**Recalculate formula 1** ..... 103  
**Recalculate formula 2** ..... 104  
**Recalculate formula 3** ..... 104  
**Recalculate induction time** ..... 102  
**Recalculate stability time** ..... 103  
**Recalculate standard time** ..... 104  
**Recalculation**  
 Definition ..... 32  
 Formulas ..... 103  
 Method ..... 102  
 Standards ..... 104  
 Window ..... 102

**Reevaluation**  
 Definition ..... 32  
 Window ..... 96  
**Reevaluation** ..... 95,96  
**Registration card**  
 (8.763.8007) ..... 15,141,143  
**Remove filter/sorting** ..... 83  
**<Remove logo>** ..... 111  
**<Rename>** ..... 55  
**Report**  
 Automatic output ..... 23  
 Settings ..... 110  
**<Reset>** ..... 96  
**Resistance** ..... 124  
**Results** ..... 63  
 Definition ..... 32  
 Extrapolation ..... 100  
 Operating tutorial ..... 24  
 Operation ..... 73  
 Settings ..... 112  
 View ..... 85  
**Results** ..... 91,100,110  
**Results window**  
 Arrangement ..... 35  
 Close ..... 36  
 Icons ..... 36  
 Menus ..... 35  
 Open ..... 35,73  
 Overview ..... 35  
**Retain zoom at determination change** ..... 114  
**Rotation ring (6.2059.000)**  
 Installation ..... 9  
 Ordering designation ..... 142  
**RS232 interface** ..... 136

**S**

**Safety notes** ..... 8  
**Safety specifications** ..... 137  
**Sample** ..... 76,77  
**Sample ID1** ..... 45,108  
**Sample identifier**  
 Delete list ..... 43  
 Entry ..... 21,70  
 Status display ..... 45  
**Sample preparation** ..... 66  
**Sample quantity** ..... 135  
**Samples**  
 Weighing in ..... 68  
**<Save>** ..... 55,56,96  
**<Save as>** ..... 55,56  
**Scaling** ..... 112  
**Sealing ring 33**  
 Figure ..... 12,51  
 Fitting ..... 10  
 Ordering designation ..... 139  
**Search** ..... 78  
**Search area** ..... 78  
**Search Fields As Formatted** ..... 79  
**Search mode** ..... 78  
**Search Only Current Field** ..... 79  
**Search template** ..... 78  
**Select all records** ..... 78  
**Select number of units** ..... 38  
**Select record** ..... 78  
**Select unit** ..... 130  
**Selection** ..... 74  
**Selection not in filter** ..... 83

Selection-based filter .....	47
Self-test .....	133
<b>Send to file</b> .....	63,75
<b>Send to printer immediately</b> .....	63
<b>Sensor name</b> .....	52,119
<b>Serial number</b> .....	44,75,85,129,130
Service .....	116,131
Service Diagnosis .....	44
Servicing .....	131
<b>Set temp.</b> .....	45
Set temperature .....	56,87
Setting up the instrument .....	9
Settings	
Company settings .....	111
Extrapolation .....	101,113
Live graph .....	114
Mains voltage .....	13
Program settings .....	110
Report settings .....	110
Results .....	112
Save .....	43,115
Scaling .....	112
<b>&lt;Settings&gt;</b> .....	100,109
<b>Shut off at end of determinations</b> .....	57,72,88
Silicone oil .....	50
Single graph	
Copy .....	30
Open .....	28
Print .....	30
Window .....	93
<b>Single graph</b> .....	109
<b>Size</b> .....	77
<b>Softswitch</b> .....	38
Softswitch (2.145.001X)	
Configuration .....	38
Ordering designation .....	142
<b>Softswitch in use</b> .....	38
<b>Softswitch port</b> .....	38
Software installation .....	15
Software license .....	143
SOP .....	116
<b>Sort</b> .....	81
Sorting	
Apply .....	83
Ascending .....	80
Descending .....	80
Remove .....	83
Special sorting .....	81
<b>Special filter/sorting</b> .....	81
Spilled chemicals .....	131
<b>Stab. time</b> .....	45,108
Stability time	
Definition .....	59
Graphical display .....	93,94
Recalculation .....	103
Set manually .....	97
<b>Stability time</b> .....	28,29,61,75,85,96,97,101,112,113
Standard equipment .....	139
<b>Standard factor</b> .....	100
Standard operating procedure .....	116
Standard solution .....	49,124
<b>Standard time</b> .....	62,75,85,112
<b>Standard width</b> .....	76
Standards	
Parameter entry .....	62
Recalculation .....	104
View .....	90
<b>&lt;Start&gt;</b> .....	49,70
<b>Start at</b> .....	41
<b>Start delay</b> .....	57,87
<b>Start heater</b> .....	41
<b>Start mode</b> .....	57,87
<b>Start options</b> .....	57,87
Static charges .....	8
<b>Status</b> .....	45
Status bar .....	33,35
Status message .....	132
Status messages .....	45
Status overview .....	45
<b>&lt;Stop&gt;</b> .....	72
<b>Stop criteria</b> .....	57,72,87
Stoppers (E.400.0010) .....	13
Switching the instrument on/off .....	14
<b>Switch-off temperature</b> .....	135
<b>T</b>	
<b>T</b> .....	52,119
Tangents	
Alter .....	98
Delete .....	98
Set manually .....	97
<b>Target temperature</b> .....	62,90,100,101,105,113
Technical data .....	135
Temperature	
GLP test .....	118
Measurement .....	135
Measuring range .....	135
Recording .....	42
Reproducibility .....	135
Resolution .....	135
Set temperature .....	56,87
Status display .....	134
<b>Temperature</b> .....	56,75,87,96,100
<b>&lt;Temperature&gt;</b> .....	128
Temperature correction .....	50,56,87
Temperature display <b>9</b>	
Figure .....	3,132
Function .....	20,69
Functioning .....	22
Status display .....	134
<b>Temperature ext. sensor</b> .....	53
Temperature sensor <b>36</b>	
Calibration data .....	119
Figure .....	51,121
Fitting .....	50,120
Technical data .....	136
Use for GLP test .....	120
<b>Test duration per channel</b> .....	119
Test resistor (6.2109.030)	
Connection .....	124
Technical data .....	136
<b>Test result</b> .....	128
<b>Test type</b> .....	118
<b>Tester</b> .....	128
<b>Time</b> .....	46,100,107
<b>Time axis</b> .....	60,91
Timer .....	41
<b>Tolerance</b> .....	117
Toolbar .....	33,35
Transport .....	9
Transport damage .....	143
<b>Type</b> .....	46
<b>U</b>	
<b>Underline</b> .....	77
<b>Undo</b> .....	83
<b>&lt;Undo&gt;</b> .....	96
<b>Unit</b> .....	45,46,75,108
<b>Unit configuration</b> .....	38
Unit information .....	44
Unit operating hours .....	44
<b>&lt;Update&gt;</b> .....	100
<b>Use test resistor</b> .....	124
User	
Add .....	39
Delete .....	40
<b>User</b> .....	31,46,75,85
User permissions .....	39
<b>Users and groups manager</b> .....	39
<b>V</b>	
Validation .....	116
<b>Value</b> .....	92
<b>W</b>	
Warnings .....	7,46,92
Warranty .....	143
<b>Weight</b> .....	138
<b>&lt;Which days&gt;</b> .....	41
Wildcards .....	78,81
<b>with offset</b> .....	111
Word .....	106
<b>Z</b>	
Zoom .....	37,95
<b>Zoom off</b> .....	93,94,95