


















Manual Operation			Normal state
 	<b>Lift position</b> Raises or lowers the lift to rest position (HOME) or work position (END).	 	<b>Lift control</b> Raises or lowers the lift.
 	<b>Initialize sample changer</b> Initializes sample changer (RESET) or resets sample changer to initial position (ENDSEQ).	 	<b>Shifting the rack</b> Shifts the rack by one position to the left (anticlockwise) or to the right (clockwise).
	<b>Stop</b> Stops all current changer functions.		
Key	Display	Entry range	Function
	<b>LIFT: 1</b> : rest mm rest, work, rinse shift, special, 0...125 mm		<b>Lift control</b> Raises or lowers the lift to the specified position.
	<b>MOVE 1</b> : sample sample, spec. 1...8, 1...999		<b>Moving a vessel</b> Shifts the rack to place the specified vessel at the needle.
	<b>SAMPLE:</b> = 1 1...999		<b>Definition of sample position</b> Defines the rack position of the first sample vessel of a series.
	<b>PUMP+</b> <b>PUMP-</b>		<b>Pump control</b> Switches on/off the peristaltic pump at the tower.
	<b>Print:</b> config config, param, usermeth, all		<b>Print a report</b> Prints out reports of the IC Sample Processor.
	<b>SCN: Rm</b> : 00000000 Interface Signal/Data Rm 8 x 1, 0 or * RS 14 ASCII characters		<b>Interface scan</b> Displays the state of the remote control lines (Input 0...7) or the received data of the serial interface (RS232).
	<b>CTL: Rm</b> : INIT Interface Signal/Data Rm 14 x 1, 0 or * RS 14 ASCII characters		<b>Interface control</b> Sends the specified signal via remote control lines (Output 0...13) in order to control (e.g. start) peripheral devices or transmits data via serial interface (RS232). <SELECT>-options see sample changer commands.

Display	
<b>Normal state</b> method name ® pump status ®	sample counter ~ <div style="border: 1px solid black; padding: 2px; display: inline-block;">                         ***** counter 1/127                          PUMP- ready                     </div> → sample changer status
<b>Series operating</b> sequence running ®	<div style="border: 1px solid black; padding: 2px; display: inline-block;">                         ***** counter 2/127                          START 03 WAIT 11 s                     </div> → parameter - current command and sequence line

## Basic Configuration




General settings  
Rack configuration  
RS232 interface settings

**Configuration menu**

Display	Default value	Entry range	Function
<b>&gt;auxiliaries</b>			
<b>dialog:</b>	<b>english</b>	english, deutsch, français, español	<b>General settings</b> Selection of dialog language
<b>display contrast</b>	<b>3</b>	0..7	Adjust display contrast
<b>beeper:</b>	<b>ON</b>	ON, OFF	Switch on/off acoustic warning signal
<b>device label</b>	<b>*****</b>	8 ASCII char.	Entry of device label
<b>program</b>	<b>5. 766. 0010</b>	read only	Program version
<b>max. lift way</b>	<b>125 mm</b>	0..125 mm	Maximum lift way; defines lowest lift position possible. <CLEAR> enters actual lift position. RESET required to effect.
<b>&gt;rack definitions</b>			
<b>rack number</b>	<b>2</b>	1..16	Number of rack
<b>code</b>	<b>010001</b>	000001..111111	Binary ID code for rack identification
<b>type:</b>	<b>M 129-2</b>	M129-2, ...	Predefined rack type
<b>work position</b>	<b>125 mm</b>	0..125 mm	Working position for needle. <CLEAR> enters actual position.
<b>rinse position</b>	<b>125 mm</b>	0..125 mm	Rinsing position for needle. <CLEAR> enters actual position.
<b>shift position</b>	<b>0 mm</b>	0 mm	Shift position for needle.
<b>special position</b>	<b>0 mm</b>	0..125 mm	Special position for needle. <CLEAR> enters actual position
<b>&gt;rack definitions 2</b>			
<b>&gt;&gt;special positions</b>			
<b>special beaker 1</b>	<b>128</b>	0..max. pos.	Rack position of special beaker 1
<b>special beaker 2</b>	<b>129</b>	0..max. pos.	Rack position of special beaker 2
<b>... up to special beaker 8</b>	<b>8</b>	...	...
<b>&gt;RS232 settings</b>			
<b>baud rate:</b>	<b>9600</b>	300, 600, 1200, 2400, 4800, 9600	<b>RS232 interface settings</b> Data transmission rate
<b>data bit:</b>	<b>8</b>	7, 8	Data bits
<b>stop bit:</b>	<b>1</b>	1, 2	Stop bits
<b>parity:</b>	<b>none</b>	none, odd, even	Parity
<b>handshake:</b>	<b>HW</b>	HWs, HWfull, SWChar, SWline, none	Handshake
<b>character set:</b>	<b>IBM</b>	IBM, Epson, Seiko, Citizen, HP	Character set for printer or PC
<b>RS control:</b>	<b>ON</b>	ON, OFF	Switch on/off data receiving via serial interface (RS232)

## Operating Sequences and Method Parameters



Number of samples  
Start sequence  
Sample sequence  
Final sequence  
Sample changer settings  
Manual stop options

**Parameter menu**  
**Method**

Display	Default value	Entry range	Function
<b>number of samples:</b>	<b>rack</b>	rack, *, 1..999	Number of samples in a series
<b>&gt;start sequence</b>			
<b>Start sequence of a series</b>			
<b>&gt;sample sequence</b>			
<b>Sample sequence of a series</b>			
<b>&gt;final sequence</b>			
<b>Final sequence of a series</b>			

# Operating Sequences and Method Parameters

continued

Display	Default value	Entry range	Function
<b>&gt;changer settings</b>			<b>Sample changer settings</b>
<b>rack number</b>	<b>0</b>	0...16	Number of rack for current method
<b>lift rate 1</b>	<b>12 mm/s</b>	3...12 mm/s	Lift rate at tower
<b>shift rate</b>	<b>20</b>	3...20	Rack shift rate (in degree/sec)
<b>&gt;manual stop</b>			<b>Manual stop options</b>
<b>CTL Rmt:</b>	*****	14 bit (1, 0 or *)	Remote output signal on manual stop command
<b>CTL RS232:</b>		&PR;\$S,14 characters	Data output via serial interface (RS232) on manual stop command



## Method Memory

	Recall a method Store a method Delete a method	<b>User method menu</b>
--	--	-------------------------



Display	Default value	Entry range	Function
<b>&gt;recall method</b>			<b>Recalling a method</b>
<b>method:</b>	<b>XXXXXXXX</b>	8 ASCII characters, *****	Selection of a method to recall (<SELECT> or enter text) ***** = blank method (<CLEAR>)
<b>&gt;store method</b>			<b>Storing a method</b>
<b>method:</b>	<b>XXXXXXXX</b>	8 ASCII characters	Name of a method to store (<SELECT> or enter text)
<b>XXXXXXXX overwrite ?</b>		<ENTER>, <QUIT>	Overwriting confirmation with <ENTER>, cancel with <QUIT>
<b>&gt;delete method</b>			<b>Deleting a method</b>
<b>method:</b>	<b>XXXXXXXX</b>	8 ASCII characters	Selection of a method to delete (<SELECT> or enter text)
<b>delete XXXXXXXX ?</b>		<ENTER>, <QUIT>	Deletion confirmation with <ENTER>, cancel with <QUIT>

## Series Operating


### Preparation




Key	Display	Entry range	Function
	<b>SAMPLE:</b>	= <b>XXX</b>  1...999	Definition of the first sample beaker of a series. (optional, if undefined SAMPLE = 1) Definitions in sequences are dominant e.g.: SAMPLE = rack (optional; if undefined in sequence, manual setting is valid)
	<b>parameters</b> <b>number of samples: rack</b>	rack, *, 1...999	Number of samples (will be saved with the method).

### Series operating / Key functions

 <ul style="list-style-type: none"> <li>starts a method</li> <li>restarts a method from 'hold' state</li> </ul>	 <ul style="list-style-type: none"> <li>stops a method (final sequence will not be executed)</li> </ul>
---	--

**Series operating**






















 <ul style="list-style-type: none"> <li>holds a running method</li> <li>may be continued with &lt;START&gt;</li> </ul>	 <ul style="list-style-type: none"> <li>quits a running command</li> <li>quits error messages</li> </ul>	 <ul style="list-style-type: none"> <li>stops a series after the current sequence (final sequence will not be executed)</li> </ul>
--	---	---

# Sample Changer Commands

Command	Default value	1 <sup>st</sup> Param.	2 <sup>nd</sup> Parameter	Function
<b>SAMPLE:</b>	= 1	<i>Function</i> =, +, -	<i>Value</i> 1...999	<b>Definition of sample position</b>
<b>MOVE 1 :</b>	<b>sample</b>	<i>Tower</i> 1	<i>Vessel</i> sample, spec. 1...8, 1...999	<b>Shift rack to place the specified vessel at the needle</b>
<b>LIFT: 1 :</b>	<b>rest mm</b>	<i>Tower</i> 1	<i>Position</i> work, rinse, shift, special, rest, 0...125 mm	<b>Setting a lift position</b>
<b>PUMP 1.1 :</b>	<b>1 s</b>	<i>Pump</i> 1.1	<i>Duration/status</i> 1...999 s, ON, OFF	<b>Pump control</b> (1 <sup>st</sup> param.: Tower.Pump)
<b>SCN: Rm :</b>	<b>Ready1</b>	<i>Interface</i> Rm  RS	<i>Input signal/data</i> Ready1, End1, End2, Wait1, Wait2, Wait*, Pump1 ?, Pump2 ?, Pump* ?, 8 bit (1, 0 or *)  14 ASCII characters (*= any char. or string)	<b>Interface scan</b> Scans the remote interface until the specified signal is received.  Scans the serial interface (RS232) until the specified data are received.
<b>CTL: Rm</b>	<b>INIT</b>	<i>Interface</i> Rm  RS	<i>Output signal/data</i> INIT, INIT 732, PROG R/S 1, PROG R/S 2, PUMP R/S 1, FILL A 1, INJECT A 1, FILL B/STEP 1, INJECT B 1, ZERO 1, PUMP 752 ON, PUMP 752 OFF, STEP MSM 753, 14 bit (1, 0 or *)  clear value: &D.S"9" 14 ASCII characters	<b>Interface control</b> Sets the specified remote output signal.  May be defined as binary pattern.  Transmits the specified data via serial interface (RS232). Use text edit mode.
<b>WAIT</b>	<b>1 s</b>	<i>Time</i> 0...9999 s		<b>Waiting time in series operation</b>
<b>ENDSEQ</b>				<b>End of a sequence</b>

## Editing

## Key functions

 	<b>Next / previous menu item</b> Displays next or previous menu item.	 	<b>First / last menu item</b> Displays first or last menu item.
 	<b>Cursor shift</b> Shifts the cursor in between the parameters of a menu line.	 	<b>Default value</b> Deletes an entry and sets its default value.
	<b>Data entry</b> Accepts data entry; opens sub-menus. Data modifications will only be accepted when confirmed with <ENTER>.	 	<b>Roll-up selector</b> Displays each entry of a given list of data to select from. Roll-up selections are marked with a colon (:).
	<b>Quit</b> Quits data entry; branches dialog to a higher menu level or normal state resp.	 	<b>Insert / delete command line</b> Inserts an empty sequence command line (NOP) or deletes the current command line.
	<b>Trace function</b> Executes the actual sequence command line and displays next line.	 	<b>Text edit mode</b> Enters text edit mode. Places text cursor at start (<) or end (>) of the text string.
 	<b>LEARN mode</b> Starts or ends the interactive LEARN mode. It allows fine tuning of the parameters of most sample changer commands in operating sequences.	<b>Adaptive sample changer commands:</b> LIFT, PUMP, SCAN, WAIT	