





















Manual Operation			Normal state
 	Lift position Raises or lowers the lift of the active tower to rest position (HOME) or work position (END).	 	Lift control Raises or lowers the lift.
 	Initialize sample processor Initializes sample processor and dosing drives (RESET) or resets the rack to initial position (RACK).	 	Shifting the rack Shifts the rack by one position to the left (anticlockwise) or to the right (clockwise).
	Changing the display Changes the display mode.		Stop Stops all current device functions. See also [PARAM] > manual stop options
Key	Display	Entry range	Function
	LIFT : 1 : rest mm rest, work, rinse, shift, special, 0...100 mm		Lift control Raises or lowers the lift to the specified position.
	MOVE 1 : sample sample, spec. 1...8, 1...999		Moving a beaker Shifts the rack to place the specified vial to the tower.
	SAMPLE: = 1 1...999		Definition of the first sample position of a series Defines the rack position of the first sample beaker of a series .
	FLOW pump : ON pump, valve ON, OFF		Flow control Switches on/off the pump or the magnetic valve.
	HEATER: INIT°C min 50...250°C, INIT, OFF 0...999		Heater control Sets the oven temperature and the heating interval to reach the specified temperature. If no heating interval is set, the maximum heating rate is applied. Cooling is possible as well.
	DOS: 1 : 1 ml Dosing unit 1...12, * Volume/Function ±0.001...999.999 ml, fill, release, prepar., empty,...		Dosimat / Dosino control Controls a connected dosing drive. Enter a volume to be dispensed or select a function to execute (only 'fill' with 685 Dosimat). With the Dosino 700, also negative volumes can be dispensed i.e. the solution can be aspirated.
	PRINT: full full, short, config, param, usermeth, all		Print a report Prints out a report.
	SCN: Rm : 00000000 Interface Signal/Data Rm 14 x 1,0 or * RS 14 ASCII-characters		Interface scan Displays the state of the remote control lines (Input0...7) or the received data of the serial interface (RS232).
	CTL: Rm : INIT Interface Signal/Data Rm 14 x 1,0 or * RS 14 ASCII-characters		Interface control Sends the specified signal via remote control lines (Output0...13) in order to control (e.g. start) peripheral devices or transmits data via serial interface (RS232). <SELECT>-options see sample processor commands.




















Manual Operation

Continued

		Definition of various settings Defines basic settings of dosing units, lifts or rack . Roll-up selection; press <DEF>-key several times.		
Display	Default value	1. Parameter	2. Parameter	Function
FILLRATE	1 160	Address 1...12	Rate 0.01...160 ml/min	Filling rate
LIFTRATE	1 12 mm/s	Tower 1,2	Rate 3...12 mm/s	Lift rate
SHIFTRATE : auto.	20	Direction auto., +, -	Rate 3...20 w/s	Shift direction and shift rate
DRIVE. PORT 1.1 : dos.		Port 1.1...12.4	Function dos., fill, rinse, prep., drain	Function assignment of Dosino ports
DOSRATE	1 160	Address 1...12	Rate 0.01...160 ml/min	Dosing rate

Editing

Key functions

 	Next / previous menu item Displays next or previous menu item.	 	First / last menu item Displays first or last menu item.
 	Cursor shift Shifts the cursor in between the parameters of a menu line.	 	Default value Deletes an entry and sets its default value.
	Data entry Accepts data entry; opens sub-menus. Data modifications will only be accepted when confirmed with <ENTER>.	 	Roll-up selector Displays each entry of a given list of data to select from. Roll-up selections are marked with a colon (:).
	Quit Quits data entry; branches dialog to a higher menu level or normal state resp.	 	Insert / delete command line Inserts an empty sequence command line (NOP) or deletes the current command line.
	Trace function Executes the actual sequence command line and displays next line.	 	Text edit mode Enters text edit mode. Places text cursor at start (<) or end (>) of the text string.
 	LEARN mode Starts or ends the interactive LEARN mode. It allows fine tuning of the parameters of most sample changer commands in operating sequences.	Adaptive sample changer commands: LIFT, DOS, SCAN, WAIT	

Basic Configuration




General settings
Heater settings
Rack configuration
Dosing unit settings
RS232 interface settings

Configuration menu

Display	Default value	Entry range	Function
>auxiliaries			General settings
dialog:	english	e/d/f/s	Selection of the dialog language
display contrast	3	0..7	Adjust display contrast
beeper:	ON	ON,OFF	Switch on/off acoustic warning signal
device label	*****	8 ASCII-chars	Entry of device label
program	5. 774. 0010	read only	Program version
mx. lift way	90 mm	0..100 mm	Maximum lift way. Defines lowest lift position possible. <CLEAR> enters actual lift position. RESET required to effect.
beaker sensor:	ON	ON,OFF	Switch on/off beaker sensor. RESET required to effect.
>oven settings			Heater and temperature settings
initial temp.:	OFF°C	50..250, OFF	Temperature to be set after switching on
mx. temperature	275°C	50..275	Maximum temperature limit for oven emergency stop..
temp. correction	0°C	-20...0...20	Correction of the temperature displayed
>rack definitions			Rack configuration
rack number	1	1..16	Number of rack
code	000001	000001..111111	Binary ID code for rack identification
type:	M2-0	M12-0 ...	Predefined rack type
work position	0 mm	0..100 mm	Working position of the lift. <CLEAR> enters actual position.
-	
shift position	0 mm	0..20 mm	Shift position of the lift. <CLEAR> enters actual position.
>rack definitions 1			Submenu for special beaker definitions
>>special positions			
special beaker 1	0	0..max.pos.	Rack position of special beaker 1
... up to special beaker 8			...
>dosing units			Dosing unit settings
see user manual			
>RS232- settings			RS 232 interface settings
baud rate:	9600	300,600,1200,2400, 4800,9600	Data transmission rate
data bit:	8	7, 8	Data bits
stop bit:	1	1, 2	Stop bits
parity:	none	none, odd, even	Parity
handshake:	HWs	HWs, HWfull, SWchar, SWline, none	Handshake
character set:	IBM	IBM, Epson, Seiko,Citizen, HP	Character set for printer or PC
RS control:	ON	ON, OFF	Switch on/off data receiving via serial Interface (RS 232)

Operating Sequences and Method Parameters



Number of samples
Start sequence
Sample sequence
Final sequence
Report definitions

Sample changer settings
Timeout settings
Gas flow settings
Dosing unit definitions
Manual stop options

**Parameter menu
Method**

<i>Display</i>	<i>Default value</i>	<i>Entry range</i>	<i>Function</i>
number of samples: rack		rack, *, 1..999	Number of samples in a series
>start sequence			Start sequence of a series
>sample sequence			Sample sequence of a series
>final sequence			Final sequence of a series
>changer settings			Sample changer settings
rack number 0		0..32	Number of rack for current method
lift rate 1 25 mm/s		3..25 mm/s	Lift rate tower 1
lift rate 2 25 mm/s		3..25 mm/s	Lift rate tower 2
shift rate 20		3..20	Rack shift rate (in degree/sec)
shift direction: auto.		auto., +, -	Shift direction of the rack
on beaker error: MOVE		MOVE, display	Sample changer reaction on missing beaker
>timeout settings			Timeout settings
HEATER timeout: off min		0..999, off	Delay time, if the temperature is not reached within heating interval
on HEAT. timeout: STOP		STOP, cont.	Reaction after HEATER timeout
SCAN timeout: off min		0..999, off	max. time interval allowed for SCAN command
on SCAN timeout: error		error, cont.	Reaction after SCAN timeout
>gas flow			Gas flow settings
unit gas flow nL/min		mL/min, L/h	Unit for the display of gas flow
min. gas flow 0.0 nL/min		0..999	Lower limit of gas flow
max. gas flow 900 nL/min		0..999	Upper limit of gas flow
gas type: air		air, others, N2	Type of gas
gas flow factor: 1		0.001..9.999	Correction factor for the flow measuring of 'other' gases
>dosing unit def.			Dosing unit definitions
<i>see user manual</i>			
>manual stop			Manual stop options
CTL Rmt: *****		STOP device1, 14 bit (1,0 or *)	Remote output signal on manual stop command
CTL RS232:		&M;\$S,14 characters	Data output via serial interface (RS232) on manual stop command
FLOW cont.		cont., off	Gas flow on manual <STOP>
HEATER init°C		init, off, 50..250	Temperature to be set on manual <STOP>

LCD Display

Normal state

sample counter ~
 method name ® ***** counter 1/36
 oven temperature ® 110°C 55 mL/min ready → sample processor status
 - gas flow

Series operating

sequence running ® ***** counter 2/36
 START 03 WAIT 11 s → parameter
 - current command and sequence line

Measured values

gas flow 55 mL/min
 oven temperature 110°C



<SELECT/DISPLAY> switches the display of the measured values on or off.

Series Operating

Preparation

Key	Display	Entry range	Function
	SAMPLE: = XXX	1...999	Definition of the first sample vial of a series. (optional, if undefined SAMPLE=1) Definitions in sequences are dominant e.g. : 01 SAMPLE = rack (optional; if undefined in sequence, manual setting is valid)
	> parameters number of samples : rack	rack, *, 1...999	Number of samples (will be saved with the method).

Series Operating / Key functions


 <ul style="list-style-type: none"> starts a method restarts a method from 'hold' state 	<ul style="list-style-type: none"> stops a method (final sequence will not be executed) 	
<p>Series operating</p> <pre> graph LR Start("<START>") -- start sequence --> SampleSequence("sample sequence") SampleSequence -- final sequence --> NormalState("normal state") SampleSequence -- "number of samples" --> SampleSequence </pre>		
 <ul style="list-style-type: none"> holds a running method may be continued with <START> 	 <ul style="list-style-type: none"> quits a running command quits error messages 	 <ul style="list-style-type: none"> stops a series after the current sequence (final sequence will not be executed)

Sample Processor Commands

Command	Default value	1. Param.	2. Parameter	Function
SAMPLE: = 1		<i>Function</i> =, +, -	<i>Value</i> 1...999	Definition of the first sample vial of a series
MOVE 1 : sample		<i>Tower</i> 1,2	<i>Beaker</i> sample, spec.1...8, 1...999	Shift the rack to place a vial to the tower
LIFT: 1 : rest mm		<i>Tower</i> 1,2, *	<i>Position</i> work, rinse, shift, special, rest, 0...325 mm	Setting a lift position
FLOW pump: on		<i>Gas flow pump, valve</i>	<i>State</i> on,off	Gas flow control <i>Switches on/off the pump or valve</i>
HEATER: init °C min		<i>Temp.</i> 30...250°C, init, off	<i>Heating interval</i> 0,1...999, blank	Heater control The temperature will be reached within the heating interval set. The temperature is remained. Max. heating rate is 15°C/min.
DOS: 1 : 1 ml		<i>Address</i> 1...12,*	<i>Function/Vol ume</i> fill, release, prepar., empty, ejct, adjust, level ±0.001...999.999 ml	Dosimat / Dosino control Dosing drives: 700 Dosinos or 685 Dosimats
SCN: Rm : ready1		<i>Interface</i> Rm, RS	<i>Input signal /data</i> Ready1, End1, Cond ok, Cond 737, no error, 8 bit (1,0 or *) 'clear' value: *R", 14 ASCII-chars (*= any character or string)	Interface scan Scans the 'remote' interface until the specified signal is received. Scans the serial interface (RS 232) until the specified data are received.
CTL: Rm INIT		<i>Interface</i> Rm, RS	<i>Output signal /data</i> START device1, START dos1, STOP device1, START 737, ENTER, INIT, 14 x 1,0 or * clear value: &M;\$G 14 ASCII-chars	Interface control Sets the specified remote output signal. May be defined as binary pattern. Transmits the specified data via serial interface (RS232). Use text edit mode.
WAIT 1 s		<i>Time</i> 0...9999 s		Waiting time in series operation
RACK				Rack initialization

More commands (DEF commands) see user manual.

Method Memory

Display	Default value	Entry range	Function
			User method menu
Recall a method Store a method Delete a method			
>recall method			Recalling a method
method: XXXXXXXX		8 ASCII-characters, *****	Selection of a method to recall (<SELECT> or enter text) ***** = blank method (<CLEAR>)
>store method			Storing a method
method: XXXXXXXX		8 ASCII-characters	Name of a method to store (<SELECT> or enter text)
XXXXXXXX overwrite ?		<ENTER>, <QUIT>	Overwriting confirmation <ENTER>, cancel with <QUIT>
>delete method			Deleting a method
method: XXXXXXXX		8 ASCII-characters	Selection of a method to delete (<SELECT> or enter text)
delete XXXXXXXX ?		<ENTER>, <QUIT>	Deletion confirmation <ENTER>, cancel with <QUIT>