





Mode selection	
	Press key <MODE> until the desired mode is displayed and confirm it with <ENTER>.
	Load method from internal method memory:  <RECALL> <X> <ENTER> X=0,1...9
<b>Modes:</b>	
<b>DOS</b>	Dosing: You may calculate a result with the dosed volume.
<b>DIS R</b>	Repetitive dispensing: The burette cylinder is refilled after each increment.
<b>DIS C</b>	Cumulative Dispensing.
<b>PIP</b>	Pipetting.
<b>DIL</b>	Diluting.
<b>CNT D</b>	Preparation of solutions with preselected content.

Mode parameters			
		For all modes: ↑ Expelling rate ↓ Filling rate OFF: Rate adjustable via potentiometer	0.001...150 mL/min, OFF 0.001...150 mL/min, OFF
		Input of volumes, different for each mode, see below	0.001...999.999 mL, OFF
<i>Display</i>	<i>Initial value</i>	<i>Meaning</i>	<i>Input range</i>
<b>Mode DOS</b>			
<b>V-LIM</b>	<b>OFF ml</b>	Security volume; stop if V-LIM is reached	0.001...999.999 mL, OFF
<b>blank</b>	<b>0 ml</b>	Blank value, key <BLANK>	0..±999.999 mL
<b>factor</b>	<b>1</b>	Factor, key <FACTOR>	0..±1E33
<b>snpl</b>	<b>1</b>	Sample size, key <SMPL>	0..±1E33
<b>unit</b>		Unit for the calculated result, key <UNIT>	ppm, %, g, mg, g/L, mg/L, mol, mol/L, mL, L, /pc, none
<b>Mode DIS R</b>			
<b>V-DIS</b>	<b>1 ml</b>	Dispensing volume	0.001...999.999 mL
<b>Mode DIS C</b>			
<b>V-DIS</b>	<b>1 ml</b>	Dispensing volume	0.001...999.999 mL
<b>V-LIM</b>	<b>OFF ml</b>	Security volume; stop if V-LIM is reached	0.001...999.999 mL, OFF
<b>Mode PIP</b>			
<b>V-PIP</b>	<b>0.1 ml</b>	Pipetting volume	0.001...49.5 mL
<b>Mode DIL</b>			
<b>V-PIP</b>	<b>0.1 ml</b>	Pipetting volume	0.001...49.5 mL
<b>V-DIL</b>	<b>1 ml</b>	Diluting volume	0.001...999.999 mL, OFF

## Sequence of mode CNT D

Sequence	Meaning	Input range
<b>CNT D 0.000 ml</b> - - <b>unit %</b> - <ENTER> - <b>cnt 1 %</b> - Input <ENTER> - - - <b>s 1 g</b> - Input <ENTER> - <b>add V X.XXX ml</b> - <GO>	<p>Mode "CNT D": Content.</p> <p>Select unit with key &lt;UNIT&gt;.</p> <p>Input of the desired content of the solution.</p> <p>Depending on the selected unit, auxiliary variables are necessary: Molar mass M, density of the solvent dens. and factor for the volume contraction f.</p> <p>Input of the weight.</p> <p>The calculated volume is displayed. The dosing can be started with &lt;GO&gt;. With &lt;CLEAR&gt; the sequence returns to the input of the content for the solution (display "cnt").</p>	<p>%, ppm, g/L, mg/L, mol/L, mmol/L, mol/kg, mmol/kg</p> <p>0..±1E33</p> <p>0..±1E33</p>