


Mode selection



Press key <MODE> until the desired mode is displayed, confirm with <ENTER>.

Modes:

KFT	Karl Fischer titration.
TITER with H2O or Std.	Titer determination with water or with methanol standard.
TITER with Na2Tart*2H2O	Titer determination with sodium tartrate.
BLANK Determination	Blank determination.

Parameters



Input of titration parameters.
Preselections for the titration sequence and the report output.

Display	Initial value	Meaning	Input range
>titration parameters		General titration parameters	
extr. time	0 s	Extraction time. Negative number: no addition of reagent during extraction time.	0...±9999 s
stop crit.:	drift	Type of stop criteria.	drift, time
stop drift	20 µl/min	Titration stops if stop drift is reached.	1...999 µL/min
t(delay)	10 s	Titration stops if there is no dosing during t(delay).	0...99 s,
stop V:	99.99 mL	Stop volume, to avoid titration vessel overflow.	0.00...99.99 mL
start V:	0.00 mL	Start volume in mode KFT. Predosing without controlling.	0.00...99.99 mL
dos.rate	max. mL/min	Dosing rate for start volume.	0.01...150 mL/min, max.
max. Rate	max. mL/min	Maximum dosing rate during titration.	0.01...150 mL/min, max.
min.volume incr.	min. µl	Minimum volume increment during titration.	0.1...9.9 µL/min, min.
>Preselections		Preselections for the sequence	
conditioning:	ON	Automatic conditioning of titration vessel.	ON, OFF
req.ident:	OFF	Requests identification after start of titration.	ON, OFF
req.smpl size:	ON	Requests sample size after start of titration.	ON, OFF
report:	OFF	Selection of result report at the end of the titration.	full, short, OFF

Configuration



KF device settings.
Settings of RS232 interface.
Settings of peripheral units.
General settings.

Display	Initial value	Meaning	Input range
>KF device settings		KF device settings	
limit KF reagent	OFF ml	Counter for KF reagent (buffer capacity, reagent supply). Actual value of counter.	0...999 mL, OFF 0...999 mL
actual KF reagent	0 ml		
polarizer:	I(pol)	Selection of polarization type.	I(pol), U(pol)
I(pol)	50 µA	Polarization current (voltametric indication).	-127...127 µA
EP	250 mV	Corresponding end point voltage.	-1500...1500 mV
U(pol)	400 mV	Polarization voltage (amperometric indication).	-1270...1270 mV
EP	25 µA	Corresponding end point current.	-150...150 µA
filling rate max.	ml/min	Filling rate after the titration.	0.01...150 mL/min, max.
>RS232 settings		RS232-settings	
baud rate:	9600	Baud rate.	300, 600, 1200, 2400, 4800, 9600
data bit:	8	Data bit.	7, 8
stop bit:	1	Stop bit.	1, 2
parity:	none	Parity.	none, odd, even
handshake:	HWS	Handshake.	HWS, HWf, SWchar, SWline, none
RS control:	ON	Receiving of commands via RS.	ON, OFF
>peripheral units		Settings of peripheral units	
send to:	IBM	Selection of printer.	Epson, Seiko, IBM
balance:	Sartorius	Selection of balance.	Sartorius, Mettler, AND, Precisa
record:	U	Selection of record type to be outputted at the analog output.	V vs. t, dV/dt vs. t, U vs. t, -U vs. t
>auxiliaries		General settings	
dialog:	english	Selection of dialog language.	english, deutsch, français, español
date	YYYY-MM-DD		
time	HH:MM		
run number	0	Current run number for result output.	0...999
electrode test:	ON	Titration sequence with or without electrode test.	ON, OFF
display KFR vol.:	ON	Result display with or without the final KFR volume.	ON, OFF
device label		Device label.	up to 8 ASCII characters
program	787.0010	Program version.	read only

Calculations




Input of the values for result and statistics calculations in the KFT mode.

- Not every request is active in the other modes.

Display	Initial value	Meaning	Input range
>calculation		Values for result calculation	
sample size	1.0 g	Sample size.	up to 6 digits: $\pm X.XXXXX$
ident.		Sample identification.	up to 8 ASCII characters
titer	5.0 mg/ml	Titer.	0.0000...99.9991mg/mL
factor	0.1	Factor.	0... \pm 1000000
divisor	1.0	Divisor, e.g. for density.	0... \pm 1000000
blank	0.0 ml	Blank.	0.0000...99.9991mL
drift corr:	OFF	Drift correction of KFR volume.	auto, man., OFF
drift value	0.0 μl/min	Input of drift value for manual drift correction.	0.0...99.9 μ L/min
>statistics		Statistics calculation	
mean	n= OFF	Statistics with n single results.	2...20, OFF
res.tab:	original	Result table for statistics calculation.	original, delete n, delete all
delete	n= 1	Delete single result number n.	1...20

Units



Choice of the units for sample size and result in mode KFT.

For the other modes the result unit is read only.

Display	Initial value	Meaning	Input range
result unit:	%	Unit of the result.	%, ppm, mg/mL, g, mg, mL, mg/pc, no unit
result unit:	%;2	Number of decimal places.	0...9
smpl size unit:	g	Unit of the sample size.	g, mg, mL, μ L, pc, no unit

Calculation formula and units

The calculation formulas are preset and depend on the selected mode.

KFT	$\text{water(content)} = \frac{(\text{volume(KFR)} - \text{blank}) \times \text{titer} \times \text{factor}}{ \text{smplsize} \times \text{divisor}}$ <p>Factor and divisor are used for the conversion of the result to different units.</p>
------------	--

<i>Unit of result</i>	<i>Sample size in...</i>	<i>Factor</i>	<i>Divisor</i>
%	g	0.1	1
%	mg	100	1
%	mL	0.1	density of sample [g/mL]
ppm	g	1000	1
ppm	mL	1000	density of sample [g/mL]
ppm	μL	1	density of sample [g/mL]
mg/mL	g	1	density of sample [g/mL]
mg/mL	mL	1	1
mg	1	1	1
mL	1	1	1000 * density H ₂ O [g/mL] ≈ 1000
mg/pc	pc	1	1

TITER	$\text{titer} = \frac{ \text{smplsize} \times \text{factor}}{\text{volume(KFR)}}$ <p>The factor is used to enter the water content of the standard.</p>
--------------	--

<i>Standard used</i>	<i>Sample size in...</i>	<i>Factor</i>
Water	g	1000
Water	μL	density H ₂ O [g/mL] ≈ 1
Methanol	mL	content of methanol [g/mL]
Methanol	μL	0.001 * content of methanol [g/mL]
Na ₂ Tart*2H ₂ O	g	156.6
Na ₂ Tart*2H ₂ O	mg	0.1566

BLANK	<p>blank = volume(KFR) × factor</p> <p>If the blank is determined with an excessive solvent volume, the result can be converted to the amount used subsequently with the aid of the factor.</p>
--------------	---