

# **789 Robotic Sample Processor XL** **778 Sample Processor**

---

Program version 5.789.0020+  
and 5.778.0020+

## **Instructions for Use**

Teachware  
Metrohm AG  
Oberdorfstr. 68  
CH-9100 Herisau

*These instructions are protected by copyright. All rights reserved.*

*Although all the information given in these instructions has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please inform the author at the address given above.*

© Metrohm AG 2007  
Printed in Switzerland

# Table of Contents

<b>1</b>	<b>Introduction</b>	<b>1</b>
1.1	<b>Instrument description</b> .....	1
1.1.1	Field of application of the Metrohm Sample Processors .....	1
1.1.2	Instrument versions .....	3
1.2	<b>Sample Processors as part of a system</b> .....	5
1.3	<b>Information about these Instructions for Use</b> .....	6
1.3.1	Please note .....	6
1.3.2	Additional documents .....	6
1.3.3	Notation and pictograms .....	6
1.4	<b>Parts and controls</b> .....	7
1.4.1	Overall view .....	7
1.4.2	Rear view .....	8
1.5	<b>Connections</b> .....	9
1.5.1	Sensors of the Sample Processor .....	9
1.6	<b>Accessories</b> .....	11
1.7	<b>The keypad</b> .....	13
1.7.1	The display .....	14
1.7.2	The keys .....	15
1.8	<b>Safety information</b> .....	23
1.8.1	General: .....	23
1.8.2	Electrical safety .....	23
1.8.3	Personal protection .....	24
<b>2</b>	<b>Installation</b>	<b>25</b>
2.1	<b>Installation flow chart</b> .....	25
2.2	<b>Instrument setup</b> .....	26
2.2.1	Setup .....	26
2.2.2	Mains connection .....	26
2.2.3	Connecting the keypad .....	27
2.2.4	Mounting a 786 Swing Head .....	27
2.2.5	Mounting the robotic arm .....	28
2.2.6	Connecting pumps .....	29
2.2.7	Connecting dosing devices and stirrers .....	30
2.3	<b>Data transmission connections</b> .....	31
2.3.1	Remote connections .....	31
2.3.2	Serial connections (RS232) .....	37
2.3.3	Connecting a printer .....	38
2.4	<b>Connecting up accessories</b> .....	40



3.7.1	Metrohm standard sample racks .....	90
3.7.2	Magnet codes .....	91
3.7.3	Rack data .....	92
<b>3.8</b>	<b>Dosing and liquid handling.....</b>	<b>95</b>
3.8.1	Dosimats and Dosinos .....	95
3.8.2	Liquid handling functions .....	100
3.8.3	The DOS command .....	101
3.8.4	Pictograms .....	101
3.8.5	Liquid handling functions in detail .....	102
3.8.6	Pipetting equipment.....	105
3.8.7	Pipetting procedure .....	105
3.8.8	Preparing the Dosing unit .....	106
3.8.9	Pipetting .....	106
<b>3.9</b>	<b>The Remote interface .....</b>	<b>109</b>
3.9.1	Output lines .....	109
3.9.2	Input lines .....	109
3.9.3	SCN command .....	110
3.9.4	CTL command .....	110
3.9.5	Manual stop options .....	111
<b>3.10</b>	<b>LEARN mode.....</b>	<b>112</b>
3.10.1	Setting lift and robotic arm positions .....	112
3.10.2	Rack adjustment .....	113
3.10.3	Parametrizing sequence commands .....	113
<b>3.11</b>	<b>TRACE function .....</b>	<b>114</b>
<b>3.12</b>	<b>Disabling keypad functions .....</b>	<b>114</b>
3.12.1	Disable whole keypad.....	114
3.12.2	Disable configuration .....	115
3.12.3	Disable parameter .....	115
3.12.4	Disable method storage functions .....	115
3.12.5	Disable display.....	115
<b>3.13</b>	<b>786 Swing Head settings .....</b>	<b>116</b>
<b>4</b>	<b>Service, maintenance, errors</b>	<b>119</b>
4.1	<b>Service.....</b>	<b>119</b>
4.1.1	Running time meter.....	119
4.2	<b>Care and maintenance .....</b>	<b>119</b>
4.3	<b>Error messages .....</b>	<b>120</b>
<b>5</b>	<b>GLP validation – diagnosis</b>	<b>122</b>
5.1	<b>Validation / GLP .....</b>	<b>122</b>
5.2	<b>Initializing the working memory .....</b>	<b>123</b>

<b>6</b>	<b>Annex</b>	<b>125</b>
6.1	<b>Technical data.....</b>	<b>125</b>
6.1.1	Keypad.....	125
6.1.2	Interfaces .....	125
6.1.3	MSB connections .....	125
6.1.4	Pumps and pump connections.....	125
6.1.5	Swing Head connection .....	126
6.1.6	Lift .....	126
6.1.7	Turntable.....	126
6.1.8	Stirrer connection (DIN socket) .....	126
6.1.9	Mains connection .....	126
6.1.10	Safety specifications.....	126
6.1.11	Electromagnetic compatibility (EMC).....	126
6.1.12	Ambient temperature.....	127
6.1.13	Dimensions and materials.....	127
6.2	<b>Standard methods .....</b>	<b>128</b>
6.2.1	Titrimo .....	129
6.2.2	PIP ext.....	130
6.2.3	KF_ext.....	132
6.2.4	pH_cal.....	134
6.2.5	Std_add .....	135
6.3	<b>Standard equipment.....</b>	<b>137</b>
6.3.1	Metrohm Sample Processor:.....	137
6.3.2	Sample racks and sample beakers .....	149
6.3.3	786 Swing Head .....	153
6.3.4	Connection cables.....	157
6.3.5	Optional accessories and additional devices.....	158
6.3.6	Electrodes for automation .....	161
6.4	<b>Warranty and conformity.....</b>	<b>1621</b>
6.4.1	Warranty.....	1621
6.4.2	Declaration of Conformity (778 Sample Processor) .....	1632
6.4.3	Declaration of Conformity (789 Robotic Sample Processor XL) .....	1643
6.4.4	Quality Management Principles .....	1654
<b>7</b>	<b>Index</b>	<b>1665</b>


## List of Illustrations

Fig. 1	System components .....	5
Fig. 2	Overall view .....	7
Fig. 3	Rear view .....	8
Fig. 4	Connection strip.....	9
Fig. 5	Magnet sensor for rack code.....	9
Fig. 6	Beaker sensor on the tower .....	10
Fig. 7	Sensor on a robotic arm .....	10
Fig. 8	Accessories .....	11
Fig. 9	Keypad.....	13
Fig. 10	Safety Shield (example shown: 6.2751.0xx for transfer robotic arm) .....	24
Fig. 11	Mounting the Swing Head .....	27
Fig. 12	Robotic arm with limiting screw .....	28
Fig. 13	Mounting the robotic arm .....	29
Fig. 13	MSB connections.....	30
Fig. 14	Remote cable.....	31
Fig. 15	RS232 connections.....	37
Fig. 16	Distributor piece.....	40
Fig. 17	Attaching the tubing.....	41
Fig. 18	Macro titration head.....	42
Fig. 19	Micro titration head .....	42
Fig. 20	Robotic arm with titration accessories .....	43
Fig. 22	Transfer robotic arm with pipetting accessories.....	44
Fig. 23	Transfer robotic arm with double-hollow needle.....	45
Fig. 24	Spray nozzles.....	45
Fig. 25	741 Magnetic stirrer .....	46
Fig. 24	Stand support .....	46
Fig. 25	Collection trough.....	47
Fig. 26	Attaching a sample rack .....	47
Fig. 27	Sample rack for XL models.....	48
Fig. 28	Mounting the safety shield.....	48
Fig. 29	Safety shield for Swing Heads .....	48
Fig. 30	Dialog arrangement .....	53
Fig. 31	Text input .....	56
Fig. 32	800 Dosino with Dosing units .....	95
Fig. 33	Dosing unit from below.....	96
Fig. 34	Dosing unit - ports .....	100
Fig. 35	Remote interface.....	125

## Tables

Table 1	Model versions (1 tower) .....	3
Table 2	Model versions (2 towers).....	4





# 1 Introduction

*This section offers you a first overview of the Metrohm Sample Processors. All the information applies to both the 789 Robotic Sample Processor XL and the 778 Sample Processor. You are informed about how you can use these versatile instruments and are introduced to the most important parts and controls.*

## 1.1 Instrument description

### 1.1.1 Field of application of the Metrohm Sample Processors

The Metrohm Sample Processors are very versatile instruments intended exclusively for use in factories and laboratories, where they cover a wide range of applications. In this way they provide an indispensable service wherever large series of samples have to be processed, no matter whether in the titration, measurement or liquid handling sectors.

As a result of the extensive communication possibilities you can work via the parallel Remote and serial RS232 interface not just with the wide range of Metrohm titrators, measuring instruments and dosing devices, but you can also work with any other instruments that have a suitable communications interface; these can be controlled by or control the Sample Processor. These abilities mean that they are predestined for all imaginable automation tasks in a modern laboratory, even within highly integrated laboratory data systems.

Despite their comprehensive range of commands and the numerous possible configurations, the Metrohm Sample Processors offer an uncomplicated type of operation that is also suitable for routine work as a result of the possibility of managing user-defined methods.

The standard methods supplied with the instruments can be used for routine tasks without any further fuss. After a short familiarization period the user can alter them to meet any particular requirements and store them in the instrument. This means that, apart from routine work, Metrohm Sample Processors can also be used for demanding special applications.

The run sequences for processing the individual samples are freely definable within wide limits. The same applies to the start sequence and

final sequence, which have to be carried out once before the start of a sample series or once after it has been completed. This offers many advantages, particularly for titrations. The electrode can be conditioned before the first titration or subjected to a special rinsing process.

A learn mode is available for creating run sequences and, with its help, command parameters can be set in manual operation.

Exchangeable standard sample racks are available for many sizes of beakers and test tubes. Freely selectable "special beaker" positions can be defined for each rack. These are then used for including rinsing or conditioning beakers, which can be addressed in any part sequence, on the rack.

By extending the system with a 786 Swing Head the number of samples to be processed on a rack can be considerably increased. The robotic arm of the 786 Swing Head allows any point on a sample rack to be addressed. This means that the number (max. 999 rack positions) and arrangement of the samples is virtually unlimited. On request we can supply customer-specific special racks for individual requirements.

Freely definable position tables can be loaded via the RS232 interface and suitable PC software for the configuration of special racks.

### 1.1.2 Instrument versions

Various equipped Metrohm Sample Processors models are available:

- All versions of the **789 Robotic Sample Processor XL** are suitable for sample racks of up to **48 cm** diameter.
- All versions of the **778 Sample Processor** are suitable for sample racks of up to **42 cm** diameter.







 <p><b>Model 2.789.0010</b> <b>Model 2.778.0010</b></p> <p>1 tower with 1 pump + 1 ext. pump connection + 1 stirrer connection + 1 Swing Head connection</p> <p>Chassis with 3 MSB sockets for dosing devices and/or stirrers + Remote connection (25-pin) + RS232 connection (9-pin) + keyboard connection</p>	 <p><b>Model 2.789.0020</b> <b>Model 2.778.0020</b></p> <p>1 tower with 2 pumps + 1 stirrer connection + 1 Swing Head connection</p> <p>Chassis with 3 MSB sockets for dosing devices and/or stirrers + Remote connection (25-pin) + RS232 connection (9-pin) + keyboard connection</p>	 <p><b>Model 2.789.0030</b> <b>Model 2.778.0030</b></p> <p>1 tower without pumps + 2 ext. pump connections + 1 stirrer connection + 1 Swing Head connection</p> <p>Chassis with 3 MSB sockets for dosing devices and/or stirrers + Remote connection (25-pin) + RS232 connection (9-pin) + keyboard connection</p>
---	---	--

Table 1 Model versions (1 tower)

 <p><b>Model 2.789.0110</b> <b>Model 2.778.0110</b></p> <p>2 towers with 2 pumps + 2 ext. pump connections + 2 stirrer connections + 2 Swing Head connections</p> <p>Chassis with 3 MSB sockets for dosing devices and/or stirrers + Remote connection (25-pin) + RS232 connection (9-pin) + keyboard connection</p>	 <p><b>Model 2.789.0120</b> <b>Model 2.778.0120</b></p> <p>2 towers with 4 pumps + 2 stirrer connections + 2 Swing Head connections</p> <p>Chassis with 3 MSB sockets for dosing devices and/or stirrers + Remote connection (25-pin) + RS232 connection (9-pin) + keyboard connection</p>	 <p><b>Model 2.789.0130</b> <b>Model 2.778.0130</b></p> <p>2 towers without pumps + 4 ext. pump connections + 2 stirrer connections + 2 Swing Head connections</p> <p>Chassis with 3 MSB sockets for dosing devices and/or stirrers + Remote connection (25-pin) + RS232 connection (9-pin) + keyboard connection</p>
---	---	--

*Table 2 Model versions (2 towers)*

## 1.2 Sample Processors as part of a system




Fig. 1 System components

Titrandos  
 (6.2148.010 Remote Box required)

## 1.3 Information about these Instructions for Use

### 1.3.1 Please note







*Please read through these Instructions for Use carefully before you start to use the Sample Processor. The instructions contain information and warnings that must be observed by the user in order to guarantee the safe use of the instrument.*

### 1.3.2 Additional documents

- Quick Reference 8.789.1013 for the Metrohm Sample Processors
- Operating Tutorial 8.789.1023 for the Metrohm Sample Processors
- Technical Reference 8.789.1033 for the Metrohm Sample Processors

### 1.3.3 Notation and pictograms

The following notation and pictograms are used in these Instructions for Use:

Location	Menu item, parameter or input value
<OK>	button, key
	<b>Danger</b> This symbol indicates a possible risk of death or injury to the user if the instructions are not followed correctly.
	<b>Warning</b> This symbol indicates a possible risk of damage to the instruments or their components if the instructions are not followed correctly.
	<b>Attention</b> This symbol indicates important information. Read the information provided before you continue.
	<b>Remarks</b> This symbol indicates additional information and tips.

# 1.4 Parts and controls

## 1.4.1 Overall view




Fig. 2 Overall view

<b>1</b>	<b>Guide chain</b>
<b>2</b>	<b>Safety shield/Splash protection</b>
<b>3</b>	<b>Sample rack</b>
<b>4</b>	<b>Stirrer rail</b>
<b>5</b>	<b>Splash protection fixing</b>

<b>6</b>	<b>Splash protection guide</b>
<b>7</b>	<b>Titration head</b>
<b>8</b>	<b>Lift</b>
<b>9</b>	<b>Tower</b>
<b>10</b>	<b>Beaker sensor</b>
<b>11</b>	<b>Chassis</b>

### 1.4.2 Rear view

This illustration shows the rear view of the standard model 2.778.0010 with one tower, one membrane pump and one connection for an external pump.




Fig. 3 Rear view

<b>12</b>	<b>Distributor block</b>	<b>16</b>	<b>Stirrer connection (Tower 1)</b> For 802 Rod Stirrer or 741 Magnetic Stirrer
<b>13</b>	<b>Membrane pump</b>	<b>17</b>	<b>Connection strip</b>
<b>14</b>	<b>PTFE tubing</b>	<b>18</b>	<b>Pump connection M8 (external)</b>
<b>15</b>	<b>Solenoid valve</b>	<b>19</b>	<b>Connection socket for the 786 Swing Head</b>
<b>a</b>	<b>Warning: Biohazard</b> See Section 1.8.3 <i>Personal protection</i>	<b>b</b>	<b>Warning: Resistance to chemicals</b> See section 2.4.1 <i>Connecting rinsing and aspiration equipment</i>

## 1.5 Connections

The electrical connections are the same for all models of the Metrohm 778/789 Sample Processor series.




Fig. 4 Connection strip

<b>20</b>	<b>Mains switch</b>	<b>24</b>	<b>Serial RS232 connection (9-pin)</b>
<b>21</b>	<b>Keyboard connection</b>	<b>25</b>	<b>Mains connection</b>
<b>22</b>	<b>MSB connections MSB1 ... MSB3</b> Metrohm Serial Bus Connection of dosing devices and stirrers	<b>26</b>	<b>Type plate</b>
<b>23</b>	<b>Remote connection (25-pin)</b>		

### 1.5.1 Sensors of the Sample Processor





Fig. 5 Magnet sensor for rack code

#### Rack sensor

The magnet sensor for recognizing the individual rack code is located beneath the turntable of the Sample Processor. The magnet code of a rack can only be read in when the rack is in the initial position. The magnet holder must be positioned directly above the sensor.


For this reason the Sample Processor should be initialized each time that a rack is changed with the <RACK> key.



*Fig. 6 Beaker sensor on the tower*

### Optical beaker sensor

Each tower of a Metrohm Sample Processor is equipped with a beaker sensor that detects the presence of a beaker in front of the tower. With this infrared sensor beakers made from different materials can be detected, provided that they are located in the correct position in front of the tower and the beaker sensor 'Tower' has been selected in the rack definition. This "Beaker test" is carried out after each MOVE command (i.e. each rotation of the rack). The beaker sensor on the tower can only be used with single-row sample racks.




*Fig. 7 Sensor on a robotic arm*

### Tactile robotic arm sensor

Robotic arms with piezo sensors may be used with multi-row racks. The sensor is activated, when the lift is run to work position and makes physical contact with a sample beaker.

## 1.6 Accessories



*Fig. 8 Accessories*

With suitable accessories a Sample Processor can be extended to form a comprehensive automated system. Depending on the tasks to be carried out, various standard components or even custom-made special parts can be used. Please consult the list of accessories on page 137ff.

### **Sample racks**

Custom-made racks for various vessel sizes with any arrangement of rack positions can also be supplied in addition to the standard racks.

### **786 Swing Head with robotic arm**

The use of multi-row sample racks or external titration cells requires the use of a 786 Swing Head. This motor drive, which is mounted on the lift of one of the Sample Processor towers, can move different types of robotic arm. Various types of standard robotic arms with accessories for titrating on the sample rack or for sample transfer to an external titration cell are available.

### **Stand support**

If an external titration cell is to be used then we recommend the use of a stand support. The stand support can accommodate a magnetic stirrer (e.g. model 801) and, with a support rod, enables you to attach a titration cell or other accessories.

**Swing head with tactile piezo sensor**

When using a multi-row rack and a 786 Swing Head, a robotic arm with a touch-sensitive piezo sensor may be used. With that kind of sensor the presence of a sample beaker can be detected reliably.

**Robotic arm for removing covers**

Whenever covered sample containers are required (e. g. with volatile samples) the lids can be removed by a special robotic arm (so called Dis-Cover) with magnetic contact before the sample treatment. Special lids are available for this purpose.

**Accessory set for sample transfer**

For pipetting samples into external titration cells an accessory set (6.5619.000) is available, comprising all necessary parts.

## 1.7 The keypad




Fig. 9 Keypad

Below the 2-line display there are three LEDs. The two LEDs '**TOWER 1**' and '**TOWER 2**' indicate the tower that is currently active. The LED '**LEARN**' lights up when the learn mode is activated.

Most keys have two functions, depending on whether the Sample Processor is in the normal operating mode or in the editing mode.

Selection menus can be accessed with the upper row of keys (<**CON-FIG**>, <**PARAM**>, <**USER METHOD**>). The other keys on the left-

hand side of the keypad are used for navigation in the menus or for altering parameters. For entering parameters the numerical block on the right-hand side of the keypad is also available.

The lowest row of keys (<HOLD>, <STOP>, <START>) are used for the direct control of a method sequence.


### 1.7.1 The display

The display consists of two lines each with 24 characters.


The first line is used as the title line, in which the current method and the sample counter count are shown. In the editing mode the menu title appears here.

The second line is used as the status line, which shows specific activities depending on the operating status. In the editing mode it is used as the input line.


#### Normal condition



#### Method sequence




#### Editing mode




If the Sample Processor is included in a computer-controlled automation system and is completely controlled via the RS232 interface then it may be advisable to switch off the display. This can be done in the setup menu of the Sample Processor, see p. 114.

## 1.7.2 The keys

### The menu keys



The **<Config>** key opens the selection menu for the configuration of the Sample Processor.




The **<PARAM>** key opens the selection menu for editing the run sequences and method parameters.



The **<User Method>** key opens the selection menu for opening, saving and deleting default or user-defined methods.

### Lift operation and sample positioning keys



With the **<↑>** and **<↓>** keys the lift of the currently active tower can be moved up and down respectively. The lowest possible lift position is defined by the configuration parameter 'max. stroke path'.


- ☞ In the editing mode the arrow keys **<↑>** and **<↓>** are used for navigation in the particular menu or submenu.



With the **<HOME>** key the lift of the currently active tower is returned to the rest position (0 mm), i.e. to the uppermost stop.


**<END>** moves the lift to the predefined working position (see p. 61).

- ☞ In the editing mode the **<HOME>** and **<END>** keys move to the first and last line of the menu or submenu respectively.



With the **<← NEXT>** and **<→ PREV>** keys the sample rack can be rotated forward or backward by one position. The beaker positions depend on the active lift. If necessary, the lift (or both lifts) are automatically raised to the shifting position. When the rack position has been reached a robotic arm which may be mounted will automatically be directed to the corresponding rack position.

- ☞ In the editing mode the arrow keys **<←>** and **<→>** are used for navigation in a menu line.




Most functions for manual operation apply specifically to a single tower. With the 2-tower models the **<SELECT/TOWER>** key can be used to switch between the towers. The currently active tower is indicated by the **TOWER 1** or **TOWER 2** LED. The following commands or keys refer to the active tower: **MOVE**, **<←>**, **<→>**, **LIFT**, **<↑>**, **<↓>**, **<HOME>**, **<END>** and **<PUMP>**.

- ☞ During data input the **<SELECT>** key is used to select a predefined entry from a selection list.

## Editing and sequence control




When editing a method sequence the **<INSERT>** and **<DELETE>** keys are used to insert or delete a command line.




The **<CLEAR/RESET>** key is used to initialize the Sample Processor and dosing devices. This corresponds to the switching-on process.

- ➔ During data input the **<CLEAR/RESET>** key is used to delete an entry or to reset the default value. In text entry mode the last character is deleted.



During a method sequence the **<QUIT>** key can be used to terminate the command which is currently being carried out. The following command is executed.

- ➔ During data input the **<QUIT>** key is used to terminate an entry. During navigation in a menu the **<QUIT>** key is used to exit the active (sub)menu and select the next highest menu level.



During data input the **<ENTER>** key is used to accept the entry.

## Command keys



The **<SAMPLE>** key is used to set the current sample position. This has to be done before a sample series is run.

At the start of a method this position is assumed to be the first sample in a series. If no sample position has been set then the Sample Processor will select rack position 1.



With **<MOVE>** a vessel or a particular rack position can be moved to the active tower or a robotic arm can be swung to an external position. The **<SELECT>** key is used to select the tower.

As well as the actual sample beaker a maximum of 16 possible special beakers can also be defined. A particular rack position can be moved to directly by entering the position number (with the numerical keys).

The direction and speed of rotation can be altered in the parameter menu or with the **<DEF>** key.

### Important:

For safety reasons it is only possible to rotate the sample rack when the lift or both lifts are located in the shifting position or above it. During a rack rotation the lift (or both lifts) are automatically first raised to the predefined shift height.






Raises or lowers the lift of the active tower. The predefined lift positions (working position, rest position, rinsing position, shifting position, special position) can be selected with the **<SELECT>** key. They can be entered and saved separately for each rack in the configuration menu.

As well as the predefined lift positions it is also possible to enter absolute lift positions in mm via the numerical keys.

In the 2-tower models the tower can be selected with **<SELECT/TOWER>**.



The **<PUMP>** key is used for switching the pumps 1 or 2 of the active tower on and off. By entering the pump number (1 or 2) the condition of the corresponding pump will be switched, i.e. if the pump is switched off it will be switched on and vice versa. With Sample Processor models that have no pump or only one built-in membrane pump the selected pump connection will be switched on or off.

The status of all pumps is shown in the display (e.g. **PUMP--++**; + means switched on, - means switched off).

Example:

**PUMP on/off No. ? <2>** Display: **PUMP -+--**

**PUMP on/off No. ? <2>** Display: **PUMP ----**

In this case pump 2 is switched on and off.

Under **<PARAM>**, **>manual stop** you can define whether the pumps should be switched off with the **<STOP>** key or not.



The **<STIR>** key is used for controlling the stirrers. A stirrer can be switched on permanently or switched on for a given period and then switched off again. The **<SELECT>** key is used to select both the stirrer and the function. The current status of the stirrer is shown directly in the display.

Example:

**STIR: T1** : **ON s** Display: **STIR +-** (+ = on - = off)

**STIR: MSB2** : **10 s** Display: **STIR 10 s**

In this case in the first line the stirrer at tower 1 is switched on. The stirrer is selected with the **<SELECT>** key. As can be seen in the second line, the duration of the stirring process can also be entered.

The stirring rate can be set for each stirrer in the parameter menu or with the **<DEF>** key.

Under **<PARAM>**, **>manual stop** you can define which stirrers can be switched off with the **<STOP>** key.



The **<DOS>** key is used to control the connected dosing devices. Both positive and negative volumes can be dosed. Negative volumes are used for aspirating liquids, e.g. during pipetting.

As well as entering the volume to be dosed (with the numerical keys),

**<SELECT>** can also be used to select additional functions:

- Filling the dosing or exchange unit (fill)
- Initializing the escange of a Dosing unit (release)
- Preparing the tubing systems and cylinder (prep.)
- Emptying the tubing system and the dosing cylinder (empty)
- Ejecting the cylinder contents (Eject)
- Driving the piston to the max. volume
- Compensating for the play between piston and spindle (compen.)
- Valve switching (port)

The first parameter of the DOS command stands for the number of the dosing instrument (1...3, \* = all) and the Dosino port (e. g. 1.1 stands for Dosino 1, port 1), the second parameter for the function or the volume to be dosed.

Example:

```
DOS: 2.1 <ENTER> 4.51 ml <ENTER>
DOS: 2.* <ENTER> <SELECT> ... fill <ENTER>
```

The dosing and filling rates can be set in the parameter menu or with the **<DEF>** key.



Shows the incoming signals or data from the Remote or the serial RS232 interfaces.

This function is used for checking the data communication with connected devices.

The first parameter shows the selected interface. The second parameter shows the signals or data that are received directly.

If the parallel Remote interface (**Rm**) is selected then the signal states of the incoming Remote lines are shown in binary form (1=line active, 0=line inactive).

If the serial RS232 interface (**RS**) is selected then the data string received via this interface will be shown (14 characters per line).

Example (Remote interface):

```
SCN: Rm :0000001
```

In this case the **Ready** line (Remote line input 0) of a connected Tritino is set.



Controls external devices via the Remote or RS232 interface.

The first parameter sets the interface (**<SELECT>**). The second parameter defines the status of the lines (Remote lines) or data (RS232 interface) to be outputted via the selected interface.

### 2nd parameter, for Remote interface

Binary pattern with 14 digits (0, 1 or \*) for the 14 output lines or prede-

defined binary pattern (<**SELECT**> selection), e.g. **START device 1, INIT** etc.


### 2nd parameter, for RS232 interface

Data string with up to 14 alphanumeric characters (any).

The default value "&M;\$G" (for starting Metrohm instruments) can be set with <**CLEAR**>.



The <**WAIT**> key has no function in the normal operating condition. It is used to insert the **WAIT** command in a run sequence.




The <**DEF**> key is used to edit various settings for manual operation. Repeated pressing of the <**DEF**> key is used to select the various settings. In order to change the entry you must first press <**ENTER**> and then enter the new value.

Alterations made in this way only apply to manual operation.

## DEF settings

---



DOSRATE

### Change dosing rate


- The dosing rate in mL/min can be set separately for each dosing device, see p. 70.
- Syntax: **DOSRATE** [dosing device] [dosing rate]



FILLRATE

### Change filling rate


- The filling rate in mL/min can be set separately for each dosing device, see p. 70.
- Syntax: **FILLRATE** [dosing device] [filling rate]



COCKMOVE

### Direction of stopcock rotation


- For each connected Dosino the direction of rotation of the stopcock switching can be defined separately, see p. 70.
- Syntax: **COCKMOVE** [Dosing device] [Direction of rotation]




LIFTRATE

### Change lift rate


- The lift speed in mm/s can be set separately for each tower (for 2-tower models), see p. 68.
  - Syntax: **LIFTRATE** [Tower] [Lift rate]
-

**SHIFTRATE****Change rack rate and direction of rotation**

- As well as the speed of rotation of the sample rack in degrees/second the direction of rotation can also be defined.
- Direction of rotation "+" means that the rack always rotates counterclockwise, i.e. in increasing rack position sequence, Direction of rotation "-" means clockwise, i.e. decreasing sequence.
- Direction of rotation **auto**: the Sample Processor itself selects the shortest possible path for a rack rotation, see p. 68.
- Syntax: **SHIFTRATE** [Direction of rotation] [Rotation rate]

**SWINGRATE****Change swing rate**


- The swing rate in degrees/s of a robotic arm can be set separately for each connected 786 Swing Head, see p. 68.
- Syntax: **SWINGRATE** [Tower] [Swing rate]

**STIRRATE****Change stirring rate**

- The stirring speed can be controlled separately for each stirrer (rod or magnetic stirrer), see p. 69.
- Syntax: **STIRRATE** [Stirrer] [Stirrer rate]

**Auxiliary commands**


---




The **<PRINT>** key is used to print out a report. The following can be selected: parameter report (method), configuration report, list of stored methods, all reports.

The type of report can be selected with **<SELECT>**.

e. g.

**Print:                    config**

The selection of the printer type and the settings of the RS232 interface must be made in the configuration menu under **>RS232 Settings**, see p. 37.



With the **<RACK>** key the sample rack can be initialized. The connected peripheral devices (e.g. Dosimats, Dosinos) are not affected by this.

The sample rack and lift (both lifts in 2-tower versions) are moved to the zero position and automatic rack recognition is carried out.

## Sequence control

---



The **<START>** key starts a method. A start is only possible when the Sample Processor is in the normal operating condition, i.e. when the display shows **ready**.


If **<START>** is activated after an interruption (**<HOLD>**, see below), then the sequence continues with the next command.

The **<START>** can also be used to carry out a single command line in a run sequence (**TRACE** function), see p. 114.



The **<STOP>** key ends a method.

If a sample series is stopped manually with **<STOP>** then the final sequence of a method will not be carried out. When the **<STOP>** key is pressed the functions listed in the parameter menu under **>manual stop** will be carried out.



The **<HOLD>** key interrupts a method sequence.

Connected peripheral devices (Titrimos, etc.) will not be stopped automatically. Only the method sequence will be interrupted. In the **HOLD** condition a method can be completely terminated with **<STOP>** or continued with **<START>**.

After an error message in the method sequence the Sample Processor switches automatically to the **HOLD** condition after **<QUIT>**.

## 1.8 Safety information



### **Warning!**

*This instrument should only be used in accordance with the information given in these Instructions for Use.*

### 1.8.1 General:

This instrument left our works in perfect condition from the point of view of its operational safety (see Technical data, safety specifications). To keep it in this condition and to continue to operate safely the following information must be carefully observed.

### 1.8.2 Electrical safety

Please observe the following guidelines:

- Only qualified Metrohm personnel should carry out service work on electronic components.
- Do not open the instrument housing as this could damage the instrument. The housing contains no components which could be serviced or exchanged by the user.

Electrical safety when handling the instrument is guaranteed within the scope of Standard IEC 61010-1.

- **Protection against electrostatic charges**



### **Warning!**

*Electronic components are sensitive to electrostatic charges and can be destroyed by a discharge. Before touching any electronic components of the Sample Processor one should ground himself and his tools by grasping a grounded object (e.g. the instrument housing or a radiator) in order to eliminate any electrostatic charges that may be present.*

- **Mains connection**



*This instrument must only be operated at the specified **mains voltages** (see rear panel of instrument).*

- **Opening a Sample Processor**



*When the Sample Processor is connected to the mains supply the instrument must not be opened, nor should any of its components be dismantled as otherwise one could come into contact with current-carrying components. Before opening the instrument it must be separated from all current sources. **The mains cable has to be removed from the mains connection socket!***


### 1.8.3 Personal protection

The different drives of the Sample Processor and the 786 Swing Head possess considerable force in order that they can freely move even under heavy loads. They are fitted with electronic overload limiters that are triggered if the mechanical resistance is too high.



When the **instrument is running** do not intrude in the working range of the Sample Processor.

For the user there is a considerable **risk of injury from manual interference** with the automatic processes of the instrument.



Running methods and individual commands can be interrupted with the **<HOLD>** key and then continued with **<START>**.

#### Emergency stop



The **<STOP>** key always terminates all running processes immediately (exception: initialization during switch-on).

#### Important!



Before you use the instrument for the first time you must install the included **6.2751.0xx Safety Shield**, see p. 48.

The Sample Processor must not be operated without the Safety Shield in position!




Fig. 10 Safety Shield  
(example shown: 6.2751.0xx for transfer robotic arm)



#### Biohazard

The Sample Processor does not provide sufficient protection when processing potentially infectious samples or reagents.


In this case install the necessary safety devices.

## 2 Installation

*This section describes the installation of all the components of a Sample Processor and shows how it is connected to other instruments. Please read through the following sections carefully and follow the instructions in detail.*

### 2.1 Installation flow chart

The following flow chart provides an overview of all the necessary installation work. Detailed information can be found in the given sections.



\* optional

## 2.2 Instrument setup

### 2.2.1 Setup

#### Packaging

A Metrohm Sample Processor and its specially packed accessories are supplied in very protective special packaging consisting of shock-absorbing plastic foam. Please store this packaging in a safe place; it is the only way in which the safe transport of the instrument can be guaranteed.

#### Checks

Please check that the delivery is complete and undamaged immediately on receipt (compare with delivery note and list of accessories given in Section 6.3, p. 137). If transport damage is evident please refer to the information given in Section 6.4.1 "Warranty", p. 162.

#### Location

The Sample Processor is a robust instrument and can therefore be used even in rough surroundings in laboratories and factories.

However, please ensure that it is not exposed to a corrosive atmosphere. The instrument should be serviced at regular intervals, particularly when operated under rough conditions.


### 2.2.2 Mains connection




*This instrument must only be operated at the specified **mains voltages** (see rear panel of instrument).*

*Protect the connection sockets against moisture.*

Mains connection **25**



### 2.2.3 Connecting the keypad




The keypad is attached to the keyboard socket provided for it on the rear panel of the instrument. To remove the plug press both sides together.


### 2.2.4 Mounting a 786 Swing Head

The procedure for the installation of the 786 Swing Head is described below. This is normally carried out by the service technician.

For Sample Processors with 2 towers Tower 1 is moved to a central lift position and Tower 2 to the parking position.

**Switch off the instrument at the mains switch.**





### Connecting to the Swing Head connection socket

Lead the connection cable of the 786 Swing Head through the guide chain of the tower and connect the mini-DIN plug of the Swing Head connection cable to the **Swing Head** socket on the rear panel of the tower.

### Configuring the Swing Head and robotic arm

As each type of robotic arm has its own dimensions, it is essential to adapt the Swing Head setup settings to suit the particular type of robotic arm.

The necessary alterations are made in the setup dialog, see Section 3.13. The individual settings affect:

- Swing offset (default: **0°**)
- Max. swing angle (default: **84°**)
- Swing radius (=length of robotic arm, default: **110 mm**)
- Swing direction (default Tower1: **-**, default Tower 2: **+**)

The setup dialog of the Sample Processor is opened by pressing the **<CONFIG>** key as the instrument is switched on. A short explanation is given in Section 3.13. Details about the individual settings and types of robotic arm can be found in the Instructions for Use of the 786 Swing Head.

### 2.2.5 Mounting the robotic arm

Certain robotic arms can, depending on the model, be mounted so that they swing to the right or the left. During the mounting process the location of the limiting screw of the robotic arm must be taken into account.




Fig. 12 Robotic arm with limiting screw



**Swing-right mounting** means: The robotic arm can swing from the zero-axis (center of the rack) to the right, as seen from the front. The swing direction has to be set to **-**.

Robotic arms that can be mounted in two different ways (e. g. 6.1462.050) can have the limiting screw positioned according to the required swing direction (see above).

**Mounting procedure**

The 786 Swing Head must already be in position and configured before the robotic arm can be mounted.




Fig. 13 Mounting the robotic arm

- Switch the instrument off and then on again. Wait until the Swing Head drive has finished the initialization movements.
- Carefully push the robotic arm over guide cam **10** from below (see drawing). Hold the robotic arm in such a way that limiting screw **9** points toward the tower of the Sample Processor while the robotic arm assumes the maximum swing angle.
- Use the three enclosed screws to attach the robotic arm from below.
- Switch the instrument off and on again.

Details about the individual robotic arm types can be found in the Instructions for Use of the 786 Swing Head.

**2.2.6 Connecting pumps**

On the Sample Processor models with one built-in membrane pump one external pump per tower can be connected and controlled. The models with no pumps even have two connections for external pumps.

**Ext.  
Pump 2**



The pump connections (for 3-pin M8 connectors) on the rear panel of the Sample Processor tower supply a feed voltage of 16 V for a max. load of 600 mA.


Suitable pump models are:

- Metrohm 823 Membrane Pump Unit (membrane pump, suitable for aqueous media without particles)
- Metrohm 772 Pump Unit (peristaltic pump, suitable for organic solvents and aqueous solutions containing particles)

### 2.2.7 Connecting dosing devices and stirrers

Metrohm Sample Processors have various connections to which dosing devices and stirrers can be connected.

#### Stirrers



The **722/802 Rod Stirrers** and **741 Magnetic Stirrer** have a DIN connector and are connected to the rear panel of a tower (Tower stirrer T1 and T2).

As with other Metrohm instruments, e. g. Titrandos, the three **MSB connectors** (MSB = Metrohm Serial Bus) are very versatile. You can connect an **801 Stirrer** or an **804 Ti Stand** to them.

#### Dosing devices

Up to three **700/800 Dosino** or **685/805 Dosimat** type dosing drives can be operated by connecting them directly to the MSB sockets of a Sample Processor. One **800 Dosino** or **805 Dosimat** can also be connected to the **MSB output connector** of another instrument connected to the MSB socket, such as the 801 Stirrer, in a so-called '*daisy chain*'. Please note that only one instrument of a particular type may be present in such a daisy chain, e. g. a 801 Stirrer and a 805 Dosimat.

A **700 Dosino** or a **685 Dosimat** must be connected directly to the **MSB socket** of the Sample Processor.




Fig. 14 MSB connections

The instrument will only recognize automatically a connected dosing device when it is switched on or after **RESET**.

## 2.3 Data transmission connections

### Connection cables

Only Metrohm cables should be used to connect a Metrohm Sample Processor to other instruments. Only these will guarantee interference-free data transmission.



**Note:**

*The Remote cables for Metrohm Sample Processors and Sample Processors have markings at each end of the cable to indicate the instrument that the particular connector is intended to be used with and the connection into which it is to be plugged.*

Example:




Fig. 15 Remote cable

The Sample Processor must be switched off before peripheral devices are connected, as otherwise the instruments could be damaged.


### 2.3.1 Remote connections




*In order to connect Remote cables to Titrandos and the 780/781 pH/Ion Meter an 8.2148.010 Remote Box is required as an additional adapter; this is connected to one of the MSB sockets (MSB1...4), see above.*

#### Sample Processor — Titrimo / Titrimo

with standard cable (and 6.2148.010 Remote Box)



Cable 6.2141.020



Cable 6.2141.020

Titrimos and Titrandos can be connected to a Sample Processor in the same way and can also be controlled.

Control commands:

```
CTL:Rm : START Device1  starts Titrimo/Titrando
CTL:Rm : *****1      "
```


Final scan:

```
SCN:Rm :          End1  waits for end of titration (EOD-impulse)
SCN:Rm :          ****1***      "
SCN:Rm :          Ready1  waits until Titrimo/Titrando is ready
SCN:Rm :          *****1      "
```

---

### Sample Processor — pH / Ion Meter / Conductometer (692/712/713)

---



Cable 6.2141.020

Control commands:

```
CTL:Rm : START device1  starts instrument
CTL:Rm : *****1      "
CTL:Rm : METER Mode pH  switches to pH measurement
CTL:Rm : *****0001*   "          (not for 712)
CTL:Rm : METER Mode T   switches to temp. measurement
CTL:Rm : *****0010*   "          (not for 712)
CTL:Rm : METER Mode U   switches to mV measurement
CTL:Rm : *****0011*   "          (not for 712)
CTL:Rm : METER Mode I   switches to Ipol (mV measurement)
CTL:Rm : *****0100*   "          (not for 712)
CTL:Rm : METER Mode C   switches to Conc measurement
CTL:Rm : *****1000*   "          (for 692/781 only)
CTL:Rm : METER Ca1 pH   switches to pH calibration
CTL:Rm : *****0101*   "          (not for 712)
CTL:Rm : METER Ca1 C    switches to Conc calibration
CTL:Rm : *****1001   "          (for 692/781 only)
```

```

CTL:Rm :    METER enter    simulates <ENTER> key
CTL:Rm :    *****1111*    "    (not for 712)
    
```

Final scan:


```

SCN:Rm :    End1    waits for end of measurement
SCN:Rm :    ****1***    "
    
```

---

**Sample Processor — 780/781 pH / Ion Meter**

---




Control commands: see above (692/712/713 pH/Ion meter/ Conductometer).

---

**Sample Processor — 781 Ion Meter — 665 Dosimat**

---

For ion measurements with calibration and standard addition.



The 781 Ion Meter automatically controls Stirrer T1 of the Sample Processor.

Start commands:

In principle these are the same commands as those given in the previous section.

Final scan:


In principle these are the same commands as those given in the previous section, with the addition of:

```
SCN:Rm :      EndMeter waits for end impulse from 781
SCN:Rm :      ***1***      "
```

---

**Sample Processor— 2 x Titrimo/Titrandos**


---

Control commands:


```
CTL:Rm :  START Device1  starts Titrimo/Titrando 1
CTL:Rm :  *****1      "
CTL:Rm :  START Device2  starts Titrimo/Titrando 2
CTL:Rm :  *****1***** "
CTL:Rm :  START Device*  starts both Titrimo/Titrandos
CTL:Rm :  *****1***** "
```

Final scan:

```
SCN:Rm :      End1      waits for end of titration Titrimo/Titrando 1
SCN:Rm :      ****1***  "
SCN:Rm :      End2      waits for end of titration Titrimo/Titrando 2
SCN:Rm :      *1*****  "
SCN:Rm :      Ready1    waits until Titrimo/Titrando 1 is ready
SCN:Rm :      *****1  "
SCN:Rm :      Ready2    waits until Titrimo/Titrando 2 is ready
SCN:Rm :      **1*****  "
SCN:Rm :      Ready*    waits until both Titrimo/Titrandos are ready
SCN:Rm :      **1*****  "
```

The combined operation of the Titrimo/Titrando/pH meter can also be carried out without any problems by the same means.

**Sample Processor — Titrino/Titrando/pH Meter — 665/725 Dosimat**



Cable 6.2141.040

Control commands:

```


CTL:Rm : START Device1  starts Titrino
CTL:Rm : *****1      "
CTL:Rm : START Dos1    starts Dosimat 1
CTL:Rm : *****1***** "
    
```

Final scan:

```

SCN:Rm : End1  waits for end of titration (impulse)
SCN:Rm : ****1*** "
SCN:Rm : Ready1  waits until Titrino is ready
SCN:Rm : *****1 "
    
```

**Sample Processor — Titrino/Titrando/pH Meter — 2x 665/725 Dosimat**



Cable 6.2141.050

Control commands:

```

CTL:Rm : START Device1  starts Titrino
CTL:Rm : *****1      "
CTL:Rm : START Dos1    starts Dosimat 1
    
```


```
CTL:Rm : *****1*****           "  
CTL:Rm :      START Dos2      starts Dosimat 2  
CTL:Rm : *****1*****           "  
CTL:Rm :      START Dos*      starts Dosimat 1 and 2  
CTL:Rm : *****1*1*****        "
```

Final scan:

```
SCN:Rm :      End1      waits for end of titration (impulse)  
SCN:Rm :      ****1***           "  
SCN:Rm :      Ready1     waits until Titrino is ready  
SCN:Rm :      *****1           "
```

**2.3.2 Serial connections (RS232)**

Many different connection to the RS232 interface are possible. As well as all Metrohm instruments which understand the Metrohm Remote language (see 'Technical Reference 8.789.1033'), a printer (requirement: serial interface or parallel/serial converter) or a personal computer can be connected. Instruments from other manufacturers that have a serial RS232 interface can also be connected.



*Fig. 16 RS232 connections*

A precondition for correct data transmission is that the transmission parameters are set correctly, these must correspond to the interface settings of the connected instrument (see next page).

The following section provides information about the settings and cables required for the connection of a printer.

### 2.3.3 Connecting a printer

A printer with a serial or parallel interface can be connected to the RS232 interface for printing out reports (e. g. parameter report).

Printers with the following printer emulation can be connected:

<b>IBM</b>	IBM Proprinter and printers with IBM emulation
<b>Epson</b>	EPSON printers and printers with EPSON emulation
<b>Seiko</b>	Seiko printer DPU-411
<b>Citizen</b>	Citizen printer IDP560 RS
<b>HP</b>	HP printers and printers with HP PCL3 emulation

If you wish to connect a different printer then you must ensure that it can emulate one of the printer modes supported by the Sample Processor.




Printers with a parallel interface require a serial/parallel converter (e.g. 2.145.0300) and the 6.2125.020 Cable.



*The Sample Processor must be switched off before a printer is connected to the RS232 interface!*


The interface parameters are entered in the configuration menu under "**>RS232 settings**".

The following table provided information about connecting some selected printers.

Printer	Cable	789/789 settings	Printer settings
Custom DP40-S4N	6.2134.110	Send to: Citizen Baud rate: 9600 Data bit: 8 Stop bit: 1 Parity: none Handshake: HWs	IDP-560 EMULATION FONT MAP = ENGLAND PRINT = REVERSE LITTLE CR CODE = VOID CR AFTER B. FULL = VOID CR ON B. EMPTY = VALID BUFFER 1K BYTE BAUD RATE = 9600 PROTOCOL = 8,N,1 FLOW CONTROL CTS-RTS
Seiko DPU-414	6.2134.110	Send to: Seiko Baud rate: 9600 Data bit: 8 Stop bit: 1 Parity: none Handshake: HWs	none
Citizen iDP562 RS	6.2134.050	Send to: Citizen Baud rate: 9600 Data bit: 8 Stop bit: 1 Parity: none Handshake: HWs	ON  SSW1
Epson LX-300+	6.2134.050	As above, or Send to: Epson	see printer manual
HP Desk Jet with serial in- terface	6.2134.050	Send to: HP Baud rate: 9600 Data bit: 8 Stop bit: 1 Parity: none Handshake: HWs	A:  A4 Papier B: 
HP Desk Jet with parallel interface	2.145.0330 RS232/ Parallel Converter	Send to: HP Baud rate: 9600 Data bit: 8 Stop bit: 1 Parity: none Handshake: HWs	see printer manual

## 2.4 Connecting up accessories

### 2.4.1 Connecting rinsing and aspiration equipment



Both the existing built-in membrane pumps and external pumps can be used for rinsing the titration head and aspirating off samples that have been processed. The necessary tubing is supplied with the corresponding Metrohm instruments and is easily mounted on the Sample Processor.

Attach the tubing for the spray nozzles and the aspiration tip according to the diagram on the left and lead them through the guide chain.

*Fig. 17 Distributor piece*

to waste canister

#### **Spray nozzles**

Spray nozzles (6.2740.030) can be used for rinsing electrodes and dosing tips.

The spray nozzles are inserted in the oblique openings of the titration head. These are connected to the distributor block on the guide chain with tubing (6.1805.420). The three distributor block openings (M6) provided for them must all be occupied. Any unoccupied openings must be sealed with an M6 threaded stopper (6.1446.040) to prevent liquid escaping during the rinsing process.

The tubing is attached to the guide chain. It may be necessary to remove the stud of the first link in the chain.

**Aspiration tip**

In the macro titration head the M8 aspiration tip (6.1543.170) is placed in a vertical opening (SGJ9). In the micro titration head a vertical opening is also used. The tip is connected with a piece of PTFE tubing (6.1805.510) and a 6.1820.030 Screw connector to Pump 2 (marking on valve '↓').

**Rinsing process**

In an ideal rinsing process (with two pumps) the sample solution is first aspirated off. The electrode is then rinsed in an empty beaker and the rinsing solution is aspirated off at the same time. A short overrun of the aspiration pump ensures that the sample vessel (or titration cell) is emptied completely.

Sequence:

<b>PUMP</b>	<b>1.2</b>	<b>:</b>	<b>10 s</b>
<b>PUMP</b>	<b>1.*</b>	<b>:</b>	<b>3 s</b>
<b>PUMP</b>	<b>1.2</b>	<b>:</b>	<b>3 s</b>

**Chemical resistance**


The pump heads of the membrane pumps, the connection nipples and the union nuts are made of PVDF. PVDF is resistant to many chemicals.

The following should **not** be used: acetone, acetic anhydride or dimethylformamide (DMF). The membrane, the valve disks and valve components are made of PTFE and are resistant to almost all chemicals.

If your sample contains solids (e.g. silver chloride) or sticky substances then, instead of the membrane pump built into the Sample Processor, you should use the 772 Pump Unit.

The polyethylene canisters (6.1621.000) used as rinsing or waste canisters are (only) suitable for aqueous solutions.

**2.4.2 Attaching the tubing**



For attaching the tubing the individual chain links can be opened with a screwdriver or similar object.

Mind the markings on the chain links.

*Fig. 18 Attaching the tubing*

### 2.4.3 Installing the titration accessories

#### Macro titration head

6.1458.010





Fig. 19 Macro titration head

- |          |  |
|----------|--|
| <b>A</b> | Electrode                              |
| <b>B</b> | Aspiration tip M8<br>6.1543.170        |
| <b>C</b> | Buret tip                              |
| <b>D</b> | 802 Stirrer<br>2.802.0020              |
| <b>E</b> | Spray nozzle M6<br>6.2740.020          |
| <b>F</b> | SGJ sleeve SGJ14<br>6.1236.020         |
| <b>G</b> | Threaded stopper<br>M8 4.658.0180      |
| <b>H</b> | Guide sleeve SGJ9<br>6.2709.070        |
| <b>I</b> | Stirrer shaft<br>6.1909.020/6.1909.040 |

#### Micro titration head

6.1458.020





Fig. 20 Micro titration head

Only special micro-electrodes can be used in the micro titration head, see list of accessories, p. 161.

Use a 741 Magnetic Stirrer instead of a Rod Stirrer.


#### Information about the macro titration head

The SGJ14 opening indicated by an arrow has been bored slightly obliquely so that a rod stirrer or electrode can be centered in a narrow titration vessel.

**Robotic arm with titration head**

6.1462.050


The titration head can be equipped with one 802 Rod Stirrer, up to two micro electrodes, three buret tips and three spray nozzles.



<b>A</b>	Micro electrode
<b>B</b>	Aspiration tip M8 6.1543.170
<b>C</b>	Buret tip (drawn back)
<b>D</b>	802 Rod Stirrer 2.802.0020
<b>E</b>	Spray nozzle M6 6.2740.030
<b>F</b>	Threaded stopper M8 4.658.0180 for micro electrode or buret tip
<b>G</b>	Guide sleeve SGJ9 6.2709.070
<b>H</b>	Stirrer shaft 6.1909.030

*Fig. 21 Robotic arm with titration accessories*


### 2.4.4 Assembling pipetting accessories



*Fig. 22 Transfer robotic arm with pipetting accessories*

Before pipetting tubing **A** is used its tip should be drawn out as fine as possible (except for high-viscosity samples).

**2.4.5 Installing a double-hollow needle**



For sample transfer from closed containers by means of a double-hollow needle the 6.1462.090 right-swing robotic arm with luer lock adapter can be installed.

- A** FEP tubing M6 (e.g. 6.1562.110 pipetting tubing 2 mL)

---

- B** M6 adapter

---

- C** Swing head with luer lock adapter (6.1462.090)

---

- D** Luer lock adapter


---

- E** Double-hollow needle with luer lock connector(6.2816.100)

---

*The parts listed above are contained in 6.5619.000 sample transfer accessory set.*

**Fig. 23 Transfer robotic arm with double-hollow needle**




**2.4.6 Spray nozzles**

The 6.2740.020 spray nozzles are installed as shown in the diagram. Their insertion height is adjustable so that an optimal rinsing effect can be achieved no matter how the head is equipped.

**Fig. 24 Spray nozzles**

### 2.4.7 Magnetic stirrer



A 741 Magnetic Stirrer can be placed in any position on the stirrer guide below the sample rack and moved about as required.


It is connected to the rear panel of the tower (MSB connector).

*Fig. 25 741 Magnetic stirrer*

The 801 Magnetic Stirrer or 804 Titration Stand is connected to an MSB socket on the rear of the Sample Processor chassis.

### 2.4.8 Mounting the stand support

An external titration or measuring cell can be mounted on the Sample Processor by using the 6.2001.070 Stand Support. The stand support is attached to stirrer guide of the Sample Processor and fixed in position with the two socket head screws.



- The height and side at which the cross strut is fixed are adjustable. It is fixed in position on the rear panel of the tower with the screw provided.
- If it is to be mounted on a 789 Robotic Sample Processor XL then the support rod must be moved back by 3 cm. A second mounting hole is provided in the base plate for this reason.

*Fig. 26 Stand support*

### 2.4.9 Installing the collection trough

Two different sizes of collection trough are available as optional accessories for the Sample Processors. They are intended to prevent any escaping liquid from spreading over the bench or penetrating the instrument housing. The two models are:

- 6.2711.060 for 778 Sample Processor with 42 cm racks
- 6.2711.070 for 789 Robotic Sample Processor XL with 48 cm racks




Fig. 27 Collection trough

Attach the supplied tubing to the collection trough outlet. Remove the sample rack and place the collection trough on the turntable so that it lies on the stirrer guide of the Sample Processor. The other end of the tubing should be located in a sink or waste container.

### 2.4.10 Sample racks





Fig. 28 Attaching a sample rack

Sample racks are attached so that the two guide pins are located in the openings provided for them in the base of the rack.


After a sample rack has been placed in position it should be initialized with **<RACK>** so that the magnet code of the rack can be read in. This is only possible when the rack is in the initial position.



Sample racks for the 789 Robotic Sample Processor XL (48 cm dia.) have a handle with fixing screws. This handle is used to fix the rack on the Sample Processor turntable by turning it clockwise.

*Fig. 29 Sample rack for XL models*

### 2.4.11 Mounting the safety shield




Attachment screw

It is essential that the safety shield is mounted before the Sample Processor is used for the first time. Not only does it provide an effective protection against splashes, but it also considerably reduces the risk of injury.

Various versions of the safety shield are available.

The safety shield is mounted on the Sample Processor tower with the screws provided. The position of the safety shield can be modified to suit the height of the sample beakers.

*Fig. 30 Mounting the safety shield*



*Fig. 31 Safety shield for Swing Heads*

## 2.5 Start-up

Some instrument settings must be made on the Sample Processor before it is used for the first time. These are:

- Dialog language
- Counter setting of running time meter
- Tower configuration
- Robotic arm settings (optional)
- Rack definitions
- Dosing unit configurations
- RS232 interface parameters

All these settings are found in the configuration menu that you can open with the **<CONFIG>** key.

Under **>auxiliaries:**

- Set the dialog language to English, see p. 57.
- Enter a warning limit under **>>running time meter**, see p. 58. Recommendation: 1000 h

Under **>tower 1** and **>tower 2:**

- Limit the max. stroke path, i.e. set the lowest possible lift position, see p. 58. Recommendation: 160...200 mm (Can be set with the **LEARN** function)
- If a 786 Swing Head is installed then, under **>>swing head 1**, set the swing position for the robotic arm, select a swing angle for at least one external position and set its working height, see p. 58ff. All this can be set using the **LEARN** function.

Under **>rack definitions:**

- Load the rack data of the attached sample rack and set the various lift positions (working position, shifting position, rinsing position, etc. for each lift) for processing the samples, see p. 61ff. This can all be done by using the **LEARN** function.
- Select a beaker sensor. For single row samples racks the optical sensor of the tower should be selected.
- Under **>>>special positions** you may want to define one or more special beaker positions. The rack position for the special beaker (e.g. rinsing beaker, conditioning beaker) must be entered and the working position must be set for each lift (with the **LEARN** function), see p. 61ff.
- The rack definitions must be saved and then the sample rack must be initialized with the **<RACK>** key.

Under **>Dosing units**:

- Select dosing drive 1...3
- For each port (1...4) of the corresponding Dosing unit set the max. dosing rate and enter the tubing dimensions, see p. 63ff. These settings are required for preparing the Dosing unit (with the **PREP** function of the DOS command). The dosing rate must be reduced for viscous liquids.

Under **>RS232 settings**:

If a PC control software (e. g. TiNet or Tiamo) is used then it may be necessary to match the transmission parameters of the serial interface to those of the software, see p. 65ff. However, the default settings are suitable for data transmission.



# 3 Operation

As well as the basic functions that can be carried out manually, this section also describes the configuration of a Sample Processor. Detailed descriptions of the way that run sequences are created and how methods are managed are given.

## 3.1 Operating principles


### 3.1.1 Display

The keypad (SC Controller) display consists of two lines each with 24 characters.


The first line is used as the title line and shows the current method together with the number of processed samples. In the editing mode the menu titles are shown here.

The second line is used as the status line and shows specific activities that depend on the operating condition. It is used as the input line in the editing mode.


#### Normal state



#### Method sequence



#### Editing mode



### 3.1.2 Instrument dialog

The instrument dialog of the Sample Processor is arranged in menu levels in which the following rules apply:

#### Main menu

The **<CONFIG>**, **<PARAM>** and **<USER METHOD>** keys of the Sample Processor open a main menu whose thematically arranged submenus are accessed by repeatedly pressing this key or with **<↓>**. The name of the main menu appears in the first line.

#### Submenu

Each submenu has its own title which is indicated by ">" and appears in the lower line of the display. From the title you can use **<ENTER>** to access the individual questions with which the most important instrument settings can be altered. The first line of the display always shows the name of the active submenu.

Navigation in the menus is by the cursor keys; with **<HOME>** you can access the first line of the menu and with **<END>** the last one.

**<QUIT>** exits the active menu and accesses either the next superior menu or the normal operating condition.

**<ENTER>** always opens a submenu or confirms the data input of the lowest menu level.

#### Input lines

For input lines without ":" the values can be entered by the numerical keys. The set value is accepted with **<ENTER>** and the next line appears.

For input lines with ":" predefined values can be selected with the **<SELECT>** key. **<ENTER>** accepts the set value and the next line appears.

Depending on the parameter, **<CLEAR>** is used to set the initial value. The **<CLEAR>** key is also used to reject incorrectly entered values.

**<QUIT>** exits the questions and returns to the submenu.

The following schematic diagram shows the instrument dialog arrangement.






Fig. 32 Dialog arrangement

### 3.1.3 Data input


#### Input line

Either one or two parameters can be entered in a menu line or a sequence. A blinking block cursor indicates where a parameter can be entered.




Switching between the parameters is carried out with the arrow keys <→> and <←>. <ENTER> moves the cursor automatically to the right, <QUIT> moves it to the left.

#### <Select> (Selection drum)



Data can usually be entered directly by the numerical keys of the keypad. For entries marked by a colon the <SELECT> key can be used to display a selection of entries. This selection is cyclically arranged like a revolving drum.


Example:



### 3.1.4 Text input

The text editor can be used wherever provision is made for the entry of a text.

Numbers can be entered directly via the keypad.



The "<" or ">" keys open the text editor. With "<" an existing character string is deleted and the text cursor moves to the left-hand end of the input field. With ">" an existing character string is retained, the text cursor moves to the last character of the existing text.

A character string is displayed that consists of all the characters that can be entered in alphabetic order. The currently selected character always blinks (text cursor).

#### Character selection


The "<" and ">" keys move the character string made up of all the selectable characters (upper and lower case letters, numbers and special characters, arranged alphabetically) in the selected direction past the text cursor. Pressing either of these keys once moves the character string in the corresponding direction by one position. Rapid character string movement is achieved by keeping the key pressed down.

#### Confirming character selection




The **<ENTER>** key appends the character currently under the text cursor to the existing text line. When the complete width of the text input field has been filled the text input mode is exited and the text line is accepted with **<ENTER>**.

#### Delete character




The **<CLEAR>** key deletes the last character of an existing text line. The text cursor automatically moves one character to the left.

#### Ending text input



**<QUIT>** exits the text input mode. The displayed text line can then be accepted with **<ENTER>** or rejected by pressing **<QUIT>** a second time.


**Arrangement:**



*Fig. 33 Text input*

Fig. 33 shows how a character string can be entered, e.g. for naming a method. Text input is concluded with **<QUIT>**. The whole of the entered character string is then shown and can be accepted with **<ENTER>** or rejected with **<QUIT>**.

## 3.2 Configuration



Main menu:

configuration >auxiliaries
configuration >tower 1
configuration >tower 2
configuration >rack definitions
configuration >Dosing units
configuration >RS232 settings

- Open submenu with **<ENTER>**
- Use **<↑>** or **<↓>** to move up and down by one menu item
- Use **<HOME>** or **<END>** to move to the first or last menu item
- Use **<QUIT>** to return to operating condition

### 3.2.1 Miscellaneous

configuration >auxiliaries	Submenu for basic settings Open the submenu with <b>&lt;ENTER&gt;</b>
-------------------------------	--

Use <b>&lt;QUIT&gt;</b> for next higher level	>auxiliaries dialog: <b>english</b>	Select dialog language
	<b>english, deutsch, français, español</b>	

>auxiliaries display contrast <b>3</b>	Set display contrast
<b>0...3...7</b>	0 = low contrast 7 = high contrast

>auxiliaries beeper: <b>ON</b>	Beep for warnings on or off
<b>ON,OFF</b>	

>auxiliaries external Start: <b>off</b>	Switch on external start via Remote line (Input 7)
<b>on, off</b>	

Normally the Sample Processor takes over the complete run control for a sample series as the "Master". This may be unwanted, particularly if an automated system has been set up that includes instruments that are not from Metrohm. This is why the Sample Processor can also be started from any external instrument via the Remote interface.

With **external Start: ON** the run sequences of the Sample Processor will

be started as soon as **input line 7** is activated at the Remote interface (low level).

The run sequences of the Sample Processor will be stopped when **input line 6** is activated (low level). The technical details can be found in the **8.789.1033 Technical Reference** of the Sample Processors.

```
>auxiliaries
>>running time meter
```

Submenu for running time meter, open with <ENTER>

Shows the elapsed time

```
>>running time meter
elapsed time      0.0 h
                  0...9999 h
```

Warning limit for running time meter

```
>>running time meter
warning          OFF h
                OFF, 0...9999 hours
```

The running time meter is used to support the regular maintenance of the Sample Processor. If a warning limit is entered then, when this period has elapsed, a message will appear indicating that maintenance is required.

```
>auxiliaries
device label
```

8 ASCII characters

Name of the instrument for identification

```
>auxiliaries
program          5.789.0011
```

read only

Program version

### 3.2.2 Tower settings

Sample Processor models with two towers can have each tower separately equipped and separately configured. For models with a single tower only the configuration menu for Tower 1 can be accessed.

configuration >tower 1	Submenu for basic settings for Tower 1 Open the submenu with <ENTER>
---------------------------	---

>tower 1 max. stroke path	235 mm
0...235 mm	

Max. lift range for Lift 1



This max. stroke path setting is **important for safety reasons**. Correct entry of this value can prevent electrode breakage, as the titration head cannot move any lower than the position entered here (0 mm = lower stop of lift).

>tower 1 min. beak. radius	* mm
*, 1...100 mm	

Sets the minimum required beaker radius for processing at Tower 1  
\* = any value



This is also a safety setting. In order to prevent a fully equipped, wide titration head from trying to enter a narrow sample beaker a limit for the minimum necessary beaker radius can be entered here. During a method sequence the value entered here will be compared with the value given for the effective beaker radius in the rack table before the lift is lowered and, if necessary, an error message will be produced.  
Entering \* prevents a comparison from being made.

>tower 1 >>swing head 1	
----------------------------	--

Submenu Swing Head 1  
Open with <ENTER>

If a 786 Swing Head equipped with a robotic arm is mounted on the Sample Processor then it is essential that the correct settings for the mounted robotic arm are entered, as otherwise the instrument will not be able to position it exactly.

Each of the following settings can be interactively defined with the **LEARN function**. Press the <LEARN> key and then move the lift or robotic arm with the arrow keys <↓> and <↑> or <←> and <→>. The set values can be accepted with <ENTER> and corrected later.

Setting the rinsing height for external positions \*)

>>swing head 1 rinse position	0 mm
0...235 mm	

Setting the swing height for external positions \*)

>>swing head 1 shift position	0 mm
0...235 mm	

Swing angle for external position 1 \*)

```
>>swing head 1
external pos.1      84.00°
```

Offset...84.00...max. angle+offset  
0° = points to rack center

The swing angle for the external positions is entered as an absolute angle. The input limits are defined by the offset of the Swing Head drive (approx. 9°) and the set max. swing angle range.

Setting the working height for external position 1 \*)

```
>>swing head 1
work position1      0 mm
```

0...235 mm

Swing angle for external position 2 \*)

```
>>Robotic arm 1
external pos.2      84.00°
```

... up to external pos. 4

\* can be set with LEARN function

```
configuration
>tower 2
```

```
Submenu for basic settings for Tower 2
Open the submenu with <ENTER>
```

... see Tower 1

### 3.2.3 Rack definitions

<b>configuration</b> <b>&gt;rack definitions</b>	Submenu for defining individual racks Open the submenu with <ENTER>
---	--

Use <QUIT>  
for next higher  
level

```
>rack definitions
>>recall rack
```

Loads rack definitions

In order to make a modification to the definition of a sample rack the rack data must first be loaded. The data of the Metrohm standard racks is stored under their ordering number.

```
>recall rack
name: 6.2041.310
```

Name of the rack to be loaded

10 ASCII characters

The <SELECT> key can be used to make a selection from the saved rack data. The rack data is loaded with <ENTER>. The first selection to appear is the rack name of the currently attached rack.

```
>recall rack
code 000001
```

Magnet code of the rack  
See table on p. 91.

000001...111111

The magnet code is used for the unambiguous identification of the rack. The magnet code is recognized during the initialization of the rack. This is why the <RACK> key should be pressed whenever a rack is changed.

```
>recall rack
work position T1 0 mm
```

Working position for sample positions at Tower 1 \*)

0...235 mm

(in mm from upper stop point)

```
>recall rack
work position T2 0 mm
```

Working position for sample positions at Tower 2 \*)

0...235 mm

(in mm from upper stop point)

```
>recall rack
rinse position T1 0 mm
```

Rinsing position for sample positions at Tower 1 \*)

0...235 mm

(in mm from upper stop point)

```
>recall rack
rinse position T2 0 mm
```

Rinsing position for sample positions at Tower 2 \*)

0...235 mm

(in mm from upper stop point)

```
>recall rack
shift position T1 0 mm
```

Shifting position for sample positions at Tower 1 \*)

0...235 mm

(in mm from upper stop point)

<pre>&gt;recall rack shift position T2      0 mm                       0...325 mm</pre>	Shifting position for sample positions at Tower 2 *) (in mm from upper stop point)
<pre>&gt;recall rack special pos. T1      0 mm                     0...325 mm</pre>	Special position for sample positions at Tower 1 *) (in mm from upper stop point)
<pre>&gt;recall rack special pos. T2      0 mm                     0...325 mm</pre>	Special position for sample positions at Tower 2 *) (in mm from upper stop point)
<pre>&gt;recall rack beaker radius        * mm                     *, 1...100 mm</pre>	Effective beaker radius for the sample positions on the rack * = any value

This setting is required for an automatic check of the beaker radius, see p. 58.

<pre>&gt;recall rack beaker sensor:      Tower                    Tower, SwingH, OFF</pre>	Selection of the beaker sensor
--	--------------------------------

The beaker sensor checks the presence of a sample beaker when executing a MOVE command.

<pre>&gt;recall rack rack offset         0.00°                    -5.00...0.00...5.00°</pre>	Correction for the rotation angle of a sample rack *)
--	---

\*) All the above lift positions and the rack offset can be set by using the **<LEARN>** function.

### Special beaker positions

<pre>&gt;&gt;recall rack &gt;&gt;&gt;special positions</pre>	Submenu special beaker positions Open with <ENTER>
--	---

*Use <QUIT>  
for next higher  
level*

Special beaker positions are predefined places on a sample rack that are not treated as sample positions. They can be occupied by rinsing beakers or conditioning beakers and selectively addressed during a run sequence. Up to 16 special beaker positions can be defined per rack. For each special beaker the work position of the lift (tower 1 and 2) and the beaker radius can be set, see above.

Selection of special beaker

<pre>&gt;&gt;&gt;special positions special beaker      1                    1...16</pre>
--

Rack position of the special beaker

```
>>special positions      1
rack position            0
```

0...number Pos.  
0 = not defined

etc. up to special beaker 16

```
>rack definitions
>>store rack
```

Save rack definitions

In order to store modifications to the definition of a sample rack the submenu **>>store rack** is selected.

```
>store rack
name:          6.2041.310
```

Name of rack

10 ASCII characters

With the **<SELECT>** key you can select one of the existing rack names. The rack data is stored with **<ENTER>**. Any rack name can be used. The input of any rack name can be made directly via the numerical keys or in the text input mode, see p. 55. Selection of the alphanumerical characters with the **<** and **>** keys or **<PRINT>** and **<RACK>**.

```
>rack definitions
>>delete rack
```

Delete rack definitions

The submenu **>>delete rack** must be selected if a rack definition is to be deleted.

```
>delete rack
name:          6.2041.310
```

Name of the rack

The **<SELECT>** key can be used to select one of the existing rack names. **<ENTER>** confirms the selection. **<QUIT>** stops the deletion process. Before the deletion a question appears.

```
>rack definitions
delete 6.2041.310 ?
```

**<ENTER>** confirms the deletion. **<QUIT>** stops the deletion process.



When the submenu **'>rack definitions'** is exited without saving modified data a request appears about saving the rack data.

```
>rack definitions
overwrite 6.2041.310 ?
```

Confirm the question (store the rack definitions) by pressing the **<ENTER>** key. Reject storage with **<QUIT>**.

### 3.2.4 Dosing units

Metrohm Exchange units (with 685 and 805 Dosimats as the dosing drive) and Dosing units for the Dosino 700 and 800 systems can be used with the Sample Processors. The following settings are used for preparing the Dosing units (**PREP** function). The tubing dimensions (length and diameter) are used for calculating the rinsing volumes.

If Exchange units are used then only the dosing rate is effective.

<b>configuration</b> <b>&gt;Dosing units</b>	Submenu for the Dosing unit settings Open the submenu with <ENTER>
---	---

Use <QUIT>  
for next higher  
level

<b>&gt;Dosing units</b>	<b>1</b>
<b>Dosing unit</b>	<b>1</b>
	<b>1...3</b>

Selection of Dosing unit or dosing device connection

The Dosing unit selection must be confirmed with <ENTER>. It will then be shown in the first menu line.

<b>&gt;Dosing units</b>	<b>1</b>
<b>max. rate 1</b>	<b>160 ml/min</b>
	<b>0.01...160 ml/min</b>

Max. dosing and filling rate for Port 1  
(depends on cylinder size)

<b>&gt;Dosing units</b>	<b>1</b>
<b>tube length 1</b>	<b>1000 mm</b>
	<b>0...1000...30000 mm</b>

Length of tubing at Dosino Port 1

<b>&gt;Dosing units</b>	<b>1</b>
<b>tube diameter 1</b>	<b>2 mm</b>
	<b>0.1...2...20 mm</b>

Diameter of tubing at Dosino Port 1

<b>&gt;Dosing units</b>	<b>1</b>
<b>max. rate 2</b>	<b>160 ml/min</b>
	<b>0.01...160 ml/min</b>

Max. dosing and filling rate for Port 2  
(depends on cylinder size)

<b>&gt;Dosing units</b>	<b>1</b>
<b>..... to port 4</b>	

Input of tubing parameters for all four ports of a Dosing unit.

### 3.2.5 RS232 settings

The correct settings of the interface parameters of the serial RS232 interface is essential for the perfect functioning of data transmission to and from the Sample Processor. This includes the control of the instrument by using a PC software such as **TiNet** or **Tiamo**, backing up of configurations or methods with the **VESUV** software from Metrohm or printing out a parameter report on a connected printer.

<b>configuration</b> <b>&gt;RS232 settings</b>	Submenu for serial interface settings Open the submenu with <ENTER>
---	--

Use <QUIT>  
for next higher  
level

<b>&gt;RS232 settings</b>	
<b>baud rate:</b>	<b>9600</b>
	300, 600, 1200, 2400, 4800, <b>9600</b> , 19200

Transmission rate  
in baud

<b>&gt;RS232 settings</b>	
<b>data bit:</b>	<b>8</b>
	7, <b>8</b>

Number of data bits

<b>&gt;RS232 settings</b>	
<b>stop bit:</b>	<b>1</b>
	1, 2

Number of stop bits

<b>&gt;RS232 settings</b>	
<b>parity:</b>	<b>none</b>
	even, odd, <b>none</b>

Parity selection

<b>&gt;RS232 settings</b>	
<b>handshake:</b>	<b>HWs</b>
	HWs, SWchar, SWline, none

Handshake selection

<b>&gt;RS232 settings</b>	
<b>character set:</b>	<b>IBM</b>
	IBM, HP, Epson, Seiko, Citizen

Character set for printer and PC  
(printer emulation)

The settings for the printers recommended by Metrohm are listed on p. 37ff. If your printer is not shown there then we recommend the setting "Epson". It may be necessary to consult the printer manual. For data transmission using personal computers you should select "IBM".

<b>&gt;RS232 settings</b>	
<b>RS control:</b>	<b>ON</b>
	ON, OFF


Switch data reception on/off

If Remote control is switched off then no data will be received, but reports can still be printed out.

## 3.3 Composition of a method

All the settings of the parameter menu form a method and can be stored as such.

### 3.3.1 Run sequences and method parameters



Main menu:

parameters	
number of samples:	rack

Number of samples to be processed

1...999,  
rack,\*

rack = one sample rack rotation  
\* = infinite

In the **rack** setting all the sample positions of the current rack will be processed (max. no. of rack positions minus number of special beakers), with only the positions occupied by sample beakers being counted. It is important that the Sample Processor can recognize the rack. This is only possible when the rack is in initial position. We recommend that each time a rack is changed the Sample Processor is initialized with **<RACK>**.

parameters >start sequence
parameters >sample sequence
parameters >final sequence
parameters >changer settings
parameters >stirring rate
parameters >Dosing unit def.
parameters >timeout settings
parameters >manual stop

Open submenu with **<ENTER>**

Use **<↑>** or **<↓>** to move up and down by one menu item

Use **<HOME>** or **<END>** to reach the first and last menu item

Use **<QUIT>** to return to normal operating condition

#### Submenus:

In each of the submenus **>start sequence**, **>sample sequence** and **>final sequence** up to 99 command lines can be entered as a run sequence. The commands can be entered directly via the keypad. The command keys on the right-hand side of the keypad are available.

parameters >start sequence	Commands for the start sequence of the sample series Open the submenu with <b>&lt;ENTER&gt;</b>
-------------------------------	--

The start sequence is **carried out once** at the start of a sample series. For example, this can be used for rinsing or conditioning an electrode.


parameters >sample sequence	Commands for the processing of each sample Open the submenu with <ENTER>
--------------------------------	---

The sample sequence is carried out when **each individual sample** in a sample series is processed.

parameters >final sequence	Commands for the final sequence of the sample series Open the submenu with <ENTER>
-------------------------------	---

The final sequence is **carried out once** at the end of a sample series. For example, it can be used to move to a rinsing or conditioning beaker.

### Series operating



### Creating methods

In principle the same input rules apply for creating run sequences as for manual operation, i. e. after selecting a command and entering the necessary parameters the entry is confirmed with **<ENTER>**. The next command line is then selected and a new command can be entered.

For more comfortable parameter input the **LEARN** function is available for certain commands. It can be used to accept direct values that are set when a single command is carried out manually. Please refer to p. 112 for details.

The **TRACE** function can be used to execute each entered command line individually in the editing mode. See p. 114.

Navigation in a sequence is the same as in other menus. In addition, the **<INSERT>** and **<DELETE>** keys are available.

**<INSERT>** inserts a new command line **above the current line** of a sequence. It is automatically occupied with the **NOP** command, which has no function. The following lines move downward by one position.

**<DELETE>** deletes the current line of a sequence. The following lines move upward by one position.

### 3.3.2 Sample Processor settings

parameters >changer settings	Submenu for setting the Sample Processor functions Open the submenu with <ENTER>
---------------------------------	---

>changer settings rack name: *	The sample rack assigned to the method
-----------------------------------	--

\*, 10 ASCII characters \* = no particular rack

With this setting the user can be forced to use a particular rack for the selected method. If this is not required then **rack name: \*** must be selected.

>changer settings lift rate T1 25 mm/s	Vertical movement rate of Lift 1
---	----------------------------------

5...25 mm/s

>changer settings lift rate T2 25 mm/s	Vertical movement rate of Lift 2 (only for 2-tower models)
---	---

5...25 mm/s

>changer settings shift rate 20°/s	Speed at which the rack rotates in degrees per second
---------------------------------------	--

3...20

>changer settings shift direction: auto.	Direction of rotation of the sample rack
---	--

+, -, auto. auto. = the Sample Processor selects the shortest path for the rotation itself.

>changer settings rotat. increment: 5.00°	Increment angle for the command MOVE +/-rotate
--	---

0.00...5.00...270.00°

>changer settings swing rate T1 55°/s	Swing rate of the robotic arm at Tower 1 in degrees/second
--	---

10...55

>changer settings swing rate T2 55°/s	Swing rate of the robotic arm at Tower 2 in degrees/second (only for 2-tower models)
--	--

10...55

>changer settings swing increment 10.00°	Increment angle for the command MOVE +/-swing
---	--

0.00...10.00...180

>changer settings on beaker error: MOVE	Definition of reaction if beaker is missing (only applicable with beaker sensor on)
--	--

MOVE, Message

**MOVE** = the last action will be carried out again. The next position according to the current SAMPLE command will be selected.  
Message = The run will be interrupted and a warning displayed.

When a rack position is accessed a check is made whether a sample beaker is present at the selected tower. If the run is not to be interrupted when a sample beaker is missing then **MOVE** can be selected. This means that, if a sample beaker is missing, the **MOVE** command moves on to the next sample position. The next sample is selected according to the last **SAMPLE** command, i.e. for a previous **SAMPLE +2** command the sample will be looked for at the next-but-one rack position, etc.

If a special beaker is missing an error message will always be shown and the run will be interrupted.

### 3.3.3 Stirrer settings

parameters >stirring rates	Submenu for the stirrer settings Open the submenu with <ENTER>
-------------------------------	---

- |                                    |          |
|------------------------------------|----------|
| >stirring rates<br>stirrer tower 1 | <b>3</b> |
|------------------------------------|----------|

Stirring speed of the stirrer at Tower 1  
(stirrer models 802, 741, 722)
- |                                    |          |
|------------------------------------|----------|
| >stirring rates<br>stirrer tower 2 | <b>3</b> |
|------------------------------------|----------|

Stirring speed of the stirrer at Tower 2  
(stirrer models 802, 741, 722)
- |                                 |          |
|---------------------------------|----------|
| >stirring rates<br>stirrer MSB1 | <b>3</b> |
|---------------------------------|----------|

Stirring speed of MSB stirrer 1  
(stirrer models 801, 804)
- |                                 |          |
|---------------------------------|----------|
| >stirring rates<br>stirrer MSB2 | <b>3</b> |
|---------------------------------|----------|

Stirring speed of MSB stirrer 2  
(stirrer models 801, 804)
- |                                 |          |
|---------------------------------|----------|
| >stirring rates<br>stirrer MSB3 | <b>3</b> |
|---------------------------------|----------|

Stirring speed of MSB stirrer 3  
(stirrer models 801, 804)

### 3.3.4 Dosing device settings

Both Dosinos (models 700 and 800, with Dosing units) and Dosimats (models 685 and 805, with exchange units) can be used as dosing drives on a Sample Processor. However, free assignment of the ports for the various functions is only possible with Dosinos.

<b>parameters</b> >Dosing unit def.	Submenu for setting the Dosing units Open the submenu with <ENTER>
--	---

>Dosing unit def. dosing drive	1	Selection of the dosing device, or the connection
	1...3	

After entering the connection of the dosing device (see p. 8) and confirming its with <ENTER> the settings for the selected dosing device should then be made.

>dosing drive dos.rate	1 max. mL/min	Set the dosing rate
	0.01...160 mL/min, max.	

The maximum possible dosing rate depends on the size of the dosing cylinder (rule of thumb: cylinder volume x 3.3).

>dosing drive fill.rate	1 max. mL/min	Set the filling rate
	0.01...160 mL/min, max.	

The maximum possible filling rate depends on the size of the dosing cylinder (rule of thumb: cylinder volume x 3.3).

The following entries only apply to 700/800 Dosinos. Details about Dosinos and Dosing units can be found on p. 95ff.

>dosing drive dosing	1 port 1	Define standard dosing port 1
	1...4	

>dosing drive dosing 2	1 port 3	Define standard dosing port 2
	1...4	

>dosing drive filling	1 port 2	Define standard filling port
	1...2...4	

>dosing drive rinsing	1 port 2	Define standard rinsing inlet (if Dosing unit is exchanged)
	1...2...4	

>dosing drive	1	Define standard outlet for the prepara-
---------------	---	---

<pre>preparation          port 1                     1...4</pre>	tion cycle
<pre>&gt;dosing drive      1 drain              port 4                     1...4</pre>	Define standard air inlet for emptying the Dosing unit
<pre>&gt;dosing drive      1 cock direction:    auto                     auto                     not over                     desc.                     asc.</pre>	Direction of rotation of Dosino stopcock = cock direction according to shortest path = never turn beyond port given below = cock rotation always in descending direction = cock rotation always in ascending direction
<pre>&gt;dosing drive      1 not over port      4                     1...4</pre>	Protected port, that it not to be moved over during a stopcock rotation (see above)

### 3.3.5 Behavior during timeout

Metrohm Sample Processors are designed to communicate with other instruments. In particular, this includes the coordination between the Sample Processor and devices connected to it. In the method sequence the **SCAN** command can be used to check whether an instrument is ready or to ask for an acknowledgement after a measurement has been made. However, it may happen sometimes that problems occur during a determination run and that the expected signal from the connected instrument is not received. This is known as a timeout. This means that it is possible to define a maximum waiting time that will always be allowed to elapse if a timeout occurs. In addition, the behavior of the Sample Processor when this waiting period has elapsed (**SCAN** timeout) can also be defined.

<pre>parameters &gt;timeout settings</pre>	Submenu for behavior if a <b>SCAN</b> timeout occurs Open the submenu with <ENTER>
--	---

<pre>&gt;timeout settings SCAN timeout:      OFF min                     off,0...999 min</pre>	Waiting time if a SCAN timeout occurs
--	---------------------------------------

When the waiting time has elapsed the function defined below will be carried out.

<pre>&gt;timeout settings on SCAN timeout:  error                     error,continue</pre>	Behavior if <b>SCAN</b> timeout occurs
--	--

If **error** is set then the run sequence will be interrupted and an error message will be shown. Otherwise the sequence will be continued.

### 3.3.6 Manual stop settings

The following settings can make a great contribution to the comfortable operation of the Sample Processor. A definition is made of what is to happen when the **<STOP>** key is pressed. This could be a normal manual cancellation of a running method, or an emergency stop of the complete system. Depending on the arrangement of the automated system and the type of application, it can be laid down in detail how the individual instrument components and connected instruments are to react when the **<STOP>** key is pressed.

<b>parameters</b> <b>&gt;manual stop options</b>	Submenu for reaction after manual stop Open the submenu with <b>&lt;ENTER&gt;</b>
---	--

<b>&gt;manual stop</b> <b>CTL Rmt:</b> *****	Set signal lines at the Remote interface
Stop device1, Stop device2, Stop device* <b>14 bit (1,0 or *)</b>	

Connected peripheral devices (e.g. a Metrohm Titrino) can be stopped automatically. The 14 signal lines of the Remote interface can be set as required, see also p. 31ff and p. 109ff.

<b>&gt;manual stop</b> <b>CTL RS232:</b>	Command or character string, for output via the RS232 interface
14 ASCII characters	default value '&M;\$S' (= stops a Titrino)

Connected peripheral devices (e.g. a Metrohm Titrino) can be stopped automatically. Any character string can be transmitted. For details of the serial RS232 interface see "**Technical Reference**".

<b>&gt;manual stop</b> <b>PUMP:</b> <span style="float: right;"><b>OFF</b></span>	Switches the pumps, or pump connections
<b>OFF, ON, cont.</b>	(cont. = retain condition)

This setting applies to all pumps and/or pump connections.

<b>&gt;manual stop</b> <b>STIR T1:</b> <span style="float: right;"><b>cont.</b></span>	Switches the stirrer at Tower 1
<b>cont., ON, OFF</b>	(cont. = retain condition)

<b>&gt;manual stop</b> <b>STIR T2:</b> <span style="float: right;"><b>cont.</b></span>	Switches the stirrer at Tower 2 (only for 2-tower models)
<b>cont., ON, OFF</b>	(cont. = retain condition)

<b>&gt;manual stop</b> <b>STIR MSB1:</b> <span style="float: right;"><b>cont.</b></span>	Switches the stirrer at MSB connection 1
<b>cont., ON, OFF</b>	(cont. = retain condition)

```
>manual stop
STIR MSB2:                cont.
                             cont., ON, OFF
```

Switches the stirrer at MSB connection 2  
(cont. = retain condition)

```
>manual stop
STIR MSB3:                cont.
                             cont., ON, OFF
```

Switches the stirrer at MSB connection 3  
(cont. = retain condition)

## 3.4 Commands

### 3.4.1 Sample Processor commands

The following commands can be programmed to form a sequence. Most of them are also available in manual operation.

The following list applies to the programming of run sequences.

#### SAMPLE



```
>start sequence
```

```
1 SAMPLE: = 1
```

```
=, +, -
```

```
1...999
```

#### Set current sample position

1<sup>st</sup> parameter: function

2<sup>nd</sup> parameter: value

With the **SAMPLE** command you can define which sample (beaker position on the rack) is to be regarded as being the current sample position ( $\text{SAMPLE} = X$ ). This is stored in a run variable and can, for example, be altered in a sample sequence ( $\text{SAMPLE} + X$  or  $\text{SAMPLE} - X$ ) in order to selectively control a sample series.

#### Examples:

```
SAMPLE: = 5
SAMPLE: + 2
SAMPLE: - 1
```

Sets **SAMPLE** variable to 5, or first sample in the series to rack position 5

Increases **SAMPLE** variable by 2

Decreases **SAMPLE** variable by 1

If in a sample sequence the **SAMPLE** command is not programmed, then each run will automatically increase the **SAMPLE** variable by 1.

The **SAMPLE** command does not have to be used for simple applications. Unless anything to the contrary is required, the first sample in a series will automatically be assumed to be in rack position 1. This is why we recommend that special beakers are not placed on the first rack positions, but are set on the highest rack positions.


Under manual operation the **<SAMPLE>** key can be used before the start of a sample series to define the position of the first sample, provided that this has not already been defined in the method itself.

If a particular arrangement of sample beakers is always to be used for a certain application then the position of the first beaker can be defined in the start sequence with **SAMPLE = X** and this setting can be stored in the particular method.



The value of the **SAMPLE** variable remains at the end of a sample series. It is only reset to 1 when the instrument is switched on, when a **RESET** command (**<RESET>** key) is executed or after a **RACK** command (**<RACK>** key or in in a method run).

**MOVE**



```
>sample sequence
2 MOVE 1 : sample
      1,2
      sample,
      ext.1...ext.4
      spec.1...16
      prev., next
      +swing, -swing
      +rotate, -rotate
      1...999
      +/-1...999
```

**Position beaker / Swing robotic arm**

1<sup>st</sup> parameter: tower  
2<sup>nd</sup> parameter: position

With the **MOVE** command rack positions are moved to, i.e. the turntable rotates until the sample rack is positioned so that the selected rack position is placed in front of Tower 1 or 2 (if present). If no 786 Swing Head is mounted then only the rack positions of a 1-row rack can be accessed.

**With 786 Swing Head**

If a 786 Swing Head is mounted then any rack positions can be moved to. The angle of rotation of the rack is then compared with the swing angle of the robotic arm and corrected accordingly.

In addition to a rack position, the robotic arm of a 786 Swing Head can also be swung to any external position. The selectable external positions **ext.1** to **ext.4** can be defined in the whole swing range from 0° up to the maximum swing angle of the robotic arm. For example, this means that it is possible to move to a titration cell located beside the sample rack.

**Examples:**

<b>MOVE 1</b>	<b>sample</b>	Moves sample beaker (defined by SAMPLE variable) in front of Tower 1
<b>MOVE 1</b>	<b>ext.1</b>	Moves robotic arm at Tower 1 to external position 1
<b>MOVE 2</b>	<b>spec.1</b>	Moves special beaker 1 in front of Tower 2
<b>MOVE 2</b>	<b>5</b>	Moves rack position 5 in front of Tower 2 (absolute positioning)
<b>MOVE 1</b>	<b>next</b>	Moves next highest rack position in front of Tower 1
<b>MOVE 1</b>	<b>+2</b>	Depending on the current sample position (SAMPLE variable), moves the next position but one in front of Tower 1 (relative positioning)

**Parameters**

**sample** – rack position corresponding to the current value of the **SAMPLE** variable, see also p. 74.

**ext.1** to **4** – predefined angle positions of the robotic arm. These are defined in the configuration of a tower, see also p. 58.

**spec.1** bis **16** – reserved special beaker positions on the sample rack. These are defined in the rack configuration, see also p. 61.

**+swing, -swing** – moves the robotic arm by a certain increment angle. The sign shows the swing direction. The value of the angle is predefined under **Parameter >changer settings**.

**+rotate, -rotate** – moves the rack by a certain increment angle. The sign shows the rotation direction. The value of the angle is predefined under **Parameter >changer settings**.

**next, prev.** – from the current rack position the next highest (next) or the next lowest (prev.) will be moved to. Special beaker positions will be skipped. If **MOVE next** is used at the highest rack position then the move will be made to position 1. If **MOVE prev.** is used at rack position 1 then the highest possible position will be moved to.

**absolute positioning** – the numbered rack position will always be moved to, even when this is a reserved special beaker position.

**relative positioning** – If a numerical rack position is given with a positive or negative sign then the selected rack position will always be relative to the value of the **SAMPLE** variable, i.e. to the current sample position.

### Remarks

In a method sequence a **MOVE** command will move the lift (or both lifts) automatically to the shifting position.



*After execution of the functions **+/-swing** or **+/-rotate** the lift can be moved without restrictions. No defined rack position is required. Thus it is possible to cause damages if the lift is lowered without precautions. These functions must be used with great care.*

The direction of rotation is normally selected automatically by the Sample Processor. In the parameter menu under **>changer settings** the direction and speed of rotation can be defined specifically for particular methods. These can also be altered in a sequence with the corresponding **DEF** command.

If no beaker is present at the selected rack position then this will be recognized by the beaker sensor of the particular tower and the appropriate reaction will take place.

The reaction of the Sample Processors to a missing beaker can be defined in the parameter menu under **>changer settings**. You can choose between interrupting the run with an error message or selecting the next rack position (see p. 68). If a special beaker is missing the run will always be interrupted.

If a robotic arm with a tactile piezo sensor is installed and in use, a 'MOVE sample' command will automatically run the lift to work position. This is for the sensing of the sample beaker.

## LIFT



```
>sample sequence
3 LIFT: 1:      rest mm
      1,2,*
           work,
           rinse, shift,
           special, rest,
           0...235 mm
```

### Lift positioning

1<sup>st</sup> parameter: tower  
2<sup>nd</sup> parameter: position

Raises or lowers one or both (\*) lifts to a defined position.

### Examples:

```
LIFT: 1:  work  Moves lift at Tower 1 to working position
LIFT: 1:  rinse Moves lift at Tower 1 to rinsing position
LIFT: 2:  shift Moves lift at Tower 2 to shifting position
LIFT: 2:  special Moves lift at Tower 2 to special position
LIFT: *:  rest  Moves both lifts to uppermost position (0 mm)
LIFT: 1:  100 mm Moves lift at Tower 1 to position 100 mm
```

Working, rinsing, shifting and special positions are defined specifically for each rack in the configuration menu under **>rack definitions** (see p. 61). The rest position is the zero position (0 mm) of the particular lift, i.e. its upper stop.

The lift speed can be set in the parameter menu under **>sample changer settings** or altered in a sequence with the appropriate **DEF** command.

Each lift can be positioned with an accuracy of one millimeter. The **LEARN** function is available for this (see p. 112).

## 3.4.2 Switching components

Pumps and stirrers can be switched as required, either together or separately. They work independently and can be operated at the same time as other functions are being carried out.

## PUMP



```
>sample sequence
4 PUMP 1.1 : 1 s
      1.1...2.2 1...999 s,
      1.*,2.*   ON,OFF
```

### Pump control

1<sup>st</sup> parameter: pump selection  
2<sup>nd</sup> parameter: status or duration

With the **PUMP** command up to 4 pumps (2 pumps per tower) can be controlled separately. The first parameter selects the pump.

Syntax: T.P (T=Tower, P=Pump)

With Sample Processor models that have external pump connections instead of built-in pumps the **PUMP** command switches the appropriate pump connections (16 V).

The pump built into a tower as Pump 1 is always used for rinsing the titration head. Pump 2 (if present) can be used for aspirating off the

sample solution. Combined rinsing and aspiration is possible by **PUMP 1.\*** or **PUMP 2.\***.

The pumps can be switched on or off selectively, or operated for a particular time. The **LEARN** mode is very useful for determining the optimal rinsing or aspiration times (see p. 113).

## STIR



```
>sample sequence
5 STIR T1 : 1 s
T1, T2, T*,
MSB1...3, MSB*, * ON,OFF
1...999 s
```

### Stirrer control

1<sup>st</sup> parameter: stirrer selection  
2<sup>nd</sup> parameter: status, duration

With the **STIR** command up to 5 stirrers can be controlled separately. The first parameter selects the stirrer. With **STIR \*** all the stirrers can be switched at the same time .

#### Selecting a stirrer

**T1, T2, T\*** Stirrer connections as Tower 1 and 2 respectively (stirrer models 802 / 722 / 741), **\*** = both connections on the towers.

**MSB1...3** Stirrer/Dosing device connectors in the chassis of the Sample Processor (stirrer models 801 / 803 / 804), **MSB\***= all stirrers connected to an MSB socket.

**\*** all stirrers

The stirrers can be switched on or off selectively, or operated for a particular time.

In the parameter menu under **>stirring rates** the speed of each individual stirrer can be defined for a particular method. This can also be altered in a sequence with the appropriate **DEF** command.

### 3.4.3 Dosing device control

Connected dosing devices can be controlled separately or all at the same time . Simultaneous addition of the same volumes or simultaneous filling of the connected dosing devices is also possible (Example: **DOS \*.2 fill** = fill all dosing devices via port 2).

## DOS



```
>sample sequence
6 DOS 1.* : 1 ml
*, *, 1.*...3.4
fill, release, prepar.
empty, eject, endVol,
compen., port,
±0.001...1...±999.999 ml
```

### Dosing device control

1<sup>st</sup> parameter: dosing device and port selection  
2<sup>nd</sup> parameter: selects function/enters volume

The **DOS** command is used for controlling Dosimats and Dosinos. Up to 3 Dosinos or Dosimats can be addressed individually or simultaneously via the **MSB** bus control.

The 1st parameter selects the dosing device and the corresponding dosing port, at which the required function is to be carried out. If '\*' is entered then the default port for the corresponding function will be used (e.g. dosing port = 1, filling port = 2, etc.)


**Please note**

After dosing the Sample Processor **will not automatically fill** the dosing cylinder. If required, this can be programmed with the command **DOS: 1.\* : fill**, see below.

It is possible to enter the volume to be added directly as the 2<sup>nd</sup> parameter or to carry out specific functions of a Dosino. Negative volumes can also be added, i.e. a certain volume can be aspirated and then ejected again (pipetting). The minus sign is entered with the **<\*>** key.

The functions listed below are used for complex liquid handling tasks such as pipetting.

**Dosing functions:**

<b>fill</b>	Fills the Dosimat and Dosino cylinder.
<b>release</b>	Prepares the Dosimat or Dosino for changing the Exchange/Dosing unit. The dosing cylinder is filled via the rinse port. The stopcock is turned to the exchange position (Port 2).
<b>prepar.</b>	Preparation cycle ( <b>PREP</b> ) for Dosinos. All tubing is rinsed and filled completely.
<b>empty</b>	Tubing system and Dosino cylinder are emptied.
<b>eject</b>	Complete Dosino cylinder contents are ejected.
<b>endVol</b>	Ejects the cylinder content to the max. volume mark.
<b>compen.</b>	Cancels the mechanical play between the dosing piston and spindle.

In the parameter menu under **>Dosing unit def.** the port assignments of the Dosinos as well as the dosing and filling rates can be defined for specific methods. This can also be carried out in a sequence with the appropriate **DEF** command.

More details about Dosino commands are given on p. 101ff.

The Sample Processor automatically recognizes whether a Dosimat or a Dosino is connected.

### 3.4.4 Communication commands

The coordination of external instruments or the explicit triggering of functions is possible via both the Remote interface and the serial RS232 interface. The interfaces can be scanned for particular signal patterns or incoming data strings. In return, individual signal lines can be set or any character strings (as control commands) can be transmitted to connected instruments.

#### SCAN



>sample sequence		<b>Scans the Remote interface</b>
7	SCN:Rm : Ready1	
Rm, RS	Ready1	= Device 1 ready
	Ready2	= Device 2 ready
	Ready*	= Devices 1+2 ready
	Cond ok	= Conditioning finished
	End1	= EOD impulse device 1
	End2	= EOD impulse device 2
	EndMeter	= End impulse 692 Ion Meter/713 pH Meter
	Continue	= Continue impulse from instrument
8 Bit (1,0 or *)		any 8-bit binary pattern

In a sequence the **SCN:Rm** command will halt the run sequence until the predefined signal pattern has been received.

Predefined binary patterns are available; these can be selected via simple names (e.g. **Ready1** or **End2**).

**Ready** describes a statically set **Ready** line (Output 0) of a connected Metrohm instrument. **End** stands for a pulsed signal, e.g. **EOD** (=End of Determination, Output line 3). When pulsed signals are being scanned it is not possible to scan several lines in parallel.

Setting special binary patterns allows the flexible monitoring of connected instruments. The following applies:

- 0** = line inactive
- 1** = line active
- \*** = any line condition

Example: **0000001** = input line 0 is active = **Device 1: Ready1**

With the **LEARN** function the binary patterns (=line conditions) can be adopted interactively (see p. 113).

Details about the Remote interface can be found in the **Technical Reference**.



>Sample sequence  
8 SCN:RS

### Scanning the RS232 interface

Rm,RS  
 default value: \*R" = scan "Ready" status message  
 14 ASCII characters any data string with 14 characters

In a sequence the **SCN:RS** command will halt the method sequence until the predefined data string (up to 14 characters) has been received via the serial RS232 interface. Incoming data is checked character by character.

Make sure that the transmission parameters of the RS232 interface are identical with those of the connected device (see configuration menu **>RS232 settings**, p. 65).

Any letters, numbers and special characters can be selected from the character set of the Sample Processor. '\*' can be set as wild card for any character string. (If '\*' is to be interpreted as an ASCII character then '\*\*' must be used.) A wild card can be set within a character string. When the first part of the character string has been recognized correctly then a search will be made for the first appearance of the character standing after '\*'. In this case the comparison will be made with the second part of the character string.

This function is primarily suitable for devices that understand the Metrohm Remote language, such as Titrinos. In this case the **AutoInfo** status messages can be scanned. The most useful of these are:

- \*.T.R" Ready, 'Ready' state has been achieved, e.g. after titration
- \*.T.F" Final, end of the determination has been achieved
- \*.T.S" Stop, instrument has been stopped manually
- \*.T.G" Go, instrument has been started
- \*.E;\* Error, error message

However, these status messages will only be transmitted when the corresponding status message has first been switched on, e.g. in the start sequence, e.g. for a Titrino with the command: CTL:RS &Se.A.T.R"ON".

More detailed information about the syntax can be found in the Instructions for Use of the instrument whose status messages are to be transmitted.

Details of the **CTL** command are given below.

With the **LEARN** function transmitted data (=character strings) can be adopted interactively (see p. 113).

## CTL



```
>sample sequence
9 CTL:Rm START Device1
```

### Setting the Remote lines

Rm,RS	START Device1	= Starts device 1
	START Device2	= Starts device 2
	START Device*	= Starts devices 1+2
	SAMPLE ready	= Continue impulse to instrument
	START Dos1	= Starts Dosimat at device 1
	START Dos2	= Starts Dosimat at device 2
	START Dos*	= Starts Dosimat at devices 1+2
	METER Mode pH	= Switches pH meter to pH measurement
	METER Mode T	= Switches pH meter to temp measurement
	METER Mode U	= Switches pH meter to mV measurement
	METER Mode I	= Switches pH meter to Ipol mode
	METER Mode C	= Switches ion meter to Conc measurement
	METER Cal pH	= Switches pH meter to pH calibration
	METER Cal C	Switches ion meter to Conc calibration
	METER enter	= <ENTER> key for pH meter simulation
	INIT	= Initializes Remote interface
14 bit	(1,0 or *)	any binary pattern with 14 bit

The **CTL:Rm** command is used for controlling external devices via the Remote interface. It sets defined line conditions or transmits impulses via the 14 Remote output lines.

Predefined binary patterns are available; these can be selected via simple names (e.g. **START device 1** or **METER Mode pH**).

**START device x** starts the set mode of a connected Metrohm device. **START Dos x** starts a Dosimat which is connected with a Metrohm titrator via the "activate" line (a special cable is required). **METER xxx** switches to a particular measuring mode with 691, 713, 780 pH Meters and 692, 781 Ion Meters.

**Sample ready** can be used as a continue impulse, e.g. for a connected Titrande.

Setting special binary patterns allows the flexible monitoring of connected instruments. The following applies

0 = line inactive  
 1 = line active  
 \* = retain line condition

Example: \*\*\*\*\*1 = Output line 0 active = start device 1

Details about the Remote interface can be found in the **Technical Reference**.



```
>sample sequence
10 CTL:RS
```

### Data transmission via serial interface

Rm,RS

default value: M;\$G = starts Metrohm instrument in current mode  
14 ASCII characters any character string with 14 characters

Data (=character strings) can be transmitted to connected devices via the serial RS232 interface.

Make sure that the transmission parameters of the RS232 interface are identical with those of the connected device (see configuration menu **>RS232 settings**, p. 65).

Any letters, numbers and special characters can be selected from the character set of the Sample Processor.

This function is primarily suitable for instruments that understand the Metrohm Remote language. These can be controlled with so-called triggers.

The most important of these are:

<b>&amp;M;\$G</b>	Go, starts instrument in current mode
<b>&amp;M;\$S</b>	Stop, stops device
<b>&amp;M;\$H</b>	Hold, interrupts determination
<b>&amp;M;\$C</b>	Continue, restarts determination

Switching on the **AutoInfo** status messages (e.g. in a start sequence) can be carried out with the following Remote commands:

<b>&amp;Se.A.T.R"ON"</b>	Status message at "Ready" condition
<b>&amp;Se.A.T.F"ON"</b>	Status message at end of a determination
<b>&amp;Se.A.T.S"ON"</b>	Status message at manual Stop
<b>&amp;Se.A.T.G"ON"</b>	Status message at start of a method
<b>&amp;Se.A.T.E"ON"</b>	Status message at error condition

Logically the corresponding **AutoInfo** messages should be switched off again in a final sequence (...**"OFF"**).

Detailed information about the syntax of the Remote language can be found in the **Technical Reference** or in the Instructions for Use of your titrator.

For communication with instruments from other manufacturers or a computer please conform with their syntax and conventions.

### 3.4.5 Auxiliary commands

#### WAIT



```
>sample sequence
11 WAIT: pause 1 s
```

**Waiting time**


pause, runtime  
0...1...9999 s

The **WAIT** command is used for setting a particular waiting time or for waiting until a particular time (run time) in a run sequence.

If **pause** is selected then the method run will be interrupted for the entered duration (in seconds).

If **runtime** is selected then the method run will be interrupted until the selected running time (seconds counter) has been reached. The running time (in seconds) always starts at the beginning of an individual sequence, i.e. start sequence, sample sequence or final sequence. If the running time of the sequence has already been reached or exceeded when the **WAIT** command is executed then the method will be continued immediately.

#### RACK



```
>Sample sequence
12 RACK
```

**Initialize rack**

The rack turns to the starting position, i.e. the lift (or both lifts) move upward to the rest position and the rack is rotated to the initial position. The rack code is read off and the **SAMPLE** variable (position of the current sample) reset to 1. The **RACK** command should only be used in a final sequence.

#### DEF



#### Redefinition of specific device settings

With the following **DEF** commands a wide range of settings can be modified during a method run. The individual entries are selected by repeatedly pressing the **<DEF>** key (selection drum).

```
>sample sequence
13 STIRRATE: T1 3
```

**Stirring rate**

T1, T2, T\*,  
MSB1, MSB2,  
MSB3, MSB\*,  
\*  
1...3...15

The stirring speed can be set individually for all stirrers. The first parameter stands for the connector of the stirrer, the second parameter allows the stirring speed to be set in 15 steps.

```
>sample sequence
13 DOSRATE 1 160
```

### Dosing rate

```
1...3
0.01...160 mL/min
```

The dosing rate can be set individually for all three dosing devices. The first parameter stands for the address (MSB1...3) of the dosing device, the second parameter allows the setting of the dosing rate in mL/min.

```
>sample sequence
13 FILLRATE 1 160
```

### Filling rate

```
1...3
0.01...160 mL/min
```

The filling rate can be set individually for all three dosing devices. The first parameter stands for the address (MSB1...3) of the dosing device, the second parameter allows the setting of the filling rate in mL/min.

```
>sample sequence
13 COCKMOVE 1 : auto
```

### Direction of rotation of the stopcock of a Dosing unit

```
1...3
auto,
not over, desc.,
asc.
```

The direction of rotation of the stopcock can be defined for each Dosing unit. Apart from the automatic selection (shortest path) a descending or ascending direction can also be selected. **Not over** means that rotation will never be carried out over the protected port (see Section 3.3.4).

```
>sample sequence
13 LIFTRATE 1 25 mm/s
```

### Lift rate

```
1,2 5...25 mm/s
```

The lift speed can be set individually for both towers. The first parameter stands for the selection of the tower, the second parameter allows the setting of the lift speed in mm/s.

```
>sample sequence
13 SHIFTRATE: auto. 20
```

### Direction and speed of rotation of the rack

```
auto., +, - 3...20 w/s
```

Direction of rotation and speed of rotation of the sample rack can be altered as required. The first parameter determines the direction of rotation.

- auto.** : The Sample Processor automatically selects the shortest path.
- +** : The rack rotates counterclockwise (to a higher rack position)
- : The rack rotates clockwise (to a lower rack position)

The second parameter determines the speed of rotation in degrees/s.

>sample sequence
13 SWINGRATE 1                      55

1,2 10...55 °/s

### Swing rate of a robotic arm

The swing rate of a robotic arm (786 Swing Head) can be set individually for both towers. The first parameter stands for the selection of the tower, the second parameter allows the swing rate to be set in degrees/s.

### 3.5 Managing methods



#### 3.5.1 User-defined methods

Main menu:

```

methods
>recall method

methods
>store method

methods
>delete method
    
```

Open submenu with **<ENTER>**

Use **<↑>** or **<↓>** up or down by one menu item

Use **<HOME>** or **<END>** to reach the first and last menu items

```

methods
>recall method
    
```

Dialog for loading methods  
Open the dialog with **<ENTER>**

Use **<QUIT>** to return to normal operating condition.

```

>recall method
name:          *****
    
```

Select method

8 ASCII characters

All the stored methods can be selected with **<SELECT>**. If an "empty" method is to be loaded then the method **\*\*\*\*\*** can be selected with **<CLEAR>**. This deletes the current working memory for methods.

```

methods
>store method
    
```

Dialog for saving methods  
Open the dialog with **<ENTER>**

```

>store method
name:          *****
    
```

Define method name

8 ASCII characters

The text input mode is activated with **'<'** or **'>'** in order to be able to enter your chosen method name (see p. 55).

```

methods
>delete method
    
```

Dialog for deleting methods  
Open the dialog with **<ENTER>**

```

>delete method
name:          *****
    
```

Select method

8 ASCII characters

```

>delete method
delete ***** ?
    
```

Confirm with **<ENTER>**  
Cancel with **<QUIT>**

### 3.5.2 **POWERUP** method

When the Sample Processor is switched on the sample rack and the lifts move to their initial positions. This means that electrodes may also be removed from the conditioning beaker. In order to re-immerses them in the conditioning beaker you can use the **POWERUP** method. This method is automatically started when the Sample Processor is switched on and the initialization of the instrument is finished.

Create a method which contains a command sequence that is to be processed when the Sample Processor is switched on. Save this method under the name **POWERUP** (see p. 87).

## 3.6 Run control



With **<START>** you can start a method from the normal operating condition. If no manual intervention is made, or if no unexpected error occurs, the sample series will be processed correctly and terminated with the final sequence. The sample sequence will be repeated several times in accordance with the entry made under **<PARAM>**, **number of samples**, starting with the sample beaker defined as **SAMPLE**.

If **External START** (see configuration 3.2.1) is switched on then the activation of the *Remote line Input 7* will also start the method.



If the sample series is halted with **<STOP>** then the Sample Processor will return directly to the normal operating condition. Unprocessed samples will not be taken into account and the final sequence will not be carried out. If under **Manual stop** settings have been defined for this case, then the corresponding actions or commands will be carried out, e.g. connected devices will be stopped.

If **external START** (see configuration 3.2.1) is switched on then the activation of *Remote line Input 6* will also stop the method.




With **<HOLD>** the method run can be interrupted. The currently active command will be terminated immediately. With **<START>** the method can be continued with the following command in the active sequence. The connected peripheral devices will **not** be interrupted with the **<HOLD>** key.



**<QUIT>** breaks off the command which is currently being executed and starts the next command line in the sequence.

If a fault occurs during a sample series an appropriate **error message** will be shown; this must be confirmed with **<QUIT>**. The Sample Processor then returns to the **HOLD** status (see above). After the fault has been remedied the run can be continued with **<START>** or, if necessary, cancelled with **<STOP>**.



**<CLEAR>** interrupts a sample series after the currently active sequence **has been carried out** (*gentle stop*). The current sample will be processed to the end.

## 3.7 Sample racks

A sample rack is a turntable that accommodates sample vessels and is placed on the Sample Processor. In order to be able to process different shapes and sizes of sample vessels various types of sample rack can be used and easily exchanged. The number of samples that a rack can accommodate depends on the diameter of the sample beakers.

- Sample racks **42 cm diameter** can be used with all Sample Processor models. With the 789 Robotic Sample Processor XL models a 786 Swing Head with robotic arm is required.
- Sample racks **48 cm diameter** are only suitable for the 789 Robotic Sample Processor XL models with 786 Swing Head and robotic arm.

### 3.7.1 Metrohm standard sample racks

#### Racks 42 cm diameter

Article no. / Rack type	No. of samples	Type of sample vessel	Vessel diameter	Magnet code (predefined)
6.2041.310	12	250 mL Metrohm titration beaker	65 mm	000001
6.2041.320	16	150 mL beaker	55 mm	000010
6.2041.340	24	75 mL Metrohm titration beaker	35 mm	001000
6.2041.350	48 *)	75 mL Metrohm titration beaker	35 mm	010000
6.2041.360	12	150 mL beaker or 200 mL disposable beaker	55 mm	100000
6.2041.370	14	200 mL disposable beaker	55 mm	000011
6.2041.380	14	8 oz disposable beaker	59 mm	000101
6.2041.400	126+2 *)	11 mL test tube 250 mL Metrohm titration beaker	16 mm 65 mm	001010
6.2041.410	141+1 *)	11 mL test tube 500 mL beaker	16 mm 71 mm	001010
6.2041.430	127+2 *)	11 mL test tube 2x 300 mL beaker	16 mm 68 mm	010001
6.2041.440	148+3 *)	11 mL test tube 3x 300 mL beaker	16 mm 68 mm	010100
6.2041.460	21	100 mL CSB tubes	43 mm	101000
6.2041.710	160	6 mL vials	22 mm	000111

\*) can only be used with 786 Swing Head

**Racks 48 cm diameter** (for 789 XL models only)

Article no. / Rack type	No. of samples	Type of sample vessel	Vessel diameter	Magnet code (predefined)
6.2041.800	100 *)	75 mL Metrohm titration beaker	35 mm	000100
6.2041.810	34 *)	150 mL beaker or 200 mL disposable beaker (Euro)	55 mm	001001
6.2041.820	28 *)	250 mL Metrohm titration beaker	65 mm	010010
6.2041.830	28 *)	200 mL PP sample beaker	55 mm	100010
6.2041.840	59 *)	120 mL PP sample beaker	42 mm	001100

\*) can only be used with 786 Swing Head

On request we can supply further user-defined racks that can be defined in the instrument via PC software. The beaker positions can be arranged as required.

**3.7.2 Magnet codes**

Each individual sample rack can be identified unambiguously by a magnet code. Magnet pins attached to the base of the rack (see Fig. 5) can be arranged to form a 6-place binary code. The Sample Processor can then automatically recognize which rack has been placed on it if the rack is located at the starting position after initialization.

When a rack is changed it should first be brought to the starting position by pressing the **<RACK>** key. In this way the rack can be positively recognized, which makes correct beaker positioning possible. Each type of rack is assigned to an internal position table, in which the angle of rotation and distance to the center of the rack is defined for each rack position.

When a sample series is started the Sample Processor automatically moves the rack to the starting position so that its magnet code can be read off. This ensures that the beaker positions coincide with the internal position table of the current rack.

The standard racks supplied by Metrohm are already provided with a predefined magnet code for each type of rack. If several racks of the same type are in use then the magnet pins can be arranged differently to allow the unambiguous identification of each individual rack.

**Magnet code format** (examples):

000001 i.e. only one magnet inserted, Bit 0

000101 i.e. two magnets inserted, Bit 0 and 2

63 different combinations are possible. Code 000000 stands for "no defined code".

### 3.7.3 Rack data

In many applications the method and beaker sizes must be observed exactly. As sample racks are intended for particular beaker sizes, the rack definitions also contain, in addition to the rack positions themselves, data about the lift positions that are directly associated with the particular beaker size.

The following data can be defined for each rack:

<b>Rack name</b>	<i>Unambiguous identification, ordering number is standard</i>
<b>Code</b>	<i>Magnet code for automatic rack recognition</i>
<b>Work position</b>	<i>Working height for Lift 1 and 2*</i>
<b>Rinse position</b>	<i>Rinsing height for Lift 1 and 2*</i>
<b>Shift position</b>	<i>Shifting height for Lift 1 and 2*</i>
<b>Special position</b>	<i>Additional height for Lift 1 and 2*</i>
<b>Beaker radius</b>	<i>Radius of sample beaker</i>
<b>Beaker sensor</b>	<i>Beaker checking mode</i>
<b>Special beaker positions</b>	<i>16 reserved rack positions with individual working heights, beaker radius and beaker sensor test</i>

*\* can be set for both towers separately*

The **rack name** is used for the unambiguous identification of a rack. The standard name stored for the rack is the ordering number. In a method, this can be assigned to a particular rack name (see Section 3.3.1). The automatic rack recognition ensures that if an incorrect rack is used this will be recognized and the user will be informed of this.

The **code** is used for automatic rack recognition. You must make sure that this 6-place binary code coincides with the actual magnet code set on the rack. Rack codes can be altered as required. However, they must always be assigned to a single rack only. Avoid the use of the predefined codes for the standard racks supplied by Metrohm, see also p. 90.

The **working position** is used for defining the lift position at which a sample is to be processed. This means that the ideal setting for a particular sample rack can be defined depending on the height of the sample beakers. In manual operation this working position can be moved to directly with the **<END>** key. In a run sequence this can be programmed with **LIFT :1 : work mm.**

The **rinsing position** is used for defining the lift position at which, for example, the electrode is to be rinsed. This means that the ideal setting for a particular sample rack can be defined depending on the height of the sample beakers. In a run sequence this can be programmed with **LIFT: 1: rinse mm.**

The **shifting position** is used for defining the lift position at which the rack can be rotated. If the lift is located below the shift position then it will be moved automatically to the shifting height before the rack is rotated. This safety feature prevents the electrode from being damaged by rotational movements of the rack to a large extent. However, a pre-requirement is that this shifting height is set correctly. In a run sequence the movement of the lift to the shift position can be programmed with **LIFT: 1 : shift mm.**

The **special position** is a further user-defined lift position. For example, when pipetting with a robotic arm it can be selected so that the pipet tip is located directly above the sample solution so that a separating bubble (air gap) is formed. In a run sequence this can be programmed with **LIFT: 1 : special mm.**

The **beaker radius** can be used to prevent a titration head from trying to enter a beaker which is too narrow for it, which could damage the electrode or sample beaker. By entering the beaker radius a Sample Processor can decide whether a titration or transfer head on the lift will “fit” into the particular sample beaker, see also Section 3.2.

The **beaker sensor** recognizes whether a beaker is present. At the moment Metrohm Sample Processors support optical infrared sensors built into each tower of the Sample Processor. In the rack definition the beaker sensor can be switched on (**Tower**) or off (**OFF**).

The setting for the beaker sensor can be made individually to suit the particular samples and the individual special beaker positions, see below.

## Special beaker

Special beakers are reserved positions on a sample rack. Up to 16 special beaker positions can be defined per rack. They can be selectively moved to during a run sequence without interrupting the sample series run or interfering with it. Special beakers can be used for rinsing the electrode during a sample sequence or for calibrating an electrode in a start sequence (buffer solutions), etc.

Reserved special beaker positions are automatically recognized in a run sequence and are skipped over when the individual sample beakers themselves are being processed.

Special beakers are moved to with **MOVE 1 : spec.1.**

The following settings can be made separately for each special beaker position on a rack:

- **Rack position**
- **Working height at Tower 1**
- **Working height at Tower 2**
- **Beaker radius**
- **Beaker sensor**

If a special beaker is required in a run sequence, but the Sample Processor cannot find a beaker at the reserved position, then an error message will always appear.

## 3.8 Dosing and liquid handling

### 3.8.1 Dosimats and Dosinos

Three 685/805 Dosimats or 700/800 Dosinos can be used as dosing instruments connected directly to an **MSB** socket. They can be controlled with the **DOS** command.

Each Dosimat or Dosino can be equipped with various exchange units or Dosing units. Before these units are exchanged the buret stopcock must always be moved to the exchange position, as otherwise when the buret is removed there is the risk of serious damage to the buret itself or to the Dosimat or Dosino drive unit.



**Before removing the Dosing unit or exchange unit always trigger the 'DOS: X.Y : release command!'**

With Dosimats and Dosinos any volumes of auxiliary solutions up to 999 mL can be added. With both types of instrument the filling function can be selectively triggered (**DOS: X.Y : fill**). When switched on the dosing or exchange unit will always be filled via Port 2 (filling port).

The Sample Processor automatically recognizes the type of dosing instrument connected to it.

Further commands are available for the 700/800 Dosinos, so that the wide range of possibilities which characterize these Dosing units can be used to the full extent.





Fig. 34 800 Dosino with Dosing units

The Dosing unit has four ports (inlets/outlets) to which various functions can be assigned. There is an additional outlet (**VENT** or **0**) that is used for venting the bottle on which the Dosing unit is mounted.



*Fig. 35 Dosing unit from below*

**VENT/0**

- is used for venting the storage bottle and is normally fitted with an absorber tube (filled with a drying agent).

**Port 1**

- is side-mounted and is defined as dosing outlet 1 as standard.

**Port 2**

- is located on the base and is defined as the filling inlet as standard; it is normally fitted with a standpipe.

**Port 3**

- is side-mounted and is defined as dosing outlet 2 as standard.

**Port 4**

- is located on the base and is defined as the air inlet opening for emptying the tubing system as standard.

The maximum dosing and filling rates, which can be entered for each port of the Dosing unit in the configuration menu under **>Dosing units**, depend on the cylinder size:

Cylinder volume	Max. dosing rate	Resolution
2 mL	6.6 mL/min	0.2 $\mu$ L
5 mL	16.6 mL/min	0.5 $\mu$ L
10 mL	33.3 mL/min	1.0 $\mu$ L
20 mL	66.6 mL/min	2.0 $\mu$ L
50 mL	160 mL/min	5.0 $\mu$ L

The following commands can be executed with Dosinos. For each command the dosing drive and the Dosing unit port at which the command is to be executed can be entered. The addition of a particular

(positive) volume and filling the cylinder are also possible with 685 and 805 Dosimats.

If \* is entered as **wild card for the drive** then the selected function will be carried out by all the connected dosing devices.

If \* is entered as **wild card for the port** then the particular function will be carried out at the port defined as standard in the parameter menu under **>def. Dosing units**. This means that this setting is method-specific, but also applies to manual operation depending on the loaded method.

## Dosing

**DOS: X.Y : xxx.xx mL** Adds a particular volume

The given volume will be ejected through the selected port. If a value with a negative sign is entered then the volume will be aspirated.

The Dosing unit will **not be refilled** after each addition. The standard dosing port (Port \*) is that defined under

```
>dosing drive      X
dosing             port Y
```

(default: Port 1).

## Filling

**DOS: X.Y : fill mL** Fills the cylinder

The Dosing unit is filled completely. The liquid is aspirated via the given port. The standard filling port (Port \*) is that defined under

```
>dosing drive      X
filling            port Y
```

(default: Port 2).

## Prepare

**DOS: X.Y : prepar. mL** Prepares = rinses and fills the connected tubing and the dosing cylinder

The tubing system of the Dosino should be freed from air bubbles at least once per day by carrying out a preparation cycle. This is a process that could take some time.

It is recommended to use this command in a start sequence.

During the preparation process the dosing cylinder as well as the connected tubing are completely filled. Several filling and dosing processes are carried out. The volumes required for this are calculated from the configuration settings of the Dosing unit, i.e. from the tubing lengths and diameters (see Section 3.2.4).

The contents of the dosing cylinder is ejected via the selected port. The standard PREP port (Port \*) is that defined under

```
>dosing drive          X
preparation          port Y
```

(default: Port 1).

## Empty

**DOS: X.Y : empty mL** Completely empties the dosing cylinder and tubing

The tubing system and Dosing unit cylinder can be completely emptied. The liquid in the dosing cylinder is ejected via the dosing port. The air required to displace the liquid from the tubing is aspirated via the given port. The standard port (Port \*) for aspirating the air is that defined under

```
>dosing drive          X
drain                  port Y
```

(default: Port 4).

The standard dosing port can be altered under

```
>dosing drive          X
dosing                 port Y
```

(default: Port 1).

## Exchange Dosing unit

**DOS: X.Y : release mL** Prepares Dosino for exchanging the Dosing unit.

Before changing the Dosing unit the **release** command must be used to fill the dosing cylinder and move the stopcock to the exchange position. The cylinder is filled by aspirating the necessary volume via the given port. If \* is entered as a wild card then the port given under

```
>dosing drive          X
rinsing                port Y
```

will be used (default: Port 2).

## Eject

**DOS: X.Y : eject mL** Ejects the complete dosing cylinder contents

The cylinder are completely ejected via the given port. The piston moves down to the stop. If \* is entered as a wild card then the port given under

```
>dosing drive          X
dosing                port Y
```

will be used (default: Port 1).

## Endvolume

**DOS: X.Y : endVol mL** Ejects the cylinder volume to the max. volume mark

The content of the dosing cylinder is ejected via the given port. The piston is run the nominal volume. If \* is entered as a wild card then the port given under

```
>dosing drive          X
dosing                port Y
```

will be used (default: Port 1).

## Compensate

**DOS: X.Y : compen. mL** Eliminates the mechanical play

The mechanical play between the dosing piston and drive spindle is eliminated after the stopcock has been rotated to the given port. If \* is entered as a wild card then the port given under


```
>dosing drive          X
dosing                port Y
```

will be used (default: Port 1).

### 3.8.2 Liquid handling functions

The Sample Processors can use the numerous capabilities of a Metrohm Dosino (700 or 800) to the full extent. The four ports of the Metrohm Dosing units for the Dosino can be used as outlet or inlet ports as defined by the user. This means that not only simple addition and filling tasks can be carried out, but complex liquid handling tasks such as pipetting or sample transfer can be carried out without any problems.

The dosing functions of the Metrohm Sample Processors are to be used in such a way that, in addition to the function, the port at which the function is to be carried out is also given. This means the Dosing unit inlet or outlet to which the stopcock first moves in order that the required function is carried out.



*Fig. 36 Dosing unit - ports*

### 3.8.3 The DOS command

The liquid handling command DOS has two parameters:

<i>General:</i>	<b>DOS:</b>	Address	Function
<i>Example:</i>	<b>DOS:</b>	<b>1.1</b>	<b>5 mL</b>
<i>Parameter:</i>		drive.port	Volume or function (negative volume permitted)


for drive: \* = all dosing drives  
 for port: \* = default port for particular function

<i>Example:</i>	<b>DOS:</b>	<b>*.*</b>	<b>fill mL</b>
<i>means:</i>	On all the connected dosing devices the cylinder is filled via the default port.		
<i>Default ports:</i>	The default ports are the ports assigned in the parameter menu under <b>&gt;Dosing unit def.</b> , see also Section 3.3.4.		

### 3.8.4 Pictograms

For complex liquid handling tasks only Dosinos (700 or 800 models) can be used as dosing drive. In order to clarify the various functions and processes the following pictograms are used on the next few pages.

The following arrangement applies:



The pictograms above show the default settings..

### 3.8.5 Liquid handling functions in detail


#### Dosing

Example:

**DOS: 1.1 : 1.000 mL**

Default port = 1

'Normal' dosing can be carried out by entering a volume. Automatic cylinder filling is triggered neither before nor after the dosing process. If the dosing piston reaches 'max. Volume' mark (10'000 impulses) during the dosing process then the cylinder will be refilled.



#### Dosing a negative volume

Example:

**DOS: 1.2 : -1.000 mL**


Default port = 1

If a negative value is set for a volume then the dosing process will take part in the opposite direction, i.e. a liquid will be aspirated through the given port. Automatic cylinder filling is triggered neither before nor after the dosing process.

If the dosing piston reaches the zero mark during the dosing process then the cylinder will be refilled.

Do not select any volume that is larger than the nominal cylinder volume. Aspiration should take place with a single piston stroke.

This function can be used for pipetting.




#### Filling

Example:

**DOS: 1.3 : fill mL**

Default port = 2

Filling the cylinder can be carried out from a freely selectable port. After the filling process the stopcock remains at the selected port.



### Exchanging Dosing unit


Example:

**DOS: 1.2 : release mL**

Default port = 2

This command can be used before exchanging a Dosing unit to fill the cylinder via the given port. In this way, for example, air can be aspirated via Port 4.

After the cylinder has been filled the stopcock will move to Port 2. The dosing drive can then be removed from the Dosing unit.



### Preparing the Dosing unit


Example:

**DOS: 1.\* : prepar. mL**

Default port = 1

In order to prepare a Dosing unit for use a complex procedure is started: first the cylinder contents are ejected via the given port, then the volume of the filling tubing is drawn in and ejected again via the given port. All the connected pieces of tubing are then filled bubble-free (tubing length > 0 mm) and finally the cylinder is filled again.

The volumes of the connected pieces of tubing (calculated from length and diameter) are taken into account during the whole procedure.



### Automatic emptying


Example:

**DOS: 1.4 : empty mL**


Default port = 1

Automatically emptying the Dosing unit takes place according to a complex procedure: first the cylinder contents are ejected via the given port, the ejection tubing is then briefly rinsed with reagent from the filling tubing and then all the pieces of tubing are emptied in sequence. Air is always aspirated via the drain port (default: Port 4) for emptying.

The volumes of the connected pieces of tubing (calculated from length and diameter) are taken into account during the whole procedure.



### Ejecting




Example:

**DOS: 1.1 : eject mL**                      Default port = 1

The whole cylinder contents are ejected via the given port. The piston is lowered beyond the max. volume mark to the stop.  
This command should be used for removing air bubbles.

### Run piston to end volume




Example:

**DOS: 1.1 : endVol mL**                      Default port = 1

The cylinder content is ejected via the given port. The piston is run to the max. volume mark.  
This command should be used for the pipetting functions.

### Compensating the mechanical play




Example:

**DOS: 1.4 : compen. mL**                      Default port = 1

As the Dosing units are exchangeable, the coupling of the Dosino connecting rod (spindle) has a slight degree of mechanical tolerance which is noticeable when the direction of the piston movement changes. This tolerance can be compensated. This is done by first making a short piston movement in the same direction as the previous piston movement, followed by an equivalent piston movement in the opposite direction.

### Stopcock rotation



Example:

**DOS: 1.3 : port mL**                      Default port = 1

The stopcock rotates to the given port. No piston movement takes place. The direction of rotation is determined by the parameter **stopcock direction** under **>Dosing unit def.** in the parameter menu.

### 3.8.6 Pipetting equipment

Together with an 700/800 Dosino as dosing drive (or "pipetting pump"), the Sample Processor is exceedingly suitable for pipetting liquids in a volume range from 0.1 to 10 mL.

#### Equipment required

- Sample Processor with 786 Swing Head
- Robotic arm with transfer head (e.g. 6.1462.030 or 6.1462.040)
- 800 Dosino with 2 mL, 5 mL, 10 or 20 mL Dosing unit
- Pipetting tubing 6.1562.020 (3 mL) or 6.1562.100 (10 mL) with drawn out tip
- Any sample rack
- Possibly Stand Support 6.2001.070 with titration cell

### 3.8.7 Pipetting procedure

#### Basic principles

Measuring out the liquid to be pipetted takes place in a pipetting tubing filled with a carrier liquid (preferably water). In order to prevent any carryover or mixing an air bubble (air gap) must be inserted between the sample and the carrier liquid. The tubing tip should be drawn out by hand to form a fine point.


A Dosino drive is used to transfer the liquids into the pipetting tubing. The pipetting tubing is connected to Port 1 of the Dosing unit. The Dosing unit is mounted on a reagent bottle containing the carrier liquid. This can be refilled via filling port 2 of the Dosing unit. The carrier liquid can also be used for rinsing and diluting the sample.

Various pipetting procedures are possible. The following conditions must always be observed:

- The sample must always be measured out in the pipetting tubing; it must never enter the Dosing unit cylinder.
- The sample should always be aspirated after the dosing cylinder has been ejected.
- The difference in level between the surface of the sample liquid and that of the carrier liquid in the cylinder should be as small as possible.
- Aqueous sample solutions can be pipetted without having to rinse the tubing. Non-aqueous samples require the pipetting tubing to be rinsed with carrier liquid (same solvent as the sample), to avoid carryover.

### 3.8.8 Preparing the Dosing unit

The Dosing unit must first be prepared, i.e. the filling tubing and pipetting tubing must be rinsed and completely filled. The cylinder contents must then be ejected.



Command sequence:

```
MOVE 1 : ext.1
LIFT: 1 : work mm
DOS: 1.1 : prepar. mL
DOS: 1.1 : EndVol mL
DOS: 1.1 : compen. mL
```



*Important!*  
Eject the cylinder contents

This step should be carried out before every sample series involving pipetting. Before each pipetting process the Dosing unit cylinder contents must be completely ejected. This is the only way to achieve reproducible, accurate pipetting.

The titration cell can then be aspirated, rinsed and prepared with fresh solvent.

### 3.8.9 Pipetting

The pipetting procedure takes place in five steps:

- Ejection of cylinder contents (see above) and formation of separating bubble
- Move to sample
- Aspirate / Measure out sample
- Move to target
- Eject sample

### Formation of separating bubble (air gap)




Command sequence:

LIFT: 1 : shift mm  
DOS: 1.1 : -0.5 mL

In order to prevent the carrier liquid and sample solution from mixing a separating bubble, must be formed that occupies a length of at least 5 mm in the pipetting tubing. However, the separating bubble must also be measured out with sufficient accuracy. This is why a Dosing unit with max. 20 mL cylinder capacity should be used for pipetting.

### Move to sample




Command sequence:

MOVE: 1 : sample  
LIFT: 1 : work mm

The (dissolved or liquid) samples can be contained in open beakers on the rack or in sealed vials. In the latter case instead of the pipetting tubing a suitable injection needle should be connected to Dosino Port 1 with a normal piece of tubing (2 mm dia.).

### Aspirate sample



Command sequence:


DOS: 1.1 : -5.0 mL  
WAIT: pause 3 s

LIFT: 1 : shift mm  
DOS: 1.1 : -0.2 mL

Aspiration of the sample should be carried out at a reduced filling rate (<10 mL/min). The appropriate setting can be made under **<PARAM> >Dosing unit def.**

With non-aqueous samples an additional small air bubble may be drawn into the pipetting tip in order to prevent the sample solution from dripping.


### Move to target



Before the pipetting tubing containing the measured out sample is moved, the corresponding Dosino port must be closed.


Command sequence:

DOS: 1.2 : port mL  
MOVE 1 : ext.1  
LIFT: 1 : work mm



The titration cell should be filled with solvent before the sample is added so that the sample can be pipetted directly into the liquid (tip immersed).

### Eject sample



Command sequence:

**DOS: 1.1 :**      **EndVol mL**

The sample should be ejected at a reduced dosing rate (<10 mL/min). The appropriate settings can be made under **<PARAM> >Dosing unit def.**

### Rinse



Command sequence:

**DOS: 1.2 :**      **fill mL**  
**DOS: 1.1 :**      **endVol mL**  
**DOS: 1.1 :**      **compen. mL**

or

**DOS: 1.2 :**      **-x.xx mL**  
**DOS: 1.1 :**      **endVol mL**  
**DOS: 1.1 :**      **compen. mL**

In order to keep sample solution carry-over as low as possible, rinsing with solvent can be carried out after the sample has been ejected. This is particularly recommended for non-aqueous samples. In preparation for the next pipetting process the cylinder should be completely ejected.

## 3.9 The Remote interface

Connected peripheral devices such as Titrinos, Titrandos, pH meters, etc. can be controlled via the Remote interface (25-pin socket).

Fourteen lines are available for the output of signals (Output 0–13). Eight lines (Input 0–7) are available for receiving signals (e.g. the "Ready" signal from a Titrino at the end of a titration).

Metrohm instruments should only be connected with the Metrohm Remote cables provided, see Section 2.3.1.

### 3.9.1 Output lines

The 14 output lines of the Remote socket can be set as required in both manual operation as well as in a run sequence with the **Control** command (**CTL**). This is done by setting a 14-place binary pattern in which each bit is assigned to an output line.

Output	13	12	11	10	9	8	7	6	5	4	3	2	1	0
Bit	13	12	11	10	9	8	7	6	5	4	3	2	1	0

(bits are always numbered from right to left)

Example: **CTL Rm \*\*\*\*\*1\***  
 sets output line 1 to active (= set), for example, with a connected Titrino this would be a Stop command. 0 resets the output line to inactive.

We recommend that non-relevant output lines are masked with an asterisk (\*) in order that the condition of these lines are not altered.

### 3.9.2 Input lines

The 8 input lines of the Remote socket can be scanned during a method sequence with the **SCAN** command (**SCN**). The method sequence will be halted until the given binary pattern coincides with the effective condition of the input lines (e.g. the status of the 'Ready' line for scanning the end of a titration from a Titrino). This is done by setting an 8-place binary pattern, in which each bit is assigned to an input line. When the patterns coincide the run sequence will continue with the next command line. In manual operation the **SCAN** command is used for displaying the status of all input lines.

Input	7	6	5	4	3	2	1	0
Bit	7	6	5	4	3	2	1	0

(bits are always numbered from right to left)

Example: **SCN Rm \*\*\*\*\*1**  
 is waiting for an active input on line 0 (1=set). For example, this line is set by a Titrino after a titration has been completed and the Titrino can again accept a start signal.

Input lines that are not of interest or for which no defined condition can be predicted should also be masked with an asterisk (\*).



If **external START** (see configuration 3.2.1) is switched on then Remote lines Input 7 and Input 6 will be reserved for external START and STOP.

In order to simplify the use of Remote commands, particularly when several devices are linked together with Metrohm cables, predefined binary patterns are available as command parameters for the commands **CTL** and **SCN** for standard conditions (1 to 2 Titrinos, possibly with auxiliary Dosimats, pH meter or ion meter). These are:

### 3.9.3 SCN command

Parameter	Signal pattern	Function
Ready1	*****1	waits for "ready" condition of device 1 (Titrino, Titrande)
Ready2	**1*****	waits for "ready" condition of device 2
Ready*	**1****1	waits for "ready" condition of devices 1 and 2
Cond ok	*****1*	waits for "cond ok" condition of device 1
End1	****1***	waits for the EOD impulse from device 1 (Titrino, Titrande)
End2	*1*****	waits for the EOD impulse from device 2
EndMeter	***11***	waits for the EOD impulse from ion meter or pH-Meter (stirrer 1 will be switched on during the waiting period)
Continue	***1****	waits for an impulse to continue e. g. "Sample ready"

With the parameter **Ready\*** it is possible to scan the condition of devices working in parallel. In this case the **Ready** line of both devices must be statically set at the end of a determination. Devices which only transmit a brief impulse at the end of e.g. a measurement cannot be controlled in parallel.

### 3.9.4 CTL command

Parameter	Signal pattern	Function
INIT	000000000000	initializes the Remote interface
START device1	*****1	starts device <sup>1</sup> (e.g. Titrino, Titrande, ...) *)
START device2	*****1*****	starts device 2 (see above, only with multiple cable)*)
START device*	*****1****1	starts devices 1 and 2 " *)
SAMPLE ready	*****1****	sends a pulse for continuing a determination
START Dos1	*****1*****	starts Dosimat at device 1 (Titrino via "activate")
START Dos2	*****1*****	starts Dosimat at device 2 "
START Dos*	*****1*1*****	starts Dosimat at devices 1 and 2 "
METER Mode pH	*****0001*	switches ion meter or pH meter to pH measurement and starts it
METER Mode T	*****0010*	switches ion meter or pH meter to temperature measurement and starts it
METER Mode U	*****0011*	switches ion meter or pH meter to mV measurement and starts it

<b>METER Mode I</b>	<b>*****0100*</b>	switches ion meter or pH meter to Ipol (mV measurement) and starts it
<b>METER Mode C</b>	<b>*****1000*</b>	switches ion meter to Conc measurement
<b>METER Cal pH</b>	<b>*****0101*</b>	switches ion meter or pH meter to pH calibration and starts it
<b>METER Cal C</b>	<b>*****1001*</b>	switches ion meter to Conc calibration
<b>METER enter</b>	<b>*****1111*</b>	simulates the <b>&lt;ENTER&gt;</b> key for ion meter or pH meter (essential for pH calibration with 691/780/481, in order to start the measurement of the second buffer)

For **START** commands the signal is transmitted as a short impulse of 200 ms.

\*) With pH meters or ion meters a result report is triggered

### 3.9.5 Manual stop options

Parameter	Signal pattern	Function
<b>STOP Device 1</b>	<b>*****1*</b>	stops device 1 (e.g. Titrino, Titrand...)
<b>STOP Device 2</b>	<b>*****1*****</b>	stops device 2 (see above, only with multiple cable)
<b>STOP Device *</b>	<b>*****1*****1*</b>	stops devices 1 and 2 "

For STOP commands the signal is transmitted as a short impulse of 200 ms.

## 3.10 LEARN mode

The **LEARN** mode is used for three types of application:

- Interactive setting of lift and robotic arm positions
- Rack adjustment
- Interactive parametrizing sequence commands

### 3.10.1 Setting lift and robotic arm positions

The fine adjustment of defined positions, such as the working height of a lift, can be made by entering the height in mm. But it is recommended to use the LEARN function to run a lift or a robotic arm to a certain position by using the arrow keys of the keypad.

#### Example: procedure for setting the 'Work height at Tower 1'

- Use manual operation to move to a rack position (at Tower 1).
- Open configuration menu (with the **<CONFIG>** key).
- Open submenu **>rack definitions** and load the data of the rack which is currently on the turntable.
- Select **work position T1**.

#### LEARN mode

The quicker an arrow key is repeatedly pressed, the higher the speed of the lift or the longer the particular movement path.

- Press **<LEARN>** key.
- Use the **<↓>** and **<↑>** keys to move the lift to the required position.
- Accept the lift position with the **<ENTER>** key.
- If necessary modify the position (in mm) by entering numbers.

#### Example: procedure for setting 'External position 1 at Tower 1'

- Use manual operation to move to a sample beaker (at Tower 1) and move the lift to a suitable lift position.
- Open configuration menu (with the **<CONFIG>**) key.
- Select submenus **>tower 1 / swing head 1** and **external pos.1**.

#### LEARN mode

The quicker an arrow key is repeatedly pressed, the higher the swing rate or the larger the particular swing angle.

- Press **<LEARN>** key.
- Use the **<←>** and **<→>** keys to move the robotic arm until it is in the required position.
- Accept the robotic arm position with **<ENTER>**.
- If necessary modify the position (in °) by entering numbers.

### 3.10.2 Rack adjustment

If necessary, each sample rack can be finely adjusted, i.e. the rack offset is determined in the direction of rotation. A requirement for this is that the working height has already been set for the particular rack and that the rack settings have been stored.

Procedure:

- Position sample rack on turntable and initialize with **<RACK>**.
- Open configuration menu (with **<CONFIG>**) and select sub-menu **>rack definitions**.
- Under **>>recall rack** call up the definitions for the current rack.
- Select **rack offset**.
- Press **<LEARN>** and confirm the question **adjust RACK ?** with **<ENTER>**.
- The rack rotates to rack position 1.
- Confirm the question **adjust tower 1 ?** with **<ENTER>**.
- The lift moves to the working position.
- Use the arrow keys to exactly align rack position 1 with the beaker sensor at Tower 1.
- Confirm the adjustment with **<ENTER>**.
- The rack offset (in °) is accepted in accordance with the adjustment carried out.

### 3.10.3 Parametrizing sequence commands

As when editing a method the parameters of a command can most easily be altered interactively, i. e. by altering them manually, particular commands are "adaptive". When editing a sequence the **LEARN** function allows certain commands to be carried out via manual operation. The resulting parameter (e.g. the lift position or the status of the input lines of the Remote interface) can then be adopted in the current command line. The **LEARN** function can be used repeatedly. Once times or volumes have been "learned" the values will always be added up. This is extremely useful for the determination of the pumping time, where the optimal duration of the rinsing process can be determined interactively in this way.

#### Procedure for editing methods:

- Enter command or select existing command line
- Press **<LEARN>** key
- The function starts, the "**LEARN**" LED lights up
- Press **<LEARN>** key
- The function is stopped, the "**LEARN**" LED blinks

- Accept the value with the **<ENTER>** key (or start the **LEARN** function again)
- The **LEARN** LED goes out; the next command line appears

The **LEARN** function is available for the following commands:

<b>Command</b>	<b>Adaptive parameter</b>	<b>Type of function</b>
LIFT	Lift position in mm	absolute
PUMP	Pumping time in s	additive
STIR	Stirring time in s	additive
WAIT	Waiting time in s	additive
DOS	Dosing volume in mL	additive
SCN Rm	Status of 8 Remote lines	"live" value
SCN RS	Received character string	"live" value

## 3.11 TRACE function

The comfortable **TRACE** function is available for testing the run of a command sequence.

The **TRACE** function is a valuable aid for testing either a complete method or parts of a method step by step for test purposes. Each command line in a sequence can be executed directly by pressing the **<START>** key. When the command has been executed the next command line is shown.

"Tracing" can be carried out directly after a command line has been entered or at any other time after the parameter menu has been opened and a sequence has been selected.

## 3.12 Disabling keypad functions

Certain user dialog areas can be made inaccessible to untrained users. Various dialog areas or keys can be disabled. For example, this can be used to prevent a method from being overwritten accidentally or parameters from being altered.

The **>keyboard options** menu for the corresponding settings is opened by keeping the **<CONFIG>** key pressed down when the Sample Processor is switched on. This menu can still be accessed even if the whole keypad has been disabled.

The individual key functions that can be disabled are:

### 3.12.1 Disable whole keypad

In routine work, if only a single particular method is used, it may be advisable to prevent manual manipulations on the Sample Processor. For this reason it is possible to disable (almost) all the keypad keys. The

**<START>**, **<STOP>** and **<CLEAR/RESET>** keys always remain accessible, so that it is always possible to start and stop methods.

**lock keyboard:** **ON** disables all the keypad keys (with the above exceptions).

### 3.12.2 Disable configuration

The configuration of the Sample Processor can be protected against being overwritten. All the settings in the configuration menu and its sub-menus can then no longer be accessed.

**lock <CONFIG>:** **ON** disables the **<CONFIG>** key.

### 3.12.3 Disable parameter

If user-defined methods are normally used then it may be advisable to prevent the stored method parameters from being overwritten. The parameter menu and its sub-menus can be made inaccessible.

**lock <PARAM>:** **ON** disables the **<PARAM>** key.

### 3.12.4 Disable method storage functions

It is extremely advisable to prevent stored methods from being accidentally deleted. It should only be possible to delete methods by deliberately switching off the disabling function.

**>>user methods + <ENTER>** opens the submenu for disabling the method storage functions.

**lock method recall:** **ON** prevents the loading of methods.

**lock method store:** **ON** prevents methods from being stored.

**lock method delete:** **ON** prevents methods from being deleted.

### 3.12.5 Disable display

If the Sample Processor is only to be operated by using external control software, the display for manual operation can be switched off.

**lock display:** **ON** switches off the display.

## 3.13 786 Swing Head settings

If a 786 Swing Head with robotic arm is to be mounted then its configuration data must be entered in the setup menu of the Sample Processor before. The most important data:

- Swing radius (=length of robotic arm)
- Swing arm offset
- Max. swing angle
- Swing direction

### Setup menu of the Sample Processor

Procedure:

- Switch off instrument
- Keep **<CONFIG>** key pressed down and switch on the instrument again
- In the menu **setup >changer setup** select the submenu **>>swing head 1** or **>>swing head 2**.
- Make the following settings:

### Axis distance

The axis distance is the horizontal distance between the axis of rotation of the rack and the swing axis of the robotic arm.

Default settings:

**axial distance 166.00 mm** (for model 778)  
**axial distance 196.00 mm** (for model 789)

### Swing offset

The swing offset is the physical angle offset of the specific robotic arm model, see the Instructions for Use for the 786 Swing Head.

Default values:

**swing offset 0.00°** (for titration robotic arm)  
**swing offset 8.00°** (for transfer robotic arm)

### Max. swing angle

The maximum swing angle is the useable swinging range (relative angle). The starting and finishing positions of this range (as absolute angle positions) are determined by the swing offset (see above) and the Swing Head drive. Because of its construction, each robotic arm model has a different

value for the max. swing range, see the Instructions for Use for the 786 Swing Head. This value can also be reduced if necessary.

Default values:

**Max. swing angle** **84.00°** (for titration robotic arm)

**Max. swing angle** **117.00°** (for transfer robotic arm)

### Swing radius

The swing radius depends on the length of the robotic arm and, together with the axial distance (see above), is the most important parameter for accurately moving to a rack position. The distance from the robotic arm axis to the middle of the processing head at the front of the robotic arm is decisive. The various robotic arm models have different swing radii, see the Instructions for Use for the 786 Swing Head.

Default values:

**swing radius** **110 mm** (for titration robotic arm)

**swing radius** **112 mm** (for transfer robotic arm)

### Rotation offset

The rotation offset does not normally need to be altered. It is only used if a Swing Head is to be mounted on the tower with a lateral displacement.

Default value:

**rot. offset** **0.00 mm** (do not alter)

### Swing direction

In principle the swing direction of the robotic arm can be selected as required. If two Swing Heads are mounted on a 2-tower model then you must take care that the two robotic arms do not come into conflict. For this reason the robotic arm at Tower 1 should always be mounted so that it swings to the right and that at Tower 2 to the left.

Right-swinging mounting means: swing direction –

Left-swinging mounting means: swing direction +

Default values:

**Swing direction:** – (at Tower 1, swing to the right)

**Swing direction:** + (at Tower 2, swing to the left)

### Adjustment rate

When the **RACK** command is carried out after switch-on (rack initialization) the robotic arm is automatically adjusted to the zero position. If the highest possible degree of positioning accuracy of the robotic arm is required the rate of the adjustment process can be reduced. However, this also lengthens the duration of the automatic adjustment process, but increases the precision.

Default value:

**Adjust rate.:**      **normal**

After the settings for the robotic arm have been entered the **<QUIT>** key must be pressed three times. The settings become active the next time that the instrument is switched on.

## 4 Service, maintenance, errors

### 4.1 Service

Servicing the Sample Processor should be carried out within the framework of an annual service visit by technicians from Metrohm or one of its local agencies. If aggressive or corrosive chemicals are used then shorter service intervals are necessary.

The Metrohm service department is always pleased to provide competent advice about the servicing and maintenance of all Metrohm instruments.

#### 4.1.1 Running time meter

The running time meter enables the service intervals to be based on the effective number of hours the instrument has been in use, See p. 57f.

### 4.2 Care and maintenance

Not just sensitive instruments, but even a Sample Processor requires a certain amount of care and attention. If the instrument becomes excessively dirty this could interfere with its functions and shorten the working life of its basically robust mechanism and electronics.

Excessive dirt on the working heads can influence the measuring results obtained. Regular cleaning of exposed parts can prevent this to a large extent.

Spilled chemicals and solvents should be cleaned up immediately. Above all, the connection strip (the mains plug in particular) must be protected against contamination.

Although this prevented to a large extent by constructive measures, if aggressive media should penetrate the interior of the instrument then the mains plug should be pulled out immediately in order to prevent massive damage to the electronics. If such damage should occur please contact the Metrohm service department.

The instrument must not be opened by untrained personnel.



* <b>invalid rack code</b>	The rack code read in by the Sample Processor cannot be found in the internal table.
* <b>missing beaker</b>	After a <b>MOVE</b> command no beaker can be recognized at the selected position.
* <b>rack data missing</b>	No sample rack in position, or no rack data can be found for the sample rack currently on the turntable.
* <b>RS232 error</b>	The transmission parameters of the RS232 interface do not coincide with those of the receiving device.
* <b>SCAN timeout</b>	The connected device has not transmitted the expected signal within the defined timeout time. The sample determination may not have been carried out properly, or the connection has been interrupted. Check the connected device.
* <b>service recommended</b>	The warning limit of the running time meter has been reached. It is time for the Sample Processor to be serviced. Please contact your local Metrohm service department.
* <b>user memory full</b>	The memory for the user-defined methods is full. Before new methods can be stored methods that are not used, or only used infrequently, must be deleted.
* <b>wrong rack</b>	The rack placed on the turntable is not the rack assigned to the method under <b>parameters</b> .
<b>trap error xxx</b>	Unexpected program error; switch instrument off and then on again.
<b>No display, LEDs Tower 1 and Tower 2 light up</b>	LCD error (system error 7). Contact Service.

# 5 GLP validation – diagnosis

## 5.1 Validation / GLP

**GLP** (Good Laboratory Practice) requires, among other things, that the precision and correctness of analytical instruments is checked at regular intervals by using SOPs (Standard Operating Procedures).

As this instrument is not a measuring instrument as such, we recommend that the Sample Processor is included in an analytical system and is validated as part of the all-embracing validation of the whole analytical system.

If the Sample Processor is primarily used for titration purposes then it is advisable to carry out the validation of the titrator with the help of the Sample Processor. In this way any possible interference (e.g. carryover of titrant or sample solutions) that could influence the results is included within the framework of the evaluation of the complete titration system.

Checking the electronic and mechanical assemblies of Metrohm instruments can and should be undertaken within the framework of regular servicing by Metrohm technicians. All Metrohm instruments are equipped with start-up check routines which check that the relevant assemblies are functioning perfectly when the instrument is switched on. If no error message appears the instrument is functioning properly. Metrohm also supplies its instruments with built-in diagnosis programs which allow the service technicians to check the functioning of particular assemblies should faults or malfunctions occur and to localize them. Diagnosis programs can also be incorporated in a validation method.

Guidelines for drawing up standard operating procedures for checking a titration system can be found in Application Bulletin 252/1 ("Validation of Metrohm titrators in accordance with GLP/ISO9001"). This can be obtained free of charge from Metrohm.

## 5.2 Initializing the working memory

With this diagnostic step the default values for instrument parameters can be entered via the keypad and the instrument returned to its original condition. This measure is important with respect to the following:



*In rare cases it could happen that massive interference signals such as current peaks, lightning, etc. could affect the contents of the data memory. If the data memory contents are undefined then this could lead to a system crash.*

The Sample Processor has different ways of initializing the working memory. Either the complete data memory (**all**) or only parts of it (**param, config, setup, assembly**) can be overwritten with the default values.




*Although the instrument number is retained, initialization should only be carried out when absolutely necessary, as this deletes the stored user data (etc.).*

- Keep key <9> pressed down when the instrument is switched on.

```
diagnosis
>RAM initialization
```

- Press the <ENTER> key to open the following diagnosis menu:



The submenus are selected in sequence by pressing the <Select> key. The individual initialization versions are accessed by the <ENTER> key and exited with the <QUIT> key.

The table shows which parts of the working memory are affected by the corresponding initialization versions. If a system crash occurs (undefined display, no reaction to pressing keys, etc.) we recommend the initialization version "**a11**".

- If necessary, press the **<SELECT>** key several times until this appears:

```
>RAM initialization
select:          a11
```

- Press **<ENTER>**.

```
diagnosis
>RAM test
```

- Press **<ENTER>**.
- Press **<QUIT>**.

The instrument finishes the diagnosis process and carries out a start-up reset.

# 6 Annex

## 6.1 Technical data

### 6.1.1 Keypad

*LC display* 2 lines with 24 characters each  
character height 5 mm

*LED display* 3 LEDs  
*Keypad* Keypad with 30 keys

### 6.1.2 Interfaces

*RS232 socket* For connecting a computers or printer, 9 pins

*Remote socket* Universal parallel interface, programmable for communication with external devices, 22 signal lines (8x input, 14x output), TTL level




Fig. 37 Remote interface

+5 V: max. load = 20 mA

### 6.1.3 MSB connections

*3 connection sockets* 9-pin mini-DIN sockets  
- for 800/700 Dosino or 805 / 685 Dosimat  
- for 801 Stirrer or 804 Ti Stand

### 6.1.4 Pumps and pump connections

Various versions depending on the model:

#### Built-in pumps with valves

*Capacity* > 450 mL/min  
head 2 m

**Pump connections**

*Output* M8 socket:  
 $U = 16 \pm 1 \text{ V}$   
 $I \leq 0.8 \text{ A}$

for 823 Membrane Pump Unit or 772 Pump Unit

**6.1.5 Swing Head connection**

*Connection socket* 9-pin mini-DIN socket  
 - for 786 Swing Head

**6.1.6 Lift**

*Max. stroke path* 235 mm  
*Max. load* approx. 30 N  
*Speed* adjustable, 5...25 mm/s

**6.1.7 Turntable**

*Speed of rotation* adjustable, 3...20 degrees/s

**6.1.8 Stirrer connection (DIN socket)**

*Stirring rate* adjustable in 15 steps  
 741 Magnetic Stirrer 180/min...2600/min  
 802/722 Rod Stirrer 180/min...3000/min

**6.1.9 Mains connection**

*Voltage* 100... 240 V ( $\pm 10\%$ )  
*Frequency* 50...60 Hz  
*Power consumption* 115 W  
*Fuses* 2.0 ATH

**6.1.10 Safety specifications**

*Construction and testing* According to EN/IEC 61010-1, UL 3101-1 protection class I

*Safety instructions* The Instructions for Use contain safety information that must be observed by the user in order to ensure the safe operation of the instrument.

**6.1.11 Electromagnetic compatibility (EMC)**

*Emission* Standards fulfilled :  
 - EN/IEC 61326  
 - EN 55022 / CISPR 22  
 - EN/IEC 61000-3-2

*Immunity* Standards fulfilled:  
 - EN/IEC 61326  
 - EN/IEC 61000-4-2  
 - EN/IEC 61000-4-3  
 - EN/IEC 61000-4-4

- EN/IEC 61000-4-5
- EN/IEC 61000-4-6
- EN/IEC 61000-4-8
- EN/IEC 61000-4-11
- EN/IEC 61000-4-14
- NAMUR

### 6.1.12 Ambient temperature

<i>Nominal working range</i>	+5...+45 °C (at max. 80 % relative humidity)
<i>Storage</i>	-20 ... +60 °C relative humidity < +40 °C < 95% relative humidity < +50 °C < 85% relative humidity < +60 °C < 50%
<i>Transport</i>	-20 ... +60 °C <b>relative humidity</b> < +40 °C < 95% relative humidity < +50 °C < 85% relative humidity < +60 °C < 50%

### 6.1.13 Dimensions and materials

<i>Height</i>	730 mm
<i>Width</i>	280 mm
<i>Depth</i>	500 mm (2.778.XXXX model) 530 mm (2.789.XXXX model)
<i>Weight</i> (without accessories)	14.5 kg (2.7XX.0010 model) 15.4 kg (2.7XX.0020 model) 13.8 kg (2.7XX.0030 model) 18.9 kg (2.7XX.0110 model) 19.9 kg (2.7XX.0120 model) 16.8 kg (2.7XX.0130 model)
<i>Materials</i>	
- <i>Sample Processor housing</i>	Metal housing, surface-treated
- <i>Keypad housing</i>	Crastin (PBTB), inside Al-metallized
- <i>Keypad film</i>	Polyester, chemical-resistant

## 6.2 Standard methods

The following pages contain a list of the user methods supplied together with explanations of important commands. A requirement for the use of these methods is that the system is properly configured. In particular, the working position, rinsing position and shifting position of Lift 1 must be defined for the sample rack being used as well as a special beaker position on the rack.

The titration or measuring methods must always be set on the corresponding instrument itself. The correct cable connections are given on pages 31ff. In the methods described as examples it is assumed that the instruments are always connected to the Remote sockets.

We recommend that each new method is worked through step by step with the **TRACE** function before it is started for the first time and adapted to suit the particular requirements. This is done by using the **<PARAM>** key to open the loaded method and selecting **>start sequence**, **>sample sequence** or **>final sequence**. The displayed command line of a sequence can be executed directly with the **<START>** key. This means that it is possible to test the complete method sequence step by step.

If a standard method is not required then it can be deleted (see p. 87). Then more memory is available for storing methods.

### 6.2.1 Titrino

This is a universal method for carrying out automated titrations directly in the sample beaker. It is extremely suitable for use as a template for further methods.

#### Required instruments and cables:

- Metrohm titrator (Titrino or Titrande with 6.2148.010 Remote Box)
- Metrohm 802 Rod Stirrer or 741 Magnetic Stirrer
- 6.2141.020 Remote cable

The method can be used with or without a 786 Swing Head on the Sample Processor.

#### Special features:

After the titration the electrode is rinsed in the sample beaker (rinsing position). The rinsing position at Lift 1 must be set in the rack configuration so that the electrode is not immersed in the sample solution. An internal or external pump is required for rinsing.

<b>789 Sample Proc. XL</b>	<b>11122 5.789.0011</b>	- Report header with instrument ID and program version
<b>parameters</b>		
<b>method</b>	<b>Titrino</b>	- Method name
<b>number of samples:</b>	<b>rack</b>	- No. of samples (here whole sample rack)
<b>&gt;start sequence</b>		- Initialize Remote interface
<b>1 CTL:Rm:</b>	<b>INIT</b>	
<b>&gt;sample sequence</b>		
<b>1 MOVE 1 :</b>	<b>sample</b>	- Move to sample
<b>2 LIFT: 1 :</b>	<b>work mm</b>	- Run lift to working height
<b>3 STIR: T1 :</b>	<b>ON s</b>	- Switch on stirrer at Tower 1
<b>4 CTL:Rm: START</b>	<b>device1</b>	- Start titrator
<b>5 SCN:Rm :</b>	<b>End1</b>	- Wait for end of titration
<b>6 STIR: T1 :</b>	<b>OFF s</b>	- Switch off stirrer at Tower 1
<b>7 LIFT: 1 :</b>	<b>rinse mm</b>	- Run lift to rinsing height
<b>8 PUMP 1.1 :</b>	<b>5 s</b>	- Rinse electrode for 5 seconds
<b>9 WAIT: pause</b>	<b>3 s</b>	- Allow electrode to drip
<b>&gt;final sequence</b>		
<b>1 MOVE 1 :</b>	<b>spec.1</b>	- Move to conditioning beaker
<b>2 LIFT: 1 :</b>	<b>work mm</b>	- Run lift to working height
<b>&gt;changer settings</b>		----- Settings for Sample Processor functions -----
<b>rack name:</b>	<b>*</b>	
<b>lift rate T1</b>	<b>25 mm/s</b>	
<b>lift rate T2</b>	<b>25 mm/s</b>	
<b>shift rate</b>	<b>20°/s</b>	
<b>shift direction:</b>	<b>auto</b>	
<b>swing rate T1</b>	<b>55°/s</b>	
<b>swing rate T2</b>	<b>55°/s</b>	
<b>on beaker error:</b>	<b>MOVE</b>	
<b>&gt;stirring rates</b>		----- Stirring rates -----
<b>stirrer tower 1</b>	<b>3</b>	
<b>stirrer tower 2</b>	<b>3</b>	
<b>stirrer MSB1</b>	<b>3</b>	
<b>stirrer MSB2</b>	<b>3</b>	
<b>stirrer MSB3</b>	<b>3</b>	
<b>&gt;Dosing unit def.</b>		----- Settings for dosing devices-----
<b>&gt;timeout settings</b>		----- Settings for interface timeout -----
<b>SCAN timeout:</b>	<b>OFF min</b>	
<b>on SCAN timeout:</b>	<b>error</b>	
<b>&gt;manual stop</b>		----- Reaction to manual stop -----
<b>CTL Rmt:</b>	<b>STOP device*</b>	- Stop titrator
<b>CTL RS232:</b>		
<b>PUMP:</b>	<b>cont.</b>	- Do not alter pump status

```

STIR T1:          cont.          - Do not alter stirrer status
STIR T2:          cont.          "
STIR MSB1:        cont.          "
STIR MSB2:        cont.          "
STIR MSB3:        cont.          "
-----

```

### 6.2.2 PIP ext

Standard method for pipetting samples into an external titration cell. Very suitable for standard titrations.

#### Required instruments and cables:

- Metrohm titrator (Titrino or Titrando with 6.2148.010 Remote Box)
- 786 Swing Head with robotic arm for sample transfer (6.1462.030 or 6.1462.040)
- 803 Magnetic Stirrer or 804 Ti Stand with 802 Rod Stirrer connected to MSB1
- 800 Dosino with 807 Dosing unit (recommended: 10 mL cylinder, 6.3032.210) and 6.1562.100 Pipetting tubing
- 6.2141.020 Remote cable

#### Special features:

Samples are pipetted directly from the sample beaker into the external titration cell. The external Position 1 with the working height for moving to the titration cell must be set in the tower configuration. No special beaker is required.

```

789 Sample Proc. XL      11122 5.789.0011 - Report header with device ID and program version
parameters
method                  PIP ext
number of samples:     rack
>start sequence
 1 CTL:Rm:              INIT
 2 MOVE 1 :             ext.1
 3 LIFT: 1 :            rinse mm
 4 PUMP 1.2 :           30 s
 5 DOS: 1.1 :           prepar. mL
 6 PUMP 1.* :           30 s
 7 LIFT: 1 :            shift mm
>sample sequence
 1 MOVE 1 :             ext.1
 2 LIFT: 1 :            work mm
 3 PUMP 1.2 :           ON s
 4 DOS: 1.1 :           endVol mL
 5 DOS: 1.1 :           compen. mL
 6 PUMP 1.2 :           OFF s
 7 LIFT: 1 :            shift mm
 8 DOS: 1.1 :           -0.5 mL
 9 MOVE 1 :             sample
10 LIFT: 1 :            work mm
11 DOS: 1.1 :           -5.0 mL
12 WAIT: pause         3 s
13 LIFT: 1 :            shift mm
14 DOS: 1.1 :           -0.2 mL
15 DOS: 1.2 :           port mL
16 MOVE 1 :             ext.1
17 STIR: MSB1 :         ON s
18 PUMP 1.1 :           10 s
19 LIFT: 1 :            work mm

```

20 DOS: 1.1 :	5.4 mL	- Eject sample (+ air-gap bubble )
21 WAIT: pause	3 s	- Wait for drips
22 DOS: 1.2 :	fill mL	- Refill solvent in dosing cylinder
23 DOS: 1.1 :	endVol mL	- Eject solvent into titration cell (rinse)
24 DOS: 1.1 :	compen. mL	- Eliminate slack
25 LIFT: 1 :	shift mm	- Run lift to shifting height
26 WAIT: pause	5 s	- Wait 5 seconds
27 CTL:Rm: START	device1	- Start titrator
28 SCN:Rm :	Ready1	- Wait for end of titration
29 PUMP 1.* :	30 s	- Rinse and aspirate titration cell
30 STIR: MSB1 :	OFF s	- Switch off stirrer
<hr/>		
>final sequence		
1 MOVE 1 :	ext.1	- Move robotic arm to titration cell
2 LIFT: 1 :	work mm	- Run lift to working height
3 DOS: 1.1 :	empty mL	- Empty dosing cylinder and tubing
4 LIFT: 1 :	shift mm	- Run lift to shifting height
5 PUMP 1.1 :	15 s	- Rinse and aspirate titration cell
>changer settings		
...		
>stirring rates		Stirring rates
>Dosing unit def.		Settings for dosing device
dosing drive	1	
dos.rate	10.0 ml/min	
fill.rate	10.0 ml/min	- Reduce dosing rate
dosing port	1	- Reduce filling rate
dosing 2 port	3	
filling port	2	
rinsing port	2	
preparation port	1	
drain port	4	
cock direction:	not over	
not over port	4	
>timeout settings		- Do not rotate to open Port 4 (in bottle)
SCAN timeout:	OFF min	- see above
on SCAN timeout:	error	
>manual stop		Actions after manual stop
CTL Rmt:	STOP device*	- Stop titrator
CTL RS232:		
PUMP:	OFF	
STIR T1:	cont.	- Switch off pump
STIR T2:	cont.	
STIR MSB1:	cont.	
STIR MSB2:	cont.	
STIR MSB3:	cont.	
-----		

### 6.2.3 KF\_ext.

Standard method for coulometric or volumetric Karl Fischer determinations. The samples are pipetted from sealed vials into an external conditioned coulometer or KF titration cell.

#### Required instruments and cables:

- Metrohm Coulometer (756 or 831) or
- Metrohm KF Titrator (e.g. 795 Titrino or 835 Titrand) with 6.2148.010 Remote Box
- 786 Swing Head with robotic arm for sample transfer (6.1462.030 or 6.1462.040)
- 803 Magnetic Stirrer or 804 Ti Stand with 802 Rod Stirrer connected to MSB1
- 800 Dosino with 807 Dosing unit (recommended: 10 mL cylinder, 6.3032.210) and 6.1562.100 Pipetting tubing
- 6.2141.020 Remote cable
- Sample transfer accessory set (6.5619.000)

#### Special features:

Samples are pipetted directly from the sample beaker into the external titration cell. Freshly titrated conditioning solution from the KF cell is used as the transport medium for pipetting.

The external Position 1 with the working height for moving to the titration cell must be set in the tower configuration.

<b>789 Sample Proc. XL</b>	<b>11122 5.789.0011</b>	- Report header with device ID and program version
<b>parameters</b>		
<b>method</b>	<b>KF_ext</b>	- Method name
<b>number of samples:</b>	<b>rack</b>	- No. of samples (here whole sample rack)
<hr/>		
<b>&gt;start sequence</b>		
<b>1 CTL:Rm:</b>	<b>INIT</b>	- Initialize Remote interface
<b>2 STIR: MSB1 :</b>	<b>ON s</b>	- Switch on stirrer
<b>3 CTL:Rm: START</b>	<b>device1</b>	- Start conditioning
<b>4 MOVE 1 :</b>	<b>ext.1</b>	- Move robotic arm to KF cell
<b>5 SCN:Rm :</b>	<b>Cond ok</b>	- Wait for conditioning
<b>6 WAIT: pause</b>	<b>60 s</b>	- Stabilization of "Cond." condition
<b>7 SCN:Rm :</b>	<b>Cond ok</b>	- Wait for conditioning
<b>8 LIFT: 1 :</b>	<b>work mm</b>	- Insert needle into KF cell
<b>9 DOS: 1.4 :</b>	<b>eject mL</b>	- Eject dosing cylinder into waste container
<b>10 DOS: 1.1 :</b>	<b>-5.0 mL</b>	- Aspirate conditioning solution into dosing cylinder
<b>11 DOS: 1.2 :</b>	<b>endVol mL</b>	- Eject dosing cylinder into waste container
<b>12 DOS: 1.1 :</b>	<b>-2.0 mL</b>	- Aspirate conditioning solution into pipetting tubing
<b>13 DOS: 1.2 :</b>	<b>eject mL</b>	- Eject dosing cylinder into waste container
<b>14 DOS: 1.2 :</b>	<b>endVol mL</b>	- Piston to max. volume mark
<b>15 MOVE 1 :</b>	<b>spec.1</b>	- Move robotic arm to conditioning beaker
<b>16 LIFT: 1 :</b>	<b>work mm</b>	- Run lift to working height
<b>17 SCN:Rm :</b>	<b>Cond ok</b>	- Wait for conditioning
<b>18 WAIT: pause</b>	<b>60 s</b>	- Stabilization of "Cond." condition
<b>19 SCN:Rm :</b>	<b>Cond ok</b>	- Wait for conditioning
<b>&gt;sample sequence</b>		
<b>1 MOVE 1 :</b>	<b>ext.1</b>	- Move robotic arm to KF cell
<b>2 DOSRATE 1</b>	<b>5</b>	
<b>3 FILLRATE 1</b>	<b>5</b>	
<b>4 DOS: 1.1 :</b>	<b>-0.1 mL</b>	- Aspirate air-gap bubble
<b>5 LIFT: 1 :</b>	<b>work mm</b>	- Pierce KF cell
<b>6 DOS: 1.1 :</b>	<b>-0.5 mL</b>	- Aspirate conditioning solution
<b>7 LIFT: 1 :</b>	<b>shift mm</b>	- Run lift to shifting height
<b>8 DOS: 1.1 :</b>	<b>-0.1 mL</b>	- Aspirate air-gap bubble

```

 9 MOVE 1 : sample
10 LIFT: 1 : work mm
11 DOS: 1.1 : -0.5 mL
12 WAIT: pause 10 s
13 DOS: 1.2 : port mL
14 MOVE 1 : ext.1
15 SCN:Rm : Cond ok
16 WAIT: pause 60 s
17 SCN:Rm : Cond ok
18 CTL:Rm: START device1
19 LIFT: 1 : work mm
20 DOS: 1.1 : endVol mL
21 DOS: 1.1 : compen mL
22 MOVE 1 : spec.1
23 LIFT: 1 : work mm
24 SCN:Rm : Cond ok
>final sequence
 1 CTL:Rm: *****1*
 2 WAIT: pause 1 s
 3 CTL:Rm: INIT
>changer settings
rack name: *
lift rate T1 25 mm/s
lift rate T2 25 mm/s
shift rate 20°/s
shift direction: auto
swing rate T1 55°/s
swing rate T2 55°/s
on beaker error: MOVE
>stirring rates
>Dosing unit def.
>timeout settings
SCAN timeout: OFF min
on SCAN timeout: error
>manual stop
CTL Rmt: STOP device*
CTL RS232:
PUMP: cont.
STIR T1: cont.
STIR T2: cont.
STIR MSB1: OFF
STIR MSB2: cont.
STIR MSB3: cont.
-----

```

- Move to sample beaker
- Pierce sample vial
- Aspirate sample volume
- Wait for pressure equilibration
- Close Port 1
- Move to KF cell
- Wait for conditioning
- Stabilization of "Cond." condition
- Wait for conditioning
- Start titration
- Pierce KF cell
- Eject sample volume
- Eliminate slack
- Move to parking position
- Move lift to work height
- Wait for end of titration

- Stop titrator
- Wait 1 second for stop impulse
- Initialize Remote interface

---

Changer settings

---

Stirring rates

---

Settings for dosing device

---

Actions after manual stop

- Stop titrator

### 6.2.4 pH\_cal

This is a universal method for automated pH measurements directly in the sample beaker. The electrode calibration is included in the start sequence.

#### Required instruments and cables:

- Metrohm pH Meter (model 691 or 780/781) with 6.2148.010 Remote box
- Metrohm 802 Propeller Stirrer or 741 Magnetic Stirrer
- 6.2141.020 Remote cable

The method can be used with or without a 786 Swing Head on the Sample Processor. In the rack definition the special beaker positions 1 to 3 must be reserved.

#### Special features:

The two buffer solutions for the electrode calibration must be in special beaker positions 1 and 2.

After each measurement the electrode is rinsed in a rinsing beaker at special beaker position 3. Internal or external pump required for rinsing.

If a 780/781 pH Meter is used, two additional command lines must be inserted in the start sequence of the method. This is to cancel the display of the calibration results:

```
17 SCN:Rm      :      end1
18 WAIT:  pause      1 s
19 CTL:Rm:     START device1
20 MOVE 1      :      spec.3
```

<b>789 Sample Proc. XL</b>	<b>11122 5.789.0011</b>	- Report header with device ID and program version
<b>parameters</b>		
<b>method</b>	<b>pH cal</b>	- Method name
<b>number of samples:</b>	<b>rack</b>	- No. of samples (here whole sample rack)
<b>&gt;start sequence</b>		
1 CTL:Rm:	INIT	- Initialize Remote interface
2 MOVE 1 :	spec.3	- Move to rinse beaker
3 LIFT: 1 :	work mm	- Lower electrode
4 PUMP 1.* :	4 s	- Rinse electrode and aspirate beaker
5 MOVE 1 :	spec.1	- Move to buffer 1
6 LIFT: 1 :	work mm	- Immerse electrode
7 STIR: T1 :	10 s	- Stir for 10 seconds
8 CTL:Rm:	METER cal pH	- Start calibration with 1st buffer
9 SCN:Rm :	End1	- Wait for calibration
10 MOVE 1 :	spec.3	- Move to rinse beaker
11 LIFT: 1 :	work mm	- Lower electrode
12 PUMP 1.* :	4 s	- Rinse electrode and aspirate beaker
13 MOVE 1 :	spec.2	- Move to 2nd buffer
14 LIFT: 1 :	work mm	- Immerse electrode
15 STIR: T1 :	10 s	- Stir for 10 seconds
16 CTL:Rm:	METER enter	- Start calibration with 2nd buffer
17 SCN:Rm :	End1	- Wait for calibration
18 MOVE 1 :	spec.3	- Move to rinse beaker
19 LIFT: 1 :	work mm	- Lower electrode
20 PUMP 1.* :	4 s	- Rinse electrode and aspirate beaker
<b>&gt;sample sequence</b>		
1 SHIFTRATE:+	20	- Rack rotation direction counterclockwise

2 MOVE 1	:	sample	- Move to sample
3 LIFT: 1	:	work mm	- Immerse electrode
4 STIR: T1	:	10 s	- Stir for 10 seconds
5 CTL:Rm:	METER mode pH		- Switch pH meter to measure
6 CTL:Rm:	START device1		- Start measurement
7 SCN:Rm	:	End1	- Wait for measurement
8 SHIFTRATE:-		20	- Rack rotation direction clockwise
9 MOVE 1	:	spec.3	- Move to rinse beaker
10 LIFT: 1	:	work mm	- Lower electrode
11 PUMP 1.*	:	4 s	- Rinse electrode and beaker
<b>&gt;final sequence</b>			
1 MOVE 1	:	spec.3	- Move to rinse beaker
2 LIFT: 1	:	work mm	- Lower electrode
<b>&gt;changer settings</b>			
rack name:		*	Changer settings
lift rate T1		25 mm/s	
lift rate T2		25 mm/s	
shift rate		20°/s	
shift direction:		auto	
swing rate T1		55°/s	
swing rate T2		55°/s	
on beaker error:		MOVE	
<b>&gt;stirring rates</b>			
<b>&gt;Dosing unit def.</b>			
<b>&gt;timeout settings</b>			
SCAN timeout:		OFF min	
on SCAN timeout:		error	Actions after manual stop
<b>&gt;manual stop</b>			
CTL Rmt:		STOP device*	- Stop pH meter
CTL RS232:			
PUMP:		OFF	- Pump off
STIR T1:		cont.	
STIR T2:		cont.	
STIR MSB1:		cont.	
STIR MSB2:		cont.	
STIR MSB3:		cont.	
-----			

### 6.2.5 Std\_add

This is a universal method for automated ion measurements directly in the sample beaker.

#### Required instruments and cables:

- Metrohm 692 Ion Meter oder 781 Ion Meter with 6.2148.010 Remote Box
- Metrohm 802 Propeller Stirrer or 741 Magnetic Stirrer
- 6.2141.020 Remote cable

The method can be used with or without a 786 Swing Head on the Sample Processor.

#### Special features:

After each measurement the electrode is rinsed in the sample beaker at the rinse position. Internal or external pump required for rinsing.

789 Sample Proc. XL parameters	11122 5.789.0011	- Report header with device ID and program version
method	std add	- Method name
number of samples:	rack	- No. of samples (here whole sample rack)
<b>&gt;start sequence</b>		
1 CTL:Rm:	INIT	- Initialize Remote interface

```

>sample sequence
  1 MOVE 1 : sample
  2 LIFT: 1 : work mm
  3 CTL:Rm: METER mode C
  4 CTL:Rm: START device1
  5 SCN:Rm : EndMeter
  6 LIFT: 1 : rinse mm
  7 PUMP 1.1 : 5 s
  8 WAIT: pause 3 s
>final sequence
  1 MOVE 1 : spec.1
  2 LIFT: 1 : work mm
>changer settings
  rack name: *
  lift rate T1 25 mm/s
  lift rate T2 25 mm/s
  shift rate 20°/s
  shift direction: auto
  swing rate T1 55°/s
  swing rate T2 55°/s
  on beaker error: MOVE
>stirring rates
>Dosing unit def.
>timeout settings
  SCAN timeout: OFF min
  on SCAN timeout: error
>manual stop
  CTL Rmt: STOP device*
  CTL RS232:
  PUMP: OFF
  STIR T1: cont.
  STIR T2: cont.
  STIR MSB1: cont.
  STIR MSB2: cont.
  STIR MSB3: cont.
  -----

```

- Move to sample
- Run lift to work height
- Switch ion meter to concentration measurement
- Start measurement
- Wait for measurement
- Runlift to rinse height
- Rinse for 5 seconds with pump
- Allow to drip for 3 seconds

- Move to rinse beaker
- runlift to work height

---

Changer settings

---

Stirring rates

---

Settings for dosing device

---

Actions after manual stop

---

- Stop ion meter







## 6.3 Standard equipment











On receipt of the instrument please check that the delivery is complete.

### 6.3.1 Metrohm Sample Processor:

**Order no. 2.789.0010 or 2.778.0010 (model versions with 1 tower and 1 pump)**

The following accessories are included as standard:





No.	Order no.	Description
1	1.789.0010 or 1.778.0010	<b>789 Robotic Sample Processor XL</b>  <b>778 Sample Processor</b> 1 working station with rinsing equipment.
2	6.1236.020	Sleeve SGJ14/12 mm 
5	6.1446.000	Plastic stopper SGJ14 
3	6.1446.010	Stopper SGJ9 
1	6.1446.040	Threaded stopper M6 
1	6.1621.000	PE canister 10 L Use as rinsing solution or waste container 
1	6.1805.110	FEP tubing with light and anti-kink protection, with 2x M6 connections L = 80 cm 










3	6.1805.420	FEP tubing with light and anti-kink protection, with 2x M6 connections L = 48 cm	
1	6.1812.000	PTFE tubing L = 400 cm	
1	6.1828.000	Connecting nipple for canister 6.1621.000	
1	6.2053.000	Cable clips 10x for mounting cables	
1	6.2142.040	Keypad for 778/789 Sample Processor	
1	6.2621.030	Hexagon key 4 mm	
1	6.2621.070	Hexagon key 5 mm	
1	6.2621.140	Hexagon key 2.5 mm	
3	6.2709.070	Guide sleeve, ETFE, for 6.1543.xxx, SGJ 9	
3	6.2740.020	Spray nozzle ETFE With valve and M6 thread	
1	6.2751.080 or	Splash protection/Safety shield (778 Model)	









1	6.2751.100	Splash protection/Safety shield (789 Model)
1	6.2122.xxx	Mains cable with coupling type IEC 320/C13  Plug as per customer's requirements: - Type SEV 12 (Switzerland...) 6.2122.020 - Type CEE(7), VII (Germany...) 6.2122.040 - Type NEMA/ASA (USA...) 6.2122.070
1	8.789.1063	Metrohm Sample Processor Instructions for Use , English
1	8.789.1013	Quick reference Metrohm Sample Processor, English
1	8.789.1023	Metrohm Sample Processor Operating Tutorial, English
1	8.789.1033	Technical Reference: Metrohm Sample Processors, English

**Order no. 2.789.0020 or 2.778.0020 (model versions with 1 tower and 2 pumps)**

The following accessories are included as standard:

No.	Order no.	Description
1	1.789.0020  1.778.0020	<b>789 Robotic Sample Processor XL</b> or <b>778 Sample Processor</b> 1 working station with rinsing and aspiration equipment.
2	6.1236.020	Sleeve SGJ14/12 mm 
5	6.1446.000	Plastic stopper SGJ14 
3	6.1446.010	Stopper SGJ9 
1	6.1446.040	Threaded stopper M6 







1	6.1446.160	<p>Stopper SGJ 14/6.6 mm</p> <p>Made of PTFE</p> <p>For positioning inserts obliquely in the titration head</p>	
1	6.1543.170	<p>Aspiration tip, PTFE</p> <p>With M8 thread</p>	
2	6.1621.000	<p>PE canister 10 L</p> <p>Use as rinsing solution or waste container</p>	
1	6.1805.110	<p>FEP tubing with light and anti-kink protection, with 2x M6 connections</p> <p>L = 80 cm</p>	
3	6.1805.420	<p>FEP tubing with light and anti-kink protection, with 2x M6 connections</p> <p>L = 48 cm</p>	
1	6.1805.510	<p>PTFE tubing with anti-kink protection</p> <p>With 2 tubing nipples M8, L = 60 cm</p> <p>Inner diameter: 3 mm</p>	
2	6.1812.000	<p>PTFE tubing</p> <p>L = 400 cm</p>	
1	6.1820.030	<p>Screw connector with M8 thread</p> <p>For 6.1812.XXX and 6.1805.XXX Tubing</p>	
2	6.1828.000	<p>Connecting nipple for canister 6.1621.000</p>	


1	6.2042.020	Buret tip clip	
1	6.2053.000	Cable clips 10x for mounting cables	
1	6.2142.040	Keypad for 778/789 Sample Processor	 1,5m
1	6.2621.030	Hexagon key 4 mm	
1	6.2621.070	Hexagon key 5 mm	
1	6.2621.140	Hexagon key 2.5 mm	
3	6.2709.070	Guide sleeve, ETFE, for 6.1543.xxx, SGJ 9	
3	6.2740.020	Spray nozzle, ETFE With valve and M6 thread	
1	6.2751.080	Splash protection/Safety shield	(778 Model)
1	6.2751.100	Splash protection/Safety shield	(789 Model)
1	6.2122.xxx	Mains cable with coupling type IEC 320/C13  Plug as per customer's requirements: - Type SEV 12 (Switzerland...) 6.2122.020	

		- Type CEE(7), VII (Germany...)	6.2122.040
		- Type NEMA/ASA (USA...)	6.2122.070
1	8.789.1063	Metrohm Sample Processor Instructions for Use , English	
1	8.789.1013	Quick reference Metrohm Sample Processor , English	
1	8.789.1023	Metrohm Sample Processor Operating Tutorial, English	
1	8.789.1033	Technical Reference: Metrohm Sample Processors, English	

**Order no. 2.789.0030 or 2.778.0030 (model versions with 1 tower, no pumps)**




The following accessories are included as standard:











No.	Order no.	Description	
1	1.789.0030 or 1.778.0030	<b>789 Robotic Sample Processor XL</b> or <b>778 Sample Processor</b> 1 working station	
2	6.1236.020	Sleeve SGJ14/12 mm	
1	6.1805.110	FEP tubing with light and anti-kink protection, with 2x M6 connections L = 80 cm	
1	6.2142.040	Keypad for 778/789 Sample Processor	
1	6.2621.030	Hexagon key 4 mm	
1	6.2621.070	Hexagon key 5 mm	
1	6.2621.140	Hexagon key 2.5 mm	




3	6.2709.070	Guide sleeve, ETFE, for 6.1543.xxx, SGJ 9	
1	6.2751.080	Splash protection/Safety shield	(778 Model)
	or		
1	6.2751.100	Splash protection/Safety shield	(789 Model)
1	6.2122.xxx	Mains cable with coupling type IEC 320/C13  Plug as per customer's requirements: - Type SEV 12 (Switzerland...) 6.2122.020 - Type CEE(7), VII (Germany...) 6.2122.040 - Type NEMA/ASA (USA...) 6.2122.070	
1	8.789.1063	Metrohm Sample Processor Instructions for Use , English	
1	8.789.1013	Quick reference Metrohm Sample Processor, English	
1	8.789.1023	Metrohm Sample Processor Operating Tutorial, English	
1	8.789.1033	Technical Reference: Metrohm Sample Processors, English	

**Order no. 2.789.0110 or 2.778.0110 (model versions with 2 towers and 2 pumps)**

The following accessories are included as standard:


No.	Order no.	Description	
1	1.789.0110 or 1.778.0110	<b>789 Robotic Sample Processor XL</b> or <b>778 Sample Processor</b> 2 working stations with rinsing equipment.	
4	6.1236.020	Sleeve SGJ14/12 mm	
10	6.1446.000	Plastic stoppers SGJ14	
6	6.1446.010	Stopper SGJ9	










2	6.1446.040	Threaded stopper M6	
2	6.1621.000	PE canister 10 L Use as rinsing solution or waste container	
2	6.1805.110	FEP tubing with light and anti-kink protection, with 2x M6 connections L = 80 cm	
6	6.1805.420	FEP tubing with light and anti-kink protection, with 2x M6 connections L = 48 cm	
2	6.1812.000	PTFE tubing L = 400 cm	
2	6.1828.000	Connection nipple for canister 6.1621.000	
2	6.2053.000	Cable clips 10x for mounting cables	
1	6.2142.040	Keypad for 778/789 Sample Processor	
1	6.2621.030	Hexagon key 4 mm	
1	6.2621.070	Hexagon key 5 mm	









1	6.2621.140	Hexagon key 2.5 mm	
6	6.2709.070	Guide sleeve, ETFE, for 6.1543.xxx, SGJ 9	
6	6.2740.020	Spray nozzle, ETFE With valve and M6 thread	
1	6.2751.090	Splash protection/Safety shield	(778 Model)
	or		
2	6.2751.100	Splash protection/Safety shield	(789 Model)
1	6.2122.xxx	Mains cable with coupling type IEC 320/C13  Plug as per customer's requirements: - Type SEV 12 (Switzerland...) 6.2122.020 - Type CEE(7), VII (Germany...) 6.2122.040 - Type NEMA/ASA (USA...) 6.2122.070	
1	8.789.1063	Metrohm Sample Processor Instructions for Use , English	
1	8.789.1013	Quick reference Metrohm Sample Processor, English	
1	8.789.1023	Metrohm Sample Processor Operating Tutorial, English	
1	8.789.1033	Technical Reference: Metrohm Sample Processors, English	



**Order no. 2.789.0120 or 2.778.0120 (model versions with 2 towers and 4 pumps)**

The following accessories are included as standard:

No.	Order no.	Description
1	1.789.0120 or 1.778.0120	<b>789 Robotic Sample Processor XL</b> or <b>778 Sample Processor</b> 2 working stations with rinsing and aspiration equipment.
4	6.1236.020	Sleeve SGJ14/12 mm  




10	6.1446.000	Plastic stopper SGJ14	
6	6.1446.010	Stopper SGJ 9	
2	6.1446.040	Threaded stopper M6	
2	6.1446.160	Stopper SGJ 14/6.6 mm Made of PTFE For positioning inserts obliquely in the titration head	
2	6.1543.170	Aspiration tip, PTFE With M8 thread	
4	6.1621.000	PE canister 10 L For Sample Processor use as rinsing solution or waste container	
2	6.1805.110	FEP tubing with light and anti-kink protection, with 2x M6 connections L = 80 cm	
6	6.1805.420	FEP tubing with light and anti-kink protection, with 2x M6 connections L = 48 cm	
2	6.1805.510	PTFE tubing with anti-kink protection With 2 tubing nipples M8, L = 60 cm Inner diameter: 3 mm	





4	6.1812.000	PTFE tubing L = 400 cm	
2	6.1820.030	Screw connector with M8 thread For 6.1812.XXX and 6.1805.XXX Tubing	
4	6.1828.000	Connection nipple for canister 6.1621.000	
2	6.2042.020	Buret tip clip	
2	6.2053.000	Cable clips 10x for mounting cables	
1	6.2142.040	Keypad for 778/789 Sample Processor	
1	6.2621.030	Hexagon key 4 mm	
1	6.2621.070	Hexagon key 5 mm	
1	6.2621.140	Hexagon key 2.5 mm	

6	6.2709.070	Guide sleeve, ETFE, for 6.1543.xxx, SGJ 9	
6	6.2740.020	Spray nozzle, ETFE With valve and M6 thread	
1	6.2751.090 or	Splash protection/Safety shield	(778 Model)
2	6.2751.100	Splash protection/Safety shield	(789 Model)
1	6.2122.xxx	Mains cable with coupling type IEC 320/C13  Plug as per customer's requirements: - Type SEV 12 (Switzerland...) 6.2122.020 - Type CEE(7), VII (Germany...) 6.2122.040 - Type NEMA/ASA (USA...) 6.2122.070	
1	8.789.1063	Metrohm Sample Processor Instructions for Use , English	
1	8.789.1013	Quick reference Metrohm Sample Processor, English	
1	8.789.1023	Metrohm Sample Processor Operating Tutorial, English	
1	8.789.1033	Technical Reference: Metrohm Sample Processors, English	


**Order no. 2.789.0130 or 2.778.0130 (model versions with 2 towers, no pumps)**












The following accessories are included as standard:











No.	Order no.	Description	
1	1.789.0130 or 1.778.0130	<b>789 Robotic Sample Processor XL</b> or <b>778 Sample Processor</b> 2 working stations	
4	6.1236.020	Sleeve SGJ14/12 mm	
2	6.1805.110	FEP tubing with light and anti-kink protection, with 2x M6 connections  L = 80 cm	
1	6.2142.040	Keypad for 778/789 Sample Processor	





1	6.2621.030	Hexagon key 4 mm	
1	6.2621.070	Hexagon key 5 mm	
1	6.2621.140	Hexagon key 2.5 mm	
6	6.2709.070	Guide sleeve, ETFE, for 6.1543.xxx, NS 9	
1	6.2751.090	Splash protection/Safety shield (778 Model)	
or			
2	6.2751.100	Splash protection/Safety shield (789 Model)	
1	6.2122.xxx	Mains cable with coupling type IEC 320/C13  Plug as per customer's requirements: - Type SEV 12 (Switzerland...) 6.2122.020 - Type CEE(7), VII (Germany...) 6.2122.040 - Type NEMA/ASA (USA...) 6.2122.070	
1	8.789.1063	Metrohm Sample Processor Instructions for Use , English	
1	8.789.1013	Quick reference Metrohm Sample Processor, English	
1	8.789.1023	Metrohm Sample Processor Operating Tutorial, English	
1	8.789.1033	Technical Reference: Metrohm Sample Processors, English	

### 6.3.2 Sample racks and sample beakers

Order no.	Description	
<b>6.2041.310</b>	<b>Sample rack 12 x 250 mL</b> (d=416 mm) for	
6.1432.320	Metrohm sample beakers, glass 250 mL	
6.1453.220	Metrohm sample beakers, PP 200 mL	
6.1453.250	Metrohm sample beakers, PP 250 mL	

6.2041.320	<b>Sample rack 16 x 150 mL</b> for Standard glass beakers (tall form)	
6.2041.340 6.1432.210 6.1458.020	<b>Sample rack 24 x 75 mL</b> (d=416 mm) for Metrohm sample beakers, glass 75 mL, Micro titration head required	
6.2041.350 6.1432.210	<b>Sample rack <sup>*)</sup> 48 x 75 mL</b> (d=416 mm) for direct titration in Metrohm sample beakers, glass 75 mL	
6.2041.360 6.1459.310	<b>Sample rack 12 x 150 mL</b> (d=416 mm) for 150 mL standard glass beakers (tall form) 200 mL disposable beakers PP (1000 pcs.)	
6.2041.370 6.1459.310	<b>Sample rack 14 x 200 mL</b> (d=416 mm) for 200 mL disposable beakers (EU) PP	
6.2041.380	<b>Sample rack 14 x 8 oz</b> (d=416 mm) for Disposable beakers (US) PP 8 oz	
6.2041.390	<b>Sample rack 16 x 8 oz</b> (d=416 mm) for Disposable beakers (US) PP 8 oz	
6.2041.400 6.1432.320 6.1453.220 6.1453.250	<b>Sample rack <sup>*)</sup> 126 x 15 mL and 2 x 250 mL</b> (d=416 mm) for test tubes 11 mL and Metrohm sample beakers, glass 250 mL or Metrohm sample beakers, PP 200 mL	
6.2041.410 6.2743.057	<b>Sample rack<sup>*)</sup> 141 x 11 mL und 1 x 500 mL</b> (d=416 mm) for 11 mL test tubes and 500 mL standard glass beakers	
6.2041.430 6.2743.057 6.1608.080	<b>Sample rack<sup>*)</sup> 127 x 11 mL und 2 x 300 mL</b> (d=416 mm) for 11 mL test tubes and 300 mL PE bottles	
6.2041.440 6.2743.057 6.1608.080	<b>Sample rack<sup>*)</sup> 148 x 11 mL und 3 x 300 mL</b> (d=416 mm) for 11 mL test tubes and 300 mL PE bottles	



6.2041.450  6.2743.057 6.2747.010	<b>Sample rack*) 56 x 11 mL + 56 x 50 mL</b> (d=476 mm) for 11 mL test tubes and 50 mL PP sample beaker	
6.2041.460  6.2058.010	<b>Sample rack 21 x 100 mL</b> (d=416 mm) for 100 mL COD tubes  Tower extension required	
6.2041.710  6.2419.000 6.1448.050	<b>Sample rack *) 160 x 6 mL</b> (d=416 mm) for 6 mL vials (1000 pcs.) with septum seals (1000 pcs.)	
6.2041.750  6.2743.057	<b>Sample rack *) 36 x 11 mL</b> (d=269 mm) for 11 mL test tubes	
6.2041.800  6.1432.210	<b>Sample rack *) 100 x 75 mL</b> (d=476 mm) For use with 789 XL and 815 XL only, for Metrohm sample beakers, glass 75 mL	
6.2041.810	<b>Sample rack *) 34 x 150 mL</b> (d=476 mm) For use with 789 XL and 815 XL only, for 150 mL standard glass beakers (tall form)	
6.2041.820  6.1432.320 6.1453.220 6.1453.250	<b>Sample rack *) 28 x 250 mL</b> (d=476 mm) For use with 789 XL and 815 XL only, for Metrohm sample beakers, glass 250 mL Metrohm sample beakers, PP 200 mL Metrohm sample beakers, PP 250 mL	
6.2041.830  6.1459.310	<b>Sample rack *) 28 x 200 mL</b> (d=476 mm) For use with 789 XL and 815 XL only, for Metrohm sample beakers, PP 200 mL	
6.2041.840  6.1459.300	<b>Sample rack *) 59 x 120 mL</b> (d=476 mm) For use with 789 XL and 815 XL only, for Metrohm sample beakers, PP 120 mL	
6.2840.000	<b>Water bath with thermostat connector *)</b> for 789 Robotic Sample Processor XL or 815 Robotic USB Sample Processor XL	

<p><b>6.2041.900</b> 6.1432.210</p>	<p><b>Rack insert for water bath 6.2840.000 *</b> <b>54 x 75 mL</b> for 75 mL Metrohm sample beakers, glass</p>	
<p><b>6.2058.010</b></p>	<p><b>Tower extension</b> for COD applications</p>	
<p><b>6.2058.020</b></p>	<p><b>Titration head extension (30 mm)</b></p>	
<p><b>6.2058.030</b></p>	<p><b>Titration head extension (9 mm) for water bath rack 6.2041.900 *</b></p>	



\*) 786 Swing Head necessary

### 6.3.3 786 Swing Head

#### Order no. 2.786.0010 (model version 'transfer' with pipetting accessories)



No.	Order no.	Description
<b>1</b>	<b>1.786.0010</b>	<b>786 Swing Head</b> High-precision drive for robotic arm 
1	6.1462.030	Robotic arm for sample transfer, for left-swinging assembly 
1	6.1562.100	Pipetting tubing 10 mL
1	6.1823.010	Guide shaft for pipetting tubing
1	6.2621.120	Hexagon key 1.5 mm
2	V.024.4012	Countersunk head screw M.I-6 KT.
1	8.786.1013	Instructions for Use for 786 Swing Head, english

#### Order no. 2.786.0020 (model version 'transfer' with pipetting accessories)


No.	Order no.	Description
<b>1</b>	<b>1.786.0010</b>	<b>786 Swing Head</b> High-precision drive for robotic arm 
1	6.1462.040	Robotic arm for sample transfer, for right-swinging assembly 
1	6.1562.100	Pipetting tubing 10 mL
1	6.1823.010	Guide shaft for pipetting tubing

1	6.2621.120	Hexagon key 1.5 mm
2	V.024.4012	Countersunk head screw M.I-6 KT.
1	8.786.1013	Instructions for Use for 786 Swing Head, english







**Order no. 2.786.0030 (model version with titration accessories)**

No.	Order no.	Description
<b>1</b>	<b>1.786.0010</b>	<b>786 Swing Head</b> High-precision drive for robotic arm
		
1	6.1462.050	Robotic arm with micro titration head
		
1	6.1909.030	Stirrer propeller PP (104 mm) KF
2	6.2042.030	Buret tip clip
1	6.2621.120	Hexagon key 1.5 mm
2	V.024.4012	Countersunk head screw M.I-6 KT.
1	8.786.1013	Instructions for Use for 786 Swing Head, english

**Order no. 2.786.0040 (model version with no robotic arm)**








No.	Order no.	Description
<b>1</b>	<b>1.786.0010</b>	<b>786 Swing Head</b> High-precision drive for robotic arm
		
1	6.2621.120	Hexagon key 1.5 mm
2	V.024.4012	Countersunk head screw M.I-6 KT.
1	8.786.1013	Instructions for Use for 786 Swing Head, english

**Robotic arms for 786 Swing Head**









Order no.	Description	
<b>6.1462.030</b>  6.1562.100	<b>Robotic arm 'transfer'</b> , L=112 mm swing-left mounting, for  10 mL pipetting tubing	
<b>6.1462.040</b>  6.1562.100	<b>Robotic arm 'transfer'</b> , L=112 mm swing-right mounting, for  10 mL pipetting tubing	
<b>6.1462.050</b>	<b>Robotic arm 'titration'</b> , L=110 mm swing-left and swing-right mounting  use with 75 mL sample beakers (6.1432.210) or greater	
<b>6.1462.060</b>  6.1458.010 6.1458.020 6.1458.030	<b>Robotic arm 'macro'</b> , L=127 mm swing-left mounting, use with  Titration head " "	
<b>6.1462.070</b>  6.1458.010 6.1458.020 6.1458.030	<b>Robotic arm 'macro'</b> , L=127 mm swing-right mounting, use with  Titration head " "	
<b>6.1462.080</b>  6.2037.050 6.2037.060	<b>Dis-Cover robotic arm</b> , L=112 mm, swing-left mounting, with magnetic contact, for decapping and covering of samples with  lid for 75 mL sample beakers or lid for 250 mL sample beakers	







<p><b>6.1462.090</b></p> <p>6.2816.100 2.1562.110</p>	<p><b>LQ robotic arm 'transfer'</b>, L=112 mm, swing-right mounting, with Luer lock adapter for</p>	
<p><b>6.1462.150</b></p>	<p><b>Robotic arm 'titration'</b>, with contact beaker sensor (piezo), L=110 mm swing-right/-left mounting, installation both ways possible</p> <p>use with 75 mL sample beakers (6.1432.210) or greater</p>	
<p><b>6.1462.160</b></p> <p>6.1458.010 6.1458.020 6.1458.030</p>	<p><b>Robotic arm 'macro'</b>, with contact beaker sensor (piezo), L=127 mm, swing-left mounting, use with</p> <p>Titration head " "</p>	
<p><b>6.1462.170</b></p> <p>6.1458.010 6.1458.020 6.1458.030</p>	<p><b>Robotic arm 'macro'</b>, with contact beaker sensor (piezo), L=127 mm, swing-right mounting, use with</p> <p>Titration head " "</p>	


**6.3.4 Connection cables**

<i>Order no.</i>	<i>Description</i>
<b>6.2141.020</b>	<b>Remote cable</b> — Titrino/Titrando/ 692/712/713/780/781 
<b>6.2141.030</b>	<b>Remote cable</b> — 2xTitrino/Titrando/692/712/713/780/781 
<b>6.2141.040</b>	<b>Remote cable</b> — Titrino/Titrando — Dosimat 765/776 
<b>6.2141.050</b>	<b>Remote cable</b> — Titrino/Titrando — Dosimat 2x 765/776 
<b>6.2141.060</b>	<b>Remote cable</b> — pH-Meter 691 
<b>6.2142.070</b>	<b>Remote cable</b> — pH-Meter 692/712/713/780/781 — Dosimat 765/776
<b>6.2125.110</b>	<b>RS232 Connection cable</b> for the serial connection Metrohm instrument — PC, 25 pin/f to 9 pin/f 
<b>6.2134.040</b>	<b>RS232 Connection cable</b> for the serial connection Metrohm instrument — PC, 9 pin/f to 9 pin/f 

### 6.3.5 Optional accessories and additional devices

Order no.	Description	
6.1458.010	Macro titration head, PTFE for Sample Processors With 6 openings SGJ 14 and 3 openings SGJ 9 Suitable for rinsing and aspiration equipment.	
6.1458.020	Micro titration head, PTFE for Sample Processors With 4xM10 thread connectors	
6.1458.030	Sample Processor titration head for Karl Fischer titrations, PTFE	
6.0340.000	Double Pt electrode for Sample Processor	
6.1543.210	3-fold stopper NS 14, with anti-diffusion micro tip	
6.1805.060	FEP tubing with 2 tubing nipples M6, inner diameter 2 mm, with light and anti-kink protection L = 60 cm	
6.1808.000	Coupling bush, ETFE 3x M6 threads	
6.1808.090	Thread adapter, PVDF Thread M8 outer / M6 inner	
6.1821.000	Aspiration tube, PTFE L = 178 mm	

6.2037.040	Foil holder for sample beaker 6.1432.210, 24 pcs. for Al foil 6.2820.000, see below	
6.2820.000	Al foil 1000 disks 80 mm dia. Thickness 0.010 mm	
E.301.0022	O-Rings 5.28/1.78 mm	
E.301.0080	O-Rings 28/5 mm	
<b>2.802.0020</b>	<b>802 Rod Stirrer</b>	
6.1909.020 6.1909.040	with PP propeller or PTFE propeller, optional	
<b>2.741.0010</b>	<b>741 Magnetic Stirrer</b>	
<b>2.772.0110</b>	<b>772 Pump Unit</b> Peristaltic pump, with standard accessories 16...24 V, M8 connector	
<b>2.772.0120</b>	<b>772 Pump Unit, with aspiration equipment</b>	
<b>2.772.0130</b>	<b>772 Pump Unit, with rinsing equipment</b>	
<b>2.823.0010</b>	<b>823 Membrane Pump Unit</b> Membrane pump, with standard accessories 16...24 V, M8 connector	
<b>2.823.0020</b>	<b>823 Membrane Pump Unit, with aspiration equipment</b>	
<b>2.823.0030</b>	<b>823 Membrane Pump Unit, with rinsing equipment</b>	

<b>6.2001.070</b>	<b>Support stand for Sample Processor</b>  Suitable for Sample Processor models 778/ 789/814/815  For mounting external titration cells and other accessories, including 801 Stirrer or 804 Ti Stand  
<b>6.6056.xx1</b>	<b>tiamo light</b> PC software for controlling a titration or dosing system. Up to two instruments can be connected.
<b>6.6056.xx2</b>	<b>tiamo full</b> PC software for controlling complex titration or dosing systems.
<b>6.6056.xx3</b>	<b>tiamo multi</b> Client/Server software for controlling complex titration or dosing systems.

**6.3.6 Electrodes for automation**

For titrations with the macro titration head we recommend the use of long-life electrodes (LL) or Titrades (without glass joint) with 6.1236.040 Sleeve made of silicone rubber.

<i>Order no.</i>	<i>Description</i>
<b>6.0229.100</b>	<b>Solvotrode</b> 12.5 cm
<b>6.0232.100</b>	<b>Ecotrode</b> 12.5 cm
<b>6.0253.100</b>	<b>Aquatrode</b> 12.5 cm
<b>6.0258.000</b>	<b>Unitrode</b> 12.5 cm
<b>6.0431.100</b>	<b>Ag-Titrode</b> 12.5 cm
<b>6.0430.100</b>	<b>Pt-Titrode</b> 12.5 cm

If the micro titration head or the 786 Swing Head is used for direct titration then the following special micro-electrodes can be used.

<i>Order no.</i>	<i>Description</i>
<b>6.0234.110</b>	<b>Comb. Micro pH electrode (LL)</b> 16 cm
<b>6.0736.110</b>	<b>Micro reference electrode Ag/AgCl</b> 16 cm
<b>6.0134.110</b>	<b>Micro glass electrode</b> 16 cm
<b>6.0433.110</b>	<b>Micro Ag-Titrode</b> 16 cm
<b>6.0434.110</b>	<b>Micro Pt-Titrode</b> 16 cm
<b>6.0435.110</b>	<b>Micro Au-Titrode</b> 16 cm
<b>6.1110.110</b>	<b>Pt 1000 Temperature sensor</b> 16 cm

## **6.4 Warranty and conformity**

### **6.4.1 Warranty (guarantee)**

Metrohm guarantees that the deliveries and services it provides are free from material, design or manufacturing errors. The warranty period is 36 months from the day of delivery; for day and night operation it is 18 months. The warranty remains valid on condition that the service is provided by an authorized Metrohm service organization.





Glass breakage is excluded from the warranty for electrodes and other glassware. The warranty for the accuracy corresponds to the technical specifications given in the 'Instructions for Use' of the instrument. For components from third parties that make up a considerable part of our instrument, the manufacturer's warranty provisions apply. Warranty claims cannot be pursued if the Customer has not complied with the obligations to make payment on time.

During the warranty period Metrohm undertakes, at its own choice, to either repair at its own premises, free of charge, any instruments that can be shown to be faulty or to replace them. Transport costs are to the Customer's account.

Faults arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc. are expressly excluded from the warranty.





**6.4.2 Declaration of Conformity (778 Sample Processor)**

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

<p><i>Name of commodity</i></p> <p><b>778 Sample Processor</b></p>	 <p>CH-9101 Herisau, Switzerland E-Mail info@metrohm.com www.metrohm.com</p>
<p><i>Description</i> Sample changer with advanced Liquid Handling abilities for the automation of batch processing of larger sample series, applying titration, dosing and measuring methods in analytical laboratories.</p>	
<p>This instrument has been built and has undergone final type testing according to the standards:</p> <p><i>Electromagnetic compatibility: Emission</i> EN/IEC 61326, EN 55022 / CISPR 22, EN/IEC 61000-3-2</p> <p><i>Electromagnetic compatibility: Immunity</i> EN/IEC 61326, EN/IEC 61000-4-2, EN/IEC 61000-4-3, EN/IEC 61000-4-4, EN/IEC 61000-4-5, EN/IEC 61000-4-6, EN/IEC 61000-4-8, EN/IEC 61000-4-11, EN/IEC 61000-4-14, NAMUR</p> <p><i>Safety specifications</i> EN/IEC 61010-1, UL 3101-1 protection class I</p> <p>It has also been certified by ElectroSuisse, a member of the International Certification Body (CB/IEC).</p>	
<p> <i>The instrument meets the requirements of the CE mark as contained in the EU directives 89/336/EEC and 73/23/EEC and fulfils the following specifications:</i></p> <p>EN 61326-1 Electrical equipment for measurement, control and laboratory use – EMC requirements EN 61010-1 Safety requirements for electrical equipment for measurement, control and laboratory use</p>	
<p>Metrohm Ltd. is holder of the SQS-certificate ISO 9001:2000 Quality management system for development, production and sales of instruments and accessories for ion analysis.</p>	
<p>The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance.</p> <p>The technical specifications are documented in the instruction manual.</p>	
<p>Herisau, March 31, 2003</p> <div style="display: flex; justify-content: space-around; align-items: flex-end;"> <div style="text-align: center;">  <p>Dr. J. Frank Vice President Head of R&amp;D</p> </div> <div style="text-align: center;">  <p>Ch. Buchmann Vice President Head of Production Responsible for Quality Assurance</p> </div> </div>	

### 6.4.3 Declaration of Conformity (789 Robotic Sample Processor XL)

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

<p>Name of commodity</p> <p><b>789 Robotic Sample Processor XL</b></p>	 <p>CH-9101 Herisau, Switzerland E-Mail info@metrohm.com www.metrohm.com</p>
<p><i>Description</i> Sample changer with advanced Liquid Handling abilities for the automation of batch processing of larger sample series, applying titration, dosing and measuring methods in analytical laboratories.</p>	
<p>This instrument has been built and has undergone final type testing according to the standards:</p> <p><i>Electromagnetic compatibility: Emission</i> EN/IEC 61326, EN 55022 / CISPR 22, EN/IEC 61000-3-2</p> <p><i>Electromagnetic compatibility: Immunity</i> EN/IEC 61326, EN/IEC 61000-4-2, EN/IEC 61000-4-3, EN/IEC 61000-4-4, EN/IEC 61000-4-5, EN/IEC 61000-4-6, EN/IEC 61000-4-8, EN/IEC 61000-4-11, EN/IEC 61000-4-14, NAMUR</p> <p><i>Safety specifications</i> EN/IEC 61010-1, UL 3101-1 protection class I</p> <p>It has also been certified by ElectroSuisse, a member of the International Certification Body (CB/IEC).</p>	
<p> <i>The instrument meets the requirements of the CE mark as contained in the EU directives 89/336/EEC and 73/23/EEC and fulfils the following specifications:</i></p> <p>EN 61326-1 Electrical equipment for measurement, control and laboratory use – EMC requirements EN 61010-1 Safety requirements for electrical equipment for measurement, control and laboratory use</p>	
<p>Metrohm Ltd. is holder of the SQS-certificate ISO 9001:2000 Quality management system for development, production and sales of instruments and accessories for ion analysis.</p>	
<p>The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance.</p> <p>The technical specifications are documented in the instruction manual.</p>	
<p>Herisau, March 31, 2003</p> <div style="display: flex; justify-content: space-around; align-items: flex-end;"> <div style="text-align: center;">  <p>Dr. J. Frank Vice President Head of R&amp;D</p> </div> <div style="text-align: center;">  <p>Ch. Buchmann Vice President Head of Production Responsible for Quality Assurance</p> </div> </div>	

#### 6.4.4 Quality Management Principles

Metrohm Ltd., CH-9101 Herisau, Switzerland

 **Metrohm**  
ion analysis  
CH-9101 Herisau/Switzerland  
E-Mail info@metrohm.com  
Internet www.metrohm.com

Metrohm Ltd. holds the ISO 9001 Certificate, registration number 10872-02, issued by SQS (Swiss Association for Quality and Management Systems). Internal and external audits are carried out periodically to assure that the standards defined by Metrohm's QM Manual are maintained.

The steps involved in the design, manufacture and servicing of instruments are fully documented and the resulting reports are archived for ten years. The development of software for PCs and instruments is also duly documented and the documents and source codes are archived. Both remain the possession of Metrohm. A non-disclosure agreement may be asked to be provided by those requiring access to them.

The implementation of the ISO 9001:2000 quality system is described in Metrohm's QM Manual, which comprises detailed instructions on the following fields of activity:

##### **Instrument development**

The organization of the instrument design, its planning and the intermediate controls are fully documented and traceable. Laboratory testing accompanies all phases of instrument development.

##### **Software development**

Software development occurs in terms of the software life cycle. Tests are performed to detect programming errors and to assess the program's functionality in a laboratory environment.

##### **Components**

All components used in the Metrohm instruments have to satisfy the quality standards that are defined and implemented for our products. Suppliers of components are audited by Metrohm as the need arises.

##### **Manufacture**

The measures put into practice in the production of our instruments guarantee a constant quality standard. Production planning and manufacturing procedures, maintenance of production means and testing of components, intermediate and finished products are prescribed.

##### **Customer support and service**

Customer support involves all phases of instrument acquisition and use by the customer, i.e. consulting to define the adequate equipment for the analytical problem at hand, delivery of the equipment, user manuals, training, after-sales service and processing of customer complaints. The Metrohm service organization is equipped to support customers in implementing standards such as GLP, GMP, ISO 900X, in performing Operational Qualification and Performance Verification of the system components or in carrying out the System Validation for the quantitative determination of a substance in a given matrix.

# 7 Index

- < ..... 55  
 << **NEXT** > ..... 15  
 <→ **PREV** > ..... 15  
 <↑>, <↓> ..... 15, 54  
 <\*> ..... 79  
 <**CLEAR/RESET**> . 16, 52, 89  
 <**CONFIG**> . 15, 114, 115, 116  
 <**CTRL**> ..... 19  
 <**DEF**> ..... 19, 84  
 <**DELETE**> ..... 16, 67  
 <**DOS**> ..... 18  
 <**END**> ..... 15, 52  
 <**ENTER**> ..... 16, 52, 54  
 <**HOLD**> ..... 22, 89  
 <**HOME**> ..... 15, 52  
 <**INSERT**> ..... 16, 67  
 <**LEARN**> ..... 59, 112, 113  
 <**LIFT**> ..... 17  
 <**MOVE**> ..... 16  
 <**PARAM**> ..... 15, 115  
 <**PRINT**> ..... 21  
 <**PUMP**> ..... 17  
 <**QUIT**> ..... 16, 52, 54, 89  
 <**RACK**> ..... 9, 21, 91  
 <**SAMPLE**> ..... 16, 74  
 <**SCAN**> ..... 18  
 <**SELECT/TOWER**> . 15, 52,  
 ..... 54, 87  
 <**START**> ..... 22, 89, 114, 129  
 <**STIR**> ..... 17  
 <**STOP**> ..... 22, 72, 89  
 <**USER METHOD**> ..... 15  
 <**WAIT**> ..... 19  
 > ..... 55  
  
 772 Pump Unit ..... 29, 41  
 778 Sample Processor ..... 138  
 786 Swing Head 27, 59, 75, 116  
 789 Robotic Sample Processor  
 XL ..... 138  
 823 Membrane Pump Unit ... 29
- A** \_\_\_\_\_  
 Absolute positioning ..... 76  
 Absorber tube ..... 96  
 Acetic anhydride ..... 41  
 Acetone ..... 41  
 Acknowledgement ..... 71  
 Active ..... 80, 82, 109  
 Adjustment ..... 113  
 Adjustment rate ..... 118
- Air bubble ..... 97, 105  
 Air inlet ..... 71  
 Air inlet opening ..... 96  
 Ambient temperature ..... 128  
 Angle ..... 116  
 Angle of rotation ..... 75  
 Angle of rotation offset ..... 117  
 Angle positions ..... 75  
 Annex ..... 125  
 Arrow keys ..... 15, 54  
 ASCII character ..... 81  
 Aspirate ..... 79  
 Aspirating air ..... 98  
 Aspiration ..... 102  
 Aspiration equipment ..... 40  
 Aspiration of the sample ..... 107  
 Aspiration pump ..... 41  
 Aspiration tip ..... 41, 42, 43  
 Assembling pipetting  
 accessories ..... 44  
 Attaching the tubing ..... 41  
 Attention ..... 6  
 AutoInfo ..... 81, 83  
 Automatic emptying ..... 103  
 Automatic rack recognition ... 92  
 Axis distance ..... 116  
 Axis of rotation ..... 116
- B** \_\_\_\_\_  
 Basic settings ..... 57, 59  
 Battery ..... 120  
 Baud ..... 65  
 Beaker ..... 90, 151, 152  
 Beaker error ..... 69  
 Beaker radius ..... 59, 62, 93  
 Beaker sensor ..... 10, 76, 93  
 Beaker test ..... 10  
 Beep ..... 57  
 Binary code ..... 91  
 Binary pattern ..... 80, 109  
 Biohazard ..... 24  
 Buret tip ..... 42, 43, 44
- C** \_\_\_\_\_  
 Cable 6.2134.050 ..... 39  
 Cable 6.2134.110 ..... 39  
 Canister ..... 41  
 Capacity ..... 125  
 Care ..... 119  
 Carrier liquid ..... 105  
 Changer settings ..... 68  
  
 Changing dosing unit ... 98, 103  
 Character selection ..... 55  
 Character set ..... 65  
 Chemical resistance ..... 41  
 Closing Dosino port ..... 107  
 Cock direction ..... 71  
 Code ..... 61  
 Collection trough ..... 47  
 Colon ..... 54  
 Command keys ..... 16  
 Command line ..... 67  
 Command sequence ..... 67  
 Commands ..... 67, 74  
 Communication commands . 80  
 Compensation ..... 104  
 Conductometer ..... 32  
 Configuration ..... 15, 57  
 Configuration menu 49, 57, 115  
 Configuration report ..... 21  
 Configure robotic arm ..... 28  
 Connecting a printer ..... 38  
 Connecting accessories ..... 40  
 Connecting dosing devices .. 30  
 Connecting pumps ..... 29  
 Connecting rod ..... 104  
 Connecting the keypad ..... 27  
 Connection cable ..... 31, 158  
 Connections sockets ..... 125  
 Continue method ..... 89  
 Contrast of display ..... 57  
 Control command ..... 109  
 Counterclockwise ..... 85  
 Coupling ..... 104  
 Creating methods ..... 67  
 CTL command ..... 109, 110  
 CTL:Rm command ..... 82  
 Current sample position ..... 16  
 Custom ..... 39
- D** \_\_\_\_\_  
 Danger ..... 6  
 Data bits ..... 65  
 Data input ..... 16, 54  
 Data reception ..... 65  
 Data transmission ..... 31  
 Data transmission via serial  
 interface ..... 83  
 DEF commands ..... 84  
 Default port ..... 79, 102  
 Default value ..... 16  
 Delete ..... 16, 55, 63

- Delete character..... 55  
Delete line ..... 67  
Delete rack definitions ..... 63  
Deleting methods ..... 87  
Depth ..... 128  
Diagnosis ..... 123  
Dialog arrangement ..... 53  
Dialog language ..... 57  
Diameter ..... 3, 64  
Diluting ..... 105  
Dimensions ..... 128  
Dimethylformamide..... 41  
Direction of rotation ..... 19, 68,  
..... 71, 76, 85  
Disable configuration..... 115  
Disable display ..... 115  
Disable keypad ..... 114  
Disable method storage ..... 115  
Disable parameter ..... 115  
Disabling keypad functions 114  
Display ..... 13, 14, 51  
Disposable beakers..... 90, 151  
Distributor block..... 8, 40  
DOS command..... 79, 101  
Dosimat... 30, 35, 64, 70, 79, 95  
Dosing..... 95, 97, 102  
Dosing cylinder ..... 101  
Dosing cylinder ejection ..... 99  
Dosing device connection 64, 70  
Dosing device control..... 78  
Dosing device settings ..... 70  
Dosing devices ..... 78  
Dosing drive..... 30, 64  
Dosing functions..... 79  
Dosing instruments..... 18, 95  
Dosing outlet..... 96  
Dosing port ..... 70, 79  
Dosing rate ..... 19, 70, 85  
Dosing unit..... 64, 96  
Dosino..... 30, 70, 79, 95, 100  
Dosino commands ..... 79  
Drawing up a method ..... 66  
Duration ..... 78
- E** \_\_\_\_\_  
Eject ..... 79, 104  
Eject cylinder contents.. 99, 106  
Ejecting ..... 104, 108  
Electrical safety ..... 23  
Electrode..... 42  
Electromagnetic compatibility  
..... 127  
Elimination of mechanical play  
..... 99  
EMC ..... 127  
Emission ..... 127  
Empty..... 98, 103
- Emptying the dosing unit..... 71  
EndVol..... 79, 104  
End volume..... 99, 104  
Entry ..... 16  
EOD ..... 110  
Epson..... 65  
EPSON Emulation..... 38  
Error message ..... 22, 71  
Exchange units ..... 64  
Ext.1 ..... 75  
External position ..... 59, 75, 112  
External pump..... 29, 77  
External Start..... 57, 89, 110
- F** \_\_\_\_\_  
Fault ..... 120  
Field of application ..... 1  
Fill ..... 97, 102  
Fill air bubble-free ..... 103  
Fill automatically ..... 79  
Fill tubing ..... 98, 103  
Filling inlet ..... 96  
Filling port ..... 70  
Filling rate ..... 19, 70, 85  
Final sequence ..... 67  
First sample ..... 16  
Frequency ..... 127  
Function selection..... 78  
Functions ..... 79  
Fuses ..... 127
- G** \_\_\_\_\_  
Gentle stop ..... 89  
GLP ..... 122  
Good Laboratory Practice .. 122  
Grounding ..... 23  
Guide sleeve ..... 42, 43
- H** \_\_\_\_\_  
Handle ..... 48  
Handshake..... 65  
Height ..... 128  
Hold method sequence ..... 22  
HOLD status ..... 22, 89, 120
- I** \_\_\_\_\_  
IBM..... 65  
IBM Emulation..... 38  
Immunity ..... 127  
Impulse ..... 82, 111  
Inactive..... 80, 82  
Infrared sensor..... 10  
Initial value ..... 52  
Initialize ..... 16  
Initialize rack ..... 21, 84  
Initializing the working memory  
..... 123
- Injection needle..... 45, 107  
Input ..... 67, 80, 109, 125  
Input line ..... 54, 57, 109  
Insert ..... 16  
Installation ..... 25  
Installing the titration  
accessories ..... 42  
Instrument dialog ..... 52  
Instrument status ..... 14  
Instrument versions ..... 3  
Interrupt method ..... 89  
Interrupt sample series ..... 89  
Interruption..... 22  
Ion meter ..... 32, 33  
ISO 9100 ..... 164, 165
- K** \_\_\_\_\_  
Keypad..... 13, 27, 125  
Keys ..... 13  
KF\_ext..... 133
- L** \_\_\_\_\_  
LEARN function..... 59, 67,  
..... 77, 80, 113  
LEARN mode ..... 78, 112  
LEDs ..... 13  
Length..... 64  
Length of robotic arm ..... 117  
Lift ..... 17, 77  
Lift operation ..... 15  
Lift position ..... 112  
Lift positioning..... 77  
Lift rate ..... 20, 68  
Lift speed ..... 85  
Limiting screw ..... 28  
Line condition ..... 80, 82, 109  
Lines ..... 66  
Liquid handling ..... 79, 95  
Liquid Handling functions ... 100  
Lower ..... 17, 77
- M** \_\_\_\_\_  
Macro titration head..... 42  
Magnet code. 47, 61, 90, 91, 92  
Magnet sensor..... 9  
Magnetic stirrer ..... 46, 127,  
..... 130, 135, 136, 160  
Mains connection ..... 23, 127  
Maintenance ..... 58, 119  
Manual stop ..... 72  
Manual stop options ..... 72, 111  
Materials ..... 128  
Max. load ..... 125, 127  
Max. swing angle..... 116  
Max. waiting time ..... 71  
Max. volume..... 99  
Measuring out..... 105

- Mechanical play..... 79, 104  
 Membrane pump..... 29, 41  
 Method name ..... 14  
 Method parameters ..... 15, 66  
 Method sequence..... 22  
 Methods..... 15, 66, 87  
 Metrohm Serial Bus ..... 9  
 Micro electrode..... 42, 43, 44  
 Micro titration head..... 42  
 Min. beaker radius ..... 59  
 Missing beaker ..... 68  
 Missing special beaker..... 69  
 Models ..... 3, 4  
 Mounting a 786 Swing Head 27  
 Mounting robotic arm ..... 28  
 MOVE..... 69  
 MOVE command ..... 75  
 Move to rack positions ..... 75  
 MSB ..... 46, 78, 79  
 MSB socket ..... 46  
 MSB stirrer ..... 69  
 Multiple cable ..... 110
- N** \_\_\_\_\_  
 Needle ..... 45  
 Negative volume..... 79, 102  
 New command line..... 67  
 Next..... 76  
 Non-aqueous samples 105, 107  
 NOP command..... 67  
 Normal condition ..... 14  
 Notation ..... 6  
 Number of samples ..... 66
- O** \_\_\_\_\_  
 Operating condition..... 51  
 Operating hours counter ..... 58  
 Operation..... 51  
 Organic Solvents ..... 29  
 Output..... 80, 109, 125  
 Output lines ..... 109  
 Overall view..... 7  
 Overview ..... 1
- P** \_\_\_\_\_  
 Parameter menu ..... 66, 115  
 Parameter report..... 21, 38  
 Parametrizing ..... 113  
 Parity ..... 65  
 Parking position..... 15, 84  
 Pause..... 84  
 Peripheral devices ..... 31, 72  
 Peristaltic pump..... 29, 160  
 pH meter..... 32, 33, 35, 130,  
 ..... 131, 133, 135, 136  
 pH\_cal ..... 135  
 Pictograms ..... 6
- Piezo sensor ..... 10  
 Pipetting..... 79, 102, 104  
 Pipetting equipment ..... 105  
 Pipetting procedure .... 105, 106  
 Pipetting tubing ..... 105  
 Port ..... 64, 96  
 Port assignments..... 79  
 Port selection..... 78  
 Position beaker ..... 75  
 Position table ..... 91  
 Positioning accuracy ..... 118  
 Power consumption..... 127  
 POWERUP method..... 88  
 Predefined values ..... 52  
 PREP function..... 64  
 Preparation ..... 64, 106  
 Preparation cycle ..... 71, 97  
 Prepare ..... 97  
 Preparing the dosing unit ... 103  
 Prevent loading..... 115  
 Prevent method deletion .... 115  
 Prevent storage ..... 115  
 Previous ..... 76  
 Printer..... 37, 65  
 Printer cable..... 37  
 Programming run sequences 74  
 Propeller stirrer .... 42, 43, 127,  
 ..... 130, 135, 136, 160  
 Protected port..... 71, 85  
 PTFE ..... 41  
 PTFE tubing ..... 41  
 Pulsed signal ..... 80  
 Pump ..... 17, 72  
 PUMP command ..... 77  
 Pump connections 17, 29, 72, 77  
 Pump control ..... 77  
 Pump outlet ..... 127  
 Pump Unit..... 160  
 PVDF..... 41
- Q** \_\_\_\_\_  
 Quit command..... 89
- R** \_\_\_\_\_  
 Rack..... 61  
 Rack adjustment..... 113  
 RACK command..... 84, 118  
 Rack data..... 61, 92  
 Rack definitions ..... 61  
 Rack name..... 61, 63, 68, 92  
 Rack offset..... 113  
 Rack position ..... 16, 63, 74, 76  
 Rack recognition..... 21, 91  
 Rack rotation..... 16  
 Raise..... 17, 77  
 RAM init..... 123  
 Ready..... 71, 110
- Rear view ..... 8  
 Redefinition..... 84  
 Relative humidity..... 128  
 Relative positioning ..... 76  
 Remarks..... 6  
 Remote Box 8.2148.010 ..... 31  
 Remote cable ..... 31, 158  
 Remote connections..... 31  
 Remote control ..... 65  
 Remote interface ... 18, 57, 72,  
 ..... 80, 82, 109  
 Remote lines..... 19  
 Removing dosing unit..... 103  
 Report ..... 38  
 Reserved position..... 94  
 Resistance ..... 41  
 Rinse..... 108  
 Rinse port ..... 70  
 Rinse position ..... 77, 93  
 Rinsing..... 40, 105  
 Rinsing equipment..... 8, 40  
 Rinsing height..... 59  
 Rinsing process..... 41  
 Rinsing volumes ..... 64  
 Robotic arm ..... 43, 75, 116  
 Robotic arm offset ..... 116  
 Robotic arm position ..... 112  
 Robotic arm sensor ..... 10  
 Rotate position..... 77, 93  
 Rotate rack ..... 93  
 Rotate robotic arm..... 75  
 Rotating height ..... 59  
 Rotating position..... 76  
 Rotating the sample rack..... 16  
 Rotation increment ..... 68  
 Rotation rate ..... 68  
 RS232 ..... 37, 65  
 RS232 interface ..... 65, 72,  
 ..... 80, 81, 83  
 RS232 settings ..... 65  
 RS232/Parallel Converter..... 39  
 Run control ..... 89  
 Run sequence ..... 15, 66  
 Running time..... 84  
 Run variable..... 74
- S** \_\_\_\_\_  
 Safety information..... 7, 23  
 Safety instructions ..... 127  
 Safety specification..... 127  
 Sample beaker 66, 69, 150, 152  
 Sample changer settings.... 116  
 SAMPLE command ..... 74  
 Sample counter ..... 14  
 Sample position..... 74  
 Sample positioning..... 15

- Sample rack..... 3, 47, 61,  
                                   68, 90, 150, 152  
 Sample series ..... 66  
 SAMPLE variable ..... 75  
 Save rack definitions..... 63  
 SCAN command..... 109  
 Scan remote interface..... 80  
 SCAN timeout ..... 71  
 Scanning RS232 interface .... 81  
 SCN command ..... 109, 110  
 SCN:Rm command ..... 80  
 Sealed vials..... 107  
 Seconds counter..... 84  
 Select method..... 87  
 Select sample position ..... 74  
 Selecting a stirrer ..... 78  
 Selection drum..... 54  
 Selection list..... 15  
 Selection of the dosing device70  
 Sensor..... 10  
 Separating bubble105, 106, 107  
 Sequence..... 67  
 Serial interface ..... 37, 65  
 Service ..... 119  
 Setting the remote lines ..... 82  
 Settings ..... 19, 49, 64  
 Setup..... 26  
 Setup menu ..... 116  
 Setup settings ..... 28  
 SGJ sleeve ..... 42  
 Signal lines..... 72  
 Signals ..... 18, 109  
 Solenoid valve..... 8  
 Solids ..... 41  
 Solvents ..... 29  
 SOP..... 122  
 Spec.1 ..... 75  
 Special beaker .... 16, 62, 74, 94  
 Special beaker positions ..... 75  
 Special height ..... 62  
 Special position ..... 62, 77, 93  
 Speed of rotation ..... 20, 85  
 Splash protection..... 7  
 Spray nozzles ..... 40, 42, 43  
 Stand support ..... 46  
 Standard dosing port..... 70  
 Standard equipment ..... 138  
 Standard methods ..... 129  
 Standard sample racks... 61, 90  
 Standards ..... 127  
 Standpipe ..... 96  
 Start..... 67  
 Start method ..... 22, 89  
 Start sequence..... 66  
 Start-up ..... 49  
 Status line ..... 14, 51  
 Status message ..... 83  
 Std\_add ..... 136  
 Step by step..... 129  
 Step by step procedure ..... 114  
 STIR command ..... 78  
 Stirrer..... 17, 72  
 Stirrer connections..... 8, 78  
 Stirrer control ..... 78  
 Stirrer guide ..... 46  
 Stirrer models 801, 804..... 69  
 Stirrer models 802, 741, 722. 69  
 Stirrer settings ..... 69  
 Stirrer shaft..... 42, 43  
 Stirring rate..... 20, 69, 84, 127  
 Stop bits..... 65  
 Stop method ..... 89  
 Stopcock..... 71, 85, 101  
 Stopcock rotation..... 104  
 Stopping a sample series..... 22  
 Submenu ..... 52  
 Swing angle ..... 60, 75  
 Swing axis ..... 116  
 Swing direction ..... 117  
 Swing increment ..... 68  
 Swing radius ..... 117  
 Swing range ..... 75, 116  
 Swing rate ..... 20, 68, 86  
 Swing to the left ..... 28, 117  
 Swing to the right..... 28, 117  
 Switching components ..... 77  
 Switching together ..... 78  
 System components..... 5
- T**
- Technical data ..... 125  
 Test command sequence... 114  
 Test tube ..... 90  
 Testing ..... 114  
 Text cursor ..... 55  
 Text editor ..... 55  
 Text input ..... 55  
 Text input mode ..... 87  
 Threaded stopper ..... 42, 43  
 Timeout ..... 71  
 Tinet ..... 65  
 Title line ..... 14, 51  
 Titrando..... 32, 35, 130,  
                                   131, 133, 135, 136  
 Titration accessories..... 42  
 Titration beaker ..... 90  
 Titration cell..... 106, 107  
 Titration head ..... 27  
 Titration stand ..... 46  
 Titrino ..... 32, 35, 130, 131, 133  
 Tolerance ..... 104  
 Tower ..... 59  
 TOWER 1 ..... 15  
 Tower settings..... 59  
 TRACE function..... 67, 114, 129  
 Transfer robotic arm ... 156, 157  
 Transmission parameters .. 37,  
                                   81, 83  
 Transmission rate ..... 65  
 Transport damage ..... 163  
 Triggers..... 83  
 Tubing..... 41, 64  
 Tubing dimensions ..... 64  
 Tubing parameters ..... 64  
 Tubing system ..... 97  
 Tubing tip ..... 105  
 Turntable ..... 9
- U**
- User dialog ..... 114
- V**
- Validation ..... 122  
 VENT ..... 96  
 Venting ..... 96  
 Versions ..... 3  
 Vertical movement rate ..... 68  
 Voltage ..... 127  
 Volume ..... 79, 97
- W**
- WAIT command ..... 19, 84  
 Waiting time ..... 71, 84  
 Warning..... 6, 57  
 Warning limit ..... 58  
 Warranty..... 163  
 Waste canister ..... 41  
 Weight..... 128  
 Width..... 128  
 Wild card..... 81, 97  
 Work pos. T1..... 112  
 Work position ..... 77, 92  
 Working height..... 60, 61  
 Working position ..... 15, 61  
 Working range ..... 128
- Z**
- Zero position ..... 77