


 **Metrohm**
Ion analysis
CH-9101 Herisau/Switzerland
Phone ++41 71 353 85 85
Fax ++41 71 353 89 01
Internet www.metrohm.ch
E-Mail info@metrohm.ch

790 Personal IC

Program «790 PC Software 1.0»



8.790.1003 Instructions for Use

14.08.2000 / dö

Table of contents

1	Introduction	1
1.1	Instrument description	1
1.2	Parts and controls	3
1.3	Information on the Instructions for Use	10
1.3.1	Organization	10
1.3.2	Notation and pictograms	11
1.4	Safety notes	12
1.4.1	Electrical safety	12
1.4.2	General precautionary rules	12
2	Installation	13
2.1	Overview	13
2.1.1	Flow chart	13
2.1.2	Connections in the 790 Personal IC	13
2.2	Setting up the instrument	15
2.2.1	Packaging	15
2.2.2	Check	15
2.2.3	Location	15
2.3	Attaching the accessories	15
2.3.1	Connection of detector block	15
2.3.2	Connection of syringe and aspirating tubing	16
2.3.3	Connection of the 6.5324.000 Bottle rack (option)	16
2.3.4	Connection of PEEK capillaries	17
2.3.5	Filter unit PEEK	18
2.4	Mains connection	19
2.4.1	Setting the mains voltage	19
2.4.2	Fuses	20
2.4.3	Mains cable and mains connection	20
2.4.4	On/off switching of the instrument	20
2.5	Connection to the PC	21
2.5.1	Connecting cable	21
2.5.2	Software installation	21
2.5.3	Basic settings	22
2.6	High-pressure pump	25
2.6.1	Removing the transport security screws	25
2.6.2	Installing the pulsation dampener (option)	25
2.6.3	Connecting the eluent bottle	27
2.6.4	Deaerating the pump and rinsing the pulsation dampener	27
2.7	Precolumns and separating columns	30
2.7.1	General information on precolumns	30
2.7.2	Precolumns with cartridge head	30
2.7.3	Precolumn glass cartridges with cartridge holder	32
2.7.4	IC anion precolumn SUPERSEP	33
2.7.5	General information on separating columns	34
2.7.6	Selection of the sample loop	34
2.7.7	Installation of the separating column without suppressor	35
2.7.8	Installation of the separating column with suppressor	36

2.8	Suppressor module	37
2.8.1	General information on suppressor module.....	37
2.8.2	Preparation of the peristaltic pump.....	37
2.8.3	Connection of supply bottles.....	40
2.8.4	Connection of the suppressor module.....	41
2.9	Putting into operation	43
2.9.1	Putting into operation without suppressor (cation system).....	43
2.9.2	Putting into operation with suppressor (anion system).....	44
2.10	Connection of external devices	46
2.10.1	Connection of the 750 Autosampler.....	46
2.10.2	Connection of the 766 IC Sample Processor.....	47
3	Operating tutorial	51
3.1	Requirements.....	51
3.2	Preparations.....	52
3.3	Calibration.....	53
3.4	Sample determination.....	60
4	Operation	63
4.1	Fundamentals of the operation	63
4.1.1	Starting/closing the program.....	63
4.1.2	Glossary.....	64
4.1.3	Overview of program windows.....	65
4.1.4	Main window elements.....	66
4.1.5	Icons of the main window.....	66
4.1.6	Overview of file types.....	67
4.1.7	Context sensitive menus.....	68
4.1.8	Keyboard and mouse functions.....	68
4.1.9	Help.....	69
4.2	Instrument and software settings	70
4.2.1	Fonts.....	70
4.2.2	Security system.....	70
4.2.3	Global settings.....	71
4.2.4	COM port.....	73
4.3	Systems	74
4.3.1	System window.....	74
4.3.2	System file handling.....	74
4.3.3	System functions.....	75
	Start/stop hardware and record baseline.....	75
	Start/stop determinations.....	75
	Print system parameters.....	76
4.3.4	PC icon.....	77
4.3.5	Watch window.....	77

4.3.6	Instrument icon.....	78
	Menu options for instrument icon	78
	System parameters for disconnected system	78
	Instrument control for connected system	79
	Time program.....	80
	Configuration.....	82
	Hardware settings	82
4.3.7	System state window	84
	Status messages	84
	Error messages.....	85
4.4	Methods.....	86
4.4.1	Method file handling.....	86
4.4.2	Passport	86
	General.....	87
	Sample	88
	Column.....	89
	Eluent	90
	Comment.....	90
	Method Log	91
	Data Log.....	92
4.4.3	Method setup	93
	General.....	93
	Processing	93
	Math	94
4.4.4	Integration	96
	Setup.....	97
	Events.....	98
4.4.5	Calibration and quantification	100
	General information.....	100
	Notations	101
	External standard calibration.....	102
	Component table	102
	Peak identification	104
	Concentration table.....	106
	Calibration curve	108
	Update calibration.....	110
	Load and save calibration data.....	111
	Import and export calibration data.....	111
	Put out calibration curves.....	111
4.4.6	Report output	112
	Report options.....	112
4.5	Chromatograms	124
4.5.1	Chromatogram window.....	124
4.5.2	Chromatogram file handling	125
	Open chromatogram.....	125
	Save chromatogram.....	126
	Close chromatogram	126
	Delete chromatogram	126
	Export chromatogram	127
4.5.3	Graphical representation.....	128
	Appearance.....	128
	Other graphical functions.....	132
4.5.4	Peak editor	133
	Switching on/off the peak editor	133
	Peak editor functions	133
	Moving the cursor	134

4.5.5	Printing	135
	Page layout for printing	135
	Printer settings	136
	Print preview	136
	Printing	136
4.5.6	Miscellaneous functions	136
	Reintegration	136
	Recalibration	136
	Subtraction of a chromatogram	137
	Data compression	137
	Invert chromatogram	137
4.6	Sample queue	138
4.6.1	Sample queue file handling	138
	Open sample queue	138
	Save sample queue	138
	Delete sample queue	138
4.6.2	Sample queue control	139
	Sample queue overview table	139
	Start sample queue	140
	Pause sample queue	141
	Cancel last run	141
	Reset sample queue	141
4.6.3	Sample queue editor	141
	Open queue editor window	141
	Sample queue editor functions	143
	Print sample queue	143
	Close sample queue editor	143
4.7	Batch reprocessing	144
4.7.1	Batch reprocessing queue file handling	144
	Open batch reprocessing queue	144
	Create new batch reprocessing queue	144
	Save batch reprocessing queue	144
4.7.2	Perform batch reprocessing	145
	Reprocess options window	145
	Merge chromatograms	148
4.7.3	Batch reprocessing queue editor	149
	Open batch reprocessing queue editor window	149
	Batch reprocessing queue editor functions	150
	Print batch reprocessing queue	150
	Close batch reprocessing queue editor	150

5	Notes – Maintenance – Faults	151
5.1	Practical notes on ion chromatography	151
5.1.1	Separating columns	151
5.1.2	High-pressure pump	152
5.1.3	Eluents.....	153
5.1.4	Peristaltic pump.....	153
5.1.5	Suppressor module.....	154
5.1.6	Connections	154
5.2	Maintenance and servicing	155
5.2.1	General information	155
5.2.2	Passivation	155
5.2.3	Recycling.....	156
5.2.4	Shutdown	156
5.2.5	Changing separating columns.....	156
5.2.6	Maintenance work at the pump head	158
5.2.7	Regeneration of the suppressor module	162
5.2.8	Cleaning the suppressor.....	163
5.2.9	Replacing the suppressor.....	165
5.2.10	Exchanging the pump tubing.....	167
5.3	Faults and malfunctions	168
5.3.1	Error messages	168
5.3.2	Malfunctions and their rectification	168
5.4	Diagnostic tests / Validation / GLP	170
6	Appendix	171
6.1	Technical data	171
6.1.1	Conductivity measurement	171
6.1.2	Conductivity detector	171
6.1.3	Injection valve.....	171
6.1.4	High-pressure pump.....	172
6.1.5	Peristaltic pump.....	173
6.1.6	Suppressor module.....	173
6.1.7	RS232 interface	173
6.1.8	Remote interface	174
6.1.9	Mains connection	174
6.1.10	Safety specifications	174
6.1.11	Electromagnetic compatibility (EMC).....	174
6.1.12	Ambient temperature.....	174
6.1.13	Housing	175
6.2	Standard equipment	176
6.3	Optional accessories	179
6.3.1	Accessories for 790 Personal IC.....	179
6.3.2	Separating columns and precolumns.....	181
6.3.3	Additional devices and cables	185
6.4	Warranty and conformity	186
6.4.1	Warranty	186
6.4.2	EU Declaration of conformity	187
6.4.3	Certificate of conformity and system validation	188
6.5	Index	189

List of figures

<u>Fig. 1:</u>	Front of the 790 Personal IC.....	3
<u>Fig. 2:</u>	Rear of the 790 Personal IC	4
<u>Fig. 3:</u>	Interior of the 2.790.0010 Personal IC.....	6
<u>Fig. 4:</u>	Interior of the 2.790.0020 Personal IC.....	8
<u>Fig. 5:</u>	Connecting diagram for 2.790.0010 Personal IC without suppressor.....	14
<u>Fig. 6:</u>	Connecting diagram for 2.790.0020 Personal IC with suppressor.....	14
<u>Fig. 7:</u>	Connectors for capillaries	17
<u>Fig. 8:</u>	6.2821.100 Filter unit PEEK.....	18
<u>Fig. 9:</u>	Setting the mains voltage.....	20
<u>Fig. 10:</u>	Connection of the pulsation dampener (option).....	26
<u>Fig. 11:</u>	Installation of precolumn cartridges with cartridge head.....	31
<u>Fig. 12:</u>	Installation of precolumn glass cartridges with cartridge holder	33
<u>Fig. 13:</u>	Installation of column without suppressor	36
<u>Fig. 14:</u>	Installing pump tubings.....	38
<u>Fig. 15:</u>	Installation of column with suppressor.....	39
<u>Fig. 16:</u>	Connections at suppressor module.....	42
<u>Fig. 17:</u>	Components of the pump head.....	159
<u>Fig. 18:</u>	Replacement of the piston seal 92	159
<u>Fig. 19:</u>	Component of inlet valve 93 and outlet valve 94	161
<u>Fig. 20:</u>	Assembling the suppressor	165

List of numbered parts and controls

1	Door to interior	3,6,8	39	Suppressor module	8,39
2	Mains pilot lamp	3	40	Tubing cartridge.....	8,38,39
3	Feedthrough for syringe tubing.....	3	41	Contact pressure lever.....	8,38,39
4	Feedthrough for aspirating tubing.....	3,6,8	42	Holding clamp	8,39
5	Opening for detector cable	4	43	Snap-action lever	8,38,39
6	Rear panel opening	4	44	Pump drive.....	8,39
7	Opening for outlet capillaries	4	45	Mounting pin.....	8,39
8	Knurled screw.....	4	46	Compression fitting.....	17,18,31,33,38
9	Detachable rear panel	4	47	Compression fitting.....	17
10	Transport security screws.....	4	48	Capillary	17,18,33
11	Mains switch.....	4,20	49	Connector with filter	18
12	Mains connection plug.....	4,20	50	Housing for filter unit PEEK.....	18
13	Fuse holder.....	4,20	51	Connector without filter	18
14	RS232 interface	4	52	Pulsation dampener.....	26
15	Connection for detector block.....	4	53	Connection to injection valve.....	26
16	Remote interface	4	54	Connection to purge valve.....	26
17	Serial number	4	55	Manufit housing	31
18	Inlet capillary for injector.....	6,8,26	56	Steel connector.....	31
19	Mounting rail.....	6,8,36,39	57	PTFE gasket.....	31
20	Column connection capillary.....	6,8,31,33,36,39	58	2 Steel meshes	31
21	Sample loop	6,8	59	Precolumn cartridge	31
22	Injection valve	6,8,26,36,39	60	Steel spacer	31
23	Aspirating tubing	6,8	61	4 Steel meshes	31
24	Coupling	6,8	62	Manufit pressure screw	31
25	Syringe tubing.....	6,8	63	Separating column.....	31,36,39
26	PEEK coupling.....	6,8,39	64	End fitting.....	33
27	Connection capillary.....	6,8,26	65	Screw cap for precolumn.....	33
28	Filter unit PEEK.....	6,8,18,26,38,39	66	Screw cap for precolumn.....	33
29	Connection capillary	6,8,26	67	Precolumn cartridge	33
30	Connection capillary	6,8	68	Connection piece.....	33
31	Purge valve	6,8,26	69	Screw cap for column.....	33
32	Aspirating capillary	6,8	70	Column holder	36,39
33	Fastening screws	6,8	71	Aspirating tubing for H ₂ O.....	38,39
34	Pump head.....	6,8,159	72	Aspirating tubing for H ₂ SO ₄	38,39
35	Connection capillary	6,8	73	Coupling	38,39
36	Connection capillary	6,8	74	Pump tubing for H ₂ SO ₄	38,39
37	Inlet capillary for detector block.....	6,8,36,39	75	Pump tubing for H ₂ O	38,39
38	Detector block	6,8,36,39	76	Stopper	38
			77	Coupling	38,39

78	Suppressor inlet capillary for eluent.....	39,42	92	Piston seal	159
79	Suppressor outlet capillary for eluent.....	39,42	93	Inlet valve	159,161
80	Suppressor inlet capillary for H ₂ O.....	38,39,42	94	Outlet valve	159,161
81	Suppressor inlet capillary for H ₂ SO ₄	38,39,42	95	Screw holder.....	159
82	Suppressor outlet capillary for H ₂ O.....	39,42	96	Special tool	159
83	Suppressor outlet capillary for H ₂ SO ₄	39,42	97	Special tool.....	159
84	Screw.....	159	98	Valve housing	161
85	Zircon piston.....	159	99	Sealing ring.....	161
86	Spring retainer	159	100	Sleeve	161
87	Spring	159	101	Sapphire sleeve	161
88	Piston cartridge	159	102	Sapphire sphere	161
89	Piston guide sleeve	159	103	Ceramic holder	161
90	Sapphire supporting ring	159	104	Seal.....	161
91	Piston guide sleeve	159	105	Screw nut.....	165
			106	Connection piece	165
			107	Suppressor rotor.....	165
			108	Suppressor holder	165

1 Introduction

1.1 Instrument description

The **Personal IC 790** is a PC-controlled system for ion chromatographic analyses. The two following versions are available:

- **2.790.0010 Personal IC without suppressor module**
- **2.790.0020 Personal IC with suppressor module**

The extremely compact housing of the 790 Personal IC contains everything needed to carry out ion chromatography at a high quality level:

- **Injection valve** – for individual injections or for use with a sample changer such as the Metrohm 766 IC Sample Processor
- **High-pressure pump** – extremely low-pulsation double piston pump with a flow range from 0.2 ... 2.5 mL/min and a maximum pressure of 25 MPa (250 bar)
- **Column chamber** – the perfect insulation of the housing provides not only thermally stable conditions for the separation column but also shields the system against electromagnetic interference
- **Columns** – whether anion columns with or without suppression, separation columns for cations or organic acids – the 790 Personal IC can accommodate them all
- **Suppressor** – the Metrohm Suppressor Module (MSM) is already integrated in the 2.790.0020 Personal IC – pressure-resistant, with fully automatic regeneration, highest performance and optimal reproducibility
- **Peristaltic pump** – integrated two-channel peristaltic pump with a flow rate of 0.5 ... 0.6 mL/min for regeneration and rinsing of the suppressor module built into the 2.790.0020 Personal IC
- **Detector** – conductivity detector with outstanding temperature stability. The detector temperature of 40 °C varies by less than 0.01°C.

All components which come into contact with eluent and sample are metal-free.

The **operation** of the 790 Personal IC takes place via a PC connected to the RS232 interface with the help of the control and evaluation program «**790 Personal IC**». This PC program can be used to open prepared systems for recording and evaluating chromatograms and change them if need be. Time programs can also be created in which a large number of instrument functions can be triggered for each program step. It is also possible to control sample changers such as the 750 Autosampler or the 766 IC Sample Processor via the remote interface.

The operating software for the 790 Personal IC meets all the requirements you could place today on a modern integration software: single or multi-point calibration, internal or external standard, selectable algorithms for non-linear calibration, various integration modes with integration parameters and integration events, different methods for peak recognition, peak editor, free scaling, superimposing several chromatograms, use of sample tables and batch reprocessing; a powerful and GLP-conform report generator with output interfaces for monitor, printer and external databases.

1.2 Parts and controls




Fig. 1: Front of the 790 Personal IC

1 Door to interior

3 Feedthrough for syringe tubing
for connection of the 6.2816.020 syringe
used for aspirating the sample

2 Mains pilot lamp
lit up when instrument switched on

4 Feedthrough for aspirating tubing




Fig. 2: Rear of the 790 Personal IC

5 Opening for detector cable	12 Mains connection plug mains connection, see <i>section 2.4</i>
6 Opening for outlet capillaries for discharge of eluent, regeneration solution, and rinsing solution from the inner compartment	13 Fuse holder changing the fuses, see <i>section 2.4</i>
7 Rear panel opening (closed with plastic stopper) for additional supply and discharge lines to and from the inner compartment	14 RS232 interface connection of the PC
8 Knurled screw for fastening the rear panel 9	15 Connection for detector block
9 Detachable rear panel access to upper part of inner compartment	16 Remote interface remote I/O lines for connection of a sample changer (750, 766)
10 Transport security screws to secure the pump head when the instrument is transported	17 Serial number
11 Mains switch to switch instrument on and off: I = ON 0 = OFF	




Fig. 3: Interior of the 2.790.0010 Personal IC (with permanently attached accessories and 1.733.0110 Detector block)

1	Door to interior	28	Filter unit PEEK (6.2821.100)
18	Inlet capillary for injector PEEK capillary, length $L = 24$ cm	29	Connection capillary PEEK capillary, length $L = 13$ cm
19	Mounting rail for 6.2027.0X0 column holder	30	Connection capillary PTFE capillary for connection of the 6.2816.020 syringe, length $L = 70$ cm
20	Column connection capillary PEEK capillary, length $L = 30$ cm	31	Purge valve
21	Sample loop 10 mL 6.1825.230 PEEK sample loop	32	Aspirating capillary Connection for 6.1834.010 aspirating tubing
22	Injection valve	33	Fastening screws for pump head 34
23	Aspirating tubing PTFE tubing for aspirating the sample; length $L = 65$ cm	34	Pump head (6.2824.100)
24	Coupling (6.2744.120) for connection of the 6.2816.020 syringe	35	Connection capillary Connection pump head – purge valve, fixed mounting
25	Syringe tubing PTFE tubing for connection of the syringe, length $L = 35$ cm	36	Connection capillary in pump head, fixed mounting
26	PEEK coupling (6.2744.040)	37	Inlet capillary for detector block PEEK capillary, fixed mounting
27	Connection capillary PEEK capillary, length $L = 13$ cm	38	Detector block (1.732.0110)




Fig. 4: Interior of the 2.790.0020 Personal IC (with permanently attached accessories and 1.733.0110 Detector block)

1 Door to interior	32 Aspirating capillary Connection for 6.1834.010 aspirating tubing
18 Inlet capillary for injector PEEK capillary, length $L = 24$ cm	33 Fastening screws for pump head 34
19 Mounting rail for 6.2027.0X0 column holder	34 Pump head (6.2824.100)
20 Column connection capillary PEEK capillary, length $L = 30$ cm	35 Connection capillary Connection pump head – purge valve, fixed mounting
21 Sample loop 20 μL 6.1825.210 PEEK sample loop	36 Connection capillary in pump head, fixed mounting
22 Injection valve	37 Inlet capillary for detector block PEEK capillary, fixed mounting
23 Aspirating tubing PTFE tubing for aspirating the sample; length $L = 65$ cm	38 Detector block (1.732.0110)
24 Coupling (6.2744.120) for connection of the 6.2816.020 syringe	39 Suppressor module (inlet and outlet capillaries are not shown)
25 Syringe tubing PTFE tubing for connection of the syringe, length $L = 35$ cm	40 Tubing cartridge (6.2755.000) for 6.1826.060 pump tubing
26 PEEK coupling (6.2744.040)	41 Contact pressure lever for adjusting the contact pressure
27 Connection capillary PEEK capillary, length $L = 13$ cm	42 Holding clamp for locking the tubing cartridge into place
28 Filter unit PEEK (6.2821.100)	43 Snap-action lever for releasing the tubing cartridge
29 Connection capillary PEEK capillary, length $L = 13$ cm	44 Pump drive roller head with contact rollers
30 Connection capillary PTFE capillary for connection of the 6.2816.020 syringe, length $L = 70$ cm	45 Mounting pin for attaching the tubing cartridges
31 Purge valve	

1.3 Information on the Instructions for Use



Please read through these Instructions for Use carefully before you put the 790 Personal IC into operation. The Instructions for Use contain information and warnings to which the user must pay attention in order to assure safe operation of the instrument.

1.3.1 Organization

These **8.790.1003 Instructions for Use** for the 790 Personal IC provide a comprehensive overview of the installation, startup procedure, operation, fault rectification and technical specifications of this instrument. The Instructions for Use are organized as follows:

- Section 1 Introduction**
General description of instrument, parts and controls and safety notes
- Section 2 Installation**
Installation of accessories, putting instrument into operation, connection of sample changers
- Section 3 Operating tutorial**
Introduction to the operation using an example
- Section 4 Operation**
Detailed description of the operation
- Section 5 Notes – Maintenance – Faults**
Notes on ion chromatography, maintenance, fault rectification, diagnostic tests, validation
- Section 6 Appendix**
Technical data, standard equipment, options, warranty, declarations of conformity, index

To find the required information on the instruments, you will find it an advantage to use either the **Table of contents** or the **Index** at the back.





As a supplement to the Instructions for Use, the **8.732.2003 Metrohm Monograph "Ion chromatography"** is also supplied. This provides an introduction to the theoretical fundamentals and general information on separating columns and sample pretreatment.

The **8.732.2013 IC Applications Collection** is also supplied; this contains all the **Application Notes** on the subject of ion chromatography.

You will find detailed information on the separating columns available from Metrohm and on special IC applications in the relevant "**Application Bulletins**", which are available on request free of charge from your Metrohm agency.

1.3.2 Notation and pictograms

The following notations and pictograms (symbols) are used in these Instructions for Use:

Range	Menu item, parameter or entry value
SYSTEM STATE	Program window
<OK>	Button
[Ctrl]	Key
35	Part or control of 790
<u>12</u>	Part or control of 750
<u>26</u>	Part or control of 766
	<p>Hazard This symbol draws attention to a possible danger to life or of injury if the associated directions are not followed correctly.</p>
	<p>Warning This symbol draws attention to possible damage to instruments or instrument parts if the associated directions are not followed correctly.</p>
	<p>Caution This symbol marks important information. First read the associated directions before you continue.</p>
	<p>Comment This symbol marks additional information and tips.</p>

1.4 Safety notes

1.4.1 Electrical safety

While electrical safety in the handling of the 790 Personal IC is assured in the context of the specifications IEC 1010-1 (protection class 1, degree of protection IP20), the following points should be noted:

- **Mains connection**



*Setting of the **mains voltage**, checking the **mains fuse** and the **mains connection** must be effected in accordance with the instructions in section 2.4.*

- **Opening the 790 Personal IC**



*If the 790 Personal IC is connected to the power supply, the instrument must not be opened nor must parts be removed from it, otherwise there is a danger of coming into contact with components which are live. Hence, always disconnect the instrument from all voltage sources before you open it and ensure that the **mains cable is disconnected from mains connection 12 !***

- **Protection against static charges**



Electronic components are sensitive to static charging and can be destroyed by discharges. Before you touch any of the components inside the 790 Personal IC, you should earth yourself and any tools you are using by touching an earthed object (e.g. housing of the instrument or a radiator) to eliminate any static charges which exist.

1.4.2 General precautionary rules

- **Handling of solvents**




Check all lines of the IC system periodically for possible leaks. Follow the relevant instructions regarding the handling of flammable and/or toxic solvents and their disposal.

2 Installation

2.1 Overview

2.1.1 Flow chart

The following flow chart provides an overview of all installation work. You will find more detailed information in the relevant sections.



2.1.2 Connections in the 790 Personal IC

The two following illustrations show the internal connections in the 790 Personal IC in schematic form. The meanings of the various numbered components are given in the detailed illustrations and descriptions in sections 2.2 – 2.10.




Fig. 5: Connecting diagram for 2.790.0010 Personal IC without suppressor




Fig. 6: Connecting diagram for 2.790.0020 Personal IC with suppressor

2.2 Setting up the instrument

2.2.1 Packaging

The 790 Personal IC is supplied together with the separately packed accessories in special packagings containing shock-absorbing foam linings designed to provide excellent protection. The instrument itself is packed in an evacuated polyethylene bag to prevent the ingress of dust. Please store all these special packagings as only they assure transport of the instrument free from damage.

2.2.2 Check

After receipt, immediately check whether the shipment is complete and has arrived without damage (compare with delivery note and list of accessories in *section 6.2*). In the case of transport damage, see instructions in *section 6.4.1 "Warranty"*.

2.2.3 Location

Position the instrument in the laboratory at a location convenient for operation, free from vibrations and protected against a corrosive atmosphere and contamination by chemicals.



To avoid disturbing temperature influences on the insulated column compartment, the instrument must be protected against direct sunlight.

2.3 Attaching the accessories

2.3.1 Connection of detector block

The metal-free **1.732.0110 Detector block** belongs to the scope of delivery of the 790 Personal IC; it must be inserted in the instrument and connected up. Proceed as follows:

1 Note the cell constant

- The cell constant **c = XX,X /cm** is printed on the rear of the detector block. Note this value; it must subsequently be entered in the software in order to ensure that an exact display of the conductivity is obtained (see *section 2.5.3*).

2 Install detector block

- Unscrew the four knurled screws **8** from the top rear panel **9** of the 790 Personal IC and remove rear panel (see *Fig. 2*).
- Position detector block **38** from the back in the space provided in the 790 Personal IC and push fully to the front (see *Fig. 3* and *Fig. 4*).

- Insert the cable permanently attached to the detector block **38** in one of the openings **5** and the outlet capillary in one of the openings **7** of the rear panel **9**.
- Replace rear panel **9** and screw to the 790 Personal IC using the four knurled screws **8**.

3 Connect detector block

- Plug the gray connecting cable permanently attached to the detector block **38** into connection **15** "Detector Block" of the 790 Personal IC and fasten to the instrument by tightening the screws in the cable connector (see *Fig. 2*).

4 Connect waste container

- Lead the outlet capillary of the detector block **38** to a sufficiently large waste container and fix in place.

2.3.2 Connection of syringe and aspirating tubing

For manual filling of the sample loop **22** mounted on the injection valve, the 6.2816.020 Syringe, the syringe tubing **25** with coupling **24** and the PTFE aspirating tubing **23** are needed. These accessories are mounted as follows:

1 Connect syringe to syringe tubing

- Pull syringe tubing **25** connected to injection valve **22** by hand out of feedthrough **3** as far as desired (see *Fig. 3* and *Fig. 4*).
- Mount a PEEK compression fitting **46** (see *section 2.3.4*) to the end of the syringe tubing **25** and screw the coupling **24** (6.2744.120) to the compression fitting **46**.
- Push 6.2816.020 Syringe (without needle) as far as it will go into the connection of coupling **24**.

2 Install aspirating tubing

- Pull aspirating tubing **23** connected to injection valve **22** by hand out of feedthrough **4** as far as desired (see *Fig. 3* and *Fig. 4*).

2.3.3 Connection of the 6.5324.000 Bottle rack (option)

The optional available 6.5324.000 Bottle rack for supply vessels can be placed on top of the 790 Personal IC. The accessories include the supply vessels for eluent (2 L), regeneration solution (1 L) and rinsing solution (1 L).

2.3.4 Connection of PEEK capillaries

For the connections between high-pressure pump and detector block **6.1831.010 PEEK capillaries** (i.d. = 0.25 mm, e.d. = 1/16") are used which are connected using either **6.2744.010 PEEK compression fittings (long)** or **6.2744.070 PEEK compression fittings (short)**. These PEEK connectors can also be used to connect 6.1822.010 PTFE microcapillaries (i.d. = 0.3 mm). Proceed as follows:



Capillaries fitted with new connectors must have a perfectly flat cut surface. To cut PEEK or PTFE capillaries it is best to use the **6.2621.080 Capillary tubing cutter**.

1 Mount compression fitting

Slide a compression fitting **46** (6.2744.010) or a compression fitting **47** (6.2744.070) over the end of the capillary **48** to be fastened as shown in *Fig. 7*.

2 Insert capillary in connection

Push capillary end in the corresponding connection as far as it will go (to avoid dead volume).

3 Tighten compression fitting

Tighten compression fitting **46** or **47** by hand (never use tools).




Fig. 7: Connectors for capillaries

46	Compression fitting (6.2744.010)
47	Compression fitting (6.2744.070)
48	Capillary 6.1831.010 PEEK capillary or 6.1822.010 PTFE microcapillary

2.3.5 Filter unit PEEK

One **6.2821.100 Filter unit PEEK** (see *Fig. 8*) is already installed between the high-pressure pump and the injection valve **22** at the 790 Personal IC. This filter unit serves to avoid contamination by abrasive particles of the piston seals.

The two other filter units PEEK supplied with the 2.790.0020 Personal IC (with suppressor) are installed between the pump tubings of the peristaltic pump and the inlet capillaries for regeneration and rinsing solution (see *section 2.8.2*). These filter units serve to protect the suppressor module from foreign particles and bacterial growth.

The Filter unit PEEK **28** consists of the housing **50** and the two connectors **49** (with filter) and **51** (without filter) screwed into the housing **50**. For the connection of capillaries **48** PEEK compression fittings **46** (6.2744.010) or **47** (6.2744.070) must be used. New connectors **49** with filter are available as an option with the ordering number 6.2821.110 (set of 10).



For the connection of the filter unit, please note the flow direction arrow printed on the housing.




Fig. 8: 6.2821.100 Filter unit PEEK

28	Filter unit PEEK (6.2821.100)	49	Connector with filter (6.2821.110) Part of 6.2821.100 Filter unit
46	Compression fitting (6.2744.010)	50	Housing for filter unit Part of 6.2821.100 Filter unit
48	Capillary 6.1831.010 PEEK capillary or 6.1822.010 PTFE microcapillary	51	Connector without filter Part of 6.2821.100 Filter unit

2.4 Mains connection



Follow the instructions below for connecting to the power supply. If the instrument is operated with a mains voltage set wrongly and/or wrong mains fuse, there is a danger of fire!

2.4.1 Setting the mains voltage

Before switching on the 790 Personal IC for the first time, check that the mains voltage set on the instrument (see Fig. 9) matches the local mains voltage. If this is not the case, you must reset the mains voltage on the instrument as follows:

1 Disconnect mains cable

Disconnect mains cable from mains connection plug **12** of the 790 Personal IC.

2 Remove fuse holder

Using a screwdriver, loosen fuse holder **13** below the mains connection plug **12** and take out completely.

3 Check and change fuse if necessary

Carefully take the fuse installed for the desired mains voltage out of fuse holder **13** and check its specifications (the position of the fuse in the fuse holder is marked by the white arrow imprinted next to the mains voltage range):

100...120 V	1.0 A (slow-blow)	Metrohm No. U.600.0016
220...240 V	0.5 A (slow-blow)	Metrohm No. U.600.0013

4 Insert fuse

Change fuse if necessary and reinsert in fuse holder **13**.

5 Install fuse holder

Depending on the desired mains voltage, insert fuse holder **13** in the 790 Personal IC so that the corresponding mains voltage range can be read normally and the adjacent white arrow points to the white bar imprinted below the fuse holder (see Fig. 9).




Fig. 9: Setting the mains voltage

2.4.2 Fuses

One of the two fuses 1 A/slow-blow for 100...120 V or 0.5 A/slow-blow for 220...240 V is installed in fuse holder **13** of the 790 Personal IC as standard.



Ensure that the instrument is never put into operation with fuses of another type, otherwise there is danger of fire!

For checking or changing fuses, process as described in section 2.4.1.

2.4.3 Mains cable and mains connection

Mains cable

The instrument is supplied with one of three mains cables

- 6.2122.020 with plug SEV 12 (Switzerland, ...)
- 6.2122.040 with plug CEE(7), VII (Germany, ...)
- 6.2133.070 with plug NEMA 5-15 (USA, ...)

which are three-cored and fitted with a plug with an earthing pin. If a different plug has to be fitted, the yellow/green lead (IEC standard) must be connected to protective earth (protection class 1).



Any break in the earthing inside or outside the instrument can make it a hazard!

Mains connection

Plug the mains cable into mains connection plug **12** of the 790 Personal IC (see Fig. 9).

2.4.4 On/off switching of the instrument

The 790 Personal IC is switched on and off using mains switch **11**. When the instrument is switched on, the mains pilot lamp **2** lights up.

2.5 Connection to the PC

2.5.1 Connecting cable



Always switch off 790 Personal IC and PC before you connect the two instruments with the 6.2134.100 Cable.

Connect the RS232 interface **14** at the 790 Personal IC to one of the serial COM ports at the PC using the 6.2134.100 Cable (9 pin/9 pin). If only a 25-pin COM interface is available on the PC then the 6.2125.110 Adapter cable or a commercially available adapter must be used.

2.5.2 Software installation

The PC program «**790 Personal IC 1.0**» is required for the operation of the 790 Personal IC; this is contained on the 6.6035.003 CD included in the accessories. This program runs under Windows 95, Windows 98, Windows NT and Windows 2000 operating systems and is installed as follows:

1 Install program

- Insert 6.6035.003 Installation CD into CD drive.
- Select <Start> and **Run**. Browse for the **setup.exe** file on the installation CD and click on <OK>. Follow the instructions given in the setup program.

The software package will be installed in the desired directory. Icons are created in the program folder and in the startup folder. In addition to the program files, the following folders are installed:

Data	Folder for storage of chromatogram files (*.chw) and batch reprocessing files (*.bar)
Devices	Folder for storage of device files (*.dev)
Methods	Folder for storage of method files (*.mtw)
Reports	Folder for storage of report files (*.txt) and graphic files (*.wmf)
Systems	Folder with subfolders with system files (*.smt) and sample queue files (*.que).

2 Registration

- Please send us your 8.790.8007 Registration card as soon as possible. Only registered users will get updated program versions at a special price.




The installed files (incl. system and method files) are generally not write-protected. To prevent these files from being deleted by mistake, switch on the write-protection or make a backup copy in another directory.

2.5.3 Basic settings


When the program is started for the first time several basic settings must be made for the 790 Personal IC. Proceed as follows:

1 Start program

- Double-click the software icon to start the program. The program window with the opening picture is opened and the Log In window appears on the screen:



- Do not enter any password here, just click on <Log In>. The following window appears:



- Enter company, division and name and click on <OK>. This window appears only one time after software installation.

2 COM port settings

*This step must only be carried out if a different **COM interface from COM1 or COM2** is used for connection to the 790 Personal IC.*


- Click on **Options / 790 Personal IC:COM1** to open the **Links** window:



- Click on **COM1** using the right mouse button and select the **Change** menu item to open a window listing all available COM ports at the PC.
- Select the desired COM port to which the 790 Personal IC has been connected and click on **<OK>**. The window is closed.
- Close the **Links** window by clicking on **<OK>**.
- If desired, repeat this procedure for **Options / 790 Personal IC:COM2**.

3 Open a system

- Click on **File / Open / System** in the main window. Open either the subfolder **790com1** (if the 790 Personal IC is connected to COM1) or the subfolder **790com2** (if the 790 Personal IC is connected to COM2). Select the system file **Anion.smt** or **Cation.smt** and click on **<Open>**. The corresponding system window is opened:




- If the connection between PC and Personal IC is working the message **Hardware settings file for 790 unit with serial number '#####' not found! Create?** will appear. Click **<Yes>** in order to create the configuration file **#####.790** for this instrument.
- If the connection between PC and Personal IC does not work then the message **Detection of hardware failed[790 Personal IC [COM#]]** appears in the **SYSTEM STATE** window. In this case check whether the instrument has been switched on, whether the connection cable is connected up properly and whether the COM interface has been set correctly (see point **2**). Then repeat point **3**.

4 Hardware settings

Of the general hardware settings only the input of the cell constant is described. Standard settings can normally be used for all other parameters.

- Click the 790 icon using the right mouse button and select the **Hardware** item. The **Hardware settings** window is opened:



- In the **Cell constant** field enter the cell constant which is printed on the 1.732.0110 Detector block (see section 2.3.1).
- Click on **<OK>** to close the window and save the settings.

2.6 High-pressure pump



In order to avoid damage to the pump it must never be operated dry. Each time that the pump is switched on always first check that the eluent supply has been connected up correctly and that sufficient eluent is present in the eluent bottle.

2.6.1 Removing the transport security screws

In order to prevent the pump drive from being damaged during transport the pump head is fitted with three transport security screws **10** (see Fig. 2). These transport security screws must be removed before the high-pressure pump is started up. Also remove the red sticker attached to the pump head **34**.



In order to avoid damage to the pump head these three security screws should be attached to the pump head each time that it is to be transported.

2.6.2 Installing the pulsation dampener (option)

To protect the column material of sensitive separating columns (e.g. methacrylate-based columns) against pressure drops caused by the injector, the use of the optional **6.2620.150 Pulsation dampener MF** is recommended. It has to be installed between the high-pressure pump and the injection valve of the 790 Personal IC as follows (see Fig. 10):

1 Install pulsation dampener

- Position the pulsation dampener **52** in the interior of the 790 Personal IC on the base.

2 Connection to the pump

- Unscrew PEEK capillary **27** of coupling **26** and attach it to connection **54** of the pulsation dampener **52**.

3 Connection to injection valve

- Unscrew PEEK capillary **18** of coupling **26** and attach it to connection **53** of the pulsation dampener **52**.



The pulsation dampener is filled with isopropanol and must be rinsed with eluent before connection to a separating column (see section 2.6.4).



The 6.2620.150 Pulsation dampener can be operated in both directions.




Fig. 10: Connection of the pulsation dampener (option)

18	Inlet capillary for injector PEEK capillary, length $L = 24$ cm	31	Purge valve
22	Injection valve	52	Pulsation dampener (6.2620.150)
27	Connection capillary PEEK capillary, length $L = 13$ cm	53	Connection to injection valve
28	Filter unit PEEK (6.2821.100)	54	Connection to purge valve
29	Connection capillary PEEK capillary, length $L = 13$ cm		

2.6.3 Connecting the eluent bottle

To connect the eluent bottle with the high-pressure pump, insert the supplied 6.1834.010 aspirating tubing (i.d. = 1.5 mm, e.d. = 2.5 mm, length = 2.5 m) into one of the openings **5** or **7** in the interior of the 790 Personal IC. Pull the aspirating tubing sufficiently far into the interior of the 790 Personal IC and push at least 5 mm of it onto aspirating capillary **32** (see *Fig. 3* and *Fig. 4*) of the high-pressure pump (it may be necessary to use emery paper). Screw on the 6.2821.090 aspirating filter at the other end of the aspirating tubing and insert it into the eluent bottle.



Only **degassed** (with N_2 , He or vacuum) and **microfiltered** (0.45 μm filter) **eluent**s should be used!

Care must be taken that the **eluent** used is **freely miscible** with any solvent remaining in the pump head (the pump head leaves the factory filled with either isopropanol or methanol/water). If this is not the case then the pump must first be rinsed with a solvent which is miscible with both the previous eluent and the following eluent (e.g. acetone).

2.6.4 Deaerating the pump and rinsing the pulsation dampener


The first time that it is started up the high-pressure pump must be deaerated. Proceed as follows:

1 Prepare for deaeration

- Open the rotary knob on purge valve **31** by approx. $\frac{1}{2}$ turn in the counterclockwise direction (see *Fig. 3* and *Fig. 4*).
- Mount a PEEK compression fitting **46** (see *section 2.3.4*) to the end of capillary **30** and screw coupling **24** (6.2744.120) to the compression fitting **46**.
- Push 6.2816.020 Syringe (without needle) into the connection on coupling **24** until the stop is reached.


2 Open system

- Start the «790 Personal IC» PC program, if it has not been already been started (see *section 2.5.3*).
- Select **File / Open / System** in the main window. Open either the subfolder **790com1** (if the 790 Personal IC is connected to COM1) or the subfolder **790com2** (if the 790 Personal IC is connected to COM2). Select the system file **Anion.smt** or **Cation.smt** and click on **<Open>**. The corresponding system window is opened:



3 Set flow rate to 2 mL/min

- Double-click the 790 icon in the system window to open the window for manual control of the 790 Personal IC (see below).
- Set the flow rate to **2 mL/min** in the **Flow** field.
- Click to **<Send to unit>** to send this value to the 790 Personal IC.



4 Deaerate pump

- Make sure that the 6.1834.010 aspirating tubing for the high-pressure pump has been immersed into the eluent.
- Click the **<On>** button for **IC pump** to switch on the high-pressure pump.
- Use the syringe connected to the purge valve **31** to aspirate air until eluent flows into the syringe.
- Click the **<Off>** button for **IC pump** to switch off the high-pressure pump.

- Close the rotary knob on purge valve **31** by turning it in a clockwise direction (see *Fig. 3* and *Fig. 4*).
- Remove the syringe from coupling **24**.

5 Rinse pulsation dampener (if present)

- Place a beaker beneath the column connection capillary **20**.
- Click the **<On>** button for **IC pump** to switch on the high-pressure pump and rinse the pulsation dampener **52** filled with isopropanol for ca. 10 min with eluent.
- Click the **<Off>** button for **IC pump** to switch off the high-pressure pump.

6 Reduce flow rate

- Reset the original flow rate under **Flow** (e.g. **0.5 mL/min**).
- Click **<Send to unit>** to send this value to the 790 Personal IC.

2.7 Precolumns and separating columns

2.7.1 General information on precolumns

The use of easily exchangeable precolumns protects the separating columns and prolongs their lifetime. The precolumns available from Metrohm (see *section 6.3.2*) are either real precolumns or precolumn cartridges, which are used together with the 6.2821.040 Cartridge head or the 6.2828.010 Precolumn cartridge holder.



New IC precolumns are normally filled with solution and sealed at both ends. Before the precolumn is installed in the system, it must be ensured that this solution is freely miscible with the eluent used (check manufacturer's specifications).

2.7.2 Precolumns with cartridge head

The 6.1005.020, 6.1005.040, 6.1005.050, 6.1007.010 and 6.1010.010 Precolumn cartridges are installed in the 6.2821.040 Cartridge head as follows (see *Fig. 11*):

1 Prepare separating column

- Remove end caps from separating column **63**.
- Unscrew fastening screw from column inlet.
- Take steel connector **56** for ferrule out of fastening screw.

2 Install cartridge

- Remove end caps for precolumn cartridge **59** (the steel mesh **58** and gaskets **57** are already installed in the cartridge).
- Mount steel spacer **60** on Manufit pressure screw **62** (the steel mesh **61** and gaskets **57** are already installed in the pressure screw).
- Mount precolumn cartridge **59** on the steel spacer **60** (comply with flow direction if specified on the precolumn).
- Screw Manufit pressure screw **62** firmly to separating column **63**.
- Mount steel connector **56** for ferrule on the inlet side of the precolumn cartridge **59**.
- Screw on Manufit housing **55** with Manufit pressure screw **62**.




Fig. 11: Installation of precolumn cartridges with cartridge head

20	Column connection capillary from injector	59	Precolumn cartridge
46	Compression fitting (6.2744.010)	60	Steel spacer (6.2821.080)
55	Manufit housing	61	4 Steel meshes (6.2821.020)
56	Steel connector (of IC separating column)	62	Manufit pressure screw
57	PTFE gasket (6.2821.010)	63	IC separating column
58	2 Steel meshes (6.2821.020)		

3 Connect precolumn

- Fit PEEK compression fitting **46** to the column connection capillary **20** mounted on the injection valve **22** (see section 2.3.4).
- Screw column connection capillary **20** to Manufit housing **55**.



The precolumn cartridge built into the cartridge head can only be rinsed together with the separating column (see section 2.7.7/8). This increases the rinsing time by approx. 20 min.

2.7.3 Precolumn glass cartridges with cartridge holder

The precolumn glass cartridge METROSEP Anion Dual 1 (6.1006.030) is inserted in the 6.2828.010 Precolumn cartridge holder as follows (see Fig. 12):

1 Insert cartridge

- Remove end fittings **64** from screw caps **65**.
- Insert precolumn cartridge **67** in sleeve **66** (mark flow direction of the column cartridge on sleeve, the column cartridge should always be operated in the same flow direction).
- Screw both screw caps **65** loosely to sleeve **66** by hand.
- Screw both end fittings **64** into screw caps **65** so that their capillary ends are seated in the PTFE seals of the cartridge.
- Tighten both screw caps **65** by hand.
- Tighten both end fittings **64** firmly using a screwdriver.

2 Connect precolumn

- Provide the column connection capillary **20** mounted on the injection valve with a PEEK compression fitting **46** (see section 2.3.4) and screw it tightly onto end fitting **64** on the inlet side of the precolumn.
- Cut off as short a piece as possible from PEEK capillary **48** and fit this with two PEEK compression fittings **46** (see section 2.3.4).
- Screw prepared capillary **48** to end fitting **64** on the other end of the precolumn.

3 Rinse precolumn

- Place a beaker beneath the outlet capillary of the precolumn.
- Open software window for manual system control.
- If necessary, modify **Flow rate** to the value suited for the inserted separating column and click on **<Send to unit>** to send this value to the 790 Personal IC.
- Switch on high-pressure pump (**IC pump**) by clicking **<On>** and rinse precolumn with eluent for ca. 10 min.
- Switch off high-pressure pump by clicking **<Off>**.



The rinsed precolumn can be connected directly to the separating column instead of by means of a PEEK capillary. Proceed as follows:

- *Screw off PEEK capillary from end fitting **64** of the precolumn.*
- *Screw off end fitting **64** from screw cap **65** of the precolumn.*
- *Screw connection piece **68** onto screw cap **65**.*
- *Screw screw cap **69** of separating column **63** onto connection piece **68**.*




Fig. 12: Installation of precolumn glass cartridges with cartridge holder

20	Column connection capillary from injector	66	Sleeve for precolumn cartridge
46	Compression fitting (6.2744.010)	67	Precolumn cartridge (6.1006.030)
48	PEEK capillary (6.1831.010)	68	Connection piece for connection precolumn – column
64	End fitting	69	Screw cap for column
65	Screw cap for precolumn		

2.7.4 IC anion precolumn SUPERSEP

The 6.1009.010 IC Anion Precolumn SUPERSEP has two connections for PEEK capillaries and is installed as follows:

1 Connect precolumn

- Remove end caps from the precolumn.
- Fit PEEK compression fitting **46** to column connection capillary **20** (see section 2.3.4).
- Screw precolumn to column connection capillary **20**.
- Cut a piece from the 6.1831.010 PEEK capillary as short as possible and fit with PEEK compression fittings **46** (see section 2.3.4).
- Fasten the prepared capillary to the other end of the precolumn.

2 Rinse precolumn

- Place a beaker beneath the outlet capillary of the precolumn.
- Open software window for manual system control.
- If necessary, modify **Flow rate** to the value suited for the inserted separating column and click on **<Send to unit>** to send this value to the 790 Personal IC.
- Switch on high-pressure pump (**IC pump**) by clicking **<On>** and rinse precolumn with eluent for ca. 10 min.
- Switch off high-pressure pump by clicking **<Off>**.

2.7.5 General information on separating columns



New IC separating columns are normally filled with solution and sealed at both ends. Before the column is installed in the system, it must be ensured that this solution is freely miscible with the eluent used (check manufacturer's specifications).

The IC separating columns and precolumns currently available from Metrohm are listed in *section 6.3.2*. A test chromatogram and an information leaflet is provided with each column. You will find additional information concerning these columns in the 8.732.2003 Metrohm Monograph «Ion chromatography» and in special «Application Bulletins», which are available on request free of charge from your local Metrohm agency.



When you install the column, always ensure that this is inserted correctly in accordance with the flow direction shown (arrow must point upwards).

2.7.6 Selection of the sample loop

Selection of the sample loop depends on the separating column used. Normally, the following sample loops are used:

Cation columns	10 µL
Anion columns with suppressor	20 µL
Anion columns without suppressor	100 µL

Depending on the instrument version, the following sample loops are installed in the 790 Personal IC:

Version	Sample loop	Volume
2.790.0010	6.1825.230 (PEEK)	10 µL
2.790.0020	6.1825.210 (PEEK)	20 µL

If desired, the built-in sample loop can be replaced by one of the sample loops available as an option (see *section 6.3.1*).

2.7.7 Installation of the separating column without suppressor

With the 2.790.0010 Personal IC without suppressor module, the IC separating column is installed as follows (see Fig. 13):

1 Connect column to injector

- Remove end caps from column **63**.
- *without precolumn:*
Screw inlet end of separating column **63** (note flow direction) to column connection capillary **20** mounted on the injector.
- *with precolumn in cartridge head:*
Install separating column **63** (note flow direction) in the cartridge head (see Fig. 11) as described in section 2.7.2.
- *with precolumn in cartridge holder:*
Screw inlet end of separating column **63** (note flow direction) to the precolumn fixed in the cartridge holder as described in section 2.7.3 (see Fig. 12).
- *with IC anion precolumn SUPERSEP:*
Screw inlet end of separating column **63** (note flow direction) to the precolumn installed as described in section 2.7.4.

2 Rinse column

- Place a beaker beneath the column outlet.
- Open software window for manual system control.
- If necessary, modify **Flow rate** to the value suited for the inserted separating column and click on **<Send to unit>** to send this value to the 790 Personal IC.
- Switch on high-pressure pump (**IC pump**) by clicking **<On>** and rinse column with eluent for ca. 10 min.
- Switch off high-pressure pump by clicking **<Off>**.

3 Connect column to detector block

- Screw outlet end of separating column **63** to the inlet capillary **37** permanently mounted on the detector block **38**.

4 Fix column

- Insert one or two column holders **70** (6.2027.030, 6.2027.040 or 6.2027.050) in the mounting rails **19** and fasten separating column **63** in the column holder **70**.




Fig. 13: Installation of column without suppressor

2.7.8 Installation of the separating column with suppressor

With the 2.790.0020 Personal IC with suppressor module, the IC separating column is first connected to the injector or precolumn. The connection to the suppressor module and the detector block is described in section 2.8.

1 Connect column to injector

- Remove end caps from column **63**.
- Screw inlet end of separating column **63** (note flow direction) to column connection capillary **20** or to the already installed precolumn (procedure see section 2.7.7).

2 Rinse column

- Place a beaker beneath the column outlet.
- Open software window for manual system control.
- If necessary, modify **Flow rate** to the value suited for the inserted separating column and click on **<Send to unit>** to send this value to the 790 Personal IC.
- Switch on high-pressure pump (**IC pump**) by clicking **<On>** and rinse column with eluent for ca. 10 min.
- Switch off high-pressure pump by clicking **<Off>**.

3 Fix column

- Insert one or two column holders **70** (6.2027.030, 6.2027.040 or 6.2027.050) in the mounting rails **19** and fasten separating column **63** in the column holder **70**.

2.8 Suppressor module

2.8.1 General information on suppressor module

The **Metrohm Suppressor Module MSM** for chemical suppression installed in the 2.790.0020 Personal IC comprises a total of 3 suppressor units which are in turn used for suppression, regenerated with sulfuric acid and rinsed with water. To record every new chromatogram under comparable conditions, work is normally carried out with freshly regenerated suppressor. Switching is either automatic together with the valve switching or manual.



The suppressor units must never be regenerated with H_2SO_4 in the same flow direction used for the eluent. You should thus always install the inlet and outlet capillaries as described in section 2.8.4 according to the scheme shown in Fig. 16.

For operation of the suppressor module, the **two-channel peristaltic pump** built into the 2.790.0020 Personal IC is used which conveys the regeneration solution (normally **20 mmol/L H_2SO_4**) and the rinsing solution (normally **dist. H_2O**) to the suppressor units (flow rate of 0.5 mL/min).

The three inlets and outlets numbered 1...3 on the suppressor module each have 2 permanently mounted PTFE capillaries, which must be connected as described in section 2.8.4 (see Fig. 15 and Fig. 16).

To avoid contamination of the suppressor module by foreign particles or bacterial growth, the two **6.2821.100 Filter units PEEK** (see section 2.3.5) supplied with the 2.790.0020 Personal IC must be installed between the peristaltic pump and the inlet capillaries of the suppressor module.



The suppressor module must never be switched in the dry state as there is a danger of blocking. Before every switching operation of the suppressor module, the three suppressor units must have been rinsed for at least $\frac{1}{2}$ h with eluent, regeneration and rinsing solution.

2.8.2 Preparation of the peristaltic pump

Before start-up the accessories for the 2-channel peristaltic pump built into the 790 Personal IC must be mounted according to Fig. 14. Proceed as follows:

1 Attach pump tubings

- Loosen both tubing cartridges **40** mounted above pump drive **44** from the holding clamp **42** by pressing down snap-action lever **43** and remove from mounting pin **45** (see Fig. 15).

- Press contact pressure lever **41** on both tubing cartridges down as far as it will go.
- Insert the pump tubings **74** and **75** (6.1826.060) into each of the tubing cartridges as shown in *Fig. 14*. The white-yellow stopper **76** must click into the corresponding holder on the left-hand side of the tubing cartridge.
- Place the tubing cartridges on mounting pin **45** and press down on the right-hand side until snap-action lever **43** clicks into position on holding clamp **42**. Take care that no kinks are formed in the pump tubing.

2 Install Filter units PEEK

- Mount a coupling **77** (6.2744.110) to the outlet end of the two pump tubings **74** and **75**.
- Attach a filter unit PEEK **28** (see *section 2.3.5*) to the other end of this coupling (note flow direction!).




Fig. 14: Installing pump tubings

28	Filter unit PEEK (6.2821.100)	73	Coupling (6.2744.030)
40	Tubing cartridge	74	Pump tubing (6.1826.060) for H ₂ SO ₄
41	Contact pressure lever	75	Pump tubing (6.1826.060) for H ₂ O
43	Snap-action lever	76	Stopper (white-yellow)
46	PEEK compression fitting (6.2744.010)	77	Coupling (6.2744.110)
71	Aspirating tubing for H ₂ O	80	Suppressor inlet capillary for H ₂ O
72	Aspirating tubing for H ₂ SO ₄	81	Suppressor inlet capillary for H ₂ SO ₄




Fig. 15: Installation of column with suppressor

19	Mounting rail for column holder 70	63	Separating column
20	Column connection capillary PEEK capillary (6.1831.010; 30 cm)	70	Column holder (6.2027.0X0)
22	Injection valve	71	Aspirating tubing for H₂O 6.1803.020 PTFE tubing
26	PEEK coupling (6.2744.040)	72	Aspirating tubing for H₂SO₄ 6.1803.020 PTFE tubing
28	Filter unit PEEK (6.2821.100)	73	Coupling (6.2744.030)
37	Inlet capillary to detector block (fixed mounting)	74	Pump tubing (6.1826.060) for H₂SO₄
38	Detector block (1.733.0110)	75	Pump tubing (6.1826.060) for H₂O
39	Suppressor module	77	Coupling (6.2744.110)

40	Tubing cartridge (6.2755.000) for pump tubings 74/75	78	Suppressor inlet capillary for eluent
41	Contact pressure lever for adjusting the contact pressure	79	Suppressor outlet capillary for eluent
42	Holding clamp for locking the tubing cartridge into place	80	Suppressor inlet capillary for H₂O
43	Snap-action lever for releasing the tubing cartridge	81	Suppressor inlet capillary for H₂SO₄
44	Pump drive Roller head with contact rollers	82	Suppressor outlet capillary for H₂O
45	Mounting pin for attaching the tubing cartridge	83	Suppressor outlet capillary for H₂SO₄

2.8.3 Connection of supply bottles

The supply lines for the regeneration and rinsing solution between the storage bottles and the peristaltic pump are installed as follows:

1 Prepare supply bottle for H₂SO₄

- Prepare regeneration solution suited for the desired application and separating column (normally 20 mmol/L H₂SO₄).
- Fill regeneration solution into supply bottle and label the bottle.

2 Connect aspirating tubing for H₂SO₄

- Prepare aspirating tubing **72**: Cut a piece of the 6.1803.020 PTFE tubing to the required length (normally ca. 120 cm).
- Insert one end of the aspirating tubing **72** into the regeneration solution storage bottle and tighten it so that the tubing is firmly held.
- Insert the free end of aspirating tubing **72** into one of the openings **5** of the 790 Personal IC (see Fig. 2) and pull it sufficiently far into the interior.
- Mount a coupling **73** (6.2744.030) to the inlet end of the rear pump tubing **74**.
- Mount a 6.2744.010 Compression fitting at the end of the aspirating tubing **72** and screw this compression fitting on to the coupling **73** (see Fig. 15).

3 Prepare supply bottle for H₂O

- Prepare rinsing solution suited for the desired application and separating column (normally dist. H₂O).
- Fill rinsing solution into a supply bottle and label the bottle.

4 Connect aspirating tubing for H₂O

- Prepare aspirating tubing **71**: Cut a piece of the 6.1803.020 PTFE tubing to the required length (normally ca. 120 cm).
- Insert one end of the aspirating tubing **71** into the rinsing solution storage bottle and tighten it so that the tubing is firmly held.
- Insert the free end of aspirating tubing **71** into one of the openings **5** of the 790 Personal IC (see *Fig. 2*) and pull it sufficiently far into the interior.
- Mount a coupling **73** (6.2744.030) to the inlet end of the front pump tubing **75**.
- Mount a 6.2744.010 Compression fitting at the end of the aspirating tubing **71** and screw this compression fitting on to the coupling **73** (see *Fig. 15*).

2.8.4 Connection of the suppressor module

The three inlets and outlets numbered 1...3 on the suppressor module **39** each have 2 permanently mounted PTFE capillaries, which must be connected as follows (see *Fig. 15* and *Fig. 16*).

1 Inlet capillary for eluent

- Screw inlet capillary **78** marked with "Eluent" of suppressor module **39** to outlet end of separating column **63** using a 6.2744.010 Compression fitting.

2 Outlet capillary for eluent


- Screw outlet capillary **79** marked with "Detector" of suppressor module **39** to coupling **26** using a 6.2744.010 Compression fitting.
- Screw inlet capillary **37** of detector block **38** to other end of coupling **26**.

3 Inlet capillary for H₂SO₄

- Attach inlet capillary **81** marked with "H₂SO₄" of suppressor module **39** using a 6.2744.010 Compression fitting to the filter unit PEEK **28** connected to the rear pump tubing **74**.

4 Outlet capillary for H₂SO₄

- Pull outlet capillary **83** marked with "Waste" of the suppressor module **39** from below through one of the openings **6** out of the inner compartment of the 790 Personal IC.
- Lead outlet capillary **83** to a sufficiently large waste container and fix it in place.



78 Suppressor inlet capillary for eluent

79 Suppressor outlet capillary for eluent

80 Suppressor inlet capillary for H₂O

81 Suppressor inlet capillary for H₂SO₄

82 Suppressor outlet capillary for H₂O

83 Suppressor outlet capillary for H₂SO₄

Fig. 16: Connections at suppressor module

5 Inlet capillary for H₂O

- Attach inlet capillary **80** marked with "H₂O" of suppressor module **39** using a 6.2744.010 Compression fitting to the filter unit PEEK **28** connected to the front pump tubing **75**.

6 Outlet capillary for H₂O

- Pull outlet capillary **82** marked with "Waste" of the suppressor module **39** from below through one of the openings **6** out of the inner compartment of the 790 Personal IC.
- Lead outlet capillary **82** to a sufficiently large waste container and fix it in place.

7 Fasten capillaries to the side walls

- If necessary, the two aspirating tubings **71** and **72** can be fixed in the required position in the interior with the help of a Y.107.0150 self-adhesive strap.
- If necessary, the two outlet capillaries **82** and **83** can be fixed in the required position in the interior with the help of a Y.107.0150 self-adhesive strap.

2.9 Putting into operation

2.9.1 Putting into operation without suppressor (cation system)


Before sample solutions can be injected at the **2.790.0010 Personal IC** (without suppressor), the entire system must be tested for leaks and then conditioned with eluent until the baseline is stable. Proceed as follows:

1 Open system

- Start the «790 Personal IC» PC program, if it has not already been started (see *section 2.5.3*).
- Select **File / Open / System** in the main window. Select the system file suited for the inserted cation column **Cation.smt** and click on <Open>.

2 Open control window

- Double-click the 790 icon in the system window. The control window for manual control of the 790 Personal IC appears, which indicates conductivity, pressure and current system parameters.



3 Start system

- Make sure, that the aspirating tubing 6.1834.010 for the high-pressure pump is immersed in the eluent.
- Select **Startup hardware (Measure baseline)** from the **Control** menu in the system window. The high-pressure pump is started, at the same time, a chromatogram window is opened where the baseline is recorded continuously.

4 Check for leaks

- Check all capillaries and their connections between the high-pressure pump and the detector block for escaping liquid. If eluent escapes anywhere, the appropriate compression fitting must be tightened further or changed.

5 Condition system

- Rinse the system with eluent until the desired stability of the baseline is reached (normally 30...60 min; if the eluent is changed, the establishment of the ion exchanger equilibrium on the separating column can take longer).
- The instrument is now ready for sample determinations using the selected system.

2.9.2 Putting into operation with suppressor (anion system)


Before sample solutions can be injected at the **2.790.0020 Personal IC** (with suppressor), the entire system must be tested for leaks and then conditioned with eluent until the baseline is stable. At the same time, the suppressor module must be conditioned. Proceed as follows:

1 Open system "Prep MSM.smt"

- Start the «790 Personal IC» PC program, if it has not already been started (see section 2.5.3).
- Select **File / Open / System** in the main window. Open the folder **790com1** or **790com2** and select the system file **Prep MSM.smt**. Click on <Open>.

2 Open control window

- Double-click the 790 icon in the system window. The control window for manual control of the 790 Personal IC appears, which indicates conductivity, pressure and current system parameters.



3 Start system

- Make sure, that the aspirating tubing 6.1834.010 for the high-pressure pump is immersed in the eluent.
- Select **Startup determination** from the **Control** menu in the system window. The high-pressure pump and the peristaltic pump are started, at the same time, a chromatogram window is opened where the baseline is recorded continuously. The suppressor module is automatically switched to the next position every 20 min and conditioned in this way.

4 Set contact pressure for pump tubings

- Press contact pressure lever **41** on both tubing cartridges **40** upwards until regeneration and rinsing solution just start to be drawn in.
- Then press contact pressure lever **41** upwards until it clicks once more to obtain optimal contact pressure.

5 Check for leaks

- Check all capillaries and their connections between the high-pressure pump and the detector block for escaping liquid. If eluent escapes anywhere, the appropriate compression fitting must be tightened further or changed.

6 Condition system

- Rinse the system with eluent until the desired stability of the baseline is reached (normally 30...60 min; if the eluent is changed, the establishment of the ion exchanger equilibrium on the separating column can take longer). After this time the suppressor module is sufficiently conditioned too.
- The instrument is now ready for sample determinations using the selected system.



Pump tubings are consumable material with a lifetime, which depends on the contact pressure. This is why the tubing cartridges should be raised completely by loosening snap-action lever on the right-hand side if the pump is to remain switched off for a considerable length of time (the set contact pressure remains unchanged).

2.10 Connection of external devices

2.10.1 Connection of the 750 Autosampler

The 750 Autosampler available from Metrohm as an option is an automatic sampler for ion chromatography. It accommodates max. 128 samples each of 730 μL , which are automatically transferred to the sample loop attached to the injection valve of the 790 Personal IC. The samples on the 750 Autosampler are changed by a remote signal that is outputted by the 790 Personal IC (790 Personal IC as "Master").

In order to connect the 750 Autosampler the 6.2128.160 Cable that is available as an option is additionally required. The 750 Autosampler is connected to the 790 Personal IC as follows:

1 Electrical connection 790 – 750

- Connect 6.2140.010 Connector plug to contact closure strip **10** of the 750 Autosampler (see *750 Instructions for Use*).
- Connect one end of the 6.2128.160 Cable to the connections 13 "EXTERNAL INJECT INPUT" and 14 "GND" on the connector plug. The two cable ends are appropriately inscribed with "750" (white) and "GND" (brown).
- Connect the other end of the 6.2128.160 Cable to the remote interface **16** of the 790 Personal IC (see *Fig. 2*).

2 Tubing connection 750 – injection valve

- Pull aspiration tubing **23** completely out from feedthrough **4** and screw it off from connection "1" of injection valve **22**.
- Cut transfer tubing **8** mounted on the 750 Autosampler (see *750 Instructions for Use*) to the shortest possible length between the needle of the 750 Autosampler and injection valve **22** of the 790 Personal IC.
- Pull the free end of transfer tubing **8** through the feedthrough **4** on the 790 Personal IC and screw it onto connection "1" of injection valve **22** using a 6.2744.010 PEEK compression fitting.

3 Tubing connection injection valve – waste

- Remove coupling **24** and 6.2744.010 Compression fitting from syringe tubing **25**.
- Lead free end of syringe tubing **25** to the waste container and fix in place.


4 Settings at the 750 Autosampler

- Set parameter "Run time" to **0.0** to activate the external start of the 750 Autosampler (see *section 3.5 of the 750 Instructions for Use*). Filling the sample loop of the injection valve only takes place when the impulse produced by the 790 Per-

sonal IC at remote output 3 (see point **5**) makes contact at the "EXTERNAL INJECT INPUT" of the 750 Autosampler. The run time will then be determined externally.

5 Settings in the «790 Personal IC» program

- A time program must be drawn up for the selected system which first switches the injection valve to the "Fill" position, then produces an impulse at remote lead 3 to start sample addition by the 750 Autosampler and finally switches the injection valve to the "Inject" position. The following example shows such a program:



2.10.2 Connection of the 766 IC Sample Processor

The 766 IC Sample Processor available from Metrohm as an option is an automatic sampler for ion chromatography. It accommodates max. 127 samples each of 2.5 mL or 11 mL, which are automatically transferred to the sample loop attached to the injection valve of the 790 Personal IC. Sample changing and filling the sample loop are started by a remote signal produced by the 790 Personal IC (790 Personal IC as "Master").

The 6.2141.110 Cable is required for connecting the 766 IC Sample Processor. The 766 IC Sample Processor is connected to the 790 Personal IC as follows:

1 Electrical connection 790 – 766

- Connect the end of the 6.2141.110 Cable marked with "766" to the remote interface **22** of the 766 IC Sample Processor (see *766 Instructions for Use*).
- Connect the end of the 6.2141.110 Cable marked with "732/1" to the remote interface **16** of the 790 Personal IC (see *Fig. 2*).

2 Tubing connection 766 – injection valve

- Pull aspiration tubing **23** completely out from feedthrough **4** and screw it off from connection "1" of injection valve **22** (see *Fig. 3* and *Fig. 4*).

- Cut PEEK capillary **18** installed at the 766 IC Sample Processor (see *766 Instructions for Use*) to the desired length.
- Pull the free end of PEEK capillary tubing **18** through the feedthrough **4** of the 790 Personal IC and screw it onto connection "1" of injection valve **22** with the help of a 6.2744.010 PEEK compression fitting.

3 Tubing connection injection valve – waste

- Remove coupling **24** and 6.2744.010 Compression fitting from syringe tubing **25**.
- Lead free end of syringe tubing **25** to the waste container and fix in place.


4 Settings at the 766 IC Sample Processor

- For operation with the 790 Personal IC it is recommended that the following program "790" is used with the 766 IC Sample Processor:

parameters		
method	790	– Method name
number of samples:	Rack	– Number of samples to be processed (entire sample rack)
>start sequence		
1 CTL:Rm:	INIT	– Initialize remote interface
>sample sequence		
1 SCN:Rm	: Wait1	– Waiting until 790 Personal IC sends signal on remote line 3
2 MOVE 1	: sample	– Move needle to sample position
3 LIFT: 1	: work mm	– Place lift with needle to working position
4 SCN:Rm	: Wait1	– Waiting until 790 Personal IC sends signal on remote line 3
5 PUMP 1.1	: 120 s	– Fill sample loop with sample during 120 s
>final sequence		
>changer settings		– Settings for changer functions
rack number	0	
lift rate 1	12 mm/s	
shift rate	20	
>manual stop		– Reaction to manual stop
CTL Rmt:	*****	
CTL RS232:	-----	

5 Settings in the «790 Personal IC» program

- A time program must be drawn up for the selected system which first produces an impulse at remote lead 3 to start sample changing on the 766 IC Sample Processor and then switches the injection valve to the "Fill" position. A further impulse is then produced to start filling the sample loop within 120 s. At the end, the injection valve is switched to the "Inject" position, which also starts data acquisition.



6 Run conditions

In order that the joint operation of 790 – 766 functions properly the program must be started first on the 766 IC Sample Processor, then the "Sample Queue" on 790 Personal IC.

3 Operating tutorial



This section introduces you to the operation of the 790 Personal IC by means of a brief operating tutorial which describes the basic operating steps needed for the recording of an ion chromatogram using one of the system files supplied.

The determination of the cationic content of a drinking water sample with the METROSEP Cation 1-2 column is used as an illustrative example. Please note that the steps and parameter settings described apply only to this example. If you use a different separating column and system file, the procedures described in the tutorial must be modified appropriately.

For further explanations of the operation, please refer to section 4.

3.1 Requirements

For the determination of cations in drinking water described in this tutorial, the following instruments, accessories and solutions are required:

- **2.790.0010 Personal IC**
(or 2.790.0020 Personal IC without connection of suppressor module)
- **6.1825.230 Sample loop (10 µL, PEEK)**
already integrated in the 2.790.0010 Personal IC
- **6.1006.100 IC cation column METROSEP Cation 1-2 (6.1010.000)**
- **Eluent**
4 mmol/L tartaric acid / 1 mmol/L dipicolinic acid in dist. H₂O
flow: 1.0 mL/min
- **Standard**
Standard solution with 1 mg/L Li⁺, 5 mg/L Na⁺, 5 mg/L NH₄⁺
10 mg/L K⁺, 10 mg/L Ca²⁺, 10 mg/L Mg²⁺ (in dist. H₂O)

3.2 Preparations

Before you start this brief tutorial, the entire IC system must be correctly installed as described in *section 2*. In what follows, the most important points for the installation are described once again (for details, see the sections mentioned).

1 Install 790 Personal IC

- ⇒ Setting up instrument *section 2.2*
- ⇒ Installing and connecting detector block *section 2.3.1*
- ⇒ Mounting syringe and aspirating tubing *section 2.3.2*
- ⇒ Mains connection *section 2.4*
- ⇒ Connection to PC *section 2.5*

2 Prepare eluent

- ⇒ Preparing eluent:
4 mmol/L tartaric acid / 1 mmol/L
dipicolinic acid in dist. H₂O
- ⇒ Microfiltering and degassing eluent *section 5.1.3*

3 Install high-pressure pump

- ⇒ Removing transport security screws *section 2.6.1*
- ⇒ Mounting pulsation dampener (option) *section 2.6.2*
- ⇒ Installing eluent supply *section 2.6.3*
- ⇒ Degassing pump *section 2.6.4*

4 Connecting separating column and suppressor

- ⇒ Connecting IC cation column *section 2.7.7*
- ⇒ Conditioning system *section 2.9.1*

3.3 Calibration

After the complete IC system has been installed as described in section 3.2, the first calibration can be started. A standard solution is required for this; it should contain the substances to be determined in approximately the same concentrations in which they can be expected in the sample.

In our example of drinking water determination using the METROSEP Cation 1-2 IC cation column a 10 µL sample loop is used; this is filled with the following standard solution:

- 1 mg/L Li⁺
- 5 mg/L Na⁺ and NH₄⁺
- 10 mg/L K⁺, Ca²⁺ and Mg²⁺

Please note that all program displays refer to the condition in which the system **Cation.smt** is loaded for the first time. If you want to work through this tutorial at a later date and the system **Cation.smt** has been altered in the meantime then differences with respect to the program display and the parameter values may occur.

In the description of the calibration procedure, it is assumed that the PC and 790 Personal IC are not in operation and that the system must first be conditioned again. If this is not the case (e.g. if you start the tutorial immediately after conditioning) then you can skip steps **1** to **5**.

1 Switch on 790 Personal IC

⇒ Switch on 790 Personal IC with mains switch **11** on the rear of the instrument. After the instrument has been switched on the mains pilot lamp **2** lights up.

2 Switch on PC

⇒ Switch on PC and start «790 Personal IC» program.

3 Open system "Cation.smt"

⇒ Select **File / Open / System** in the main window. Select the system file **Cation.smt** in the folder **790com1** or **790com2** and click on **<Open>**.

4 Start system "Cation.smt"

⇒ Select **Start determination** of the **Control** menu in the system window. The high-pressure pump is started, at the same time, a chromatogram window opens where the baseline is recorded continuously.

5 Condition system "Cation.smt"


⇒ Rinse IC system with eluent until the desired stability of the baseline is achieved (at least 1 h).

6 Start determination

⇒ Click on **<Start>** in the system window or select **Start determination** of the **Control** menu. A chromatogram window opens where the baseline is recorded continuously. The status bar of this window shows the message **Measure(Baseline)**, beside this running time, analysis time, absolute conductivity and number of measuring points per second are displayed. The message **Waiting for INJECT [790COM#]** appears in the **SYSTEM STATE** window

7 Enter information for determination

⇒ Enter the desired information concerning the sample in the automatically opened **Edit sample description** window. Set **Calibration level** to **0** and confirm with **<OK>**.



8 Condition system


⇒ Let the system run for several minutes. If the system has already been conditioned then a stable conductivity value of approx. 700 $\mu\text{S}/\text{cm}$ will be obtained after a few minutes.

⇒ Double-click the chromatogram or select **View all** of the **View** menu. The sensitivity is set automatically so that all measuring points are visible.




*This function is only available if the integration has been started. In the present example the integration only starts after a delay period **Delay = 2.2 min.***

⇒ Select the desired sensitivity using the cursor keys **<↑>** or **<↓>**. In the chromatogram window the following display is seen, for example:




9 Switch injection valve to "FILL" position

⇒ Click on the  button of the 790 icon in the system window to switch the injection valve to the "FILL" position. At the same time, the suppressor module is switched to the next position.


10 Fill sample loop


⇒ Immerse the aspirating tubing **23** in the standard solution.
 ⇒ Using the syringe fixed to syringe tubing **25**, siphon in ca. 1 mL standard solution.

11 Switch injection valve to "INJECT" position


⇒ Click on the  button of the 790 icon in the system window to switch the injection valve to the "INJECT" position. At the same time, the data recording is started automatically. The status bar of the chromatogram window shows the message **Measure**, the **SYSTEM STATE** window shows the message **INJECT done [790COM#]**.
 ⇒ When the analysis time of 14 min set in the method has expired, the recorded chromatogram is automatically evaluated and stored. The chromatogram window then closes.


12 Open chromatogram

⇒ Click on  or **File / Open / Chromatogram** in the main window. Select the chromatogram file ***.chw** just recorded and click on **<OK>**. The chromatogram window is opened where the found peaks are numbered and the baselines are drawn in.




13 Modify integration parameters


- ⇒ Click on  or select **Integration** from the **Method** menu in the main window to open the integration parameters window.
- ⇒ If required (e.g. if the lithium peak is not detected), change the time delay **Delay** before starting peak integration.
- ⇒ Click on **<Apply>**. The **Integration parameters** window is opened and the chromatogram is reintegrated.



- ⇒ Repeat this procedure for all other integration parameters until the result satisfies your expectations. Click on **<Apply>** for each parameter change.
- ⇒ Close the **Integration parameters** window with **<OK>**.

14 Modify peak allocation

⇒ Click on  or select **Calibration / Components** from the **Method** menu. The prepared table for the components lithium, sodium, ammonium, potassium, calcium, and magnesium is displayed below the chromatogram. For all peaks which could be allocated unambiguously to a component this table contains the peak number in the **Peak** column and the retention time in the **Time** column. In the chromatogram itself, this retention time is indicated by the cursor (vertical line).




- ⇒ If the **Peak** column contains a **0** the peak must be manually allocated to a component. In this case, insert the peak number from the chromatogram into the **Peak** column and click the **Time** field on this row with the mouse. The corresponding retention time is then inserted automatically as the new time for this component.
- ⇒ For each correctly allocated peak, click on the retention time in the **Time** column. This time is automatically shown in the chromatogram by the cursor. If the retention time is not located in the area of the middle of the peak, then the value in the **Time** column should be adapted accordingly. This is done by clicking on the chromatogram and moving the cursor with the aid of the cursor keys <←> and <→> to the middle of the peak and reading off the corresponding retention time in the status line. Enter this value (rounded off if necessary) into the **Time** column. The optimized table could, for example, then appear as follows:


	<i>Peak</i>	<i>Time</i>	<i>Wind. %</i>	<i>Ref.</i>	<i>Name</i>	<i>Group</i>	<i>Index</i>	<i>RF</i>	<i>min C</i>	<i>max C</i>
1	1	2.6	5	No	lithium	0	0	0.696929	0	0
2	2	3.1	5	No	sodium	0	0	2.67996	0	0
3	3	3.4	5	No	ammonium	0	0	2.40064	0	0
4	4	4.2	5	No	potassium	0	0	5.3861	0	0
5	5	6.3	5	No	calcium	0	0	3.17997	0	0
6	6	9.3	5	No	magnesium	0	0	1.48385	0	0

15 Start calibration

⇒ Click the **<Concentrations>** button in the component window.



⇒ Click on **<Calibrate>**. The following window appears:




⇒ Confirm **Level 1** with **<OK>**.

⇒ Click on **<OK>** in the **Concentrations** window.

16 Display calibration curves

- ⇒ Click on <Graphs> in the **Components** window.
- ⇒ Select the desired component in the **Component** field for which the calibration curve should be displayed (e.g. for **sodium**).



- ⇒ Close the **Component - sodium** window with <OK>.
- ⇒ Close the **Components** window below the chromatogram by clicking on <OK>.

17 Save chromatogram and method

- ⇒ Close the chromatogram window. A window appears with the message **Changes in *.chw. Modified: Calibration. Save changes?**
- ⇒ Click on <Yes>. A window appears with the message ***.chw already exists. Overwrite?**
- ⇒ Click on <Yes>. A window appears with the message **Method Cation1.mtw was modified. Save changes?**
- ⇒ Click on <Yes>.

3.4 Sample determination


Following the calibration of the IC system as described in *section 3.3* the first sample solution can be injected.

1 Start determination

⇒ Click on **<Start>** or select **Start determination** of the **Control** menu in the system window. An empty chromatogram window opens where the baseline is recorded continuously. The message **Initialisation [790 COM#]** appears in the **SYSTEM STATE** window. Then the **Edit sample description** window is opened automatically.

2 Enter information for determination

⇒ Enter the desired information concerning the sample in the **Edit sample description** window and confirm with **<OK>**. Set **Calibration level** to **0** and confirm with **<OK>**.



Edit sample description: 790COM1

Ident: Tap water Calibration level: 0

Info 1:

Info 2:

Volume: 10 µL Dilution: 1 Vial number: 1

Amount: 1 Internal standard amount: 100

Date/time when sample was collected (if different from injection time):
 / / : :

Ok Cancel

3 Switch injection valve to "FILL" position


⇒ Click on the **Fill** button of the 790 icon in the system window to switch the injection valve to the "FILL" position. At the same time, the suppressor module is switched to the next position.


4 Fill sample loop

⇒ Immerse the aspirating tubing **23** in the vessel containing the drinking water sample.

⇒ Using the syringe fixed to the syringe tubing **25**, siphon in ca. 1 mL drinking water.


5 Switch injection valve to "INJECT" position

- ⇒ Click on the  button of the 790 icon in the system window to switch the injection valve to the "INJECT" position. At the same time, the data recording is started automatically. The status bar of the chromatogram window shows the message **Measure**, the **SYSTEM STATE** window shows the message **INJECT done [790COM#]**.
- ⇒ When the analysis time of 14 min set in the method has expired data recording is automatically stopped and the chromatogram is integrated and evaluated. In the chromatogram window the peaks which have been found are numbered and the baseline drawn. In the **SYSTEM STATE** window the message **Finished [790COM#]** appears. It is now possible to record further samples with the same system.





- ⇒ At the end of the determination, the chromatogram window is closed and the chromatogram is saved automatically. The file name is generated automatically and contains date and time in a coded form:
 - 1st numeral: Alphabetical code for year
(e.g. j = 1999, k = 2000, l = 2001, etc.)
 - 2nd numeral: Code for month
(1...9 = Jan....Sept., a...c = Oct....Dec.)
 - 3rd+4th numeral: Day (01 ... 31)
 - 5th - 8th numeral: Time (hh:mm)

6 Open chromatogram

- ⇒ Click on  or **File / Open / Chromatogram** in the main window. Select the chromatogram file ***.chw** just recorded and click on **<OK>**. The chromatogram window is opened where the found peaks are numbered and the baselines are drawn in.


7 Report output

⇒ Click on  or select **Make report** of the **Process** menu in the main window. The **Report Options** window appears.



- ⇒ Under **Report destination**, switch on the desired target option for report output.
- ⇒ Under **Items to report** and **More items to report**, click all desired elements for the report output.
- ⇒ Click on **<Report>** to start the report output to the selected target.

8 Print report

⇒ Click on  or select **Print** of the **File** menu in the main window. The standard printing window is opened where printer, printing range and number of copies can be selected. After confirmation with **<OK>** the results including the chromatogram are printed out.

4 Operation



This section describes the most important points concerning the operation of the 790 Personal IC. For further details, please refer to the on-line help in the PC program, which can provide you with the required information rapidly and conveniently from any place in the program.

4.1 Fundamentals of the operation


4.1.1 Starting/closing the program

Start the «790 Personal IC» program



Start the program

Double-click this icon or the **Metro790.exe** file to start the «790 Personal IC 1.0» program. The Login window appears:



Enter your password and click on <Log In>.



*After software installation, the program can be started without entering a **Password**. For the definition of users, see section 4.2.2.*

Close the «790 Personal IC» program

790 PERSONAL IC / File / Exit

Exit the «790 Personal IC» program.

The program is also quit by clicking on in the upper right part of the main window **790 PERSONAL IC**.

4.1.2 Glossary

System

The term “system“ is used to describe the combination of **instrument settings, time program** and **processing method** which have been optimized for the specific separating column and the determination to be carried out with it. A system is used to start single determinations or determinations with the help of a sample queue.

Systems are stored as **system files (*.smt)** in the **Systems** directory.

Method

A method contains all information necessary for **data acquisition, integration, peak evaluation** and **quantification**. It can be considered as the chromatogram template, i.e. chromatogram without raw data.

Methods are stored as **method files (*.mtw)** in the **Methods** directory.

Each system is linked to a method. This method is called **processing method** and is opened automatically at the start of a new determination.

Chromatogram

A chromatogram is a graphic plot of the elution curve (signal vs. time) recorded following a chromatographic separation on a separating column.

Chromatograms are stored as **chromatogram files (*.chw)** in the **Data** directory. As well as the measuring data, the chromatogram files also contain the method parameters and system settings, which have been used for data recording, data processing and remote control.

Determination

In order to carry out a determination a suitable **system** must be selected for the separating problem. The result of the determination is a **chromatogram**, in which the measuring data and results of the determination are stored.

Calibration

Calibration is used to describe the method of determining the relationship between the peak height or peak area found for one component and its concentration in the sample. The result of the calibration is a **calibration function** (calibration curve), which shows the relationship between the amount of sample and the evaluated quantity.

The determination of the calibration function with reference solutions can be carried out as a **one-point** or as a **multiple-point calibration**. The calibration method that is mainly used in ion chromatography is the **external standard calibration** (absolute calibration); calibrations with an **internal standard** (relative calibration) or **tabulated calibrations** are also possible.

Integration

Integration is to be understood as being the method for determining the peak area and peak height with the aid of approximate baselines. The integration algorithm included in the program is influenced by the **integration parameters** and the optionally programmable **integration events**, which are defined in the method. In addition, the integration can be manually corrected later with the aid of the **peak editor**.

Sample queue

A sample queue is used for the automated processing of series of samples, particularly in combination with a sample changer.

Batch reprocessing

Batch reprocessing is understood to be the subsequent reprocessing of a series of chromatograms, which have been loaded in a batch reprocessing queue. During reprocessing with a selected method the settings for calibration, integration, passport, appearance and report can be altered at will.


4.1.3 Overview of program windows

The «790 Personal IC» program consists of different windows whose functionality is linked together. The different windows are:

790 PERSONAL IC	Main program window for file administration, printing, opening of systems, methods and chromatograms, login and user rights, optional settings and window handling.
CHROMATOGRAM	Window for graphic plot of running or recorded chromatograms.
SYSTEM	Window for loaded system with possibility for manual control of the 790 Personal IC.
SYSTEM STATE	Window for display of status messages for the connected system.
WATCH WINDOW	Window for display of conductivity and pressure.
QUEUE EDITOR	Window for edition of sample queues and batch reprocessing queues.















4.1.4 Main window elements






The elements of the main window **790 PERSONAL IC** are the menu bar, the toolbar and the status bar, indicating prompts and logged-in user.



4.1.5 Icons of the main window

The following icons are displayed in the **790 PERSONAL IC** main window:

-  Open chromatogram
-  Save chromatogram
-  Open last batch reprocessing file
-  Report settings
-  Print preview
-  Print report
-  Passport
-  General method settings
-  Integration parameters
-  Component table
-  Cascade all opened chromatogram windows
-  Vertical tiling of open chromatogram windows
-  Horizontal tiling of open chromatogram windows
-  Appearance of the chromatogram window

-  Enable/disable peak editor mode
-  View whole chromatogram
-  Help
-  Lock system
-  Opened system

4.1.6 Overview of file types

The following file types are produced by the «790 Personal IC» software:

- *.bar** **Batch reprocessing file**
This binary file contains data of the batch reprocessing queue. The ***.bar** file is stored automatically in the **Data** folder.
- *.cal** **Calibration file**
This binary file contains calibration data, which can be exported with **790 PERSONAL IC / Method / Calibration / Export calibration**. The ***.cal** file is stored automatically in the **Methods** folder.
- *.chw** **Chromatogram file**
This binary file contains chromatogram, system and method data of a determination. The ***.chw** file is stored automatically in the **Data** folder.
- *.mtw** **Method**
This binary file contains the data acquisition method linked to a system. The ***.mtw** file is stored automatically in the **Methods** folder.
- *.que** **Sample queue file**
This binary file contains sample table data. The ***.que** file is stored automatically in the **Methods** folder.
- *.rtt** **Report template**
This ASCII file contains a report template. The ***.rtt** file is stored in the program folder.
- *.smt** **System file**
This ASCII file contains the system settings. The ***.stm** file is stored automatically in the **Systems** folder.
- *.dev** **Device file**
This ASCII file contains drivers for devices. The ***.dev** file is stored in the **Devices** folder.

4.1.7 Context sensitive menus

Some of the menu functions of the program windows are also accessible by clicking on the desired window or item and pressing the **right mouse button**. The pop up windows have different contents and functions depending on the selected active window or item type.

4.1.8 Keyboard and mouse functions

The **mouse** can be used to carry out the normal program operating functions such as the selection of menu items and fields. It can additionally be used for magnifying a section of a chromatogram (**zooming**). To **zoom** a portion of the plot it is necessary to place the mouse cursor to the upper left corner of the square to zoom, press the left mouse button and drag the cursor to the lower right corner of the rectangle. After releasing of the left mouse button the selected region will be zoomed full-screen. If the cursor is active in the peak editor mode then it can be moved by pressing down the right-hand mouse key.

The **keyboard** can also be used to scale a chromatogram in the window, as described below.

Keyboard quick reference

Cursor is inactive:

[up]	Increases sensitivity on the Y axis.
[down]	Reduces sensitivity on the Y axis.
[right]	Expands a chromatogram on the X axis.
[left]	Shrinks a chromatogram on the X axis.
[Ctrl] + [Home]	Autoscale procedure on the X axis (shows all on X).
[Ctrl] + [End]	Autoscale procedure on the Y axis (shows all on Y).
[PageUp]	Shifts a chromatogram on $1/10$ part of a screen upwards.
[PageDown]	Shifts a chromatogram on $1/10$ part of a screen downwards.
[Shift] + [up]	Increases a distance between channels of a chromatogram.
[Shift] + [down]	Reduces a distance between channels of a chromatogram.
[0 (Zero)]	Adjusts a zero on the last point of a chromatogram (running chromatogram) or its lowest level (finished run).



Only part of the chromatogram is on screen:

[Ctrl] + [right]	Moves one window right (without change of scale on X and Y axes).
[Ctrl] + [left]	Moves one window left (without change of scale on X and Y axes).
[Home]	Shows the beginning of a chromatogram (without change of scale on X and Y).
[End]	Shows the end of a chromatogram (without change of scale on X and Y).
[0 (Zero)]	Adjusts a zero on the lowest level in the window.

Cursor is active:

[0 (Zero)]	Adjust a zero in site of the cursor.
[right]	Moves cursor left to right.
[Shift] + [right]	Quickly moves cursor left to right.
[left]	Moves cursor right to left.
[Shift] + [left]	Quickly moves cursor right to left.
[Home]	Moves cursor to beginning of a window.
[End]	Moves cursor to end of a window.
[Shift] + [End]	Sets the beginning of a window in site of the cursor.
[Shift] + [Home]	Sets the end of a window in site of the cursor.

4.1.9 Help

By clicking on  , by clicking on  , by selecting the **Help / Contents** menu item, or by pressing the [F1] key you can get on-line help on the current topic anywhere in the program.

Green texts can be clicked to jump to a different Help topic.

Violet texts identify the dialog item, parameter or button in the corresponding window.

Blue texts identify titles and important information.

4.2 Instrument and software settings

4.2.1 Fonts

790 PERSONAL IC / Options / Fonts

This option allows the selection of fonts used by the system.


Font for dialog...	Selection of font used for dialog boxes. Default setting: MS Sans Serif / Standard / 8 pt.
Font for reports...	Selection of font used for report output to the screen or printer. Default setting: Courier New / Standard / 10 pt.
Font for tables...	Selection of font used for data presentation in tables on the screen. Default setting: MS Sans Serif / Bold / 8 pt.
Font for plots...	Selection of font used for labels on chromatogram plots and calibration curves. Default setting: Times New Roman / Bold / 10 pt.
Save fonts configuration	Save chosen font configuration.

4.2.2 Security system

The «790 Personal IC» program has a security system based on the list of users. Every user has his unique password and one of the following access levels:

Novice	Restricted access to program functions. Allows only start and stop of determinations using existing system and method files and manual control of the 790 Personal IC. Modifications of system, method and data files are not allowed.
Master	Access to all program functions with few exceptions: the user cannot set Global preferences , Hardware settings and security system.
Administrator	Access to all program functions.

It is recommended to make user lists and enter passwords as a first action after system installation. Therefore, select the menu item **790 PERSONAL IC / Options / Security** and click on **<Log In>** in the Log In window without entering a password. The **Security options** window appears:




Enter user name, password and access level for all users. Don't forget to make one of the users an **Administrator**, otherwise this window will not open again. At the end, click on **<OK>**.

After configuration of the security system the program prompts for the password every time the system starts. This user name stamps all methods, chromatograms and reports. It is possible at any time to change the user with the menu item **790 PERSONAL IC / Options / Lock system**.

4.2.3 Global settings

790 PERSONAL IC / Options / Global preferences

This window is used for **global program settings**.



*This window is only accessible for users with **Administrator** access level.*


Data file overwrite	Overwriting of chromatogram files:
Never	Chromatogram files cannot be overwritten. A modified chromatogram is saved as a new file with the file name number raised by 1.
Ask	The user is asked if the chromatogram should be overwritten.
Always	Chromatogram files are always overwritten without confirmation.
<hr/>	
If method changed	Saving of method files:
Don't ask to save	The method file is not saved automatically. It can be saved only with File / Save / Method .
Ask to save	The user is asked if the method should be saved.
Overwrite without asking	Method files are overwritten without confirmation.
<hr/>	
If disk method is newer	Overwriting of method files:
Ask to overwrite	The user is asked if the method should be overwritten.
Overwrite without asking	Method files are overwritten without confirmation.
<hr/>	
Opening chromatograms	Opening of chromatograms:
Ignore last data directory	If this option is enabled, the default data directory Data is opened.
<hr/>	
GLP On	If this option is enabled, the following parameters are automatically set: Data file overwrite = Never If method changed = Don't ask to save If disk method is newer = Ask to overwrite
<hr/>	
OK, Apply buttons mean	
Save object to file	If this option is enabled, the <Apply> button in system settings windows is replaced by the <Save> button. The system settings are saved if the <Save> or <OK> button is clicked.
<hr/>	
Chromatogram units	Units for chromatograms:
Flow	Unit for flow rate: µL/min, mL/min
Pressure	Unit for pressure: MPa, psi, bar, atm
<hr/>	
Print via print spooler	Switch on/off printing via print spooler. Switch off this option if you use a GDI printer.

Measuring priority	Setting the priority of program execution. real-time means that the «790 Personal IC» program has the highest priority, normal means that all active programs have the same priority.
<Default colors>	Default colors for chromatographic windows (details see <i>section 4.5.3</i>).

4.2.4 COM port


790 PERSONAL IC / Options / 790 Personal IC:COM1

This menu item opens the **Links** window for COM port (serial RS232 interface) selection and settings.



If **COM1** is clicked with the right mouse button, the following menu items appear:

- Put on desktop** Possibility for setting COM port parameters for recording data transfer (details see on-line help).
- Change** Possibility for changing the COM port (default setting: **COM1**). The following window is opened, where the interface can be changed by clicking on the desired COM port.




4.3 Systems

The term “System“ describes the combination of **instrument settings**, **time program** and **processing method** which has been optimized for a specific separating column and the determination to be carried out with it. A system is used to start single determinations or determinations with the help of a sample queue. Systems are stored as **system files** (*.smt) in the **Systems** directory.

4.3.1 System window

A system window is opened with **790 PERSONAL IC / File / Open / System** and the selection of the desired system file. It contains icons for **PC**, **Watch window** (screen) and **790 Personal IC**, and two buttons to **start** and **stop** a determination.



4.3.2 System file handling

The following menu items are used for opening, changing, saving and closing of systems:

790 PERSONAL IC / File / Open / System

Load an existing system file (*.smt) from the **Systems** directory and open the corresponding system window.

The directory and name of the opened system file are displayed in the title bar of the system window. A star (*) at the end of the name indicates that the system settings have been changed since the last saving.

SYSTEM / System / Change

Load a new system and open automatically the system window.

SYSTEM / System / Save

Save the current system settings in a system file (*.smt) in the **Systems** directory.

SYSTEM / System / Close


Close the selected system.

4.3.3 System functions

Start/stop hardware and record baseline

SYSTEM / Control / Startup hardware (Measure Baseline)

Starting the hardware at the 790 Personal IC includes sending of **System startup values**, starting of the high-pressure pump and (if present) starting of the peristaltic pump.

At the same time, the **recording of the measurement signal** using the method defined for the system is started. Independently of the set chromatogram **Duration**, the measurement signal is recorded until a new determination is started. Alternatively the baseline recording can be stopped by clicking the  icon of the chromatogram window. In this case, the user is asked if the recorded baseline should be saved or not.

SYSTEM / Control / Shutdown hardware

High-pressure pump and (if present) peristaltic pump at the 790 Personal IC are immediately stopped. A running determination and an active sample queue are also stopped.

Start/stop determinations

SYSTEM / Control / Start determination

Start determination using the settings of the selected system. At this start command, the **System startup values** are set at the 790 Personal IC. The high-pressure pump and the peristaltic pump are started if they are not already running. The time program is started immediately; the data recording is started after switching the injection valve to the "Inject" position.

At each start of a determination, the **Edit sample description** window is opened automatically. This window can be used for entering information concerning the sample:

Edit sample description: 790COM1

Ident: Tap water Calibration level: 0

Info 1:

Info 2:

Volume: 10 μL Dilution: 1 Vial number: 1

Amount: 1 Internal standard amount: 100

Date/time when sample was collected (if different from injection time):

/ / : :


Ok Cancel

Ident	User defined identifier (title) for the chromatogram to be displayed in the title bar of the chromatogram window and in the Chromatogram open window.
Calibration level	Calibration level (0 = sample; 1...n = calibration solutions).
Info 1 / Info 2	Sample description.
Volume	Injected volume in μL .
Dilution	Dilution of the sample.
Vial number	Autosampler vial position to take sample from.
Amount	Sample amount. If this value is different for the calibration run (c) and the sample run (s), the component concentrations of the sample are calculated as follows: $C_s = C_c \cdot \text{Amount}_s / \text{Amount}_c$
Internal standard amount	Concentration of the internal standard component for relative concentration calculations.
Date/time when...	Date and time of sample collection (the default values are equal to the date and time when the chromatogram starts).

The **Edit sample description** window is not opened if determinations are made with an active sample queue.

SYSTEM / Control / Stop determination

Stop running determination. Data acquisition and time program are terminated immediately. The recorded chromatogram is saved automatically if the **Save chromatogram after the run** option on the **Passport / Processing** tab is enabled.


Alternatively the determination can be stopped by clicking the  icon of the chromatogram window. In this case, the user is asked always if the determination should be saved or not.

Print system parameters

SYSTEM / System / Print parameter

A report of the system parameters including the time program is created and opened using the «Wordpad» program. The ***.txt** file opened can be printed, saved and exported into other programs. The system report includes the name of the method linked to the system, the configuration settings for the 790 Personal IC, the system startup values and the time program (incl. **Remote configuration**), if such a program exists and is switched on (**ENABLED**).


4.3.4 PC icon



The PC icon is one of the components of the system window. If the system is connected and this icon is clicked with the right mouse button, the following menu item appears:


- Open** Load the processing method linked to the system and open an empty chromatogram window.

4.3.5 Watch window



The watch window icon is one of the components of the system window. If the system is connected and this icon is clicked with the right mouse button, the following menu item appears:

- Open** Open the **WATCH WINDOW** for live display of conductivity and pressure:



The colors of the watch window fields can be changed by clicking the fields with the right mouse button and selecting the menu item **Choose color / ...**.

4.3.6 Instrument icon

Menu options for instrument icon



System disconnected



System connected
Injection valve in
"INJECT" position



System connected
Injection valve in
"FILL" position

The instrument icon for the 790 Personal IC is one of the components of the **SYSTEM** window. If the system is connected, the 790 icon contains two buttons for manual control of the injection valve:


- <Inject> Switch injection valve to "INJECT" position.
- <Fill> Switch injection valve to "FILL" position.

If the system is connected and the 790 icon is clicked with the right mouse button, the following menu appears:

- Open** Open the **system settings** window.
- Hardware** Open the **hardware settings** window.
- Diagnostics** Open the **diagnostics** window.

System parameters for disconnected system

A double-click on the instrument icon or selection of the **Open** menu item using the right mouse button opens the system settings window. For a disconnected system, the **Initial** tab is displayed.




Flow, mL/min

Startup value for flow rate of the high-pressure pump.
Entry range: **0.20 ... 2.50 mL/min**

Pressure max, MPa	Startup value for maximum pressure limit for high-pressure pump. Entry range: 0.0 ... 25.0 MPa
Pressure min, MPa	Startup value for minimum pressure limit for high-pressure pump. Entry range: 0.0 ... 25.0 MPa

Instrument control for connected system

A double-click on the instrument icon or selection of the **Open** menu item using the right mouse button opens the system settings window. For a connected system, the **Control** tab is displayed. It allows manual control of the 790 Personal IC functions and setting of startup values to be sent to the instrument. This tab shows the current measurement values for conductivity and pressure.




Conductivity, µS/cm	Live display of measured conductivity.
Pressure, MPa	Live display of measured pressure. The color settings for this two fields can be changed by clicking the fields with the right mouse button and selecting the menu item Choose color / ...

Actual values	Display of actual values
Flow, mL/min	Display of flow rate of the high-pressure pump.
Pressure max, MPa	Display of maximum pressure limit for high-pressure pump.
Pressure min, MPa	Display of minimum pressure limit for high-pressure pump.

Valve	Injection valve
<Inject>	Switch injection valve to "INJECT" position.
<Fill>	Switch injection valve to "FILL" position
Suppressor	Suppressor module
<Step>	Switch the suppressor module to the next position. The time since the last switching of the suppressor module is displayed in the field beside the <Step> button.
IC pump	High-pressure pump
<On>	Start high-pressure pump.
<Off>	Stop high-pressure pump.
Peristaltic	Peristaltic pump
<On>	Start peristaltic pump.
<Off>	Stop peristaltic pump.
System startup values	The system startup values are sent and applied to the 790 Personal IC each time the system is connected, a determination is started, or the values are sent manually with <Send to unit>.
Flow, mL/min	Startup value for flow rate of the high-pressure pump. Entry range: 0.20 ... 2.50 mL/min
Pressure max, MPa	Startup value for maximum pressure limit for high-pressure pump. This limit is controlled even without connection to the PC. Entry range: 0.0 ... 25.0 MPa
Pressure min, MPa	Startup value for minimum pressure limit for high-pressure pump. This limit is controlled even without connection to the PC. Entry range: 0.0 ... 25.0 MPa

Time program

On the **Program** tab of the system settings window a user-defined time program for instrument control can be entered. This program is started automatically at the moment the determination is started. The time program contains program steps including time, program instruction and parameter.



Time (1st column)	Time at which program instruction is applied. Entry range: 0.0 ... 999.9 min If no time is entered, the program instruction is applied together with the last instruction with time entry.
Command (2nd column)	Program instruction (see <i>List of program instructions</i>).
Parameter (3rd column)	Parameter for program instruction (see <i>List of program instructions</i>).
ENABLED	Enable program start (a disabled program is not started).
<Add>	Add new program instruction.
<Delete>	Delete selected program instruction.
<Verify>	Test the time program (error messages are displayed if program is faulty).

List of program instructions


The following program instructions can be added to the time program on the **Program** subpage:

Instruction	Parameter entry	Meaning
Valve	Inject, Fill	Switch injection valve to "INJECT" or "FILL" position.
ICPump	on, off	Switch on or off the high-pressure pump .
Flow	0.2 ... 2.5 mL/min	Set flow rate of the high-pressure pump to the desired value.
Pmax	0.0 ... 25.0 MPa	Set maximum pressure limit for the high-pressure pump to the desired value.
Pmin	0.0 ... 25.0 MPa	Set minimum pressure limit for the high-pressure pump to the desired value.
Remote		Pulse output (pulse length = 200 ms) on remote output line 3 for sample changer control.
Program	END, RESET	The END flag can be used to end a program, especially if the program time should be longer than the chromatogram duration. Additional steps after this flag are not allowed. The RESET flag is used to reset the parameters to the system startup values.

Suppressor	Switch suppressor module to the next position.
Peristaltic on, off	Switch on or off the peristaltic pump .

Configuration

The **Configuration** tab in the system settings window contains configuration settings for the 790 Personal IC.



790 unit version	Selection of instrument version: 790.0010 790 Personal IC without suppressor 790.0020 790 Personal IC with suppressor
-------------------------	---

Suppressor Autostep with	Suppressor module: Automatic switching to the next position each time the injection valve is switched either to the Fill or to the Inject position.
---	--

Use peristaltic pump	If this option is disabled, the peristaltic pump is not switched on automatically at the start of a determination or with Startup hardware .
-----------------------------	---


Polarity	Selection of the polarity of the output signal: + positive polarity (for anions) - negative polarity (for cations)
-----------------	--

Hardware settings

Selection of the **Hardware** menu item using the right mouse button opens the hardware settings window consisting of the two tabs **Hardware** and **Alarm stops**.

Hardware

This tab of the hardware settings window defines general parameters which are set automatically at power on of the instrument.



Flow correction

Factor for correction of the difference between displayed and actual flow rate of the high-pressure pump.

Range: **0.9 ... 1.09**

The flow correction is determined by measurement of the actual flow rate with the aid of a measuring cylinder as follows:

$$\text{Flow correction} = \frac{\text{Displayed flow rate}}{\text{Measured flow rate}}$$

Cell constant

Cell constant of the conductivity cell for correct display of the absolute conductivity. Enter the cell constant printed on the detector block into this field.

Range: **0.1 ... 1000 /cm**

For a precise determination of the cell constant, pump a calibration solution of known conductivity through the IC system, observe the displayed conductivity and change the cell constant until the correct conductivity value is displayed.

<Break times>


Possibility for changing the break times for injection valve **Valve** and suppressor module **Suppressor**.



These values should only be changed after consulting the Metrohm Service.


Alarm stops

The **Alarm stops** tab of the hardware settings window defines the events for which the 790 Personal IC is stopped immediately. At an alarm stop, high-pressure pump and peristaltic pump are stopped immediately; the running determination and the active sample queue are also stopped.



Event	Events for alarm stop:
Pump stop	Pump stopped because pressure limits are exceeded.
Hardware error	Hardware error detected at the 790 Personal IC (high-pressure pump, injection valve, or suppressor not working correctly).
Power on	Temporary power failure at the 790 Personal IC hardware.

4.3.7 System state window



The **SYSTEM STATE** window is automatically opened if a system is opened. It shows status and error messages for this system. Messages concerning the 790 Personal IC hardware are followed by **[790 Personal IC]** or **[790 Personal IC 2]**, messages concerning the loaded system are followed by **[790COM#]** (name of subfolder containing the system file).

Status messages

Checking on-line	Checking connection between PC and 790 Personal IC.
On-line	Connection between PC and 790 Personal IC ok.
UploadStartupValues	System startup values have been loaded to the 790 Personal IC.
Initialization	Hardware or system initialization.
Ready	System is ready for starting a new determination.
Starting	Starting program or chromatogram data acquisition.
Running	Running program or chromatogram data acquisition.
Running program (xxx min left)	Running time program (in brackets: time left).
Waiting for INJECT	Waiting for "INJECT" to start program and/or chromatogram data acquisition.

INJECT done	Injection valve has been switched to the "INJECT" position.
Finished	Determination has been finished.
Executing *.smt [*.que]	The sample queue *.que is processed using the system file *.smt.
SHUTDOWN	790 Personal IC hardware is shutdown.

Error messages

Detection of hardware failed	Bad connection between PC and 790 Personal IC or instrument switched off (check connecting cable or switch on instrument).
E1	Program checksum wrong (call Metrohm service).
E2	RAM faulty (call Metrohm service).
E200	Invalid instrument adjustment (call Metrohm service).
E237	Storage of configuration values failed (repeat last action; if error reappears, call Metrohm service).
E238	Storage of instrument number failed (repeat last action; if error reappears, call Metrohm service).
E240	EEPROM faulty (call Metrohm service).
E258	Storage of setup values failed (repeat last action; if error reappears, call Metrohm service).
E295	Storage of memory values failed (repeat last action; if error reappears, call Metrohm service).
E296	Instrument stopped (restart instrument; if error reappears, call Metrohm service).
E298	Storage of flow correction value failed (repeat last action; if error reappears, call Metrohm service).
E299	Storage of break time values failed (repeat last action; if error reappears, call Metrohm service).
E300	High-pressure pump faulty (restart pump; if error reappears, call Metrohm service).
E301	Injection valve blocked (check injection valve; if error reappears, call Metrohm service).
E302	Suppressor module blocked (check suppressor; if error reappears, call Metrohm service).
E303	Storage of maintenance information failed (repeat last action; if error reappears, call Metrohm service).

4.4 Methods

A method contains all information necessary for **data acquisition, integration, peak evaluation** and **quantification**. It can be considered as the chromatogram template, i.e. chromatogram without raw data.

Six prepared methods are stored as **method files (*.mtw)** in the **Methods** directory. These methods cannot be renamed or copied, but only changed and saved under the same name.

Each system is linked to the method having the same name as the system. This method is called **processing method** and is opened automatically at the start of a new determination.

4.4.1 Method file handling

The following menu items are used for opening and saving of methods:

SYSTEM / PC icon / Open

By clicking this menu item or by double-clicking the PC icon in the system window, the method file (***.mtw**) belonging to the system is loaded from the **Methods** directory. At the same time, an empty chromatogram window is opened.

The name of the method is displayed in the title bar of the system window. A star (*****) at the end of the name indicates that the method has been changed since the last saving.

790 PERSONAL IC / File / Save / Method

Save the method of the current chromatogram in a method file (***.mtw**) in the **Methods** directory.

4.4.2 Passport



790 PERSONAL IC / Method / Passport

This menu item opens the **Passport** window that includes a detailed textual description of the chromatographic run and consists of the following tabs:

General	General information on determination.
Sample	Sample information.
Column	Column information.
Eluent	Eluent information.
Comment	User comments on chromatogram.
Method Log	GLP log with date/time stamp for method.
Data Log	GLP log with date/time stamp for chromatogram.

General

Tab **General** of the passport with general description of the method and determination.



Ident	User defined identifier (title) for the chromatogram to be displayed in the title bar of the chromatogram window and in the Chromatogram open window.
Duration	Duration of the chromatogram in minutes. You can modify this value when the chromatogram is running.
METHOD	Path and file name of the method used for data acquisition (read-only).
DATA	Path and file name of the current chromatogram (read-only).
Date/time	Date and time of the chromatogram start (read-only).
Last update	Date and time of the last chromatogram modification and saving (read-only).
Calibration level	If the current run is used for calibration, the calibration level is a positive integer, otherwise it equals zero. You can modify this value when the chromatogram is running.
User	Name of the current user. It is taken from the list of users, in accordance with the password entered at the program starting.
Detector	Name of the detector (read-only).
Run	Number of the current run starting from the very first one. All runs are automatically numbered by the system.



To reset the **Run** number press [Ctrl] [F8]. The **Analysis number** window is opened where the run number for the next run can be modified to the desired value.

Batch

X/Y: total number of started sample queues / current analysis number in the sample queue (read-only).

Sample

Tab **Sample** of the passport with sample information.



Info 1	First sample description (max. 256 characters).
Info 2	Second sample description (max. 256 characters).
Volume	Injected volume in µL.
Dilution	Dilution of the sample.
Vial number	Autosampler vial position to take sample from.
Amount	Sample amount. If this value is different for the calibration run (c) and the sample run (s), the component concentrations of the sample are calculated as follows: $C_s = C_c \cdot \text{Amount}_s / \text{Amount}_c$
Internal standard amount	Concentration of the internal standard component for relative concentration calculations.
Date/time when...	Date and time of sample collection (the default values are equal to the date and time when the chromatogram starts).

Column

Tab **Column** of passport with column information.

The screenshot shows a software dialog box titled "Eigenschaften von Passport" with a "Column" tab selected. The fields are as follows:

- Number: []
- ID: 4 mm
- Length: 125 mm
- Packing material: 6.1010.000 Metrosep Cation 1-2
- Particle size: 7 µm
- Void volume: 0 %
- Precolumn (set length = 0 if absent): ID: 2 mm, Length: 0 mm

Buttons at the bottom: OK, Abbrechen, Übernehmen, Hilfe.

Number Serial number of column (max. 256 characters).
ID Internal diameter of the column in mm.
Length Length of the column in mm. This parameter is used to calculate **Linear flow**.

Packing material

Column Column description (max. 256 characters).
Particle size Particle size of the column in µm. This value is used for the calculation of reduced theoretical plate height **Reduced TP height**.
Void volume Void volume for the column in %. Used to calculate **Logarithmic index**, **Capacity factor** and **Linear flow**. It is calculated by the system in accordance with the settings defined on **Method setup / Math**.

Precolumn

ID Internal diameter of the precolumn in mm.
Length Length of the precolumn in mm (set length to 0 if no precolumn is used).

Eluent

Tab **Eluent** of passport with eluent information.



Eluent A	Description of eluent composition (max. 256 characters).
Eluent B	Description of eluent composition (max. 256 characters).
Eluent C	Description of eluent composition (max. 256 characters).

Flow	Flow rate of the high-pressure pump in mL/min. The flow rate is used to recalculate the time axis into volumetric units. The system startup value for the flow rate is entered automatically into this field at the start of a determination.
Pressure	Pressure in MPa. The measured value for the pressure is entered automatically into this field at the start of a determination.
Temp.	Temperature in °C. Here you can enter the temperature of the column's thermostat or the ambient temperature.

Comment

Tab **Comment** of passport for entry of free-text user comments into the chromatogram description. Use this feature to enter any additional information about the chromatogram not included in other sections of the method.



Method Log

Tab **Method Log** of passport with GLP Log for the method. The **Method Log** displays all automatically produced GLP messages concerning modifications made on a method. In addition to these messages, free-text user comments can be entered. All this entries are stamped automatically with date and time of the entry.



<Add new>

Add a new GLP comment. Text can be entered in the **New comment** window, date/time and logged-in user are automatically added.

<Clear>

Clear all GLP messages and add a new GLP comment. Text can be entered in the **New comment** window, date/time and logged-in user are automatically added.

Data Log

Tab **Data Log** of passport with GLP Log for the recorded chromatogram. Here a GLP message is automatically produced for every modification made on a recorded chromatogram. In addition to these messages, free-text user comments can be entered. All this entries are stamped automatically with date and time of the entry and cannot be cleared.



<Add new>

Add a new GLP comment. Text can be entered in the **New comment** window, date/time and logged-in user are automatically added.

4.4.3 Method setup



790 PERSONAL IC / Method / Method setup

This menu item opens the **Method setup** window that includes the most common parameters for data acquisition and evaluation of the method and consists of the following tabs:

General	General information on determination.
Processing	Set actions that are performed when the chromatogram finishes.
Math	Parameters that are used for various types of calculations.

General

Tab **General** of **Method setup** window with general description of the method and determination. This tab is identical to the **General** tab called with **Method / Passport** (see section 4.4.2).

Processing

Tab **Processing** of **Method setup** window for definition of actions that are performed during the run or when the chromatogram finishes. Automatic post-run data processing is especially useful for chromatographic systems that include an autosampler.



Save chromatogram after the run

If checked, autosaving of the chromatogram will be performed on data acquisition ending.

while running every ... min

Save the chromatogram during the run for data safety at the time interval set (if 0 is set, the chromatogram is not saved during the run).


Show all

If checked, scales on X and Y axes will be automatically chosen so that the recorded chromatogram fits the window.

Make report	If checked, report auto-printing will be performed on data acquisition ending.
Close window	If checked, the window is closed automatically after the run (or sample queue) is finished.
Chromatogram directory	Directory where the current chromatogram is saved. Use the <Browse> button to select a new directory.

Math

Tab **Math** of **Method setup** window with parameters that are used for various types of calculations.



Parameter	Selection of calculation parameter :
Formula set	Formula setting by customer (Custom formulae) or according either to European Pharmacopoeia or US Pharmacopoeia .
Theoretical plates	Selection of calculation formula for number of theoretical plates if Formula set = Custom formulae is set.
Resolution	Selection of calculation formula for resolution if Formula set = Custom formulae is set.
Asymmetry	Selection of calculation formula for asymmetry if Formula set = Custom formulae is set.
Formula	Selection of calculation formula for the calculation parameter selected for Parameter . With Parameter = Formula set the following settings are possible:
Custom formulae	Formula setting by customer.

European Pharmacopoeia	Formula setting according to European Pharmacopoeia.
US Pharmacopoeia	Formula setting according to US Pharmacopoeia.
<hr/>	
Void time/volume	Calculation of void time/void volume :
Calculation method	Selection of the method for void time calculation:
None	Manual entry in the Void time field.
1st component	The peak corresponding to the first component is selected as a void time marker and its retention time replaces the previous value of the void time. If the first component is not identified, the expected retention time is used.
1st peak	The first detected peak is used as a void time marker for the run and its retention time is stored in the Void time field.
From void volume %	For the calculation of the void time, the % value entered in the Void volume field is multiplied by the ratio of empty column volume and eluent flow rate.
Void volume	Void volume in mL and % of column volume. The % value can be entered manually if Calculation method = From void volume % is set.
Void time	Void time in s. This value can be entered manually if Calculation method = None is set.
<hr/>	
Index	Calculation of Retention index :
Interpolation	Use of linear or logarithmic interpolation scale for retention index calculation.
Internal	The index scale is constructed on the basis of the current chromatogram. All components used for index scale calibration should be present in the current sample.
External	The index scale is constructed based on another standard chromatogram.

4.4.4 Integration



790 PERSONAL IC / Method / Integration


This menu item opens the **Integration parameters** window that contains parameters and events for peak integration and consists of the following two tabs:

Setup	Integration parameters.
Events	Integration events.

The «790 Personal IC» software includes a built-in automatic integration algorithm to detect peaks on the chromatographic curve and to evaluate them using calculated baselines. The integration procedure is tuned by the integration parameters on the **Setup** tab and the integration events on the **Events** tab. Integration events have higher priority than integration parameters. The integration can be modified manually using the peak editor (see section 4.5.4).

The built-in integrator algorithm is based on the use of the first derivative (slope). In order to decide whether the slope at some point is significant, the first derivative value is divided by the baseline noise and the result is compared with a threshold value called “**Slope**”. The threshold values for the upslope and down slope differ.

The baseline and the start and end points are determined for each peak. For separated peaks, the baseline consists of a straight line between the start point and the end point. For overlapping neighboring peaks, a common baseline is normally constructed; this connects the start point of the first peak with the end point of the last peak. The peaks are separated from one another by dropping a perpendicular line from the lowest point (valley) between the peaks to the common baseline. Alternatively, the baseline can also be drawn directly between the valleys of the peaks with the event **Enable valley-to-valley** (see integration events for baseline).




Setup

Tab **Setup** of **Integration parameters** window with integration parameters. The two most important of these parameters for peak recognition are **Width** and **Slope**.



If the signal-to-noise ratio of the chromatogram is high enough (e.g. 100:1 or higher), the integration algorithm is not too sensitive to changes in integration parameters. On the contrary, if the baseline noise is high, careful tuning of these parameters may be needed.




Number of peaks	Number of peaks detected (read-only).
Delay	Time delay before starting peak detection (in min). Entry range: 0...1440 min
Width	Peak width (at the baseline, in s). This parameter is used for setting baseline start and end points. If this value is too low, excessive narrow peaks are detected from the baseline noise. If it is too large, several adjacent peaks or baseline drift can be integrated as a single peak. For an optimal evaluation, it is recommended to enter the width of the narrowest peak on the chromatogram (normally 2...10 s). Entry range: 0.1...480 s
Slope	Threshold for peak recognition. The value of the first derivation (slope) of the curve is divided by the baseline noise (which is estimated using a special algorithm) and the result is compared with the Slope threshold value. Reasonable range of Slope parameter is 0.5...25. Entry range: 0.1...400

Min height	Minimum height of peaks to be detected.
Negative peaks	If this option is enabled, negative peaks are also detected.

Events

Tab **Events** of **Integration parameters** window for definition of integration events. Integration events have a higher priority than integration parameters and are used for fine-tuning of the integration process. They should be used only in the case when problems cannot be solved by tuning a set of parameters from the integration parameters window.



Number of events Number of integration events (read-only).


<Add> Add event to the list.

<Modify> Edit the selected event.

<Delete> Delete the selected event.

<Apply> Accept integration parameters and perform re-integration.


If an integration event is added or modified, the **Edit integration event** window is opened where the following parameters can be set:



Time	Start time for integration event in min.
Event	List box with integration events to be selected (see below).

Integration events

Disable detection	Set this mode to stop detection of new peaks. If a peak is available at the moment of the event it is either terminated (down slope peaks) or rejected (upslope peaks).
Enable detection	Clear Disable detection mode.
Set peak start	Force the beginning of a new peak. If a peak is available at the moment of the event it is either rejected (upslope) or terminated (down slope).
Set peak end	Force peak end at the time of the event. Upslope peaks are rejected (except those born by Set peak start event), down slope peaks are terminated.



Enable valley-to-valley Disable valley-to-valley

Enable valley-to-valley Disable perpendicular drop peak separation. All peaks are considered to be baseline-separated. The bottom of the valley becomes the baseline point.


Disable valley-to-valley Enable perpendicular drop peak separation (default setting).

4.4.5 Calibration and quantification

General information

The aim of any chromatographic analysis is to answer the question "What components are present in the sample and what are their concentrations?". Two procedures are used to achieve this goal: the first step is called **calibration**, the second step includes **quantification**.

Calibration has two aims: to get retention characteristics for all components of interest (these data are stored in the component table) and to establish a relation between injected amounts and corresponding instrumental responses for all components of interest (stored in the concentration table). Calibration is performed by running one or several chromatograms of samples with known composition and known concentration of components (standards). For each calibrated component a **calibration curve** is constructed as a result.



With the 790 Personal IC three different procedures can be used for the construction of the calibration curve. By far the most important method for ion chromatography is the **external standard calibration** (absolute calibration) that is described in detail in this section. The other methods of **internal standard calibration** (relative calibration) and **tabulated calibration** (relative gradient factor, a modified method for external standard calibration) are of lesser importance and are not described in detail here (for details please refer to on-line help).

Identification is a procedure that enables to decide what peaks on the chromatogram correspond to what components. The identification is performed based on the **Component table** created for calibration.

Quantification is a calculation procedure that determines components concentrations, based on instrumental response (peak height or area), using the calibration curves obtained earlier for each component.

Notations

R	Stands for response value, either area or height , depending on setting selected in the Calibration graphs window.
V	Sample Volume injected.
D	Dilution coefficient, shows number of times to which the initial solution is dissolved before injection.
$V' = V / D$	Adjusted volume of injected sample. A correction is made for the dilution coefficient.
C	Concentration of the component in the initial solution (before dilution).
$Q = C \cdot V'$	Quantity of component, used for calibration curve construction.
t	Retention time . Time needed by the mobile phase to flow through the separation system.
t_0	Void time . Dead time needed by the mobile phase to flow through the separation system.
$t' = t - t_0$	Corrected retention time , called also net retention time.
L	Column length.
$v = L / t_0$	Linear Flow rate .
$W(R) = k_2R^2 + k_1R + k_0$	Calibration function (component quantity W vs. detector response R). In the case of the most common linear calibration curve $Q = W(R) = k_1R$ it comes through the origin. The concentration of the component in the analyzed mixture is calculated by the formula $C = W(R) / V'$.
RSD(Q, R)	Procedure, used for computation of regression coefficients (k_0 , k_1 and k_2) of the calibration function W(R) using RSD (Residual Standard Deviation). The procedure gets input as a set of calibration points (quantity Q vs. response R) and outputs the calibration function W(R) used for prediction of the component quantity $Q_i = W(R_i)$.

Subscript values used:

j	Stands for j-th calibration level run.
s	Stands for standard component.
i	Stands for component number.

External standard calibration

External standard calibration (absolute calibration) is a key way of calibration. Basically this procedure calculates the dependence of injected component's quantity versus area (or height) of the corresponding peak. This dependence is showed as a calibration plot with injected quantity along the Y axis and peak area (or height) along the X axis. The injected quantity Q_i is calculated as product of component's concentration C_i to the reduced injected sample volume V' :

$$Q_i = C_i \cdot V'$$

The calibration function $W_i(R)$ for each component is calculated using the RSD method.

$$W_i(R) = RSD(Q_{ij}, R_{ij}) = RSD((C_i \cdot V'), R_{ij})$$

When quantification procedure is performed, a concentration C_i (absolute concentration of a component in the sample) is determined as a ratio of the computed quantity $W_i(R_i)$ to the corrected volume of sample injected V' :


$$C_i = W_i(R_i) / V'$$

Component table



790 PERSONAL IC / Method / Calibration / Components

If the component table is opened, the chromatogram window is split into two parts. The upper one shows the chromatogram, the component table appears in the bottom part of the screen. When moving into the component table, a special cursor in the upper part of the window jumps to the peak corresponding to the current component.



The **Components table** stores the following information on the components to be analyzed:

Number	Row number (read-only).
Peak	Number of the peak in the chromatogram that corresponds to the given component. If this field contains a 0 , the peak number must be entered manually.
Time	Expected retention time of the component in the calibration run. After a new entry of the peak number the corresponding retention time is entered automatically as soon as this field is clicked on with the mouse.
Wind.%	Identification window of the component. This value is the maximum allowed difference of the actual retention time of the component and its expected retention time, measured as % of expected retention time. The component will be identified within its identification window only.
Ref.	Reference component. The entry Yes defines components as reference components; these ensure improved peak identification of the other Ordinary components (details see on-line help).
Name	Component's name (cannot include spaces). Rows with empty names are excluded from the table of components when you leave the component table editor.
Group	Number of the group to which the components are allocated. If group numbers ≥ 1 are entered here then a group report can be outputted for each group; this also contains the intermediate totals for this group. This group report is produced after the main report table.
Index	Retention index for components with known index. In the case of an unknown index this value should be equal to 0 . For index calculation the user must define index values for at least two peaks, and all other values will be calculated by the software using linear or logarithmic approximation (details see on-line help).
RF	Response factor. This value corresponds to the slope of the calibration curve (coefficient k₁ of the calibration formula).
min C (max C)	Minimum (Maximum) concentration value for the component. Components whose concentration leaves the range min C...max C are marked in the peak table by the sign " ! ".

Additionally, the **Components table** contains the following buttons:

<Add>	Add a new component (an empty column) to the component table.
	Delete the current component from the component table.
<Identification>	Open the Peak identification window (see <i>Identification</i>).
<Concentrations>	Open the Concentration table (see <i>Concentration table</i>).
<Graphs>	Open the component window with display of calibration curve and calibration parameters.
<OK>	Accept all changes and close the Component table .
<Cancel>	Reject all changes and close the Component table .

The standard methods supplied with the «790 Personal IC» program all contain a component table with pre-defined components with entries for the names and retention times.




*If the output of all peaks is required in the result report then a so-called universal component can be defined in the component table for this purpose; this receives the value **0** in the **Time** field and no peak is allocated to it. The name entered under **Name** (e.g. **unknown**) is then used for all peaks for which no component can be allocated.*

Peak identification

790 PERSONAL IC / Method / Calibration / Identification

Selecting this menu item or clicking on <Identification> in the component table opens the **Peak identification** window with parameters for tuning the peak identification procedure.

Number of components	Number of components in the component table (read-only).
Scheme	Quick choice of identification parameters scheme.
Standard	Set default meanings of identification parameters: Height for reference components, Time for all other components (ordinary components).
Nonstandard	Set custom meanings of identification parameters.




Identification	Identification parameters.
Reference peaks	Identification parameter for reference components. Default is Height .
Other peaks	Identification parameter for ordinary components. Default is Time .
Retention units	Choice of retention units. Default is min .
Retention time	Retention time.
<Update>	Updates expected retention times of components according to the current chromatogram.
Worst case...	Information on the component with the worst (largest) deviation of actual and expected retention time. Deviation is given as a part of component's identification window.
Average relative deviation	Relative deviation averaged for all components.

Concentration table

790 PERSONAL IC / Method / Calibration / Concentrations

Selecting this menu item or clicking on **<Concentrations>** in the component table opens the **Concentrations** window with the concentration table containing concentrations of all components for all **Calibration levels**. Each calibration level corresponds to a sample used for calibration and to a point on the calibration curve.



The screenshot shows the 'Concentrations' window with the following data:

	<i>Name</i>	<i>This run</i>	<i>Level 1</i>
1	fluoride	0.00409721	0.5
2	chloride	0.0283832	5
3	nitrate	0.0298117	10
4	sulphate	0.039948	10

Below the table, there are buttons for 'OK', 'Cancel', 'Add', 'Delete', 'Calibrate', and 'Help'. A 'Levels' group box is also visible.

Concentration units User defined units of concentration that will appear in the report. Nevertheless, entry of new concentration units does not cause recalculation of concentrations.

Data type Choice of the data type that will be shown in the concentration table: **concentrations** (default setting), peak **heights**, peak **areas**, **volumes**.


Number Number of the component (read-only).

Name Name of the component taken from the component table (read-only).

This run Contains concentrations (or other chosen values) obtained in the current run (read-only, except for universal component).

Level 1...Level N Contains concentrations (or other chosen values) of components for corresponding **Calibration level**.

<Add> Add a new **Calibration level** to the concentration table. The following window appears:



Create calibration level

Number of calibration level to add.

Fill level with concentration

Concentration to be filled in for all components.

Copy concentrations from level


Number of calibration level from which concentrations are copied.

Calibrate immediately

Recalibration with new calibration levels.

<Delete> Delete the current calibration level (where cursor is situated) from the concentration table.

<Calibrate> Calibrate the current chromatogram with the selected calibration level. The following window appears:



Level

Calibration level to be used for recalibration.

Calibration curve

790 PERSONAL IC / Method / Calibration / Graphs

Selecting this menu item or clicking on **<Graphs>** in the component table opens the **Component** window.

Component - sodium

Copy to clipboard Print/preview...

Q = 2.67996·A

RSD = 0.000 %; Corr. = 0.00000

Concentration

6.25

1

Area

5 10 15 E+00

k2= 0 k1= 2.67996 k0= 0

k3= 0 Exclude point

Level	Conc.	Area	File	Used
1	5	18.7	k4110840.chw	Yes

Component: sodium

Retention time: 3.109

Concentration: 5.000

Calibration method: External standard

Local

Base: Area

Reference channel: Cond

Formula: Y=K1·X

Statistical weight: 1

Standard component:

Concentration of standard: 100

OK Cancel Help

In this window the results of the calibration are shown for each component together with the calibration curve. The parameters for calculating the calibration curve can also be entered here.

Calibration results and curve

Analytical expression A line above the calibration curve looking like

$$Q = k_3 \cdot A^3 + k_2 \cdot A^2 + k_1 \cdot A + k_0.$$

This formula is used for approximation of the calibration curve.

RSD RSD (Residual Standard Deviation) value to evaluate the error of calibration curve approximation.

Corr. Correlation coefficient value is available only in the case of linear calibrations without weighting of data points.

Calibration curve Plot of measured calibration points and calculated calibration curve.

k0, k1, k2, k3 Calibration coefficients **k0**, **k1**, **k2** and **k3** (coefficients of the calibration formula).

Calibration points table

This table contains the basic information used to construct the calibration curve:

Level	Number of Calibration level .
Conc.	Concentration of the current component in the calibration sample. It is taken from the concentration table.
Area (or Height)	Peak height or area of the current component, depending on calibration base.
File	File name of calibration chromatogram that stores data on the given calibration level.
Used	Information whether the calibration level is used for calculation of the calibration curve or not.
<Exclude point>	Exclude the selected calibration point from the list and recalculate calibration coefficients for the current component. Repeated pressing includes the point again. You can exclude points that drop out of the calibration curve watching for RSD values.

Component information

Component	Allows to select the current component from the list. It is possible also to scroll the list of components by mouse, using the special arrow buttons on the right.
Retention time	Retention time of the selected component (read-only).
Concentration	Concentration of the selected component (read-only).

Calibration parameters


Calibration method	Method that is used for calibration procedure. There are three basic methods of calibration:
External standard calibration	Absolute calibration. It is the basic calibration procedure in ion chromatography.
Internal standard calibration	Relative calibration.
Tabulated calibration	Relative response factor calibration. It is a simplified method of external standard calibration.
Local	If this checkbox is marked then the corresponding parameter is valid for the current component only. Otherwise, this parameter is global (i.e. default for all other components).

Base	Base for calculations (Area or Height) in quantification and calibration procedures. Indicates, which parameter is to be used as a peak response.
Reference channel	Channel which is used to measure area or height (read-only).
Formula	Formula for calibration function. There are six possible calibration dependencies for linear and nonlinear calibration curves.
Statistical weight	Weighting parameter used to weight calibration points deviations on calculating coefficients of the calibration formula. Apart from the uniform weighting 1 the values 1/x , 1/x² , 1/y and 1/y² are also possible.
Standard component	Name of the standard component in the case of Internal standard calibration or Tabulated calibration method (details see on-line help).
Concentration of standard	Concentration of the standard component in the case of Internal standard calibration or Tabulated calibration method (details see on-line help).

Update calibration

790 PERSONAL IC / Method / Calibration / Update

Selection of this menu item opens the **Update** window.



Retention time	Replace the expected retention times of the components with the retention times of the corresponding peaks obtained in the current run.
Indexes	Fill Index field for components with zero index.
Coefficients	Recalculate all calibration coefficients for all components (performs re-calibration).
Coef. (balance)	Recalculate coefficient of the universal component. The aim of this action is to adjust the re-

sponse factor of this component so, that the total amount of the injected components becomes equal to the declared value. This declared summary amount is entered into the **Concentration** field of the universal component for **This run** level.

Load and save calibration data

790 PERSONAL IC / Method / Calibration / Load from method

Load the calibration from the method. This option is used to update the calibration for the current run by the one taken from the method used for acquisition.

790 PERSONAL IC / Method / Calibration / Save to method

Save the current calibration to the method. This option is used to transfer the updated calibration from the run to the current method (i.e. the method associated with the current chromatogram window).

Import and export calibration data

790 PERSONAL IC / Method / Calibration / Import calibration

This option imports a calibration from a calibration file ***.cal** to the current method or chromatogram. It is used to transfer the updated calibration from chromatogram-to-chromatogram or from method-to-method. Different methods can be processed in this case using the same calibration.

790 PERSONAL IC / Method / Calibration / Export calibration

This option exports a calibration to a in calibration file ***.cal**. It is useful to transfer a calibration from chromatogram-to-chromatogram for chromatograms that were obtained by different methods.

Put out calibration curves

COMPONENT / Copy to clipboard

This option allows to copy the calibration curve of the selected component to the clipboard so that it is available for other Windows applications, such as WinWord, Excel, etc.

COMPONENT / Print/Preview / Preview this

Display calibration curve print preview of the selected component.

COMPONENT / Print/Preview / Preview all

Display calibration curve print preview of all components.

COMPONENT / Print/Preview / Print this

Print selected calibration curve of the selected component.

COMPONENT / Print/Preview / Print all

Print calibration curves of all components.

4.4.6 Report output


Report options



790 PERSONAL IC / Method / Report options

790 PERSONAL IC / Process / Make report

The **Report options** window is divided into several regions that combine parameters for report output on their functionality and includes several different areas. The report output is started by clicking on **<Report>**.



Items to report

This part of the **Report options** window contains a list of important report items that can be included into the report on the user's choice:

General

General description on the analysis from the Passport. The **Report date** and the logged-in user **Printed by** are also put out.

Report date:	05/04/2000 12:34:00
Printed by:	Roland Dörig
Ident:	tap water
Analysis from:	11/04/2000 17:45:46
File:	k4111745.chw Last save: 13/04/2000 06:33:12
Method:	Anion1.mtw Last save: 11/04/2000 15:13:10
Run operator:	Urs Waldburger
Analysis number:	29

Sample Sample information from the **Sample** tab of the passport.

SAMPLE Info 1:	Tap water from Application Laboratory
Info 2:	Herisau, Switzerland
Vial number:	2
Volume:	20.0 µl
Dilution:	1.00
Amount:	1.0000


Column Column information from the **Column** tab of the passport.

COLUMN:	6.1005.320 Metrosep A Supp3
Size:	4.6 x 250 mm
Number:	A106
Part.size:	9.0 µm

Eluent Eluent information from the **Eluent** tab of the passport.

ELUENT A:	1.7 mM NaHCO3 / 1.8 mM Na2CO3
B:	
C:	
Flow:	1.00 mL/min
Temperature:	20.0°C
Pressure:	8.0 MPa

Chromatogram plot Plot of the chromatogram (only for **Printer** or **File** destination).



Peak table Peak table with results. The structure of the peak table depends on the selected **Quantification method**. The following report is put out as default:

PEAK TABLE					
Quantification method: Custom					
No	Retention min	Height uS/cm	Area uS/cm*sec	Conc. ppm	Name
1	4.15	0.06	0.688	0.057	fluoride
2	5.13	5.10	70.561	9.092	chloride
3	7.84	1.53	37.506	8.743	nitrate
4	13.59	1.09	38.682	6.584	sulfate
4	16.01	7.77	147.437	24.477	

Comment User defined comment from the **Comment** tab of the passport.

```
COMMENT
Beispielmethode zur Wasserbestimmung.
```

GLP Log Method and Data GLP messages from the **Method Log** and **Data Log** tabs of the passport.

```
METHOD GLP LOG
13/04/2000 08:32:49 W. Terzer

04/05/2000 10:45:23
790 Personal IC 1.0 Apr 27 2000 14:05:01

DATA GLP LOG
Analysing data of channel :Cond
Time constant, ms: 100.00
Sigma: 0.0
Rejection difference: 1.0

11/04/2000 18:01:45
Status: Processing
Data saved to disk file c:\programme\ic 790\data\K4111745.CHW
```

More items to report

This part of the **Report options** window contains a list of rarely used report items that can be included into the report on the user's choice:

Acquisition Parameters used for data acquisition.

```
ACQUISITION PARAMETERS
Method duration: 16.00 min
Run duration: 16.01 min
Measurements (method): 9606
Measurements (run): 9606
Freq.divisor: 1
Sampling: 10.00 pts/sec
Device: 790 Personal IC
```

Integration Integration parameters and integration events.

```
INTEGRATION DEFAULTS
Delay: 3.60 min
Width: 5.00 sec
Slope: 20.00
MinHeight: 0.01

No. min
1 0.00 Enable valley-to-valley
2 6.00 Disable valley-to-valley
3 9.50 Set peak end
```

Calibration defaults Calibration parameters defined in the component window and on the **Math** tab of the **Method setup** window.

CALIBRATION	
Method:	External standard
Response:	Area
Standard:	No
IDENTIFICATION	
Reference peaks:	Height
Other peaks:	Time
Retention units:	min

Component table Component table.

COMPONENT TABLE							
No	Retention	Window%	RF	Conc.	Index	Type	Group Name
1	4.05	5.0	1.660e+00	0.06	0.000		0 fluoride
2	5.01	5.0	2.735e+00	9.09	0.000	S	0 chloride
3	5.82	5.0	4.137e+00	0.00	0.000		0 nitrite
4	6.63	5.0	6.027e+00	0.00	0.000		0 bromide
5	7.74	5.0	4.763e+00	8.74	0.000		0 nitrate
6	11.67	5.0	9.164e+00	0.00	0.000		0 phosphate
7	13.20	5.0	3.449e+00	6.58	0.000		0 sulfate

Calibration results Calibration results of the component window. For each component a new page with calibration results and calibration curve is put out.

CALIBRATION OF COMPONENT chloride

Method: Anion1.mtw
 Equation: $Q = -0.0022434 \cdot A^2 + 2.73542 \cdot A$
 RSD: 0.072 %
 Correlation coefficient: 0.999709

K3 = 0 K2 = -0.0022434 K1 = 2.73542 K0 = 0
 Base: Area
 Ref.channel: Cond
 ISTD:
 Formula: $Y = K2 \cdot X^2 + K1 \cdot X$
 Weight: 1

Level	Height	Area	Conc.	Vol/Dil	Retention	Used File
1	0.6004	7.311	1	20	5.013	Yes k4101634.chw
2	3.092	37.74	5	20	5.013	Yes k4101653.chw
3	14.64	169.9	20	20	5.013	Yes k4101712.chw

Channel table

Parameters of the data acquisition channel used.

CHANNELS TABLE										
No	Name	Units	Input	Minimum	Zero	Maximum	Range	Coefficient	Noise	Shift
1	Cond	uS/cm	1	-8388607	0	8388607	8.389e+01	1.000e-05	1.52	0
#	Name/Units	Noise		RMS	PeakToPeak	Drift/hour				
1	Cond	1.5		52871.6	512493.0	4232.464				
	uS/cm	1.52e-05		0.529	5.12	0.0423				

Name	Channel's name. This text appears as Channel label in the chromatogram.
Units	Units of detector response ($\mu\text{S}/\text{cm}$).
Input	Number of data acquisition channel.
Minimum	Minimum value of the linear range of the AD converter (in ADC conversion units). This parameter is used to detect an underflow condition.
Zero	A signal level on a baseline (in ADC conversion units). This parameter is a result of the ADC calibration.
Maximum	Maximum value of the linear range of the AD interface (in ADC conversion units). This parameter is used to detect overflow condition.
Range	An input signal value in maximum in $\mu\text{S}/\text{cm}$. Range = (Maximum - Zero) • Coef
Coef	ADC sensitivity coefficient (weight of ADC bit in $\mu\text{S}/\text{cm}$).
Noise	Estimated baseline noise value of the channel in ADC conversion units (bits).
Shift	<i>Not used for 790 Personal IC.</i>

Report destination

The following output devices can be selected simultaneously in any combination as targets for report output:

Screen	Output report (without curves) to screen.
Printer	Output report to printer.
File	Output report to file (see <i>File output options</i>).

Peak table

The following parameters can be selected for peak table output with determination results:

Quantification method Selection of the method to calculate concentrations of components:

Response normalization

The method normalizes the sum of responses for all peaks of the chromatogram to the **NORM** factor entered in the **Total % for normalization** field:

$$R_i\% = \text{NORM} \cdot R_i / \sum R_i$$

All peaks are calculated and reported no matter whether they have been recognized or not.

The peak table includes the following columns: **number, retention time, area + area% or height + height%, name**

Normalized concentration

The method normalizes the sum of absolute concentrations for all peaks of the chromatogram to the **NORM** factor entered in the **Total % for normalization** field:

$$C_i\% = \text{NORM} \cdot W_i(R_i) / \sum W_i(R_i)$$

As long as no universal component is defined, only peaks with non-zero concentrations are listed.

The peak table includes the following columns: **number, retention time, height, area, response factor, concentration%, name**

Absolute concentration

This quantification method reports raw quantity (absolute concentration) of the component that is calculated directly from the calibration formula:

$$Q_i = W_i(R_i) / V'$$

As long as no universal component is defined, only peaks with non-zero concentrations are listed.

The peak table includes the following columns: **number, retention time, height, area, response factor, concentration, concentration%, name**

Relative concentration

This quantification method uses the internal standard method to calculate the relative concentrations of the components. For this procedure, a **Standard component** must be selected

	<p>and the concentration of this internal standard must be entered in the Concentration of std field.</p> <p>As long as no universal component is defined, only peaks with non-zero concentrations are listed.</p> <p>The peak table includes the following columns: number, retention time, height, area, response factor, relative concentration, relative concentration%, name</p>
Index	<p>This column reports retention indexes of identified components.</p> <p>The peak table includes the following columns: number, retention time, width/2, height, index, name</p>
Column test	<p>Column test quantification method calculates multiple values that are necessary to evaluate column performance.</p> <p>The peak table includes the following columns: number, retention time, K' (capacity factor), TP (number of theoretical plates per column), TP/m (number of theoretical plates per m), HETP/dp (reduced height equivalent to theoretical plate), asym. (peak asymmetry), name</p>
Custom	<p>This quantification method enables to customize the peak table on the user's choice using the <Customize> button.</p> <p>The peak table for the methods supplied contains the following columns as standard: number, retention time, height, area, response factor, raw concentration, name</p>
Standard component	Internal standard component for quantification procedure.
Concentration of std	Concentration of the internal standard component for relative concentration calculations. This value is stored in This run column of the concentration table.
Total % for normalization	A value to which the sum of concentrations is normalized. It is used for Response normalization and Normalized concentrations quantification methods. The default value is 100.
Printing order	Determines the order of the components in the peak table.
By peaks	Lists run results for all peaks detected. Unidentified peaks (corresponding to the universal component) are included, but missing components are not listed.

By components	Missing components are always included into the report, even those with concentration equal to zero. In the case of presence of universal component the summary for all unidentified peaks is presented as a single line.
<Customize>	Calls up a box with a list of columns that can be included into the report (see <i>Report elements</i>). It is available only if Quantification method = Custom has been selected.
Report all peaks	Checkbox that includes all peaks irrespective of the concentration of the substance. Otherwise the report generator excludes lines with zero concentration of the component. It is available only if Quantification method = Custom has been selected.
Groups	If this checkbox is checked, a separate peak table report is put out for every group specified in the component table.

Template options

For report output using a report template the following parameters can be modified:

Template	Selection of the report template file *.rtt . Several templates for different languages and applications are available.
Separator	Defines which character separates report columns in tables. A separator other than Space is useful when the report is printed to a file and this file is imported into a spreadsheet software.
Tab size	Sets tab stop value for the screen window. Important only when Separator = Tabulation has been selected or there are tab characters in the template file *.rtt .

File output options

If **File** has been selected for the report destination, the following parameters can be modified:

Directory	Directory for report output. For creation of a new directory, use the <Browse> button.
Name	File name to save report to. The report is saved in text form in the ANSI or ASCII format. So add an extension like *.txt to the file name. If the chromatogram plot checkbox is checked, the plot is saved in a separate file *.wmf in the WMF format under the same name.

Mode	The two options Overwrite (overwrite the file) or Append (append to existing file) are available.
Character set	The two options Windows (ANSI) or DOS (ASCII) are available. These settings are essential for printing of non-English symbols (e.g. ö, ä, ü).
Custom program	Path and name of the program to be started after report output. This option enables to transfer the report to a database, an electronic table or another application for further processing.

Report elements

If **Quantification method = Custom** has been selected, the peak table can be customized on the user's choice. After clicking the **<Customize>** button, the following columns can be included in the peak table.

number	Peak number.
retention time	Retention time of the component (in minutes, irrespective of the chosen retention units on the chromatogram graph axes). The total value in the column is equal to the chromatogram duration.
halfwidth	Width of the peak at half height (in minutes).
height	Height of the peak (in $\mu\text{S}/\text{cm}$). The total value in this column is the sum of heights for all identified peaks.
height%	Normalization of peak heights for all peaks using the NORM value entered in the Total % for normalization field (default setting 100%): $H_i\% = \text{NORM} \cdot H_i / \sum H_i$
area	Area of the peak. Depends on the units on the X and Y axes of the chromatogram. The total value in this column is the sum of areas for all identified peaks (including universal component).
area%	Normalization of peak areas for all peaks using the NORM value entered in the Total % for normalization field (default setting 100%): $A_i\% = \text{NORM} \cdot A_i / \sum A_i$
capacity factor	The capacity factor k'_i of the component is equal to the ratio of its corrected retention time ($t - t_0$) to the void time of the system t_0 : $k'_i = (t_i - t_0) / t_0$ <p>Total for this column is equal to the capacity factor of the last peak of the chromatogram.</p>

resolution	<p>Resolution R for two neighboring peaks is calculated as:</p> $R = (t_{i+1} - t_i) / (w_{0.607i} + w_{0.607(i+1)})$ <p>where i and i+1 indexes refer to the neighboring peaks, and w_{0.607} stands for the peak width at 60.7 % of the peak height.</p>
effectivity, TP	<p>Effectivity for the peak in number of theoretical plates. The number of theoretical plates N_i per column for a chosen peak is calculated for a chromatographic peak by one of two formulas:</p> $N_i = 2 PI (t_i \cdot H_i / A_i)^2,$ <p>where PI = 3.1415926..., t_i = retention time, H_i = height, A_i = area of the peak. The more commonly used formula is:</p> $N_i = 5.54 (t_i / w_i)^2,$ <p>where w_i is the width on the half-height of the peak. The first formula offers better estimates for fused or unresolved peaks, because the half-width errors for those peaks are much greater than height or area errors.</p> <p>Total for this column includes average value for the peaks listed.</p>
effectivity, TP/m	<p>Effectivity for the peak in number of theoretical plates per meter. The number of theoretical plates per meter N' for the given component is calculated as:</p> $N' = N_i \cdot 1000 / L,$ <p>where L is length of the column in mm and N_i is effectivity of the column for i-th component.</p> <p>Total for this column includes average value for the peaks listed.</p>
reduced TP height, HETP/dp	<p>The height of theoretical plate divided by particle size, called also reduced height, is calculated by formula:</p> $H_i = 1000 \cdot L / (N_i dp).$ <p>where L is length of the column in mm, dp is particle diameter in μm.</p>
asymmetry	<p>Peak asymmetry A_s is calculated at ¹/₁₀ of the peak height as a ratio of width after the top of the peak w₂ to the width before the top w₁.</p> $A_s = w_2 / w_1$
response factor	<p>Coefficient k₁ of the calibration curve.</p>


raw concentration	<p>Absolute concentrations (raw quantity) of the components, calculated as</p> $C_i = W_i(R_i) / V'$ <p>Total value for the column is the sum of concentrations for all components.</p>
concentration%	<p>Normalization of concentrations for all peaks using the NORM value entered in the Total % for normalization field (default setting 100%):</p> $C_i\% = \text{NORM} \cdot W_i(R_i) / \sum W_i(R_i)$ <p>Total value for the column equals NORM for the main table, and the sum of concentrations in the group for the group table.</p>
rel. concentration	<p>Relative concentration of the component relative to the standard component, assumes concentration of one of the components to be known in advance. The concentration of the component is calculated by the formula:</p> $C'_i = W_i(R_i) / V' = W_i(R_i) \cdot C_s / W_s(R_s),$ <p>where $V' = W_s(R_s) / C_s$ is the effective (reduced) volume of injected sample. The standard component concentration C_s is entered by the user.</p> <p>Total value for relative concentration column does not include the concentration of the standard component.</p>
rel. concentration%	<p>Relative concentration percent differs from concentration% by exclusion of standard component(s) from summation, so that sum of all concentrations (excluding standard) equals NORM.</p> <p>Total value for the column equals NORM value for main peak table, and the sum of concentrations in the group for Group table.</p>
index	<p>Linear or logarithmic retention index for identified components.</p> <p>Total value for the column is a weighted index average with absolute concentration as weight.</p> $I = \sum (I_i C_i) / \sum C_i.$
type	<p>The types of components are designated by a one-letter code:</p> <ul style="list-style-type: none"> R Reference component (used for peak identification). S Standard component. C Calibration standard in the case it is different from the quantification standard.

?	Some points of corresponding peak are out of ADC or detector range – suspicious result.
!	Component concentration is outside of minimum (min C) and maximum limits (max C), set in the component table.
p	Special component (it has one of the special flags checked on in the COMPONENT window).
N	The component response is outside of calibrated region.
	In the case of peaks additional letters precede the component type information, indicating, how peaks are separated from other peaks, for example:
BD_	Peak that starts on the baseline (B) and ends on the drop line (D) that separates it from another adjacent peak.
BBR	Special case of a rider (R) peak that is tangentially separated from the main one. A main peak will have a third letter H (horse) in this case. The complete component type may look like: BBD : !R.
group	Number of the group for the component.
spectral ratio	<i>Not used for 790 Personal IC.</i>
name	Name of the component.
file name	File name of chromatogram. This column is very useful for processing of exported data.
ident	Title (identifier) for the chromatogram. This column is useful for the processing of exported data.

4.5 Chromatograms

4.5.1 Chromatogram window

The **CHROMATOGRAM** window is used to show a running or recorded chromatogram.



The **title bar** of the chromatogram window contains buttons for minimizing, maximizing and closing the window. The window name consists of the elements "**Ident (method name) chromatogram name**". A star (*) at the end of the name indicates that the chromatogram has been changed since the last saving.

The **status bar** contains two fields. The first field indicates the current measurement status, one of the listed below:

Ready	Chromatogram is ready to start.
Waiting	Chromatogram is waiting for the first measuring points.
Measure	Chromatogram is being measured.
Measure(Baseline)	Recording of baseline.
Finished	Measurement finished, but the chromatogram is not processed.
Processing	Chromatogram is being processed after finishing.
Failure	Failure, e.g. unexpected pump stop, etc.
COM error	Error on COM port interface.

During active data acquisition the elapsed time, the duration of the chromatogram or program and the current X and Y values are shown in the second field of the status bar. If the peak editor is switched on then the position of the cursor will be shown.

A chromatogram can be scaled with the help of keyboard or mouse functions or through the **Chromatogram axes** window opened with **790 PERSONAL IC / View / Appearance / Chromatogram axes**. Some of the window control functions are collected in the **Window** menu.


4.5.2 Chromatogram file handling

Open chromatogram



790 PERSONAL IC / File / Open / Chromatogram

Load an existing chromatogram file (*.chw) from the **Data** directory and open the chromatogram window. The following window appears:



Filename Filename (common wildcards as * and ? are possible).

Selected Number of selected chromatograms and their total size.

Directory Path name of the working directory.

File window Box with the list of files residing in the working directory sorted by time (the last recorded chromatogram is at the top of the list). As additional short information on the chromatogram the two parameters **Ident** and **Method** are also displayed. Files in this box can be selected by moving the selection bar to the desired chromatogram and pressing the [Space] key or by a single mouse click.

It is possible to select several chromatograms at once. Press [Shift] and the left mouse button to select all chromatograms up to last item clicked, press [Ctrl] and the left mouse button

to add a single chromatogram to the selection already made. All selected items are painted.

Directories	Box with the list of directories. Allows to change the working directory defined in the method, if necessary.
Chromatogram description	A set of fields with the description of the current chromatogram. Selected fields from the passport are listed here.
<OK>	Load all selected chromatogram files (*.chw) and open each in its own window.
<Cancel>	Close this window without any action.
<To Batch>	Add selected chromatograms to the specified batch reprocessing file.
<Copy>	Copy selected chromatograms to the specified location.
<Move>	Move selected chromatograms to the specified location.
<Delete>	Delete selected chromatograms and move them to the Windows wastebasket.

Save chromatogram




790 PERSONAL IC / File / Save / Chromatogram

Save the selected chromatogram in a chromatogram file (*.chw) in the working directory. If this chromatogram has been already stored in this directory, the message **File ... exists. Overwrite?** appears. If it is preferable to overwrite the previous copy (for example, after the reprocessing), press the **<OK>** button. Pressing **<No>** saves the chromatogram into the new file (a new copy of the chromatogram is created).

Close chromatogram

790 PERSONAL IC / File / Close

Close the selected chromatogram window. If the data were not saved or the method was modified, a warning appears. Usually it is more convenient to close the window by clicking the  button at the top corner area of the window.

Delete chromatogram

790 PERSONAL IC / File / Delete

Delete the selected chromatogram window. A warning will be generated by the software in any case.


Export chromatogram

790 PERSONAL IC / File / Export / AIA file

Export the selected chromatogram in the AIA format (Analytical Instrument Association) as CDF file (*.cdf).

790 PERSONAL IC / File / Export / Raw data to txt

Export chromatographic raw data into an ASCII text file (*.txt). The following window appears:



Represent raw data as Output of raw data:

Default units (float) Y values are exported in default units for the channel ($\mu\text{S}/\text{cm}$).

ADC response (integer) Y values are exported in bits.

Include retention time

In seconds (float) X values are exported in seconds.

As point number X values are exported as point number.

4.5.3 Graphical representation

Appearance



790 PERSONAL IC / View / Appearance

This menu item opens the **Appearance** window. It determines the appearance of the chromatogram and consists of four tabs:


Chromatogram axes	Scaling of chromatogram axes.
Labels	Settings for peak labels and baseline drawing.
Select channel	Select channels to be displayed (only available for multi-channel chromatograms).
Colors	Color settings for chromatogram elements.



*The appearance settings for a chromatogram are not saved automatically on closing the chromatogram. In order to have the same appearance after reopening the chromatogram, it must be saved explicitly after modification of these parameters with **790 PERSONAL IC / File / Save / Chromatogram**.*

Chromatogram axes

Tab **Chromatogram axes** of the **Appearance** window with parameters for scaling of the chromatogram axes.




Time axis

X from	Beginning of the window on X axis.
X to	End of the window on X axis.
X units	Selection of unit (Retention unit) for X axis.

Response axis	
Y from	Beginning of the window on Y axis.
Y to	End of the window on Y axis.
Tic marks	Plot axes labels and tic marks on Y axis with no , relative or absolute scaling.
<hr/>	
View all	Set scales on X and Y axes so that the whole chromatogram is visible.
<hr/>	
Drift compensation	Turns the chromatogram screen image so that the last and first points of the chromatogram are on the same level. This option is useful for gradient chromatograms. Drift compensation has no effect while a chromatogram is running.
<hr/>	
Grid	Plot dotted grid lines in the selected chromatogram window.
<hr/>	
Set all	Set scales on X and Y axes of all opened chromatograms automatically to the settings in the selected chromatogram window.

Labels

Tab **Labels** of the **Appearance** window with parameters for peak labeling and baseline drawing.




Peak labels

None	No peak labels.
Peak number	Peak number.
Retention time	Retention time.
Component name	Name of the component.
Component name + Quantity	Name and quantity of the component.

Show baseline and peaks	Baselines are shown.
Baseline marker	Start and end of baselines are marked.
Label always visible	Label is shown always with zooming.
How to draw chromatogram point	
every	Normal chromatogram drawing with all raw data points.
middle	Chromatogram drawing with additional smoothing of raw data points.
random	Fast chromatogram draw for slow PCs.
Channel labels	
At start	Channel label is shown at the start of the chromatogram.
In corner	Channel label is shown in the upper right corner of the chromatogram.
Set all	Set peak labeling of all opened chromatograms automatically to the settings in the selected chromatogram window.

Select channel


Tab **Select channel** of the **Appearance** window for selection of channels for multi-channel chromatograms.



Selection window	The selection window shows all chromatograms of the batch reprocessing queue which can be selected for display. If total is selected, all chromatograms are added and the resulting addition curve is shown.
<Show all>	Select all chromatograms in the selection window.
<Hide all>	Deselect all chromatograms in the selection window.

Colors

Tab **Colors** of the **Appearance** window with color settings parameters for chromatograms.



this color is default	Reset the selected element to the default color.
<Choose>	Choose a new color for the selected element. An example of the chosen color is shown beside the button.
Hair	Cursor's color.
Background	Background color.
Axes	Axes and axes labels color.
Baseline	Baseline color.
AlterGround	Alternative background color for Measure baseline option.
Channel 1...8	Color of the selected channel. Each channel can be of a unique color.
Line width	Line width for the selected element (Axes or Channel 1...8) in % of selected plot font height. Range: 0 ... 232
Set all	Set colors of all opened chromatograms automatically to the settings in the selected chromatogram window.
<Load defaults>	Set colors in the selected chromatogram window to the set defaults colors.
<Save defaults>	Save color settings in the selected chromatogram window as default colors.

Other graphical functions

790 PERSONAL IC / View / X full scale

By selecting this menu item or pressing [Ctrl] + [Home] the X axis is scaled so that the chromatogram perfectly fits the window horizontally.

790 PERSONAL IC / View / Y full scale

By selecting this menu item or pressing [Ctrl] + [End] the Y axis is scaled so that the chromatogram perfectly fits the window vertically.



790 PERSONAL IC / View / View all

By selecting this menu item or pressing [Alt] + [V] the X and Y axes are scaled so that the chromatogram perfectly fits the window horizontally and vertically.

790 PERSONAL IC / View / Recorder autoscale

This option allows to look at the chromatogram so that the last point is always visible on the screen while data are acquired.

If the **Recorder autoscale** option is switched on, the following autoscaling rules apply:

- If the last acquired point comes out of the window down, autozero is performed.
- If the last point is too high, the recorder scale shrinks twice until the point comes to screen.
- If the last point comes outside of the plotting area to the right, the window is shifted half-screen right.

If the **Recorder autoscale** option is switched off, the window scale does not change automatically during data acquisition.

4.5.4 Peak editor


Switching on/off the peak editor



790 PERSONAL IC / Process / Peak editor

The peak editor mode is used for subsequent manual correction of the automatic integration for chromatogram peaks (see *section 4.4.4*). It can be used to select the most important points for the peak evaluation (start, end and top of peak, valley between peaks) and move them to the required position.

If the peak editor mode is enabled, the vertical bar cursor appears in the chromatogram window. At the same time, the **Peak** menu appears in the menu bar and the peak editor icons appear in the icon bar of the main window.

The peak editor mode is switched on or off by clicking on  or by pressing [Alt] + [C]. It is also possible to click the right mouse button somewhere in the chromatogram window and to select the **Peak editor** item.

You cannot edit the peak pattern while the component table is active and vice versa, when the peak editor mode is switched on the component table option is disabled.

Peak editor functions

The functions of the peak editor can be triggered with the corresponding menu items of the **Peak** menu, by the peak editor symbols in the symbol bar or with key combinations.



790 PERSONAL IC / Peak / Undo

Undo the last action.



790 PERSONAL IC / Peak / Insert peak

[Insert]

Add a peak to the chromatogram.



790 PERSONAL IC / Peak / Delete peak

[Delete]

Delete selected peak.

790 PERSONAL IC / Peak / Select nearest point

Move the cursor to the nearest start, top, end, or valley point and select the peak.



790 PERSONAL IC / Peak / Select start point

Move the cursor to the beginning of the nearest peak and select the peak.



790 PERSONAL IC / Peak / Select top point

Move the cursor to the top of the nearest peak and select the peak.

**790 PERSONAL IC / Peak / Select end point**

Move the cursor to the end of the nearest peak and select the peak.

**790 PERSONAL IC / Peak / Select valley point**

Move the cursor to the valley of the two nearest peaks.

**790 PERSONAL IC / Peak / Unselect peak**

Delete the selection.

**790 PERSONAL IC / Peak / Move selected point** [-]

Move a selected point (start, top, end, valley) to the cursor position.

**790 PERSONAL IC / Peak / Merge peaks** [+]

Merge two neighboring peaks into a single peak.

**790 PERSONAL IC / Peak / Fuse peaks** [*]

Fuse the beginning of the previous and the end of the next peak at the cursor position.

**790 PERSONAL IC / Peak / Split peaks** [/]

Split a single peak into two peaks at the cursor position.

**790 PERSONAL IC / Peak / Delete all left**

Delete all peaks to the left of the selected point.

**790 Personal IC / PEAK / DELETE ALL RIGHT**

Delete all peaks to the right of the selected point.

Moving the cursor

The cursor can be dragged by the **mouse** when the right mouse button is pressed. Releasing the button will leave the cursor at the new position. The position of the cursor is displayed in the status bar of the chromatogram window.

The cursor can also be moved using the **keyboard**:


[←]	Move cursor left.
[Shift] + [←]	Move cursor left quickly.
[→]	Move cursor right.
[Shift] + [→]	Move cursor right quickly.
[Home]	Move cursor to the beginning of the window.
[End]	Move cursor to the end of window.

4.5.5 Printing

Page layout for printing

790 PERSONAL IC / File / Page layout

The selection of this menu item or clicking on <Page...> in the **Report options** window opens the **Page layout** window for page layout parameters.



Length units	Selection of length units.
inches	Inches.
centimeters	Centimeters.
Page margins	Selection of page margins.
Left	Left margin.
Right	Right margin.
Top	Top margin.
Bottom	Bottom margin.
Chromatogram plot size	Selection of plot area for chromatograms.
Width	Width of plot area.
Height	Height of plot area.
Calibration graph size	Selection of plot area for calibration curves.
Width	Width of plot area.
Height	Height of plot area.
<Load default>	Load default values for page layout parameters.
<Save as default>	Save set page layout parameters as default values.

Printer settings

790 PERSONAL IC / File / Printer setup

By selecting this menu item the **PRINT SETUP** window is opened where printer, paper size and format can be defined.

Print preview



790 PERSONAL IC / File / Preview

By selecting this menu item the **Preview** window appears on the screen where the report is shown in the appearance formatted for the desired printer as defined in the **Report options** window.

Printing



790 PERSONAL IC / File / Print

By selecting this menu item the **Printing** window is opened where printer, printing range and number of copies can be defined. All settings defined in the **Report options** window will be active except the destination of the printout.

In order for the chromatogram to be printed, the option **Chromatogram plot** must be switched on in the **Report options** window under **Items to report**.



790 PERSONAL IC / Process / Make report

Selection of this menu item opens the **Report options** window for report output to screen, printer or file (details see *section 4.4.6*).

4.5.6 Miscellaneous functions

Reintegration



790 PERSONAL IC / Process / Reintegrate

This menu item opens the **Integration parameters** window where the reintegration of the chromatogram can be started (details see *section 4.4.4*).

Recalibration


790 PERSONAL IC / Process / Calibrate

This menu item opens the **Recalibration** window for entry of the calibration level and recalibration of the chromatogram using this level (details see *section 4.4.5*).

Subtraction of a chromatogram

790 PERSONAL IC / Process / More / Subtract

This menu item opens the **Subtract chromatogram** window allowing the subtraction of any opened chromatogram from the active, selected chromatogram.




Select the chromatogram that should be subtracted, click on **<Subtract>** and then on **790 PERSONAL IC / View / View all**. The result is shown in the active chromatogram window, which is considered as a new chromatogram and will be stored under a new name to avoid overwriting of old data.

Data compression

790 PERSONAL IC / Process / More / Compress

Selection of this menu item opens the **Raw data compression** window for compression of the raw data of the active chromatogram by summing up several neighboring points.



Divisor Coefficient of compression (the number of data points is reduced by this value). The default value of the compression coefficient is calculated so that the width at half-height of the narrowest peak will be digitized by at least 30 points. If the user enters a coefficient of compression that exceeds the value offered by the software, integration precision may decrease.

Invert chromatogram

790 PERSONAL IC / Process / More / Invert!

This menu item inverts the response curves for all channels of the chromatogram so that negative peaks become positive and vice versa (useful for chromatograms with wrong input polarity).

4.6 Sample queue

A **sample queue** is a table containing sample-specific data which is used to facilitate work with autosamplers and multiple analyses. The number of rows of the sample queue table defines the number of determinations that are run automatically. Once the queue has been started, the sample-specific data are transferred line-by-line to the running determination overwriting the corresponding fields of the method.

A sample queue is stored in a sample queue file ***.que** in a subfolder of the **Methods** directory.

4.6.1 Sample queue file handling

Open sample queue

SYSTEM / System / Sample queue

This menu item opens the **OPEN** window for opening an existing sample queue file ***.que** or creation of a new sample queue file ***.que** by entering a new name. After confirmation with **<OK>**, the sample queue overview table is opened (see *section 4.6.2*).

Save sample queue



QUEUE EDITOR / File / Save

Save the sample queue table in a sample queue file (***.que**) in the working directory. The sample queue editor window remains open.



QUEUE EDITOR / File / Save & exit

Save the sample queue table in a sample queue file (***.que**) in the working directory and close the sample queue editor window.

SAMPLE QUEUE OVERVIEW / File / Save as

Save a copy of the current sample queue into a new sample queue file (***.que**) in the working directory.

Delete sample queue

SAMPLE QUEUE OVERVIEW / File / Delete

Close the sample queue overview window and delete the sample queue file (***.que**).

4.6.2 Sample queue control

Sample queue overview table

SYSTEM / System / Sample queue

This menu item opens the **OPEN** window for opening an existing sample queue file *.que or creation of a new sample queue file *.que by entering a new name. After confirmation with <OK> the sample queue overview table is opened showing the current status of the sample queue:

No	System	Ident	Vial	Inj-s	Started	Level	Volume	Amount
1	Anion.smt	Probe 1	1	1	0	0	20.00	1.00
2	Anion.smt	Probe 2	2	1	0	0	20.00	1.00
3	Anion.smt	Probe 3	3	1	0	0	20.00	1.00
4	Anion.smt	Probe 4	4	1	0	0	20.00	1.00
5	Anion.smt	Probe 5	5	1	0	0	20.00	1.00
6	Anion.smt	Probe 6	6	1	0	0	20.00	1.00

No Row number.

System System file to be used for the determination.

Ident User defined identifier for chromatogram. Will be placed into appropriate passport field when starting the chromatogram.

Vial Autosampler vial position to take sample from.



This value will not be transferred automatically to the autosampler. Make sure that the vial positions entered in this column are identical to the vial positions used at the autosampler.

Inj-s Number of injections for the same vial position (only 1 injection is permissible for **Level ≥ 1**).




This value will not be transferred automatically to the autosampler. Make sure that the numbers of injections entered in this column are identical to the numbers of injections used at the autosampler.

Started Indication whether the sample has been started or not:

- 0** sample not started
- 1** sample started
- 2...n** number of injections done

Level	Calibration level (see <i>section 4.4.5</i>) for the sample: Level 0 stays for normal analysis, levels 1 and greater for calibration runs. Correct filling of this column enables to calculate calibration coefficients automatically.
Volume	Injected volume of sample in μL .
Amount	Sample amount. If this value is different for the calibration run (c) and the sample run (s), the component concentrations of the sample are calculated as follows: $C_s = C_c \cdot \text{Amount}_s / \text{Amount}_c$
Shut down system after the queue finishes	If this option is enabled, the system hardware is shut down after the sample queue has been finished (all pumps are switched off), but it remains connected.
Close queue dialog window after the queue finishes	If this option is checked, the sample queue window is closed automatically after the sample queue has been finished.
<Start>	Start execution of the sample queue from the first line with Started = 0 .
<Pause>	Stop execution of the sample queue after the current determination has been finished.
<Edit>	Open the sample queue editor program for editing the sample queue table (see <i>section 4.6.3</i>).
<Close>	Close the sample queue overview table.

Start sample queue




SAMPLE QUEUE OVERVIEW / Control / Start

Start execution of the sample queue from the first row with **Started = 0**. For each row, the sample-specific data are transferred to the method of the current determination and the parameter **Started** is set to **1** for this row.

If processing the sample queue is interrupted with **<Pause>** then it can be continued again with **<Start>** (with a running chromatogram only when data acquisition is finished).

Pause sample queue



SAMPLE QUEUE OVERVIEW / Control / Pause

Interrupt execution of the sample queue after the current determination has been finished.

If processing the sample queue is interrupted with <Pause> then it can be continued again with <Start> (with a running chromatogram only when data acquisition is finished).

Cancel last run

SAMPLE QUEUE OVERVIEW / Control / Cancel last run

Cancel the last run (set **Started = 0** for this row). If the queue is restarted, this row is executed again.

This function is only available if the sample queue has been interrupted by pressing the <Pause> button.

Reset sample queue


SAMPLE QUEUE OVERVIEW / Control / Reset

Reset sample queue table (set **Started = 0** for all rows).

This function is only available if the sample queue has been interrupted by pressing the <Pause> button.

4.6.3 Sample queue editor



Open queue editor window



SAMPLE QUEUE OVERVIEW / Control / Edit

This menu item opens the **Queue editor** window for editing the sample queue table.

	System	Ident	Vial	Volume	Dilution	Amount	Internal Standard Amount	Level	Injections	Done	Sample Info 1	Sample Info 2
1	Anion.smt	Probe 1	1	20.0	1.0	1.0	100.0	0	1	0		
2	Anion.smt	Probe 2	2	20.0	1.0	1.0	100.0	0	1	0		
3	Anion.smt	Probe 3	3	20.0	1.0	1.0	100.0	0	1	0		
4	Anion.smt	Probe 4	4	20.0	1.0	1.0	100.0	0	1	0		
5	Anion.smt	Probe 5	5	20.0	1.0	1.0	100.0	0	1	0		
6	Anion.smt	Probe 6	6	20.0	1.0	1.0	100.0	0	1	0		
7												
8												
9												
10												

No	Row number.
System	System file to be used for the determination. The system is automatically opened and connected.
Ident	User defined identifier for chromatogram. Will be placed into appropriate passport field when starting the chromatogram.
Vial	Autosampler vial position to take sample from. <div style="display: flex; align-items: center; margin-top: 10px;">  <div style="border: 1px solid black; padding: 5px; font-style: italic;"> <p>This value will not be transferred automatically to the autosampler. Make sure that the vial positions entered in this column are identical to the vial positions used at the autosampler.</p> </div> </div>
Volume	Injected volume of sample in μL .
Dilution	Sample dilution prior to injection.
Amount	Sample amount. If this value is different for the calibration run (c) and the sample run (s), the component concentrations of the sample are calculated as follows: $C_s = C_c \cdot \text{Amount}_s / \text{Amount}_c$
Internal standard amount	Concentration of the internal standard component for relative concentration calculations.
Level	Calibration level (see section 4.4.5) for the sample: Level 0 stays for normal analysis, levels 1 and greater for calibration runs. Correct filling of this column enables to calculate calibration coefficients automatically.
Injections	Number of injections for the same vial position (only 1 injection is permissible for Level ≥ 1). <div style="display: flex; align-items: center; margin-top: 10px;">  <div style="border: 1px solid black; padding: 5px; font-style: italic;"> <p>This value will not be transferred automatically to the autosampler. Make sure that the numbers of injections entered in this column are identical to the numbers of injections used at the autosampler.</p> </div> </div>
Done	Indication whether the sample determination has been done or not: 0 determination not yet done 1...n determination or injection done
Sample Info 1	First sample description.
Sample Info 2	Second sample description.

Sample queue editor functions

The functions in the editor window for sample queues can be triggered with the corresponding menu items of the **Edit** menu or with the corresponding symbols in the symbol bar.



QUEUE EDITOR / Edit / Undo

Undo the last modification of the sample queue table.



QUEUE EDITOR / Edit / Cut row(s)

Cut the selected rows of the sample queue table and copy them into the clipboard.



QUEUE EDITOR / Edit / Copy row(s)

Copy the selected rows of the sample queue table into the clipboard.



QUEUE EDITOR / Edit / Paste row(s)

Paste the rows from the clipboard into the sample queue table.



QUEUE EDITOR / Edit / Delete row(s)

Delete the selected rows from the sample queue table.



QUEUE EDITOR / Edit / Duplicate row(s)

Duplicate the selected rows of the sample queue table.



QUEUE EDITOR / Edit / Increment

Fill selected column fields with values automatically incremented by 1. The last character of the first selected field must be a number. This function is available for the fields **Ident**, **Vial**, **Level**, **Sample Info 1** and **Sample Info 2**.



QUEUE EDITOR / Edit / Propagate

Fill selected column fields with the same value as the first field.



QUEUE EDITOR / Edit / Reset

Reset sample queue table (set **Done = 0** for all rows).



QUEUE EDITOR / Edit / Change system

Change system file (*.smt) to be used for the selected rows.

Print sample queue



QUEUE EDITOR / File / Print

Print the sample queue table in landscape format.

Close sample queue editor



QUEUE EDITOR / File / Save & Exit

Save the modified sample queue table and close the sample queue editor window.

QUEUE EDITOR / File / Exit

Close the sample queue editor window. An inquiry appears asking whether or not the sample queue should be saved.

4.7 Batch reprocessing

Reprocessing is understood to be the subsequent reprocessing of a series of chromatograms which have been loaded in a **queue (Batch reprocessing queue)**. For reprocessing according to a selected method the settings for calibration, integration, passport, appearance and report can be altered at will.

A batch reprocessing queue is stored in a batch reprocessing queue file ***.bar** in the **Data** directory.

4.7.1 Batch reprocessing queue file handling

Open batch reprocessing queue

790 PERSONAL IC / File / Open / Batch reprocessing

Load an existing batch reprocessing file (***.bar**) from the **Data** directory and open the **Reprocess** window.



790 PERSONAL IC / File / Open / Last batch

Load the last opened batch reprocessing file (***.bar**) from the **Data** directory and open the **Reprocess** window.

Create new batch reprocessing queue



790 PERSONAL IC / File / Open / Chromatogram

For creation of a new batch reprocessing queue, open the **Chromatogram open** window. Select the desired chromatograms (***.chw**) and click the **<To Batch>** button. Enter the name of the new batch reprocessing file and press **<OK>**. The selected chromatograms are then added to this batch reprocessing queue.



If the reprocessing includes reintegration and/or recalibration, only chromatograms recorded with the same method should be loaded into the batch reprocessing queue.

Save batch reprocessing queue



QUEUE EDITOR / File / Save

Save the batch reprocessing queue in a batch reprocessing queue file (***.bar**) in the working directory. The batch reprocessing queue editor window remains open.



QUEUE EDITOR / File / Save & exit

Save the batch reprocessing queue in a batch reprocessing queue file (***.bar**) in the working directory and close the batch reprocessing editor window.


4.7.2 Perform batch reprocessing

Reprocess options window

790 PERSONAL IC / File / Open / Batch reprocessing

790 PERSONAL IC / File / Open / Last batch

These two menu items are used to open the **Reprocess** window, in which the various options for reprocessing are set and can then be triggered.



Use method from file for reprocessing

Selection of the desired chromatogram file whose method should be used for reprocessing.



If the reprocessing includes reintegration and/or recalibration, only chromatograms recorded with the same method should be loaded into the batch reprocessing queue.

<Open example>

Open the chromatogram selected above in the **Use method...** field.

<Open all files>

Open all chromatograms of the batch reprocessing queue.

<Edit sample table>	Open the batch reprocessing editor program for editing the batch reprocessing table (see <i>section 4.7.3</i>).
----------------------------------	--

Reprocess sample runs	Reprocess all sample chromatograms (calibration level = 0).
------------------------------	--

Reprocess calibration runs	Reprocess all calibration chromatograms (calibration level > 0).
-----------------------------------	--

Update method file in <METHODS> directory after reprocessing	Save the method file *.mtw after reprocessing if the method is changed.
---	--

Reintegrate	Reintegrate chromatograms according to the current settings of the integration parameters and integration events.
--------------------	---



<Edit integration parameters>	Open the Integration parameters window for modification of integration parameters and integration events.
--	--

Recalibrate	Reprocess all calibration chromatograms (if Reprocess calibration runs is switched on) and apply new calibration to all sample chromatograms (if Reprocess sample runs is switched on) by updating the concentration table.
--------------------	---

Default scheme	Default setting for recalibration reprocessing. The two options Apply final calibration... and Forget calibration points... are switched on. A new calibration is performed with the calibration runs and the resulting new calibration parameters (component table, concentration table and calibration curve) are applied to all sample runs.
-----------------------	---

Apply final calibration to all reprocessed files	Apply the updated calibration to all calibration and sample runs. If this option is switched off, the calibration stored in the first chromatogram is used for all other chromatograms.
---	---


Forget calibration points before reprocessing	Forget old calibration points of the calibration curve and perform a new calibration using all calibration runs of the batch reprocessing queue. In this case, each calibration chromatogram adds a new point to the calibration curve. If this option is switched off, the calibration curve stored in the first chromatogram remains active. Each further calibration chromatogram in the batch reprocessing queue adds a new point to this calibration curve.
--	--

Recalculate only	<p>Enable recalculation of chromatograms with values for Volume, Dilution, Amount and Internal standard amount entered in the batch reprocessing table.</p>
	<p><i>The recalculation is done automatically if the Reintegrate and/or Recalculate options are enabled. If the Recalculate only option is enabled, the Reintegrate and Recalculate options are disabled automatically.</i></p>
Change passport	<p>If this option is enabled, those parameters of the passport which have been changed after clicking the <Edit passport> button are applied to all chromatograms.</p>
<Edit passport>	<p>Open the Passport window for modification of the passport parameters.</p>
	<p><i>Only some of the passport parameters can be modified. The passport parameters Ident, Sample Info 1 and Sample Info 2 entered in the batch reprocessing table are overwritten if these values are modified in the Passport window.</i></p>
Modify chromatogram appearance	<p>If this option is enabled, the chromatogram Appearance parameters which have been changed after clicking the <Edit appearance> button are applied to all chromatograms.</p>
<Edit appearance>	<p>Open the Appearance window for modification of the settings for chromatogram axes, labels and colors.</p>
Make report	<p>Print report for all chromatograms using the current report settings of the selected chromatogram in the Use method... field.</p>
<Edit report options>	<p>Open the Report options window for modification of the report settings.</p>
<Reprocess>	<p>Start reprocessing.</p>
<Merge>	<p>Combine all chromatograms of the batch reprocessing queue into a single multi-channel chromatogram.</p>
<Close>	<p>Close the Reprocess window.</p>

Merge chromatograms

REPROCESS / <Merge>

Combine all chromatograms of the batch reprocessing queue into a single multi-channel chromatogram.



The chromatograms are displayed in the same order as in the batch reprocessing table slightly displaced one upon the other. The distance between the curves can be increased by pressing [Shift] + [↑] and decreased by pressing [Shift] + [↓].

The chromatogram axes, labels and colors can be set in the **Appearance** window.

The multi-channel chromatogram can be saved with **File / Save / Chromatogram**.

4.7.3 Batch reprocessing queue editor

Open batch reprocessing queue editor window

REPROCESS / <Edit sample table>

This menu item opens the batch reprocessing editor program for editing the batch reprocessing table.

	File Name	Method	Ident	Vial	Volume	Dilution	Amount	Internal Standard Amount	Calibration Level	Sample Info 1	Sample Info 2
1	k4101653.chw	Anion1.mtw	std 2	1	20.0	1.0	1.0	100.0	2	in ultra pure water	
2	k4110840.chw	Cation1.mtw	standard	1	10.0	1.0	1.0	100.0	1	in 2mmol/L HNO3	
3	k4111022.chw	Cation1.mtw	tap water	1	10.0	1.0	1.0	100.0	0	in 2mmol/L HNO3	
4	k4111745.chw	Anion1.mtw	tap water	3	20.0	1.0	1.0	100.0	0		

No	Row number.
File name	Name of the chromatogram file (read-only).
Method	Name of the method file *.mtw used for recording the chromatogram (read-only).
Ident	User defined identifier for chromatogram. Will be placed into appropriate passport field when starting the chromatogram.
Vial	Autosampler vial position to take sample from. Will be placed into appropriate passport field.
Volume	Injected volume of sample in μL . Will be placed into appropriate passport field.
Dilution	Sample dilution prior to injection. Will be placed into appropriate passport field.
Amount	Sample amount. If this value is different for the calibration run (c) and the sample run (s), the component concentrations of the sample are calculated as follows: $C_s = C_c \cdot \text{Amount}_s / \text{Amount}_c$ Will be placed into appropriate passport field.
Internal standard amount	Concentration of the internal standard component for relative concentration calculations. Will be placed into appropriate passport field.
Calibration level	Calibration level (see section 4.4.5) for the sample: Level 0 stays for normal analysis, levels 1 and greater for calibration runs. Correct filling of this column enables to calculate calibration coefficients automatically.

Sample Info 1	First sample description. Will be placed into appropriate passport field.
Sample Info 2	Second sample description. Will be placed into appropriate passport field.

Batch reprocessing queue editor functions

The functions in the editor window for the batch reprocessing queues can be triggered with the corresponding menu items of the **Edit** menu or with the corresponding symbols in the symbol bar.



QUEUE EDITOR / Edit / Undo

Undo the last modification of the batch reprocessing queue.



QUEUE EDITOR / Edit / Delete row(s)

Delete the selected rows from the batch reprocessing queue.



QUEUE EDITOR / Edit / Increment

Fill selected column fields with values automatically incremented by 1. The last character of the first selected field must be a number. This function is available for the fields **Ident**, **Vial**, **Calibration level**, **Sample Info 1** and **Sample Info 2**.



QUEUE EDITOR / Edit / Propagate

Fill selected column fields with the same value as the first field.



QUEUE EDITOR / Edit / Rotate rows

Rotate selected rows by one position (the bottom line is moved to the first position; all other lines are moved down by 1 position).

Print batch reprocessing queue



QUEUE EDITOR / File / Print

Print the batch reprocessing table in landscape format.

Close batch reprocessing queue editor



QUEUE EDITOR / File / Save & Exit

Save the modified batch reprocessing queue and close the **QUEUE EDITOR** window.

QUEUE EDITOR / File / Exit

Close the batch reprocessing **QUEUE EDITOR** window without saving the modified batch reprocessing table.

5 Notes – Maintenance – Faults

5.1 Practical notes on ion chromatography

5.1.1 Separating columns

Separation efficiency

The attainable quality of analyses with the 790 Personal IC depends to a large extent on the separation efficiency of the column used. When purchasing an IC column you should ensure that the separation efficiency suffices for the analysis problems at hand. Ascertain the **characteristic data of the IC column** on the standard chromatogram enclosed with the column such as capacity factors, selectivity, plate number and resolution and check these data with your own measurements. If any difficulties arise, you should always first check the quality of the column by recording a **standard chromatogram**.

You will find additional general tips on handling IC separating columns in the **8.732.2003 Metrohm Monograph "Ion chromatography"**, as well as detailed information on the separating columns available from Metrohm (see *section 6.3.2*) in the leaflets supplied and in the special **Application Bulletins** available free of charge from your nearest Metrohm agency.

Protection

To protect the column against foreign particles which could have an adverse influence on the separation efficiency, we advise you to subject both the eluents and all samples to **microfiltration** (0.45 µm filter) and to siphon the eluent through the **6.2821.090 Aspirating filter**.

To avoid contamination by abrasive particles arising from piston seals of the high-pressure pump, it is advantageous to install an **in-line filter** between the pump and the injection valve. In the 790 Personal IC a **6.2821.100 Filter unit PEEK** is already mounted for this purpose (see *section 2.6.2*).

The use of readily interchangeable **precolumns** serves to protect the actual separating columns and increase their service life appreciably. The precolumns available from Metrohm (see *section 6.3.2*) are either actual precolumns or so-called precolumn cartridges that are used together with the 6.2821.040 Cartridge head or the 6.2828.010 Precolumn cartridge holder (see *section 2.7*).

Storage

Always store the separating columns closed when not in use and filled in accordance with the manufacturer's specifications.

Dead volume

Dead volume at the end of a column can be the cause of extreme peak broadening or splitting (appearance of double peaks). Filling the column with glass beads ($\varnothing \leq 100 \mu\text{m}$) frequently improves the separation efficiency.

Regeneration

If the separation properties of the column have deteriorated, it can be regenerated in accordance with the column manufacturer's specifications. With the separating columns available from Metrohm (see section 6.3.2), the instructions for regeneration can be found on the leaflet enclosed with every column.



*In the case of separating columns with carrier material based on silica, **only solutions with pH 2...7** may be used for regeneration, otherwise the columns could be damaged.*

5.1.2 High-pressure pump

Pulsation dampener

The **6.2620.150 Pulsation dampener MF** belongs to the standard accessories of the 790 Personal IC; its installation is described in section 2.6.2. It is used to reduce interfering pulsation in highly sensitive measurements and also protects the column material against pressure shocks caused by the injection.

Maintenance

To protect the pump against foreign particles, we advise you to subject the eluent to **microfiltration** (0.45 μm filter) and siphon the eluent through the **6.2821.090 Aspirating filter**.

In many cases, an unstable baseline (pulsation, flow fluctuations) can be traced to contaminated valves or faulty, leaky piston seals.

Contaminated valves are cleaned by rinsing with water, RBS solution or acetone (see section 5.2.6). When the cleaned valves are reinstalled, you must ensure that the flow direction is correct.

The **replacement of piston seals** is described in section 5.2.6.

Salt crystals between the piston and the seal are the cause of abrasive particles, which can enter the eluent. These lead to contaminated valves, pressure rise and in extreme cases to scratched pistons. It is thus essential to ensure that **no precipitates** can appear (see also section 5.1.3).

5.1.3 Eluents

Treatment

For the preparation of the eluents, only chemicals of a purity degree of at least "p.a." should be used. For diluting, please use only **high purity water**.

Fresh eluents should always be **microfiltered** (0.45 µm filter) and **de-gassed** (with N₂, He or vacuum). For alkaline eluents and eluents with low buffering capacity one should preferably use a **CO₂ absorber**.

The supply vessel containing the eluent must be closed as tightly as possible to avoid excessive evaporation. This is primarily important with eluents containing organic solvents (e.g. acetone), the evaporation of which can lead to drifts in the long term. If work is performed in a very sensitive range, even if one drop of condensate falls back in the eluent this can cause a noticeable change in the background conductivity.

Influence of various parameters on anion columns

- *Concentration:* An increase in the concentration usually leads to shorter retention times and quicker separation, but also to a higher background conductivity.
- *pH:* pH alterations lead to shifts in the dissociation equilibrium and thus to changes in the retention times.
- *Organic modifiers:* Addition of an organic solvent (e.g. methanol, acetone, acetonitrile) to aqueous eluents generally accelerates lipophilic ions.

Eluent change

When the eluent is changed, it must be ensured that **no precipitates** can be formed. Solutions used in direct succession must therefore be miscible. If the system has to be rinsed with an organic solution, several solvents with increasing or decreasing lipophilic character may possibly have to be used (e.g. water ↔ acetone ↔ chloroform).

5.1.4 Peristaltic pump

The 6.1826.060 pump tubings used by the peristaltic pump are consumable material with a limited lifetime and should be exchanged at regular intervals (approx. every 4 weeks under continuous use; see *section 5.2.10*).

The working life of pump tubing depends to a considerable extent on the contact pressure. This is why the contact pressure must be correctly set as described in *section 2.9.2*. If the pump is to remain switched off for a lengthy period of time the tubing cartridges **40** should be raised completely by loosening the snap-action lever **43** on the right-hand side (the set contact pressure remains unchanged).

5.1.5 Suppressor module

Protection

To avoid contamination of the suppressor module by foreign particles or bacterial growth, the two **6.2821.100 Filter units PEEK** (see section 2.3.5) supplied must be installed between peristaltic pump and inlet capillaries of the suppressor module (see section 2.8.2).

Operation

The **Metrohm Suppressor Module MSM** comprises a total of 3 suppressor units which are in turn used for suppression, regenerated with sulfuric acid and rinsed with water. To record every new chromatogram under comparable conditions, work is normally carried out with freshly regenerated suppressor. Switching takes place automatically together with the valve switching.

For operation of the suppressor module, a **two-channel peristaltic pump** is needed which conveys the regeneration solution (normally **20 mmol/L H₂SO₄**) and the rinsing solution (normally **dist. H₂O**) to the suppressor units (flow rate of 0.5 mL/min).



The suppressor units must never be regenerated with H₂SO₄ in the same flow direction used for the eluent. You should thus always install the inlet and outlet capillaries as described in section 2.8.4 according to the scheme shown in Fig. 16.



The suppressor module must never be switched in the dry state as there is a danger of blocking. Before every switching operation of the suppressor module, the three suppressor units must have been rinsed for at least ½ h with eluent, regeneration and rinsing solution.

The capacity of the suppressor units is exhausted after approx. 2 h. With longer interruptions between the individual measurements, it is recommended that the **Prep MSM.smt** system is started during the pauses; this automatically switches the suppressor module further every 20 min.

Maintenance

In the event of reduced capacity or high counter-pressure the suppressor module must be regenerated (section 5.2.7), cleaned (section 5.2.8) or exchanged (section 5.2.9).

5.1.6 Connections

All connections between injector, column and detector must be as short as possible, have a low dead volume and be absolutely tight. The PEEK capillary after the detector block must be free from constriction (the measuring cell is tested to 5 MPa = 50 bar back pressure).

5.2 Maintenance and servicing

5.2.1 General information

Care

The 790 Personal IC requires proper care and attention. Excessive contamination of the instrument could possibly lead to malfunctions and a shorter service life of the inherently rugged mechanical and electronic parts.

Spilled chemicals and solvents should be wiped up immediately. It is especially important to protect the plug connections at the rear of the instrument (particular the mains plug) against contamination.



Although constructional measures have been designed to virtually eliminate such a situation, should corrosive media penetrate the interior of the instruments the mains plug of the 790 Personal IC must be immediately disconnected to prevent extensive damage to the instrument electronics. Inform Metrohm service if your instrument(s) have been damaged in such a way.



The instrument must not be opened by untrained personnel. Please comply with the safety notes in section 1.4.1.

Maintenance by Metrohm service

Maintenance of the 790 Personal IC is best done as part of an annual service performed by specialists from the Metrohm company. If work is frequently performed with caustic and corrosive chemicals, it may be necessary to shorten the interval between servicing.

The Metrohm service department is always willing to offer expert advice on the maintenance and servicing of all Metrohm instruments.

5.2.2 Passivation

Passivation of the complete IC System (without column) by rinsing with 20...50 mL 0.2 mol/L HNO₃ is only required when unusual alterations to the measuring properties of the cell are observed. For passivation, the separating column **63** is removed from the 790 Personal IC. The two capillaries **20** and **37** (see Fig. 13 and Fig. 15) are directly connected to each other using the Coupling **26** (6.2620.060).

5.2.3 Recycling

To keep the eluent consumption between injections to a minimum when the system is at rest (e.g. overnight), the so-called recycling procedure can be used. In recycling the eluent exiting the outlet capillary of the detector block is led back directly to the eluent supply vessel. The IC system is thus quickly ready for new injections without a long conditioning period.



The recycling procedure must **not** be used

- in operation with the suppressor module,
- with alkaline eluents,
- with the 6.1010.000 IC Cation column METROSEP Cation 1-2.

5.2.4 Shutdown

If the 790 Personal IC is shut down for a considerable length of time, the entire IC system (**without** column and suppressor) must be **rinsed free from salt** with methanol/water (1:4) to avoid crystallization of eluent salts with the corresponding subsequent damage.

For rinsing, the connections to the separating column and the suppressor module are removed; the two capillaries **20** and **37** (see *Fig. 13* and *Fig. 15*) are directly connected to each other using coupling **26** (6.2620.060). Rinse with methanol/water (1:4) until the conductivity drops below 10 $\mu\text{S}/\text{cm}$.

5.2.5 Changing separating columns

Identical separation system

If you wish to replace an IC separating column by a column of the same type, proceed as follows (see *Fig. 13* and *Fig. 15*):

1 Remove old column

- Switch off high-pressure pump and wait for pressure to drop.
- Unscrew column **63** from inlet capillary **37** of the detector block or from suppressor inlet capillary **78**.
- Unscrew column **63** from column connection capillary **20** or the precolumn.

2 Connect new column to injector

- Remove end caps from column **63**.
- Screw inlet end of separating column **63** (note flow direction) to column connection capillary **20** or to the precolumn (see *section 2.7.7/section 2.7.8*).

3 Rinse column

- Place beaker beneath the column outlet.
- Switch on high-pressure pump and rinse column with eluent for ca. 10 min, then switch off pump.

4 Connect column to detector block

- Screw outlet end of separating column **63** to inlet capillary **37** or suppressor inlet capillary **78**.

Changing the separation system

If you wish to replace an IC separating column by a column of a different type, proceed as follows (see *Fig. 13* and *Fig. 15*):

1 Remove old column

- Switch off high-pressure pump and wait for pressure to drop.
- Unscrew column **63** from inlet capillary **37** of the detector block or from the suppressor inlet capillary **78**.
- Unscrew column **63** from column connection capillary **20** or the precolumn.

2 Rinse with eluent

- Place beaker beneath the column connection capillary **20**.
- Rinse IC system with eluent used for the separating column (flow rate 1 mL/min) for ca. 15 min.

3 Connect new column to injector

- Remove end caps from column **63**.
- Screw inlet end of separating column **63** (note flow direction) to column connection capillary **20** or to the precolumn (see *section 2.7.7/section 2.7.8*).

4 Rinse column

- Place beaker beneath the column outlet.
- Set flow rate for new separating column.
- Switch on high-pressure pump and rinse column with eluent for ca. 10 min, then switch off pump.

5 Connect column to detector block

- Screw outlet end of separating column **63** to inlet capillary **37** or suppressor inlet capillary **78**.

5.2.6 Maintenance work at the pump head

In many cases, an unstable baseline (pulsation, flow fluctuations) can be traced to contaminated valves or faulty, leaky piston seals at the high-pressure pump. For cleaning contaminated valves and/or replacement of wear parts such as pistons, piston seals and valves, proceed as follows:

1 Detach pump head

- Disconnect aspirating tubing from aspirating capillary **32** at the pump head **34** (see *Fig. 3* and *Fig. 4*).
- Unscrew connection capillary **35** from the pump head **34**.
- Remove pump head **34** by loosening the 4 hexagon screws **33** using the 6.2621.030 Hexagon key. The main piston is on the left (when viewed from front), the auxiliary piston on the right.

2 Disassemble pump head

- Strip down pump head **34** in accordance with *Fig. 17*. Main and auxiliary pistons are identical with the following exceptions:
 - The spring **87** of the auxiliary piston (right piston) is more powerful (longer) than that of the main piston (left piston).
 - Inlet and outlet valve are not present in the secondary cylinder.



*To prevent the piston **85** suddenly jumping out of the piston cartridge **87**, the screw **84** must be undone very carefully by hand.*

3 Cleaning/replacement of piston **85**

- Pistons contaminated by abrasive particles or deposits are cleaned with scouring powder and rinsed free of any particles with dist. water.
- Relatively badly contaminated or scratched pistons must be replaced (spare part: 6.2824.070 Zircon piston).

4 Replacement of piston seal **92**

- To remove damaged piston seals **92** the special tool **96** is used. This is screwed into the seal **92**, which can then be pulled out (see *Fig. 18A*).



*When the special tool **96** is screwed into the piston seal **92** the latter is completely destroyed!*



If you use only aqueous eluents, the 6.2741.000 Piston seal can be replaced by the 6.2741.010 PE Piston seal available as an option.




Fig. 17: Components of the pump head




Fig. 18: Replacement of the piston seal 92

34	Pump head (6.2824.100)	91	Piston guide sleeve (4.709.4370)
84	Screw (3.709.1100) for piston cartridge 88	92	Piston seal (6.2741.000/6.2741.010)
85	Zircon piston (6.2824.070) with piston shaft	93	Inlet valve (6.2824.090)
86	Spring retainer 4.709.0730	94	Outlet valve (6.2824.080)
87	Spring (6.2824.050) for main piston or Spring (6.2824.060) for auxiliary piston	95	Screw holder for valve
88	Piston cartridge (4.709.0760)	96	Special tool (6.2617.010) to remove the piston seal 92
89	Piston guide sleeve (4.709.4380)	97	Special tool (6.2617.010) to install the piston seal 92
90	Sapphire supporting ring (6.2824.030)		

- To install a new piston seal **92** the special tool **97** is used.
- First the new seal is inserted firmly in the recess of tool **97** by hand (see *Fig. 18B*). The seal spring must be located on the outside.
- The tool **97** together with the seal is then inserted in the pump head **34** and the seal pressed into the pump head recess with the aid of tool **96** (see *Fig. 18C*).



*The seal surface in the pump head **34** must not be damaged (avoid contact with tool)!*

5 Cleaning/replacement of inlet valve **93** and outlet valve **94**

- Contaminated or blocked valves are cleaned by rinsing with dist. water, RBS solution or acetone. The rinsing effect can be reinforced by brief treatment in an ultrasonic bath (max. 20 s; if longer the sapphire sphere of the valve can be damaged).
- If this does not have the desired effect, the valves can be disassembled as shown in *Fig. 19*. The valve components are pushed out with the aid of a syringe needle inserted through the upper opening in the valve housing **98**. The individual components are rinsed with dist. water and/or acetone, and the sapphire sphere cleaned with a paper towel. The valve is then reassembled in accordance with *Fig. 19*. The components of the inlet and outlet valves are identical, they are distinguished only by the positioning of the sapphire sleeve **101** and the ceramic holder **103** (see *Fig. 19*).
- Valves that fail to function faultlessly after such cleaning must be replaced.
- In the reinstallation of the inlet valve **93** or the outlet valve **94** on no account must the two outwardly identical valves be interchanged. To determine which valve is which, note that the liquid flows through the pump head from the bottom up. The flow direction of the valves can be checked simply by blowing through the clean valve. Both valves are installed with the black face in the direction of the pump head (see *Fig. 19*).



*If by mistake an inlet valve **93** is installed instead of the outlet valve **94**, an extreme pressure buildup occurs within the working cylinder, which is not detected by the pressure transducer and will destroy the piston seal **92**!*

5.2.7 Regeneration of the suppressor module

Regenerating a suppressor operating at reduced capacity

If the suppressor units are exposed to certain heavy metals (e.g. iron) or organic contaminants for long periods, these can no longer be completely removed by the regeneration solution normally used (20 mmol/L H_2SO_4). This diminishes the capacity of the suppressor units, which, in milder cases, results in a reduced sensitivity to phosphates and, in severe cases, in a strong increase in the baseline. If such capacity problems occur at one or several positions, the suppressor units must be treated as follows:

1 Disconnect suppressor from IC system

- Disconnect suppressor from separating column and detector.

2 Regenerate suppressor

- Rinse each suppressor unit for about 15 min with one of the following solutions:

Contamination with heavy metals

1 mol/L H_2SO_4

Contamination with organic cationic complexing agents

0.1 mol/L H_2SO_4 / 0.1 mol/L oxalic acid / acetone 5%

Severe contamination with organic substances

0.2 mol/L H_2SO_4 / acetone \geq 20%



The 6.1826.060 pump tubing is made of PP and must not be used for rinsing with solutions that contain organic solvents. In such cases, rinse with different pump tubing or a different pump.

3 Connect suppressor to IC system

- Reconnect suppressor to the IC system. If capacity problems persist, replace the suppressor rotor (see section 5.2.9).

Regenerating the suppressor in case of high counter pressure

If excessive counter pressure is observed in one or several suppressor units, treat the units as follows:

1 Disconnect suppressor from IC system

- Disconnect suppressor from separating column and detector.

2 Regenerate suppressor

- Connect the inlet capillary **81** marked " H_2SO_4 " to the inlet capillary **18** using coupling **26** (see Fig. 10 and Fig. 15). This connects the suppressor module directly to the high-pressure

pump.

- Set the flow at the high-pressure pump to 0.5 mL/min and rinse the suppressor unit with 1 mol/L H₂SO₄ for 5 to 10 min.
- As the pressure falls, slowly increase the flow at the high-pressure pump to 2 mL/min. Do not exceed a maximum pressure of 2 MPa (20 bar).
- Switch off high-pressure pump.
- Switch suppressor to the next position using the **Step** button.
- Connect the inlet capillary **80** marked "H₂O" to the inlet capillary **18** using coupling **26** (see *Fig. 10* and *Fig. 15*).
- Set the flow at the high-pressure pump to 0.5 mL/min and rinse the suppressor unit with 1 mol/L H₂SO₄ for 5 to 10 min.
- As the pressure falls, slowly increase the flow at the high-pressure pump to 2 mL/min. Do not exceed a maximum pressure of 2 MPa (20 bar).
- Switch off high-pressure pump.

3 Connect suppressor to IC system

- Connect inlet capillaries **80** and **81** to the peristaltic pump (see *section 2.8.4*).
- If the pressure problems persist, replace the suppressor rotor (see *section 5.2.9*).

5.2.8 Cleaning the suppressor

It may be necessary to clean the suppressor in the following cases:

- High counter pressure on the suppressor connection tubing
- Irremediable blockage of the suppressor (the suppressor can no longer deliver solutions)
- Irremediable obstruction of the suppressor (the suppressor can no longer be switched to next position)

To clean the connection piece and the rotor, proceed as follows (see *Fig. 20*):

1 Disconnect suppressor from IC system

- Disconnect input capillary **78** of the suppressor module **39** from the separating column **63** (see *Fig. 15*).
- Disconnect output capillary **79** from inlet capillary **37**.
- Disconnect inlet capillaries **80** and **81** from the filter units **28** (supply from peristaltic pump).

2 Dismantle suppressor

- Unscrew nut **105** from suppressor holder **108**.
- Pull out connection piece **106** and suppressor rotor **107** from suppressor holder **108** (the connection piece and the rotor normally stick together).
- Loosen connection piece **106** from suppressor rotor **107**.

3 Clean input and output leads

- Connect each of the 6 capillary tubings attached to connection piece **106** to the high-pressure pump one after another, and pump through ultrapure water.
- Check whether solution emerges from connection piece **106**. If one of the input or output leads remains blocked, replace the connection piece **106** (order number 6.2832.010).

4 Clean suppressor rotor

- Clean the sealing surface of suppressor rotor **107** using a lint-free cloth and ethanol.

5 Insert suppressor rotor

- Insert suppressor rotor **107** in suppressor holder **108** in such a way that the tubing connections at the rear of the rotor fit in the corresponding openings inside the rotor, and that one of the three holes in the rotor can be seen from below in one of the openings of the holder.
- If the rotor has been inserted correctly, its sealing surface will be about 4 mm inside the holder. If this is not the case, bring the rotor into the correct position from below with the aid of a sharp object (e.g. a screwdriver).

6 Clean connection piece

- Clean the sealing surface of connection piece **106** using a lint-free cloth and ethanol.

7 Insert connection piece

- Insert connection piece **106** in suppressor holder **108** in such a way that connection "1" is at the top, and that the three lugs on the connection piece fit in the corresponding openings of the holder.
- Screw nut **105** onto the thread of suppressor holder **108** manually (do not use tools).

8 Connect and condition the suppressor

- Reconnect the suppressor to the IC system.
- Before switching the suppressor to the next position for the first time, rinse all 3 suppressor units with solution for 5 min.




Fig. 20: Assembling the suppressor

105	Screw nut	107	Suppressor rotor (6.2832.000)
106	Connection piece (6.2832.010) with input and output leads	108	Suppressor holder

5.2.9 Replacing the suppressor

The suppressor in the suppressor block may have to be replaced in the following cases:

- Irremediable loss of suppressor capacity (reduced phosphate sensitivity and/or strong rise in baseline)
- Irremediable blockage of the suppressor (the suppressor can no longer deliver solutions)

Both the 6.2832.000 Suppressor rotor and the 6.2832.010 Connection piece with the input and output leads can be replaced. To replace these components proceed as follows (see Fig. 20):

1 Disconnect suppressor from IC system

- Disconnect all input and output leads of the suppressor from IC system and peristaltic pump.

2 Dismantle suppressor

- Unscrew nut **105** from suppressor holder **108**.
- Pull out connection piece **106** and suppressor rotor **107** from suppressor holder **108** (the connection piece and the rotor normally stick together).
- Loosen connection piece **106** from suppressor rotor **107**.

3 Clean suppressor rotor

- Clean the sealing surface of new suppressor rotor **107** (6.2832.000) using a lint-free cloth and ethanol.

4 Insert suppressor rotor

- Insert new suppressor rotor **107** in suppressor holder **108** in such a way that the tubing connections at the rear of the rotor fit in the corresponding openings inside the rotor and that one of the three holes in the rotor can be seen from below in one of the openings of the holder.
- If the rotor has been inserted correctly, its sealing surface will be about 4 mm inside the holder. If this is not the case, bring the rotor into the correct position from below with the aid of a sharp object (e.g. a screwdriver).

5 Clean connection piece

- Clean the sealing surface of new connection piece **106** (6.2832.010) with the aid of a lint-free cloth and ethanol.

6 Insert connection piece

- Insert new connection piece **106** in suppressor holder **108** in such a way that connection "1" is at the top and that the three lugs on the connection piece fit in the corresponding openings of the holder.
- Screw nut **105** onto the thread of suppressor holder **108** manually (do not use tools).

7 Connect and condition the suppressor

- Reconnect the suppressor to the IC system.
- Before switching the suppressor to the next position for the first time, rinse all 3 suppressor units with solution for 5 min.

5.2.10 Exchanging the pump tubing

The pump tubings used by the peristaltic pump are consumable material with a limited lifetime and should be exchanged at regular intervals (approx. every 4 weeks under continuous use).

The working life of pump tubing depends to a considerable extent on the contact pressure. This is why the contact pressure must be correctly set as described in *section 2.9.2*. If the pump is to remain switched off for a lengthy period of time the tubing cartridges **40** should be raised completely by loosening the snap-action lever **43** on the right-hand side (the set contact pressure remains unchanged).

As the pump is always operated on the same side, the 6.1826.060 Pump tubings supplied can be used on both sides. To exchange a pump tubing proceed as follows:

1 Remove old pump tubing

- Press contact pressure lever **41** on the tubing cartridge **40** down as far as it will go.
- Release tubing cartridge **40** from holding clamp **42** by pressing down snap-action lever **43** and remove from mounting pin **45** (see *Fig. 15*).
- Remove old pump tubing **74** or **75**.

2 Insert new pump tubing

- Insert the new pump tubing **74** or **75** in the tubing cartridge **40** as shown in *Fig. 14*. The white-yellow stopper **76** must click into the corresponding holder on the left-hand side of the tubing cartridge.
- Place the tubing cartridge **40** on mounting pin **45** and press down on the right-hand side until snap-action lever **43** clicks into position on holding clamp **42**. Take care that no kinks are formed in the pump tubing.

3 Set contact pressure

- Switch on peristaltic pump.
- Press contact pressure lever **41** upwards until the solution just starts to be drawn in. Then press contact pressure lever upwards until it clicks once more to obtain optimal contact pressure.
- Switch off peristaltic pump.

5.3 Faults and malfunctions

5.3.1 Error messages

If any type of malfunction occurs during operation of the 790 Personal IC, this is shown by error messages in the PC program, which appear either in an **error window** or in the **SYSTEM STATE** window.

Follow the instructions given in the **error window** and close this window with <OK>.

You will find further details of the error messages of the **SYSTEM STATE** window, their possible causes and the procedure for their rectification in *section 4.3.7*.

5.3.2 Malfunctions and their rectification

If difficulties appear with the 790 Personal IC during analyses, their causes are best investigated in the order **separating column → high-pressure pump → eluent → connections**. Several of the malfunctions which may appear are listed in the following table with details of possible causes and countermeasures.

Malfunction	Cause	Rectification
Baseline with high noise level, pulsation	<ul style="list-style-type: none"> Contaminated pump valves Faulty piston seals 	<ul style="list-style-type: none"> Clean the valves (see <i>section 5.1.5</i>) Replace the piston seals (see <i>section 5.1.5</i>)
Drift of the baseline	<ul style="list-style-type: none"> Thermal equilibrium not yet reached Leak in system Evaporation of organic solvent in eluent 	<ul style="list-style-type: none"> Condition system with heating switched on Check connections and make leakproof Ensure better closure of eluent supply vessel
Considerable pressure drop	<ul style="list-style-type: none"> Leak in system 	<ul style="list-style-type: none"> Check connections and make leakproof
Considerable pressure rise	<ul style="list-style-type: none"> Contamination of the filter in the 6.2821.100 Filter unit PEEK Contamination of the column inlet filter Change of column packing by injection of contaminated samples 	<ul style="list-style-type: none"> Replace the 6.2821.110 Filter (see <i>section 2.3.5</i>) Clean or replace 6.2821.020 Steel mesh(es) Regenerate column (see <i>section 5.1.1</i>) or replace column <p><i>Note:</i> Samples should always be microfiltered.</p>

Malfunction	Cause	Rectification
Chromatograms with poor resolution, change in the retention times	<ul style="list-style-type: none"> • Deterioration in separation efficiency of the IC column 	<ul style="list-style-type: none"> • Regenerate column (see <i>section 5.1.1</i>) or replace column
Extreme peak broadening, splitting (double peaks)	<ul style="list-style-type: none"> • Dead volume at the column ends 	<ul style="list-style-type: none"> • Fill column with glass beads ($\varnothing \leq 100 \mu\text{m}$) or replace column
No feed of regeneration or rinsing solution for suppressor	<ul style="list-style-type: none"> • Contact pressure too low • Leak in system • Defective pump tubing • Contamination of the filter in the 6.2821.100 Filter unit PEEK • Counterpressure at suppressor module too high 	<ul style="list-style-type: none"> • Adjust contact pressure (see <i>section 2.9.2</i>) • Check connections and make leakproof • Replace pump tubing (see <i>section 5.2.9</i>) • Replace the 6.2821.110 Filter (see <i>section 2.3.5</i>) • Clean or replace suppressor (see <i>section 5.2.6...5.2.8</i>)

5.4 Diagnostic tests / Validation / GLP

The requirements of **GLP** (**Good Laboratory Practice**) include a periodic check of analytical measuring instruments with regard to their reproducibility and accuracy using **Standard Operating Procedures, SOP**. Under the title «**Application Bulletin No. 277 – Validation of Metrohm Ion Chromatography Systems by using Standard Operating Procedures (SOP)**» an example of such a standard operating procedure is available from Metrohm; it can be adapted and used with the 790 Personal IC.

Further information on the subjects of QA, GLP and validation can also be found in the brochure «**Quality management with Metrohm**», which is available from your local Metrohm agency.

Testing of the electronic and mechanical function groups of Metrohm instruments can and should be performed as part of a regular service by trained personnel of the manufacturing company (see *section 5.2.1*). All Metrohm instruments are equipped with start-up-test routines that check for perfect functioning of the relevant assemblies when the instrument is switched on. If no error message is displayed, it may be assumed the instrument is operating without faults.

The Metrohm company also supplies its instruments with an integrated diagnostic program which, in the case of possible malfunctions or faulty behavior, allows the service technician to check the functioning of certain assemblies and localize the fault.

6 Appendix


6.1 Technical data

6.1.1 Conductivity measurement

<i>Measurement range</i>	0...1000 $\mu\text{S/cm}$ (resolution: 0.56 nS/cm)
<i>Maximum error</i>	$\pm 10 \mu\text{S/cm}$ and $\pm 1 \%$ of measurement value ($k = 16.7/\text{cm}$)
<i>Linearity</i>	Deviations $< \pm 5 \mu\text{S/cm}$
<i>Noise</i>	typ. 3 nS/cm
<i>Drift (electronic)</i>	typ. $< 10 \text{ nS/cm / h}$
<i>Temperature dependence</i>	typ. $< 40 \text{ nS/cm / }^\circ\text{C}$
<i>Reserve range</i>	$> 33 \%$ ($k = 16.7/\text{cm}$)
<i>Sampling rate</i>	10 measurements/s (fixed)

6.1.2 Conductivity detector

<i>Construction</i>	Thermostatted conductivity detector with 2 ring-shaped steel electrodes
<i>Measurement principle</i>	Alternating current measurement with 1 kHz frequency and ca. 1.7 V amplitude (peak to peak).
<i>Cell volume</i>	1.5 μL
<i>Cell constant</i>	approx. 17 /cm (the exact value is printed on the detector)
<i>Maximum back pressure for measuring cell</i>	5.0 MPa (50 bar)
<i>Thermostatting</i>	Connectable dynamic control to operating temperature
<i>Operating temperature</i>	40 $^\circ\text{C}$
<i>Max. temperature deviation</i>	$\pm 2.5^\circ\text{C}$
<i>Heating time</i>	$\geq 30 \text{ min}$
<i>Temperature stability</i>	$\leq 0.01^\circ\text{C}$ at constant ambient temperature
<i>Connection for detector block</i>	Dsub 15 pin (female)



6.1.3 Injection valve

<i>Actuator switching duration</i>	100...150 ms
<i>Pressure resistance</i>	25 MPa (250 bar)

6.1.4 High-pressure pump

<i>Type</i>	Serial dual piston pump with two valves
<i>Pump capacity</i>	
<i>Flow range</i>	0.20...2.5 mL/min
<i>Maximum error</i>	< ± 2 % of set value
<i>Flow constancy</i>	< 0.5 % of set value
<i>Reproducibility of eluent flow</i>	typ. better than ± 0.1 %
<i>Pressure measurement</i>	
<i>Pressure range</i>	0...25.0 MPa (0...250 bar)
<i>Residual pulsation</i>	< 10 % (at 1 mL/min water and 10 MPa pressure, without pulsation dampener)
<i>Measurement principle</i>	Piezoresistive measurement Response time: 3 ms Measurement volume: ca. 50 µL
<i>Maximum error</i>	± 3 % of set value
<i>Resolution</i>	0.1 MPa (conductivity measurements)
<i>Sampling rate</i>	1 measurement/piston stroke (pump running) 1 measurement/s (pump not running)
<i>Safety shutdown</i>	
<i>Function</i>	Automatic shutdown when upper and lower pressure limits violated
<i>Maximum pressure limit</i>	adjustable between 0.1...25.0 MPa (1...250 bar) Response time: 1 pump cycle
<i>Minimum pressure limit</i>	adjustable between 0.1 ... 25.0 MPa (1...250 bar), inactive at 0 MPa Response time: 5 pump cycles
<i>Pump head</i>	
<i>Pump head volumes</i>	Main piston: 40 µL Priming piston: 20 µL
<i>Pump displacement volumes</i>	Main piston: 28.5 µL Priming piston: 14.25 µL
<i>Length of stroke</i>	Main piston: 3.6 mm Priming piston: 1.8 mm

6.1.5 Peristaltic pump


<i>Type</i>	2-channel peristaltic pump
<i>Pump capacity</i>	
<i>Rotational speed</i>	20 U/min at 50 Hz 24 U/min at 60 Hz
<i>Flow range</i>	0.5...0.6 mL/min with 6.1826.060 Pump tubing
<i>Maximum error</i>	± 5 %
<i>Maximum pressure</i>	0.4 MPa (4 bar)
<i>Pumpable liquids</i>	Clear liquids with no solid contents
<i>Pump tubing material</i>	PP (polypropylene)

6.1.6 Suppressor module

<i>Switching duration</i>	140 ms
<i>Pressure resistance</i>	2.5 MPa (25 bar)

6.1.7 RS232 interface


Connector Dsub 9 pin (male)



Function TxD and RxD signal for connection with software handshake


Default settings 9600 baud, 8 bit, 1 stop bit, no parity, XON/XOFF

Pin assignment




6.1.8 Remote interface


Connector Dsub 25 pin (female)



Circuit diagram for output line 3



Potentials



6.1.9 Mains connection

<i>Voltage</i>	115 V: 100...120 V \pm 10 % 230 V: 220...240 V \pm 10 %
<i>Frequency</i>	50...60 Hz
<i>Power consumption</i>	100 VA
<i>Fuse</i>	5 mm dia., 20 mm length 100...120 V: 1.0 A (slow-blow) 220...240 V: 0.5 A (slow-blow)

6.1.10 Safety specifications

<i>Construction/testing</i>	According to IEC 1010 / EN 61010 / UL 3101-1, protection class 1, degree of protection IP20
<i>Safety directions</i>	The Instructions for Use include information and warnings that must be heeded by the user to assure safe operation of the instrument.

6.1.11 Electromagnetic compatibility (EMC)

<i>Emitted interference</i>	Standards met: EN55011 (class B), EN55022 (class B), EN50081-1, IEC61326 (class B)
<i>Immunity to interference</i>	Standards met: IEC801-2/IEC1000-4-2 (class 3), IEC801-3/ IEC1000-4-3 (class 2), IEC801-4/IEC1000-4-4 (class 4), IEC801-5/IEC1000-4-5 (class 2/3), IEC801-6/IEC1000-4-6 (class 2), EN50082-1, EN61000-3-2/IEC1000-3-2, EN61000-3-3/ IEC1000-3-3, EN61000-4-11/IEC1000-4-11, IEC61326

6.1.12 Ambient temperature

<i>Nominal operating range</i>	+5...+45°C (at 20...80 % atmospheric humidity)
<i>Storage</i>	-20...+70°C
<i>Transport</i>	-40...+70°C

6.1.13 Housing

<i>Material of cover</i>	Polyurethane rigid foam (PUR) with fire protection for fire class UL94VO, CFC-free
<i>Material of base</i>	Steel, enameled
<i>Width</i>	255 mm
<i>Height</i>	385 mm
<i>Depth</i>	343 mm
<i>Weight (incl. accessories)</i>	2.790.0010: 14.7 kg 2.790.0020: 16.0 kg

6.2 Standard equipment
















*Subject to changes !
All dimensions are given in mm.*




The 790 Personal IC is available in two versions:

- **2.790.0010** Personal IC without suppressor module
- **2.790.0020** Personal IC with suppressor module

These instruments include the following parts:

Quant.		Order No.	Description													
2.790.0010	2.790.0020															
1	1	1.732.0110	Detector block (metal-free) with permanently attached connecting cable to 790 Personal IC													
-	1	6.1803.020	PTFE capillary Length = 5 m													
-	2	6.1826.060	Pump tubing made of PP (polypropylene) with 2 firmly attached white-yellow stoppers; in.d. = 0.51 mm, ex.d. = 2.31 mm													
1	1	6.1834.010	Aspirating tubing made of PTFE, with connector for 6.2821.090 aspirating filter Length = 2.5 m For the connection high-pressure pump – eluent bottle													
1	1	6.2027.030	Column holder Diameter d = 8.5 mm													
1	1	6.2027.040	Column holder Diameter d = 11.3 mm													
1	1	6.2122.0X0	Mains cable to customer's specifications: <table style="width: 100%; border: none;"> <tr> <td style="width: 50%;"><u>Cable socket</u></td> <td style="width: 50%;"><u>Cable connector</u></td> <td></td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type SEV 12 (CH...)</td> <td>6.2122.020</td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type CEE (7), VII (D...)</td> <td>6.2122.040</td> </tr> <tr> <td>Type CEE (22), V</td> <td>Type NEMA 5-15 (USA...)</td> <td>6.2122.070</td> </tr> </table>	<u>Cable socket</u>	<u>Cable connector</u>		Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020	Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040	Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070	
<u>Cable socket</u>	<u>Cable connector</u>															
Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020														
Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040														
Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070														







Quant.		Order No.	Description
2.790.0010	2.790.0020		
1	1	6.2134.100	<p>Connecting cable Connecting cable 790 Personal IC (RS232) – PC</p> 
1	1	6.2617.010	<p>Special tool For removing the piston seal of the pump head</p> 
1	1	6.2621.030	<p>Hexagon key 4 mm For mounting the pump head of the high-pressure pump.</p> 
1	3	6.2744.010	<p>PEEK compression fitting For the connection of 6.1831.010 PEEK capillaries or 6.1822.010 PTFE microcapillaries, set of 5</p> 
-	1	6.2744.030	<p>PEEK coupling Connection between 6.2744.010 PEEK compression fitting and 6.1826.060 Pump tubing; Set of 4</p> 
-	1	6.2744.040	<p>PEEK coupling For connection of 1/16" capillaries</p> 
-	2	6.2744.110	<p>PEEK coupling Connection between 6.2821.100 Filter unit PEEK and 6.1826.060 Pump tubing</p> 
1	1	6.2744.120	<p>Coupling 1/16" – Luer Connection between 6.2744.010 PEEK compression fitting and 6.2816.020 Syringe</p> 

Quant.		Order No.	Description
2.790.0010	2.790.0020		
1	1	6.2816.020	Syringe made of PP, volume = 10 mL; for manual filling of the sample loop 
1	1	6.2821.090	Aspirating filter Pore dimension 20 µm For 6.1834.000 Aspirating tubing. Set of 5 
-	2	6.2821.100	Filter unit PEEK 2 µm To avoid contamination due to abrasive particles of piston seals. Spare part: 6.2821.110 Filter 
1	1	6.6035.003	Software CD «790 PC Software 1.0»
1	3	Y.107.0150	Cable strap
1	1	8.732.2003	Metrohm Monograph «Ion chromatography» (English)
1	1	8.732.2013	«IC Application Notes» (English)
1	1	8.790.1003	Instructions for Use (English) for 790 Personal IC
1	1	8.790.8007	Registration card (German/English) for PC program «Metrodata 790 PC Software 1.0»









6.3 Optional accessories










6.3.1 Accessories for 790 Personal IC









Order No.	Description	
6.1825.XXX	Sample loop, made of PEEK For injection valve; incl. two 6.2744.010 PEEK compression fittings 6.1825.230: Volume = 10 μ L 6.1825.210: Volume = 20 μ L 6.1825.220: Volume = 100 μ L	
6.1831.010	PEEK capillary Length = 3 m	
6.2620.150	Pulsation dampener MF Metal-free pulsation dampener to reduce pulsation and prolong the life of separating columns.	
6.2621.080	Capillary tubing cutter for plastic capillaries For 6.1831.010 PEEK capillaries and 6.1822.010 PTFE microcapillaries incl. 5 additional blades	
6.2741.000	Piston seal Spare part for 6.2824.100 Pump head	
6.2741.010	PE piston seal Spare part for 6.2824.100 Pump head (only suitable for aqueous eluents)	
6.2744.070	PEEK compression fitting short Spare part for 6.2824.100 Pump head Set of 5	
6.2821.110	Filter for filter unit PEEK 2 μm Spare part for 6.2821.100 Filter unit PEEK. Set of 10	











Order No.	Description	
6.2824.030	Sapphire supporting ring Spare part for 6.2824.100 Pump head	
6.2824.050	Spring for main piston Spare part for 6.2824.100 Pump head	
6.2824.060	Spring for auxiliary piston Spare part for 6.2824.100 Pump head	
6.2824.070	Zircon piston Spare part for 6.2824.100 Pump head	
6.2824.080	Outlet valve (metal-free) Spare part for 6.2824.100 Pump head	
6.2824.090	Inlet valve (metal-free) Spare part for 6.2824.100 Pump head	
6.2824.100	Pump head (metal-free) Complete, with fixations screws	
6.5324.000	Bottle rack Complete, incl. supply vessels for eluent (2 L), regeneration solution (1 L) rinsing solution (1 L)	



6.3.2 Separating columns and precolumns

Order No.	Description	
6.1005.000	IC anion column PRP-X100 (125 mm) For the determination of anions without chemical suppression. Column dimensions: 125 × 4.0 mm Precolumn: 6.1005.020	
6.1005.010	IC anion column PRP-X100 (250 mm) For the determination of anions without chemical suppression. Column dimensions: 250 × 4.0 mm Precolumn: 6.1005.020	
6.1005.020	IC precolumn cartridge PRP-X100 To prolong the service life of 6.1005.000 and 6.1005.010 IC anion columns PRP-X100. Column dimensions: 20 × 4.0 mm Installation using 6.2821.040 Cartridge head.	
6.1005.030	IC exclusion column PRP-X300 For the determination of organic acids without chemical suppression. Column dimensions: 250 × 4.0 mm Precolumn: 6.1005.040	
6.1005.040	IC precolumn cartridge PRP-X300 To prolong the service life of the 6.1005.030 IC exclusion column PRP-X300. Column dimensions: 20 × 4.0 mm Installation using 6.2821.040 Cartridge head.	
6.1005.050	IC precolumn cartridge PRP-1 To prolong the service life of the 6.1009.000 IC anion column SUPERSEP and the 6.1006.100 IC anion column METROSEP Anion Dual 2. Column dimensions: 20 × 4.0 mm Installation using 6.2821.040 Cartridge head.	
6.1005.100	IC anion column Star-Ion A300 For the determination of anions with chemical suppression. Column dimensions: 100 × 4.6 mm	
6.1005.200	IC anion column Organic Acids For the determination of organic acids. Column dimensions: 250 × 7.5 mm	





Order No.	Description	
6.1006.020	<p>IC column cartridge METROSEP Anion Dual 1 For the determination of anions with and without chemical suppression. Column dimensions: 150 × 3.0 mm Installation using 6.2828.000 Glass cartridge holder.</p>	
6.1006.030	<p>IC precolumn cartridge METROSEP Anion Dual 1 Set of 3 To prolong the service life of the 6.1006.020 IC column cartridge METROSEP Anion Dual 1. Column dimensions: 30 × 3.0 mm Installation using 6.2828.010 Precolumn cartridge holder.</p>	
6.1006.040	<p>IC column cartridge METROSEP Anion Dual 1 For the determination of anions with and without chemical suppression. Column dimensions: 70 × 3.0 mm Installation using 6.2828.020 Glass cartridge holder.</p>	
6.1006.100	<p>IC anion column METROSEP Anion Dual 2 For the determination of anions with and without chemical suppression. Column dimensions: 75 × 4.6 mm Precolumn: 6.1005.050 (installation with 6.2821.050 Twin cartridge holder)</p>	
6.1006.200	<p>Preconcentration cartridge METROSEP Anion For anion preconcentration. Column dimensions: 30 × 3.0 mm Installation using 6.2828.010 Precolumn cartridge holder.</p>	
6.1007.000	<p>IC cation column Nucleosil 5SA For the determination of divalent cations without chemical suppression. Column dimensions: 125 × 4.0 mm Precolumn: 6.1007.010</p>	
6.1007.010	<p>IC precolumn cartridge Nucleosil 5SA To prolong the service life of the 6.1007.000 IC cation column Nucleosil 5SA. Column dimensions: 20 × 4.0 mm Installation using 6.2821.040 Cartridge head.</p>	
6.1008.010	<p>IC cation column Hyperrez Monovalent For the determination of monovalent cations without chemical suppression. Column dimensions: 50 × 4.6 mm</p>	

Order No.	Description	
6.1009.000	IC anion column SUPERSEP For the determination of anions without chemical suppression. Column dimensions: 100 × 4.6 mm Precolumns: 6.1009.010 IC anion column SUPERSEP or 6.1005.010 IC precolumn cartridge	
6.1009.010	IC anion precolumn SUPERSEP To prolong the service life of the 6.1009.000 IC anion column SUPERSEP.	
6.1010.000	IC cation column METROSEP Cation 1-2 For the determination of monovalent and divalent cations without chemical suppression. Column dimensions: 125 × 4.0 mm Precolumn: 6.1010.010	
6.1010.010	IC precolumn cartridge METROSEP Cation 1-2 To prolong the service life of the 6.1010.000 IC cation column METROSEP Cation 1-2. Column dimensions: 20 × 4.0 mm Installation using 6.2821.040 Cartridge head.	
6.1012.X00	Sample pretreatment cartridge IC-RP For non-polar solid phase extraction. Removes organic substances; for the enrichment of heavy metals. With Luer connection. 6.1012.000: set of 50 6.1012.100: set of 10	
6.1012.X10	Sample pretreatment cartridge IC-H Cation exchanger in H ⁺ form. Removes interfering cations, CO ₃ ²⁻ , HCO ₃ ⁻ or for alkaline samples. With Luer connection. 6.1012.010: set of 50 6.1012.110: set of 10	
6.1012.X20	Sample pretreatment cartridge IC-Ag Cation exchanger in Ag ⁺ form. Removes halides. With Luer connection. 6.1012.020: set of 50 6.1012.120: set of 10	
6.1012.X30	Sample pretreatment cartridge IC-OH Cation exchanger in OH ⁻ form. For highly acidic samples. With Luer connection. 6.1012.030: set of 50 6.1012.130: set of 10	

Order No.	Description	
6.1012.200	Sample pretreatment cartridge Chromafix C18 Removes organic substances (not suitable for fluoride determinations). With Luer connection. 6.1012.200: set of 50	
6.2027.050	Column holder Diameter d = 15.0 mm	
6.2821.010	PTFE gasket Spare part for 6.1005.000, 6.1005.010, 6.1005.030, 6.1007.000 and 6.1010.000 separating columns; set of 4.	
6.2821.020	Steel mesh Spare part for 6.1005.000, 6.1005.010, 6.1005.030, 6.1007.000 and 6.1010.000 separating columns; set of 4.	
6.2821.040	Cartridge head To hold precolumn cartridges; installed in the inlet capillary of the separating column.	
6.2821.080	Steel spacer Spare part for 6.2821.040 Cartridge head.	
6.2828.000	Glass cartridge holder To hold the 6.1006.0020 Column cartridge METROSEP Anion Dual 1.	
6.2828.010	Precolumn cartridge holder To hold the 6.1006.0030 Precolumn cartridge METROSEP Anion Dual 1.	
6.2828.020	Glass cartridge holder To hold the 6.1006.0040 Column cartridge METROSEP Anion Dual 1.	

Order No.	Description	
6.2832.000	Suppressor rotor Replacement cartridge for Metrohm suppressor module	
6.2832.010	Connection piece for suppressor rotor with input and output leads	

6.3.3 Additional devices and cables

Order No.	Description	
2.750.0020	750 Autosampler Sampler for the automation of sample injection. Capacity: 128 sample vessels each with an effective volume of ca. 700 µL; incl. accessories. Accessories: 6.2413.000 Glass sample vessels, set of 1000 6.2743.000 PP sample vessels, set of 1000 6.2743.010 Polyethylene stoppers, transparent, set of 1000 6.2743.020 Polyethylene stoppers, red, set of 1000 6.2743.030 Filter stoppers, set of 100	
2.766.0010	766 IC Sample Processor Sampler for the automation of sample injection. Capacity: 127 sample vessels each with an effective volume of ca. 11 mL; incl. accessories. Accessories: 6.2743.050 PP sample vessels, set of 2000 6.2743.060 Polyethylene stoppers, transparent, set of 1000	
6.2125.010	Cable Connecting cable 790 Personal IC (RS232 interface) – PC Adapter cable 25-pin to 9-pin.	
6.2125.020	Cable RS232 extension cable	
6.2125.060	Cable Connecting cable 790 Personal IC (RS232 interface) – PC	
6.2125.110	Cable Connecting cable 790 Personal IC (RS232 interface) – PC	

6.4 Warranty and conformity

6.4.1 Warranty

The warranty on our products is limited to defects that are traceable to material, construction or manufacturing error that occur within 12 months from the day of delivery. In this case, the defects will be rectified in our workshops free of charge. Transport costs are to be paid by the customer.

For day and night operation, the warranty is limited to 6 months.

Glass breakage in the case of electrodes or other parts is not covered by the warranty. Checks, which are not a result of material or manufacturing faults, are also charged during the warranty period. For parts of outside manufacture insofar as these constitute an appreciable part of our instrument, the warranty stipulations of the manufacturer in question apply.




With the regard to the guarantee of accuracy, the technical specifications in the instruction manual are authoritative.

Concerning defects in material, construction or design as well as the absence of guaranteed features, the orderer has no rights or claims except those mentioned above.



If damage of the packaging is evident on receipt of a consignment or if the goods show signs of transport damage after unpacking, the carrier must be informed immediately and a written damage report demanded. Lack of an official damage report releases Metrohm from any liability to pay compensation.

If any instruments and parts have to be returned, the original packaging should be used if at all possible. This applies above all to instruments, electrodes, burette cylinders and PTFE pistons. Before embedment in wood shavings or similar material, the parts must be packed in a dust-proof package (for instruments, use of a plastic bag is imperative). If open assemblies are enclosed in the scope of delivery that are sensitive to electromagnetic voltages (e.g. data interfaces etc.) these must be returned in the associated original protective packaging (e.g. conductive protective bag). (Exception: assemblies with built-in voltage source belong in a non-conductive protective packaging). For damage, which arises as a result of non-compliance with these instructions, no warranty responsibility whatsoever will be accepted by Metrohm.

6.4.2 EU Declaration of conformity

 <h2 style="margin: 0;">EU Declaration of Conformity</h2>							
<p>The METROHM AG company, Herisau, Switzerland hereby certifies, that the instrument:</p> <p style="text-align: center;">790 Personal IC</p> <p>meets the requirements of EC Directives 89/336/EEC and 73/23/EEC.</p> <p>Source of the specifications:</p> <table style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 20%; padding: 5px;">EN 50081-1</td> <td style="padding: 5px;">Electromagnetic compatibility, basic specification; Emitted Interference</td> </tr> <tr> <td style="padding: 5px;">EN 50082-1</td> <td style="padding: 5px;">Electromagnetic compatibility, basic specification; Interference Immunity</td> </tr> <tr> <td style="padding: 5px;">EN 61010</td> <td style="padding: 5px;">Safety requirements for electrical laboratory measurement and control equipment</td> </tr> </table> <p>Description of the instrument:</p> <p style="padding: 5px;">Instrument for recording ion chromatograms with electronic or chemical suppression.</p> <p style="padding: 10px;">Herisau, June 27, 2000</p> <div style="display: flex; justify-content: space-around; align-items: flex-end; margin-top: 20px;"> <div style="text-align: center;">  <p>Dr. J. Frank Development Manager</p> </div> <div style="text-align: center;">  <p>Ch. Buchmann Production and Quality Assurance Manager</p> </div> </div>		EN 50081-1	Electromagnetic compatibility, basic specification; Emitted Interference	EN 50082-1	Electromagnetic compatibility, basic specification; Interference Immunity	EN 61010	Safety requirements for electrical laboratory measurement and control equipment
EN 50081-1	Electromagnetic compatibility, basic specification; Emitted Interference						
EN 50082-1	Electromagnetic compatibility, basic specification; Interference Immunity						
EN 61010	Safety requirements for electrical laboratory measurement and control equipment						

6.4.3 Certificate of conformity and system validation

Certificate of Conformity and System Validation	
This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.	
Name of commodity:	790 Personal IC
System software:	Stored in ROMs
Name of manufacturer:	Metrohm Ltd., Herisau, Switzerland
Principal technical information:	Voltages: 100...120, 220...240 V Frequency: 50...60 Hz
<p>This Metrohm instrument has been built and has undergone final type testing according to the standards:</p> <p style="text-align: center;">IEC801-2/IEC1000-4-2 (class 3), IEC801-3/IEC1000-4-3 (class 2), IEC801-4/IEC1000-4-4 (class 4), IEC801-5/IEC1000-4-5 (class 2/3), IEC801-6/IEC1000-4-6 (class 2), EN50082-1, EN61000-3-2/3/IEC1000-3-2/3, EN61000-4-11/IEC1000-4-11, IEC61326 (class B), EN55011 (class B), EN55022 (class B), EN50081-1 — <i>Electromagnetic compatibility</i></p> <p style="text-align: center;">IEC1010, EN61010, UL3101-1 — <i>Security specifications</i></p> <p>It has also been certified by the Swiss Electrotechnical Association (SEV), which is member of the International Certification Body (CB/IEC).</p> <p>The technical specifications are documented in the instruction manual.</p> <p>The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance. The features of the system software are documented in the instruction manual.</p>	
Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.	
Herisau, June 27, 2000	
 	
Dr. J. Frank	Ch. Buchmann
Development Manager	Production and Quality Assurance Manager

6.5 Index

1	Ordering designation..... 178	<Break times> 83
1st component 95		By components 119
1st peak 95		By peaks 118
7		C
750..... <i>see Autosampler 750</i>	Aspirating tubing 23	Cable (6.2125.010) 185
766..... <i>see IC Sample Processor 766</i>	Aspirate sample 60	Cable (6.2125.020) 185
790 PERSONAL IC 65,66	Figure 7,9	Cable (6.2125.060) 185
790 unit version 82	Installation 16	Cable (6.2125.110) 21
	Aspirating tubing 71	Cable (6.2128.160) 46
	Figure 38,39	Cable (6.2134.100) 21,177
	Fixing 42	Cable (6.2141.110) 47
	Installation 41	Cable strap (Y.107.0150) 178
	Ordering designation..... 176	Calculation formulae 94
	Aspirating tubing 72	Calculation method 95
	Figure 38,39	Calculations 94
	Fixing 42	Calibrate 136
	Installation 40	Calibrate immediately 107
	Ordering designation..... 176	Calibration
	Aspirating tubing (6.1834.010)	Definition 64
	Installation 27	File 111
	Ordering designation..... 176	File type 67
	Assembling the suppressor 165	Formula 108
	asymmetry 121	General information..... 100
	Asymmetry 94	Method 109
	At start 130	Operating tutorial 53
	Autosampler 750	Recalibration 136
	Connection 46	Results..... 108
	Ordering designation..... 185	Update 110
	Settings 46	Calibration / Load from method 111
	Average relative deviation 105	Calibration / Save to method 111
	Axes 131	Calibration coefficients..... 108
	Axes labels 129	Calibration curve
		Display 108
	B	Formula 110
	Back pressure 171	General information..... 100
	Background 131	Graphical representation..... 108
	Background color..... 131	Output 111
	Base 110	Plot area 135
	Baseline..... 131	Preview 111
	Color settings 131	Print 111
	Display..... 130	Report..... 115
	Evaluation..... 96	Results..... 108
	Marks..... 130	Weighting 110
	Baseline marker 130	Calibration data
	Baseline noise 96,116	Export 111
	Baseline with high noise level..... 168	Import 111
	Basic settings 22	Load 111
	Batch 88	Save 111
	Batch reprocessing 144,145	Calibration defaults 115
	Definition..... 65,144	Calibration function 101,110
	Perform..... 145	Calibration graph size 135
	Start..... 147	Calibration level 76,87,106,107,109, 136,140,142,146,149
	Batch reprocessing queue	Calibration method 109
	Close 150	Calibration parameters..... 109
	Create 144	Calibration point..... 106,109
	Edit 146	Calibration points table 109
	Editor 149	Calibration results 115
	File handling..... 144	Cancel last run 141
	File type 67	capacity factor 120
	Functions in editor window..... 150	Capacity factor 89,151
	Print 150	Capacity problems..... 162
	Reprocess options 145	Capillaries 17
	Save..... 144,150	Capillary 48
	Blockage 163,165	Connection..... 17
	Bottle rack (6.5324.000)	
	Connection 16	
	Ordering designation..... 180	

Figure.....	17,18,33	Chromatogram units	72	of detector block.....	15
Ordering designation.....	179	Cleaning the suppressor.....	163	of eluent bottle.....	27
Capillary tubing cutter (6.2621.080).....	17,179	Close	126	of waste container.....	16
Care.....	155	<Close>	140,147	to the PC.....	21
Cartridge head (6.2821.040) Figure.....	31	Close queue dialog window after the queue finishes	140	Connection 15 Connection of detector block.....	16
Installation.....	30	Close window	94	Figure.....	5
Ordering designation.....	184	Coef	116	Connection 53 Connect capillary 18	25
Cautions.....	11	Coef. (balance)	111	Figure.....	26
Cell constant.....	83,171	Coefficient of compression.....	137	Connection 54 Connect capillary 27	25
Entry.....	24	Coefficients	110	Figure.....	26
Printed.....	15	Color settings.....	131	Connection capillary 27 Connection of pulsation dampener 52	25
Cell volume.....	171	Colors	131	Figure.....	7,9,26
Ceramic holder 103 Figure.....	161	Column..... <i>see Separating column</i>		Connection capillary 29 Figure.....	26
Certificate of conformity and system validation.....	188	Column	89,113	Connection capillary 30 Figure.....	7,9
Change passport	147	Column connection capillary 20 Figure.....	7,9,31,33,36,39	Connection capillary 35 Detaching.....	158
Change system	143	Installation.....	31,33,35,36	Figure.....	7,9
Changing of separating column.....	156	Mount precolumn.....	32	Connection capillary 36 Figure.....	7,9
Changing the separation system.....	157	Column holder 70 Figure.....	36,39	Connection capillary 37 Figure.....	7,9
Channel.....	110	Fixing column.....	35,36	Connection of 750 Autosampler	46
Channel 1...8	131	Ordering designation.....	176,184	Connection of 766 IC Sample Processor.....	47
Channel labels	130	Column test	118	Connection of syringe and aspirating tubing.....	16
Channel table	116	COM error	124	Connection piece 68 Figure.....	33
Character set	120	COM port settings.....	22,73	Installation.....	32
Check.....	15	Comment	90,114	Connection piece 106 Cleaning.....	164,166
Checking on-line	84	Comment.....	11	Figure.....	165
Chemical suppression.....	37	Component Information.....	109	Ordering designation.....	185
Choose color	77,79	Selection.....	109	Replacement.....	165
Chromatogram Appearance.....	128	Component	108,109	Connection scheme.....	13
Automatic closing of the window.....	94	Component name	129	Connections.....	154
Autosaving.....	93	Component name + Quantity	129	Connections at suppressor module.....	42
Batch reprocessing.....	144	Component table.....	102,115	Connector 49 Figure.....	18
Chromatogram window.....	124	Components	102	Installation.....	18
Close.....	126	Compress	137	Ordering designation.....	179
Color settings.....	131	Compression fitting 46 Figure.....	17,18,31,33,38	Connector 51 Figure.....	18
Data compression.....	137	Installation.....	17	Installation.....	18
Definition.....	64	Ordering designation.....	177	Ordering designation.....	179
Delete.....	126	Compression fitting 47 Figure.....	17	Connector plug (6.2140.010).....	46
Directory for saving.....	94	Installation.....	17	Connectors for capillaries.....	17
Duration.....	87	Ordering designation.....	179	Contact closure strip 10 Insert connector plug.....	46
Export.....	127	Conc.	109	Contact pressure.....	45,167
File handling.....	125	Concentration.....	101,106,109, 110,117,122,140,142,149	Contact pressure lever 41 Figure.....	9,38,40
File type.....	67	Concentration	109	Install pump tubings.....	37
GLP Log.....	92	Concentration of internal standard.....	118,142,149	Install tubing cartridge.....	167
Identifier.....	76,87	Concentration of std	110	Set contact pressure.....	45
Intermediate saving.....	93	Concentration table	106	Contaminated valves.....	152
Invert.....	137	Concentration units	106	Contamination with heavy metals.....	162
Merge.....	147,148	concentration%	122	Contamination with organic cationic complexing agents.....	162
Open.....	72,125	Concentrations	106	Contamination with organic substances.....	162
Overwrite.....	72	<Concentrations>	104,106	Context sensitive menus.....	68
Path.....	87	Conductivity.....	77,79	Control	79
Plot.....	113	Conductivity detector.....	171	Copy concentrations from level	107
Plot area.....	135	Conductivity measurement.....	171	Copy row(s)	143
Print.....	136	Configuration	82		
Recalibration.....	136	Configuration settings.....	82		
Reintegration.....	136	Conformity.....	187		
Save.....	126	Connection Bottle rack (6.5324.000).....	16		
Scaling.....	128,132	of accessories.....	15		
Start time.....	87				
Title.....	123,139,142,149				
Units.....	72				
CHROMATOGRAM	65				
Chromatogram axes	128				
Chromatogram description	126				
Chromatogram directory	94				
Chromatogram plot	113				
Chromatogram plot size	135				

Copy to clipboard 111
Corr...... 108
Corrected retention time 101
 Correlation coefficient..... 108
Coupling 24
 Figure 7,9
 Ordering designation..... 177
 Remove 46,48
Coupling 26
 Figure 7,9,39
 Installation..... 41
 Ordering designation..... 177
 Passivation 155
 Shutdown..... 156
Coupling 73
 Figure 38,39
 Installation..... 40,41
 Ordering designation..... 177
Coupling 77
 Connect pump tubing..... 38
 Figure 38,39
 Ordering designation..... 177
Create calibration level..... 107
 Cursor..... 131,133,134
Custom..... 118
Custom formulae 94
Custom program..... 120
<Customize>..... 119
Cut row(s)..... 143

D

DATA..... 87
 Data acquisition
 Channel 131
 Parameter report..... 114
 Start 75
 Data acquisition channel 130
 Data compression 137
Data file overwrite 72
Data Log 92,114
Data type 106
 Data type for concentration table 106
 Date 76,87,88
Date/time 87
Date/time when... 76,88
 Dead volume 152,169
 Declaration of conformity 187
<Default colors> 73
Default scheme 146
Default units 127
 Degassing of eluents..... 153
 Degree of protection..... 12
Delay 97
Delete 126,138
<Delete> 81
Delete all left 134
Delete all right..... 134
Delete peak..... 133
Delete row(s) 143,150
Detection of hardware failed..... 85
Detector 87
Detector block 38
 Connection 15
 Connection of column 35
 Figure 7,9,36,39
 Ordering designation..... 176

Detector name..... 87
 Determination
 Definition..... 64
 Example..... 51
 Run number..... 87
 Start 75
 Stop 76
 Diagnostic tests..... 170
Diagnostics 78
Dilution..... 76,88,101,142,149
 Dimensions 175
Directories..... 126
Directory 119,125
Disable detection 99
Disable valley-to-valley 99
 Disposal 12
Divisor..... 137
Done 142
Door 1
 Figure 3,7,9
 Double peaks 152,169
 Drift..... 153,168,171
Drift compensation 129
 Drivers 67
Duplicate row(s)..... 143
Duration 75,87

E

E### 85
 Earthing 12
<Edit>..... 140
<Edit appearance>..... 147
Edit integration event..... 98
<Edit integration parameters>..... 146
<Edit passport> 147
<Edit report options> 147
Edit sample description 75
<Edit sample table> 146,149
 Editor window 150
effectivity, TP..... 121
effectivity, TP/m..... 121
 Electrical safety 12
 Electromagnetic compatibility 174
 Eluent
 Composition 90
 Information 90
 Report..... 113
Eluent..... 90,113
Eluent A 90
Eluent B 90
 Eluent bottle
 Connection 27
Eluent C 90
 Eluent change 153
 Eluent vessel
 Recycling..... 156
 Eluents
 Notes 153
 EMC 174
 Emitted interference 174
Enable detection..... 99
Enable valley-to-valley 99
ENABLED 81
END..... 81

End fitting **64**
 Figure 33
 Installation 32
 Error messages..... 85,168
 Table of messages..... 168
 Error window 168
 EU Declaration of conformity 187
European Pharmacopoeia 95
 Evaporation 153,168
Event..... 84,99
Events..... 98
every 130
 Exchanging the pump tubing..... 167
<Exclude point>..... 109
Exit 63,143,150
Export / AIA file 127
Export / Raw data to txt..... 127
Export calibration 111
External..... 95
 External devices 46
 External standard
 calibration..... 100,102,109

F

Failure..... 124
 Fastening screws **33**
 Detach pump head **34**..... 158
 Figure 7,9
 Faults..... 151,168
Feedthrough 3
 Figure 3
 Install syringe tubing **25** 16
Feedthrough 4
 Figure 3
 Install aspirating tubing **23** 16
File..... 109,116
file name 123
File name 149
File output options..... 119
 File types 67
File window 125
Filename 125
<Fill>..... 78,80
Fill level with concentration 107
Filter unit PEEK 28
 Figure 7,9,18,26,38,39
 Information 151
 Installation 18,38
 Notes..... 154
 Ordering designation 178
Finished 85,124
Flow 72,78,79,80,81,90
 Flow chart
 Installation 13
 Flow constancy 172
Flow correction 83
 Flow range 172
 Flow rate..... 101
 Correction factor 83
 Display 79,90
 Program instruction..... 81
 System startup value..... 78,80
 Unit..... 72
 Fonts 70
Forget calibration points before reprocessing..... 146

Format	136
Formula	94,110
Formula for calibration function.....	110
Formula set	94
From void volume %	95
Front	3
Fundamentals of the operation.....	63
Fuse holder 13	
Figure.....	5,20
Replacement of fuses.....	19
Fuse peaks	134
Fuses	20,174
G	
GDI printer	72
General	87,93,112
General information	
on maintenance.....	155
on precolumns.....	30
on separating columns.....	34
on suppressor module	37
General precautionary rules	12
Glass beads	152
Glass cartridge holder (6.2828.000) ..	184
Global preferences	71
Global settings.....	71
Glossary.....	64
GLP.....	72,91,92,114,170
GLP comment.....	91,92
GLP Log.....	91,92,114
GLP messages	114
GLP On	72
Graphical representation	128
Graphs	108
<Graphs>	104,108
Grid	129
group	123
Group	103
Groups	119
H	
Hair	131
halfwidth	120
Handling of solvents	12
Hardware	
Error.....	84
Settings.....	24,82
Start	75
Hardware error	84
Hardware?	78,82
Hazard	11
Heating time	171
height	120
Height	109
Height of theoretical plate.....	121
height%	120
Help	69
HETP/dp	121
Hexagon key (6.2621.030)	158,177
<Hide all>	130
High-pressure pump	
Alarm stop	83
Connection of eluent bottle.....	27
Deaerating.....	27
Installation	25
Maintenance.....	152
Optional accessories.....	179
Practical notes	152
Program instruction	81
Start	75,80
Stop	75,80
Technical data	172
Holding clamp 42	
Figure	9,40
Install tubing cartridge.....	37
Housing	175
Housing 50 for filter unit PEEK	
Figure	18
Installation	18
How to draw chromatogram point	130
I	
IC anion column METROSEP	
Anion Dual 2	182
IC anion column Organic Acids	181
IC anion column PRP-X100	181
IC anion column Star-Ion A300	181
IC anion column SUPERSEP	183
IC anion precolumn SUPERSEP	
Installation	33
Ordering designation.....	183
IC Applications Collection	
(8.732.2013)	10
IC cation column	
Hyperrez Monovalent.....	182
IC cation column METROSEP	
Cation 1-2	156,183
IC cation column Nucleosil 5SA	182
IC column cartridge	
METROSEP Anion Dual 1	182
IC exclusion column PRP-X300	181
IC precolumn cartridge	
METROSEP Anion Dual 1	182
IC precolumn cartridge	
METROSEP Cation 1-2.....	183
IC precolumn cartridge	
Nucleosil 5SA	182
IC precolumn cartridge PRP-1	181
IC precolumn cartridge PRP-X100	181
IC precolumn cartridge PRP-X300	181
IC pump	80
IC Sample Processor 766	
Connection	47
Ordering designation.....	185
Settings.....	48
IC separating columns	<i>see Separating column</i>
Icons	
790 Personal IC	74,78
Main window.....	66
PC.....	74,77
Watch window	74,77
ICPump	81
ID	89
ident	123
Ident	76,87,124,125,139,142,149
Identical separation system.....	156
Identification	100,104,105
<Identification>	104
Identification parameters.....	105
Identification window	103
Identifier.....	76,87,123,139,142,149
If disk method is newer	72
If method changed	72
Ignore last data directory	72
Immunity to interference	174
Import calibration	111
In corner	130
In seconds	127
Include retention time	127
Increase in the baseline	162
Increment	143,150
index	122
Index	95,103,118
Indexes	110
Info 1	76,88
Info 2	76,88
Information on the	
Instructions for Use	10
Initial	78
Initialization	84
<Inject>	78,80
INJECT done	85
Injected sample volume	101
Injected volume	76,88,140,142,149
Injection valve 22	
Break time	83
Connect 750 Autosampler.....	46,47
Figure	7,9,26,36,39
Program instruction	81
Switch to "FILL"	55,60,78,80
Switch to "INJECT"	55,61,78,80
Technical data.....	171
Injections	142
Inj-s	139
Inlet capillary 18	
Connection of	
pulsation dampener 52	25
Figure	7,9,26
Inlet capillary 37	
Connection to	
suppressor module 39	41
Figure	7,9,36,39
Installation	35
Inlet valve 93	
Cleaning	160
Components.....	161
Figure	159,161
Ordering designations.....	180
Replacement	160
Input	116
Insert peak	133
Installation	13,52
of pump tubings.....	38
PEEK Filter unit	18
Instructions for Use (8.790.1003) .	10,178
Instrument control	80
Instrument description.....	1
Instrument icon.....	78
Instrument settings.....	70
Instrument version.....	82
Integration	96,114
Algorithm	96
Definition	65
Events.....	98
Parameters	56,97
Reintegration	136
Integration parameters	96,146
Interior	

of 2.790.0010 Personal IC 6
of 2.790.0020 Personal IC 8
Internal 95
Internal standard
 amount 76,88,142,149
Internal standard
 calibration 100,109,117
Interpolation 95
Introduction 1
Invert! 137
Ion chromatography 151
Items to report 112

K

Keyboard 68,134
Knurled screw **8**
 Figure 5
 Open rear panel 15

L

Label always visible 130
Labels 129
Last batch 144,145
Last update 87
Leaks 12,44,45,168
Length 89
Length units 135
Level 106,107,109,140,142
Line width 131
Linear flow 89
Linearity 171
Links 73
List of program instructions 81
List of users 70
<Load default> 135
<Load defaults> 131
Local 109
Location 15
Lock system 71
Log In window 22
Logarithmic index 89
Login 63
Loss of suppressor capacity 165

M

Main window
 elements 66
 Icons 66
Mains cable
 Connection 20
 Ordering designation 176
Mains connection
 Procedure 19,20
 Safety notes 12
 Technical data 174
Mains connection plug **12**
 Connection of mains cable 20
 Figure 5,20
Mains frequency 174
Mains pilot lamp **2**
 Figure 3

Function 20
Mains switch **11**
 Figure 5,20
 On/off switching
 of the instruments 20
Mains voltage
 Setting 19
 Technical data 174
Maintenance 151,155
Make report 94,112,136,147
Malfunctions 168
Manufit housing **55**
 Figure 31
 Installation 30
Manufit pressure screw **62**
 Figure 31
 Installation 30
Master 70
Math 94
Maximum 116
Maximum back pressure 171
Maximum concentration 103
Maximum pressure limit
 Display 79
 Program instruction 81
 System startup value 79,80
 Technical data 172
Measure 124
Measure Baseline 75
Measure(Baseline) 124
Measurement principle 172
Measuring priority 73
Menu bar 66
<Merge> 147,148
Merge peaks 134
Method 125,149
 Definition 64,86
 File handling 86
 File type 67
 GLP Log 91
 Open 86
 Overwrite 72
 Path 87
 Save 72,86
METHOD 87
Method Log 91,114
Method setup 93
Metrohm Monograph
 "Ion chromatography" 10,151
Metrohm service 155
Metrohm suppressor module
 *see* *Suppressor module*
Microfiltration 151,152,153
middle 130
min C (max C) 103
Min height 98
Minimum 116
Minimum concentration 103
Minimum pressure limit
 Display 79
 Program instruction 81
 System startup value 79,80
 Technical data 172
Mode 120
Modify chromatogram
 appearance 147
 Monograph
 "Ion chromatography" 10,151
More items to report 114

Mounting pin **45**
 Figure 9,40
 Install tubing cartridge 37
Mounting rail **19**
 Figure 7,9,36,39
 Installation of column holder 35,36
Mouse 68,134
Move selected point 134
MSM *see* *Suppressor module*
Multi-channel
 chromatogram 130,147,148

N

name 123
Name 103,106,116,119
Name and quantity
 of the component 129
Name of component 103,106,123,129
Name of user 87
Negative peaks 98,137
No 139,142,149
Noise 116,171
None 95,129
Normalized concentration 117
Notation 11
Notations 101
Notes 151
Novice 70
number 120
Number 89,103,106
Number of components 104
Number of copies 136
Number of events 98
Number of injections 139,142
Number of peaks 97
Number of the component 106
Number of theoretical plates 94,121

O

Obstruction 163
OK, Apply buttons mean 72
On/off switching of the instrument 20
On-line 84
Online help 69
Open 77,78
Open / Chromatogram 125,144
Open / System 74
<Open all files> 145
Open control window 43,44
<Open example> 145
Opening **5**
 Figure 5
 Insert aspirating tubing 27,40,41
Opening **6**
 Figure 5
 Insert outlet capillary **82** 42
 Insert outlet capillary **83** 41
Opening **7**
 Insert aspirating tubing 27
Opening chromatograms 72
Opening the 790 Personal IC 12
Operating temperature 171

Operating tutorial.....	51
Operation.....	63
Optional accessories.....	179
Ordinary components	103
Organic modifiers	153
Other peaks	105
Outlet valve 94	
Cleaning.....	160
Components.....	161
Figure.....	159,161
Ordering designations.....	180
Replacement.....	160
Overview	
File types.....	67
Program windows.....	65
P	
Packaging.....	15
Packing material	89
Page layout	135
Page margins	135
Paper size.....	136
Parameter	94
Parameter for program instruction.....	81
Particle size	89
Parts and controls.....	3
Passivation.....	155
Passport.....	86,147
Password	63,70,87
Paste row(s)	143
Pause	141
<Pause>.....	140,141
PC icon.....	77
Peak.....	103
Allocation.....	57
area.....	106
Area.....	101,109,120
Delete.....	133
Drop line.....	96
End.....	96
Fuse.....	134
Height.....	101,106,109,120
Identification.....	104
Insert.....	133
Integration.....	96
Integration events.....	99
Label.....	129
Merged peaks.....	96
Minimum height.....	98
Negative peaks.....	98
Number.....	57,103
Number of peaks detected.....	97
Peak editor.....	133
Peak number.....	120,129
Peak table.....	113,117,120
Recognition.....	96
Recognition threshold.....	97
Single peak.....	96
Split.....	134
Start.....	96
Top.....	96
Valley.....	96
Width.....	97
Width at half height.....	120
Peak asymmetry.....	121
Peak broadening.....	152,169
Peak editor	
Functions.....	133
Switching on/off.....	133
Peak editor	133
Peak labels	129
Peak number	129
Peak table	113,117
PEEK Capillary tubing 18	
Mount.....	47
Peristaltic	80,82
Peristaltic pump	
Alarm stop.....	83
Deactivation.....	82
Notes.....	153
Preparation.....	37
Program instruction.....	82
Start.....	75,80
Stop.....	75,80
Technical data.....	172
pH.....	153
Pictograms.....	11
Piston cartridge 88	
Figure.....	159
Piston guide sleeve 89	
Figure.....	159
Piston guide sleeve 91	
Figure.....	159
Piston seal 92	
Figure.....	159
Ordering designations.....	179
Replacement.....	158
Plate number.....	151
Plot area.....	135
Pmax	81
Pmin	81
Polarity	82
Position on	
Autosampler.....	76,88,139,142,149
Power consumption.....	174
Power failure.....	84
Power on	84
Practical notes.....	151
Precautionary rules.....	12
Precipitates.....	152,153
Precolumn	89
Precolumn cartridge 59	
Figure.....	31
Installation.....	30
Precolumn cartridge 67	
Figure.....	33
Installation.....	32
Precolumn cartridge holder (6.2828.010)	
Figure.....	33
Installation.....	32
Ordering designation.....	184
Precolumns	
General information.....	30
Installation.....	30
Internal diameter.....	89
Length.....	89
Notes.....	151
Ordering designations.....	181
Precolumn glass cartridge with cartridge holder.....	32
with cartridge head.....	30
Preconcentration cartridge	
METROSEP Anion.....	182
Preparations for operating tutorial.....	52
Pressure	
Display.....	77,79
drop.....	168
measurement.....	172
range.....	172
rise.....	168
Unit.....	72
Pressure	72,79,90
Pressure max	79,80
Pressure min	79,80
Preview	136
Preview all	111
Preview this	111
Print	136,143,150
Print all	111
Print preview.....	136
Print spooler.....	72
Print this	111
Print via print spooler	72
Printed by	112
Printer	116
Printer settings.....	136
Printer setup	136
Printing.....	135
Printing order	118
Printing range.....	136
Processing	93,124
Processing method.....	64
Program.....	80,81
Close.....	63
Installation.....	21
Program flag.....	81
Program for 750 Autosampler.....	47
Program for 766	
IC Sample Processor.....	48
Program icon.....	22
Program instructions.....	81
Program step.....	80
Program windows.....	65
Start.....	63
Prompts.....	66
Propagate	143,150
Protection class.....	12,20
Protective earth.....	20
PTFE gasket 57	
Figure.....	31
Installation.....	30
Ordering designation.....	184
PTFE microcapillary (6.1822.010).....	17
Pulsation.....	168
Pulsation dampener 52	
Connection.....	25
Figure.....	26
Ordering designation.....	179
Practical notes.....	152
Rinsing.....	29
Pump capacity.....	172
Pump drive 44	
Figure.....	9,40
Pump head 34	
Components.....	159
Disassemble.....	158
Figure.....	7,9,159
Maintenance.....	158
Mounting.....	161
Ordering designation.....	180
Technical data.....	172
Transport security.....	25
Pump stop	84
Pump tubing 74	
Connect aspirating tubing 72	40
Connection to suppressor module 39	41
Exchange.....	153,167
Figure.....	38,39
Installation.....	37

Scaling chromatogram axes.....	128	Set all	129,130	Figure	31
Scaling of X axis	132	Set peak end	99	Installation	30
Scaling of Y axis	132	Set peak start	99	Ordering designation	184
Scheme	104	Setting up the instrument	15	<Step>	80
Screen	116	Setup	97	Stop determination	76
Screw 84		Show all	93	Stopper 76	
Figure.....	159	<Show all>	130	Figure	38
Undo	158	Show baseline and peaks	130	Install pump tubing.....	167
Screw cap 65		Show shutdown system after the queue finishes 140		Install pump tubings	37
Figure.....	33	Shutdown	156	Strap (Y.107.0150)	42
Installation.....	32	SHUTDOWN	85	Subtract	137
Screw cap 69		Shutdown hardware	75	Subtraction of a chromatogram	137
Figure.....	33	Sleeve 66		Sunlight	15
Installation.....	32	Figure	33	Supply bottle	40
Screw holder 95		Installation	32	Suppression	37
Figure.....	159	Sleeve 100		Suppressor	80,82
Screw nut 105		Figure	161	Suppressor holder 108	
Figure.....	165	Slope	96,97	Figure	165
Mounting.....	166	Snap-action lever 43		Mounting	166
Seal 104		Figure	9,38,40	Suppressor inlet capillary 78	
Figure.....	161	Install tubing cartridge.....	37	Figure	40,42
Sealing ring 99		Raise tubing cartridge	153,167	Installation	41
Figure.....	161	Software installation.....	21	Suppressor inlet capillary 80	
Security options	70	Software settings	70	Figure	38,40,42
Security system	70	SOP	170	Installation	42
Select channel	130	Special tool 96		Suppressor inlet capillary 81	
Select end point	134	Figure	159	Figure	38,40,42
Select nearest point	133	Ordering designation.....	177	Installation	41
Select start point	133	Remove piston seal 92	158	Suppressor module 39	
Select top point	133	Special tool 97		Autostep	82
Select valley point	134	Figure	159	Break time	83
Selected	125	Install piston seal 92	160	Capacity	154
Selection		Ordering designation.....	177	Cleaning	163
of channels	130	Spilled chemicals	155	Connection	41
of printer	136	Split peaks	134	Figure	9,39
of sample loop.....	34	Splitting.....	152,169	General information.....	37,154
Selectivity	151	Spring 87		Notes	154
Separating column 63		Figure	159	Program instruction	82
Changing	156	Ordering designations.....	180	Protection	154
Characteristic data	151	Spring retainer 86		Regeneration	162
Connection at detector block 38	35	Figure	159	Replacement	165
Connection to		Standard chromatogram	151	Switching to next position.....	80
suppressor module 39	41	Standard component	110,118	Technical data.....	173
Dead volume	152	Standard equipment.....	176	Suppressor outlet capillary 79	
Figure.....	31,36,39	Standard operating procedures	170	Figure	40,42
Fixing	35,36	Start	140	Installation	41
Information.....	89	<Start>	140	Suppressor outlet capillary 82	
Install precolumn	30	Start determination	75	Figure	40,42
Installation with suppressor	36	Start time for event	99	Fixing	42
Installation without suppressor	35	Started	139	Installation	41
Internal diameter	89	Starting	84	Suppressor outlet capillary 83	
Length.....	89,101	Startup hardware	75	Figure	40,42
Ordering designations	181	Static charges	12	Fixing	42
Packing material	89	Statistical weight	110	Installation	41
Particle size.....	89	Status bar	66,124	Suppressor rotor 107	
Practical notes	151	Status messages.....	84	Cleaning	166
Protection.....	151	Steel connector 56		Figure	165
Regeneration	152	Figure	31	Inserting	166
Report	113	Installation	30	Ordering designation	185
Rinsing.....	35,36	Steel meshes 58		Replacement	165
Separation efficiency	151	Figure	31	Syringe (6.2816.020)	
Serial number	89	Installation	30	Aspirate sample	60
Storage	152	Ordering designation.....	184	Connection of syringe	16
Test results	118	Steel meshes 61		Deaerate pump	27
Void volume	89	Figure	31	Ordering designation	178
Separation efficiency	151	Installation	30	Syringe tubing 25	
Separator	119	Ordering designation.....	184	Connection to waste container ...	46,48
Serial number 17		Steel meshes 61		Figure	7,9
Figure.....	5	Figure	31	Installation	16
Service.....	155,170	Installation	30	System	
Servicing	155	Ordering designation.....	184	Change.....	74
		Steel spacer 60		Close	74
				Conditioning	44,45

Definition..... 64,74
 File handling 74
 File type 67
 Open..... 74
 Print system parameters..... 76
 Save..... 74
 System functions 75
 System settings 78,79
 System startup values 80
 System state window 84
 System window..... 74
SYSTEM 65
System 139,142
System / Change..... 74
System / Close..... 74
System / Save..... 74
System startup values..... 75,80
System state..... 84
SYSTEM STATE 65,168

T

Tab size 119
 Tabulated calibration..... 100,109
 Technical data 171
Temp. 90
 Temperature 90
 Temperature dependence 171
 Temperature deviation..... 171
 Temperature stability 171
Template..... 119
 Test chromatogram 34
Theoretical plates 94,121
 Thermostatting..... 171
this color is default 131
This run 106
 Threshold for peak recognition..... 97
Tic marks..... 129
 Time..... 76,87,88,99,103
Time axis 128
 Time program
 Enable 81
 Entry 80
 Report 76
 Verify 81
 Title bar 124
<To Batch> 126,144
 Tool bar 66
Total % for normalization 118
 Transfer tubing **8**
 Installation..... 46
 Transport 15
 Transport damage..... 186
 Transport security screws **10**
 Figure 5
 Remove 25
 Treatment of eluents..... 153
 Tubing cartridge **40**
 Figure 9,38,40
 Install pump tubings 37
 Raise..... 153,167
type 122

U

Undo 133,143,150

Unit for X axis 128
Units 116
 Universal component 104,111
Unselect peak 134
Update 110
<Update> 105
**Update method file in <METHODS>
 directory after reprocessing** 146
UploadStartupValues 84
US Pharmacopoeia 95
**Use method from file
 for reprocessing** 145
Use peristaltic pump 82
Used 109
User 87

V

Validation..... 170
Valve 80,81
 Valve housing **98**
 Figure 161
<Verify> 81
Vial 139,142,149
Vial number 76,88
View all 129,132
Void time 95,101,120
Void time/volume 95
Void volume 89,95
Volume 76,88,101,106,140,142,149

W

Waiting 124
Waiting for INJECT 84
 Warning 11
 Warranty 186
 Waste container..... 16,41,42,46,48
 Watch window 77
 Menu options..... 77
WATCH WINDOW 65
Weight..... 175
while running every ... min 93
Width..... 97
Wind.% 103
 Working directory 126
Worst case..... 105

X

X from 128
X full scale..... 132
X to 128
X units 128

Y

Y from 129
Y full scale..... 132
Y to 129

Z

Zero 116
 Zircon piston **85**
 Cleaning 158
 Figure 159
 Ordering designation 180
 Replacement 158

