



Mode selection



Press key <MODE> until the desired mode is displayed, press <ENTER>, select desired measured quantity with <SELECT>, and confirm with <ENTER>.




Recall method from internal method memory with <USER METH>. Select method with <SELECT> or by entering its name.

Modes:

DET	pH, U, Ipo1, Upo1	D ynamic E quivalence Point T itration.
MET	pH, U, Ipo1, Upo1	M onotonic E quivalence Point T itration.
SET	pH, U, Ipo1, Upo1	S et E ndpoint T itration.
MEAS	pH, U, Ipo1, Upo1, T	M EASuring.
CAL		C ALibration.
TIP		T itration- P rocedure. Linking of various commands and methods to a titration procedure.

Sample data





Input for sample data:

- Sample identification. Can be used as sample specific calculation values
- Sample size

Inquiries with silo = OFF (LED "Silo" is OFF).

Display	Initial value	Meaning	Input range
id#1 or C21 id#2 or C22 id#3 or C23		Sample identification. Can be used as sample specific calculation values C21...C23.	up to 8 ASCII characters
smpl size smpl unit:	1.0 g g	Sample size, value C00. Unit of sample size.	0...±999 999 g, mg, mL, µL, pc or up to 5 ASCII characters

Calculations and definitions for data output



Input of formula.
Preselection of automatic report output at the end of a determination.
Allocations of values for statistics, common variables, silo calculations and temporary variables (TIP calculations).

Display	Initial value	Meaning	Input range
>formula		Input of formulas	
RS? RS1=EP1*C01/C00		Enter formula number. Enter formula by means of 3 rd functions of keyboard and confirm with <ENTER>. CXX are calculations variables, see below.	1...9
RS1 text	RS1	Text for result output.	RS1 or up to 8 ASCII characters
RS1 decimal places	2	Number of decimal places for result.	0...5
RS1 unit:	%	Selection of result unit. Select the unit with <SELECT> or enter a unit. Enter the values of calculation variables C01..C19 with <C-FMLA>.	%, ppm, g/L, mg/mL, mol/L, g, mg, mL, mg/pc, no unit or up to 6 ASCII characters
>silo calculations		Allocations for silo calculations	
>common variables		Allocations for common variables	
>report		Selection of report blocks for data output	
report:		Press <SELECT> for selection. If you wish several reports, use ";" as separator.	full, short, mplist, curve, derive, comb, scalc full, scalc srt, calc, param, calib, ff
>mean		Allocations for statistics calculations	
>temporary variables		Allocations for calculations in TIP	

Meaning of the calculation variables CXX

Variable	Meaning
C00	Sample size, key <SMPL DATA>.
C01...C19	Method specific calculation values, such as molecular mass, factors, key <C-FMLA>.
C21...C23	Sample specific calculation values, such as dilution factors, key <SMPL DATA>.
C24, C25	Variables for storing determination results in the silo memory.
C26, C27	Means from silo calculations.
C30...C39	Common variables, i.e. for titer.
C40	Initial measured value of the sample, in MEAS last measured value.
C41	End volume.
C42	Determination time.
C43	Volume drift for SET with conditioning.
C44	Temperature.
C45	Dispensed start volume.
C46	Asymmetry-pH (calibration).
C47	Electrode slope (calibration).
C51...C59	Fix EP for DET and MET.
C61...C69	pK/HNP values for DET and MET.
C70...C79	Temporary variables for calculations in TIP.

Configuration



Settings of peripheral units.
 General settings.
 Settings of RS232 interface, values of the common variables.

<i>Display</i>	<i>Initial value</i>	<i>Meaning</i>	<i>Input range</i>
>peripheral units		Settings of peripheral units	
send to:	IBM	Selection of printer.	Epson, Seiko, Citizen, HP, IBM
balance:	Sartorius	Selection of balance.	Sartorius, Mettler, Mettler AT, AND, Precisa
record:	U	Selection of record type to be outputted at the analog output.	U, dU/dt, V, dV/dt, U(rel), T
>auxiliaries		General settings	
dialog:	english	Selection of dialog language.	english, deutsch, français, español, italiano, portugese, svenska
date time	YYYY-MM-DD HH:MM		
run number	0	Current run number for result output.	0...9999
auto start	OFF	Automatic starts of titrations.	1...9999, OFF
start delay	0 s	Waiting time before start of titration.	0...999999 s
device label		Device label.	8 ASCII characters
program	794.0010	Program version.	read only
>RS232 settings		RS232-settings	
baud rate:	9600	Baud rate.	300, 600, 1200, 2400, 4800, 9600
data bit:	8	Data bit.	7, 8
stop bit:	1	Stop bit.	1, 2
parity:	none	Parity.	none, odd, even
handshake:	HWS	Handshake.	HWS, HWf, SWchar, SWline, none
RS control:	ON	Receiving of commands via RS	ON, OFF
>common Variables		Values of common variables	

Parameters for DET and MET

Display	Initial value	Meaning	Input range
>titration parameters		General titration parameters	
meas.pt.density	4	Parameter for DET: Measuring point density. 0 is highest.	0...9
min.incr.	10.0 µl	Parameter for DET: Minimum increment.	0...999.9 µL
[V step	0.10 ml]	Parameter for MET: Volume increment.	0...9.999 mL
titr.rate	max. ml/min	Dosing rate for increments.	0.01...150 mL/min, max.
signal drift	50 mV/min	Drift criterion for measured value acquisition.	pH, U: 0.5...999 mV/min, OFF I: 0.1...50...99.9 µA, OFF
equilibr.time	26 s	Waiting time for measured value acquisition.	0...9999 s, OFF
start V:	OFF	Type of start volume.	abs., rel. OFF
start V	0.00 ml	Volume for <i>absolute</i> start volume.	0...999.99 mL
factor	0.0	Factor for <i>relative</i> start volume: factor * smpl size.	0... ±999999
dos.rate	max. ml/min	Dosing rate for start volume.	0.01...150 mL/min, OFF
pause	0 s	Waiting time, i.e. after start volume.	0...999999 s
meas.input:	1	Measuring input for pH and U.	1, 2, diff.
temperature	25.0 °C	Temperature.	-170.0...500.0 °C
>stop conditions		Stop conditions for titration	
stop V:	abs.	Type of stop volume.	abs., rel., OFF
stop V	99.99 ml	Volume for <i>absolute</i> stop volume.	0.00...9999.99 mL, OFF
factor	999999	Factor for <i>relative</i> stop volume: factor * smpl size.	0... ±999999
stop pH	OFF	Stop at measured value pH, U, I.	pH: 0.00...±20.00, OFF U: 0... ±2000 mV, OFF I: 0.0... ±200 µA, OFF
stop EP	9	Stop after a number of EP's have been found.	1...9, OFF
filling rate	max. ml/min	Filling rate.	0.01...150 mL/min, max.
>statistics		Statistics calculation	
status:	OFF	Status of statistics calculation.	ON, OFF
mean	n= 2	Number n of single values for statistics.	2...20
res.tab:	original	Result table for statistics calculation.	original, delete n, delete all
delete	n= 1	Delete data from sample number n.	1...20
>evaluation		EP evaluation/recognition	
EPC	5	Endpoint criterion.	DET 0...200 MET pH: 0.10...0.5...9.99 U: 1...30...999 mV I: 0.1...2...99.9 µA
EP recognition:	all	Recognition of EP's which fulfill the EP criterion. If "window" is selected, lower and upper limits of windows are inquired.	all, greatest, last, window, OFF
Fix EP1 at pH	OFF	Interpolation of volume at a given measured value pH, U or I. Up to 9 fix EP's (→ C5X).	pH: 0.00... ±20.00, OFF U: 0... ±2000 mV, OFF I: 0.0 ... ±200 µA, OFF
pK/HNP:	OFF	Evaluation of pK or half neutralization potential (HNP → C6X).	ON, OFF
>preselections		Preselections for the sequence	
req.ident:	OFF	Request of identifications after start of titration.	id1, id1 & 2, all, OFF
req.smpl size:	OFF	Request of sample size after start of titration.	value, unit, all, OFF
activate pulse:	OFF	Pulse output on I/O line L6.	ON, OFF

Parameters for SET

Display	Initial value	Meaning	Input range
>SET1			
Control parameters for EP1			
EP1 at pH	OFF	Preset EP1 at pH, U or I.	pH: 0.00... ±20.00, OFF U: 0... ±2000 mV, OFF I: 0.0 ... ±200 µA, OFF
dynamics	OFF	Distance from EP where constant dosing should stop and controlling begins.	pH: 0.01... ±20.00, OFF U: 1... ±2000 mV, OFF I: 0.1 ... ±200 µA, OFF
max.rate	10.0 ml/min	Maximum dosing rate.	0.01...150 mL/min, OFF
min.rate	25.0 µl/min	Minimum dosing rate.	0.01...9999 µL/min
stop crit: stop drift t(delay) stop time	drift 20 µl/min 10 s OFF s	Type of stop criteria. Titration stops if stop drift is reached. Titration stops if there is no dosing during t(delay). If t(delay) is "INF" stop after at time.	drift, time 1...999 µL/min 0...999 s, inf. 0...999999 s, OFF
>SET2			
Parameters for EP2. Identical as for EP1.			
>titration parameters			
General titration parameters			
titr.direction:	auto	auto: Direction is set automatically. +: Titration to higher pH, voltage or current.	+, -, auto
start V: start V factor dos.rate. max.	OFF 0.00 ml 0.0 .mL/min	Type of start volume: absolute or relative. Volume for <i>absolute</i> start volume. Factor for <i>relative</i> start volume: factor * smpl size Dosing rate for start volume.	abs., rel., OFF 0...999.99 ml 0... ±999999 0.01...150 mL/min, OFF
pause	0 s	Waiting time <i>after</i> start volume.	0...999999 s
meas. input:	1	Measuring input for pH and U or polarization current I _{pol} or polarization voltage for U _{pol} in steps of 10 mV and test for polarized electrodes.	1, 2, diff.
temperature	25.0 °C	Titration temperature.	-170.0...500.0 °C
>stop conditions			
Stop conditions for titration			
stop V: stop V: factor	abs. 99.99 ml 999999	Type of stop volume. Volume for <i>absolute</i> stop volume. Factor for <i>relative</i> stop volume: factor * smpl size.	abs., rel., OFF 0.00...9999.99 mL, OFF 0... ±999999
filling rate max.	ml/min	Filling rate.	0.01...150 mL/min, max.
>statistics			
Statistics calculation			
status: mean res.tab: delete	OFF n= 2 original n= 1	Status of statistics calculation. Number n of single values for statistics calculation. Result table for statistics calculation. Delete data from sample number n.	ON, OFF 2...20 original, delete n, delete all 1...20
>preselections			
Preselections for the sequence			
conditioning: display drift:	OFF ON	Automatic conditioning of titration vessel. Display of drift during conditioning.	ON, OFF ON, OFF
req.ident:	OFF	Request of identifications after start of titration.	id1, id1 & 2, all, OFF
req.smpl size:	OFF	Request of sample size after start of titration.	value, unit, all, OFF
activate pulse:	OFF	Pulse output on I/O line L6.	first, all, cond., OFF

Parameters for MEAS

Display	Initial value	Meaning	Input range
>measuring parameters		Measuring parameters	
signal drift	OFF mV/min	Drift criterion for measured value acquisition.	pH, U, Ipol: 0.5...999 mV/min, OFF Upol: 0.05...99.9 µA/min, OFF T: 0.5...999 °C/min, OFF
equibr.time	OFF s	Waiting time for measured value acquisition.	0...9999s, OFF
meas.input:	1	Measuring input for pH and U.	1, 2, diff.
I(pol)	1 µA	Polarization current.	-127...127 µA
U(pol)	400 mV	Polarization potential.	-1270...1270 mV
electrode test:	OFF	Electrode test.	ON, OFF
temperature	25.0 °C	Temperature.	-170...500.0 °C
>statistics		Statistics calculation	
status:	OFF	Status of statistics calculation.	ON, OFF
mean	n= 2	Number n of single values for statistics calculation.	2...20
res.tab:	original	Result table for statistics calculation.	original, delete n, delete all
delete	n= 1	Delete data from sample number n.	1...20
>preselections		Preselections for the sequence	
req.ident:	OFF	Request of identifications after start of titration.	id1, id1 & 2, all, OFF
req.smpl size:	OFF	Request of sample size after start of titration.	value, unit, all, OFF
activate pulse:	OFF	Pulse output on I/O line L6.	ON, OFF

Parameters for CAL

Display	Initial value	Meaning	Input range
>calibration parameters		Calibration parameters	
meas.input:	1	Measuring input.	1, 2, diff.
cal.temp.	25 °C	Calibration temperature.	-20...120.0 °C
Puffer #1 pH	7.00	pH value of first buffer (up to 9 buffers). <CLEAR> sets OFF. Buffers are requested until "OFF" is set.	0... ±20, OFF
Puffer #2 pH	4.00		
Puffer #3 pH	OFF		
signal drift	2.0 mV/min	Signal drift for measured value acquisition.	0.5...999 mV/min, OFF
equibr.time	110 s	Equilibration time.	0...9999 s, OFF
electr.id		Electrode identification.	up to 8 characters
sample changer cal:	OFF	Calibration with sample changer.	ON, OFF
activate pulse:	OFF	Pulse output on the I/O line L6.	all, first, OFF
>statistics		Statistics calculation	
status:	OFF	Status of statistics calculation.	ON, OFF
mean	n= 2	Number n of single values for statistics calculation.	2...20
res.tab:	original	Result table for statistics calculation.	original, delete n, delete all
delete	n= 1	Delete data from sample number n.	1...20