---

# 799 GPT Titrino

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# Contenido

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## 3 Operation via RS232 Interface

### 3.1 General rules

The Titrino has an extensive remote control facility that allows full control of the Titrino via the RS 232 interface, i.e. the Titrino can receive data from an external controller or send data to an external controller.  $C_R$  and  $L_F$  are used as terminators for the data transfer. The Titrino sends  $2xC_R$  and  $L_F$  as termination of a data block, to differentiate between a data line which has  $C_R$  and  $L_F$  as terminators. The controller terminates its commands with  $C_R$  and  $L_F$ . If more than one command per line is sent by the controller, “;” is used as a separator between the individual commands.

The data are grouped logically and easy to understand. Thus e.g., for the selection of the dialog language, the following must be sent

**&Config.Aux.Language "english"**

whereby it is sufficient to only transmit the boldface characters, thus:

**&C.A.L "english"**

The quantities of the commands above are:

<b>Config</b>	configuration data
<b>Aux</b>	auxiliaries, various data
<b>Language</b>	setting the dialog language

The data are hierarchically structured (tree form). The quantities that occur in this tree are called **objects** in the following. The dialog language is an object which can be called up with the

**&Config.Aux.Language**

command.

If one is in the desired location in the tree, the value of the object can be queried.

**&Config.Aux.Language \$Q** Q means Query


The query command \$Q initiates the issuing of the value on the instrument and the value emission is triggered. Entries which start with \$, trigger something. They are thus called **triggers**.

Values of objects can not only be queried, they can also be modified. Values are always entered in quotes, for example:

**&Config.Aux.Language "english"**

### 3.1.1 Call up of objects

An excerpt from the object tree is represented below:

Rules	Example
The root of the tree is designated by &.	
The branches (levels) of a tree are marked with a dot (.) when calling up an object.	
When calling up an object, it is sufficient to give only as many letters as necessary to uniquely assign the object. If the call is not unequivocal, the first object in the series will be recognized.	Calling up the dialog language <b>&amp;Config.Aux.Language</b> or <b>&amp;C.A.L</b>
Upper- or lowercase letters may be used.	<b>&amp;C.A.L</b> or <b>&amp;c.a.l</b>
An object can be assigned a value. Values are signified at the beginning and end by quotes (""). They may contain up to 24 ASCII characters. Numerical values can contain up to 6 digits, a negative sign, and a decimal point. Numbers with more than 6 characters are not accepted; more than 4 decimal places are rounded off. For numbers <1, it is necessary to enter leading zeros.	Entering the dialog language: <b>&amp;C.A.L"english"</b>  correct entry of numbers: <b>"0.1"</b>  incorrect entry of numbers "1,5" or "+3" or ".1"
The current object remains until a new object is called.	entry of another dialog language: <b>"deutsch"</b>
New objects can be addressed relative to the old object: <b>A preceding dot</b> leads <b>forwards</b> to the next level in the tree.	From the root to node 'Aux': <b>&amp;C.A</b> Forward from node 'Aux' to 'Prog': <b>.P</b>
<b>More than one preceding dot</b> leads one level <b>backwards</b> in the tree. n node backwards require n+1 preceding dots.	Jump from node 'Prog' to node 'Aux' and select a new object 'Language' at this level: <b>..L</b>
If you must jump back to the root, enter a preceding &.	Change from node 'Language' via the root to node 'Mode': <b>&amp;M</b>

### 3.1.2 Triggers

Triggers initiate an action on the Titrimo, for example, starting a process or sending data. Triggers are marked by the introductory symbol \$.

The following triggers are possible:

<b>\$G</b>	<b>Go</b>	Starts processes, for ex. starting the mode run or setting the RS 232 interface parameters
<b>\$S</b>	<b>Stop</b>	Stops processes
<b>\$H</b>	<b>Hold</b>	Holds processes
<b>\$C</b>	<b>Continue</b>	Continues processes after Hold
<b>\$Q</b>	<b>Query</b>	Queries all information from the current node in the tree forward up to and including the values
<b>\$Q.P</b>	<b>Path</b>	Queries the path from the root of the tree up to the current node
<b>\$Q.H</b>	<b>Highest Index</b>	Queries the number of son nodes of the current node
<b>\$Q.N"i"</b>	<b>Name</b>	Queries the name of the son node with index i, $i = 1 - n$
<b>\$D</b>	<b>Detail-Info</b>	Queries the detailed status information
<b>\$U</b>	<b>qUit</b>	Aborts the data flow of the instrument, for example, after \$Q

The triggers \$G and \$S are linked to particular objects, see the summary table page 116ff.

All other triggers can be used at any time and at all locations on the object tree.

#### Examples:

Querying the value of the baud rate: **&Config.RSSet.Baud \$Q**

Querying all values of the node "RSSet": **&Config.RSSet \$Q**

Querying the path of the node "RSSet": **&Config.RSSet \$Q.P**

Start mode: **&Mode \$G**

Querying the detailed status: **\$D**

### 3.1.3 Status messages

In order to have an efficient control by an external control device, it must also be possible to query status conditions; they provide information on the status of the Titrino. The trigger \$D initiates output of the status. Status messages consist of the global status, the detailed status and eventual error messages, e.g. \$S.Mode.SET;E26. The global status informs on the activity of the process, while the detailed status conditions show the exact activity within the process.

The following **global status conditions** are possible:

\$G	Go:	The Titrino is executing the last command.
\$H	Hold:	The Titrino has been held (\$H, key <meas/hold> or by an error which effects the hold status)
\$C	Continue:	The Titrino has been restarted actively after hold
\$R	Ready:	The Titrino has executed the last command and is ready
\$S	Stop:	A process has been aborted in an "unnatural manner". e.g. stopped or aborted because there was an error.

#### Detailed status conditions

##### Status conditions of the global \$G:

\$G	.Mode.DET	.Inac:	Instrument at the beginning or at the end of a titration.
		.Req .Id1:	Instrument in the DET mode, requesting Id1 after titration start.
		.Id2:	Instrument in the DET mode, requesting Id2 after titration start.
		.Id3:	Instrument in the DET mode, requesting Id3 after titration start.
		.Smp1:	Instrument in the DET mode, requesting sample size after titration start.
		.Unit:	Instrument in the DET mode, requesting unit of sample size after titration start.
		.Start:	Instrument in the DET mode, processing the start conditions.
		.Titr:	Instrument in the DET mode, titrating.
\$G	.Mode.MET...		As DET.
\$G	.Mode.SET	.Inac:	Instrument at the beginning or at the end of a titration.
		.Req .Id1:	Instrument in the SET mode, requesting Id1 after start.
		.Id2:	Instrument in the SET mode, requesting Id2 after start.
		.Id3:	Instrument in the SET mode, requesting Id3 after start.
		.Smp1:	Instrument in the SET mode, requesting sample size after start.
		.Unit:	Instrument in the SET mode, requesting unit of sample size after start.
		.Start:	Instrument in the SET mode, processing the start conditions.
		.SET1:	Instrument in the SET mode, titrating to the first endpoint.
		.SET2:	Instrument in the SET mode, titrating to the second endpoint.
		.Cond.Ok:	Instrument in the SET, conditioning, endpoint reached (after the first startup from the standby mode).
		.Cond.Prog:	Instrument in the SET mode, conditioning, endpoint not reached (Conditioning progressing).
\$G	.Mode.KFT...		As SET.
\$G	.Mode.STAT	.Inac:	Instrument at the beginning or at the end of a determination.
		.Req .Id1:	Instrument in the STAT mode, requesting Id1 after start.
		.Id2:	Instrument in the STAT mode, requesting Id2 after start.
		.Id3:	Instrument in the STAT mode, requesting Id3 after start.
		.Smp1:	Instrument in the STAT mode, requesting sample size after start.
		.Unit:	Instrument in the STAT mode, requesting unit of sample size after start.
		.Start:	Instrument in the STAT mode, processing the start conditions.
		.Titr:	Instrument in the STAT mode, controlling.
		.Wait:	Instrument in the STAT mode, waiting because a monitored limit has been violated.
\$G	.Mode.DOS...		As STAT.

\$G .Mode.DOC...	As STAT.
\$G .Mode.MEAS .Inac:	Instrument at the beginning or at the end of a titration.
.Req .Id1:	Instrument in the MEAS mode, requesting Id1 after start.
.Id2:	Instrument in the MEAS mode, requesting Id2 after start.
.Id3:	Instrument in the MEAS mode, requesting Id3 after start.
.Smp1:	Instrument in the MEAS mode, requesting sample size after start.
.Unit:	Instrument in the MEAS mode, requesting unit of sample size after start.
.Meas:	Instrument in the MEAS mode, measuring.
\$G .Mode.CAL .Inac:	Instrument at the beginning or at the end of a calibration
.Req .Temp:	Instrument in the CAL mode, requesting calibration temperature.
.Meas .Temp:	Instrument in the CAL mode, measuring calibration temperature.
.Req .Buf1:	Instrument in the CAL mode, requesting pH of buffer 1.
.Meas .Buf1:	Instrument in the CAL mode, measures buffer 1.
.Req .Buf2:	Instrument in the CAL mode, requesting pH of buffer 2.
.Meas .Buf2:	Instrument in the CAL mode, measures buffer 2.
etc.	
\$G .Assembly.Bur .Fill:	Buret in filling process
.ModeDis:	Buret in DIS mode
\$G .Prep.X.Active:	Preparing buret X, X=0, 1, 2.
.Empty.X.Active:	Emptying buret X, X=1, 2.

In TIP, its global status as well as the step number (X) is available.

\$G .TIP.X .Inac:	Instrument at the beginning or at the end of a TIP.
.Req .Id1:	Instrument in the TIP mode, requesting Id1 after start.
.Id2:	Instrument in the TIP mode, requesting Id2 after start.
.Id3:	Instrument in the TIP mode, requesting Id3 after start.
.Smp1:	Instrument in the TIP mode, requesting sample size after start.
.Unit:	Instrument in the TIP mode, requesting unit of sample size after start.
.Pause:	Instrument in the TIP mode, in pause.
.Info:	Instrument in the TIP mode, in info.
.Mode...:	Instrument in the TIP mode, working off a submethod. The detailed status messages of the submethod appear, see above.

#### Status conditions of the global \$H:

The status message of the action which has been held appears.  
If the process is held because a monitored limit has been violated, its status message is \$H.Mode.XXX.Titr.

#### Status conditions of the global \$C:

The status conditions of the global \$C are identical with the ones of the global status \$G. They appear when the process has been restarted actively from the status "Hold" (\$C, key <meas/hold> or automatically after elimination of an error).

**Status conditions of the global \$R:**

**\$R .Mode.XXXX.QuickMeas:** Quick manual measurement from the initial status in mode XXXX.

**\$R .Mode.DET .Inac:** Instrument in the DET mode, inactive.  
**\$R .Mode.MET .Inac:** Instrument in the MET mode, inactive.  
**\$R .Mode.SET .Inac:** Instrument in the SET mode, inactive.  
     **.Cond.Ok:** Instrument in the SET mode, conditioning, endpoint reached.  
     **.Cond.Prog:** Instrument in the SET mode, conditioning, endpoint not reached.  
**\$R .Mode.KFT...** As SET.  
**\$R .Mode.STAT .Inac:** Instrument in the STAT mode, inactive.  
**\$R .Mode.DOS .Inac:** Instrument in the DOS mode, inactive.  
**\$R .Mode.DOC .Inac:** Instrument in the DOC mode, inactive.  
**\$R .Mode.MEAS .Inac:** Instrument in the MEAS mode, inactive.  
**\$R .Mode.CAL .Inac:** Instrument in the CAL mode, inactive.  
**\$R .Assembly.Bur.ModeDis:** Buret in the DIS mode, inactive.  
**\$R .TIP .Inac:** Instrument in TIP, inactive.

**Status conditions of the global \$\$:**

**\$\$ .Mode.XXXX.QuickMeas:** Quick manual measurement from the initial status in mode XXXX.

The instrument gives the status from which it has been stopped. The detailed status information is therefore identical to for the global status \$G.  
 Violation of monitored limits with action "end" give the status message  
 \$\$ .Mode.XXX.Inac;EYYY.

### 3.1.4 Error messages

Error messages are added to the status messages and separated from them by the sign ";".

<b>E8</b>	Card read/write error. Exit: Send new command.
<b>E9</b>	Wrong card, a card has been removed/inserted during the inquiry. Exit: Send new command.
<b>E10</b>	The card has lost data. Exit: Send new command.
<b>E18</b>	Card battery low (it is between 2.37...2.64 V). Exit: Send new command.
<b>E20</b>	Check exchange unit. Exit: Mount Exchange Unit (properly) or &m \$\$.
<b>E21</b>	Check electrode, short circuit. Exit: Rectify fault or &m \$\$.
<b>E22</b>	Check electrode, break. Exit: Rectify fault or &m \$\$.
<b>E23</b>	Division by zero. Exit: The error message disappears on next startup or on recalculation.
<b>E24</b>	Check drive unit. Exit: Connect drive unit (correctly) or &m \$\$.
<b>E26</b>	Manual stop. Exit: The error message disappears on next startup.
<b>E27</b>	Stop V reached in SET, STAT, DOS or DOC. Exit: The error message disappears on next startup.
<b>E28</b>	Wrong object call up Exit: Send correct path for object. Start path at root.
<b>E29</b>	Wrong value or no value allowed. Exit: Send correct value or call up new object.
<b>E30</b>	Wrong trigger, this trigger is not allowed or carrying-out of action not possible. Exit: Send correct trigger (exception: \$D) or call up new object.
<b>E31</b>	Command is not possible in active status. Repeat command in inactive status. Exit: Send new command.
<b>E32</b>	Command is not possible during titration. Repeat command during the conditioning phase or in inactive status. Exit: Send new command.
<b>E33</b>	Value has been corrected automatically. Exit: Send new command.
<b>E34</b>	Instrument at the end of the titration and sample data is edited; the instrument at rest or editing during filling. Exit: &m \$\$.

**RS receive errors:**

- E36** Parity  
Exit: <QUIT> and ensure settings of appropriate parameters at both devices are the same.
- E37** Framing error  
Exit: <QUIT> and ensure settings of appropriate parameters at both devices are the same.
- E38** Overrun error. At least 1 character could not be read.  
Exit: <QUIT>
- E39** The internal working-off buffer of the Titrino is full (>82 characters).  
Exit: <QUIT>

**RS send errors:**

- E42** CTS=OFF No proper handshake for more than 1 s.  
Exit: <QUIT> Is the receiver switched on and ready to receive?
- E43** The transmission of the Titrino has been interrupted with XOFF for at least 6 s.  
Exit: Send XON or <QUIT>
- E45** The receive buffer of the Titrino contains an incomplete command ( $L_f$  missing). Sending from the Titrino is therefore blocked.  
Exit: Send  $L_f$  or <QUIT>.
- E120** Overrange of the primary measured value (pH, U,  $I_{pol}$ ,  $U_{pol}$  or T with MEAS T). The secondary measured value (temperature) may be instable as well.  
Exit: Correct error or &m \$\$.
- E121** Measuring point list overflow (more than 500 measuring points).  
Exit: The error message disappears on next startup.
- E122** EP overflow.  
Exit: The error message disappears on next startup or on recalculation.
- E123** Missing EP for calculation.  
Exit: The error message disappears on next startup or on recalculation.
- E124** Number of EP does not correspond with the set windows.  
Exit: The error message disappears on next startup or on recalculation.
- E125** Missing fix EP for calculation, has not been defined.  
Exit: The error message disappears on next startup or on recalculation.
- E126** Fix-EP outside of measuring point list.  
Exit: The error message disappears on next startup or on recalculation.
- E128** No new mean.  
Exit: The error message disappears on next startup or on recalculation.
- E129** No new common variable, old value remains.  
Exit: The error message disappears on next startup or on recalculation.

- E130** Wrong sample. For SET, KFT or DOC with preset titration direction the first measured value lies behind the endpoint.  
Exit: The error message disappears on next startup.
- E131** No EP set for SET, STAT.  
Exit: The error message disappears on next startup.
- E132** Silo empty and it has been started with open silo or empty silo has been opened.  
Exit: Send a silo entry.
- E133** Silo full.  
Exit: Send new command.
- E134** No method. A method, which is required from the silo memory or in TIP, does not exist.  
Exit: The error message disappears on next startup.
- E135** Check temp.sensor in MEAS T or with activated temperature monitoring.  
Exit: Correct error or &m \$\$.
- E136** Same buffer in CAL. Measured value of the second buffer differs less than 6 mV from the measured value of the first buffer.  
Exit: Correct error or &m \$\$.
- E137** XXX Bytes are missing so that the method, the silo line could not be stored or not enough RAM for running TIP.  
Exit: Send new command.
- E149** Rate missing for calculation.  
Exit: The error message disappears on next start or on recalculation.
- E150** Number of rates does not agree with the set windows or the mean rate C80 could not be calculated.  
Exit: The error message disappears on next start or on recalculation.
- E151** Fix volume missing for calculation.  
Exit: The error message disappears on next startup or on recalculation.
- E155** No new silo result (C24 or C25).  
Exit: The error message disappears on next start or on recalculation.
- E157** No sequence defined in TIP.  
Exit: The error message disappears on next start.
- E158** A second TIP has been called up in TIP.  
Exit: The error message disappears on next start.
- E160** No new temporary variable.  
Exit: The error message disappears on next start.
- E161** Measurement range of the secondary measured value (temperature) exceeded. The primary measured value (pH, U, lpol, Upol) can also be unstable.  
Exit: Rectify error or &m \$\$.
- E162** No ramp end set in DOC.  
Exit: The error message disappears on next start or &m \$\$.
- E166** Save lines is "OFF" although a submethod of TIP includes an assignment to C24 or C25.  
Exit: The error message disappears on next start. Attention: The data of this sample will not be stored.
- E167** Rate too high in DOS. No dispensing possible with the Exchange Unit currently mounted.  
Exit: The error message disappears on next start or &m \$\$.
- E170** Fix time missing for calculation.  
Exit: The error message disappears on next start or on recalculation.

<b>E171</b>	Rate too low in DOS. No dispensing possible with the Exchange Unit currently mounted. Exit: The error message disappears on next start or &m \$\$.
<b>E172</b>	In TIP or DOS a QuickMeas was started, without defining a measuring quantity. Exit: The error message disappears on next start or &Mode.QuickMeas \$\$.
<b>E173</b>	The warning interval of the internal buret D0 called. Exit: Execute prep &a.b.p \$G or start next method.
<b>E174</b>	The warning interval of the external buret D1 called. Exit: Execute prep &a.b.p \$G or start next method.
<b>E175</b>	The warning interval of the external buret D2 called. Exit: Execute prep &a.b.p \$G or start next method.
<b>E176</b>	The function &Assembly.Buret.Prepare or &Assembly.Buret.Empty was interrupted manually. Exit: The error message disappears on next start.
<b>E177</b>	Accessing to the memory card, the card was not (properly) inserted. Exit: The error message disappears on next start.
<b>E178</b>	The date of changing the battery of the card is expired. Exit: The error message disappears on next start.
<b>E180</b>	Memory card write-protected. Exit: Send new command.
<b>E181</b>	Memory card not formatted. Exit: Send new command.
<b>E182</b>	Memory card not accessible. Exit: Send new command.
<b>E183</b>	A directory with the same name exists already on the memory card. Exit: Send new command.
<b>E184</b>	Measured value below lower limit. Exit: The error message disappears when again within the limits or &m \$\$.
<b>E185</b>	Measured value above upper limit. Exit: The error message disappears when again within the limits or &m \$\$.
<b>E186</b>	Temperature below lower limit. Exit: The error message disappears when again within the limits or &m \$\$.
<b>E187</b>	Temperature above upper limit. Exit: The error message disappears when again within the limits or &m \$\$.
<b>E188</b>	Rate below lower limit. Exit: The error message disappears when again within the limits or &m \$\$.
<b>E189</b>	Rate above upper limit. Exit: The error message disappears when again within the limits or &m \$\$.
<b>E196</b>	Result is out of limits. Exit: The error message disappears on next start or on recalculation.
<b>E197</b>	Sample size is out of limits. Exit: The error message disappears on next start or on introduction of new sample size.
<b>E198</b>	Validation interval is expired. Exit: The error message disappears on next start or clear counter with &Config.Monitoring.Validation.ClearCount \$G.
<b>E199</b>	Service date is reached. Exit: The error message disappears on next start or change date in &Config.Monitoring.Service.Date.

- E203** No Oven parameters: Oven not (correctly) connected.  
Exit: The error message disappears on next start. If you don't wish oven parameters in your report, select &Mode.Parameter.Presel.Oven "no" in your method(s).
- E205** Calibration interval is expired.  
Exit: The error message disappears on next calibration or if you delete the calibration.
- E212** Transmission error from Remote Box. Unknown characters.  
Exit: Rectify error and switch Titrino off and on again.
- E213** Time-out error from PC keyboard (Remote Box)  
Exit: Rectify error and switch Titrino off and on again.
- E214** Check Remote Box. Remote Box not (properly) connected but activated in &Config.Periph.RemoteBox.  
Exit: Rectify error and switch Titrino off and on again.
- E270** Overload in dosing element: piston is blocked.  
Exit: <STOP>, <QUIT>. Rectify fault at the wet part. Mount it again onto the dos driver and press <CLEAR>. The dosing element will be initialized.
- E282** Overload in dosing element: cock is blocked.  
Exit: <STOP>, <QUIT>. Rectify fault at the wet part. Mount it again onto the dos driver and press <CLEAR>. The dosing element will be initialized.

## 3.2 Remote control commands

### 3.2.1 Overview

The internal object tree can be divided into the following branches:

<b>&amp;</b>	<b>Root</b>
Mode	Method parameters
UserMeth	Administration of the internal user-memory for methods
MemoryCard	Administration of the memory card
Config	Instrument configuration
SmplData	Sample specific data
Hotkey	Keys with direct access
Info	Current Data
Assembly	Component data
Setup	Setting the operating mode
Diagnosis	Diagnostics program

## &Mode

Object	Description	Input range	Reference
<b>&amp; Root</b>			
<b>Mode</b>	<b>Mode</b>	\$G, \$\$, \$H, \$C	3.2.2.1.
. <b>QuickMeas</b>	Rapid meas. in basic mode	\$G, \$\$	3.2.2.2.
. <b>Select</b>	Mode selection	<b>DET</b> , MET, SET, KFT, STAT, DOS, DOC, MEAS, CAL, TIP	3.2.2.3.
. <b>DETQuantity</b>	Measured quantity for DET	<b>pH</b> , U, lpol, Upol	ditto
. <b>METQuantity</b>	Measured quantity for MET	<b>pH</b> , U, lpol, Upol	ditto
. <b>SETQuantity</b>	Measured quantity for SET	<b>pH</b> , U, lpol, Upol	ditto
. <b>MEASQuantity</b>	Measured quantity for MEAS	<b>pH</b> , U, lpol, Upol, T	ditto
. <b>KFTQuantity</b>	Measured quantity for KFT	<b>lpol</b> , Upol	ditto
. <b>STATQuantity</b>	Measured quantity for STAT	<b>pH</b> , U, lpol, Upol	ditto
. <b>DOCQuantity</b>	Measured quantity for DOC	<b>pH</b> , U, lpol, Upol	ditto
. <b>Name</b>	Name of current method	read only/read + write	3.2.2.4.
. <b>Parameter*</b>	Parameter of current mode, see below		
. <b>Def</b>	<b>Definitions for data output</b>		
. <b>Formulas</b>	<b>Calculation formulas</b>		
. <b>1</b>	for result 1		
. <b>Formula</b>	Calculation formula	special	3.2.2.5.
. <b>TextRS</b>	Text for result output	up to 8 ASCII char	ditto
. <b>Decimal</b>	Number of decimal places	0...2...5	ditto
. <b>Unit</b>	Unit for result output	up to 6 ASCII char	ditto
. <b>Limits</b>	Limits for result	<b>ON</b> , <b>OFF</b>	ditto
. <b>LoLim</b>	Lower limit	0...±999 999	ditto
. <b>UpLim</b>	Upper limit	0...±999 999	ditto
. <b>Output</b>	Output on L13	active, pulse, <b>OFF</b>	ditto
:	up to 9 results		
. <b>SiloCalc</b>	<b>Silo calculations</b>		
. <b>Assign</b>	Assignment		
. <b>C24</b>	Store as variable C24	RSX, EPX, CXX	3.2.2.6.
. <b>C25</b>	Store as variable C25	RSX, EPX, CXX	
. <b>MatchId</b>	Matching of Id's	id1, id1&2, all, <b>OFF</b>	
. <b>ComVar</b>	<b>Assignment of common variables</b>		
. <b>C30</b>	for C30	RSX, EPX, CXX, MNX	3.2.2.7.
up to C39			
. <b>Report</b>	<b>Reports at the end of determination</b>		
. <b>Assign1</b>	Output to COM1	special	3.2.2.8.
. <b>Assign2</b>	Output to COM 2	as COM1	
. <b>Mean</b>	<b>Assignment for mean calculation</b>		
. <b>1</b>	MN1		
. <b>Assign</b>	Input of variable	RSX, EPX, CXX	3.2.2.9.
:			
. <b>TempVar</b>	<b>Assignment of temporary variables</b>		
. <b>C70</b>	for C70	RSX, EPX, CXX	3.2.2.10.
up to C79			
. <b>CFmla</b>	<b>Calculation constants</b>		
. <b>1</b>	Calculation constant C01		
. <b>Value</b>	Input of value	0...±999 999	3.2.2.11.
up to C19			

*Parameter	Tree part "Parameters for DET"		
<b>.TitrPara</b>			
<b>Titration parameters</b>			
.MptDensity	Measuring point density	0... <b>4</b> ...9	3.2.2.12.
.MinIncr	Minimum increment	0... <b>10.0</b> ...999.9	ditto
.DosRate	Dispensing rate for increments	0.01...150.0, <b>max.</b>	3.2.2.13.
.SignalDrift	Drift for meas. value acquisition	depends on meas.quant.	3.2.2.14.
.UnitSigDrift	Unit of measured value drift	read only	ditto
.EquTime	Equilibrium time	0... <b>26</b> ...9999, OFF	ditto
<b>.StartV</b>			
<b>Start volume</b>			
.Type	Type of start volume	abs., rel., <b>OFF</b>	3.2.2.15.
.V	Volume for absolute start volume	<b>0</b> ...999.99	ditto
.Factor	Factor for relative start volume	<b>0</b> ...±999 999	ditto
.Rate	Dispensing rate for start volume	0.01...150.0, <b>max.</b>	ditto
.Pause	Waiting time	<b>0</b> ...999 999	3.2.2.16.
.DosUnit	Selection of the dosing unit	<b>internal D0</b> , external D1, external D2	3.2.2.17.
.MeasInput	Measuring input	<b>1</b> , 2, diff.	3.2.2.18.
.Ipol	Polarization current	0... <b>1</b> ...±127	ditto
.Upol	Polarization voltage	0... <b>400</b> ...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON, <b>OFF</b>	ditto
.Temp	Titration temperature	-170.0... <b>25.0</b> ...500.0	3.2.2.19.
<b>.StopCond</b>			
<b>Stop conditions</b>			
<b>.VStop</b>			
<b>Stop volume</b>			
.Type	Type of stop volume	<b>abs.</b> , rel., OFF	3.2.2.20.
.V	Volume for absolute stop volume	0... <b>99.99</b> ...9999.99	ditto
.Factor	Factor for relative stop volume	0...± <b>999 999</b>	ditto
.MeasStop	Stop measured value pH, U, I	depends on meas.quant.	3.2.2.21.
.UnitMStop	Unit of stop measured value	read only	ditto
.EPStop	Stop after a number of EP's	1... <b>9</b> , OFF	3.2.2.22.
.FillRate	Filling rate	0.01...150.0, <b>max.</b>	3.2.2.23.
<b>.Statistics</b>			
<b>Statistics</b>			
.Status	Status of statistics calculation	ON, <b>OFF</b>	3.2.2.24.
.MeanN	No. of individual determinations	<b>2</b> ...20	ditto
<b>.ResTab</b>			
<b>Result table</b>			
.Select		<b>original</b> , delete n, delete all	ditto
.DelN	Deletion of individual results	<b>1</b> ...20	ditto
<b>.Evaluation</b>			
<b>Evaluation</b>			
.EPC	EP criterion	0... <b>5</b> ...200	3.2.2.25.
<b>.Recognition</b>			
<b>EP recognition</b>			
.Select	Type of EP recognition	<b>all</b> , greatest, last, window, OFF	ditto
<b>.Window</b>			
<b>Window</b>			
.1	up to 9 windows		
.LowLim	Lower limit window 1	depends on meas.quant.	ditto
.UpLim	Upper limit window 1	depends on meas.quant.	ditto
.FixEP	Fix endpoints		
.1	up to 9 fix EP's		
.Value	Measured value for fix EP1	depends on meas.quant.	ditto
.pK	pK or HNP evaluation	ON, <b>OFF</b>	ditto
<b>.Presel</b>			
<b>Preselections</b>			
.IReq	Request of Id's after start	id1, id1&2, all, <b>OFF</b>	3.2.2.26.
.Sreq	Request of smpl size after start	value, unit, all, <b>OFF</b>	ditto
<b>.LimSmplSize</b>			
<b>Limits for sample size</b>			
.Status	Status of limit control	ON, <b>OFF</b>	3.2.2.27.
.LoLim	Lower limit	<b>0.0</b> ...999 999	ditto
.UpLim	Upper limit	0.0... <b>999 999</b>	ditto
.ActPulse	Output of a pulse	ON, <b>OFF</b>	3.2.2.28.

*Parameter	Tree part "Parameters for MET"		
<b>.TitrPara</b>			
<b>Titration parameters</b>			
.VStep	Volume increment	0... <b>0.10</b> ...999.9	3.2.2.12.
.DosRate	Dispensing rate for increments	0.01...150.0, <b>max.</b>	3.2.2.13.
.SignalDrift	Drift for meas. value acquisition	depends on meas.quant.	3.2.2.14.
.UnitSigDrift	Unit of measured value drift	read only	ditto
.EquTime	Equilibrium time	0... <b>26</b> ...9999, OFF	ditto
<b>.StartV</b>			
<b>Start volume</b>			
.Type	Type of start volume	abs., rel., <b>OFF</b>	3.2.2.15.
.V	Volume for absolute start volume	0...999.99	ditto
.Factor	Factor for relative start volume	0...±999 999	ditto
.Rate	Dispensing rate for start volume	0.01...150.0, <b>max.</b>	ditto
.Pause	Waiting time	0...999 999	3.2.2.1.
.DosUnit	Selection of the dosing unit	<b>internal D0</b> , external D1, external D2	3.2.2.17.
.MeasInput	Measuring input	1, 2, diff.	3.2.2.18.
.Ipol	Polarization current	0...1...±127	ditto
.Upol	Polarization voltage	0... <b>400</b> ...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON, <b>OFF</b>	ditto
.Temp	Titration temperature	-170.0... <b>25.0</b> ...500.0	3.2.2.19.
<b>.StopCond</b>			
<b>Stop conditions</b>			
<b>.VStop</b>			
<b>Stop volume</b>			
.Type	Type of stop volume	<b>abs.</b> , rel., OFF	3.2.2.20.
.V	Volume for absolute stop volume	0... <b>99.99</b> ...9999.99	ditto
.Factor	Factor for relative stop volume	0...± <b>999 999</b>	ditto
.MeasStop	Stop measured value pH, U, I	depends on meas.quant.	3.2.2.21.
.UnitMStop	Unit of stop measured value	read only	ditto
.EPStop	Stop after a number of EP's	1... <b>9</b> , OFF	3.2.2.22.
.FillRate	Filling rate	0.01...150.0, <b>max.</b>	3.2.2.23.
<b>.Statistics</b>			
<b>Statistics</b>			
.Status	Status of statistics calculation	ON, <b>OFF</b>	3.2.2.24.
.MeanN	No. of individual determinations	<b>2</b> ...20	ditto
<b>.ResTab</b>			
<b>Result table</b>			
.Select		<b>original</b> , delete n, delete all	ditto
.DelN	Deletion of individual results	<b>1</b> ...20	ditto
<b>.Evaluation</b>			
<b>Evaluation</b>			
.EPC	EP criterion	depends on meas.quant.	3.2.2.25.
<b>.Recognition</b>			
<b>EP recognition</b>			
.Select	Type of EP recognition	<b>all</b> , greatest, last, window, OFF	ditto
<b>.Window</b>			
<b>Window</b>			
.1	up to 9 windows		
.LowLim	Lower limit window 1	depends on meas.quant.	ditto
.UpLim	Upper limit window 1	depends on meas.quant.	ditto
<b>.FixEP</b>			
<b>Fix endpoints</b>			
.1	up to 9 fix EP's		
.Value	Measured value for fix EP1	depends on meas.quant.	ditto
.pK	pK or HNP evaluation	ON, <b>OFF</b>	ditto
<b>.Presel</b>			
<b>Preselections</b>			
.IReq	Request of Id's after start	id1, id1&2, all, <b>OFF</b>	3.2.2.26.
.SReq	Request of sample size after start	value, unit, all, <b>OFF</b>	ditto
<b>.LimSmplSize</b>			
<b>Limits for sample size</b>			
.Status	Status of limit control	ON, <b>OFF</b>	ditto
.LoLim	Lower limit	<b>0.0</b> ...999 999	ditto
.UpLim	Upper limit	0.0... <b>999 999</b>	ditto
.ActPulse	Output of a pulse	ON, <b>OFF</b>	3.2.2.28.

*Parameter	Tree part "Parameters for SET"		
<b>.SET1</b>	<b>Control parameters for EP1</b>		
.EP Endpoint 1	depends on meas.quant.	3.2.2.29.	
.UnitEp	Unit of endpoint	read only	ditto
.Dyn	Dynamics	depends on meas.quant.	3.2.2.30.
.UnitDyn	Unit of dynamics	read only	ditto
.MaxRate	Maximum dosing rate	0.01... <b>10</b> ...150, max.	ditto
.MinRate	Minimum dosing rate	0.01... <b>25.0</b> ...9999	ditto
.Stop	Titration stop		
.Type	Type of stop criterion	<b>drift</b> , time	3.2.2.31.
.Drift	Stop drift	1... <b>20</b> ...999	ditto
.Time	Switch-off delay time	0... <b>10</b> ...999, inf	ditto
.StopT	Stop time	0...999 999, <b>OFF</b>	ditto
<b>.SET2</b>	<b>Control parameters for EP2, as for EP1</b>		
<b>.TitrPara</b>	<b>Titration parameters</b>		
.Direction	Titration direction	+, -, <b>auto</b>	3.2.2.32.
.XPause	Waiting time before start volume	<b>0</b> ...999 999	3.2.2.33.
.Start V	Start volume		
.Type	Type of start volume	abs., rel., <b>OFF</b>	3.2.2.15.
.V	Volume for absolute start volume	<b>0</b> ...999.99	ditto
.Factor	Factor for relative start volume	<b>0</b> ...±999 999	ditto
.Rate	Dispensing rate for start volume	0.01...150.0, <b>max.</b>	ditto
.Pause	Waiting time after start volume	<b>0</b> ...999 999	3.2.2.16.
.ExtrT	Extraction time	<b>0</b> ...999 999	3.2.2.34.
.DosUnit	Selection of the dosing unit	<b>internal D0</b> , external D1, external D2	3.2.2.17.
.MeasInput	Measuring input	<b>1</b> , <b>2</b> , diff.	3.2.2.18.
.Ipol	Polarization current	0... <b>1</b> ...±127	ditto
.Upol	Polarization voltage	0... <b>400</b> ...±1270	ditto
.PolElectrTest	Test for polarized electrodes	<b>ON</b> , <b>OFF</b>	ditto
.Temp	Titration temperature	-170.0... <b>25.0</b> ...500.0	3.2.2.19.
.TDelta	Time interv. for meas.acquisition	1... <b>2</b> ...999 999	3.2.2.35.
<b>.StopCond</b>	<b>Stop conditions</b>		
.VStop	Stop volume		
.Type	Type of stop volume	<b>abs.</b> , rel., <b>OFF</b>	3.2.2.20.
.V	Volume for absolute stop volume	0... <b>99.99</b> ...9999.99	ditto
.Factor	Factor for relative stop volume	0...± <b>999 999</b>	ditto
.FillRate	Filling rate	0.01...150.0, <b>max.</b>	3.2.2.23.
<b>.Statistics</b>	<b>Statistics</b>		
.Status	Status of statistics calculation	<b>ON</b> , <b>OFF</b>	3.2.2.24.
.MeanN	No. of individual determinations	<b>2</b> ...20	ditto
.ResTab	Result table		
.Select		<b>original</b> , delete n, delete all	ditto
.DelN	Deletion of individual results	<b>1</b> ...20	ditto
<b>.Presel</b>	<b>Preselections</b>		
.Cond	Conditioning	<b>ON</b> , <b>OFF</b>	3.2.2.36.
.DriftDisp	Display of drift during cond.	<b>ON</b> , <b>OFF</b>	ditto
.DCor	Drift correction		
.Type	Type of drift acquisition	auto, man., <b>OFF</b>	ditto
.Value	Drift value for manual drift corr.	<b>0.0</b> ...99.9	ditto
.IReq	Request of Id's after start	id1, id1&2, all, <b>OFF</b>	3.2.2.26.
.SReq	Request of smpl size after start	value, unit, all, <b>OFF</b>	ditto
.LimSmplSize	Limits for sample size		3.2.2.27
.Status	Status of limit control	<b>ON</b> , <b>OFF</b>	ditto
.LoLim	Lower limit	<b>0.0</b> ...999 999	ditto
.UpLim	Upper limit	0.0... <b>999 999</b>	ditto
.ActPulse	Output of a pulse	first, all, cond., <b>OFF</b>	3.2.2.28.

*Parameter	Tree part "Parameters for KFT"		
<b>.CtrlPara</b>			
.EP	Endpoint	depends on meas.quant.	3.2.2.37.
.UnitEp	Unit of endpoint		read only ditto
.Dyn	Dynamics	depends on meas.quant.	3.2.2.38.
.UnitDyn	Unit of dynamics		read only ditto
.MaxRate	Maximum dosing rate		0.01...150, <b>max.</b> ditto
.MinIncr	Minimum increment		0.1...9.9, <b>min.</b> ditto
<b>.Stop</b>			
.Type	Type of stop criterion		<b>drift</b> , time 3.2.2.30.
.Drift	Stop drift		1...20...999 ditto
.Time	Switch-off delay time		0...10...999, inf ditto
.StopT	Stop time		0...999 999, <b>OFF</b> ditto
<b>.TitrPara</b>			
<b>Titration parameters</b>			
.Direction	Titration direction		+, -, auto 3.2.2.31.
.XPause	Waiting time before start volume		0...999 999 3.2.2.32.
<b>.StartV</b>			
.Type	Type of start volume		abs., rel., <b>OFF</b> 3.2.2.15.
.V	Volume for absolute start volume		0...999.99 ditto
.Factor	Factor for relative start volume		0...±999 999 ditto
.Rate	Dispensing rate for start volume		0.01...150.0, <b>max.</b> ditto
.Pause	Waiting time after start volume		0...999 999 3.2.2.16.
.ExtrT	Extraction time		0...999 999 3.2.2.33.
.DosUnit	Selection of the dosing unit		<b>internal D0</b> , external D1, external D2 3.2.2.17.
.MeasInput		without meaning	
.Ipol	Polarization current		0...50...±127 3.2.2.18.
.Upol	Polarization voltage		0...400...±1270 ditto
.PolElectrTest	Test for polarized electrodes		ON, <b>OFF</b> ditto
.Temp	Titration temperature		-170.0...25.0...500.0 3.2.2.19.
.TDelta	Time interv. for meas.acquisition		1...2...999 999 3.2.2.34.
<b>.StopCond</b>			
<b>Stop conditions</b>			
.VStop	Stop volume		
.Type	Type of stop volume		<b>abs.</b> , rel., <b>OFF</b> 3.2.2.20.
.V	Volume for absolute stop volume		0...99.99...9999.99 ditto
.Factor	Factor for relative stop volume		0...±999 999 ditto
.FillRate	Filling rate		0.01...150.0, <b>max.</b> 3.2.2.23.
<b>.Statistics</b>			
<b>Statistics</b>			
.Status	Status of statistics calculation		ON, <b>OFF</b> 3.2.2.24.
.MeanN	No. of individual determinations		2...20 ditto
<b>.ResTab</b>			
.Select			<b>original</b> , delete n, delete all ditto
.DelN	Deletion of individual results		1...20 ditto
<b>.Presel</b>			
<b>Preselections</b>			
.Cond	Conditioning		<b>ON</b> , <b>OFF</b> 3.2.2.35.
.DriftDisp	Display of drift during cond.		<b>ON</b> , <b>OFF</b> ditto
<b>.DCor</b>			
.Type	Type of drift acquisition		auto, man., <b>OFF</b> ditto
.Value	Drift value for manual drift corr.		0.0...99.9 ditto
.IReq	Request of Id's after start		id1, id1&2, all, <b>OFF</b> 3.2.2.26.
.SReq	Request of smpl size after start		value, unit, all, <b>OFF</b> ditto
<b>.LimSmplSize</b>			
.Status	Status of limit control		ON, <b>OFF</b> 3.2.2.27
.LoLim	Lower limit		0.0...999 999 ditto
.UpLim	Upper limit		0.0...999 999 ditto
.Oven	KF Oven connected		COM1, COM2, <b>no</b> 3.2.2.39.
.ActPulse	Output of a pulse		first, all, cond., <b>OFF</b> 3.2.2.28.

*Parameter	Tree part "Parameters for STAT"		
<b>.CtrlPara</b>			
<b>Control parameters</b>			
.EP	Measuring point	depends on meas.quant.	3.2.2.36.
.UnitEp	Unit of endpoint	read only	ditto
.Dyn	Dynamics	depends on meas.quant.	3.2.2.37.
.UnitDyn	Unit of dynamics	read only	ditto
.MaxRate	Maximum dosing rate	0.01... <b>10</b> ...150, max.	ditto
.MinRate	Minimum dosing rate	0.01... <b>25.0</b> ...9999	3.2.2.40.
<b>.TitrPara</b>			
<b>Titration parameters</b>			
.StartV	Start volume		
.Type	Type of start volume	abs., rel., <b>OFF</b>	3.2.2.15.
.V	Volume for absolute start volume	<b>0</b> ...999.99	ditto
.Factor	Factor for relative start volume	<b>0</b> ...±999 999	ditto
.Rate	Dispensing rate for start volume	0.01...150.0, <b>max.</b>	ditto
.Pause	Waiting time after start volume	<b>0</b> ...999 999	3.2.2.16.
.Tstart	Start time	<b>0</b> ...999 999	3.2.2.41.
.StartMeas	Start measuring value	depends on meas.quant.	ditto
.RateStart	Start rate	0.01...150, <b>OFF</b>	ditto
.TDelta	Time interv.for meas.acquisition	1... <b>2</b> ...999 999	3.2.2.34.
.Direction	Titration direction	+, -, <b>auto</b>	3.2.2.31.
.DosUnit	Selection of the dosing unit	<b>internal D0</b> , external D1, external D2	3.2.2.17.
.MeasInput	Measuring input	<b>1</b> , 2, diff.	3.2.2.18.
.Ipol	Polarization current	0... <b>1</b> ...±127	ditto
.Upol	Polarization voltage	0... <b>400</b> ...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON, <b>OFF</b>	ditto
.Temp	Titration temperature	-170.0... <b>25.0</b> ...500.0	3.2.2.19.
<b>.StopCond</b>			
<b>Stop conditions</b>			
.TStop	Stop time		
.Type	Type of stop time	abs.,rel.,delta,delay, <b>OFF</b>	3.2.2.42.
.Time	Time for absolute stop time	0... <b>999 999</b>	ditto
.Factor	Factor for relative stop time	0...± <b>999 999</b>	ditto
.Delta	Time after EP is once reached	0... <b>999 999</b>	ditto
.Delay	Time after last dosing step	0... <b>999 999</b>	ditto
.VStop	Stop volume		
.Type	Type of stop volume	<b>abs.</b> , rel., <b>OFF</b>	3.2.2.20.
.V	Volume for absolute stop volume	0... <b>99.99</b> ...9999.99	ditto
.Factor	Factor for relative stop volume	0...± <b>999 999</b>	ditto
.RateStop	Stop rate	0.01...150.0, <b>OFF</b>	3.2.2.43.
.FillRate	Filling rate	0.01...150.0, <b>max.</b>	3.2.2.23.
<b>.Statistics</b>			
<b>Statistics</b>			
.Status	Status of statistics calculation	ON, <b>OFF</b>	3.2.2.24.
.MeanN	No. of individual determinations	<b>2</b> ...20	ditto
.ResTab	Result table		
.Select		<b>original</b> ,delete n,delete all	ditto
.DelN	Deletion of individual results	<b>1</b> ...20	ditto

*Parameter	Tree part "Parameters for STAT", continuation		
<b>Evaluation</b>			
.TimeWin	Time windows for rate evaluation		
.1	up to 9 windows		
.LowLim	Lower limit window 1	0...999 999, <b>OFF</b>	3.2.2.44.
.UpLim	Upper limit window 1	0...999 999, <b>OFF</b>	ditto
.FixVol	Fix volumes		
.1	up to 9 fix volumes		
.Value	Value for fix volume 1	0...999 999, <b>OFF</b>	3.2.2.45.
.FixTime	Fix times		
.1	up to 9 fix times		
.Value	Value for fix time 1	0.01...1, <b>OFF</b>	3.2.2.46.
<b>Monitoring</b>			
.MeasVal	Monitoring of measured values		
.Status	Status	<b>ON, OFF</b>	3.2.2.47.
.LowLim	Lower limit	depends on meas.quant.	ditto
.UpLim	Upper limit	depends on meas.quant.	ditto
.Action	Action if out of limits	end, hold, wait, <b>none</b>	ditto
.Rate	Monitoring of rates		
.Status	Status	<b>ON, OFF</b>	3.2.2.48.
.LowLim	Lower limit	<b>0.000...150</b>	ditto
.UpLim	Upper limit	<b>0.000...150</b>	ditto
.Action	Action if out of limits	end, hold, wait, <b>none</b>	ditto
.Temp	Monitoring of temperatures		
.Status	Status	<b>ON, OFF</b>	3.2.2.49.
.LowLim	Lower limit	<b>-170.0...500.0</b>	ditto
.UpLim	Upper limit	<b>-170.0...500.0</b>	ditto
.Action	Action if out of limits	end, hold, wait, <b>none</b>	ditto
.L10Output	Assignment of output L10		
.AssignOutput	To type of limit violation	meas,temp,rate,all, <b>none</b>	3.2.2.50.
.Limit	Assign output to violated limit	upper, lower, <b>any</b>	ditto
.Output	Signal to be set on output L10	active, <b>pulse</b>	ditto
.L11Output	as for output L10		
.L12Output	as for output L10		
.L13Output	as for output L10		
<b>Preselections</b>			
.IReq	Request of Id's after start	id1, id1&2, all, <b>OFF</b>	3.2.2.26.
.SReq	Request of sample size after start	value, unit, all, <b>OFF</b>	ditto
.LimSmplSize	Limits for sample size		3.2.2.27
.Status	Status of limit control	<b>ON, OFF</b>	ditto
.LoLim	Lower limit	<b>0.0...999 999</b>	ditto
.UpLim	Upper limit	<b>0.0...999 999</b>	ditto
.RateDisp	Display of drift during cond.	<b>ON, OFF</b>	3.2.2.51.
.ActPulse	Output of a pulse	<b>ON, OFF</b>	3.2.2.28.

*Parameter	Tree part "Parameters for DOS"		
<b>.DosPara</b>	<b>Dosing parameters</b>		
.Type	Type of dosing	<b>volume</b> , time, rate	3.2.2.52.
.Volume	Volume dosing		
.Volume	Volume	0... <b>10</b> ...99 999.99	ditto
.DisType	Second dosing criterion	time, <b>rate</b>	ditto
.Rate	Rate	0.001...150, <b>max.</b>	ditto
.Time	Dosing time	1... <b>100</b> ...999 999	ditto
.Time	Time dosing		
.Time	Dosing time	1... <b>100</b> ...999 999	ditto
.DisType	Second dosing criterion	<b>volume</b> , rate	ditto
.Rate	Rate	0.001...150, <b>max.</b>	ditto
.Volume	Volume	0... <b>10</b> ...99 999.99	ditto
.Rate	Rate dosing		
.Rate	Rate	0.001...150, <b>max.</b>	ditto
.StopType	Second dosing criterion	<b>volume</b> , time	ditto
.Time	Dosing time	1... <b>100</b> ...999 999	ditto
.Volume	Volume	0... <b>10</b> ...99 999.99	ditto
.Pause	Waiting time after start volume	<b>0</b> ...999 999	3.2.2.53.
.TDelta	Time interv.for meas.acquisition	1... <b>10</b> ...999 999	3.2.2.54.
.DosUnit	Selection of the dosing unit	<b>internal D0</b> , external D1, external D2	3.2.2.17.
.Temp	Titration temperature	-170.0... <b>25.0</b> ...500.0	3.2.2.55.
<b>.StopCond</b>	<b>Stop conditions</b>		
.VStop	Stop volume		
.Type	Type of stop volume	abs., rel., <b>OFF</b>	3.2.2.20.
.V	Volume for absolute stop volume	0... <b>99.99</b> ...9999.99	ditto
.Factor	Factor for relative stop volume	0...± <b>999 999</b>	ditto
.FillRate	Filling rate	0.01...150.0, <b>max.</b>	3.2.2.23.
<b>.Statistics</b>	<b>Statistics</b>		
.Status	Status of statistics calculation	ON, <b>OFF</b>	3.2.2.24.
.MeanN	No.of individual determinations	<b>2</b> ...20	ditto
.ResTab	Result table		
.Select		<b>original</b> , delete n, delete all	ditto
.DelN	Deletion of individual results	<b>1</b> ...20	ditto
<b>.Monitoring</b>	<b>Monitoring</b>		
.MeasVal	Monitoring of measured values		
.MeasMode	Measured value	pH, U, Ipol, Upol, <b>OFF</b>	3.2.2.56.
.MeasInput	Measuring input	<b>1</b> , 2, diff.	ditto
.PCurrent	Polarization current	0... <b>1</b> ...±127	ditto
.PVoltage	Polarization voltage	0... <b>400</b> ...±1270	ditto
.PoIElectrTest	Test for polarized electrodes	ON, <b>OFF</b>	ditto
.pH	Limits for monitoring of pH values		
.LowLim	Lower limit	<b>-20.00</b> ...20.00	ditto
.UpLim	Upper limit	<b>-20.00</b> ...20.00	ditto
.U	Limits for monitoring of U values		
.LowLim	Lower limit	<b>-2000</b> ...2000	ditto
.UpLim	Upper limit	<b>-2000</b> ...2000	ditto
.Ipol	Limits for monitoring of Ipol values		
:	as for U		
.U	Limits for monitoring of Upol values		
:	as above		
.Action	Action if out of limits	end, hold, wait, <b>none</b>	ditto

*Parameter	Tree part "Parameters for DOS, monitoring", continuation		
-.Monitoring	<b>Monitoring</b>		
-.Temp	Monitoring of temperatures		
-.Status	Status	ON, <b>OFF</b>	3.2.2.47.
-.LowLim	Lower limit	<b>-170.0</b> ...500.0	ditto
-.UpLim	Upper limit	-170.0... <b>500.0</b>	ditto
-.Action	Action if out of limits	end, hold, wait, <b>none</b>	ditto
-.L10Output	Assignment of output L10		
-.AssignOutput	To type of limit violation	meas,temp,all, <b>none</b>	3.2.2.48.
-.Limit	Assign output to violated limit	upper, lower, <b>any</b>	ditto
-.Output	Signal to be set on output L10	active, <b>pulse</b>	ditto
-.L11Output	as for output L10		
-.L12Output	as for output L10		
-.L13Output	as for output L10		
-.Presel	<b>Preselections</b>		
-.IReq	Request of Id's after start	id1, id1&2, all, <b>OFF</b>	3.2.2.26.
-.SReq	Request of smpl size after start	value, unit, all, <b>OFF</b>	ditto
-.LimSmplSize	Limits for sample size		3.2.2.27
-.Status	Status of limit control	ON, <b>OFF</b>	ditto
-.LoLim	Lower limit	<b>0.0</b> ...999 999	ditto
-.UpLim	Upper limit	0.0... <b>999 999</b>	ditto
-.ActPulse	Output of a pulse	ON, <b>OFF</b>	3.2.2.28.

*Parameter	Tree part "Parameters for DOC"		
<b>.DosPara</b>			
<b>Dosing parameters</b>			
.BeginMeas	Begin of measuring	depends on meas.quant.	3.2.2.57.
.EndMeas	End of measuring	depends on meas.quant.	ditto
.UnitMeas	Unit of measured quantity	read only	ditto
.SweepTime	Sweep time	0... <b>300</b> ...999 999	3.2.2.58.
.Dyn	Dynamics	depends on meas.quant.	3.2.2.59.
.UnitDyn	Unit of dynamics	read only	ditto
.MaxRate	Maximum dosing rate	0.01... 150, <b>max.</b>	ditto
.MinRate	Minimum dosing rate	0.01... <b>5.0</b> ...9999	ditto
.Direction	Direction	+, -, <b>auto</b>	3.2.2.60.
.StartV	Start volume		
.Type	Type of start volume	abs., rel., <b>OFF</b>	3.2.2.61.
.V	Volume for absolute start volume	0...999.99	ditto
.Factor	Factor for relative start volume	0...±999 999	ditto
.Rate	Dispensing rate for start volume	0.01...150.0, <b>max.</b>	ditto
.Pause	Waiting time after start volume	0...999 999	ditto
.TDelta	Time interv. for meas.acquisition	1... <b>10</b> ...999 999	3.2.2.52.
.DosUnit	Selection of the dosing unit	<b>internal D0</b> , external D1, external D2	3.2.2.17.
.MeasInput	Measuring input	1, 2, diff.	3.2.2.62.
.Ipol	Polarization current	0... <b>1</b> ...±127	ditto
.Upol	Polarization voltage	0... <b>400</b> ...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON, <b>OFF</b>	ditto
.Temp	Titration temperature	-170.0... <b>25.0</b> ...500.0	3.2.2.53.
<b>.StopCond</b>			
<b>Stop conditions</b>			
.VStop	Stop volume		
.Type	Type of stop volume	<b>abs.</b> , rel., OFF	3.2.2.20.
.V	Volume for absolute stop volume	0... <b>99.99</b> ...9999.99	ditto
.Factor	Factor for relative stop volume	0...± <b>999 999</b>	ditto
.FillRate	Filling rate	0.01...150.0, <b>max.</b>	3.2.2.23.
<b>.Statistics</b>			
<b>Statistics</b>			
.Status	Status of statistics calculation	ON, <b>OFF</b>	3.2.2.24.
.MeanN	No. of individual determinations	<b>2</b> ...20	ditto
.ResTab	Result table		
.Select		<b>original</b> , delete n, delete all	ditto
.DeIN	Deletion of individual results	<b>1</b> ...20	ditto
<b>.Monitoring</b>			
<b>Monitoring</b>			
.MeasVal	Monitoring of measured values		
.Status	Status	ON, <b>OFF</b>	3.2.2.45.
.LowLim	Lower limit	depends on meas.quant.	ditto
.UpLim	Upper limit	depends on meas.quant.	ditto
.Action	Action if out of limits	end, hold, wait, <b>none</b>	ditto
.Temp	Monitoring of temperatures		
.Status	Status	ON, <b>OFF</b>	3.2.2.47.
.LowLim	Lower limit	<b>-170.0</b> ...500.0	ditto
.UpLim	Upper limit	-170.0... <b>500.0</b>	ditto
.Action	Action if out of limits	end, hold, wait, <b>none</b>	ditto

*Parameter	Tree part "Parameters for DOC, monitoring", continuation		
<b>Monitoring</b>			
.L10Output	Assignment of output L10		
.AssignOutput	To type of limit violation	meas,temp,all, <b>none</b>	3.2.2.48.
.Limit	Assign output to violated limit	upper, lower, <b>any</b>	ditto
.Output	Signal to be set on output L10	active, <b>pulse</b>	ditto
.L11Output	as for output L10		
.L12Output	as for output L10		
.L13Output	as for output L10		
<b>Preselections</b>			
.IReq	Request of Id's after start	id1, id1&2, all, <b>OFF</b>	3.2.2.26.
.SReq	Request of sample size after start	value, unit, all, <b>OFF</b>	ditto
.LimSmplSize	Limits for sample size		3.2.2.27
.Status	Status of limit control	ON, <b>OFF</b>	ditto
.LoLim	Lower limit	<b>0.0</b> ...999 999	ditto
.UpLim	Upper limit	0.0... <b>999 999</b>	ditto
.RateDisp	Display of drift during cond.	ON, <b>OFF</b>	3.2.2.49.
.ActPulse	Output of a pulse	ON, <b>OFF</b>	3.2.2.28.

*Parameter	Tree part "Parameters for MEAS"		
<b>.Measuring</b>	<b>Measuring parameters</b>		
.SignalDrift	Drift for meas.value acquisition	depends on meas.quant.	3.2.2.63.
.UnitSigDrift	Unit of measured value drift	read only	ditto
.EquTime	Equilibrium time	0...9999, <b>OFF</b>	ditto
.MeasInput	Measuring input	1, 2, diff.	3.2.2.64.
.Ipol	Polarization current	0...1...±127	ditto
.Upol	Polarization voltage	0...400...±1270	ditto
.PolElectrTest	Test for polarized electrodes	ON, <b>OFF</b>	ditto
.Temp	Titration temperature	-170.0...25.0...500.0	3.2.2.65.
.TDelta	Time interv.for meas.acquisition	1...2...999 999	3.2.2.66.
<b>.Statistics</b>	<b>Statistics</b>		
.Status	Status of statistics calculation	ON, <b>OFF</b>	3.2.2.24.
.MeanN	No. of individual determinations	2...20	ditto
.ResTab	Result table		
.Select		<b>original</b> , delete n, delete all	ditto
.DelN	Deletion of individual results	1...20	ditto
<b>.Presel</b>	<b>Preselections</b>		
.IReq	Request of Id's after start	id1, id1&2, all, <b>OFF</b>	3.2.2.26.
.SReq	Request of sample size after start	value, unit, all, <b>OFF</b>	ditto
.LimSmplSize	Limits for sample size		3.2.2.27
.Status	Status of limit control	ON, <b>OFF</b>	ditto
.LoLim	Lower limit	0.0...999 999	ditto
.UpLim	Upper limit	0.0...999 999	ditto
.ActPulse	Output of a pulse	ON, <b>OFF</b>	3.2.2.28.

*Parameter	Tree part "Parameters for CAL"		
<b>.Calibration</b>	<b>Calibration parameters</b>		
.MeasInput	Measuring input	1, 2, diff.	3.2.2.67.
.CalTemp	Calibration temperature	-20.0...25.0...120.0	3.2.2.68.
.Buffer			
.1			
.Value	pH value of buffer 1	0...7.00...±20.00	3.2.2.69.
.2			
.Value	pH value of buffer 2	0...4.00...±20.00, OFF	ditto
.	up to 9 buffers		
.SignalDrift	Drift for meas.value acquisition	depends on meas.quant.	3.2.2.70.
.EquTime	Equilibrium time	0...26...9999, OFF	ditto
.Electrodel	Electrode identification	8 ASCII char.	3.2.2.71.
.SmplChanger	Calibration on a Titrimo	ON, <b>OFF</b>	3.2.2.72.
.ActPulse	Output of a pulse	first, all, <b>OFF</b>	3.2.2.73.
<b>.Statistics</b>	<b>Statistics</b>		
.Status	Status of statistics calculation	ON, <b>OFF</b>	3.2.2.24.
.MeanN	No. of individual determinations	2...20	ditto
.ResTab	Result table		
.Select		<b>original</b> , delete n, delete all	ditto
.DelN	Deletion of individual results	1...20	ditto

*Parameter	Tree part "Parameters for TIP"		
<b>.Sequence</b>	<b>Sequence</b>		
<b>.1</b>	Step 1		
<b>.Select</b>	Step selection	method,pause,L4 output,L6 output, info,prep,stirrer, <b>OFF</b>	3.2.2.74.
<b>.Method</b>	Method from mem.or card	special	3.2.2.75.
<b>.Pause</b>	Waiting time	<b>0</b> ...999 999, INF	ditto
<b>.L4Output</b>	Line L4	active,inactive,pulse, <b>OFF</b>	ditto
<b>.L6Output</b>	Line L6	active,inactive,pulse, <b>OFF</b>	ditto
<b>.Info</b>	Display information	up to 16 ASCII char.	ditto
<b>.Prep</b>	Preparation of titrating buret	<b>internal D0</b> , external D1, external D2	ditto
<b>.Stirrer</b>	Stirrer	<b>ON</b> , <b>OFF</b>	ditto
<b>:</b>	up to 30 steps		
<b>.Statistics</b>	<b>Statistics</b>		
<b>.Status</b>	Status of statistics calculation	<b>ON</b> , <b>OFF</b>	3.2.2.24.
<b>.MeanN</b>	No. of individual determinations	<b>2</b> ...20	ditto
<b>.ResTab</b>	Result table		
<b>.Select</b>		<b>original</b> ,delete n,delete all	ditto
<b>.DelN</b>	Deletion of individual results	<b>1</b> ...20	ditto
<b>.Presel</b>	<b>Preselections</b>		
<b>.IReq</b>	Request of Id's after start	id1, id1&2, all, <b>OFF</b>	3.2.2.26.
<b>.SReq</b>	Request of sample size after start	value, unit, all, <b>OFF</b>	ditto
<b>.LimSmplSize</b>	Limits for sample size		3.2.2.27
<b>.Status</b>	Status of limit control	<b>ON</b> , <b>OFF</b>	ditto
<b>.LoLim</b>	Lower limit	<b>0.0</b> ...999 999	ditto
<b>.UpLim</b>	Upper limit	<b>0.0</b> ... <b>999 999</b>	ditto
<b>.MeasMode</b>	Measuring mode for man.meas.	pH,U,lpol,Upol,T, <b>OFF</b>	3.2.2.76.
<b>.MeasInput</b>	Measuring input	<b>1</b> , <b>2</b> , diff.	ditto
<b>.Ipol</b>	Polarization current	<b>0</b> ... <b>1</b> ...±127	ditto
<b>.Upol</b>	Polarization voltage	<b>0</b> ... <b>400</b> ...±1270	ditto
<b>.PolElectrTest</b>	Test for polarized electrodes	<b>ON</b> , <b>OFF</b>	ditto
<b>.Temp</b>	Titration temperature	-170.0... <b>25.0</b> ...500.0	ditto

<b>.RSSet1</b>			
<b>Settings RS232, 1</b>			
-	<b>.Baud</b>	Baud rate	\$G 3.2.2.108. 300,600,1200,2400,4800, <b>9600</b> ,19200,38400,57600, 115200 ditto
-	<b>.DataBit</b>	Number of data bits	7, <b>8</b> ditto
-	<b>.StopBit</b>	Number of stop bits	<b>1</b> , 2 ditto
-	<b>.Parity</b>	Parity	even, odd, <b>none</b> ditto
-	<b>.Handsh</b>	Handshake	<b>HWs</b> , SWchar, SWline, none ditto
<b>.RSSet2</b> as for RS1			
<b>.ComVar</b> <b>Values of common variables</b>			
-	<b>.C30</b>	C30	<b>0... ±999 999</b> 3.2.2.109.
-	up to C39	<b>0... ±999 999</b>	
<b>.DosPrep</b> <b>Parameters for the preparation of burets</b>			
-	<b>.PowerOnPrep</b>	Warning after power on	ON, <b>OFF</b> 3.2.2.110.
-	<b>.Report</b>	Preparation report	ON, <b>OFF</b> 3.2.2.111.
-	<b>.Select</b>	Selection of dosing unit	<b>internal D0</b> , external D1, external D2 3.2.2.112.
<b>.D0</b> Internal Titrino buret			
-	<b>.WarnInterval</b>	Warning interval for preparation	5...9999, <b>OFF</b> 3.2.2.113.
-	<b>.V</b>	Volume	0... <b>3.5</b> ...99999.99 3.2.2.114.
-	<b>.Repeat</b>	Number of cycles	1... <b>2</b> ...9 3.2.2.115.
-	<b>.DosRate</b>	Dosing rate	0.01...150, <b>max.</b> 3.2.2.116.
-	<b>.FillRate</b>	Filling rate	0.01...150, <b>max.</b> ditto
<b>.D1</b> Buret D1			
-	<b>.WarnInterval</b>	Warning interval for preparation	5...9999, <b>OFF</b> 3.2.2.104.
-	<b>.Select</b>	Type of dosing unit	<b>Dosimat</b> , Dosino 3.2.2.117.
-	<b>.Dosimat</b>	Parameters for Dosimats	
-	<b>.V</b>	Volume	0... <b>3.5</b> ...99999.99 3.2.2.107.
-	<b>.Repeat</b>	Number of cycles	1... <b>2</b> ...9 3.2.2.108.
-	<b>.DosRate</b>	Dosing rate	0.01...150, <b>max.</b> 3.2.2.109.
-	<b>.FillRate</b>	Filling rate	0.01...150, <b>max.</b> ditto
-	<b>.Dosino</b>	Parameters for Dosinos	
-	<b>.Outlet</b>	Expelling outlet	<b>tip</b> , flask 3.2.2.118.
-	<b>.DosTubing</b>	Size of dosing tubes	
-	<b>.Length</b>	Length	1... <b>40.0</b> ...999.9 3.2.2.119.
-	<b>.Diam</b>	Diameter	0.1... <b>2.0</b> ...9.9 ditto
-	<b>.AspirTubing</b>	Size of aspirating tubes	
-	<b>.Length</b>	Length	1... <b>25.0</b> ...999.9 ditto
-	<b>.Diam</b>	Diameter	0.1... <b>2.0</b> ...9.9 ditto
-	<b>.DosRate</b>	Dosing rate	0.01...150, <b>max.</b> 3.2.2.120.
-	<b>.FillRate</b>	Filling rate	0.01...150, <b>max.</b> ditto
<b>.D2</b> Buret D2 as for D1			

## &Info

Object	Description	Input range	Reference
<b>&amp; Root</b>			
Info	<b>Current data</b>		
- .Report	<b>Transmission of formatted reports</b>	\$G	3.2.2.130.
- .Select	Report type	configuration, parameters, smpl data, statistics, silo, calib C-fmla, def, user method, <b>full</b> , short, mplist, curve, derive, comb, meas crv, temp crv, adj para, scalc full, scalc srt, prep, calc, act dir, mem card, all, ff	ditto
- .CalibrationData	<b>pH calibration data</b>	\$G	3.2.2.131.
- .Inp1	For measuring input 1		
- .pHas	Asymmetry pH	0... <b>7.00</b> ...±20.00	ditto
- .Slope	Slope of electrode	0... <b>1.000</b> ...±9.999	ditto
- .Temp	Calibration temperature	-170.0... <b>25.0</b> ...500.0	ditto
- .Date	Date of calibration	read only	ditto
- .ElectrodeId	Id of the calibrated electrode	read only	ditto
- .Inp2	For measuring input 2, as for input 1		
- .Diff	For differential input, as for input 1		
- .PrepData	<b>Preparation data</b>		
- .D0	For the internal Titrino buret		
- .Date	Preparation date	read only	3.2.2.132.
- .Time	Preparation time	read only	ditto
- .D1	For buret D1		
- .Type	Buret type	read only	ditto
- .Date	Preparation date	read only	ditto
- .Time	Preparation time	read only	ditto
- .D2	For buret D2, as for D1		
- .Checksums	<b>Checksums</b>	\$G	3.2.2.133.
- .MPList	Checksum of meas.point list	read only	ditto
- .ActualMethod	Checksum of current method	read only	ditto
- .DetermData	<b>Determination data</b>	\$G	3.2.2.134.
- .Write	Read/write for several nods	<b>ON, OFF</b>	
- .ExV	Volume of Exchange/Dosing unit	read only/read + write	ditto
- .MPList	Measuring point list		
- .1	Measuring point 1		
- .Attribute	Attribute	read only/read + write	ditto
- .X	X coordinate	read only/read + write	ditto
- .Y	Y coordinate	read only/read + write	ditto
- .Z1	Z1 coordinate	read only/read + write	ditto
- .Z2	Z2 coordinate	read only/read + write	ditto
	for each measuring point		
- .TitrResults	<b>Titration results</b>		
- .RS	Calculated results		
- .1	1 <sup>st</sup> result		
- .Value	Value	read only	3.2.2.135.
	up to 9 results		

<b>"Info", continuation</b>			
.EP	Endpoint		
.1	1 <sup>st</sup> result		
.V	Value	read only	ditto
.Meas	Measured value	read only	
.Mark	Mark if more than 1 EP per window	read only	
up to 9 EP's			
.Var	Variables C4X		
.C40	Start measured value	read only/read + write	ditto
.C41	Titration end volume	read only/read + write	
.C42	Titration time	read only/read + write	
.C43	Volume drift in SET/KFT	read only/read + write	
.C44	Titration temperature	read only/read + write	
.C45	Start volume	read only/read + write	
.C46	Asymmetry pH	read only	
.C47	Slope of electrode	read only	
.C48	Volume at maximum voltage	read only/read + write	
.C49	Volume at minimum voltage	read only/read + write	
.DTime	Time for drift corr. or dosing time	read only/read + write	
.FixEP	Fix EP or fix volume		
.51	C51		
.Value	Value	read only	3.2.2.136.
up to 59			
.pK	pK/HNP or fix times		
.61	C61		
.Value	Value	read only	ditto
up to 69			
.TempVar	Temporary variables C7X		
.C70	up to C79	read only/read + write	ditto
.TimeWin	Rates in time windows		
.81	C81		
.Mean	Rate	read only	3.2.2.137.
.Dev	Standard deviation of the rate	read only	ditto
up to 89			
.MeanRateC80	Mean rate C80		
.Mean	Rate	read only	ditto
.Dev	Standard deviation of the rate	read only	ditto
<b>Statistics values</b>			
.StatisticsVal			
.ActN	Number of results in chart	read only	3.2.2.138.
.1	1 <sup>st</sup> mean		
.Mean	Mean	read only	ditto
.Std	Absolute standard deviation	read only	ditto
.RelStd	Relative standard deviation	read only	ditto
up to 9 mean values			
<b>Values of silo calculations</b>			
.SiloCalc			
.C24	Values of variable C24		
.Name	Name	read only	3.2.2.139.
.Value	Value	read only	ditto
.Unit	Unit	read only	ditto
.C25	as for C24		
.C26	Values of variable C26		
.ActN	Number of single values	read only	ditto
.Mean	Mean value	read only	ditto
.Std	Absolute standard deviation	read only	ditto
.RelStd	Relative standard deviation	read only	ditto
.C27	as for C26		

-	<b>.ActualInfo</b>	<b>"Info", continuation</b>		
-	<b>.Inputs</b>	<b>Current data</b>		
-	<b>.Status</b>	I/O Inputs		
-	<b>.Status</b>	Line status	read only	3.2.2.140.
-	<b>.Change</b>	Change of line status	read only	ditto
-	<b>.Clear</b>	Clear change	\$G	ditto
-	<b>.Outputs</b>	as for I/O Inputs		
-	<b>.Assembly</b>	From Assembly		
-	<b>.CyclNo</b>	Cycle number	read only	3.2.2.141.
-	<b>.Counter</b>	Assembly counter	read only	3.2.2.142.
-	<b>.V</b>	Volume counter	read only	ditto
-	<b>.Clear</b>	Clears counter	\$G	ditto
-	<b>.Meas</b>	Measured value	read only	3.2.2.143.
-	<b>.Titrator</b>	From Titrator		
-	<b>.CyclNo</b>	Cycle number	read only	3.2.2.144.
-	<b>.V</b>	Volume	read only	ditto
-	<b>.Meas</b>	Measured indicator voltage	read only	ditto
-	<b>.dVdt</b>	Volume drift dV/dt	read only	ditto
-	<b>.dMeasdt</b>	Measured value drift	read only	ditto
-	<b>.dMeasdV</b>	1st deviation of titration curve	read only	ditto
-	<b>.ERC</b>	ERC from DET	read only	ditto
-	<b>.T</b>	Temp.as secondary meas.	read only	ditto
-	<b>.MeasPt</b>	Entry in measuring point list		
-	<b>.Index</b>	Index of entry	read only	3.2.2.145.
-	<b>.X</b>	X coordinate	read only	ditto
-	<b>.Y</b>	Y coordinate	read only	ditto
-	<b>.Z1</b>	Z1 coordinate	read only	ditto
-	<b>.Z2</b>	Z2 coordinate	read only	ditto
-	<b>.EP</b>	EP entry		
-	<b>.Index</b>	Index of entry	read only	ditto
-	<b>.X</b>	X coordinate	read only	ditto
-	<b>.Y</b>	Y coordinate	read only	ditto
-	<b>.Oven</b>	<b>Oven data</b>		
-	<b>.HeatTime</b>	Heating time	read only	3.2.2.146.
-	<b>.SampleTemp</b>	Sample temperature	read only	ditto
-	<b>.LowTemp</b>	Lowest temperature	read only	ditto
-	<b>.HighTemp</b>	Highest temperature	read only	ditto
-	<b>.GasFlow</b>	Gas flow	read only	ditto
-	<b>.UnitFlow</b>	Unit of gas flow	read only	ditto
-	<b>.Display</b>	<b>Display</b>		
-	<b>.L1</b>	Text line 1	up to 32 ASCII char	3.2.2.147.
-	up to line 8			
-	<b>.DelAll</b>	Delete display	\$G	ditto
-	<b>.Comport</b>	Comport		
-	<b>.Number</b>	COM where PC is connected	read only	3.2.2.148.
-	<b>.Assembly</b>	<b>Assembly</b>		
-	<b>.CycleTime</b>	Cycle time	read only	3.2.2.149.
-	<b>.ExV</b>	Volume of Exchange/Dosing unit	read only	ditto

## &Assembly

Object	Description	Input range	Reference
<b>&amp; Root</b>			
-			
- Assembly	<b>Assembly control</b>		
- .Bur	<b>Buret</b>		
- .Select	Selection of buret	<b>internal D0</b> , external D1 external D2	3.2.2.150.
- .Empty	Empties the buret	\$G,\$S,\$H,\$C	3.2.2.151.
- .Prep	Prepares the buret	\$G,\$S,\$H,\$C	3.2.2.152.
- .Rates	Rates		
- .Forward	Forward rate		
- .Select	Type of rate control	digital, <b>analog</b>	3.2.2.153.
- .Digital	Digital rate	0...150, <b>max.</b>	ditto
- .Reverse	as for forward rate		
- .Select	Type of rate control	digital, <b>analog</b>	ditto
- .Digital	Digital rate	0...150, <b>max.</b>	ditto
- .Fill	Fill	\$G,\$H,\$C	3.2.2.154.
- .ModeDis	Dispensing	\$G,\$S,\$H,\$C	3.2.2.155.
- .Select	Type of dispensing control	<b>volume</b> , time	ditto
- .V	Volume to be dispensed	0.0001... <b>0.1</b> ...9999	ditto
- .Time	Time to dispense	0.25... <b>1</b> ...86 400	ditto
- .VStop	Limit volume	0.0001...9999, <b>OFF</b>	ditto
- .AutoFill	Filling after each increment	ON, <b>OFF</b>	ditto
- .Meas	<b>Measuring</b>		
- .Status	Measuring ON/OFF	ON, <b>OFF</b>	3.2.2.156.
- .MeasInput	Selection of measuring input	<b>1</b> , 2, diff., lpol, Upol, Temp	ditto
- .Ipol	Polarization current	0... <b>1</b> ...±127	ditto
- .Upol	Polarization voltage	0... <b>400</b> ...±1270	ditto
- .Outputs	<b>I/O outputs</b>		
- .AutoEOD	Automatic output of EOD	<b>ON</b> , OFF	3.2.2.157.
- .SetLines	Set I/O lines	\$G	ditto
- .LO	Signal on LO	active,inactive,pulse, <b>OFF</b>	ditto
- up to L13			
- .ResetLines	Reset I/O lines	\$G	ditto
- .Stirrer	<b>Stirrer</b>	ON, <b>OFF</b>	3.2.2.158.

## &Setup

Object	Description	Input range	Reference
<b>&amp; Root</b>			
.			
└ Setup			
	<b>Settings for the operating mode</b>		
└ .Comport	Output of automatic info	1,2,1&2	3.2.2.159.
└ .Keycode	Send key code	ON, OFF	3.2.2.160.
└ .Tree	Sending format of path info		
└ .Short	Short format of path	ON, OFF	3.2.2.161.
└ .ChangedOnly	Paths of modified nodes only	ON, OFF	ditto
└ .Trace	Message on changed values	ON, OFF	3.2.2.162.
	<b>Lock key functions</b>		
└ .Lock			
└ .Keyboard	Lock all keyboard keys	ON, OFF	3.2.2.163.
└ .Config	Lock <CONFIG> key	ON, OFF	ditto
└ .Parameter	Lock <PARAM> key	ON, OFF	ditto
└ .SmplData	Lock <SMPL DATA> key	ON, OFF	ditto
└ .UserMeth	Lock functions		
└ .Recall	Lock "loading"	ON, OFF	ditto
└ .Store	Lock "saving"	ON, OFF	ditto
└ .Delete	Lock "deletion"	ON, OFF	ditto
└ .Display	Lock display function	ON, OFF	ditto
	<b>Setting waiting intervals</b>		
└ .Mode			
└ .StartWait	Waiting time after start	ON, OFF	3.2.2.164.
└ .FinWait	Waiting time after run	ON, OFF	ditto
	<b>Automatic sending of measured values</b>		
└ .SendMeas			
└ .SendStatus	Connect/disconnect sending	ON, OFF	3.2.2.165.
└ .Interval	Time interval	0.08...4...16200, MPList	ditto
└ .Select	Selection	Assembly, <b>Titrator</b>	3.2.2.166.
└ .Assembly	From assembly		
└ .CyclNo	Cycle number	ON, OFF	3.2.2.167.
└ .V	Volume	ON, OFF	ditto
└ .Meas	Measured indicator voltage	ON, OFF	ditto
└ .Titrator	From Titrator		
└ .CyclNo	Cycle number	ON, OFF	3.2.2.168.
└ .V	Volume	ON, OFF	ditto
└ .Meas	Measured indicator voltage	ON, OFF	ditto
└ .dVdt	Volume drift dV/dt	ON, OFF	ditto
└ .dMeasdt	Measured value drift	ON, OFF	ditto
└ .dMeasdV	1st deviation of titration curve	ON, OFF	ditto
└ .ERC	ERC from DET	ON, OFF	ditto
└ .T	Temp.as secondary meas.	ON, OFF	ditto

	<b>"Setup", continuation</b>		
.AutolInfo	<b>Automatic message for changes</b>		3.2.2.169.
.Status	Switch AutolInfo on/off	ON, <b>OFF</b>	ditto
.P	When mains is switched on	ON, <b>OFF</b>	ditto
.T	<b>Titration infos</b>		
.R	When "ready"	ON, <b>OFF</b>	ditto
.G	When method started	ON, <b>OFF</b>	ditto
.GC	When start is initiated	ON, <b>OFF</b>	ditto
.S	When stopped	ON, <b>OFF</b>	ditto
.B	Begin of method	ON, <b>OFF</b>	ditto
.F	End of process	ON, <b>OFF</b>	ditto
.E	Error	ON, <b>OFF</b>	ditto
.H	When "hold"	ON, <b>OFF</b>	ditto
.C	Continue after "hold"	ON, <b>OFF</b>	ditto
.O	Conditioning OK	ON, <b>OFF</b>	ditto
.N	Conditioning not OK	ON, <b>OFF</b>	ditto
.Re	Request after start	ON, <b>OFF</b>	ditto
.Si	Silo empty	ON, <b>OFF</b>	ditto
.M	Entry in measuring point list	ON, <b>OFF</b>	ditto
.EP	Entry in EP list	ON, <b>OFF</b>	ditto
.RC	Recalculation of results done	ON, <b>OFF</b>	ditto
.C	<b>Comport infos</b>		
.B1	When COM1 sends a report	ON, <b>OFF</b>	ditto
.R1	When COM1 is ready again	ON, <b>OFF</b>	ditto
.B2	When COM2 sends a report	ON, <b>OFF</b>	ditto
.R2	When COM2 is ready again	ON, <b>OFF</b>	ditto
.I	Changing an I/O input	ON, <b>OFF</b>	ditto
.O	Changing an I/O output	ON, <b>OFF</b>	ditto
.Graphics	Changing the curve output		
.Grid	Grid on curve	<b>ON</b> , OFF	3.2.2.170.
.Frame	Frame on curve	<b>ON</b> , OFF	ditto
.Scale	Type of depending axis	<b>Full</b> , Auto	ditto
.Recorder	Length of axes		
.Right	Length of meas value axis	0.2... <b>0.5</b> ...1.00	ditto
.Feed	Length of paper drive axis	0.01... <b>0.05</b> ...1.00	ditto
.PowerOn	RESET (power on)	\$G	3.2.2.171.
.Initialise	Set default values	\$G	3.2.2.172.
.Select	Selection of branch	<b>ActMeth</b> , Config, Silo, Calib Assembly, Setup, All	ditto
.RamInit	Initialization of working mem.	\$G	3.2.2.173.
.InstrNo	Device Identification	\$G	3.2.2.174.
.Value	Input of device identification	8 ASCII characters	ditto

**&Diagnose**

Object	Description	Input range	Reference
<b>&amp; Root</b> . ├ Diagnose └ .Report	<b>Diagnose</b> Output of adjustment parameters	\$G	3.2.2.175.

## 3.2.2 Description of the remote control commands

**3.2.2.1. Mode** \$G, \$S, \$H, \$C  
 Start and stop (\$G, \$S) or hold of the current method (3.2.2.3) with \$H and continue with \$C.

\$G also serves to continue after inquiries of identifications and sample size after the start (see 3.2.2.26) as well as after inquiries of calibration temperature and pH values of buffers (see 3.2.2.68 and 3.2.2.69).

**3.2.2.2. Mode.QuickMeas** \$G, \$S  
 Start and stop of a measurement in the basic mode with the parameters (measured quantity, measuring input) of the current method. Corresponds to the <meas/hold> key. In TIP and DOS, the measured quantity is selected with &Mode.Parameter, see 3.2.2.76 and 3.2.2.56.  
 With an ongoing measurement, the current mode can be started. This stops the measurement automatically.

**3.2.2.3. Mode.Select** DET, MET, SET, KFT, STAT, DOS, DOC, MEAS, CAL, TIP

<b>Mode.DETQuantity</b>	<b>pH, U, Ipol, Upol</b>
<b>Mode.METQuantity</b>	<b>pH, U, Ipol, Upol</b>
<b>Mode.SETQuantity</b>	<b>pH, U, Ipol, Upol</b>
<b>Mode.MEASQuantity</b>	<b>pH, U, Ipol, Upol, T</b>
<b>Mode.KFTQuantity</b>	<b>Ipol, Upol</b>
<b>Mode.STATQuantity</b>	<b>pH, U, Ipol, Upol</b>
<b>Mode.DOCQuantity</b>	<b>pH, U, Ipol, Upol</b>

Selection of the standard mode. Mode and the measured quantity belong to the complete selection.

If a method is selected from the method memory, the nodes &Mode.Select and &Mode.XXXQuantity are overwritten with mode and measured quantity of the corresponding user method.

**3.2.2.4. Mode.Name** read only  
 Name of the current method in the working memory. \$Q sends 8 ASCII characters. Standard methods carry the name \*\*\*\*\*. The node can be set read + write, see 3.2.2.134.

**3.2.2.5. Mode.Def.Formulas.1.Formula** EPX, CXX, RSX, +, -, \*, /, (, )  
**Mode.Def.Formulas.1.TextRS** up to 8 ASCII characters  
**Mode.Def.Formulas.1.Decimal** 0...2...5  
**Mode.Def.Formulas.1.Unit** up to 6 ASCII characters  
**Mode.Def.Formulas.1.Limits** ON, OFF  
**Mode.Def.Formulas.1.LoLim** 0...±999 999  
**Mode.Def.Formulas.1.UpLim** 0...±999 999  
**Mode.Def.Formulas.1.Output** active, pulse, OFF  
**Mode.Def.Formulas.2.Formula**  
 etc. up to .9

Entry of formulas. Rules for formula entry, see page 74.  
 Example: "(EP2-EP1)\*C01/C00"

In addition to the formula, a text for result output, the number of decimal places and a unit for the result output can be selected. "No unit" is selected with the blank string.

In place of "RSX", a result name may be entered (.TextRS). This name is outputted in the report full, short, scalc full and scalc srt. It is used for the result and the corresponding mean value.

The limit control for results can also be activated. If a result is out of limit, a message appears in the result report, E196 is sent, and output line L13 can be set.

**3.2.2.6. Mode.Def.SiloCalc.Assign.C24** RSX, EPX, CXX  
**Mode.Def.SiloCalc.Assign.C25** RSX, EPX, CXX  
**Mode.Def.SiloCalc.MatchId** id1, id1&2, all, **OFF**

.Assign.C2X: Assignment to store results in the silo as C2X.

.MatchId: Indication which sample identification(s) have to match so that the results can be combined.

**3.2.2.7. Mode.Def.ComVar.C30** RSX, MNX, EPX, CXX  
**Mode.Def.ComVar.C31**  
 etc., up to **.C39**

Assignment of common variables.

The values of the common variables are to be found in &Config.ComVar. They can be viewed and entered there, see 3.2.2.109.

**3.2.2.8. Mode.Def.Report.Assign1**  
**Mode.Def.Report.Assign2**  
 DET: full, short, mplist, curve, derive comb, scalc full, scalc srt, calc, param, calib, ff  
 MET, SET, MEAS: full, short, mplist, curve, scalc full, scalc srt, calc, param, calib, ff  
 KFT: full, short, mplist, curve, scalc full, scalc srt, calc, param, ff  
 CAL: full, short, scalc full, scalc srt, calc, param, calib, ff  
 STAT, DOC: full, short, mplist, curve, meas crv, temp crv, scalc full, scalc srt, calc, param, calib, ff  
 DOS: full, short, mplist, curve, meas crv, temp crv, scalc full, scalc srt, calc, param, ff  
 TIP: full, short, scalc full, scalc srt, calc, param, ff

Definition of the report sequence, which is outputted automatically at the end of the determination. Entries of more than one block have to be separated with ",".

.Assign1: Output to COM1 of the Titrimo. Identical for COM2.

**3.2.2.9. Mode.Def.Mean.1.Assign** RS1, RSX, EPX, CXX  
**Mode.Def.Mean.2.Assign**  
 etc., up to **.9**

Assignment of the statistics calculations. Valid assignments are a requirement for statistics calculations. In addition, the statistics calculation must be switched on, see 3.2.2.24. Rules for statistics calculations see page 77.

**3.2.2.10. Mode.Def.TempVar.C70** RSX, EPX, CXX  
etc. up to **.C79**

Assignment of temporary variables in a submethod for calculations in TIP.

**3.2.2.11. Mode.CFmla**  
**Mode.CFmla.1.Value** 0...±999 999  
**Mode.CFmla.2.Value**  
etc., up to **.19**

Calculation constants specific to a method. Stored in the method memory of the Titrino. Operands specific to the sample (3.2.2.117 and 3.2.2.135) and values of common variables (3.2.2.109) on the other hand are not stored with the methods.

**3.2.2.12. Mode.Parameter.TitrPara.MptDensity** 0...4...9  
**Mode.Parameter.TitrPara.MinIncr** 0...10.0...999.9  
**Mode.Parameter.TitrPara.VStep** 0...0.10...9.999

.MptDensity: Parameter for DET: Measuring point density.

.MinIncr: Parameter for DET: Minimum increment in uL. If the minimum increment is set to 0, measured values are stored vs. time.

.VStep: Parameter for MET: Volume increment in mL. With "0", there is no dispensing and measured values vs. time are entered in the measuring point list.

**3.2.2.13. Mode.Parameter.TitrPara.DosRate** 0.01...150, **max.**  
Parameters for DET and MET: Dispensing rate for the volume increments in mL/min. Max. means maximum possible dispensing rate with the Exchange Unit in current use.

**3.2.2.14. Mode.Parameter.TitrPara.SignalDrift** pH,U,Ipol:0.5...50...999, OFF  
Upol: 0.05...50...99.9, OFF  
**Mode.Parameter.TitrPara.UnitSigDrift** read only  
**Mode.Parameter.TitrPara.EquTime** 0...26...9999, OFF

Parameters for DET and MET: Criteria for the measured value acquisition.

Measured value drift in mV/min (with pH, U, Ipol) or uA/min (with Upol), equilibration time in s. OFF means that the corresponding criterion is switched off. If both criteria are OFF, the measured values are acquired immediately after dispensing.

If the equilibration time has never been edited, it is automatically calculated by the instrument to match the drift, see page 20. After it has been edited once, it remains in force with the set value.

**3.2.2.15. Mode.Parameter.TitrPara.StartV.Type** abs., rel., **OFF**  
**Mode.Parameter.TitrPara.StartV.V** 0...999.99  
**Mode.Parameter.TitrPara.StartV.Factor** 0...±999 999  
**Mode.Parameter.TitrPara.StartV.Rate** 0.01...150, **max.**

Parameters for DET, MET, SET, KFT, STAT: Start volume.

If an absolute start volume (abs.) has been selected, the volume in mL is valid.

A relative start volume (rel.) is dispensed as a function of the sample size:

Start volume in mL =  $\text{smpl size} \cdot \text{factor}$

The factor is valid.

The dispensing rate in mL/min applies to both cases. Max. means maximum possible dispensing rate with the Exchange Unit in current use.

**3.2.2.16. Mode.Parameter.TitrPara.Pause** 0...999 999  
Parameters for DET, MET, SET, KFT, STAT: Pause time in s. Is waited off after the dispensing of the start volume.

**3.2.2.17. Mode.Parameter.TitrPara.DosUnit** internal D0, external D1/D2  
**Mode.Parameter.DosPara.DosUnit** internal D0, external D1/D2  
.TitrPara.DosUnit: Parameter for DET, MET, SET, KFT, STAT. Selection of dosing unit.  
.DosPara.DosUnit: Parameter for DOS, DOC. Selection of dosing unit.

**3.2.2.18. Mode.Parameter.TitrPara.MeasInput** 1, 2, diff.  
**Mode.Parameter.TitrPara.Ipol** -127...1...+127  
**Mode.Parameter.TitrPara.Upol** -1270...400...+1270  
**Mode.Parameter.TitrPara.PolElectrTest** ON, OFF

Parameters for DET, MET, SET, KFT STAT:

Selection of the measuring input; valid with measured quantities pH and U.

"diff." means differential amplifier, see page 216.

With Ipol, the inquiries for the polarization current in uA (Ipol) and

.PolElectrTest are valid.

With Upol, the inquiry for the polarization voltage in mV (Upol) is valid. Entry in steps of 10 mV.

Besides .PolElectrTest is valid.

If the test for polarized electrodes is switched on, it is performed on change-over from the inactive state to an active state (titration or conditioning).

**3.2.2.19. Mode.Parameter.TitrPara.Temp** -170.0...25.0...500.0  
Parameters for DET, MET, SET, KFT, STAT: Titration temperature in °C. If a Pt100 or Pt1000 is connected, the temperature is measured continuously and the parameter .Temp is updated.  
The temperature is used for the temperature correction in pH measurements.

**3.2.2.20. Mode.Parameter.StopCond.VStop.Type** abs., rel., OFF  
**Mode.Parameter.StopCond.VStop.V** 0...99.99...9999.99  
**Mode.Parameter.StopCond.VStop.Factor** 0...±999 999

Parameters for DET, MET, SET, KFT, STAT, DOS, DOC: Stop volume. (With DOS the default value for .Type is "OFF" and the input range for V is 0...99999.99.)

If an absolute stop volume (abs.) has been selected, the volume in mL is valid.

A relative stop volume (rel.) is dispensed as a function of the sample size:

**3.2.2.25. Mode.Parameter.Evaluation.EPC**      DET: 0...**5**...200  
 MET pH: 0.1...**0.50**...9.99  
 U, Ipol: 1...**30**...999  
 Upol: 0.1...**2**...99.9

**Mode.Parameter.Evaluation.Recognition.Selected**      **all**, greatest,  
 last, window, OFF

**Mode.Parameter.Evaluation.Recognition.Window.1.LowLim**  
 pH: 0...±20.00, **OFF**  
 U, Ipol: 0...±2000, **OFF**  
 Upol: 0...±200.0, **OFF**

**Mode.Parameter.Evaluation.Recognition.Window.1.UpLim**  
 Input range as LowLim

etc. up to 9 windows

**Mode.Parameter.Evaluation.FixEP.1.Value** pH: 0...±20.00, **OFF**  
 U; Ipol: 0...±2000, **OFF**  
 Upol: 0...±200.0, **OFF**

etc. up to 9 fix EP's

**Mode.Parameter.Evaluation.pK**      ON, **OFF**

Parameters for DET and MET: Evaluation of the EP's, see page 21.

.EPC: EP criterion in pH, in mV (with U and Ipol) or in uA (with Upol).

.Recognition.Selected: EP recognition.

all: All endpoints found are recognized.

great: Only the largest EP is recognized.

last: Only the last EP is recognized.

window: Only EP's that lie within set windows are recognized.

OFF: The EP evaluation is switched off.

.Recognition.Window.1.LowLim: Lower limit for window in pH, mV (with U and Ipol) or uA (with Upol).

.Recognition.Window.1.UpLim: Upper limit for window in pH, mV (with U and Ipol) or uA (with Upol).

Windows are opened until the lower limit is set to OFF. For every expected EP, an individual window must be set, see page 23.

.FixEP.1.Val: Fix-EP's in pH, mV (for U, Ipol) resp. uA (for Upol). Fix EP's are evaluated until the setting OFF is found.

.pK: pK or HNP evaluation. Possible only in pH and U titrations.

**3.2.2.26. Mode.Parameter.Presel.IReq**      id1, id1&2, all, **OFF**  
**Mode.Parameter.Presel.SReq**      value, unit, all, **OFF**

Parameters for DET, MET, SET, KFT, STAT, DOS, DOC, MEAS: Automatic inquiry after the start of the determination. From such an inquiry, the determination continues if the requested entry/entries is/are made, e.g.

&SmplData.OFFSilo.Id1 (see 3.2.2.122) or with &M \$G, see 3.2.2.1.

\$H is not possible in requests.

**3.2.2.27. Mode.Parameter.Presel.LimSmplSize.Status**      ON, **OFF**  
**Mode.Parameter.Presel.LimSmplSize.LoLim**      **0.0**...999 999  
**Mode.Parameter.Presel.LimSmplSize.UpLim**      **0.0**...**999 999**

Limit control for the sample size.

**3.2.2.28. Mode.Parameter.Presel.ActPuls** ON, OFF  
 for SET, KFT: first, all, cond., OFF  
 Output of a pulse on the I/O line "Activate", see page 224.

**3.2.2.29. Mode.Parameter.SET1.EP** pH: 0...±20.00, OFF  
 U, Ipol: 0...±2000, OFF  
 Upol: 0...±200.0, OFF  
**Mode.Parameter.SET1.UnitEp** read only  
 Parameters for SET: Setting the 1st endpoint as pH value, in mV (with U and Ipol) resp. uA (with Upol). The corresponding unit can be read with .UnitEP. If the value is on "OFF", no further nodes will appear from SET1.

**3.2.2.30. Mode.Parameter.SET1.Dyn** pH: 0.01...20.00, OFF  
 U, Ipol: 1...2000, OFF  
 Upol: 0.1...200.0, OFF  
**Mode.Parameter.SET1.UnitDyn** read only  
**Mode.Parameter.SET1.MaxRate** 0.01...10...150, max.  
**Mode.Parameter.SET1.MinRate** 0.01...25.0...999.9  
 Parameters for SET: Control parameters, see page 31.  
 .Dyn: Dynamics, control range in pH, mV (with U and Ipol) or uA (with Upol). The corresponding unit can be read with .UnitDyn.  
 .MaxRate: Maximum allowed titration rate in mL/min. Max. means maximum possible rate with the Exchange Unit in current use.  
 .MinRate: Minimum titration rate in ul/min.

**3.2.2.31. Mode.Parameter.SET1.Stop.Type** drift, time  
**Mode.Parameter.SET1.Stop.Drift** 1...20...999  
**Mode.Parameter.SET1.Stop.Time** 0...10...999, inf  
**Mode.Parameter.SET1.Stop.StopT** 0...99 999, OFF  
 Parameters for SET: Type and size of the stop criterion of the titration.  
 .Type: Type of stop criterion after stop drift or switch-off delay time.  
 .Drift: Stop drift in ul/min. Applies when "drift" has been selected.  
 .Time: Switch-off delay time in s. Applies when "time" has been selected. "inf" means infinite.  
 .StopT: Stop time in s. Applies when "time" has been selected and the value of .Time is set to "inf".

**3.2.2.32. Mode.Parameter.TitrPara.Direction** +, -, auto  
 Parameters for SET, KFT, STAT: Titration direction.  
 "auto" means the titration direction is determined automatically by the instrument. If 2 EP's have been set in a SET titration, the titration direction is given by the two EP's. The entry of the titration direction is then invalid.

**3.2.2.33. Mode.Parameter.TitrPara.XPause** 0...999 999  
 Parameter for SET, KFT: Pause time in s. Runs before dosing the start volume.

**3.2.2.34. Mode.Parameter.TitrPara.ExtrT** 0...999 999  
Parameter for SET, KFT: Extraction time in s.

**3.2.2.35. Mode.Parameter.TitrPara.TDelta** 1...2...999 999  
Parameter for SET, KFT, STAT: Time interval in s for the entry of a measurement point in the list of measured points.

**3.2.2.36. Mode.Parameter.Presel.Cond** ON, OFF  
**Mode.Parameter.Presel.DriftDisp** ON, OFF  
**Mode.Parameter.Presel.DCor.Type** auto, man., OFF  
**Mode.Parameter.Presel.DCor.Value** 0.0...99.9

Parameters for SET, KFT:

.Cond: Conditioning ON/OFF

.DriftDisp: Drift display during conditioning ON/OFF.

.DCor.Type: Type of drift take-over for the drift correction. auto: Take-over of the drift value at start.

.DCor.Value: Drift value for the manual drift correction.

**3.2.2.37. Mode.Parameter.CtrlPara.EP** pH: 0...±20.00, OFF  
U, Ipol: 0...±2000, OFF  
Upol: 0...±200.0, OFF  
**Mode.Parameter.CtrlPara.UnitEp** read only

Parameters for KFT, STAT: Setting of the EP, resp. Control point as a pH value, in mV (with U and Ipol) or uA (with Upol). The corresponding unit can be read with .UnitEP.

With KFT the standard value for Ipol is 250 mV, for Upol 25 uA.

**3.2.2.38. Mode.Parameter.CtrlPara.Dyn** pH: 0.01...20.00, OFF  
U, Ipol: 1...2000, OFF  
Upol: 0.1...200.0, OFF  
**Mode.Parameter.CtrlPara.UnitDyn** read only  
**Mode.Parameter.CtrlPara.MaxRate** 0.01...10...150, max.  
**Mode.Parameter.CtrlPara.MinIncr** 0.1...9.9, min.

Parameters for KFT, STAT: Control parameters.

.Dyn: Dynamics (control range) in pH, mV (with U and Ipol) or uA (with Upol). The corresponding unit can be read with .UnitDyn. Standard with KFT Ipol 100, with KFT Upol 10.

.MaxRate: Maximum allowed titration rate in ml/min. Max. means maximum possible rate with mounted Exchange Unit.

.MinIncr: Parameter for KFT. Minimum volume increment in ul.

**3.2.2.39. Mode.Parameter.Presel.Oven** COM1, COM2, no  
Parameter for KFT.

If an Oven is connected, its results will be incorporated into the result report of the Titrino.

If there is no Oven connected via RS232, this parameter has to be on "no".

**3.2.2.40. Mode.Parameter.CtrlPara.MinRate** 0.01...**25.0**...9999  
 Parameter for STAT: Control parameters.  
 .MinRate: Minimum titration rate in ul/min.

**3.2.2.41. Mode.Parameter.TitrPara.TStart** 0...999 999  
**Mode.Parameter.TitrPara.StartMeas** pH: 0.00...±20.00, **OFF**  
 U, Ipol: 0...±2000, **OFF**  
 Upol: 0.0...±200.0, **OFF**  
**Mode.Parameter.TitrPara.RateStart** 0.01...150, **OFF**  
 Parameters for STAT: Start of the measured point acquisition after the start condition has been met.  
 .TStart: Start time in s.  
 .StartMeas: Start measured value.  
 .RateStart: Start rate in mL/min. The current rate must be less than the start rate. This criterion does not become effective until 10 s after the start.

**3.2.2.42. Mode.Parameter.StopCond.TStop.Type** abs.,rel.,delta,delay,**OFF**  
**Mode.Parameter.StopCond.TStop.Time** 0...**999 999**  
**Mode.Parameter.StopCond.TStop.Factor** 0...±**999 99**  
**Mode.Parameter.StopCond.TStop.Delta** 0...**999 999**  
**Mode.Parameter.StopCond.TStop.Delay** 0...**999 999**  
 Parameters for STAT: Stop time.  
 If an absolute stop time (abs.) has been selected, the time is in s (Time).  
 A relative stop time (rel.) is calculated as a function of the sample size:  
 Stop time in s = sample size \* factor  
 The factor selected (Factor) applies.  
 delta: Time in s after control point is reached for the first time.  
 delay: Time in s after last dosing step.  
 OFF means that the criterion is not monitored.

**3.2.2.43. Mode.Parameter.StopCond.RateStop** 0.01...150, **OFF**  
 Parameters for STAT: Stop rate in ml/min. The current rate must be less than the start rate. This criterion does not become effective until 10 s after the start.

**3.2.2.44. Mode.Parameter.Evaluation.TimeWin.1.LowLim** 0...999 999, **OFF**  
**Mode.Parameter.Evaluation.TimeWin.1.UpLim** 0...999 999, **OFF**  
 Parameters for STAT: Evaluation of rates in time windows. Entry of the times in s. The evaluated rates are available as variables C8X (X = 1...9).

**3.2.2.45. Mode.Parameter.Evaluation.FixVol.1.Value** 0...999 999, **OFF**  
 Parameters for STAT: Interpolation of the volume at preset times. Entry of the times in s. The interpolated volumes are available as variables C5X.

**3.2.2.46. Mode.Parameter.Evaluation.FixTime.1.Value** 0.01...1.00, **OFF**  
 Parameters for STAT: Interpolation of the time at preset fraction of the final volume. Entry as part of V(tot). The interpolated volumes are available as variable C6X.

**3.2.2.47. Mode.Parameter.Monitoring.MeasVal.Status** ON, **OFF**  
**Mode.Parameter.Monitoring.MeasVal.LowLim** pH: 0...±20.00  
**Mode.Parameter.Monitoring.MeasVal.UpLim** U, Ipol: 0...±2000  
 Upol: 0...±200.0  
**Mode.Parameter.Monitoring.MeasVal.Action** end, hold, wait, **none**

Parameters for STAT, DOC: Monitoring of limit values. Measured values (pH, U or I) are entered in the list of measured points only if the monitoring is active.

On limit value violation, an action can be selected:

end: Abort of the determination.

hold: Interrupt dosing until the sequence is continued by manual switching (&m \$C).

wait: Interrupt dosing until the limit value is no longer violated, then continue automatically.

**3.2.2.48. Mode.Parameter.Monitoring.Rate.Status** ON, **OFF**  
**Mode.Parameter.Monitoring.Rate.LowLim** **0.000...150**  
**Mode.Parameter.Monitoring.Rate.UpLim** **0.000...150**  
**Mode.Parameter.Monitoring.Rate.Action** end, hold, wait, **none**

Parameters for STAT: Monitoring of limit values of the rate. Entry in ml/min.

On limit value violation, an action can be selected:

end: Abort of the determination.

hold: Interrupt dosing until the sequence is continued by manual switching (& m \$C).

wait: Interrupt dosing until the limit value is no longer violated, then continue automatically. On violation of the lower limit, the sequence can no longer be continued as the rate would become even less than that when dosing was interrupted.

**3.2.2.49. Mode.Parameter.Monitoring.Temp.Status** ON, **OFF**  
**Mode.Parameter.Monitoring.Temp.LowLim** **-170.0...500.0**  
**Mode.Parameter.Monitoring.Temp.UpLim** **-170.0...500.0**  
**Mode.Parameter.Monitoring.Temp.Action** end, hold, wait, **none**

Parameters for STAT, DOS, DOC: Monitoring of the limit values of the temperature in °C. Temperature values are entered in the list of measured points only if the monitoring is active.

On limit value violation, an action can be selected:

end: Abort of the determination.

hold: Interrupt dosing until the sequence is continued by manual switching (&m \$C).

wait: Interrupt dosing until the limit value is no longer violated, then continue automatically.

**3.2.2.50. Mode.Parameter.Monitoring.L10Output.AssignOutput**

meas, temp, rate, all, **none**  
**Mode.Parameter.Monitoring.L10Output.Limit** upper, lower, **any**  
**Mode.Parameter.Monitoring.L10Output.Output** active, **pulse**  
 Parameters for STAT, DOS, DOC: Output of a signal on L10 output (pin 8) of the remote socket on limit value violation.  
 The output can be assigned a definite or all monitoring value(s) as well as the violated limit and the type of signal to be outputted determined.

**3.2.2.51. Mode.Parameter.Presel.RateDisp** ON, **OFF**  
 Parameters for STAT, DOC: Display of the current rate during the determination.

**3.2.2.52. Mode.Parameter.DosPara.Type** **volume**, time, rate  
**Mode.Parameter.DosPara.Volume.Volume** 0.000...**10**...99999.99  
**Mode.Parameter.DosPara.Volume.DisType** time, **rate**  
**Mode.Parameter.DosPara.Volume.Rate** 0.001...150, **max.**  
**Mode.Parameter.DosPara.Volume.Time** 1...**100**...999 999  
**Mode.Parameter.DosPara.Time.Time** 1...**100**...999 999  
**Mode.Parameter.DosPara.Time.DisType** **volume**, rate  
**Mode.Parameter.DosPara.Time.Rate** 0.001...150, **max.**  
**Mode.Parameter.DosPara.Time.Volume** 0.000...**10**...99999.99  
**Mode.Parameter.DosPara.Rate.Rate** 0.001...150, **max.**  
**Mode.Parameter.DosPara.Rate.StopType** **volume**, time  
**Mode.Parameter.DosPara.Rate.Time** 1...**100**...999 999  
**Mode.Parameter.DosPara.Rate.Volume** 0.000...**10**...99999.99

Parameters for DOS: In the equation Rate = Volume/Time, two parameters are preset, the third is calculated. Depending on the entry of the type, the following information is required: .Volume, .Time or .Rate.  
 Entries: .Volume in ml, .Time in s, .Rate in ml/min.

**3.2.2.53. Mode.Parameter.DosPara.Pause** **0**...999 999  
 Parameters for DOS: Pause time in s.

**3.2.2.54. Mode.Parameter.DosPara.TDelta** 1...**10**...999 999  
 Parameters for DOS, DOC: Time interval for the entry of a measured point in the list of measured points.

**3.2.2.55. Mode.Parameter.DosPara.Temp** -170.0...**25.0**...500.0  
 Parameters for DOS, DOC: Temperature in °C. The temperature is measured continuously and updated if a T sensor is attached.  
 The temperature is used for temperature compensation in pH measurements.

**3.2.2.56. Mode.Parameter.Monitoring.MeasVal.MeasModepH,U,Ipol,Upol**, **OFF**  
**Mode.Parameter.Monitoring.MeasVal.MeasInput** **1**, 2, diff.  
**Mode.Parameter.Monitoring.MeasVal.PCurrent** 0...**1**...±127  
**Mode.Parameter.Monitoring.MeasVal.PVoltage** 0...**400**...±1270  
**Mode.Parameter.Monitoring.MeasVal.PolElectrTest** ON, **OFF**

**Mode.Parameter.Monitoring.MeasVal.pH.LowLim** 0...±20.00  
**Mode.Parameter.Monitoring.MeasVal.pH.UpLim** 0...±20.00  
**Mode.Parameter.Monitoring.MeasVal.U.LowLim** 0...±2000  
**Mode.Parameter.Monitoring.MeasVal.U.UpLim** 0...±2000  
**Mode.Parameter.Monitoring.MeasVal.Ipol.LowLim** 0...±2000  
**Mode.Parameter.Monitoring.MeasVal.Ipol.UpLim** 0...±2000  
**Mode.Parameter.Monitoring.MeasVal.Upol.LowLim** 0...±200.0  
**Mode.Parameter.Monitoring.MeasVal.Upol.UpLim** 0...±200.0  
**Mode.Parameter.Monitoring.MeasVal.Action** end, hold, wait, **none**

Parameters for DOS: Selection of the measured quantity for measured value monitoring. Measured values (pH, U and I) are entered in the list of measured points only when the monitoring is active.

Selection of the measurement input (MeasInput) applies with measured quantities pH and U. diff. means differential amplifier, see page 216.

With Ipol, the requests for the polarization current in uA (Ipol) and .PolElectrTest apply.

With Upol, the request for the polarization voltage in mV (Upol) applies. Entry in steps of 10 mV. .PolElectrTest also applies. If the test for polarized electrodes is switched on (ON), it is performed on the change from the inactive basic mode to an active mode.

On limit value violation, an action can be selected:

stop: Abort of the determination.

hold: Interrupt dosing until the sequence is continued by manual switching (&m \$C).

wait: Interrupt dosing until the limit value is no longer violated, then continue automatically.

**3.2.2.57. Mode.Parameter.DosPara.BeginMeas** pH: 0...±20.00, **init**  
 U, Ipol: 0...±2000, **init**  
 Upol: 0...±200.0, **init**

**Mode.Parameter.DosPara.EndMeas**

as with .Begin, but OFF instead of init

**Mode.Parameter.DosPara.UnitMeas**

read only

Parameters for DOS: Setting of ramp start and ramp end. init means ramp start with default measured value.

**3.2.2.58. Mode.Parameter.DosPara.SweepTime** 0...**300**...999 999  
 Sweep time for the ramp in s.

**3.2.2.59. Mode.Parameter.DosPara.Dyn** pH: 0.01...**0.25**...20.00, OFF  
 U: 1...**15**...2000, OFF  
 Ipol: 1...**6**...2000, OFF  
 Upol: 0.1...**10**...200.0, OFF

**Mode.Parameter.DosPara.UnitDyn**

read only

**Mode.Parameter.DosPara.MaxRate**

0.01...150, **max.**

**Mode.Parameter.DosPara.MinRate**

0.01...**5.0**...9999

Parameters for DOC: Control parameters, see page 68.

.Dyn: Dynamics in pH, mV (with U and Ipol) or uA (with Upol). The corresponding unit can be read with .UnitDyn.

.MaxRate: Maximum allowed titration rate in ml/min. Max. means maximum possible rate with the mounted Exchange Unit.

.MinRate: Minimum titration rate in ul/min.

**3.2.2.60. Mode.Parameter.DosPara.Direction** +, -, **auto**

Parameters for DOC: Titration direction.

auto means the titration direction is determined automatically by the instrument. The direction is used only when the ramp start is the same as init.

**3.2.2.61. Mode.Parameter.DosPara.StartV.Type** abs., rel., **OFF**  
**Mode.Parameter.DosPara.StartV.V** **0...999.99**  
**Mode.Parameter.DosPara.StartV.Factor** **0...±999.999**  
**Mode.Parameter.DosPara.StartV.Rate** 0.01...150, **max.**  
**Mode.Parameter.DosPara.Pause** **0...999.999**

Parameters for DOC: Start volume and pause time in s. Runs after dispensing of the start volume.

If an absolute start volume (abs.) has been selected, the volume is in ml (Volume).

A relative start volume (rel.) is dispensed as a function of the sample size:

Start volume in ml = sample size \* factor (the factor selected (Factor))

The dosing rate in ml/min holds for both cases. Max. means the maximum possible dosing rate with the Exchange Unit currently mounted.

**3.2.2.62. Mode.Parameter.DosPara.MeasInput** **1, 2, diff.**  
**Mode.Parameter.DosPara.Ipol** **0...1...±127**  
**Mode.Parameter.DosPara.Upol** **0...400...±1270**  
**Mode.Parameter.DosPara.PolElectrTest** **ON, OFF**

Parameters for DOC:

Selection of the measurement input (MeasInput) applies to the measured quantities pH and U. diff. means differential amplifier, see page 216.

With Ipol, the requests for the polarization current in uA (Ipol) and

.PolElectrTest apply.

With Upol, the request for the polarization voltage in mV (Upol) applies. Entry in steps of 10 mV. .PolElectrTest also applies.

If the test for polarized electrodes is switched on (ON), it will be performed on the change from the inactive basic mode to an active mode.

**3.2.2.63. Mode.Parameter.Measuring.SignalDrift**  
 pH, U, Ipol, T: **0.5...999, OFF**  
 Upol: **0.05...99.9, OFF**  
**Mode.Parameter.Measuring.UnitSigDrift** **read only**  
**Mode.Parameter.Measuring.EquTime** **0...9999, OFF**

Parameters for MEAS: Criteria for the measured value acquisition. Measured value drift in mV/min (with pH, U, Ipol, T), uA/min (with Upol), resp. °C/min (with T). Equilibration time in s. OFF means that the corresponding criterion is switched off. If both criteria are OFF, the measurement continues indefinitely. If the equilibration time has never been edited, it is automatically calculated by the instrument to match the drift, see page 69. After it has been edited once, it remains in force with the set value.

<b>3.2.2.64.</b>	<b>Mode.Parameter.Measuring.MeasInput</b>	<b>1, 2, diff.</b>
	<b>Mode.Parameter.Measuring.Ipol</b>	<b>±127...1...+127</b>
	<b>Mode.Parameter.Measuring.Upol</b>	<b>±1270...400...+1270</b>
	<b>Mode.Parameter.Measuring.PolElectrTest</b>	<b>ON, OFF</b>

Parameters for MEAS:

Selection of the measuring input; valid with measured quantities pH and U.

"diff." means differential amplifier, see page 216.

With Ipol, the inquiries for the polarization current in uA (Ipol) and .PolElectrTest are valid.

With Upol the inquiry for the polarization voltage in mV (Upol) is valid. Entry in steps on 10 mV.

Besides .PolElectrTest is valid.

If the test for polarized electrodes is switched on, it is performed on change-over from the inactive state to the measurement.

<b>3.2.2.65.</b>	<b>Mode.Parameter.Measuring.Temp</b>	<b>-170.0...25.0...500.0</b>
------------------	--------------------------------------	------------------------------

Parameters for MEAS: Measurement temperature in °C. If a Pt100 or Pt1000 is connected, the temperature is measured.

The temperature is used for the temperature correction in pH measurements.

<b>3.2.2.66.</b>	<b>Mode.Parameter.Measuring.TDelta</b>	<b>1...2...999 999</b>
------------------	--	------------------------

Parameters for MEAS: Time interval in s for the acquisition of a measured value into the measuring point list.

<b>3.2.2.67.</b>	<b>Mode.Parameter.Calibration.MeasInput</b>	<b>1, 2, diff.</b>
------------------	---	--------------------

Parameters for CAL: Selection of the measuring input. "diff." means differential amplifier, see page 215.

<b>3.2.2.68.</b>	<b>Mode.Parameter.Calibration.CalTemp</b>	<b>-20.0...25.0...120.0</b>
------------------	---	-----------------------------

Parameters for CAL: Calibration temperature in °C. If a Pt 100 or Pt1000 is connected, the temperature is measured.

<b>3.2.2.69.</b>	<b>Mode.Parameter.Calibration.Buffer.1.Value</b>	<b>0...7.00...±20.00</b>
------------------	--	--------------------------

**Mode.Parameter.Calibration.Buffer.2.Value** 0...4.00...±20.00, OFF etc. up to 9 buffers

Parameters for CAL: pH of buffers. The first buffer which is set to "OFF" determines the number of buffers in the calibration.

<b>3.2.2.70.</b>	<b>Mode.Parameter.Calibration.SignalDrift</b>	<b>0.5...2...999, OFF</b>
------------------	---	---------------------------

**Mode.Parameter.Calibration.EquTime** 0...110...9999, OFF

Parameters for CAL: Criteria for measured value acquisition. Measured value drift in mV/min, equilibration time in s. OFF means that the corresponding criterion is switched off. If both criterions are on OFF, the measured value is acquired immediately.

If the equilibration time has never been edited, it is automatically calculated by the instrument to match the drift, see page 20. After it has been edited once, it remains in force with the set value.

**3.2.2.71. Mode.Parameter.Calibration.Electrodeld** up to 8 ASCII char  
Parameters for CAL: Electrode identification. It is classified under calibration data, see 3.2.2.131.

**3.2.2.72. Mode.Parameter.Calibration.SmplChanger** ON, OFF  
Parameters for CAL: Calibration at Titrino.  
With "ON", there are no hold points in the calibration sequence for entries, the first buffer is measured directly.

**3.2.2.73. Mode.Parameter.Calibration.ActPulse** first, all, OFF  
Parameters for CAL: Output of a pulse on the I/O line "Activate", see page 224.

**3.2.2.74. Mode.Parameter.Sequence.X.Select** method, pause, L4 output, L6 output, info, prep, stirrer, OFF  
Parameters for TIP: Selection of an element for step X (X = 1...30). For the parameters of the elements see 3.2.2.75.

**3.2.2.75. Mode.Parameter.Sequence.X.Method** Method name  
**Mode.Parameter.Sequence.X.Pause** 0...999 999, INF  
**Mode.Parameter.Sequence.X.L4Output** active, inactive, pulse, OFF  
**Mode.Parameter.Sequence.X.L6Output** as for L4  
**Mode.Parameter.Sequence.X.Info** up to 16 ASCII characters  
**Mode.Parameter.Sequence.X.Prep** internal D0, external D1, external D2  
**Mode.Parameter.Sequence.X.Stirrer** ON, OFF

Parameters for TIP: Parameters of the elements of TIP.

.Method: Method name of a method available in the user memory or on the card. Up to 8 ASCII characters.

.Pause: Pause time in s. INF means infinite. Continue the sequence with &m \$G.

.L4 Output: Warning: A pulse triggered by the limit value monitoring at L4 (pin 3) in a submethod sets an output set to active in TIP to inactive.

.L6 Output: Warning: An activate pulse at L6 output (pin 1) in a submethod sets an output set to active in TIP to inactive.

.Info: Entry of a message which is written into the display. The sequence remains in the display with the corresponding message. Continue with &m \$G.

.Prep: Preparation of titrating buret.

.Stirrer: Switching stirrer on/off.

**3.2.2.76. Mode.Parameter.Presel.MeasMode** pH, U, Ipol, Upol, OFF  
**Mode.Parameter.Presel.MeasInput** 1, 2, diff.  
**Mode.Parameter.Presel.Ipol** 0...1...±127  
**Mode.Parameter.Presel.Upol** 0...400...±1270  
**Mode.Parameter.Presel.PolElectrTest** ON, OFF

**Mode.Parameter.Presel.Temp** -170....**25.0**...500.0

Parameters for TIP: Selection of the measured quantity for manual measurements in the inactive state, see 3.2.2.2. Selection of the measuring input (MeasInput) applies to measured quantities pH and U. "diff." means differential amplifier, see page 216. With Ipol the requests for the polarization current in uA (Ipol) and .PolElectrTest apply. With Upol the request for the polarization voltage in mV (Upol) applies. Entry in steps of 10 mV. .PolElectrTest also applies. If the test for polarized electrodes is switched on (ON), it will be performed on the change from the inactive state to an active state. The temperature applies to pH measurements.

**3.2.2.77. UserMeth.FreeMem** read only  
Memory space, available for user methods or silo lines. \$Q sends the number of free bytes, e.g.  
"4928".

**3.2.2.78. UserMeth.Recall** \$G  
**UserMeth.Recall.Name** up to 8 ASCII characters  
**UserMeth.Store** \$G  
**UserMeth.Store.Name** up to 8 ASCII characters  
**UserMeth.Delete** \$G  
**UserMeth.Delete.Name** up to 8 ASCII characters  
**UserMeth.DeIAI** \$G

Management of the internal method memory: Load, store and delete methods. An action is performed if "\$G" is sent to the corresponding node just after entering the name.

Do not use blank characters before and after method name!

.DelAll: Deletes all methods in the user memory.

**3.2.2.79. UserMeth.List.1.Name** read only  
**UserMeth.List.1.Mode** read only  
**UserMeth.List.1.Quantity** read only  
**UserMeth.List.1.DosUnit** read only  
**UserMeth.List.1.Bytes** read only  
**UserMeth.List.1.Checksum** read only  
for each method

List of the methods in the user method memory with the following characteristics:

.Name: Name of the method  
.Mode: Mode  
.Quantity: Measured quantity  
.DosUnit: Buret of the method  
.Bytes: Number of bytes of the user memory used by the method  
.Checksum: Checksum of the method, see 3.2.2.133.

**3.2.2.80. MemoryCard.Recall** \$G  
**MemoryCard.Recall.Name** up to 8 ASCII characters  
**MemoryCard.Store** \$G  
**MemoryCard.Store.Name** up to 8 ASCII characters  
**MemoryCard.Delete** \$G

**MemoryCard.Delete.Name** up to 8 ASCII characters  
 Administration of the method memory of the current directory of the memory card: load, store and delete methods. The action is carried out, if "\$G" is transmitted to the corresponding node.  
 Do not use blank characters before and after method name!

**3.2.2.81. MemoryCard.ChangeDir** \$G  
**MemoryCard.ChangeDir.Name** up to 10 ASCII characters  
**MemoryCard.ChangeDir.Checksum** \$G  
**MemoryCard.ChangeDir.Checksum.Value** read only  
 Changing the current directory. The action is carried out if "\$G" is transmitted to the corresponding node.  
 Do not use blank characters before and after name!  
 The checksum characterizes the content of the directory.

**3.2.2.82. MemoryCard.CreateDir** \$G  
**MemoryCard.CreateDir.Name** up to 10 ASCII characters  
 Delete directory. The action is carried out if "\$G" is transmitted to the corresponding node.  
 Do not use blank characters before and after name!

**3.2.2.83. MemoryCard.DeIDir** \$G  
**MemoryCard.DeIDir.Name** up to 10 ASCII characters  
 Delete directory. The action is carried out if "\$G" is transmitted to the corresponding node.  
 Do not use blank characters before and after name!

**3.2.2.84. MemoryCard.Backup** \$G  
**MemoryCard.Backup.Name** up to 10 ASCII characters  
 Backup of the internal memory onto the card. The action is carried out if "\$G" is transmitted to the corresponding node.  
 Do not use blank characters before and after name!

**3.2.2.85. MemoryCard.Reload** \$G  
**MemoryCard.Reload.Name** up to 10 ASCII characters  
 Reload a backup from the memory card into the internal memory. The action is carried out if "\$G" is transmitted to the corresponding node.  
 Do not use blank characters before and after name!

**3.2.2.86. MemoryCard.Format** \$G  
**MemoryCard.CardLabel.Name** up to 8 ASCII characters  
 Format the memory card. The action is carried out if "\$G" is transmitted to the corresponding node.  
 Do not use blank characters before and after name!

**3.2.2.87. MemoryCard.FreeMem** read only  
Free memory on the card. \$Q sends number of free bytes, e.g. "4928".

**3.2.2.88. MemoryCard.BatteryChange** \$G  
**MemoryCard.BatteryChange.Date** YYYY-MM-DD  
Date for changing battery. The date will be set with "\$G".

**3.2.2.89. MemoryCard.List.Card.1.Name** read only  
**MemoryCard.List.Card.1.Bytes** read only  
for each directory

List of all directories on the memory card with the following information:

.Name: Name of the directory  
.Bytes: Number of bytes used by the directory

**3.2.2.90. MemoryCard.List.ActDir.1.Name** read only  
**MemoryCard.List.ActDir.1.Mode** read only  
**MemoryCard.List.ActDir.1.Quantity** read only  
**MemoryCard.List.ActDir.1.DosUnit** read only  
**MemoryCard.List.ActDir.1.Bytes** read only  
**MemoryCard.List.ActDir.1.Checksum** read only  
for each method

List of all methods of the current card directory with the following information:

.Name: Name of the method  
.Mode: Mode  
.Quantity: Measured quantity  
.DosUnit: Buret of the method  
.Bytes: Number of bytes used by the method  
.Checksum: Checksum of the method, see 3.2.2.133.

**3.2.2.91. Config.Monitoring.Validation.Status** ON, OFF  
**Config.Monitoring.Validation.Interval** 1...365...9999  
**Config.Monitoring.Validation.Counter** 0...9999  
**Config.Monitoring.Validation.ClearCount** \$G

Monitoring of validation.

.Interval: Time interval in days for validation.  
.Counter: Time counter in days since last validation.  
.ClearCount: Clears the above counter.

**3.2.2.92. Config.Monitoring.Calibration.Status** ON, OFF  
**Config.Monitoring.Calibration.MeasInput** 1, 2, diff  
**Config.Monitoring.Calibration.Interval** 1...7...9999  
**Config.Monitoring.Calibration.Counter** 0...9999

Monitoring of pH calibration.

.MeasInput: Measuring input.  
.Interval: Time interval in days for calibration from the last calibration date.  
.Counter: Time counter in days since last calibration. The counter is reset to zero if a new calibration is carried out or if the calibration for the corresponding measuring input is entered manually.

- 3.2.2.93.** **Config.Monitoring.Service.Status** ON, OFF  
**Config.Monitoring.Service.Date** XXXX-XX-XX  
 Monitoring of service interval.
- 3.2.2.94.** **Config.Monitoring.DiagRep** ON, OFF  
 Printing of system test report after each switching on of the Titrino.
- 3.2.2.95.** **Config.PeriphUnit.CharSet1** Epson, Seiko, Citizen, HP, IBM  
**Config.PeriphUnit.CharSet2**  
 Selection of the character set and the graphics control characters for COM1 resp. COM2 of the Titrino.  
 IBM means the IBM character set following character set table 437 and IBM graphics control characters. Select 'IBM' for work with the computer.
- 3.2.2.96.** **Config.PeriphUnit.RepToComport** 1, 2, 1&2  
 Selection of COM of the Titrino where manually triggered reports should be outputted.
- 3.2.2.97.** **Config.PeriphUnit.Balance** Sartorius, Mettler, Mettler AT, AND, Precisa  
 Selection of the balance type.
- 3.2.2.98.** **Config.PeriphUnit.Stirrer** ON, OFF  
 Automatic stirrer control. With "ON" the stirrer will be switched on after starting a method. At the end of the method it is switched off again. Valid for DET, MET, SET, KFT, STAT, DOC.
- 3.2.2.99.** **Config.PeriphUnit.RemoteBox.Status** ON, OFF  
**Config.PeriphUnit.RemoteBox.Keyboard** US, deutsch, francais, español, schweiz.  
**Config.PeriphUnit.RemoteBox.Barcode** input, method, id1, id2, id3, smpl size  
 Connections via Remote Box.  
 .Status: Select if a Remote Box is connected.  
 .Keyboard: Type of keyboard which is connected to the Remote Box.  
 .Barcode: Select target in Titrino where you wish to have the string from the barcode reader. "input" means that the string comes into the field where the cursor is currently placed.
- 3.2.2.100.** **Config.Aux.Language** english, deutsch, francais, espanol, italiano, portugese, svenska  
 Selection of the dialog language.
- 3.2.2.101.** **Config.Aux.Set** \$G  
**Config.Aux.Set.Date** YYYY-MM-DD

- Config.Aux.Set.Time** HH:MM  
 Date and time.  
 Input format of the date: Year-month-day, two-digit, enter leading zeros.  
 Input format for the time: Hours:minutes, two-digit, enter leading zeros.  
 Date and time have to be set with &Config.Aux.Set \$G just after entry of the value.
- 3.2.2.102. Config.Aux.RunNo** 0...9999  
 Current sample number.  
 Set to 0 on power on and initialization. After 9999, counting starts again at 0.
- 3.2.2.103. Config.Aux.AutoStart** 1...9999, **OFF**  
 Number of automatic, internal starts.
- 3.2.2.104. Config.Aux.StartDelay** 0...999 999  
 Start delay time in s. During this time, the data of the preceding determination are retained.
- 3.2.2.105. Config.Aux.ResDisplay** **bold**, standard  
 Character set for the result display at the end of the determination.
- 3.2.2.106. Config.Aux.DevName** up to 8 ASCII characters  
 Name of the instrument for connections with several units. It is advisable to use only the letters A...Z (ASCII No. 65...90), a...z (ASCII No. 97...122) and the numbers 0...9 (ASCII No. 48...57) when the function Setup.AutoInfo (3.2.2.169) is used at the same time.  
 If a name has been entered, it will be printed out in the result report (full, short).
- 3.2.2.107. Config.Aux.Prog** read only  
 Output of the program version.  
 The Titrimo sends "799.0010" on requests with \$Q.
- 3.2.2.108. Config.RSSet1** \$G  
**Config.RSSet1.Baud** 300, 600, 1200, 2400, 4800, **9600**,  
19200, 38400, 57600, 115200  
**Config.RSSet1.DataBit** 7, **8**  
**Config.RSSet1.StopBit** **1**, 2  
**Config.RSSet1.Parity** even, odd, **none**  
**Config.RSSet1.Handsh** **HWs**, SWchar, SWline, none  
 \$G sets all RS settings. The changes are performed only if the instrument is inactive. After the setting of the interface parameters, wait at least 2 s to allow the components to equilibrate.  
 Settings of the values for the data transmission via the RS interface: baud rate, data bit, stop bit, parity and type of handshake, see also page 180 ff.

Baud rates >9600 need a PC which is equipped accordingly (e.g. with 16550 component).

The setting of the values must be initiated with \$G immediately after entry of the values.

**3.2.2.109. Config.ComVar.C30** 0... ±999 999  
with up to **.C39**, etc.

Values of the common variables from C30 up to C39. Insert the common variables directly or describe the determination results directly from the method, see 3.2.2.7

**3.2.2.110. Config.DosPrep.PowerOnPrep** ON, OFF  
Warning after power on, that the preparation of the buret has to be executed.

**3.2.2.111. Config.DosPrep.Report** ON, OFF  
Report output after preparation of burets.

**3.2.2.112. Config.DosPrep.Select** internal D0, external D1, external D2  
Selection of buret.

**3.2.2.113. Config.DosPrep.DX.WarnInterval** 5...9999, OFF  
Time interval in min for automatic warning that a preparation has to be executed.

**3.2.2.114. Config.DosPrep.DX(.Dosimat).V** 0...3.5...99999.99  
Volume in ml, which will be expelled at the preparation. Start of sequence see 3.2.2.152.

**3.2.2.115. Config.DosPrep.DX(.Dosimat).Repeat** 1...2...9  
Number of cycles for expelling the volume at the preparation. Start of sequence see 3.2.2.152.

**3.2.2.116. Config.DosPrep.DX(.Dosimat).DosRate** 0.01...150, max.  
**Config.DosPrep.DX(.Dosimat).FillRate** 0.01...150, max.  
Dosing and filling rate in ml/in for the preparation. Start of sequence see 3.2.2.152.

**3.2.2.117. Config.DosPrep.DX.Select** Dosimat, Dosino  
Selection of buret type. Start of sequence see 3.2.2.152.

**3.2.2.118. Config.DosPrep.DX.Dosino.Outlet** tip, flask  
Place where the liquid of the preparation is expelled. Start of sequence see 3.2.2.152.

tip: To the tip.  
flask: Back to the flask.

<b>3.2.2.119.</b>	<b>Config.DosPrep.DX.Dosino.DosTubing.Length</b>	1... <b>40</b> ...999.9
	<b>Config.DosPrep.DX.Dosino.DosTubing.Diameter</b>	0.1... <b>2</b> ...9.9
	<b>Config.DosPrep.DX.Dosino.AspirTubing.Length</b>	1... <b>25</b> ...999.9
	<b>Config.DosPrep.DX.Dosino.AspirTubing.Diameter</b>	0.1... <b>2</b> ...9.9

.DosTubing: Dosing tube. Length in cm, diameter in mm.

.AspirTubing: Aspiration tube. Length in cm, diameter in mm.

Start of sequence see 3.2.2.152.

<b>3.2.2.120.</b>	<b>Config.DosPrep.DX.Dosino.DosRate</b>	0.01...150, <b>max.</b>
	<b>Config.DosPrep.DX.Dosino.FillRate</b>	0.01...150, <b>max.</b>

Dosing and filling rate in ml/min for preparation. Start of sequence see 3.2.2.152.

<b>3.2.2.121.</b>	<b>SmplData.Status</b>	<b>ON, OFF</b>
-------------------	------------------------	----------------

On/off switching of silo memory. When the silo memory is switched on, the sample data are fetched from the lowest valid silo line.

<b>3.2.2.122.</b>	<b>SmplData.OFFSilo.Id1</b>	up to 8 ASCII characters
	<b>SmplData.OFFSilo.Id2</b>	up to 8 ASCII characters
	<b>SmplData.OFFSilo.Id3</b>	up to 8 ASCII characters
	<b>SmplData.OFFSilo.ValSmpl</b>	6-digits, sign and decimal point
	<b>SmplData.OFFSilo.UnitSmpl</b>	up to 5 ASCII characters

Current sample data.

The identifications Id1...Id3 can be used in formulas as sample-specific calculation constants C21...C23.

If "no unit" is desired for the unit of the sample size, the blank string must be entered.

<b>3.2.2.123.</b>	<b>SmplData.ONSilo.Counter.MaxLines</b>	read only
	<b>SmplData.ONSilo.Counter.FirstLine</b>	read only
	<b>SmplData.ONSilo.Counter.LastLine</b>	read only

Information on silo memory.

.MaxLines: Maximum possible number of silo lines.

.FirstLine: Lowest valid silo line.

.LastLine: Last occupied silo line.

<b>3.2.2.124.</b>	<b>SmplData.ONSilo.EditLine.1.Method</b>	up to 8 ASCII characters
	<b>SmplData.ONSilo.EditLine.1.Id1</b>	up to 8 ASCII characters
	<b>SmplData.ONSilo.EditLine.1.Id2</b>	up to 8 ASCII characters
	<b>SmplData.ONSilo.EditLine.1.Id3</b>	up to 8 ASCII characters
	<b>SmplData.ONSilo.EditLine.1.ValSmpl</b>	6-digits, sign and dec.point
	<b>SmplData.ONSilo.EditLine.1.UnitSmpl</b>	up to 5 ASCII characters
	<b>SmplData.ONSilo.EditLine.1.C24</b>	read only
	<b>SmplData.ONSilo.EditLine.1.C25</b>	read only

**SmplData.ONSilO.EditLine.1.Mark** read only  
 etc., up to **.255**  
 Contents of a silo line.  
 .Method: Method used to process the sample, from the method memory or from the card.  
 .Id: The identifications Id1...Id3 can also be used as sample-specific calculation constants C21...C23 in formulas.  
 .UnitSmpl: If "no unit" is desired for the sample size, the blank string must be entered.  
 .C24, .C25: Results which have been assigned to C24 and C25.  
 .Mark: Mark of the silo line: "\*" = deleted line, "+" = line which is worked off, "-" = line which is worked off and not valid for silo calculations (deleted), "/" = last worked-off line, where recalculation can still be done. Silo lines which have been worked off are "read only".

**3.2.2.125. SmplData.ONSilO.DeLine** \$G  
**SmplData.ONSilO.DeLine.LineNum** 1...255, **OFF**  
 Deletion of a silo line. The line # is deleted with &SmplData.ONSilO.DeLine \$G. If a formerly deleted line is edited again, it becomes valid (function "undelete").

**3.2.2.126. SmplData.ONSilO.DeAll** \$G  
 Deletes the entire silo memory. Must be triggered with \$G.

**3.2.2.127. SmplData.ONSilO.CycleLines** **ON, OFF**  
 Silo data cycling.  
 With "ON", executed lines are copied to the next free silo lines, see page 96.  
 Exercise caution if you edit the silo memory during the determinations!

**3.2.2.128. SmplData.ONSilO.SaveLines** **ON, OFF**  
 Silo lines are not deleted when they are worked off. Assigned results are stored as C24 and C25. "Save lines" can only be set to "ON" if the silo is completely empty. Delete the silo, see 3.2.2.126.

**3.2.2.129. HotKey.User.Name** up to 10 ASCII characters  
**HotKey.User.Delete** \$G  
**HotKey.User.Delete.Name** up to 10 ASCII characters  
**HotKey.User.DeAll** \$G  
**HotKey.User.List.1.Name** read only  
 Management of user names.  
 .Name: Input of user names.  
 .Delete.Name: Deletes selected user name with &HotKey.User.Delete \$G.  
 .List: List of all user names.

**3.2.2.130. Info.Report** \$G  
**Info.Report.Select** configuration, parameters, smpl data,

statistics, silo, calib, C-fmla, def, user method, **full**, short, mplist, curve, deriv, comb, meas crv, temp crv, adj para, scalc full, scalc srt, prep, calc, act dir, mem card, all, ff

\$G sends the selected report to the COM which is set in &Config.PeriphUnit.RepToComport:

configuration: Configuration report. Is not accessible during a running determination.

parameters: Parameter report of the current method. During a running determination only "live"-parameters are accessible.

smpl data: Current sample data.

statistics: Statistics table with the individual results.

silo: Contents of the silo memory.

calib: Calibration data of the measuring input in the current method.

C-fmla: Contents of the <C-fmla> key.

def: Contents of the <def> key.

user method: Contents of the method memory.

full: Full result report of the last completed determination.

short: Short result report of the last completed determination.

mplist: Measuring point list of the running determination.

curve: Titration curve of the last determination.

derive: 1st derive of titration curve of the last determination (with DET).

comb: 1st derive combined with the titration curve of the last determination (with DET).

meas crv: Curve measured value vs. time (with STAT, DOS, DOC and activated meas monitoring).

temp crv: Curve temperature vs. time (with STAT, DOS, DOC and activated temperature monitoring).

adj para: Adjustment parameters.

scalc full: Full report of the silo calculations.

scalc srt: Short report of the silo calculations.

prep: Preparation report.

calc: Calculation report of the current method.

act dir: Methods of the current directory of the memory card.

mem card: All methods of the memory card.

all: All reports.

ff: Form feed on printer.

Reports which are sent from the Titrino are marked with space (ASCII 32) and ' at the beginning. Then an individual identifier for each report follows. Reports which are triggered by RS232 (\$G) have the same introducer but without preceding space, i.e. they start with '.

**3.2.2.131. Info.CalibrationData** \$G

Info.CalibrationData.Inp1.pHas	-20.00... <b>7.00</b> ... +20.00
Info.CalibrationData.Inp1.Slope	-9.999... <b>1.000</b> ... +9.999
Info.CalibrationData.Inp1.Temp	-170.0... <b>25.0</b> ... +500.0
Info.CalibrationData.Inp1.Date	read only
Info.CalibrationData.Inp1.ElectrodeId	read only
identical for .Inp2 and .Diff	

pH calibration data for measuring input 1. After the calibration, the data are entered automatically together with the date of the calibration and the electrode identification, see 3.2.2.71.

Calibration data can be entered. They are accepted with &Info.CalibrationData \$G. If calibration data are entered, the calibration date is deleted.

**3.2.2.132.** Info.PrepData.D0.Date read only  
 Info.PrepData.D0.Time read only  
 Info.PrepData.D1.Type read only  
 Info.PrepData.D1.Date read only  
 Info.PrepData.D1.Time read only  
 identical to .D2

Data from the last correctly executed preparation.

**3.2.2.133.** Info.Checksums \$G  
 Info.Checksums.MPList read only  
 Info.Checksums.ActualMethod read only

The checksums can be used to identify the content of a file unequivocally, e.g. files with identical content

have identical results of the checksums. An empty file has checksum "0". The calculation of the checksums is triggered with \$G.

.MPList: Result of the checksum of the current measuring point list.

.ActualMethod: Result of the checksum of the current method in the working memory. Identical methods with different method names have the same results of the checksum.

**3.2.2.134.** Info.DetermData \$G  
 Info.DetermData.Write ON, OFF  
 Info.DetermData.MPList.1.Attribute read only/read + write  
 Info.DetermData.ExV read only/read + write  
 Info.DetermData.MPList.1.X read only/read + write  
 Info.DetermData.MPList.1.Y read only/read + write  
 Info.DetermData.MPList.1.Z1 read only/read + write  
 Info.DetermData.MPList.1.Z2 read only/read + write  
 for every measuring point

Determination data in hexadecimal format. A measuring point list is available in mode DET, MET, STAT, DOS, DOC, SET, KFT, and MEAS.

Recalculation of the measuring data is triggered with \$G.

.Write: With "ON", the following nodes can be overwritten:  
 &Info.DetermData.MP.List, &Info.TitrResults.Var.C4X (X = 0...5),  
 &Info.TitrResults.TempVar.C7X (X = 0...9), and &Mode.Name.

.ExV: Volume of the exchange unit, with which the determination was executed

.MPList.1.Attribute: Attribute

.MPList.X: X coordinate, time

.MPList.Y: Y coordinate, volume

.MPList.Z1: Z1 coordinate, measuring value

.MPList.Z2: Z2 coordinate, temperature

**3.2.2.135.** Info.TitrResults.RS.1.Value read only  
 etc., up to .9  
 Info.TitrResults.EP.1.V read only

**Info.TitrResults.EP.1.Meas** read only  
 etc., up to **.2**  
**Info.TitrResults.Var.C40** read only/read+ write  
 etc., up to **.C47**  
**Info.TitrResults.Var.DTime** read only/read+ write  
**.RS:** Values of the calculated results.  
**.EP:** Endpoints with DET, MET, SET, KFT:  
 Volume coordinate in mL, e.g. "1.2340"  
 Measured value coordinate in pH "5.12", mV (with U and Ipol) "-  
 241" or uA (with Upol) "43.7".  
**.Var:** Various variables. You may overwrite the variables C40...C45, see  
 3.2.2.134.  
 C40: Initial measured value in pH "5.12", mV (with U and Ipol)  
 "41", uA (with Upol) "43.7" or °C (with T) "25.0". In MEAS final  
 measured value.  
 C41: End volume with SET, STAT, DOS, DOC in ml, "12.5360".  
 C42: Time from start of titration to end in s, "62".  
 C43: Volume drift on start of a SET/KFT titration from the condi-  
 tioning in ul/min, "3.5".  
 C44: (Last measured) temperature in °C. Used for the temperature  
 compensation in pH measurements.  
 C45: Start volume with SET, STAT, DOC in ml, "2.800".  
 C46: Asymmetry pH of CAL, "6.89".  
 C47: Relative electrode slope of CAL, "0.9950".  
 C48: Voltage at maximum voltage (not valid in CAL and TIP)  
 C49: Voltage at minimum voltage (not valid in CAL and TIP)  
 DTime: Dosing time in s in DOS/DOC or time for the drift correc-  
 tion in SET/KFT with conditioning, resp.

**3.2.2.136. Info.TitrResults.FixEP.51.Value** read only  
 etc. up to **.59**  
**Info.TitrResults.pK.61.Value** read only  
 etc. up to **.69**  
**Info.TitrResults.TempVar.C70** read only/read+ write  
 etc. up to **.C79**  
**.FixEP:** Fix EP with DET, MET resp. Fix V with STAT. C5X corresponds to the  
 fix volumes X, X = 1...9.  
**.pK:** With DET, MET resp. time with given part of the end volume in s in STAT.  
 C6X corresponds to X = 1...9.  
**.TempVar:** Temporary variables in TIP corresponding to the assignments in  
 the submethods.

**3.2.2.137. Info.TitrResults.TimeWin.81.Mean** read only  
**Info.TitrResults.TimeWin.81.Dev** read only  
 etc. up to **.89**  
**Info.TitrResults.MeanRateC80.Mean** read only  
**Info.TitrResults.MeanRateC80.Dev** read only

For STAT: Mean rates and standard deviations in mL/min which have been  
 calculated with linear regression. C80 is the rate calculated over all points of  
 the measuring point list, C8X are the rates in the time windows X, X = 1...9.

<b>3.2.2.138.</b>	<b>Info.StatisticsVal.ActN</b>	read only
	<b>Info.Statistics.1.Mean</b>	read only
	<b>Info.Statistics.1.Std</b>	read only
	<b>Info.Statistics.1.RelStd</b>	read only
	etc. up to <b>.9</b>	

The current values of the statistics calculation.

\$Q sends, e.g.

ActN: Current value of the individual results "3"

Data for MN1:

Mean: Mean value (decimal places as in result) "3.421"

Std: Standard deviation (1 decimal place more than in result) "0.0231"

RelStd: Relative standard deviation (in %, 2 decimal places) "0.14"

<b>3.2.2.139.</b>	<b>Info.SiloCalc.C24.Name</b>	read only
	<b>Info.SiloCalc.C24.Value</b>	read only
	<b>Info.SiloCalc.C24.Unit</b>	read only
	for <b>.C25</b> as for <b>.C24</b>	
	<b>Info.SiloCalc.C26.ActN</b>	read only
	<b>Info.SiloCalc.C26.Mean</b>	read only
	<b>Info.SiloCalc.C26.Std</b>	read only
	<b>Info.SiloCalc.C26.RelStd</b>	read only
	for <b>.C27</b> as for <b>.C26</b>	

The current values from the silo calculations. C26 is the mean value out of the C24 variables; C27 comes from C25.

\$Q sends:

C24.Name: Name of the assigned value "RS1"

C24.Value: Value "2.222"

C24.Unit: Unit of the assigned value "%"

C26.ActN: Number of single results "3"

C26.Mean: Mean (decimal places as for the result itself) "3.421"

C26.Std: Standard deviation (decimal places as for the result + 1) "0.0231"

C26.RelStd: Relative standard deviation (in %, 2 decimal places) "0.14"

<b>3.2.2.140.</b>	<b>Info.ActualInfo.Inputs.Status</b>	read only
	<b>Info.ActualInfo.Inputs.Change</b>	read only
	<b>Info.ActualInfo.Inputs.Clear</b>	\$G
	<b>Info.ActualInfo.Outputs.Status</b>	read only
	<b>Info.ActualInfo.Outputs.Change</b>	read only
	<b>Info.ActualInfo.Outputs.Clear</b>	\$G

Status sends the current status of the I/O lines, Change sends the information regarding whether a change in status of a line has taken place since the last clearing, Clear clears the change information. For the output, there is a conversion from binary to decimal, e.g.

	0   0   0   0   0   0   0   0   0   0   1   0   1   0
Line No.	13  12  11  10  9   8   7   6   5   4   3   2   1   0
Output:	$2^1 + 2^3 = "10"$

1 means ON or change; 0 means OFF or no change.

The lines are assigned as follows (see also pages 222ff):

<b>Inputs:</b>	<b>Outputs:</b>
0 Start (pin 21)	0 Ready (pin 5)
1 Stop (pin 9)	1 Cond. ok (pin 18)
2 Enter (pin 22)	2 Titration (pin 4)

3	Clear (pin 10)	3	EOD (pin 17)
4	Smpl Ready (pin 23)	4	Monitoring, line L4 (pin 3)
5	pin 11	5	Error (pin 16)
6	pin 24	6	Activate, line L6 (pin 1)
7	pin 12	7	Pulse for recorder (pin 2)
		8	not used (pin 6)
		9	not used (pin 7)
		10	Monitoring (pin 8)
		11	Monitoring (pin 13)
		12	Monitoring (pin 19)
		13	Monitoring (pin 20)

**3.2.2.141. Info.ActualInfo.Assembly.CyclNo** read only  
 \$Q sends the current cycle number of the voltage measurement cycle, e.g. "127". From the cycle number and the cycle time (see 3.2.2.149), a time frame can be set up.  
 The cycle number is set to 0 on switching on the instrument, on every start and for QuickMeas. It is incremented as long as the instrument remains switched on.

**3.2.2.142. Info.ActualInfo.Assembly.Counter.V** read only  
**Info.ActualInfo.Assembly.Counter.Clear** \$G  
 \$Q sends the volume. With the function &Info.Assembly.Counter.Clear \$G, the volume counter is set to zero.

**3.2.2.143. Info.ActualInfo.Assembly.Meas** read only  
 \$Q sends the current measured value from the assembly.

**3.2.2.144. Info.ActualInfo.Titrator.CyclNo** read only  
**Info.ActualInfo.Titrator.V** read only  
**Info.ActualInfo.Titrator.Meas** read only  
**Info.ActualInfo.Titrator.dVdt** read only  
**Info.ActualInfo.Titrator.dMeasdt** read only  
**Info.ActualInfo.Titrator.dMeasdV** read only  
**Info.ActualInfo.Titrator.ERC** read only  
**Info.ActualInfo.Titrator.T** read only

\$Q sends the current values in the following formats:

	DET	MET	SET	KFT	STAT, DOC	DOS	MEAS	CAL
CyclNo	127	127	127	127	127	127	127	127
V(ml)	1.2345	1.2345	1.2345	1.2345	1.2345	1.2345	-	-
Meas:								
pH	3.345	3.345	3.6(mV)	-	3.6(mV)	3.6(mV)	3.345	3.345
U, Ipol (mV)	-345.6	-345.6	-345.6	-345.6	-345.6	-345.6	-345.6	-
Upol (uA)	-12.5	-12.5	-12.5	-12.5	-12.5	-12.5	-12.5	-
T (°C)	-	-	-	-	-	-	25.0	-
dVdt (ul/s)	-	-	2.5142	2.5142	2.5142	-	-	-
dMeasdt								

pH,U,lpol mV/s	0.7957	0.7957	0.7957	0.7957	0.7957	-	0.7957	0.7957
Upol (uA/s)	0.7957	0.7957	0.7957	0.7957	0.7957	-	0.7957	-
T (°C/s)	-	-	-	-	-	-	0.7957	-
dMeasdV (mV/ul)	-	-	10.6326	10.6326	10.6326	-	-	-
ERC	34	-	-	-	-	-	-	-
T (°C)	25.9	25.9	-	-	25.9	25.9	25.9	-

NV: Not Valid. If in the signal drift is OFF in modes MEAS and CAL, the signal drift is NV.

OV will be sent for "overrange".

A time frame can be set up from the cycle number and the cycle time (see 3.2.2.149). The cycle number is set to 0 at the start of a method and it is incremented until the end of the method.

<b>3.2.2.145.</b>	<b>Info.ActuallInfo.MeasPt.Index</b>	read only
	<b>Info.ActuallInfo.MeasPt.X</b>	read only
	<b>Info.ActuallInfo.MeasPt.Y</b>	read only
	<b>Info.ActuallInfo.MeasPt.Z1</b>	read only
	<b>Info.ActuallInfo.MeasPt.Z2</b>	read only
	<b>Info.ActuallInfo.EP.Index</b>	read only
	<b>Info.ActuallInfo.EP.X</b>	read only
	<b>Info.ActuallInfo.EP.Y</b>	read only

\$Q sends the last entry into the measuring point list (.MeasPt) or the last entry into the list of EP's with DET, MET, SET, KFT.

.MeasPt.X"165" Volume (DET, MET), time (STAT, DOS, DOC, SET, KFT, MEAS), resp. of the MPList in s

.MeasPt.Y"3.654" Measured value (DET, MET), volume (STAT, DOS, DOC, SET, KFT), resp. of the MPList in mL

.MeasPt.Z1"6.34" Measured value (STAT, DOS, DOC, SET, KFT, MEAS) of the MPList, format depends on the measured quantity

.MeasPt.Z2"25.8" Temperature of the MPList in °C

.EP.X"1.234" Volume coordinate of the EP

.EP.Y"5.34" Measured value coordinate of the EP

<b>3.2.2.146.</b>	<b>Info.ActuallInfo.Oven.HeatTime</b>	read only
	<b>Info.ActuallInfo.Oven.SampleTemp</b>	read only
	<b>Info.ActuallInfo.Oven.LowTemp</b>	read only
	<b>Info.ActuallInfo.Oven.HighTemp</b>	read only
	<b>Info.ActuallInfo.Oven.GasFlow</b>	read only
	<b>Info.ActuallInfo.Oven.UnitFlow</b>	read only

\$Q sends the current values from a connected KF Oven. If no Oven is connected, the values are empty.

.HeatTime: Heating time of sample in s.

.SampleTemp: Nominal sample temperature in °C.

.LowTemp: Lowest temperature during the sample heating time in °C.

.HighTemp: Highest temperature during the sample heating time in °C.

.GasFlow: Average gas flow during sample heating time.

.UnitFlow: Unit of gas flow.

<b>3.2.2.147.</b>	<b>Info.ActuallInfo.Display.L1</b>	up to 32 ASCII characters
	<b>Info.ActuallInfo.Display.L8</b>	up to 32 ASCII characters
	<b>Info.ActuallInfo.Display.DelAll</b>	\$G

Lines of the display. The display can be written to from the computer. Proceed as follows:

1. Lock the display, see 3.2.2.163.
  2. Delete the whole display (.DelAll).
  3. For writing onto the display, the standard character set will be used.
  4. Unlock the display, see 3.2.2.163
  5. Delete the whole display (.DelAll).
  6. Send a value to nod &Config.Aux.ResDisplay (see 3.2.2.105) to refresh the display.
- \$Q sends the contents of the corresponding display line.

**3.2.2.148. Info.ActualInfo.Comport.Number** read only  
\$Q sends the comport number of the Titrino where the PC is connected.

**3.2.2.149. Info.Assembly.CycleTime** read only  
**Info.Assembly.ExV** read only  
Inquiries regarding basic variables of the assembly: Cycle time in s, volume of the active Exchange Unit in mL.

**3.2.2.150. Info.Assembly.Bur.Select** internal D0, external D1, external D2  
Selection of buret for assembly functions.

**3.2.2.151. Info.Assembly.Bur.Empty** \$G, \$S, \$H, \$C  
Starts the function "empty". Only possible with the buret type "Dosino". The parameters are under the function &Config.DosPrep, see 3.2.2.118 up to 3.2.2.120.

**3.2.2.152. Info.Assembly.Bur.Prep** \$G, \$S, \$H, \$C  
Starts the function "preparation". The parameters are under the function &Config.DosPrep, see 3.2.2.110 up to 3.2.2.120.

**3.2.2.153. Assembly.Bur.Rates.Forward.Selected** digital, analog  
**Assembly.Bur.Rates.Forward.Digital** 0...150, max.  
**Assembly.Bur.Rates.Reverse.Selected** digital, analog  
**Assembly.Bur.Rates.Reverse.Digital** 0...150, max.

Expel and aspirating rate.

Digital or analog control. With digital control, the inputted value applies (in mL/min). "max." means maximum possible rate with the Exchange Unit in current use.

Analog means rate control with the analog potentiometer on Titrino.

**3.2.2.154. Assembly.Bur.Fill** \$G, \$H, \$C  
\$G starts the 'FILL' mode of the buret function.

**3.2.2.155. Assembly.Bur.ModeDis** \$G, \$S, \$H, \$C  
**Assembly.Bur.ModeDis.Selected** volume, time  
**Assembly.Bur.ModeDis.V** 0.0001...0.1...9999

<b>Assembly.Bur.ModeDis.Time</b>	0.25...1...86400
<b>Assembly.Bur.ModeDis.VStop</b>	0.0001...9999, <b>OFF</b>
<b>Assembly.Bur.ModeDis.AutoFill</b>	<b>ON, OFF</b>

Dispensing mode with parameters. The dispensing mode can only be started and stopped via the RS Control. During a running dosification, no method can be started at the Titrino.

.Selected: Dispensing of volume increments or during a preset time.

.Volume, .Time: Size of the volume increments or entry of time.

.VStop: Limit volume for the dispensing.

.AutoFill: ON means automatic filling after every dispensing.

<b>3.2.2.156. Assembly.Meas.Status</b>	<b>ON, OFF</b>
<b>Assembly.Meas.MeasInput</b>	1, 2, Diff., Ipol, Upol, Temp
<b>Assembly.Meas.Ipol</b>	±127...1...+127
<b>Assembly.Meas.Upol</b>	±1270...400...+1270

Measurement in assembly. The measuring function can only be started via RS Control. When the measuring function is switched on, no method can be started at the Titrino.

.Input: Selection of the potentiometric measuring input 1, 2, diff., polarized electrodes or temperature.

.Ipol: Polarization current in uA.

.Upol: Polarization potential in mV, entry in steps of 10 mV.

<b>3.2.2.157. Assembly.Outputs.AutoEOD</b>	<b>ON, OFF</b>
<b>Assembly.Outputs.SetLines</b>	<b>\$G</b>
<b>Assembly.Outputs.SetLines.L0</b>	active, inactive, pulse, <b>OFF</b>
up to .L13	
<b>Assembly.Outputs.ResetLines</b>	<b>\$G</b>

Setting the I/O output lines.

.AutoEOD: The automatic output of the EOD (End of Determination) at the end of the determination can be switched off. Thus, for example, in conjunction with a Titrino several determinations can be performed in the same beaker. Before AutoEOD is switched on, line 3 must be set to "OFF".

.SetLines: With \$G, all lines are set.

.SetLines.LX: Set the line LX. "active" means setting of a static signal, "inactive" means resetting of the signal, "pulse" means output of a pulse of app. 150 ms, "OFF" means the line is not operated, see also page 224.

Warnings:

- If you have "AutoEOD" to "ON", an active line 3 is set to "inactive" by the EOD pulse.
- L6 is the line of the activate pulse. An active line 6 is set to "inactive" by the activate pulse.
- L5 is the error line. It is continuously controlled by the Titrino program and can therefore not be set freely.

Line assignments in Titrino program:

L0	Ready, inactive state
L1	Conditioning OK
L2	Titration in progress

L3	EOD (End Of Determination)
L4	Can be set in TIP
L5	Error
L6	Activate pulse + can be set in TIP
L7	Buret volume pulses
L8, 9	---
L10-13	Monitoring in modes STAT, DOS, DOC

.ResetLines: Lines are set to the inactive status (= high).

**3.2.2.158. Assembly.Stirrer** ON, OFF  
Switching stirrer ON/OFF.

**3.2.2.159. Setup.Comport** 1, 2, 1&2  
Selects the Titrino COM for the output of automatic info:  
&Setup.Keycode  
&Setup.Trace  
&Setup.SendMeas  
&Setup.AutoInfo

**3.2.2.160. Setup.Keycode** ON, OFF  
ON means the key code of a key pressed on the Titrino is outputted. The key code comprises 2 ASCII characters; table of the keys with their code, see page 195. A keystroke of key 11 is sent as follows:  
#11  
The beginning of the message is marked by a space (ASCII 32).

**3.2.2.161. Setup.Tree.Short** ON, OFF  
**Setup.Tree.ChangedOnly** ON, OFF  
Definition of the type of answer to \$Q.  
.Short: With "ON", each path is sent with only the necessary amount of characters in order to be unequivocal (printed in bold in this manual). A combination of .Short and .ChangedOnly is not possible.  
.ChangedOnly: Sends only the changed values, i.e. values which have been edited. All paths are sent absolute, i.e. from the root.

**3.2.2.162. Setup.Trace** ON, OFF  
The Titrino automatically reports when a value has been confirmed with <enter> at the Titrino. Message, e.g.:  
&SmplData.OFFSilo.Id1"Trace"  
The beginning of the message is marked by a space (ASCII 32).

**3.2.2.163. Setup.Lock.Keyboard** ON, OFF  
**Setup.Lock.Config** ON, OFF  
**Setup.Lock.Parameter** ON, OFF  
**Setup.Lock.SmplData** ON, OFF

<b>Setup.Lock.UserMeth.Recall</b>	<b>ON, OFF</b>
<b>Setup.Lock.UserMeth.Store</b>	<b>ON, OFF</b>
<b>Setup.Lock.UserMeth.Delete</b>	<b>ON, OFF</b>
<b>Setup.Lock.Display</b>	<b>ON, OFF</b>

ON means disable the corresponding function:

- .Keyboard: Disable all keys of the Titrinos
- .Config: Disable the <configuration> key
- .Parameter: Disable the <parameter> key
- .SmplData: Disable the <smpl data> key
- .UserMeth.Recall: Disable "recall" in <user meth> key
- .UserMeth.Store: Disable "store" in <user meth> key
- .UserMeth.Delete: Disable "delete" in <user meth> key
- .Display: Disable the display, i.e. it will not be written to by the device program of the Titrino and can be operated from the computer.

<b>3.2.2.164. Setup.Mode.StartWait</b>	<b>ON, OFF</b>
<b>Setup.Mode.FinWait</b>	<b>ON, OFF</b>

Holding points in the method sequence. If they are "ON", the sequence stops until "OFF" is sent. Switching the instrument on sets both nodes to OFF:

- .StartWait: Holding point right after starting a method or submethod in TIP (holding point after AutoInfo !"T.GC").
- .FinWait: Holding point at the end a method or submethod in TIP (holding point after AutoInfo !"T.F").

<b>3.2.2.165. Setup.SendMeas.SendStatus</b>	<b>ON, OFF</b>
<b>Setup.SendMeas.Interval</b>	<b>0.08...4...16200, MPList</b>

- .SendStatus: ON means the automatic transmission of measured values (see 3.2.2.167 and 3.2.2.168) in the inputted interval is active.
- .Interval: Time interval (in s) for the automatic transmission of associated measured values defined under points 3.2.2.167 and 3.2.2.168. The inputted value is rounded off to a multiple of 0.08. The smallest possible time interval depends on the number of measured values which have to be sent, on the baud rate, on the load on the interface and on the type of device connection. With "MPList" the measured values are sent at the time of their entry into the measured point list.

The automatic transmission is switched on/off with 'SendStatus'.

<b>3.2.2.166. Setup.SendMeas.Select</b>	<b>Assembly, Titrator</b>
---	---------------------------

Selection of the unit of which the measured values should be sent (3.2.2.167 or 3.2.2.168).

<b>3.2.2.167. Setup.SendMeas.Assembly.CyclNo</b>	<b>ON, OFF</b>
<b>Setup.SendMeas.Assembly.V</b>	<b>ON, OFF</b>
<b>Setup.SendMeas.Assembly.Meas</b>	<b>ON, OFF</b>

Selection of the values from Assembly for the output in the set time interval (see 3.2.2.165):

- .CyclNo: Cycle number of the potential measurement. Together with the cycle time (3.2.2.149), a time frame can be set up.

The cycle number is set to 0 on switching on the instrument and it is always incremented as long as the instrument remains switched on.

.V: Volume  
 .Meas: Measured value associated to the cycle number.  
 The unit "assembly" must be preset (see 3.2.2.166).

**3.2.2.168. Setup.SendMeas.Titrator.CyclNo** ON, OFF  
**Setup.SendMeas.Titrator.V** ON, OFF  
**Setup.SendMeas.Titrator.Meas** ON, OFF  
**Setup.SendMeas.Titrator.dVdt** ON, OFF  
**Setup.SendMeas.Titrator.dMeasdt** ON, OFF  
**Setup.SendMeas.Titrator.dMeasdV** ON, OFF  
**Setup.SendMeas.Titrator.ERC** ON, OFF  
**Setup.SendMeas.Titrator.T** ON, OFF

Selection of the values from the titrator which are sent in the set time interval (see 3.2.2.165, formats see 3.2.2.144):

.CyclNo: Cycle number. Together with the cycle time (3.2.2.149), a time frame can be set up. The other data belong to the corresponding cycle number. The cycle number is set to 0 at the start of a method and it is incremented until the end of the method.

.V: Volume.  
 .Meas: Measuring value with STAT, DOS, DOC with activated measured value monitoring  
 .dVdt: associated volume drift.  
 .dMeasdt: associated measured value drift.  
 .dMeasdV: associated 1st derivative of the titration curve.  
 .ERC: ERC in DET.  
 .T: In STAT, DOS, DOC, with activated temperature monitoring.  
 The unit "titrator" must be preset (see 3.2.2.166).

**3.2.2.169. Setup.AutoInfo.Status** ON, OFF  
**Setup.AutoInfo.P** ON, OFF  
**Setup.AutoInfo.T.R** ON, OFF  
**Setup.AutoInfo.T.G** ON, OFF  
**Setup.AutoInfo.T.GC** ON, OFF  
**Setup.AutoInfo.T.S** ON, OFF  
**Setup.AutoInfo.T.B** ON, OFF  
**Setup.AutoInfo.T.F** ON, OFF  
**Setup.AutoInfo.T.E** ON, OFF  
**Setup.AutoInfo.T.H** ON, OFF  
**Setup.AutoInfo.T.C** ON, OFF  
**Setup.AutoInfo.T.O** ON, OFF  
**Setup.AutoInfo.T.N** ON, OFF  
**Setup.AutoInfo.T.Re** ON, OFF  
**Setup.AutoInfo.T.Si** ON, OFF  
**Setup.AutoInfo.T.M** ON, OFF  
**Setup.AutoInfo.T.EP** ON, OFF  
**Setup.AutoInfo.T.RC** ON, OFF  
**Setup.AutoInfo.C.B1** ON, OFF  
**Setup.AutoInfo.C.R1** ON, OFF

<b>Setup.AutoInfo.C.B2</b>	<b>ON, OFF</b>
<b>Setup.AutoInfo.C.R2</b>	<b>ON, OFF</b>
<b>Setup.AutoInfo.I</b>	<b>ON, OFF</b>
<b>Setup.AutoInfo.O</b>	<b>ON, OFF</b>

ON means that the Titrino reports automatically the moment the corresponding change occurs.

.Status: Global switch for all set AutoInfo.

.P PowerOn: Simulation of power on (3.2.2.171). Not from mains.

Messages from node .T, Titrator:

.T.R Ready: Status 'Ready' has been reached.

.T.G Go: Instrument has been started.

.T.GC GoCommand: Instrument (or submethod in TIP) has received a go command.

.T.S Stop: Status 'Stop' has been reached.

.T.B Begin of sequence (or submethod).

.T.F Final: End of determination (or submethod), the final steps will be carried out.

.T.E Error. Message together with error number, see page 111ff.

.T.H Hold: Status 'Hold' has been reached.

.T.C Continue: Continue after hold.

.T.O Conditioning OK: EP reached (in SET, KFT with conditioning).

.T.N Conditioning Not OK: EP not reached (in SET, KFT with conditioning).

.T.Re Request: In the inquiry of an identification or the sample size after start of titration.

.T.Si SiloEmpty: Silo empty, i.e. the last line has been removed from the silo memory.

.T.M MeasList: Entry in the measuring point list (with DET, MET, STAT, DOS, DOC, SET, KFT, MEAS).

.T.EP EPList: Entry into EP list (with DET, MET, SET, KFT)

.T.RC Results have been recalculated.

Messages from node .C, Comport:

.C.B1 COM1: A report is outputted on COM1. During this time, COM2 will be blocked. COM2 is generally blocked, if COM1 is busy.

.C.R1 COM1 is ready again. (Comes also when you <QUIT> an error.)

.C.B2, .R2 Identical for COM2.

Messages for changings in the I/O lines. If the changings are made simultaneously, there is 1 message. Pulses receive 2 messages: one message each for line active and inactive.

.I Input: Change of an input line.

.O Output: Change of an output line (except 7, pin 2, for recorder pulses).

If a change occurs that requires a message, the Titrino sends space (ASCII 32) and ! as an introducer. This is followed by the name of the device (see 3.2.2.106). Special ASCII characters in the device name are ignored. If no device name has been entered, only ! is sent. Finally the Titrino sends the information which node has triggered the message.

Example: !John".T.Si": The message was triggered from instrument "John", node .T.Si

<b>3.2.2.170. Setup.Graphics.Grid</b>	<b>ON, OFF</b>
<b>Setup.Graphics.Frame</b>	<b>ON, OFF</b>
<b>Setup.Graphics.Scale</b>	<b>Full, Auto</b>
<b>Setup.Graphics.Recorder.Right</b>	<b>0.2...0.5...1.00</b>

**Setup.Graphics.Recorder.Feed** 0.01...**0.05**...1.00

Change in the appearance and the format of the curve for the output. The settings are valid for both Titrimo COM ports.

.Grid: On/off switching of grid over curve.

.Frame: On/off switching of frame surrounding the curve. If grid and frame are switched off, the curve is printed faster as the print head does not have to move to the end of the paper.

.Scale: Type of scaling of the measured value axis: Full means that the scale runs from the smallest up to the greatest measured point. With auto, the smallest measured value is taken and the next smaller tick defines the beginning of the scale; the next greater tick to the greatest measured value is the end of the scale.

.Right: Relative specification of the width of the output medium (e.g. paper width) for the length of the measured value axis. 1 means the measured value axis is plotted over the entire width of the paper (largest possible width). In extreme cases, the writing of the right tick may lie outside.

.Feed: Length of the volume axis referred to the buret cylinder volume, V(B) per cm (0.1 means, e.g. 1 mL/cm with a 10 mL Exchange Unit). Depending on the printer, the measure in cm may not always be correct.

**3.2.2.171. Setup.PowerOn** \$G

Simulation of 'power on'. The device has the same status as after power on: The cylinder is filled, error messages deleted and the current sample number set to 0. The method last used is ready for operation.

**3.2.2.172. Setup.Initialise** \$G

**Setup.Initialise.Select** **ActMeth, Silo, Calib, Config, Assembly, Setup, All**

Setting of default values for the following areas:

ActMeth: Current method. Parameters, calculations, and assignments for the data output, operands C01...C19.

Silo: The silo memory is deleted. Same function as delete entire silo.

Calib: pH calibration data for all measuring inputs.

Config: All values under &Config.

Assembly: All values under &Assembly.

Setup: All values under &Setup.

All: Values of the entire tree (except silo and method memory).

The action must be triggered with &Setup.Initalise \$G.

**3.2.2.173. Setup.RamInit** \$G

Initializes instrument, see page 205. All parameters are set to their default value and error messages are cleared. The user and silo memories will be deleted. The user memory contains the default user methods from Metrohm.

**3.2.2.174. Setup.InstrNo** \$G

**Setup.InstrNo.Value** **serial number, 8 ASCII characters**

Instrument identification for report output.  
Set the value with &Setup.InstrNo \$G .

**3.2.2.175. Diagnose.Report** \$G  
Output of the report containing the adjustment parameters. The Titrino has to be in its inactive basic state.

### 3.3 Properties of the RS 232 Interface

Data Transfer Protocol

The Titrino is configured as DTE (Data Terminal Equipment).

The RS 232 interface has the following technical specifications:

- Data interface according to the RS 232C standard, adjustable transfer parameters, see page 10.
- Max. line length: 512 characters
- Control characters:  $C_R$  (ASCII DEC 13)  
 $L_F$  (ASCII DEC 10)  
XON (ASCII DEC 17)  
XOFF (ASCII DEC 19)
- Cable length: max. approx. 15 m

Start	7 or 8 Data Bit	Parity Bit	1 or 2 Stop Bit
-------	-----------------	------------	-----------------

Only a shielded data cable (for example, METROHM D.104.0201) may be used to couple the Titrino with foreign devices. The cable shield must be properly grounded on both instruments (pay attention to current loops; always ground in a star-head formation). Only plugs with sufficient shielding may be used (for example, METROHM K.210.0381 with K.210.9045).

#### 3.3.1 Handshake


##### Software-Handshake, SWchar

Handshake inputs on the Titrino (CTS) are not checked.


Handshake outputs (DTR, RTS) are set by the Titrino.

The Titrino sends XOFF when its input buffer contains 384 characters. After this it can receive 128 extra characters (including  $L_F$ ).

Titirino as Receiver :

Titrimo as Sender :


max. characters: 2 characters at 300...9600 baud  
16 characters at ≥ 19200 baud


**Software-Handshake, SWline**

Handshake input ports on the Titrimo (CTS) are not checked.  
Handshake output ports (DTR, RTS) are set by the Titrimo.  
The Titrimo has an input buffer which can accept up to 512 characters.

Titrimo as Receiver :


Titrimo as Sender:


Titrimo transmission can be stopped by external instruments with XOFF. After XOFF is received the Titrimo completes sending the line already started. If data output is disabled for more than 6 s by XOFF, E43 appears in the display.

**Hardware-Handshake, HWs**

Titrimo as Receiver :


Titrimo as Sender:

The data flow can be interrupted by deactivating the CTS line.

### 3.3.2 Pin Assignment

#### RS232C Interface

Protective earthing

Direct connection from cable plug to the protective ground of the instrument.

Polarity allocation of the signals


- Data lines (TxD, RxD)
  - voltage negative (< -3 V): signal state "ON"
  - voltage positive (> +3 V): signal state "ZERO"
- control or message lines (CTS, RTS, DTR)
  - voltage negative (< -3 V): OFF state
  - voltage positive (> +3 V): ON state

In the transitional range from +3 V to -3 V the signal state is undefined.

Driver 14C88      according to EIA RS 232C specification

Receiver 14C89      "      "

**Contact arrangement at plug (female) for RS 232C socket (male)**

View of soldered side of plug

Ordering numbers:  
K.210.0381 and K.210.9045

No liability whatsoever will be accepted for damage or injury caused by improper interconnection of instruments.

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
- 1: In CAL, EOD is sent after every buffer. Automatic output of EOD can be switched off via RS232, see page 174.
  - 2: The error line is reset when the error is rectified.
  - 3: Line also used in TIP as L6.
  - 4: Line is active if a remote box is registered, see page 9.
  - 5: Lines for monitoring L10, L11, L12, L13 in modes STAT, DOS, DOC. Line L12 is active if sample size is out; L13 can be set if the result is out of limits, see page 75
- The output L9 is free.

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<p><b>Certificate of Conformity and System Validation</b></p> <p>This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.</p>	
Name of commodity:	799 GPT Titrino
System software:	Stored in ROMs
Name of manufacturer:	Metrohm Ltd., Herisau, Switzerland
<p>This Metrohm instrument has been built and has undergone final type testing according to the standards:</p> <p><i>Electromagnetic compatibility: Emission</i> EN50081-1, EN50081-2, EN55022 (class B)</p> <p><i>Electromagnetic compatibility: Immunity</i> EN50082-1, IEC61000-6-2, Namur, IEC61000-4-2, IEC61000-4-3, ENV50204, IEC61000-4-4, IEC61000-4-5, IEC61000-4-6, IEC61000-4-11</p> <p><i>Safety specifications</i> IEC61010-1, EN61010-1, UL3101-1</p> <p>It has also been certified by the Swiss Electrotechnical Association (SEV), which is member of the International Certification Body (CB/IEC).</p> <p>The technical specifications are documented in the instruction manual.</p> <p>The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance. The features of the system software are documented in the instruction manual.</p> <p>Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.</p>	
<p>Herisau, October 30, 2001</p> <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p>Dr. J. Frank Development Manager</p> </div> <div style="text-align: center;">  <p>Ch. Buchmann Production and Quality Assurance Manager</p> </div> </div>	

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# Index

Keys are marked with < >, **display texts** are in bold characters, and pages concerning the green part are *printed in italic*.

## A

Accessories .....	230
<b>activate pulse:</b>	19, 29, 38, 49, 58, 66, 70, 72
<b>action:</b> .....	47, 57, 65
Analog potentiometer .....	2
<b>assign output L1X:</b> .....	48, 57, 65
Automatic reports .....	80
<b>auto start</b> .....	9
<b>auxiliaries</b> .....	9

## B

<b>backup</b> .....	92
<b>balance:</b> .....	8
Balance connection .....	210
Barcode reader .....	213
<b>barcode:</b> .....	9
<b>baud rate:</b> .....	10
<b>buffer 1 pH</b> .....	71
Buret	
connection .....	208
preparation .....	103
preselection .....	17, 55
<b>bytes missing</b> .....	187

## C

Cables .....	231
CAL .....	71
Calculations .....	74
<CAL.DATA> .....	94
Calibration .....	72
<b>cal.date</b> .....	94
<b>calibrate electrode</b> .....	187
<b>calibration:</b> .....	7
<b>cal.temp.</b> .....	94
<CARD> .....	91
<b>card...</b> .....	187
Card .....	91
Certificates .....	228f
CE .....	228
<C-FMLA> .....	76
<b>change card battery</b> .....	187
Characters input .....	6
<b>check ...</b> .....	187, 188
<CLEAR> .....	4
<b>common variables</b> .....	10, 79
Computer connection .....	212

## COM X

characteristics .....	180ff
configuration .....	10
pin assignment .....	183
<b>conditioning:</b> .....	28, 37
<CONFIG> .....	7
Configuration .....	7
Conformity .....	227
Connection	
balance .....	210
barcode reader .....	213
buret .....	208
computer .....	212
electrode .....	216
PC keyboard .....	214
printer .....	209
remote box .....	213
Sample Changer .....	211
stirrer .....	207
Ti Stand .....	207
Connection cables .....	231f
Contrast of display .....	2
Control parameters	
DOC .....	68
KFT .....	40
SET .....	31
STAT .....	51
<b>curve</b> .....	13
Curve	
changing output .....	13
display .....	83
examples .....	61, 81
printing .....	80
<CURVE> .....	83
<b>cycle lines:</b> .....	98
<b>cycles DX</b> .....	11

## D

<b>DX overload</b> .....	188
<b>data bit:</b> .....	10
<b>data lost</b> .....	188
Data	
input .....	5, 6
output .....	80
reproduction .....	81
transfer (RS232) .....	105ff
<b>data set reevaluation</b> .....	188
<b>date</b> .....	9
Decimal places .....	75

<DEF> .....	74ff	Entry of	
<b>delete all:</b> .....	98	data .....	5
<b>delete all silo lines</b> .....	97	text .....	6
<b>delete line n</b> .....	96	<b>EP at</b> .....	25, 34 63
<b>delete n</b> .....	77, 97	<b>EPC</b> .....	18
Delete		<b>EP overflow</b> .....	188
common variables .....	79	<b>EP recognition:</b> .....	18
formulas .....	74	EP titration .....	25
methods .....	89, 91	<b>equibr.time</b> .....	16, 69, 71
statistic values .....	77	Error messages .....	111, 187
text .....	6	<b>error XX</b> .....	192
user name .....	84	<b>evaluation</b> .....	18, 47
<b>delete silo lines</b> .....	97	Exchange Unit .....	230
<b>delta time</b> .....	52	<b>external DX</b> .....	103
DET .....	15ff	<b>extr.time</b> .....	27, 36
<b>dev.label</b> .....	10	<b>F</b>	
Diagnosis .....	193	Fabrication number .....	3
Diagnostic report .....	226	<b>factor</b> .....	16, 17, 46
<b>dialog:</b> .....	9	Filling rate	
<b>diam...</b> .....	11	after determination .....	17
Differential potentiometry .....	216	after manual dosing .....	2
<b>direction:</b> .....	63	Filling times .....	60
<b>directory...</b> .....	91	<b>fix EP</b> .....	19
<b>directory already exists</b> .....	188	<b>fix time</b> .....	47
<b>disp.crit:</b> .....	55	<b>fix-V</b> .....	47
<b>dispensing type:</b> .....	55	<b>formula</b> .....	74
Display, contrast .....	2	<b>G</b>	
<b>display drift:</b> .....	28, 37	General rules for RS232 .....	105ff
<b>display rate:</b> .....	49, 66	Graphics	
<b>disp.time</b> .....	55	changing output .....	13
<b>division by zero</b> .....	188	display .....	83
DOC .....	62ff	examples .....	61, 81
DOS .....	55ff	printing .....	80
<DOS> .....	2	<b>H</b>	
<b>dos.drive:</b> .....	11	<b>handshake:</b> .....	10
<b>dos.element:</b> .....	11, 17, 55	Handshake .....	180
Dosimat		Hardware handshake .....	82
connection .....	208	<b>I</b>	
preselection .....	17, 55	<b>id#1 or C21</b> .....	95
Dosing .....	2, 55, 103	Identification .....	95
Dosino		request .....	18
connection .....	208	<b>info</b> .....	83
preselection .....	17, 55	Initialize RAM .....	205
<b>dos.rate</b> .....	16	Input of	
Drift .....	20, 69	values .....	5
<b>drift corr:</b> .....	28, 37	text .....	6
<b>drift value</b> .....	28 37	Inquiries .....	5
<b>dynamics</b> .....	25, 34, 43, 62	<b>internal DO</b> .....	103
<b>E</b>		I/O lines .....	222f
<b>edit silo lines</b> .....	96	<b>I(pol)</b> .....	17
<b>electr.id</b> .....	71	ISO .....	229
Electrode connection .....	216		
<b>electrode test:</b> .....	17		
<b>end at</b> .....	62		
<ENTER> .....	4		

**K**

Keyboard, PC .....	214
<b>keyboard:</b> .....	9
Keypad .....	4
Keys	
<ABC> .....	6
<CAL.DATA> .....	94
<CARD> .....	91
<C-FMLA> .....	76
<CLEAR> .....	4
<CONFIG> .....	7
<CURVE> .....	83
<DEF> .....	74ff
<DOS> .....	2
<ENTER> .....	4
<MEAS/HOLD> .....	4
<MODE> .....	14
<PARAM> .....	15ff
<PREP> .....	103
<PRINT> .....	82
<QUIT> .....	4
<REPORTS> .....	82
<SILO> .....	96
<SMPL DATA> .....	95ff
<START> .....	4
<STATISTICS> .....	77
<STOP> .....	4
<USER> .....	84
<USER METH> .....	89
KF Oven .....	38, 232
KFT .....	34ff

**L**

<b>LX output:</b> .....	86
LED .....	2
<b>length...</b> .....	11
<b>limit</b> .....	188
<b>limit smpl size:</b> .....	19
Lock functions .....	12
<b>low lim.</b> .....	18, 19, 47, 75

**M**

Mains .....	3, 220
Malfunctions .....	111, 185, 186ff
<b>man.reports to COM:</b> .....	8
Manual operation .....	4ff
<b>manual stop</b> .....	188
<b>match id:</b> .....	100
Maximum/minimum evaluation .....	24
<b>max.rate</b> .....	25, 34, 43, 62
<b>mean n</b> .....	77
Mean values	
calculation .....	77
delete results .....	78
report .....	81
silo .....	100

MEAS .....	69
<MEAS/HOLD> .....	4
<b>meas.input:</b> .....	8, 17, 56, 64, 69, 71, 73
<b>meas</b> .....	188
<b>meas.mode:</b> .....	56, 73
<b>meas out of limit</b> .....	188
<b>meas.pt.density</b> .....	15
<b>meas.pt list overflow</b> .....	188
<b>measured quantity:</b> .....	14
<b>measuring parameters</b> .....	69
Measuring point list .....	61
printing .....	80
<b>meas.val:</b> .....	47, 57, 65
Memory card .....	91
MET .....	15ff
<b>method:</b> .....	85, 89, 91
Method name .....	89
Minimum/maximum evaluation .....	24
<b>min.incr.</b> .....	15
<b>min.rate</b> .....	25, 44, 62
<b>min.volume incr.</b> .....	34
<b>missing...</b> .....	188
<MODE> .....	14
Mode	
CAL .....	71
DET .....	25
DOC .....	62
DOS .....	55
KFT .....	34
MEAS .....	69
MET .....	15
SET .....	25
STAT .....	43
TIP .....	73
Mode selection .....	14
<b>monitoring</b> .....	47, 57, 65
Monitoring	
calibration interval .....	7
measured values .....	47, 57, 65
rates .....	47, 57
result limits .....	75
sample size limits .....	19
service interval .....	8
temperature .....	47, 57, 65
validation .....	9

**N**

<b>name:</b> .....	84
<b>next service</b> .....	8
<b>no...</b> .....	188, 189
Non-aqueous titrations .....	216
<b>not possible</b> .....	189
<b>not valid</b> .....	189

**O**

Object tree .....	116
-------------------	-----

- Ordering designations ..... 230
- outlet:** ..... 11
- Outputs..... 222f
- output LIX:** ..... 48, 57, 65, 75
- outside** ..... 189
- Oven ..... 38, 232
- Oven:** ..... 38
- overrange** ..... 189
- P**
- Parameters
- CAL ..... 71
- DET ..... 25
- DOC ..... 62
- DOS ..... 55
- KFT ..... 34
- MEAS ..... 69
- MET ..... 15
- Präp..... 10
- SET ..... 25
- STAT..... 43
- TIP ..... 73
- <PARAM> ..... 15ff
- parity:** ..... 10
- pause (X)** ..... 17, 26, 73
- PC keyboard ..... 214
- peripheral units** ..... 8
- pH(as)** ..... 94
- pH(init) ..... 75
- Pin assignment
- "Remote" socket..... 222f
- RS232 ..... 183
- pK/HNP:** ..... 19
- Power connection ..... 3, 220
- power on prep:** ..... 10
- prep...** ..... 190
- <PREP> ..... 4, 103
- prep.dosing elements** ..... 10
- preselections** . 19, 28, 37, 49, 58, 66, 70, 73
- <PRINT> ..... 82
- Printer
- connection ..... 209
- problems ..... 185
- selection..... 8
- Problems
- KF titrations ..... 42
- printer ..... 185
- SET titrations ..... 33
- STAT..... 51
- program** ..... 10
- Q**
- <QUIT> ..... 4
- R**
- RAM, initialize ..... 205
- rate:** ..... 48
- rate...** ..... 189, 190
- recall method**..... 89, 91
- reload** ..... 92
- Remote box..... 213
- remote box:** ..... 89
- Remote control
- via "Remote" lines ..... 222f
- via RS232..... 105ff
- report COMX:** ..... 80
- Report
- output ..... 80
- printing..... 81
- reproduction ..... 81
- selection ..... 80
- <REPORTS> ..... 81
- req.ident:** ..... 18
- req.smpl size:** ..... 19
- res.tab:** ..... 77
- Result
- calculation ..... 74
- delete ..... 77
- display ..... 10
- report ..... 80
- store..... 99
- text ..... 75
- result display:** ..... 10
- result out of limits** ..... 190
- RS1 decimal places**..... 75
- RS1 limit control:**..... 75
- RS1 text**..... 75
- RS1 unit:** ..... 75
- RS232 settings COMX**..... 10
- RS232 interface
- characteristics ..... 180ff
- configuration..... 10
- pin assignment..... 183
- run number** ..... 9
- S**
- same buffer** ..... 190
- sample size out**..... 190
- save lines:** ..... 98
- save lines off**..... 190
- Sample changer ..... 211
- sample changer:** ..... 71
- Sample identification ..... 95
- Sample size ..... 95
- Scope of delivery ..... 230
- second TIP call**..... 190
- Selection of
- balance..... 8
- printer..... 8
- send to COMX:**..... 8
- Sequence
- CAL..... 72

DET.....	20	<b>stop V</b> .....	17, 28, 37, 46, 56, 64
DOC.....	67	Store	
DOS.....	59	methods.....	89, 91
Sequence		results.....	99
KFT.....	39	Submethods in TIP.....	87
MET.....	20	Survey.....	2
SET.....	30	<b>sweep time</b> .....	62
STAT.....	50	<b>system error XX</b> .....	190, 191
TIP.....	87	System test report.....	226
Serial number.....	3	<b>system test report:</b> .....	8
<b>service:</b> .....	8		
<b>service is due</b> .....	190		
SET.....	25ff		
<b>SET X</b> .....	25		
Setting up.....	207ff		
<b>signal drift</b> .....	16, 69, 71		
<SILO>.....	96		
<b>silos calculations</b> .....	100		
<b>silos empty</b> .....	190		
<b>silos full</b> .....	190		
<b>silos lines</b> .....	97		
Silos reports.....	82, 100		
<b>slope</b> .....	94		
<SMPL DATA>.....	95ff		
<b>smpl unit:</b> .....	95		
Software handshake.....	180		
Special messages.....	187		
<START>.....	4		
<b>start at...</b> .....	62		
<b>start delay</b> .....	9		
<b>start pH</b> .....	44		
<b>start rate</b> .....	45		
<b>start time</b> .....	44		
<b>start V</b> .....	16, 26, 35, 44, 63		
STAT.....	43ff		
<STATISTICS>.....	4, 77		
<b>statistics</b> .....	77		
Statistics calculations.....	77		
Statistics values			
delete results.....	78		
report.....	81		
silos.....	100		
Stirrer connection.....	207		
<b>stirrer:</b> .....	86		
<b>stirrer control:</b> .....	9		
<STOP>.....	4		
<b>stop bit:</b> .....	10		
<b>stop conditions</b> .....	17, 28, 37, 46, 56, 64		
<b>stop crit:</b> .....	26, 35		
Stop criterion ...	17, 26, 35, 26, 35, 46, 56, 64		
<b>stop drift</b> .....	26, 35		
<b>stop EP</b> .....	18		
<b>stop pH</b> .....	18		
<b>stop... reached</b> .....	190		
<b>stop rate</b> .....	46		
<b>stop time</b> .....	46		
Stop times (STAT).....	52		
		<b>T</b>	
		Technical specifications.....	219
		<b>t(delay)</b> .....	26, 35, 46
		<b>t(delta)</b> .....	46
		<b>temp.</b> .....	191
		<b>temperature</b> .....	17, 27, 36, 45, 56, 64
		<b>temperature:</b> .....	48, 57, 65
		<b>temporary variables</b> .....	88
		<b>temp.out of limit</b> .....	191
		Text input.....	6
		<b>time</b> .....	9
		<b>time counter</b> .....	7, 8
		<b>time interval</b> .....	7, 8, 27, 36, 45, 55, 63, 71
		<b>time-out PC keyboard</b> .....	191
		TIP.....	73ff
		<b>TIP terminated</b> .....	191
		Titration	
		problems.....	33, 42, 51
		vessel.....	218
		Titration curve	
		changing output.....	13
		display.....	83
		examples.....	61, 81
		printing.....	80
		Titration modes.....	15
		Titration parameters	
		DET.....	25
		DOC.....	62
		KFT.....	34
		MET.....	25
		SET.....	25
		STAT.....	43
		Titration sequence	
		DET.....	20
		DOC.....	67
		KFT.....	39
		MET.....	20
		SET.....	30
		STAT.....	50
		TIP.....	87
		<b>titr.direction:</b> .....	26, 35, 45
		<b>transmission error</b> .....	191
		Tree.....	116

**U**

## Unit

result .....	75
sample .....	95
<b>up lim.</b> .....	18, 19, 47, 75
<b>U(pol)</b> .....	17
<USER> .....	84
<USER METH> .....	89
User	
methods .....	89, 91
name .....	84

**V**

<b>validate instrument</b> .....	191
<b>validation:</b> .....	9
Validation .....	226
Values input .....	5, 6
<b>violated limit:</b> .....	48, 57, 65
<b>V step</b> .....	15
<b>volume</b> .....	55
<b>volume DX</b> .....	11

**W**

<b>warn.interv.DX</b> .....	11
Warranty .....	227
<b>wrong card</b> .....	191
<b>wrong sample</b> .....	191