



Ion analysis

CH-9101 Herisau/Switzerland

Phone +41 71 353 85 85

Fax +41 71 353 89 01

E-Mail info@metrohm.ch

<http://www.metrohm.com>

813 Compact Autosampler

Program version 5.813.0010

Instructions for use

D. Möckli
Teachware
Metrohm Ltd.
Oberdorfstr. 68
CH-9100 Herisau

1st Edition 2001

Copyright © 2001 by Metrohm Ltd. All rights reserved.

This manual was carefully produced. Nevertheless the author does not warrant the information contained herein to be free of errors. The author appreciates your comments and feedback.


Table of content

1	Introduction	1
1.1	Instrument description	1
1.2	Parts and controls	2
1.3	Rear view	4
1.4	The keyboard	5
2	Installation	6
2.1	Setting up the instrument.....	6
2.1.1	Packaging	6
2.1.2	Delivery check	6
2.1.3	Location.....	6
2.2	Mains connection	7
2.2.1	Setting the mains voltage.....	7
2.2.2	Fuses.....	8
2.2.3	Mains cable and mains connection	8
2.2.4	Switching the instrument on/off.....	9
2.3	Attaching the accessories	9
2.3.1	Connecting the keyboard.....	9
2.3.2	Installing the plug cover	9
2.3.3	Installing the needle	10
2.3.4	Installing the splash protection	11
2.3.5	Placing the sample rack.....	12
2.3.6	Adjusting the sample rack	12
2.3.7	Tubing connections.....	14
2.3.8	Tubing connection to Metrohm 761 / 790 IC-Systems	15
2.3.9	Tubing connection to a Metrohm 747 VA Stand.....	16
2.3.10	Tubing connection to a Metrohm Computrace 757.....	16
2.4	Instrument connections.....	17
2.4.1	Remote interface	17
2.4.2	Interconnection with Metrohm IC systems.....	17
2.4.3	Interconnection with a Metrohm 746 VA Trace Analyzer	18
2.4.4	Interconnection with a Metrohm 757 VA Computrace	19
2.4.5	Interconnection with other instruments	19

3	Operation	20
3.1	Keyboard functions	20
3.2	Configuration	21
3.2.1	Course of configuration	23
3.3	RAM initialization	24
3.4	The methods	25
3.5	The Metrohm IC control software	27
3.5.1	Operational conditions	27
3.6	Modular Metrohm IC System 732/733	28
3.7	Metrohm 746 VA Trace Analyzer	28
3.7.1	Monitoring	28
3.7.2	Command sequence	28
3.7.3	Conditions to be met	29
4	Appendix.....	30
4.1	Error Messages	30
4.2	Sequences	31
4.3	Technical specifications	35
4.4	Maintenance and servicing	37
4.4.1	Maintenance by Metrohm service	37
4.4.2	Care of the unit	37
4.4.3	Replacing the pump tubing	37
4.5	Warranty and conformity	39
4.5.1	Warranty	39
4.5.2	Certificate of Conformity and System Validation	40
4.6	Standard equipment	42
4.7	Optional accessories	46
5	Index.....	48

Table of Figures

Fig. 1	Side view	2
Fig. 2	Rear view	4
Fig. 3	Keyboard	5
Fig. 4	Setting the mains voltage	8
Fig. 5	Keyboard connection	9
Fig. 7	Installing the splash protection	11
Fig. 8	Sample rack placing	12
Fig. 9	Installing the pump tubing	15
Fig. 10	Connection of a 747 VA Stand	16
Fig. 11	Measuring head of the 747 VA Stand	16
Fig. 12	Interconnection with 761 or 790 IC System	17
Fig. 13	Interconnection with a modular Metrohm IC System	18
Fig. 14	Interconnection with a Metrohm 746 VA Trace Analyzer	18
Fig. 15	Interconnection with 757 VA Computrace	19
Fig. 16	Configuration	23
Fig. 17	Time program of the IC control software	27



1 Introduction

1.1 Instrument description

The **813 Compact Autosampler** can be used for automating different tasks, e.g. ion chromatographic determinations with the Metrohm IC system instruments or voltammetric determinations with the Metrohm 746 Trace Analyzer or 757 VA Computrace. The sample rack provided may hold up to 36 vials of 2.5 or 11 mL volume respectively. Sample tubes made of polypropylene are standard. To protect the samples from external contamination, the tubes can be hermetically sealed.

Sample introduction from the 813 Compact Autosampler is achieved by means of its integrated peristaltic pump. The sample is conveyed by the pump through the capillary into e. g. the sample loop of the injector located within an ion chromatograph. A PEEK needle for sealed or a PEEK tube for open sample tubes can be used alternatively.

Preprogrammed methods for the most common modes of operation allow to use the 813 Compact Autosampler directly, without programming effort.

1.2 Parts and controls




Fig. 1 Side view

Safety note:



Never operate the instrument without splash protection and plug cover being mounted. The plug cover prevents any contamination of the connectors, caused by spilled solvents or chemicals.

1 Splash protection 6.2751.040 Must be installed always in operation!	14 Holding clamp For locking the tubing cartridge into place
2 Needle PEEK needle 6.1835.020 or PEEK tube 6.1835.010 or 6.1835.030	15 Snap-action lever For releasing the tubing cartridge
3 Lift	16 Screw thread for splash protection
4 Steel holder 4.766.4330 attached (part of 6.2833.000 Needle holder)	17 PEEK capillary 6.1831.060 or 6.1822.410 (1 m or 1.2 m) For conveying the sample to e.g. an injection valve
5 PEEK compression fitting 4.766.4320 for connection of PEEK capillary 7 (part of 6.2833.000 Needle holder)	18 Tower
6 PEEK compression fitting 6.2744.010	19 Remote connection
7 PEEK capillary 6.1831.050 (40 cm) connection needle – pump tubing	20 Plug cover 6.2752.010
8 Pump tubing 6.1826.040 For conveying the sample	21 Fuse holder Changing the fuses, see <i>section 2.2.2</i>
9 Screw thread for splash protection	22 Mains connection plug Mains connection, see <i>section 2.2.3</i>
10 Mounting pin For attaching the tubing cartridge	23 Mains switch For switching the instrument on/off: I = ON 0 = OFF
11 Tubing cartridge 6.2755.000 For 6.1826.0X0 pump tubing 6.1826.0X0	24 Sample rack 6.2041.7500
12 Pump drive Roller head with contact rollers	25 PP sample tube 6.2743.050 (can be sealed with 6.2743.060 PE caps)
13 Contact pressure lever For adjusting the contact pressure	26 Needle guide

1.3 Rear view




Fig. 2 Rear view

9	Mounting screws for splash protection	24	Sample rack 6.2041.750
11	Tubing cartridge 6.2755.000 For 6.1826.0X0 pump tubing	25	PP sample tubes 6.2743.050 (can be sealed with 6.2743.060 PE caps)
16	Mounting screws for splash protection	27	Serial number
18	Tower	28	Keyboard connection
20	Plug cover 6.2752.010		

1.4 The keyboard




Fig. 3 Keyboard

- The <START> key starts a sample series or continues an interrupted series.


- The <STOP> key terminates a sample series or resets (if pressed twice) the 813 Compact Autosampler (RESET).

- The <HOLD> key interrupts the course of a sample series and establishes the 'HOLD' state. The green LED blinks. With the <START> key the interrupted sample series can be continued.

- The green LED indicates the status of the Autosampler.
 - The LED is lit in the 'ready' state.
 - It blinks steady during the course of a method.
 - It blinks in fast frequency when a series is interrupted ('HOLD' state).

- The red LED indicates error states or parameter settings.
 - If an error occurred the type of error is indicated by the number of blinks, see *Appendix, Error messages*.

2 Installation



2.1 Setting up the instrument

2.1.1 Packaging

The 813 Compact Autosampler is supplied together with the separately packed accessories in special packagings containing shock-absorbing foam linings designed to provide excellent protection. The actual instrument is packed in an evacuated polyethylene bag to prevent the ingress of dust. Please store all these special packagings as only they can assure damage-free transport of the instrument.

2.1.2 Delivery check

After receipt, immediately check whether the shipment is complete and undamaged (compare with delivery note and list of accessories in *section 4.6*). In the case of transport damage, see instructions in *section 4.5.1 "Warranty"*.

2.1.3 Location

Position the instrument in the laboratory at a location convenient for operation, free from vibrations and protected against a corrosive atmosphere and contamination by chemicals.



*Do not operate the 813 Compact Autosampler without splash protection **1** and plug cover **20** being installed!*



Take precautions to ensure that any leaks from pump tubings or connections cannot cause more damage.



Take precautions to ensure that any leaks from pump tubings or connections cannot cause more damage.

2.2 Mains connection



Follow the instructions below for connecting to the power supply. If the instrument is operated with the mains voltage set wrongly and/or wrong mains fuse there is a danger of fire!

2.2.1 Setting the mains voltage

Before switching on the 813 Compact Autosampler for the first time, check that the mains voltage set on the instrument (see Fig. 4) matches the local mains voltage. If not, reset the mains voltage on the instrument as follows:

1 Disconnect mains cable

Disconnect mains cable from mains connection plug **22** of the 813 Compact Autosampler.

2 Remove fuse holder

Using a screwdriver, loosen fuse holder **21** beside the mains connection **22** and take out completely.

3 Check fuse

Carefully take the fuse installed for the desired mains voltage out of fuse holder **21** and check its specifications (the position of the fuse in the fuse holder is marked by the white arrow imprinted next to the mains voltage range):

100...120 V 0.5 A (slow-blow)

Metrohm-Nr. U.600.0013

220...240 V 0.25 A (slow-blow)

Metrohm-Nr. U.600.0010

4 Insert fuse

Change fuse if necessary and reinsert in fuse holder **21**.

5 Install fuse holder

Depending on the desired mains voltage, insert fuse holder **21** in the 813 Compact Autosampler so that the corresponding mains voltage range can be read normally and the adjacent white arrow points to the white bar imprinted on the right side of the fuse holder (see Fig. 4).




Fig. 4 Setting the mains voltage

2.2.2 Fuses

One of the two fuses 0.5 A/slow-blow for 100...120 V or 0.25 A/slow-blow for 220...240 V is installed in fuse holder **21** of the 813 Compact Autosampler as standard



Ensure that the instrument is never put into operation with fuses of another type, otherwise there is danger of fire!

For checking or changing fuses, proceed as described in section 2.2.2.

2.2.3 Mains cable and mains connection

Mains cable

The instrument is supplied with one of three mains cables

- 6.2122.020 with plug SEV 12 (Switzerland, ...)
- 6.2122.040 with plug CEE(7), VII (Germany, ...)
- 6.2133.070 with plug NEMA 5-15 (USA, ...)

which are three-cored and fitted with a plug with an earthing pin. If a different plug has to be fitted, the yellow/green lead (IEC standard) must be connected to protective earth (protection class 1).



Any break in the earthing inside or outside the instrument can make it a hazard!

Mains connection

Plug the mains cable into mains connection plug **22** of the 813 Compact Autosampler (see Fig. 4).

2.2.4 Switching the instrument on/off

The 813 Compact Autosampler is switched on and off using the mains switch **23**.

2.3 Attaching the accessories



For attaching the accessories at the 813 Compact Autosampler, proceed in the order described below.

2.3.1 Connecting the keyboard

1 Connecting the keyboard

Connect the 6.2142.110 keyboard to the keyboard connection **28** "Keyboard". For disconnection, press the plug together slightly on both sides.

2 Switch on the instrument

Switch on the 813 Compact Autosampler with mains switch **23**. The instrument is initialized and the lift is raised completely. The green LED lights up.




Fig. 5 Keyboard connection

2.3.2 Installing the plug cover



To prevent any contamination of the mains and remote connection by spilled solvents or chemicals, the 6.2752.010 plug cover must always be installed when operating the 813 Compact Autosampler!

Install the plug cover **20** in the corresponding guide groove above mains connection plug **22** and remote connection **19** (see Fig. 1 and Fig. 2).

2.3.3 Installing the needle

1 Remove PEEK compression fitting 5

Remove PEEK compression fitting **5** screwed onto the needle holder **4**.

2 Insert needle

- Needle **2** (6.1835.010 PEEK needle or 6.1835.020 or 6.1835.010 PEEK tube) is inserted a short distance into the opening of steel holder **4**.
- Push PEEK round seal **29** over needle **2** from above. The smaller end of the seal must face upwards.

3 Fix needle

- Screw PEEK compression fitting **5** into steel holder **4** while pressing needle **2** gently upwards from below.
- Tighten compression fitting **5** in steel holder **4** by hand (never use tools!).

2 PEEK needle (6.1835.010) or PEEK tube (6.1835.020 or 6.1835.030)

4 Steel holder (4.766.4330)
(part of 6.2833.000 Needle holder)

5 PEEK compression fitting (4.766.4320)
(part of 6.2833.000 Needle holder)

6 PEEK compression fitting (6.2744.010)

7 PEEK capillary (6.1831.050)

29 PEEK ferrule




Fig. 6 Needle installation



If you are using the 6.1835.020 or 6.1835.030 PEEK tube as needle, the sample tubes must **not be sealed with caps** because they cannot be pierced by the PEEK tube and this could damage the tube!

2.3.4 Installing the splash protection



To avoid any danger of injury by the needle, the 6.2751.070 splash protection must always be installed when operating the 813 Compact Autosampler!

1 Remove holding screws

Remove the holding screws **16** and the washer mounted on the screw threads **9** at tower **18** using the 6.2621.100 Allen key.

2 Remove protective film from splash protection

Remove the plastic film glued on both sides of the splash protection **1**.

3 Install splash protection

Attach splash protection **1** onto the screw threads **9** and **16** at the tower **18** and fix it with the holding screws **9** and the washer using the 6.2621.100 Allen key.




Fig. 7 Installing the splash protection

2.3.5 Placing the sample rack

1 Place sample rack

Place sample rack **24** on the turntable of the 813 Compact Autosampler acc. to *Fig. 8*.

2 Initialize rack

Press <RESET> to move the rack to the home position, in which the magnetic code can be read (coding magnet **30**).




Fig. 8 Sample rack placing

2.3.6 Adjusting the sample rack

If a new sample rack is placed on the 813 Compact Autosampler for the first time, it must be adjusted. Make sure you have mounted a steel needle or a PEEK needle before. Proceed as follows:

1 Place sample rack

Place empty sample rack **24** on the turntable of the 813 Compact Autosampler and switch off the instrument.

2 Switch on the instrument

Keep the <HOLD> key pressed and simultaneously press the mains switch on the right hand side of the 813 Compact Autosampler. The instrument is initialized and the needle is placed above rack position 1. Press the <START> key. The green LED is now blinking twice.

3 Check needle position

- Press <HOLD> until the needle **2** is about 1 cm above the upper level of the sample rack **24**.
- Check needle position: If the needle **2** cannot be lowered unhindered through the upper hole of the rack position 1, continue directly with point **4**.
- Press <HOLD> once and continue lowering the needle by pressing <HOLD> again until the needle **2** stops.
- Check needle position: If the needle **2** points to the middle of lower hole of the rack position 1, the sample rack must not be adjusted (continue directly with point **5**).

4 Adjusting the sample rack

- Loosen the four adjusting screws **33** on the lower level of the sample rack using the 6.2621.100 Allen key .
- Carefully turn the two upper levels of the sample rack **27** by hand until the lowered needle **2** is exactly in the middle of the lower hole at the rack position 1.
- Tighten the adjusting screws **33**.

5 Move to rest position

Press <STOP> twice to move the sample rack to the initial position.

2.3.7 Tubing connections



*Pump tubings are consumable material with a lifetime which depends on the contact pressure (see section 4.4 Maintenance and servicing). This is why the tubing cartridges should be raised completely by loosening snap-action lever **15** on the right-hand side if the pump is to remain switched off for a considerable length of time (the set contact pressure remains unchanged).*



The 6.1826.0X0 pump tubing is made of PVC and must not be used for rinsing with solutions which contain acetone. In such cases, rinse with different pump tubing or a different pump.

For transferring the sample from the 813 Compact Autosampler to an injection valve the following tubing connections must be made:

1 Install pump tubing

- Release tubing cartridge **11** from holding clamp **14** by pressing down snap-action lever **15** and remove it from mounting pin **10** on the 813 Compact Autosampler (see Fig. 9).
- Press contact pressure lever **13** on the tubing cartridge down as far as it will go.
- Insert pump tubing **8** in the tubing cartridge **11** as shown in Fig. 9. The black stopper **32** must click into the corresponding holder on the left-hand side of the tubing cartridge.
- Place the tubing cartridge **11** on mounting pin **10** and press down on the right-hand side until snap-action lever **15** clicks into position on holding clamp **14**. Take care that no kinks are formed in the pump tubing.

2 Connection needle – pump tubing

- Mount a PEEK compression fitting **6** on both ends of the PEEK capillary **7**.
- Screw the PEEK capillary **7** with the PEEK compression fitting **6** on to the PEEK compression fitting **5** already mounted on needle holder **4** (see Fig. 9).
- Screw a coupling **31** on to the PEEK compression fitting at the other end of PEEK capillary **7**.
- Push coupling **31** on to the inlet end of the pump tubing **8** (see Fig. 9).




Fig. 9 Installing the pump tubing

6	PEEK compression fitting 6.2744.010	15	Snap-action lever
7	PEEK capillary 6.1831.050	17	PEEK capillary 6.1831.060 or 6.1822.410
8	Pump tubing 6.1826.040	31	Coupling 6.2744.030
11	Tubing cartridge 6.2755.000	32	Stopper (black)
13	Contact pressure lever		

2.3.8 Tubing connection to Metrohm 761 / 790 IC-Systems

3 Connection pump tubing – injection valve

- Use a 6.2744.010 PEEK compression fitting and a 6.2744.030 coupling to tighten the 6.1831.060 PEEK capillary **17** to the pump tubing **8**.
- Loosen the screw thread of the outer connector of the suction tubing inside the Ion chromatograph (see the Ion chromatograph's user manual).
- Replace the suction tubing of the Ion chromatograph with PEEK capillary **17**. Pass it through the outer connector of the front side of the housing.
- Use a 6.2744.010 PEEK compression fitting to tighten the capillary to the connector "1" of the injection valve.

4 Tubing connection injection valve – waste

- Remove the coupling and the 6.2744.010 compression fitting from syringe tubing.
- Lead the free end of the syringe tubing to the waste container and fix in place.

2.3.9 Tubing connection to a Metrohm 747 VA Stand

3 Connection Pumpschlauch – VA Stand 747

- Use a 6.2744.010 PEEK compression fitting and a 6.2744.030 coupling to tighten the 6.1831.060 PEEK capillary **17** to the pump tubing **8**.

To keep the dead volume as small as possible shorten the pump tubing to a suitable length.

4 Measuring head of the 747 VA Standes

- Insert the other end of the PEEK capillary **17** into a threaded opening of the Measuring head of the VA Stand and fix it with a 4.420.2580 fitting screw.




Fig. 10 Connection of a 747 VA Stand




Fig. 11 Measuring head of the 747 VA Stand

2.3.10 Tubing connection to a Metrohm Computrace 757

Install tubing connections 813 Compact Autosampler – 757 VA Computrace as shown above for the 747 VA Stand. See the 757 VA Computrace user manual for more details.

2.4 Instrument connections

2.4.1 Remote interface

External devices can be connected to the 25 pin remote connection (see Fig. 1). The remote line Input 2 (pin 22, see section 4.3 Technical description) allows the 813 Compact Autosampler to get synchronized with an other instrument.

Connecting instruments to the 813 Compact Autosampler requires Metrohm cables. Otherwise a safe signal transmission may not be guaranteed.



*Before an external device is connected to remote connection **22** the 813 Compact Autosampler must always be switched off using mains switch **26**!*

2.4.2 Interconnection with Metrohm IC systems

The 813 Compact Autosampler is connected to a Metrohm IC system (761 Compact IC or 790 Personal IC) using the 25 pin 6.2141.130 Remote cable.




Fig. 12 Interconnection with 761 or 790 IC System

The interconnection of a 813 Compact Autosampler with a modular Metrohm IC System (IC Detector 732, IC Separation Center 733 ...) can be done in different ways. Depending on the individual instrument combinations you need a 6.2141.130 or 6.2125.090 Remote cable to connect the 813 Compact Autosampler to a 732 IC Detector. A 6.2125.120 Remote adapter may be required as well. See the user manual of your IC System or contact your Metrohm distributor for more details.

Example for the interconnection of a modular Metrohm IC System:




Fig. 13 Interconnection with a modular Metrohm IC System

2.4.3 Interconnection with a Metrohm 746 VA Trace Analyzer

The 813 Compact Autosampler is connected to a Metrohm 746 VA Trace Analyzer with 747 VA Stand using the 25 pin 6.2141.020 Remote cable. Connect the Remote cable to the 'Control lines' socket of the Trace analyzer.




Fig. 14 Interconnection with a Metrohm 746 VA Trace Analyzer

2.4.4 Interconnection with a Metrohm 757 VA Computrace

You need a 6.2141.150 Remote cable to interconnect a 813 Compact Autosampler with a Metrohm VA Computrace 757, see below.





Fig. 15 Interconnection with 757 VA Computrace

2.4.5 Interconnection with other instruments

The Metrohm 813 Compact Autosampler can be connected to other Metrohm instruments or instruments of other equipment manufacturers. It is strongly recommended to use a parallel Metrohm 25-pin Remote cable (6.2141.130 or 6.2141.020) for the interconnection to function properly. The parallel interface of the other instruments has to fulfill the electrical requirements of the Metrohm Remote interface, see 4.3 *Technical specifications*.

The method sequences of the Metrohm 813 Compact Autosampler require trigger signals (pulses) via the Remote connection. A trigger signal will be accepted, whenever a negative slope (TTL voltage) or a negative pulse (min. duration 200 ms) on Remote input line 2 (**pin 22**) is recognized, see section 4.3 *Technical specifications*.

3 Operation



3.1 Keyboard functions


The 813 Compact Autosampler provides three keys and two LEDs, which may offer different functions depending on the current instrument status.

After switching on the instrument the **readiness** of the 813 Compact Autosampler is displayed by lighting up of the green LED.


Function of the keys



Starts processing of a sample series.



Stops processing of a sample series or initializes the 813 Compact Autosampler (by pressing twice).



Interrupts processing of a sample series, green LED blinks. An interrupted process can be continued with <START>.



*If using different proceedings (single injections, double or triple injections), you have to select the appropriate **method before starting a sample series.***



*Depending on the size of the sample tubes and the desired lift height for the suction of the sample solution the **work position** for the needle has to be set **before starting a sample series.***

Read the following pages to learn how to configure the 813 Compact Autosampler.

3.2 Configuration


To select a method or to define the work position switch on the 813 Compact Autosampler while holding down the <HOLD> key. The green LED blinks once a second while you can enter the desired settings.

Keyboard functions for configuration:


(after **Power on** while <HOLD> depressed, green **LED blinks**)



Accepts the setting of a value and moves on to the next value.



Selection of method or manual lift movement for the definition of the work position.



Cancels a setting or closes the configuration mode after a settings is accepted with <START>.

Four different **Autosampler methods** are available. They apply to the following tasks:

- *Single injection with air gap in after each determination, the needle will be drawn out of the sample.*
- *Single injection without air gap, the needle rests in the sample during the determination.*
- *Double injection of each sample, no air gap, the needle remains in the sample during the determination.*
- *Triple injection of each sample, no air gap, the needle remains in the sample during the determination.*















Before running the 813 Compact Autosampler for the first time, the **work position** of the lift has to be set correctly. The work position defines the position of the suction needle while aspirating the sample. It depends mainly from the size of the sample tubes.

For setting the work position place a filled sample tube on the sample rack of the Autosampler on rack position 1.

The needle should remain dipped in the sample solution during the aspiration interval or reach the bottom of the sample vial, depending on the application.

Press <START> to accept a setting. To cancel a setting press <STOP>. Press <STOP> twice to abort the configuration procedure.

The configuration procedure

Key	LEDs	Explanation
 + "Power ON"	 green 1*	<p>In order to start the configuration procedure the <HOLD> key must be pressed while the Autosampler is switched on.</p> <p>The green LED lights up once to show that the configuration procedure is active. The method can now be selected.</p>
	 red 1-4*  green 1*	<p>When the <HOLD> key is pressed the red LED blinks once for the first method. Each time <HOLD> is pressed again the next method is selected up to method 4, then method 1 again..</p>
	 green 2*	<p><START> accepts the setting. The green LED blinks twice. The working position of the lift can now be set.</p>
	 red blinks  green 2*	<p>When the <HOLD> key is pressed the lift can be moved manually. Each further time the <HOLD> key is pressed the direction in which the lift moves is altered. In this way the working position can be set accurately.</p>
	 blinks	<p><START> accepts the setting. The working position is stored. The configuration procedure is now complete. It can be restarted by pressing <HOLD> in order to select the method.</p>
	 lights up	<p>The configuration procedure is terminated with <STOP>. When the green LED is illuminated the 813 Compact Autosampler is in the basic state and is ready to process a series of samples.</p>

See next page for a figure of the whole configuration process.

3.2.1 Course of configuration




Fig. 16 Configuration

If you just want to select a method:

1. Switch off Autosampler.
2. Press <HOLD> and simultaneously switch on Autosampler. The green LED blinks.
3. Press <HOLD> again, to select method 1. The red LED blinks slowly. To select method 2 press <HOLD> again. The red LED blinks double times. To select method 3 or 4 press <HOLD> again. The red LED will display the selected method by the blinking frequency.

4. Accept the selected method with <START>.
5. To return to the basic state, press the <STOP> key.

If you just want to define the work position:


1. Switch off Autosampler .
2. Press <HOLD> and simultaneously switch on Autosampler. The green LED blinks.
3. Press <START>. The green LED blinks twice.
4. Press <HOLD> to run the lift downward. The red LED slowly. Press <HOLD> again and again to set the work position exactly. Each time the <HOLD> key is pressed, the direction of the lift movement is reversed. The PEEK needle should be immersed deep enough to assure a trouble-free aspiration of the sample or. Depending on your application you may want to aspirate the whole sample. In the latter case drive the lift down until the needle reaches the bottom of the sample vial.
5. The lift position set can be accepted with the <START> key.
6. Reestablish basic state of the 813 Compact Autosampler by pressing the <STOP> key.

3.3 RAM initialization

The internal software of the 813 Compact Autosampler may be initialized to reset all instrument settings after serious malfunctions, which should be very scarce of course. All settings (including the not accessible ones) will be reset to factory defaults.

Initialization can be done by pressing all keys simultaneously while switching on the instrument.

It may be appropriate to initialize the 813 Compact Autosampler before using it the first time.



Method selection and the definition of the working position has to be redone after RAM initialization.

3.4 The methods

All methods of the 813 Compact Autosampler make use of the built-in peristaltic pump to aspirate the sample and to fill the sample loop of the Ion chromatograph connected. The pump duration is controlled by the time program of the Ion chromatograph's PC software by sending Remote signals to the Autosampler.

The number of samples to be processed is to be set in the sample queue of the control software. The Autosampler should be stopped manually after processing the last sample.

Method 1 (detailed command sequence, see *section 4.2 Sequences*)

Single determination with air gap after each run, needle lifted after sample transfer.

Sequence:

1. Await first Remote signal of the instrument.
2. Get sample tube and immerse needle.
3. On receiving second Remote signal, switch on the peristaltic pump. Sample is transferred.
4. Await third Remote signal. The instrument starts its determination.
5. Peristaltic pump is switched off after 12 s.
6. Raise needle out of sample and aspirate air for 6 s.

Method 2 (detailed command sequence, see *section 4.2 Sequences*)

Single determination, no air gap, needle remains in sample after sample transfer.

Sequence:

1. Await first Remote signal of the instrument.
2. Get sample tube and immerse needle.
3. On receiving second Remote signal, switch on the peristaltic pump. Sample is transferred.
4. Await third Remote signal. The instrument starts its determination.
5. Peristaltic pump is switched off after 12 s.

Method 3 (detailed command sequence, see 4.2 Sequences)

Double determination of each sample, no air gap, needle remains in sample after sample transfer.

Sequence:

1. Await first Remote signal of the instrument.
2. Get sample tube and immerse needle.
3. On receiving second Remote signal, switch on the peristaltic pump. Sample is transferred.
4. Await third Remote signal. The instrument starts its determination.
5. Peristaltic pump is switched off after 12 s.
6. Await first Remote signal of the instrument again.
7. On receiving second Remote signal, switch on the peristaltic pump. Sample is transferred.
8. Await third Remote signal. The starts the second determination.
9. Peristaltic pump is switched off after 12 s.

Method 4 (detailed command sequence, see 4.2 Sequences)

Triple determination of each sample, no air gap, needle remains in sample after sample transfer.

Sequence:

1. Await first Remote signal of the instrument.
2. Get sample tube and immerse needle.
3. On receiving second Remote signal, switch on the peristaltic pump. Sample loop is transferred.
4. Await third Remote signal. The instrument starts its determination.
5. Peristaltic pump is switched off after 12 s.
6. Await first Remote signal of the instrument again.
7. On receiving second Remote signal, switch on the peristaltic pump. Sample is transferred.
8. Await third Remote signal. The instrument starts the second determination.
9. Peristaltic pump is switched off after 12 s.
10. Await first Remote signal of the IC again.
11. On receiving second Remote signal, switch on the peristaltic pump. Sample is transferred.
12. Await third Remote signal. The instrument starts the third determination.
13. Peristaltic pump is switched off after 12 s.

3.5 The Metrohm IC control software

To use the 813 Compact Autosampler with a 761 Compact IC or a 790 Personal IC one has to create a time program. First of all the IC has to send a Remote signal (Remote line 3) to synchronize the command sequence (first Remote signal) of the 813 Compact Autosampler. After that the injection valve is switched to 'Fill' position. Then a second Remote signal starts the filling of the sample loop. The pump duration for filling the sample loop is defined by the third Remote signal. After sending this Remote signal the injection valve is switched to 'Inject' position immediately which starts the data acquisition. During sample injection the peristaltic pump of the 813 Compact Autosampler is still running. 12 s after receiving the Remote signal the Autosampler will stop the pump.




Fig. 17 Time program of the IC control software

3.5.1 Operational conditions

To ensure a correct interconnection (Autosampler – Ion chromatograph), the following conditions must be met:

- When the Ion chromatograph is operational the 813 Compact Autosampler has to be started before the 'Sample Queue' of the IC.
- When using the 761 Compact IC the initial state of the Remote line 3 has always to be set to 0 (set **System startup values: Remote line 3 = 0**).

3.6 Modular Metrohm IC System 732/733

The 813 Compact Autosampler can be controlled directly by a Metrohm 732 IC Detector or by a PC control software, such as the Metrohm IC Net 2.0 using the remote control lines of the 732 IC Detector. See the 732 IC Detector user manual for details or read the IC Net 2.0 documentation.

3.7 Metrohm 746 VA Trace Analyzer

Using a 813 Compact Autosamplers with a 746 VA Trace Analyzer require modified command sequences and monitoring settings.

Between each sample determination a rinsing procedure for the tubings has to be executed. For this purpose place a sample tube containing rinsing solution on each second (even) rack position. The samples to be determined have to be placed between those rinsing tubes. This reduces the amount of a sample series to 18 samples max.

3.7.1 Monitoring

The number of samples of a series is defined on the monitoring page of the 746 VA Trace Analyzer, see below.

```
Auto.samples on 18 of 160
Auto.batch on 0 of 18
Start interval 0.0 h
```

3.7.2 Command sequence

A command sequence of the 746 VA Trace Analyzer has to contain explicit trigger signals for the 813 Compact Autosampler that are to be sent via the Remote output line 2. Use the 'SETCTRL' command to send the required pulses, see line 4 and 5 of the example method on the next page.

Sample method for the 746 VA Trace Analyzer:

	Instructions	t/s	Main parameters	Auxiliary parameters
1	STIR		Rot.speed 2000 /min	
2	PURGE			
3	SMPL/M		V.fraction mL	V.total L
4	SETCTRL	1.0	Code *****1**	Message
5	SETCTRL		Code *****0**	Message
6	REM	30.0	move to sample	
7	SETCTRL	1.0	Code *****1**	Message
8	SETCTRL		Code *****0**	Message
9	REM	300.0	transfer sample	
10	SETCTRL	1.0	Code *****1**	Message
11	SETCTRL		Code *****0**	Message
12	REM		stop pump	
13	DOS>M		Soln.name Buffer	V.add 2.000 mL
14	PURGE	90.0		
15	(ADD			
16	STIR		Rot.speed 2000 /min	
17	PURGE	30.0		
18	(REP			
19	SEGMENT		Segm.name NiCoAdSV	
20	REP) 1			
21	ADD>M		Soln.name Standard	V.add 0.010 mL
22	ADD) 2			
23	RINSE		Cycles 1	V.rinse 20.000 mL
24	SETCTRL	1.0	Code *****1**	Message
25	SETCTRL		Code *****0**	Message
26	REM	20.0	move to rinsing position	
27	SETCTRL	1.0	Code *****1**	Message
28	SETCTRL		Code *****0**	Message
29	REM	300.0	rinse tubes	
30	SETCTRL	1.0	Code *****1**	Message
31	SETCTRL		Code *****0**	Message
32	REM		stop pump	
33	RINSE		Cycles 1	V.rinse 40.000 mL
34	RINSE		Cycles 2	V.rinse 15.000 mL
35	END			


3.7.3 Conditions to be met

For the correct intercommunication of the 813 Compact Autosampler and the 746 VA Trace Analyzer the following conditions must be met:

- Method 2 of the 813 Compact Autosampler has to be selected.
- Always start the method of the 813 Compact Autosampler first, then you can start the 746 VA Trace Analyzer.
- Ensure that the control line 2 (Remote output 2) of the 746 VA Trace Analyzer is reset (inactive state) before you start the sample series.

4 Appendix

4.1 Error Messages



An error occurring during a sample series is displayed by a blinking red LED. The number of flashes indicates the error number. If an error occurs during processing of a sample series, the changer will then be switched into the 'HOLD' state. The current method has to be halted by pressing <STOP>.

The list of possible error messages and their causes:

1* RAM defect	Call Metrohm-Service.
3* battery low	The battery for the permanent storage of the settings must be replaced.
6* changer low power	The power supply cannot deliver enough power for the simultaneous operation of all components currently in use (pump and lift).
7* rack data missing	No sample rack is in position or no rack data can be found for the sample rack that is in place.
8* invalid rack code	The rack code read by the changer could not be found in the internal position tables.
10* raise lift first	Turning of a rack could not be carried out because the lift was below the rest (shift) position.
12* changer overload	Load or resistance too large to turn the rack.

4.2 Sequences

The available Autosampler methods are described briefly in section 3.4 The following listings explain the sequences in detail.

Method 1

Single determination with air gap after each run, needle lifted after sample transfer.

```

813 Compact Autosampler      813.0010 ← Program version
parameters
  method          813_1      ← Method name
  number of samples:  *      ← Number of samples unlimited
>start sequence
  1 SAMPLE:  =          1      ← First rack position
>sample sequence
  1 SCN:Rm:          *****1** ← Await 1st Remote signal
  2 MOVE  1      :      sample ← Move next sample to tower
  3 LIFT:  1      :          work mm ← Lift to work position / immerse needle
  4 SCN:Rm:          *****1** ← Await 2nd Remote signal
  5 PUMP  1.1      :          on ← Switch on pump
  6 WAIT                          5 s ← Waiting time
  7 SCN:Rm          :      *****1** ← Await 3rd Remote signal
  8 WAIT                          12 s ← Waiting time for sample transfer
  9 PUMP  1.1      :          off ← Switch off pump
 10 LIFT:  1      :          0 mm ← Lift to rest position / raise needle
 11 PUMP  1.1      :          6 s ← Switch on pump for 6 seconds
>final sequence
>changer settings          ----- Settings for changer functions -----
  rack number          0
  lift rate            12 mm/s
  shift rate           20
  shift direction:     auto.
>manual stop
  -----

```

Method 2

Single determination, no air gap, needle remains in sample after sample transfer.

```

813 Compact Autosampler      813.0010 ← Program version
parameters
  method          813_2      ← Method name
  number of samples: *      ← Number of samples unlimited
>start sequence
  1 SAMPLE: =          1      ← First rack position
>sample sequence
  1 SCN:Rm:          *****1** ← Await 1st Remote signal
  2 MOVE 1          :      sample ← Move next sample to tower
  3 LIFT: 1          :      work mm ← Lift to work position / immerse needle
  4 SCN:Rm:          *****1** ← Await 2nd Remote signal
  5 PUMP 1.1        :      on      ← Switch on pump
  6 WAIT           :      5 s      ← Waiting time
  7 SCN:Rm          :      *****1** ← Await 3rd Remote signal
  8 WAIT           :      12 s     ← Waiting time for sample transfer
  9 PUMP 1.1        :      off     ← Switch off pump
>final sequence
>changer settings          ----- Settings for changer functions -----
  rack number          0
  lift rate            12 mm/s
  shift rate           20
  shift direction:     auto.
>manual stop
  -----

```

Method 3

Double determination of each sample, no air gap, needle remains in sample after sample transfer.

```

813 Compact Autosampler      813.0010 ← Program version
parameters
method                        813_3      ← Method name
number of samples:           *          ← Number of samples unlimited
>start sequence
  1 SAMPLE: =                1          ← First rack position
>sample sequence
  1 SCN:Rm:                   *****1** ← Await 1st Remote signal
  2 MOVE 1 :                   sample    ← Move next sample to tower
  3 LIFT: 1 :                   work mm  ← Lift to work position / immerse needle
  4 SCN:Rm:                   *****1** ← Await 2nd Remote signal
  5 PUMP 1.1 :                   on      ← Switch on pump
  6 WAIT                       5 s      ← Waiting time
  7 SCN:Rm :                   *****1** ← Await 3rd Remote signal
  8 WAIT                       12 s     ← Waiting time for 1st sample transfer
  9 PUMP 1.1 :                   off     ← Switch off pump
 10 SCN:Rm:                   *****1** ← Await 1st Remote signal again
 11 WAIT                       5 s      ← Waiting time
 12 SCN:Rm:                   *****1** ← Await 2nd Remote signal
 13 PUMP 1.1 :                   on      ← Switch on pump
 14 WAIT                       5 s      ← Waiting time
 15 SCN:Rm :                   *****1** ← Await 3rd Remote signal
 16 WAIT                       12 s     ← Waiting time for 2nd sample transfer
 17 PUMP 1.1 :                   off     ← Switch off pump
>final sequence
>changer settings             ----- Settings for changer functions -----
rack number                   0
lift rate                     12 mm/s
shift rate                    20
shift direction:              auto.
>manual stop
-----

```

Method 4

Triple determination of each sample, no air gap, needle remains in sample after sample transfer.

```


813 Compact Autosampler      813.0010 ← Program version
parameters
  method          813_4      ← Method name
  number of samples:  *      ← Number of samples unlimited
>start sequence
  1 SAMPLE:  =          1      ← First rack position
>sample sequence
  1 SCN:Rm:      *****1** ← Await 1st Remote signal
  2 MOVE 1      :      sample ← Move next sample to tower
  3 LIFT: 1      :      work mm ← Lift to work position / immerse needle
  4 SCN:Rm:      *****1** ← Await 2nd Remote signal
  5 PUMP 1.1    :      on      ← Switch on pump
  6 WAIT                5 s    ← Waiting time
  7 SCN:Rm        :      *****1** ← Await 3rd Remote signal
  8 WAIT                12 s    ← Waiting time for 1st sample transfer
  9 PUMP 1.1      :      off    ← Switch off pump
 10 SCN:Rm:        *****1** ← Await 1st Remote signal again
 11 WAIT                5 s    ← Waiting time
 12 SCN:Rm:        *****1** ← Await 2nd Remote signal
 13 PUMP 1.1      :      on      ← Switch on pump
 14 WAIT                5 s    ← Waiting time
 15 SCN:Rm        :      *****1** ← Await 3rd Remote signal
 16 WAIT                12 s    ← Waiting time for 2nd sample transfer
 17 PUMP 1.1      :      off    ← Switch off pump
 18 SCN:Rm:        *****1** ← Await 1st Remote signal again
 19 WAIT                5 s    ← Waiting time
 20 SCN:Rm:        *****1** ← Await 2nd Remote signal
 21 PUMP 1.1      :      on      ← Switch on pump
 22 WAIT                5 s    ← Waiting time
 23 SCN:Rm        :      *****1** ← Await 3rd Remote signal
 24 WAIT                12 s    ← Waiting time for 3rd sample transfer
 25 PUMP 1.1      :      off    ← Switch off pump
>final sequence
>changer settings          ----- Settings for changer functions -----
  rack number              0
  lift rate                 12 mm/s
  shift rate                20
  shift direction:         auto.
>manual stop
  -----

```


4.3 Technical specifications

Dimensions	Height: 0.50 m, Width: 0.28 m, Depth: 0.49 m	
Weight	12.3 kg (incl. sample rack)	
Material	Housing:	Metal case, multiple enameling
	Splash protection:	Polymethylmethacrylate (PMMA)
	Sample rack:	Polypropylene (PP)
Lift path	ca. 125 mm	
Lift	Load:	ca. 30 N
	Stroke speed:	12 mm/s
Turntable	Rotational speed:	20 angular degrees/s
Pump	<i>Pump type</i>	1-channel peristaltic pump with rotational speed of 20 /min
	<i>Pump capacity</i>	(with water, without counterpressure) with 6.1826.040 pump tubing: typ. 0.9...1.1 mL/min (depends on contact pressure)
	<i>Pressure</i>	max. 1.5 bar (0.15 MPa)
	<i>Pump tubing material</i>	PVC (Tygon®)

**Remote-
RS232 interface** Parallel interface for remote control by external devices



The +5 V supply line may be charged with 20 mA maximally.



Temperatures	Nominal operating range 5...40 °C		
	at 20...80 % atmospheric humidity		
	Storage, transport -20...+60 °C		
	60 °C	at atmospheric humidity	<50%
	50 °C	" "	<85%
40 °C	" "	<95%	

Mains connection	Voltage 100...120 V, 220...240 V		
	Frequency 50...60 Hz		
	Power consumption 40 VA		
	Fuse 0.5 AT (110 V), 0.25 AT (220 V)		

Safety specifications

Construction / Testing According to IEC 1010 / EN 61010 /
UL 3101-1, protection class 1

Degree of protection IP 22

The Instructions for Use include information and warnings which must be heeded by the user to assure safe operation of the instrument.

Electromagnetic compatibility (EMC)

Emitted interference

Standards met:

EN55011 (class B), EN55022 (class B), EN50081-1 01.92

Immunity to interference

Standards met:

IEC801-2/IEC1000-4-2 (class 4), IEC801-3/IEC1000-4-3 (class 2),
IEC801-4/IEC1000-4-4 (class 3), IEC801-5/IEC1000-4-5 (class 2/3),
IEC801-6/IEC1000-4-6 (class 2), EN55011 (class B), EN55022 (class
B), EN50081-1/2 01.92, EN50082-1 01.92, EN61000-3, EN61316-1
03.97

4.4 Maintenance and servicing

4.4.1 Maintenance by Metrohm service

is best done as part of an annual service performed by specialists from the Metrohm company. If work is frequently performed with caustic and corrosive chemicals, it may be necessary to shorten the interval between servicing.

The Metrohm service department is always willing to offer expert advice on the maintenance and servicing of all Metrohm instruments

4.4.2 Care of the unit

The 813 Compact Autosampler requires proper care and attention. Excessive contamination of the instrument could possibly lead to malfunctions and a shorter service life of the inherently rugged mechanical and electronic parts.

Wipe up spilled chemicals and solvents immediately. The connectors (in particular the power supply) should be protected from contamination. The 813 Compact Autosampler should never be operated without plug cover



The unit has been constructed in such a way as to virtually eliminate the possibility of penetration of corrosive media into the interior of the instruments. If such a situation does occur, disconnect the mains plug of the 766 Compact Autosampler immediately to prevent extensive damage to the instrument electronics. Inform Metrohm service if your instrument has been damaged in such a way.



The instrument must not be opened by untrained personnel. Please comply with the safety notes in section 2.2.

4.4.3 Replacing the pump tubing

Pump tubings are consumable material with a limited lifetime and should be replaced at regular intervals (approx. every 2 weeks under continuous use).

The working life of pump tubing depends to a considerable extent on the contact pressure. This is why the contact pressure must be correctly set as described in *section 2.3.7*. If the pump is to remain switched off for a lengthy period of time the tubing cartridges should be raised completely by loosening snap-action lever **15** on the right-hand side (the pre-set contact pressure remains unchanged).

To replace a pump tubing proceed as follows:

1 Remove old pump tubing

- Press contact pressure lever **11** on the tubing cartridge down as far as it will go.
- Release tubing cartridge **11** from holding clamp **14** by pressing down snap-action lever **15** and remove from mounting pin **10** at the 813 Compact Autosampler (see Fig. 1).
- Remove old pump tubing.

2 Insert new pump tubing

- Insert the new pump tubing **8** (6.1826.040) in the tubing cartridge as shown in Fig. 9. The Stopper **32** must click into the corresponding holder on the left-hand side of the tubing cartridge.
- Place the tubing cartridge on mounting pin **10** and press down on the right-hand side until snap-action lever **15** clicks into position on holding clamp **14**. Take care that no kinks are formed in the pump tubing.

3 Set contact pressure

- Press contact pressure lever **13** upwards until the solution just starts to be drawn in. Then press contact pressure lever upwards until it clicks once more to obtain optimal contact pressure.



The 6.1826.0X0 pump tubing is made of PVC and must not be used for rinsing with solutions which contain acetone. In such cases, rinse with different pump tubing or a different pump.

4.5 Warranty and conformity

4.5.1 Warranty

The warranty on our products is limited to defects that are traceable to material, construction or manufacturing error which occur within 12 months from the day of delivery. In this case, the defects will be rectified in our workshops free of charge. Transport costs are to be paid by the customer.

For day and night operation, the warranty is limited to 6 months.

Glass breakage in the case of electrodes or other parts is not covered by the warranty. Checks which are not a result of material or manufacturing faults are also charged during the warranty period. For parts of outside manufacture insofar as these constitute an appreciable part of our instrument, the warranty stipulations of the manufacturer in question apply.

With the regard to the guarantee of accuracy, the technical specifications in the instruction manual are authoritative.

Concerning defects in material, construction or design as well as the absence of guaranteed features, the orderer has no rights or claims except those mentioned above.

If damage of the packaging is evident on receipt of a consignment or if the goods show signs of transport damage after unpacking, the carrier must be informed immediately and a written damage report demanded. Lack of an official damage report releases Metrohm from any liability to pay compensation.

If any instruments and parts have to be returned, the original packaging should be used if at all possible. This applies above all to instruments, electrodes, burette cylinders and PTFE pistons. Before embedment in wood shavings or similar material, the parts must be packed in a dustproof package (for instruments, use of a plastic bag is imperative). If open assemblies are enclosed in the scope of delivery that are sensitive to electromagnetic voltages (e.g. data interfaces etc.) these must be returned in the associated original protective packaging (e.g. conductive protective bag). (Exception: assemblies with built-in voltage source belong in a non-conductive protective packaging).

No warranty responsibility whatsoever will be accepted by Metrohm for damage which arises as a result of non-compliance with these instructions.

4.5.2 Certificate of Conformity and System Validation

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

Name of commodity:	813 Compact Autosampler
System software:	stored in ROMs
Name of manufacturer:	Metrohm Ltd., Herisau, Switzerland
Technical information:	Voltages: 100...120, 220...240 V Frequency: 50...60 Hz

This Metrohm instrument has been built and has undergone final type testing according to the standards:

Electromagnetic compatibility

IEC 801-2 / level 3, IEC 801-3 / level 2, IEC 801-4 / level 3,
EN 55011 / class B, EN 55022 / class B, EN 50081-1/2 1992,
EN 50082-1 1997

Security specifications

IEC 1010, EN 61010, UL 3101-1

It has also been certified by the Swiss Electrotechnical Association (SEV), which is member of the International Certification Body (CB/IEC).

The technical specifications are documented in the instruction manual.

The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance. The features of the system software are documented in the instruction manual.

Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.

Herisau, March 5, 2001



Dr. J. Frank

Ch. Buchmann

Development Manager

Production and
Quality Assurance Manager

Ionenanalytik • Analyse des ions • Ion analysis • Análisis iónico
813 Compact Autosampler



EU Declaration of conformity

The METROHM AG company, Herisau, Switzerland hereby certifies that the instrument:

813 Compact Autosampler

meets the requirements of EC Directives 89/336/EWG and 73/23/EWG.

Source of the specifications:

EN 50081	Electromagnetic compatibility, basic specification Emitted Interference
EN 50082-1	Electromagnetic compatibility, basic specification Interference Immunity
EN 61010	Safety requirements for electrical laboratory measurement and control equipment

Description of the instrument:

Sample changer for automatic processing of sample series with ion chromatographic methods

Herisau, March 5, 2001



Dr. J. Frank

Ch. Buchmann

Development Manager

Production and
Quality Assurance Manager






4.6 Standard equipment







Subject to changes! All dimensions are given in mm.








The following parts are included:



2.813.0010 Compact Autosampler (for IC applications)

Quant.	Order No.	Description												
2	6.1826.040	Pump tubing made of PVC (Tygon®); with 2 firmly attached black-black stoppers; i.d. = 0.76 mm, e.d. = 2.27 mm 												
1	6.1831.050	PEEK Capillary Length = 40 cm 												
2	6.1831.060	PEEK Capillary Length = 100 cm 												
1	6.2041.750	Sample rack for 36 sample tubes 6.2743.050 (11 mL) or 6.2743.040 (2,5 mL)												
1	6.2122.0X0	Mains cable to customers specifications: <table border="0"> <tr> <td><u>Cable socket</u></td> <td><u>Cable plug</u></td> <td></td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type SEV 12 (CH...)</td> <td>6.2122.020</td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type CEE (7), VII (D...)</td> <td>6.2122.040</td> </tr> <tr> <td>Type CEE (22), V</td> <td>Type NEMA 5-15 (USA...)</td> <td>6.2122.070</td> </tr> </table>	<u>Cable socket</u>	<u>Cable plug</u>		Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020	Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040	Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070
<u>Cable socket</u>	<u>Cable plug</u>													
Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020												
Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040												
Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070												
1	6.2141.130	Connection cable (Remote) Connection cable 813 Compact Autosampler – Compact IC 761, Personal IC 790 or 732 IC Detector												
1	6.2142.110	Keyboard for 813 Compact Autosampler												
1	6.2621.060	Open-end spanner 5/16" 												
2	6.2621.090	Open-end spanner 1/2" 												







Quant.	Order No.	Description
1	6.2621.100	Allen key 3 mm For Allen screws on sample rack and for splash protection 
1	6.1835.010	PEEK needle For aspiration of solutions from sealed sample tubes 
1	6.2743.057	PP Sample tubes (11 mL) For 6.2041.750 sample rack set of 200
1	6.2743.077	PE Caps For sealing the 6.2743.057 sample tubes set of 200
1	6.2744.010	PEEK Compression fitting For the connection of 6.1831.0X0 PEEK capillaries or 6.1803.0X0 PTFE capillaries; set of 5 
1	6.2744.030	PEEK Coupling Connection between 6.2744.010 PEEK compression fitting and 6.1826.0X0 pump tubing set of 4 
1	6.2751.070	Splash protection Must be installed at the tower of the 813 Compact Autosampler.
1	6.2752.010	Plug cover Must be installed at the 813 Compact Autosampler.
1	6.2755.000	Tubing cartridge For peristaltic pump at 813 Compact Autosampler
1	8.813.1003	Instructions for Use (English) for 813 Compact Autosampler

2.813.0020 Compact Autosampler (for VA applications)


Quant.	Order No.	Description												
2	6.1826.020	Pump tubing made of PVC (Tygon®); with 2 firmly attached black-black stoppers; i.d. = 0.76 mm, e.d. = 2.27 mm 												
1	6.1831.050	PEEK Capillary Length = 40 cm 												
1	6.1822.410	PEEK Capillary 1x M6 Length = 120 cm 												
1	6.2041.750	Sample rack for 36 sample tubes 6.2743.050 (11 mL) or 6.2743.040 (2,5 mL)												
1	6.2122.0X0	Mains cable to customers specifications: <table border="0"> <tr> <td><u>Cable socket</u></td> <td><u>Cable plug</u></td> <td></td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type SEV 12 (CH...)</td> <td>6.2122.020</td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type CEE (7), VII (D...)</td> <td>6.2122.040</td> </tr> <tr> <td>Type CEE (22), V</td> <td>Type NEMA 5-15 (USA...)</td> <td>6.2122.070</td> </tr> </table>	<u>Cable socket</u>	<u>Cable plug</u>		Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020	Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040	Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070
<u>Cable socket</u>	<u>Cable plug</u>													
Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020												
Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040												
Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070												
1	6.2142.110	Keyboard for 813 Compact Autosampler												
1	6.2621.060	Open-end spanner 5/16" 												
2	6.2621.090	Open-end spanner 1/2" 												
1	6.2621.100	Allen key 3 mm For Allen screws on sample rack and for splash protection 												
1	6.1835.030	PEEK tube For aspiration of solutions from open sample tubes 												

Quant.	Order No.	Description	
1	6.2743.057	PP Sample tubes (11 mL) For 6.2041.750 sample rack set of 200	
1	6.2744.010	PEEK Compression fitting For the connection of 6.1831.0X0 PEEK capillaries or 6.1803.0X0 PTFE capillaries; set of 5	
1	6.2744.030	PEEK Coupling Connection between 6.2744.010 PEEK compression fitting and 6.1826.0X0 pump tubing set of 4	
1	6.2751.070	Splash protection Must be installed at the tower of the 813 Compact Autosampler.	
1	6.2752.010	Plug cover Must be installed at the 813 Compact Autosampler.	
1	6.2755.000	Tubing cartridge For peristaltic pump at 813 Compact Autosampler	
1	8.813.1003	Instructions for Use (English) for 813 Compact Autosampler	

4.7 Optional accessories

Order No.	Description	
6.1835.020	PEEK needle For aspiration of solutions from open sample tubes	
6.1835.000	PEEK needle For aspiration of solutions from open sample tubes	
6.2624.000	Steel needle For aspiration of solutions from sealed sample tubes	
6.1831.040	PEEK capillary Length = 15 cm	
6.2833.000	Needle holder (1/8") for PEEK needle 6.1835.010 and 6.1835.020	
6.2833.010	Needle holder (1/16") For steel needle 6.2624.000 and PEEK needle 6.1835.000	
6.2743.040	PP Sample tubes (2,5 mL) For sample rack 6.2041.750 set of 2000	
6.2743.047	PP Sample tube (2,5 mL) For sample rack 6.2041.750 set of 200	

6.2743.050	PP Sample tube (11 mL) For sample rack 6.2041.750 set of 2000
6.2743.077	PE Caps For sealing the 6.2743.057 sample tubes set of 200
6.2743.060	PE Caps For sealing the 6.2743.0x0 sample tubes set of 1000
6.2125.120	Adaptor for remote connection with 5 connections for remote cables
6.2141.150	Connection cable (Remote) Connection cable 813 Compact Autosampler – 757 Computrace/ 731 Relay Box/Dosimat
6.1826.020	Pump tubing made of PVC (Tygon®); with two firmly attached blue-blue stoppers; i.d. = 1.6 mm, e.d. = 2.27 mm. For use with preconcentration columns.



5 Index

- <HOLD> key 5
 <START> key 5
 <STOP> key 5
- 6.1835.010 PEEK needle
 Insertion 10
 6.1835.020 PEEK tube
 Insertion 10
 6.2141.020 Remote cable
 Connection 18
 6.2141.130 Remote cable
 Connection 17
 746 VA Trace Analyzer
 Interconnection 18
- A**
-
- Abort the configuration
 procedure 22
 Accept setting 22
 Accessories
 Attaching 9
 Adapter 6.2125.120
 Ordering designation 47
 Adjusting screws 13
 Adjusting the sample rack 13
 Allen key 6.2621.100
 Installing the splash
 protection 11
 Ordering designation 43, 44
 Appendix 30
 Attaching the accessories .. 9
- B**
-
- Battery 30
- C**
-
- Cable 6.2141.130
 Ordering designation 42, 44
 Cable 6.2141.020
 Interconnection 19
 Changer overload 30
 Coding magnet 12
 Configuration 21
 Connecting
 Keyboard 9
 Connection plug
 Mains connection 8
 Contact pressure 14, 38
 Contact pressure lever
 Install tubing cartridge ... 38
 Contact pressure lever
 Figure 3, 15
 Coupling
 Figure 15
 Installation 14
 Ordering designation 43, 45
 Remove 15
- Course of configuration 23
- D**
-
- Degree of protection 36
 Delivery check 6
 Dimensions 35
 Double determination 26
- E**
-
- Earthing 8
 Electromagnetic compatibility
 36, 40, 41
 Emitted interference 36
 Error messages 30
 Error state 5
 EU Declaration of
 conformity 41
- F**
-
- Factory defaults 24
 Fuse 7
 Fuse holder
 Figure 3
 Fuse holder
 Replace fuses 8
 Fuse holder
 Setting the mains voltage ..
 7
 Fuses 8
- H**
-
- HOLD 5
 HOLD key 20
 Holding clamp
 Figure 3
 Holding clamp
 Install tubing cartridge ... 14
 Holding screw
 Figure 11
 Holding screws
 Installing the splash
 protection 11
- I**
-
- IC Software 27
 Immunity to interference ... 36
 Initialization 24
 Initialize 20
 Injection valve 27
 Input 35
 Installation 6
 Installing
 Splash protection 11
 Installing the pump tubing 15
 Instructions for Use
 8.813.1003 43, 45
- Instrument description 1
 Instrument status 20
 Interconnection
 746 VA Trace Analyzer ... 18
 IC systems 17
 other instruments 19
 Interrupt 20
 Introduction 1
 Invalid rack code 30
- K**
-
- Keyboard
 Connecting 9
 Figure 5
 Ordering designation 42, 44
 Keyboard connection 4
 Figure 4
 Keyboard functions 20, 21
- L**
-
- Leaks 6
 LED
 configuration 22
 green 5, 21
 red 5, 22
 Lifetime of pump tubings.. 14
 Lift 24, 35
 Figure 3
 Location 6
 Low power 30
- M**
-
- Mains cable
 Attaching 8
 Ordering designation 42, 44
 Mains connection 7, 36
 Mains connection plug
 Figure 3
 Mains switch
 Figure 3
 Switching on/off 9
 Maintenance 37
 Material 35
 Measuring head 16
 Method
 Selection 22, 23
 Methods 21, 25
 Metrohm service department
 37
 Mounting pin
 Install tubing cartridge ... 14
 Mounting pin
 Figure 3
 Mounting screws
 Figure 4
 Needle
 Figure 3, 10

Installation 10
 Needle guide
 Figure 3
 Needle holder
 Figure 10
 Needle installation 10
 Ordering designation 46
 Needle position 13

O

Open-end spanner
 6.2621.060 42, 44
 Optional accessories 46
 Output 35

P

Packaging 6
 PE Caps 6.2743.060
 Ordering designation 47
 PE Caps 6.2743.077
 Ordering designation 43, 45
 PEEK capillary
 Figure 10
 PEEK capillary
 Connection to injection
 valve 15
 Connection to 747 VA Stand
 16
 Figure 3, 15
 Installation 14
 Ordering designation 42, 44
 PEEK capillary 6.1831.040
 Ordering designation 46
 PEEK compression fitting
 Connect PEEK capillary . 14
 Figure 3, 10, 15
 Fix needle 10
 Installation 14
 Needle installation 10
 Ordering designation 43, 45
 PEEK needle 6.1835.020
 Ordering designation 46
 PEEK needle 6.2624.000
 Ordering designation 43, 44
 Plug cover
 Mounting 9
 Ordering designation 43, 45
 Plug cover
 Figure 3, 4
 Pressure 35
 Protection class 36
 Protective earth 8
 Pump
 Capacity 35
 Technical specifications 35
 Tubing material 35
 Type 35
 Pump drive
 Figure 3

Pump tubing
 Connection to needle 14
 Figure 3, 15
 Installation 14
 Life time 37
 Ordering designation 42, 44
 Replacement 37
 Pump tubing 6.1826.020 .. 47

R

Rack data missing 30
 Raise lift first 30
 RAM defect 30
 RAM initialization 24
 ready 5
 Rear view 4
 Remote connection
 Figure 3
 Remote interface 35
 Instrument connections . 17
 Remote-Adapter 6.2125.120
 17
 Reset 24
 RESET 5

S

Safety note 2
 Safety specifications 36
 Sample rack
 Adjustment 12
 Figure 3, 4
 Ordering designation 42, 44
 Placing 12
 Sample tube
 Figure 3, 4
 Sample tubes
 Ordering designation
 43, 45, 46, 47
 Screw thread for splash
 protection
 Installation 11
 Security specifications 40
 Sequence 25, 31
 Sequences 31
 Serial number
 Figure 4
 Servicing 37
 Setting the mains voltage ... 7
 Setting up 6
 Side view 2
 Single determination 25
 Snap-action lever
 Figure 3, 15
 Install tubing cartridge ... 14
 Software
 IC Chromatograph 27
 Splash protection
 Figure 3, 4, 11

Installation 11
 Ordering designation 43, 45
 Remove protection film .. 11
 Standard equipment 42
 Start configuration 22
 Start data acquisition 27
 Start key 20
 Status 5, 20
 Steel holder
 Figure 3
 Steel needle 6.2624.000
 Ordering designation 46
 Stop key 20
 Stopper
 Figure 15
 Install pump tubing 38
 Pump tubing installation 14
 Switching on 24
 Switching on/off 9
 Synchronize 27
 Syringe tubing
 Connection to waste
 container 15

T

Technical specifications ... 35
 Temperatures 36
 Terminate
 Configuration 22
 Tower
 Figure 3, 4
 Transport 6
 Transport damage 39
 Triple determination 26
 Tubing cartridge
 Figure 3, 4, 15
 Installation 38
 Install pump tubing 14
 Ordering designation 43, 45
 Tubing connections 14
 Turntable 35

V

VA Computrace 757 19
 VA Stand 747
 Connection of pump tubing
 16

W

Warranty 39
 Waste container 15
 Weight 35
 Work position 21, 24
 working position 22