

816 IC Eluent Selector





Ionanalytik

CH-9101 Herisau/Switzerland

E-Mail info@metrohm.com

Internet www.metrohm.com

816 IC Eluent Selector

Instructions for Use

Teachware
Metrohm AG
Oberdorfstrasse 68
CH-9101 Herisau
teachware@metrohm.com

1st Edition 2002

These instructions are protected by copyright. All rights reserved.

Although all the information given in these instructions has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please inform the author at the address given above.

Table of Contents

1	Introduction	1
1.1	Description of the instrument	1
1.2	Parts and controls	2
1.3	Information about the Instructions for Use	4
1.3.1	Organization	4
1.3.2	Notation and pictograms	5
1.4	Safety information	6
1.4.1	Electrical safety	6
1.4.2	General safety rules	6
2	Installation	7
2.1	Setting up the instrument	7
2.1.1	Packaging	7
2.1.2	Checks	7
2.1.3	Location	7
2.2	Mains connection	7
2.2.1	Setting the mains voltage	7
2.2.2	Fuses	9
2.2.3	Mains cable and mains connection	9
2.2.4	Switching the instrument on/off	9
2.3	Electrical connection	10
2.3.1	Operation in a modular IC system	10
2.3.2	Operation with a 761 Compact IC	10
2.4	Valve connections	11
2.5	Software installation	12
3	Operation	13
3.1	Manual operation	13
3.2	Operation via «IC Net»	13
3.2.1	Symbol for 816 IC Eluent Selector	14
3.2.2	Instrument parameters in the "816 IC Eluent Selector" window	14
4	Annex	17
4.1	Technical data	17
4.1.1	Parts and controls	17
4.1.2	RS 232 Interface	17
4.1.3	Mixer speed	17
4.1.4	Mains connection	18
4.1.5	Safety specifications	18
4.1.6	Electromagnetic compatibility (EMC)	18
4.1.7	Ambient temperature	18
4.1.8	Housing	19
4.2	Standard equipment	19
4.3	Optional accessories	20
4.4	Validation / GLP	21
4.5	Warranty and conformity	22
4.5.1	Warranty	22
4.5.2	EU Declaration of Conformity	23
4.5.3	Certificate of conformity and system validation	24
4.6	Index	25

List of illustrations


Figure 1:	Front panel of the 816 IC Eluent Selector	2
Figure 2:	Rear panel of the 816 IC Eluent Selector	3
Figure 3:	Setting the mains voltage.....	8
Figure 4:	Connecting the 816 IC Eluent Selector to 762 IC Interface.....	10
Figure 5:	Connecting the 816 IC Eluent Selector to a PC	10
Figure 6:	Flow diagram of the IC Eluent Selector.....	11

1 Introduction

1.1 Description of the instrument

The **816 IC Eluent Selector** has two 3-way valves and can be integrated in an IC system via an RS232 interface. When used in a Metrohm IC system it can be operated completely by remote control.

Up to three different eluents can be used with the 816 IC Eluent Selector. In this way the 816 IC Eluent Selector makes it easy to use step gradients in order to optimize the separating performance of your system with difficult samples. A further frequent application is its use in a 2-channel system for investigating samples which require different separating systems for their complete analysis, e.g. the determination of cations and anions in a single sample. In combination with the 2-channel version of the 812 IC Valve Unit and the 761 Compact IC with suppressor, the 816 IC Eluent Selector forms a compact IC system with which the cations and anions in a single sample can be determined.



1.2 Parts and controls




Figure 1: Front panel of the 816 IC Eluent Selector

Keys

for manual operation of the 816 IC Eluent Selector.

LEDs

The LED shows the instrument status.

1 A, B, C

These keys are only active in the manual mode. They can be used for the manual selection of the required eluent. The LED **2** for the selected eluent lights up.

2 A, B, C

The LEDs show which eluent has been selected.
LED lights up: valve is open; the particular eluent can pass through the IC Eluent Selector.
LED is off: valve is closed.

3 EXT.

Switches between **external** control and **manual** control of the instrument.

4 EXT.

The LED shows whether the instrument is under **external (LED lights up)** or **manual (LED is off)** control.

5 Valve I

Switches eluents **A** and **B**.

6 Mixing chamber

Has no function in the 816 IC Eluent Selector; the eluents simply pass through this component.

7 Valve II

Switches eluent **C**.




Figure 2: Rear panel of the 816 IC Eluent Selector

<p>8 Mains switch for switching the instrument on and off: I = ON 0 = OFF Mains connection: see Section 2.2.</p>	<p>9 Mains connection Mains connection: see Section 2.2.</p>
<p>10 Model and serial number</p>	<p>11 Ground connection</p>
<p>12 Fuse holder Changing the fuses: see Section 2.2.1.</p>	<p>13 Mixer speed Dip switch for setting the stirring rate in the mixing chamber: see Section 4.1.3. We recommend that the mixing chamber is kept <u>switched off</u> (all switches down).</p>
<p>14 RS 232 RS 232 interface for remote control of the 816 IC Eluent Selector: see Section 4.1.2.</p>	

1.3 Information about the Instructions for Use



Please study the instructions carefully before you start to use the 816 IC Eluent Selector. The instructions contain information and warnings that must be observed by the user in order to guarantee the safe use of the instrument. Please keep these instructions near the instrument so that they are always to hand when required.

1.3.1 Organization

These **8.816.1003 Instructions for Use** for the 816 IC Eluent Selector provide you with a comprehensive overview of the installation, start up, operation, troubleshooting and technical specifications of this instrument. The instructions are arranged as follows:

Section 1 Introduction

General description of the instrument, parts and controls, safety information

Section 2 Installation

Setup, mains connection, connection to the IC system, flow diagram

Section 3 Operation

Manual operation, remote control





Section 4 Annex

Technical data, standard equipment, options, warranty, declaration of conformity, index

In order to find the information you require about the 816 IC Eluent Selector you should either use the **Contents** or the **Index**.

1.3.2 Notation and pictograms

The following notation and pictograms (symbols) are used in these Instructions:

Program	Menu item, parameter or input value in the «IC Net» program
<OK>	Button in the «IC Net» program
15	Numbers of the parts and controls: 816
	Danger/Warning This symbol indicates a possible risk of injury to the user and possible damage to the instrument or its components from electricity.
	Danger/Warning This symbol indicates a possible risk of injury to the user and possible damage to the instrument or its components.
	Attention This symbol indicates important information that you should read before continuing.
	Information This symbol indicates additional information and tips which may be of particular use to you.

1.4 Safety information

**Warning!**

This instrument should only be used in accordance with the information given in these installation instructions.

1.4.1 Electrical safety

Electrical safety when handling the 816 IC Eluent Selector is guaranteed within the framework of the IEC 61010 Standard. The following point must be observed:

- **Mains connection**



***Mains connection** and checking the **fuses** must be carried out in accordance with the instructions given in Section 2.2.*

- **Opening the 816 IC Eluent Selector**



*When the 816 IC Eluent Selector is connected to the mains it must not be opened, nor should any components be removed, as otherwise you run the risk of contacting parts which are “live”. Always disconnect the instrument from all voltage sources before opening it and make sure that the **mains cable has been removed from mains connection 9!***

- **Protection against electrostatic charges**



Electronic components are sensitive to electrostatic charges and can be destroyed by a discharge. Before you touch any component inside the 816 IC Eluent Selector you should ground yourself and your tool by touching a grounded object (e.g. instrument housing or radiator) to eliminate any existing electrostatic charges.

1.4.2 General safety rules

- **Handling solutions**



Check all inlets and outlets for leaks at regular intervals. Observe the appropriate regulations for the handling and disposal of flammable and/or toxic solutions.

2 Installation

2.1 Setting up the instrument

2.1.1 Packaging

The 816 IC Eluent Selector and the separately packed accessories are delivered in special packaging that provides excellent protection. It contains impact-absorbing plastic foam. The instrument itself is packed in a dustproof evacuated polyethylene bag. Please store this packaging in a safe place; it is the only way in which the safe transport of the instrument can be guaranteed.

2.1.2 Checks

Please check that the delivery is complete and undamaged immediately on receipt (compare with delivery note and list of accessories given in Section 4.2). If transport damage is evident please refer to the information given in Section 4.5.1 Warranty.

2.1.3 Location

Place the instrument on a suitable vibration-free laboratory bench, protected as much as possible from corrosive atmospheres and contact with chemicals.

Choose a location where the temperature is usually between +5 °C and +45 °C. The instrument should be protected against excessive variations in temperature and direct sunlight.

2.2 Mains connection



Please observe the following rules when connecting the instrument to the electricity supply. If the instrument is operated with an incorrectly set mains voltage and/or an incorrect mains fuse then it represents a fire hazard!

2.2.1 Setting the mains voltage

Before you switch on the 816 IC Eluent Selector for the first time please check that the mains voltage set on the instrument (see Figure 2) corresponds to your local mains voltage. If this is not the case then you must alter the mains voltage as follows:

1 Pull out mains cable

Remove the mains cable from mains supply connection **9** of the 816 IC Eluent Selector.

2 Remove the fuse holder

Use a screwdriver to loosen fuse holder **12** below the mains supply connection and remove it completely.

3 Check the fuse

Carefully remove the built-in fuse for the intended voltage from the fuse holder and check its specifications (the position of the fuse in the fuse holder is indicated by the white arrow beside the voltage range):

100...120 V 1.0 A (slow blow) Metrohm No. U.600.0016

220...240 V 0.5 A (slow blow) Metrohm No. U.600.0013

4 Insert fuse

Exchange the fuse if necessary and replace it in the fuse holder.

5 Insert fuse holder

Depending on the required mains voltage, insert the fuse holder in the 816 IC Eluent Selector so that the correct voltage range can be read normally and the white arrow beside it points to the white bar printed below the fuse holder (see Figure 3).




Figure 3: **Setting the mains voltage**

2.2.2 Fuses

The fuse holder of the IC Eluent Selector is fitted with one of the following two fuses as standard: 1.0 A/slow-blow for 100...120 V or 0.5 A/slow-blow for 220...240 V.



Never use different fuses in the instrument from those mentioned above as otherwise there is a fire risk!

Checking and exchanging fuses is described in Section 2.2.1.

2.2.3 Mains cable and mains connection

Mains

The instrument is supplied with one of the following mains cables:

- 6.2122.020 with SEV 12 plug (Switzerland, ...)
- 6.2122.040 with CEE(7), VII plug (Germany, ...)
- 6.2133.070 with NEMA 5-15 plug (USA, ...)

which has three wires and is fitted with a plug with a grounding pin. If a different plug has to be used then the yellow/green wire (IEC standard) must be connected to the grounding pin (protection class I).



Any interruption to the grounding inside or outside the instrument can represent a hazard!

Mains connection

Insert the mains cable into **mains connector 9** of the 816 IC Eluent Selector.

2.2.4 Switching the instrument on/off

The 816 IC Eluent Selector is switched on and off with **mains switch 8**.

2.3 Electrical connection

2.3.1 Operation in a modular IC system



Always switch off the 762 IC Interface and 816 IC Eluent Selector before connecting the two instruments.

Use a 6.2134.080 or 6.2134.090 Cable to connect RS 232 interface **14** of the 816 IC Eluent Selector with a free RS 232 interface on the 762 IC Interface according to the following diagram.




Figure 4: Connecting the 816 IC Eluent Selector to 762 IC Interface

2.3.2 Operation with a 761 Compact IC



Always switch off the 816 IC Eluent Selector and PC before connecting the two instruments.



Figure 5: Connecting the 816 IC Eluent Selector to a PC

The 816 IC Eluent Selector can also be operated with the 761 Compact IC without a 762 IC Interface. The RS 232 interface of the two instruments must then be connected directly to the RS 232 interfaces of the PC as shown in Figure 5.

2.4 Valve connections

The following diagram shows the flow path in the 816 IC Eluent Selector according to the selected eluent.




Figure 6: Flow diagram of the IC Eluent Selector

The 816 IC Eluent Selector is installed between the eluent storage containers and the high-pressure pump of the IC system. The following connections must be made:

1 Attaching the eluent inlet tubing

- Screw the PEEK fittings of the three pieces of **6.1834.020 Aspiration tubing** supplied loosely by hand into connections "A" and "B" of **Valve I 5** and "C" of **Valve II 7** of the IC Eluent Selector. Tighten the screw fittings by using the **6.2621.050 ¼" Wrench** (also included in the standard equipment) to rotate them a further ¼ turn.
- Make sure that the other ends of the aspiration tubing are immersed in the eluents.

2 Attaching the outlet to the IC pump

- Screw the **6.1834.030 Tubing connection 816-709** by hand into the "PUMP" connection of **Mixing chamber 6**. Tighten up the connection by using the **6.2621.050 ¼" Wrench** to rotate it a further ¼ turn.
- Slide the tubing end onto the capillary at the aspiration connection of the IC pump.



Make sure that the tubing between the outlet of the IC Eluent Selector and the IC pump is kept as short as possible by shortening the 6.1834.030 Tubing connection 816-709 accordingly. In this way you minimize the dead volume and therefore the reaction time after switching to a different eluent.

2.5 Software installation

The PC program «**IC Net 2.1**» is required for controlling the 816 IC Eluent Selector from a PC. This program runs under the operating systems Windows 95, Windows 98, Windows NT and Windows 2000 and is installed as described in Section 1.4.2 of the **8.110.8223 Instructions for Use** for the Metrodata «**IC Net 2.1**» software.

The installation of the 816 IC Eluent Selector is described in Section 6.12 of the **8.110.8223 Instructions for Use** for the Metrodata «**IC Net 2.1**» software.

Make sure that the RS 232 settings of the 816 IC Eluent Selector and the interface of the connected PC/IC interface are identical. The settings for the 816 IC Eluent Selector are given in Section 4.1.2.


3 Operation

3.1 Manual operation

Switching the instrument on/off

- The 816 IC Eluent Selector is switched on and off with mains switch **8** (see Section 2.2.4).


Switching to manual operation



LED "EXT." **4** shows whether the instrument is controlled **externally (LED lights up)** or **manually (LED is out)**.

- Press the "EXT." **3** key to switch to manual operation of the 816 IC Eluent Selector.

Select eluent




LEDs **2** show which eluent is currently selected.

- The eluent is selected by pressing key "A", "B" or "C" **1**.



The settings are retained when the instrument is switched off.



3.2 Operation via «IC Net»

Switching the instrument on/off

- The 816 IC Eluent Selector is switched on and off with mains switch **8** (see Section 2.2.4).

Switching to external operation



LED "EXT." **4** shows whether the instrument is controlled **externally (LED lights up)** or **manually (LED is out)**.

- Press the "EXT." **3** key to switch to the external mode for remote control of the 816 IC Eluent Selector by «IC Net».



Only the most important functions and settings for operating the 816 IC Eluent Selector are described in this section. Further information can be found in the Instructions for Use for «IC Net» and in the online help for the program.

3.2.1 Symbol for 816 IC Eluent Selector




2.816.0010 816 IC Eluent Selector

3.2.2 Instrument parameters in the "816 IC Eluent Selector" window

The **816 IC Eluent Selector** window for parameter settings is opened with a double-click with the left-hand mouse key on the 816 symbol in the system window or by clicking on this symbol with the right-hand mouse key and then selecting the menu option **Open**. It consists of four register cards **Manual**, **Aliases**, **Program** and **Links**.

Manual

The **Manual** register card of the **816 IC Eluent Selector** window is used for the manual selection of the eluents.



Current eluent


Selection of **Eluent A**, **Eluent B** or **Eluent C**.

<Set>

The valves are switched so that the required eluent is delivered.

Aliases

In the **Aliases** register card of the **816 IC Eluent Selector** window you can freely define names for connections "A", "B" **5** and "C" **7** of the 816 IC Eluent Selector.




- Flow A** You can rename **Eluent A**.
- Flow B** You can rename **Eluent B**.
- Flow C** You can rename **Eluent C**.

The names of the eluents entered here replace the default names **Eluent A – C** in all the register cards of the **816 IC Eluent Selector** window.

Program

On the register card **Program** in the **816 IC Eluent Selector** window you can enter a user-specific time program. Depending on the settings made in the **Start mode** window (for details see Section 4.4.3 **8.110.8223 Instructions for Use** for the Metrodata «**IC Net 2.1**» software) the program will be started automatically at the start of the determination (**Start with determination**) or when the sample is injected (**Start with inject**, in this case the **Enabled** checkbox must be activated).



Time (Column 1)	Time for command to be executed. Range: 0.0 ... 999.9 min If no time is entered then the command will be carried out together with the last command for which a time has been entered.
Command (Column 2)	Switches the IC Eluent Selector to the required eluents. The alias names entered above will be shown for selection.
ENABLED	Program for program start is activated (if a program is not activated it will not start).
<Add>	Adds new program command.
<Delete>	Deletes selected program command.
<Verify>	Checks time program (any fault will produce an error message).

Links

The register card **Links** of the **816 IC Eluent Selector** window is used to select the serial RS 232 interface (for details see **8.110.8223 Instructions for Use** for the Metrodata «**IC Net 2.1**» software).

4 Annex

4.1 Technical data

4.1.1 Parts and controls

Keypad

Chemical-resistant touch-sensitive polyester keypad with function keys


Indicators

LEDs for indicating valve status

4.1.2 RS 232 Interface

Socket

D-sub socket 25-pin (male)



Standards

EIA RS232C, CCITT V.24, ISO 2110, DIN 66020

Parameters

Baud rate: 9600
 Handshake: no handshake
 Parity: none
 Stop bits: 1
 Data bits: 8

4.1.3 Mixer speed

Switch

Dip switch, 8-position



Function

Not required for the 816. We recommend that the mixing chamber is kept switched off (all switches down).

The left-hand six switches control the stirring speed in mixing chamber **6**, the right-hand two have no function.

Only one of the six left-hand switches can be in the ON position, the stirring speed increases from left to right.

4.1.4 Mains connection

<i>Voltage</i>	115 V: 100...120 V \pm 10 % 230 V: 220...240 V \pm 10 % Switched with mains voltage selector in the fuse holder; see Section 2.2.1
<i>Frequency</i>	50...60 Hz
<i>Power consumption</i>	40 VA
<i>Fuses</i>	5 mm dia., 20 mm long 100...120 V: 1.0 A (slow blow) 220...240 V: 0.5 A (slow blow) replace only with the same type

4.1.5 Safety specifications

<i>Construction and testing</i>	According to EN/IEC 61010-1 / UL 3101-1, protection class I
<i>Safety information</i>	The Instructions for Use contain safety information that must be observed by the user in order to ensure the safe operation of the instrument.

4.1.6 Electromagnetic compatibility (EMC)

<i>Emission</i>	Standards fulfilled: - EN/IEC 61326-1 - EN 55022
<i>Immunity</i>	Standards fulfilled: - EN/IEC 61326-1 - EN/IEC 61000-4-2 - EN/IEC 61000-4-3 - EN/IEC 61000-4-4 - EN/IEC 61000-4-5 - EN/IEC 61000-4-6 - EN/IEC 61000-4-11 - NAMUR

4.1.7 Ambient temperature

<i>Nominal working range</i>	+5...+45 °C (at 20...80% relative humidity)
<i>Transport</i>	-40...+70 °C




4.1.8 Housing


<i>Cover material</i>	Rigid polyurethane foam(PUR) with flammability proofing for UL94VO, CFC-free
<i>Base material</i>	Lacquered steel
<i>Width</i>	260 mm
<i>Height</i>	129 mm
<i>Depth</i>	366 mm
<i>Weight</i>	4.5 kg

4.2 Standard equipment






*We reserve the right to make alterations !
All dimensions given in mm.*

The 2.816.0010 IC Eluent Selector includes the following accessories:

No.	Order no.	Description												
1	1.816.0010	IC Eluent Selector												
1	6.21220X0	Mains cable to customer's requirements: <table border="0" style="width: 100%;"> <tr> <td style="text-align: left;"><u>Cable socket</u></td> <td style="text-align: left;"><u>Cable plug</u></td> <td></td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type SEV 12 (CH...)</td> <td>6.2122.020</td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type CEE (7), VII (D...)</td> <td>6.2122.040</td> </tr> <tr> <td>Type CEE (22), V</td> <td>Type NEMA 5-15 (USA)</td> <td>6.2122.070</td> </tr> </table>	<u>Cable socket</u>	<u>Cable plug</u>		Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020	Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040	Type CEE (22), V	Type NEMA 5-15 (USA)	6.2122.070
<u>Cable socket</u>	<u>Cable plug</u>													
Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020												
Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040												
Type CEE (22), V	Type NEMA 5-15 (USA)	6.2122.070												
3	6.1834.020	Aspiration tubing for 816 made of PTFE, with connection piece for 6.2821.090 Suction filter  Length = 1.5 m, i.d. = 1.5 mm, o.d. = 2.5mm For the connection 816 – eluent container												
1	6.1834.030	Tubing connection 816 – 709 made of PTFE  Length = 0.5 m, i.d. = 1.5 mm, o.d. = 2.5mm For the connection 816 – high-pressure pump												
1	6.2621.050	Wrench 1/4" 												

1	6.2821.090	Aspiration filter pore size 20 µm For 6.1834.020 Aspiration tubing. Set of 5 pieces.	
1	8.816.1003	Instructions for Use (English) for 816 IC Eluent Selector	

4.3 Optional accessories

No.	Order no.	Description	
1	6.2134.080	RS232 connection cable 816 - 762, length = 2 m	
1	6.2125.110	RS232 connection cable 816 - PC, length = 3 m	
1	U.600.0013	Fuse 0.5 AT slow blow	
1	U.600.0016	Fuse 1.0 AT slow blow	

4.4 Validation / GLP

GLP (Good Laboratory Practice) requires, among other things, that the precision and correctness of analytical instruments is checked at regular intervals by using SOPs (Standard Operating Procedures, **SOP**). An example of such a standard operating procedure is available from Metrohm under the title «**Application Bulletin No. 277 – Validation of Metrohm Ion Chromatographs**». This SOP can be adapted for your ion chromatography system and used for its validation.

The 816 IC Eluent Selector must be included as a part of the whole ion chromatography system, whose most important components include the pumps, separating columns, detector and evaluation system, in the all-embracing validation of the whole system.

Please contact your local Metrohm agency in order to receive support in validating your 816 IC Eluent Selector. It can also provide you with validation documentation which will help you to carry out your installation qualification (IQ) and operational qualification (OQ).

Further information about QA, GLP and validation can also be found in the brochure «**Quality Management with Metrohm**» which is also obtainable from your local Metrohm agency.

Checking the electronic and mechanical assemblies of Metrohm instruments can and should be undertaken within the framework of regular servicing by Metrohm technicians.

4.5 Warranty and conformity

4.5.1 Warranty

The warranty on our products is limited to defects that are traceable to material, construction or manufacturing error which occur within 12 months from the day of delivery. In this case the defects will be rectified in our workshops free of charge. Transport costs are to be paid by the customer.

For day and night operation the warranty is limited to 6 months.

Glass breakage in the case of electrodes or other parts is not covered by the warranty. Checks which are not a result of material or manufacturing faults are also charged during the warranty period. For parts from outside manufacturers, insofar as these constitute an appreciable part of our instrument, the warranty stipulations of the manufacturer in question apply.

With the regard to the guarantee of accuracy the technical specifications in the instruction manual are authoritative.

Concerning defects in materials, construction or design as well as the absence of guaranteed features the purchaser has no rights or claims except those mentioned above.

If damage of the packaging is evident on receipt of a consignment or if the goods show signs of transport damage after unpacking, the carrier must be informed immediately and a written damage report demanded. Lack of an official damage report releases Metrohm from any liability to pay compensation.

If any instruments and parts have to be returned then the original packaging should be used if at all possible. This applies above all to instruments, electrodes, buret cylinders and PTFE pistons. Before embedment in wood shavings or similar material the parts must be packed in a dustproof package (for instruments the use of a plastic bag is essential). If open assemblies are included that are sensitive to electromagnetic voltages (e.g. data interfaces, etc.) then these must be returned in the associated original protective packaging (e.g. conductive protective bag). (Exception: assemblies with a built-in voltage source belong in non-conductive protective packaging).

For damage which arises as a result of non-compliance with these instructions no warranty responsibility whatsoever will be accepted by Metrohm.

4.5.2 EU Declaration of Conformity

 <p>EU Declaration of Conformity</p>						
<p>The company Metrohm AG, Herisau, Switzerland, certifies herewith, that the following instrument:</p> <p style="text-align: center;">816 IC Eluent Selector</p> <p>meets the CE mark requirements of EU Directives 89/336/EEC and 73/23/EEC.</p>						
<p>Source of specifications:</p> <table style="width: 100%; border: none;"> <tr> <td style="padding-left: 20px;">EN 61326-1</td> <td>Electrical equipment for measurement, control and laboratory use – EMC requirements</td> </tr> <tr> <td style="padding-left: 20px;">EN 61010-1</td> <td>Safety requirements for electrical equipment for measurement, control and laboratory use</td> </tr> </table>	EN 61326-1	Electrical equipment for measurement, control and laboratory use – EMC requirements	EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use		
EN 61326-1	Electrical equipment for measurement, control and laboratory use – EMC requirements					
EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use					
<p>Description of apparatus:</p> <p>External controllable device with two 3-way valves and mixing chamber.</p>						
<p>Herisau, March 9th, 2001</p> <table style="width: 100%; border: none; margin-top: 20px;"> <tr> <td style="text-align: center; width: 50%;">  </td> <td style="text-align: center; width: 50%;">  </td> </tr> <tr> <td style="text-align: center;">Dr. J. Frank</td> <td style="text-align: center;">Ch. Buchmann</td> </tr> <tr> <td style="text-align: center;">Development Manager</td> <td style="text-align: center;">Production and Quality Assurance Manager</td> </tr> </table>			Dr. J. Frank	Ch. Buchmann	Development Manager	Production and Quality Assurance Manager
						
Dr. J. Frank	Ch. Buchmann					
Development Manager	Production and Quality Assurance Manager					

4.5.3 Certificate of conformity and system validation

Certificate of Conformity and System Validation

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

Name of commodity:	816 IC Eluent Selector
System software:	Stored in ROMs
Name of manufacturer:	Metrohm Ltd., Herisau, Switzerland

This Metrohm instrument has been built and has undergone final type testing according to the standards:

Electromagnetic compatibility: Emission
IEC 61326-1, EN 55022

Electromagnetic compatibility: Immunity
IEC 61326-1, IEC 61000-4-2, IEC 61000-4-3, IEC 61000-4-4, IEC 61000-4-5,
IEC 61000-4-6, IEC 61000-4-11, NAMUR

Safety specifications
IEC 61010-1, UL3101-1


It has also been certified by the Swiss Electrotechnical Association (SEV), which is member of the International Certification Body (CB/IEC).

The technical specifications are documented in the instruction manual.

The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance.

Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.

Herisau, March 9th, 2001



Dr. J. Frank
Development Manager



Ch. Buchmann
Production and
Quality Assurance Manager

4.6 Index

A

A	
Connection 5	2
Key 1	2, 13
LED 2	2
Accessories	20
Aliases	15
Ambient temperature	18
An/Cat	1
Annex	17
Aspiration filter	
Order number	20
Aspiration tubing	12
Order number	19
Attention	5

B

B	
Connection 5	2
Key 1	2, 13
LED 2	2

C

C	
Connection 7	2
Key 1	2, 13
LED 2	2
Cable	
6.2125.110	20
6.2134.080	20
816-762	20
816-PC	20
CE sign	23
Certificate of conformity and system validation	24
Checks	7
Conformity	22
Connection	
762 IC Interface	10
A 5	2
B 5	2
C 7	2
Eluent	11
Modular system	10
To PC	10
Valve	11
Contents	1
Current eluent	14

D

Danger	5
Dead volume	12
Declaration of Conformity	23
Description	
Eluent Selector	1
of the instrument	1

E

Electrical connection	10
Electrical safety	6
Electromagnetic compatibility	18

Electrostatic charges	6
Eluent Selector	
Description	1
EU Declaration of Conformity	23
EXT.	
Key 3	2
LED 4	2

F

Flow A	15
Flow B	15
Flow C	15
Flow diagram	
Illustration	11
Front panel	2
Fuse	8, 9
Order number	20
Fuse holder 12	
Illustration	3, 9
Fuses	18

G

General safety rules	6
GLP	21
Ground connection 11	3
Grounding	6, 9
Guarantee	22

H

Handling solutions	6
Housing	19

I

IC Eluent Selector	
Description	1
Order number	19
IC Net	
816 operation	13
Illustrations	
List	11
Information	5
Installation	7
Instructions for Use	
Order number	4, 20
Organization	4
Instrument on/off	13
Interface	
RS 232 14	3
Introduction	1
IQ	21
ISO 9001	24

K

Key	
A 1	2, 13
B 1	2, 13
C 1	2, 13
EXT. 3	2
Keys	2

L

Leaks	6
LED	2
A 2	2
B 2	2
C 2	2
EXT. 4	2
Links	16
List of illustrations	11
Location	7

M

Mains cable	9
Order number	19
Mains connection	6, 7
Technical data	18
Mains connection 9	3
Illustration	9
Mains connector 9	9
Mains frequency	18
Mains switch 8	3
Illustration	9
Switching the instrument on/off	9
Mains voltage	18
Manual	14
Mixer speed 13	3
Technical data	17
Mixing chamber 6	2

N

Notation	5
----------	---

O

Operation	13
In modular system	10
Manual	13
Via IC Net	13
With 761 Compact IC	10
Optional accessories	20
OQ	21
Order numbers	19
Organization	4

P

Packaging	7
part numbers	19
Parts and controls	5
Pictograms	5
power consumption	18
Program	15
Protection class	9

Q

QA	21
Quality assurance	21

R

Reaction time.....	12
Rear panel	3
Room temperature	7
RS 232	
Connection	3
Interface 14	3
Technical data.....	17
Settings.....	12
RS232 connection cable	
816-762.....	20
816-PC.....	20
Order number.....	20

S

Safety information.....	6
Safety rules	6

Safety specifications.....	18
Serial number 10	3
Setting the mains voltage.....	7, 8
Setup.....	7
Software installation.....	12
SOP.....	21
Standard equipment.....	19
Standard operating procedures..	21
Step gradients	1
Switching on/off	13
Switching the instrument on/off9,	13

T

Table of contents	I
Technical data	17
Time program	15
Transport.....	7

Transport damage	22
Tubing connection 816-709.....	12
Order number	19

V

Validation.....	21
Valve	
Connection.....	11
I 5	2
II 7	2

W

Warning	5
Warranty	22
Wrench ¼"	
6.2621.050	12
Order number	19