

819 IC Detector

820 IC Separation Center

5.819.0010 Program



819 IC Detector 820 IC Separation Center

5.819.0010 Program

Instructions for Use

Teachware
Metrohm AG
Oberdorfstrasse 68
CH-9101 Herisau
teachware@metrohm.com

2nd Edition 2004

These instructions are protected by copyright. All rights reserved.

Although all the information given in these instructions has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please inform the author at the address given above.

Table of contents

1	Introduction	1
1.1	Instrument description	1
1.1.1	819 IC Detector	1
1.1.2	820 IC Separation Center	2
1.1.3	Schematic arrangement of the IC systems	3
1.2	Parts and controls	4
1.2.1	819 IC Detector	4
1.2.2	820 IC Separation Center	6
1.3	Information on the Instructions for Use	10
1.3.1	Organization	10
1.3.2	Notation and pictograms	11
1.4	Safety notes	12
1.4.1	General precautionary rules	12
2	Installation	13
2.1	Setting up the instrument	13
2.1.1	Packaging	13
2.1.2	Check	13
2.1.3	Location	13
2.1.4	Arrangement of the instruments	13
2.2	Installation of accessories	15
2.2.1	Insert detector block	15
2.2.2	Connection of column heating	18
2.2.3	Connection of syringe and aspirating tubing	19
2.2.4	Connection of the drain tube	19
2.2.5	Connection of the 6.5324.000 Bottle rack (option)	19
2.3	Electrical connection	20
2.3.1	Connection to 830 IC Interface	20
2.3.2	Connection to PC	23
2.4	Mains connection	25
2.4.1	819 IC Detector	25
2.4.2	820 IC Separation Center	27
2.4.3	Column heating	27
2.5	Software-installation	27
2.6	Capillary connections	28
2.7	Connection of 818 IC Pump	29
2.7.1	Electrical connection	29
2.7.2	Pulsation dampener	29
2.7.3	Filter unit PEEK	30
2.7.4	Connection to injection valve	31
2.8	Precolumns	33
2.9	Separating columns and suppressor module	34
2.9.1	General information on separating columns	34
2.9.2	Installing a separating column in the column heating	35
2.9.3	Selection of the sample loop	37
2.9.4	General information on suppressor module	37
2.9.5	One-channel system without suppressor module	39
2.9.6	Two-channel system without suppressor module	41
2.9.7	One-channel system with suppressor module	42
2.9.8	Two-channel system with suppressor module	47
2.9.9	Leak testing and conditioning	48

3	Operation	50
3.1	819 IC Detector	50
3.1.1	819 IC Detector icon	50
3.1.2	Settings in the "819 IC Detector" window.....	50
3.2	820 IC Separation Center	60
3.2.1	820 IC Separation Center icon.....	60
3.2.2	"820 IC Separation Center" window	60
4	Notes – Maintenance – Faults	66
4.1	Practical notes on ion chromatography	66
4.1.1	Separating columns	66
4.1.2	Pumps	67
4.1.3	Eluents.....	68
4.1.4	Suppressor module.....	69
4.1.5	Connections	69
4.2	Maintenance and servicing	69
4.2.1	General information.....	69
4.2.2	Passivation	70
4.2.3	Recycling.....	70
4.2.4	Shutdown	70
4.2.5	Changing separating columns	71
4.2.6	Regeneration of suppressor	72
4.2.7	Cleaning the suppressor.....	73
4.2.8	Replacement of suppressor module	75
4.3	Faults and malfunctions	77
4.3.1	Malfunctions and their rectification	77
5	Interfaces	79
5.1	RS 232 interface.....	79
5.1.1	Data transmission protocol.....	79
5.1.2	Pin assignment.....	80
5.2	Remote interfaces.....	81
5.2.1	"Remote" interface	81
5.2.2	"IC Separation Center" interface.....	83
5.3	Analog output.....	85
5.4	External power supply for 820 IC Separation Center	85
5.5	Valve interfaces.....	86
6	Appendix	87
6.1	Technical data	87
6.1.1	819 IC Detector	87
6.1.2	820 IC Separation Center.....	90
6.2	Scope of delivery	92
6.2.1	819 IC Detector	92
6.2.2	820 IC Separation Center.....	94
6.3	Optional accessories	103
6.3.1	Accessories for 820 IC Separation Center	103
6.3.2	Column heating.....	104
6.3.3	Cable	105
6.3.4	Literature.....	105
6.4	Validation / GLP	106

6.5	Warranty and Conformity	107
6.5.1	Warranty.....	107
6.5.2	Declaration of Conformity	108
6.5.3	Declaration of Conformity	109
6.5.4	Quality Management Principles	110
6.6	Index	111

List of figures

Fig. 1:	Block diagram of the ion chromatography systems	3
Fig. 2:	Front of the 819 IC Detector	4
Fig. 3:	Rear of the 819 IC Detector	5
Fig. 4:	Front of the 820 IC Separation Center.....	6
Fig. 5:	Rear of the 820 IC Separation Center	7
Fig. 6:	Column heating closed	9
Fig. 7:	Insert detector block in the Separation Center	15
Fig. 8:	Connection Detector 819 – 2.820.0X10/2.820.0X30	16
Fig. 9:	Connection 2 x 819 Detector – 2.820.0X20	17
Fig. 10:	Install column heating in Separation Center	18
Fig. 11:	Connection 820 – column heating	19
Fig. 12:	Connection of 820 via 819 to 830.....	20
Fig. 13:	Connection of 819 and 820 to 830	22
Fig. 14:	Connection of 819 and 820 to a PC	24
Fig. 15:	Setting the mains voltage	26
Fig. 16:	Connect capillaries	28
Fig. 17:	Connection of 818 IC Pump to 830 IC Interface.....	29
Fig. 18:	Filter unit PEEK (6.2821.120).....	30
Fig. 19:	Connection to injection valve	31
Fig. 20:	Insert column in column heating	35
Fig. 21:	Interior of the 2.820.0310 Separation Center	39
Fig. 22:	Interior of the 2.820.0320 Separation Center	41
Fig. 23:	Interior of 2.820.0330 Separation Center	42
Fig. 24:	Suppressor module connections	44
Fig. 25:	Interior of 2.820.0320 Separation Center with suppressor module.....	47
Fig. 26:	Operating principle of Auto-zero	52
Fig. 27:	Diagram of measuring and full-scale range	54
Fig. 28:	Flow path of the injection valve	62
Fig. 29:	Assembling the suppressor	74

1 Introduction


1.1 Instrument description

1.1.1 819 IC Detector

The **819 IC Detector** is a conductivity detector especially designed for ion chromatography with an extensive operating range and high sensitivity for the recording of chromatograms with and without chemical suppression. The associated thermostatable detector block is normally installed in the 819 IC Separation Center, but can also be used as a separate detector. The two following versions are available:

- **2.819.0010** IC Detector with metal detector block
- **2.819.0110** IC Detector with metal-free detector block

The **819 IC Detector** can be fully remotely controlled via the Metrohm «IC Net» software. All functions can be quickly and easily accessed by the software and can be clearly shown on the screen. Apart from setting the measuring parameters, it is possible to draw up any time programs required. These can be used to control a large number of instrument functions and freely programmable events.



The 819 IC Detector is connected via the 830 IC Interface or directly to a PC. The control of the instrument and the recording and evaluation of the chromatograms take place under «IC Net». It is also possible to control external devices or start functions in the IC system from these devices via a “Remote” interface using programmable signals.

1.1.2 820 IC Separation Center

The **820 IC Separation Center** is a thermally and electronically insulated wet-chemistry component that accommodates injectors, columns, column heating, detector blocks, suppressor module, pulsation dampeners and various optional sample preparation modules. The 820 IC Separation Center is also remotely controlled via the Metrohm «IC Net» software; it can either be connected to the 819 IC Detector and controlled through it, or connected directly to the 830 IC Interface and controlled directly from «IC Net».



The following versions are available:

- **2.820.0210** IC Separation Center with 1 injector for a one-channel system with column heating, metal-free
- **2.820.0220** IC Separation Center with 2 injectors for a two-channel system with column heating, metal-free
- **2.820.0230** IC Separation Center with 1 injector and 1 Metrohm Suppressor Module MSM for a one-channel system with column heating, metal-free
- **2.820.0310** IC Separation Center with 1 injector for a one-channel system, metal-free
- **2.820.0320** IC Separation Center with 2 injectors for a two-channel system, metal-free
- **2.820.0330** IC Separation Center with 1 injector and 1 Metrohm Suppressor Module MSM for a one-channel system; metal-free

1.1.3 Schematic arrangement of the IC systems

The 819 IC Detector and 820 IC Separation Center are the main components of a modular **ion chromatography system** that can be expanded to meet the wishes of the individual user (see Fig. 1). The minimum configuration of the one-channel system also includes a 818 IC Pump, a separating column, a 830 IC Interface and an PC. The two-channel system requires at least a second 819 IC Detector and a second 818 IC Pump. Both systems can be extended virtually without any limits by further modular Metrohm IC devices such as sample changers for automation, sample preparation devices, components for post-column derivation and other detectors. Further, practically all HPLC peripherals and parts available on the market such as precolumns, additional separating columns, additional detectors and other injection systems can be seamlessly integrated in the system.

However, the individual IC units can also be freely combined with common HPLC instruments. This offers the possibility of expanding your system to a stand-alone ion chromatograph.




Fig. 1: Block diagram of the ion chromatography systems

C	Separating Column	M	Suppressor module	818	IC Pump
D	Detector	PC	PC	819	IC Detector
E	Eluent	PR	Printer	820	IC Separation Center
I	Injector	S	Probe		
IF	Interface	W	Waste		

1.2 Parts and controls



In this Section you will find the numbers and designations of the parts and controls of the 819 IC Detector and 820 IC Separation Center. The numbering applies throughout the instructions for use, i.e. bold numbers in the text (e.g. **3**) refer to the parts and controls illustrated here.

1.2.1 819 IC Detector




Fig. 2: Front of the 819 IC Detector

1 Mains pilot lamp

Lights up when instrument is switched on




Fig. 3: Rear of the 819 IC Detector

2 Mains switch

For switching instrument on/off:

I = ON 0 = OFF

3 Mains connection plug

Mains connection see *Section 2.4.1*

4 Fuse holder

Changing the fuses, see *Section 2.4.1*

5 Analog output

0...1 V

6 Connection for detector block

7 Model plate with serial number

8 Remote interface

remote I/O lines for connection of external devices

9 RS 232 interface for connection of IC Interface 830 or PC

10 Connection for 820 IC Separation Center

1.2.2 820 IC Separation Center




Fig. 4: Front of the 820 IC Separation Center

11 Mains pilot lamp

Lights up when instrument is switched on

12 Thermostat lamp

Lights up when column heating is connected, see *Section 2.2.2* and *Section 2.4.3*

13 Connection for syringe 6.2816.020

14 Feedthrough

for sample, eluent, suppressor reagent, ...




Fig. 5: Rear of the 820 IC Separation Center

14 Feedthrough

for sample, eluent, suppressor reagent,...

15 Knurled screw

for fastening rear panel **21**, resp. **28**.

16 Opening for detector cable B

opening for connecting cable:
detector block B – 819

17 Opening for outlet capillary B

discharge of the eluent of column B to waste

18 Rear panel opening

(closed with plastic stopper) for additional supply and discharge lines to and from the inner compartment

19 Opening for outlet capillary A

discharge of the eluent of column A to waste

20 Opening for detector cable A and for connecting cable to the column heating

opening for connecting cable:
detector block B – 819
column heating – 820 (6.2108.120)

21 Detachable rear panel

access to top part of the inner compartment

22 RS 232 interface for column heating

to connect to 830 IC Interface or PC with cable 6.2134.040

23 Connection for column heating

with cable 6.2108.120

24 Connection for external power supply of column heating

with power supply unit 6.2152010

25 Connection for external power supply of the Separation Center

connection of power supply unit
6.2152.000

(5 V, 0.5 A / 24 V, 2 A) in operation without 819 IC Detector

26 Terminal block for valve B

Ground, Fill, Inject:

inputs for control of the valve

RUN Pos. / COM Fill:

output signal on switching of the valve to position "FILL"

RUN Integr. / COM Start:

output signal on switching of the valve to position "INJECT"

27 Connection for 819 IC Detector**28 Detachable rear panel**

access to bottom part of the inner compartment

29 Model plate with serial number**30 Connection for 819 IC Detector A****31 Terminal block for Valve A**

Ground, Fill, Inject:

inputs for control of the valve

RUN Pos. / COM Fill:

output signal on switching of the valve to position "FILL"

RUN Integr. / COM Start:

Output signal on switching of the valve to position "INJECT"

32 Opening for inlet capillary B

supply of the eluent for column B

33 Connection for drain tube

for discharge of spilled liquid from the inner compartment

34 Opening for inlet capillary A

supply of the eluent for column A

Column heating




Fig. 6: Column heating closed

35 Feedthrough
for the capillary from the separating column to detector

36 Knurled screw
for fastening the heater cover **37**

37 Cover of column heating

38 Connection for Separation Center
with cable 6.2108.120

39 Magnetic holder
for fastening the column heater in the inner compartment of the Separation Center

40 Model plate with serial number

41 Feedthrough
for column connection capillary **50**

1.3 Information on the Instructions for Use



Please read through these Instructions for Use carefully before you put the 819 IC Detector and 820 IC Separation Center into operation. The Instructions for Use contain information and warnings which must be heeded the user to assure safe operation of the instruments.

1.3.1 Organization

These **8.819.1013 Instructions for Use** for the 819 IC Detector and 820 IC Separation Center provide a comprehensive overview of the installation, startup procedure, operation, fault rectification and technical specifications of these instruments. The Instructions for Use are organized as follows:

Section 1 Introduction

General description of instruments, parts and controls and safety notes

Section 2 Installation

Installation of 819 IC Detector / 820 IC Separation Center, electrical connections, tubing connections, mains connection

Section 3 Operation

Operation via «IC Net»

Section 4 Notes - Maintenance - Faults

Practical notes, maintenance, fault rectification

Section 5 Interfaces

Description of the interfaces

Section 6 Appendix





Technical data, standard equipment, options, validation, warranty, declaration of conformity, index

To find the information you require about the instrument please use either the **Table of contents** or the **Index** at the back.

As an addition to the Instructions for Use the **Metrohm Monographs** listed in *Section 6.3.4* can be requested on the Internet under <http://www.metrohm.com> in the literature section, or free of charge from your local Metrohm agency. Detailed information about separating columns available from Metrohm can be found in the **Metrohm IC Column Catalog** or on the Internet under <http://www.metrohm.com> in the ion chromatography range of products. Information about special IC applications is given in the relevant "**Application Bulletins**" or "**Application Notes**"; these can be found on the Internet under <http://www.metrohm.com> in the applications section or can be requested free of charge from your local Metrohm agency.

1.3.2 Notation and pictograms

The following notations and pictograms (symbols) are used in these Instructions for Use:

Range	Menu item, parameter or entry value in «IC Net» program
SYSTEM STATE	Program window in «IC Net» program
<OK>	Button in «IC Net» program
7	Part or control of 819 / 820
<u>3</u>	Part or control of 818
<u>21</u>	Part or control of 833
	Danger/Warning This symbol indicates a possible risk of death or injury to the user and possible damage to the instrument or its components by electricity.
	Danger/Warning This symbol indicates a possible risk of death or injury to the user and possible damage to the instrument or its components.
	Attention This symbol indicates important information that you should read before continuing.
	Information This symbol indicates additional information and tips which may be of particular use to you.

1.4 Safety notes

While electrical safety in the handling of the 819 IC Detector and 820 IC Separation Center is assured in the context of the specifications EN/IEC 61010-1 / UL 3101-1 (protection class I), the following points should be noted:

- **Mains connection**



Setting of the **mains voltage**, checking the **mains fuse** and the **mains connection** must be effected in accordance with the instructions in Section 2.4.

- **Opening the 819 IC Detector**

Inside the instrument there are no parts which must be set or adjusted by the user.



If the 819 IC Detector is connected to the mains then the instrument must not be opened, nor should any components be removed from it as otherwise you run the risk of coming into contact with current-carrying components. For this reason you should always separate the instrument from sources of electricity before opening it and also ensure that the **mains plug is removed from mains connector 3!**

- **Opening the 820 IC Separation Center**

The instrument contains no components that can be set or adjusted by the user.



Disconnect **all connecting cables on the rear of the 820 IC Separation Center** before you remove the middle housing panel with connectors.

- **Protection against static charges**



Electronic components are sensitive to static charging and can be destroyed by discharges. Before you touch any of the components inside the 819 IC Detector or 820 IC Separation Center, you should earth yourself and any tools you are using by touching an earthed object (e.g. housing of the instrument or a radiator) to eliminate any static charges which exist.

1.4.1 General precautionary rules

- **Handling of solvents**



Check all lines of the IC system periodically for possible leaks. Follow the relevant instructions regarding the handling of flammable and/or toxic solvents and their disposal.

2 Installation

2.1 Setting up the instrument

2.1.1 Packaging

819 IC Detector and 820 IC Separation Center are supplied together with the separately packed accessories in special packagings containing shock-absorbing foam linings designed to provide excellent protection. The actual instrument is packed in an evacuated polyethylene bag to prevent the ingress of dust. Please store all these special packagings as only they assure transport of the instruments free from damage.

2.1.2 Check

After receipt, immediately check whether the shipment is complete and has arrived without damage (compare with delivery note and list of accessories in *Section 6.2*). In the case of transport damage, see instructions in *Section 6.5.1 "Warranty"*.

2.1.3 Location

Position the instrument the laboratory at a location convenient for operation, free from vibrations and protected against a corrosive atmosphere and contamination by chemicals.



To avoid disturbing temperature influences on the insulated column compartment, the pump and eluent reservoir must be protected against direct sunlight.

2.1.4 Arrangement of the instruments

Modular IC instruments can be stacked on top of each other in any sequence.



In a modular system any instruments containing liquid-transporting assemblies should always be located as low down in the stack as possible, so that any leaks which may occur in tubing or connections will cause no damage to the other instruments by escaping liquids such as acids.

In one-channel operation, the 818 IC Pump, 820 IC Separation Center and 819 IC Detector are best stacked on top of one another in this order.

In two-channel operation (2.820.0X20 IC Separation Center), the optimum arrangement (1, 2 or 3 towers) depends on the laboratory space available. However, the 818 IC Pumps should be set up at the very bottom and the 819 IC Detectors at the very top.



To ensure that the arrangement of pumps and detectors for the two channels A and B is clearly apparent in two-channel operation, it is advantageous to mark the instruments. The 6.2248.000 Magnetic plate is enclosed with the 819 IC Detector for this purpose. It can be cut to the desired size, labeled (e.g. with "A" or "B") and affixed to the appropriate instrument.

2.2 Installation of accessories

2.2.1 Insert detector block

In order to operate an 820 IC Separation Center with an 819 IC Detector the detector block contained in the standard equipment of the detector (**1.732.0110 Detector block metal-free** or **1.732.0100 Metal**) must be inserted in the upper part of the Separation Center.




Fig. 7: Insert detector block in the Separation Center

2.820.0X10/2.820.0X30 IC Separation Center

The instrument versions 2.820.0X10 and 2.820.0X30 of the IC Separation Center are operated with one 819 IC Detector. It is best to proceed as follows when connecting the two instruments and the detector block:

1 Install detector block

- Unscrew the four knurled screws **15** from the top rear panel **21** of the 820 IC Separation Center and remove rear panel (see *Fig. 5*).
- Position detector block from the back in the space provided in the 820 IC Separation Center on the right and push fully to the front (see *Fig. 7*, resp. *Fig. 21* and *Fig. 23*).
- Insert the cable permanently attached to the detector block in opening **20** and the outlet capillary in opening **19** "Waste A" of the rear panel **21**.
- Replace rear panel **21** and screw to the 820 IC Separation Center using the four knurled screws **15**.




Fig. 8: Connection Detector 819 – 2.820.0X10/2.820.0X30

2 Connect detector block

- Plug the gray connecting cable permanently attached to the detector block into connection **6** „Detector Block“ of the 819 IC Detector and fasten to the instrument by tightening the screws in the cable connector (see *Fig. 8*).

3 Connect waste container

- Lead the outlet capillary of the detector block to a sufficiently large waste container and fix in place.

2.820.0X20 IC Separation Center

The two-channel versions 2.820.0X20 of the IC Separation Center can be operated with two 819 IC Detectors. It is best to proceed as follows when connecting the instruments and the two detector blocks:

1 Install detector blocks

- Unscrew the four knurled screws **15** from the top rear panel **21** of the 820 IC Separation Center and remove rear panel (see *Fig. 5*).
- Position first detector block from the back in the space provided in the 820 IC Separation Center on the right and push fully to the front (see *Fig. 7*, resp. *Fig. 22* and *Fig. 25*).
- Position second detector block from the back in the space provided in the 820 IC Separation Center on the left and push fully to the front (see *Fig. 7*, resp. *Fig. 22* and *Fig. 25*).
- Insert the cable permanently attached to the detector block A in opening **20** and the outlet capillary in opening **19** "Waste A" of the rear panel **21**.
- Insert the cable permanently attached to the detector block B in opening **16** and the outlet capillary in opening **17** "Waste B" of the rear panel **21**.

- Replace rear panel **21** and screw to the 820 IC Separation Center using the four knurled screws **15**.

2 Connect detector blocks

- Plug the gray connecting cable permanently attached to the detector block A into connection **6** „Detector Block“ of the **first 819 IC Detector** and fasten to the instrument by tightening the screws in the cable connector (see Fig. 9).
- Plug the gray connecting cable permanently attached to the detector block B into connection **6** „Detector Block“ of the **second 819 IC Detector** and fasten to the instrument by tightening the screws in the cable connector (see Fig. 9).

3 Connect waste container

- Lead the outlet capillary of both detector blocks to a sufficiently large waste container and fix in place.




Fig. 9: Connection 2 x 819 Detector – 2.820.0X20

2.2.2 Connection of column heating

The instrument versions 2.820.02X0 of the IC Separation Center contain a column heating for thermostating the separating column as an option. If you wish to use this column heating then it should be installed as follows:




Fig. 10: Install column heating in Separation Center

1 Laying the connection cable

- Screw off the 4 knurled screws **15** from upper rear panel **21** of the 820 IC Separation Center and remove the rear panel (see *Fig. 5*).
- Lead 6.2108.120 Cable through opening **20**.
- Replace rear panel **21** and screw it onto 820 IC Separation Center using the 4 knurled screws **15**.

2 Insert column heating

- Carefully place the column heating inside the Separation Center. If the column has already been inserted in the heating then be careful with the capillaries at feedthroughs **35** and **41** (see *Section 2.9.2*).
- Use magnetic holder **39** and the corresponding counter-pieces to attach the column heating to the left-hand inside wall of the Separation Center.

3 Connect column heating

- Insert 6.2108.120 Cable at connection **38** of the column heating (see *Fig. 6*) and at connection **23** on the rear panel of the Separation Center (see *Fig. 5*).




Fig. 11: Connection 820 – column heating

2.2.3 Connection of syringe and aspirating tubing

For manual filling of the sample loops mounted on the injection valves, the 6.2816.020 Syringe and the PTFE aspirating tubing already screwed to the valve are needed. These accessories are mounted or adjusted as follows:

1 Connect syringe

- Push 6.2816.020 Syringe (without needle) as far as it will go into connection socket **13** (left for valve A, right for valve B) (see Fig. 4).

2 Lead aspirating tubing to the outside

- Lead the PTFE aspirating tubing **63** (see Fig. 21, Fig. 22 resp. Fig. 23) through one of the feedthroughs **14** at the front or side of the separation center (see Fig. 4 and Fig. 5) to the outside.

2.2.4 Connection of the drain tube

The 820 IC Separation Center has a connection at the rear to which a drain tube for discharged liquids can be attached. Proceed as follows:

1 Connect drain tube

- Mount 6.1816.00 Silicone tubing on connection nipple **33** (see Fig. 5).

2 Lead drain tube to collecting vessel

- Lead the other end of the drain tube to a suitable collecting vessel and fix in place.

2.2.5 Connection of the 6.5324.000 Bottle rack (option)

The optional available 6.5324.000 Bottle rack for supply vessels can be placed on top of the IC system tower. The accessories include the supply vessels for eluent (2 L), regeneration solution (1 L) and rinsing solution (1 L). For the connection of the supply capillaries leading to the 818 IC Pump and the suppressor module, see the instructions given on the enclosed leaflet.

2.3 Electrical connection

2.3.1 Connection to 830 IC Interface

The instruments are controlled via the Metrohm «**IC Net**» software. They are normally connected to a **PC** via an **830 IC Interface**.

The **820 IC Separation Center** can be operated in two different ways. It can be operated via the **819 IC Detector**, or included as an independent instrument in the «**IC Net**» via the **830 IC Interface**. Both the cable connections and the installation under «**IC Net**» depend on the type of control.

Control of the Separation Center via 819

Normally the Separation Center is controlled with «**IC Net**» via the time program of the detector. When a system is being set up either **820.02X0 to IC Detector** or **820.03X0 to IC Detector** must be selected for the Separation Center in «**IC Net**».

In a 2-channel system each channel of the Separation Center is controlled via the time program of the assigned detector.




Fig. 12: Connection of 820 via 819 to 830

1 Connect 819 to 820

- Use 6.2125.090 Cable to connect Connection **10** “IC Separation Center” of the **819 IC Detector** to Connection **30** “IC Detector” of the **820 IC Separation Center** (see Fig. 12).

For the 2-channel versions 820.0X20 additionally:

- Use a further 6.2125.090 Cable to connect Connection **10** “IC Separation Center” of the **second 819 IC Detector** to Connection **27** “IC Detector” of the **820 IC Separation Center** (see Fig. 12).

2 Connect 819 to 830

- Use 6.2134.090 Cable (830 accessory) to connect Connection **9** “RS 232” of the **819 IC Detector** to a free RS 232 interface of **System 1** of the **830 IC Interface** (see Fig. 12).
- Use 6.2134.130 Cable to connect Connection **5** “Output” of the **819 IC Detector** to analog connection “Channels 1/2” of **System 1** of the **830 IC Interface** (see Fig. 12).

For 2-channel versions 820.0X20 additionally:

- Use a further 6.2134.090 Cable (830 accessory) to connect Connection **9** “RS 232” of the **second 819 IC Detector** to a free RS 232 interface of **System 2** of the **830 IC Interface** (see Fig. 12).
- Use 6.2134.130 Cable to connect Connection **5** “Output” of the **819 IC Detector** to analog connection “Channels 3/4” of **System 2** of the **830 IC Interface** (see Fig. 12).

3 Connect 820 to 830

This cable connection starts data recording (chromatogram window turns blue) when the valve in the 820 is switched to INJECT via a program command; not necessary if valve is operated manually:

- Use 6.2115.070 Cable to connect the positions **Integr.** and **Start** of terminal block **31** of the **820 IC Separation Center** to the event line **Start** of **System 1** of the **830 IC Interface** (see Fig. 12).

For 2-channel versions 820.0X20 additionally:

- Use a further 6.2115.070 cable to connect the positions **Integr.** and **Start** of terminal block **26** of the **820 IC Separation Centers** to the event line **Start** of **System 2** of the **830 IC Interface** (see Fig. 12).



Take care that each of the **RUN** connections at 830 and 820 and the **COM** connections at 830 and 820 are connected to each other.

4 Connect column heating to 830

For 820.02X0 only:

- Use 6.2134.040 Cable to connect RS 232 interface **22** of the **820 IC Separation Center** to a free RS 232 interface of the **830 IC Interface** (see Fig. 12).

Separate control of the Separation Center

Alternatively you have the possibility of controlling the Separation Center with its own time program from «IC Net». This requires you to select the appropriate **820.0210 to Interface**, **820.0220 to Interface**, **820.0230 to Interface**, **820.0310 to Interface**, **820.0320 to Interface** or **820.0330 to Interface** for the Separation Center when setting up the systems in «IC Net».




Fig. 13: Connection of 819 and 820 to 830

1 Connect 819 to 830

- Use 6.2134.090 Cable (830 accessory) to connect Connection **9** "RS 232" of the **819 IC Detector** to a free RS 232 interface of **System 1** of the **830 IC Interface** (see Fig. 13).
- Use 6.2134.130 Cable to connect Connection **5** "Output" of the **819 IC Detector** to the analog connection "Channels 1/2" of **System 1** of the **830 IC Interface** (see Fig. 13).

2 Connect 820 to 830

Control lines to switch the valve to FILL or INJECT:

- Use 6.2128.100 Cable to connect positions **Ground**, **Fill** and **Inject** of terminal block **31** of the **820 IC Separation Center** to two free event lines of **System 1** of the **830 IC Interface** (see *Fig. 13*).

Start data recording:

- Use 6.2115.070 Cable to connect positions **Integr.** and **Start** of terminal block **31** of the **820 IC Separation Center** to the event line **Start** of **System 1** of the **830 IC Interface** (see *Fig. 13*).

This cable connection starts data recording (chromatogram window turns blue) when the valve in the 820 is switched to INJECT via a program command; not necessary if valve is operated manually.

3 Connect column heating to 830

For 820.02X0 only:

- Use 6.2134.040 Cable to connect RS 232 interface **22** of the **820 IC Separation Center** to a free RS 232 interface of the **830 IC Interface** (see *Fig. 13*).

If a 2-channel system is used then the cable connections for the second detector and the second valve are made in a similar way.

2.3.2 Connection to PC

The 819 IC Detector can also be connected directly to a PC. In this case the 820 IC Separation Center must then be connected to the detector and controlled via its time program. In addition a 771 IC Compact Interface is required; this is included between the analog output of the detector and the RS 232 interface of the PC and carries out the A/D-conversion of the detector signal. If present, the column heating is also connected directly to an RS 232 interface of the PC.

This means that three free RS 232 interfaces are required at the PC for controlling the detector and Separation Center. If a complete IC system is to be operated directly from a PC then at least four free RS 232 interfaces are required (2 for detector, 1 for column heating and 1 for the high-pressure pump) at the PC.



Always switch off the 819 IC Detector, column heating and PC before connecting up the instruments with cables.

1 Connect 819 to 820

- Use 6.2125.090 Cable to connect Connection **10** "IC Separation Center" of the **819 IC Detector** to Connection **30** "IC Detector" of the **820 IC Separation Center** (see *Fig. 14*).

2 Connect 819 to PC

- Use 6.2125.110 Cable to connect Connection **9** “RS 232” of the **819 IC Detector** to a free RS 232 interface at the **PC** (see Fig. 14).

3 Connect 819 and 820 to 771

- Use 6.2134.140 Cable to connect Connection **5** “Output” of the **819 IC Detector** to Connection “Channel 1” of the **771 IC Compact Interface** (see Fig. 14). The cable end with the marking “Analog” is inserted into the detector.
- Attach the two cable ends “RUN” and “COM” of 6.2134.140 Cable to the positions **Integr.** and **Start** of terminal block **27** of the **820 IC Separation Center** (see Fig. 14).

3 Connect 771 to PC

- Use 6.2134.100 Cable to connect Connection “RS 232, PC” of the **771 IC Compact Interface** to a free RS 232 interface at the **PC** (see Fig. 14).

4 Connect column heating to PC

For 820.02X0 only:

- Use 6.2134.040 Cable to connect RS 232 interface **22** of the **820 IC Separation Center** to a free RS 232 interface at the **PC** (see Fig. 14).




Fig. 14: Connection of 819 and 820 to a PC

If a 2-channel system is used then the cable connections for the second detector and second valve are made in a similar way.

2.4 Mains connection



Follow the instructions below for connecting to the power supply. If the instrument is operated with a mains voltage set wrongly and/or wrong mains fuse, there is a danger of fire!

2.4.1 819 IC Detector

Setting the mains voltage

Before switching on the 819 IC Detector for the first time, check whether the mains voltage set on the instrument (see *Fig. 15*) matches the local mains voltage. If this is not the case, you must reset the mains voltage on the instrument as follows:

1 Disconnect mains cable

Disconnect mains cable from mains connection plug **3** of the 819 IC Detector.

2 Remove fuse holder

Using a screwdriver, loosen fuse holder **4** below the mains connection plug **3** and take out completely.

3 Change and check fuse if necessary

Carefully take the fuse installed for the desired mains voltage out of fuse holder **4** and check its specifications (the position of the fuse in the fuse holder is marked by the white arrow imprinted next to the mains voltage range):

100...120 V 0.63 A (slow blow) Metrohm-No. U.600.0014

220...240 V 0.315 A (slow blow) Metrohm-No. U.600.0011

4 Insert fuse

Change fuse if necessary and reinsert in fuse holder **4**.

5 Install fuse holder

Depending on the desired mains voltage, insert fuse holder **4** in the 819 IC Detector so that the corresponding mains voltage range can be read normally and the adjacent white arrow points to the white bar imprinted below the fuse holder (see *Fig. 15*).




Fig. 15: Setting the mains voltage

Fuses

One of the two fuses 0.63 A/slow-blow for 100...120 V or 0.315 A/slow-blow for 220...240 V is installed in fuse holder **4** of the 819 IC Detector as standard.



Ensure that the instrument is never put into operation with fuses of another type, otherwise there is danger of fire!

For checking or changing fuses, process as described above.

Mains cable

The instrument is supplied with one of three mains cables

- 6.2122.020 with plug SEV 12 (Switzerland, ...)
- 6.2122.040 with plug CEE(7), VII (Germany, ...)
- 6.2133.070 with plug NEMA 5-15 (USA, ...)

which are three-cored and fitted with a plug with an earthing pin. If a different plug has to be fitted, the yellow/green lead (IEC standard) must be connected to protective earth (protection class 1).



Any break in the earthing inside or outside the instrument can make it a hazard!

Plug the mains cable into mains connection plug **3** of the 819 IC Detector (see Fig. 3).

On/off switching of the instrument

The 819 IC Detector is switched on and off using mains switch **2** (see Fig. 3). When the instrument is switched on the mains pilot lamp **1** lights up.

2.4.2 820 IC Separation Center

If an 820 IC Separation Center is included in a modular IC system via an 819 IC Detector then the power supply for the Separation Center is provided by the 6.2125.090 Cable (see *Fig. 12* and *Fig. 14*).

If the 820 IC Separation Center is controlled separately (see *Fig. 13*) then the optionally available 6.2152.000 Power Supply must be used and inserted in Connection **25** (see *Fig. 5*).

2.4.3 Column heating

The column heating contained in the 2.820.02X0 versions has its own power supply. You must connect the 6.2152.010 Power Supply contained in the standard equipment to Connection **24** on the rear panel of the Separation Center (see *Fig. 5*).

2.5 Software-installation

The PC program «**IC Net 2.3**» is required for the operation of 819 IC Detector and 820 IC Separation Center. This program runs under Windows 2000 and Windows XP operating systems and is installed according to *Section 1.5.2* of the **Instructions for Use 8.110.8283** from the Metrodata software «**IC Net 2.3**».

The installation of the instruments is described in *Section 6.1.1 / Section 6.1.2* of the **Instructions for Use 8.110.8283** from the Metrodata software «**IC Net 2.3**». For the settings of the 819 IC Detector see *Section 6.4*, for the settings of 820 Separation Center see *Section 6.8*.

If the Separation Center is to be controlled via the time program of the detector (see *Fig. 12* and *Fig. 14*) then during installation the Separation Center must be installed appropriately as **820.02X0 to IC Detector** or **820.03X0 to IC Detector**. If the 820 IC Separation Center is to be controlled separately (see *Fig. 13*), then it must be installed as **820.0210 to Interface**, **820.0220 to Interface**, **820.0230 to Interface**, **820.0310 to Interface**, **820.0320 to Interface** or **820.0330 to Interface**.

2.6 Capillary connections

Some of the connections under high pressure between the feed pump and the detector block must be set up by the user. For that the **6.1831.010 PEEK capillary** (i.d. = 0.25 mm, o.d. = 1/16", length = 3 m) and the **6.2744.010 Compression fittings** supplied in the accessories of the 820 IC Separation Center can be used. These capillaries can be used in the pressure range from 0...30 MPa (0...300 bar).



Capillaries fitted with new connectors must have a perfectly flat cut surface. For PEEK capillaries it is best to use the **6.2621.080 Capillary tubing cutter**




Fig. 16: Connect capillaries

42 PEEK capillary 6.1831.010

43 Compression fitting 6.2744.010

For the connection of 6.1831.010 PEEK capillaries or of the 6.1834.000 PTFE aspirating tubing, the supplied **6.2744.010 PEEK Compression fittings** are used. It is best to proceed as follows:

1 Mount compression fitting

Slide a compression fitting **43** (6.2744.010) over the end of the capillary **42** to be fastened as shown in Fig. 16.

2 Insert capillary in connection

Push capillary end in the corresponding connection as far as it will go (to avoid dead volume).

3 Tighten compression fitting

Tighten compression fitting **43** by hand (never use tools).

2.7 Connection of 818 IC Pump

2.7.1 Electrical connection

For operation of the 819 IC Detector and 820 IC Separation Center you can use any commercial HPLC pump. However, as the attainable sensitivity depends to a large extent on the quality of the pump, Metrohm advises use of the **818 IC Pump**, which has been specially developed for the demands of ion chromatography and has minimal pump pulsation and an outstanding flow constancy.

Startup and operation of the 818 IC Pump are described in the *818 Instructions for Use*. The eluent, which must be degassed and filtered (cf. *Section 4.1.3*), is selected on the basis of the separating column installed in the 820 IC Separation Center and the current separation problem.

The **818 IC Pump** is normally integrated into the modular IC system via 830 IC Interface. Connect the RS 232 interface **17** of the 818 IC Pump with a free RS 232 interface at the 830 IC Interface using the 6.2134.080 or 6.2134.090 cable according to the figure below. To ensure proper functioning of the communication between the 830 IC Interface and 818 IC Pump, the sliding switch **16** on the IC pump must be set to "RS 232".




Fig. 17: Connection of 818 IC Pump to 830 IC Interface

2.7.2 Pulsation dampener

To protect the column material against pressure drops caused by the injector, the use of a pulsation dampener connected between the pump and the injection valve of the 820 IC Separation Center is recommended. The optional **6.2620.150 Pulsation dampener** (see *Section 6.3*) is very well suited to this purpose.

The metal-free 6.2620.150 Pulsation Dampener is supplied fully assembled and has two connections for capillaries, for which either the con-

nectors supplied or two 6.2744.010 PEEK compression fittings can be used. The flow direction is arbitrary. The pulsation dampener is positioned in the interior of the 820 IC Separation Center on the base below the injection valve (see Fig. 21, Fig. 22, Fig. 23, resp. Fig. 25).



The pulsation dampener is filled with isopropanol when new. Rinse your IC system carefully after the first installation of a new pulsation dampener.

2.7.3 Filter unit PEEK

The **6.2821.120 Filter unit PEEK 47** supplied (see Fig. 18) serves to avoid contamination of the piston seals by abrasive particles and can be used in the pressure range 0...30 MPa (0...300 bar). The filter unit consists of the housing **44**, the filter **45** and the connector **46** to be screwed into the housing **44**. For the connection of capillaries **42**, PEEK compression fittings **43** (6.2744.010) must be used.

New filter **45** are available as an option with the ordering number 6.2821.110 (10 pieces).



For the connection of the filter unit, please note the flow direction arrow printed on the housing.

The filter unit is filled with isopropanol when new. Rinse your IC system carefully after the first installation of a new filter unit.




Fig. 18: Filter unit PEEK (6.2821.120)

42 PEEK capillary (6.1831.010)

43 Compression fitting (6.2744.010)

44 Filter housing

**45 Filter
6.2821.130**

46 Filter connector

2.7.4 Connection to injection valve

It is recommended to use 6.1831.010 PEEK capillaries, a 6.2620.150 Pulsation dampener (see Section 2.7.2) and a 6.2821.120 Filter unit PEEK (see Section 2.7.3) to connect the 818 IC Pump to the injection valve of the 820 IC Separation Center. Proceed as follows:




Fig. 19: Connection to injection valve

3 Connection capillary to 818 IC Pump PEEK capillary (6.1831.010)	50 Column connection capillary PEEK capillary (6.1831.010)
47 Filter unit PEEK (6.2821.120)	51 Injection valve
48 PEEK capillary (6.1831.010)	
49 Pulsation dampener 6.2620.150	

- 1 Connection to 818 IC Pump**
 - Cut connection capillary **3** (6.1831.010 PEEK capillary) to the required length and equip it with compression fittings.
 - Attach connection capillary **3** to the upper connection of the purge valve (see Fig. 1 818 Instructions for Use).
- 2 Connection of filter unit PEEK**
 - Attach the other end of connection capillary **3** at the filter housing **44** of the filter unit **47** (see Fig. 18).
 - Connect a PEEK capillary 6.1831.010 cut to the required length and equipped with connectors to the filter connector **46** of the filter unit **47**.

3 Installation of the capillary in the IC Separation Center

- Unscrew the four knurled screws **15** of the bottom rear panel **28** of the 820 IC Separation Center and remove rear panel (see *Fig. 5*).
- Insert PEEK capillary **48** from the back into the inner compartment of the IC Separation Center.
- Install rear panel **28** so that the capillary is positioned in opening **34** „Inlet A“ or **32** „Inlet B“ and screw on with the four knurled screws **15**.

4 Connection of the pulsation dampener

- Procedure see section 2.7.2.

5 Connection to the injection valve

- Connect PEEK capillary **48** to pulsation dampener **49** (see *Section 2.7.2*). Using another PEEK capillary **48**, connect pulsation dampener to connection „**5**“ (see *Fig. 28*) of injection valve **51**.

6 Mount column connection capillary

This step is only necessary, if no column heating is used.

- Connect column connection capillary **50** (ca. 20 cm of 6.1831.010 PEEK capillary) to connection „**4**“ (see *Fig. 28*) of injection valve **51**.

2.8 Precolumns

The use of easily exchangeable precolumns protects the separating columns and prolongs their lifetime. The precolumns available from Metrohm are either real precolumns or precolumn cartridges, which are used together with a cartridge holder. For the installation of a precolumn cartridge into the accompanying cartridge holder see the attached leaflet.



The precolumn that is suitable for your separating column can be found in the **Metrohm IC Column Catalog** which can be obtained from your local Metrohm agency, the data sheet accompanying your separating column, the product information about separating columns that can be found under <http://www.metrohm.com>, ion chromatography products, or let your agency advise you directly.



New IC precolumns are normally filled with solution and sealed at both ends. Before the precolumn is installed in the system, it must be ensured that this solution is freely miscible with the eluent used (check manufacturer's specifications).

1 Connect precolumn

- Remove end caps from the precolumn.
- Fit compression fitting to the connection capillary **50** mounted to the injection valve (see Section 2.6).
- Screw column connection capillary **50** to precolumn.
- Cut off a small piece, ca. 5 cm, from the PEEK capillary 6.1831.010 and fit compression fittings to both ends (see Section 2.6).
- Mount the prepared capillary to the outlet of the precolumn.



When you install the column, always ensure that this is inserted correctly in accordance with the flow direction (if existent) shown.

2 Rinse the precolumn

- Place a beaker beneath outlet capillary.
- Switch on 818 IC Pump in «IC Net» and rinse precolumn for ca. 10 min with eluent.
- Switch off 818 IC Pump.

2.9 Separating columns and suppressor module

2.9.1 General information on separating columns



The precolumn that is suitable for your separating column can be found in the **Metrohm IC Column Catalog**, the product information about separating columns that can be found under <http://www.metrohm.com>, ion chromatography products, or let your agency advise you directly.



New IC separating columns are normally filled with solution and sealed at both ends. Before the column is installed in the system, it must be ensured that this solution is freely miscible with the eluent used (check manufacturer's specifications).

The separating columns currently available from Metrohm and precolumns can be found in the **Metrohm IC Column Catalog**, or on the Internet under <http://www.metrohm.com> in the ion chromatography product section. Each column is supplied with a test chromatogram and a data sheet. More detailed information about special IC applications can be found in the relevant "**Application Bulletins**" or "**Application Notes**"; these are available on the Internet under <http://www.metrohm.com> in the applications sector, or can be requested free of charge from your local Metrohm agency.



When you install the column, always ensure that this is inserted correctly in accordance with the flow direction shown.

2.9.2 Installing a separating column in the column heating




Fig. 20: Insert column in column heating

35 Feedthrough
for the capillary from separating column to detector

36 Knurled screw
for fastening the heater cover **37**

37 Cover of the column heating

41 Feedthrough
for column connection capillary **50**.

50 Column connection capillary
PEEK capillary 6.1831.010

52 Capillary to detector block
PEEK capillary 6.1831.010

53 Separating column

54 Holder plate
for separating column

55 Capillary prewarming
for column connection capillary **50**.

1 Open column heating

- Open the two knurled screws **36** and remove cover **37** of the column heating.

2 Prepare column heating

- Unscrew the six screws of holder plate **54** and remove the holder plate from the column heating.
- Lead column connection capillary **50** (approx. 1 m PEEK-Capillary 6.1831.010) through feedthrough **41** in the column heating and place in the recesses in capillary prewarmer **55** as shown in *Fig. 20*.
- Then screw down holder plate **54** again with the six screws.
- Cut off a further piece of 6.1831.010 PEEK Capillary, lead this Capillary **52** out through feedthrough **35** and provide its inner end with a 6.2744.010 Compression fitting **43**.



*If not much space is available then the **6.2744.070 PEEK compression fitting short** contained in the standard equipment can also be used for attaching the PEEK capillary.*

3 Insert separating column

- Press separating column **53** into the holding clips and fasten the column connection capillary **50** to the separating column inlet (bottom) with a 6.2744.010 Compression fitting **43**.
- Fasten capillary **52** from the detector to the outlet of the separating column (top) with a 6.2744.010 Compression fitting **43**.

4 Close column heating

- Screw down cover **37** on to the column heating again with the two knurled screws **36**.

2.9.3 Selection of the sample loop

Selection of the sample loop depends on the separating column used. Normally, the following sample loops are used:

Anion columns	100 µL
Cation columns	10 µL
Columns for suppressor technique	20 µL

Depending on the instrument version, the following sample loops are installed in the 820 IC Separation Center:

Version	Valve	Sample loop	Volume
2.820.0X10	A	6.1825.210 (PEEK)	20 µL
2.820.0X20	A	6.1825.220 (PEEK)	100 µL
	B	6.1825.210 (PEEK)	20 µL
2.820.0X30	A	6.1825.210 (PEEK)	20 µL

If desired, the built-in sample loop can be replaced by one of the sample loops available as an option (see Section 6.3.1).

2.9.4 General information on suppressor module

The **Metrohm Suppressor Module MSM** for chemical suppression installed in the 2.820.0X30 IC Separation Center comprises a total of 3 suppressor units which are in turn used for suppression, regenerated with sulfuric acid and rinsed with water. To record every new chromatogram under comparable conditions, work is normally carried out with freshly regenerated suppressor. Switching is either automatic together with the valve switching or manual.



The suppressor units must never be regenerated with H₂SO₄ in the same flow direction used for the eluent. You should thus always install the inlet and outlet capillaries as described in Section 2.9.7 according to the scheme shown in Fig. 23.

For operation of the suppressor module, a **two-channel peristaltic pump** is needed which conveys the regeneration solution (normally **50 mmol/L H₂SO₄**) and the rinsing solution (normally **dist. H₂O**) to the suppressor units. We advise working with a **flow rate of 0.5 mL/min**.



*The **2.833.0010 IC Liquid Handling Pump Unit** is available from Metrohm as an option. Two lengths of pump tubing (6.1826.050) are enclosed with this pump (flow rate 0.5 mL/min). Startup and operation of the 833 Pump Unit is described in the Instructions for Use enclosed with the pump.*

The three inlets and outlets numbered 1...3 on the suppressor module each have 2 permanently mounted PTFE capillaries, which must be connected as described in Section 2.9.7 (see Fig. 23 und Fig. 24).

To avoid contamination of the suppressor module by foreign particles or bacterial growth, it is advantageous to install an **in-line filter** between the 833 IC Liquid Handling Pump Unit and the suppressor module. For this the two **6.2821.120 Filter units PEEK** (mounting see *Fig. 18 Section 2.7.3*) supplied with the 2.833.0010 IC Liquid Handling Pump Unit should be used.



The suppressor module must never be switched in the dry state as there is a danger of blocking. Before every switching operation of the suppressor module, the three suppressor units must have been rinsed for at least ½ h with eluent, regeneration and rinsing solution.

2.9.5 One-channel system without suppressor module



Make sure that the connection separating column outlet – detector inlet is as short as possible in order to prevent unnecessary peak widening after the separation.




Fig. 21: Interior of the 2.820.0310 Separation Center

48 PEEK capillary (6.1831.010)
to 818 IC Pump

49 Pulsation dampener (6.2620.150)

50 Column connection capillary
PEEK capillary 6.1831.010

51 Injection valve

53 Separating column	60 Capillary for syringe 13 PEEK capillary, fixed mounting
56 Detector block	61 Sample loop (6.1825.210) 20 µL, PEEK, fixed mounting
57 Inlet capillary for detector block fixed mounting	62 Inlet capillary for injector PEEK-Capillary 6.1831.010
58 Mounting rail	63 PTFE aspirating tube fixed mounting
59 Column holder (6.2027.0X0)	

With the one-channel system without suppressor module, the IC separating column is installed in the 2.820.0X10 IC Separation Center as follows (see *Fig. 21*):

1 Connect column to injector

- Remove end caps from column **53**.
- *without column heating:*
Screw inlet end of separating column **53** (note flow direction) to column connection capillary **50** mounted on the injector.
- *with column heating:*
Prepare column heating according to *Section 2.9.2* and screw column connection capillary **50** (see *Fig. 20*) with a compression fitting to injection valve **51**.
- *With precolumn:*
Install precolumn according to the supplied leaflet between inlet of the separating column and the injection valve.

2 Rinse column

- Place a beaker beneath the column outlet.
- Start 818 IC in «IC Net» with suitable flow (see *leaflet of the column*) and rinse column for ca. 10 min with eluent.
- Stop 818 IC Pump.

3 Connect column to detector block

- *without column heating*
Screw outlet end of separating column **53** to the inlet capillary **57** permanently mounted on the detector block **56**.
- *with column heating*
Connect outlet end of separating column **52** (see *Fig. 20*) with coupling (6.2744.040) to the inlet capillary **57** permanently mounted on the detector block **56**.

4 Fix column

- *without column heating*
Insert one or two column holders **59** (6.2027.030, 6.2027.040 or 6.2027.050) in the mounting rails **58** and fasten separating column in the column holder.
- *with column heating*
Insert column heating according to *Fig. 10* into the Separation Center.

2.9.6 Two-channel system without suppressor module



Make sure that the connection separating column outlet – detector inlet is as short as possible in order to prevent unnecessary peak widening after the separation.

With the two-channel system without suppressor module (2.820.0X20 IC Separation Center), the first IC separating column is connected on the left side to injection valve A and detector block A as with the one-channel system (see Section 2.9.5 und Fig. 21). The second column is connected on the right side to injection valve B and detector block B in an analogous manner.




Fig. 22: Interior of the 2.820.0320 Separation Center

2.9.7 One-channel system with suppressor module



Make sure that the connection separating column outlet – detector inlet is as short as possible in order to prevent unnecessary peak widening after the separation.

In the case of the one-channel system with suppressor module, first the IC separating column is installed in the 2.820.0X30 IC Separation Center (see Fig. 23) and then the suppressor module is connected to the 833 IC Liquid Handling Pump Unit needed for operation (see Fig. 24). Proceed as follows:




Fig. 23: Interior of 2.820.0330 Separation Center

48 PEEK capillary (6.1831.010) to 818 IC Pump	62 Inlet capillary for injector PEEK-Capillary 6.1831.010
49 Pulsation dampener (6.2620.150)	63 PTFE aspirating tube fixed mounting
50 Column connection capillary PEEK capillary 6.1831.010	64 Suppressor inlet capillary for eluent
51 Injection valve	65 Suppressor inlet capillary for H₂SO₄
53 Separating column	66 Suppressor outlet capillary for H₂SO₄
56 Detector block	67 Suppressor outlet capillary for H₂O
57 Inlet capillary for detector block fixed mounting	68 Suppressor inlet capillary for H₂O
58 Mounting rail	69 Suppressor outlet capillary for eluent
59 Column holder (6.2027.0X0)	70 Suppressor module
60 Capillary for syringe 13 PEEK capillary, fixed mounting	71 Coupling (6.2744.040; PEEK)
61 Sample loop (6.1825.210) 20 µL, PEEK, fixed mounting	

1 Connect column to injector

- Remove end caps from column **53**
- *without column heating:*
Screw inlet end of separating column **53** (note flow direction) to column connection capillary **50** mounted on the injector.
- *with column heating:*
Prepare column heating according to *Section 2.9.2* and screw column connection capillary **50** (see *Fig. 20*) with a compression fitting to injection valve **51**.
- *With precolumn:*
Install precolumn according to the supplied leaflet between inlet of the separating column and the injection valve.

2 Rinse column


- Place a beaker beneath the column outlet.
- Start 818 IC Pump in «IC Net» with suitable flow (see *leaflet of the column*) and rinse column for ca. 10 min with eluent.
- Stop 818 IC Pump.

3 Connect column to suppressor module (Suppressor connection 1)

- Cut inlet capillary **64** (marked with "Eluent") of suppressor module **70** to the required length. Use the 6.2126.080 Capillary cutter available as an option.
- *without column heating*
Screw inlet capillary **64** on to the outlet end of separating column **53** using a 6.2744.010 Compression fitting.
- *with column heating*
Connect capillary **52** at the outlet of the separating column (see Fig. 20) to the inlet capillary **64** using the 6.2744.040 Coupling.

4 Fix column

- *without column heating*
Insert one or two column holders **59** (6.2027.030, 6.2027.040 or 6.2027.050) in the mounting rails **58** and fasten separating column in the column holder.
- *with column heating*
Insert column heating according to Fig. 10 into the Separation Center.



64 Suppressor inlet capillary for eluent

65 Suppressor inlet capillary for H₂SO₄

66 Suppressor outlet capillary for H₂SO₄

67 Suppressor outlet capillary for H₂O

68 Suppressor inlet capillary for H₂O

69 Suppressor outlet capillary for eluent

Fig. 24: Suppressor module connections

5 Connect suppressor module to detector block (Suppressor connection 1)

- Cut outlet capillary **69** (marked with "Detector") of suppressor module **70** to the required length. Use the 6.2126.080 Capillary cutter available as an option.
- Screw outlet capillary 69 on to coupling 71 by using a 6.2744.010 Pressure screw.
- Screw inlet capillary **56** of detector block **57** on to the other end of coupling **71**.

6 Prepare 833 IC Liquid Handling Pump Unit

- Mount the two lengths of **6.1826.050** Pump tubing (**white-yellow** stopper) with aspirating tube **16** and filter unit PEEK **22** 6.2821.120 as described in Section 3.2, Fig. 7 of the 833 Instructions for Use.
For mounting of the filter unit PEEK see also section 2.7.3.



Pump tubing is a consumable whose working life depends on the contact pressure. This is why when the pump is switched off for a long time, the tubing cassette should be lifted by fully loosening spring lever on the right-hand side (in this way the correctly set optimal contact pressure is retained).

7 Suppressor connection 2: H₂SO₄

- Pull inlet capillary **65** (marked with "H₂SO₄") by hand as far as required through one of the feedthroughs **14** to the outside (see Fig. 23 und Fig. 24).
- Screw the filter unit PEEK on to PEEK coupling **21** (6.2744.160) as described in Section 3.2 of 833 Instructions for Use at the outlet end of the **first** piece of pump tubing **19**.
- Fasten inlet capillary **65** (**23**) to the outlet of filter unit PEEK (6.2821.120) **47** by using a compression fitting screw (see Fig. 7 of 833 Instructions for Use).
- At the inlet end of the **first** piece of tubing **19** fasten a suitably long piece of PTFE tubing **16** (6.1803.020) to coupling **18** (6.2744.034) by using a compression fitting (6.2744.010) (see Fig. 7 of 833 Instructions for Use).
- Immerse the other end of the aspiration tubing in a container with regeneration solution (normally 50 mmol/L H₂SO₄) and fix it in position.
- Pull outlet capillary **66** of the suppressor module (marked with "Waste") through opening **18** on the rear panel, lead it into a sufficiently large waste bottle and fix it in position.

8 Suppressor connection 3: H₂O

- Pull inlet capillary **68** (marked with "H₂SO₄") by hand as far as required through one of the feedthroughs **14** to the outside (see *Fig. 23* und *Fig. 24*).
- Screw the filter unit PEEK on to PEEK coupling **21** (6.2744.160) as described in *Section 3.2* of *833 Instructions for Use* at the outlet end of the **second** piece of pump tubing **19**.
- Fasten inlet capillary **68** (**23**) to the outlet of filter unit PEEK (6.2821.120) **47** by using a compression fitting screw (see *Fig. 7* of *833 Instructions for Use*).
- At the inlet end of the **second** piece of tubing **19** fasten a suitably long piece of PTFE tubing **16** (6.1803.020) to coupling **18** (6.2744.034) by using a compression fitting (6.2744.010) (see *Fig. 7* of *833 Instructions for Use*).
- Immerse the other end of the aspiration tubing in a container with regeneration solution (normally 50 mmol/L H₂SO₄) and fix it in position.
- Pull outlet capillary **67** of the suppressor module (marked with "Waste") through opening **18** on the rear panel, lead it into a sufficiently large waste bottle and fix it in position.

2.9.8 Two-channel system with suppressor module

Make sure that the connection separating column outlet – detector inlet is as short as possible in order to prevent unnecessary peak widening after the separation.




Fig. 25: Interior of 2.820.032 Separation Center with suppressor module

In order to operate a 2-channel Separation Center with chemical suppression a 2.833.0020 IC Liquid Handling Suppressor Unit is additionally required. For this reason the 1.753.0100 Suppressor Module of the 833 IC Liquid Handling Suppressor Unit is placed on the floor of the inner compartment of the Separation Center (see *Fig. 25*).

In a 2-channel system with suppressor module the first IC separating column, just like in a 1-channel system, is connected on the left-hand side to injection valve A and detector block A (similar to *Section 2.9.5*, *Fig. 21*). Components of this channel are shown in dark gray in *Fig. 25*.

The second column, just like in a 1-channel system with suppressor module (see *Section 2.9.7* and *Fig. 23*), is connected on the right-hand side to injection valve B and detector block B. The installation of the 833 IC Liquid Handling Suppressor Unit is described in the *833 Instructions for Use*; connection to the 2.820.0X20 Separation Center is described there in detail in *Section 5*. Components of this channel are shown in light gray in *Fig. 25*.

2.9.9 Leak testing and conditioning

Before sample solutions can be injected in the IC system, the entire system must be tested for leaks and then conditioned with eluent until the baseline is stable. Proceed as follows:

1 Switch on instruments

- Switch on 818 IC Pump at mains switch.
- Switch on 819 IC Detector at mains switch **2**.
- The 820 IC Separation Center is normally supplied via the 819 IC Detector, see also *Section 2.4.2*.
- If present, switch on 833 IC Liquid Handling Pump Unit at the mains switch.
- If present, switch on 830 IC Interface at the mains switch.
- Switch on PC and start «IC Net», load a system corresponding to the device configuration or create a new one (see «*IC Net*» *Instructions for Use*).
- Start the system with **Control/Startup hardware**. The 818 IC Pump and, if present, the 833 IC Liquid Handling Unit will start to run, and the baseline will be shown in the chromatogram window.

2 Check 818 IC Pump

- Check that the aspiration capillary of 818 IC Pump is immersed in the eluent.
- In the loaded system open the **818 IC Pump** window and check/set the values for the flow rate and the maximum switch-off pressure under **Settings**.
- The recommended flow rate for the separating column used can be found in the corresponding *column data sheet* (normally 0.5...2 mL/min).
- The maximum switch-off pressure should normally be approx. 3 MPa above the pressure observed with the used column, for details see *818 Instructions for Use*.
- Start 818 IC Pump with **<Start>**.

3 Check for leaks

- Check all capillaries and their connections between the 818 IC Pump and the detector block for escaping liquid. If eluent escapes anywhere, the appropriate compression fitting must be tightened further or changed..
- Check also all tubing's and their connections between the 833 IC Liquid Handling Unit and the suppressor module (if present) for escaping liquid and tighten or change the fittings if necessary.

4 Check 819 IC Detector

- Open the window **819 IC Detector** of the loaded system and check/set the values of operating temperature, cell constant, measuring range and full-scale range on the **Method Parameter** tab (see sec. 3.1.2).
- Set **measuring range Range** so that the displayed absolute conductivity value of the eluent lies within the selected range (default value = 1 mS/cm).
- Set **full-scale range, Full scale** (default value = range/1). To start with, it is advisable not to set too narrow a full-scale range as the conductivity of the eluent in the conditioning phase can change greatly until attainment of a stable temperature.

5 Condition IC system

- Rinse the IC system with eluent until the desired stability of the baseline is reached (normally 30...60 min).

6 Condition suppressor (if present)

- Open window **820 IC Separation Center** and switch suppressor module at least three times to the next position every 10 min.

3 Operation


The operation of 819 IC Detector and 820 IC Separation Center is made completely via the Metrohm software «**IC Net**».



This Section describes only the most important points concerning the operation of the 819 IC Detector and 820 IC Separation Center. For further details please refer to the «IC Net» Instructions for Use and to the on-line help in the PC program.

3.1 819 IC Detector

3.1.1 819 IC Detector icon




3.1.2 Settings in the "819 IC Detector" window

By double-clicking the 819 icon or by clicking this icon with the right mouse button and selecting the **Open** menu item the **819 IC Detector** window for parameter settings is opened. It consists of the four tabs **Display**, **Method parameter**, **Program** and **Links**.

Display

The **Display** tab of the **819 IC Detector** window is only available for a **connected system**. It displays the current measurement values and parameters settings of the 819 IC Detector.



- Conductivity** Live display of measured conductivity (auto-zero value) in μS/cm. This field is available for the **WATCH WINDOW**. The watch window color is changed to the **Out of range** color if the measured value is more than 80% outside the **Full Scale** range.

- abs. Cond.** Live display of absolute conductivity value in μS/cm. This field is available for the **WATCH WINDOW**. The watch window color is changed to the **Out of range** color if the measured value is bigger or equal to 110% of **Range**.

- Ion** Display of ion settings:
Anion: polarity = +, temperature coeff. = 2.5
Cation: polarity = -, temperature coeff. = 1.5

- Range** Display of measuring range in μS/cm.
- Full scale** Display of full scale range in μS/cm.
- Cell const.** Display of cell constant of the conductivity cell in /cm.

- Thermostat** Display of temperature within in the detector block.

- Offset** Display of offset of the zero point.
- Damping** Display of electronic damping of the analog output.

- Remote lines** Display value of remote line settings 1...8.

- Zero <ON>** Triggers the auto-zero function.
- Zero <OFF>** Switch off the auto-zero function.

"Auto-zero" is the name used for the automatic electronic background compensation, i.e. the current measured **Conductivity** value is set to zero and hence lies in the middle of the selected full-scale range.

The auto-zero function is disabled, if the actual value of the absolute conductivity **abs. Cond.** out of the limit set by the measuring **Range** (see Fig. 27), the diagram below shows the operating principle of the auto-zero function.




Fig. 26: Operating principle of Auto-zero

Before the auto-zero function is activated for the first time or after it has been switched off, the measuring signal lies outside the full-scale range; it is automatically set to 0 $\mu\text{S}/\text{cm}$ each time the **Zero <ON>** key is pressed.

The auto-zero function can be triggered at any time using the **Zero <ON>** key.




The auto-zero function works properly only when the measured value remains relatively stable (e.g. not during the appearance of the injection peak at the start of the chromatogram).

If the measured conductivity value exceeds $\pm 180\%$ of the set full-scale range, the live display of **Conductivity** in the **WATCH WINDOW** changes its color to the settings for **Out of range**. (default: red). At the same time the absolute conductivity is displayed instead of the auto-zero value.

Method parameter

The **Method parameter** tab of the **819 IC Detector** window contains the method startup values for the instrument. These values are sent and applied to the 819 IC Detector by **Start determination** and **Startup hardware**, or if the values are sent manually with **<SET>**.



Range

Measuring range in $\mu\text{S}/\text{cm}$
 Selection: **100, 200, 500 $\mu\text{S}/\text{cm}$** and
1, 2, 5, 10 mS/cm

Select the measuring range so that the conductivity value of the eluent used is always within the selected range. If the measured conductivity value exceeds the upper range limit by more than 10%, the live display of **abs. Cond.** in the **WATCH WINDOW** changes its color to the settings for **Out of range**. (default: red). The selected measuring range limits the possible settings of the full-scale range, which can be selected to be maximum 2000 times more sensitive than the measuring range.

Full scale

The full scale range (operating range) sets the desired sensitivity for the analog output. The possible values of the full scale range in $\mu\text{S}/\text{cm}$ depend on the preset measuring **Range**.

This value can be reduced in a total of 11 stages by a factor 1 (lowest sensitivity) up to maximum a factor of 2000 (highest sensitivity).

The full-scale range is primarily used to increase the sensitivity for the recording of the chromatogram following the electronic background compensation of the eluent sensitivity (auto-zero function, see also description of the auto-zero function above).

For the selected full-scale range (FS), the linearity of the conductivity measurement is assured in the range -1.8 Full Scale... $+1.8$ Full Scale. If the measured value violates these limits, the live display of **Conductivity** in the **WATCH WINDOW** changes its color to the settings for **Out of range**. (default: red). At the same time the absolute conductivity is displayed instead of the auto-zero value.

The full-scale range determines the limits for the output of the analog signal at analog output **5** (0...1 V). The polarity of the output signal can be changed at any time with the parameter **Polarity** to + (0...+Full Scale) or - (0...-Full Scale).





Fig. 27: Diagram of measuring and full-scale range

Ion	Combination of polarity and temperature coefficient:
Anion	polarity = +, temperature coeff. = 2.5
Cation	polarity = -, temperature coeff. = 1.5
Custom	independent selection of polarity and temperature coeff.
Polarity	Polarity of the analog output signal: + positive polarity (e.g. for anions) - negative polarity (e.g. for cations)
Temperature coeff.	Converts automatically the conductivity from the operating temperature of the measuring cell to the reference temperature of 20°C. Selection: 1.5, 2.5 %/°C
Cell constant	Cell constant of the conductivity cell in /cm. Entry range: 13 ... 21 /cm In ion chromatography, interest is normally centered only on relative changes in the conductivity and not on the absolute value. With the cell constant of 16.7/cm set in the factory, the error in the display of the absolute conductivity can be up to ca. ±10%. If you wish to have a more accurate display of the absolute conductivity, the cell constant must be determined by means of a calibration solution. To do this, pump a solution of known conductivity through the IC system, observe the display of the absolute conductivity and change the cell constant until the displayed value matches the actual value.
Thermostat	Operating temperature of the conductivity cell within in the detector block. Selection: 25, 30, 35, 40, 45 °C and off At constant ambient temperature the heating built into the 819 IC Detector regulates the selected measuring cell temperature to an accuracy of ± 0.01 °C and thus establishes the precondition for highly sensitive determinations. It normally takes 30...60 min after the instrument has been switched on until this temperature stability is attained.
Offset	Offset of the zero point of the conductivity. Selection: 0, 10, 50 %full-scale range

Damping	Electronic damping of the analog output Selection: on, off Damping to reduce disturbing pulsations or a high noise level. Standard solutions and samples must both be recorded with damping as this can influence the peak heights.
Remote lines	Startup value for remote line settings 1...8. Selection: 0, 1
<SET>	Send current parameters immediately to the 819 IC Detector. Parameters are not stored in the system file (*.smt) as long as the file is not saved.

Time program

On the **Program** tab of the **819 IC Detector** window a user-defined time program for the remote output lines of the 819 IC Detector can be entered. This program is started automatically as defined in the **Start mode** window (see Section 4.4.3 of the «IC Net» *Instructions for Use*) either at the moment the determination is started (**Start with determination**) or at the moment the sample is injected (**Start with inject**).



The **Program** tab contains the two following subpages:

Program	Main time program with all program steps.
Remote configuration	Possibility for creation of user-defined remote commands.

Program

On the **Program** subpage, program steps including time, program instruction and parameter can be entered.

First column	Time at which program instruction is applied. Entry range: 0.0 ... 999.9 min If no time is entered, the program instruction is applied together with the last instruction with time entry.
Second column	Program instruction (see below). In addition to these predefined instructions, user-defined remote commands can be entered if activated on the Remote configuration tab.
Third column	Parameter for program instruction (see below).
ENABLED	Enable program start (a disabled program is not started).
<Add>	Add new program instruction.
<Delete>	Delete selected program instruction.
<Verify>	Test the time program (error messages are displayed if program is wrong).

List of program instructions


The following program instructions can be added to the time program on the **Program** subpage:

Instruction	Parameter entry	Meaning
Remote	0, 1, *	Set remote output lines 1...8 to the desired values. For entry of the first value, enter 1, 0 or * . For entry of the other values, move the cursor in front of the value to be changed and enter 1, 0 or * .
Range	100, 200, 5000 µS/cm 1, 2, 5, 10 mS/cm	Set measuring range to the selected value.
FS	0.05 µS/cm ... 10 mS/cm	Set full scale range to the selected value (depends on Range setting).
Zero	on, off	Switch on or off the auto-zero function.
Polarity	+, -	Reset the polarity of the output signal.
Mark		Trigger marking signal.

<i>Instruction</i>	<i>Parameter entry</i>	<i>Meaning</i>
ValveA	fill, inject	Switch injection valve A to "inject" or "fill" position.
ValveB	fill, inject	Switch injection valve B to "inject" or "fill" position.
Suppressor		Switch suppressor module to the next position.

Remote configuration

On the **Remote configuration** subtab user-defined remote commands can be defined, which can be inserted into a time program after being activated with **<Activate>**.



Name (1st column) User-definable name of the remote command (e.g. **Start_766**).

Remote command (2nd column)

Setting the remote output lines 1...8.

0 (line off, inactive, open)

1 (line on, active, 0 V)

***** (leave line in current status)

For entry of the first value, enter **1**, **0** or *****. For entry of the other values, move the cursor in front of the value to be changed and enter **1**, **0** or *****.


<Add>	Add new remote command.
<Delete>	Delete selected remote command.
<Activate>	Activate the defined remote commands for insertion into the time program.

Links

The **Links** tab of the **819 IC Detector** is used for COM port selection and settings (details see *Section 5.3.4 Links* of the «*IC Net*» *Instructions for Use*).

3.2 820 IC Separation Center

3.2.1 820 IC Separation Center icon



3.2.2 "820 IC Separation Center" window

By double-clicking the 819 icon or by clicking this icon with the right mouse button and selecting the **Open** menu item the **820 IC Separation Center** window for parameter settings is opened. It consists of the four tabs **Manual**, **Program**, **Thermostat** and **Links**.

The title bar of the window displays **820 IC Separation Center**, (if connected to IC Detector) or **820.01X0 IC Separation Center** (if connected to IC Interface).




Pay attention when adding the Separation Center to the system! The kind of Separation Center you have to choose from the device list depends on the wiring of the instruments in your system.

*Choose the **820.02X0 to IC Detector OR 820.03X0 to IC Detector** from the device list, if the **Separation Center** is connected to the **819 IC Detector**. The type of Separation Center is recognized automatically.*


*Choose the appropriate **820.0210 to Interface, 820.0220 to Interface, 820.0230 to Interface, 820.0310 to Interface, 820.0320 to Interface OR 820.0330 to Interface** from the device list, if the **Separation Center** is connected to the **830 IC Interface**.*

Manual

The **Manual** tab of the **820 IC Separation Center** window is only available for a **connected system**. The parameters and buttons for manual operation depend on the type of instrument connected.



820 IC Separation Center connected to **819**



820.0230 Separation Center connected to **830**

Configuration

- valve A**
- valve A + valve B**
- valve A + suppressor**

Display of current instrument configuration.*

820 IC Separation Center with 1 injector.

820 IC Separation Center with 2 injectors.

820 IC Separation Center with 1 injector and 1 suppressor module.

Valve A

- fill/inject**
- <FILL>**
- <INJECT>**

Display of current valve position.*

Switch valve A to "fill" position.

Switch valve A to "inject" position.

Valve B

- fill/inject**
- <FILL>**
- <INJECT>**

Display of current valve position.*

Switch valve B to "fill" position.

Switch valve B to "inject" position.

Suppressor

- in position**
- <STEP>**


Display of current suppressor status.*

Switch suppressor module to next position.

Time program

The **Program** tab is only available for a **820.0XX0 Separation Center connected to a 830 IC Interface**. The **820 IC Separation Center** that is **connected to a 819 IC Detector** is controlled by the time program of the detector.

On the **Program** tab of the **820.0XX0 Separation Center** window a user-defined time program can be entered. This program is started automatically as defined in the **Start mode** window (see *Section 4.4.3* of the «IC Net» *Instructions for Use*) either at the moment the determination is started (**Start with determination**) or at the moment the sample is injected (**Start with inject**).



First column	Time at which program instruction is applied. Entry range: 0.0 ... 999.9 min If no time is entered, the program instruction is applied together with the last instruction with time entry.
Second column	Program instruction (see below).
Third column	Parameter for program instruction (see below).
ENABLED	Enable program start (a disabled program is not started).
<Add>	Add new program instruction.
<Delete>	Delete selected program instruction.
<Verify>	Test the time program (error messages are displayed if program is wrong).

List of program instructions

The following program instructions can be added to the time program on the **Program** page:


Instruction	Parameter entry	Meaning
ValveA	fill, inject	Switch injection valve A to “fill” or “inject” position.
ValveB	fill, inject	Switch injection valve B to “fill” or “inject” position.
Step		Rotate suppressor to the next position.

Thermostat

The **Thermostat** tab of the **820 IC Separation Center** or **820.02X0 IC Separation Center** window contains the temperature control of the column heating.



*Only available for a 2.820.02X0 Separation Center that is connected either to the 819 IC Detector and installed as **820.02X0 to IC Detector**, or to the 830 IC Interface and installed as **820.02X0 to Interface**.*



Actual temperature

Display of **actual temperature** within the column heating in °C.

Set point

Operating temperature of the thermostat within the column heating in °C.
Entry Range: **10.0 ... 80.0**.

<On>	Only active if the thermostat is off . Switch the thermostat of the column heating on and apply the operating temperature specified under Set point , the button turns to <Set> .
<Set>	Only active if the thermostat is on and a new operating temperature is entered at Set point . Sets the new operating temperature.
<Off>	Only active if the thermostat is on . Switch the thermostat of the column heating off .

Start with startup hardware

Start the **Thermostat** with the startup of the system hardware.

Enabled

The temperature of the column heating **is monitored** by «IC Net», an injection is triggered, if the set operating temperature of the column heating is reached.



The operating temperature is reached, when the actual temperature is within $\pm 1\text{ }^{\circ}\text{C}$ of the temperature set above.

Disabled

The temperature of the column heating **is not monitored** by «IC Net». The proceeding of a determination or a queue is not affected if the set value is not reached.

Links

The **Links** tab of the **820 IC Separation Center** is used for COM port selection and settings (details see *Section 5.2.4 Links* of the «IC Net» *Instructions for Use*), resp. for Event line selection.

4 Notes - Maintenance - Faults

4.1 Practical notes on ion chromatography

4.1.1 Separating columns

The attainable quality of analyses with the 819/820 IC system depends to a large extent on the separation efficiency of the column used. When purchasing an IC column you should ensure that the separation efficiency suffices for the analysis problems at hand. Ascertain the **characteristic data of the IC column** on the standard chromatogram enclosed with the column such as capacity factors, selectivity, plate number and resolution and check these data with your own measurements. If any difficulties arise, you should always first check the quality of the column by recording a **standard chromatogram**.

Detailed information about the separating columns available from Metrohm can be found in the data sheet accompanying your separating column, in the **Metrohm IC Column Catalog**, which is available from your local Metrohm agency, or on the Internet under <http://www.metrohm.com> in the ion chromatography product sector. Information about special IC applications can be found in the relevant "**Application Bulletins**" or "**Application Notes**"; these are available on the Internet under <http://www.metrohm.com> in the applications sector, or can be requested free of charge from your local Metrohm agency.

Protection

To protect the column against foreign particles which could have an adverse influence on the separation efficiency, we advise you to subject both the eluents and all samples to **microfiltration** (0.45 µm filter) and to siphon the eluent through the **6.2821.090 Aspirating Filter**.

To avoid contamination by abrasive particles arising from piston seals of the 818 IC Pump, it is advantageous to install an **in-line filter** between the pump and the 820 IC Separation Center. **The 6.2821.120 Filter unit PEEK 47** (see *Section 2.7.3*) supplied with the 818 IC Pump for operation with PEEK capillaries in the pressure range 0...30 MPa.

The use of readily interchangeable **precolumns** serves to protect the actual separating columns and increase their service life appreciably. The precolumn that is suitable for your separating column can be found in the **Metrohm IC Column Catalog** which is available from your local Metrohm agency, the data sheet accompanying your separating column, the product information about separating columns that can be found under <http://www.metrohm.com>, ion chromatography products, or let your agency advise you directly.

Storage

Always store the separating columns closed when not in use and filled in accordance with the manufacturer's specifications.

Dead volume

Dead volume at the end of a column can be the cause of extreme peak broadening or splitting (appearance of double peaks). Filling the column with glass beads ($\varnothing \leq 100 \mu\text{m}$) frequently improves the separation efficiency.

Regeneration

If the separation properties of the column have deteriorated, it can be regenerated in accordance with the column manufacturer's specifications. With the separating columns available from Metrohm (see **Metrohm's "Range of IC columns"**), the instructions for regeneration can be found on the leaflet enclosed with every column.



*In the case of separating columns with carrier material based on silica, **only solutions with pH 2...7** may be used for regeneration, otherwise the columns could be damaged.*

4.1.2 Pumps

Pulsation dampener

For sensitive measurements with the 819/820 IC System, high pressure pumps as free from pulsation as possible and with a very constant capacity are required. If the pulsations that appear are too large, a pulsation dampener between the pump and the 820 IC Separation Center may possibly help. The optional **6.2620.150 Pulsation dampener** (see *Section 6.3.1*), whose installation is described in *Section 2.7.2* is eminently suitable for this purpose. A pulsation dampener also offers protection against pressure shocks on the column material caused by injection.

Maintenance

To protect the pump against foreign particles, we advise you to subject the eluent to **microfiltration** (0.45 μm filter) and siphon the eluent through the **6.2821.090 Aspirating Filter**.

In many cases, an unstable baseline (pulsation, flow fluctuations) can be traced to contaminated valves or faulty, leaky piston seals.

Contaminated valves are cleaned by rinsing with water, RBS solution or organic solvents. The rinsing effect is improved by brief treatment in an ultrasonic bath. When the cleaned valves are reinstalled, you must ensure that the flow direction is correct.

The replacement of piston seals has to be done in accordance with the pump manufacturer's directions. The corresponding maintenance work for the 818 IC Pump is described in *Section 4.2* of the *818 Instructions for Use*.

Salt crystals between the piston and the seal are the cause of abrasive particles, which can enter the eluent. These lead to contaminated

valves, pressure rise and in extreme cases to scratched pistons. It is thus essential to ensure that **no precipitates** can appear (see also Section 4.1.4).

4.1.3 Eluents

Treatment

For the preparation of the eluents only chemicals of a purity degree of at least "**p.a.**" should be used. For diluting please use only **high purity water**.

Fresh eluents should always be **microfiltered** (0.45 µm filter) and **de-gassed** (with N₂, He or vacuum). For high sensitive measurements, the eluent should be continuously **stirred** with a magnetic stirrer, particularly when the recycling procedure is employed or when alkaline eluents are used. For alkaline eluents and eluents with low buffering capacity one should preferably use a **CO₂ absorber** (e.g. the absorber supplied with the optional 6.5324.000 Bottle rack).

The supply vessel containing the eluent must be closed as tightly as possible to avoid excessive evaporation. This is primarily important with eluents containing organic solvents (e.g. acetone), the evaporation of which can lead to drifts in the long term. If work is performed in a very sensitive range, even if one drop of condensate falls back in the eluent this can cause a noticeable change in the background conductivity.

Influence of various parameters on anion columns

- *Concentration:* An increase in the concentration usually leads to shorter retention times and quicker separation, but also to a higher background conductivity.
- *pH:* pH alterations lead to shifts in the dissociation equilibrium and thus to changes in the retention times.
- *Organic modifiers:* Addition of an organic solvent (e.g. methanol, acetone, acetonitrile) to aqueous eluents generally accelerates lipophilic ions.

Eluent change

When the eluent is changed, it must be ensured that **no precipitates** can be formed. Solutions used in direct succession must therefore be miscible. If the system has to be rinsed with an organic solution, several solvents with increasing or decreasing lipophilic character may possibly have to be used (e.g. water ↔ acetone ↔ chloroform).

4.1.4 Suppressor module

To avoid contamination of the suppressor module by foreign particles or bacterial growth, it is advantageous to install an **in-line filter** between the 833 IC Liquid Handling Pump Unit and the suppressor module. The two **6.2821.120 Filter units PEEK** (see *Section 2.7.3* und *Section 2.9.4*) supplied with the 2.833.0010 IC Liquid Handling Pump Unit are eminently suitable for this purpose.

4.1.5 Connections

All connections between injector, column and detector must be as short as possible, have a low dead volume and be absolutely tight. The PEEK capillary after the detector block must be free from constriction (the measuring cell is tested to 5 MPa = 50 bar back pressure).

4.2 Maintenance and servicing

4.2.1 General information
Care

The 819 IC Detector and 820 IC Separation Center require proper care and attention. Excessive contamination of the instruments could possibly lead to malfunctions and a shorter service life of the inherently rugged mechanical and electronic parts.

Spilled chemicals and solvents should be wiped up immediately. It is especially important to protect the plug connections at the rear of the instrument (particular the mains plug) against contamination.



Although constructional measures have been designed to virtually eliminate such a situation, should corrosive media penetrate the interior of the instruments the mains plug of the 819 IC Detector must be immediately disconnected to prevent extensive damage to the instrument electronics. Inform Metrohm service if your instrument(s) have been damaged in such a way.



The instrument must not be opened by untrained personnel. Please comply with the safety notes in Section 1.4.1.

Maintenance by Metrohm service

Maintenance of the 819 IC Detector and 820 IC Separation Center is best done as part of an annual service performed by specialists from the Metrohm company. If work is frequently performed with caustic and corrosive chemicals, it may be necessary to shorten the interval between servicing.

The Metrohm service department is always willing to offer expert advice on the maintenance and servicing of all Metrohm instruments.

4.2.2 Passivation

Passivation of the entire IC system (without column) by rinsing with 20...50 mL 0.2 mol/L HNO₃ is appropriate only in the following case:

- if exceptional changes are observed in the measurement properties of the cell

For passivation the separating column is removed from the 820 IC Separation Center. The two capillaries **50** and **57** (see Fig. 21 and Fig. 23) are connected directly to each other with the **6.2744.040 Coupling 71** included in the accessories.

4.2.3 Recycling

To keep the eluent consumption between injections to a minimum when the system is at rest (e.g. overnight), the so-called recycling procedure can be used. In recycling the eluent exiting the outlet capillary of the detector block is led back directly to the eluent supply vessel. The IC system is thus quickly ready for new injections without a long conditioning period.



*The recycling procedure must **not** be used*

- *in operation with the suppressor module,*
 - *with alkaline eluents,*
- with the 6.1010.000 IC Cation column METROSEP Cation 1-2.*

4.2.4 Shutdown

If the 819/820 IC-System is shut down for a considerable length of time, the entire IC system (**without** column and suppressor) must be **rinsed free from salt** with methanol/water (1:4) to avoid crystallization of eluent salts with the corresponding subsequent damage.

For rinsing the separating column is removed; the two capillaries **50** and **57** (see Fig. 21 and Fig. 23) are directly connected to each other with the 6.2744.040 Coupling enclosed with the accessories. Rinse with methanol/water (1:4) until the conductivity drops below 10 µS/cm.

4.2.5 Changing separating columns

Identical separation system

If you wish to replace an IC separating column by a column of the same type, proceed as follows (see *Fig. 21* and *Fig. 23*):

1 Remove old column

- Switch off pump drive of the 818 IC Pump.
- Unscrew column from inlet capillary **57** of the detector block or from suppressor inlet capillary **64**.
- Unscrew column from column connection capillary **50** or the precolumn.

2 Connect new column to injector

- Remove end caps from column **53**.
- Screw inlet end of separating column **53** (note flow direction) to column connection capillary **50** or to the precolumn (see *Section 2.9.5/Section 2.9.7*).

3 Rinse column

- Place beaker beneath the column outlet.
- Switch on 818 IC Pump and rinse column with eluent for ca. 10 min, then switch off pump.

4 Connect column to detector block

- Screw outlet end of separating column **53** to inlet capillary **57** or suppressor inlet capillary **64**.

Changing the separation system

If you wish to replace an IC separating column by a column of a different type, proceed as follows (see *Fig. 21* and *Fig. 23*):

1 Remove old column

- Switch off pump drive of the 818 IC Pump.
- Unscrew column from inlet capillary **57** of the detector block or from suppressor inlet capillary **64**.
- Unscrew column from column connection capillary **50** or the precolumn.

2 Connect detector block to injection valve

- Connect column connection capillary **50** using a **6.2744.040 Coupling 71** directly to inlet capillary **57** of the detector block (see *Fig. 21*).

3 Rinse with dist. H₂O

- Rinse IC system with dist. or deionized water (flow rate 1 mL/min) for approx. 10 min.

4 Rinse with eluent

- Rinse IC system with the eluent needed for the separating column you intend to install until the absolute conductivity displayed on the 819 IC Detector is stable.

5 Remove coupling

- Remove **6.2744.040 Coupling 71** between column connection capillary **50** and inlet capillary **57** (see *Fig. 21*).

6 Connect new column to injector

- Remove end caps from column **53**.
- Screw inlet end of separating column **53** (note flow direction) to column connection capillary **50** or to the precolumn (see *Section 2.9.5/Section 2.9.7*).

7 Rinse column

- Place beaker beneath the column outlet.
- Switch on 818 IC Pump and rinse column with eluent for ca. 10 min, then switch off pump.

8 Connect column to detector block

- Screw outlet end of separating column **53** to inlet capillary **57** or suppressor inlet capillary **64**.

4.2.6 Regeneration of suppressor

Regenerating a suppressor with reduced capacity

If the suppressor units are exposed to certain heavy metals (e.g. iron) or organic contaminants for long periods of time, these can no longer be completely removed by the regeneration solution normally used (20 mmol/L H₂SO₄). This diminishes the capacity of the suppressor units, which, in milder cases, results in a reduced sensitivity to phosphates and, in severe cases, in a strong increase in the baseline. If such capacity problems occur at one or several positions, the suppressor units must be treated as follows:

1 Disconnect suppressor from IC system

- Disconnect suppressor from separating column and detector.

2 Regenerate suppressor

- Rinse each suppressor unit for about 10 min with one of the following solutions:

Contamination with heavy metals

0.2 mol/L H₂SO₄
+ 0.1 mol/L oxalic acid (with Fe²⁺/Fe⁺³)

Contamination with organic substances

0.2 mol/L H₂SO₄ / acetone ≥ 20%



The 6.1826.050 pump tubing is made of PVC and must not be used for rinsing with solutions which contain acetone. In such cases, rinse with different pump tubing or a different pump.

3 Connect suppressor to IC system

- Reconnect suppressor to the IC system. If capacity problems persist, replace the suppressor rotor (see *Section 4.2.8*).

4.2.7 Cleaning the suppressor

It may be necessary to clean the suppressor in the following cases:

- High counterpressure on the suppressor connection tubing
- Irremediable blockage of the suppressor (the suppressor can no longer deliver solutions)
- Irremediable obstruction of the suppressor (the suppressor can no longer be switched to next position)

To clean the connection piece and the rotor, proceed as follows (see *Fig. 29*):

1 Disconnect suppressor from IC system

- Disconnect input lead of suppressor from the separating column and output lead to the detector.

2 Dismantle suppressor

- Unscrew nut **72** from suppressor holder **75**.
- Pull out connection piece **73** and suppressor rotor **74** from suppressor holder **75** (the connection piece and the rotor normally stick together).
- Loosen connection piece **73** from suppressor rotor **74**.

3 Clean input and output leads

- Connect each of the 6 capillary tubings attached to connection piece **73** to the pump one after another, and pump through ultrapure water.
- Check whether solution emerges from connection piece **73**. If one of the input or output leads remains blocked, replace the connection piece **73** (order number 6.2832.010).

4 Clean suppressor rotor

- Clean the sealing surface of suppressor rotor **74** using a lint-free cloth and ethanol.

5 Insert suppressor rotor

- Insert suppressor rotor **74** in suppressor holder **75** in such a way that the tubing connections at the rear of the rotor fit in the corresponding openings inside the rotor, and that one of the three holes in the rotor can be seen from below in one of the openings of the holder.
- If the rotor has been inserted correctly, its sealing surface will be about 4 mm inside the holder. If this is not the case, bring the rotor into the correct position from below with the aid of a sharp object (e.g. a screwdriver).

6 Clean connection piece

- Clean the sealing surface of connection piece **73** using a lint-free cloth and ethanol.




Fig. 29: Assembling the suppressor

72 Screw nut

73 Connection piece with input and output leads
(6.2832.010) for suppressor module

74 Suppressor rotor(6.2832.000)

75 Suppressor holder

7 Insert connection piece

- Insert connection piece **73** in suppressor holder **75** in such a way that connection "1" is at the top, and that the three lugs on the connection piece fit in the corresponding openings of the holder.
- Screw nut **72** onto the thread of suppressor holder **75** manually (do not use tools).

8 Connect and condition the suppressor

- Reconnect the suppressor to the IC system.
- Before switching the suppressor to the next position for the first time, rinse all 3 suppressor units with solution for 5 min.

4.2.8 Replacement of suppressor module

The suppressor in the suppressor block may have to be replaced in the following cases:

- Irremediable loss of suppressor capacity (reduced phosphate sensitivity and/or strong rise in baseline)
- Irremediable blockage of the suppressor (the suppressor can no longer deliver solutions)

Both the 6.2832.000 Suppressor rotor and the 6.2832.010 Connection piece with the input and output leads can be replaced. To replace these components proceed as follows (see *Fig. 29*):

1 Disconnect suppressor from IC system

- Disconnect all input and output leads of the suppressor from IC system and the 833 IC Liquid Handling Pump Unit.

2 Dismantle suppressor

- Unscrew nut **72** from suppressor holder **75**.
- Pull out connection piece **73** and suppressor rotor **74** from suppressor holder **75** (the connection piece and the rotor normally stick together).
- Loosen connection piece **73** from suppressor rotor **74**.

3 Clean suppressor rotor

- Clean the sealing surface of new suppressor rotor **74** (6.2832.000) using a lint-free cloth and ethanol.

4 Insert suppressor rotor

- Insert suppressor rotor **74** in suppressor holder **75** in such a way that the tubing connections at the rear of the rotor fit in the corresponding openings inside the rotor, and that one of the three holes in the rotor can be seen from below in one of the openings of the holder.
- If the rotor has been inserted correctly, its sealing surface will be about 4 mm inside the holder. If this is not the case, bring the rotor into the correct position from below with the aid of a sharp object (e.g. a screwdriver).

5 Clean connection piece

- Clean the sealing surface of connection piece **73** using a lint-free cloth and ethanol.

6 Insert connection piece

- Insert connection piece **73** in suppressor holder **75** in such a way that connection "1" is at the top, and that the three lugs on the connection piece fit in the corresponding openings of the holder.
- Screw nut **72** onto the thread of suppressor holder **75** manually (do not use tools).

7 Connect and condition the suppressor

- Reconnect the suppressor to the IC system.
- Before switching the suppressor to the next position for the first time, rinse all 3 suppressor units with solution for 5 min.

4.3 Faults and malfunctions

4.3.1 Malfunctions and their rectification

If difficulties appear with the IC system during analyses, their causes are best investigated in the order **separating column → pump → eluent → 819/820 IC System**. Several of the malfunctions which may appear are listed in the following table with details of possible causes and countermeasures.

Malfunction	Cause	Rectification
Baseline with high noise level, pulsation	<ul style="list-style-type: none"> Contaminated pump valves Faulty piston seals Quality of the pump does not suffice for the selected sensitivity 	<ul style="list-style-type: none"> Clean the valves (see Section 4.1.3) Replace the piston seals (see Section 4.1.3) Use pulsation dampener, use more powerful pump or lower the sensitivity
Drift of the baseline	<ul style="list-style-type: none"> Thermal equilibrium not yet reached Leak in system Evaporation of organic solvent in eluent 	<ul style="list-style-type: none"> Condition system with heating switched on Check connections and make leakproof Ensure better closure of eluent supply vessel
Considerable pressure drop	<ul style="list-style-type: none"> Leak in system 	<ul style="list-style-type: none"> Check connections and make leakproof
Considerable pressure rise	<ul style="list-style-type: none"> Contamination of the filter in the 6.2821.120 Filter unit PEEK Change of column packing by injection of contaminated samples 	<ul style="list-style-type: none"> Replace the 6.2821.130 Filter (see Section 2.7.3) Regenerate the column (see Section 4.1.2) or replace column <i>Note:</i> <i>Samples should always be microfiltered.</i>
Chromatograms with poor resolution, change in the retention times	<ul style="list-style-type: none"> Deterioration in separation efficiency of the IC column 	<ul style="list-style-type: none"> Regenerate the column (see Section 4.1.2) or replace column
Extreme peak broadening, splitting (double peaks)	<ul style="list-style-type: none"> Dead volume at the column ends 	<ul style="list-style-type: none"> Fill column with glass beads ($\varnothing \leq 100 \mu\text{m}$) or replace column

No feed of regeneration or rinsing solution for suppressor	<ul style="list-style-type: none">• Contamination of the filter in the 6.2821.120 Filter unit PEEK• Counterpressure too high in suppressor module	<ul style="list-style-type: none">• Replace the 6.2821.130 filter (see <i>Section 2.7.3</i>)• Clean or replace suppressor module (see <i>Section 4.2.6...4.2.8</i>)
Detector is not recognized by «IC Net»	<ul style="list-style-type: none">• No RS connection	<ul style="list-style-type: none">• Switch instrument on• Check wiring• Check that the connection specified in the 819 IC Detector window (see <i>Section 3.1.2 Links</i>) matches the actual wiring

5 Interfaces

5.1 RS 232 interface

The RS 232 interface **9** of the 819 IC Detector is used for the connection of the 830 IC Interface or a PC (see *Section 2.3.1 – 2.3.2*) and allows remote operation of 819 IC Detector and 820 IC Separation Center.

5.1.1 Data transmission protocol

The RS 232 interface **9** of the 819 IC Detector is configured as DTE (Data Terminal Equipment) with the following technical specifications:

- *Standard* Data interface in accordance with EIA standard RS 232C (DIN 66020, page 1)
- *Baud rate* 9600
- *Data Bits* 8
- *Stop Bits* 1
- *XON/XOFF* on
- *Control characters*

C _R	DEC 13	HEX 0D
L _F	DEC 10	HEX 0A
XON	DEC 17	HEX 11
XOFF	DEC 19	HEX 13
- *Max. line length* 80 characters + C_R L_F
- *Mode* Full duplex (simultaneous transmission and receiving)

Restriction:
If data are first received by the interface, transmission is not started until the receipt is at an end.
- *Cable* For interconnection of the 819 IC Detector with third-party devices, only a shielded data cable may be used (e.g. Metrohm D.104.0201). The cable shield must be perfectly earthed at both devices (pay attention to current loops: always use star earthing). Use only connectors with sufficient shielding (e.g. Metrohm K.210.0001 with K.210.9004).
- *Cable length* max. ca. 15 m

5.1.2 Pin assignment









RS 232 C interface	external
<p>Transmitted Data (TxD) If no data transmission takes place, the line is maintained in the "ON" condition. Data are sent only when CTS and DSR are in the "ON" condition and DCD is in the "OFF" condition.</p> <p>Received Data (RxD) Data are received only when DCD is "ON".</p> <p>Request to Send (RTS) ON condition: 819 IC Detector is ready to send data.</p> <p>Clear to Send (CTS) ON condition: Remote station is ready to receive data.</p> <p>Data Set Ready (DSR) ON condition: The transmission line is connected.</p> <p>Signal Ground (GND)</p> <p>Data Carrier Detect (DCD) ON condition: The received signal level is within the tolerance range (remote station is ready to send data).</p> <p>Data Terminal Ready (DTR) ON condition: 819 IC Detector is ready to receive data.</p>	<p>Pin 2 Transmitted Data</p> <p>Pin 3 Received Data</p> <p>Pin 4 Request to Send</p> <p>Pin 5 Clear to Send</p> <p>Pin 6 Data Set Ready</p> <p>Pin 7 Signal Ground</p> <p>Pin 8 Data Carrier Detect</p> <p>Pin 20 Data Terminal Ready</p>
<p>Protective earth Direct connection from cable connector to protective earth of the device.</p> <p>Polarity assignment of the signals</p> <ul style="list-style-type: none"> Data lines (TxD, RxD) <ul style="list-style-type: none"> Voltage negative (<-3 V): signal status "ONE" Voltage positive (> +3 V): signal status "ZERO" Control or message lines (CTS, DSR, DCD, RTS, DTR) <ul style="list-style-type: none"> Voltage negative (<-3 V): OFF status Voltage positive (> +3 V): ON status <p>In the transition region from +3 V to -3 V, the signal status is undefined.</p> <p>Driver 14C88 to EIA RS 232C specification</p> <p>Receiver 14C89 to EIA RS 232C specification</p>	<p>Contact arrangement at the socket "RS 232" (male)</p> <p>Connection cables to external devices must have an appropriate 25-pin connector (female). Metrohm offers connector sockets (K.210.9004) and the associated housing (K.210.0001) as an option for the preparation of such cables.</p>
<p><i>No liability whatsoever will be accepted for damage arising from the improper connection of devices.</i></p>	

5.2 Remote interfaces

The 819 IC Detector has the two remote interfaces **8** "Remote" and **10** "IC Separation Center" (see Fig. 3) which each have 8 input and 8 output lines.

5.2.1 "Remote" interface

The remote interface **8** "Remote" is used for the connection of any type of external device and has the following pin assignment:

819 IC Detector	Pin	Function																																																																																																																														
Inputs																																																																																																																																
	21	Triggering the print command  Print																																																																																																																														
	9																																																																																																																															
	22	2⁰ Triggering of functions (remote control) Assignment of the input lines <table border="1" data-bbox="746 1012 1396 1556"> <thead> <tr> <th></th> <th colspan="4">Input</th> <th>Decimal code</th> <th>Function (default settings)</th> </tr> <tr> <th></th> <th>23</th> <th>10</th> <th>22</th> <th>9</th> <th></th> <th></th> </tr> </thead> <tbody> <tr><td>H</td><td>H</td><td>H</td><td>H</td><td>H</td><td>0</td><td>Inactive</td></tr> <tr><td>H</td><td>H</td><td>H</td><td>L</td><td>L</td><td>1</td><td>PROG R/S</td></tr> <tr><td>H</td><td>H</td><td>L</td><td>H</td><td>H</td><td>2</td><td>PUMP R/S</td></tr> <tr><td>H</td><td>H</td><td>L</td><td>L</td><td>L</td><td>3</td><td>FILL B/STEP</td></tr> <tr><td>H</td><td>L</td><td>H</td><td>H</td><td>H</td><td>4</td><td>FILL A</td></tr> <tr><td>H</td><td>L</td><td>H</td><td>L</td><td>L</td><td>5</td><td>REPORT</td></tr> <tr><td>H</td><td>L</td><td>L</td><td>H</td><td>H</td><td>6</td><td>ZERO</td></tr> <tr><td>H</td><td>L</td><td>L</td><td>L</td><td>L</td><td>7</td><td>MARK</td></tr> <tr><td>L</td><td>H</td><td>H</td><td>H</td><td>H</td><td>8</td><td>INJECT A</td></tr> <tr><td>L</td><td>H</td><td>H</td><td>L</td><td>L</td><td>9</td><td>ZERO OFF</td></tr> <tr><td>L</td><td>H</td><td>L</td><td>H</td><td>H</td><td>10</td><td>PLOT</td></tr> <tr><td>L</td><td>H</td><td>L</td><td>L</td><td>L</td><td>11</td><td>Switch polarity</td></tr> <tr><td>L</td><td>L</td><td>H</td><td>H</td><td>H</td><td>12</td><td>INJECT B</td></tr> <tr><td>L</td><td>L</td><td>H</td><td>L</td><td>L</td><td>13</td><td>SELECT</td></tr> <tr><td>L</td><td>L</td><td>L</td><td>H</td><td>H</td><td>14</td><td>QUIT</td></tr> <tr><td>L</td><td>L</td><td>L</td><td>L</td><td>L</td><td>15</td><td>ENTER</td></tr> </tbody> </table>		Input				Decimal code	Function (default settings)		23	10	22	9			H	H	H	H	H	0	Inactive	H	H	H	L	L	1	PROG R/S	H	H	L	H	H	2	PUMP R/S	H	H	L	L	L	3	FILL B/STEP	H	L	H	H	H	4	FILL A	H	L	H	L	L	5	REPORT	H	L	L	H	H	6	ZERO	H	L	L	L	L	7	MARK	L	H	H	H	H	8	INJECT A	L	H	H	L	L	9	ZERO OFF	L	H	L	H	H	10	PLOT	L	H	L	L	L	11	Switch polarity	L	L	H	H	H	12	INJECT B	L	L	H	L	L	13	SELECT	L	L	L	H	H	14	QUIT	L	L	L	L	L	15	ENTER
	Input				Decimal code	Function (default settings)																																																																																																																										
	23	10	22	9																																																																																																																												
H	H	H	H	H	0	Inactive																																																																																																																										
H	H	H	L	L	1	PROG R/S																																																																																																																										
H	H	L	H	H	2	PUMP R/S																																																																																																																										
H	H	L	L	L	3	FILL B/STEP																																																																																																																										
H	L	H	H	H	4	FILL A																																																																																																																										
H	L	H	L	L	5	REPORT																																																																																																																										
H	L	L	H	H	6	ZERO																																																																																																																										
H	L	L	L	L	7	MARK																																																																																																																										
L	H	H	H	H	8	INJECT A																																																																																																																										
L	H	H	L	L	9	ZERO OFF																																																																																																																										
L	H	L	H	H	10	PLOT																																																																																																																										
L	H	L	L	L	11	Switch polarity																																																																																																																										
L	L	H	H	H	12	INJECT B																																																																																																																										
L	L	H	L	L	13	SELECT																																																																																																																										
L	L	L	H	H	14	QUIT																																																																																																																										
L	L	L	L	L	15	ENTER																																																																																																																										
	10																																																																																																																															
	23	2³ 																																																																																																																														
	11 24 12																																																																																																																															
		not assigned																																																																																																																														

819 IC Detector	Pin	Function
Outputs		
	5	<p>Ready L = ready (program inactive or not started) H = not ready (program running)</p> <p>$V_{CE0} = 40\text{ V}$ $IC = 20\text{ mA}$</p>
	18	<p>Remote 1 Programmable remote output line 1 (command "remote" in program or event)</p> <p>$V_{CE0} = 40\text{ V}$ $IC = 20\text{ mA}$</p>
	4	<p>Remote 2 Programmable remote output line 2 (command "remote" in program or event)</p> <p>$V_{CE0} = 40\text{ V}$ $IC = 20\text{ mA}$</p>
	17	<p>Advance pulse (EOD) The EOD signal is outputted in the following cases:</p> <ul style="list-style-type: none"> • <u>without</u> program: on every "INJECT A" • <u>with</u> program: each time a program is ended and on every return to the loop start with loop programs <p></p>
	3	<p>Remote 3 Programmable remote output line 3 (command "remote" in program or event)</p> <p>$V_{CE0} = 40\text{ V}$ $IC = 20\text{ mA}$</p>
	16	<p>Error L = instrument error (as long as it exists) H = instrument ok</p> <p>$V_{CE0} = 40\text{ V}$ $IC = 20\text{ mA}$</p>
	1	<p>Remote 4 Programmable remote output line 4 (command "remote" in program or event)</p> <p>$V_{CE0} = 40\text{ V}$ $IC = 20\text{ mA}$</p>
	2	<p>Remote 5 Programmable remote output line 5 (command "remote" in program or event)</p> <p>$V_{CE0} = 40\text{ V}$ $IC = 20\text{ mA}$</p>
Voltages		
	15	<p>5 V: inactive High (H) $I \leq 40\text{ mA}$ 0 V: active Low (L) $R_i \cong 12\ \Omega$</p>
	14	
	25	
Contact arrangement at socket "Remote" (female)		
		<p>Connecting cables to external devices must have an appropriate 25-pin connector (male). Metrohm offers connector sockets (K.210.9060) and the associated housing (K.210.0002) as an option for the preparation of such cables.</p>
<p><i>No liability whatsoever will be accepted for damage arising from the improper connection of devices.</i></p>		

5.2.2 "IC Separation Center" interface


The 820 IC Separation Center is normally attached to remote interface **10** "IC Separation Center" (see *Section 2.3*). If the 819 IC Detector is operated without a Separation Center, this remote interface can also be used for the connection of external devices. The remote interface "IC Separation Center" has the following pin assignment:

819 IC Detector	Pin	Function															
Inputs																	
	21	<p>Position "FILL" at injection valve A L = Valve A in position "FILL" H = Valve A not in position "FILL"</p>															
	9	<p>Position "INJECT" at injection valve A L = Valve A in position "INJECT" H = Valve A not in position "INJECT"</p>															
	22	<p>Position "FILL" at injection valve B or "in position" at suppressor L = Valve B in position "FILL" or suppressor set to "in position" H = Valve B not in position "FILL" or suppressor not set to "in position"</p>															
	10	<p>Position "INJECT" at injection valve B L = Valve B in position "INJECT" H = Valve B not in position "INJECT"</p>															
Configuration of the 820 IC Separation Center																	
	2^0	<table border="1"> <thead> <tr> <th>Input</th> <th>Decimal code</th> <th>Status</th> </tr> </thead> <tbody> <tr> <td>H H</td> <td>0</td> <td>undefined</td> </tr> <tr> <td>H L</td> <td>1</td> <td>valve A</td> </tr> <tr> <td>L H</td> <td>2</td> <td>valve A + suppressor</td> </tr> <tr> <td>L L</td> <td>3</td> <td>valve A + valve B</td> </tr> </tbody> </table>	Input	Decimal code	Status	H H	0	undefined	H L	1	valve A	L H	2	valve A + suppressor	L L	3	valve A + valve B
Input	Decimal code	Status															
H H	0	undefined															
H L	1	valve A															
L H	2	valve A + suppressor															
L L	3	valve A + valve B															
	2^1																
	11																
	24	<p>Connection 820 IC Separation Center L = 820 IC Separation Center connected H = 820 IC Separation Center not accessible</p>															
	12	not assigned															

819 IC Detector	Pin	Function
Outputs		
	5	<p>Valve A → "FILL" Valve A is switched to position "FILL".</p>
	18	<p>Valve A → "INJECT" Valve A is switched to position "INJECT".</p>
	4	<p>Lock Terminal block 26, 31 L = functions <FILL> and <INJECT> are locked H = functions <FILL> and <INJECT> are not locked</p> <p>$V_{CEO} = 40\text{ V}$ $I_C = 20\text{ mA}$</p>
	17	<p>Valve B → "FILL" or suppressor → "STEP" Valve B is switched to position "FILL" or suppressor module is switched to next position.</p>
	3	<p>Valve B → "INJECT" Valve B is switched to position "INJECT".</p>
	2	<p>Remote 6 Programmable remote output line 6 (command "remote" in program or event)</p> <p>$V_{CEO} = 40\text{ V}$ $I_C = 20\text{ mA}$</p>
	1	<p>Remote 7 Programmable remote output line 7 (command "remote" in program or event)</p> <p>$V_{CEO} = 40\text{ V}$ $I_C = 20\text{ mA}$</p>
	16	<p>Remote 8 Programmable remote output line 8 (command "remote" in program or event)</p> <p>$V_{CEO} = 40\text{ V}$ $I_C = 20\text{ mA}$</p>
Voltages		
<p>+5V → 15</p> <p>0V → 14</p> <p>24VAKT → 13</p> <p>0VAKT → 19</p> <p>→ 25</p> <p>→ 6</p> <p>→ 8</p> <p>→ 13</p> <p>→ 19</p> <p>→ 20</p>		<p>5 V: inactive High (H) $I \leq 40\text{ mA}$</p> <p>0 V: active Low (L) $R_i \cong 12\ \Omega$</p>
Contact arrangement at socket "IC Separation Center" (female)		
<p>Connecting cables to external devices must have an appropriate 25-pin connector (male). Metrohm offers connector sockets (K.210.9060) and the associated housing (K.210.0002) as an option for the preparation of such cables.</p>		
<p><i>No liability whatsoever will be accepted for damage arising from the improper connection of devices.</i></p>		


5.3 Analog output

Analog output **5** "Output" (0...1 V) of the 819 IC Detector (see Fig. 3) is used for connecting a recording channel to 830 IC Interface or 771 IC Compact Interface or from an external recorder or integration system. The connection from 830 IC Interface and 771 IC Compact Interface is described in Section 2.3. The polarity of the output signal at the analog output connectors can be changed at any time under **Polarity** in the window **819 IC Detector** in «IC Net» (see Section 3.1.2). The switching of the analog output is as follows:

819 IC Detector	Pin	Voltage
Contact arrangement at socket "Output" (female) 	1	0 V
	6	Signal 0...1 V The other pins are not assigned.






5.4 External power supply for 820 IC Separation Center

If the 820 IC Separation Center is operated without the 819 IC Detector, an external power supply must be attached to connection **25** "External Power". Power supply and connecting cable must comply with the following conditions:

"External Power" socket	Function	
	1 +5 V DC / 0.5 A	
	2 not assigned	
	3 +24 V DC / 2 A	
	4 0 V (digital) *	* Pin 4 and pin 5 must be connected in the vicinity of the power source ("zero point of star")
	5 0 V (analog) *	
	6 Ground	
	7 Ground	

5.5 Valve interfaces

The two valve interfaces **26** and **31** can be used to connect the 820 IC Separation Center directly to the 830 IC Interface (see *Section 2.3.1*) or external devices (e.g. 766 IC Sample Processor or 788 IC Filtration/Dialysis Sample Processor). The two valve interfaces for the injection valves A and B (or the suppressor module) are identical and have the following connections:

820 IC Sep. Center	Connection	Function
Inputs		
	Ground	Ground
	Fill	Valve → "FILL" Valve is switched to the position "FILL" (or suppressor module is switched further)
	Inject	Valve → "INJECT" Valve is switched to the position "INJECT"
Outputs		
	COM RUN	Position "Fill" A pulse is outputted when the valve is switched to the position "FILL".
	COM RUN	Integrator Start A pulse is outputted when the valve is switched to the position "INJECT".
<i>No liability whatsoever will be accepted for damage arising from the improper connection of devices.</i>		

6 Appendix

6.1 Technical data



Provided that nothing to the contrary is mentioned, the published data are typical values at an ambient temperature of 25°C.

6.1.1 819 IC Detector

Conductivity measurement

<i>Measurement ranges</i>	100, 200, 500, 1000, 2000, 5000, 10'000 µS/cm
<i>Full Scale</i>	0.05...10'000 µS/cm
<i>Temperature correction</i>	Automatic correction of the measured conductivity to reference temperature 20°C with the adjustable temperature coefficients (1.5 %/°C or 2.5 %/°C)
<i>Accuracy, absolute value</i>	< 4 % of full-scale value for all ranges under the following reference conditions: Adjustment with reference resistance to 1.000 V (range 1 mS/cm, cell constant 16.7 /cm, temp. coeff. 2.5 %/°C, frequency 5 kHz, room temperature) to accuracy < ± 2 mV, measurement accuracy < 0.5 % after 30 min warm-up time
<i>Linearity</i>	Up to ± 150 % of each full-scale range Deviations < 0.5 % of full-scale range
<i>Temperature dependence</i>	Typically 25 ppm/°C
<i>Dependence of cell const..</i>	< ± 1 % over entire adjustment range
<i>Drift (electronic)</i>	< 0.0013 % of selected measurement range /h/°C
<i>Noise (electronic)</i>	Typically < 0.0003 % of selected meas. range
<i>Damping</i>	2-stage damping (Bessel 4 th order) "Damping off" 0.25 s (10...90 %) "Damping on" 2.00 s (10...90 %)

Autozero

<i>Function</i>	Automatic zero setting (electronic background compensation) over entire selectable measurement range
<i>Maximum error</i>	± 0.6 % of full-scale range
<i>Linear range</i>	± 150 % of full-scale range

Conductivity detector

Construction Thermostatted conductivity detector with 2 ring-shaped steel electrodes

Measurement principle Alternating current measurement with following cell supply:

<i>Meas. Range</i>	<i>Amplitude</i>	<i>Frequency</i>
100 $\mu\text{S/cm}$	1 V peak	714 Hz
200 $\mu\text{S/cm}$	1 V peak	5000 Hz
500 $\mu\text{S/cm}$	1 V peak	5000 Hz
1 mS/cm	1 V peak	5000 Hz
2 mS/cm	0.5 V peak	5000 Hz
5 mS/cm	0.2 V peak	5000 Hz
10 mS/cm	0.1 V peak	5000 Hz

Effective cell volume 0.8 μL

Cell constant 16.7 /cm; adjustable 13.0...21.0 /cm

Maximum back pressure for measuring cell 5.0 MPa (50 bar)

Thermostating Connectable dynamic control to adjustable operating temperature

Operating temperature Adjustable in steps of 5°C from 25...45°C

Max. temperature deviation $\pm 1.5^\circ\text{C}$

Temperature stability $\leq 0.01^\circ\text{C}$ at constant ambient temperature

Parts and controls

Indicator Power LED **1**

ON/Off switch At rear of instrument

Mains connection

Voltage 115 V: 100...120 V $\pm 10\%$
230 V: 220...240 V $\pm 10\%$
Switching with mains voltage selector in fuse holder (see Section 2.4.1)

Frequency 50...60 Hz

Power consumption 70 VA

Fuse 5 mm \varnothing , 20 mm length
100...120 V: 0.63 A (slow-blow)
220...240 V: 0.315 A (slow-blow)

Interfaces

RS 232 interface Specifications, see Section 5.1

Remote interfaces Specifications, see Section 5.2

Analog output Specifications, see Section 5.3

Safety specifications

<i>Construction/testing</i>	According to EN/IEC 61010-1 / UL 3101-1, protection class I
<i>Safety directions</i>	The Instructions for Use include information and warnings which must be heeded by the user to assure safe operation of the instrument.

Electromagnetic compatibility (EMC)

<i>Emitted interference</i>	Standards met: - EN/IEC 61326 - EN 55022 - CISPR 22 - EN/IEC 61000-3-2
<i>Immunity to interference</i>	Standards met: - EN/IEC 61326 - EN/IEC 61000-4-2 - EN/IEC 61000-4-3 - EN/IEC 61000-4-4 - EN/IEC 61000-4-5 - EN/IEC 61000-4-6 - EN/IEC 61000-4-8 - EN/IEC 61000-4-11 - EN/IEC 61000-4-14 - NAMUR

Ambient temperature

<i>Nominal operating range</i>	+5...+45°C (at 20...80 % atmospheric humidity)
<i>Storage, transport</i>	-20...+70°C

Reference conditions

<i>Ambient temperature</i>	+25°C (± 3 °C)
<i>Rel. humidity</i>	$\leq 60\%$
<i>Instrument condition</i>	Instrument in operation for at least 30 min

Diagnostic test

<i>Self-diagnostic test</i>	Automatic self-diagnosis when instrument switched on
-----------------------------	--

Housing

<i>Material of cover</i>	Polyurethane rigid foam (PUR) with fire protection for fire class UL94VO, FCH-free
<i>Material of base</i>	Steel, enameled
<i>Width</i>	255 mm

<i>Height</i>	130 mm
<i>Depth</i>	343 mm
<i>Weight</i>	8.2 kg (with detector block, without accessories)

6.1.2 820 IC Separation Center

Temperature control of the column heating

<i>Temperature range</i>	Ambient temperature +5°C...80°C, adjustable in steps of 0.1 °C.
<i>Accuracy</i>	± 0.5 °C
<i>Stability</i>	better than 0.05 °C

Parts and controls

<i>Indicators</i>	Power LED 11 for Separation Center Power LED 12 for Column heating
-------------------	---

Power supply of separation center

<i>Internal</i>	Internal supply of 819 IC Detector via connecting cable
<i>External</i>	External supply via DIN socket with 6.2152.000 power supply unit: 5 V / 0.5 A 24 V / 2 A (transient, 200...300 ms)

Power supply of column heating

Via jack bush with 6.2152.010 power supply unit:
24 Vdc / 2 A

Interfaces

<i>Valve interfaces</i>	Specifications, see <i>Section 5.5</i>
-------------------------	--

Safety specifications

<i>Construction/testing</i>	According to EN/IEC 61010-1 / UL 3101-1, Protection class III
<i>Safety notes</i>	The Instructions for Use include information and warnings which must be heeded by the user to assure safe operation of the instrument.

Electromagnetic compatibility (EMC)

<i>Emitted interference</i>	Standards met: - EN/IEC 61326 - EN 55022 - CISPR 22
-----------------------------	--

<i>Immunity to interference</i>	Standards met: - EN/IEC 61326 - EN/IEC 61000-4-2 - EN/IEC 61000-4-3 - EN/IEC 61000-4-4 - EN/IEC 61000-4-5 - EN/IEC 61000-4-6 - EN/IEC 61000-4-8 - EN/IEC 61000-4-11 - EN/IEC 61000-4-14 - NAMUR
---------------------------------	---

Ambient temperature

<i>Nominal operating range</i>	+5...+45°C (at 20...80 % atmospheric humidity)
<i>Storage</i>	-20...+65°C
<i>Transport</i>	-20...+70°C

Reference conditions

<i>Ambient temperature</i>	+25°C (± 3 °C)
<i>Rel. humidity</i>	≤ 60%
<i>Instrument condition</i>	Instrument in operation for at least 30 min

Housing

<i>Material</i>	Polyurethane rigid foam (PUR) with fire protection for fire class UL94VO, FCH-free
<i>Width</i>	255 mm
<i>Height</i>	385 mm
<i>Depth</i>	343 mm
<i>Weight</i>	6.4 kg (without accessories)

6.2 Scope of delivery



Subject to changes !
All dimensions are given in mm.




6.2.1 819 IC Detector

The 819 IC Detector is available in two versions:




- **2.819.0010** IC Detector with metal detector block
- **2.819.0110** IC Detector with metal-free detector block

These instruments include the following parts:

2.819.0010 IC Detector

Quant.	Order-No.	Description												
1	1.819.0010	819 IC Detector												
1	1.732.0100	Detector block with permanently attached connecting cable to 819 IC Detector												
1	6.2125.110	Connecting cable Connecting cable 819 IC Detector (RS 232) – PC 												
1	6.2125.090	Connecting cable Connecting cable 819 IC Detector – 820 IC Separation Center 												
1	6.2248.0000	Magnetic plate for identification of the 819 IC Detec- tor 												
1	6.2122.0X0	Mains cable to customer's specifications: <table style="width: 100%; border: none;"> <tr> <td style="border: none;"><u>Cable socket</u></td> <td style="border: none;"><u>Cable connector</u></td> <td style="border: none;"></td> </tr> <tr> <td style="border: none;">Type IEC 320/C 13</td> <td style="border: none;">Type SEV 12 (CH...)</td> <td style="border: none; text-align: right;">6.2122.020</td> </tr> <tr> <td style="border: none;">Type IEC 320/C 13</td> <td style="border: none;">Type CEE (7), VII (D...)</td> <td style="border: none; text-align: right;">6.2122.040</td> </tr> <tr> <td style="border: none;">Type CEE (22), V</td> <td style="border: none;">Type NEMA 5-15 (USA...)</td> <td style="border: none; text-align: right;">6.2122.070</td> </tr> </table>	<u>Cable socket</u>	<u>Cable connector</u>		Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020	Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040	Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070
<u>Cable socket</u>	<u>Cable connector</u>													
Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020												
Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040												
Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070												
1	8.819.1013	Instructions for Use (English) for 819 IC Detector and 820 IC Separation Center												

2.819.0110 IC Detector

Quant.	Order-No.	Description												
1	1.819.0010	819 IC Detector												
1	1.732.0110	Detector block with permanently attached connecting cable to 819 IC Detector												
1	6.2125.110	Connecting cable Connecting cable 819 IC Detector (RS 232) – PC 												
1	6.2125.090	Connecting cable Connecting cable 819 IC Detector – 820 IC Separation Center 												
1	6.2248.0000	Magnetic plate for identification of the 819 IC Detec- tor 												
1	6.2122.0X0	Mains cable to customer's specifications: <table border="0"> <tr> <td><u>Cable socket</u></td> <td><u>Cable connector</u></td> <td></td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type SEV 12 (CH...)</td> <td>6.2122.020</td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type CEE (7), VII (D...)</td> <td>6.2122.040</td> </tr> <tr> <td>Type CEE (22), V</td> <td>Type NEMA 5-15 (USA...)</td> <td>6.2122.070</td> </tr> </table>	<u>Cable socket</u>	<u>Cable connector</u>		Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020	Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040	Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070
<u>Cable socket</u>	<u>Cable connector</u>													
Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020												
Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040												
Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070												
1	8.819.1013	Instructions for Use (English) for 819 IC Detector and 820 IC Separation Center												






6.2.2 820 IC Separation Center









The 820 IC Separation Center is available in the following versions:

- **2.820.0210** IC Separation Center with 1 injector for a one-channel system with column heating, metal-free
- **2.820.0220** IC Separation Center with 2 injectors for a two-channel system with column heating, metal-free
- **2.820.0230** IC Separation Center with 1 injector and 1 Metrohm Suppressor Module MSM for a one-channel system with column heating, metal-free
- **2.820.0310** IC Separation Center with 1 injector for a one-channel system, metal-free
- **2.820.0320** IC Separation Center with 2 injectors for a two-channel system, metal-free
- **2.820.0330** IC Separation Center with 1 injector and 1 Metrohm Suppressor Module MSM for a one-channel system; metal-free










These instruments include the following parts:





2.820.0210 IC Separation Center

Quant.	Order-No.	Description	
1	1.820.0210	820 IC Separation Center	
1	1.820.0500	IC Column heating	
1	6.2108.120	Connecting cable Connecting cable Column heating – 820 IC Separation Center	 0.8 m
1	6.2134.040	Connecting cable RS 232-connecting cable 820 IC Separation Center (Connection 22) – 830 IC Interface	 3m 9 p./f 9 p./f
1	6.1816.000	Silicone tubing Drain tube for inner compartment, length = 0.5 m	
1	6.1831.010	PEEK capillary Length = 3 m	 0.25 1/16"
1	6.2027.030	Column holder diameter d = 8.5 mm	 42





1	6.2027.040	Column holder diameter d = 11.3 mm													
1	6.2027.050	Column holder diameter d = 15.0 mm													
1	6.2744.010	PEEK compression fitting For the connection of 6.1831.010 PEEK capillaries, set of 5													
1	6.2744.020	Coupling 1/16" – Luer Coupling for connection of a 6.1803.000 PTFE capillary to connection 13 of the 820 IC Separation Center when a sample changer (766, 788, 813) is used													
2	6.2744.040	PEEK coupling for the connection of 1/16" capillaries													
1	6.2744.070	PEEK compression fitting-(short) for the connection of 6.1831.010 PEEK capillaries in cramped conditions, e.g. in the column heating, set of 5													
1	6.2816.020	Syringe made of PP, volume = 10 mL; for manual filling of the sample loop													
1	6.2152.010	Power supply unit for column heating													
1	6.2122.0X0	Mains cable to customer's specifications: <table border="0" style="width: 100%; border-collapse: collapse;"> <tr> <td style="width: 50%;"><u>Cable socket</u></td> <td style="width: 50%;"><u>Cable connector</u></td> <td></td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type SEV 12 (CH...)</td> <td>6.2122.020</td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type CEE (7), VII (D...)</td> <td>6.2122.040</td> </tr> <tr> <td>Type CEE (22), V</td> <td>Type NEMA 5-15 (USA...)</td> <td>6.2122.070</td> </tr> </table>	<u>Cable socket</u>	<u>Cable connector</u>		Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020	Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040	Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070	
<u>Cable socket</u>	<u>Cable connector</u>														
Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020													
Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040													
Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070													
1	8.819.1013	Instructions for Use (English) for 820 IC Separation Center and 819 IC Detector													










2.820.0220 IC Separation Center

Quant.	Order-No.	Description	
1	1.820.0220	820 IC Separation Center	
1	1.820.0500	IC Column heating	
1	6.2108.120	Connecting cable Connecting cable Column heating – 820 IC Separation Center	 0.8 m
1	6.2134.040	Connecting cable RS 232-connecting cable 820 IC Separation Center (Connection 22) – 830 IC Interface	 3m 9 p.f 9 p.f
1	6.1816.000	Silicone tubing Drain tube for inner compartment, length = 0.5 m	
1	6.1831.010	PEEK capillary Length = 3 m	
2	6.2027.030	Column holder diameter d = 8.5 mm	
2	6.2027.040	Column holder diameter d = 11.3 mm	
2	6.2027.050	Column holder diameter d = 15.0 mm	
2	6.2744.010	PEEK compression fitting For the connection of 6.1831.010 PEEK capillaries, set of 5	
2	6.2744.020	Coupling 1/16" – Luer Coupling for connection of a 6.1803.000 PTFE capillary to connec- tion 13 of the 820 IC Separation Cen- ter when a sample changer (766, 788, 813) is used	

4	6.2744.040	PEEK coupling for the connection of 1/16" capillaries	
1	6.2744.070	PEEK compression fitting-(short) for the connection of 6.1831.010 PEEK capillaries in cramped conditions, e.g. in the column heating, set of 5	
2	6.2816.020	Syringe made of PP, volume = 10 mL; for manual filling of the sample loop	
1	6.2152.010	Power supply unit for column heating	
1	6.2122.0X0	Mains cable to customer's specifications: <u>Cable socket</u> <u>Cable connector</u> Type IEC 320/C 13 Type SEV 12 (CH...)..... 6.2122.020 Type IEC 320/C 13 Type CEE (7), VII (D...)..... 6.2122.040 Type CEE (22), V Type NEMA 5-15 (USA...)..... 6.2122.070	
1	8.819.1013	Instructions for Use (English) for 820 IC Separation Center and 819 IC Detector	









2.820.0230 IC Separation Center



Quant.	Order-No.	Description	
1	1.820.0230	820 IC Separation Center	
1	1.820.0500	IC Column heating	
1	6.2108.120	Connecting cable Connecting cable Column heating – 820 IC Separation Center	
1	6.2134.040	Connecting cable RS 232-connecting cable 820 IC Separation Center (Connection 22) – 830 IC Interface	
1	6.1816.000	Silicone tubing Drain tube for inner compartment, length = 0.5 m	
1	6.1831.010	PEEK capillary Length = 3 m	

1	6.2027.030	Column holder diameter d = 8.5 mm	
1	6.2027.040	Column holder diameter d = 11.3 mm	
1	6.2027.050	Column holder diameter d = 15.0 mm	
2	6.2744.010	PEEK compression fitting For the connection of 6.1831.010 PEEK capillaries, set of 5	
1	6.2744.020	Coupling 1/16" - Luer Coupling for connection of a 6.1803.000 PTFE capillary to connection 13 of the 820 IC Separation Center when a sample changer (766, 788, 813) is used	
2	6.2744.040	PEEK coupling for the connection of 1/16" capillaries	
1	6.2744.070	PEEK compression fitting-(short) for the connection of 6.1831.010 PEEK capillaries in cramped conditions, e.g. in the column heating, set of 5	
1	6.2816.020	Syringe made of PP, volume = 10 mL; for manual filling of the sample loop	
1	6.2152.010	Power supply unit for column heating	
1	6.2122.0X0	Mains cable to customer's specifications: <u>Cable socket</u> Type IEC 320/C 13 Type IEC 320/C 13 Type CEE (22), V <u>Cable connector</u> Type SEV 12 (CH...)..... 6.2122.020 Type CEE (7), VII (D...)..... 6.2122.040 Type NEMA 5-15 (USA...)..... 6.2122.070	







1	8.819.1013	Instructions for Use (English) for 820 IC Separation Center and 819 IC Detector
---	------------	---





2.820.0310 IC Separation Center

Quant.	Order-No.	Description	
1	1.820.0310	820 IC Separation Center	
1	6.1816.000	Silicone tubing Drain tube for inner compartment, length = 0.5 m	
1	6.1831.010	PEEK capillary Length = 3 m	
1	6.2027.030	Column holder diameter d = 8.5 mm	
1	6.2027.040	Column holder diameter d = 11.3 mm	
1	6.2027.050	Column holder diameter d = 15.0 mm	
1	6.2744.010	PEEK compression fitting For the connection of 6.1831.010 PEEK capillaries, set of 5	
1	6.2744.020	Coupling 1/16" - Luer Coupling for connection of a 6.1803.000 PTFE capillary to connec- tion 13 of the 820 IC Separation Cen- ter when a sample changer (766, 788, 813) is used	
2	6.2744.040	PEEK coupling for the connection of 1/16" capillaries	





1	6.2744.070	PEEK compression fitting-(short) for the connection of 6.1831.010 PEEK capillaries in cramped conditions, e.g. in the column heating, set of 5	
1	6.2816.020	Syringe made of PP, volume = 10 mL; for manual filling of the sample loop	
1	8.819.1013	Instructions for Use (English) for 820 IC Separation Center and 819 IC Detector	







2.820.0320 IC Separation Center

Quant.	Order-No.	Description	
1	1.820.0320	820 IC Separation Center	
1	6.1816.000	Silicone tubing Drain tube for inner compartment, length = 0.5 m	
1	6.1831.010	PEEK capillary Length = 3 m	
2	6.2027.030	Column holder diameter d = 8.5 mm	
2	6.2027.040	Column holder diameter d = 11.3 mm	
2	6.2027.050	Column holder diameter d = 15.0 mm	
2	6.2744.010	PEEK compression fitting For the connection of 6.1831.010 PEEK capillaries, set of 5	

2	6.2744.020	Coupling 1/16" – Luer Coupling for connection of a 6.1803.000 PTFE capillary to connection 13 of the 820 IC Separation Center when a sample changer (766, 788, 813) is used	
4	6.2744.040	PEEK coupling for the connection of 1/16" capillaries	
1	6.2744.070	PEEK compression fitting-(short) for the connection of 6.1831.010 PEEK capillaries in cramped conditions, e.g. in the column heating, set of 5	
2	6.2816.020	Syringe made of PP, volume = 10 mL; for manual filling of the sample loop	
1	8.819.1013	Instructions for Use (English) for 820 IC Separation Center and 819 IC Detector	








2.820.0330 IC Separation Center



Quant.	Order-No.	Description	
1	1.820.0330	820 IC Separation Center	
1	6.1816.000	Silicone tubing Drain tube for inner compartment, length = 0.5 m	
1	6.1831.010	PEEK capillary Length = 3 m	
1	6.2027.030	Column holder diameter d = 8.5 mm	
1	6.2027.040	Column holder diameter d = 11.3 mm	

1	6.2027.050	Column holder diameter d = 15.0 mm	
2	6.2744.010	PEEK compression fitting For the connection of 6.1831.010 PEEK capillaries, set of 5	
1	6.2744.020	Coupling 1/16" - Luer Coupling for connection of a 6.1803.000 PTFE capillary to connection 13 of the 820 IC Separation Center when a sample changer (766, 788, 813) is used	
2	6.2744.040	PEEK coupling for the connection of 1/16" capillaries	
1	6.2744.070	PEEK compression fitting-(short) for the connection of 6.1831.010 PEEK capillaries in cramped conditions, e.g. in the column heating, set of 5	
1	6.2816.020	Syringe made of PP, volume = 10 mL; for manual filling of the sample loop	
1	8.819.1013	Instructions for Use (English) for 820 IC Separation Center and 819 IC Detector	

6.3 Optional accessories

6.3.1 Accessories for 820 IC Separation Center

Order No.	Description	
6.1825.XXX	<p>Sample loop PEEK For injection valve; incl. 2 PEEK compression fittings 6.2744.010</p> <p>6.1825.230: Volume = 10 µL 6.1825.210: Volume = 20 µL 6.1825.220: Volume = 100 µL</p>	
6.2821.120	<p>Filter unit PEEK 2 µm To avoid contamination.</p> <p>Spare part: 6.2821.130 Filter</p>	
6.2821.130	<p>Filter for filter unit PEEK 2 µm Spare part for 6.2821.120 Filter unit PEEK Set of 10</p>	
6.2832.000	<p>Suppressor rotor Replacement cartridge for Metrohm suppressor module</p>	
6.2832.010	<p>Connection piece for suppressor rotor with input and output leads</p>	
6.2620.150	<p>Pulsation dampener Metal-free pulsation dampener for reduction of pulsation and preservation of separating columns.</p>	
6.2621.080	<p>Capillary tubing cutter for 6.1831.010 PEEK capillaries and 6.1822.010 PTFE micro capillaries incl. 5 additional blades</p>	
6.2744.160	<p>PEEK coupling with tubing security device for connection to pressure side of 6.1826.0X0 Pump tubing</p>	







6.2821.090	Aspirating filter PE 20 µm For 6.1834.000 Aspirating tubing (supplied with the 818 IC Pump); set of 5.													
6.5324.000	Bottle rack For holding 3 supply bottles with eluent, regeneration solution and rinsing solution, incl. accessories (bottles, siphon, etc.)													
6.2152.000	Power supply unit for 820 Separation Center (without mains cable, please order the cable suitable for your country separately)													
<table border="0"> <tr> <td><u>Cable socket</u></td> <td><u>Cable connector</u></td> <td></td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type SEV 12 (CH...)</td> <td>6.2122.020</td> </tr> <tr> <td>Type IEC 320/C 13</td> <td>Type CEE (7), VII (D...)</td> <td>6.2122.040</td> </tr> <tr> <td>Type CEE (22), V</td> <td>Type NEMA 5-15 (USA...)</td> <td>6.2122.070</td> </tr> </table>		<u>Cable socket</u>	<u>Cable connector</u>		Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020	Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040	Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070	
<u>Cable socket</u>	<u>Cable connector</u>													
Type IEC 320/C 13	Type SEV 12 (CH...)	6.2122.020												
Type IEC 320/C 13	Type CEE (7), VII (D...)	6.2122.040												
Type CEE (22), V	Type NEMA 5-15 (USA...)	6.2122.070												

6.3.2 Column heating

Order No.	Description
2.820.0500	Add-on kit Column Heating including power supply unit and all connection cables
6.5858.000	Conversion set for column heating

To upgrade an IC Separation Center without column heating (2.820.03X0) to one with column heating (2.820.02X0) you require the 2.820.0500 Add-on kit Column Heating and the 6.5858.000 conversion set. This conversion can only be carried out by Metrohm Service technicians.

6.3.3 Cable

Order No.	Description
6.2115.070	Remote Connection cable Connection cable 820 IC Separation Center – 830 IC Interface 
6.2128.100	Remote Connection cable Connection cable 820 IC Separation Center – 830 IC Interface for direct control of 820 
6.2134.130	Connection cable Connection cable 830 IC Interface – 819 IC Detector 
6.2134.140	Connection cable Connection cable 771 IC Compact Interface – 819 IC Detector and 820 IC Separation Center 
6.2134.080	RS 232 Connection cable Connection cable 830 IC Interface – 1 external device (818, 819, 766, etc.) 
6.2134.090	RS 232-Connection cable Connection cable 830 IC Interface – 2 external devices (818, 819, 766, etc.) 

6.3.4 Literature

Order No.	Description
8.732.2003	Monograph Ion Chromatography (English) Fundamentals of Ion Chromatography
8.792.5003	Monograph Practical Ion Chromatography (English)
8.025.5003	Monograph Sample Preparation Techniques for Ion Chromatography (English)

6.4 Validation / GLP

GLP (Good Laboratory Practice) requires, among other things, that the precision and correctness of analytical instruments is checked at regular intervals by using SOPs (**Standard Operating Procedures, SOP**). An example of such a standard operating procedure is available from Metrohm under the title «**Application Bulletin No. 277 – Validation of Metrohm Ion Chromatography Systems by using Standard Operating Procedures (SOP)**». This SOP can be adapted for your ion chromatography system and used for its validation.

819 IC Detector and 820 IC Separation Center must be included as a part of the whole ion chromatography system, whose most important components include the pumps, separation columns and evaluation system, in the all-embracing validation of the whole system.

Please contact your local Metrohm agency in order to receive support in validating your 819 IC Detector and 820 IC Separation Center. It can also provide you with validation documentation which will help you to carry out your installation qualification (IQ) and operational qualification (OQ).

Further information about QA, GLP and validation can also be found in the brochure «**Quality management with Metrohm**» which is also obtainable from your local Metrohm agency.

Checking the electronic and mechanical assemblies of Metrohm instruments can and should be undertaken within the framework of regular servicing by Metrohm technicians. All Metrohm instruments are equipped with start-up check routines which check that the relevant assemblies are functioning perfectly when the instrument is switched on. If no error message appears it can be assumed that the instrument is functioning properly.

819 IC Detector and 820 IC Separation Center also contains a built-in diagnosis program which allows the service technicians to check the functioning of particular assemblies should faults or malfunctions occur and to localize them.

6.5 Warranty and Conformity

6.5.1 Warranty

The warranty on our products is limited to defects that are traceable to material, construction or manufacturing error which occur within 12 months from the day of delivery. In this case, the defects will be rectified in our workshops free of charge. Transport costs are to be paid by the customer.

For day and night operation, the warranty is limited to 6 months.

Glass breakage in the case of electrodes or other parts is not covered by the warranty. Checks which are not a result of material or manufacturing faults are also charged during the warranty period. For parts of outside manufacture insofar as these constitute an appreciable part of our instrument, the warranty stipulations of the manufacturer in question apply.

With the regard to the guarantee of accuracy, the technical specifications in the instruction manual are authoritative.

Concerning defects in material, construction or design as well as the absence of guaranteed features, the orderer has no rights or claims except those mentioned above.

If damage of the packaging is evident on receipt of a consignment or if the goods show signs of transport damage after unpacking, the carrier must be informed immediately and a written damage report demanded. Lack of an official damage report releases Metrohm from any liability to pay compensation.

If any instruments and parts have to be returned, the original packaging should be used if at all possible. This applies above all to instruments, electrodes, burette cylinders and PTFE pistons. Before embedment in wood shavings or similar material, the parts must be packed in a dust-proof package (for instruments, use of a plastic bag is imperative). If open assemblies are enclosed in the scope of delivery that are sensitive to electromagnetic voltages (e.g. data interfaces etc.) these must be returned in the associated original protective packaging (e.g. conductive protective bag). (Exception: assemblies with built-in voltage source belong in a non-conductive protective packaging).



No warranty responsibility whatsoever will be accepted by Metrohm for damage which arises as a result of non-compliance with these instructions.

6.5.2 Declaration of Conformity



This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

CH-9101 Herisau, Switzerland
 Tel. +41 71 353 85 85
 Fax +41 71 353 89 01
 www.metrohm.com



<i>Name of commodity</i>	819 IC Detector
<i>System software</i>	Stored in ROMs
<i>Name of manufacturer</i>	Metrohm Ltd., Herisau, Switzerland
<i>Description</i>	Thermostatted conductivity detector for ion chromatographic analyses with electronic or chemical suppression.
<p>This Metrohm instrument has been built and has undergone final type testing according to the standards:</p> <p><i>Electromagnetic compatibility: Emission</i> EN/IEC 61326, EN 55022 / CISPR 22, EN/IEC 61000-3-2</p> <p><i>Electromagnetic compatibility: Immunity</i> EN/IEC 61326, EN/IEC 61000-4-2, EN/IEC 61000-4-3, EN/IEC 61000-4-4, EN/IEC 61000-4-5, EN/IEC 61000-4-6, EN/IEC 61000-4-8, EN/IEC 61000-4-11, EN/IEC 61000-4-14, NAMUR</p> <p><i>Safety specifications</i> EN/IEC 61010-1, UL 3101-1 protection class I</p> <p>It has also been certified by ElectroSuisse, which is member of the International Certification Body (CB/IEC).</p> <p><i>The instrument meets the requirements of the CE mark as contained in the EU directives 89/336/EEC and 73/23/EEC and fulfils the following specifications:</i></p>	
EN 61326	Electrical equipment for measurement, control and laboratory use – EMC requirements
EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use
<p>Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.</p>	
<p>The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance.</p> <p>The technical specifications are documented in the instruction manual.</p>	
<p>Herisau, May 28, 2003</p> <div style="display: flex; justify-content: space-around; align-items: flex-end;"> <div style="text-align: center;">  <p>Dr. J. Frank Vice President Head of R&D</p> </div> <div style="text-align: center;">  <p>Ch. Buchmann Vice President Head of Production Responsible for Quality Assurance</p> </div> </div>	

6.5.3 Declaration of Conformity



This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

CH-9101 Herisau, Switzerland
 Tel. +41 71 353 85 85
 Fax +41 71 353 89 01
 www.metrohm.com

<i>Name of commodity</i>	820 IC Separation Center
<i>System software</i>	Stored in ROMs
<i>Name of manufacturer</i>	Metrohm Ltd., Herisau, Switzerland
<i>Description</i>	Thermally and electronically isolated wet room for ion chromatographic analyses, optionally with thermostatted column heating.
<p>This Metrohm instrument has been built and has undergone final type testing according to the standards:</p> <p><i>Electromagnetic compatibility: Emission</i> EN/IEC 61326, EN 55022 / CISPR 22</p> <p><i>Electromagnetic compatibility: Immunity</i> EN/IEC 61326, EN/IEC 61000-4-2, EN/IEC 61000-4-3, EN/IEC 61000-4-4, EN/IEC 61000-4-5, EN/IEC 61000-4-6, EN/IEC 61000-4-8, EN/IEC 61000-4-11, EN/IEC 61000-4-14, NAMUR</p> <p><i>Safety specifications</i> EN/IEC 61010-1, UL 3101-1 protection class III</p> <p>It has also been certified by ElectroSuisse, which is member of the International Certification Body (CB/IEC).</p> <p><i>The instrument meets the requirements of the CE mark as contained in the EU directives 89/336/EEC and 73/23/EEC and fulfils the following specifications:</i></p>	
EN 61326	Electrical equipment for measurement, control and laboratory use – EMC requirements
EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use
<p>Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.</p>	
<p>The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance.</p> <p>The technical specifications are documented in the instruction manual.</p>	
<p>Herisau, May 28, 2003</p> <div style="display: flex; justify-content: space-around; align-items: flex-end;"> <div style="text-align: center;">  <p>Dr. J. Frank Vice President Head of R&D</p> </div> <div style="text-align: center;">  <p>Ch. Buchmann Vice President Head of Production Responsible for Quality Assurance</p> </div> </div>	

6.5.4 Quality Management Principles

Metrohm Ltd., CH-9101 Herisau, Switzerland



Metrohm
Ion analysis
CH-9101 Herisau/Switzerland
E-Mail info@metrohm.com
Internet www.metrohm.com

Metrohm Ltd. holds the ISO 9001 Certificate, registration number 10872-02, issued by SQS (Swiss Association for Quality and Management Systems). Internal and external audits are carried out periodically to assure that the standards defined by Metrohm's QM Manual are maintained.

The steps involved in the design, manufacture and servicing of instruments are fully documented and the resulting reports are archived for ten years. The development of software for PCs and instruments is also duly documented and the documents and source codes are archived. Both remain the possession of Metrohm. A non-disclosure agreement may be asked to be provided by those requiring access to them.

The implementation of the ISO 9001 quality system is described in Metrohm's QM Manual, which comprises detailed instructions on the following fields of activity:

Instrument development

The organisation of the instrument design, its planning and the intermediate controls are fully documented and traceable. Laboratory testing accompanies all phases of instrument development.

Software development

Software development occurs in terms of the software life cycle. Tests are performed to detect programming errors and to assess the program's functionality in a laboratory environment.

Components

All components used in the Metrohm instruments have to satisfy the quality standards that are defined and implemented for our products. Suppliers of components are audited by Metrohm as the need arises.

Manufacture

The measures put into practice in the production of our instruments guarantee a constant quality standard. Production planning and manufacturing procedures, maintenance of production means and testing of components, intermediate and finished products are prescribed.

Customer support and service

Customer support involves all phases of instrument acquisition and use by the customer, i.e. consulting to define the adequate equipment for the analytical problem at hand, delivery of the equipment, user manuals, training, after-sales service and processing of customer complaints. The Metrohm service organisation is equipped to support customers in implementing standards such as GLP, GMP, ISO 900X, in performing Operational Qualification and Performance Verification of the system components or in carrying out the System Validation for the quantitative determination of a substance in a given matrix.

6.6 Index

8

818 IC Pump	
Connection at 819	29
Electrical connection	29
Maintenance	67
Practical notes	67
819 IC Detector	
Connection to 830	20
Control via IC Net	50
Description	1
Icon	50
Links	59
Mains connection	25
Method parameter	53
Opening	12
Ordering designation	92, 93
Parameter display	50
Parts and controls	4
program	56
program instructions	57
Scope of delivery	92
Technical data	87
Window	50
820 IC Separation Center	
Connection to 818	29
Connection to 830	20
Connection to detector block	15
Control via IC Net	60
Description	2
External power supply	85
Icon	60
Interior	39, 41, 42, 47
links	65
Mains connection	27
Manual	61
Opening	12
Optional accessories	103
Ordering designation	94, 96, 97, 99, 100, 101
Parts and controls	6
program	63
program instructions	64
Scope of delivery	94
Technical data	90
Thermostat	64
Valve interfaces	86
via 819	60
via 830	60
Window	60
820 IC Separation Center	
Scope of delivery (2.820.0210)	94
Scope of delivery (2.820.0220)	96
Scope of delivery (2.820.0230)	97
Scope of delivery (2.820.0310)	99
Scope of delivery (2.820.0320)	100
Scope of delivery (2.820.0330)	101
830 IC Interface	
Connection to 819	20
Connection to 820	20
833 IC Liquid Handling Pump Unit	
Preparation	45
833 IC Liquid Handling Unit	
General information	37
A	
Abrasive particles	67

abs. Cond.	51, 52, 53
Accuracy	87, 90
<Activate>	59
Actual temperature	64
<Add>	57, 59, 63
Advance pulse	82
Alternating current measurement	88
Ambient temperature	89, 91
Analog output 5	
Figure	5
Signal output	54
Switching	85
Appendix	87
Application Bulletins	10, 34, 66
Application Notes	10, 34, 66
Arrangement of the instruments	13
Aspirating tube 63	
Figure	40, 43
Asprating filter	
6.2821.090	66, 67, 104
Aspirating tubing	19
Assembly	45, 46
Attention	11
Auto Step	62
Autozero	87
Auto-zero function	52, 54
B	
Back pressure	88
Background compensation	see <i>Auto-zero</i>
Baseline with high noise level	77
Baud rate	79
Block diagram	3
Blockage	73, 75
Bottle rack (6.5324.000)	
Connection	19, 104
C	
Cable	105
6.2108.120	94, 96, 97
6.2125.070	105
6.2125.090	92, 93
6.2125.110	92, 93
6.2128.100	105
6.2134.000	105
6.2134.040	94, 96, 97
6.2134.080	29, 105
6.2134.090	29
6.2134.100	23
6.2134.130	105
6.2134.140	105
Capacity factor	66
Capacity problems	72
Capillary 52	
Figure	35
Capillary connections	28
Capillary for syringe 60	
Figure	40, 43
Capillary prewarming 55	
Figure	35

Capillary tubing cutter	
(6.2621.080)	28, 103
Care	69
Cell const.	51
Cell constant	55, 88
Cell supply	88
Cell volume	88
Changing of separating column	71
Changing the separation system	71
Check	13
Chemical suppression	37
Cleaning the suppressor	73
Column	see <i>Separating column</i>
Column connection capillary 50	
Changing separating column	71
Figure	31, 35, 39, 43
Mounting	32, 33, 36, 40, 43
Column heating	64, 104
Add-on kit	104
Connection power supply 24	8
Connection socket 23	8
Connection socket 38	9
Conversion set	104
Figure	9
Install column	35
Install in Separation Center	18
Mains connection	27
Parts and controls	9
RS 232 interface 22	8
Technical data	90
Column holder 59	
Figure	40, 43
Ordering designation	94, 95, 96, 98, 99, 100, 101, 102
Communication	1
Compression fitting	
6.2744.070	36, 95, 97, 98, 100, 101, 102
6.2744.010	28, 30
Mounting	28
Ordering designation	95, 96, 98, 99, 100, 102
Conditioning	48
Conductivity	51, 52, 53, 54, 55
Conductivity detector	88
Conductivity measurement	87
Configuration	61
Connect suppressor module	45
Connection	
818 IC Pump	29
Column heating	18
Suppressor module	44
Connection 10	
Figure	5
Connection 13	
Connection of syringe	19
Figure	6
Connection 22	
Figure	8
Connection 23	
Figure	8
Connection 24	
Figure	8

Connection 25		Detector block 56		Ordering designation	103
Figure	8	Connect column	40	Filter unit PEEK 6.2821.120	
Pin assignment	85	Connection at 819	16, 17	47	31
Connection 26		Figure	40, 43	Connection to 833 Pump Unit....	45, 46
Figure	8	Insert	15	Connector	30
Connection 27		Ordering designation.....	92, 93	Figure	30
Figure	8	Detector block connection 6		Filter	30
Connection 30		Figure	5	Housing	30
Figure	8	Dimensions	89, 91	Mounting.....	30
Connection 31		Disposal	12	Front	
Figure	8	Double peaks	67, 77	819	4
Connection 33		Drain tube	19	820	6
Connection of drain tube	19	Connection 33	8	FS	57
Figure	8	Drift	68, 77, 87	Full duplex	79
Connection 38		E		Full scale	49, 51, 54
Figure	9	Earthing	12	Full-scale range	49
Connection 5		Electrical connection	20	Fuse	25, 88
Figure	5	Electromagnetic compatibility	89, 90	Fuse holder 4	
Connection 6		Eluent change	69	Figure	5, 26
Figure	5	Eluents	68	Replacement of fuses	25
Connection 820		EMC	89, 90	G	
Figure	7	Emitted interference	89, 90	General information	
Connection 9		ENABLED	57, 63	on maintenance	69
Figure	5	EOD	82	General information	
Connection capillary 3		Error signal	82	on separating columns	34
Connection 818 – 820	31	Evaporation	68, 77	General information	
Mounting	31	External device	81	on suppressor module	37
Connection of 818 IC Pump	29	External devices	83, 86	General precautionary rules	12
Connection of syringe and aspirating		External power supply		Glass beads	67
tubing	19	Connection 25	8	GLP	106
Connection of waste container	16, 17	Connection 25	8	H	
Connection piece 73		External power supply for 820	85	Handling of solvents	12
Cleaning	73	F		Heater cover 37	
Figure	74	Faults	66, 77	Figure	9
Ordering designation	103	Feedthrough 14		Figure	35
Connections	69	Figure	6, 8	Heating	55
Connections 819		Mounting of		Holder plate 55	
Figure	5	suppressor inlet capillary	45, 46	Figure	35
Connections 820		Feedthrough 16		Housing	89, 91
Figure	6	for detector cable	8	I	
Connectors for capillaries	28	Feedthrough 17		IC cation column	
Contaminated valves	67	Waste	8	METROSEP Cation 1-2	70
Contamination with heavy metals	72	Feedthrough 18	8	IC Column heating	94, 96, 97
Contamination with organic substances		Feedthrough 19		IC Net	
Contents	1	Waste	8	no connection	78
Control characters	79	Feedthrough 20		IC separating columns	see
Conversion kit column heating	104	for connecting cable to column		<i>Separating columns</i>	
Coupling		heating	8	IC system	
6.2744.020	95, 96, 98, 99, 101, 102	for detector cable	8	Arrangement	3
6.2744.040	95, 97, 98, 99, 101, 102	Feedthrough 32		Identical separation system	71
Ordering designation	103	Inlet B	8	Immunity to interference	89, 91
Coupling 71		Feedthrough 34		Increase in the baseline	72
6.2744.040	70, 71, 72	Inlet A	8	Information	11
Figure	43	Feedthrough 35		Information on the	
Mounting	45	Figure	9, 35	Instructions for Use	10
D		Feedthrough 41		Injection	62
Damping	51, 56, 87	Figure	9, 35	Injection valve	
Danger	11	Figures		Flow path	62
Data transmission protocol	79	list	IV	Injection valve 51	
Dead volume	77	Filling of sample loop	62	Connection of 818	32
Declaration of Conformity		Filter 45		Connection to 818	32
819	108	Ordering designation	103	Figure	31, 39, 43
820	109	Filter unit PEEK 47		Inlet capillary 57	
Degassing of eluents	68	Connection	31	Changing separating column	71
<Delete>	57, 59, 63	Connection to Pump Unit 833....	45, 46		
Detector		Mounting	30		
Connection 27	8	Notes	66, 69		
Connection 30	8				

Figure..... 40, 43
 Mounting..... 40
Inlet capillary 62
 Figure..... 40, 43
 Installation 13
 Accessories 15
 Installing the filter unit PEEK 30
 Instructions for Use
 819 IC Detector 92, 93
 820 IC Separation Center 95, 97,
 99, 100, 101, 102
 Instructions for Use (8.819.1013)..... 10
 Instrument description 1
 Integrator start 86
 Interfaces 79
 Interior of Separation Center
 2.820.0310..... 39
 2.820.0320..... 41
 2.820.0330..... 42
 Interior Separation Center
 2.820.0320..... 47
 Introduction 1
Ion..... 51, 55
 Ion chromatography..... 66
 Ion chromatography system..... 3
 IQ 106
 ISO 9001..... 108, 109

K

Knurled screw 15
 Figure..... 8
 Open rear panel 15, 16, 18, 32
Knurled screw 36
 Figure..... 9, 35

L

Leak testing 48
 Leaks 12, 13, 77
 Linearity 54, 87
Links..... 59, 65
 List of figures IV
 List of program instructions 57, 64
 Literature..... 105
 Location 13
 Loss of suppressor capacity..... 75

M

Magnetic holder 39
 Figure..... 9
Magnetic plate (6.2248.000)
 Installation 14
 Ordering designation 92, 93
Mains cable 26
 Connection 26
 Ordering designation 92, 93,
 95, 97, 98
Mains connection 12, 25
 Technical data..... 88
Mains connection plug 3
 Figure..... 5, 26
 Mains connection 26
Mains frequency..... 88
Mains pilot lamp 1
 Figure..... 4
 Operating display 26
Mains pilot lamp 10
 Figure..... 6

Mains switch 2
 Figure 5, 26
 On/off switching of the instrument..... 26
Mains voltage
 Setting 25
 Technical data..... 88
Maintenance..... 66, 69
Malfunctions..... 77
Mark..... 57
Maximum back pressure 88
Measuring cell temperature 55
Measuring range 49
Metrohm service 70
Metrohm suppressor module see
Suppressor module
Microfiltration..... 66, 67, 68
Model plate 34
 Figure 8
Model plate 40
 Figure 9
Model plate 7
 Figure 5
Monograph..... 10
 Ion Chromatography
 8.732.2003 105
 Practical Ion Chromatography
 8.792.5003 105
 Sample Preparation Techniques
 for Ion Chromatography
 8.025.5003 105
Mounting
 Suppressor..... 74
Mounting rail 58
 Figure 40, 43
 Installation of column holder 40, 44
MSM see *Suppressor module*

N

No connection 78
 Noise 56, 87
 Notation 11
 Notes 66

O

Obstruction..... 73
<Off> 65
Offset 51, 55
 On/off switching of the instrument..... 26
<On> 65
 One-channel system..... 3, 39, 42
Opening 16
 Figure 8
Opening 17
 Figure 8
Opening 18
 Figure 8
Opening 19
 Figure 8
Opening 20
 Figure 8
Opening 32
 Insert capillary 32
Opening 32
 Figure 8
Opening 34
 Insert capillary 32
Opening 34
 Figure 8

Opening the instrument 12
 Operating temperature 88
 Operation
 Via IC Net 50
 Optional accessories 103
 OQ 106
 Organic modifiers 68

P

Packaging 13
Parts and controls..... 4, 88
 820 6
Parts and Controls
 819 4
Passivation..... 70
PC connection 9
 Figure 5
Peak broadening 67, 77
PEEK capillary
 6.1831.010 94, 96, 97, 99, 100, 101
PEEK capillary 3
 6.1831.010 31
PEEK capillary 42
 6.1831.010 28, 30, 31, 39, 43
PEEK capillary 48
 6.1831.010 31, 39, 43
PEEK capillary 42
 6.1831.010 28
 pH..... 68
Pictograms 11
Pin assignment 80
Plate number 66
Polarity..... 54, 55, 57, 85
Power consumption 88
Power Supply 6.2152.000
 Separation Center 27
Power supply 6.2152.010
 Column heating 27
Power supply unit
 6.2152.000 104
 6.2152.010 95, 97, 98
Power supply unit 6.2152.000
 of separation center..... 8
Power supply unit 6.2152.010
 of column heating..... 8
Practical notes 66
Precautionary rules 12
Precipitates 68, 69
Precolumns
 Installation 33
 Notes..... 66
Pressure drop 77
Pressure rise 77
Program
 819 time program 56
 820 time program 63
Program# 56, 63
Protection class 12, 26, 89, 90
Protective earth 26
PTFE aspirating tubing 63
 Figure 40, 43
 Lead to the outside 19
PTFE tubing 16
 Mounting 45, 46
Pulsation dampener 49
 6.2620.150 31, 39, 43
 Connection..... 29, 32
 Ordering designation..... 103

Practical notes	67	Screw nut 71		Connection 1	45
Pulsations.....	56, 77	Figure.....	74	Connection 2	45
Pump.....	<i>see 818 IC Pump</i>	Selection of the sample loop.....	37	Connection 3	46
Pump tubing (6.1826.050)		Selectivity	66	Connection of column.....	42
General information	37	Separating column		Connections	44
Q		Changing	71	Figure.....	43
QA	106	Characteristic data	66	General information.....	37
Quality assurance	106	Dead volume	67	Mounting.....	44
R		Fixing.....	40, 44	Protection	69
Range	49, 51, 52, 53, 57	General information.....	34	Regeneration	72
Ready signal	82	Install in column heating.....	35	Replacement	75
Rear		Installation.....	40	Suppressor outlet capillary 66	
819	5	Mounting.....	42	Figure.....	43, 44
820	7	Practical notes	66	Mounting.....	45, 46
Rear panel 21		Protection.....	66	Suppressor outlet capillary 67	
Close rear panel	15, 16	Regeneration	67	Figure.....	43, 44
Figure	8	Rinsing.....	40, 43	Suppressor outlet capillary 69	
Replacing rear panel	18	Separation efficiency.....	66	Figure.....	43, 44
Rear panel 28		Storage	67	Mounting.....	45
Figure	8	Separating column 53		Suppressor rotor 74	
Open rear panel.....	32	Changing	71	Figure.....	74
Rectification of malfunctions	77	Figure.....	35, 40, 43	Ordering designation	103
Recycling	68, 70	Mounting.....	40, 43	Syringe	
Reference conditions.....	89, 91	Separation Center connection 10		Connection	6
Reference temperature.....	87	Figure.....	5	Syringe (6.2816.020)	
Regeneration.....	67, 72	Serial number 34		Connection of syringe	19
Regeneration solution.....	37, 45, 46, 72	Figure.....	8	Ordering designation	95, 97, 98, 100, 101, 102
Remote	57	Serial number 40		T	
Remote configuration	58	Figure.....	9	Technical data	87
Remote input lines		Serial number 7		Temperature	
Functions	81	Figure.....	5	of column heating	64
Remote interface 10		Service	70	of thermostat	64
Description.....	83	Servicing	69	Temperature coeff.	55
Figure	5	Set point	64	Temperature coefficient	87
Pin assignment.....	83	<Set>	65	Temperature control	
Remote interface 8		<SET>	56	Column heating.....	90
Description.....	81	Setting up the instrument	13	Temperature correction.....	87
Pin assignment	81	Shutdown.....	70	Temperature dependence	87
Remote lines	51, 56	Silicone tubing		Temperature deviation	88
Remote output lines		6.1816.000.....	94, 96, 97, 99,	Temperature range.....	90
Functions	82	100, 101		Temperature stability.....	55, 88
Replacing the suppressor	75	Single-channel system	3, 40	Terminal block 26	
Resolution	66, 77	Software installation.....	27	Pin assignment.....	86
Rinse with dist. H ₂ O.....	71	SOP.....	106	Terminal block 31	
Rinse with eluent.....	71	Spilled chemicals	69	Pin assignment.....	86
Rinsing of IC system.....	70	Splitting	67, 77	Test chromatogram.....	34
Rinsing solution	37	Stability.....	90	Thermostat	51, 55, 64
RS 232 interface 22		Standard chromatogram.....	66	actual temperature.....	64
Figure	8	Standard operating procedures.....	106	operating temperature	64
RS 232 interface 9		Start with startup hardware		Thermostat lamp 12	
Baud rate	79	Thermostat.....	65	Figure.....	6
Data transmission protocol	79	Static charges.....	12	Thermostatting	88
Description.....	79	Step	64	Transport	13
Figure	5	Stirring of the eluent.....	68	Transport damage.....	107
Pin assignment	80	Sunlight	13	Treatment of eluents	68
S		Suppression.....	37	Triggering of functions	81
Safety notes	12	Suppressor	58, 61	Tubing connections	
Safety specifications.....	89, 90	Suppressor holder 75		Connection to injection valve.....	31
Sample loop 61		Figure.....	74	Two-channel system	3, 41, 47
Figure	40, 43	Suppressor inlet capillary 64		V	
Fill.....	62	Figure.....	43, 44	Validation	106
Ordering designation.....	103	Mounting.....	44	Valve A.....	61
Selection	37	Suppressor inlet capillary 65		Terminal block 31	8
Scope of delivery	92	Figure.....	43, 44	Valve B.....	61
		Mounting.....	45	Terminal block 26	8
		Suppressor inlet capillary 68		Valve interfaces	86
		Figure.....	43, 44		
		Mounting.....	46		
		Suppressor module 70			
		Cleaning.....	73		

ValveA	58, 64
ValveB	58, 64
<Verify>	57, 63

W

Warning	11
Warranty	107
Waste container	16, 17

Z

Zero	57
Zero <OFF>	52
Zero <ON>	52