
822 Titration Curve Simulator

Mode d'emploi
8.822.1012

Metrohm AG

CH-9100 Herisau

Suisse

Téléphone : +41 71 353 85 85

Fax : +41 71 353 89 01

info@metrohm.com

www.metrohm.com

822 Titration Curve Simulator

Mode d'emploi

Teachware
Metrohm AG
CH-9100 Herisau
teachware@metrohm.com

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SER_OUT - bit 0:	BUSY	Status message for active cock switching
SER_OUT - bit 1:	STOP	Status message if difference counter = zero
SER_OUT - bit 2:	ESU	End switch down: cylinder at minimal stroke
SER_OUT - bit 3:	EXPO	Bit 0 for coding the cylinder volume
SER_OUT - bit 4:	EXP1	Bit 1 for coding the cylinder volume
SER_OUT - bit 5:	EXP2	Bit 2 for coding the cylinder volume
SER_OUT - bit 6:	LOCK	Exchange unit present
SER_OUT - bit 7:	DOS_TYP1	Bit 1 definition of the chip function (formerly "GTYP")

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SER_IN - bit 0:	POS0	Bit 0 for the default of the desired cock position
SER_IN - bit 1:	POS1	Bit 1 for the default of the desired cock position
SER_IN - bit 2:	DOWN	Dosing impulses for spindle motor
SER_IN - bit 3:	UP	Filling impulses for spindle motor
SER_IN - bit 4:	NRES	Reset for the Dosimat controller (L activ)
SER_IN - bit 5:	POS2	Bit 2 for the default of the desired cock position
SER_IN - bit 6:	R_NRW	Use data or set STROBE impulse
SER_IN - bit 7:	DIR	Definition of the direction of rotation of the cock motor

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SER_OUT - bit 0:	STOP	Status message if difference counter = zero
SER_OUT - bit 1:	ESU	End switch down: cylinder at minimal stroke
SER_OUT - bit 2:	HOTLOCK	Bit for prompting the "Hot Plug & Play" of the exchange unit
SER_OUT - bit 3:	ESO	End switch up: cylinder at maximum stroke
SER_OUT - bit 4:	LOCK	Exchange unit present
SER_OUT - bit 5:	MOT_DO / IN0	Bit 0 D/A converter or digital input IN0
SER_OUT - bit 6:	MOT_D1 / IN1	Bit 1 D/A converter or digital input IN1
SER_OUT - bit 7:	MOT_D2 / IN2	Bit 2 D/A converter or digital input IN2
SER_OUT - bit 8:	MOT_D3 / NHWR	Bit 3 D/A converter or status message of the reset signal
SER_OUT - bit 9:	MOT_D4 / NPOR	Bit 4 D/A converter or of the power-on-reset
SER_OUT - bit 10:	MOT_SN / OVERL	Algebraic sign of D/A converter or motor overload
SER_OUT - bit 11:	MOT_INFO	Definition of the signals SER_OUT bit 5 ... 10, 15 ... 17
	MOT_INFO = L:	Information about current speed of the spindle motor (MOT_D0...SN) and cock position (ISTPOS0...2)
	MOT_INFO = H:	Information of the digital inputs IN0...2, NHWR, NPOR, OVERL, EXP0...2
SER_OUT - bit 12:	BUSY	Status message for active cock switching
SER_OUT - bit 13:	STOP_0	Driver stage spindle or cock motor inactive
SER_OUT - bit 14:	ESM	End switch middle
SER_OUT - bit 15:	ISTPOS0 / EXP0	Bit 0 for cock position or coding the cylinder volume
SER_OUT - bit 16:	ISTPOS1 / EXP1	Bit 1 for cock position or coding the cylinder volume
SER_OUT - bit 17:	ISTPOS2 / EXP2	Bit 2 for cock position or coding the cylinder volume
SER_OUT - bit 18:	DOS_TYPO	Bit 0 definition of the chip function
SER_OUT - bit 19:	DOS_TYP1	Bit 1 definition of the chip function (formerly "GTYP")

SER_OUT - bit 20:	DOS_TYP2	Bit 2 definition of the chip function
SER_OUT - bit 21:	DOS_TYP3	Bit 3 definition of the chip function
SER_OUT - bit 22:	P10000	Bit 0 for stating the number of pulses for the total volume
SER_OUT - bit 23:	P20000	Bit 1 for stating the number of pulses for the total volume
SER_OUT - bit 24:	PWR	Enables "Hot Plug & Play" of the device and recognize temporarily unusual supply and control data

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SER_IN - bit 0:	UP	Filling impulses for spindle motor
SER_IN - bit 1:	DOWN	Dosing impulses for spindle motor
SER_IN - bit 2:	R_NRW	Use data or set STROBE impulse
SER_IN - bit 3:	MOT_INFO	Definition of the digital signals SER_OUT bit 5 ... 10, 15 ... 17
	MOT_INFO = L:	Information about current speed of the spindle motor (MOT_D0...SN) and cock position (ISTPOS0...2)
	MOT_INFO = H:	Information of the digital inputs IN0...2, NHWR, NPOR, OVERL, EXP0...2
SER_IN - bit 4:	NRES	Reset for dosimat controller (L activ)
SER_IN - bit 5:	OUT0	Bit 0 digital output
SER_IN - bit 6:	OUT1	Bit 1 digital output
SER_IN - bit 7:	OUT2	Bit 2 digital output
SER_IN - bit 8:	POS0	Bit 0 for the default of the desired cock position
SER_IN - bit 9:	POS1	Bit 1 for the default of the desired cock position
SER_IN - bit 10:	POS2	Bit 2 for the default of the desired cock position
SER_IN - bit 11:	DIR	Definition of the direction of rotation of the cock motor
SER_IN - bit 12:	STOP_HM	Signal to switch off the driver stage of the cock motor
SER_IN - bit 13:	STOP_SM	Signal to switch off the driver stage of the spindle motor
SER_IN - bit 14:	P10000	Bit 0 for stating the number of pulses for the total volume
SER_IN - bit 15:	P20000	Bit 1 for stating the number of pulses for the total volume

SER_IN - bit 16:	TITR_NLIQ	Definition of the mechanism of the motor controller
SER_IN - bit 17:	SET_HOTLOCK	Enables "Hot Plug & Play" of the exchange unit
SER_IN - bit 18:	EXT_D0	Bit 0 for the theoretical speed of the spindle motor
SER_IN - bit 19:	EXT_D1	Bit 1 for the theoretical speed of the spindle motor
SER_IN - bit 20:	EXT_D2	Bit 2 for the theoretical speed of the spindle motor
SER_IN - bit 21:	EXT_D3	Bit 3 for the theoretical speed of the spindle motor
SER_IN - bit 22:	EXT_D4	Bit 4 for the theoretical speed of the spindle motor
SER_IN - bit 23:	EXT_SN	SIGN for the specification of the direction of rotation of the spindle motor
SER_IN - bit 24:	SET_PWR	Enables "Hot Plug & Play" of the device and recognize temporally unusual supply and control data

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
4.5.2 Declaration of Conformity

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

<p><i>Name of commodity</i></p> <p>822 Titration Curve Simulator</p>	 <p>CH-9101 Herisau/Switzerland E-Mail info@metrohm.com www.metrohm.com</p>
<p><i>Description</i> GLP test device for system check. Different types of curves can be simulated for verification of dosing output, electrode input and titration software.</p>	
<p>This instrument has been built and has undergone final type testing according to the standards:</p> <p><i>Electromagnetic compatibility: Emission</i> EN/IEC 61326, EN 55022 / CISPR 22</p> <p><i>Electromagnetic compatibility: Immunity</i> EN/IEC 61326, EN/IEC 61000-4-2, EN/IEC 61000-4-3, EN/IEC 61000-4-4, EN/IEC 61000-4-5, EN/IEC 61000-4-6</p> <p><i>Safety specifications</i> EN/IEC 61010-1, UL3101-1</p>	
<div style="display: flex; align-items: flex-start;"> <div style="font-size: 2em; margin-right: 10px;">  </div> <div> <p><i>The instrument meets the requirements of the CE mark as contained in the EU directives 89/336/EEC and 73/23/EEC and fulfils the following specifications:</i></p> <p>EN 61326 Electrical equipment for measurement, control and laboratory use – EMC requirements</p> <p>EN 61010-1 Safety requirements for electrical equipment for measurement, control and laboratory use</p> </div> </div>	
<p>Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.</p>	
<p>The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance.</p> <p>The technical specifications are documented in the instruction manual.</p>	
<p>Herisau, February 28, 2002</p> <div style="display: flex; justify-content: space-around; margin-top: 20px;"> <div style="text-align: center;">  <p>Dr. J. Frank Vice President Head of R&D</p> </div> <div style="text-align: center;">  <p>Ch. Buchmann Vice President Head of Production Responsible for Quality Assurance</p> </div> </div>	

4.5.3 Quality Management Principles

Metrohm Ltd., CH-9101 Herisau, Switzerland

**Metrohm**
Ion analysis
CH-9101 Herisau/Switzerland
E-Mail info@metrohm.com
Internet www.metrohm.com

Metrohm Ltd. holds the ISO 9001 Certificate, registration number 10872-02, issued by SQS (Swiss Association for Quality and Management Systems). Internal and external audits are carried out periodically to assure that the standards defined by Metrohm's QM Manual are maintained.

The steps involved in the design, manufacture and servicing of instruments are fully documented and the resulting reports are archived for ten years. The development of software for PCs and instruments is also duly documented and the documents and source codes are archived. Both remain the possession of Metrohm. A non-disclosure agreement may be asked to be provided by those requiring access to them.

The implementation of the ISO 9001 quality system is described in Metrohm's QM Manual, which comprises detailed instructions on the following fields of activity:

Instrument development

The organization of the instrument design, its planning and the intermediate controls are fully documented and traceable. Laboratory testing accompanies all phases of instrument development.

Software development

Software development occurs in terms of the software life cycle. Tests are performed to detect programming errors and to assess the program's functionality in a laboratory environment.

Components

All components used in the Metrohm instruments have to satisfy the quality standards that are defined and implemented for our products. Suppliers of components are audited by Metrohm as the need arises.

Manufacture

The measures put into practice in the production of our instruments guarantee a constant quality standard. Production planning and manufacturing procedures, maintenance of production means and testing of components, intermediate and finished products are prescribed.

Customer support and service

Customer support involves all phases of instrument acquisition and use by the customer, i.e. consulting to define the adequate equipment for the analytical problem at hand, delivery of the equipment, user manuals, training, after-sales service and processing of customer complaints. The Metrohm service organization is equipped to support customers in implementing standards such as GLP, GMP, ISO 900X, in performing Operational Qualification and Performance Verification of the system components or in carrying out the System Validation for the quantitative determination of a substance in a given matrix.

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