

838 Advanced Sample Processor

Program version 5.838.0012

Installation Instructions

Teachware
Metrohm AG
Oberdorfstr. 68
CH-9100 Herisau
teachware@metrohm.com

These instructions are protected by copyright. All rights reserved.

Although all the information given in these instructions has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please inform the author at the address given above.

Table of contents

| | | |
|----------|---|-----------|
| 1 | Introduction | 1 |
| 1.1 | Instrument description | 1 |
| 1.1.1 | Scope of application of the 838 Advanced Sample Processor | 1 |
| 1.1.2 | Instrument models and available versions | 2 |
| 1.2 | Information about this Installation Instructions | 4 |
| 1.2.1 | Note | 4 |
| 1.2.2 | Notation and pictograms | 4 |
| 1.3 | Parts and controls | 5 |
| 1.3.1 | Overall view | 5 |
| 1.3.2 | Rear panel | 6 |
| 1.3.3 | Connections | 7 |
| 1.3.4 | Injection valve | 7 |
| 1.3.5 | Peristaltic pump | 8 |
| 1.3.6 | Sensors of the Advanced Sample Processor | 8 |
| 2 | Safety notes | 9 |
| 2.1 | General | 9 |
| 2.2 | Electrical safety | 9 |
| 2.3 | Personal protection | 10 |
| 3 | Installation | 11 |
| 3.1 | Installation flowchart | 11 |
| 3.2 | Instrument setup | 13 |
| 3.2.1 | Setup | 13 |
| 3.2.2 | Mains connection | 13 |
| 3.3 | Fitting the accessories | 14 |
| 3.3.1 | Connecting the keypad | 14 |
| 3.3.2 | Connecting the Swing Head | 14 |
| 3.3.3 | Fitting the retention plate | 15 |
| 3.3.4 | Fitting the needle or tubing | 15 |
| 3.3.5 | Fitting the filtration cell holder or Dosino holder | 16 |
| 3.3.6 | Fitting the ultrafiltration cell | 17 |
| 3.3.7 | Fitting the dialysis cell | 20 |
| 3.3.8 | Fitting the Dosino | 22 |
| 3.3.9 | Tubing connections | 23 |
| 3.3.10 | Pump tubings | 25 |
| 3.3.11 | Connection of the ultrafiltration cell | 27 |
| 3.3.12 | Connecting the dialysis cell | 29 |
| 3.3.13 | Inserting capillaries or tubings | 34 |
| 3.3.14 | Fitting the stand bracket | 34 |

| | | |
|------------|---|-----------|
| 3.3.15 | Connecting pumps | 35 |
| 3.3.16 | Connector dosing drives and stirrers | 35 |
| 3.3.17 | Installing the catch pan | 36 |
| 3.3.18 | Sample racks | 36 |
| 3.3.19 | Fitting the safety guard/splash guard | 37 |
| 3.4 | Connecting external instruments | 38 |
| 3.4.1 | General information on interface connections | 38 |
| 3.4.2 | Serial connections (RS232) | 39 |
| 3.4.3 | Remote interface..... | 39 |
| 3.5 | Connection to Metrohm IC systems..... | 40 |
| 3.5.1 | IC system remotely controlled by PC | 40 |
| 3.5.2 | Metrohm IC systems..... | 40 |
| 3.5.3 | MIC-2 modular anionic system with chemical suppression | 41 |
| 3.5.4 | MIC-2 modular anionic system with chemical suppression and sample dilution | 42 |
| 3.5.5 | Capillary connections for sample dilutions | 44 |
| 3.5.6 | Equipping the work head of the mixing container | 45 |
| 3.5.7 | MIC-2 modular anionic system with chemical suppression and ultrafiltration of the samples | 46 |
| 3.5.8 | MIC-3 modular anionic and cationic system with chemical suppression..... | 47 |
| 3.5.9 | MIC-7 anionic system with chemical suppression and dialysis | 48 |
| 3.5.10 | MIC-8 modular system with 817 Bioscan for the sugar analysis..... | 52 |
| 3.5.11 | Connection to 761 Compact IC or 790 Personal IC..... | 53 |
| 3.6 | Connection to Metrohm VA systems..... | 56 |
| 3.6.1 | Connection to a Metrohm 797 VA Computrace | 56 |
| 3.6.2 | MVA-13 modular system for automatic analysis of electroplating baths..... | 57 |
| 3.6.3 | Accessories for sample volumes greater than 10 ml..... | 57 |
| 3.7 | Sample racks | 59 |
| 3.7.1 | Metrohm standard sample racks | 59 |
| 3.7.2 | Magnetic codes | 60 |
| 4 | Service, maintenance, faults | 61 |
| 4.1 | Servicing | 61 |
| 4.2 | Care and maintenance | 61 |
| 4.2.1 | Exchanging pump tubings | 61 |
| 4.2.2 | Replacing the filtration membrane | 63 |
| 4.2.3 | Replacing the dialysis membrane | 65 |
| 5 | GLP validation – diagnosis | 66 |
| 5.1 | Validation / GLP | 66 |
| 6 | Annex | 67 |
| 6.1 | Filtration | 67 |
| 6.1.1 | Selection of possible sample types..... | 67 |

| | | |
|------------|---|-----------|
| 6.1.2 | Filter working life | 68 |
| 6.1.3 | Filter membrane selection | 68 |
| 6.2 | Technical data | 70 |
| 6.2.1 | 2 Channel peristaltic pump | 70 |
| 6.2.2 | Pump connections | 70 |
| 6.2.3 | Injection valve..... | 70 |
| 6.2.4 | Keypad with LCD display..... | 70 |
| 6.2.5 | Interfaces..... | 71 |
| 6.2.6 | MSB connections | 71 |
| 6.2.7 | Swing Head connection..... | 71 |
| 6.2.8 | Lift..... | 71 |
| 6.2.9 | Turntable | 71 |
| 6.2.10 | Stirrer connection (DIN socket)..... | 71 |
| 6.2.11 | Mains connection..... | 71 |
| 6.2.12 | Safety specification | 72 |
| 6.2.13 | Electromagnetic compatibility (EMC) | 72 |
| 6.2.14 | Ambient temperature | 72 |
| 6.2.15 | Dimensions and material | 73 |
| 6.3 | Standard equipment..... | 74 |
| 6.3.1 | Metrohm 838 Advanced Sample Processor | 74 |
| 6.3.2 | Keyboard (optional) | 91 |
| 6.3.3 | Dosing drive and units | 91 |
| 6.3.4 | Rührer..... | 91 |
| 6.3.5 | Sample racks and vessels | 92 |
| 6.3.6 | 786 Swing Head..... | 93 |
| 6.3.7 | Connecting cables | 95 |
| 6.3.8 | Filtration and dialysis accessories..... | 95 |
| 6.3.9 | Pump tubing..... | 96 |
| 6.3.10 | Needles and tubings..... | 97 |
| 6.4 | Warranty..... | 98 |
| 6.5 | Declaration of Conformity..... | 99 |
| 6.5.1 | Quality Management Principles..... | 100 |

7 Index


101

Tables

| | | |
|---------|---|----|
| Table 1 | Instrument models..... | 2 |
| Table 2 | Available versions..... | 3 |
| Table 3 | Pump tubings | 26 |
| Table 4 | MIC systems..... | 40 |
| Table 5 | Sample racks..... | 59 |
| Table 6 | Filtering various samples..... | 67 |
| Table 7 | Selection of the filtration membrane..... | 69 |

List of illustrations

| | | |
|-----------|--|----|
| Figure 1 | Models 2.838.0x10 without injection valve | 2 |
| Figure 2 | Models 2.838.0x20 with injection valve..... | 2 |
| Figure 3 | Overall view..... | 5 |
| Figure 4 | Rear panel..... | 6 |
| Figure 5 | Connection panel..... | 7 |
| Figure 6 | Injection valve | 7 |
| Figure 7 | Peristaltic pump with 2 tubing cartridges | 8 |
| Figure 8 | Magnet sensor for rack code..... | 8 |
| Figure 9 | Safety Guard (example shown. 6.2751.0xx for transfer robotic arm) | 10 |
| Figure 10 | Mains connection socket..... | 13 |
| Figure 11 | Connecting the keypad..... | 14 |
| Figure 12 | Connecting the Swing Head..... | 14 |
| Figure 13 | Fitting the retention plate | 15 |
| Figure 14 | Fitting the needle | 15 |
| Figure 15 | Filtration cell holder..... | 16 |
| Figure 16 | Dosino holder | 16 |
| Figure 17 | Fitting the ultrafiltration cell | 17 |
| Figure 18 | Fitting the dialysis cell..... | 20 |
| Figure 19 | Fitting the dialysis cell into the 6.2057.010 holder | 21 |
| Figure 20 | Fitting the Dosino..... | 22 |
| Figure 21 | Threaded adapter GL 45 (6.1618.020) | 22 |
| Figure 22 | Fitting pump tubing..... | 23 |
| Figure 23 | Connecting the ultrafiltration cell | 27 |
| Figure 24 | Capillary connections for dialysis | 29 |
| Figure 25 | Inserting capillaries or tubings | 34 |
| Figure 26 | Stand bracket..... | 34 |
| Figure 27 | MSB connections..... | 35 |
| Figure 28 | Catch pan | 36 |
| Figure 29 | Fitting the sample rack | 36 |
| Figure 30 | Fitting the safety guard | 37 |
| Figure 31 | MIC-2 interconnection with anionic system with chemical suppression | 41 |
| Figure 32 | MIC-2 interconnection with anionic system with chemical suppression and sample dilution..... | 42 |
| Figure 33 | Capillary connections for sample dilutions | 44 |
| Figure 34 | Equipping the mixing container | 45 |
| Figure 35 | MIC-2 interconnection with anionic system with chemical suppression | 46 |
| Figure 36 | MIC-3 interconnection with modular anionic and cationic system with chemical suppression | 47 |
| Figure 37 | MIC-7 interconnection with anionic system with chemical suppression and dialysis | 48 |
| Figure 38 | MIC-8 interconnection with modular system with 817 Bioscan for sugar analysis..... | 52 |
| Figure 39 | Interconnection with the 761 Compact IC via PC (RS232 only) | 53 |
| Figure 40 | Interconnection with compact systems via remote connection (PC-controlled) | 54 |
| Figure 41 | Interconnection with 797 VA Computrace | 56 |
| Figure 42 | MVA-13 connection to an analysis system for electroplating baths..... | 57 |
| Figure 43 | Accessories for larger sample volumes..... | 57 |
| Figure 44 | Remote interface..... | 71 |



1 Introduction

This Chapter offers you an initial overview of the Metrohm 838 Advanced Sample Processor and its instrument versions. It tells you how you can use this versatile instrument and familiarises you with the most important operating controls.

1.1 Instrument description

1.1.1 Scope of application of the 838 Advanced Sample Processor

The Metrohm 838 Advanced Sample Processor is an instrument with diverse applications and is designed specifically for sample preparation in ion chromatography. Thanks to its flexibility however, it can also be used for other tasks. Online filtration, dialysis, preparation of dilutions, sample enrichment, inline calibration and partial loopfill etc. are techniques which can be used easily with the 838 Advanced Sample Processor. The entire range of tasks in sample preparation can be automated to save time.

Thanks to the very extensive communication options, the Sample Processor operates not only together with the broad range of Metrohm instruments via the parallel remote interface and serial RS232 interface. Any instruments featuring a suitable communication interface can be coordinated with this instrument.

Interchangeable standard sample racks are available for many vessel sizes. Freely selectable "special beaker" positions can be used on each rack.


Equipping with a 786 Swing Head with robotic arm allows any point on a sample rack to be approached. This means that the number of samples (max. 999 rack positions) and arrangement of the samples on the sample rack can be selected virtually at will. Special customer-specific racks for individual needs can be manufactured on request.

1.1.2 Instrument models and available versions

There are various models and available versions of the Metrohm 838 Advanced Sample Processor with various accessories.

Instrument models


- 838 Advanced Sample Processor with Swing Head, peristaltic pump and 2 pump connections
- 838 Advanced Sample Processor with Swing Head, peristaltic pump, injection valve, and 2 pump connections



*Figure 1 Models 2.838.0x10
without injection valve*

Tower with 786 Swing Head
+ peristaltic pump
+ 2 external pump connections
+ 1 stirrer connection

Chassis with 3 MSB sockets for dosing drive
and/or stirrer
+ remote socket (25-pin)
+ RS232 connection (9-pin)
+ keypad connection



*Figure 2 Models 2.838.0x20
with injection valve*

Tower with 786 Swing Head
+ peristaltic pump
+ injection valve
+ 2 external pump connections
+ 1 stirrer connection

Chassis with 3 MSB sockets for dosing drive
and/or stirrer
+ remote socket (25-pin)
+ RS232 connection (9-pin)
+ keypad connection

Table 1 Instrument models

Available versions

| | |
|-------------------|---|
| 2.838.0010 | Advanced IC Sample Processor Model without injection valve, without keypad, with standard accessories |
| 2.838.0020 | Advanced IC Injection Sample Processor Model with injection valve, without keypad, with standard accessories |
| 2.838.0110 | Advanced IC Dialysis Sample Processor Model without injection valve, without keypad, with dialysis accessories |
| 2.838.0120 | Advanced IC Dilution Sample Processor Model with injection valve, without keypad, with accessories for sample dilution |
| 2.838.0210 | Advanced IC Ultrafiltration Sample Processor Model without injection valve, without keypad, with filtration accessories |
| 2.838.0310 | Advanced VA Sample Processor Model without injection valve, with keypad, with VA accessories |

Table 2 Available versions

The 2.838.00xx and 2.838.02xx versions (no keypad) are intended as ion chromatography systems, controlled by «IC Net».

The 2.838.0310 version (with keypad) is the stand-alone system version for voltammetric applications. Sample treatment is controlled by instrument methods of the Sample Processor.



«IC Net 2.3» (Service release 3 or later) sets a parameter in the 838 Advanced Sample Processor, which forces the Sample Processor to reset all of its device settings and delete all stored methods on start-up (RAM init). This behavior is intended and required for controlling the Sample processor with «IC Net».

If a 838 Advanced Sample Processor is to be used as a stand-alone system with keypad control and instrument methods, which previously was controlled by «IC Net 2.3», the setting mentioned above has to be switched off beforehand. In «IC Net» uncheck the check box "Autoinit" which is accessible after right-clicking the Sample Processor icon in the working system.

1.2 Information about this Installation Instructions

1.2.1 Note







Attention!

Please study these Installation Instructions carefully before you start to use the Instrument. The instructions contain information and warnings that must be observed by the user in order to guarantee the safe use of the instrument. Please keep these instructions near the instrument so that they are always to hand when required.

1.2.2 Notation and pictograms

The following notation and pictograms are used in these Instructions:

| | |
|---|---|
| 9 | Control element, instrument element |
| <OK> | Button, Key |
|  | Danger/Warning This symbol indicates a possible risk of death or injury to the user and possible damage to the instrument or its components by electricity. |
|  | Danger/Warning This symbol indicates a possible risk of death or injury to the user and possible damage to the instrument or its components. |
|  | Attention This symbol indicates important information that you should read before continuing. |
|  | Information This symbol indicates additional information and tips which may be of particular use to you. |

1.3 Parts and controls

1.3.1 Overall view

Safety note

The safety guard **2** must always be fitted before operating an 838 Advanced Sample Processor.




Figure 3 Overall view

| | | | |
|----------|---------------------------------------|-----------|--|
| 1 | Guide chain | 6 | Swing Head |
| 2 | Safety guard/splash protection | 7 | Injection valve * |
| 3 | Sample rack | 8 | Robotic arm** with needle adapter |
| 4 | Stirrer rail | 9 | Retention plate * |
| 5 | Peristaltic pump | 10 | Holder for ultrafiltration cell *** |
| | | 11 | Chassis |

* optional, depending on available version

** robotic arm model 6.1462.030

*** depending on available version, holder for ultrafiltration cell/dialysis cell or Dosino holder

1.3.2 Rear panel




Figure 4 Rear panel

12 Pump connections M8 (external)
for 772 or 823 pump models

13 Stirrer connection
for 802 rod stirrer or
741 magnetic stirrer

**14 Connection socket for the
786 Swing Head**

15 Connection panel

a Warning:
Biohazard
See Chapter 2.3 Personal protection

1.3.3 Connections




Figure 5 Connection panel

| | | | |
|-----------|--|-----------|--|
| 16 | Power switch | 20 | Serial RS232 connection (9-pin) |
| 17 | Keypad connection | 21 | Mains connection |
| 18 | MSB connections MSB1 ... MSB3 Metrohm Serial Bus Connection of dosing drives and stirrers | 22 | Instrument number |
| 19 | Remote connection (25-pin) | | |

1.3.4 Injection valve





Figure 6 Injection valve

The injection valve on the right-hand side of the tower, with its connections 1 to 6 (see Figure at left), can be used for complex sample preparation steps. The schematic below shows the assignments of the connections in the switch positions **FILL** and **INJECT**.



If not in use, the injection valve should be protected against dust contamination with the red protective cap.

1.3.5 Peristaltic pump

The peristaltic pump of the 838 Advanced Sample Processor can be used as a single-channel or two-channel pump. One or two tubing cartridges 6.2755.000 may be fitted.





Figure 7 Peristaltic pump with 2 tubing cartridges

| | | | |
|-----------|--|-----------|--|
| 48 | Tubing cartridge (6.2755.000) For pump tubing (6.1826.0X0) | 51 | Snap-action lever For detaching the tubing cartridge |
| 49 | Pressure lever For regulating the contact pressure | 52 | Pump drive Roller head with contact rollers |
| 50 | Retaining bracket | 53 | Retaining cam For engaging the tubing cartridges |

1.3.6 Sensors of the Advanced Sample Processor



The magnet sensor for detecting the individual rack codes is mounted beneath the turntable of the Advanced Sample Processor. The magnetic code of a rack can be read only if the rack is in initial position, i.e. if the magnet holder is precisely above the sensor. For this reason, the Sample Processor should be switched off and back on again with each rack change.

Figure 8 Magnet sensor for rack code

2 Safety notes

**Warning!**

This instrument should only be used in accordance with the information given in these installation instructions.

2.1 General

This instrument left our works in perfect condition from the point of view of its operational safety (see Technical data, safety specifications). To keep it in this condition and to continue to operate safely the following information must be carefully observed.

2.2 Electrical safety

Please observe the following guidelines:

- Only qualified Metrohm personnel should carry out service work on electronic components.
- Do not open the instrument housing as this could damage the instrument. The housing contains no components which could be serviced or exchanged by the user.

Electrical safety when handling the instrument is guaranteed within the scope of Standard IEC 61010-1.

- **Protection against electrostatic charges**

**Warning!**

Electronic components are sensitive to electrostatic charges and can be destroyed by a discharge. Before you touch any electronic components of the Sample Processor you should ground you and your tools by grasping a grounded object (e.g. the instrument housing or a radiator) in order to eliminate any electrostatic charges that may be present.

- **Opening a Sample Processor**



*When the Sample Processor is connected to the mains supply the instrument must not be opened, nor should any of its components be dismantled as otherwise you could come into contact with current-carrying components. Before opening the instrument separate it from all current sources and make sure that **mains cable has been removed from the mains connection socket!***

- **Mains connection**



*This instrument must only be operated at the specified **mains voltages** (see rear panel of instrument).*

2.3 Personal protection

The different drives of the Advanced Sample Processor and the Swing Head possess considerable force in order that they can freely move even under heavy loads. They are fitted with electronic overload limiters that are triggered if the mechanical resistance is too high.



*When the **instrument is running** do not intrude in the working range of the Sample Processor.*


*For the user there is a considerable **risk of injury from manual interference** with the automatic processes of the instrument.*

Important!



*Before you use the instrument for the first time you must install the included **6.2751.0xx Safety Guard**, see p. 37.*

The Sample Processor must not be operated without the Safety Guard in position!



*Figure 9 Safety Guard
(example shown. 6.2751.0xx for transfer robotic arm)*



Biohazard

The Advanced Sample Processor does not provide sufficient protection when processing potentially infectious samples or reagents.


In this case install the necessary safety devices.


3 Installation

This Chapter describes installation of all components of a Metrohm Sample Processor and illustrates how to connect it to other instruments. Please read the following sections carefully and follow the instructions in detail. This is the only way of ensuring that a Sample Processor will function correctly.

3.1 Installation flowchart

The following installation flowchart provides an overview of the required installation work. Please refer to the specified Chapters for further information.





* optional

3.2 Instrument setup

3.2.1 Setup

Packaging

An 838 Metrohm Advanced Sample Processor and its specially packed accessories are supplied in very protective special packaging consisting of shock-absorbing plastic foam. Please store this packaging in a safe place; it is the only way in which the safe transport of the instrument can be guaranteed.

Checks

Please check that the delivery is complete and undamaged immediately on receipt (compare with delivery note and list of accessories given in Section 6.3). If transport damage is evident please refer to the information given in Section 6.4 "Warranty".

Location

The 838 Advanced Sample Processor is a robust instrument and can therefore be used even in rough surroundings in laboratories and factories.

However, please ensure that it is not exposed to a corrosive atmosphere. The instrument should be serviced at regular intervals, particularly when operated under rough conditions.

3.2.2 Mains connection



*This instrument must only be operated at the specified **mains voltages** (see rear panel of instrument).*

Protect the connection sockets against moisture.

Mains connection **21**




Figure 10 Mains connection socket

3.3 Fitting the accessories

3.3.1 Connecting the keypad




Figure 11
Connecting the keypad

(Only required if Sample Processor methods are to be used and edited).

The 6.2142.050 keypad must be connected to the instrument rear panel at the keypad socket **17** provided. The instrument must be switched off before doing this. Squeeze the connector together at both sides in order to disconnect the connector.

After connecting the keypad, you can switch on the 838 Advanced Sample Processor at the power switch **16**. The display on the keypad starts to light. The instrument is initialised and the lift is raised fully.

3.3.2 Connecting the Swing Head




Figure 12
Connecting the Swing Head

Connect to the Swing Head connection socket **14**

Route the connection cable of the Swing Head through the guide chain of the tower and connect the Mini-DIN connector of the Swing Head connection cable to the **Swing Head** socket on the rear panel of the tower.

If a robotic arm other than the Type 6.1462.030 supplied as standard is to be used, it will be necessary to make configuration settings. You will require a 6.2142.050 keypad to do this.

Configuring the Swing Head and robotic arm

Since each type of robotic arm has different geometrical dimensions, it is absolutely essential to adapt the Setup settings for the Swing Head to the relevant type of robotic arm. No settings need to be made for the robotic arm supplied as standard.

Changes are made in the Setup dialog box. The individual settings relate to:

- swing offset
- maximum swing angle
- swing radius
- swing direction

You can open the Sample Processor's Setup dialog by pressing the **<CONFIG>** key on the keypad during switch-on of the instrument.

Please refer to the 786 Swing Head Instructions for Use, which you can obtain from Metrohm by quoting order No. 8.786.1001, for details of the individual settings and robotic arm types.

3.3.3 Fitting the retention plate

When sampling from covered vials with the needle, the retention plate **9** serves to wipe the vials when the lift is rising. It is fixed in position on the front side of the tower with the three supplied Allen screws.





Figure 13 Fitting the retention plate

3.3.4 Fitting the needle or tubing



- | | |
|----|--|
| 22 | PEEK needle (6.1835.010) or PEEK tubing (6.1835.020 / 6.1835.040 / 6.1835.050) |
| 23 | Needle holder |
| 24 | PEEK pressure screw (4.766.4320 or 6.2833.020 + 6.2744.080) |
| 25 | PEEK pressure screw (6.2744.010) |
| 26 | PTFE capillary (6.1803.070 or 6.1831.050 / 6.1831.060 / 6.1831.080) |



When using a 6.1835.xxx PEEK tubing, **do not use any stoppers** on the sample vessels since these cannot be pierced by the tubing and this would damage the PEEK tubing!
With PEEK needles **perforated stoppers** can be used.

Figure 14 Fitting the needle

1 Remove the PEEK pressure screw 24

Undo the PEEK pressure screw **24** screwed onto the needle holder **23** and remove it.

2 Insert the needle

- Insert the needle **22** (6.1835.010 PEEK needle or 6.1835.020 PEEK tubing) from the top a little into the opening of the needle holder **23**.
- Slide the PEEK ring wedge from the top over the needle **22**. The narrow side of the seal must point upwards.

3 Secure the needle

- Screw the PEEK pressure screw **24** into the needle holder **23**. When doing this, press the needle **22** upwards lightly from the bottom.
- Firmly screw the PEEK pressure screw **24** tight by hand in the needle holder **23** (do not use tools!).

3.3.5 Fitting the filtration cell holder or Dosino holder




Figure 15 Filtration cell holder




Figure 16 Dosino holder

The 6.2057.030 filtration cell holder **10** or the 6.2057.040 Dosino holder **27** can be fitted on the side panel of the tower of the 838 Advanced Sample Processor, see above.

First remove the second and third screw from the bottom from the side panel. Then fix the filtration cell holder or Dosino holder in place with the two supplied screws.

3.3.6 Fitting the ultrafiltration cell

This section describes initial fitting of the ultrafiltration cell including filtration membrane. Please note *Chapter 6.1* when assessing the need to exchange the membrane and when selecting membranes. The procedure for exchanging a filtration membrane is also described in full in *Chapter 4.2.2*.




Figure 17 Fitting the ultrafiltration cell

| | | | |
|-----------|--|-----------|---|
| 28 | Ultra filtration cell (6.2729.110) | 33 | Outlet for filtrate |
| 29 | Upper section of the filtration cell | 34 | Filtration membrane (6.2714.020) |
| 30 | Sealing ring (E.301.0111) | 35 | Inlet for sample solution |
| 31 | Lower section of the filtration cell | 36 | Outlet for sample solution |
| 32 | Screws (V.022.6030) including washer (4.754.4090) | | |

1 Preparing the ultrafiltration cell

- Remove the ultrafiltration cell **28** from its packaging and remove the three screwed-on blind stoppers.
- Using the 6.2621.070 Allen key (5 mm), fully undo the 5 screws **32**, disconnect the upper section **29** from the lower section **31** and remove the sealing ring **30**.
- Thoroughly rinse off the sealing ring, lower section and upper section of the ultrafiltration cell with ultra-pure water.



Only ultra-pure water or ethanol may be used to clean the ultrafiltration cell. Using organic solvents (e.g. acetone) will damage the Plexiglass cell!

2 Preparing the filtration membrane

- Remove a new 6.2714.020 filtration membrane **34** from the packaging using the 6.2831.010 tweezers and immerse it for approx. 2 minutes in a petri dish filled with ultra-pure water until the membrane is fully soaked with water.



It is imperative that the filtration membrane be soaked with water before inserting it into the filtration cell since this causes it to expand. If it is inserted dry, this will then lead to creasing and folding in the filtration cell which may result in blockage of the cell.

3 Inserting the filtration membrane

- **Place** the upper section **29** of the ultrafiltration cell with the inside pointing upwards on a clean paper towel.
- Insert the sealing ring **30** into the recess provided on the upper section.
- Place the wet filtration membrane, using the tweezers, centred inside the sealing ring, onto the upper section.

4 Closing the ultrafiltration cell

- Fit the upper section onto the lower section **31** so that the two parts fit perfectly.
- Using the 6.2621.070 Allen key, screw in the 5 screws **32** with washers fully and then tighten firmly.

5 Fitting the ultrafiltration cell

- After screwing the ultrafiltration cell together, fit it, as shown in *Figure 17*, in the holder **10** so that the heads of the screws **32** are located in the holes provided in the holder.



*Please use only the specified **6.2744.000 PVDF pressure screws for connection to the ultrafiltration cell**. If you use the 6.2744.010 PEEK pressure screws, this may result in stress cracking on the ultrafiltration cell!*

6 Rinsing the ultrafiltration cell

- First make the required tubing connections (see *section 3.3.11*) and then proceed as follows:
- After fitting a new filtration membrane, it is necessary to remove any air still in the filtration cell and in the lines. In order to do this, you should flush all lines with ultra-pure water from one of the special beakers for instance:
- Position a 300 ml PE bottle containing ultra-pure water as the rins-

ing solution at a special beaker position on the sample rack.

- Use the IC NET software to control the 838 Advanced Sample Processor directly. Choose the corresponding system and open the **Manual control** dialog box. You can run the individual functions of the Sampler Processor under **Autosampler**. Please refer to the IC NET software Online Help documentation for details of operation.
- Move the beaker with the rinsing solution beneath the needle with **MOVE**.
- You can lower the needle into the rinsing position with the **LIFT** function.
- Switch on the peristaltic pump with **PUMP**, and then all connected lines are rinsed and any air is removed.
- After allowing an adequate time (e.g. 5 minutes), switch the pump back off again with **PUMP** and move the needle back to home position.

3.3.7 Fitting the dialysis cell

Before fitting the dialysis cell, you must first insert a dialysis membrane.




Figure 18 Fitting the dialysis cell

| | | | |
|-----------|--|-----------|---|
| 37 | Dialysis cell (6.2729.100) | 43 | Slot for engaging the dialysis cell in the cell holder |
| 38 | Upper section of the dialysis cell | 44 | Outlet for acceptor solution |
| 39 | Sealing ring (E.301.0111) | 45 | Inlet for acceptor solution |
| 40 | Dialysis membrane (e.g. 6.2714.010) | 46 | Inlet for sample solution |
| 41 | Lower section of the dialysis cell | 47 | Outlet for sample solution |
| 42 | Screws (V.022.6030) including washer (4.754.4090) | | |

Proceed as follows:

1 Preparing the dialysis cell

- Remove the 6.2729.100 dialysis cell **37** from its packaging and remove the four screwed-on 6.2744.060 blind stoppers.
- Using the 6.2621.070 Allen key, fully undo the 5 screws **42**, disconnect the upper section **38** from the lower section **41** and remove the sealing ring **39**.

2 Cleaning the dialysis cell

- Thoroughly rinse the sealing ring **39**, upper section **38** and lower section **41** of the dialysis cell with ultra-pure water and dry them with N₂ or with a lint-free cloth.



Only ultra-pure water or ethanol may be used to clean the dialysis cell since other organic solvents (e.g. acetone) would damage the Plexiglass cell!

3 Preparing the dialysis membrane

- Remove a new dialysis membrane **40** (e.g. 6.2714.010) from the packaging using the 6.2831.010 tweezers and immerse it for approx. 2 minutes in a petri dish filled with ultra-pure water until the membrane is fully soaked with water.

4 Inserting the dialysis membrane

- Place the lower section **41** with the inside pointing upwards onto a paper towel.
- Insert the sealing ring **39** into the recess provided on the lower section **41**.
- Insert the wet dialysis membrane **40** (e.g. 6.2714.010) using the 6.2831.010 tweezers inside the sealing ring **39** on the lower section **41**.



Ensure that the water-soaked dialysis membrane does not dry out before fitting it, otherwise, it can no longer be used!

5 Closing the dialysis cell

- Fit the upper section **38** onto the lower section **41** so that the two parts fit perfectly.
- Using the 6.2621.070 Allen key, fully screw in and tighten the 5 screws **42**.




Figure 19 Fitting the dialysis cell into the 6.2057.010 holder

6 Fitting the dialysis cell

- Slide the 6.2057.010 cell holder into a free mounting rail of the 820 IC Separation Centre or the 761 Compact IC and, now that the dialysis cell has been screwed together, insert the dialysis cell **37** into the holder.

3.3.8 Fitting the Dosino

If a Dosino is required, it can be installed in the 6.2057.040 Dosino holder.

Fit the Dosino holder as specified in 3.3.5.




Figure 20 Fitting the Dosino

- Insert the 6.1618.020 threaded adapter (included in the scope of delivery of the 2.838.310) from below into the Dosino holder.
- Fit the dosing unit onto the Dosino holder and screw the threaded adapter tight.
- Connect the tubing.




Figure 21 Threaded adapter GL 45 (6.1618.020)

3.3.9 Tubing connections



Pump tubing is consumable material whose service life depends on the contact pressure. You should thus fully raise the tubing cartridges by releasing the snap-action lever **51** at the right-hand side if the pump is to be switched off for a long time (this will maintain the set optimum contact pressure).



The 6.1826.0X0 pump tubing consists of PVC or PP and may thus not be used for rinsing with solutions containing acetone. In such cases, please use different pump tubing or use a different pump for rinsing.




Figure 22 Fitting pump tubing

| | | | |
|-----------|--|-----------|--|
| 48 | Tubing cartridge | 56 | Pump tubing (6.1826.0X0) |
| 49 | Contact pressure lever | | Type, depending on intended application |
| 51 | Snap-action lever | 57 | Stopper |
| 26 | Aspirator tubing/ capillary (6.1803.0X0) PTFE tubing/capillary or PEEK capillary | | The colour indicates the material and the dimensions of the pump tubing. |
| 54 | PEEK pressure screw (6.2744.010) | 58 | PEEK coupling (6.2744.160) with tubing guard |
| 55 | PEEK coupling (6.2744.034) | 59 | PEEK filter unit (6.2821.120) optional |
| | | 60 | Delivery capillary Delivery end for connection to: dialysis cell ultrafiltration cell etc. |

1 Detaching the tubing cartridge(s)

- Detach the tubing cartridge(s) **48** by pressing in the snap-action lever **51** on the retaining bracket and disengage it (them) from the retaining cam **53**.

2 Fitting the pump tubing


- Press the contact pressure lever **49** of the tubing cartridge(s) fully down.
- Insert one length of pump tubing into the tubing cartridge(s) as shown above.
- The stopper **57** must engage in the corresponding fixture at the top of the tubing cartridge.

3 Fitting the tubing cartridge(s)

- Engage the tubing cartridge(s) in the retaining cam **53** and press it (them) down until the snap-action lever **51** on the retaining bracket **50** engages, see *Figure 7*. Ensure that the pump tubing is not kinked.

4 Fitting the coupling on the pump tubing

- Connect the PEEK coupling **55** to the **aspiration end** of the pump tubing **56** (see illustration below).
- Fit the PEEK coupling **58** with 6.2744.160 tubing guard onto the **delivery end** of the pump tubing **56**. In order to do this, dismantle the tubing guard and slide the union nut first, followed by the thrust piece, onto the tubing.



Connect the tubing to the PEEK coupling and screw the union nut onto the coupling to lock the tubing.

5 Fitting the filter unit

- Using two 6.2744.010 PEEK pressure screws **54** and a section of aspiratory tubing/capillary **26** (6.1803.0X0), screw the optionally available 6.2821.120 PEEK filter unit **59** onto the 6.2744.160 PEEK coupling **58** at the outlet end of the pump tubing **56**.

Correctly set the contact pressure:







Push the contact pressure lever **49** up until the solutions are just being aspirated. Then push the contact pressure lever one more detent position up to achieve an optimum contact pressure.

The delivery rate of the peristaltic pump primarily depends on the inside diameter of the pump tubing **56**, besides the correct contact pressure. Different pump tubings are used depending on application.

The table on the following page provides information on the characteristics and use of the available pump tubings.

3.3.10 Pump tubings

| <i>Used with</i> | <i>Order No.</i> | <i>Description</i> |
|---|-------------------|---|
| Advanced IC Ultrafiltration Sample Processor 2.838.0210: both channels | 6.1826.010 | Pump tubing made of PVC (Tygon® ST) with 2 permanently mounted white-white stoppers; Inner diameter = 1.02 mm ± 0.05 mm, Delivery 1.41 ml/min (20 rpm) 1.69 ml/min (24 rpm) |
| Advanced VA Sample Processor 2.838.0310 | 6.1826.020 | Pump tubing made of PVC (Tygon® ST) with 2 permanently mounted blue-blue stoppers; Inner diameter = 1.65 mm ± 0.05 mm, Delivery 3.75 ml/min (20 rpm) 4.50 ml/min (24 rpm) |
| Advanced IC Dialysis Sample Processor 2.838.0110: channel for pumping the acceptor solution Advanced IC Ultrafiltration Sample Processor 2.838.0210: channel for pumping the filtered solution | 6.1826.030 | Pump tubing made of PVC (Tygon® ST) with 2 permanently mounted orange-yellow stoppers; Inner diameter = 0.51 mm ± 0.05 mm, Delivery 0.40 ml/min (20 rpm) 0.48 ml/min (24 rpm) |
| Advanced IC Dialysis Sample Processor 2.838.0110: channel for pumping the sample solution | 6.1826.040 | Pump tubing made of PVC (Tygon® ST) with 2 permanently mounted black-black stoppers; Inner diameter = 0.76 mm ± 0.05 mm, Delivery 0.75 ml/min (20 rpm) 0.90 ml/min (24 rpm) |
| Optional | 6.1826.050 | Pump tubing made of PVC (Tygon® ST) with 2 permanently mounted white-yellow stoppers; Inner diameter = 0.57 mm ± 0.05 mm, Delivery 0.43 ml/min (20 rpm) 0.52 ml/min (24 rpm) |
| Optional | 6.1826.060 | Pump tubing made of PP (PharMed®) with 2 permanently mounted orange-yellow stoppers; Inner diameter = 0.51 mm ± 0.05 mm, Delivery 0.47 ml/min (20 rpm) 0.56 ml/min (24 rpm) |









| | | | |
|---|--------------------------|---|---|
| <p>Advanced IC Ultrafiltration Sample Processor 2.838.0210: channel for pumping the sample solution</p> | <p>6.1826.070</p> | <p>Pump tubing made of PVC (Tygon® ST) with 2 permanently mounted yellow-yellow stoppers; Inner diameter = 1.42 mm ± 0.05 mm, Delivery 2.55 ml/min (20 rpm) 3.06 ml/min (24 rpm)</p> |  |
| <p>Optional</p> | <p>6.1826.110</p> | <p>Long-life pump tubing made of PVC (Tygon® LFL) with 2 permanently mounted orange-yellow stoppers; Inner diameter = 0.51 mm ± 0.0102 mm, Delivery 0.40 ml/min (20 rpm) 0.48 ml/min (24 rpm)</p> |  |
| <p>Optional</p> | <p>6.1826.120</p> | <p>Long-life pump tubing made of PVC (Tygon® LFL) with 2 permanently mounted orange-white stoppers; Inner diameter = 0.59 mm ± 0.05 mm, Delivery 0.44 ml/min (20 rpm) 0.53 ml/min (24 rpm)</p> |  |
| <p>Advanced IC Sample Processor 2.838.0010: both channels Advanced IC Injection Sample Processor 2.838.0020: both channels</p> | <p>6.1826.130</p> | <p>Long-life pump tubing made of PVC (Tygon® LFL) with 2 permanently mounted white-white stoppers; Inner diameter = 1.02 mm ± 0.0127 mm, Delivery 1.41 ml/min (20 rpm) 1.69 ml/min (24 rpm)</p> |  |
| <p>Advanced IC Dilution Sample Processor 2.838.0120: both channels</p> | <p>6.1826.140</p> | <p>Long-life pump tubing made of PVC (Tygon® LFL) with 2 permanently mounted grey-grey stoppers; Inner diameter = 1.25 mm ± 0.0127 mm, Delivery 2.18 ml/min (20 rpm) 9.7 ml/min (stage 15)</p> |  |
| <p>Advanced VA Sample Processor 2.838.0310</p> | <p>6.1826.150</p> | <p>Pump tubing made of PVC (Tygon® ST) with 2 permanently mounted violet-orange stoppers; Inner diameter = 2.54 mm, Delivery 23.8 ml/min (stage 10)</p> |  |

Table 3 Pump tubings

3.3.11 Connection of the ultrafiltration cell

| | |
|-----------|---|
| 60 | PEEK pressure screw (6.2744.010) |
| 56 | Pump tubing (6.1826.070) |
| 61 | Pump tubing (6.1826.030) |
| 62 | PEEK coupling (6.2744.030) |
| 63 | PTFE capillary (6.1803.050) |
| 64 | PVDF pressure screw (6.2744.000) |
| 65 | PTFE capillary (6.1803.060) |
| 66 | PTFE capillary (6.1803.080) |
| 67 | PEEK capillary (6.1831.060) |
| 29 | Upper section of the filtration cell |
| 31 | Lower section of the filtration cell |




Figure 23 Connecting the ultrafiltration cell



Use only the specified **6.2744.000 PVDF pressure screws** to connect to the ultrafiltration cell. If the 6.2744.010 PEEK pressure screws are used, this may result in stress cracking on the ultrafiltration cell!

1 Pump tubing - filtration cell connection

- Fit a PEEK pressure screw **60** to one end of the PTFE capillary **65**.
- Screw a coupling **62** onto this pressure screw.
- Fit a PVDF pressure screw **64** onto the other end of the PTFE capillary.
- Fit this end of the PTFE capillary with the PVDF pressure screw to the inlet opening **35** on the underside of the ultrafiltration cell (see Figure 9).
- Fit the coupling **62** onto the outlet end of the pump tubing **56** (see Figure 10).

2 Filtration cell - waste tubing connection

- Fit a PVDF pressure screw **64** to one end of the PTFE capillary **66**.
- Fit this end of the PTFE capillary with the PVDF pressure screw to the outlet opening **36** on the underside of the ultrafiltration cell (see Fig-

ure 9).

- Insert the second end of the capillary into the waste container (e.g. optional: 6.1608.070 bottle with 6.1602.150 bottle top GL45).

3 Filtration cell – pump tubing connection

- Fit a PEEK pressure screw **60** to one end of the PTFE capillary **63**.
- Screw a coupling **62** onto this pressure screw.
- Fit a PVDF pressure screw **64** onto the other end of the PTFE capillary.
- Fit this end of the PTFE capillary with the PVDF pressure screw to the outlet opening for the filtrate **36** on the upper side of the ultrafiltration cell (see Figure 9).
- Fit a coupling **62** onto the inlet end of the pump tubing **61**.

4 Pump tubing – injection valve connection

- Fit a PEEK pressure screw **60** to one end of the PEEK capillary **67**.
- Screw a coupling **62** onto this pressure screw.
- Fit a coupling **62** onto the outlet end of the pump tubing **61**.
- Pull the PEEK capillary **67** through one of the feedthroughs on the front side of the 820 IC Separation Centre and, using a PEEK pressure screw **60**, screw it onto connection "1" of the injection valve in the 820 IC Separation Centre.
- Close the rotary nipple on the inner compartment side of the corresponding feedthrough and thus fix the PEEK capillary in position.

3.3.12 Connecting the dialysis cell




Figure 24 Capillary connections for dialysis

68 PTFE capillary (6.1803.040, 1 m)
supply of the sample/acceptor solution

37 Dialysis cell (6.2729.100)

61 Pump tubing (6.1826.030)
with **orange-yellow** stoppers for delivery of the **acceptor solution**

69 Pump tubing (6.1826.040)
with **black-black** stoppers for delivery of the **sample**

70 PEEK sample loop (6.1825.210, 20 µl)

71 Injection valve A

72 PTFE capillary (6.1803.040, 1 m)
peristaltic pump connection (**838 Advanced Sample Processor channel 1** for supply of fresh sample) –dialysis cell inlet (sample end)

73 PTFE capillary (6.1803.040, 1 m)
Dialysis cell outlet (sample end) – peristaltic pump connection (**838 Advanced Sample Processor channel 2** for disposal of the dialysed sample solution)

74 PTFE capillary (6.1803.040, 1 m)
Peristaltic pump (**838 Advanced Sample Processor channel 2** for disposal of the dialysed sample solution) – waste container connection

75 PTFE capillary (6.1803.040, 1 m)
Peristaltic pump (**833 Pump Unit channel 1** for supply of fresh acceptor solution) – dialysis cell inlet (acceptor end) connection

76 PEEK capillary (6.1831.050, 40 cm)
Outlet dialysis cell (acceptor end) – connection 1 valve A connection

77 PEEK capillary (6.1831.060, 1 m)
Connection 2 valve A – peristaltic pump (**833 Pump Unit channel 2** for disposal of the acceptor solution) connection

78 PTFE capillary (6.1803.040, 1 m)
peristaltic pump (**833 Pump Unit channel 2** for disposal of the acceptor solution) – waste container

64 PVDF pressure screw (6.2744.000)

Proceed as follows in order to connect the dialysis cell **37** to the IC system:



Capillaries provided with new connectors must feature a flawless, flat cut edge. In order to ensure this, it is best to use the **6.2621.080 capillary cutter**.

Connections for sample solution

1 Fitting the pump tubing to the 838 Advanced Sample Processor

- Fit the two **6.1826.040** pump tubings (**black-black** stoppers), as described in *Chapter 3.3.8*, to the peristaltic pump of the 838 Advanced Sample Processor.

2 Connecting the aspiratory capillary for sample

- Attach one end of the **6.1803.040** aspiratory capillary **68** (L =1 m) with a **6.2744.010** PEEK pressure screw to the PEEK coupling on the aspiration end of the **black-black** pump tubing of **channel 1** of the **838 Advanced Sample Processor** (see *Chapter 3.3.12* and *Figure 22*).
- Connect the other end of the aspiratory capillary **68** to the PEEK needle **22** (see *Chapter 3.3.4* and *Figure 14*) and fix the capillary in position in the guide chain of the 838 Advanced Sample Processor.

3 Connecting the dialysis cell

- Attach the **6.1803.040** PTFE capillary **72** (L =1 m) with a **6.2744.010** PEEK pressure screw to the **6.2744.034** PEEK coupling on the delivery end of the **black-black** pump tubing of **channel 1** of the peristaltic pump of the **838 Advanced Sample Processor** (see *Chapter 3.3.12* and *Figure 22*). Guide the other end of the capillary through one of the openings at the front or on the side of the 820 Separation Centre into the column chamber and screw it on to the inlet opening **46** of the cell lower section **41** of the dialysis cell using a **6.2744.000** PVDF pressure screw **64** (see *Chapter 3.3.7* and *Figure 18*).
- Screw the **6.1803.040** PTFE capillary **73** (L =1 m) with a **6.2744.000** PVDF pressure screw **64** to the outlet opening **47** of the cell lower section **41** of the dialysis cell (see *Chapter 3.3.7* and *Figure 18*) and guide it through one of the openings at the front or at the side of the 820 Separation Centre out of the column chamber. Attach the other end of the PTFE capillary **73** with a **6.2744.010** PEEK pressure screw **64** to the PEEK coupling at the aspiration end of the **black-black** pump tubing of **channel 2** of the peristaltic pump of the **838 Advanced Sample Processor** (see *Chapter 3.3.12* and *Figure 22*).



Use only the specified **6.2744.000 PVDF pressure screws** for connection to the dialysis cell. If the 6.2744.010 PEEK pressure screws are used, this could cause stress cracking on the dialysis cell!

4 Connecting the waste container

- Attach the 6.1803.040 PTFE capillary **74** (L =1 m) with a 6.2744.010 PEEK pressure screw to the 6.2744.034 PEEK coupling on the delivery end of the **black-black** pump tubing of **channel 2** of the peristaltic pump of the **838 Advanced Sample Processor** (see *Chapter 3.3.12* and *Figure 22*). Insert the other end of the capillary into an adequately large waste container and secure it there.

Connections for acceptor solution

1 Fitting the pump tubing to the 833 Pump Unit

- Fit the two pump tubings **6.1826.030** (**orange-yellow** stoppers) as described in *Chapter 3.3.9*, to the **833 Pump Unit**.

2 Connecting the aspiratory capillary for acceptor solution

- Attach one end of the 6.1803.040 aspiratory capillary **68** (L =1 m) with a 6.2744.010 PEEK pressure screw to the PEEK coupling on the aspiration end of the **orange-yellow** pump tubing of **channel 1** of the **833 Pump Unit** (see *Chapter 3.3.9* and *Figure 22*).
- Insert the other end of the aspiratory capillary **68** into the vessel with **acceptor solution** and secure it there. **Ultra-pure water** which has been **degassed** beforehand (with N₂, He or vacuum) is normally used as the acceptor solution.

3 Connecting the dialysis cell

- Attach the 6.1803.040 connection capillary **75** (L =1 m) with a 6.2744.010 PEEK pressure screw to the 6.2744.034 PEEK coupling on the delivery end of the **orange-yellow** pump tubing of **channel 1** of the **833 Pump Unit** (see *Chapter 3.3.9* and *Figure 22*). Guide the other end of the capillary through one of the openings at the front or on the side of the 820 Separation Centre into the column chamber and screw it onto the inlet opening **45** of the cell upper section **38** of the dialysis cell using a 6.2744.000 PVDF pressure screw **64** (see *Chapter 3.3.7* and *Figure 18*).
- Screw the 6.1831.050 PEEK capillary **76** (L =40 cm) with a 6.2744.000 PVDF pressure screw **64** to the outlet opening **44** of the cell upper section **38** of the dialysis cell (see *Chapter 3.3.7* and *Figure 18*).



Use only the specified **6.2744.000 PVDF pressure screws** to connect to the dialysis cell. If the 6.2744.010 PEEK pressure screws are used, this could cause stress cracking on the dialysis cell!

4 Connecting the injection valve

- Screw the free end of the 6.1831.050 PEEK capillary **76** (L = 40 cm) with a 6.2744.010 PEEK pressure screw to connection "1" of the injection valve A **71**.
- Screw the 6.1831.060 PEEK capillary **77** (L = 1 m) using a 6.2744.010 PEEK pressure screw to connection "2" of the injection valve A **71** and guide the other end of the capillary through one of the openings at the front or at the side of the 820 Separation Centre to the outside. Attach the other end of the PEEK capillary **77** with a 6.2744.010 PEEK pressure screw to the PEEK coupling on the aspiration end of the **orange-yellow** pump tubing of **channel 2** of the **833 Pump Unit** (see Chapter 3.3.9 and Figure 22).

5 Connecting the waste container

- Attach the 6.1803.040 PTFE capillary **78** (L = 1 m) with a 6.2744.010 PEEK pressure screw to the 6.2744.034 PEEK coupling on the delivery end of the **orange-yellow** pump tubing of **channel 2** of the **833 Pump Unit** (see Chapter 3.3.9 and Figure 22). Insert the other end of the capillary into an adequately large waste container and secure it there.

Conditioning the dialysis system

Before the first analysis, the dialysis cell with the fitted dialysis membrane and all tubing connections must be rinsed with ultra-pure water. Please proceed as follows:

1 Settings on the 820 IC Separation Centre

- In the «IC Net» software, switch injection valve A in the 820 IC Separation Centre to position "FILL".


2 Commissioning the peristaltic pumps

- Immerse the two aspiratory tubings **68** and **72** (see *Figure 24*) in the acceptor solution (degassed ultra-pure water).
- Switch on the 838 Advanced Sample Processor and the 833 Pump Unit with their power switches.
- Switch on the peristaltic pumps of both instruments in «IC Net» in the corresponding 'Manual Control' dialog box. Please refer to the «IC Net» software Online Help documentation for details of operation.
- Set the contact pressure on both tubing cartridges: Push the pressure lever **49** (see *Chapter 1.3.5*) up until the solutions are just being aspirated. Then push the pressure lever 1 further detent position upwards to achieve an optimum contact pressure.
- Flush the dialysis system for approx. 10 minutes with ultra-pure water and check whether solution is emerging uniformly on the two inlet capillaries leading into the waste container.

3 Conditioning the dialysis membrane

- Flush the dialysis system for approx. 20 minutes with ultra-pure water. Check whether solution is emerging uniformly on the two inlet capillaries leading into the waste container.
- Check all tubings from the supply vessels through the tubing cartridges and the dialysis cell through to the waste vessels for the presence of emerging fluid. If fluid is escaping at any point, the corresponding connection must be tightened or exchanged.
- If air bubbles remain in the dialysis cell, unscrew the PEEK capillary **76** (acceptor solution) and PTFE capillary **73** (sample) from the outlets **44 resp. 47** of the dialysis cell and wait until the air bubbles have disappeared. Then screw the tubings back onto the dialysis cell.

3.3.13 Inserting capillaries or tubings




The individual chain links of the guide chain can be opened using a screwdriver or a similar implement in order to insert capillaries or tubings.

Please note the marking on the chain links.

Figure 25 Inserting capillaries or tubings

3.3.14 Fitting the stand bracket

An external dilution or measuring cell can be fitted with the aid of a 6.2001.070 stand bracket on the Advanced Sample Processor. The stand bracket must be engaged on the stirrer rail and fixed in position with the two Allen screws.



- The height and lateral alignment of the cross-brace can be adjusted. It is fixed in position with the supplied screw on the rear panel of the tower.

Figure 26 Stand bracket

3.3.15 Connecting pumps

Two external pumps can be connected and controlled on the 838 Advanced Sample Processor.

Ext.
Pump 2



The pump connections (for 3-pin M8 connectors) on the rear side of the tower of the Advanced Sample Processor supply a voltage of 16 V and may be loaded with max. 600 mA.


The appropriate pump models are:

- Metrohm 823 Membrane Pump Unit (membrane pump, suitable for aqueous media without precipitation)
- Metrohm 772 Pump Unit (peristaltic pump, suitable for organic solvents and aqueous solutions containing precipitation)

3.3.16 Connector dosing drives and stirrers

Metrohm Sample Processors feature various connections to which stirrers or dosing drives can be connected.

Stirrers




The **722/802 Rod stirrer** and the **741 Magnetic stirrers** feature a DIN connector and are connected on the rear panel of the tower.

As is the case on other Metrohm instruments, the three **MSB sockets** (MSB = Metrohm Serial Bus) can be used in diverse ways. Connect the **801 Stirrer** at this point.

Dosing drives

Up to three dosing drives of Type **700/800 Dosino** or **685/805 Dosimat** can be connected directly to the MSB sockets of a Sample Processor and operated. A **800 Dosino** or **805 Dosimat** can also be connected via the **MSB output socket** of another instrument connected to the MSB, such as the 801 Stirrer in a so-called 'daisy chain'. It must be noted that only one single instrument of a specific type may occur in such a chain, e.g. one 801 Stirrer and one 800 Dosino.

The **700 Dosino** and the **685 Dosimat** must be connected directly to an **MSB socket** of the Sample Processor.



A connected dosing drive is detected automatically by the instrument only when switching on.

3.3.17 Installing the catch pan

The 6.2711.060 catch pan is available as an optional accessory. It prevents escaping fluid spreading over the laboratory bench or penetrating the interior of the instrument.




Figure 28 Catch pan

Attach the supplied tubing to the drain connector of the catch pan. Detach the sample rack and move the catch pan over the turntable so that it lies on the Sample Processor's stirrer rail. Insert the other end of the tubing into a drain or waste canister.

3.3.18 Sample racks




Figure 29 Fitting the sample rack

Sample racks must be fitted so that the two guide cams of the turntable engage in the guide openings provided in the base of the rack.

After fitting a sample rack, the Sample Processor should be switched off and then back on again so that the magnetic code of the rack can be read in. This is possible only if the rack is in initial position.

3.3.19 Fitting the safety guard/splash guard

It is essential to fit the safety guard before using the 838 Advanced Sample Processor for the first time. It not only acts as an effective splash guard but also prevents the risk of injury. The safety guard must be fitted with the securing screw on the side of the tower of the 838 Advanced Sample Processor. The height of the safety guard can be adapted to the sample vessels.




Figure 30 Fitting the safety guard

3.4 Connecting external instruments

There are numerous options for connecting the 838 Advanced Sample Processor to Metrohm systems. The 838 Advanced Sample Processor features a serial RS-232 interface via which it is possible to access the instrument functions directly using suitable control software. If working with the sequence methods of the instrument, Sample Processor and connected instruments can be coordinated with control signals via the 25-pin remote interface.

3.4.1 General information on interface connections



*The 838 Advanced Sample Processor must always be switched off with the power switch **34** before connecting an external instrument to the remote connection **30** or to the RS232 interface **40**!*

Remote interface

Any external instruments can be connected to the 25-pin remote connection **30**. The 838 Advanced Sample Processor can be controlled remotely via the 8 input leads. External instruments can be controlled via the 14 output lines.

The connector pin assignment of the remote interface, its functions, electrical conditions and statuses and the assignment of the remote input lines are described in *Chapter 5.1*.

RS232 interface

The connection options at the serial RS232 interface **40** are diverse. Besides Metrohm instruments (in this case the 830 IC Interface) featuring Metrohm remote control language, it is also possible to connect a PC (see *Chapter 2.7*). Third-party instruments featuring a serial RS232 interface can also be connected.


Connecting cables

Only Metrohm cables should be used for interconnection of the 838 Advanced Sample Processor with other devices. Only these cables guarantee interference-free data transfer.



The Metrohm remote cables have a designation at the ends of the cable, indicating for what instrument the relevant connector is scheduled and at what connector position it must be connected.


Example:



3.4.2 Serial connections (RS232)

The connection options at the serial RS232 interface are diverse. Besides all Metrohm instruments featuring the Metrohm remote control language, it is possible to connect a personal computer. Any third-party instruments featuring a serial RS232 interface may also be connected.

Correct data transfer requires correct setting of the transmission parameters. They must correspond to the settings of the interface on the connected instrument.



3.4.3 Remote interface

Connected peripheral instruments, such as the 819 Advanced IC Detector etc., can be controlled via the remote interface (25-pin socket).

14 lines (output 0–13) are available for output of signals.
8 lines (input 0–7) are available for reception of signals.

Only the designated Metrohm remote cables should be used to connect Metrohm instruments.

3.5 Connection to Metrohm IC systems

3.5.1 IC system remotely controlled by PC

Data acquisition and control of the IC components can be performed from the PC using the «IC Net 2.x» PC software and the 830 IC Interface. In this case, the 838 Advanced Sample Processor is connected to the IC system both via the serial RS232 interface and via the remote interface. The preprogrammed instrument methods of the 838 Advanced Sample Processor are not required in this case since the corresponding sequences are programmed and controlled by the PC software. This Chapter describes the various combination options of the 838 Advanced Sample Processor with Metrohm IC systems on the basis of this technology.

3.5.2 Metrohm IC systems

Metrohm IC systems can be subdivided in to modular IC systems ("MIC") and compact IC systems (e. g. **761 Compact IC**). The following table provides an overview of preconfigured modular Metrohm IC systems.

| Designation | Function of the modular IC system |
|---------------|--|
| MIC-1 | Anionic or cationic system with electronic suppression |
| MIC-2 | Anionic system with chemical suppression |
| MIC-3 | Anionic und cationic system with chemical suppression |
| MIC-4 | Anionic or cationic system with electrical suppression and enrichment |
| MIC-5 | Anionic system with chemical suppression and enrichment |
| MIC-6 | Anionic system with chemical suppression, enrichment and matrix elimination |
| MIC-7 | Anionic system with chemical suppression and dialysis |
| MIC-8 | Modular system with 817 Bioscan for sugar analysis |
| MIC-9 | Combined system for anionic and sugar analysis |
| MIC-10 | Binary high pressure-gradient system for anion determination with chemical suppression |

Table 4 MIC systems

The electrical interconnection of some of these systems with the 838 Advanced Sample Processor with complete control of the entire IC system by the «**IC Net 2.3**» software is described below. Since certain versions differ essentially as the result of the additional enrichment technology (see MIC-1 / MIC-4 and MIC-2 / MIC-5), we shall describe the **MIC-2** and **MIC-3** systems at this point. Systems **MIC-7** and **MIC-8** are also listed as further typical versions.

3.5.3 MIC-2 modular anionic system with chemical suppression

This system requires the **model 2.838.0010** of the 838 Advanced Sample Processor.

The 838 Advanced Sample Processor is connected to an anionic IC system with chemical suppression with 819 IC detector, 820 IC Separation Centre (2.820.0230), 818 IC pump, 833 IC Liquid Handling Unit and 830 IC Interface in accordance with *Figure 31* using the 6.2141.190 remote cable and the 6.2134.100 RS232 cable.




Figure 31 MIC-2 interconnection with anionic system with chemical suppression

838-specific settings in the «IC Net 2.3» program

Program for the 838 Advanced Sample Processor:

| | | | |
|-----|-------|------------|--|
| 001 | Ctrl | INIT 819 | – Initialise remote lines on 819 |
| 002 | Move | sample | – Move needle to sample position |
| 003 | Lift | work | – Move lift with needle to working height |
| 004 | Ctrl | ZERO 1 | – Trigger autozero on 819/1 IC Detector |
| 005 | Ctrl | FILL A 1 | – Switch over injection valve on 820 "Fill" |
| 006 | Pump | 240 s | – Fill sample loop for 240 s with sample |
| 007 | Ctrl | INJECT A 1 | – Switch over injection valve on 820 to "Inject" |
| 008 | Start | DA | – Start data recording |

| | | |
|----------------------|---------------|---|
| 001 Valve | fill | – Switch over injection valve on 838 to "Fill" |
| 002 Move | spec1 | – Move needle to rinsing beaker * |
| 003 Lift | special | – Position lift with needle to special position * |
| 004 Dosino MSB1 rate | 5 | – Set dosing rate |
| 005 Dosino MSB1 | 3.5 | – Eject 3.5 ml of ultra-pure water (prepare pipet.) |
| 006 Move | sample | – Move needle to sample |
| 007 Dosino MSB1 rate | -0.1 | – Form air bubble |
| 008 Lift | work | – Move lift with needle to working position * |
| 009 Dosino MSB1 rate | 5 | – Set dosing rate |
| 010 Dosino MSB1 | -3 | – Aspirate 3 ml of sample ** |
| 011 Move | spec1 | – Move needle to rinsing beaker |
| 012 Lift | special | – Move lift with needle to special position |
| 013 Dosino MSB1 | compen. | – Compensate for mech. play of Dosino coupling |
| 014 Wait | 5 | – Waiting time |
| 015 Valve | inject | – Move injection valve on 838 to "Inject" position |
| 016 Wait | 1 | – Waiting time |
| 017 Dosino MSB1 rate | 2 | – Set dosing rate |
| 018 Dosino MSB1 | Concentr.vol. | – Sample volume is dosed |
| 019 Wait | 5 | – Waiting time |
| 020 Valve | fill | – Switch over injection valve on 838 to "Fill" |
| 021 Dosino MSB1 rate | 5 | – Set dosing rate |
| 022 Dosino MSB1 | 3.5 | – Eject rest of sample into rinsing beaker |
| 023 Valve | inject | – Set injection valve on 838 to "Inject" position |
| 024 Wait | 1 | – Waiting time |
| 025 Dosino MSB1 rate | 2 | – Set dosing rate |
| 026 Dosino MSB1 | Dilution vol. | – Dilution volume (ultra-pure water) is dosed |
| 027 Stirrer MSB2 | 15 | – Set stirring speed |
| 028 Valve | fill | – Switch over injection valve on 838 to "Fill" |
| 029 Ctrl | Fill A1 | – Switch over injection valve on 820 to "Fill" |
| 030 Pump speed | 4 | – Set delivery rate |
| 031 Pump | 65 | – Transfer diluted sample to sample loop |
| 032 Ctrl | Zero 1 | – Trigger autozero on 819/1 IC Detector |
| 033 Ctrl | Inject A1 | – Switch over injection valve on 820 to "Inject" |
| 034 Pump speed | 10 | – Set delivery rate |
| 035 Pump | 30 | – Mixing vessel is emptied |
| 036 Stirrer MSB2 | on | – Switch on stirrer |
| 037 Valve | inject | – Set injection valve on 838 to "Inject" position |
| 038 Dosino MSB1 rate | 20 | – Set dosing rate |
| 039 Dosino MSB1 | 5 | – Rinse mixing vessel with ultra-pure water |
| 040 Pump | on | – Switch on pump |
| 041 Valve | fill | – Switch over injection valve on 838 to "Fill" |
| 042 Wait | 60 | – Empty mixing vessel again |
| 043 Pump | off | – Switch off pump |
| 044 Valve | inject | – Set injection valve on 838 to "Inject" position |
| 045 Dosino MSB1 rate | 20 | – Set dosing rate |
| 046 Dosino MSB1 | 5 | – Rinse mixing vessel again |
| 047 Valve | fill | – Switch over injection valve on 838 to "Fill" |
| 048 Pump speed | 10 | – Set delivery rate |
| 049 Pump | on | – Switch on pump |
| 050 Wait | 60 | – Empty mixing vessel again |
| 051 Stirrer MSB2 | off | – Switch off stirrer |
| 052 Wait | 10 | – Waiting time |
| 053 Pump | off | – Switch off pump |
| 054 Dosino MSB1 | fill | – Fill Dosino for next dilution |

* must be predefined in «IC Net»

** max. volume must be taken into account. The sample must not reach the Dosino cylinder.

3.5.5 Capillary connections for sample dilutions




Figure 33 Capillary connections for sample dilutions

| | | | |
|----------|--|----------|--|
| A | Capillary (6.1803.120) IØ = 0.97 mm, L = 3.0 m | E | Capillary (6.1803.060) IØ = 0.97 mm, L = 0.2 m |
| B | Capillary (6.1803.070) IØ = 0.97 mm, L = 0.4 m | F | Capillary (6.1803.080) IØ = 0.97 mm, L = 1.0 m |
| C | Capillary (6.1803.070) * IØ = 0.97 mm, L = 0.4 m | G | Capillary (6.1803.080) IØ = 0.97 mm, L = 1.0 m |
| D | Capillary (6.1803.080) IØ = 0.97 mm, L = 1.0 m | H | Injection loop |

* The **capillary C** serves to aspirate. It must extend down to the bottom of the mixing container down to the lowest point of the inclined container bottom.

| | | | |
|----------|--|----------|--|
| 1 | 800 Dosino with 10 ml dosing unit, filled with water | 5 | PE bottle (6.1608.080) as rinsing vessel, on sample rack (6.2041.440) |
| 2 | Switchover valve on the 838 Note tubing connection | 6 | Peristaltic pump on 838 with pump tubing LFL grey/grey (6.1826.140) |
| 3 | Mixing container (6.1456.210) with 801 magnetic stirrer and titration head 6.1414.070 | 7 | Injection valve in 820 Separation Centre Note tubing connection |
| 4 | Sample vessels (6.2743.050) with PP stopper (6.2743.070), on sample rack (6.2041.440) | 8 | IC column |

3.5.6 Equipping the work head of the mixing container

The 6.1414.070 work head of the mixing container must be provided with a capillary for supply of sample and water and a capillary for aspiration of diluted solution. The other openings must be fitted with blind stoppers.




Figure 34 Equipping the mixing container

* The PEEK capillary for aspiration must extend down to the bottom of the mixing container!

3.5.7 MIC-2 modular anionic system with chemical suppression and ultrafiltration of the samples

This system requires the **model 2.838.0210** of the 838 Advanced Sample Processor. Install the ultrafiltration cell, see *Chapter 3.3.6*.

The 838 Advanced Sample Processor is connected to an anionic IC system with chemical suppression with 819 IC detector, 820 IC Separation Centre, 818 IC pump, 833 IC Liquid Handling Unit and 830 IC Interface in accordance with *Figure 35* using the 6.2141.190.remote cable and the 6.2134.100 RS232 cable.




Figure 35 MIC-2 interconnection with anionic system with chemical suppression

838-specific settings in the «IC Net 2.3» program

Program for the 838 Advanced Sample Processor:

| | | | |
|-----|-------|------------|--|
| 001 | Ctrl | INIT 819 | – Initialise remote lines on 819 |
| 002 | Move | sample | – Move needle to sample position |
| 003 | Lift | work | – Move lift with needle to working height |
| 004 | Ctrl | ZERO 1 | – Trigger autozero on 819/1 IC Detector |
| 005 | Ctrl | FILL A 1 | – Switch over injection valve on 820 to "Fill" |
| 006 | Pump | 240 s | – Fill sample loop for 240 s with sample |
| 007 | Ctrl | INJECT A 1 | – Switch over injection valve on 820 to "Inject" |
| 008 | Start | DA | – Start data recording |

3.5.8 MIC-3 modular anionic and cationic system with chemical suppression

This system requires the **model 2.838.0010** of the 838 Advanced Sample Processor.

The combined anionic/cationic IC system can operate either with chemical suppression or with electronic suppression. The 838 Advanced Sample Processor is connected to the System MIC-3 with two 819 IC detectors, 820 IC Separation Centre (2 injectors), two 818 IC pumps, 833 IC Suppressor Module and 830 IC Interface in accordance with *Figure 36* using the 6.2125.120 remote adapter, the 6.2125.090 remote cable and the 6.2134.100 RS232 cable.




Figure 36 MIC-3 interconnection with modular anionic and cationic system with chemical suppression

838-specific settings in the «IC Net program 2.3»

Configuration of additional remote lines for the 838 Advanced Sample Processor:

```

Fill A 2    ***0*010***
Inject A 2  ***1*000***
Reset 2    ***0*000***
ZERO 2     ***0*011***
    
```

Program for the 838 Advanced Sample Processor:

| | | | |
|-----|------|--------------|--|
| 001 | Move | sample | – Move needle to sample position |
| 002 | Lift | work | – Move lift with needle to working height |
| 003 | Ctrl | ZERO 1 | – Trigger autozero on 819/1 IC Detector |
| 004 | Ctrl | FILL A 1 | – Switch over injection valve 1 on 820 to "Fill" |
| 005 | Ctrl | STEP MSM 833 | – Trigger step in the Suppressor Unit |
| 006 | Pump | 120 s | – Fill sample loop for 120 s with sample |
| 007 | Ctrl | INJECT A 1 | – Switch over injection valve 1 on 820 to "Inject" |
| 008 | Move | Sample+1 | – Move needle to next sample |
| 009 | Lift | work | – Move lift with needle to working height |

| | | | |
|-----|------|------------|--|
| 010 | Ctrl | ZERO 2 | – Trigger autozero on 819/2 IC Detector |
| 012 | Scan | Wait1 | – Wait for Remote signal |
| 013 | Ctrl | FILL A 2 | – Switch over injection valve 2 on 820 to "Fill" |
| 015 | Pump | 120 s | – Fill sample loop for 120 s with sample |
| 016 | Ctrl | INJECT A 2 | – Switch over injection valve 2 on 820 to "Inject" |

3.5.9 MIC-7 anionic system with chemical suppression and dialysis

This system requires the **model 2.838.0110** of the 838 Advanced Sample Processor.

The 838 Advanced Sample Processor is connected to an anionic IC system with chemical suppression with 819 IC detector, 820 IC Separation Centre, 818 IC pump, 833 IC Liquid Handling Unit and 830 IC Interface in accordance with *Figure 37* using the 6.2141.190 remote cable and the 6.2134.090 RS232 cable.




Figure 37 MIC-7 interconnection with anionic system with chemical suppression and dialysis

838-specific settings in the «IC Net 2.3» program

The sequences of an automated IC analysis with prior dialysis can be programmed in two different ways:

- A) A sample is first prepared by dialysis for each analysis and this sample is then chromatographically separated.

⇒ The total analysis time per sample comprises the times for dialysis and for the chromatogram.

B) Dialysis of the next sample is performed as early as the phase in which one sample is being chromatographically separated.

⇒ The total analysis time per sample is shortened to the time required for one chromatogram.

A) Dialysis and chromatogram of a sample in one sequence

Entry of the program for the 838 Advanced Sample Processor:

- Dialysis by stop flow procedure, dialysis time 10 min.

| | | | |
|-----|------|---------------|---|
| 001 | Move | sample | – Move needle to sample position |
| 002 | Lift | work | – Move lift with needle to working height |
| 003 | Ctrl | Pump 833 on | – Switch on pump 833 (acceptor solution) |
| 004 | Pump | on | – Pump sample solution |
| 005 | Wait | 120 s | – Wait for rinsing time |
| 006 | Ctrl | Pump 833 off | – Switch off pump 833, dialysis starts |
| 007 | Wait | 420 s | – Dialysis time + 180 s |
| 008 | Ctrl | FILL B/STEP 1 | – Trigger step of suppressor unit |
| 009 | Wait | 180 s | – Wait until baseline is stable (= 10 min. dialysis time) |
| 010 | Ctrl | ZERO 1 | – Trigger autozero on 819/1 IC Detector |
| 011 | Ctrl | FILL A 1 | – Switch over injection valve in 820 to "Fill" |
| 012 | Ctrl | Pump 833 on | – Pump acceptor solution |
| 013 | Wait | 30 s | – Transfer time (must be determined) |
| 014 | Ctrl | INJECT A 1 | – Switch over injection valve in 820 to "Inject" |
| 015 | Ctrl | Pump 833 off | – Stop pumping acceptor solution |
| 016 | Pump | off | – Stop pumping sample solution |

B) Chromatogram with simultaneous dialysis of the next sample

The following program descriptions prerrequisite chromatograms with a duration of e. g. 20 minutes. Dialysis of the next sample is started 10 minutes prior to the end of each chromatogram. A waiting time of 600 seconds (see below) is entered for this purpose in line 008 of the 838 program in System **Dialysis.smt**. This time information may need to be adapted to specific requirements.

A program which controls both injection of the dialysed sample and also dialysis of the next sample ('Dialysis.smt') is written for the 838 Advanced Sample Processor in the existing system.

In order to also process the first and last sample of a Sample Queue correctly, it is necessary to create programs in separate system files for these samples. The 838 program in system file 'Start-Dialysis.smt' only starts dialysis of the first sample. 'End-Dialysis.smt' is required for injection of the last dialysed sample.

Dialysis.smt

Entry of the program for the 838 Advanced Sample Processor

Dialysis program:

| | | | |
|-----|------|---------------|---|
| 001 | Pump | on | – Pump sample solution |
| 002 | Ctrl | FILL A 1 | – Switch over injection valve in 820 to "Fill" |
| 003 | Ctrl | Pump 833 on | – Pump acceptor solution |
| 004 | Wait | 30 s | – Transfer time (must be determined) |
| 005 | Ctrl | INJECT A 1 | – Switch over injection valve in 820 to "Inject" |
| 006 | Ctrl | Pump 833 off | – Stop pumping acceptor solution |
| 007 | Pump | off | – Stop pumping sample solution |
| 008 | Wait | 600 | – Waiting time for synchronisation |
| 009 | Move | Sample+1 | – Approach next sample |
| 010 | Lift | work | – Move lift with needle to working height |
| 011 | Pump | on | – Pump sample solution |
| 012 | Ctrl | Pump 833 on | – Pump acceptor solution |
| 013 | Wait | 120 | – Rinsing time |
| 014 | Ctrl | Pump 833 off | – Stop pumping acceptor solution, start of dialysis |
| 015 | Wait | 480 s | – Dialysis time (120 s + 480 s = 10 min.) |
| 016 | Ctrl | FILL B/STEP 1 | – Trigger step of suppressor unit |
| 017 | Wait | 120 | – Wait until baseline is stable |
| 018 | Ctrl | ZERO 1 | – Trigger autozero on 819/1 IC Detector |
| 019 | Pump | off | – Stop pumping sample solution |

Start-Dialysis.smt

- Remove the data recorder by clicking on the recorder icon with the right mouse button and selecting '*unlink*'.

Start program:

| | | | |
|-----|------|---------------|---|
| 001 | Move | sample | – Move needle to sample position |
| 002 | Lift | work | – Move lift with needle to working height |
| 003 | Ctrl | Pump 833 on | – Pump acceptor solution |
| 004 | Pump | on | – Pump sample solution |
| 005 | Wait | 120 s | – Rinsing time |
| 006 | Ctrl | Pump 833 off | – Stop pumping acceptor solution, start of dialysis |
| 007 | Wait | 480 s | – Dialysis time (780 s + 120 s = 10 min) |
| 008 | Ctrl | FILL B/STEP 1 | – Trigger step of suppressor unit |
| 009 | Wait | 180 s | – Wait until baseline is stable |
| 010 | Ctrl | ZERO 1 | – Trigger autozero on 819/1 IC Detector |
| 011 | Pump | off | – Switch over injection valve the 820 to "Inject" |

End-Dialysis.smt

End program:

| | | | | |
|-----|------|--------------|--|--|
| 001 | Pump | on | | – Pump sample solution |
| 002 | Ctrl | FILL A 1 | | – Switch over injection valve in 820 to "Fill" |
| 003 | Ctrl | Pump 833 on | | – Pump acceptor solution |
| 004 | Wait | 30 s | | – Transfer time (must be determined) |
| 005 | Ctrl | INJECT A 1 | | – Switch over injection valve in 820 to "Inject" |
| 006 | Ctrl | Pump 833 off | | – Stop pumping acceptor solution |
| 007 | Pump | off | | – Stop pumping sample solution |

The Sample Queue to be processed then contains the following entries for instance:

| | System | Ident. | Vial | Chrom. No. | Dialysis No. |
|---|---------------------------|-----------------|-------------|-------------------|---------------------|
| 1 | Start-Dialysis.smt | Sample 0 | 1 | - | 1 |
| 2 | Dialysis.smt | Sample 1 | 1 | 1 | 2 |
| 3 | Dialysis.smt | Sample 2 | 2 | 2 | 3 |
| 4 | Dialysis.smt | Sample 3 | 3 | 3 | 4 |
| 5 | Dialysis.smt | Sample 4 | 4 | 4 | 5 |
| 6 | Dialysis.smt | Sample 5 | 5 | 5 | 6 |
| 7 | End-Dialysis.smt | Sample 6 | 6 | 6 | - |

This example shows processing of a sample series of 6 samples. Columns 'Chrom. No.' and 'Dialysis No.' serve only for information purposes in this case, Please note that the sample position specified here is not approached in the 838 program of system 'Dialysis.smt' for dialysis but that this position is incremented by one, thus processing the next sample. Consequently, all sample vessels must be located in a closed sequence on the sample rack.

3.5.10 MIC-8 modular system with 817 Bioscan for the sugar analysis

This system requires the **model 2.838.0020** of the 838 Advanced Sample Processor.

The 838 Advanced Sample Processor is connected to a modular Bioscan system with 819 IC detector and 817 Bioscan in accordance with *Figure 38* via PC (RS232 connections).





Figure 38 MIC-8 interconnection with modular system with 817 Bioscan for sugar analysis


Capillary connections of the injection valve of the 838 Advanced Sample Processors:



Notice:

* cut a short fragment of the 6.1831.060 PEEK capillary ($\varnothing=0.5$ mm).

** cut a fragment (approx 1 m) of the 6.1831.010 PEEK Capillary ($\varnothing=0.25$ mm).



838-specific settings in the «IC Net 2.3» program

Program for the 838 Advanced Sample Processor:

| | |
|---|---|
| <p>001 Move Sample</p> <p>002 Lift work</p> <p>003 Valve Fill</p> <p>004 Pump 240 s</p> <p>005 Valve Inject</p> <p>006 Start DA</p> | <ul style="list-style-type: none"> – Move needle to sample position – Move lift with needle to working height – Switch over injection valve on 838 to "Fill" – Fill sample loop with sample – Switch over injection valve on 838 to "Inject" – Start data recording |
|---|---|

3.5.11 Connection to 761 Compact IC or 790 Personal IC

The compact 761 Compact IC is also controlled from the PC. The 838 Advanced Sample Processor is controlled via an RS232 connection. If necessary, Compact IC and Sample Processor may also be synchronised via a remote connection.





Figure 39 Interconnection with the 761 Compact IC via PC (RS232 only)

838-specific settings in the «IC Net 2.3» program

Program for the 838 Advanced Sample Processor:

| | | | |
|-----|------|--------|---|
| 001 | Move | Sample | – Move needle to sample position – Move lift with needle to working height – Fill sample loop with sample |
| 002 | Lift | work | |
| 003 | Pump | 240 s | |

Time program for 761



Switchover of the injection valve is programmed in the time program of the 761 Compact IC. The times may need to be adapted for the application.

The functions of the 838 Advanced Sample Processor may be triggered time-controlled via remote connection for more complex sample preparation steps. A simple example is shown below.





Figure 40 Interconnection with compact systems via remote connection (PC-controlled)

838-specific settings in the «IC Net 2.3» program

Program for the 838 Advanced Sample Processor:


| | | | | |
|-----|------|--------|--|---|
| 001 | Scan | Wait1 | | - Wait for pulse on remote lead |
| 002 | Move | Sample | | - Move needle to sample position |
| 003 | Lift | work | | - Move lift with needle to working height |
| 004 | Scan | Wait1 | | - Wait for pulse on remote lead |
| 005 | Pump | 240 s | | - Fill sample loop with sample |

Time program for 761



Remote pulses are programmed time-controlled in the time program of the 761 Compact IC. The times may need to be adapted to the application.

3.6 Connection to Metrohm VA systems



If the 838 Advanced Sample Processor is operated with a Metrohm VA system, the predefined sequence methods of the 838 Advanced Sample Processor are used. It is absolutely essential to connect the 6.2142.050 keypad in order to adapt these methods.

The possible instrument combinations depend greatly on the individual determination methods in voltammetry. Two systems are shown by way of example below.

3.6.1 Connection to a Metrohm 797 VA Computrace

The 838 Advanced Sample Processor is connected to a Metrohm 797 VA Computrace with the 6.2141.180 multiple remote cable, see below. Connection of pumps to the 731 Relay Box requires one 6.2160.010 adapter cable.




Figure 41 Interconnection with 797 VA Computrace

3.6.2 MVA-13 modular system for automatic analysis of electroplating baths

The 838 Advanced Sample Processor is connected to a Metrohm 797 VA Computrace with the multiple 6.2141.180 remote cable, see below. Connection of pumps to the 731 Relay Box requires one 6.2160.010 adapter cable in each case.




Figure 42 MVA-13 connection to an analysis system for electroplating baths

3.6.3 Accessories for sample volumes greater than 10 ml




Figure 43 Accessories for larger sample volumes

The 6.1835.040 PEEK tubing, the 6.2833.000 needle holder and the 6.1831.070 PEEK capillary can be used for sample transfer of small volumes. Please refer to the 797 VA Computrace Instructions for Use or the documentation of the MVA-13 system for further information.

3.7 Sample racks

A sample rack is a turntable for accommodating sample vessels which is fitted on the Sample Processor. Various types of sample rack can be used and exchanged easily in order to support different types of sample vessel. The rack offers space for different numbers of samples depending on the diameter of the sample vessels.

3.7.1 Metrohm standard sample racks

| Article No. / Rack type | Number of samples | Type of sample vessel | Vessel diameter | Magnetic code predefined |
|-------------------------|-------------------|--|-----------------|--------------------------|
| 6.2041.310 | 12 | 250 ml Metrohm titration beaker | 65 mm | 000001 |
| 6.2041.320 | 16 | 150 ml beaker | 55 mm | 000010 |
| 6.2041.340 | 24 | 75 ml Metrohm titration beaker | 35 mm | 001000 |
| 6.2041.350 | 48 | 75 ml Metrohm titration beaker | 35 mm | 010000 |
| 6.2041.360 | 12 | 150 ml beaker or 200 ml disposable beaker | 55 mm | 100000 |
| 6.2041.370 | 14 | 200 ml disposable beaker | 55 mm | 000011 |
| 6.2041.380 | 14 | 8 oz disposable beaker | 59 mm | 000101 |
| 6.2041.400 * | 126 + 2 | 11 ml sample tube 250 ml Metrohm rinsing beaker | 16 mm 65 mm | 001010 |
| 6.2041.410 * | 141 + 1 | 11 ml sample tube 500 ml beaker | 16 mm 71 mm | 001010 |
| 6.2041.430 * | 127 + 2 | 11 ml sample tube 2 x 300 ml PE bottles | 16 mm 68 mm | 010001 |
| 6.2041.440 * | 148 + 3 | 11 ml sample tube 3 x 300 ml PE bottles | 16 mm 68 mm | 010100 |
| 6.2041.450** | 56 + 56 | 50 ml PP sample vessel 11 ml sample tube | 30 mm 16 mm | 100100 |
| 6.2041.750 * | 36 | 11 ml sample tube | 16 mm | 011000 |

* recommended for IC applications

** recommended for VA applications

Table 5 Sample racks

On request, other user-defined racks can be supplied and defined via PC software in the instrument. Any arrangements of the sample positions are possible.

3.7.2 Magnetic codes

Each individual sample rack can be uniquely identified by means of a magnetic code. Magnetic pins on the underside of the rack (see Figure 8) can be combined to form a binary, six-digit code. The Sample Processor can thus automatically detect what rack is fitted.

When changing the rack, it should first be moved to initial position by pressing the **<RACK>** key (if a keypad is connected) or the instrument should be switched off and then on again. This allows unique detection of the rack and, thus, correct beaker positioning. An internal position table defining the angle of rotation and the distance from the centre of the rack is assigned to each rack type.

When a sample series is started, the Sample Processor automatically moves the rack first to initial position in order to read out the magnetic code. This always ensures that the beaker positions correspond to the internal position table of the fitted rack.

The standard racks supplied by Metrohm already feature a predefined magnetic code for each type. If several racks of the same type are used, the magnetic pins may be arranged differently so as to allow unique identification of a sample rack.

Magnetic code format (example):

000001, i.e. only one magnet is fitted, bit 0

000101, i.e. two magnets are fitted, bits 0 and 2

63 different combinations are possible. Code 000000 means "no code defined".

4 Service, maintenance, faults

4.1 Servicing

Servicing the Advanced Sample Processor should be carried out within the framework of an annual service visit by technicians from Metrohm or one of its local agencies. If aggressive or corrosive chemicals are used then shorter service intervals are necessary.

The Metrohm service department is always pleased to provide competent advice about the servicing and maintenance of all Metrohm instruments.

4.2 Care and maintenance

Not just sensitive instruments, but even a Sample Processor requires a certain amount of care and attention. If the instrument becomes excessively dirty this could interfere with its functions and shorten the working life of its basically robust mechanism and electronics.

Excessive dirt on the working heads can influence the measuring results obtained. Regular cleaning of exposed parts can prevent this to a large extent.

Spilled chemicals and solvents should be cleaned up immediately. Above all, the connection strip (the mains plug in particular) must be protected against contamination.

Although this prevented to a large extent by constructive measures, if aggressive media should penetrate the interior of the instrument then the mains plug should be pulled out immediately in order to prevent massive damage to the electronics. If such damage should occur please contact the Metrohm service department.

The instrument must not be opened by untrained personnel.

4.2.1 Exchanging pump tubings

Pump tubings are consumable material with a restricted service life. You should thus periodically replace pump tubings (approx. every 2 weeks if used continuously).

The service life of pump tubings depends greatly on the contact pressure. You should thus set the contact pressure correctly in accordance with *Chapter 3.3.9* and fully raise the tubing cartridges by releasing the snap-action lever **51** if the pump is to be switched off for a long period (this maintains the optimum contact pressure once set).

Proceed as follows in order to exchange pump tubing:

1 Removing the old pump tubing

- Press the contact pressure lever **49** on the tubing cartridge down fully.
- Detach the tubing cartridge **48** from the retaining bracket **50** by pushing the snap-action lever **51** in and disengage the tubing cartridge from the retaining cam **53** on the 838 Advanced Sample Processor.
- Remove the old pump tubing.

2 Fitting the new pump tubing

- Insert the new pump tubing **56** (6.1826.xxx) into the tubing cartridge. The stopper **57** must engage in the corresponding fixture on the upper side of the tubing cartridge.
- Engage the tubing cartridge in the retaining cam **53** and push it at the rear to the front until the snap-action lever **51** engages on the retaining bracket **50**. Ensure that the pump tubing is not kinked when doing this.

3 Setting the contact pressure

- Push the pressure lever **49** upwards until the solution is just being aspirated. Then push the pressure lever 1 further detent position up to achieve an optimum contact pressure.



Pump tubings made of PVC or PP may not be used to rinse with solutions containing acetone. Use different pump tubings or use another pump for rinsing in this case.

4.2.2 Replacing the filtration membrane

Perfect condition of the filtration membrane used is absolutely essential in order to achieve constant quality of the analysis results. Consequently, this membrane must be exchanged at regular intervals.

Please follow the information in *Chapter 6.1* when assessing the need to exchange a membrane.

Please proceed as follows (see also *Chapter 3.3.6*) when replacing a filtration membrane:

1 Preparing the ultrafiltration cell

- Detach all capillary connections on the ultrafiltration cell **28** by unscrewing the PVDF pressure screws **22**.
- Using the 6.2621.070 Allen key (5 mm) fully undo the 5 screws **46**, separate the upper section **43** from the lower section **45** and remove the sealing ring **44** and the spent filtration membrane.
- Thoroughly rinse off the sealing ring, lower section and upper section of the ultrafiltration cell with ultra-pure water.



Only ultra-pure water or ethanol may be used to clean the ultrafiltration cell. Organic solvents (e.g. acetone) will result in damage to the Plexiglass cell!

2 Preparing the filtration membrane

- Remove a new 6.2714.020 filtration membrane **50** using the 6.2831.010 tweezers from the packaging and immerse it for approx. 2 minutes in a petri dish filled with ultra-pure water until the membrane is fully soaked with water.

3 Fitting the filtration membrane

- Place the upper section **43** of the ultrafiltration cell with the inside pointing upwards onto a clean paper towel.
- Fit the sealing ring **44** into the recess provided on the upper section.
- Place the wet filtration membrane using the tweezers, centred inside the sealing ring, onto the upper section.



It is imperative that the filtration membrane be soaked with water before inserting it into the filtration cell since this causes it to expand. If it is inserted dry, this will then lead to creasing and folding in the filtration cell which may result in blockage of the cell.



Please use only the specified **6.2744.000 PVDF pressure screws** for connection to the ultrafiltration cell. If you use the 6.2744.010 PEEK pressure screws, this may result in stress cracking on the ultrafiltration cell!

4 Closing the ultrafiltration cell

- Fit the lower section **45** onto the upper section **43** so that the two parts fit perfectly.
- Using the 6.2621.070 Allen key, fully screw in and firmly tighten the 5 screws **46**.

5 Fitting the ultrafiltration cell

- After screwing the ultrafiltration cell together, insert it into the holder **26** as shown in *Figure 17* so that the heads of the screws **46** are located in the holes provided in the holder.
- Fit the three capillary connections which were detached previously back on the inlet and outlet openings provided on the ultrafiltration cell using the PVDF pressure screws **22**.

6 Rinsing the ultrafiltration cell

- Each time the ultrafiltration membrane is replaced, any air still in the filtration cell and in the lines must be removed. In order to do this, rinse all lines with ultra-pure water from one of the special beakers for instance:
- Place a 300 ml PE bottle with ultra-pure water as the rinsing solution onto a special beaker position on the sample rack.
- Use the IC Net software to control the 838 Advanced Sample Processor directly. Choose the corresponding system and open the **Manual control** dialog box. You can run the individual functions of the Sample Processor under **Autosampler**. Please refer to the IC NET software Online Help documentation for details of operation.
- Move the beaker with the rinsing solution above the lift with **MOVE**.
- You can use the **LIFT** function to lower the needle into rinsing position.
- Switch on the peristaltic pump with **PUMP**, and all connected lines are then rinsed and any air they contain is removed.
- After an adequate time (e.g. 5 min), switch the pump back off again with **PUMP** and move the needle back to home position.

4.2.3 Replacing the dialysis membrane

It may be necessary to replace the dialysis membrane in the following cases:

- Reduced yield in dialysis
- Dried-out membrane or membrane damage as the result of deposits or bacterial growth
- Irremediable blockage of the sample channel (sample is no longer able to be pumped through the dialysis cell)

Proceed as follows in order to exchange the membrane (see also 3.3.7):

1 Detaching the dialysis cell

- Unscrew the four inlet resp. outlet capillaries from the dialysis cell **37** and remove the cell from the cell holder.
- Using the 6.2621.070 Allen key, fully undo the 5 screws **42** and separate the lower section **41** from the upper section **38**
- Remove the old dialysis membrane **40**.

2 Cleaning the dialysis cell, fitting the dialysis membrane and fitting the cell back on

- See *Chapter 3.3.7* Fitting the dialysis cell.

3 Connecting the dialysis cell

- See *Chapter 3.3.12* Connecting the dialysis cell.

4 Conditioning the dialysis membrane

See also *Chapter 3.3.12*.

- Immerse the two aspiratory tubings **68** and **72** (see Figure 24) into the acceptor solution (degassed ultra-pure water).
- Switch on the peristaltic pump in «IC Net».
- Rinse the dialysis system for approx. 20 minutes with ultra-pure water. Check whether solution is emerging uniformly at the two inlet capillaries into the waste vessel.
- Check all tubings from the supply vessels through the tubing cartridges and the dialysis cell as far as the waste vessels for escaping fluid. If fluid is escaping, at any point, the corresponding connection must be tightened or replaced.
- If air bubbles remain in the dialysis cell, unscrew the PEEK capillary **76** (acceptor solution) and PTFE capillary **73** (sample) from the outlets **44 resp. 47 of the** dialysis cell and wait until the air bubbles have disappeared. Then screw the tubings back onto the dialysis cell.
- Switch off the peristaltic pump.

5 GLP validation – diagnosis

5.1 Validation / GLP

GLP (Good Laboratory Practice) requires, among other things, that the precision and correctness of analytical instruments is checked at regular intervals by using SOPs (Standard Operating Procedures).

As this instrument is not a measuring instrument as such, we recommend that the Sample Processor is included in an analytical system and is validated as part of the all-embracing validation of the whole analytical system.

Checking the electronic and mechanical assemblies of Metrohm instruments can and should be undertaken within the framework of regular servicing by Metrohm technicians. All Metrohm instruments are equipped with start-up check routines which check that the relevant assemblies are functioning perfectly when the instrument is switched on. If no error message appears it can be assumed that the instrument is functioning properly. Metrohm also supplies its instruments with built-in diagnosis programs which allow the service technicians to check the functioning of particular assemblies should faults or malfunctions occur and to localize them. Diagnosis programs can also be incorporated in a validation method.

6 Annex

6.1 Filtration

6.1.1 Selection of possible sample types

Each filtration process using a filtration membrane with a small pore size could be subject to a membrane blockage.

The following table lists some types of sample that have been filtered with the Metrohm 6.2729.110 Ultrafiltration Cell with the standard 6.2714.020 Filtration membrane (0.15 μm) or a second filter membrane (0.2 μm) and then analyzed on a Metrohm IC system. The concentrations of the following 7 anions were determined: F, Cl, NO₂⁻, Br, NO₃⁻, HPO₄²⁻, SO₄²⁻.

| Sample type | Pore size of membrane | No. of samples per filter |
|------------------------------|-----------------------|---------------------------|
| Orange juice with fruit pulp | 0.15 | 40 |
| Surface water | 0.15 | 500 |
| Drinking water | 0.15 | 1000 |
| Ground water | 0.15 | 500 |
| Wastewater 1 | 0.15 | 1000 |
| Wastewater 2 | 0.15 | 130 |
| Wastewater 3 | 0.15 | 40 |
| Wastewater 4 | 0.15 | 80 |
| NaCl solution (1%) | 0.2 | 5000 |
| Schöninger digestion soln. | 0.2 | 100 |
| Acidic soil extracts | 0.2 | 1000 |
| Aqueous soil extracts | 0.2 | 200 |

Table 6 Filtering various samples

The given number of samples that can be filtered on a filter membrane without any loss of quality being observed are empirical values. They have been determined at Metrohm AG and by various customers and should be used as an orientation guide for estimating the application of the ultrafiltration cell for sample preparation. These values must be determined individually for each new application.

6.1.2 Filter working life

A reduction in the recovery rate when analyzing standard solutions can be used as a possible indicator for the early recognition of an impending blockage. Ideally these solutions should be made up in the sample matrix being analyzed.

This means that if a large number of samples are to be analyzed then it is advisable to measure standard solutions at regular intervals, for example after every 5th or 10th sample for sample with a high particle load. However, it is not possible to make any general prediction about the number of possible filtration processes. The development of the recovery rate with the number of samples may also be subject to large variations. While with one sample matrix the recovery rate may remain stable rate for many samples and then drop sharply, with a different sample matrix it may diminish slowly and continuously.

The time when the filter membrane should be replaced ultimately depends on the sample matrix and the specifications of the analytical method used. In our experience very fine and suspended particles in the sample matrix will block the filter membrane more rapidly than coarse particles, which tend to be carried past the membrane in the sample flow.

Instructions for replacing the filter membrane are given in Section 4.2.2.

6.1.3 Filter membrane selection

You can apply existing sample preparation rules to filtration with the Metrohm 6.2729.110 Ultrafiltration Cell. If you wish to use a different filter membrane from that supplied as standard then please note that, even if the particle size is known, the selection of a membrane with an appropriate pore size does not automatically produce the required results.

Some investigations have shown that the retention capability of normal filter membranes does not always correspond to their specified pore size. The following table shows the qualitative filtering effect of filter membranes with different nominal pore sizes. The tests were carried out on aqueous solutions containing silica particles with particles sizes 1.5 μm and 5 μm .

| Test solutions: silica particles in water | Pore size of filter membrane 1 | Effect |
|--|---|---------------------------|
| 0.5 %, 5 μm | 0.15 μm | no breakthrough |
| 0.5 %, 5 μm | 3 μm | no breakthrough |
| 0.5 %, 5 μm | 8 μm | no breakthrough |
| 0.5 %, 5 μm | 10 μm | breakthrough ² |
| 0.5 %, 5 μm | 12 μm | no breakthrough |
| 0.5 %, 1.5 μm | 0.15 μm | no breakthrough |
| 0.5 %, 1.5 μm | 3 μm | breakthrough |

⁽¹⁾ Nominal pore size according to information supplied by manufacturer.

⁽²⁾ All the membranes came from a single manufacturer except this one.

Table 7 Selection of the filtration membrane

Please also note that, because of their lack of thickness, the retention capability of filter membranes may be lower than that of thicker filters with the same nominal pore size. This should be taken into account when choosing a suitable filter membrane.

6.2 Technical data

6.2.1 2 Channel peristaltic pump

| | | |
|---|--|-----------------|
| <i>Speed</i> | ± 6.7/min ... ± 100/min adjustable in ±15 steps | |
| <i>Delivery rate (typical)</i> | <i>Tubing</i> | <i>Delivery</i> |
| <i>(with water, no back-pressure, dependent on contact pressure and tubing type, at 20/min apart from 6.1826.150 in the case of 67/min)</i> | 6.1826.010: | 1.4 ml/min |
| | 6.1826.020: | 3.8 ml/min |
| | 6.1826.030: | 0.4 ml/min |
| | 6.1826.040: | 0.8 ml/min |
| | 6.1826.070: | 2.6 ml/min |
| | 6.1826.130: | 1.4 ml/min |
| | 6.1826.150: | 23.8 ml/min |
| <i>Pressure</i> | max. 4 bar (0.4MPa) | |
| <i>Pumpable fluids</i> | Clear fluids, no solids | |
| <i>Pump tubing material</i> | PVC (Tygon® ST) PVC (Tygon® LFL) PP | |

6.2.2 Pump connections

| | |
|---|---------------------------|
| <i>Outlets</i> | 2 x M8 sockets: |
| <i>for 823 Membrane Pump or 772 Pump Unit</i> | U = 16 ± 1 V I ≤ 0.8 A |

6.2.3 Injection valve

| | |
|------------------|-----------------|
| <i>Positions</i> | Fill ... Inject |
|------------------|-----------------|

6.2.4 Keypad with LCD display

| | |
|-----------------------------|--------|
| <i>Number of characters</i> | 2 x 24 |
| <i>Character height</i> | 5 mm |
| <i>LED indicators</i> | 3 |
| <i>Number of keys</i> | 30 |

6.2.5 Interfaces

RS232 connector

For connection of a computer, 9-pin

Remote socket

Universal parallel interface, programmable for communication with external instruments, 22 signal lines (8 x Input, 14 x Output), TTL level








Figure 44 Remote interface



$t_p > 20 \text{ ms}$
approx. 50 k Ω
Pull-up



$t_p > 200 \text{ ms}$
 $V_{CE0} = 40 \text{ V}$
 $I_C = 20 \text{ mA}$
Open Collector



+5 V: maximum load = 20 mA

6.2.6 MSB connections

3 connection sockets

9-pin Mini-DIN sockets
- for 800/700 Dosino or 805 / 685 Dosimat
- for 801 Stirrer or 804 Titration Stand

6.2.7 Swing Head connection

Connection socket

9-pin Mini-DIN socket
- for 786 Swing Head

6.2.8 Lift

Max. lift travel
Max. load
Lift rate

235 mm
approx. 30 N
adjustable, 5...25 mm/s

6.2.9 Turntable

Shifting rate

adjustable, 3...20 angular degrees/s

6.2.10 Stirrer connection (DIN socket)

Stirring speed

adjustable in 15 steps
741 magnetic stirrer 180/min...2600/min
802 / 722 rod stirrer 180/min...3000/min

6.2.11 Mains connection

Voltage
Frequency
Power consumption
Fuse

100... 240 V ($\pm 10\%$)
50...60 Hz
115 W
2.0 A slow-blow

6.2.12 Safety specification

| | |
|---------------------------|---|
| <i>Design and testing</i> | in accordance with EN/IEC 61010-1, EN/IEC 61010-2-081, UL 3101-1 safety class I |
| <i>Safety notes</i> | The Instructions for Use contain information and warnings which must be followed by the user in order to ensure safe operation of the instrument. |

6.2.13 Electromagnetic compatibility (EMC)

| | |
|-----------------------------|---|
| <i>Emitted interference</i> | Standards met: - EN/IEC 61326 - EN 55022 / CISPR 22 - EN/IEC 61000-3-2 |
|-----------------------------|---|

| | |
|------------------------------|---|
| <i>Interference immunity</i> | Standards met: - EN/IEC 61326 - EN/IEC 61000-4-2 - EN/IEC 61000-4-3 - EN/IEC 61000-4-4 - EN/IEC 61000-4-5 - EN/IEC 61000-4-6 - EN/IEC 61000-4-8 - EN/IEC 61000-4-11 - EN/IEC 61000-4-14 - NAMUR |
|------------------------------|---|

6.2.14 Ambient temperature

| | |
|--------------------------------|---|
| <i>Nominal operating range</i> | +5...+45 °C (at max. 80 % relative humidity) |
| <i>Storage</i> | -20 ... +60 °C Relative humidity < +40 °C < 95 % Relative humidity < +50 °C < 85 % Relative humidity < +60 °C < 50 % |
| <i>Transport</i> | -20 ... +60 °C Relative humidity < +40 °C < 95 % Relative humidity < +50 °C < 85 % Relative humidity < +60 °C < 50 % |

6.2.15 Dimensions and material

| | |
|---------------|--------|
| <i>Height</i> | 0.73 m |
| <i>Width</i> | 0.28 m |
| <i>Depth</i> | 0.50 m |

| | |
|------------------------------------|-----------------------|
| <i>Weight</i> | 15.50 kg (1.838.0010) |
| <i>(not including accessories)</i> | 15.90 kg (1.838.0020) |

Materials

| | |
|---------------------------------|---|
| - <i>Sample changer housing</i> | Metal housing, surface-treated |
| - <i>Keypad housing</i> | Crastin (PBTB), aluminium-sputtered on inside |
| - <i>Keypad film</i> | Polyester, chemical-resistant |

6.3 Standard equipment






On receipt of the instrument please check that the delivery is complete.









6.3.1 Metrohm 838 Advanced Sample Processor

Order no. 2.838.0010

Advanced IC Sample Processor

The following accessories are included as standard:








| No. | Order no. | Description |
|-----|--|--|
| 1 | 1.838.0010 1.786.0010 6.1462.090 | 838 Advanced Sample Processor For Metrohm IC systems <ul style="list-style-type: none"> • 1 tower • Swing Head • Robotic arm with needle adapter • Peristaltic pump • no injection valve • no sample rack Sample racks see 6.3.2.  |
| 3 | 6.1608.080 | PE bottle 300 ml (40 mm opening)  |
| 2 | 6.1826.130 | Pump tubing LFL White/White |
| 1 | 6.1831.050 | PEEK capillary 1/16 " d = 0.5 mm, L = 0.40 m  |
| 2 | 6.1831.060 | PEEK capillary 1/16 " d = 0.5 mm, L = 1 m  |
| 1 | 6.1835.010 | PEEK needle 1/8"  |










| | | | |
|---|------------|--|---|
| 1 | 6.2064.000 | Retention plate |  |
| 1 | 6.2134.100 | RS-232 cable DB9/f — DB9/f (5 m) |  |
| 1 | 6.2141.190 | Remote cable 732/819 — 766/788/838 | |
| 1 | 6.2621.300 | Tool kit for 838 |  |
| 1 | 6.2739.000 | Spanner for nipple |  |
| 1 | 6.2744.010 | Pressure screw made of PEEK with UNF 10/32 connection, 5 ea. |  |
| 1 | 6.2744.034 | Coupling, olive/UNF 10/32 2 ea. |  |
| 1 | 6.2744.160 | Coupling with lock Olive/UNF 10/32 |  |
| 1 | 6.2751.110 | Safety guard/splash guard |  |
| 1 | Y.1070.150 | Cable strap, self-adhesive | |
| 1 | 6.2122.xxx | Mains cable with coupling, Type IEC 320/C13 Cable connector to customer specifications: - Type SEV 12 (Switzerland...) 6.2122.020 - Type CEE(7), VII (Germany...) 6.2122.040 - Type NEMA/ASA (USA...) 6.2122.070 | |
| 1 | 8.838.1313 | Installation Instructions, 838 Advanced Sample Processor, English | |

Order No. 2.838.0020

Advanced IC Injection Sample Processor

The following accessories are included in the scope of delivery:







| Qty. | Order No. | Description |
|------|--|--|
| 1 | 1.838.0020 1.786.0010 6.1462.090 | 838 Advanced Sample Processor For Metrohm IC systems <ul style="list-style-type: none"> • 1 tower • Swing Head • Robotic arm with needle adapter • Peristaltic pump • Injection valve • Without sample rack Sample racks, see 6.3.2. |
| | |  |
| 3 | 6.1608.080 | PE bottle 300 ml (40 mm opening) |
| | |  |
| 1 | 6.1825.210 | Sample loop 20 μ L made of PEEK |
| | |  |
| 2 | 6.1826.130 | Pump tubing LFL White/White |
| 1 | 6.1831.010 | PEEK capillary 1/16 " d = 0.25 mm, L = 3.0 m |
| | |  |
| 1 | 6.1831.050 | PEEK capillary 1/16 " d = 0.5 mm, L = 0.40 m |
| | |  |
| 2 | 6.1831.060 | PEEK capillary 1/16 " d = 0.5 mm, L = 1 m |
| | |  |
| 1 | 6.1835.010 | PEEK needle 1/8" |
| | |  |











| | | | |
|---|------------|--|---|
| 1 | 6.2064.000 | Retention plate |  |
| 1 | 6.2134.100 | RS-232 cable DB9/f — DB9/f (5 m) |  |
| 1 | 6.2141.190 | Remote cable 732/819 — 766/788/838 | |
| 1 | 6.2621.300 | Tool kit for 838 |  |
| 1 | 6.2739.000 | Spanner for nipple |  |
| 1 | 6.2744.010 | Pressure screw made of PEEK with UNF 10/32 connection, 5 ea. |  |
| 1 | 6.2744.034 | Coupling, olive/UNF 10/32 2 ea. |  |
| 1 | 6.2744.110 | Coupling UNF 10/32 external / olive |  |
| 1 | 6.2744.160 | Coupling with lock Olive/UNF 10/32 |  |
| 1 | 6.2751.110 | Safety guard/splash guard |  |
| 1 | Y.1070.150 | Cable strap, self-adhesive | |
| 1 | 6.2122.xxx | Mains cable with coupling, Type IEC 320/C13 Cable connector to customer specifications: - Type SEV 12 (Switzerland...) 6.2122.020 - Type CEE(7), VII (Germany...) 6.2122.040 - Type NEMA/ASA (USA...) 6.2122.070 | |
| 1 | 8.838.1313 | Installation Instructions, 838 Advanced Sample Processor, English | |







Order No. 2.838.0110

Advanced IC Dialysis Sample Processor

The following accessories are included in the scope of delivery:

| Qty. | Order No. | Description |
|------|--|---|
| 1 | 1.838.0010 1.786.0010 6.1462.090 | 838 Advanced Sample Processor For Metrohm IC systems <ul style="list-style-type: none"> • 1 tower • Swing Head • Robotic arm with needle adapter • Peristaltic pump • Without injection valve • Without sample rack Sample rack, see 6.3.2. |
| | |  |
| 1 | 6.1602.150 | Bottletop GL 45-3x10/32 |
| | |  |
| 3 | 6.1608.070 | Transparent glass bottle 2 L GL45 |
| | |  |
| 3 | 6.1608.080 | PE bottle 300 ml (40 mm opening) |
| | |  |
| 6 | 6.1803.040 | PTFE capillary 1/16 " d = 0.5 mm, L = 1.0 m |
| | |  |
| 1 | 6.1803.100 | PTFE capillary 1/16 " d = 0.5 mm, L = 0.40 m |
| | |  |
| 2 | 6.1826.030 | Pump tubing AME-05 Orange/Yellow |
| 2 | 6.1826.040 | Pump tubing AME-08, Black/Black |





| | | | |
|---|------------|--|---|
| 1 | 6.1831.050 | PEEK capillary 1/16 " d = 0.5 mm, L = 0.40 m |  |
| 1 | 6.1831.060 | PEEK capillary 1/16 " d = 0.5 mm, L = 1 m |  |
| 1 | 6.1835.010 | PEEK needle 1/8" |  |
| 1 | 6.2057.010 | Holder for dialysis cell 6.2729.100 |  |
| 1 | 6.2064.000 | Retention plate |  |
| 1 | 6.2134.100 | RS-232 cable DB9/f — DB9/f (5 m) |  |
| 1 | 6.2141.190 | Remote cable 732/819 — 766/788/838 | |
| 1 | 6.2621.300 | Tool kit for 838 |  |
| 1 | 6.2714.010 | Dialysis membrane made of cellulose acetate 0.2 μm, 50 ea. |  |
| 1 | 6.2729.100 | Dialysis cell made of PMMA for the 754 Pump Unit with 4 x stopper 6.2744.060 |  |
| 1 | 6.2739.000 | Spanner for nipple |  |











| | | | |
|---|------------|--|---|
| 1 | 6.2744.000 | Pressure screw made of PVDF, UNF 10/32 5 ea. |  |
| 1 | 6.2744.010 | Pressure screw made of PEEK with UNF 10/32 connection, 5 ea. |  |
| 1 | 6.2744.034 | Coupling, olive/UNF 10/32 2 ea. |  |
| 1 | 6.2744.160 | Coupling with lock Olive/UNF 10/32 |  |
| 1 | 6.2751.110 | Safety guard/splash guard |  |
| 1 | 6.2831.010 | PBTP tweezers |  |
| 1 | Y.1070.150 | Cable strap, self-adhesive | |
| 1 | 6.2122.xxx | Mains cable with coupling, Type IEC 320/C13 Cable connector to customer specifications: - Type SEV 12 (Switzerland...) 6.2122.020 - Type CEE(7), VII (Germany...) 6.2122.040 - Type NEMA/ASA (USA...) 6.2122.070 | |
| 1 | 8.838.1313 | Installation Instructions, 838 Advanced Sample Processor, English | |










Order No. 2.838.0120




Advanced IC Dilution Sample Processor

The following accessories are included in the scope of delivery:

| Qty. | Order No. | Description |
|------|--|--|
| 1 | 1.838.0020 1.786.0010 6.1462.090 | <p>838 Advanced Sample Processor For Metrohm IC systems</p> <ul style="list-style-type: none"> • 1 tower • Swing Head • Robotic arm with needle adapter • Peristaltic pump • With injection valve • Without sample rack <p>Sample racks, see 6.3.2.</p>  |
| 1 | 6.1414.070 | <p>Titration vessel upper section</p>  |
| 1 | 6.1456.210 | <p>Mixing vessel 10 ... 90 ml</p>  |
| 3 | 6.1608.070 | <p>Transparent glass bottle 2 l GL45</p>  |

| | | | |
|---|------------|---|---|
| 3 | 6.1608.080 | PE bottle 300 ml (40 mm opening) |  |
| 1 | 6.1619.010 | Absorber |  |
| 1 | 6.1803.060 | PTFE capillary 1/16 " d = 0.97 mm, L = 0.20 m |  |
| 1 | 6.1803.070 | PTFE capillary 1/16 " d = 0.97 mm, L = 0.40 m |  |
| 1 | 6.1803.080 | PTFE capillary 1/16 " d = 0.97 mm, L = 1.0 m |  |
| 1 | 6.1803.110 | PTFE capillary 1/16 " d = 0.97 mm, L = 0.50 m |  |
| 1 | 6.1803.120 | PTFE capillary 1/16 " d = 0.97 mm, L = 3.0 m |  |
| 2 | 6.1826.140 | Pump tubing LFL Grey/Grey | |
| 1 | 6.1835.010 | PEEK needle 1/8" |  |
| 1 | 6.2001.070 | Stand bracket for the sample changers 778/789/814/815/838 |  |
| 1 | 6.2013.010 | Adjusting ring for stand rods 10 mm |  |







| | | | |
|---|------------|---|---|
| 1 | 6.2036.000 | Retaining ring for titration vessel |  |
| 1 | 6.2064.000 | Retention plate |  |
| 1 | 6.2134.100 | RS-232 cable DB9/f — DB9/f (5 m) |  |
| 1 | 6.2141.190 | Remote cable 732/819 — 766/788/838 | |
| 1 | 6.2621.300 | Tool kit for 838 |  |
| 1 | 6.2730.030 | Stopper with nipple M10 |  |
| 1 | 6.2739.000 | Spanner for nipple |  |
| 1 | 6.2744.010 | Pressure screw made of PEEK with UNF 10/32 connection, 5 ea. |  |
| 1 | 6.2744.034 | Coupling, olive/UNF 10/32 2 ea. |  |
| 1 | 6.2744.060 | Threaded stopper POM UNF 10/32 |  |









| | | | |
|---|------------|--|---|
| 1 | 6.2744.080 | PEEK coupling, M6 external/UNF |  |
| 1 | 6.2744.160 | Coupling with lock Olive/UNF 10/32 |  |
| 1 | 6.2751.110 | Safety guard/splash guard |  |
| 1 | Y.1070.150 | Cable strap, self-adhesive | |
| 1 | 6.2122.xxx | Mains cable with coupling, Type IEC 320/C13 Cable connector to customer specifications: - Type SEV 12 (Switzerland...) 6.2122.020 - Type CEE(7), VII (Germany...) 6.2122.040 - Type NEMA/ASA (USA...) 6.2122.070 | |
| 1 | 8.838.1313 | Installation Instructions, 838 Advanced Sample Processor, English | |







Order No. 2.838.0210

Advanced IC Ultrafiltration Sample Processor

The following accessories are included in the scope of delivery:

| Qty. | Order No. | Description |
|------|--|---|
| 1 | 1.838.0010 1.786.0010 6.1462.090 | <p>838 Advanced Sample Processor For Metrohm IC systems</p> <ul style="list-style-type: none"> • 1 tower • Swing Head • Robotic arm with needle adapter • Peristaltic pump • Without injection valve • Without sample rack <p>Sample racks, see 6.3.2.</p>  |
| 3 | 6.1608.080 | <p>PE bottle 300 ml (40 mm opening)</p>  |
| 1 | 6.1803.050 | <p>PTFE capillary 1/16 " d = 0.5 mm, L = 0.20 m</p>  |
| 1 | 6.1803.060 | <p>PTFE capillary 1/16 " d = 0.97 mm, L = 0.20 m</p>  |
| 1 | 6.1803.070 | <p>PTFE capillary 1/16 " d = 0.97 mm, L = 0.40 m</p>  |
| 2 | 6.1803.080 | <p>PTFE capillary 1/16 " d = 0.97 mm, L = 1.0 m</p>  |
| 2 | 6.1826.010 | <p>Pump tubing ENE-11 White/White</p> |
| 2 | 6.1826.030 | <p>Pump tubing AME-05 Orange/Yellow</p> |
| 2 | 6.1826.070 | <p>Pump tubing ENE-16 Yellow/Yellow</p> |









| | | | |
|---|------------|--|---|
| 2 | 6.1831.060 | PEEK capillary 1/16 " d = 0.5 mm, L = 1 m |  |
| 1 | 6.1835.010 | PEEK needle 1/8" |  |
| 1 | 6.2057.030 | Filtration cell holder |  |
| 1 | 6.2064.000 | Retention plate |  |
| 1 | 6.2134.100 | RS-232 cable DB9/f — DB9/f (5 m) |  |
| 1 | 6.2141.190 | Remote cable 732/819 — 766/788/838 | |
| 1 | 6.2621.300 | Tool kit for 838 |  |
| 1 | 6.2714.020 | Filtration membranes 0.2 μm, 50 ea. | |
| 1 | 6.2729.110 | Ultrafiltration cell |  |
| 1 | 6.2739.000 | Spanner for nipple |  |









| | | | |
|---|------------|--|---|
| 1 | 6.2744.000 | Pressure screw made of PVDF, UNF 10/32 5 ea. |  |
| 1 | 6.2744.010 | Pressure screw made of PEEK with UNF 10/32 connection, 5 ea. |  |
| 1 | 6.2744.030 | Coupling, olive/UNF 10/32 4 ea. |  |
| 1 | 6.2744.160 | Coupling with lock Olive/UNF 10/32 |  |
| 1 | 6.2751.110 | Safety guard/splash guard |  |
| 1 | 6.2831.010 | PBTP tweezers |  |
| 1 | Y.1070.150 | Cable strap, self-adhesive | |
| 1 | 6.2122.xxx | Mains cable with coupling, Type IEC 320/C13 Cable connector to customer specifications: - Type SEV 12 (Switzerland...) 6.2122.020 - Type CEE(7), VII (Germany...) 6.2122.040 - Type NEMA/ASA (USA...) 6.2122.070 | |
| 1 | 8.838.1313 | Installation Instructions, 838 Advanced Sample Processor, English | |









Order No. 2.838.0310

Advanced VA Sample Processor


The following accessories are included in the scope of delivery:

| Qty. | Order No. | Description | |
|------|--|--|---|
| 1 | 1.838.0010 1.786.0010 6.1462.090 | 838 Advanced Sample Processor For Metrohm VA systems <ul style="list-style-type: none"> • 1 tower • Swing Head • Robotic arm with needle adapter • Peristaltic pump • Without injection valve |  |
| 1 | 6.1618.020 | Threaded adapter S40 to GL45 |  |
| 1 | 6.1805.020 | FEP tubing connection M6 d = 2 mm, L = 0.52 m |  |
| 1 | 6.1805.060 | FEP tubing connection M6 d = 2 mm, L = 0.60 m |  |
| 1 | 6.1805.120 | FEP tubing connection M6 d = 2 mm, L = 1.0m |  |
| 1 | 6.1808.000 | Coupling sleeve 2 x M6 |  |
| 2 | 6.1808.020 | Tubing olive, external/M6 internal |  |
| 1 | 6.1819.010 | PTFE tubing |  |



| | | | |
|---|------------|---|---|
| 1 | 6.1822.420 | PEEK capillary 1x M6, 1/16" d = 0.75 mm | |
| 2 | 6.1826.020 | Pump tubing ENE-18 Blue/Blue | |
| 2 | 6.1826.150 | Pump tubing ENE-23 Violet/Orange | |
| 1 | 6.1831.070 | PEEK capillary 1/16" d = 0.75 mm, L = 0.70 m |  |
| 1 | 6.1831.080 | PEEK capillary 1/16" d = 0.75 mm, L = 1.10 m |  |
| 1 | 6.1835.040 | PEEK needle 1/8" | |
| 1 | 6.1835.050 | PEEK needle 1/8" d = 1.58 mm | |
| 1 | 6.2041.450 | Sample rack 56 x 15 ml + 56 x 50 ml |  |
| 1 | 6.2057.040 | Dosino holder |  |
| 1 | 6.2142.050 | Keypad for 838 |  |
| 1 | 6.2621.090 | Fork spanner 1/2" |  |
| 1 | 6.2621.100 | Allen key 3 mm |  |
| 1 | 6.2621.130 | Allen key 2 mm |  |

| | | | |
|---|------------|--|---|
| 1 | 6.2739.000 | Spanner for nipple |  |
| 1 | 6.2743.057 | PP sample vessels 11 ml 200 ea. |  |
| 1 | 6.2744.010 | Pressure screw made of PEEK with UNF 10/32 connection, 5 ea. |  |
| 2 | 6.2744.080 | Coupling M6 external/UNF 10/32 internal |  |
| 1 | 6.2744.160 | Coupling with lock Olive/UNF 10/32 |  |
| 4 | 6.2747.010 | PP sample vessels 50 ml 25 ea. |  |
| 1 | 6.2751.110 | Safety guard/splash guard |  |
| 1 | 6.2833.020 | Needle holder 1/8" to M6 |  |
| 2 | Y.1070.150 | Cable strap, self-adhesive | |
| 1 | 6.2122.xxx | Mains cable with coupling, Type IEC 320/C13 Cable connector to customer specifications: - Type SEV 12 (Switzerland...) 6.2122.020 - Type CEE(7), VII (Germany...) 6.2122.040 - Type NEMA/ASA (USA...) 6.2122.070 | |
| 1 | 8.838.1013 | Instructions for Use Metrohm 838 Advanced Sample Processor | |
| 1 | 8.838.1313 | Installation Instructions, 838 Advanced Sample Processor, English | |



6.3.2 Keyboard (optional)

| Order No. | Description | |
|------------|---|---|
| 6.4142.050 | Keyboard for the 838 Advanced Sample Processor |  1,5m |








6.3.3 Dosing drive and units









| Order No. | Description | |
|-------------------|---------------------------------------|---|
| 2.800.0010 | 800 Dosino |  |
| 6.3032.120 | 807 Dosing unit, 2 mL glass cylinder | |
| 6.3032.150 | 807 Dosing unit, 5 mL glass cylinder | |
| 6.3032.210 | 807 Dosing unit, 10 mL glass cylinder | |
| 6.3032.220 | 807 Dosing unit, 20 mL glass cylinder |  |
| 6.3032.250 | 807 Dosing unit, 50 mL glass cylinder | |

6.3.4 Stirrer

| Order No. | Description | |
|-------------------|---|---|
| 2.801.0040 | 801 Magnetic stirrer with support rod |  |
| 2.802.0020 | 802 Rod stirrer |  |
| 6.1909.020 | Incl. stirring propeller PP | |
| 6.1909.040 | or PTFE optional | |


6.3.5 Sample racks and vessels


| Order No. | Description | |
|---|---|---|
| 6.2041.310 6.1432.320 6.1453.220 6.1453.250 | Sample rack 12 x 250 ml (d = 416 mm) for Metrohm sample beakers, glass, 250 ml Metrohm sample beaker PP 200 ml Metrohm sample beaker PP 250 ml |  |
| 6.2041.320 | Sample rack 16 x 150 ml for Standard beakers (high type) | |
| 6.2041.340 6.1432.210 | Sample rack 24 x 75 ml (d = 416 mm) for Metrohm sample beakers, glass 75 ml, |  |
| 6.2041.350 6.1432.210 | Sample rack ^{*)} 48 x 75 ml (d = 416 mm) for direct titration to Metrohm sample beakers, glass, 75 ml |  |
| 6.2041.360 6.1459.310 | Sample rack 12 x 150 ml (d=416 mm) for 150 ml standard beakers (high type) 200 ml disposable beakers (EU) PP (1000 ea.) |  |
| 6.2041.370 6.1459.310 | Sample rack 14 x 200 ml (d = 416 mm) for 200 ml disposable beakers (EU) PP | |
| 6.2041.380 | Sample rack 14 x 8 oz (d = 416 mm) for disposable beakers (US) PP 8 oz | |
| 6.2041.400 6.2743.057 6.1432.320 6.1453.220 6.1453.250 | Sample rack 126 x 11 ml and 2 x 250 ml (d = 416 mm) for sample tubes 11 ml and Metrohm sample beakers, glass, 250 ml or Metrohm sample beakers, PP 200 ml or Metrohm sample beakers, PP 250 ml |  |
| 6.2041.430 6.2743.057 6.1608.080 | Sample rack 127 x 11 ml and 2 x 300 ml (d = 416 mm) for sample tubes 11 ml and 300 ml PE bottles |  |
| 6.2041.440 6.2743.057 6.1608.080 | Sample rack 148 x 11 ml and 3 x 300 ml (d = 416 mm) for sample tubes 11 ml and 300 ml PE bottles |  |

| | | |
|--|--|---|
| <p>6.2041.450 6.2743.057 6.2747.010</p> | <p>Sample rack 56x 11 ml + 56x 50 ml (d = 476 mm,) for sample tubes 11 ml and PP sample vessel 50 ml</p> |  |
| <p>6.2041.760 6.2743.057 6.1432.320 6.1453.220 6.1453.250</p> | <p>Sample rack 54x 11 ml + 1x 250 ml (d = 476 mm,) for sample tubes 11 ml and Metrohm sample beakers, glass, 250 ml or Metrohm sample beakers, PP 200 ml or Metrohm sample beakers, PP 250 ml</p> | |
| <p>6.2743.040</p> | <p>PP sample vessels 2.5 mL 2000 ea., not for use with ultra filtration cell</p> |  |
| <p>6.2743.047</p> | <p>PP sample vessels 2.5 mL 200 ea., not for use with ultra filtration cell</p> |  |
| <p>6.2743.050</p> | <p>PP sample vessels 11 mL 2000 ea.</p> |  |
| <p>6.2743.057</p> | <p>PP sample vessels 11 mL 200 ea.</p> |  |
| <p>6.2743.070</p> | <p>PP stopper with perforation 2000 ea. for 6.2743.xxx sample vessels</p> |  |
| <p>6.2743.075</p> | <p>PP stopper with perforation 200 ea. for 6.2743.xxx sample vessels</p> |  |
| <p>6.2747.010</p> | <p>PP sample vessels 50 ml 25 each</p> |  |










6.3.6 786 Swing Head

Order No. 2.786.0010 (model version "Transfer" with pipetting accessories)

| No. | Order no. | Description |
|-----------------|--------------------------|---|
| <p>1</p> | <p>1.786.0010</p> | <p>786 Swing Head Precision drive for robotic arms</p>  |

| | | | |
|---|------------|---|---|
| 1 | 6.1462.030 | Robotic arm for sample transfer, for left-swinging installation |  |
| 1 | 6.1562.100 | Pipetting tubing 10 ml | |
| 1 | 6.1823.010 | Guide shank for pipetting tubing | |
| 1 | 6.2621.120 | Allen key 1.5 mm | |
| 1 | 6.2751.100 | Splash guard | |
| 2 | V.024.4012 | Countersunk screws M.I hexagon | |










6.3.7 Connecting cables





| <i>Order No.</i> | <i>Description</i> | |
|-------------------|---|---|
| 6.2141.020 | Standard remote cable 25-pin/25-pin |  |
| 6.2134.090 | Remote cable 762/830 IC Interface — Device 1...4 or 5... |  |
| 6.2141.110 | Remote cable 766/788/813 — 709 or 762 and 732 |  |
| 6.2141.140 | Remote cable 766/788/813 — 817 |  |
| 6.2141.180 | Remote cable 797— 766/813/838 — 731, 25-pin/25-pin/25-pin |  |
| 6.2125.090 | Connecting cable RS232 Wired 1:1, 1m, 25-pin/m to 25-pin/m |  |
| 6.2125.110 | Connecting cable RS232 for connecting Metrohm instrument to PC 25-pin/f to 9-pin/f |  |
| 6.2134.040 | Connecting cable RS232 for connecting Metrohm instrument to PC 9-pin/f to 9-pin/f |  |
| 6.2125.120 | Remote adapter Distributor from remote lines to 5 x 25-pin sockets. |  |

6.3.8 Filtration and dialysis accessories




| <i>Order No.</i> | <i>Description</i> |
|-------------------|---|
| 6.2714.010 | Dialysis membrane for inline dialysis made of cellulose acetate, nominal pore size 0.2 μm , d = 47 mm, 50 ea. |
| 6.2714.020 | Filtration membrane made of cellulose acetate, nominal pore size 0.15 μm , d = 47 mm, 50 ea. |
| 6.2714.030 | Dialysis membrane for inline dialysis, made of poly- amide, nominal pore size 0.2 μm , d = 47 mm, 50 ea. |

6.3.9 Pump tubing

| Order No. | Description | |
|------------|---|---|
| 6.1826.010 | Pump tubing made of PVC (Tygon® ST), White/White Inner diam. = 1.02 mm ± 0.05 mm |  |
| 6.1826.020 | Pump tubing made of PVC (Tygon® ST), Blue/Blue Inner diam. = 1.65 mm ± 0.05 mm |  |
| 6.1826.020 | Pump tubing made of PVC (Tygon® ST), Blue/Blue Inner diam. = 1.65 mm ± 0.05 mm |  |
| 6.1826.030 | Pump tubing made of PVC (Tygon® ST), Or- ange/Yellow Inner diam. = 0.51 mm ± 0.05 mm |  |
| 6.1826.040 | Pump tubing made of PVC (Tygon® ST), Black/Black Inner diam. = 0.76 mm ± 0.05 mm |  |
| 6.1826.050 | Pump tubing made of PVC (Tygon® ST), White/Yellow Inner diam. = 0.57 mm ± 0.05 mm |  |
| 6.1826.060 | Pump tubing made of PP (PharMed®) Or- ange/Yellow Inner diam. = 0.51 mm ± 0.05 mm |  |
| 6.1826.070 | Pump tubing made of PVC (Tygon® ST), Yel- low/Yellow Inner diam. = 1.42 mm ± 0.05 mm |  |
| 6.1826.110 | Long-life pump tubing made of PVC (Tygon® LFL), Or- ange/Yellow Inner diam. = 0.51 mm ± 0.0102 mm |  |

| | | |
|-------------------|--|---|
| 6.1826.120 | Long-life pump tubing made of PVC (Tygon® LFL), Orange/White Inner diam. = 0.59 mm ± 0.05 mm |  |
| 6.1826.130 | Long-life pump tubing made of PVC (Tygon® LFL), White/White Inner diam. = 1.02 mm ± 0.0127 mm |  |
| 6.1826.140 | Long-life pump tubing made of PVC (Tygon® LFL), Grey/Grey Inner diam. = 0.51 mm ± 0.0102 mm |  |
| 6.1826.150 | Pump tubing made of PVC (Tygon® ST), Violet/Orange Inner diam. = 2.54 mm |  |

6.3.10 Needles and tubings

| <i>Order No.</i> | <i>Description</i> | |
|------------------|--|---|
| 6.1835.010 | PEEK needle 1/8" For aspirating solutions from closed sample vessels with perforated stopper |  L = 150 mm |
| 6.1835.020 | PEEK tube 1/8" For aspirating solutions from open sample vessels, for filtration and dialysis applications |  L = 150 mm |
| 6.1835.040 | PEEK sample aspiration tube 1/8" For aspirating solutions from open sample vessels, inner diameter 0.75 mm | |
| 6.1835.050 | PEEK sample aspiration tube 1/8" For aspirating solutions from open sample vessels, inner diameter 1.58 mm, sphenoidal tip | |
| 6.2624.000 | Steel needle, inner diameter 0.66 mm For aspirating solutions from open sample vessels |  L = 150 mm |

6.4 Warranty

The warranty on our products is limited to defects that are traceable to material, construction or manufacturing error which occur within 12 months from the day of delivery. In this case the defects will be rectified in our workshops free of charge. Transport costs are to be paid by the customer.

For day and night operation the warranty is limited to 6 months.

Glass breakage in the case of electrodes or other parts is not covered by the warranty. Checks which are not a result of material or manufacturing faults are also charged during the warranty period. For parts from outside manufacturers, insofar as these constitute an appreciable part of our instrument, the warranty stipulations of the manufacturer in question apply.

With the regard to the guarantee of accuracy the technical specifications in the instruction manual are authoritative.

Concerning defects in materials, construction or design as well as the absence of guaranteed features the purchaser has no rights or claims except those mentioned above.





If damage of the packaging is evident on receipt of a consignment or if the goods show signs of transport damage after unpacking, the carrier must be informed immediately and a written damage report demanded. Lack of an official damage report releases Metrohm Ltd. from any liability to pay compensation.

If any instruments and parts have to be returned then the original packaging should be used if at all possible. This applies above all to instruments, electrodes, burette cylinders and PTFE pistons. Before embedment in wood shavings or similar material the parts must be packed in a dustproof package (for instruments the use of a plastic bag is essential). If open assemblies are included that are sensitive to electromagnetic voltages (e.g. data interfaces, etc.) then these must be returned in the associated original protective packaging (e.g. conductive protective bag). (Exception: assemblies with a built-in voltage source belong in non-conductive protective packaging).

For damage which arises as a result of non-compliance with these instructions no warranty responsibility whatsoever will be accepted by Metrohm.

6.5 Declaration of Conformity

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

| | |
|--|--|
| <p>Name of commodity</p> <p>838 Advanced Sample Processor</p> |  <p>Metrohm l o n a n a l y s i s</p> <p>CH-9101 Herisau/Switzerland E-Mail info@metrohm.com www.metrohm.com</p> |
| <p><i>Description</i> Sample changer with advanced Liquid Handling abilities for the automation of batch processing of sample preparation in analytical laboratories.</p> | |
| <p>This instrument has been built and has undergone final type testing according to the standards:</p> <p><i>Electromagnetic compatibility: Emission</i> EN/IEC 61326, EN 55022 / CISPR 22, EN/IEC 61000-3-2</p> <p><i>Electromagnetic compatibility: Immunity</i> EN/IEC 61326, EN/IEC 61000-4-2, EN/IEC 61000-4-3, EN/IEC 61000-4-4, EN/IEC 61000-4-5, EN/IEC 61000-4-6, EN/IEC 61000-4-8, EN/IEC 61000-4-11, EN/IEC 61000-4-14, NAMUR</p> <p><i>Safety specifications</i> EN/IEC 61010-1, EN/IEC 61010-2-081, UL 3101-1 protection class I</p> <p>It has also been certified by ElectroSuisse, a member of the International Certification Body (CB/IEC).</p> | |
| <p> <i>The instrument meets the requirements of the CE mark as contained in the EU directives 89/336/EEC and 73/23/EEC and fulfils the following specifications:</i></p> <p>EN 61326 Electrical equipment for measurement, control and laboratory use – EMC requirements EN 61010-1 Safety requirements for electrical equipment for measurement, control and laboratory use EN 61010-2-081 Part 2-081: Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes.</p> | |
| <p>Metrohm Ltd. is holder of the SQS-certificate of the quality system ISO 9001 for quality assurance in design/development, production, installation and servicing.</p> <p>The system software, stored in Read Only Memories (ROMs) has been validated in connection with standard operating procedures in respect to functionality and performance. The technical specifications are documented in the instruction manual.</p> | |
| <p>Herisau, February 16, 2004</p> <div style="display: flex; justify-content: space-around;"> <div data-bbox="414 1702 630 1937">  <p>D. Strohm Vice President Head of R&D</p> </div> <div data-bbox="853 1713 1316 1971">  <p>Ch. Buchmann Vice President Head of Production Responsible for Quality Assurance</p> </div> </div> | |

6.5.1 Quality Management Principles

Metrohm AG, CH-9101 Herisau, Switzerland

 **Metrohm**
l o n a n a l y s i s
CH-9101 Herisau/Switzerland
E-Mail info@metrohm.com
Internet www.metrohm.com

Metrohm Ltd. holds the ISO 9001 Certificate, registration number 10872-02, issued by SQS (Swiss Association for Quality and Management Systems). Internal and external audits are carried out periodically to assure that the standards defined by Metrohm's QM Manual are maintained.

The steps involved in the design, manufacture and servicing of instruments are fully documented and the resulting reports are archived for ten years. The development of software for PCs and instruments is also duly documented and the documents and source codes are archived. Both remain the possession of Metrohm. A non-disclosure agreement may be asked to be provided by those requiring access to them.

The implementation of the ISO 9001 quality system is described in Metrohm's QM Manual, which comprises detailed instructions on the following fields of activity:

Instrument development

The organisation of the instrument design, its planning and the intermediate controls are fully documented and traceable. Laboratory testing accompanies all phases of instrument development.

Software development

Software development occurs in terms of the software life cycle. Tests are performed to detect programming errors and to assess the program's functionality in a laboratory environment.

Components

All components used in the Metrohm instruments have to satisfy the quality standards that are defined and implemented for our products. Suppliers of components are audited by Metrohm as the need arises.

Manufacture

The measures put into practice in the production of our instruments guarantee a constant quality standard. Production planning and manufacturing procedures, maintenance of production means and testing of components, intermediate and finished products are prescribed.

Customer support and service

Customer support involves all phases of instrument acquisition and use by the customer, i.e. consulting to define the adequate equipment for the analytical problem at hand, delivery of the equipment, user manuals, training, after-sales service and processing of customer complaints. The Metrohm service organisation is equipped to support customers in implementing standards such as GLP, GMP, ISO 900X, in performing Operational Qualification and Performance Verification of the system components or in carrying out the System Validation for the quantitative determination of a substance in a given matrix.

7 Index

<RACK> 60

731 Relay Box 56, 57

761 Compact IC 21, 40, 53

772 Pump Unit 35

797 VA Computrace 56, 57

820 IC Separation Centre 21, 33

823 Membrane Pump Unit 35

833 Dialysis Unit 31, 33

A

Acceptor solution 29, 31, 33

Accessories .. 76, 78, 81, 85, 88

Acetone 23, 62

Adapter cable (6.2160.010)

..... 56, 57

Air bubbles 33, 65

Allen key 17, 65

Ambient temperature 72

Aspirate 44

Aspirator tubing 23

Aspiratory capillary 23

Available versions 2

B

Bacterial growth 65

Beaker 59, 92

Binary code 60

Biohazard 10

Blind stoppers 45

Blockage of the membrane .. 68

Blockages 65

C

Cable (6.2141.020) 56

Cable (6.2141.170) 57

Cable (6.2141.190) 41, 42,

..... 46, 48

Capillary connections 44

Capillary cutter 30

Care 61

Catch pan 36

Cleaning the dialysis cell 20

Commissioning the peristaltic

pumps 33

Conditioning the dialysis

membrane 33

Conditioning the dialysis system

..... 33

Configuring the robotic arm.. 14

Connecting cable 38, 94

Connecting dosers 35

Connecting stirrers 35

Connecting the dialysis cell

..... 29, 31

Connecting the injection valve 32

Connecting the keypad 14

Connecting the Swing Head. 14

Connecting the ultrafiltration cell

..... 27

Connecting the waste container

..... 31, 32

Connection of a dialysis cell ...

..... 29, 30

Connection panel 7

Connection to IC systems 38

Connections for acceptor

solution 31

Contact pressure 23, 24, 33, 62

Coupling 23, 24, 31

D

Danger 4

Declaration of Conformity 98

Delivery 25

Delivery capillary 23

Delivery rate 24, 70

Depth 73

Dialysis 48

Dialysis accessories 93

Dialysis cell 20, 29

connection to the IC system

..... 65

detaching 65

Dialysis flow pathway 29

Dialysis membrane 20, 21

conditioning 65

replacement 65

Dialysis membrane 21

Dialysis program 50

Dimensions 73

Disposable beaker 59, 92

Dosimat 35

Dosing drives 35

Dosino 35, 44

Dosino holder 22

E

Electrical safety 9

Electromagnetic compatibility 72

Electrostatic charges 9

EMC 72

Emitted interference 72

Equipping the work head 45

Exchanging a pump tubing .. 62

External instruments 38

External pump 35

F

FILL 7, 33, 53

Filter unit (PEEK) 23

Filter useful life 68

Filtration

pore size 67

sample type 67

Filtration accessories 94

Filtration cell 27

Filtration cell – pump tubing

connection 28

Filtration cell - waste connection

..... 27

Filtration membrane.. 17, 18, 68

(6.2714.020) 67

Fitting pump tubing 23

Fitting the aspiratory capillary.

..... 30, 31

Fitting the dialysis cell 20

Fitting the Dosino 22

Fitting the Dosino holder 16

Fitting the filtration cell holder 16

Fitting the ultrafiltration cell ... 17

Frequency 71

Fuse 71

G

GLP 66

Good Laboratory Practice 66

Grounding 9

Guide chain 5, 34

H

Height 73

Holder 5

I

IC Interface 762 40

IC Net ... 18, 33, 41, 40, 42, 46,

..... 47, 48, 53, 54, 55, 64

Inadequate yield 65

INJECT 7

Injection loop 44

Injection valve 2, 5, 7, 29, 52

Inlet for acceptor solution 20

Inlet for sample solution.. 17, 20

Inlet opening 30

Input 39, 71

Input lines 38

- Inserting a dialysis membrane 20
 Inserting capillaries 34
 Inserting the filtration membrane 18
 Installation 11, 13
 Instrument models 2
 Instrument number 7
 Interconnection 38
 797 VA Computrace 56
 analysis system for electroplating baths (MVA-13) 57
 anionic system with chemical suppression (MIC-2) 41, 42, 46
 anionic system with chemical suppression and dialysis (MIC-7) 48
 modular anionic/cationic system with chemical suppression (MIC-3) 47
 modular system with 817 Bioscan for sugar analysis (MIC-8) 52
 Interface connections 38
 Interfaces 71
 Interference immunity 72
 ISO 9001 98
- K** _____
- Keypad 14, 56
 Keypad connection 7
- L** _____
- Leakage 33, 65
 Lower section 17, 20
- M** _____
- Magnet sensor 8
 Magnetic code 36, 59, 60
 Magnetic stirrer 71
 Mains connection 7, 10, 71
 Mains voltages 13
 Maintenance 61
 Manual control 18, 64
 Materials 73
 Maximum load 71
 Membrane pump 35
 Metrohm IC systems (MIC) ... 40
 Metrohm Serial Bus 7
 MIC – Metrohm IC systems .. 40
 Mini-DIN connector 14
 Mixing container 44, 45
 Model
 2.838.0010 41, 47, 48
 2.838.0020 52
 2.838.0120 42
 2.838.0210 46
 Models 2
- MSB connections 7, 71
 MSB sockets 35
- N** _____
- Needle holder 15, 16
 Needles 96
- O** _____
- Operating range 72
 Organic solvents 35
 Outlet for acceptor solution .. 20
 Outlet for filtrate 17
 Outlet for sample solution 17, 20
 Outlet opening 30, 31
 Output 39, 71
 Output lines 38
 Overall view 5
 Overview 1
- P** _____
- PEEK capillary 29, 32, 58
 PEEK coupling 23, 27, 30
 PEEK filter unit 24
 PEEK needle/tubing 15
 PEEK pressure screw .. 15, 16, 27, 30
 PEEK sample loop 29
 PEEK tubing 58
 Peristaltic pump .. 2, 5, 8, 35, 70
 Personal protection 10
 PharMed® 25
 Position table 60
 Power consumption 71
 Power switch 7
 PP 25, 62
 Preparing the dialysis cell 20
 Preparing the ultrafiltration cell 17
 Pressure lever 8, 23, 24, 33
 Pressure screw 15, 23, 31
 Printer 39
 Printer cable 39
 Program 41, 42, 46, 47, 49, 53, 54, 55
 PTFE capillary 15, 27, 29, 30
 Pump connections 6, 35, 70
 Pump drive 8
 Pump tubing 23, 24, 27, 29, 30, 31, 61, 94
 Pump tubing - filtration cell connection 27
 Pump tubing – injection valve connection 28
 Pump tubings
 service life 61
 Pumps 56, 57
 PVC 25, 62, 95, 96
 PVDF pressure screw . 18, 27, 29, 30, 63
- PVDF pressure screws 32
- Q** _____
- Quality Management Principles 99
- R** _____
- Rack 36
 Rack change 8
 Rack code 8
 Rack detection 60
 Rear panel 6
 Recovery rate 68
 Relative humidity 72
 Remote adapter (6.2125.120) 47
 Remote cable 41, 42, 46, 47, 57, 94
 Remote connection 7, 38
 Remote interface 39
 Replacing the dialysis membrane 65
 Replacing the filtration membrane 63
 Retaining bracket 8, 24
 Retaining cam 8
 Retention plate 5, 15
 Ring wedge 15
 Rinsing the ultrafiltration cell 18, 64
 Risk of injury 37
 Robotic arm 5, 14
 Rod stirrer 71
 RS232 39
 RS232 cable . 41, 42, 46, 47, 48
 RS232 connection 7
 RS232 interface 38
- S** _____
- Safety guard 5, 10, 37
 Safety notes 72
 Safety specification 72
 Sample beakers 92
 Sample dilution 42, 44
 Sample queue 51
 Sample rack 5, 36, 59, 92
 Sample transfer 58
 Sample tube 59, 92
 Scope of application 1
 Screw
 detaching the dialysis cell 65
 Sealing ring 17, 20
 Selection of the filter membrane 68
 Sequence methods 56
 Service 61
 Service life of pump tubing ... 23
 Setup 13
 Setup settings 14

Signals 39
 Simultaneous dialysis 49
 Snap-action lever 8, 23, 24
 Solvents 35
 SOP 66
 Splash guard 37
 Stand bracket 34
 Standard equipment 74
 Standard sample racks 59
 Standards 72
 Steel needle 96
 Stirrer connection 6
 Stirrer rail 5, 34
 Stirring speed 71
 Stopper 23, 24, 25
 Swing Head 5, 93
 Switch positions 7
 Switchover valve 44

T

Third-party instruments 38
 Threaded adapter 22
 Time program 54, 55
 Time-controlled 54
 Titration beaker 59
 Transmission parameters 39
 Transport damage 97
 Tubing cartridge 8, 23, 24
 Tubing connections 23
 Tubings 96
 Turntable 8
 Tweezers 18, 21
 Tygon® 25, 94, 96

U

Ultrafiltration cell 17, 27

Ultra-pure water 31, 33, 63
 Upper section 17, 20

V

VA systems 56
 Validation 66
 Vials 15
 Voltage 71
 Voltammetry 56

W

Warning 4, 6
 Warranty 97
 Waste container 27
 Weight 73
 Width 73