

Robotic Acid/Base Analyzer



Applications

8.855.2003

 **Metrohm**
Ion analysis

Metrohm AG
CH-9101 Herisau
Switzerland

Applications

Robotic Acid/Base Analyzer

 **Metrohm**
Ion analysis
CH-9101 Herisau/Switzerland
E-Mail info@metrohm.ch
www.metrohm.com



Symbols used in this document:

$c(X)$	Molar concentration of substance X in mol/L, often written as $[X]$
$M(X)$	Molar mass of substance or atomic mass (relative mass) of substance X in g/mol
$w(X)$	Mass fraction of substance X, e.g. $w(\text{NaOH}) = 40\%$
$\beta(X)$	Mass concentration of substance X, e.g. $\beta(\text{NaCl}) = 20 \text{ g/L}$
$a(X)$	Activity of substance X (only corresponds to the molar concentration in very dilute solutions; as the concentration increases, dissolved particles mutually influence each other so that their activity is lower than would be expected from the concentration).
$\text{p}K_p$	Autoprotolysis constant of a solvent
{ }	In this document braces such as used in the formula $\{\text{AgCl}\}$ indicate solid substances (precipitates) that do not consist of single molecules (in this case AgCl), but whose ionic components form an extended ionic lattice.

***tiamo* methods**

The methods described in this document are part of the *tiamo* installation on your PC system. To use these methods you can import them into your method groups with the built-in method manager in the *tiamo* software. Please consult the *tiamo* user manual or online-help for further instructions.

The methods can be found in the following path:

C:\Program Files\Metrohm\tiamo\examples\methods\english\855 Robotic Analyzer\Acid Base

(The *tiamo* installation path of your system may be different.)

The *tiamo* user manual (in pdf file format) can be found in the following path:

C:\Program Files\Metrohm\tiamo\doc\english

Accessories of the Robotic Acid/Base Analyzer

1.855.0020	855 Robotic Titrator
1.800.0010	Dosino
1.802.0010	Rod stirrer
6.0253.100	LL-Aquatrode plus OK
6.1236.020	Sleeve with SGJ 14/12 mm 2x
6.1458.040	Titration head insert 3x SGJ 14
6.1459.300	Sample beakers
6.1462.170	Robotic arm with sensor
6.1543.170	Aspiration tip M8
6.1546.030	Piston rod
6.1608.023	Bottle 1 L
6.1621.000	PE container 10 L 2x
6.1805.060	FEP tubing / M6 / 60 cm 3x
6.1805.120	FEP tubing / M6 / 100 cm
6.1805.510	PTFE tubing M8, 60 cm
6.1812.000	PTFE tubing 4/6 mm, 4m
6.1828.000	PVDF connection nipple 2x
6.1909.020	Stirring propeller 104 mm
6.2001.120	Bottle holder base
6.2041.840	Sample rack 59 x 120 mL
6.2053.000	Cable clip 10x
6.2061.010	Bottle holder
6.2104.020	Electrode cable 1 m
6.2151.000	Cable USB A - mini-DIN 8P
6.2308.050	KCl solution 3 mol/L (50 mL)
6.2323.000	Storage solution
6.2621.030	Hexagon key 4 mm
6.2621.070	Hexagon key 5 mm
6.2621.130	Hexagon key 2 mm
6.2621.140	Hexagon key 2.5 mm
6.2740.020	Spray nozzle 3x
6.2751.100	Safety shield
6.3032.220	Dosing Unit 20 mL
6.6056.112	tiamo 1.1 full version CD
T.2400.102	Ferrite cores 4x

Table of contents

	Page
1 Application Information	1
1.1 Introduction: titration means counting!	1
1.2 Acid/Base reactions	2
Protolysis	2
Acid-base reactions with water	3
Protolytes	4
1.3 Acid-base titrations	5
Titration of strong protolytes (strong acid with strong base and vice versa)	5
Titration of a weak acid with a strong base	5
Titration of a weak base with a strong acid	6
Titration of a weak acid with a weak base	6
Titration curve	6
1.4 Calibrating the pH glass electrode	6
Recommended accessory	6
Reagents	6
General	6
Electrode calibration and electrode handling	9
Handling	10
1.5 <i>tiamo</i> method: Calibrating the pH glass electrode	10
Application note	10
Procedure	10
Remark	10
Calibration report (example)	11
2 Titrants	12
2.1 Preparing the most important alkaline and acidic titrants	12
General	12
2.2 Titer determination of alkaline titrants	13
– $c(\text{NaOH}) = 0.1 \text{ mol/L}$	13
Recommended accessories	13
Titer determination	13
Calculation	13
2.3 <i>tiamo</i> method: Titer determination of alkaline titrants	14
Application note	14
Sample preparation	14
Procedure	14
Remarks	14
Result report (example)	15
Statistical results (example)	16
2.4 Titer determination of acidic titrants	17
– $c(\text{HCl}) = 0.1 \text{ mol/L}$ in H_2O	17
– $c(\text{H}_2\text{SO}_4) = 0.05 \text{ mol/L}$	17
Recommended accessories	17
Titer determination	17
Calculation	17
2.5 <i>tiamo</i> method: Titer determination of acidic titrants	18
Application note	18
Sample preparation	18
Procedure	18
Remarks	18
Result report (example)	19
Statistical results (example)	20

3	Typical applications for the Robotic Acid/Base Analyzer	21
3.1	pH value and acid capacity in water	21
	General	21
	Recommended accessories	21
	Reagents	21
	Analysis	21
	Calculations	21
	Remarks	21
3.2	<i>tiamo</i> method: pH value and acid capacity in water	22
	Application note	22
	Sample preparation	22
	Procedure	22
	Remarks	22
	Result report (example)	23
	Statistical results (example)	24
3.3	Determination of the total titratable acidity in wine	25
	Recommended accessories	25
	Reagents	25
	General	25
	General procedure	25
	Detailed procedures	25
	Calculation	26
	Remarks	26
3.4	<i>tiamo</i> method: Determination of the total titratable acidity in wine	27
	Application note	27
	Sample preparation and procedure	27
	Remark	27
	Result report	28
	Statistical results (example)	29
3.5	Determination of different pharmaceutical substances	30
	Recommended accessories	30
	Reagents	30
	Citric acid (USP)	30
	Calculation	30
	N-Acetyltirosine (Pharm.Europe)	30
	Allantoin (Pharm.Europe)	30
	Biotin (USP)	31
	Busulfan (USP)	31
	Cilazapril (Pharm.Europe)	31
	Enalapril maleate (Pharm.Europe)	31
	Glutamic acid (Pharm.Europe)	31
	Histidine hydrochloride monohydrate (Pharm.Europe)	31
	Lisinopril dihydrate (Pharm.Europe)	32
	Mandelic acid (USP)	32
	Methimazole (USP)	32
	Phenylmercuric borate (Pharm.Europe)	32
	Potassium hydrogen tartrate (Pharm.Europe)	32
	Procarbazin hydrochloride (USP)	32
	Ramipril (Pharm.Europe)	32
	Saccharin calcium (USP)	32
	Saccharin sodium (USP)	33
	Tiaprofenic acid (Pharm.Europe)	33
3.6	<i>tiamo</i> method (example): Titration of citric acid	34
	Application note	34
	Sample preparation	34
	Procedure	34

Remarks	34
Result report.....	35
Statistical results (example).....	36
3.7 Determination of phosphoric acid in cola drinks	37
Recommended accessories.....	37
Reagents.....	37
Sample preparation	37
Analysis	37
Calculations	37
Remarks	37
3.8 <i>tiamo</i> method: Determination of phosphoric acid in cola drinks	38
Application note	38
Sample preparation	38
Procedure	38
Remark.....	38
Result report.....	39
Statistical results (example).....	40
4 Troubleshooting pH glass electrodes	41
4.1 Troubleshooting glass membrane.....	41
4.2 Cleaning and care of pH-electrodes with fixed-sleeve diaphragm	42
4.3 Storage	42
5 Method reports	43
5.1 Calibrating the pH glass electrode.....	43
5.2 Titer determination of alkaline titrants	50
5.3 Titer determination of acidic titrants	58
5.4 pH value and acid capacity in water.....	67
5.5 Determination of the total titratable acidity in wine.....	76
5.6 Titration of citric acid.....	83
5.7 Determination of Phosphoric acid in cola drinks	91

1 Application Information

1.1 Introduction: titration means counting!

Together with gravimetry, titration is one of the oldest analytical methods. Both belong to a group of analytical methods that is based on chemical reaction.

In a titration one determines the volume of a standard solution (titrant) that is necessary for complete chemical reaction with the analyte. The titrant contains a known amount of a particular substance.

Since Loschmidt and Avogadro we know that one gram molecule of a substance contains a defined number of particles. A standard solution is produced by dissolving a particular weight of a substance in a solvent. Each volume fraction of this standard solution contains a defined number of particles of the dissolved substance. This means that measuring the volume of a standard solution is a method of counting particles: Titration means counting!

Despite many new, mainly physical instrumental analytical methods, titrimetry as a «wet-chemistry method» still remains a standard procedure for quantitative analysis today. This is because it has a number of specific advantages:

- Titration is one of the absolute content determination methods, i.e. the result of the analysis provides direct information about the amount of substance to be determined, without instrument or method-specific factors having to be calibrated (such as is normal in relative methods, for example HPLC, atomic spectroscopy or UV/VIS photometry).
- Titrations are easy to carry out: The equipment and the procedures to be performed are simple. They are easy to understand – the fundamentals of titrimetric methods are widely known or can be learned in a short time.
- Titrations are carried out rapidly: If the total time for setting up the workplace to obtaining the analytical result is taken into consideration, then titrimetric determinations require much less time than other methods.
- Titration is a versatile method: Numerous titration methods have been drawn up, these range from the determination of inorganic ions up to the determination of complex organic compounds. The analyte concentrations can range from 100% or virtually 100% (analysis of ultra pure substances, purity determinations) down to the ppm range. Sample amounts of a few micrograms are adequate, amounts in the gram range are also possible.
- Titration supplies highly reproducible and correct results. A typical reproducibility is <1%. In high-precision titrations values of 0.1% are demanded and also achieved. For such demands the accuracy should lie within the limits of the standard deviation.
- Titrations can be automated: Titrimetric determinations can be automated to a high degree. This means that, for example, they are suitable for analyzing the content in routine operation (e.g. active ingredient content in tablets).

1.2 Acid/Base reactions


The term acid-base has been revised several times since the 18th century in accordance with the state of knowledge:

- Acids contain oxygen (Lavoisier)
- Acids contain hydrogen that can be replaced by a metal (Liebig)
- Acids dissociate in aqueous solution to release protons and bases release hydroxide ions (Arrhenius, Ostwald)

These theories can only partly explain chemical reactions. They refer to aqueous solutions; ampholytes are not included. Examples are aqueous solutions of salts: KHSO_4 reacts acidic, Na_2HPO_4 alkaline.

Brønsted defined the acid-base reaction as the exchange of protons. The *acid* acts as *proton donor* and the *base* as *proton acceptor*. A deprotonated acid is called the conjugate base and a protonated base is called the conjugate acid. In this way the theory can also be applied without any problem to ampholytes, salts and prototropes of non-aqueous solvents. The equilibrium constants of the reactions allow a quantitative reactivity series of acids and bases to be drawn up (see below).

Lewis extended the Brønsted theory (Brønsted theory). A Lewis base has a free electron pair, a Lewis acid an electron **pair** gap. Lewis acid and Lewis base react to form a covalent or coordinative compound. Examples:




Typical Lewis acids are molecules with an incomplete valency shell, cations as the central atoms of complexes and molecules with multiple polar bonds (e.g. acid anhydrides). Typical Lewis bases are molecules and ions with free electron pairs, anions as complex ligands and molecules with multiple bonds. The Lewis model can also be used to describe precipitation and complex-forming reactions. As no universal reference acid or base exists, there is no practical quantitative reactivity series for Lewis acids and bases; however, the Lewis conception allows numerous chemical reactions to be systemized qualitatively.

Protolysis

Polar (prototropic, protic) solvents undergo self-dissociation, i.e. they form an acid-base pair to a slight extent.

Example water:



Neither free protons nor free electrons are found in a solution. For protons this means that H^+ should not be used as an expression, but rather H_3O^+ to represent the clusters of protons and water effectively present in an aqueous solution.

According to the law of mass action, the ionic product of water K_w is calculated as follows:

$$a(\text{X}^+) \times a(\text{Y}^-) / a^2(\text{XY})$$

for H_2O this means:

$$a(\text{H}_3\text{O}^+) \times a(\text{OH}^-) / a^2(\text{H}_2\text{O}) = K_w = 10^{-14}$$

K is temperature-dependent. In tables it is usually given for 25 °C (the values for H_2O are 14.9 at 0 °C, 14.2 at 20 °C, 14.0 at 25 °C, 13.5 at 40 °C, etc.).

The relative «strength» of an acid or base in an aqueous solution is described by the equilibrium constants. These provide information about the extent to which an acid or base dissolved in water at equilibrium can accept or donate protons. Water is the reference base for each acid and the reference acid for each base. These acidity and basicity constants are

known as K_a and K_b respectively. We have simplified their formulation here by using the molar concentration instead of the activities of the participating species. The molar concentration of water in a mixture is regarded as being constant for diluted solutions (55.5 mol/L) and not included.

$$K_a = [\text{H}_3\text{O}^+] \times [\text{A}^-] / [\text{HA}]$$

$$K_b = [\text{BH}^+] \times [\text{OH}^-] / [\text{B}]$$

The negative common logarithms of these values are, analogous to the pH definition, known as $\text{p}K_a$ and $\text{p}K_b$ respectively. The higher the pK value, the weaker the acid or base.

The neutral point of a non-aqueous solvent is calculated in a similar manner to that for water. Clearly, the pH term only applies to purely aqueous solutions!

$$\text{H}_2\text{O, neutral pH} = -\log \sqrt{K_w} = -\log \sqrt{10^{-14.0}} = +7.0$$


The following table shows a small selection of autoprotolysis constants of prototropic solvents (at 25 °C).

Solvent	Autoprotolysis	$\text{p}K_p$
Ethanol	$2 \text{C}_2\text{H}_5\text{OH} \rightleftharpoons \text{C}_2\text{H}_5\text{OH}_2^+ + \text{C}_2\text{H}_5\text{O}^-$	19.1
Methanol	$2 \text{CH}_3\text{OH} \rightleftharpoons \text{CH}_3\text{OH}_2^+ + \text{CH}_3\text{O}^-$	16.7
Acetic acid	$2 \text{CH}_3\text{COOH} \rightleftharpoons \text{CH}_3\text{COOH}_2^+ + \text{CH}_3\text{COO}^-$	14.5
Water	$2 \text{H}_2\text{O} \rightleftharpoons \text{H}_3\text{O}^+ + \text{OH}^-$	14.0
Formic acid	$2 \text{HCOOH} \rightleftharpoons \text{HCOOH}_2^+ + \text{HCOO}^-$	6.2
Sulfuric acid	$2 \text{H}_2\text{SO}_4 \rightleftharpoons \text{H}_3\text{SO}_4^+ + \text{HSO}_4^-$	3.6


Acid-base reactions with water

In aqueous solutions the protolysis reaction with water normally predominates (as an ampholyte H_2O can act as a Brønsted acid or a Brønsted base):

acidic reaction of solvent



basic reaction of solvent



Just as with the ionic product of water, $a(\text{solvent})$ is set = 1 for dilute solutions. The resulting equilibrium constants are known as the *acidity constant* K_a and the *basicity constant* K_b .

$$K_a \times K_b = K_w$$

$$\text{p}K_a + \text{p}K_b = \text{p}K_w = 14 \text{ (25 °C)}$$

The stronger the acid or base, the weaker the conjugate base or acid.

At the same time the pK values of the aqueous system limit the acidic or basic strength of protolytes in aqueous solution (differentiation).

Stronger acids and bases ($\text{p}K < 0$ or $\text{p}K > 14$ for the conjugate protolyte) are converted to oxonium ions (H_3O^+) and hydroxide ions (OH^-) in aqueous solution. *Leveling* occurs – the separate determination of mixed very strong acids or mixed very strong bases in an aqueous solution is not possible by titration. (Only one endpoint is found in the titration of mixtures of e.g. KOH/NaOH or $\text{HCl/H}_2\text{SO}_4$).


The known corresponding acid-base pairs are arranged according to increasing pK_a (or decreasing pK_b value) and normally divided into the following categories:

- very strong $pK_a < 0$
- strong $pK_a 0 \dots 4$
- weak $pK_a 4 \dots 10$
- very weak $pK_a 10 \dots 14$
- extremely weak $pK_a > 14$

Protolytes

As a result of the leveling effect of water, differences in the protolysis constants of strong protolytes can only be determined in non-aqueous solvents. In such solvents the relationship $pK_a + pK_b = pK_{\text{solvent}}$ also applies for the conjugate acid-base pair. This means that the solvent cation (*lyonium ion*) and the solvent anion (*lyate ion*) are always the strongest acids and bases in the affected system.

Whereas acetic acid behaves in water like a weak acid, in HClO_4 it behaves like a base:




In order to be able to achieve the separate determination of the components in an aqueous solution by titration, the differences in their pK values should be approx. 5 units. In contrast, in suitable non-aqueous solvents a difference of 2 to 3 pK units is normally sufficient.

pKa values of some selected acids		pKb values of some selected bases	
Acids	pKa	Bases	pKb
Acetic acid	4.73	2-Picoline	7.52
Acrylic acid	4.26	Acridine	9.89
Benzoic acid	4.20	Ammonia	4.75
Boric acid	9.24	Aniline	9.42
Chloroacetic acid	2.81	Benzylamine	4.62
Hydrobromic acid	≈ -6	Calcium hydroxide	1.30
Hydrochloric acid	≈ -3	Coffeine	13.39
Hydrocyanic acid	9.40	Cyclohexylamine	3.36
Hydrofluoric acid	3.14	Diethanolamine	5.12
Hydrogen sulfide acid 1 st stage	6.90	Ethanolamine	4.56
Hydrogen sulfide acid 2 nd stage	12.9	Ethylamine	3.33
Lactic acid	3.86	Imidazole	7.00
Nitric acid	1.32	Lithium hydroxide	-0.10
Oxalic acid 1st stage	4.31	Magnesium hydroxide	3.36
Oxalic acid 2nd stage	1.42	Naphthylamine	10.08
Perchloric acid	≈ -9	o-aminobenzoic acid	10.92
Phenol	9.95	o-toluidine	9.61
Phosphoric acid 1st stage	1.96	Piperidine	2.80
Phosphoric acid 2nd stage	7.12	Pyridine	8.81
Phosphoric acid 3rd stage	12.36	Triethanolamine	6.23
Salicylic acid	2.98	Triethylamine	3.28
Sulfuric acid 1st stage	≈ -3	Trimethylamine	4.20
Sulfuric acid 2nd stage	1.92	Urea	13.8


The strongest acid in the list is perchloric acid, the strongest base is lithium hydroxide.

The following empirical rules apply for the *relative acid strength* of inorganic acids:

- *Hydracids*: the acidity increases as the atomic number increases, both within the group and within the period. (This means that HI is the strongest acid.):



- *Oxyacids* are stronger the fewer H atoms and more O atoms they contain. (The strongest mono-oxyacid is HClO_4).



1.3 Acid-base titrations

Acid-base reactions occur very rapidly (proton transfer takes far less than one millionth of a second).

In mixtures of acids the titration is always carried out in the sequence of the relative acid strengths. This means that the strongest acid is always titrated first, the weakest acid last.

Example:

The titration curve of a mixture of HCl and CH_3COOH shows two equivalence points/endpoints. The consumption up to EP1 corresponds to the molar concentration of HCl and the consumption for the difference EP2 – EP1 to the molar concentration of acetic acid.

Titration of strong protolytes (strong acid with strong base and vice versa)

The salts of strong acids with strong bases do not undergo protolysis with water, which means that an equivalent amount of these acids and bases react neutrally. The equivalence point in the titration of a strong acid with a strong base and vice versa is therefore identical with the neutral point ($\text{pH} = 7.0$ in H_2O).

The shape of the titration curve results from the fact that in order to increase the pH from 1 to 2 it is necessary to reduce the H_3O^+ concentration to 1/10 of the original value (or that 90% of the equivalent base has already been added). From $\text{pH} = 2$ to $\text{pH} = 3$ it is then 9%, from $\text{pH} = 3$ to $\text{pH} = 4$ then 0.9%, etc. With the excess base it is exactly the opposite; from $\text{pH} = 7$ to $\text{pH} = 8$ 0.0009% is necessary, from $\text{pH} = 8$ to $\text{pH} = 9$ requires 0.009%, from $\text{pH} = 9$ to $\text{pH} = 10$ it is 0.09% excess base, etc.

This is why typical symmetric titration curves with a steep «jump region» are obtained, whose equivalence point (titration endpoint, corresponds to the point of inflection of the curve) lies at $\text{pH} = 7.0$.

Titration of a weak acid with a strong base

Before the start of the titration it can be seen that the pH is acidic, which corresponds to the weak acid not being completely dissociated in water. After the start of the titration a buffering range (10...90% of neutralization) can be observed in which the weak acid and its salt or acid anion are both present. In this buffering range the addition of the base only alters the pH slightly. The point of inflection of the buffering range corresponds to the semi-equivalence point of the titration (half reaction) and therefore, according to the Henderson-Hasselbalch equation, the pK_a value of the weak acid is given by: $\text{pH} = \text{pK}_a + \log \frac{[\text{A}^-]}{[\text{HA}]}$, where at half reaction the molar concentrations of acid anion and acid are equal – resulting in $\log 1$, which is equal to 0.

In the above case the equivalence point no longer lies at $\text{pH} = 7$, because the added base competes with the conjugate base of the acid. At the equivalence point the *salt solution* present has an *alkaline* reaction, as the acid anions of a weak acid undergo protolysis with water (hydrolysis). The equivalence point of such a titration lies – depending on the relative strength of the acid – in the pH range 7.5...10.

The higher the concentration of acid and base in the solution, the wider the buffering range on the volume axis. The pH value achieved by adding additional base beyond the equivalence point finally approaches that of the base used as the titrant.

Titration of a weak base with a strong acid

Before the start of the titration an alkaline pH is measured that corresponds to that of the not fully dissociated weak base in water. After the start of the titration a buffering range can also be observed (10...90% of the neutralization), in which the weak base and its conjugated protonated form are both present. Addition of acid hardly alters the pH in this buffering range at all. The point of inflection of the buffering range corresponds to the semi-equivalence point of the titration (half reaction) and lies at that pH value which corresponds to the pK_b value of the base. The equivalence point is again no longer at $pH = 7$, because the added acid competes with the conjugate acid of the base. The corresponding salt solution is acidic as the protonated base undergoes protolysis with water (hydrolysis). The equivalence point lies – depending on the relative strength of the base – in the pH range 6.5...4.

Titration of a weak acid with a weak base

This titration is unsuitable for quantitative determinations! There is no marked jump in the titration curve, and the slope of the curve in the equivalence region does not achieve a maximum value. No suitable titration endpoint for a quantitative determination is reached. This means that, when selecting a titrant for the titration of weak acids and bases, weak bases and acids should not be chosen.

Titration curve

The following applies in general: The «jump» of a titration curve is more marked the higher the K_a or K_b values of the two reaction partners and the higher the concentrations of both species.

1.4 Calibrating the pH glass electrode

Recommended accessory


- Comb. pH glass electrode for aqueous titrations 6.0257.000 Aquatrode Plus with 6.2104.020 electrode cable

Reagents

- Buffer solutions $pH = 4.00$, 7.00 and 9.00 , e.g. 6.2307.100, 6.2307.110 and 6.2307.120

General

Free protons (H^+ ions) occur in solutions just as little as free electrons. They combine with water to form oxonium ions:



The pH value is defined as the negative logarithm of the oxonium ion activity, i.e. of the concentration of free, dissociated oxonium ions in mol/L: $pH = -\log [H_3O^+]$

Strictly speaking, the term pH only applies to purely aqueous solutions.

The pH scale ranges from 0 to 14 with the neutral point at $pH = 7.0$, where the H_3O^+ and OH^- ions are present in equilibrium. pH values below 7 result from an H_3O^+ excess, pH values above 7 from an OH^- ion excess. The more acidic a solution the higher its H_3O^+ ion concentration and the lower its pH value.

Weak acids, e.g. tartaric acid, do not dissociate completely. This means that only a small fraction (approx. 2...3%) of their acid ions are released. This also means that only on very rare occasions can the pH value be used as a measure of the concentration of acids or bases.


As the pH scale is logarithmic this means that small differences in pH correspond to large differences in the concentration of H_3O^+ ions. For example, at pH = 3.0 there are ten times more H_3O^+ ions present than at pH = 4.0, and at pH = 3.1 there are twice as many H_3O^+ ions present than at pH = 3.4.

The measuring setup for potentiometric measurements always consists of two electrodes – a *measuring or indicator electrode* and a *reference electrode*. For practical reasons these two electrodes are usually contained in a single *combined* electrode.

The indicator electrode (in this case the pH glass electrode) produces a potential that is dependent on the composition of the sample solution.

The reference electrode (usually Ag/AgCl) has the task of providing a potential that is as independent as possible of the sample solution (reference potential).

The potential measurement itself takes place virtually current-free by using a «voltmeter» with a high-impedance measuring input (this is necessary to avoid unwanted potential drops). The measured potential U is made up from the individual potentials produced by the indicator and reference electrodes. The following illustration shows a schematic diagram with a separate pH glass electrode (left) and a reference electrode (right):



U_1 : Galvani potential between measuring electrode and measuring solution

U_2 : Galvani potential between internal buffer and glass membrane

U_3 : Galvani potential between internal reference electrode and internal buffer

U_4 : Galvani potential of reference electrode

U_5 : Diffusion potential at the diaphragm

The individual potentials U_2 , U_3 and U_4 are determined by the construction of the electrodes and are therefore constant for a given electrode pair. The diffusion potential U_5 should be kept relatively constant and low by taking suitable measures. These measures include an optimal and clean diaphragm, constant stirrer speed during the measurements as well as a suitable reference electrolyte solution whose anions and cations have similar ionic mobilities – e.g. KCl. In this way the potential U_1 measured between the electrodes depends only on the sample solution. In the pH measurement this potential is again dependent on the activity a_i of the measuring ion (H_3O^+ ion / OH^- ion). This relationship is described by the *Nernst equation*:

$$U = U_0 + \frac{2.303 \times R \times T}{z_i \times F} \times \log a_i = U_0 + U_N \times \log a_i$$

where:

U measured difference in potential between indicator and reference electrode

U_0 standard potential of the combined electrode (depends on its construction)

R gas constant (8.31441 J / (K mol))

T absolute temperature in K (273.15 + t / °C)

z_i charge on the measuring ion i including its sign (+1 for H_3O^+ and –1 for OH^-)

a_i activity of measuring ion

- U_N Nernst slope (59.16 mV at 25 °C and $z = 1$)
- 2.303 conversion factor from natural to common logarithm

The Nernst slope U_N describes the theoretical electrode slope and corresponds to the change in potential produced by altering a_i by a factor of ten. It depends on the temperature and charge z of the measuring ion. **Please note: The instrument compensates the effect of temperature on U_N but not on the pH value of the solution!**

The following table shows values of U_N as a function of $t / ^\circ\text{C}$ for $z = 1$:

Temperature $t / ^\circ\text{C}$	U / mV	Temperature $t / ^\circ\text{C}$	U / mV
0	54.20	40	62.14
10	56.18	50	64.12
20	58.17	60	66.10
25	59.16	70	68.09
30	60.15	80	70.07
38	61.74	90	72.06

An ideal pH glass electrode has a slope of 1 (100% of the Nernst slope) and an electrode zero point pH_{as} of 7.0, the latter corresponding to $U_{\text{as}} = 0$ mV.

Things are different in practice. The electrode zero point should have a value for U_{as} of ± 15 mV (corresponds to $\text{pH}_{\text{as}} = 6.75 \dots 7.25$) and the slope should be > 0.95 (> 56.2 mV / pH at 25 °C).

In order to «inform» the instrument of the true electrode data it is necessary to calibrate the electrode. Buffer solutions have a defined pH value which, however, is temperature-dependent. The relevant information about the buffer solution is entered during the calibration procedure. The following table shows the pH values of Metrohm buffer solutions as a function of the temperature:

Temperature $t / ^\circ\text{C}$	pH = 4.00 ± 0.02	pH = 7.00 ± 0.02	pH = 9.00 ± 0.02
10	3.99	7.06	9.13
20	3.99	7.02	9.04
25	4.00	7.00	9.00
30	4.00	6.99	8.96
38	4.02	6.98	8.91
40	4.02	6.98	8.90
50	4.04	6.97	8.84
60	4.07	6.97	8.79
70	4.11	6.98	8.74
80	4.15	7.00	8.71
90	4.20	7.01	8.68


Buffer solutions are not stable!! They can be decomposed by bacteria and/or molds or – this applies to alkaline buffer solutions – alter their pH value by absorbing CO_2 from the atmosphere. This is why you should always use only fresh buffer solutions and reject them after use, i.e. not pour them back into the storage bottle.

Electrode calibration and electrode handling

We recommend the following procedure for **calibrating** a pH glass electrode:

- Remove the electrode from its storage vessel, attach a cable if necessary and connect it to the instrument.
- Open the electrolyte filling opening and, if necessary, top up the electrolyte solution.
- Rinse the electrode thoroughly with dist. H_2O and dab dry with a soft paper tissue (do not rub).
- Fill pH = 7.0 buffer solution into a beaker and add a stirrer bar.
- Immerse the electrode into the buffer solution and stir for approx. 1 min.
- On the **tiamo** enter the pH value of the buffer solution (at the corresponding temperature) and start the calibration with buffer 1 under stirring.
- When the measured value has been accepted, remove the electrode from the solution, rinse it thoroughly with dist. H_2O and dab dry with a soft paper tissue.
- Add pH = 4.0 or 9.0 buffer solution to a second beaker, add a stirrer bar, immerse the electrode and stir for approx. 1 min (the second buffer solution must have the same temperature as the first one).
- On the **tiamo** enter the pH value of the second buffer solution (at the corresponding temperature) and continue the calibration under stirring.
- After the measured value has been accepted, end the calibration. Remove the electrode from the solution, rinse the electrode thoroughly with dist. H_2O and dab dry with a soft paper tissue.

What happens in the instrument during calibration can be seen in the following plot:



Handling

Laboratory electrodes should have a long lifetime. Their characteristics (slope, response behavior, $\text{pH}_{\text{as}} / U_{\text{as}}$) must lie within the given criteria. In order to ensure this a few basic rules must be observed:

- After use rinse the electrode thoroughly with dist. H_2O and dab dry with a soft paper tissue. Close off the electrolyte filling opening and store the electrode by immersing it in storage solution to an adequate depth. Dry storage leads to delayed and poor response behavior. The electrolyte solution may become concentrated and $\text{pH}_{\text{as}} / U_{\text{as}}$ could alter. Storage in dist. H_2O could result in the diaphragm being blocked by AgCl .
- If the electrode responds sluggishly and/or the slope is unsatisfactory then the electrode membrane must be etched. This is done by immersing the membrane in a 10% solution of ammonium difluoride (NH_4HF_2 , use plastic beaker) for 1 min, then swirling it for approx. 10 s in $c(\text{HCl}) = 5 \text{ mol/L}$, rinsing it thoroughly with dist. H_2O and then wiping off the silicate residue with a moist tissue. In order to build up a new gel layer the electrode is placed in $c(\text{KCl}) = 3 \text{ mol/L}$ for 24 h (or for 5 h in the same solution at 50°C).
- If the diaphragm becomes blocked please refer to the electrode data sheet that accompanies each electrode. Removing such a blockage is complicated and time-consuming – it is better to send the electrode to your local Metrohm distributor for this.
- Contamination by fats, oils, lacquers, paints, etc.: remove the contamination with an organic solvent (acetone, petroleum benzene, toluene), rinse thoroughly with ethanol and dist. H_2O , dab dry and place in electrolyte solution.
- Contamination by proteins: immerse the electrode in a solution of 5% pepsin in $c(\text{HCl}) = 0.1 \text{ mol/L}$ for a few hours. Then rinse thoroughly with dist. H_2O , dab dry and place in electrolyte solution.

1.5 *tiamo* method: Calibrating the pH glass electrode

Application note

This is a method to calibrate a pH glass electrode. It consists of three different tracks.

In the series start track the rack is initialized and the system prepared for the measurement. Additionally the electrode is removed from the storage beaker and rinsed thoroughly.

The main track deals the calibration measurements followed by the rinsing of the electrode. After the measurements the results are calculated, a report is printed and the data saved in the previously defined database 'Robotic Acid Base Analyzer' (can be modified).

Procedure

Fill three beakers with Metrohm buffer solutions pH 4, 7 and 9. Make sure that the electrode is properly immersed in the buffer solution. The special beaker for storing the electrode has to be filled with Metrohm storage solution.

Remark

To run this method the settings of the 855 Robotic Titrator have to be adjusted.

The lift positions have to be defined according to the rack and beakers used. As the robotic swing arm is equipped with a sensor the head has to touch the beaker brim when moving to the work position. If this requirement can not be fulfilled no beaker will be recognized. Two special beakers have to be defined including their own work positions. Special beaker 2 is used for rinsing the electrode while special beaker 1 contains storage solution for a proper electrode treatment between the determination series.

Calibration report (example)



Robotic Acid / Base Analyzer
Calibration

Programm version tiamo 1.1

2005-09-13 11:12:24 UTC+2

Results report

Determination

Method Calibrating the pH glass electrode
 Method saving date 2005-08-19 11:02:38 UTC+2
 Method version 1
 Method state original
 Determination ID -4f2e45c7:105cd5e9155:-7e1c
 Determination start 2005-08-19 11:05:00 UTC+2
 Determination state original
 Determination version 1
 Run number 3
 User (full name) Metrohm

Sample data


Sample Identification 1 Calibration

Results

slope **99.30 %**
zero point **6.89**

Calibration curve

calibration loop pH.1 - CAL LOOP pH



2 Titrants

2.1 Preparing the most important alkaline and acidic titrants

General

Titrants are standard solutions, i.e. solutions that contain a defined content of a reactant. This content is given as the molar concentration c in mol/L. The «normality», which was frequently used previously, is no longer valid today and should therefore not be used.

Examples:

- 0.1 N HCl $\Rightarrow c(\text{HCl}) = 0.1 \text{ mol/L}$
- 0.1 N H_2SO_4 $\Rightarrow c(\text{H}_2\text{SO}_4) = 0.05 \text{ mol/L}$

Not all titrants have a stable titer, i.e. their concentration can vary with time.

Examples:

- Hydroxides absorb CO_2 from the atmosphere to form carbonates

This all means that the titer of the titrant may alter as time passes. In order to know the true titer concentration the titer must be determined at regular intervals.

The so-called standard titrimetric substances are used for determining the titer. Their content hardly changes, they are available with a defined degree of purity, can be dried and can be traced back directly to standard reference materials (e.g. National Institute of Standards and Technology – NIST, USA).

Such standard titrimetric substances /secondary standards are:

- For bases potassium hydrogen phthalate, $M = 204.23 \text{ g/mol}$
- For acids tris(hydroxymethyl)-aminomethane, $M = 121.14 \text{ g/mol}$

Most standard solutions/titrants are commercially available as ready-to-use solutions with a titer adjusted by the manufacturer at 20 °C to 1.000. We recommend that you purchase such ready-to-use solutions and do not prepare them yourselves.

In principle titer determinations should always be carried out at the same temperature at which the analyses are later to be carried out. Please note that solutions expand as their temperature increases.


For aqueous solutions a temperature difference of 5 °C for a theoretical consumption of 10.00 mL results in a difference in volume of 12.5 μL . This means that a titer of 1.0000 at 20 °C becomes a titer of 0.9988 at 25 °C – and for non-aqueous solutions this difference is even larger. In such a case a titer of 1.0000 at 20 °C becomes a titer of 0.9950 at 25 °C (possible error 0.5%)!

Normally the titer determination is carried out three times and the mean value is used. The mean value of the titer** is saved as «TITER» in the titrant list of tiamo.

** In Europe a dimensionless factor, with at least 4 decimal places (e.g. 1.0015).
In the USA already multiplied by c , with at least 4 decimal places (e.g. 1.0015 mol/L).

2.2 Titer determination of alkaline titrants

Alkaline titrants do not have a stable titer. They may absorb CO₂ from the atmosphere to form carbonates, e.g.:



This not only reduces the titer. The different strengths of the bases NaOH and Na₂CO₃ may negatively influence the titration curves and therefore the results – high-bias results are simulated, particularly for weak acids. In order to reduce CO₂ absorption as much as possible, soda lime (e.g. Merck no. 106839) is placed in the drying/absorber tube of the Exchange Unit.

– c(NaOH) = 0.1 mol/L

If possible this titrant should be bought ready for use. Otherwise dissolve 4.0 g NaOH in CO₂-free dist. H₂O, make up to 1 liter and mix.

Recommended accessories

- Comb. pH glass electrode for aqueous titrations 6.0257.000 Aquatrode Plus with 6.2104.020 electrode cable
- 20 mL Dosing Unit: 6.3032.220

Titer determination

Potassium hydrogen phthalate is dried overnight in a drying oven at 105 °C and allowed to cool down in a desiccator for at least 1 h.

Approx. 200 mg potassium hydrogen phthalate is weighed out into the titration beaker with an accuracy of 0.1 mg and dissolved in approx. 50 mL dist. H₂O. It is then immediately titrated to after the first endpoint with c(base) = 0.1 mol/L.

Calculation

$$\text{Titer} = \text{'MV.Sample size'} * 1000 / (\text{'DET pH.EP\{1\}.VOL'} * \text{'DET pH.CONC'} * 204.23)$$

where:

'MV.Sample size'	weight of potassium hydrogen phthalate [g]
1000	conversion factor for milliliter
'DET pH.EP{1}.VOL'	mL base up to endpoint
'DET pH.CONC'	concentration of the NaOH [mol/L]
204.23	molecular mass of potassium hydrogen phthalate [g/mol]

2.3 *tiamo* method: Titer determination of alkaline titrants

Application note

This is a method to determine the titer of an alkaline titrant. Alkaline titrants tend to absorb CO_2 from the atmosphere which results in a reduced titer value. With the titration of a primary standard the actual titer is determined. The method consists of five different tracks. In the start series track the rack is initialized and the system prepared for the measurements. The main track handles the determination followed by the rinsing of the electrode, which is done after every determination. In the exit track the results are calculated, a report is printed and the data is saved in the predefined database Robotic Acid Base Analyzer (can be modified). In the series end track the electrode is moved to a storage beaker where it is stored in Metrohm storage solution for a proper electrode treatment between the determination series. In case of an error, the error track guarantees that the electrode is moved to the storage beaker.

Sample preparation

The potassium hydrogen phthalate is dried over night in a drying oven at 105°C and allowed to cool down in a desiccator for at least 1 hour. Approximately 200 mg potassium hydrogen phthalate are weighed out in a beaker with an accuracy of 0.1 mg.

Procedure

The beakers are placed on the rack and ca. 60 mL distilled water are added automatically for dissolving. The solution is titrated with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ and the titration stopped after the first endpoint. Then the solution is removed and the electrode is rinsed thoroughly before the next determination.

Remarks

To run this method the settings of the 855 Robotic Titrosampler have to be adjusted. The lift positions have to be defined according to the rack and beakers used. As the robotic swing arm is equipped with a beaker sensor the head has to touch the beaker brim when moving to the work position. If this requirement can not be fulfilled no beaker will be recognized. One special beaker has to be defined including its work position. Special beaker 1 is used to ensure a proper electrode treatment between the determination series.

The value '204.23' used in the calculation formula is the molecular weight of potassium hydrogen phthalate in [g/mol].

Result report (example)



Robotic Acid / Base Analyzer
Titer determination

Programm version tiamo 1.1

2005-09-13 10:49:57 UTC+2

Results report

Determination

Method Titer determination of alkaline titrants
 Method saving date 2005-08-19 13:54:25 UTC+2
 Method version 1
 Method state original
 Determination ID -4f2e45c7:105cd5e9155:-7d4a
 Determination start 2005-08-19 14:06:00 UTC+2
 Determination state original
 Determination version 1
 Run number 7
 User (full name) Metrohm

Sample data


Sample Identification 1 Titer NaOH
 Sample size 154.8 mg

End points

DET pH DET pH.1
 EP1 8.302 pH 7.7528 mL

Results

Titer 0.9777
 mean value titer 0.9776



Statistical data (short)



Method titer determination of alkaline titrants

Number of single determinations 3

Result name	n	Mean value	s +/-	s rel
Titer	3	0.9776	0.00410	0.42 %

Statistical results (example)

Titer determination NaOH (0.1M)

Sample size [g]	Titer
0.1797	0.9735
0.2079	0.9817
0.1548	0.9777
0.2078	0.9821
0.1798	0.9728
0.1732	0.9743
0.2113	0.9772
0.1815	0.9773
0.1604	0.9787
0.1739	0.9781
Mean value	0.9773
abs. std. dev.	0.0031
rel. std. dev. %	0.32

2.4 Titer determination of acidic titrants

– $c(\text{HCl}) = 0.1 \text{ mol/L}$ in H_2O

If possible this titrant should be bought ready for use. Otherwise place approx. 800 mL dist. H_2O in a 1000 mL volumetric flask, add 9.8 mL $w(\text{HCl}) = 32\%$, make up to the mark with dist. H_2O or ethanol and mix.

– $c(\text{H}_2\text{SO}_4) = 0.05 \text{ mol/L}$

If possible this titrant should be bought ready for use. Otherwise place approx. 800 mL dist. H_2O in a 1000 mL volumetric flask and carefully add 2.7 mL $w(\text{H}_2\text{SO}_4) = 96\%$ while swirling the flask. Make up to the mark with dist. H_2O and mix.

Recommended accessories

- Comb. pH glass electrode for aqueous titrations 6.0257.000 Aquatrode Plus
- 20 mL Dosing Unit: 6.3032.220

Titer determination

Tris(hydroxymethyl)-aminomethane (TRIS) is dried overnight in a drying oven at 105 °C and allowed to cool down in a desiccator for at least 1 h. Approx. 100 mg TRIS is weighed out into the titration beaker with an accuracy of 0.1 mg and dissolved in approx. 50 mL dist. H_2O . It is then immediately titrated with HCl or H_2SO_4 to after the first endpoint.

Calculation

$$\text{Titer} = \text{'MV.Sample size'} / (\text{'DET pH.EP\{1\}.VOL'} * \text{'DET pH.CONC'} * 121.14)$$

where:

'MV.Sample size'	weight of TRIS [mg]
'DET pH.EP{1}.VOL'	mL acid up to endpoint
'DET pH.CONC'	concentration of the used acid [mol/L]
121.14	molecular weight of the TRIS [g/mol]

2.5 *tiamo* method: Titer determination of acidic titrants

Application note

This is a method to determine the titer of an acidic titrant. Once the titrant bottle is open the titer value changes. With the titration of a primary standard the actual titer is determined. The method consists of five different tracks. In the start series track the rack is initialized and the system prepared for the measurements. The main track handles the titration followed by the rinsing of the electrode, which is done after every determination. In the exit track the results are calculated, a report is printed and the data is saved in the predefined database Robotic 'Acid Base Analyzer' (can be modified). In the series end track, the electrode is moved to a storage beaker where it is stored in Metrohm storage solution for a proper electrode treatment between the determination series. In case of an error, the error track is carried out, which guarantees that the electrode is moved to the storage beaker.

Sample preparation

The tris(hydroxymethyl)-aminoethane (TRIS) is dried over night in a drying oven at 105°C and allowed to cool down in a desiccator for at least 1 hour. Approximately 100 mg TRIS are weighed out in a beaker with an accuracy of 0.1 mg.

Procedure

The beakers are placed on the rack and ca. 60 mL distilled water are added automatically for dissolving. The solution is then titrated with $c(\text{HCl}) = 0.1 \text{ mol/L}$ and the titration stopped after the first endpoint. Then the solution is removed and the electrode is rinsed thoroughly before the next determination.

Remarks

To run this method the settings of the 855 Robotic Titrosampler have to be adjusted. The lift positions have to be defined according to the rack and beakers used. As the robotic swing arm is equipped with a beaker sensor the head has to touch the beaker brim when moving to the work position. If this requirement can not be fulfilled no beaker will be recognized. One special beaker has to be defined including its work position. Special beaker 1 is used for a proper electrode treatment between the determination series.

The value '121.14' in the calculation formula is the molecular weight of TRIS in [g/mol].

Result report (example)



Robotic Acid / Base Analyzer
Titer determination

Programm version tiamo 1.1

2005-09-13 11:02:29 UTC+2

Results report

Determination

Method Titer determination of acidic titrants
 Method saving date 2005-08-24 10:35:22 UTC+2
 Method version 1
 Method state original
 Determination ID 6de54dcc:105e79b0d60:-7f6b
 Determination start 2005-08-24 10:52:20 UTC+2
 Determination state original
 Determination version 1
 Run number 6
 User (full name) Metrohm

Sample data


Sample Identification 1 Titer HCl
 Sample size 0.0872 g

End points

DET pH DET pH.1
 EP1 4.745 pH 7.2152 mL

Results

Titer 0.9977
 mean value titer 0.9986



Statistical data (short)



Method Titer determination of acidic titrants

Number of single determinations 10

Result name	n	Mean value	s +/-	s rel
Titer	4	0.9986	0.00194	0.19 %

Statistical results (example)

Titer determination HCl (0.1M)

Sample size [g]	Titer
0.1040	0.9970
0.1035	0.9830
0.1200	1.0140
0.0872	0.9977
0.0769	1.0021
0.0939	1.0041
0.0911	0.9988
0.0841	1.0020
0.0997	1.0002
0.1019	1.0050
Mean value	1.0007
abs. std. dev.	0.0027
rel. std. dev. %	0.27

3 Typical applications for the Robotic Acid/Base Analyzer

3.1 pH value and acid capacity in water

General

The acid capacity of a water (K_A) is the molar amount of oxonium ions that the water can accept before reaching certain predefined pH values. It is given in mmol/L.

For the acid capacity up to pH 8.2 ($K_{A_{8.2}}$, previously p value) HCl is used to titrate to pH = 8.2; for the acid capacity up to pH 4.3 ($K_{A_{4.3}}$, previously m value) the titration is to pH = 4.3.

Recommended accessories

- Comb. pH glass electrode for aqueous titrations 6.0257.000 Aquatrode Plus with 6.2104.020 electrode cable
- 20 mL Dosing Unit: 6.3032.220

Reagents

- Titrant: $c(\text{HCl}) = 0.1 \text{ mol/L}$. For very soft water possibly use $c(\text{HCl}) = 0.02 \text{ mol/L}$.

Analysis

The Aquatrode has already been calibrated.

50 mL water sample is added to the titration beakers which are placed on their positions on the sample rack. The pH value is first measured under stirring. The samples are then titrated with $c(\text{HCl}) = 0.1 \text{ mol/L}$ using pH = 3.5 as the stop criterion.

Calculations

(Results given to two decimal places)

$$K_{A_{8.2}} = \text{'SET pH.EP\{1\}.VOL'} * \text{'SET pH.CONC'} * \text{'SET pH.TITER'} * 1000 / \text{'MV.Sample size'}$$

$$K_{A_{4.3}} = \text{'SET pH.EP\{2\}.VOL'} * \text{'SET pH.CONC'} * \text{'SET pH.TITER'} * 1000 / \text{'MV.Sample size'}$$

where

'SET pH.EP{1}.VOL' mL HCl until first endpoint is reached or a fixed EP at pH = 8.2

'SET pH.EP{2}.VOL' mL HCl until second EP is reached or a fixed EP at pH = 4.3

'SET pH.CONC' concentration of HCl solution [mol/L] (0.1 or 0.02)

'SET pH.TITER' titer of HCl solution

1000 1000 (conversion factor for mmol/L)

'MV.Sample size' volume of the sample [mL]

Remarks

Many waters have an initial pH value below 8.2 – $K_{A_{8.2}}$ is not found and therefore cannot be calculated (only EP1 at pH = 4.3).

- Because the pH value depends on the temperature of the sample it is recommended to measure the temperature (temperature sensor 6.1110.100) during determinations

3.2 *tiamo* method: pH value and acid capacity in water

Application note

This is a method to measure the pH value and to determine the acid capacity of water samples. The acid capacity of water (KA) is the molar amount of oxonium ions that water can accept before reaching certain predefined pH values. There are existing two KA values, one at pH 8.2 and one at a pH of 4.3. Many waters have an initial pH value below 8.2. Therefore it is impossible to determine the KA at a pH of 8.2. The method consists of six different tracks. In the start series track the rack is initialized and the system is prepared for the determination runs. The main track handles the measurement of the pH and the titration of the sample. The electrode is rinsed after every determination according to the procedure in the rinsing track. In the exit track the results are calculated, a report is printed and the data is saved in the previously defined database 'Robotic Acid Base Analyzer' (can be modified). When the determination series is finished the series end track is executed in which the electrode is moved to a storage beaker where it is stored in a Metrohm storage solution for a proper treatment between the determination series. In case of an error, the error track is carried out, which guarantees that the electrode is moved to the storing beaker.

Sample preparation

50 mL water sample are filled in a titration beaker.

Procedure

The beakers are placed on the rack and the determination series started. The solution is titrated with $c(\text{HCl}) = 0.1 \text{ mol/L}$ and the titration stopped after the endpoint at pH 4.3 is reached. Then the beaker is emptied and the electrode is thoroughly rinsed before the next determination.

Remarks

To run this method the settings of the 855 Robotic Titrosampler have to be adjusted. The lift positions have to be defined according to the rack and beakers used. As the robotic swing arm is equipped with a beaker sensor the head has to touch the beaker brim when moving to work position. If this requirement can not be fulfilled no beaker will be recognized. One special beaker has to be defined including its work position. Special beaker 1 is used for a proper electrode treatment between the determination series.

Result report (example)Robotic Acid / Base Analyzer
pH value and acid capacity in water

Programm version tiamo 1.1

2005-09-13 11:16:05 UTC+2

Results report**Determination**

Method pH value and acid capacity in water
Method saving date 2005-08-25 13:56:24 UTC+2
Method version 1
Method state original
Determination ID 5ba76f4:105ec3853fe:-7d0b
Determination start 2005-08-25 14:22:23 UTC+2
Determination state original
Determination version 1
Run number 9
User (full name) Metrohm

Sample data

Sample Identification 1 Tab water
Sample size 50 mL

End points**MEAS pH MEAS pH.1**

EME 7.979 pH

SET pH SET pH.1

EP1 0.0000 mL 6.31 s
EP2 2.8180 mL 158.41 s

Results

pH 7.98
Temperature 25.33 °C
K (s 8,2) 0.00 mmol/L
K (s 4,3) 11.18 mmol/L

Statistical results (example)*pH value and acid capacity in water*

Sample size	pH value	Temperature	Acid capacity
[mL]		[°C]	[mmol/L]
50	7.95	25.62	11.20
50	7.96	25.50	11.20
50	7.97	25.52	11.20
50	7.97	25.46	11.20
50	7.97	25.42	11.20
50	7.98	25.41	11.19
50	7.98	25.33	11.18
50	7.99	25.29	11.19
50	8.00	25.21	11.18
50	8.00	25.18	11.18
Mean value	7.98	25.39	11.19
abs. std. dev.	0.02	0.14	0.01
rel. std. dev. %	0.21	0.55	0.08

3.3 Determination of the total titratable acidity in wine

Recommended accessories

- 20 mL Dosing Unit: 6.3032.220
- Comb. pH glass electrode for aqueous titrations 6.0257.000 Aquatrode Plus with 6.2104.020 electrode cable

Reagents

- Titrant: $c(\text{NaOH}) = 0.1 \text{ mol/L}$
- Possibly nitrogen from a pressure cylinder

General

The total titratable acidity is understood to be that fraction of acids contained in the must and/or wine (with the exception of carbonic acid) that are determined when the must is neutralized with NaOH to an agreed or predefined pH value. These titratable acids are mainly weak acids (tartaric acid, malic acid, etc.). Their neutral point is above $\text{pH} = 7.0$ (at approx. $\text{pH} = 8.2$). It is therefore obvious that lower values will be obtained by titrating to $\text{pH} = 7.0$ than when titrating to $\text{pH} = 8.2$. To be able to compare the analytical values it is essential that the titration is carried out to the agreed pH value. As this latter is defined differently, it is best to record the whole titration curve and to allow the titrant consumption up to the given pH value to be recalculated (by the titrator).

- CH, EU, Israel and RSA Titration to $\text{pH} = 7.0$
- Au and USA Titration to $\text{pH} = 8.2$ or to the point of inflection

General procedure

The sample is degassed by passing nitrogen through it, by briefly boiling and cooling it down rapidly or under vacuum (CO_2 removal). The given sample volume is diluted with CO_2 -free dist. H_2O and titrated immediately with $c(\text{NaOH}) = 0.1 \text{ mol/L}$.

Detailed procedures

A) CH, EU, Israel, RSA

Manual titration: Titrate 10 mL degassed sample plus 10 mL dist. H_2O (CO_2 -free) with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ up to $\text{pH} = 7.0$.

Automated titration: It is recommended to use 50 mL dist. H_2O (CO_2 -free) to make sure that the electrode is immersed in the sample.

B) Au

Titrate 10 mL degassed sample plus 50 mL dist. H_2O (CO_2 -free) with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ up to $\text{pH} = 8.2$ or, preferably, to the inflection point of the titration curve.

C) USA

Titrate 5 mL degassed sample plus 100 mL dist. H_2O (CO_2 -free) with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ up to $\text{pH} = 8.2$ or, preferably, to the inflection point of the titration curve.

D) SA

Titrate 20 mL degassed sample with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ up to $\text{pH} = 7.0$. Result also given in $\text{g/L H}_2\text{SO}_4$.

Calculation

The total titratable acidity is usually given in g/L tartaric.

total acidity at pH 7 = 'DET pH.FP{1}.VOL' * 'DET pH.TITER' * 'DET pH.CONC' * 150.1 / 'MV.Sample size'

total acidity at pH 8.2 = 'DET pH.FP{2}.VOL' * 'DET pH.TITER' * 'DET pH.CONC' * 150.1 / 'MV.Sample size'

where:

'DET pH.FP{1}.VOL'	mL NaOH up to pH = 7.0
'DET pH.FP{2}.VOL'	mL NaOH up to pH = 8.2
'DET pH.TITER'	titer of the NaOH solution
'DET pH.CONC'	concentration of NaOH solution [mol/L]
150.1	molecular weight of tartaric acid [g/mol]
'MV.Sample size'	sample volume [mL]

Remarks

The determination of the total titratable acidity is important for:

- musts, to be able to adjust the acidity and add the correct amount of SO₂;
- wines, to be able to follow the change in pH and the tartaric acid concentration.

After acid degradation the values for the **total titratable acidity** normally lie between **4.0 and 6.5 g/L** (red wines 4.0...5.5 g/L, «dry» wines 6.0...9.0 g/L).

3.4 *tiamo* method: Determination of the total titratable acidity in wine

Application note

This is a method to determine the acidity in wine. The total titratable acidity is understood to be the fraction of acids contained in the wine that are determined when the wine is neutralized with sodium hydroxide to a predefined pH value. These acids are mainly weak acids and their neutral point is above $\text{pH} = 7.0$. It is therefore obvious that lower values will be obtained by titrating to $\text{pH} = 7.0$ than when titrating to $\text{pH} = 8.2$. To be able to compare the analytical values it is essential that the titration is carried out to the agreed pH value. As there exist different defined values it is best to record the whole titration curve and to allow the titrant consumption up to the given pH value to be recalculated.

CH, EU, Israel and RSA: Titration to $\text{pH} = 7$

Au and USA: Titration to $\text{pH} = 8.2$ or to the inflection point

The method consists of six different tracks. In the start track the rack is initialized and the system prepared for the determination series. The main track handles the moving of the electrode to the sample beaker and the titration. The electrode is rinsed after every determination according to the procedure in the rinsing track. In the exit track the results are calculated, a report is printed and the data is stored in the predefined database 'Robotic Acid Base Analyzer' (can be modified). When the determination is finished the series end track is executed in which the electrode is moved to a storage beaker where it is stored in Metrohm storage solution for a proper electrode treatment between the determination series. In case of an error the error track is carried out which guarantees that the electrode is moved to the storage beaker.

Sample preparation and procedure

The sample has to be degassed before the determination either by passing nitrogen through it, by briefly boiling and cooling it down rapidly or under vacuum.

Depending on the country the sample preparation has to be adapted.

CH, EU, Israel and RSA:

Manual titration: 10 mL of the degassed sample plus 10 mL of dist. water are titrated with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ up to $\text{pH} = 7.0$.

Automated titration: It is recommended to use 50 mL of dist. water to make sure that the electrode is immersed in the sample.

Au:

10 mL of the degassed sample plus 50 mL of dist. water are titrated with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ up to $\text{pH} = 8.2$ or to the inflection point of the titration curve.

USA:

5 mL of the degassed sample plus 100 mL of dist. water are titrated with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ up to $\text{pH} = 8.2$ or to the inflection point of the titration curve.

The value 150.1 in the calculation formula is the molecular weight of tartaric acid in $[\text{g/mol}]$. Two values for the total acidity will be obtained, one at a pH value of 7.0 and one at 8.2.

Remark

To run this method the settings of the 855 Robotic Titrosampler have to be adjusted. The lift positions have to be defined according to the rack and beakers used. As the robotic swing arm is equipped with a beaker sensor the head has to touch the beaker brim when moving down to work position. If that requirement can not be fulfilled no beaker will be recognized. One special beaker has to be defined including its own work position. Special beaker 1 is used for a proper electrode treatment between the determination series.

Result report



Robotic Acid / Base Analyzer
Total acidity of wine

Programm version tiamo 1.1

2005-09-13 10:55:05 UTC+2

Results report

Determination

Method Total acidity in wine
 Method saving date 2005-08-23 10:59:21 UTC+2
 Method version 1
 Method state original
 Determination ID -1bd7dce7:105e1e3c5c4:-7d79
 Determination start 2005-08-23 11:50:49 UTC+2
 Determination state original
 Determination version 1
 Run number 24
 User (full name) Metrohm

Sample data

Sample Identification 1 acidity in wine
 Sample size 10 mL


End points

DET pH **DET pH.1**

FP1	7.000	pH	7.4442	mL
FP2	8.200	pH	7.8674	mL
FP3	invalid	pH	invalid	mL

Results

total acidity at pH 7 10.93 g/L
total acidity at pH 8,2 11.55 g/L



Statistical results (example)
Total acidity of wine

Sample size	Content tartaric acid at pH 7	Content tartaric acid at pH 8.2
[mL]	[g/L]	[g/L]
10	10.94	11.57
10	10.91	11.52
10	10.92	11.54
10	10.90	11.53
10	10.91	11.52
10	10.93	11.55
10	10.92	11.54
10	10.91	11.54
10	10.94	11.60
10	10.91	11.56
Mean value	10.92	11.55
abs. std. dev.	0.01	0.02
rel. std. dev. %	0.13	0.21

3.5 Determination of different pharmaceutical substances

Example titration with citric acid

Recommended accessories

- 20 mL Dosing Unit: 6.3032.220
- Comb. pH glass electrode for aqueous titrations 6.0257.000 Aquatrode Plus with 6.2104.020 electrode cable

Reagents

- $c(\text{NaOH}) = 0.1 \text{ mol/L}$
- $c(\text{NaOH}) = 1.0 \text{ mol/L}$
- Dist. H_2O , CO_2 -free (or demineralized H_2O)
- Ethanol and methanol, analytical grade
- Sorbitol

Citric acid (USP)

Weigh approx. 80 mg sample into the titration beaker and dissolve in 50 mL dist. H_2O . Titrate with $c(\text{NaOH}) = 1 \text{ mol/L}$ past the last endpoint.

Molecular weight ($\text{C}_6\text{H}_8\text{O}_7$) = 192.12 g/mol

Calculation

Mass fraction w in % = 'DET pH.EP{3}.VOL' * 'DET pH.TITER' * 'DET pH.CONC' * 192.12 * 100 / ('MV.Sample size' 3)

where:

'DET pH.EP{3}.VOL'	mL NaOH until third endpoint is reached
'DET pH.TITER'	titer of NaOH solution
'DET pH.CONC'	concentration of NaOH solution [mol/L]
192.12	molecular weight of citric acid [g/mol]
100	factor for conversion to %
'MV.Sample size'	sample weight [mg]

The following substances can also be titrated according to this method, though the calculation formula has to be adjusted:

N-Acetyltyrosine (Pharm. Europe)

Weigh approx. 0.180 g sample into the titration beaker and dissolve in 50 mL dist. H_2O . Titrate with $c(\text{NaOH}) = 0.1 \text{ mol/L}$.

Molecular weight ($\text{C}_{11}\text{H}_{13}\text{NO}_4$) = 223.22 g/mol

Allantoin (Pharm. Europe)

Weigh approx. 0.120 g sample into the titration beaker and dissolve in 40 mL dist. H_2O . Titrate with $c(\text{NaOH}) = 0.1 \text{ mol/L}$.

Molecular weight ($\text{C}_4\text{H}_6\text{N}_4\text{O}_3$) = 158.11 g/mol

Biotin (USP)

Weigh approx. 0.500 g (0.250 g) sample into the titration beaker, add 100 mL dist. H₂O and heat to 50 °C under stirring. Titrate slowly with c(NaOH) = 0.1 mol/L past the first endpoint (suspension).

Molecular weight (C₁₀H₁₆N₂O₃S) = 244.31 g/mol

Busulfan (USP)

Weigh approx. 0.160 g sample into an Erlenmeyer flask and dissolve in 50 mL dist. H₂O. Titrate with c(NaOH) = 0.1 mol/L past the first endpoint. Mount a reflux condenser and boil for 30 min. After the mixture has cooled down, titrate with c(NaOH) = 0.1 mol/L past the first endpoint. Use the total NaOH consumption for the calculation.

Molecular weight (C₆H₁₄O₆S₂) = 246.30 g/mol

Cilazapril (Pharm. Europe)

Weigh approx. 0.300 g sample into the titration beaker, dissolve in 25 mL ethanol, add 50 mL dist. H₂O and titrate with c(NaOH) = 0.1 mol/L.

Molecular weight (C₂₂H₃₁N₃O₅) = 417.50 g/mol

Enalapril maleate (Pharm. Europe)

Weigh approx. 0.100 g sample into the titration beaker and dissolve in 30 mL dist. H₂O. Titrate with c(NaOH) = 0.1 mol/L.

Molecular weight (C₂₄H₃₂N₂O₉) = 492.51 g/mol

Glutamic acid (Pharm. Europe)

Weigh approx. 0.130 g sample into the titration beaker and dissolve in 50 mL dist. H₂O under slight warming. Titrate with c(NaOH) = 0.1 mol/L.

Molecular weight (C₅H₉NO₄) = 147.13 g/mol

Histidine hydrochloride monohydrate (Pharm. Europe)

Weigh approx. 0.160 g sample into the titration beaker and dissolve in 50 mL dist. H₂O. Titrate with c(NaOH) = 0.1 mol/L.

Molecular weight (C₆H₁₀ClN₃O₂) = 191.61 g/mol

Lisinopril dihydrate (Pharm. Europe)

Weigh approx. 0.350 g sample into the titration beaker and dissolve in 50 mL dist. H₂O. Titrate with $c(\text{NaOH}) = 0.1 \text{ mol/L}$.

Molecular weight ($\text{C}_{21}\text{H}_{31}\text{N}_3\text{O}_5$) = 405.48 g/mol

Mandelic acid (USP)

Weigh approx. 0.500 g (0.200 g) previously dried sample into the titration beaker and dissolve in 100 mL dist. H₂O. Titrate with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ past the first endpoint.

Molecular weight ($\text{C}_8\text{H}_8\text{O}_3$) = 152.12 g/mol (C01)

Methimazole (USP)

Weigh approx. 0.250 g sample into the titration beaker and dissolve in 75 mL dist. H₂O. Under stirring add 15.00 mL $c(\text{NaOH}) = 0.1 \text{ mol/L}$ and 30 mL AgNO₃ solution. Titrate with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ past the first endpoint. Use the total NaOH consumption for the calculation, i.e. 15 mL + mL EP1.

Molecular weight ($\text{C}_4\text{H}_6\text{N}_2\text{S}$) = 114.17 g/mol

Phenylmercuric borate (Pharm. Europe)

Weigh approx. 0.600 g sample into the titration beaker and dissolve in 25 mL dist. H₂O under slight warming. Dissolve 10 g sorbitol in the warm solution, leave to cool and titrate with $c(\text{NaOH}) = 0.1 \text{ mol/L}$.

Molecular weight (H_3BO_3) = 61.83g/mol

Potassium hydrogen tartrate (Pharm. Europe)

Weigh approx. 0.170 g sample into the titration beaker and dissolve in 100 mL dist. H₂O. Heat and titrate the warm solution with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ past the first endpoint.

Molecular weight ($\text{C}_4\text{H}_5\text{KO}_6$) = 188.17g/mol

Procarbazine hydrochloride (USP)

Weigh approx. 0.750 g (0.250 g) sample into the titration beaker and dissolve in 100 mL dist. H₂O. Titrate with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ past the first endpoint.

Molecular weight = ($\text{C}_{12}\text{H}_{20}\text{ClN}_3\text{O}$) = 257.76 g/mol

Ramipril (Pharm. Europe)

Weigh approx. 0.300 g sample into the titration beaker and dissolve in 25 mL methanol. Add 25 mL dist. H₂O and titrate with $c(\text{NaOH}) = 0.1 \text{ mol/L}$.

Molecular weight ($\text{C}_{23}\text{H}_{32}\text{N}_2\text{O}_5$) = 416.15 g/mol

Saccharin calcium (USP)

Weigh 0.500 g (0.200 g) sample into a separating funnel, add 10 mL dist. H₂O and 6 mL $c(\text{HCl}) = 1 \text{ mol/L}$ and shake/extract first with 30 mL and subsequently five times with 20 mL

chloroform/ethanol 9:1 (vol/vol). Evaporate the combined extracts to dryness and dissolve in 40 mL ethanol. Add 40 mL dist. H₂O and titrate with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ past the first end-point.

Molecular weight ($\text{C}_{14}\text{H}_8\text{CaN}_2\text{O}_6\text{S}_2$) = 404.43 g/mol

Saccharin sodium (USP)

Same sample preparation and titration as under saccharin calcium.

Molecular weight ($\text{C}_7\text{H}_4\text{NNaO}_3\text{S}$) = 205.17 g/mol

Tiaprofenic acid (Pharm. Europe)

Weigh approx. 0.500 g sample into the titration beaker and dissolve in 25 mL ethanol. Add 25 mL dist. H₂O and titrate with $c(\text{NaOH}) = 0.1 \text{ mol/L}$.

Molecular weight ($\text{C}_{14}\text{H}_{12}\text{O}_3\text{S}$) = 260.31 g/mol

3.6 *tiamo* method (example): Titration of citric acid

Application note

This is a method to determine the content of citric acid. The method consists of five different tracks. In the start series track the rack is initialized and the system prepared for the determination series. The main track handles the dissolving of the sample and the titration followed by the rinsing of the electrode, which is done after every determination. In the exit track the results are calculated, a report is printed and the data is saved in the predefined database Robotic Acid Base Analyzer (can be modified). In the series end track, the electrode is moved to a storing beaker where it is stored in Metrohm storage solution for a proper electrode treatment between the determination series. In case of an error, the error track is carried out, which guarantees that the electrode is moved to the storage beaker.

Sample preparation

Approximately 0.08 g of the sample are weighed out in a titration beaker.

Procedure

The beakers are placed on the rack and the determination series started. Approximately 60 mL dist. water are added automatically to the sample for dissolving. The addition is followed by a waiting time of 30 seconds while the solution is mixed thoroughly. In case the sample can not be dissolved completely within these 30 seconds the waiting time should be adjusted. The solution is then titrated with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ and the titration stopped after 20 mL titrant have been consumed. Then the electrode is rinsed before the next determination.

The value '192.12' in the calculation formula is the molecular weight of citric acid in [g/mol].

Remarks

To run this method the settings of the 855 Robotic Titrosampler have to be adjusted. The lift positions have to be defined according to the rack and beakers used. As the robotic swing arm is equipped with a beaker sensor the head has to touch the beaker brim when moving to the work position. If this requirement can not be fulfilled no beaker will be recognized.

One special beaker has to be defined including its own work position. Special beaker 1 is used for a proper electrode treatment between determination series.

Result report



Robotic Acid / Base Analyzer
Citric acid

Programm version tiamo 1.1

2005-08-23 16:11:31 UTC+2

Results report

Determination

Method Titration of citric acid
 Method saving date 2005-08-23 08:24:37 UTC+2
 Method version 1
 Method state original
 Determination ID -1bd7dce7:105e1e3c5c4:-7f0b
 Determination start 2005-08-23 08:46:24 UTC+2
 Determination state original
 Determination version 1
 Run number 8
 User (full name) Metrohm

Sample data


Sample Identification 1 citric acid
 Sample size 0.0493 g

End points

DET pH **DET pH.1**
 EP1 8.802 pH 7.8387 mL

Results

percent by weight 99.59 %



Statistical results (example)*Determination of citric acid*

Sample size	Percent by weight
[g]	[%]
0.0758	99.32
0.0589	99.49
0.0554	99.58
0.0719	98.81
0.0630	99.36
0.0493	99.59
0.0580	99.13
0.0690	99.24
0.0488	99.30
0.0733	99.07
Mean value	99.29
abs. std. dev.	0.24
rel. std. dev. %	0.24

3.7 Determination of phosphoric acid in cola drinks

Recommended accessories

- Comb. pH glass electrode for aqueous titrations 6.0257.000 Aquatrode Plus with 6.2104.020 electrode cable
- 20 mL Dosing Unit: 6.3032.220

Reagents

- Titrant: $c(\text{NaOH}) = 0.1 \text{ mol/L}$

Sample preparation

Cola drink is degassed in an ultrasonic bath for 5 min (degassing can also be carried out under vacuum).

Analysis

50 mL of the degassed sample solution is titrated with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ to after the first endpoint (first phosphoric acid step at $\text{pH} = \text{approx. } 4.6$).

Calculations

$$\text{g/L H}_3\text{PO}_4 = \text{'DET pH.EP\{1\}.VOL'} * \text{'DET pH.TITER'} * \text{'DET pH.CONC'} * 97.993 / \text{'MV.Sample size'}$$

where:

'DET pH.EP{1}.VOL'	mL NaOH up to endpoint
'DET pH.TITER'	titer of NaOH solution
'DET pH.CONC'	concentration of NaOH solution [mol/L]
97.993	molecular weight of phosphoric acid [g/mol]
'MV.Sample size'	sample volume [mL]

Remarks

The determination only works in cola drinks that contain no citric acid, as the citric acid is also determined. Such samples must be analyzed according to a different method, in which the sample is ashed (after previous neutralization to $\text{pH} = \text{approx. } 8$ with NaOH).

The 761 SD Compact IC is an ion chromatography system that has been developed by Metrohm specially for the rapid and reliable determination of phosphoric acid (and citric acid) in these drinks. This instrument is recommended worldwide by leading manufacturers of cola drinks.

3.8 *tiamo* method: Determination of phosphoric acid in cola drinks

Application note

This is a method to determine the content of phosphoric acid in cola drinks. It only works in cola drinks that contain no citric acid, because the citric acid is also determined.

The method is built up of six different tracks. In the start track the rack is initialized and the system prepared for the measurement. The main track handles the measurement as the electrode is moved to the sample beaker and immersed in the solution. After the measurement the electrode is moved to a special beaker, where the rinsing takes place. For this step a special beaker has to be defined. The electrode is thoroughly rinsed with water after every determination. In the exit track the results are calculated, a report is printed and the data is stored in the predefined database Robotic Acid Base Analyzer. The end track contains the command to move the electrode to second special beaker (has to be defined). This beaker is filled with storing solution in order to store the electrode in an adequate way.

In case of an error the error track is carried out which guarantees that the electrode is moved to the storing beaker.

Sample preparation

The cola drink has to be degassed either in an ultrasonic bath (for 5 minutes) or under vacuum.

Procedure

50 mL of the degassed sample are titrated with $c(\text{NaOH}) = 0.1 \text{ mol/L}$ and the titration stopped after the first endpoint.

Remark

To run this method the settings of the 855 Robotic Titrosampler have to be adjusted. The lift positions have to be defined according to the rack and beakers used. As the robotic swing arm is equipped with a beaker sensor the head has to touch the beaker brim when moving to work position. If this requirement can not be fulfilled no beaker will be recognized.

For the rinsing as well as for the storing procedure a special beaker has to be defined. Each special beaker requires an own work position.

The value '97.993' in the calculation formula is the molecular weight of phosphoric acid in [g/mol].

Result report



Robotic Acid / Base Analyzer
Phosphoric acid in cola drinks

Programm version tiamo 1.1

2005-08-23 16:03:31 UTC+2

Results report

Determination

Method Determination of phosphoric acid in cola drinks
 Method saving date 2005-08-23 14:55:12 UTC+2
 Method version 1
 Method state original
 Determination ID -1bd7dce7:105e1e3c5c4:-7b42
 Determination start 2005-08-23 14:59:54 UTC+2
 Determination state original
 Determination version 1
 Run number 47
 User (full name) Metrohm

Sample data


Sample Identification 1 phosphoric acid in cola drinks
 Sample size 50 mL

End points

DET pH **DET pH.1**
 EP1 4.372 pH 2.7134 mL

Results

content phosphoric acid 0.520 g/L



Statistical results (example)*Phosphoric acid in cola drinks*

Sample size	Phosphoric acid
[mL]	[g/L]
50	0.516
50	0.518
50	0.518
50	0.517
50	0.518
50	0.520
50	0.520
50	0.519
50	0.520
50	0.518
Mean value	0.518
abs. std. dev.	0.0013
rel. std. dev. %	0.26

4 Troubleshooting pH glass electrodes

4.1 Troubleshooting glass membrane

Source of error	Effects	Cleaning	Alternatives
HF-containing solutions	Etching and dissolution of the glass membrane → corrosion potential during the measurement/short working life		Use of the Sb electrode
High pH value and high alkali content	Increased alkali error → pH too low		Use of electrodes with U glass
High temperatures	Rapid rise in membrane resistance by aging → increased polarizability and drift		Use of electrodes with U glass. Careful treatment with etching salt and HCl
Measurements at low temperatures	High membrane resistance → polarization effects		Use of electrodes with T glass and Idrolyte as reference electrolyte
Dry storage	Zero point drift	Storage in 6.2323.000 storage solution	
Non-aqueous media	Reduced sensitivity	Store in water	T glass / non-aqueous electrolyte solution
Deposition of solids on membrane surface	Slow response, zero point shift, slope reduction	Solvent or strong acids	
Deposition of proteins on membrane surface	Slow response, zero point shift, slope reduction	5% pepsin in 0.1 mol/L HCl	

4.2 Cleaning and care of pH-electrodes with fixed-sleeve diaphragm

Diaphragm type	Type of contamination	Cleaning
General	Precipitates of silver halides and silver sulfides	Immerse diaphragm for several hours in a solution of 7% thiourea in 0.1 mol/L HCl.
	Proteins, polypeptides	Immerse diaphragm for several hours in a solution of 5% pepsin in 0.1 mol/L HCl.
	Suspensions, solids, resins, glues, oils, fats	Clean electrode with suitable solvent
Fixed ground joint	All types of contamination	Aspirate off reference electrolyte and immerse electrode in the corresponding cleaning solution.
Separable ground joint	All types of contamination	Loosen the ground-joint sleeve (using hot water if necessary) and clean according to the type of contamination.

4.3 Storage

Aquatrode Plus: Store in 6.2323.000 storage solution only! The filling plug for the bridge electrolyte must not be opened.

5 Method reports

5.1 Calibrating the pH glass electrode



License ID 124049905 Program version tiamo 1.1 - 31
 Client name TITRATION14
 User Metrohm 2005-09-13 09:17:19 UTC+2

Method parameters

Method Calibrating the pH glass electrode
 Method saving date 2005-09-13 09:17:07 UTC+2
 Method version 1
 Method group Robotic Acid Base Analyzer
 Method status original
 Method saved by (full name) Metrohm
 Method saved by (short name) Metrohm

START

Main track

General

Workplace view
 Current view on
 Track view for live window
 Live display 1 Main track
 Live display 2 Main track
 Statistics off
 Conditioning
 Automatic conditioning off

Application note

See attached documents

Method variables

Name	Type	Assignment	Value	Comment	Monitoring
Sample size	Number	Sample size		Sample size	off
Sample size unit	Text	Sample size unit		Sample size unit	off
Sample position	Number	Sample position		Sample position number	off
ID1	Text	ID1		Sample identification 1	off
ID2	Text	ID2		Sample identification 2	off
ID3	Text	ID3		Sample identification 3	off

Name **Sample size**
 Type Number
 Assignment on, Sample size
 Value off,
 Check at start on
 Comment Sample size
 Variable monitoring off
 Lower limit
 Upper limit
 Message
 Display message on
 Record message on
 Message by e-mail off



License ID 124049905
Client name TITRATION14
User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 09:17:19 UTC+2

Mail to
Subject Message from tiamo - Method 'New method 1' - Command 'Main track'
User
Mail from
SMTP Server
POP3 Server
Acoustic signal off
Action off
Stop determination on
Stop determination and series off

Name **Sample position**
Type Number
Assignment on, Sample position
Value off,
Check at start on
Comment Sample position number
Variable monitoring off
Lower limit
Upper limit
Message
Display message on
Record message on
Message by e-mail off
Mail to
Subject Message from tiamo - Method 'New method 1' - Command 'Main track'
User
Mail from
SMTP Server
POP3 Server
Acoustic signal off
Action off
Stop determination on
Stop determination and series off

Name **Sample size unit**
Type Text
Assignment on, Sample size unit
Value off,
Check at start on
Comment Sample size unit

Name **ID1**
Type Text
Assignment on, ID1



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 09:17:19 UTC+2

Value off
 Check at start on
 Comment Sample identification 1

Name ID2
 Type Text
 Assignment on. ID2
 Value off.
 Check at start on
 Comment Sample identification 2

Name ID3
 Type Text
 Assignment on. ID3
 Value off.
 Check at start on
 Comment Sample identification 3

**CAL LOOP
pH**

calibration loop pH
 Buffers
 Number of buffers 3
 Buffer type Metrohm
 Request for buffer exchange off
 Subsequent command calculation

MOVE

to sample
 Device
 Device name 855_1
 Target
 Tower 1
 Move Rack position
 Number =0 + 'calibration loop pH.LCO'
 Beaker test
 Display message off
 Stop determination off
 Stop determination and series on
 Parameters
 Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

**CAL MEAS
pH**

CAL MEAS pH
General/Hardware
 Device
 Device name 855_1
 Sensor
 Measuring input 1



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31
 2005-09-13 09:17:19 UTC+2

Sensor Aquatrode Plus
 Temperature measurement automatic
 Stirrer
 Stirrer 1
 Stirring rate 8
 Switch off automatically on

Measuring parameters

Measurement with drift control
 Signal drift 2.0 mV/min
 Min. waiting time 10 s
 Max. waiting time 110 s
 Measuring interval 2.0 s
 Temperature
 Temperature 25.0 °C

MOVE to rinsing beaker

Device
 Device name 855_1
 Target
 Tower 1
 Move Special beaker
 Number 2
 Beaker test
 Display message off
 Stop determination on
 Stop determination and series off
 Parameters
 Shift rate 20 °/s
 Shift direction auto
 Swing rate 55 °/s

PUMP rinse and aspirate

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 1+2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 10.0 s

PUMP aspirate

Device
 Device name 855_1
 Pumps
 Tower 1



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 09:17:19 UTC+2

Pump(s) 2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 5 s

CALC calculation

Result name	Formula	Unit	Decimal places	Assignment	Statistics
slope	= 'calibration loop pH.SLO'	%	2	RS01	off
zero point	= 'calibration loop pH.ENP'		2	RS02	off

Result name **slope**
 Formula = 'calibration loop pH.SLO'
 Unit %
 Decimal places 2
 Assignment RS01
 Statistics off
 Description RS.'Result name'[.VAL]Result value.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer off
 Solution name

Result name **zero point**
 Formula = 'calibration loop pH.ENP'
 Unit
 Decimal places 2
 Assignment RS02
 Statistics off
 Description RS.'Result name'[.VAL]Result value.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer off
 Solution name

REPORT report

Report template
 Report template calibration
 Report output
 Printer off
 PDF file on
 PDF file C:\Programme\Metrohm\tiamo\Data\Calibrating the pH Glass Electrode.pdf



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31
 2005-09-13 09:17:19 UTC+2

Send e-mail off

DATABASE database

Database
 Robotic Acid Base Analyzer

SERIES START Series start track

RACK initialize rack

Device
 Device name 855_1
 Rack test off

MOVE to rinse beaker

Device
 Device name 855_1
 Target
 Tower 1
 Move Special beaker
 Number 2
 Beaker test
 Display message off
 Stop determination on
 Stop determination and series off
 Parameters
 Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

PUMP aspirate and rinse

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 1+2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 10.0 s

PUMP empty beaker

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2



License ID 124049905
Client name TITRATION14
User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 09:17:19 UTC+2

Action

Switch on off
Switch off off
Duration on
Time 5 s

SERIES
END

Series end track

MOVE

to storage beaker

Device

Device name 855_1

Target

Tower 1

Move Special beaker

Number 1

Beaker test

Display message off

Stop determination on

Stop determination and series off

Parameters

Shift rate 20 °/s

Shift direction auto

Swing rate 55 °/s

5.2 Titer determination of alkaline titrants



License ID 124049905 Program version tiamo 1.1 - 31
 Client name TITRATION14
 User Metrohm 2005-09-13 09:26:24 UTC+2

Method parameters

Method Titer determination of alkaline titrants
 Method saving date 2005-09-13 09:25:27 UTC+2
 Method version 1
 Method group Robotic Acid Base Analyzer
 Method status original
 Method saved by (full name) Metrohm
 Method saved by (short name) Metrohm

START Main track

General

Workplace view
 Current view on
 Track view for live window
 Live display 1 Main track
 Live display 2 Main track
 Statistics on
 Number of single determinations not defined
 Conditioning
 Automatic conditioning off

Application note

See attached documents

Method variables

Name	Type	Assignment	Value	Comment	Monitoring
Sample size	Number	Sample size		Sample size	off
Sample size unit	Text	Sample size unit		Sample size unit	off
Sample position	Number	Sample position		Sample position number	off
ID1	Text	ID1		Sample identification 1	off
ID2	Text	ID2		Sample identification 2	off
ID3	Text	ID3		Sample identification 3	off

Name **Sample size**
 Type Number
 Assignment on. Sample size
 Value off.
 Check at start on
 Comment Sample size
 Variable monitoring off
 Lower limit
 Upper limit
 Message
 Display message on
 Record message on



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 09:26:24 UTC+2

Message by e-mail off
 Mail to
 Subject Message from tiamo - Method 'New method 5' - Command 'Main track'
 User
 Mail from
 SMTP Server
 POP3 Server
 Acoustic signal off
 Action off
 Stop determination on
 Stop determination and series off

Name **Sample position**
 Type Number
 Assignment on. Sample position
 Value off.
 Check at start on
 Comment Sample position number
 Variable monitoring off
 Lower limit
 Upper limit
 Message
 Display message on
 Record message on
 Message by e-mail off
 Mail to
 Subject Message from tiamo - Method 'New method 5' - Command 'Main track'
 User
 Mail from
 SMTP Server
 POP3 Server
 Acoustic signal off
 Action off
 Stop determination on
 Stop determination and series off

Name **Sample size unit**
 Type Text
 Assignment on. Sample size unit
 Value off.
 Check at start on
 Comment Sample size unit

Name **ID1**
 Type Text



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 09:26:24 UTC+2

Assignment on. ID1
 Value off.
 Check at start on
 Comment Sample identification 1

Name **ID2**
 Type Text
 Assignment on. ID2
 Value off.
 Check at start on
 Comment Sample identification 2

Name **ID3**
 Type Text
 Assignment on. ID3
 Value off.
 Check at start on
 Comment Sample identification 3

MOVE to sample position

Device
 Device name 855_1
 Target
 Tower 1
 Move Sample position
 Beaker test
 Display message off
 Stop determination off
 Stop determination and series on
 Parameters
 Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

PUMP add water

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 1
 Action
 Switch on off
 Switch off off
 Duration on
 Time 9 s

STIR stirrer on



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 09:26:24 UTC+2

Device

Device name 855_1

Stirrer

Stirrer 1
 Stirrer type unknown
 Stirring rate 10

Action

Switch on on
 Switch off off
 Duration off

WAIT

dissolving

Wait

Stop track and waiting for [Continue] off
 Stop all tracks and waiting for [Continue] off
 Waiting time on
 Time 90
 Unit s

Message

Record message on
 Message by e-mail off
 Acoustic signal off

DET pH

DET pH

General/Hardware

Device

Device name 855_1

Dosing device

Dosing device 1
 Solution NaOH

Sensor

Measuring input 1
 Sensor Aquatrode Plus
 Temperature measurement automatic

Stirrer

Stirrer off

Start conditions

Initial measured value

Signal drift off mV/min
 Min. waiting time 0 s
 Max. waiting time 1 s

Start volume

Start volume 0 mL
 Dosing rate maximum mL/min

Start measured value

Start measured value pH off
 Dosing rate 5 mL/min



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31
 2005-09-13 09:26:24 UTC+2

Start slope
 Start slope off pH/mL
 Dosing rate 5 mL/min
 Pause
 Pause 0 s

Titration parameters

Titration rate
 Titration rate user
 Measured value acceptance
 Signal drift 50.0 mV/min
 Min. waiting time 0 s
 Max. waiting time 38 s
 Dosing of increments
 Measuring point density 4
 Min. increment 50.0 µL
 Max. increment off µL
 Dosing rate maximum mL/min
 Temperature
 Temperature 25.0 °C

Stop conditions

Stop volume 20 mL
 Stop measured value pH off
 Stop EP 1
 Volume after EP 0.5 mL
 Stop time off s
 Filling rate maximum mL/min

Potentiometric evaluation

Evaluation without window on
 EP criterion 5
 EP recognition all
 Evaluation with measured value window (pH) off
 Evaluation with volume window (mL) off

Additional evaluations

Fix end point evaluation off
 pK/HNP evaluation off
 Minimum evaluation off
 Maximum evaluation off
 Break point evaluation off

Additional measured values

Additional calculated measured values off
 Additional external measured values off

STIR

stirrer off

Device
 Device name 855_1
 Stirrer
 Stirrer 1



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 09:26:24 UTC+2

Stirrer type unknown
 Stirring rate 10
 Action
 Switch on off
 Switch off on
 Duration off

PUMP aspirate sample solution

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 5 s

PUMP rinse and aspirate

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 1+2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 5 s

PUMP aspirate

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 3 s

SERIES START Series start track

RACK initialize rack

Device
 Device name 855_1



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 09:26:24 UTC+2

Rack test off

EXIT Exit track

CALC calculation

Result name	Formula	Unit	Decimal places	Assignment	Statistics
Titer	= 'MV.Sample size' * 1000 / ('DET pH.EP {1}.VOL' * 'DET pH.CONC' * 204.23)		4	RS01	on
mean value titer	= 'RS.Titer.MNV'		4	RS02	off

Result name **Titer**
 Formula = 'MV.Sample size' * 1000 / ('DET pH.EP {1}.VOL' * 'DET pH.CONC' * 204.23)
 Unit
 Decimal places 4
 Assignment RS01
 Statistics on
 Description RS.'Result name' [.VAL] Result value.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer on
 Solution name NaOH

Result name **mean value titer**
 Formula = 'RS.Titer.MNV'
 Unit
 Decimal places 4
 Assignment RS02
 Statistics off
 Description RS.'Result name' [.VAL] Result value.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer on
 Solution name NaOH

REPORT report

Report template
 Report template titer determination
 Report output
 Printer off
 PDF file on
 PDF file C:\Programme\Metrohm\tiamo\Data\Titer determination of alkaline titrants.pdf
 Send e-mail off



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 09:26:24 UTC+2

DATABASE database

Database
 Robotic Acid Base Analyzer

**SERIES Series end track
 END**

MOVE to storing beaker

Device
 Device name 855_1
 Target
 Tower 1
 Move Special beaker
 Number 1
 Beaker test
 Display message off
 Stop determination on
 Stop determination and series off
 Parameters
 Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

ERROR Error track

MOVE to storage beaker

Device
 Device name 855_1
 Target
 Tower 1
 Move Special beaker
 Number 1
 Beaker test
 Display message off
 Stop determination on
 Stop determination and series off
 Parameters
 Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

5.3 Titer determination of acidic titrants



License ID 124049905 Program version tiamo 1.1 - 31
 Client name TITRATION14
 User Metrohm 2005-09-13 10:06:48 UTC+2

Method parameters

Method Titer determination of acidic titrants
 Method saving date 2005-09-13 10:06:43 UTC+2
 Method version 1
 Method group Robotic Acid Base Analyzer
 Method status original
 Method saved by (full name) Metrohm
 Method saved by (short name) Metrohm

START Main track

General

Workplace view
 Current view on
 Track view for live window
 Live display 1 Main track
 Live display 2 Main track
 Statistics on
 Number of single determinations not defined
 Conditioning
 Automatic conditioning off

Application note

See attached documents

Method variables

Name	Type	Assignment	Value	Comment	Monitoring
Sample size	Number	Sample size		Sample size	off
Sample size unit	Text	Sample size unit		Sample size unit	off
Sample position	Number	Sample position		Sample position number	off
ID1	Text	ID1		Sample identification 1	off
ID2	Text	ID2		Sample identification 2	off
ID3	Text	ID3		Sample identification 3	off

Name **Sample size**
 Type Number
 Assignment on. Sample size
 Value off.
 Check at start on
 Comment Sample size
 Variable monitoring off
 Lower limit off
 Upper limit off
 Message
 Display message on
 Record message on



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:06:48 UTC+2

Message by e-mail off
 Mail to
 Subject Message from tiamo - Method 'New method 4' - Command 'Main track'
 User
 Mail from
 SMTP Server
 POP3 Server
 Acoustic signal off
 Action off
 Stop determination on
 Stop determination and series off

Name **Sample position**
 Type Number
 Assignment on. Sample position
 Value off.
 Check at start on
 Comment Sample position number
 Variable monitoring off
 Lower limit
 Upper limit
 Message
 Display message on
 Record message on
 Message by e-mail off
 Mail to
 Subject Message from tiamo - Method 'New method 4' - Command 'Main track'
 User
 Mail from
 SMTP Server
 POP3 Server
 Acoustic signal off
 Action off
 Stop determination on
 Stop determination and series off

Name **Sample size unit**
 Type Text
 Assignment on. Sample size unit
 Value off.
 Check at start on
 Comment Sample size unit

Name **ID1**
 Type Text



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:06:48 UTC+2

Assignment on. ID1
 Value off.
 Check at start on
 Comment Sample identification 1

Name ID2
 Type Text
 Assignment on. ID2
 Value off.
 Check at start on
 Comment Sample identification 2

Name ID3
 Type Text
 Assignment on. ID3
 Value off.
 Check at start on
 Comment Sample identification 3

MOVE to sample position

Device
 Device name 855_1

Target
 Tower 1
 Move Sample position

Beaker test
 Display message off
 Stop determination off
 Stop determination and series on

Parameters
 Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

PUMP add water

Device
 Device name 855_1

Pumps
 Tower 1
 Pump(s) 1

Action
 Switch on off
 Switch off off
 Duration on
 Time 20 s

STIR stirrer on



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:06:48 UTC+2

Device

Device name 855_1

Stirrer

Stirrer 1
 Stirrer type unknown
 Stirring rate 8

Action

Switch on on
 Switch off off
 Duration off

WAIT

dissolving

Wait

Stop track and waiting for [Continue] off
 Stop all tracks and waiting for [Continue] off
 Waiting time on
 Time 90
 Unit s

Message

Record message on
 Message by e-mail off
 Acoustic signal off

DET pH

DET pH

General/Hardware

Device

Device name 855_1

Dosing device

Dosing device 1
 Solution HCl

Sensor

Measuring input 1
 Sensor Aquatrode Plus
 Temperature measurement automatic

Stirrer

Stirrer off

Start conditions

Initial measured value

Signal drift off mV/min
 Min. waiting time 0 s
 Max. waiting time 1 s

Start volume

Start volume 1 mL
 Dosing rate maximum mL/min

Start measured value

Start measured value pH off
 Dosing rate 5 mL/min



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31
 2005-09-13 10:06:48 UTC+2

Start slope
 Start slope off pH/mL
 Dosing rate 5 mL/min
 Pause
 Pause 5 s

Titration parameters

Titration rate
 Titration rate user
 Measured value acceptance
 Signal drift 50.0 mV/min
 Min. waiting time 0 s
 Max. waiting time 26 s
 Dosing of increments
 Measuring point density 4
 Min. increment 100 µL
 Max. increment off µL
 Dosing rate maximum mL/min
 Temperature
 Temperature 25.0 °C

Stop conditions

Stop volume 20 mL
 Stop measured value pH off
 Stop EP 1
 Volume after EP 0.5 mL
 Stop time off s
 Filling rate maximum mL/min

Potentiometric evaluation

Evaluation without window on
 EP criterion 5
 EP recognition all
 Evaluation with measured value window (pH) off
 Evaluation with volume window (mL) off

Additional evaluations

Fix end point evaluation off
 pK/HNP evaluation off
 Minimum evaluation off
 Maximum evaluation off
 Break point evaluation off

Additional measured values

Additional calculated measured values off
 Additional external measured values off

STIR

stirrer off

Device
 Device name 855_1
 Stirrer
 Stirrer 1



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:06:48 UTC+2

Start slope
 Start slope off pH/mL
 Dosing rate 5 mL/min
 Pause
 Pause 5 s

Titration parameters

Titration rate
 Titration rate user
 Measured value acceptance
 Signal drift 50.0 mV/min
 Min. waiting time 0 s
 Max. waiting time 26 s
 Dosing of increments
 Measuring point density 4
 Min. increment 100 µL
 Max. increment off µL
 Dosing rate maximum mL/min
 Temperature
 Temperature 25.0 °C

Stop conditions

Stop volume 20 mL
 Stop measured value pH off
 Stop EP 1
 Volume after EP 0.5 mL
 Stop time off s
 Filling rate maximum mL/min

Potentiometric evaluation

Evaluation without window on
 EP criterion 5
 EP recognition all
 Evaluation with measured value window (pH) off
 Evaluation with volume window (mL) off

Additional evaluations

Fix end point evaluation off
 pK/HNP evaluation off
 Minimum evaluation off
 Maximum evaluation off
 Break point evaluation off

Additional measured values

Additional calculated measured values off
 Additional external measured values off

STIR

stirrer off

Device
 Device name 855_1
 Stirrer
 Stirrer 1



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31
 2005-09-13 10:06:48 UTC+2

Stirrer type unknown
 Stirring rate 8
 Action
 Switch on off
 Switch off on
 Duration off

PUMP aspirate sample solution

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 7 s

PUMP rinse and aspirate

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 1+2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 5 s

PUMP aspirate

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 3 s

SERIES START Series start track

RACK initialize rack

Device
 Device name 855_1



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:06:48 UTC+2

Rack test off

EXIT Exit track

CALC calculation

Result name	Formula	Unit	Decimal places	Assignment	Statistics
Titer	= 'MV.Sample size' * 1000 / ('DET pH.EP {1}.VOL' * 'DET pH.CONC' * 121.14)		4	RS01	on
mean value titer	= 'RS.Titer.MNV'		4	RS02	off

Result name **Titer**
 Formula = 'MV.Sample size' * 1000 / ('DET pH.EP {1}.VOL' * 'DET pH.CONC' * 121.14)
 Unit
 Decimal places 4
 Assignment RS01
 Statistics on
 Description RS.'Result name'[.VAL]Result value.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer off
 Solution name

Result name **mean value titer**
 Formula = 'RS.Titer.MNV'
 Unit
 Decimal places 4
 Assignment RS02
 Statistics off
 Description RS.'Result name'[.VAL]Result value.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer on
 Solution name HCl

REPORT report

Report template
 Report template titer determination
 Report output
 Printer off
 PDF file on
 PDF file C:\Programme\Metrohm\tiamo\Data\titer determination acid titrants.pdf
 Send e-mail off



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:06:48 UTC+2

DATABASE database

Database
 Robotic Acid Base Analyzer

**SERIES Series end track
 END**

MOVE to storing beaker

Device
 Device name 855_1
 Target
 Tower 1
 Move Special beaker
 Number 1
 Beaker test
 Display message off
 Stop determination on
 Stop determination and series off
 Parameters
 Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

ERROR Error track

MOVE to storage beaker

Device
 Device name 855_1
 Target
 Tower 1
 Move Special beaker
 Number 1
 Beaker test
 Display message off
 Stop determination on
 Stop determination and series off
 Parameters
 Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

5.4 pH value and acid capacity in water



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:11:14 UTC+2

Method parameters

Method pH value and acid capacity in water
 Method saving date 2005-09-13 10:11:10 UTC+2
 Method version 1
 Method group Robotic Acid Base Analyzer
 Method status original
 Method saved by (full name) Metrohm
 Method saved by (short name) Metrohm

START

Main track

General

Workplace view
 Current view on
 Track view for live window
 Live display 1 Main track
 Live display 2 Main track
 Statistics off
 Conditioning
 Automatic conditioning off

Application note

See attached documents

Method variables

Name	Type	Assignment	Value	Comment	Monitoring
Sample size	Number	Sample size		Sample size	off
Sample size unit	Text	Sample size unit		Sample size unit	off
Sample position	Number	Sample position		Sample position number	off
ID1	Text	ID1		Sample identification 1	off
ID2	Text	ID2		Sample identification 2	off
ID3	Text	ID3		Sample identification 3	off

Name **Sample size**
 Type Number
 Assignment on Sample size
 Value off
 Check at start on
 Comment Sample size
 Variable monitoring off
 Lower limit
 Upper limit
 Message
 Display message on
 Record message on
 Message by e-mail off



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:11:14 UTC+2

Mail to
 Subject Message from tiamo - Method 'New method 2' - Command 'Main track'
 User
 Mail from
 SMTP Server
 POP3 Server
 Acoustic signal off
 Action off
 Stop determination on
 Stop determination and series off

Name **Sample position**
 Type Number
 Assignment on, Sample position
 Value off,
 Check at start on
 Comment Sample position number
 Variable monitoring off
 Lower limit
 Upper limit
 Message
 Display message on
 Record message on
 Message by e-mail off
 Mail to
 Subject Message from tiamo - Method 'New method 2' - Command 'Main track'
 User
 Mail from
 SMTP Server
 POP3 Server
 Acoustic signal off
 Action off
 Stop determination on
 Stop determination and series off

Name **Sample size unit**
 Type Text
 Assignment on, Sample size unit
 Value off,
 Check at start on
 Comment Sample size unit

Name **ID1**
 Type Text
 Assignment on, ID1



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:11:14 UTC+2

Value off
 Check at start on
 Comment Sample identification 1

Name **ID2**
 Type Text
 Assignment on ID2
 Value off
 Check at start on
 Comment Sample identification 2

Name **ID3**
 Type Text
 Assignment on ID3
 Value off
 Check at start on
 Comment Sample identification 3

MOVE

to sample

Device
 Device name 855_1
 Target
 Tower 1
 Move Sample position
 Beaker test
 Display message off
 Stop determination off
 Stop determination and series on
 Parameters
 Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

MEAS pH

MEAS pH

General/Hardware

Device
 Device name 855_1
 Sensor
 Measuring input 1
 Sensor Aquatrode Plus
 Temperature measurement automatic
 Stirrer
 Stirrer 1
 Stirring rate 8
 Switch off automatically on

Measuring parameters

Measurement



License ID 124049905 Program version tiamo 1.1 - 31
 Client name TITRATION14
 User Metrohm 2005-09-13 10:11:14 UTC+2

Measurement with drift control on
 Signal drift 5 mV/min
 Min. waiting time 30 s
 Max. waiting time 120 s
 Measuring interval 1 s
 Stop measured value pH off
 Measurement without drift control off

Temperature
 Temperature 25.0 °C

Evaluations

Fix end point evaluation off
 Minimum evaluation off
 Maximum evaluation off
 Break point evaluation off

Additional measured values

Additional calculated measured values off
 Additional external measured values off

WAIT

wait

Wait

Stop track and waiting for [Continue] off
 Stop all tracks and waiting for [Continue] off
 Waiting time on
 Time 2
 Unit s

Message

Record message off
 Message by e-mail off
 Acoustic signal off

SET pH

SET pH

General/Hardware

Device

Device name 855_1

Dosing device

Dosing device 1
 Solution HCl

Sensor

Measuring input 1
 Sensor Aquatrode Plus
 Temperature measurement automatic

Stirrer

Stirrer 1
 Stirring rate 8
 Switch off automatically on

Start conditions

Initial measured value



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31
 2005-09-13 10:11:14 UTC+2

Signal drift off mV/min
 Min. waiting time 0 s
 Max. waiting time 1 s
 Pause 1
 Pause 1 0 s
 Start volume
 Start volume 0 mL
 Dosing rate maximum mL/min
 Pause 2
 Pause 2 0 s

Control parameters

End point 1
 EP1 at pH 8.2
 Titration rate
 Titration rate user
 Control
 Dynamics pH 2
 Max. rate 10 mL/min
 Min. rate 5 µL/min
 Stop criterion
 Stop criterion drift
 Stop drift 20 µL/min
 End point 2
 EP2 at pH 4.3
 Titration rate
 Titration rate user
 Control
 Dynamics pH 3
 Max. rate 10.00 mL/min
 Min. rate 5 µL/min
 Stop criterion
 Stop criterion drift
 Stop drift 20 µL/min

Titration parameters

Titration direction auto
 Extraction time 0 s
 Temperature 25.0 °C
 Time interval measuring point 2.0 s

Stop conditions

Stop volume 20 mL
 Stop time off s
 Filling rate maximum mL/min

Conditioning

Conditioning off

Additional evaluations

Fix end point evaluation off



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31
 2005-09-13 10:11:14 UTC+2

Minimum evaluation off
 Maximum evaluation off
Additional measured values
 Additional calculated measured values off
 Additional external measured values off

PUMP aspirate sample solution

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 5 s

PUMP rinse and aspirate

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 1+2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 5 s

PUMP empty beaker

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 3 s

SERIES START Series start track

RACK initialize rack

Device
 Device name 855_1
 Rack test off



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:11:14 UTC+2

EXIT Exit track

CALC calculation

Result name	Formula	Unit	Decimal places	Assignment	Statistics
pH	= 'MEAS pH.LP.MEA'		2	RS01	off
Temperature	= 'MEAS pH.LP.TEM'	°C	2	RS02	off
K (s 8,2)	= 'SET pH.EP{1}.VOL' * 'SET pH.CONC' * 2 * 'SET pH.TITER' * 1000 / 'MV.Sample size'	mmol/L	2	RS03	off
K (s 4,3)	= 'SET pH.EP{2}.VOL' * 'SET pH.CONC' * 2 * 'SET pH.TITER' * 1000 / 'MV.Sample size'	mmol/L	2	RS04	off

Result name **pH**
 Formula = 'MEAS pH.LP.MEA'
 Unit
 Decimal places 2
 Assignment RS01
 Statistics off
 Description RS.'Result name'[.VAL]Result value.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer off
 Solution name

Result name **Temperature**
 Formula = 'MEAS pH.LP.TEM'
 Unit °C
 Decimal places 2
 Assignment RS02
 Statistics off
 Description
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer off
 Solution name

Result name **K (s 8,2)**
 Formula = 'SET pH.EP{1}.VOL' * 'SET pH.CONC' * 2 * 'SET pH.TITER' * 1000 / 'MV.Sample size'
 Unit mmol/L
 Decimal places 2
 Assignment RS03
 Statistics off
 Description RS.'Result name'[.VAL]Result value.



License ID 124049905
Client name TITRATION14
User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:11:14 UTC+2

Stop determination on
Stop determination and series off

Parameters

Shift rate 20 %/s
Shift direction auto
Swing rate 55 %/s

ERROR Error track

MOVE to storage beaker

Device

Device name 855_1

Target

Tower 1
Move Special beaker
Number 1

Beaker test

Display message off
Stop determination on
Stop determination and series off

Parameters

Shift rate 20 %/s
Shift direction auto
Swing rate 55 %/s

5.5 Determination of the total titratable acidity in wine



License ID 124049905 Program version tiamo 1.1 - 31
 Client name TITRATION14
 User Metrohm 2005-09-13 10:14:07 UTC+2

Method parameters

Method Total acidity of wine
 Method saving date 2005-09-13 10:14:02 UTC+2
 Method version 1
 Method group Robotic Acid Base Analyzer
 Method status original
 Method saved by (full name) Metrohm
 Method saved by (short name) Metrohm

START Main track

General

Workplace view
 Current view on
 Track view for live window
 Live display 1 Main track
 Live display 2 Main track
 Statistics off
 Conditioning
 Automatic conditioning off

Application note

See attached documents

Method variables

Name	Type	Assignment	Value	Comment	Monitoring
Sample size	Number	Sample size		Sample size	off
Sample size unit	Text	Sample size unit		Sample size unit	off
Sample position	Number	Sample position		Sample position number	off
ID1	Text	ID1		Sample identification 1	off
ID2	Text	ID2		Sample identification 2	off
ID3	Text	ID3		Sample identification 3	off

Name **Sample size**
 Type Number
 Assignment on Sample size
 Value off
 Check at start on
 Comment Sample size
 Variable monitoring off
 Lower limit
 Upper limit
 Message
 Display message on
 Record message on
 Message by e-mail off
 Mail to



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:14:07 UTC+2

Subject Message from tiamo - Method 'New method 7' - Command 'Main track'

User

Mail from

SMTP Server

POP3 Server

Acoustic signal off

Action off

Stop determination on

Stop determination and series off

Name **Sample position**

Type Number

Assignment on. Sample position

Value off.

Check at start on

Comment Sample position number

Variable monitoring off

Lower limit

Upper limit

Message

Display message on

Record message on

Message by e-mail off

 Mail to

 Subject Message from tiamo - Method 'New method 7' - Command 'Main track'

 User

 Mail from

 SMTP Server

 POP3 Server

Acoustic signal off

Action off

Stop determination on

Stop determination and series off

Name **Sample size unit**

Type Text

Assignment on. Sample size unit

Value off.

Check at start on

Comment Sample size unit

Name **ID1**

Type Text

Assignment on. ID1

Value off.



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:14:07 UTC+2

Check at start on
 Comment Sample identification 1

Name **ID2**
 Type Text
 Assignment on. ID2
 Value off.
 Check at start on
 Comment Sample identification 2

Name **ID3**
 Type Text
 Assignment on. ID3
 Value off.
 Check at start on
 Comment Sample identification 3

MOVE to sample position

Device
 Device name 855_1

Target
 Tower 1
 Move Sample position

Beaker test
 Display message off
 Stop determination off
 Stop determination and series on

Parameters
 Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

DET pH DET pH

General/Hardware

Device
 Device name 855_1

Dosing device
 Dosing device 1
 Solution NaOH

Sensor
 Measuring input 1
 Sensor Aquatrode Plus
 Temperature measurement automatic

Stirrer
 Stirrer 1
 Stirring rate 8
 Switch off automatically on



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:14:07 UTC+2

Start conditions

Initial measured value
 Signal drift off mV/min
 Min. waiting time 0 s
 Max. waiting time 1 s
 Start volume
 Start volume 0 mL
 Dosing rate maximum mL/min
 Start measured value
 Start measured value pH off
 Dosing rate 5 mL/min
 Start slope
 Start slope off pH/mL
 Dosing rate 5 mL/min
 Pause
 Pause 0 s

Titration parameters

Titration rate
 Titration rate user
 Measured value acceptance
 Signal drift 50.0 mV/min
 Min. waiting time 0 s
 Max. waiting time 26 s
 Dosing of increments
 Measuring point density 4
 Min. increment 10.0 µL
 Max. increment off µL
 Dosing rate maximum mL/min
 Temperature
 Temperature 25.0 °C

Stop conditions

Stop volume 20 mL
 Stop measured value pH 10
 Stop EP off
 Volume after EP off mL
 Stop time off s
 Filling rate maximum mL/min

Potentiometric evaluation

Evaluation without window on
 EP criterion 5
 EP recognition off
 Evaluation with measured value window (pH) off
 Evaluation with volume window (mL) off

Additional evaluations

Fix end point evaluation on
 Fixed quantity Measured value



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31
 2005-09-13 10:14:07 UTC+2

Fix EP1 at pH 7
 Fix EP2 at pH 8.2
 Fix EP3 at pH off
 pK/HNP evaluation off
 Minimum evaluation off
 Maximum evaluation off
 Break point evaluation off

Additional measured values

Additional calculated measured values off
 Additional external measured values off

PUMP

aspirate

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 5 s

PUMP

rinse and aspirate

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 1+2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 5 s

PUMP

aspirate rinsing solution

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 3 s

**SERIES
 START**

Series start track



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:14:07 UTC+2

RACK initialize rack

Device
 Device name 855_1
 Rack test off

EXIT Exit track

CALC calculations

Result name	Formula	Unit	Decimal places	Assignment	Statistics
total acidity at pH 7	= 'DET pH.FP{1}.VOL' * 'DET pH.TITER' * 'DET pH.CONC' * 150.1 / 'MV.Sample size'	g/L	2	RS01	off
total acidity at pH 8,2	= 'DET pH.FP{2}.VOL' * 'DET pH.CONC' * 'DET pH.TITER' * 150.1 / 'MV.Sample size'	g/L	2	RS02	off

Result name **total acidity at pH 7**
 Formula = 'DET pH.FP{1}.VOL' * 'DET pH.TITER' * 'DET pH.CONC' * 150.1 / 'MV.Sample size'
 Unit g/L
 Decimal places 2
 Assignment RS01
 Statistics off
 Description RS.'Result name'[.VAL]Result value.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer off
 Solution name

Result name **total acidity at pH 8,2**
 Formula = 'DET pH.FP{2}.VOL' * 'DET pH.CONC' * 'DET pH.TITER' * 150.1 / 'MV.Sample size'
 Unit g/L
 Decimal places 2
 Assignment RS02
 Statistics off
 Description RS.'Result name'[.VAL]Result value.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer off
 Solution name

REPORT report

Report template
 Report template report acidity of wine



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:14:07 UTC+2

Report output

Printer off
 PDF file on
 PDF file C:\Programme\Metrohm\tiamo\Data\total acidity of wine.pdf
 Send e-mail off

DATABASE database

Database

Robotic Acid Base Analyzer

**SERIES Series end track
 END**

MOVE to storing beaker

Device

Device name 855_1

Target

Tower 1
 Move Special beaker
 Number 1

Beaker test

Display message off
 Stop determination on
 Stop determination and series off

Parameters

Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

ERROR Error track

MOVE to storage beaker

Device

Device name 855_1

Target

Tower 1
 Move Special beaker
 Number 1

Beaker test

Display message off
 Stop determination on
 Stop determination and series off

Parameters

Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

5.6 Titration of citric acid



License ID 124049905 Program version tiamo 1.1 - 31
 Client name TITRATION14
 User Metrohm 2005-09-13 10:16:25 UTC+2

Method parameters

Method Titration of citric acid
 Method saving date 2005-09-13 10:16:21 UTC+2
 Method version 1
 Method group Robotic Acid Base Analyzer
 Method status original
 Method saved by (full name) Metrohm
 Method saved by (short name) Metrohm

START

Main track

General

Workplace view
 Current view on
 Track view for live window
 Live display 1 Main track
 Live display 2 Main track
 Statistics off
 Conditioning
 Automatic conditioning off

Application note

See attached documents

Method variables

Name	Type	Assignment	Value	Comment	Monitoring
Sample size	Number	Sample size		Sample size	off
Sample size unit	Text	Sample size unit		Sample size unit	off
Sample position	Number	Sample position		Sample position number	off
ID1	Text	ID1		Sample identification 1	off
ID2	Text	ID2		Sample identification 2	off
ID3	Text	ID3		Sample identification 3	off

Name **Sample size**
 Type Number
 Assignment on Sample size
 Value off
 Check at start on
 Comment Sample size
 Variable monitoring off
 Lower limit
 Upper limit
 Message
 Display message on
 Record message on
 Message by e-mail off
 Mail to



License ID 124049905
Client name TITRATION14
User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:16:25 UTC+2

Subject Message from tiamo - Method 'New method 6' - Command 'Main track'

User

Mail from

SMTP Server

POP3 Server

Acoustic signal off

Action off

Stop determination on

Stop determination and series off

Name **Sample position**

Type Number

Assignment on. Sample position

Value off.

Check at start on

Comment Sample position number

Variable monitoring off

Lower limit

Upper limit

Message

Display message on

Record message on

Message by e-mail off

Mail to

Subject Message from tiamo - Method 'New method 6' - Command 'Main track'

User

Mail from

SMTP Server

POP3 Server

Acoustic signal off

Action off

Stop determination on

Stop determination and series off

Name **Sample size unit**

Type Text

Assignment on. Sample size unit

Value off.

Check at start on

Comment Sample size unit

Name **ID1**

Type Text

Assignment on. ID1

Value off.



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:16:25 UTC+2

Check at start on
 Comment Sample identification 1

Name ID2
 Type Text
 Assignment on. ID2
 Value off.
 Check at start on
 Comment Sample identification 2

Name ID3
 Type Text
 Assignment on. ID3
 Value off.
 Check at start on
 Comment Sample identification 3

MOVE to sample position

Device
 Device name 855_1
 Target
 Tower 1
 Move Sample position
 Beaker test
 Display message off
 Stop determination off
 Stop determination and series on
 Parameters
 Shift rate 20 °/s
 Shift direction auto
 Swing rate 55 °/s

PUMP add water

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 1
 Action
 Switch on off
 Switch off off
 Duration on
 Time 10.0 s

STIR stirrer on

Device
 Device name 855_1



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31
 2005-09-13 10:16:25 UTC+2

Stirrer

Stirrer 1
 Stirrer type unknown
 Stirring rate 8

Action

Switch on on
 Switch off off
 Duration off

WAIT

dissolving sample

Wait

Stop track and waiting for [Continue] off
 Stop all tracks and waiting for [Continue] off
 Waiting time on
 Time 60
 Unit s

Message

Record message on
 Message by e-mail off
 Acoustic signal off

DET pH

DET pH

General/Hardware

Device

Device name 855_1

Dosing device

Dosing device 1
 Solution NaOH

Sensor

Measuring input 1
 Sensor Aquatrode Plus
 Temperature measurement automatic

Stirrer

Stirrer off

Start conditions

Initial measured value

Signal drift off mV/min
 Min. waiting time 0 s
 Max. waiting time 1 s

Start volume

Start volume 0 mL
 Dosing rate maximum mL/min

Start measured value

Start measured value pH off
 Dosing rate 5 mL/min

Start slope

Start slope off pH/mL



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:16:25 UTC+2

Dosing rate 5 mL/min
 Pause
 Pause 0 s

Titration parameters

Titration rate
 Titration rate user
 Measured value acceptance
 Signal drift 50.0 mV/min
 Min. waiting time 0 s
 Max. waiting time 26 s
 Dosing of increments
 Measuring point density 4
 Min. increment 100 µL
 Max. increment off µL
 Dosing rate maximum mL/min
 Temperature
 Temperature 25.0 °C

Stop conditions

Stop volume 20 mL
 Stop measured value pH off
 Stop EP 1
 Volume after EP 1 mL
 Stop time off s
 Filling rate maximum mL/min

Potentiometric evaluation

Evaluation without window on
 EP criterion 20
 EP recognition all
 Evaluation with measured value window (pH) off
 Evaluation with volume window (mL) off

Additional evaluations

Fix end point evaluation off
 pK/HNP evaluation off
 Minimum evaluation off
 Maximum evaluation off
 Break point evaluation off

Additional measured values

Additional calculated measured values off
 Additional external measured values off

STIR

stirrer off

Device
 Device name 855_1
 Stirrer
 Stirrer 1
 Stirrer type unknown
 Stirring rate 8



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:16:25 UTC+2

Action
 Switch on off
 Switch off on
 Duration off

PUMP aspirate sample solution

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 5 s

PUMP rinse and aspirate

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 1+2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 5 s

PUMP aspirate

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 3 s

SERIES START Series start track

RACK initialize rack

Device
 Device name 855_1
 Rack test off

SERIES END Series end track



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:16:25 UTC+2

MOVE to storing beaker

Device
 Device name 855_1
 Target
 Tower 1
 Move Special beaker
 Number 1
 Beaker test
 Display message off
 Stop determination on
 Stop determination and series off
 Parameters
 Shift rate 20 °/s
 Shift direction auto
 Swing rate 55 °/s

EXIT Exit track

CALC calculation

Result name	Formula	Unit	Decimal places	Assignment	Statistics
percent by weight	= 'DET pH.EP{1}.VOL' * 192.12 * 'DET pH.CONC' * 'DET pH.TITER' * 100 / ('MV.Sample size' * 3 * 1000)	%	2	RS01	off

Result name **percent by weight**
 Formula = 'DET pH.EP{1}.VOL' * 192.12 * 'DET pH.CONC' * 'DET pH.TITER' * 100 / ('MV.Sample size' * 3 * 1000)
 Unit %
 Decimal places 2
 Assignment RS01
 Statistics off
 Description RS.'Result name'[.VAL]Result value.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer off
 Solution name

REPORT report

Report template
 Report template citric acid
 Report output
 Printer off
 PDF file on
 PDF file C:\Programme\Metrohm\tiamo\Data\titration of citric acid.pdf



License ID 124049905
Client name TITRATION14
User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:16:25 UTC+2

Send e-mail off

DATABASE database

Database

Robotic Acid Base Analyzer

ERROR Error track

MOVE to storage beaker

Device

Device name 855_1

Target

Tower 1

Move Special beaker

Number 1

Beaker test

Display message off

Stop determination on

Stop determination and series off

Parameters

Shift rate 20 %/s

Shift direction auto

Swing rate 55 %/s

5.7 Determination of Phosphoric acid in cola drinks



License ID 124049905 Program version tiamo 1.1 - 31
 Client name TITRATION14
 User Metrohm 2005-09-13 10:20:41 UTC+2

Method parameters

Method Determination of phosphoric acid
 in cola drinks
 Method saving date 2005-09-13 10:20:36 UTC+2
 Method version 1
 Method group Robotic Acid Base Analyzer
 Method status original
 Method saved by (full name) Metrohm
 Method saved by (short name) Metrohm

START

Main track

General

Workplace view
 Current view on
 Track view for live window
 Live display 1 Main track
 Live display 2 Main track
 Statistics off
 Conditioning
 Automatic conditioning off

Application note

See attached documents

Method variables

Name	Type	Assignment	Value	Comment	Monitoring
Sample size	Number	Sample size		Sample size	off
Sample size unit	Text	Sample size unit		Sample size unit	off
Sample position	Number	Sample position		Sample position number	off
ID1	Text	ID1		Sample identification 1	off
ID2	Text	ID2		Sample identification 2	off
ID3	Text	ID3		Sample identification 3	off

Name **Sample size**
 Type Number
 Assignment on Sample size
 Value off
 Check at start on
 Comment Sample size
 Variable monitoring off
 Lower limit
 Upper limit
 Message
 Display message on
 Record message on
 Message by e-mail off



License ID 124049905
Client name TITRATION14
User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:20:41 UTC+2

Mail to
Subject Message from tiamo - Method 'New method 3' - Command 'Main track'
User
Mail from
SMTP Server
POP3 Server
Acoustic signal off
Action off
Stop determination on
Stop determination and series off

Name **Sample position**
Type Number
Assignment on, Sample position
Value off,
Check at start on
Comment Sample position number
Variable monitoring off
Lower limit
Upper limit
Message
Display message on
Record message on
Message by e-mail off
Mail to
Subject Message from tiamo - Method 'New method 3' - Command 'Main track'
User
Mail from
SMTP Server
POP3 Server
Acoustic signal off
Action off
Stop determination on
Stop determination and series off

Name **Sample size unit**
Type Text
Assignment on, Sample size unit
Value off,
Check at start on
Comment Sample size unit

Name **ID1**
Type Text
Assignment on, ID1



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:20:41 UTC+2

Value off
 Check at start on
 Comment Sample identification 1

Name **ID2**
 Type Text
 Assignment on ID2
 Value off
 Check at start on
 Comment Sample identification 2

Name **ID3**
 Type Text
 Assignment on ID3
 Value off
 Check at start on
 Comment Sample identification 3

MOVE to sample position

Device
 Device name 855_1
 Target
 Tower 1
 Move Sample position
 Beaker test
 Display message off
 Stop determination off
 Stop determination and series on
 Parameters
 Shift rate 20 %/s
 Shift direction auto
 Swing rate 55 %/s

DET pH DET pH

General/Hardware

Device
 Device name 855_1
 Dosing device
 Dosing device 1
 Solution NaOH
 Sensor
 Measuring input 1
 Sensor Aquatrode Plus
 Temperature measurement automatic
 Stirrer
 Stirrer 1
 Stirring rate 8



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31
 2005-09-13 10:20:41 UTC+2

Switch off automatically on

Start conditions

Initial measured value
 Signal drift off mV/min
 Min. waiting time 0 s
 Max. waiting time 1 s
 Start volume
 Start volume 0 mL
 Dosing rate maximum mL/min
 Start measured value
 Start measured value pH off
 Dosing rate 5 mL/min
 Start slope
 Start slope off pH/mL
 Dosing rate 5 mL/min
 Pause
 Pause 0 s

Titration parameters

Titration rate
 Titration rate user
 Measured value acceptance
 Signal drift 50.0 mV/min
 Min. waiting time 0 s
 Max. waiting time 26 s
 Dosing of increments
 Measuring point density 4
 Min. increment 100 µL
 Max. increment off µL
 Dosing rate maximum mL/min
 Temperature
 Temperature 25.0 °C

Stop conditions

Stop volume 20 mL
 Stop measured value pH off
 Stop EP 1
 Volume after EP 1 mL
 Stop time off s
 Filling rate maximum mL/min

Potentiometric evaluation

Evaluation without window on
 EP criterion 5
 EP recognition all
 Evaluation with measured value window (pH) off
 Evaluation with volume window (mL) off

Additional evaluations

Fix end point evaluation off



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:20:41 UTC+2

pK/HNP evaluation off
 Minimum evaluation off
 Maximum evaluation off
 Break point evaluation off

Additional measured values

Additional calculated measured values off
 Additional external measured values off

PUMP aspirate sample solution

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 7 s

PUMP rinse and aspirate

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 1+2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 5 s

PUMP aspirate

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 3.0 s

SERIES START Series start track

RACK initialize rack

Device



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31
 2005-09-13 10:20:41 UTC+2

Device name 855_1
 Rack test off

EXIT Exit track

CALC calculation

Result name	Formula	Unit	Decimal places	Assignment	Statistics
content phosphoric acid	'DET pH.EP{1}.VOL' * 97.993 * 'DET pH.TITER' * 'DET pH.CONC' / 'MV.Sample size'	g/L	3	RS01	off

Result name **content phosphoric acid**
 Formula ='DET pH.EP{1}.VOL' * 97.993 * 'DET pH.TITER' * 'DET pH.CONC' / 'MV.Sample size'
 Unit g/L
 Decimal places 3
 Assignment RS01
 Statistics off
 Description RS.'Result name'[_VAL]Result value.
 Result monitoring off
 Save result as common variable off
 Name
 Save result as titer off
 Solution name

REPORT report

Report template
 Report template phosphoric acid in cola drinks
 Report output
 Printer off
 PDF file on
 PDF file C:\Programme\Metrohm\tiamo\Data\Determination of phosphoric acid in cola drinks.pdf
 Send e-mail off

DATABASE database

Database
 Robotic Acid Base Analyzer

SERIES END Series end track

MOVE to storing beaker

Device
 Device name 855_1
 Target



License ID 124049905
 Client name TITRATION14
 User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:20:41 UTC+2

Tower 1
 Move Special beaker
 Number 1
 Beaker test
 Display message off
 Stop determination off
 Stop determination and series on
 Parameters
 Shift rate 20 °/s
 Shift direction auto
 Swing rate 55 °/s

ERROR Error track

MOVE to rinsing beaker

Device
 Device name 855_1
 Target
 Tower 1
 Move Special beaker
 Number 2
 Beaker test
 Display message off
 Stop determination off
 Stop determination and series on
 Parameters
 Shift rate 20 °/s
 Shift direction auto
 Swing rate 55 °/s

PUMP aspirate and rinse

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 1+2
 Action
 Switch on off
 Switch off off
 Duration on
 Time 10 s

PUMP empty rinsing beaker

Device
 Device name 855_1
 Pumps
 Tower 1
 Pump(s) 2



License ID 124049905
Client name TITRATION14
User Metrohm

Program version tiamo 1.1 - 31

2005-09-13 10:20:41 UTC+2

Action

Switch on off
Switch off off
Duration on
Time 10 s

MOVE to storage beaker

Device

Device name 855_1

Target

Tower 1
Move Special beaker
Number 1

Beaker test

Display message off
Stop determination off
Stop determination and series on

Parameters

Shift rate 20 %/s
Shift direction auto
Swing rate 55 %/s