

872 Extension Module



IC Module – 2.872.0030

Manual
8.872.8004EN



Metrohm AG
CH-9100 Herisau
Switzerland
Phone +41 71 353 85 85
Fax +41 71 353 89 01
info@metrohm.com
www.metrohm.com

872 Extension Module

IC Module – 2.872.0030

Manual

Teachware
Metrohm AG
CH-9100 Herisau
teachware@metrohm.com

This documentation is protected by copyright. All rights reserved.

Although all the information given in this documentation has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please send us your comments using the address given above.

Documentation in additional languages can be found on <http://products.metrohm.com> under **Literature/Technical documentation**.

Table of contents

1	Introduction	1
1.1	Instrument description	1
1.2	Intended use	2
1.3	About the documentation	2
1.3.1	Content and scope	2
1.3.2	Symbols and conventions	3
1.4	Safety instructions	4
1.4.1	General notes on safety	4
1.4.2	Electrical safety	4
1.4.3	Tubing and capillary connections	5
1.4.4	Flammable solvents and chemicals	5
1.4.5	Recycling and disposal	5
2	Overview of the instrument	7
2.1	Front	7
2.2	Rear	8
3	Assembly	9
3.1	General	9
3.2	Mounting the extension module onto the IC instrument	10
3.3	Mounting the extension module below the IC instrument	13
3.4	Setting up the extension module next to the IC instrument	16
3.5	Transport locking screws	20
4	Installation	21
4.1	About this chapter	21
4.2	Installation overview	21
4.3	Installation diagram	22
4.4	Eluent	23
4.4.1	Connecting eluent bottle	23
4.5	Eluent degasser	27
4.6	High pressure pump	29
4.6.1	Capillary connections high pressure pump/purge valve	29
4.6.2	Deaerating the high pressure pump	31



4.7	Inline filter	33
4.8	Pulsation damper	34
4.9	Injection valve	35
4.9.1	Connecting the injection valve	35
4.9.2	Mode of operation of the injection valve	37
4.9.3	Selecting the sample loop	38
5	Set to work	39
6	Operation and maintenance	40
6.1	General information	40
6.1.1	Care	40
6.1.2	Maintenance by Metrohm Service	40
6.1.3	Operation	40
6.1.4	Shutting down	41
6.2	Door	41
6.3	Eluent	41
6.3.1	Production	41
6.3.2	Operation	42
6.4	High pressure pump	42
6.4.1	Protection	42
6.4.2	Maintenance	43
6.5	Inline filter	53
6.5.1	Maintenance	53
6.6	Inline sample preparation	55
6.7	Rinsing the sample path	55
6.8	Injection valve	57
6.8.1	Protection	57
6.9	Quality Management and validation with Metrohm	57
7	Troubleshooting	58
7.1	Problems and their solutions	58
8	Technical specifications	60
8.1	Reference conditions	60
8.2	Instrument	60
8.3	Ambient conditions	60
8.4	Housing	61
8.5	Eluent degasser	61
8.6	High pressure pump	61
8.7	Injection valve	62



8.8	Interfaces	62
8.9	Safety specification	62
8.10	Electromagnetic compatibility (EMC)	63
8.11	Weight	63
9	Conformity and warranty	64
9.1	Declaration of Conformity	64
9.2	Quality Management Principles	65
9.3	Warranty (Guarantee)	66
10	Accessories	68
10.1	Scope of delivery	68
	Index	75



Table of figures

Figure 1	Front 872 Extension Module IC Module	7
Figure 2	Rear 872 Extension Module IC Module	8
Figure 3	Setup versions	10
Figure 4	Dismounting the bottle holder	11
Figure 5	Mounting the bottle holder	12
Figure 6	Removing the base tray	14
Figure 7	Mounting the base tray	15
Figure 8	Mounting the base tray	17
Figure 9	Mounting the bottle holder	18
Figure 10	Connecting the drainage tubings	19
Figure 11	Installation diagram	23
Figure 12	Installing eluent bottle attachment	24
Figure 13	Mounting aspiration filter	24
Figure 14	Install tubing weighting and aspiration filter	25
Figure 15	Eluent aspiration tubing fully equipped.	25
Figure 16	Eluent bottle – connected	26
Figure 17	Eluent degasser	28
Figure 18	Capillary connections high pressure pump/purge valve	29
Figure 19	High pressure pump – Connect inlet	30
Figure 20	Deaerating the high pressure pump	32
Figure 21	Connecting inline filter	34
Figure 22	Pulsation damper – Connection	35
Figure 23	Injection valve – connected	36
Figure 24	Injection valve – Positions	37
Figure 25	Removing piston	44
Figure 26	Components of the piston cartridge	45
Figure 27	Tool for piston seal 6.2617.010	46
Figure 28	Removing the piston seal	47
Figure 29	Insert the piston seal into the tool	47
Figure 30	Inserting the piston seal into the pump head	48
Figure 31	Removing valves	49
Figure 32	Dismantling valve	50
Figure 33	Components of the inlet valve and outlet valve	51
Figure 34	Changing the filter	53

1 Introduction

1.1 Instrument description

Existing 850 Professional IC instruments can be expanded to include additional functions by means of 872 Extension Modules. Every 850 Professional IC instrument can be supplemented with up to 3 extension modules.

The **872 Extension Module IC Module** allows to extend an 850 Professional IC instrument by another IC Module complete with high pressure pump, purge valve, inline filter, pulsation damper and injection valve.

The **872 Extension Module IC Module** can be used wherever an additional IC Module is necessary.

With the 872 Extension Module IC Module, a one-channel low pressure gradient instrument can this way be extended to a two-channel AnCat instrument with low pressure gradient. All instruments with sample preparation can be extended to two-channel systems as well.

In a system with photometric detection with a 886 Professional Thermostat/Reactor as column holder and thermostat and an 887 Professional UV/VIS Detector, it is also possible to use the 872 Extension Module IC Module for conveying reagents.

The extension module is operated with **MagIC Net** software, just like the IC instrument. When it is connected to an 850 Professional IC instrument, **MagIC Net** recognizes the extension module automatically and checks its functional capability. It controls and monitors the unit IC instrument – extension module, evaluates the measured data and administers it in a database.

The **872 Extension Module – IC Module** comprises the following components:

Eluent degasser

The eluent degasser removes gas bubbles and dissolved gases from the eluent. For degassing, the eluent flows into a vacuum chamber through a special fluoropolymer capillary.

High pressure pump

The intelligent and low pulsation high pressure pump pumps the eluent through the system. It is equipped with a chip on which its technical specifications and "life history" (operating hours, service data, ...) are saved.



Inline filter

Inline filters protect the separation column securely against possible contamination from the eluent. Inline filters can however also just as well be used for the purpose of protecting other sensitive components against contaminations in the solutions used. The fine 2 µm material of the readily and easily replaceable filter platelets removes particles such as bacteria and algae from the solutions.

Pulsation damper

The pulsation damper protects the separation column from damage caused by pressure fluctuations when switching the injection valve, and reduces interfering pulsations during highly sensitive measurements.

Injection valve

The injection valve connects the eluent and sample path through rapid and precise valve switchover. A precisely measured amount of sample solution is injected and rinsed with eluent onto the separation column.

1.2 Intended use

The present instrument is suitable for processing chemicals and flammable samples. The usage of the 872 Extension Module therefore requires that the user has basic knowledge and experience in the handling of toxic and caustic substances. Knowledge with respect to the application of the fire prevention measures prescribed for laboratories is also mandatory.

1.3 About the documentation

1.3.1 Content and scope







This document describes the **872 Extension Module – IC Module**, its assembly and connection to the IC instrument, as well as the installation, operation and maintenance of the individual components. Technical specifications, troubleshooting and information concerning scope of delivery and optional accessories makes up the rest of the manual.

This document does not on the other hand describe the functions of the IC instrument - IC extension module unit, nor does it describe the capillary connections that proceed from the extension module. For this purpose, please refer to the manual for the IC instrument and that for the sample processor.

Additional information concerning the configuration of MagIC Net can be found on the online help for MagIC Net.

1.3.2 Symbols and conventions

The following symbols and styles are used in this documentation:

(5-12)	<p>Cross-reference to figure legend</p> <p>The first number refers to the figure number, the second to the instrument part in the figure.</p>
1	<p>Instruction step</p> <p>Carry out these steps in the sequence shown.</p>
	<p>Warning</p> <p>This symbol draws attention to a possible life hazard or risk of injury.</p>
	<p>Warning</p> <p>This symbol draws attention to a possible hazard due to electrical current.</p>
	<p>Warning</p> <p>This symbol draws attention to a possible hazard due to heat or hot instrument parts.</p>
	<p>Warning</p> <p>This symbol draws attention to a possible biological hazard.</p>
	<p>Caution</p> <p>This symbol draws attention to a possible damage of instruments or instrument parts.</p>
	<p>Note</p> <p>This symbol marks additional information and tips.</p>

Protection against electrostatic charges



Warning

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Always pull the mains cable out of the mains connection socket before connecting or disconnecting electrical appliances on the rear panel of the instrument.

1.4.3 Tubing and capillary connections



Caution

Leaks in tubing and capillary connections are a safety risk. Tighten all connections well by hand. Avoid applying excessive force to tubing connections. Damaged tubing ends lead to leakage. Appropriate tools can be used to loosen connections.

Check the connections regularly for leakage. If the instrument is used mainly in unattended operation, then weekly inspections are mandatory.

1.4.4 Flammable solvents and chemicals




Warning

All relevant safety measures are to be observed when working with flammable solvents and chemicals.

- Set up the instrument in a well-ventilated location (e.g. laboratory flue).
- Keep all sources of flame far from the workplace.
- Clean up spilled fluids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.

1.4.5 Recycling and disposal



This product is covered by European Directive 2002/96/EC, WEEE – Waste from Electrical and Electronic Equipment.

The correct disposal of your old equipment will help to prevent negative effects on the environment and public health.



More details about the disposal of your old equipment can be obtained from your local authorities, from waste disposal companies or from your local dealer.

2 Overview of the instrument

2.1 Front




Figure 1 Front 872 Extension Module IC Module

1	Standby indicator	2	High pressure pump <i>See Chapter 4.6</i>
3	Coupling 6.2744.230 For connecting the eluent aspiration tubing.	4	Purge valve <i>See Chapter 4.6</i>
5	Pulsation damper <i>See Chapter 4.8</i>	6	Injection valve <i>See Chapter 4.9</i>
7	Inline filter <i>See Chapter 4.7</i>	8	Eluent degasser <i>See Chapter 4.5</i>



2.2 Rear




Figure 2 Rear 872 Extension Module IC Module

1 Connector In

To connect the extension module to the IC instrument or to a previous extension module.

2 Connector Out

To connect an additional extension module.

3 Connector Vacuum

For connecting further degassing chambers in extension modules (labeled with *Vacuum*).

4 Waste air opening

For extracting the air from the vacuum chamber (labeled with *Exhaust*).

5 Transport locking screws

For securing the high pressure pump and the vacuum pump when transporting the instrument.

3 Assembly

3.1 General

The extension modules are fitted directly to the 850 Professional IC instrument and connected with it via 6.2156.060 connection cable. Extension modules have no power supply of their own, but rather draw the electricity they require from the instrument with which they are connected.

Up to three extension modules can be connected to an 850 Professional IC instrument. The following restrictions are to be taken into account:

Restrictions

The 850 Professional IC instruments and their extension modules must not have more than 4 identical components in common, i.e.:

- a maximum of 4 high pressure pumps,
- a maximum of 4 peristaltic pumps,
- a maximum of 4 injection valves,
- a maximum of 4 suppressors (MSM, SPM incl.),

BUT

- only a maximum of 3 degassers
- and a maximum of 3 CO₂ suppressors (MCS)



Note

If all 4 high pressure pumps are being used at once, then not all of them are permitted to run at maximum flow for longer periods of time.

Extension modules can be mounted in the following setup versions:

- above, between instrument and bottle holder (**3-A**), or
- below, between instrument and base tray (**3-B**), or
- next to the instrument (**3-C**) with a separate 6.2061.110 base tray and a 6.2061.100 bottle holder (to be ordered additionally). For this, the longer 6.2156.070 connection cable (to be ordered additionally), is necessary, too.




Figure 3 Setup versions

A Extension module above the IC instrument
Between the 850 Professional IC and the bottle holder.

B Extension module on the below the IC instrument
Between the base plate and the 850 Profession IC.

C Extension module next to the IC instrument
With its own base plate and its own bottle holder next to the 850 bottle holder.

Position the extension module in such a way that the capillary connections can be kept as short as possible. If several extension modules are used, then they should all be installed in the same location if possible – either above, below or next to the IC instrument. If this is not possible, then the extension modules that are located at greater distances from one another must be connected with one another by means of the longer 6.2156.070 connection cable (available as optional accessory).

3.2 Mounting the extension module onto the IC instrument

1 Switching off the IC instrument

Switch the IC instrument off and disconnect the mains cable.

2 Clearing the bottle holder

If there are bottles and other items on the bottle holder, remove them.

3 Removing drainage tubings

Loosen the drainage tubing from the drainage tubing connector on the bottle holder.

4 Dismounting the bottle holder




Figure 4 Dismounting the bottle holder

1 Cover stoppers

2 Cylinder screws

- Remove covering stoppers (4-1).
- Loosen the cylinder screws with a 6.2621.100 3 mm hexagon key.
- Remove the bottle holder

5 Attaching the extension module(s)

Place the extension module(s) on the IC instrument.



6 Mounting the bottle holder




Figure 5 Mounting the bottle holder

1 Cover stoppers

2 Cylinder screws

- Attach the bottle holder on the extension module.
- Tighten the cylinder screws (4-**2**) with an 6.2621.100 3 mm hexagon key
- Insert covering stoppers (4-**1**).

7 Connecting the extension module

- Plug a 6.2156.060 cable into the connector **In** of the extension module and screw it tight.
- Plug the other end of the cable into the **Extension module** connector of the IC instrument and screw it tight.

8 Optional: Connecting a further extension module

- Plug the cable 6.2156.060 or a longer 6.2156.070 cable (optional accessory) into the connector **In** of the second extension module and screw it tight.
- Plug the other end of the cable into the **Out** connector of the first extension module and screw it tight.

9 Mounting the drainage tubing

Reconnect the drainage tubing to the drainage tubing connector of the bottle holder.

Possibly, a longer section of silicone tubing 6.186.020 must be cut to fit and mounted (*see also the manual for the IC instrument*).

3.3 Mounting the extension module below the IC instrument

1 Switching off the IC instrument

Switch the IC instrument off and disconnect the mains cable.

2 Clearing the bottle holder

If there are bottles and other things on the bottle holder, remove them.

3 Disconnecting all connections on the rear of the instrument

- Disconnect the mains cable,
- Disconnect the MSB cable,
- Disconnect the USB cable,
- Disconnect the leak sensor,
- Remove the drainage tubings.

4 Removing the detector(s)

Disconnect the detector cable(s) and remove the detector(s) from the IC instrument (*see manual for the IC instrument*).

5 Removing the base tray

- Tilt the IC instrument sideways and lay it down flat.
- Loosen the cylinder screws with a 6.2621.100 3 mm hexagon key.
- Remove base tray.




Figure 6 Removing the base tray

1 Base tray

2 Cylinder screws
With washer.

6 Mounting the base tray

- Tilt the extension module sideways and lay it down flat.
- Attach base tray.
- Slide the washers onto the cylinder screws (7-2) and tighten these with a 3 mm hexagon key (6.2621.100).




Figure 7 Mounting the base tray

1 Base tray

2 Cylinder screws

With washer.

- Set up the extension module.
- Optional: Set up further extension modules.
- Set up the IC instrument.

7 Connecting the extension module

- Plug a 6.2156.060 cable into the connector **In** of the extension module and screw it tight.
- Plug the other end of the cable into the **Extension module** connector of the IC instrument and screw it tight.

8 Optional: Connecting a further extension module

- Plug the cable 6.2156.060 or a longer 6.2156.070 cable (optional accessory) into the connector **In** of the second extension module and screw it tight.
- Plug the other end of the cable into the **Out** connector of the first extension module and screw it tight.

**9 Inserting the detector(s) again and connecting it (them)**

See manual for the IC instrument.

10 Restoring the loosened connections

- Reconnect the drainage tubings,
Possibly, a longer section of silicone tubing 6.186.020 must be cut to fit and mounted (*see also the manual for the IC instrument*).
- Connect the leak sensor (*see manual for the IC instrument*),
- Connect the USB cable,
- Connect the MSB cable,
- Plug in the mains cable.

3.4 Setting up the extension module next to the IC instrument

1 Switching off the IC instrument

Switch the IC instrument off and disconnect the mains cable.

2 Mounting the base tray

- Tilt the extension module sideways and lay it down flat.
- Attach base tray.
- Tighten the washers and the cylinder screws (8-2) with a 3 mm hexagon key (6.2621.100).




Figure 8 Mounting the base tray

1 Base tray

2 Cylinder screws

With washer.

- Set up the extension module.
- Optional: Set up further extension modules.



3 Mounting the bottle holder




Figure 9 Mounting the bottle holder

1 Cover stoppers

2 Cylinder screws

- Place the bottle holder on the extension module.
- Tighten the cylinder screws (9-2) with an 6.2621.100 3 mm hexagon key.
- Insert covering stoppers (9-1).

4 Connecting the extension module

- Plug a 6.2156.060 cable into the connector **In** of the extension module and screw it tight.
- Plug the other end of the cable into the **Extension module** connector of the IC instrument and screw it tight.

5 Optional: Connecting a further extension module

- Plug the cable 6.2156.060 or a longer 6.2156.070 cable (optional accessory) into the connector **In** of the second extension module and screw it tight.
- Plug the other end of the cable into the **Out** connector of the first extension module and screw it tight.

6 Connecting the leak sensor

- Plug the 6.2103.170 adapter into the leak sensor connector of the IC instrument.
- Connect the leak sensor cable of the IC instrument to the adapter.
- Connect the leak sensor cable of the extension module to the adapter.

7 Connecting the drainage tubings




Figure 10 Connecting the drainage tubings

1 Drainage tubing connector

For draining escaped liquid from the bottle holder.

2 Drainage tubing

Section of the 6.1816.020 silicon tubing. For draining escaped liquid from the bottle holder.

**3 Drainage tubing connection**

For supplying escaped liquid through the connected drainage tubing to the leak sensor.

4 Drainage tubing connector

For draining escaped liquid from the base tray through the connected drainage tubing.

5 Drainage tubing

Section of the 6.1816.020 silicon tubing. Guides escaped liquid into a waste container.

- Plug the drainage tubing (10-2) onto the drainage tubing connector (10-1) of the bottle holder and shorten to required length.
- Plug the other end of the drainage tubing (10-2) onto the drainage tubing connector (10-3) of the base tray.
- Plug the drainage tubing (10-5) onto the drainage tubing connector (10-4) and guide the other end into a waste container.

3.5 Transport locking screws

To avoid damage to the high pressure pump and vacuum pump during transport, the pumps are secured with transport locking screws .

Remove these transport locking screws before the initial start-up.

Removing transport locking screws

- 1 Remove all of the transport locking screws with the 6.2621.030 4 mm hexagon key and keep them in a safe place.



Warning

In order to avoid damage to the pump, the transport locking screws must be remounted each time the instrument undergoes major transport.

4 Installation

4.1 About this chapter

The *Installation* chapter contains

- this overview
- a brief set of instructions for the installation of the 872 Extension Module IC Module (see *Chapter 4.2, page 21*). At each step you will find cross-references to more detailed installation instructions for individual components, should you require such aids.
- an installation diagram (see *Chapter 4.3, page 22*), showing a completely installed instrument.
- several chapters with detailed installation instructions for all components, including those that are already installed at the time the instrument is delivered.

4.2 Installation overview



Note

A number of the capillary connections are already connected at the time the instrument is delivered.

The following work steps must still be carried out:

Installing 872 Extension Module IC Module

1 Installing the eluent path

- Tighten the clamping screw of the 6.1834.080 eluent aspiration tubing to the eluent degasser input.
Equip the free end of the eluent aspiration tubing (**11-1**) and connect it to the eluent bottle (see *Chapter 4.4, page 23*).
- Tighten one end of the 6.1831.010 PEEK capillary with a 6.2744.014 PEEK pressure screw to the connector 4 of the injection valve (see *Chapter 4.9.1, page 35*).
Shorten to required length with the capillary cutter and with the aid of a 6.2744.040 coupling and two 6.2744.010 pressure screws connect to the column outlet capillary in the IC instrument.

**Note**

The separation column must not be connected to the system before having it started up for the first time (see manual for the IC instrument).

2 Installing the sample path

- Connect one end of the 6.1803.040 PTFE capillary (11-3) with a 6.2744.014 PEEK pressure screw to connector 1 of the injection valve (see Chapter 4.9.1, page 35).
Guide the other end through a suitable capillary feed-through out of the extension module to the Sample Processor and connect it there (see Sample Processor manual).
- Connect one end of the second 6.1803.040 PTFE capillary (11-4) with a 6.2744.014 PEEK pressure screw to connector 2 of the injection valve.
Guide the other end through a capillary feed-through out of the extension module to a waste container and fasten it there.

3 Putting the extension module into operation

See Chapter 5, Page 39

- Putting the extension module into operation together with the IC instrument (see chapter *Start-up* in the manual for the IC instrument).
- Deaerate the high pressure pump at the same time (see Chapter 4.6.2, page 31).

4.3 Installation diagram

Figure 11 Installation diagram shows the capillary connections of the 872 Extension Module IC Module as an additional IC Module which can be integrated into any existing system.

The arrangement of the modules in the diagram corresponds to the front view of the extension module.

A few of the capillaries are already pre-installed at the time the instrument is delivered. Capillaries on which nothing needs to be done at the time of initial installation are not numbered in the diagram.




Figure 11 Installation diagram

<p>1 Eluent aspiration tubing 6.1834.080 Connected to the eluent degasser. Connect the other end to the eluent bottle.</p>	<p>2 Eluent connection capillary 6.1831.010</p>
<p>3 Sample aspirating capillary 6.1803.040</p>	<p>4 Sample outlet capillary 6.1803.040</p>

The following chapters describe the individual installation steps in detail.

4.4 Eluent

4.4.1 Connecting eluent bottle

The eluent is aspirated out of the eluent bottle via the eluent aspiration tubing (12-1).

The eluent aspiration tubing is connected to the eluent degasser (see Chapter 4.5, page 27). The tubing must be threaded through a suitable capillary feed-through of the instrument before the other end can be equipped.

You will require the parts from the following accessories for equipping the eluent aspiration tubing:

- 6.1602.160 Eluent bottle attachment GL 45
- 6.2744.210 tubing adapter for aspiration filter
- 6.2821.090 aspiration filter

To equip the eluent aspiration tubing proceed as follows:

Assembling eluent aspiration tubing

- 1** Guide the free end of the eluent aspiration tubing (12-1) out of the instrument through a suitable capillary feed-through.



2 Installing eluent bottle attachment 6.1602.160

- Slide tubing nipple (12-2) and O-ring (12-3) onto the eluent aspiration tubing (12-1).
- Push eluent aspiration tubing (12-1) through the bottle attachment (12-4) and screw tight.




Figure 12 Installing eluent bottle attachment

<p>1 Eluent aspiration tubing 6.1834.080</p>	<p>2 Tubing nipple From accessories set 6.1602.160.</p>
<p>3 O-ring From accessories set 6.1602.160.</p>	<p>4 Bottle attachment From accessories set 6.1602.160.</p>

3 Mounting aspiration filter

- Insert filter holder (13-1) into the aspiration filter (13-2) and screw tight.




Figure 13 Mounting aspiration filter

<p>1 Filter holder From accessories set 6.2744.210.</p>	<p>2 Aspiration filter 6.2821.090</p>
--	--




Figure 14 Install tubing weighting and aspiration filter

1 Eluent aspiration tubing 6.1834.080

2 Eluent bottle attachment 6.1602.160

3 Tubing weighting
From accessories set 6.2744.210.

4 Clamping screw
From accessories set 6.2744.210.

5 Aspiration filter 6.2821.090
With filter holder from accessories set 6.2744.210.

- Slide the tubing weighting (14-3) onto the eluent aspiration tubing (14-1).
- Slide the clamping screw (14-4) onto the eluent aspiration tubing (14-1).
- Insert the eluent aspiration tubing (14-1) into the aspiration filter (14-5). The end of the tubing should approximately reach to the center of the aspiration filter.
- Screw together clamping screw (14-4) and filter holder (13-1).




Figure 15 Eluent aspiration tubing fully equipped.

5 Mounting eluent aspiration tubing to the eluent bottle

- Insert the eluent aspiration tubing into the eluent bottle (16-10).



- Fasten the completely equipped bottle attachment (16-10) on the eluent bottle. The aspiration filter (16-6) must rest on the base of the eluent bottle.
- Close the remaining small opening on the bottle attachment with a threaded stopper from the accessories set.

6 Mounting the adsorber tube



Note

In the case of alkaline eluents and eluents with lower buffer capacity, the eluent bottle must be equipped with a CO₂ adsorber (16-4).

- First, place a piece of cotton (16-3), then the CO₂ adsorber (16-4) in the large opening of the adsorber tube (16-2) and close with the plastic cover.
- Fasten the adsorber tube (16-2) using the SGJ clip (16-12) onto the bottle attachment (16-11).




Figure 16 Eluent bottle – connected

1 Eluent aspiration tubing 6.1834.080
For aspirating the eluent. Pre-installed.

2 Adsorber tube 6.1609.000

3	Wadding	4	CO₂ adsorber Adsorbs CO ₂ from the air (e.g. Merck soda lime with indicator, no. 6839.10).
5	Eluent	6	Aspiration filter 6.2821.090
7	Filter holder From accessories set 6.2744.210.	8	Clamping screw From accessories set 6.2744.210.
9	Tubing weighting From accessories set 6.2744.210.	10	Eluent bottle 6.1608.070
11	Bottle attachment 6.1602.160	12	SGJ clip 6.2023.020
13	Tubing nipple	14	Threaded stopper

4.5 Eluent degasser

Gas bubbles in the eluent lead to an unstable baseline, as high pressure pumps can transport liquids, but not gases. The eluent therefore has to be degassed, before it reaches the high pressure pump.

The eluent degasser removes gas bubbles and dissolved gases from the eluent. For degassing, the eluent flows into a vacuum chamber through a special fluoropolymer capillary.



Note

The eluent degasser is already installed in the newly delivered instrument. The following installation instructions must only be followed, if the connections to the degasser had to be disconnected for maintenance.



Connecting the eluent degasser




Figure 17 Eluent degasser

1	Eluent degasser input	2	Eluent degasser output
3	Tubing flare With tubing nipple.	4	Clamping screw
5	Eluent aspiration tubing (6.1834.080) For aspirating the eluent. The clamping screw (17-4) is firmly mounted.	6	Connection tubing (6.1834.090) Connection from the eluent degasser to the high pressure pump (see Chapter 4.6, page 29). The clamping screw (17-4) is firmly mounted.

1



Caution

The clamping screws (17-4) must be tightened carefully. Use the wrench (6.2621.050) to do this.

- Insert the eluent aspiration tubing (17-5) into the eluent degasser input (17-1).
- Carefully tighten the clamping screw (17-4).

2

- Insert connection capillary (17-6) (the end with the longer clamping screw (17-4)) into the eluent degasser output (17-2).
- Carefully tighten the clamping screw (17-4).
- Connect the other end of the connection capillary (17-6) (with the shorter clamping screw) to the high pressure pump (18-9) (see "Connecting inlet to the high pressure pump", page 30).

4.6 High pressure pump

The intelligent and low pulsation high pressure pump pumps the eluent through the system. It is equipped with a chip on which its technical specifications and "life history" (operating hours, service data, ...) are saved.

The purge valve is used for deaerating (see Chapter 4.6.2, page 31) the high pressure pump.

4.6.1 Capillary connections high pressure pump/purge valve



Note

All of the capillary connections of the high pressure pump and the purge valve are already installed in the newly delivered instrument.




Figure 18 Capillary connections high pressure pump/purge valve

1 Connection capillary
PEEK capillary, connects main piston and auxiliary piston.

3 Outlet valve holder

2 PEEK pressure screw, short 6.2744.070

4 Pump head 6.2824.110

5 Fastening screws

For fastening the pump head.

7 Pump head input capillary

PEEK capillary at the input of the pump head.

9 Coupling

For the connection of the eluent path at the input of the high pressure pump. Can be ordered together with the pressure screw (18-8) under the number 6.2744.230.

11 Purge valve

For deaerating the high pressure pump. With rotary knob in the center and pressure sensor.

13 Connection capillary

Connects the output of the pump head with the purge valve.

6 Inlet valve holder**8 Pressure screw**

For connecting a PEEK capillary to the coupling (18-9).

10 Deaerating capillary

For aspirating the eluent when deaerating the high pressure pump (see Chapter 4.6.2, page 31).

12 Connection capillary

For connecting the inline filter (see Chapter 4.7, page 33)

**Note**

The eluent aspiration tubing is already installed in the newly delivered instrument. The following installation instructions need **not** be carried out at the time of initial installation.


Connecting inlet to the high pressure pump

Figure 19 High pressure pump – Connect inlet

1 Pressure screw

For connecting the coupling (19-2) to the pump head input capillary (18-7). Can be ordered together with the coupling under the number 6.2744.230.

2 Coupling 6.2744.230

For connecting the eluent aspiration tubing (19-4) to the input of the high pressure pump.

3 Clamping screw**4 Eluent aspiration tubing**

Eluent aspiration tubing 6.1834.080 or
6.1834.090.

5 Backup ring**1 Connecting coupling**

Fasten the coupling (19-2) with a pressure screw (19-1) on the pump head input capillary (18-7).

2 Connecting eluent aspiration tubing**Caution**

The clamping screws must be tightened carefully. To tighten, grip the coupling (19-2) with the 6.2739.000 key and grip the clamping screw (19-3) with the 6.2621.050 wrench.

- Plug the eluent aspiration tubing (19-4) into the coupling (19-2).
- Tighten clamping screw (19-3).

4.6.2 Deaerating the high pressure pump

The high pressure pump will only operate perfectly if the pump head contains no more air bubbles. Therefore it must be deaerated during initial start-up and after every change of eluent.

**Caution**

The high pressure pump must **not** be deaerated before the initial start-up.

Deaerate the high pressure pump as follows (see Figure 20, page 32):

Deaerating the high pressure pump

The instrument must be connected to the PC and switched on to deaerate the high pressure pump.




Figure 20 Deaerating the high pressure pump

1	Syringe 10 mL 6.2816.020 For aspirating the eluent.	2	Luer connector On purging needle.
3	Purging needle 6.2816.040	4	Deaerating capillary
5	PEEK pressure screws, short 6.2744.070	6	Purge valve
7	Purge valve rotary knob		

1 Connecting the purging needle

- Push the end of the purging needle (20-3) over the end of the deaerating capillary (20-4) on the purge valve.

2 Connecting the syringe

- Insert syringe (20-1) in the Luer connector (20-2) of the purging needle (see Figure 20, page 32).

3 Opening purge valve

- Open the rotary knob (20-7) by approx. $\frac{1}{2}$ rotation counterclockwise.

4 Setting the flow rate

- Start MagIC Net (if not yet started).
- Ensure that the eluent aspiration tubing is immersed sufficiently in the eluent.
- Let the high pressure pump run.

5 Aspirating eluent

- Aspirate with the syringe (20-1) until bubble-free eluent flows into the syringe.

6 Completing deaerating

- Turn off high pressure pump.
- Close rotary knob (20-7).
- Remove syringe (20-1) from the Luer connector (20-2).
- Pull the purging needle (20-3) out of the deaerating capillary (20-4).

4.7 Inline filter

Between the purge valve and the pulsation damper the 6.2821.120 inline filter is installed as protection against particles.

Inline filters protect the separation column securely against possible contamination from the eluent. Inline filters can however also just as well be used for the purpose of protecting the suppressor against contaminations in the regeneration or rinsing solutions. The fine 2 µm material of the readily and easily replaceable filter platelets removes particles such as bacteria and algae from the solutions.



Note

The inline filter is already installed in the newly delivered instrument. The following installation instructions need **not** be carried out at the time of initial installation.

Installing the inline filter



Caution

Observe the flow direction marked on the filter housing for the connection of the inline filter.




Figure 21 Connecting inline filter

1 Connection capillary Connects the purge valve with the inline filter	2 PEEK pressure screws, short 6.2744.070
3 Inline filter 6.2821.120 Protects against particles.	4 Connection capillary Connects the inline filter with the pulsation damper.

- 1 Screw on the connection capillary running from the purge valve with a 6.2744.070 pressure screw to the input side of the inline filters.
- 2 Screw on the connection capillary running to the pulsation damper with a 6.2744.070 pressure screw to the output side of the inline filter.

4.8 Pulsation damper



Note

The pulsation damper is already installed in the newly delivered instrument.



Caution

The pulsation damper is maintenance-free and may not be opened.

The pulsation damper protects the separation column from damage caused by pressure fluctuations when switching the injection valve, and reduces interfering pulsations during highly sensitive measurements. In order to ensure these functionalities, it must be connected between the high pressure pump (see Chapter 4.6, page 29) and injection valve (see Chapter 4.9, page 35).

The pulsation damper can be operated in both directions.




Figure 22 Pulsation damper – Connection

1	Connection capillary Connection to the inline filter.	2	Fastening screws
3	PEEK pressure screws, short 6.2744.070	4	Holder for pulsation damper
5	Pulsation damper 6.2620.150	6	Connection capillary Connection to the injection valve.

4.9 Injection valve

The injection valve connects the eluent and sample path. Through rapid and precise valve switchover a precise amount of sample solution defined by the size of the sample loop is injected and rinsed with eluent onto the separation column.

4.9.1 Connecting the injection valve

The injection valve has six connectors: two for the sample path (connectors 1 and 2), two for the eluent path (connectors 4 and 5) and two for the sample loop (connectors 3 and 6).



Note

The capillaries of the eluent path and the sample path and the sample loop are already installed in the newly delivered instrument.




Figure 23 Injection valve – connected

1 Injection valve	2 Sample loop Connected to connectors 3 and 6.
3 Connection capillary Connected to connector 4. Carries eluent to the injection valve.	4 Connection capillary (column inlet capillary) Connected to connector 5. Carries eluent to the separation column.
5 Connection capillary Connected to connector 1. Carries sample to the injection valve.	6 Connection capillary Connected to connector 2. Carries sample to the waste container.
7 PEEK pressure screw 6.2744.010	

Replacing the sample loop

The sample loop can be replaced, depending on requirements. For additional information concerning selection of the appropriate sample loop, see *Chapter 4.9.3, page 38*.



Note

Use only 6.2744.010 PEEK pressure screws for connecting capillaries and sample loop to the injection valve.

1 Removing existing sample loop

- Loosen 6.2744.010 pressure screws at connector 3 and connector 6.
- Remove sample loop.

2 Mounting new sample loop

- Fasten one end of the sample loop (23-2) with a 6.2744.010 PEEK pressure screw (23-7) to connector 3.
- Fasten the other end of the sample loop (23-2) with a second 6.2744.010 PEEK pressure screw (23-7) to connector 6.

4.9.2 Mode of operation of the injection valve

The injection valve (see Figure 24, page 37) can adopt two valve positions - **FILL** and **INJECT**. Switching back and forth between the two valve positions determines whether the sample path or the eluent path is guided through the sample loop. The following figure provides a schematic display of the flow paths of the two valve positions.




Figure 24 Injection valve – Positions

A	Position FILL	B	Position INJECT
1	Eluent input Capillary coming from the high pressure pump.	2	Eluent output Capillary to the column.
3	Sample input Sample aspirating capillary.	4	Sample output Capillary to waste container.
5	Sample loop		

Position A

In the position **FILL**, the sample solution flows through the sample loop to the waste container. The eluent flows directly to the separation column at the same time.

Position B

In the position **INJECT**, the eluent flows through the sample loop to the separation column. If sample solution is to be found in the sample loop at the time of the valve switchover, then this will be conveyed along with the eluent, thus making its way to the separation column. The flow in the sample path is either stopped or the sample flows directly to the waste container.



4.9.3 Selecting the sample loop

The amount of sample solution injected depends on the volume of the sample loop. The choice is made on the basis of the application. The following sample loops are normally used:

Cation determination	10 µL
Anion determination with suppression	20 µL
Anion determination without suppression	100 µL

5 Set to work

The extension module is put into operation together with the IC instrument.

Putting IC instrument into operation with extension modules

- 1 Start MagIC Net.
- 2 Connect the 872 Extension Module to the 850 Professional IC.
- 3 Connect the 850 Professional IC instrument to the PC and switch on.
The extension module is recognized automatically by MagIC Net.
- 4 Follow the other instructions in the chapter *Initial start-up* of the manual for the 850 Professional IC.
- 5 Deaerate the high pressure pump at the same time (*see Chapter 4.6.2, page 31*).

6.1.4 Shutting down

If the instrument is shut down for a longer period of time, the entire IC system must be rinsed as follows to rid it of salts in order to prevent eluent salts from forming crystals which may cause subsequent damage.

- rinse all capillaries and the Dosino with methanol/ultra pure water (1:4),
- rinse all pump tubings of the peristaltic pump with ultra pure water.

6.2 Door



Caution

The door is made of PMMA (polymethylmetacrylate). It must never be cleaned with abrasive media or solvents.



Caution

Never use the door as a handle.

6.3 Eluent

6.3.1 Production

The chemicals used for the production of eluents should have a degree of purity of at least "p.a.". Only ultra pure water (resistance > 18.2 MΩ *cm) may be used for dilution (this generally applies for reagents which are used in ion chromatography).

Newly produced eluents should always be microfiltered (filter 0.45 µm).



Caution

Only microfiltered (filter 0.45 µm) eluents may be used.

The composition of the eluent has a crucial effect on the chromatographic analysis:

Concentration

An increase in the concentration generally leads to shorter retention times and faster separation, but also to higher background conductivity.



pH	pH changes result in shifts in the dissociation equilibria and hence changes in the retention times.
Organic solvents	The addition of an organic solvent (e.g. methanol, acetone, acetonitrile) to aqueous eluents generally accelerates lipophilic ions.

6.3.2 Operation

6.3.2.1 Supply bottle

The supply bottle with the eluent must be connected as indicated in *chapter 4.4.1, page 23*. This is above all important for eluents with volatile solvents (e.g. acetone).

Moreover, condensation must also be prevented in the eluent bottle. Drop formation can change the concentration ratio in the eluent.

6.3.2.2 Aspiration filter

To protect the IC system against foreign particles, we recommend aspirating the eluents via a 6.2821.090 aspiration filter (**13-2**). This aspiration filter must be replaced should it show signs of yellow discoloration (but no later than every 3 months).

In the case of very sensitive measurements, the eluent should be stirred constantly with a magnetic stirrer.

6.3.2.3 Changing the eluent

When changing the eluent, it must be ensured that no precipitates can occur. Solutions following one another in direct succession must therefore be miscible. If the system has to be rinsed organically, several solvents with rising or falling lipophilia must be used.

6.4 High pressure pump

6.4.1 Protection



Caution

The pump head is filled ex works with methanol/ultra pure water. It must be ensured that the eluent used is freely miscible with the solvent remaining in the pump head.

To protect the high pressure pump against **foreign particles**, we recommend that the eluent undergoes a **microfiltration** (filter 0.45 µm) before being aspirated via a 6.2821.090 aspiration filter (*see "Assembling eluent aspiration tubing", page 23*).

Salt crystals between the piston and seal cause abrasion particles which can find their way into the eluent. These lead to contaminated valves, a rise in pressure and in extreme cases scratched pistons. It is therefore essential to ensure that **no precipitates** can occur (see Chapter 6.3.2.3, page 42).



Caution

In order to spare the pump seals, the pump should not be operated dry. Therefore ensure that the eluent supply is correctly connected and that there is enough eluent in the eluent bottle each time before turning on the pump.

6.4.2 Maintenance



Caution

Maintenance work on the high pressure pump may not be carried out unless the **instrument is switched off**.

Pump head maintenance

An unstable baseline (pulsation, flow fluctuations) is in many cases the result of contaminated valves (31-2), (31-3) or defective, leaking piston seals on the high pressure pump. Proceed as follows for cleaning contaminated valves and/or replacing worn parts such as pistons, piston seal and valves:

This maintenance work should be carried out at least once a year.

Removing the pump head

- 1 Turn off high pressure pump and wait until pressure is released.
- 2 Loosen the pressure screw on the inlet valve holder (18-2) and unscrew the coupling (18-9), the pump head input capillary (18-7) and the eluent aspiration tubing from the pump head.
- 3 Unscrew the pump head output capillary (18-13) from the pump head.
- 4 Remove pump head from the pump housing by loosening the 4 fastening screws (18-5) using the 6.2621.030 hexagon key. The main



piston is on the left (viewed from the front), and the auxiliary piston is on the right.

Cleaning/replacing the zirconium oxide piston

Clean one piston after the other as follows:

1 Removing the piston cartridge from the pump head

Loosen the piston cartridge with a wrench and unscrew from the pump head by hand.




Figure 25 Removing piston

- 1 Pump head**
- 2 Piston**

2 Dismantle the piston



Caution

On the inside of the piston cartridge there is a taut spring than can jump out of the piston cartridge if suddenly loosing tension.

When opening the piston cartridge, hold pressure towards the spring and unscrew carefully.

- Loosen the screw of the piston cartridge with a wrench and unscrew carefully by hand and by holding pressure towards the taut spring.
- Remove the zirconium oxide piston and lay on a tissue.

- Remove the spring retainer, spring and the inner plastic sleeve from the piston cartridge and lay by.
- Remove the backup ring from the pump head and lay to the other parts.




Figure 26 Components of the piston cartridge

1	Piston cartridge screw	2	Retaining washer
3	Zirconium oxide piston with piston shaft Order number: 6.2824.070.	4	Spring retainer
5	Spring Order number: 6.2824.060.	6	Inner plastic sleeve Protects from metallic abrasion.
7	Piston cartridge	8	Backup ring

3 Cleaning the components of the piston

- Clean zirconium oxide pistons contaminated by abrasion or deposits with pure abrasive cleaning powder, rinse particle free with ultra pure water and dry.
Replace highly contaminated or scratched zirconium oxide pistons (spare part: 6.2824.070 zirconium oxide piston).
- Rinse the other parts of the piston and dry with a lint-free cloth.

4 Assembling the piston

- Insert the inner plastic sleeve, spring and spring retainer into the piston cartridge.
- Slide the zirconium oxide piston carefully into the piston cartridge until its tip emerges from the small opening of the piston cartridge.
- Attach screw and tighten by hand.



Replacing the piston seal

The 6.2617.010 special tool (see Figure 27, page 46) is necessary in order to remove the piston seal from the pump head. It consists of two parts: a tip for removing the old piston seal and a sleeve for inserting the new piston seal.




Figure 27 Tool for piston seal 6.2617.010

1 Pin
Pin for removing the old piston seal.

2 Sleeve
Sleeve for inserting the new piston seal.



Caution

Screwing the 6.2617.010 special tool for the piston seal into the piston seal destroys this completely!

1 Removing the piston seal



Caution

Avoid touching the sealing surface in the pump head (18-4) with the tool.

Screw the special tool for the piston seal (27-1) with the narrow side just as far into the piston seal as the same can be removed.




Figure 28 Removing the piston seal

1 Piston seal

2 Tool for piston seal
Pin of the tool.

2 Inserting the new piston seal into the tool

Insert the new piston seal tightly by hand into the recess of the sleeve of the tool for the piston seal (27-2). The sealing springs must be visible from the outside.




Figure 29 Insert the piston seal into the tool

1 Tool for piston seal 6.2617.010
Sleeve for inserting the new piston seal.

2 Piston seal
Order number: 6.2741.020

3 Inserting the new piston seal into the pump head

Guide the sleeve of the tool for the piston seal (27-2) with inserted piston seal into the pump head and press the seal with the wide end of the tool for the piston seal (27-1) into the pump head recess.




Figure 30 Inserting the piston seal into the pump head

4 Replacing the piston cartridge

Screw the assembled piston cartridge back into the pump head and tighten, first by hand, then additionally by approx. 15° with a wrench.

Cleaning the inlet valve and outlet valve

1 Removing valves

- Unscrew the connection capillary for the auxiliary piston (18-1) from the outlet valve holder.
- Unscrew the holders for the inlet and outlet valves and remove valves.




Figure 31 Removing valves

1 Outlet valve holder**2 Outlet valve**

Order number: 6.2824.160

3 Inlet valve

Order number: 6.2824.170

4 Inlet valve holder**2 Cleaning undissected valve**

Clean contaminated or blocked valves initially **without** dismantling them completely.

- Rinse the valve in eluent flow and counterflow direction using a spray bottle filled with ultra pure water, RBS solution or acetone.
- The rinsing effect is further increased through a short treatment (lasting for a maximum of 20 s) in an ultrasonic bath.

**Note**

Longer lasting ultrasonic baths can damage the ruby ball of the valve.

Only if this cleaning is useless, dismantle the valves separately and clean the components.

3 Dismantling valve

Dismantle every valve separately.

**Note**

For dismantling the valve the 6.2617.020 tool for valve cartridges is required.



- Place the valve with the seal faced downwards above the recess in the holder.
- Push the valve components out of the valve housing using the needle of the tool.




Figure 32 Dismantling valve

1 Needle

For pushing the valve components out of the valve housing.

2 Valve

3 Holder

4 Recess

For collecting the valve components.

The components of the valve are collected in the recess of the holder.



Note

The components of the valve are very small. In order not to lose them, put the components into a dish.

- The inlet valve and the outlet valve consist of the same, just differently arranged components (see Figure 33, page 51).




Figure 33 Components of the inlet valve and outlet valve

1	Inlet valve 6.2824.170	2	Outlet valve 6.2824.160
3	Inlet valve housing	4	Outlet valve housing
5	Sealing ring (black)	6	Sleeve
7	Sapphire sleeve The shiny side must point to ruby ball.	8	Ruby ball
9	Ceramic holder for ruby ball	10	Seal The larger opening must point outwards.

4 Clean the components of the valve

Rinse the valve components with ultra pure water and/or acetone and dry with a lint-free cloth.

5 Reassemble the valve

Reassemble valve components *according to figure 33, page 51*.

- Insert the seal with the larger opening faced downwards into the recess of the tool.
- Lay the other valve components above another in the correct sequence (*see Figure 33, page 51*).



- Place the valve housing over the stacked components and hold it tightly.
- By tilting the tool, the valve components slide into the valve housing.
- Press the seal by hand well on the valve housing.

6 Checking the flow direction

Rinse the valve in the direction of the arrow on the valve housing and check whether liquid is escaping on the other end.

If this is not the case, the valve has to be dismantled again and be assembled correctly (see *Figure 33, page 51*).

7 Inserting the valves back into the pump head



Caution

If by mistake, the inlet valve is mounted instead of the outlet valve, an extreme pressure builds up within the working cylinder, which can destroy the piston seal!

When inserting the valves, please take into account that the liquid is being pumped through the pump head from bottom to top.

- Insert the inlet valve into the inlet valve holder the way the seal is visible.
- Screw the inlet valve holder into the bottom of the pump head and tighten with a wrench (31-4).
- Insert the outlet valve into the outlet valve holder the way the seal is visible.
- Screw the outlet valve holder into the top of the pump head and tighten with a wrench (31-1).

Mounting the pump head



Note

To prevent the pump head from being positioned the wrong way, it is provided with different bore hole depths for the fastening bolts, i. e. a fastening bolt is longer than all others. The bore hole with the greatest depth must therefore be assigned to the longest bolt. If this is not the case, the pump will not function perfectly.

- 1** Mount the pump head on the pump again using the four fastening screws (18-5). Firmly tighten the screws with the 6.2621.030 hexagon key.
- 2** Screw connection capillaries (18-1), (18-7) and (18-13) onto the pump head again.

6.5 Inline filter

6.5.1 Maintenance

The 6.2821.120 inline filters comprise the filter housing (34-2), the filter screw (34-4) and the filter (34-3). New filters (34-3) are available under the order number 6.2821.130 (10 items).

The 6.2821.130 filters (21-3) should be changed every 3 months (more frequently at higher backpressure).




Figure 34 Changing the filter

1 PEEK pressure screws, short
6.2744.070

3 6.2821.130 filter
Packaging contains 10 items.

5 Connection capillaries

2 Filter housing
Housing of the inline filter. Part of the 6.2821.120 accessories.

4 Filter screw
Screw of the inline filter. Part of the 6.2821.120 accessories.



Changing the filter

The flow must be stopped before changing the filter.

1 Removing the inline filter

- Unscrew the pressure screws (34-1) from the inline filter.

2 Unscrewing the filter screw

- Screw the filter screw (34-4) out of the filter housing (34-2) with the aid of two 6.2621.000 adjustable wrenches.

3 Inserting the filter

- Remove the old filter (34-3) with tweezers.
- Place the new filter (34-3) flat in the filter housing with tweezers (34-2).

4 Mounting filter screw

- Screw the filter screw (34-4) back into the filter housing (34-2) and tighten by hand. Then additionally tighten slightly with two 6.2621.000 adjustable wrenches.

5 Remounting the inline filter

- Screw the pressure screws (34-1) back onto the inline filter.

6 Rinsing the inline filter

- Dismantle the guard column (if present) and the separation column and replace with a 6.2744.040 coupling.
- Rinse the instrument with eluent.

6.6 Inline sample preparation

To protect the separation column against foreign particles which can affect the separating efficiency, we recommend that all samples undergo a microfiltration (filter 0.45 µm). The ultrafiltration cell can be used for **filtration** (see manual of the *IC Equipment for Ultrafiltration*).

Matrix-loaded samples (e.g. blood, oil) should be prepared for the measurement by means of dialysis (see manual on the *IC Equipment for Dialysis*).

If the concentration of the sample is too high, the sample should be **diluted** before feeding (see documentation on the *IC Equipment for Sample Dilution*).

For the sample preparation methods **Neutralization** (replacement of e.g. Na⁺ with H⁺) and **cation exchange** (replacement of e.g. heavy metals with H⁺), a sample preparation module (SPM) is used.

For an overview of all Metrohm inline sample preparation methods go to the following website: <http://misp.metrohm.com>

6.7 Rinsing the sample path

Before a new sample can be measured, the sample path must be rinsed with it so that the measuring result is not falsified by the previous sample (**Sample carry-over**).

In the case of automated sample feeding, the rinsing time should be at least 3 times the **transfer time**.

The transfer time is the time required by the sample to flow from the sample vessel to the end of the sample loop. **##NO_MATCH##**.

Ascertaining the transfer time

To ascertain the transfer time, proceed as follows:

1 Emptying the sample path

Pump air through the sample path (pump tubing, tubing connections, sample loop) for several minutes until all liquid is displaced by the air.



2 Aspirating the sample and measuring time

Aspirate a sample typical for the later application and use a stop watch to measure the time required by the sample to travel from the sample vessel to the end of the sample loop.

The time measured corresponds to the "transfer time". The rinsing time should be at least 3 times the transfer time.

Checking the rinsing time

It is possible to determine whether the rinsing time is adequate via a direct measurement of the sample carry-over. Proceed as follows:

1 Preparing two samples

- **Sample A:** A typical sample for the application.
- **Sample B:** Ultrapure water.

2 Determining "Sample A"

Let "Sample A" pass through the sample path for the duration of the rinsing time, then inject and measure.

3 Determining "Sample B"

Let "Sample B" pass through the sample path for the duration of the rinsing time, then inject and measure.

4 Calculating the sample carry-over

The degree of the sample carry-over corresponds to the ratio of the peak areas of the measurement for sample B to the measurement for sample A. The lower the ratio, the lower the sample carry-over. This ratio can be modified by varying the rinsing time – thus allowing the rinsing time required for the application to be ascertained.

6.8 Injection valve

6.8.1 Protection

To prevent contamination of the injection valve, a 6.2821.120 inline filter (see Chapter 4.7, page 33) should be mounted between the high pressure pump and the pulsation damper.

6.9 Quality Management and validation with Metrohm

Quality Management

Metrohm offers you comprehensive support in implementing quality management measures for instruments and software. Further information on this can be found in the brochure «**Quality Management with Metrohm**» available from your local Metrohm agent.

Validation

Please contact your local Metrohm agent for support in validating instruments and software. Here you can also obtain validation documentation to provide help for carrying out the **Installation Qualification** (IQ) and the **Operational Qualification** (OQ). IQ and OQ are also offered as a service by the Metrohm agents. In addition, various application bulletins are also available on the subject, which also contain **Standard Operating Procedures** (SOP) for testing analytical measuring instruments for reproducibility and correctness.

Maintenance

Electronic and mechanical functional groups in Metrohm instruments can and should be checked as part of regular maintenance by specialist personnel from Metrohm. Please ask your local Metrohm agent regarding the precise terms and conditions involved in concluding a corresponding maintenance agreement.



Note

You can find information on the subjects of quality management, validation and maintenance as well as an overview of the documents currently available at www.metrohm.com/com/ under **Support**.



7 Troubleshooting

7.1 Problems and their solutions

Problem	Cause	Remedy
Marked rise in pressure	<i>6.2821.120 inline filter blocked.</i>	Replace 6.2821.130 filter (<i>see Chapter 6.5, page 53</i>).
	<i>Injection valve – valve blocked.</i>	Have the valve cleaned (by Metrohm service technicians).
Drift of the baseline	<i>Eluent – Evaporation of the organic solvent in the eluent.</i>	Check the eluent bottle attachment (<i>see Figure 14, page 25</i>).
Peak areas lower than expected	<i>Sample – leak in the sample path.</i>	Check the sample path.
	<i>Sample – blockage in the sample path.</i>	Check the sample path.
	<i>Sample – sample loop not (completely) filled.</i>	Prolong the sample transfer time.
Very noisy baseline	<i>High pressure pump – contaminated pump valves.</i>	Clean pump valves (<i>see Chapter 6.4.2, page 43</i>).
	<i>Eluent – Leak in the eluent path.</i>	Check the eluent path
	<i>Eluent – Blockage in the eluent path.</i>	Check the eluent path
	<i>High pressure pump – defective piston seals.</i>	Replace piston seals (<i>see Chapter 6.4.2, page 43</i>).
	<i>Pulsation damper not connected.</i>	Connect pulsation damper (<i>see Chapter 4.8, page 34</i>).
Background conductivity too high	<i>Incorrect eluent.</i>	Change eluent (<i>see Chapter 6.3.2.3, page 42</i>).
Individual peaks greater than expected	<i>Sample – carry-over of the samples from previous measurement.</i>	Rinse system longer between two samples.
Poor reproducibility of the retention times	<i>Eluent – Leak in the eluent path.</i>	Check the eluent path

Problem	Cause	Remedy
	<i>Eluent – Blockage in the eluent path.</i>	Check the eluent path
Precision problems - significant scattering of the measured values	<i>Injection valve – sample loop.</i>	Check installation of the sample loop (see Chapter 4.9.1, page 35).
	<i>Sample – rinsing volume too low.</i>	Increase rinsing time (see Chapter 6.7, page 55).
	<i>Injection valve – defective.</i>	Request Metrohm Service.
Unexpected change to the retention times in the chromatograms	<i>Eluent – Gas bubbles in the eluent.</i>	Check connections of the eluent degasser (see Chapter 4.5, page 27).
	<i>High pressure pump – defective.</i>	Request Metrohm Service.
Vacuum is not being built	<i>Eluent Degasser – Connector Vacuum on the rear of the instrument not (tightly) sealed.</i>	<ul style="list-style-type: none"> ▪ Seal the connector Vacuum tightly with a 6.1446.040 threaded stopper.



8 Technical specifications

8.1 Reference conditions

The technical specifications listed in this chapter refers to the following reference conditions:

<i>Ambient temperature</i>	+25 °C (±3 °C)
<i>Instrument status</i>	> 40 minutes in operation (equilibrated)

8.2 Instrument

<i>IC system</i>	Metal-free IC system
<i>Material</i>	Painted polyurethane hard foam without CFCs, fire class V0
<i>Intelligent components</i>	MagIC Net

8.3 Ambient conditions

<i>Operation</i>	
<i>Ambient temperature</i>	+5...+45 °C
<i>Humidity</i>	20...80 % relative humidity
<i>Storage</i>	
<i>Ambient temperature</i>	-20...+70 °C
<i>Transport</i>	
<i>Ambient temperature</i>	-40...+70 °C

8.4 Housing

Dimensions

<i>Width</i>	365 mm
<i>Height</i>	131 mm
<i>Depth</i>	380 mm

<i>Material of base tray, housing and bottle holder</i>	Polyurethane hard foam (PUR) with flame retardation for fire class UL94V0, CFC-free, coated
---	---

8.5 Eluent degasser

<i>Material</i>	fluoropolymer
<i>Resistance to solvents</i>	No restriction (apart from PFC)
<i>Build-up time for the vacuum</i>	< 60 s

8.6 High pressure pump

<i>Type</i>	<ul style="list-style-type: none"> ▪ Serial dual-piston pump ▪ Intelligent pump head recognition ▪ Chemically inert ▪ Metal-free pump heads ▪ Materials in contact with eluent: PEEK, ZrO₂, PTFE/PE ▪ Self-optimizing flow and pressure
-------------	--

Flow rate

<i>Adjustable flow range</i>	0.001...20.0 mL/min
<i>Flow increment</i>	1 µL/min
<i>Reproducibility of the eluent flow</i>	< 0.1 % deviation

Pressure range

<i>Pump</i>	0...50.0 MPa (0...500 bar)
<i>Pump head</i>	0...35.0 MPa (0...350 bar) (applies for the standard PEEK pump head)
<i>Residual pulsation</i>	< 1 %

*Safety shutdown*

<i>Function</i>	Automatic shutdown upon reaching the pressure limit values
<i>Maximum pressure limit</i>	<ul style="list-style-type: none"> ▪ Adjustable from 0.1...50 MPa (1...500 bar) ▪ The pump is automatically shut down at the first piston stroke above the maximum limit value
<i>Minimum pressure limit</i>	<ul style="list-style-type: none"> ▪ Adjustable from 0...49 MPa (0...490 bar) ▪ The shutdown mechanism is inactive at 0 MPa ▪ The shutdown mechanism only becomes active 2 minutes after system start ▪ The pump is automatically shut down after 3 piston strokes below the minimum pressure limit
<i>Gradient capacity</i>	Isocratic or gradient (extendable to quaternary)
<i>Profile</i>	Step, linear, convex and concave
<i>Resolution</i>	< 1 nL/min flow increments

8.7 Injection valve

<i>Actuator time</i>	typ. 100 ms
<i>Max. operating pressure</i>	35 MPa (350 bar)
<i>Material</i>	PEEK

8.8 Interfaces

<i>Auxiliary</i>	1 DSUB plug 15-pin (female)
<i>Analog Output</i>	Analog Output (optional)

8.9 Safety specification

<i>Design / testing</i>	<ul style="list-style-type: none"> ▪ EN/IEC 61010-1 ▪ UL 61010-1 ▪ CSA-C22.2 No. 61010-1 ▪ Protection class III
-------------------------	---

8.10 Electromagnetic compatibility (EMC)

- Emission*
- EN/IEC 61326-1
 - EN/IEC 61000-6-3
 - EN 55011 / CISPR 11

- Immunity*
- EN/IEC 61326-1
 - EN/IEC 61000-6-1
 - EN/IEC 61000-4-2
 - EN/IEC 61000-4-3

8.11 Weight

1.872.0030 7.7 kg (without accessories)



9 Conformity and warranty

9.1 Declaration of Conformity

This is to certify the conformity to the standard specifications for electrical appliances and accessories, as well as to the standard specifications for security and to system validation issued by the manufacturing company.

Name of commodity

872 Extension Module

The 872 Extension Module is an expansion tool for upgrading all 850 Professional IC instruments.

This instrument has been built and has undergone final type testing according to the standards:

Electromagnetic compatibility

Emission: EN/IEC 61326-1: 2006, EN/IEC 61000-6-3: 2004, EN 55011 / CISPR 11: 2003

Immunity: EN/IEC 61326-1: 2006, EN/IEC 61000-6-1: 2007, EN/IEC 61000-4-2: 2001, EN/IEC 61000-4-3: 2002.

Safety specifications

EN/IEC 61010-1: 2001, UL 61010-1: 2004, CSA-C22.2 No. 61010-1: 2004, protection class III



This instrument meets the requirements of the CE mark as contained in the EU directives 2006/95/EC (LVD), 2004/108/EC (EMC). It fulfils the following specifications:

EN 61326-1: 2006 Electrical equipment for measurement, control and laboratory use – EMC requirements

EN 61010-1: 2001 Safety requirements for electrical equipment for measurement, control and laboratory use

Manufacturer

Metrohm Ltd., CH-9101 Herisau/Switzerland

Metrohm Ltd. is holder of the SQS certificate ISO 9001:2000 Quality management system for development, production and sales of instruments and accessories for ion analysis.

Herisau, 31 March, 2008



D. Strohm

Vice President, Head of R&D



Ch. Buchmann

Vice President, Head of Production

Responsible for Quality Assurance

9.2 Quality Management Principles

Metrohm Ltd. holds the ISO 9001:2000 Certificate, registration number 10872-02, issued by SQS (Swiss Association for Quality and Management Systems). Internal and external audits are carried out periodically to assure that the standards defined by Metrohm's QM Manual are maintained.

The steps involved in the design, manufacture and servicing of instruments are fully documented and the resulting reports are archived for ten years. The development of software for PCs and instruments is also duly documented and the documents and source codes are archived. Both remain the possession of Metrohm. A non-disclosure agreement may be asked to be provided by those requiring access to them.

The implementation of the ISO 9001:2000 quality management system is described in Metrohm's QM Manual, which comprises detailed instructions on the following fields of activity:

Instrument development

The organization of the instrument design, its planning and the intermediate controls are fully documented and traceable. Laboratory testing accompanies all phases of instrument development.

Software development

Software development occurs in terms of the software life cycle. Tests are performed to detect programming errors and to assess the program's functionality in a laboratory environment.

Components

All components used in the Metrohm instruments have to satisfy the quality standards that are defined and implemented for our products. Suppliers of components are audited by Metrohm as the need arises.

Manufacture

The measures put into practice in the production of our instruments guarantee a constant quality standard. Production planning and manufacturing procedures, maintenance of production means and testing of components, intermediate and finished products are prescribed.



Customer support and service

Customer support involves all phases of instrument acquisition and use by the customer, i.e. consulting to define the adequate equipment for the analytical problem at hand, delivery of the equipment, user manuals, training, after-sales service and processing of customer complaints. The Metrohm service organization is equipped to support customers in implementing standards such as GLP, GMP, ISO 900X, in performing Operational Qualification and Performance Verification of the system components or in carrying out the System Validation for the quantitative determination of a substance in a given matrix.

9.3 Warranty (Guarantee)

Metrohm guarantees that the deliveries and services it provides are free of errors in materials, design or manufacturing.

The general warranty period is 36 months (exclusions below) from the date of delivery or 18 months in the event of continuous operation. The warranty remains valid on the condition that the service is provided by an authorized Metrohm Service Organization at defined intervals and with a defined scope.

The warranty period for the suppressors "MSM II" and "MSM-HC" is 120 months from the date of delivery or 60 months in the event of continuous operation.

The warranty period for IC separation columns is 12 months from the date of delivery or 6 months in the event of continuous operation. The technical specifications contained in the manual are authoritative for warranty of accuracy.

For third-party components that are recognizable as such, the manufacturer's warranty regulations apply.

Consumables and materials with limited storage life and glass breakage in the case of electrodes or other glass parts are excluded from the warranty.

Warranty claims cannot be asserted if the customer has failed to meet his payment obligations according to schedule.

During the warranty period, Metrohm undertakes either to replace free of charge or to credit the purchaser for any instruments, assemblies or components that can be shown to be faulty. Any transport or customs fees that may apply are the ordering party's responsibility.

The precondition for this is that the ordering party must use the Return Material Authorization (RMA) to report the faulty part, along with specification of the article number, the article designation, an adequate error description, the delivery date and (if applicable) the serial number. In addi-

tion, the ordering party undertakes to store the faulty part for at least 2 years in accordance with current storage directives (in compliance with ESD guidelines) and to hold it in readiness for onsite inspection or for return shipment to Metrohm. Metrohm reserves the right to invoice the ordering party for these articles, including retroactively, in the event of noncompliance with these pre-conditions.

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc. are expressly excluded from the warranty.

Metrohm also offers a 120-month spare parts warranty and a 5-year PC software support warranty, calculated from the date on which the product is withdrawn from the market. The content of this warranty is the ability of the customer to obtain functioning spare parts or appropriate software support at market prices during the time of the warranty period.

If Metrohm AG is unable to meet this obligation due to circumstances beyond the control of Metrohm AG, then the ordering party shall be offered alternative solutions at preferential conditions.



10 Accessories








Note

Subject to change without notice.





10.1 Scope of delivery

2.872.0030 872 Extension Module – IC Module

Qty.	Order no.	Description	
1	1.872.0030	872 Extension Module – IC Module	
1	6.1602.160	Eluent bottle attachment GL 45 For eluent bottles; with connections for drying tube and aspiration tubing. Opening ground joint: A-14/15	
1	6.1608.070	Eluent bottle / 2 L / GL 45 Material: Clear glass Height (mm): 262 Volume (mL): 2000	

Qty.	Order no.	Description	
1	6.1609.000	Adsorbing tube, large and bent For filling with adsorber material. Material: Glass Height (mm): 129 Inner diameter (mm): 32 SGJ size: B-14/15	
2	6.1803.040	PTFE capillary 0.5 mm i.d. / 1 m Capillary for sample handling in IC. Material: PTFE Outer diameter (inches): 1/16 Inner diameter (mm): 0.5 Length (m): 1	
1	6.1825.230	PEEK sample loop 10 µL For injection valve, with 2 PEEK pressure screws Material: PEEK (metal-free) Outer diameter (inches): 1/16 Volume (mL): 0.01	







Qty.	Order no.	Description	
1	6.1831.010	PEEK capillary 0.25 mm i.d. / 3 m	
		For all IC components.	
		Material: PEEK	
		Outer diameter (inches): 1/16	
		Inner diameter (mm): 0.25	
		Length (m): 3	
1	6.1834.080	Aspiration tubing, 2 m	
		Aspiration tubing for Professional IC instruments	
		Material: PTFE	
		Outer diameter (mm): 2.5	
		Inner diameter (mm): 1.5	
		Length (m): 2	
1	6.2023.020	Clip for SGJ 14/15	
		Clip for SGJ 14/16	
		Material: POM	
1	6.2156.060	Cable Extension Module - Professional IC, 40 cm	
		Cable connecting the extension module with a Professional IC instrument	
		Length (m): 0.4	





Qty.	Order no.	Description	
1	6.2617.010	Tool for piston seal For removing and assembling the piston seal for all standard pump heads	
1	6.2621.000	Adjustable wrench Max. opening: 20 mm. For IC instruments Length (mm): 150	
1	6.2621.030	Hexagon key 4 mm Length (mm): 73	
1	6.2621.050	1/4 in. wrench For 1/4 in. screws. For IC instruments Length (mm): 73	



Qty.	Order no.	Description	
1	6.2621.080	Capillary cutter For plastic capillaries. Used with IC instruments Length (mm): 118	
1	6.2621.100	Hexagon key 3 mm Hexagon key 3 mm for IC Sample Processors Length (mm): 73	
2	6.2739.000	Wrench For tightening connectors Length (mm): 68	
1	6.2744.010	Pressure screw 5x With UNF 10/32 connection. For the connection of PEEK capillaries Material: PEEK Length (mm): 26	

Qty.	Order no.	Description
1	6.2744.014	Pressure screw 2x With UNF 10/32 connection. For the connection of PEEK capillaries Material: PEEK Length (mm): 26
		
1	6.2744.070	Pressure screw short Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries Material: PEEK Length (mm): 21
		
1	6.2744.120	Coupling pressure screw / Luer-F Connection pressure screw and syringe. Material: PEEK
		
1	6.2744.210	Tubing adaptor for aspiration filter (ProfIC) For Professional IC instruments
		



Qty.	Order no.	Description	
1	6.2816.020	Syringe 10 mL with Luer connection	
		For various applications in IC and VA	
		Material:	PP
		Length (mm):	102
		Volume (mL):	10
1	6.2816.040	Purging needle	
		With PTFE tubing and Luer connection. For syringes. For aspirating eluents.	
1	6.2821.090	Aspiration filter	
		Pore size 20 µm, set of 5 pieces. For 6.1834.000 Aspiration tubing and 6.1821.040 and 6.1821.050 Filter tubes.	
		Material:	PE
		Outer diameter (mm):	9.5
		Length (mm):	35.5
1	6.2821.130	Spare filter for inline filter	
		Spare filters for inline filter.	
1	8.872.8004EN	Manual for 872 Extension Module, 2.872.0030 IC Module, English	

Index

Numbers/Symbols

- 6.2821.090 aspiration filter 42
6.2821.130 filter 53

A

- Accessories
 Scope of delivery 68
Ambient conditions 60
Aspiration filter 6.2821.090 42
Aspiration tubing for eluent 23

B

- Baseline
 Unstable 43
Blood 55

C

- Carry-over 55
Changing
 Eluent 42
Cleaning
 Valves of the high pressure
 pump 48
Contamination
 High pressure pump 43
 Valves in the high pressure
 pump 43
Crystallization
 High pressure pump 43

D

- Deaerating
 High pressure pump 31
 Purge valve 29
Degasser
 Eluent degasser 27
Degassing
 Eluent 27
Dilution 55
Dimensions 61
Door 41

E

- Electrostatic charge 5
Eluent
 Aspiration 23
 Production 41
Eluent bottle
 Figure 26
 Installation 23

- Operation 42
Eluent degasser
 Installation 27
 Technical specifications 61

F

- Fill
 Injection valve 37
Filter
 see also "inline filter" 33
Filter 6.2821.090
 Aspiration filter 42
Flow fluctuations 43
Flow increment 61
Flow range 61
Flow rate 61
Front 7

G

- Gas 27
GLP 57
Guarantee 66

H

- High pressure pump
 Installation 29
 Maintenance 42
 Protection 20, 42
 Technical specifications 61
 Tubing connection 29
 Valves 51
Housing 61
Humidity 60

I

- Inject
 Injection valve 37
Injection valve 2
 Fill 37
 Inject 37
 Installation 35, 62
 Maintenance 57
 Protection 57
Inline filter 33
Inline sample preparation 55
Installation
 Eluent bottle 23
 Eluent degasser 27
 High pressure pump 29
 Injection valve 35, 62

- Overview 21
 Pulsation damper 34
Instrument
 Front 7
 Rear 8

L

- Leak 43
Leaking piston seals 43
Loop
 see also "Sample loop" 38

M

- Mains voltage 4
Maintenance
 High pressure pump 42
 Injection valve 57
 Pump head 43
Material 61

O

- Oil 55
Overview of the instrument 7

P

- Piston seal 43
Pistons of the high pressure pump
 43
Precipitates 43
Pressure limit 62
Pressure range 61
Protection
 Injection valve 57
 Inline filter 33
Pulsation 43
Pulsation damper
 Installation 34
Pump head
 Maintenance 43
Purge valve 29

Q

- Quality Management 57

R

- Rear 8
Reference conditions 60
Regeneration 40
Rinsing
 Sample path 55

Index

Rinsing time 56
Rise in pressure 43

S

Safety instructions 4
Safety shutdown 62
Sample
 Carry-over 55
 Sample loop 38
 Transfer time 55
Sample loop 38
Sample path
 Rinsing 55
Sample preparation 55
Scope of delivery 68

Separation column
 Protection 2, 34
Service 4, 40
Service Agreement 57
Shutting down 41
Storage 60

T

Technical specifications
 Eluent degasser 61
 High pressure pump 61
 Reference conditions 60
Temperature 60
Transfer time 55
Transport 60

Transport locking screws 20

V

Vacuum pump
 Protection 20
Validation 57
Valve
 see also "Injection valve" 35
Valves of the high pressure pump
..... 51

W

Warranty 66