

894 Professional CVS



Manual – Short Instructions

8.894.8002EN / v3 / 2025-09-30



Metrohm AG
CH-9100 Herisau
Switzerland
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

894 Professional CVS

Manual – Short Instructions

Technical Communication
Metrohm AG
CH-9100 Herisau

This documentation is protected by copyright. All rights reserved.

This documentation is an original document.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Disclaimer

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g., conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

Table of contents

1	About these Short Instructions	1
2	Introduction	2
2.1	Device description	2
2.2	Displaying accessories	3
2.3	About the documentation	4
2.3.1	Further information and literature	4
2.3.2	Symbols and conventions	4
3	Overview of the instrument	6
3.1	Front	6
3.2	Rear	7
3.3	RDE measuring head	9
3.4	Measuring head connector plate and measuring head insert	10
3.5	Tubing connector (measuring head arm)	13
4	Installation	15
4.1	Setting up the device	15
4.1.1	Packaging	15
4.1.2	Checks	15
4.1.3	Location	15
4.2	Equipping the RDE measuring head	15
4.2.1	Preparing the RDE measuring head	16
4.2.2	Preparing electrodes and inserting them in the RDE measuring head	19
4.2.3	Inserting the RDE measuring head	25
4.2.4	Connecting the inert gas supply	28
4.3	Establishing the tubing connections	30
4.3.1	Removing the measuring head cover	30
4.3.2	Installing the 4-way micro dosing tip	31
4.3.3	Installing capillaries	34
4.3.4	Installing FEP tubing	36
4.4	Connecting devices electrically	43
4.4.1	Connecting the instrument to the power grid	43
4.4.2	Connecting the 894 Professional CVS	44
4.4.3	Connecting an 800 Dosino	46
4.4.4	Connecting a Sample Processor	49
4.4.5	Connecting an external pump	51
4.4.6	Connecting USB devices directly to the measuring device	53



- 5 Start-up** **54**
- 6 Troubleshooting** **58**
 - 6.1 894 Professional CVS **58**
 - 6.2 Peripheral devices **60**
- 7 Displaying accessories** **62**
- Index** **63**

Table of figures

Figure 1	894 Professional CVS front	6
Figure 2	894 Professional CVS rear	7
Figure 3	RDE measuring head - Overview	9
Figure 4	Measuring head connector plate	10
Figure 5	Measuring head insert	12
Figure 6	Tubing connector (measuring head arm)	13
Figure 7	Removing the stopper from the pipetting opening	16
Figure 8	Removing the measuring head cover	17
Figure 9	Inserting the gas inlet	18
Figure 10	Connecting the gas inlet	19
Figure 11	Removing the protective cap from the electrode tip	20
Figure 12	Tightening the electrode tip to the driving axle	20
Figure 13	Working electrode, installed	20
Figure 14	Fastening the drive belt	21
Figure 15	Connecting the working electrode	22
Figure 16	Assembling the reference electrode with the electrolyte vessel	23
Figure 17	Connecting a reference electrode	24
Figure 18	Connecting the auxiliary electrode	25
Figure 19	Installing the measuring head cover	26
Figure 20	Inserting the stopper in the pipetting opening	26
Figure 21	Do not touch the drive disk	27
Figure 22	Inserting the measuring head	28
Figure 23	Gas washing glass	29
Figure 24	Nipple for inert gas supply	30
Figure 25	Removing the stopper from the pipetting opening	31
Figure 26	Removing the measuring head cover	31
Figure 27	Inserting the 4-way micro dosing tip	32
Figure 28	Screwing a PTFE capillary to a dosing unit	33
Figure 29	Sealing the PTFE capillary of a 4-way micro dosing tip	33
Figure 30	Inserting a PEEK capillary into the measuring head	34
Figure 31	Screwing a capillary to a dosing unit	36
Figure 32	Inserting the FEP tubing into the measuring head	38
Figure 33	Installing the measuring head cover	39
Figure 34	Inserting the stopper in the pipetting opening	40
Figure 35	Installing a bottle cap with pieces of tubing on a rinsing canister	41
Figure 36	Connecting the five-way tubing connector to a waste canister	42
Figure 37	Connecting FEP tubing to the dosing unit (auxiliary solution)	43
Figure 38	Connecting a dosing device to an MSB socket	46
Figure 39	Connecting an external pump using the remote cable	52
Figure 40	Connecting an external pump to the tower	52

1 About these Short Instructions

This short instruction manual contains important chapters from the comprehensive manual. In addition to an introduction, safety instructions and an overview of the instrument, you will also find information about installing and starting up the 894 Professional CVS as well as a chapter on troubleshooting. You can download the detailed manual as a PDF file from the Internet.

Downloading the manual


You can find the detailed manual on the Internet under <https://www.metrohm.com/>:

1. Enter the order number of your instrument as the search term (e. g. **2.894.X210**)
2. Click on **>More information**.
3. Click on **Documents**.
All available documents for the instrument will be displayed.
4. Click on the PDF link to download the desired manual.

2.2 Displaying accessories

Up-to-date information on the scope of delivery and on optional accessories can be found on the Metrohm website.

1 Searching for a product on the website

- Go to <https://www.metrohm.com>.
- Click on .
- Enter the article number of the product (e.g. **2.1001.0010**) into the search field and press **[Enter]**.

The search result is displayed.

2 Displaying product information

- To display the products matching the search term, click on **Product models**.
- Click on the desired product.

Detailed information regarding the product is displayed.

3 Displaying accessories and downloading the accessories list

- To display the accessories, scroll down to **Accessories and more**.
 - The **scope of delivery** is displayed.
 - Click on **[Optional parts]** for the optional accessories.
- To download the accessories list, click on **[Download accessories PDF]** under **Accessories and more**.



NOTE

Metrohm recommends keeping the accessories list for reference purposes.

2.3 About the documentation



CAUTION

Please read through this documentation carefully before putting the device into operation. The documentation contains information and warnings which the user must follow in order to ensure safe operation of the device.




2.3.1 Further information and literature

You can find additional information regarding the 894 Professional CVS in the following publications:

- CVS tutorial
- Multimedia guide "Electrodes in Voltammetry"
- Online help for the **viva** computer software
- Monograph "Voltammetric analysis methods in electroplating"

2.3.2 Symbols and conventions

The following symbols and formatting may appear in this documentation:

<i>(5-12)</i>	Cross-reference to figure legend The first number refers to the figure number, the second to the device part in the figure.
1	Instruction step Perform the steps one after the other.
Method	Dialog text, parameter in the software
File ► New	Menu or menu item
[Continue]	Button or key
	WARNING This symbol draws attention to a possible life-threatening hazard or risk of injury.
	WARNING This symbol draws attention to a possible hazard due to electrical current.
	WARNING This symbol draws attention to a possible hazard due to heat or hot instrument parts.

**WARNING**

This symbol draws attention to a possible biological hazard.

**WARNING**

Warning of optical radiation

**CAUTION**

This symbol draws attention to possible damage to devices or device parts.

**NOTICE**

This symbol highlights additional information and tips.

3 Overview of the instrument

The following figures provide a detailed overview of the parts of the 894 Professional CVS. Some of these parts are not specifically relevant for CVS analyses, but for the sake of completeness they are listed nevertheless; however, they are grayed out and labeled "Not relevant for CVS analyses".

3.1 Front




Figure 1 894 Professional CVS front

1 "Status" LED

Continuously on: instrument ready for operation. Blinking regularly: instrument operating. Blinking pattern "LED on a long time - off a short time - on a long time - off a short time, etc.": Standby potential is being applied to the electrodes. Do not remove the electrode cables!

2 "On" LED

Illuminated if the 894 Professional CVS is connected to the power grid.

3 Holder for measuring vessel For inserting the measuring vessel.	4 Connector plate measuring head arm For inserting the measuring head.
5 Measuring head arm (tiltable)	6 Gas washing glass (6.2405.030) <i>Not relevant for CVS analyses.</i>
7 Decanting glass (6.2405.030) For the deposition of solid materials that may be present.	8 Holder for drip pan For positioning the drip pan.

3.2 Rear




Figure 2 894 Professional CVS rear

1 Temperature sensor connector (Temp.) For connecting a temperature sensor of the type Pt1000. Two B sockets, 2 mm.	2 Type plates With serial number.
--	--

**3 Calibrator**

For service procedures carried out by Metrohm.

5 Tubing connector (measuring head arm)

See Chapter 3.5, page 13.

7 "Controller" connector

For connecting to a PC with the **viva** computer software installed. Mini DIN, 8-pin.

9 Power socket**4 Fan**

Runs during operation.

6 USB connectors (USB 1 and USB 2), type A

For connecting barcode readers, keyboards, etc.

8 MSB connectors (MSB 1 to 4)

Metrohm Serial Bus. For connecting dosing devices (800 Dosino) and Remote Boxes. Mini DIN, 8-pin.

3.3 RDE measuring head




Figure 3 RDE measuring head - Overview

1 Measuring head insert

With openings for inserting electrodes and tubing connections (see figure 5, page 12).

2 Measuring head connector plate

For connecting the RDE measuring head to the connector plate of the measuring head arm (1-4).

For connecting the electrodes and tubing (see figure 4, page 10).



3 Measuring head cover
For shielding against electromagnetic interference.

4 Slide lock
For measuring head cover.

5 Stopper (6.2709.100)
For closing the pipetting opening (5-22).

3.4 Measuring head connector plate and measuring head insert




Figure 4 Measuring head connector plate

1 M6 threaded opening (3)
For adding or aspirating solutions. Can be connected to one of the openings **24 - 27** (FEP tubing from 6.1829.070).

2 M6 threaded opening (2)
For adding or aspirating solutions. Can be connected to one of the openings **24 - 27** (FEP tubing from 6.1829.070).

<p>3 M6 threaded opening (N₂) With preinstalled stopper. <i>Not relevant for CVS analyses.</i></p>	<p>4 M6 threaded opening (WASTE) For aspirating the measuring solution. Can be connected to one of the openings 24 - 27 (FEP tubing from 6.1829.070).</p>
<p>5 UNF 10/32 threaded opening (SMP) For the automated sample addition. Must be connected to opening 28 (PEEK capillary 6.1831.020).</p>	<p>6 M6 threaded opening (1) For adding or aspirating solutions. Can be connected to one of the openings 24 - 27 (FEP tubing from 6.1829.070).</p>
<p>7 Electrode connector (AE, RE, WE) With electrode cables, for connecting electrodes.</p>	<p>8 Threaded opening (PURGE) For tubing connection to opening 19 - gas inlet. <i>Not relevant for CVS analyses.</i></p>
<p>9 Drive shaft for rotating disk electrode (RDE)</p>	<p>10 M6 threaded opening (N₂) With preinstalled stopper. <i>Not relevant for CVS analyses.</i></p>
<p>11 M6 threaded opening (OUT) With preinstalled tubing connection to opening 18 - gas outlet. <i>Not relevant for CVS analyses.</i></p>	<p>12 Guide roller Transfers the rotary movement of the motor to the driving axle of the working electrode.</p>
<p>13 M6 threaded opening (TAP) With preinstalled tubing connection to threaded opening 17 - gas inlet. <i>Not relevant for CVS analyses.</i></p>	<p>14 - 30: see next figure</p>




Figure 5 Measuring head insert

<p>14 Opening For feeding through a four-way micro dosing tip (6.1824.000) from below.</p>	<p>15 Opening For positioning the driving axle.</p>
<p>16 Threaded opening With preinstalled screw nipple and stopper. Can be equipped with a four-way micro dosing tip (6.1824.000).</p>	<p>17 M6 threaded opening With preinstalled tubing connection to the threaded opening 13 (TAP) - gas inlet. <i>Not relevant for CVS analyses.</i></p>
<p>18 Opening With preinstalled tubing connection to threaded opening 11 (OUT) - gas outlet. <i>Not relevant for CVS analyses.</i></p>	<p>19 Opening For tubing connection to threaded opening 8 (PURGE) - gas inlet. <i>Not relevant for CVS analyses.</i></p>

20 Opening for electrode
For inserting the reference electrode (RE).

22 Pipetting opening
For manually dosing solutions. Is closed with a 6.2709.100 stopper (3-5).

24 Opening
For adding or aspirating solutions. Can be connected with threaded opening **1, 2, 3** or **WASTE** (FEP tubing from 6.1829.070).

26 Opening
For adding or aspirating solutions. Can be connected with threaded opening **1, 2, 3** or **WASTE** (FEP tubing from 6.1829.070).

28 Opening
For the automated sample addition. Must be connected to threaded opening **5** (SMP) (PEEK capillary 6.1831.020).

30 Opening
For feeding through the temperature sensor cables from above.

21 Opening for electrode
For inserting the driving axle for the RDE (working electrode - WE).

23 Opening for electrode
For inserting the auxiliary electrode (AE).

25 Opening
For adding or aspirating solutions. Can be connected with threaded opening **1, 2, 3** or **WASTE** (FEP tubing from 6.1829.070).

27 Opening
For adding or aspirating solutions. Can be connected with threaded opening **1, 2, 3** or **WASTE** (FEP tubing from 6.1829.070).

29 Opening for sensor
For inserting a temperature sensor (Pt1000).

3.5 Tubing connector (measuring head arm)




Figure 6 Tubing connector (measuring head arm)

1 Nipple (N₂)
Is connected to threaded openings **TAP**, **PURGE** and **N₂** via the measuring head arm. *Not relevant for CVS analyses.*

2 M6 threaded opening (3)
For connecting tubing for aspirating or adding solutions.

3.5 Tubing connector (measuring head arm)



3 M6 threaded opening (2)

For connecting tubing for aspirating or adding solutions.

5 M6 threaded opening (1)

For connecting tubing for aspirating or adding solutions.

7 M6 threaded opening (OUT)

Not relevant for CVS analyses.

4 UNF 10/32 threaded opening (SMP)

For connecting a PEEK capillary for the automated sample addition.

6 M6 threaded opening (WASTE)

For connecting tubing for aspirating the measuring solution.

4 Installation

4.1 Setting up the device

4.1.1 Packaging

The product and accessories are supplied in protective special packaging. Keep this packaging to ensure safe transportation of the product. If a transport lock is present, keep this as well for future reuse.

4.1.2 Checks

Inspect the delivery immediately upon receipt:

- Check the delivery against the delivery note to ensure completeness.
- Check the product for damage.
- If the delivery is incomplete or damaged, contact your regional Metrohm representative.

4.1.3 Location

The instrument has been developed for operation indoors and may not be used in explosive environments.

Place the instrument in a location of the laboratory which is suitable for operation and free of vibrations and which provides protection against corrosive atmosphere and contamination by chemicals.

The instrument should be protected against excessive temperature fluctuations and direct sunlight. Do not place the instrument in the vicinity of an air-conditioning unit outlet.

4.2 Equipping the RDE measuring head

The RDE measuring head can be used for both VA trace analysis and CVS analysis. If the measuring head is used for the determination of organic additives with CVS, then chapter 3.3.1, step 3 and chapter 3.3.4 can be omitted.



CAUTION

The measuring head insert (3-1) is made of PTFE. Do not use sharp objects around it to ensure the material is not damaged.

4.2.1 Preparing the RDE measuring head



NOTE

For equipping, Metrohm recommends placing the RDE measuring head in the measuring head holder and only then inserting it onto the measuring head arm.

- 1 Remove the stopper (3-5) from the pipetting opening to remove the measuring head cover.




Figure 7 Removing the stopper from the pipetting opening

- 2 Pull the slide lock (3-4) on the top of the measuring head cover towards you and, at the same time, tilt the measuring head cover to an angle of approx. 45° and remove it.




Figure 8 Removing the measuring head cover

3 Connecting the gas inlet



NOTE

Only carry out this step if the RDE measuring head is used for VA trace analysis. If the RDE measuring head is used for CVS analysis, then no gas inlet is necessary.

- Insert the PTFE tubing for adding gas to the solution (6.1829.030) through the opening (5-19).
- Pull the transparent inner tubing through as far as it will go.
- Ensure that the green kink protection is protecting the entire piece of tubing.




Figure 9 Inserting the gas inlet

- Connect the tubing to the PURGE threaded opening (4-**8**) and tighten it hand-tight.
- Finally, tighten the tubing nipple using the wrench provided (6.2739.000).




Figure 10 Connecting the gas inlet

4.2.2 Preparing electrodes and inserting them in the RDE measuring head

The 894 Professional CVS uses the potentiostatic three-electrode principle. The following electrodes are used:

- Rotating disk electrode (RDE) as working electrode (WE)
- Reference electrode (RE)
- Auxiliary electrode (AE)



NOTE

Also observe the notes contained in the electrode leaflets that can be downloaded from the [Metrohm website](#) using the corresponding article number. In addition, you can learn how to best handle the electrodes from the multimedia guide (A.717.0003).

4.2.2.1 Working electrode (WE)

The working electrode consists of the following 2 articles:

- Electrode tip (e.g. 6.1204.610)
- Driving axle for rotating disk electrode (RDE) (e.g. 6.1204.510)

Preparing and inserting the working electrode

Proceed as follows:

1 Preparing the electrode tip

Remove the protective cap from the electrode tip.




Figure 11 Removing the protective cap from the electrode tip

2 Assembling the working electrode

Hold the driving wheel by the driving axle and tighten the electrode tip to the driving axle.




Figure 12 Tightening the electrode tip to the driving axle

3 Inserting the working electrode into the measuring head insert

Insert the working electrode into the opening (5-21) of the measuring head insert.

Make sure that the pin on the lower part of the driving axle is positioned in the opening (5-15) of the measuring head insert.




Figure 13 Working electrode, installed

4 Fastening the drive belt

- Slide the drive belt (6.1244.050) over the drive shaft (4-9),
- guide it on both sides over the guide roller (4-12) from below,
- pull it over the working electrode and fasten it in the driving wheel of the driving axle.




Figure 14 Fastening the drive belt



NOTE

Make sure that the drive belt does not rub against the driving axle or surrounding components (tubing, cables, etc.).

5 Connecting the working electrode

Plug the electrode cable (4-7) with the **WE** marking on the plug onto the metal contact of the driving axle.




Figure 15 Connecting the working electrode

4.2.2.2 Reference electrode (RE)

The reference electrode consists of the following 2 articles:

- Reference electrode filled with reference electrolyte (e.g. 6.0728.130)
- Electrolyte vessel filled with bridge electrolyte (e.g. 6.1245.010)

Preparing and inserting the reference electrode

Proceed as follows:

- 1 Remove the reference electrode from the storage vessel.
The reference electrode that is part of the accessories is already filled with reference electrolyte ($c(\text{KCl}) = 3 \text{ mol/L}$).
- 2 Fill the electrolyte vessel with bridge electrolyte (e.g. $c(\text{KNO}_3) = 1 \text{ mol/L}$) in accordance with the information in the electrode leaflet.

- 3 Allow the bridge electrolyte to react in the electrolyte vessel until the diaphragm is soaked with bridge electrolyte.
- 4 Place the reference electrode in the filled electrolyte vessel and screw it in place.

The electrolyte solution that is displaced in the electrolyte vessel is forced out of the deaeration openings.




Figure 16 Assembling the reference electrode with the electrolyte vessel

- 5 Rinse the installed reference electrode with ultrapure water.
- 6 Insert the installed reference electrode into the opening (5-20) of the measuring head insert.

7



CAUTION

The electrode cables for the reference and the auxiliary electrode look identical. Observe the markings on the plugs, because the 2 cables must not be mixed up.



Plug the electrode cable with the **RE** marking on the plug onto the metal contact of the reference electrode.




Figure 17 Connecting a reference electrode

4.2.2.3 Auxiliary electrode (AE)

The auxiliary electrode (e.g. 6.0343.100) can be placed directly in the measuring head.

Inserting the auxiliary electrode

Proceed as follows:

- 1 Insert the auxiliary electrode into the opening (5-23) of the measuring head insert.

2

**CAUTION**

The electrode cables for the reference and the auxiliary electrode look identical. Observe the markings on the plugs, because the 2 cables must not be mixed up.

Plug the electrode cable (4-7) with the **AE** marking on the plug onto the metal contact of the auxiliary electrode.




Figure 18 Connecting the auxiliary electrode

4.2.3 Inserting the RDE measuring head

Once the RDE measuring head has been fully equipped, it can be inserted on the measuring head arm. Proceed as follows:

1 Installing the measuring head cover

Place the measuring head cover at an angle of approx. 45° in the guide bolt in the front part of the measuring head insert.



Fold back the measuring head cover and gently push it in place.




Figure 19 Installing the measuring head cover

The measuring head cover must snap into place with an audible click.

2 Inserting the stopper

Insert the stopper (3-5) into the pipetting opening.




Figure 20 Inserting the stopper in the pipetting opening

3 Inserting the measuring head



WARNING

If the measuring head arm is folded down without due care, this may result in injuries to the hands.

Make sure that your fingers do not get caught between the measuring head arm and the device housing.



CAUTION

Do not apply pressure to the drive disk on the connector plate of the measuring head arm. Otherwise, the stirrer motor may be damaged.




Figure 21 Do not touch the drive disk

Use one hand to hold the measuring head arm on the rear and use the other hand to insert the measuring head into the connector plate of the measuring head arm.

The measuring head must snap into place with an audible click.




Figure 22 Inserting the measuring head

4.2.4 Connecting the inert gas supply



NOTE

Only connect the inert gas supply if the RDE measuring head is used for VA trace analysis.

If the RDE measuring head is used for CVS analysis, then no inert gas supply is necessary.

Generally, nitrogen (N_2) is used as inert gas in VA trace analysis for purging the measuring solution and for operating the MME pro, the SPE and the RDE. Only nitrogen of sufficient purity may be used for this.

For general polarography/voltammetry:

- 4.5 (w(N_2) = 99.995%)

For analyses in organic solvents; for determinations that result in very high current strengths (such as for determining the smallest concentrations without preceding deposition)

- 5.0 (w(N_2) = 99.999%)

1 Filling the gas washing glass

- Unscrew the gas washing glass (1-6) from the measuring head arm.

- Fill the gas washing glass as follows:
 - Standard: Fill the gas washing glass halfway with distilled H₂O.
 - For long-term measurement with base electrolytes like acetic acid/acetate buffer solution or ammonia/ammonium chloride buffer solution, add the base electrolyte.
 - For measurements in organic solvents, fill with the solvent being used.
- Screw the gas washing glass back onto the measuring head arm.




Figure 23 Gas washing glass

2 Connecting the inert gas inlet

- Connect one end of the PVC tubing (6.1801.080) on the **N₂** nipple of the 894 Professional CVS.
- Connect the other end of the PVC tubing (6.1801.080) to the inert gas bottle connection.
- Set the inert gas pressure on the gas bottle to $p = 1.0\text{--}1.2$ bar (or 14.5–17.4 psi or 0.1–0.12 MPa) using the reducing valve.
- Open the gas inlet on the gas bottle.




Figure 24 Nipple for inert gas supply

4.3 Establishing the tubing connections

Depending on the peripheral devices you are using with the measuring instrument (manual, semiautomated or automated determinations), different tubing connections have to be established. All tubing connections that can be used are described below.

4.3.1 Removing the measuring head cover

Remove the measuring head cover to establish the tubing connections.



CAUTION

The measuring head insert (3-1) is made of PTFE. Do not use sharp objects around it to ensure the material is not damaged.

Removing the measuring head cover

- 1 Remove the stopper (3-5) from the pipetting opening to remove the measuring head cover.




Figure 25 Removing the stopper from the pipetting opening

- 2 Pull the slide lock (3-4) on the top of the measuring head cover towards you and, at the same time, tilt the measuring head cover to an angle of approx. 45° and remove it.




Figure 26 Removing the measuring head cover

4.3.2 Installing the 4-way micro dosing tip

The 4-way micro dosing tip (6.1824.000) can be used to connect the 894 Professional CVS to dosing devices of the type 800 Dosino and to add auxiliary solutions and standard solutions automatically. Information on the electrical connection of dosing devices can be found in the manual.



Installing a 4-way micro dosing tip in the measuring head

Proceed as follows:

- 1** Remove the stopper from the screw nipple of the threaded opening (5-16).
- 2** Loosen the screw nipple in the threaded opening a little.
This slightly loosens the O-ring located on the bottom of the screw nipple.
- 3** Feed the 4-way micro dosing tip through the opening from below (5-14).
- 4** Insert the 4-way micro dosing tip into the screw nipple of the threaded opening (5-16) until it stops.




Figure 27 Inserting the 4-way micro dosing tip

- 5** Tighten the screw nipple in the threaded opening again.

Connecting the PTFE capillaries of the 4-way micro dosing tip to a dosing unit

Proceed as follows:

- 1 Screw the PTFE capillaries to the dosing units (port 1).




Figure 28 Screwing a PTFE capillary to a dosing unit

Sealing unused PTFE capillaries of the 4-way micro dosing tip

To prevent accidental aspiration of solution from the measuring vessel, unused PTFE capillaries must be sealed. Proceed as follows:

- 1 Screw a coupling (6.1808.000) to each unused PTFE capillary.
- 2 Screw a threaded stopper (6.1446.040) to each coupling (6.1808.000).




Figure 29 Sealing the PTFE capillary of a 4-way micro dosing tip

- 2 Slide the pressure screw over the capillary. Ensure that the capillary protrudes 1 to 2 mm from the tip of the pressure screw.
- 3 Push the capillary into the threaded opening **SMP (4-5)** of the measuring head connector plate until it stops.
- 4 Only then start turning the pressure screw, while holding the capillary firmly in place.
- 5 Cut the capillary to the desired length using the capillary cutter.
In order to prevent diffusion between the solution in the capillary and the solution in the measuring vessel, make sure that the end of the capillary is positioned above the measuring solution.

Installing the capillary on the tubing connector of the measuring head arm

We recommend using PTFE capillaries (e.g. 6.1803.020) between the tubing connector of the measuring head arm (*see chapter 3.5, page 13*) and the dosing unit and/or peristaltic pump. Proceed as follows:

- 1 Slide the pressure screw over the capillary. Ensure that the capillary protrudes 1 to 2 mm from the tip of the pressure screw.
- 2 Push the capillary into the threaded opening **SMP (4-5)** of the tubing connector of the measuring head arm as far as it will go.
- 3 Only then start turning the pressure screw, while holding the capillary firmly in place.

Connecting a capillary to a dosing unit

To attach the capillary to the dosing unit using a pressure screw, you need an adapter for the thread on the dosing unit. Proceed as follows:

- 1 Screw the adapter (6.2744.080) to port 2 of the dosing unit.
- 2 Slide the pressure screw (e.g. 6.2744.010) over the capillary (e.g. 6.1803.020). Ensure that the capillary protrudes 1 to 2 mm from the tip of the pressure screw.



- 3 Push the capillary into the threaded opening of the adapter until it stops.
- 4 Only then start turning the pressure screw, while holding the capillary firmly in place.




Figure 31 Screwing a capillary to a dosing unit

Connecting a capillary to the peristaltic pump

Proceed as follows:

- 1 Install the peristaltic pump with all tubing connections (including capillaries).

Proceed in accordance with the information contained in the manual of the Sample Processor.

4.3.4 Installing FEP tubing

With FEP tubing, you can add or aspirate the following solutions:

- Rinsing solution – adding via 843 Pump Station
- Waste solution – aspirating via 843 Pump Station
- Auxiliary solutions (e.g. buffer, electrolyte, VMS, etc.) – adding via 800 Dosino with dosing unit

You can find several illustrations of the tubing setup with FEP tubing in the chapter *Automated operation* in the detailed manual.

Installing FEP tubing in the measuring head



NOTE

Use the appropriate FEP tubing from the tubing set (6.1829.070): For threaded openings **1** and **WASTE**, use the tubing with the long kink protection; and for threaded openings **2** and **3**, use the tubing with the short kink protection.

Proceed as follows:

- 1** Insert the conical end of the FEP tubing into one of the openings (**24** - **27**) in the measuring head insert.
- 2** Carefully pull the conical end of the tubing protruding on the bottom side of the measuring head insert through the opening using the sandpaper supplied.

Pull the tubing through the opening until the kink protection of the tubing sits flush on the opening of the measuring head insert.
- 3** Screw the tubing nipple of the FEP tubing into the required threaded opening (**1**, **2**, **3** or **WASTE**) of the measuring head connector plate.
- 4** Cut the FEP tubing on the underside of the measuring head insert to the required length.

In order to prevent diffusion between the solution in the tubing and the solution in the measuring vessel, make sure that the end of the tubing is positioned above the measuring solution. This does not apply to the tubing for draining the measuring solution; this tubing must touch the base of the measuring vessel.




Figure 32 Inserting the FEP tubing into the measuring head



NOTE

If more than 1 piece of tubing is to be installed, we recommend pulling in all of the pieces of tubing into the measuring head insert first and shortening them to the required length.

In order to simplify fastening the tubing nipples to the threaded openings, remove the measuring head from the device and screw the tubing nipples in place from bottom to top (in the order **3**, **2**, **WASTE**, **1**).

Installing the measuring head cover

The measuring head cover must be on the measuring head during determinations. Proceed as follows:




Figure 33 Installing the measuring head cover

- 1 Place the measuring head cover at an angle of approx. 45° in the guide bolt in the front part of the measuring head insert.

2



WARNING

If the measuring head cover is installed without due care, this may result in injury to the hand.

Make sure that your fingers do not get caught between the measuring head cover and the measuring head.

Fold back the measuring head cover and gently push it in place.

The measuring head cover must snap into place with an audible click.

- 3 Insert the stopper (3-5) into the pipetting opening.




Figure 34 Inserting the stopper in the pipetting opening

Installing a piece of FEP tubing on the tubing connector of the measuring head arm

Proceed as follows:

- 1 Screw the tubing nipple of the FEP tubing (e.g. 6.1805.530) into the required threaded opening (**1**, **2**, **3** or **WASTE**) of the tubing connector of the measuring head arm.

An illustration of this tubing connection can be found in the chapter *Automated operation* in the detailed manual.

Connecting FEP tubing to an 843 Pump Station

Proceed as follows:

- 1 Install the inlet and outlet tubing on the two pumps or the 843 Pump Station [\[Link target not found in publication context!\]](#).

Connecting FEP tubing to a rinsing canister

A bottle cap (6.1602.115) is required to connect a piece of FEP tubing to a rinsing canister. A figure illustrating the tubing setup between the rinsing canister, the 843 Pump Station and the 894 Professional CVS can be found in the chapter *Automated operation* in the detailed manual.




Figure 35 Installing a bottle cap with pieces of tubing on a rinsing canister

Proceed as follows:

- 1 Insert a piece of FEP tubing (e.g. 6.1819.110) from above into the smallest threaded opening of the bottle cap (6.1602.115) until the tubing flare of the FEP tubing sits flush in the threaded opening.
- 2 Screw the bottle cap with inserted FEP tubing onto the larger of the two openings of the rinsing canister.
- 3 Screw the tubing nipple of the FEP tubing (e.g. 6.1805.530) into the smallest threaded opening of the bottle cap.

Connecting FEP tubing to the waste canister

A 5-way tubing connector (6.1828.020) is required to connect FEP tubing to a waste canister. A figure illustrating the tubing setup between the waste canister, the 843 Pump Station and the 894 Professional CVS can be found in the chapter *Automated operation* of the detailed manual.




Figure 36 Connecting the five-way tubing connector to a waste canister

Proceed as follows:

- 1 Unscrew the cap with the hole in the center from the waste canister.
- 2 Feed each of the 5 tubings from the 5-way tubing connector with the couplings already installed from below through the hole of the unscrewed cap.
- 3 Screw the cap with the 5-way tubing connector inserted back onto the waste canister.
- 4 Screw the tubing nipples of the pieces of FEP tubing (e.g. 6.1805.530) to the couplings of the pieces of tubing.



NOTE

In order to ensure that liquid is transported smoothly from and to the canisters, the canisters must not be air-tight. If necessary, loosen the screw caps a little.

Connecting FEP tubing to a dosing unit with auxiliary solution

Proceed as follows:




Figure 37 Connecting FEP tubing to the dosing unit (auxiliary solution)

- 1 Screw the tubing nipple of the FEP tubing (e.g. 6.1805.120) to port 1 of the dosing unit containing the auxiliary solution (e.g. buffer, electrolyte, VMS, etc.).

4.4 Connecting devices electrically

4.4.1 Connecting the instrument to the power grid



WARNING

Electric shock from electrical potential

Risk of injury by touching live components or through moisture on live parts.

- Never open the housing of the instrument while the power cord is still connected.
- Protect live parts (e.g. power supply unit, power cord, connection sockets) against moisture.
- Unplug the power plug immediately if you suspect that moisture has gotten inside the instrument.
- Only personnel who have been issued Metrohm qualifications may perform service and repair work on electrical and electronic parts.

Connecting the power cord

Accessories

Power cord with the following specifications:

- Length: max. 2 m
- Number of cores: 3, with protective conductor

- Instrument plug: IEC 60320 type C13
- Conductor cross-section 3x min. 1.0 mm² / 18 AWG
- Power plug:
 - according to customer requirement (6.2122.XX0)
 - min. 10 A

**NOTE**

Do not use a not permitted power cord!

1 Plugging in the power cord

- Plug the power cord into the instrument's power socket.
- Connect the power cord to the power grid.

4.4.2 Connecting the 894 Professional CVS

The 894 Professional CVS is connected to the computer with the supplied controller cable.

Connecting the computer

- 1** Connect the controller cable (6.2151.000) to the "Controller" connector of the 894 Professional CVS.

**NOTE**


The plug on the controller cable is protected against accidental disconnection by means of a pull-out protection feature. If you wish to pull out the plug, you first need to pull back the outer plug sleeve.

- 2** Connect the USB plug of the controller cable to an available USB connector on the computer.

Initializing the 894 Professional CVS in viva


- 1** Start **viva**.

The following dialog window is displayed:



2 Click on **Yes**.

The following dialog window is displayed:



3 Change the suggested device name if required.

4 Confirm with **OK**.

The instrument will be automatically listed in the device table of the **Configuration** program part.

4.4.3 Connecting an 800 Dosino

Up to 4 dosing devices of the type 800 Dosino can be connected directly to the 894 Professional CVS. As an alternative, the dosing devices can be connected via an 846 Dosing Interface, a sample changer or any other supported instrument that also has MSB outputs.



NOTE

An 800 Dosino is connected to Metrohm instruments via an MSB connector. Make sure that the flat side of the plug matches the marking on the socket.




Figure 38 Connecting a dosing device to an MSB socket




NOTE

Information regarding assembling a Dosino with the dosing unit as well as operation and maintenance is available in the manual of the 800 Dosino.

Connecting an 800 Dosino with dosing unit directly to the 894 Professional CVS

- 1 Connect the connection cable of the 800 Dosino to one of the 4 MSB connectors (2-8) of the 894 Professional CVS.

The following dialog window is displayed:




- 2 Confirm with **OK**.

Initializing a dosing unit in viva

- 1 Select the 894 Professional CVS in the device table of the **Configuration** section of the program.


- 2 In the device table, click on the **Edit** button and select **Initialize**.

If a brand-new dosing unit is connected, the following dialog window is displayed:



Or:

If a dosing unit is connected that has been previously configured, then the following dialog window is displayed:



- 3 Click on **OK** if a brand-new dosing unit is used.

The following dialog window is displayed:

Dosing unit - GLP

Hardware

Name

Comment

Device name / dosing device

Order number

Serial number

Cylinder volume mL

Cylinder serial number

Parameters for preparation

Dosing port Prep/Empty

Dosing rate Dosing port 1 mL/min

Dosing rate Dosing port 2 mL/min

Dosing rate Fill port mL/min

Dosing rate Special port mL/min

Tubing parameters

	Port	Length	Diameter
Dosing port 1	<input type="text" value="Port 1"/>	<input type="text" value="80.0"/> cm	<input type="text" value="0.3"/> mm
Dosing port 2	<input type="text" value="Port 3"/>	<input type="text" value="0.0"/> cm	<input type="text" value="2.0"/> mm
Fill port	<input type="text" value="Port 2"/>	<input type="text" value="25.0"/> cm	<input type="text" value="2.0"/> mm
Special port	<input type="text" value="Port 4"/>	<input type="text" value="0.0"/> cm	<input type="text" value="2.0"/> mm

Valve disk

Rotating direction

Not over

OK Cancel

- 4 Configure the new dosing unit in this dialog window.



NOTE

Under **Tubing parameters**, adapt the parameters **Length** and **Diameter** to the actual installation. This ensures that functions such as **Prepare** or **Empty** work properly in **viva**. The lengths and diameters of the tubing connections incorporated in the measuring head arm can be found in the appendix of the detailed manual.

- 5 Click on **Yes** if you are using the already configured dosing unit.

The dosing unit will be automatically displayed in the **Dosing units** subwindow of the **Configuration** program part.



NOTE

Check the **Tubing parameters** and adapt them to the actual installation, if necessary.

4.4.4 Connecting a Sample Processor

The following Sample Processors can be used in conjunction with the 894 Professional CVS:

- 858 Professional Sample Processor
- 919 IC Autosampler plus
- 814 USB Sample Processor
- 815 Robotic USB Sample Processor XL

Sample Processors are connected to the power supply with a power supply cable and to the PC with a controller cable.



WARNING

An incorrect supply voltage can damage the instrument. Operate the instrument only with the supply voltage specified for it.

In order to ensure a constant power supply, we recommend using an uninterruptible power supply (UPS).

Connecting to power supply and PC

Proceed as follows:

- 1 Connect the Sample Processor to the power supply with the power supply cable (6.2122.0x0).
- 2 Connect the controller cable (6.2151.000) to the "Controller" connector of the Sample Processor.

**NOTE**


The plug on the controller cable is protected against accidental disconnection by means of a pull-out protection feature. If you wish to pull out the plug, you first need to pull back the outer plug sleeve.

- 3 Connect the USB plug of the controller cable to an available USB connector on the PC.

**NOTE**

We recommend connecting the Sample Processor directly to the PC and not to one of the two USB connections on the rear of the 894 Professional CVS.

The following dialog window is displayed:



- 4 Click on **Yes**.

The following dialog window is displayed:

Properties - 858.0020 Professional Sample Processor - 858 Professional Sa...

General | GLP

Device name: 858 Professional Sample Processor 1

Device type: 858.0020 Professional Sample Processor

Program version: 5.858.0011

Device serial number: 03693

Set to work: 2013-03-28 14:23:57 UTC+1

Remarks:

OK Cancel

5 Change the suggested instrument name if required.

6 Confirm with **OK**.

The instrument will be automatically listed in the device table of the **Configuration** section of the program.

4.4.5 Connecting an external pump

If solutions are to be added or aspirated via an external pump, we recommend using the 843 Pump Station with membrane pump. The 843 Pump Station is connected to the tower of the Sample Processor with the 6.2141.300 cable.

Connecting an 843 Pump Station

Proceed as follows:

- 1** Connect the 843 Pump Station to the power grid with the power cord (6.2122.0x0).
- 2** Insert the remote plug of the cable (6.2141.300) into the "Remote 2" socket of the 843 Pump Station.




Figure 39 Connecting an external pump using the remote cable

- 3 Connect the 2 plugs of the cable (6.2141.300) for the pump connectors to the tower of the Sample Processor.

Make sure that the markings on the cables match the markings on the tower ("Ext. Pump 1" and "Ext. Pump 2").




Figure 40 Connecting an external pump to the tower

4.4.6 Connecting USB devices directly to the measuring device

Additional USB devices can be incorporated using the USB interfaces on the measuring device. You can, for example, connect the following USB devices:

- Barcode reader
- Keyboard
- Mouse
- USB hub



NOTE

Use only USB connectors on the PC to connect Metrohm devices with a USB connector.

5 Start-up

The 894 Professional CVS is operated exclusively via the **viva** computer software. You can find information on operating **viva** in the online help and in the CVS tutorial.

Proceed as follows for the initial start-up of the 894 Professional CVS:

Preparing the system for start-up



WARNING

Uncontrolled splashing of reagents

Splashing reagents may result in injuries.

Operate the 894 Professional CVS only with the measuring head in place and the measuring head arm lowered.



WARNING

Drive belt in operation

Hair may become entangled or caught in the running drive belt, for example.

Operate the 894 Professional CVS only with the measuring head cover in place.

- 1** Insert the measuring head as described in *"Equipping the RDE measuring head"*, page 15.
- 2** Equip the measuring head with the electrodes (see chapter 4.2.2, page 19).
- 3** Establish the tubing connections in the measuring head and those between the 894 Professional CVS and the peripheral devices (see chapter 4.3, page 30).
- 4** Attach the measuring head cover and the stopper (3-5) to the measuring head as described in *"Installing the measuring head cover"*, page 38.

5 Insert the measuring vessel into the holder (1-3).

6



CAUTION

If the measuring head arm is folded down without due care, this may result in injuries to the hands.

Make sure that your fingers do not get caught between the measuring head arm and the device housing.

Lower the measuring head arm.

7 Place the drip pan (6.2711.090) into the drip pan holder (1-8).

8 Connect the devices electrically (see chapter 4.4, page 43).

Activating the calibrator


The calibrator built into the 894 Professional CVS guarantees a very high measuring accuracy for 3 years after the initial start-up. You can find more detailed information in the chapter *Calibrator* of the detailed manual.

In **viva**, proceed as follows:

1 Select the 894 Professional CVS in the device table of the **Configuration** section of the program.

2 In the device table, click on the **Edit** button and select **Properties...**

3 Open the **Calibrator** tab.



- 4 Click on **Activate calibrator** and then on **OK**.

The set-to-work date of the calibrator and the calibrator's serial number are displayed on the **Calibrator** tab after the dialog window has been closed and opened again.

Properties - 894 Professional CVS - 894_1

General MSB 1 MSB 2 MSB 3 MSB 4 **Calibrator** GLP

Calibrator serial number 01126

Set to work 2013-02-13 13:49:14 UTC+1

Remarks

Monitor calibrator

Exchange interval 1095 days

Next check 2016-02-13 ...

Message

Message by e-mail E-mail...

Acoustic signal

Action

Record message

Display message

Cancel determination

OK Cancel



NOTE

The calibration certificate is available online. Go to <http://www.metrohm.com/com/Support/zertifikate/instrument-Certificate.html> and enter the serial number of the calibrator and download the certificate.

- 5 If required, adjust the settings on the **Calibrator** tab (e.g. **Message** or **Action**).

Problem	Cause	Remedy
	<i>The measuring vessel contains too little solution or is empty.</i>	Check the fill level in the measuring vessel and add solution if necessary.
The electrode tip of the working electrode and/or the auxiliary electrode shows copper deposits.	<i>An incorrect potential has been applied because the reference electrode is defective.</i>	<ul style="list-style-type: none"> ▪ Check the reference electrode; perform maintenance procedure in accordance with the electrode leaflet, if necessary. ▪ Dip the electrode tip of the working electrode and/or the auxiliary electrode into concentrated nitric acid for 1 to 2 seconds and then thoroughly rinse with distilled water.
	<i>An incorrect potential was applied because the reference electrode is not connected.</i>	<ul style="list-style-type: none"> ▪ Check the electrode connectors. ▪ Dip the electrode tip of the working electrode and/or the auxiliary electrode into concentrated nitric acid for 1 to 2 seconds and then thoroughly rinse with distilled water.
	<i>An incorrect potential was applied because the measuring vessel contains too little measuring solution.</i>	<ul style="list-style-type: none"> ▪ All three electrodes must be immersed in the measuring solution. ▪ Dip the electrode tip of the working electrode and/or the auxiliary electrode into concentrated nitric acid for 1 to 2 seconds and then thoroughly rinse with distilled water.
The measured values are widely scattered.	<i>The solutions are pipetted manually.</i>	Use an 800 Dosino with dosing unit.
	<i>The solutions are not added via the pipetting opening.</i>	Add solutions only via the pipetting opening (5-22).
The measurement curves are noisy.	<i>Contact problem at the driving axle (6.1204.510 or 6.1204.520).</i>	<ul style="list-style-type: none"> ▪ Remove the abrasion residue from the driving axle. ▪ Replace the driving axle.
The measuring signal for Cu VMS fluctuates.	<i>The ambient and/or solution temperature is not constant.</i>	Keep the ambient and solution temperature stable during measurements (± 4 °C).
	<i>The reference potential drifts off.</i>	Perform the maintenance procedure for the reference electrode in accordance with the electrode leaflet. However, for CVS analyses,



Problem	Cause	Remedy
		replace the reference electrolyte every other day and the bridge electrolyte every day.
	<i>The measuring vessel and the electrodes have been contaminated with organic additives.</i>	Thoroughly rinse the measuring vessel and the electrodes and use fresh VMS.
The measuring vessel overflows.	<i>Incorrect pump times have been defined in viva.</i>	Adjust the pump times.
	<i>The volumes defined in the dosing commands in viva are too large.</i>	Reduce the volumes.
The peak is no longer recognized.	<i>The reference potential has shifted.</i>	Perform the maintenance procedure for the reference electrode in accordance with the electrode leaflet.
The signal does not decrease in spite of suppressor addition.	<i>No chloride is contained in the Cu VMS.</i>	Check the preparation of the Cu VMS and modify it if required.

6.2 Peripheral devices


Problem	Cause	Remedy
The 800 Dosino cannot be actuated by the 894 Professional CVS.	<i>The connection between the 800 Dosino and the 894 Professional CVS is either interrupted or an error has occurred on the 800 Dosino.</i>	<ul style="list-style-type: none"> ▪ Check the cable connections. ▪ Disconnect the 894 Professional CVS from the power grid and connect it again. ▪ Check the dosing and filling rate. ▪ Contact your regional Metrohm service representative if necessary.
The data of the dosing unit cannot be read.	<i>The memory chip of the dosing unit is mechanically damaged or impaired by chemicals.</i>	<ul style="list-style-type: none"> ▪ Remove the dosing drive and attach it again. ▪ Clean the memory chip and the contact surfaces. ▪ Have the memory chip replaced by your regional Metrohm service representative.
The dosing unit is blocked and/or leaking.	<i>Crystals have formed (in the cylinder, on the valve disk or in the capillary).</i>	<ul style="list-style-type: none"> ▪ Check the flow path. ▪ Rinse the dosing unit and the connected tubing and capillaries (Prepare function) when the measuring system is not in use.

Problem	Cause	Remedy
		<ul style="list-style-type: none"> ▪ Clean the dosing unit at least every two weeks.
The dosing unit is recognized either not at all or incorrectly.	<i>The dosing drive was not attached correctly.</i>	<ul style="list-style-type: none"> ▪ Remove the dosing drive and attach it again. ▪ Check whether the dosing drive is correctly seated. ▪ Disconnect the 894 Professional CVS from the power grid and connect it again. ▪ Contact your regional Metrohm service representative if necessary.
The membrane pumps of the 843 Pump Station do not operate at full pump capacity.	<i>The cables are not or not correctly connected.</i>	Connect the cables in accordance with the instructions in the comprehensive manual.
	<i>The tubing connections are leaking.</i>	Check the tubing connections and tighten, if necessary.
	<i>The rinse and/or waste canister are sealed airtight.</i>	Loosen the lids on the canisters a little or remove them.
The pump time of the peristaltic pump increases.	<i>The pump tubing of the peristaltic pump has aged or is defective.</i>	Replace the pump tubing.
The sample is not completely transferred from the Sample Processor to the measuring vessel via the peristaltic pump.	<i>The PEEK sample needle on the Sample Processor is positioned more than 0.5 mm from the base of the sample vial.</i>	Position the PEEK sample needle in accordance with the instructions in the detailed manual.
	<i>The contact pressure set for the tubing cartridge is insufficient.</i>	Set the contact pressure of the tubing cartridge in accordance with the information in the Sample Processor manual.
	<i>The selected pump times are too short.</i>	Prolong the pump times.

7 Displaying accessories

Up-to-date information on the scope of delivery and on optional accessories can be found on the Metrohm website.

1 Searching for a product on the website

- Go to <https://www.metrohm.com>.
- Click on .
- Enter the article number of the product (e.g. **2.1001.0010**) into the search field and press **[Enter]**.

The search result is displayed.

2 Displaying product information

- To display the products matching the search term, click on **Product models**.
- Click on the desired product.

Detailed information regarding the product is displayed.

3 Displaying accessories and downloading the accessories list

- To display the accessories, scroll down to **Accessories and more**.
 - The **scope of delivery** is displayed.
 - Click on **[Optional parts]** for the optional accessories.
- To download the accessories list, click on **[Download accessories PDF]** under **Accessories and more**.



NOTE

Metrohm recommends keeping the accessories list for reference purposes.

Index

- 4-way micro dosing tip
 Connect 33
 Install 32
- 800 Dosino
 Connect (electrically) 46
 Troubleshooting 60
- 843 Pump Station
 Connect (electrically) 51
 Troubleshooting 61
 Tubing setup 40
- 894 Professional CVS
 Start-up 54
 Troubleshooting 58
- A**
 Auxiliary electrode
 Connect 25
- C**
 Calibration certificate 57
 Calibrator
 Activate 55
 Certificate 57
 Capillary
 Install 34
 Connect
 Connect to power grid 43
- D**
 Dosing unit
 Connect 4-way micro dosing
 tip 33
 Connect capillary 35
 Connect tubing 42
 Initialize 47
 Troubleshooting 60, 61, 60
- Drive belt
 Fasten 21
- E**
 Electrical connection
 800 Dosino 46
- 843 Pump Station 51
 Instrument 44
 Sample Processor 49
- F**
 FEP tubing
 See also "Tubing" 36
 Front of the instrument 6
- I**
 Initialization
 Dosing unit 47
 Instrument 44
 Instrument
 Connect (electrically) 44
 Initialize 44
- M**
 Measuring head
 Equip 15
 Insert 27
 Parts 9
 Measuring head arm
 Drive disk 27
 Tubing connector 13
 Measuring head connector plate
 Detail view 10
 Measuring head cover
 Install 38
 Remove 16, 31
 Measuring head insert
 Detail view 12
 MSB device
 Connect 46
- P**
 Peristaltic pump
 Troubleshooting 61
 Power connection 43, 44
- R**
 RDE
 See also "Working electrode"
 20
 RDE measuring head
 See also "Measuring head" ... 9
 Rear of the instrument 7
 Reference electrode
 Connect 24
 Prepare 22
 Rinsing canister
 Tubing setup 40
- S**
 Sample changer
 See also "Sample Processor"
 49
 Sample Processor
 Connect (electrically) 49
 Standby potential
 Blinking pattern 6
 Start-up
 Prepare 54
- T**
 Troubleshooting 58
 Tubing
 Install 36
 Tubing setup 36
 843 Pump Station 40
 Rinsing canister 40
 Waste canister 41
- W**
 Waste canister
 Tubing setup 41
 Working electrode
 Connect 21
 Insert 20
 Prepare 20