

901 Titrando



Manual

8.901.8003EN / v9 / 2026-01-09



Metrohm AG
CH-9100 Herisau
Switzerland
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

901 Titrando

Manual

Technical Communication
Metrohm AG
CH-9100 Herisau

This documentation is protected by copyright. All rights reserved.

This documentation is an original document.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Disclaimer

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g., conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

Table of contents

1	Introduction	1
1.1	The Titrand system	1
1.2	Device description	2
1.3	Titration modes – Measuring modes – Dosing commands	3
1.4	Displaying accessories	4
1.5	Symbols and conventions	5
2	Safety	6
2.1	Intended use	6
2.2	Responsibility of the operator	6
2.3	Personnel requirement	7
2.4	Safety instructions	7
2.4.1	Electrical safety	7
2.4.2	Tubing and capillary connections	8
2.4.3	Flammable solvents and chemicals	8
3	Overview of the instrument	9
4	Installation	11
4.1	Setting up the device	11
4.1.1	Packaging	11
4.1.2	Checks	11
4.1.3	Location	11
4.2	Connecting a controller	11
4.2.1	Operation	11
4.3	Connecting MSB devices	16
4.3.1	Connecting a dosing device	17
4.3.2	Connecting a stirrer or titration stand	18
4.3.3	Connecting a Remote Box	19
4.4	Connecting USB devices	20
4.4.1	General	20
4.4.2	Connecting a USB hub	21
4.4.3	Connecting a printer	21
4.4.4	Connecting a balance	22
4.4.5	Connecting a keyboard (only for operation with Touch Control)	23
4.4.6	Connecting a barcode reader	24



4.5	Setting up the titration vessel	25
4.5.1	General	25
4.5.2	Titration vessel for volumetric KF titration	26
4.6	Connecting sensors	30
4.6.1	Connecting a pH, metal or ion-selective electrode	30
4.6.2	Connecting a reference electrode	30
4.6.3	Connecting a polarizable electrode	31
4.6.4	Connecting a temperature sensor or an electrode with integrated temperature sensor	31
4.6.5	Connecting the iConnect	32
4.6.6	Differential potentiometry	35
5	Karl Fischer titration	36
5.1	Volumetric titration	36
5.1.1	Principle of the volumetric Karl Fischer titration	36
5.1.2	Endpoint determination	36
5.1.3	Karl Fischer reagents	37
5.1.4	Application of the Karl Fischer titration	37
5.1.5	Working with water standards	37
5.1.6	Sample addition	39
5.1.7	Optimum working conditions	41
6	Operation and maintenance	43
6.1	General notes	43
6.1.1	Care	43
6.1.2	Maintenance by Metrohm Service	43
7	Troubleshooting	44
7.1	General	44
7.2	Karl Fischer titration	44
7.2.1		44
7.3	SET titration	46
7.3.1		46
8	Appendix	47
8.1	Remote interface	47
8.1.1	Pin assignment of the remote interface	47
9	Technical specifications	51
9.1	Measuring interface	51
9.1.1	Potentiometry	51
9.1.2	Temperature	51
9.1.3	Polarizer	52
9.2	Power connection	53
9.3	Ambient conditions	53



9.4	Reference conditions	53
9.5	Dimensions	54
9.6	Interfaces	54
	Index	55

1 Introduction

1.1 The Titrando system

The Titrando is the heart of the modular Titrando system. Operation is carried out either by a Touch Control with a touch-sensitive screen ("stand-alone titrator") or by a computer with a corresponding software.

A Titrando system can contain numerous, various kinds of devices. The following figure provides an overview of the peripheral devices you can connect to the 901 Titrando.




Figure 1 The Titrando system



Up to 3 control instruments (Titrandos, Dosing Interface, USB Sample Processor, etc.) can be controlled via USB connection during operation with the 900 Touch Control.

You can request information on special applications in the "Application Bulletins" and "Application Notes", available free of charge via the regional Metrohm representative. Various monographs on the subjects of titration techniques and electrodes are also available.

Updating the device software is described in the help of the corresponding PC software.

1.2 Device description

The 901 Titrandos has the following characteristics:

- **Operation**
Operation is carried out by means of a touch-sensitive Touch Control or with a high-performance PC software.
- **MSB connectors**
4 MSB connectors (Metrohm Serial Bus) for connecting dosing devices (Dosimat or Dosino), stirrers, titration stands, and Remote Boxes.
- **USB connectors**
Two USB connectors, through which devices such as printers, PC keyboards, barcode readers or additional control devices (USB Sample Processor, Titrandos, Dosing Interface, etc.) can be connected.
- **Measuring interface**
One measuring input each for:
 - a potentiometric electrode (pH, metal or ion-selective electrode)
 - a separate reference electrode
 - a temperature sensor (Pt1000 or NTC)
 - a polarizable electrode
 - an iConnect (measuring interface for electrodes with integrated data chip, so-called iTrodes)








1.3 Titration modes – Measuring modes – Dosing commands

The 901 Titrande supports the following titration modes, measuring modes and dosing commands:

- **SET**
Endpoint titration at one or two specified endpoints.
Measuring modes:
 - **pH** (pH measurement)
 - **U** (potentiometric voltage measurement)
 - **Ip_{ol}** (voltametric measurement with selectable polarization current)
 - **Up_{ol}** (amperometric measurement with selectable polarization voltage)
- **KFT**
Volumetric water content determination according to Karl Fischer.
Measuring modes:
 - **Ip_{ol}** (voltametric measurement with selectable polarization current)
 - **Up_{ol}** (amperometric measurement with selectable polarization voltage)
- **MEAS**
The following measuring modes can be selected for measurements:
 - **pH** (pH measurement)
 - **U** (potentiometric voltage measurement)
 - **Ip_{ol}** (voltametric measurement with selectable polarization current)
 - **Up_{ol}** (amperometric measurement with selectable polarization voltage)
 - **T** (temperature measurement)
- **CAL**
Electrode calibration.
Measuring mode:
 - **pH** (calibration of pH electrodes)
- **ELT**
Electrode test for pH electrodes.
This mode is listed separately only in *tiamo*[™]. In Touch Control, the electrode test is a component part of the CAL calibration mode.

1.5 Symbols and conventions

The following symbols and formatting may appear in this documentation:

(5-12)	Cross-reference to figure legend
	The first number refers to the figure number, the second to the device part in the figure.
1	Instruction step
	Perform the steps one after the other.
Method	Dialog text, parameter in the software
File ► New	Menu or menu item
[Continue]	Button or key
	WARNING
	This symbol draws attention to a possible life-threatening hazard or risk of injury.
	WARNING
	This symbol draws attention to a possible hazard due to electrical current.
	WARNING
	This symbol draws attention to a possible hazard due to heat or hot instrument parts.
	WARNING
	This symbol draws attention to a possible biological hazard.
	WARNING
	Warning of optical radiation
	CAUTION
	This symbol draws attention to possible damage to devices or device parts.
	NOTICE
	This symbol highlights additional information and tips.

2.3 Personnel requirement

Only qualified personnel may operate the product. Qualified personnel are persons who meet the following requirements:

- Knowledge of basic regulations on occupational safety and accident prevention.
- Knowledge regarding the application of fire prevention measures for laboratories.
- In-depth knowledge of handling hazardous chemicals.
- Personnel have been instructed and are able to operate the product safely and independently recognize and avoid potential dangers.
- The user documentation has been read and understood. The personnel operate the product according to the instructions in the user documentation.

2.4 Safety instructions

2.4.1 Electrical safety

Electrical safety when working with the device is ensured as part of the international standard IEC 61010.



WARNING

Only personnel qualified by Metrohm are authorized to carry out service work on electronic components.



WARNING

Never open the housing of the device. The device could become damaged. There is a considerable risk of injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.

Supply voltage



WARNING

An incorrect supply voltage can damage the device.

Operate this device only with a supply voltage specified for it (refer to the rear of the device).

Protection against electrostatic charges



WARNING

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Do not fail to pull the power cord out of the power socket before setting up or disconnecting electrical plug connections at the rear of the device.

2.4.2 Tubing and capillary connections



CAUTION

Leaks in tubing connections and capillary connections are a safety risk. Tighten all connections well by hand. Avoid applying excessive force to tubing connections. Damaged tubing ends lead to leakage. Suitable tools can be used for disconnecting connections.

The leak-tightness of the connections must be checked regularly. If the device is used mainly in unattended operation, then weekly inspections are mandatory.

2.4.3 Flammable solvents and chemicals



WARNING

All relevant safety measures are to be observed when working with flammable solvents and chemicals.

- Set up the device in a well-ventilated location (e.g., fume cupboard).
- Keep all sources of ignition far from the workplace.
- Clean up spilled liquids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.

3 Overview of the instrument




Figure 2 Front 901 Titrando

1 Bottle holder

With holding clamps, for two reagent bottles.

2 "On" LED

Lights up when the Titrando is ready for operation.




Figure 3 Rear 901 Titrando

<p>1 Type plate Contains specifications concerning supply voltage, instrument type and serial number.</p>	<p>2 USB connector (USB 1 and USB 2) USB ports (type A) for connecting printer, keyboard, barcode reader, additional Titrandos, USB Sample Processor, etc.</p>
<p>3 Connector (Controller) For connecting a Touch Control or a PC with installed PC software. Mini DIN, 9-pin.</p>	<p>4 MSB connector (MSB 1 to MSB 4) Metrohm Serial Bus. For connecting external dosing devices, stirrers or Remote Boxes. Mini DIN, 9-pin.</p>
<p>5 Power socket</p>	<p>6 Measuring interface 1 (Input 1)</p>
<p>7 Electrode connector (Ind.) For connecting pH, metal or ion-selective electrodes with integrated or separated reference electrode. Socket F.</p>	<p>8 Electrode connector (Ref.) For connecting reference electrodes, e.g. Ag/AgCl reference electrode. Socket B, 4 mm.</p>
<p>9 Temperature sensor connector (Temp.) For connecting temperature sensors (Pt1000 or NTC). Two B sockets, 2 mm.</p>	<p>10 Electrode connector (Pol.) For connecting polarizable electrodes, e.g. double Pt wire electrodes. Socket F.</p>
<p>11 Electrode connector (iConnect) For connecting electrodes with integrated data chip (iTrodes).</p>	

4 Installation

4.1 Setting up the device

4.1.1 Packaging

The product and accessories are supplied in protective special packaging. Keep this packaging to ensure safe transportation of the product. If a transport lock is present, keep this as well for future reuse.

4.1.2 Checks

Inspect the delivery immediately upon receipt:

- Check the delivery against the delivery note to ensure completeness.
- Check the product for damage.
- If the delivery is incomplete or damaged, contact your regional Metrohm representative.

4.1.3 Location

The instrument has been developed for operation indoors and may not be used in explosive environments.

Place the instrument in a location of the laboratory which is suitable for operation and free of vibrations and which provides protection against corrosive atmosphere and contamination by chemicals.

The instrument should be protected against excessive temperature fluctuations and direct sunlight.

4.2 Connecting a controller

4.2.1 Operation

Two different versions are available for operating the 901 Titrando:

- A Touch Control with touch-sensitive screen. It forms a "stand-alone instrument" together with the 901 Titrando.
- A computer enables operation of the 901 Titrando with the help of a PC software, e.g. *tiamo*.



CAUTION

Take care to ensure that the power cord is pulled out of the power socket before either setting up or disconnecting connections between the instruments.

If you would prefer not to position the Touch Control directly next to the Titrando, then you can lengthen the connection with the 6.2151.010 cable. The maximum connection length permitted is 5 m.

4.2.1.2 Connecting the instrument to the power grid



WARNING

Electric shock from electrical potential

Risk of injury by touching live components or through moisture on live parts.

- Never open the housing of the instrument while the power cord is still connected.
- Protect live parts (e.g. power supply unit, power cord, connection sockets) against moisture.
- Unplug the power plug immediately if you suspect that moisture has gotten inside the instrument.
- Only personnel who have been issued Metrohm qualifications may perform service and repair work on electrical and electronic parts.

Connecting the power cord

Accessories

Power cord with the following specifications:

- Length: max. 2 m
- Number of cores: 3, with protective conductor
- Instrument plug: IEC 60320 type C13
- Conductor cross-section 3x min. 0.75 mm² / 18 AWG
- Power plug:
 - according to customer requirement (6.2122.XX0)
 - min. 10 A



NOTE

Do not use a not permitted power cord!

1 Plugging in the power cord

- Plug the power cord into the instrument's power socket.
- Connect the power cord to the power grid.

4.2.1.3 Connecting a computer

The 901 Titrando requires a USB connection to a computer in order to be able to be controlled by a computer software. Using a 6.2151.000 controller cable, the instrument can be connected directly, either to a USB port on a computer, to a connected USB hub or to a different Metrohm control instrument.

You need administrator rights for the installation of driver software and control software on your computer.

Cable connection and driver installation

A driver installation is required in order to ensure that the 901 Titrando is recognized by the computer software. To accomplish this, you must comply with the procedures specified. The following steps are necessary:

1 Installing the software

- Insert the installation CD of the computer software and carry out the instructions of the installation program.
- Exit the program if you have started it after the installation.

2 Establishing the cable connections

- Connect all peripheral devices to the instrument, *see chapter 4.3, page 16*.
- Connect the instrument to the power grid if you have not already done this (*see chapter 4.2.1.2, page 13*).
The "On" LED on the 901 Titrando is not yet illuminated!
- Connect the instrument to a USB connector (type A) of your computer (see manual of your computer). The 6.2151.000 cable is used for this purpose.




Figure 5 Connecting the computer

The instrument is recognized. Depending on the version of the Windows operating system used, the driver installation proceeds dif-

ferently afterwards. Either the necessary driver software is installed automatically or an installation wizard is started.

3 Follow the instructions of the installation wizard.


The "On" LED on the 901 Titrande lights up when the driver installation has been completed and the instrument is ready for operation.

If problems should occur during installation, contact your company's IT support team.



NOTE

The plug on the instrument end of the 6.2151.000 controller cable is protected against accidental disconnection by means of a pull-out protection feature. If you wish to pull out the plug, you first need to pull back the outer plug sleeve marked with arrows.



Registering and configuring the instrument in the computer software

The instrument must be registered in the configuration of your computer software. Once that has been done, you can then configure the instrument according to your requirements.

1 Setting up the instrument

- Start up the computer software.
The instrument is automatically recognized. The configuration dialog for the instrument is displayed.
- Make configuration settings for the instrument and its connectors.

More detailed information concerning the instrument configuration can be found in the documentation for the respective computer software.

4.3 Connecting MSB devices

In order to connect MSB devices, e.g., stirrers or dosing devices, Metrohm devices are equipped with up to a maximum of 4 connectors on what is referred to as the *Metrohm Serial Bus* (MSB). Various kinds of peripheral devices can be connected in sequence (in series, as a "daisy chain") at a single MSB connector (8-pin Mini DIN socket) and controlled simultaneously by the respective control instrument. In addition to the connection cable, stirrers and the Remote Box are each equipped with their own MSB socket for this purpose.

The following figure provides an overview of the instruments that can be connected to an MSB socket, along with a number of different cabling variations.




Figure 6 MSB connections

The control instrument determines which peripheral devices are supported.



NOTE

When connecting MSB devices together, the following must be observed:

- Only one device of the same type can be used at a single MSB connector at one time.
- Dosing devices of the 700 Dosino and 685 Dosimat plus type cannot be connected together with other MSB devices on a shared connector. These dosing devices must be connected separately.

**CAUTION**

Exit the control software before you plug in MSB devices. When it is switched on, the control instrument automatically recognizes which device is connected to which MSB connector. The operating unit or the control software enters the connected MSB devices into the system configuration (device manager).

MSB connectors can be extended with the 6.2151.010 cable. The maximum connection length permitted is 15 m.

4.3.1 Connecting a dosing device

Four dosing devices can be connected to the instrument (**MSB 1 to MSB 4**).

The types of dosing devices that are supported are:

- 800 Dosino
- 700 Dosino
- 805 Dosimat
- 685 Dosimat plus

Connecting a dosing device

1 Connecting a dosing device

- Exit the control software.
- Connect the connection cable of the dosing device to 1 of **MSB** sockets on the rear of the control instrument.
- Start the control software.




Figure 7 Connecting a dosing device

4.3.2 Connecting a stirrer or titration stand

You can use the following devices:

These devices have a built-in magnetic stirrer (stirring "from below"):

- 801 Stirrer
- 803 Ti Stand

This device has no built-in magnetic stirrer (stirring "from above"):

- 804 Ti Stand with propeller stirrer 802 Stirrer

Connecting the stirrer or titration stand

- 1 Exit the control software.
- 2 Connect the connection cable of the magnetic stirrer or titration stand to an **MSB** socket on the rear of the control instrument.
804 Ti Stand only: Connect the propeller stirrer to the stirrer connector (socket with stirrer symbol) of the titration stand.
- 3 Start the control software.




Figure 8 Connecting an MSB stirrer




Figure 9 Connecting the propeller stirrer to the titration stand

4.3.3 Connecting a Remote Box

Instruments that are controlled via remote lines and/or that send control signals via remote lines can be connected via the 6.2148.010 Remote Box. In addition to Metrohm, other instrument manufacturers also use similar connectors that make it possible to connect different instruments together. These interfaces are also frequently given the designations "TTL Logic", "I/O Control" or "Relay Control" and they generally have a signal level of 5 volts.

You can connect the following devices to the remote connector, among others:

- 849 Level Control (fill level monitoring in a canister)
- 731 Relay Box (switch box for 230/110 volt alternating current sockets and low-voltage direct current outlets)
- 843 Pump Station (for complex sample preparations or for rinsing external titration vessels)

The Remote Box also has an MSB socket at which a further MSB device, e.g., a dosing device or a stirrer, can be connected. You will find pre-

cise information concerning the pin assignment of the interface on the Remote Box in the *Appendix*.

Control signals are understood to be electrical line statuses or electrical pulses (> 200 ms) which display the operating status of an instrument or which trigger or report an event. Runs on a variety of instruments and can thus be coordinated in a single complex automation system. However, no data exchange is possible.

Connecting the Remote Box

- 1 Exit the control software.
- 2 Connect the connection cable of the Remote Box to an **MSB** socket on the rear of the control instrument.
- 3 Start the control software.




Figure 10 Connecting the Remote Box

4.4 Connecting USB devices

4.4.1 General

The 901 Titrando has 2 USB connectors (type A sockets) for peripheral devices with USB interfaces. The Titrando functions as a USB hub (distributor), no matter how it is operated. If you wish to connect more than 2 devices to the USB, you can also use an additional, commercially available USB hub.

**CAUTION**

If you operate the 901 Titrandò with the aid of the Touch Control, then take care to ensure that the Touch Control is switched off when you set up or disconnect connections between the various devices. If you use a computer software to control the 901 Titrandò, then you should exit the program before you set up or disconnect the USB connections.

4.4.2 Connecting a USB hub

If you wish to connect more than 2 devices to the USB connector of the 901 Titrandò, then you can also use an additional commercially available USB hub (distributor). If you operate the 901 Titrandò with the help of the Touch Control, then you should use a USB hub with its own energy supply.

- 1 Switch off the Touch Control and/or exit the computer software.
- 2 With the aid of the 6.2151.020 cable, connect the USB connector of the 901 Titrandò (type A) with the USB connector of the hub (type B, see manual for the hub).
- 3 Switch on the Touch Control.

The USB hub is recognized automatically.

4.4.3 Connecting a printer

Printers that are connected to the 901 Titrandò with Touch Control must meet the following requirements:

- Printer languages: HP-PCL (PCL 3 to 5, PCL 3GUI), Canon BJJ Commands or Epson ESC P/2
- Printer resolution: 300 dots/inch or 360 dots/inch (Epson)
- Paper size: A4 or Letter, single-sheet feed.

- 1 Switch off the Touch Control.
- 2 With the aid of the 6.2151.020 cable, connect the USB connector of the 901 Titrandò (type A) with the USB connector of the printer (type B, see manual for the printer).
- 3 Switch on the printer first, then the Touch Control.



- 4 Configure the printer in the device manager of the Touch Control (see Touch Control manual).




Figure 11 Connecting a printer

4.4.4 Connecting a balance

- Operation with a computer software:
 - Connect the balance directly to the serial connector (COM) of the computer. This is usually 9-pin and marked with the symbol **IOIOI**.
- Operation with Touch Control:
 - You will need the 6.2148.050 USB/RS-232 adapter to connect a balance.

The following table offers an overview of the balances that you can use together with the 901 Titrando and of which cable you will need for connection to the RS-232 interface:

Balance	Cables
AND ER, FR, FX with RS-232 interface (OP-03)	6.2125.020 + 6.2125.010
Mettler AB, AG, PR (LC-RS9)	In the scope of delivery for the balance
Mettler AM, PM, PE with interface option 016 or Mettler AJ, PJ with interface option 018	6.2146.020 + 6.2125.010 Also from Mettler: ME 47473 adapter and either ME 42500 hand switch or ME 46278 foot switch
Mettler AT	6.2146.020 + 6.2125.010 Also from Mettler: ME 42500 hand switch or ME 46278 foot switch
Mettler AX, MX, UMX, PG, AB-S, PB-S, XP, XS	6.2134.120

Balance	Cables
Mettler AE with interface option 011 or 012	6.2125.020 + 6.2125.010 Also from Mettler: ME 42500 hand switch or ME 46278 foot switch
Ohaus Voyager, Explorer, Analytical Plus	Cable AS017-09 from Ohaus
Precisa balances with RS-232-C interface	6.2125.080 + 6.2125.010
Sartorius MP8, MC, LA, Genius, Cubis	6.2134.060
Shimadzu BX, BW	6.2125.080 + 6.2125.010

Operation with Touch Control

- 1 Plug in the USB plug of the USB/RS-232 adapter at the USB connector of the 901 Titrando.
- 2 Connect the RS-232 interface of the USB/RS-232 adapter with the RS-232 interface of the balance (see table for cable).
- 3 Switch on the Touch Control.
- 4 Switch on the balance.
- 5 Activate the RS-232 interface of the balance if necessary.
- 6 Configure the RS-232 interface of the USB/RS-232 adapter in the device manager of the Touch Control (see Touch Control manual).

4.4.5 Connecting a keyboard (only for operation with Touch Control)

The computer keyboard is used as an aid for text and numerical input.

- 1 Insert the USB plug of the keyboard into 1 of the USB port of the 901 Titrando.
- 2 Switch on the Touch Control.

The keyboard is recognized automatically and entered in the device manager.



- 3 Configure the keyboard in the device manager of the Touch Control (see Touch Control manual).

4.4.6 Connecting a barcode reader

The barcode reader is used as an aid for text and numerical input. You can connect a barcode reader with USB interface.

Operation with Touch Control

- 1 Insert the USB plug of the barcode reader into 1 of the USB ports of the 901 Titrando.
- 2 Switch on the Touch Control.
The barcode reader is recognized automatically and entered in the device manager.
- 3 Configure the barcode reader in the device manager of the Touch Control (see Touch Control manual).

Settings on the barcode reader:

- 1 Switch the barcode reader to programming mode.
- 2 Specify the desired layout for the keyboard (USA, Germany, France, Spain, German-speaking Switzerland).
This setting must match the setting in the device manager (see the Touch Control manual).
- 3 Make sure that the barcode reader is set in such a way that Ctrl characters (ASCII 00 to 31) can be sent.
- 4 Program the barcode reader in such a way that the ASCII character 02 (STX or Ctrl B) is sent as the first character. This first character is normally referred to as the "Preamble" or "Prefix Code".
- 5 Program the barcode reader in such a way that the ASCII character 04 (EOT or Ctrl D) is sent as the last character. This last character is normally referred to as the "Postamble", "Record Suffix" or "Postfix Code".
- 6 Exit the programming mode.

4.5 Setting up the titration vessel

4.5.1 General

During the titration, it is important that the solution is mixed well. The stirring rate should be high enough to form a small vortex. If the stirring rate is too high, then air bubbles will be aspirated. This results in incorrect measured values. If the stirring rate is too low, then the solution at the electrode will not be mixed correctly. In order for the measurement to be taken in a well-mixed solution after the addition of the titrant, the tip should be positioned where turbulence is high. Furthermore, the distance between the addition of the titrant and the electrode should be as large as possible. Therefore, take into account the stirring direction (counterclockwise or clockwise) when positioning the electrode and tip (see figure below).




Figure 12 Diagrams showing magnetic stirrer, electrode, and tip during a titration. a) clockwise stirring direction, b) counterclockwise stirring direction.


1 Magnetic stirrer

2 Electrode


3 Tip

4.5.2 Titration vessel for volumetric KF titration

Mounting the KF titration cell



- 1** Screw the 6.2013.010 clamping ring tightly to the support rod.
- 2** Fix the 6.1414.030 vessel lid of the KF titration cell (with correctly inserted sealing ring from the 6.1244.040 sealing set) to the support rod. Keep the locking lever pressed down until it can be released at the desired position.
- 3** Fasten the 6.1415.220 (or 6.1415.250) titration vessel with a 6.1903.020 (or 6.1903.030) stirring bar inside on the vessel lid. Fold back the support bracket upwards while doing so. The markings on the vessel lid and on the plastic ring must be aligned above one another. Afterwards, press the support bracket downwards in order to fasten the titration vessel in place. The levers of the support bracket must enclose the pins of the plastic ring on the titration vessel in order to ensure a secure hold.




- 4** Adjust the height of the KF titration cell by pressing the locking lever. It should almost touch the surface of the stirrer. The position can now be fixed by readjusting the clamping ring.

If the height of the KF titration cell has been adjusted correctly, then the entire cell can now be raised and swiveled as required by pressing the locking lever.


Filling the adsorber tube

Before insertion, the 6.1403.040 adsorber tube must be filled with the 6.2811.000 / 6.2811.010 molecular sieve.




- 1** Insert a small cotton plug into the bottom of the adsorber tube. Do not pack the cotton too tightly.
- 2** Fill the molecular sieve up to $\frac{3}{4}$ of the height.
- 3** Place a small cotton plug on the molecular sieve. Do not pack the cotton too tightly.
- 4** Seal the adsorber tube with the associated lid.

Inserting the tips



- 1** Place the 3 screw nipples from 6.2730.030 (including O-rings, but without stoppers) in the rear openings of the vessel lid.
- 2** Insert the 6.1543.110 dosing tip through the screw nipple in the middle rear opening.
- 3** Insert the 6.1543.120 aspiration tip through the screw nipple in the right rear opening.
If solvent is aspirated, then the end of the aspiration tip must touch the vessel base, but it must not inhibit the action of the stirring bar.
The aspiration tip can, if needed, be pulled out of the solvent.
- 4** Insert the 6.1543.060 tip through the screw nipple in the left rear opening.
- 5** Screw the 6.1805.200 PTFE M8 tubing of the aspiration bottle onto the aspiration tip.
- 6** Screw the 6.1805.200 PTFE M8 tubing of the solvent bottle onto the dosing tip.

Inserting electrode, adsorber tube and septum stoppers



- 1** Introduce the 6.2730.020 septum stopper (with septum inserted) into the front opening of the vessel lid.
- 2** Insert the O-rings of the electrode and of the adsorber tube into the middle openings of the vessel lid.
- 3** Screw the 2 6.2730.010 screw nipples into the openings with the O-rings. Do not screw too tightly.
- 4** Introduce the 6.0338.100 double Pt electrode into the left-hand opening and then tighten the screw nipple until it seals.
- 5** Screw the 6.2104.020 electrode cable tightly onto the electrode.
- 6** Insert the filled 6.1403.040 adsorber tube to the right of the electrode into the remaining opening and then tighten the screw nipple until it seals.




Figure 14 Connecting a reference electrode

4.6.3 Connecting a polarizable electrode

- 1 Plug the electrode plug into the **Pol.** socket of the 901 Titrande.




Figure 15 Connecting a polarizable electrode



NOTE

The electrode cable is protected against accidental disconnection of the cable by means of a pull-out protection. If you wish to pull out the plug again, you first need to pull back the outer plug sleeve.

4.6.4 Connecting a temperature sensor or an electrode with integrated temperature sensor

A temperature sensor of the Pt1000 or NTC type can be connected to the **Temp.** connector.

- 1 Plug the temperature sensor plug into the **Temp.** sockets of the Titrande.




Figure 16 Connecting a temperature sensor or an electrode with integrated temperature sensor



NOTE

Always insert the red plug into the red socket. This is the only way that shielding against electrical interference can be ensured.

4.6.5 Connecting the iConnect

The external measuring interface 854 iConnect can be connected to the **iConnect** connector.

Connecting the 6.2168.000 mini USB adapter cable to an instrument

- 1 Plug the mini USB adapter cable (2) into the **iConnect** connector of the instrument (1). Observe the correct orientation (markings).




Figure 17 Plugging in the mini USB adapter cable

Depending on the instrument's model version, there may be several **iConnect** connectors.

- 2 Leave the adapter cable plugged in to protect the connector inside the instrument (1) against mechanical impact.

Connecting the 854 iConnect to the adapter cable

Check whether the mini USB adapter cable is connected to the instrument.

The 854 iConnect can also be connected while the instrument is switched on.

- 1 Plug the plug of the 854 iConnect (3) into the socket of the mini USB adapter cable (2). Observe the correct orientation (markings).




Figure 18 Plugging in the 854 iConnect

As soon as the instrument is switched on, the 854 iConnect is detected automatically and entered as measuring input into the device properties.

Connecting the electrode

The 854 iConnect is used as measuring input for iTrodes (electrodes with integrated memory chip).

- 1 Remove the protective cap of the 854 iConnect.




Figure 19 Removing the protective cap

- 2 Align the guide pin (5) of the electrode with the recess in the 854 iConnect (4).




Figure 20 *Aligning the guide pin*

- 3** Attach the electrode to the 854 iConnect.




Figure 21 *Attaching the electrode*

The guide pin guarantees correct connection in such a way that the contact pins cannot be damaged.

- 4** Tighten the screw cap by hand.

If there is an electrode in the sensor list of the firmware or software, the electrode is detected automatically when connecting it.



NOTE

Mount the protective cap again as soon as the 854 iConnect is not in use anymore and no electrode is connected.

Removing the 854 iConnect

The 854 iConnect can also be removed while the instrument is switched on.

- 1** Unplug the 854 iConnect (3) from the socket of the mini USB adapter cable (2).




Figure 22 Unplugging the 854 iConnect

- 2** Leave the mini USB adapter cable (2) plugged into the socket of the instrument (1).



NOTE

Position the mini USB adapter cable in such a way that it cannot be removed inadvertently.



NOTE

For additional information on the 854 iConnect, see *Manual iConnect 854*, article number 8.854.8002.

4.6.6 Differential potentiometry

Potentiometric measurements with high-ohm measuring chains can be disrupted by electrostatic and electromagnetic fields in media with low conductivity. Use our 6.0229.100 Solvotrode or other special electrodes for pH measurements in organic solvents. If no reliable measurements are possible with these, then a 6.5104.030 (230 V) or 6.5104.040 (115 V) differential amplifier can be used. The differential amplifier is connected to the high-ohm measuring input (**Ind.**).

5.1.3 Karl Fischer reagents

1-component reagents

They contain all the reactive parts in a single solution – iodine, sulfur dioxide, and a base, dissolved in a suitable alcohol.

Two-component reagents

The reactive parts are distributed among 2 separated solutions. The titration reagent contains iodine in methanol. The KF solvent is a solution of sulfur dioxide and a base in methanol. It is used as a working medium in the KF titration cell.

5.1.4 Application of the Karl Fischer titration

The volumetric Karl Fischer titration is preferably used for the determination of an amount of water between 0.1 and 100%. It has the advantage, that even solid and pasty samples can directly be added to the titration vessel. In addition, various organic solvents can be used which fit the corresponding samples.

5.1.5 Working with water standards

5.1.5.1 Certified water standards

Commercially available, certified water standards with a water content of 10.0 ± 0.1 mg/g should be used for validating the instrument as a whole, integrated system.

5.1.5.2 Practical recommendations

For validation, it is essential to work very accurately. In order to minimize any measurement inaccuracies that could occur, the sample preparation and the sample processing should proceed in accordance with a defined scheme:

- 1 Put on gloves (always for Karl Fischer titration).
- 2 Use a clean syringe.
- 3 Take a new ampoule of water standard and shake it briefly.
- 4 With a folded paper towel held between thumb and index finger, break open the ampoule at the marking.
- 5 Draw approx. 1 mL of the water standard into the syringe.



- 6** Pull the plunger of the syringe up to the end and shake the syringe back and forth somewhat.

The inside of the syringe is rinsed by water standard and freed of water contamination.

- 7** Dispose of the used water standard in a waste bottle.

- 8** Draw the rest of the water standard into the syringe, aspirating as little air as possible.

- 9** Push out any air bubbles that may be present in the syringe.

- 10** Wipe off the needle with a lint-free paper towel and cover it with the appropriate cap.

- 11** Place the syringe on the balance and press **[TARA]**.

- 12** As soon as the drift on the 901 Titrande is stable, take the syringe in your hand, press **[START]** and inject approx. 1 mL of the water standard through the septum.

There are two possibilities:

- Version 1:
Inject the water standard without immersing the needle in the reagent liquid. If a little drop remains on the end of the needle, then it must be aspirated back before pulling the needle out of the septum.
The water standard should not be sprayed from the syringe onto the electrode nor onto the wall of the titration vessel.
- Version 2:
Inject the water standard directly under the surface of the reagent liquid.
Take care to ensure that you do not aspirate any liquid at the time you withdraw the syringe from the reagent liquid.

- 13** Close the syringe with the same cap and place it back on the balance.

- 14** Read off the value displayed by the balance and enter it as sample size on the Touch Control or in the computer software (e.g. *tiamo*).

- 15** The next determination can be started as soon as the determination has been finished and the titration cell has been conditioned (drift stable) again.

5.1.6 Sample addition

This chapter contains a few notes concerning sample addition. Further notes can be found in the publications of the reagent manufacturers and in the Karl Fischer monograph published by Metrohm.

5.1.6.1 Size of the sample size

The sample weight should be small in order to be able to titrate as many samples as possible in the same electrolyte solution and thus to keep the titration time short. However, ensure that the sample contains at least 50 µg of H₂O. The following tables provide clues for the sample size.

Table 1 Approximate sample size in grams (5 mL cylinder volume)

Water content of the sample	KF reagent 1	KF reagent 2	KF reagent 5
0.5%	0.1–0.9	0.2–1.8	0.5–4.5
1.0%	0.05–0.45	0.1–0.9	0.25–2.25
5.0%		0.02–0.18	0.05–0.45
10.0%			0.03–0.22
25.0%			
50.0%			

Table 2 Approximate sample size in grams (10 mL cylinder volume)

Water content of the sample	KF reagent 1	KF reagent 2	KF reagent 5
0.5%	0.2–1.8	0.4–3.6	
1.0%	0.1–0.9	0.2–1.8	0.5–4.5
5.0%	0.02–0.18	0.04–0.36	0.1–0.9
10.0%		0.02–0.18	0.05–0.45
25.0%			0.02–0.18
50.0%			0.02–0.09

Table 3 Approximate sample size in grams (20 mL cylinder volume)

Water content of the sample	KF reagent 1	KF reagent 2	KF reagent 5
0.5%	0.4–3.6		
1.0%	0.2–1.8	0.4–3.6	
5.0%	0.04–0.36	0.08–0.72	0.2–1.8
10.0%	0.02–0.18	0.04–0.36	0.1–0.9
25.0%		0.02–0.14	0.04–0.36
50.0%			0.02–0.18

KF reagent 1: 1 mL KF reagent reacts with around 1 mg H₂O

KF reagent 2: 1 mL KF reagent reacts with around 2 mg H₂O

KF reagent 5: 1 mL KF reagent reacts with around 5 mg H₂O

5.1.6.2 Working with liquid samples

Liquid samples are added with a syringe. The samples can be injected two different ways:

- Use a syringe with a long needle that you immerse in the reagent during the injection.
- Use a syringe with a short needle and aspirate the last drop back into the needle.

The best way for you to determine the injected sample amount is to weigh the sample by difference.

Glass syringes should be used for the **determination of traces and validations**. We recommend obtaining these from a specialized syringe manufacturer.

Highly volatile samples and samples of low viscosity should be cooled before sampling. Doing so avoids losses while working. The syringe must, however, not be cooled directly, as condensation could be formed. For the same reason, no air may be aspirated into a syringe into which a cooled sample has been aspirated beforehand.

Samples of high viscosity can be thinned by heating. The syringe must be heated as well. The same target can be reached by diluting with a suitable solvent. In this case, the water content of the solvent has to be determined and subtracted as a blank value.

Viscous samples can be added to the measuring cell with a syringe without needle. You can use the ground-joint opening for this. The best way for you to determine the added sample amount is by weighing by difference the sample.

5.1.6.3 Working with solid samples

If possible, solid samples are to be extracted or dissolved in a suitable solvent. The resulting solution is injected, during which a blank-value correction for the solvent must be carried out.

If no suitable solvent can be found for a solid sample, or if the sample reacts with the Karl Fischer reagent, then a Karl Fischer oven should be used.

If solid samples have to be directly added to the titration cell, then they can be inserted through the ground-joint opening. While doing so, take care to ensure that

- the sample releases its moisture completely.
- no side reaction with the Karl Fischer reagent takes place.
- the surface of the electrode is not covered by the sample substance (incomplete KF reaction!).
- the Pt wires of the indicator electrode do not become damaged.

5.1.7 Optimum working conditions

5.1.7.1 Drift

A constant drift in the range of $\leq 10 \mu\text{L}/\text{min}$ is all right. Lower values are, however, quite possible. Higher but stable values will still produce good results because it is possible to compensate for the drift.

A constantly high drift can be caused by water-containing deposits in inaccessible parts of the titration vessel. In these cases, shaking the titration vessel can reduce the value. Make sure that there are no drops above the liquid level in the titration vessel.

When you work with a Karl Fischer oven, a drift $\leq 10 \mu\text{L}/\text{min}$ is all right. The drift depends on the gas flow (the smaller the gas flow, the lower the drift).

5.1.7.2 Reagent replacement

The electrolyte solution must be replaced in the following cases:

- The titration vessel is too full.
- The drift is too high, and cannot be reduced by shaking the titration vessel.

Exhausted electrolyte solution is best disposed of by aspiration. To do this, you can use, for example, an *803 Ti Stand* with built-in membrane pump. An advantage is that the titration vessel does not have to be disassembled.

In the event of severe contamination, the titration vessel can be rinsed with a suitable solvent which is also aspirated.



5.1.7.3 Indicator electrode

A new indicator electrode may require a certain warm-up time for surface formation. During this time, unexpectedly long titration times and excessively high measurement results can occur. This phenomenon will, however, disappear after a short time of use. In order to accelerate the setting of a new indicator electrode, the 901 Titrande can be conditioned (e.g. over night).

A contaminated indicator electrode can be carefully cleaned with an abrasive cleaner (6.2802.000 polishing set or toothpaste). After the cleaning, rinse with ethanol.

The two Pt wires of the indicator electrode should run as parallel as possible to one another. Check the Pt wires before inserting the electrode.

6 Operation and maintenance

6.1 General notes

6.1.1 Care

The 901 Titrande requires appropriate care. Excess contamination of the instrument may result in functional disruptions and a reduction in the service life of the otherwise sturdy mechanics and electronics.

Spilled chemicals and solvents should be removed immediately. Above all, the plug connections on the rear of the instrument (in particular the power socket) should be protected from contamination.



CAUTION

Although this is extensively prevented by design measures, the power plug should be unplugged immediately if aggressive media have penetrated the inside of the instrument, so as to avoid serious damage to the instrument electronics. In such cases, Metrohm Service must be informed.

6.1.2 Maintenance by Metrohm Service

Maintenance of the 901 Titrande is best carried out as part of an annual service, which is performed by specialist personnel of the Metrohm company. A shorter maintenance interval may be necessary if you frequently work with caustic and corrosive chemicals.

Metrohm Service offers every form of technical advice for maintenance and service of all Metrohm instruments.

7 Troubleshooting

7.1 General

Problem	Cause	Remedy
The "On" LED is not illuminated, even though the instrument is connected to the power supply.	<i>The Touch Control or the computer has not been switched on yet or the plugs are not correctly plugged in.</i>	<ol style="list-style-type: none"> 1. Check the plug connections. 2. Switch on the Touch Control or the computer.

7.2 Karl Fischer titration

Problem	Cause	Remedy
The drift is very high during conditioning.	<i>The titration cell is leaking.</i>	<ul style="list-style-type: none"> ▪ Check the seals and the septum. Replace if necessary. ▪ Replace the molecular sieve.
The drift becomes greater after each titration.	<i>The sample releases water very slowly.</i>	<ul style="list-style-type: none"> ▪ Adjust the method. ▪ Add solubility promoter. ▪ Increase the temperature (possibly using a KF oven). ▪ See technical literature.
	<i>A side reaction is taking place.</i>	<ul style="list-style-type: none"> ▪ Use special reagents. ▪ Adjust the method (increase/decrease the temperature, external extraction). ▪ See technical literature.
	<i>The pH value is no longer in the optimum range.</i>	Add buffer (see technical literature).
The titration will not be finished.	<i>The titration cell is leaking.</i>	<ul style="list-style-type: none"> ▪ Check the seals and the septum. Replace if necessary. ▪ Replace the molecular sieve.
	<i>The minimum increment is too low.</i>	Select the user-defined titration rate and increase the minimum volume increment (see manual/help of the software used).
	<i>The stop criterion is unsuitable.</i>	Adjust the control parameters (see manual/help of the software used): <ul style="list-style-type: none"> ▪ Increase the stop drift.

Problem	Cause	Remedy
		<ul style="list-style-type: none"> Select a short delay time.
	<i>See also: The drift becomes greater after each titration.</i>	
The sample is over-titrated.	<i>The increments at the end of the titration are too high.</i>	<ul style="list-style-type: none"> Select the user-defined titration rate and reduce the dosing rate (see manual/help of the software used). The following experiment provides a clue for the optimum dosing rate: During conditioning, display the drift and add sample without starting the titration. Select a value below the highest drift as dosing rate. Stir faster.
	<i>The amount of methanol in the working medium is too low.</i>	<ul style="list-style-type: none"> Replace the working medium. Reduce the amount of solubility promoter, if working with solvent mixtures, see technical literature.
	<i>The electrode may be covered.</i>	Wipe off the electrode with ethanol or a suitable solvent.
The solution becomes darker after each titration.		Replace the working medium.
	<i>The electrode may be covered.</i>	Wipe off the electrode with ethanol or a suitable solvent.
	<i>The electrode has a short circuit.</i>	<ol style="list-style-type: none"> Check the Pt wires. Activate the electrode check.
The endpoint is reached too quickly.	<i>The dosing rate outside the control range is too high.</i>	Select the user-defined titration rate and reduce the dosing rate (see manual/help of the software used).
The titration times with volumetric titration are constantly longer.	<i>The buffer capacity of the solvent may be exhausted with two-component reagents.</i>	Replace the working medium.



7.3 SET titration

Problem	Cause	Remedy
The titration will not be finished.	<i>The minimum dosing rate is too low.</i>	Select the user-defined titration rate and increase the minimum rate (see manual/help of the software used).
	<i>The stop criterion is unsuitable.</i>	Adjust the control parameters (see manual/help of the software used): <ul style="list-style-type: none"> ▪ Increase the stop drift. ▪ Select a short delay time.
The sample is over-titrated.	<i>The control parameters are unsuitable.</i>	Adjust the control parameters (see manual/help of the software used): <ul style="list-style-type: none"> ▪ Select Titration rate = slow. ▪ Select the user-defined titration rate and increase the control range. ▪ Select the user-defined titration rate and reduce the maximum rate. ▪ Select the user-defined titration rate and reduce the minimum rate. ▪ Stir faster. ▪ Arrange the electrode and tip to an optimum.
	<i>The electrode responds too slowly.</i>	Replace the electrode.
	<i>The control parameters are unsuitable.</i>	Adjust the control parameters (see manual/help of the software used): <ul style="list-style-type: none"> ▪ Select Titration rate = optimal or fast. ▪ Select the user-defined titration rate and decrease the control range. ▪ Select the user-defined titration rate and increase the maximum rate. ▪ Select the user-defined titration rate and increase the minimum rate.
The titration time is too long.	<i>The control parameters are unsuitable.</i>	Adjust the control parameters (see manual/help of the software used): <ul style="list-style-type: none"> ▪ Select Titration rate = optimal or fast. ▪ Select the user-defined titration rate and decrease the control range. ▪ Select the user-defined titration rate and increase the maximum rate. ▪ Select the user-defined titration rate and increase the minimum rate.
	<i>The electrode responds too slowly.</i>	Replace the electrode.
The results are spread widely.	<i>The minimum dosing rate is too high.</i>	Select user-defined titration rate and decrease the minimum rate (see manual/help of the software used).
	<i>The electrode responds too slowly.</i>	Replace the electrode.

8 Appendix

8.1 Remote interface

The 6.2148.010 Remote Box allows devices to be controlled which cannot be connected directly to the MSB interface of the Titrande.




Figure 23 Connectors of the Remote Box

1 Cable

For connecting to the Titrande.

2 MSB connector

Metrohm Serial Bus. For connecting external dosing devices or stirrers.

3 Remote connector

For connecting instruments with a remote interface.

8.1.1 Pin assignment of the remote interface





Figure 24 Pin assignment of remote socket and remote plug

The above figure of the pin assignment applies for all Metrohm instruments with 25-pin D-Sub remote connector.

Inputs




approx. 50 k Ω Pull-up

$t_p > 20$ ms

active = low, inactive = high

Outputs



Open Collector

$t_p > 200$ ms

active = low, inactive = high

$I_C = 20$ mA, $V_{CE0} = 40$ V

+5 V: maximum load = 20 mA

The following tables offer information concerning the assignment of the individual pins and their function:

Table 4 Inputs and outputs of the remote interface

Assignment	Pin No.	Function*
Input 0	21	Start
Input 1	9	Stop
Input 2	22	
Input 3	10	Quit
Input 4	23	–
Input 5	11	
Input 6	24	
Input 7	12	
Output 0	5	Ready
Output 1	18	Conditioning OK
Output 2	4	Determination
Output 3	17	EOD
Output 4	3	
Output 5	16	Error
Output 6	1	
Output 7	2	Warning

Assignment	Pin No.	Function*
Output 8	6	
Output 9	7	
Output 10	8	
Output 11	13	
Output 12	19	
Output 13	20	
0 volts / GND	14	
+5 volts	15	
0 volts / GND	25	

* Signal activated only for operation with Touch Control.

Table 5 Explanation of the individual functions

Function	Explanation
Start	The current method is started at the time of activation. $t_{\text{pulse}} > 100 \text{ ms}$
Stop	The current method is canceled (Stop) at the time of activation. $t_{\text{pulse}} > 100 \text{ ms}$
Quit	The current command in the determination run will be canceled at the time of activation. $t_{\text{pulse}} > 100 \text{ ms}$
Ready	The instrument is ready to receive a start signal.
Conditioning OK	The line is set when Conditioning with SET titration and KFT titration is at OK. The line remains set until the determination is started with [START] .
Determination	The instrument performs a data-generating determination.
EOD	End of Determination. Pulse ($t_{\text{pulse}} = 200 \text{ ms}$) after a determination or after a buffer/standard solution during calibration using a Sample Processor.
Error	The line is set for error message display.



Function	Explanation
Warning	The line is set for warning message display.

9 Technical specifications

9.1 Measuring interface

The 901 Titrande has one galvanically isolated measuring interface.

The measuring cycle is 100 ms for all measuring modes.

9.1.1 Potentiometry

One high-ohm measuring input (**Ind.**) for pH, metal or ion-selective electrodes and one measuring input (**Ref.**) for separate reference electrodes.

Input resistance > $1 \cdot 10^{12} \Omega$

Offset current < $1 \cdot 10^{-12} \text{ A}$ (under reference conditions)

Measuring mode
pH

Measuring range -13 to +20 pH

Resolution 0.001 pH

Measuring accuracy $\pm 0.003 \text{ pH}$
(± 1 digit, without sensor error, under reference conditions)

Measuring mode
U

Measuring range -1,200 to +1,200 mV

Resolution 0.1 mV

Measuring accuracy $\pm 0.2 \text{ mV}$
(± 1 digit, without sensor error, under reference conditions)

9.1.2 Temperature

A measuring input (**Temp.**) for temperature sensors of the Pt1000 or NTC type with automatic temperature compensation.

R (25 °C) and B value can be configured for NTC sensors.

Measuring range

Pt1000 -150 to +250 °C

NTC -20 to +250 °C

(R (25 °C) = 30,000 Ω and B (25/50) = 3,964 K)

Resolution

Pt1000 0.1 °C



NTC 0.1 °C

Measuring accuracy

Pt1000 ±0.2 °C
(Applies for measuring range –20 to +150 °C; ±1 digit; without sensor error, under reference conditions)

NTC ±0.6 °C
(Applies for measuring range +10 to +40 °C; ±1 digit; without sensor error, under reference conditions)

9.1.3 Polarizer

One measuring input (**Pol.**) for polarizable electrodes.

Measuring mode
Ipol Determination with adjustable polarization current

Polarization current –122.5 to +122.5 µA (increment: 0.5 µA)
–125.0 to +125.0 µA: non-guaranteed values, dependent on reference voltage +2.5 V

Measuring range –1,200 to +1,200 mV

Resolution 0.1 mV

Measuring accuracy ±0.2 mV
(±1 digit, without sensor error, under reference conditions)

Measuring mode
Upol Determination with adjustable polarization voltage

Polarization voltage –1,225 to +1,225 mV (increment: 25 mV)
–1,250 to +1,250 mV: non-guaranteed values, dependent on reference voltage +2.5 V

Measuring range –120 to +120 µA

Resolution 0.1 µA

9.2 Power connection

<i>Supply voltage</i>	100–240 V ($\pm 10\%$)
<i>Frequency</i>	50–60 Hz ($\pm 3\%$)
<i>Power consumption</i>	max. 45 W
<i>Fuse</i>	Electronic overload protection

9.3 Ambient conditions

Automatic interior temperature monitoring

<i>Nominal function range</i>	+5 to +45 °C at max. 80% relative humidity, non-condensing
<i>Storage</i>	+5 to +45 °C at max. 80% relative humidity, non-condensing
<i>Altitude / Pressure range</i>	max. 2,000 m. above sea level / min. 780 mbar
<i>Overvoltage category</i>	II
<i>Pollution degree</i>	2

9.4 Reference conditions

<i>Ambient temperature</i>	+25 °C (± 3 °C)
<i>Relative humidity</i>	$\leq 60\%$
<i>Device status</i>	The device is in operation for at least 30 min.



9.5 Dimensions

<i>Width</i>	142 mm
<i>Height</i>	227 mm
<i>Depth</i>	231 mm
<i>Weight</i>	2.6 kg (without accessories)
<i>Material (housing)</i>	Polybutylene terephthalate (PBT)

9.6 Interfaces

USB connectors

<i>USB ports</i>	2 USB downstream ports (type A sockets), 500 mA, for connecting peripheral devices such as printers, keyboards, barcode readers or RS-232/USB boxes (6.2148.020).
------------------	---

"Controller" connector

<i>Controller port</i>	USB upstream port with auxiliary power supply (Mini DIN socket) for connecting Touch Control or computer for controlling the 901 Titrande.
<i>Touch Control</i>	With integrated Touch Control cable.
<i>Computer</i>	With 6.2151.000 cable

MSB connectors (Metrohm Serial Bus)

<i>Dosing device</i>	Connection of a max. of 4 external dosing devices of the Dosimat or Dosino (MSB 1 to MSB 4) types.
<i>Stirrer</i>	Connection of a max. of 4 stirrers. Stirrer control: Switching on/off manually or coordinated with the titration sequence. Speed in 15 steps and shift direction can be selected.
<i>Remote Box</i>	Connection of a max. of 4 Remote Boxes. Remote Boxes can be used to actuate and monitor external devices.

Index

6.2151.000 controller cable	14
685 Dosimat plus	17
700 Dosino	17
800 Dosino	17
801 Stirrer	18
803 Ti Stand	18
804 Ti Stand	18
805 Dosimat	17

A

ADD	4
Ambient conditions	53

B

Balance	22
Barcode reader	
Connect	24

C

Calibration mode	
CAL	3
Computer	
Connect	14
Computer keyboard	
Connect	23
Connect	
Balance	22
Barcode reader	24
Computer	14
Computer keyboard	23
Dosing device	17
MSB devices	16
Power grid	13
Printer	21
Remote Box	19
Stirrer	18
Titration stand	18
Touch Control	12
USB hub	21

Connector

MSB	2
USB	2

Controller

Connector	10
-----------------	----

D

Device description	2
Device software	
Update	2
Differential potentiometry	35

Dosing command	3
ADD	4
EMPTY	4
LQH	4
PREP	4

Dosing device

Connect	17
---------------	----

Driver software

Install	14
---------------	----

E

Electrode	
Connect	30
Electrode test	
ELT	3
Electrostatic charge	8
EMPTY	4

H

Humidity	53
----------------	----

I

Indicator electrode	
Warm-up time	42
Installation	
Driver software	14
Instrument type	10

K

Keyboard	
Connect	23
KFT	3

L

LED	
On	44
LQH	4

M

Maintenance	43
MEAS	3
Measuring interface	2, 10
Measuring mode	3
MEAS	3
Metrohm Serial Bus MSB, see also	
"MSB"	16
MSB	
Connect devices	16
Connector	10
MSB connector	2

O

Operation	53
Overview of the instrument	9
Overvoltage category	53

P

Pin assignment	47
Power connection	10, 13
PREP	4
Printer	21

R

Reagent replacement	41
Remote	
Interface	47
Pin assignment	47
Remote Box	
Connect	19

S

Safety instructions	6
Sample	
Liquid	40
Solid	41
Sample size	
Size	39
Sea level	53
Sensor	
Connect	30
Serial number	10
Service	7
SET	3
Stirrer	
Connect	18
Stirring direction	25
Storage	53
Supply voltage	8
System test	12

T

Temperature	53
Temperature sensor	
Connect	30
Titrand system	1
Titration mode	3
KFT	3
SET	3
Titration stand	
Connect	18

Index

Titration vessel	
Equip	25
Touch Control	
Connect	12

U

Update	
Device software	2
USB	
Connector	10
USB connector	2

USB hub	
Connect	21

V

Volumetric KF titration	
1-component reagents	37
2-component reagents	37
Application	37
Endpoint determination, biamperometric (Upol)	36



Endpoint determination, biamperometric (Ipol)	36
Principle	36
Sample size	39
Tips for daily practice	37
Water standard	37
Working conditions	41

W

Water standard	
Certified	37