

909 UV Digester



Manual
8.909.8001EN



Metrohm AG

CH-9100 Herisau

Switzerland

Phone +41 71 353 85 85

Fax +41 71 353 89 01

info@metrohm.com

www.metrohm.com

909 UV Digester

Manual

Teachware
Metrohm AG
CH-9100 Herisau
teachware@metrohm.com

This documentation is protected by copyright. All rights reserved.

Although all the information given in this documentation has been checked with great care, errors cannot be entirely excluded. Should you notice any mistakes please send us your comments using the address given above.

Documentation in additional languages can be found on
<http://documents.metrohm.com>.

Table of contents

| | | |
|----------|--|-----------|
| 1 | Introduction | 1 |
| 1.1 | Instrument description | 1 |
| 1.1.1 | Supply voltage | 1 |
| 1.1.2 | Intended use | 1 |
| 1.2 | Function description | 2 |
| 1.3 | About the documentation | 3 |
| 1.3.1 | Symbols and conventions | 3 |
| 1.4 | Safety instructions | 4 |
| 1.4.1 | General notes on safety | 4 |
| 1.4.2 | Electrical safety | 4 |
| 1.4.3 | Personnel safety | 5 |
| 1.4.4 | Flammable solvents and chemicals | 6 |
| 1.4.5 | Recycling and disposal | 6 |
| 2 | Overview of the instrument | 7 |
| 3 | Installation | 9 |
| 3.1 | Setting up the instrument | 9 |
| 3.1.1 | Packaging | 9 |
| 3.1.2 | Checks | 9 |
| 3.1.3 | Location | 9 |
| 3.2 | Inserting the UV lamp | 9 |
| 3.3 | Mounting the UV protective shield | 11 |
| 3.4 | Preparing and connecting the temperature sensor | 14 |
| 3.5 | Inserting the sample holder | 17 |
| 3.6 | Connecting the power supply cable | 18 |
| 4 | Operation | 20 |
| 4.1 | Switching the instrument on and off | 20 |
| 4.2 | Dialog window overview | 20 |
| 4.2.1 | Main dialog | 20 |
| 4.2.2 | "Parameters" dialog | 20 |
| 4.2.3 | "Lamp info" dialog | 21 |
| 4.2.4 | "Lamp setup" dialog | 21 |
| 4.3 | The keypad | 23 |
| 4.3.1 | The key functions | 23 |
| 4.4 | The display | 24 |
| 4.4.1 | Display elements | 24 |
| 4.4.2 | Status display | 25 |



| | | |
|-----------|---|-----------|
| 4.5 | Changing numerical values | 25 |
| 5 | Carrying out sample digestion | 26 |
| 6 | Application examples | 31 |
| 7 | Operation and maintenance | 33 |
| 7.1 | General | 33 |
| 7.2 | Care | 34 |
| 7.2.1 | Maintaining and replacing the UV lamp | 34 |
| 7.2.2 | Clearing a jammed fan | 35 |
| 7.2.3 | Cleaning sample glasses | 36 |
| 7.3 | Quality Management and qualification with Metrohm | 36 |
| 8 | Troubleshooting | 37 |
| 8.1 | Error messages | 37 |
| 8.2 | Additional faults and their remedy | 38 |
| 9 | Technical specifications | 39 |
| 9.1 | Radiating | 39 |
| 9.2 | Mercury vapor lamp | 39 |
| 9.3 | Temperature measuring input | 39 |
| 9.4 | Power connection | 40 |
| 9.4.1 | 220 - 240 V power connection | 40 |
| 9.4.2 | 100 - 127 V power connection with a transformer | 40 |
| 9.5 | Safety specifications | 40 |
| 9.6 | Electromagnetic compatibility (EMC) | 41 |
| 9.7 | Ambient temperature | 41 |
| 9.8 | Reference conditions | 42 |
| 9.9 | Dimensions | 42 |
| 10 | Warranty (guarantee) | 43 |
| 11 | Accessories | 45 |
| 11.1 | Scope of delivery | 45 |
| 11.2 | Optional accessories | 49 |
| | Index | 50 |

Table of figures

| | | |
|-----------|---|----|
| Figure 1 | UV lamp spectrum | 2 |
| Figure 2 | Front 909 UV Digester | 7 |
| Figure 3 | Rear 909 UV Digester | 8 |
| Figure 4 | Guiding the UV lamp downward into the wet end | 10 |
| Figure 5 | Installing the UV lamp in the wet end | 11 |
| Figure 6 | Removing the stop | 11 |
| Figure 7 | Positioning the UV protective shield over the wet end | 12 |
| Figure 8 | Positioning the stop over the threaded hole | 13 |
| Figure 9 | Screwing the stop in place | 13 |
| Figure 10 | Removing the temperature sensor from the storage vessel | 14 |
| Figure 11 | Removing the spacer from the temperature sensor | 14 |
| Figure 12 | Removing the SGJ sleeve from the temperature sensor | 14 |
| Figure 13 | Attaching the thermometer plug on the temperature sensor | 15 |
| Figure 14 | Attaching the spacer to the temperature sensor | 15 |
| Figure 15 | Inserting the storage vessel for the temperature sensor into the holder | 16 |
| Figure 16 | Connecting a temperature sensor | 16 |
| Figure 17 | Guiding the sample holder into the wet end | 17 |
| Figure 18 | Positioning the cylindrical guide for the sample holder | 18 |
| Figure 19 | Inserting the temperature sensor at sample position 6 | 18 |
| Figure 20 | Keypad for the operating unit | 23 |
| Figure 21 | Main dialog - Display elements | 24 |

1 Introduction

1.1 Instrument description

The 909 UV Digester is used for the digestion of aqueous samples moderately loaded with organic material. Its primary use is sample preparation for analyses in voltammetry, ion chromatography, AAS and ICP-AES. The digestion of samples is based on the photolytic generation of hydroxyl radicals, whose formation is initiated by a hydrogen peroxide radical initiator. The hydroxyl radicals immediately react with the organic compounds and work to break them down. The temperature and time for digestion can be entered conveniently using the operating unit. The temperature is controlled automatically using air cooling and the UV lamp's emitter power (medium pressure mercury lamp), which is regulated between 100% and 50% as needed.

The 909 UV Digester combines the wet end and control unit in one housing. The UV lamp is located in the center of the wet end. The sample vessels are inserted in the sample holder encircling the UV lamp. The fan for air cooling is located in the lower portion of the wet end. The control unit consists of the operating unit and all of the electronics including the electronic ballast for regulating the UV lamp. The 909 UV Digester is used as a standalone instrument.

1.1.1 Supply voltage

The 909 UV Digester can be operated only with a supply voltage of 220 to 240 V. Therefore, for a supply voltage of 100 to 120 V, a transformer has to be connected between the power and the instrument. You can find details on the transformer's specifications in Chapter *Technical specifications*, page 39.

1.1.2 Intended use

The 909 UV Digester is a sample preparation instrument that can be used in analytical laboratories for preparing samples for the following analysis processes:

- Voltammetry
- Ion chromatography
- AAS
- ICP-AES



Chemicals and flammable substances can be used as a sample additive or auxiliary reagent in this instrument. Usage of the 909 UV Digester therefore requires the user to have basic knowledge and experience in handling toxic and caustic substances. Knowledge with respect to the application



of the fire prevention measures prescribed for laboratories is also mandatory.

1.2 Function description

Through the effects of UV radiation, dissolved organic matter that disrupts the trace analysis of heavy metals is decomposed. However, the UV radiation itself does not cause the digestion of the organic sample components; rather, hydroxyl radicals that form break down the organic matter. Hydroxyl radicals are initially formed from the added hydrogen peroxide that functions as the initiator for the radical chain reaction. Subsequently, additional hydroxyl radicals are generated from the water if the temperature is high enough. Depending on the sample matrix, it may be necessary to add more H_2O_2 during digestion.



Spectrum of the UV lamp

The following figure shows a typical spectrum for the UV lamp. Certain deviations from this spectrum are possible; this can be traced back to production-related factors.




Figure 1 UV lamp spectrum

1.3 About the documentation



Caution

Please read through this documentation carefully before putting the instrument into operation. The documentation contains information and warnings which the user must follow in order to ensure safe operation of the instrument.

1.3.1 Symbols and conventions

The following symbols and formatting may appear in this documentation:

(5-12)

Cross-reference to figure legend

The first number refers to the figure number, the second to the instrument part in the figure.

1

Instruction step

Carry out these steps in the sequence shown.

Method

Dialog text, parameter in the software

File ▶ New

Menu or menu item

[Continue]

Button or key



Warning

This symbol draws attention to a possible life-threatening hazard or risk of injury.



Warning

This symbol draws attention to a possible hazard due to electrical current.



Warning

This symbol draws attention to a possible hazard due to heat or hot instrument parts.



Warning

This symbol draws attention to hazardous UV radiation.



Warning

This symbol draws attention to a possible biological hazard.

**Caution**

This symbol draws attention to possible damage to instruments or instrument parts.

**Note**

This symbol marks additional information and tips.

1.4 Safety instructions

1.4.1 General notes on safety

**Warning**

This instrument may only be operated in accordance with the specifications in this documentation.

This instrument has left the factory in a flawless state in terms of technical safety. To maintain this state and ensure non-hazardous operation of the instrument, the following instructions must be observed carefully.

1.4.2 Electrical safety

The electrical safety when working with the instrument is ensured as part of the international standard IEC 61010.

**Warning**

Only personnel qualified by Metrohm are authorized to carry out service work on electronic components.

**Warning**

Never open the housing of the instrument. The instrument could be damaged by this. There is also a risk of serious injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.

Supply voltage



Warning

An incorrect supply voltage can damage the instrument.

Operate this instrument only with a supply voltage specified for it (refer to the rear of the instrument). Use a transformer if necessary (for specifications for the transformer see *Chapter 9.4.2, page 40*).

Protection against electrostatic charges



Warning

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Do not fail to pull the mains cable out of the mains connection socket before you set up or disconnect electrical plug connections at the rear of the instrument.

1.4.3 Personnel safety



Warning

Hazardous chemical substances

Chemical substances can cause chemical burns and skin damage.

Correctly install the UV protective shield and wear protective glasses, protective gloves and work clothes appropriate for laboratory work when operating the 909 UV Digester.



Warning

Harmful UV radiation

UV radiation damages the eyes.

Protect your eyes from UV light during operation. Correctly install the UV protective shield and wear UV protective glasses.

The 909 UV Digester has to be operated in a closed fume cupboard or in a closed laboratory and may be operated only by trained personnel.

Attach the supplied UV warning shield included in the accessories to the fume cupboard.



Warning

Escaping ozone

Ozone can damage mucous membranes.

The 909 UV Digester has to be operated in a closed fume cupboard or in a closed laboratory and may be operated only by trained personnel.



Warning

Hot surfaces

The UV lamp, sample vessels and sample holder become very hot during operation. The UV lamp can reach a temperature of 900 °C.

Skin contact during and immediately after a radiating procedure will cause burns.

Avoid skin contact until the instrument is completely cool or wear heat-insulated gloves if necessary.

1.4.4 Flammable solvents and chemicals




Warning

All relevant safety measures are to be observed when working with flammable solvents and chemicals.

- Set the instrument up in a closed fume cupboard.
- Keep all sources of flame far from the workplace.
- Clean up spilled liquids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.

1.4.5 Recycling and disposal



This product is covered by European Directive 2002/96/EC, WEEE – Waste from Electrical and Electronic Equipment.

The correct disposal of your old equipment will help to prevent negative effects on the environment and public health.

More details about the disposal of your old equipment can be obtained from your local authorities, from waste disposal companies or from your local dealer.

2 Overview of the instrument




Figure 2 Front 909 UV Digester

| | |
|--|---|
| 1 Wet end | 2 Warning symbols "UV radiation" and "Hot surfaces". |
| 3 Stops For positioning the UV protective shield. | 4 Lamp holder With upper lamp socket. |
| 5 Holder For inserting the sample holder. | 6 Stop With hex socket head screw, for mounting the UV protective shield. |
| 7 "On" LED Lights up if the instrument is switched on. | 8 Operating unit With keyboard and display. |
| 9 Control unit | |




Figure 3 Rear 909 UV Digester

1 Holder

For storage vessel for temperature sensor.

2 Temperature sensor connection (Temp.)

For connecting the Pt1000 temperature sensor. Two B sockets, 2 mm.

3 Power switch

For switching the instrument on and off.

4 Power socket

5 Type plate

Contains the serial number.

3 Installation

3.1 Setting up the instrument

3.1.1 Packaging

The instrument is supplied in highly protective special packaging together with the separately packed accessories. Keep this packaging, as only this ensures safe transportation of the instrument.

3.1.2 Checks

Immediately after receipt, check whether the shipment has arrived complete and without damage by comparing it with the delivery note.

3.1.3 Location

The instrument has been developed for operation indoors and may not be used in explosive environments.

In the wet end of the 909 UV Digester, ozone from the air is released due to the UV radiation. Set up the instrument in a fume cupboard or closed laboratory at a location in the laboratory that is suitable for operation and free of vibrations. The location should be protected from corrosive atmospheres and contamination from chemicals.

The instrument should be protected against excessive temperature fluctuations and direct sunlight.

3.2 Inserting the UV lamp

The UV lamp (6.2804.090) is included in the scope of delivery for the 909 UV Digester.



Caution

Spots on the quartz tube

Spots will be burned into the quartz glass if you touch the UV lamp's quartz tube with bare hands and put it into operation without cleaning it first. This lowers the power of the UV lamp.

Always use the cotton glove included with the UV lamp to install the lamp and only touch the UV lamp at the two flattened ends. Remove any grease spots or dust particles stuck to the quartz tube before starting the UV lamp.



Proceed as follows to insert the UV lamp:

- 1** Ensure that the power for the 909 UV Digester has been disconnected.
- 2** Remove the sample holder from the wet end if it has already been inserted.
- 3** Put on the cotton glove included with the UV lamp.
- 4** Wipe off any spots on the quartz tube using a clean towel wetted with pure alcohol.
- 5** Hold the UV lamp by one of the two flattened ends and guide it downward into the wet end.




Figure 4 Guiding the UV lamp downward into the wet end

- 6** Insert the metal contact projecting on the lower end of the UV lamp into the lower lamp socket and press it down gently.
- 7** While holding the UV lamp down, click the metal contact projecting upwards into the upper lamp socket.




Figure 5 Installing the UV lamp in the wet end

3.3 Mounting the UV protective shield

A UV protective shield (6.2745.100) pre-installed on a metal ring is included in the scope of delivery for the 909 UV Digester.

Proceed as follows:

- 1 Using the hex key (6.2621.100), unscrew the stop's hex socket head screw (2-6) on the back of the housing and remove the stop.

Keep the stop with the hex socket head screw in reach. It will have to be screwed back in place on the housing afterwards.




Figure 6 Removing the stop

- 2 Place the UV protective shield over the wet end's opening so that the metal ring rests on the two front stops (2-3).



Ensure that the UV protective shield is positioned towards the front.




Figure 7 Positioning the UV protective shield over the wet end

- 3** Place the stop removed earlier over the threaded hole in the back on the housing.

Ensure that the metal ring for the UV protective shield is located in the slot between the lower and upper part of the stop.



3.4 Preparing and connecting the temperature sensor

The temperature sensor (6.1110.010) is included in the scope of delivery for the 909 UV Digester.

Preparing

Proceed as follows:

- 1 Remove the temperature sensor from the storage vessel.




Figure 10 Removing the temperature sensor from the storage vessel

- 2 Remove the spacer in the lower part of the glass tube.




Figure 11 Removing the spacer from the temperature sensor

Keep the spacer in reach. It will have to be attached back on the glass tube afterwards.

- 3 Remove the SGJ sleeve.




Figure 12 Removing the SGJ sleeve from the temperature sensor

- 4 Guide the thermometer plug (6.1446.190) over the temperature sensor from the bottom until it is flush with the head of the temperature sensor.




Figure 15 Inserting the storage vessel for the temperature sensor into the holder

For storing the temperature sensor between digestion processes.

- 2 Insert the temperature sensor equipped with the temperature plug and spacer into the storage vessel.

Connecting

Proceed as follows:

- 1 Insert the plugs for the temperature sensor into the **Temp.** sockets.




Figure 16 Connecting a temperature sensor

**Note**

Always insert the red plug into the red socket. This is the only way that shielding against electrical interference can be ensured.

3.5 Inserting the sample holder

The sample holder has to be fully filled to ensure uniform temperature distribution in the samples (see "Preparing samples", page 26).

Proceed as follows:

- 1** Guide the sample holder on the mounting bracket into the wet end so that the slot between sample positions **1** and **12** ends up located at the front at the lamp holder (2-4).




Figure 17 Guiding the sample holder into the wet end

Ensure that the cylindrical guide on the bottom of the sample holder's mounting bracket clicks in place in the holder (2-5).




Figure 18 Positioning the cylindrical guide for the sample holder

- 2 Always insert the temperature sensor into the sample vessel filled with sample at **Position 6** of the sample holder.




Figure 19 Inserting the temperature sensor at sample position 6

3.6 Connecting the power supply cable



Warning

An incorrect supply voltage can damage the instrument.

The 909 UV Digester is designed for 220 to 240 V. You need to use a transformer if your power supply provides 100 to 120 V. You can find the specification for the transformer in *Chapter 9.4.2, page 40*.

Proceed as follows to connect the instrument to the power supply.



Without a transformer

- 1** Connect the UV Digester to the power supply using the power supply cable included in the scope of delivery.

With a transformer

- 1** Connect the transformer to the power supply.
- 2** Connect the UV Digester to the transformer using the power supply cable included in the scope of delivery.



4 Operation

4.1 Switching the instrument on and off

The toggle switch for switching on and off is located on the rear of the instrument (3-3).


After switching on the instrument, the main dialog is displayed with the most recently configured parameters.

When starting for the first time, the standard parameters **Sample Temp.** = 90 °C and **Rad. Time** = 90 min appear.

4.2 Dialog window overview

The following four dialogs are available for operating the 909 UV Digester:

4.2.1 Main dialog




The main dialog is displayed after switching on the 909 UV Digester.

The following information is displayed:

- The temperature currently being measured in the sample (**Curr. Temp.**).
- The sample setpoint temperature configured in the **Parameters** dialog (**Set Temp.**).

4.2.2 "Parameters" dialog



You can open the **Parameters** dialog from the main dialog by pressing the **[SET]** key.

You must configure the following parameters in this dialog:

- The sample setpoint temperature (**Sample Temp.**).
- The radiating time (**Rad. Time**).

Sample Temp.

Setpoint temperature for the samples.

| | |
|---------------|---|
| Input range | 80 - 100 °C (Increment: 1) |
| Default value | 90 °C |

Rad. Time

Radiating time for the samples.

| | |
|---------------|---|
| Input range | 1 - 999 min (Increment: 1) |
| Default value | 90 min |



Note

- The configured parameters can be modified during digestion.
- The configured parameters remain stored in the instrument after it is switched off.

4.2.3 "Lamp info" dialog

```
Lamp info
>Counters
Hours           75 h
Ignitions       80
```

You can open the **Lamp info** dialog from the main dialog by pressing the **[OK]** key twice.

The following information is displayed:

- The cumulative number of operating hours for the current UV lamp in use since the last replacement (**Hours**).
- The number of times the current UV lamp in use has been switched on since the last replacement (**Ignitions**).

4.2.4 "Lamp setup" dialog

```
Lamp setup
Replace after 750 h
Reset Counters
```



You can open the **Lamp setup** dialog by pressing and holding the **[SET]** key while switching on the 909 UV Digester.

You have to edit the following parameters if you have inserted a new UV lamp:

- Reset the counters back to zero for the cumulative operating hours and number of times the lamp has been switched on (**Reset Counters**).
- Configure the number of operating hours for the new UV lamp until the next replacement (**Replace after**).

Replace after


Maximum number of operating hours for the UV lamp until the next replacement.

| | |
|---------------|---|
| Input range | 0 - 1000 h (Increment: 1) |
| Default value | 750 h |
| | We recommend replacing the UV lamp after 750 h of operating time. |

Reset counters

For resetting both counters for cumulative operating hours and the number of times the UV lamp has been switched on (ignitions) to 0.

Confirming with **[OK]** triggers the following confirmation prompt:



Selection **[OK] | [QUIT]**

[OK]

Both counters are reset to 0.

[QUIT]

Both counters remain at their current values.



Caution

The counters cannot be reset during the ongoing service life of a UV lamp currently in use.

4.3 The keypad

Below you will find an overview of the operating unit's keypad and the functions of the keys in the various dialog windows.




Figure 20 Keypad for the operating unit

4.3.1 The key functions

| Key | Main dialog | Parameters | Lamp info | Lamp setup |
|--------------|---|--|---|--|
| [START STOP] | <ul style="list-style-type: none"> Start lamp radiating. Stop lamp radiating. | – | – | – |
| [SET] | Open the Parameters dialog. | – | – | – |
| [QUIT] | Acknowledge an error message. | Depending on the cursor position: <ul style="list-style-type: none"> Discard a change. Leave the dialog. | Leave the dialog. | Depending on the cursor position: <ul style="list-style-type: none"> Discard a change. Leave the dialog. |
| [OK] | Open the Lamp info dialog. | Depending on the cursor position: <ul style="list-style-type: none"> Save a change and switch to the next parameter. Save a change and leave the dialog. | Display the number of cumulative operating hours and times the lamp has been switched on (ignitions). | Depending on the cursor position: <ul style="list-style-type: none"> Save a change. Trigger the confirmation prompt for resetting the counters. Reset the counters to zero. |



| Key | Main dialog | Parameters | Lamp info | Lamp setup |
|-----|----------------------------------|---|-------------------|---|
| [▲] | Increase the display's contrast. | Depending on the cursor position: <ul style="list-style-type: none"> Go up one line. Page up by one value list. | Go up one line. | Depending on the cursor position: <ul style="list-style-type: none"> Go up one line. Page up by one value list. |
| [▼] | Decrease the display's contrast. | Depending on the cursor position: <ul style="list-style-type: none"> Go down one line. Page down by one value list. | Go down one line. | Depending on the cursor position: <ul style="list-style-type: none"> Go down one line. Page down by one value list. |
| [▶] | – | Depending on the cursor position: <ul style="list-style-type: none"> Move the cursor to the value range. Move the cursor one position to the right. | – | Depending on the cursor position: <ul style="list-style-type: none"> Move the cursor to the value range. Move the cursor one position to the right. |
| [◀] | – | Depending on the cursor position: <ul style="list-style-type: none"> Move the cursor to the value range. Move the cursor one position to the left. | – | Depending on the cursor position: <ul style="list-style-type: none"> Move the cursor to the value range. Move the cursor one position to the left. |

4.4 The display

4.4.1 Display elements




Figure 21 Main dialog - Display elements

1 Radiating time

Shows the remaining radiating time since the UV lamp was started.

2 Temperature display

Shows the current temperature measured in the sample.

3 Temperature display





Shows the configured sample setpoint temperature.

4 Status line

Contains the current status or an error message along with the respective symbol.

4.4.2 Status display


The following symbols and the respective text can be displayed in the status line:

| Symbol | Text | Meaning |
|---|------------------------|--|
|  | Heating | The samples are being heated. |
|  | Temp. reached | The configured sample setpoint temperature has been reached. |
|  | Lamp cooling | The samples are being cooled. |
|  | Various error messages | see Chapter 8.1, page 37. |

4.5 Changing numerical values

The various parameters can be selected using the arrow keys [▼] and [▲]. The colors on the current line are reversed.

Change parameters as follows:

- 1 Use the arrow keys [▶] or [◀] to select a numerical value.
The selected digit is displayed in black on white. The colors for the rest of the line remain reversed.
- 
- 2 Change the digit using the arrow keys [▲] or [▼].
 - 3 Select the next digit and change it, etc.
 - 4 Apply the modified value with the [OK] key.
Press [QUIT] if the previous value is to be restored.



Note

Use the provided sample vessels

Only use the sample vessels made of quartz glass (6.2414.000) that have been provided. Normal reagent glasses do not let UV radiation through and do not ensure complete digestion. Clean the sample glasses before use as described in *Chapter 7.2.3, page 36*.

- 2 Add nitric acid ($w(\text{HNO}_3) = 65\%$, for trace analysis) and hydrogen peroxide ($w(\text{H}_2\text{O}_2) = 30\%$, for trace analysis).

Samples with a low organic load are offset with the following reagent concentrations:

- 0.01 mL HNO_3 ($w(\text{HNO}_3) = 65\%$, for trace analysis)
- 0.05 mL H_2O_2 ($w(\text{H}_2\text{O}_2) = 30\%$, for trace analysis)

The sample volume, the volume of the reagents and the digestion parameters must be adjusted accordingly for samples with increased organic load.

- 3 Place the PTFE stopper (6.1446.100) on the sample vessels. Leave the sample vessel at **Position 6** free for measuring the temperature.
- 4 Insert the sample holder into the wet end (*see Chapter 3.5, page 17*).



Note

Fully filled sample holder

Only a fully filled sample holder enables even temperature distribution to the samples. If you use fewer than twelve samples, fill the remaining sample vessels with distilled water.

- 5 Insert the temperature sensor (6.1110.010) with attached thermometer plug (6.1446.190) and spacer into the sample vessel at **Position 6** and press it down.

**Note**

Leave the PTFE stoppers on the sample vessels until you are ready to continue using the samples.

- 3 Switch off the 909 UV Digester or carry out additional sample digestion processes.

Canceling sample digestion

Proceed as follows if you would like to stop a sample digestion process before the radiating time has elapsed:

- 1 Press the **[START / STOP]** key.
The UV lamp is switched off.
- 2 Wait for the samples to cool (the fan switches off after 10 min).
- 3 Remove the sample holder with the samples from the wet end.
- 4 Switch off the 909 UV Digester or carry out additional sample digestion processes.

**Note****Opened samples smell like ammonia or fish (amine)**

If an open sample smells like ammonia or fish, you need to increase the amount of acid and/or H_2O_2 being added.

6 Application examples

You can find useful examples of analysis methods in the Application Bulletins and Application Notes below.

Table 1 Application Bulletins

| Application Bulletin No. | Analytes |
|--------------------------|----------------------------|
| 74 | Tl, Sb, Bi, Fe, Cu, V |
| 96 | Hg |
| 116 | Cr |
| 117 | Se |
| 123 | Fe, Mn |
| 131 | Al |
| 146 | Mo |
| 176 | Sn, Pb |
| 186 | Al |
| 207 | Ag |
| 220 | Pt, Rh |
| 226 | As |
| 231 | Zn, Cd, Pb, Cu, Tl, Ni, Co |
| 242 | W |
| 243 | Cr |
| 266 | Ti, U |

Table 2 Application Notes

| Application Note No. | Title |
|----------------------|--|
| V-83 | Zinc, cadmium, lead and copper in waste water after UV digestion |
| V-84 | Total chromium in waste water after UV digestion |
| V-89 | Mercury in waste water |
| V-92 | Nickel in white wine |



| Application Note No. | Title |
|-----------------------------|---|
| V-93 | Zinc, cadmium, lead and copper in red wine |
| V-94 | Rhodium and platinum in red wine |
| V-96 | Platinum in urine |
| V-106 | Nickel and cobalt in waste water |
| V-107 | Tin in waste water |
| V-108 | Thallium in waste water |
| V-109 | Selenium in waste water |
| V-110 | Chromium in waste water |
| V-120 | Nickel in glycol after UV digestion |
| V-121 | Determination of total iron in waste water of UV digestion |
| V-176 | Total selenium in tap water after reduction of selenium(VI) to selenium(IV) by means of UV Digester |
| V-201 | Nickel and cobalt in red wine |

7 Operation and maintenance



Warning

Never open the housing of the instrument. The instrument could be damaged by this. There is also a risk of serious injury if live components are touched.

There are no parts inside the housing which can be serviced or replaced by the user.

7.1 General

The 909 UV Digester requires appropriate care. Excess contamination of the instrument may result in functional disruptions and a reduction in the lifetime of the sturdy mechanics and electronics of the instrument.

Severe contamination can also have an influence on sample digestion. Regular cleaning of exposed parts can prevent this to a large extent.

Spilled chemicals and solvents must be removed immediately. In particular, the power plug should be protected from contamination.



Warning

Hot surfaces

The UV lamp, sample vessels and sample holder become very hot during operation. The UV lamp can reach a temperature of 900 °C.

Skin contact during and immediately after sample digestion will cause burns.

Only clean the wet end once it has cooled completely.



7.2 Care

7.2.1 Maintaining and replacing the UV lamp

Keep the following information in mind so that the lifetime of the UV lamp is not shortened unnecessarily:

- Clean the UV lamp as needed using a clean towel wetted with pure alcohol. This lets you remove particles from the quartz glass that have not yet been burned in.
Do not touch any part of the UV lamp other than the two flattened ends during cleaning and wear the provided cotton glove when removing or installing the lamp.
- Avoid switching the UV lamp on and off unnecessarily. Each time the lamp is switched on shortens its lifetime.

Removing the old UV lamp



Warning

Risk of electric shock due to defective electronic components

Coming into contact with potentially defective electronic components in the instrument could result in an electric shock.

Ensure that the instrument is unplugged before replacing the UV lamp.

Proceed as follows:

- 1 Ensure that the wet end has cooled completely.
- 2 Remove the sample holder from the wet end if it has been inserted.
- 3 Put on the glove included with the UV lamp.
- 4 Hold the UV lamp by the upper flattened end and gently press it down.
- 5 While holding the UV lamp down, release the upper end from the lamp socket.
- 6 Take the UV lamp out of the lower lamp socket and remove it from the wet end.

Inserting the new UV lamp

- 1 Insert the new UV lamp in accordance with *Chapter 3.2, page 9*.

Configuring parameters for the new UV lamp

- 1 Reset the counters to zero for the number of operating hours and the number of times the lamp has been switched on under **Lamp setup ▶ Reset Counters** (see "*Reset counters*", page 22).
- 2 Set the number of operating hours for the new UV lamp under **Lamp setup ▶ Replace after** (see "*Replace after*", page 22).

We recommend **750 h** of operating time.

7.2.2 Clearing a jammed fan

The wet end with inserted sample holder is designed so that no parts such as PTFE stoppers can fall down. If the fan happens to get jammed anyway, keep the following safety instruction in mind:



Warning

Jammed fan

Improper work on the fan may result in injuries.

Ensure that the instrument is unplugged before clearing the jammed fan.

Proceed as follows to remove PTFE stoppers that have fallen into the wet end of the 909 UV Digester:

- 1 Remove the sample holder.
- 2 Turn the instrument upside down.
- 3 Gently shake the overturned instrument until the PTFE stoppers come out.



7.2.3 Cleaning sample glasses

Clean the sample glasses after each sample digestion process:

Proceed as follows:

- 1 Thoroughly rinse the sample glasses with double-distilled water.
- 2 Place the sample glasses in dilute nitric acid overnight or until their next use.
- 3 Rinse the sample glasses again using double-distilled water.

7.3 Quality Management and qualification with Metrohm

Quality management

Metrohm offers you comprehensive support in implementing quality management measures for instruments and software. Further information on this can be found in the brochure "**Metrohm Quality Management**" available from your local Metrohm representative.

Qualification

Please contact your local Metrohm representative for support in qualification of instruments and software. The **Installation Qualification (IQ)** and **Operational Qualification (OQ)** are offered by Metrohm representatives as a service. They are carried out by trained employees using standardized qualification documents and in accordance with the currently applicable requirements of the regulated industry. Further information on this can be found in the brochure "**Analytical Instrument Qualification – Confidence in quality with IQ/OQ**".

Maintenance

The electronic and mechanical functional groups of Metrohm instruments can and should be checked by specialist personnel from Metrohm as part of a regular preventive maintenance schedule. Please ask your local Metrohm representative regarding the precise terms and conditions involved in concluding a corresponding maintenance agreement. Further information on this can be found in the brochure "**Metrohm Care Contracts – Protect your investment the smart way**" available from your local Metrohm representative.

8 Troubleshooting

8.1 Error messages

The following error messages can occur when operating the 909 UV Digester. Follow the instructions in the "Remedy" column and acknowledge each of the error messages by pressing the **[QUIT]** button.

| Error message | Cause | Remedy |
|---------------------------|--|--|
| Max. temp. reached | The temperature in the wet end is above 105 °C. The UV lamp is switched off automatically. | <ul style="list-style-type: none"> ▪ Check the sample matrix and adjust it if necessary. ▪ Check whether sample fluid may have evaporated and refill the sample if necessary. ▪ Check the fan and clear it if jammed (<i>see Chapter 7.2.2, page 35</i>). ▪ Contact your local Metrohm Service department. |
| Adj. data missing | The instrument is not ready for operation since the checksum for the calibration data was incorrect when the instrument was switched on. | <ul style="list-style-type: none"> ▪ Switch the instrument off and back on. ▪ Contact your local Metrohm Service department. |
| Lamp exhausted | The maximum lifetime configured for the UV lamp in the Lamp setup dialog has been reached. | Replace the UV lamp (<i>see Chapter 7.2.1, page 34</i>). |
| No sensor connectd | The temperature sensor is either not connected or incorrectly connected. | <ul style="list-style-type: none"> ▪ Check the temperature sensor's connection. ▪ Connect the temperature sensor. |



| Error message | Cause | Remedy |
|---------------------------|--|--|
| Check Sensor/Lamp | <p>No increase in temperature is being recorded. Possible causes:</p> <ul style="list-style-type: none"> ▪ The temperature sensor has not been inserted in the sample. ▪ The sample holder including the temperature sensor has not been inserted into the wet end. ▪ The UV lamp is defective. ▪ The electronics are defective. | <ul style="list-style-type: none"> ▪ Insert the temperature sensor into the sample. ▪ Insert the sample holder including the temperature sensor into the instrument. ▪ Replace the UV lamp (<i>see Chapter 7.2.1, page 34</i>). ▪ Contact your local Metrohm Service department. |
| Dev. busy cooling | The UV lamp is too hot to be switched on again. | Wait for the end of the cooling process (Lamp cooling status message disappears). |
| Thermo switch open | The instrument is defective. Overtemperature protection is switching off the electronic ballast. | Contact your local Metrohm Service department. |

8.2 Additional faults and their remedy

| Problem | Cause | Remedy |
|--|---|---|
| The configured digestion temperature cannot be saved. | <i>The entered temperature is outside the input range of 80 - 100 °C.</i> | Set a temperature within the input range. |

9 Technical specifications

9.1 Radiating

| | |
|------------------------------|---------------|
| <i>Digestion temperature</i> | +80 - +100 °C |
| <i>Ambient temperature</i> | +5 - +35 °C |
| <i>Duration</i> | 1 - 999 min |
| <i>Power</i> | 600 - 300 W |
| <i>Control accuracy</i> | ±3 °C |

9.2 Mercury vapor lamp

| | |
|-----------------------|----------------------|
| <i>UV-A radiation</i> | 315 - 400 nm 35 W |
| <i>UV-B radiation</i> | 280 - 315 nm 40 W |
| <i>UV-C radiation</i> | 200 - 280 nm 70 W |
| <i>Dimensions</i> | |
| <i>Length</i> | 137 - 139 mm |
| <i>Outer diameter</i> | approx. 16 mm |

9.3 Temperature measuring input

| | |
|------------------------|---|
| | A measuring input for the Pt1000 type temperature sensor. |
| <i>Accuracy</i> | ±0.2 °C (under reference conditions) |
| <i>Measuring range</i> | -20.0 - +150.0 °C |



9.4 Power connection

9.4.1 220 - 240 V power connection

| | |
|--------------------------|--------------------------|
| <i>Supply voltage</i> | 220 - 240 V (\pm 10%) |
| <i>Frequency</i> | 50 - 60 Hz (\pm 3%) |
| <i>Power consumption</i> | 630 W |
| <i>Fuse</i> | 6.3 ATH |

9.4.2 100 - 127 V power connection with a transformer

Specification for the transformer

| | |
|-----------------------|---------------|
| <i>Input voltage</i> | 110 V typical |
| <i>Output voltage</i> | 230 V typical |
| <i>Frequency</i> | 50 - 60 Hz |
| <i>Power</i> | > 700 VA |


9.5 Safety specifications

This instrument fulfills the following electrical safety requirements:



CE marking in accordance with the EU directives:

- 2006/95/EC (Low Voltage Directive, LVD)
- 2004/108/EC (EMC Directive, EMC)



Federal Inspectorate for Heavy Current Installations ESTI (Accreditation Number SCESp 033)

- Safety mark for certification type 2 in accordance with NEV (type testing with market monitoring, EMC conformity)

Design and testing

According to EN/IEC/UL 61010-1, CSA-C22.2 No. 61010-1, protection class I, EN/IEC 61010-2-010, EN/IEC 60529.

Safety instructions

This document contains safety instructions which have to be followed by the user in order to ensure safe operation of the instrument.

9.6 Electromagnetic compatibility (EMC)

Emission

- | | |
|----------------------------|---|
| <i>Standards fulfilled</i> | <ul style="list-style-type: none"> ▪ EN/IEC 61326-1 ▪ EN/IEC 61000-6-3 ▪ EN 55011 / CISPR 11 ▪ EN/IEC 61000-3-2 ▪ EN/IEC 61000-3-3 |
|----------------------------|---|

Immunity

- | | |
|----------------------------|---|
| <i>Standards fulfilled</i> | <ul style="list-style-type: none"> ▪ EN/IEC 61326-1 ▪ EN/IEC 61000-6-2 ▪ EN/IEC 61000-4-2 ▪ EN/IEC 61000-4-3 ▪ EN/IEC 61000-4-4 ▪ EN/IEC 61000-4-5 ▪ EN/IEC 61000-4-6 ▪ EN/IEC 61000-4-11 ▪ EN/IEC 61000-4-14 ▪ EN/IEC 61000-4-28 |
|----------------------------|---|

9.7 Ambient temperature

- | | |
|-------------------------------|----------------------------------|
| <i>Nominal function range</i> | +10 - +35 °C (Humidity < 80%) |
| <i>Storage</i> | -20 - +70 °C (Humidity < 95%) |
| <i>Transport</i> | -40 - +70 °C (Humidity < 95%) |



9.8 Reference conditions

| | |
|-------------------------------------|---|
| <i>Ambient temperature</i> | +25 °C (± 3 °C) |
| <i>Supply voltage</i> | 230 V |
| <i>Relative humidity</i> | $\leq 60\%$ |
| <i>Operating temperature status</i> | Instrument in operation at least 30 min |
| <i>Validity</i> | After adjustment |

9.9 Dimensions

| | |
|-------------------------------------|------------------------|
| <i>Width</i> | 330 mm |
| <i>Height</i> | 290 mm |
| <i>With UV protective shield</i> | 500 mm |
| <i>Depth</i> | 310 mm |
| <i>Weight (without accessories)</i> | 10.90 kg |
| <i>Material of housing</i> | Metal, surface-treated |

10 Warranty (guarantee)

Metrohm guarantees that the deliveries and services it provides are free of defects in materials, design or manufacturing.

The general warranty period is 36 months (exclusions below) from the date of delivery or 18 months in the event of continuous operation. The warranty remains valid on the condition that the servicing is provided by a service organization authorized by Metrohm at defined intervals and with a defined scope.

The warranty period for anion suppressors of the type "MSM" is 120 months from the date of delivery or 60 months in the case of continuous operation.

The warranty period for IC separation columns is 90 days after start-up.

For third-party components that are recognizable as such, the manufacturer's warranty regulations apply.

Consumables and materials with limited storage life and glass breakage in the case of electrodes or other glass parts are excluded from the warranty.

Warranty claims cannot be asserted if the ordering party has failed to meet its payment obligations according to schedule.

During the warranty period, Metrohm undertakes either to replace free of charge or to credit the purchaser for any modules or components that can be shown to be faulty. Any transport or customs fees that may apply are the ordering party's responsibility.

The precondition for this is that the ordering party has to specify the article number, the article designation, an adequate error description, the delivery date and (if applicable) the serial number or chip data in the Support Tracker. Metrohm then decides whether a replacement or a credit note is to be issued or whether the faulty part has to be returned using the Return Material Authorization (RMA). If a replacement or credit note is issued, the ordering party undertakes to store the faulty part for at least 24 months in accordance with the current storage directives (in compliance with ESD guidelines) and to hold it in readiness for onsite inspection or for return shipment to Metrohm. Metrohm reserves the right to invoice the ordering party for these articles, including retroactively, in the event of noncompliance with these preconditions.

The same warranty periods that are specified for a corresponding new part apply to parts that are replaced or repaired within the above-mentioned warranty periods. However, replacement or repair of a part does not extend the warranty period of the entire system.



Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty.

Metrohm also offers a 120 month spare parts availability guarantee and a 60 month PC software support warranty, calculated from the date on which the product is withdrawn from the market. The content of this warranty is the ability of the customer to obtain functioning spare parts or appropriate software support at market prices during the time of the warranty period.

If Metrohm AG is unable to meet this obligation due to circumstances beyond the control of Metrohm AG, then the ordering party shall be offered alternative solutions at preferential conditions.

11 Accessories




Note

Subject to change without notice.




11.1 Scope of delivery




2.909.0014 909 UV Digester

| Qty. | Order no. | Description | | | | | | | | | | | | | | | | | | |
|--------------------------------|--|---|-----------------|-------|------------------|-------------|-----------------|----|--------------------------------|-----|-------------|---------------------|--------------------------|----|-----------------------------|---|----------------------------|----|-------------------------|--|
| 1 | 1.909.0014 | <p>909 UV Digester 220-240 V/50-60 Hz</p> <p>Digestion instrument for UV photolysis of water samples with low to medium organic load. For sample preparation in trace element determination by means of voltammetry, ion chromatography and spectroscopy (AAS, ICP). Integrated instrument with operating unit and wet end. With air cooling and automatic control of digestion temperature and time. For 12 samples with a maximum of 12 mL sample volume each. Instrument for 220 - 240 V and 50 - 60 Hz.</p> | | | | | | | | | | | | | | | | | | |
| 1 | 6.1110.010 | <p>Pt1000 temperature sensor</p> <p>Pt1000 temperature sensor (Class B), 1 mm plug for 909 UV Digester</p> <table> <tr> <td>Shaft material:</td> <td>Glass</td> </tr> <tr> <td>Measuring range:</td> <td>-50 ... 180</td> </tr> <tr> <td>Measuring unit:</td> <td>°C</td> </tr> <tr> <td>Max. installation length (mm):</td> <td>120</td> </tr> <tr> <td>SGJ sleeve:</td> <td>Flexible SGJ sleeve</td> </tr> <tr> <td>Shaft diameter top (mm):</td> <td>12</td> </tr> <tr> <td>Shaft diameter bottom (mm):</td> <td>5</td> </tr> <tr> <td>Min. immersion depth (mm):</td> <td>20</td> </tr> <tr> <td>Electrode plug-in head:</td> <td>Fixed cable (l = 0.5 m) with plug 2 x 1 mm</td> </tr> </table> | Shaft material: | Glass | Measuring range: | -50 ... 180 | Measuring unit: | °C | Max. installation length (mm): | 120 | SGJ sleeve: | Flexible SGJ sleeve | Shaft diameter top (mm): | 12 | Shaft diameter bottom (mm): | 5 | Min. immersion depth (mm): | 20 | Electrode plug-in head: | Fixed cable (l = 0.5 m) with plug 2 x 1 mm |
| Shaft material: | Glass | | | | | | | | | | | | | | | | | | | |
| Measuring range: | -50 ... 180 | | | | | | | | | | | | | | | | | | | |
| Measuring unit: | °C | | | | | | | | | | | | | | | | | | | |
| Max. installation length (mm): | 120 | | | | | | | | | | | | | | | | | | | |
| SGJ sleeve: | Flexible SGJ sleeve | | | | | | | | | | | | | | | | | | | |
| Shaft diameter top (mm): | 12 | | | | | | | | | | | | | | | | | | | |
| Shaft diameter bottom (mm): | 5 | | | | | | | | | | | | | | | | | | | |
| Min. immersion depth (mm): | 20 | | | | | | | | | | | | | | | | | | | |
| Electrode plug-in head: | Fixed cable (l = 0.5 m) with plug 2 x 1 mm | | | | | | | | | | | | | | | | | | | |








| Qty. | Order no. | Description | |
|------|------------|---|---|
| 12 | 6.1446.100 | PTFE stopper for UV quartz sample vessel 12 mL |  |
| | | For quartz sample vessel 6.2414.000 | |
| | | Material: PTFE black | |
| | | Height (mm): 35 Outer diameter (mm): 19 | |
| 1 | 6.1446.190 | Thermometer stopper for Pt1000 tempera- ture sensor |  |
| | | For 6.1110.010 Pt1000 temperature sensor and 6.2414.000 12 mL quartz sample vessel. | |
| | | Material: PP | |
| | | Material remark: Stopper | |
| | | Material 2: Nitrile rubber | |
| | | Material remark 2: O-ring Height (mm): 22 Outer diameter (mm): 18 Inner diameter (mm): 12.2 | |
| 1 | 6.2008.040 | Storage vessel |  |
| | | Together with 6.2008.050 storage vessel holder provides a support for the electrode on 6.3032.XXX dosing units. | |

| Qty. | Order no. | Description | | |
|------|------------|---|---|----------------------------|
| 1 | 6.2041.240 | Sample holder for 12 quartz sample vessels 12 mL |  | |
| | | For 909 UV Digester. | | |
| | | Material: | | Stainless steel (AISI 304) |
| | | Height (mm): | | 185 |
| | | Outer diameter (mm): | 100 | |
| 12 | 6.2414.000 | Quartz sample vessel 12 mL |  | |
| | | For UV Digester | | |
| | | Material: | | Quartz |
| | | Outer diameter (mm): | | 15.6 |
| | | Length (mm): | | 125 |
| | | Volume (mL): | 12 | |
| 1 | 6.2621.100 | Hexagon key 3 mm |  | |
| | | Hexagon key 3 mm for IC Sample Processors | | |
| | | Length (mm): | | 73 |




| Qty. | Order no. | Description | |
|----------|---------------------|--|---|
| 1 | 6.2745.100 | UV protective shield |  |
| | | For 909 UV Digester. | |
| | | Material: | Plexiglas (PMMA) |
| | | Material remark: | Protective shield |
| | | Material 2: | Aluminum |
| | | Material remark 2: | Mounting ring |
| | | Height (mm): | 212 |
| | | Outer diameter (mm): | 160 |
| | | Inner diameter (mm): | 140 |
| 1 | 6.2804.090 | UV mercury vapor lamp |  |
| | | Medium pressure mercury vapor lamp for 909 UV Digester | |
| | | Material: | Quartz |
| | | Outer diameter (mm): | 16 ... 18 |
| | | Length (mm): | 137 ... 139 |
| 1 | 6.2225.000 | UV warning plate |  |
| | | For 909 UV Digester | |
| | | Material: | PVC |
| | | Width (mm): | 100 |
| | | Height (mm): | 92 |
| 1 | 6.2122.0x0 | Power supply cable with IEC 60320 line socket, type C13 | |
| | | Cable plug according to customer requirements. | |
| | | Switzerland: | SEV 1011, Type 12 6.2122.020 |
| | | Germany, ...: | CEE 7, Type VII 6.2122.040 |
| | | USA, ...: | NEMA 5-15, Type 498 6.2122.070 |
| | | Length: | 1.5 m |
| 1 | 8.909.8001EN | 909 UV Digester Manual | |

11.2 Optional accessories

2.909.0014 909 UV Digester

| Qty. | Order no. | Description |
|------|----------------------|--------------------------------------|
| 1 | 6.2042.050 | Spacer for Pt1000 temperature sensor |
| | Material: | PTFE |
| | Height (mm): | 1.5 |
| | Outer diameter (mm): | 11.5 |
| | Inner diameter (mm): | 4.8 |





Index

A

| | |
|-------------|----|
| Accessories | |
| Optional | 49 |
| Application | |
| Examples | 31 |

C

| | |
|--------------|---|
| Control end | 1 |
| Control unit | 7 |

D

| | |
|-----------------------|----|
| Dialog window | |
| "Lamp info" | 21 |
| "Lamp setup" | 21 |
| "Parameters" | 20 |
| Display elements | 24 |
| Main dialog | 20 |
| Symbols | 25 |
| Digestion | |
| Cancel | 29 |
| End | 28 |
| Interference | 30 |
| Start | 28 |
| Digestion temperature | |
| Cannot be saved | 38 |

E

| | |
|----------------------|----|
| Electrostatic charge | 5 |
| Error message | 37 |

F

| | |
|-------------------------|----|
| Fan | |
| Jammed | 35 |
| Front of the instrument | 7 |
| Function description | 2 |

G

| | |
|-----------|----|
| GLP | 36 |
| Guarantee | 43 |

I

| | |
|----------------------|----|
| Installation | |
| Storage vessel | 15 |
| UV lamp | 9 |
| UV protective shield | 11 |

K

| | |
|----------------|----|
| Keypad | |
| Operating unit | 23 |

Keys

| | |
|----------|----|
| Function | 23 |
|----------|----|

L

| | |
|-----------------------|----|
| Lamp ignitions | |
| Reset | 22 |
| Lamp use | |
| Number | 21 |
| Literature | |
| Application Bulletins | 31 |
| Application Notes | 31 |

M

| | |
|-----------------------|----|
| Maintenance Agreement | 36 |
| Mercury vapor lamp | |
| Radiation data | 39 |

O

| | |
|-------------------------|--------|
| Operating hours | |
| UV lamp | 21, 22 |
| Operating hours counter | |
| Reset | 22 |
| Operation | 2 |
| Keypad | 23 |
| Organic load | |
| High | 26 |
| Ozone emission | |
| Protective measures | 6 |

P

| | |
|------------------|-------|
| Parameter | |
| Change | 25 |
| Configure | 20 |
| Power connection | |
| 100 to 120 volts | 19 |
| 220 to 240 volts | 19 |
| Transformer | 1, 19 |

Q

| | |
|--------------------|----|
| Quality Management | 36 |
|--------------------|----|

R

| | |
|------------------------|----|
| Radiating time | |
| Configure | 21 |
| Reaction equation | 2 |
| Reagents | |
| Add | 27 |
| Rear of the instrument | 8 |
| Reset | |
| Lamp ignitions | 22 |

| | |
|-------------------------|----|
| Operating hours counter | 22 |
|-------------------------|----|

S

| | |
|---------------------|-------|
| Safety instructions | 4 |
| Personnel safety | 5 |
| Sample digestion | |
| Procedure | 26 |
| Sample glasses | |
| Clean | 36 |
| Sample holder | |
| Insert | 17 |
| Sample temperature | |
| Configure | 21 |
| Samples | |
| Prepare | 26 |
| Scope of delivery | 45 |
| Service | 4 |
| Setup location | 9 |
| Supply voltage | 5, 18 |
| Symbols | |
| Meaning | 25 |

T

| | |
|--------------------|----|
| Temperature | |
| Configure | 21 |
| Temperature sensor | |
| Insert | 27 |
| Temperature sensor | |
| Connect | 16 |
| Prepare | 14 |
| Storage | 15 |
| Transformer | |
| Power connection | 1 |

U

| | |
|----------------------|--------|
| UV digester | |
| Prepare | 26 |
| UV lamp | |
| Clean | 34 |
| Handling | 9 |
| Ignitions | 21 |
| Install | 9 |
| New | 22 |
| Operating hours | 21, 22 |
| Radiation data | 39 |
| Replace | 34 |
| Spectrum | 2 |
| UV protective shield | |
| Mount | 11 |



UV radiation
Protective measures 28
Protective measures 5

V
Validation 36

W
Warranty 43
Wet end 1, 7