

920 Absorber Module



Manual

8.920.8004EN / v6 / 2023-12-31



Metrohm AG
CH-9100 Herisau
Switzerland
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

920 Absorber Module

2.920.0010

Manual

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Disclaimer

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g. conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

Table of contents

1	Introduction	1
1.1	Instrument description	1
1.2	Intended use	2
1.3	About the documentation	2
1.3.1	Content and scope	2
1.3.2	Symbols and conventions	3
1.4	Safety instructions	4
1.4.1	General notes on safety	4
1.4.2	Electrical safety	4
1.4.3	Tubing and capillary connections	5
1.4.4	Flammable solvents and chemicals	5
1.4.5	Recycling and disposal	6
2	The Combustion IC system	7
2.1	Overview	7
2.2	System description	8
2.3	Sequence of a determination	8
3	Overview of the instrument	10
3.1	Front	10
3.2	Rear	14
4	Installation	16
4.1	Setting up the instrument	16
4.1.1	Packaging	16
4.1.2	Checks	16
4.1.3	Location	16
4.2	Overview	16
4.3	Dosino for ultrapure water	19
4.3.1	Assembling the Dosino	22
4.3.2	Feeding the capillary into the instrument for Combustion Module (AJ)	23
4.3.3	Feeding the capillary into the instrument for Combustion Oven (TEI)	24
4.3.4	Connecting the ports – Overview	25
4.3.5	Inserting the Dosino in the instrument	33
4.3.6	Connecting the Metrosep I Trap 1	34
4.3.7	Connecting the ultrapure water bottle	37
4.4	Dosino for the absorber solution	39
4.4.1	Assembling the Dosino	40

4.4.2	Connecting the ports – Overview	41
4.4.3	Inserting the Dosino in the instrument	46
4.4.4	Connecting the T connector	47
4.4.5	Connecting the Metrosep A Trap 1	49
4.4.6	Connecting the bottle with the absorber solution	53
4.5	10-port-valve	55
4.5.1	Connecting the ports	56
4.5.2	Connecting the bottle with the standard solution	59
4.5.3	Connecting the bottle with the check standard	60
4.5.4	Function of the 10-port valve	61
4.6	6-port valve	62
4.6.1	Connecting the ports	63
4.6.2	Mode of operation of the 6-port valve	66
4.7	Absorber tube	68
4.8	Connecting the drainage tubing and leak sensor	72
4.8.1	Installing the drainage tubing	72
4.8.2	Connecting the leak sensor	73
4.9	Connecting the instrument to a computer	74
4.10	Connecting the instrument to the power grid	75
4.11	Connecting the 920 Absorber Module to the remaining components	76
4.11.1	Connecting the Combustion Module (AJ) inlet: water infeed ..	76
4.11.2	Connecting the Combustion Module (AJ) outlet: absorber solution	79
4.11.3	Connecting the Combustion Oven (TEI) inlet	81
4.11.4	Connecting the Combustion Oven (TEI) outlet	82
4.11.5	Connecting the 930 Compact IC Flex	82
5	Start-up	85
6	Operation and maintenance	88
6.1	General	88
6.2	Door	88
6.3	Dosino	88
6.4	Dosing unit	88
6.5	Metrosep I Trap 1	89
6.6	Metrosep A Trap 1	89
6.7	10-port valve	89
6.8	6-port valve	89
6.9	Absorber tube	89

Index 97

1 Introduction

1.1 Instrument description

The 920 Absorber Module forms part of the Combustion-IC-System, which includes a Combustion Module (AJ) / Combustion Oven (TEI) and an ion chromatograph in addition to the 920 Absorber Module. The Combustion-IC-System can be used for digesting solid, liquid and gaseous combustible samples. The Combustion-IC-System's prime application is the determination of halogenides and sulfur.

The 920 Absorber Module is the connecting link between combustion and IC analysis in the Combustion-IC-System. The 920 Absorber Module absorbs the combustion gases in an aqueous medium and transfers it into the ion chromatography system. The 920 Absorber Module performs the necessary Liquid Handling and enables the application of various sample preparation techniques, such as matrix elimination, preconcentration and automated calibration.

The 920 Absorber Module is comprised of the following modules:

Dosino for ultrapure water

A Dosino with a 5 mL dosing unit is inserted in the back in the 920 Absorber Module. This Dosino serves to dose water for the water infeed during combustion and for the transfer of the absorption solution to the IC instrument. Furthermore, it is used for liquid handling.

Dosino for the absorber solution

A Dosino with a 10 mL dosing unit is inserted in the front of the 920 Absorber Module. This Dosino serves to dose the absorber solution into the absorber tube.

10-port valve

The 10-port valve is the centerpiece of the Liquid Handling; it serves as a switch-over point between the various liquid flows.

6-port valve

The 6-port valve serves as a switch-over point between the two liquid flows. The sample loop allows a precisely defined sample amount to be measured before it is transferred to the analyzer.

Absorber tube

The combustion gases are absorbed by a solution in the absorber tube.

Leak sensor

The leak sensor detects leaking liquid that collects in the instrument's base tray. Liquid that leaks in the instrument is routed to the base tray using drainage tubing and detected there.

1.2 Intended use

The 920 Absorber Module represents the connecting link between combustion and IC analysis in a Combustion-IC-System. It absorbs the combustion gases in an absorber solution and transfers the resulting sample to an ion chromatograph. Furthermore, the 920 Absorber Module allows additional sample preparation techniques to be applied.

This instrument is suitable for processing chemicals and flammable samples. Therefore, the use of the 920 Absorber Module requires the user to have basic knowledge and experience in handling toxic and caustic substances. Knowledge regarding the application of fire prevention measures prescribed for laboratories is also mandatory.

1.3 About the documentation



CAUTION

Please read through this documentation carefully before putting the instrument into operation. The documentation contains information and warnings which the user must follow in order to ensure safe operation of the instrument.

1.3.1 Content and scope

This manual focuses on the 920 Absorber Module instrument, which forms part of the Combustion-IC-System.

Separate manuals exist for the other modules of the Combustion-IC-System:

- Manuals for the Combustion Module including Auto Boat Drive (ABD) and optional autosampler (MMS 5000) from Analytik Jena (PDFs available on the enclosed CD)
- Manuals for the Combustion Oven and the optional samplers Liquid Autosampler CIC, Solid Autosampler CIC or GLS Sampler CIC from Trace Elemental Instruments (PDFs available on the enclosed USB flash drive)
- Manual for 930 Compact IC Flex
- Manual for 800 Dosino
- Manual for 807 Dosing Unit

1.3.2 Symbols and conventions

The following symbols and formatting may appear in this documentation:

(5- 12)	Cross-reference to figure legend The first number refers to the figure number, the second to the instrument part in the figure.
1	Instruction step Perform the steps one after the other.
Method	Dialog text, parameter in the software
File ► New	Menu or menu item
[Continue]	Button or key
	WARNING This symbol draws attention to a possible life-threatening hazard or risk of injury.
	WARNING This symbol draws attention to a possible hazard due to electrical current.
	WARNING This symbol draws attention to a possible hazard due to heat or hot instrument parts.
	WARNING This symbol draws attention to a possible biological hazard.
	CAUTION This symbol draws attention to possible damage to instruments or instrument parts.
	NOTE This symbol highlights additional information and tips.

Protection against electrostatic charges



WARNING

Electronic components are sensitive to electrostatic charges and can be destroyed by discharges.

Do not fail to pull the power cord out of the power socket before you set up or disconnect electrical plug connections at the rear of the instrument.

The device is to be operated only with the door closed.

1.4.3 Tubing and capillary connections



CAUTION

Leaks in tubing and capillary connections are a safety risk. Tighten all connections well by hand. Avoid applying excessive force to tubing connections. Damaged tubing ends lead to leakage. Appropriate tools can be used to loosen connections.

Check the connections regularly for leakage. If the instrument is used mainly in unattended operation, then weekly inspections are mandatory.

1.4.4 Flammable solvents and chemicals




WARNING

All relevant safety measures are to be observed when working with flammable solvents and chemicals.

- Set up the instrument in a well-ventilated location (e.g. fume cupboard).
- Keep all sources of flame far from the workplace.
- Clean up spilled liquids and solids immediately.
- Follow the safety instructions of the chemical manufacturer.



1.4.5 Recycling and disposal



Properly dispose of chemicals and of the product to reduce negative effects on the environment and public health. Local authorities, waste disposal companies or dealers provide more detailed information on disposal. Observe the WEEE EU directive (WEEE = Waste Electrical and Electronic Equipment) for the proper disposal of waste electronic equipment within the European Union.

2 The Combustion IC system

2.1 Overview

Combustion is a sample preparation technique that uses burning to digest a sample. Solid, liquid and gaseous samples that are flammable can be digested using combustion. Combustion IC is particularly suited for determining the following elements:

	1	2	13	14	15	16	17	18
	Ia							VIIa
1 (K)	001H 1.00794	IIa	IIIa	IVa	Va	VIa	VIIa	002He 4.002602
2 (L)	003Li 6.941	004Be 9.012182	005B 10.811	006C 12.0107	007N 14.00674	008O 15.9994	009F 18.998403	010Ne 20.1797
3 (M)	011Na 22.98977	012Mg 24.305	013Al 26.981538	014Si 28.0855	015P 30.973761	016S 32.066	017Cl 35.4527	018Ar 39.948
4 (N)	019K 39.0983	020Ca 40.078	031Ga 69.723	032Ge 72.61	033As 74.9216	034Se 78.96	035Br 79.904	036Kr 83.8
5 (O)	037Rb 85.4678	038Sr 87.62	049In 114.818	050Sn 118.71	051Sb 121.76	052Te 127.6	053I 126.90447	054Xe 131.29
6 (P)	055Cs 132.90545	056Ba 137.327	081Tl 204.3833	082Pb 207.2	083Bi 208.98038	084Po [209]	085At [210]	086Rn [222]
7 (Q)	087Fr [223]	088Ra [226]	113Uut [287]	114Uuq [289]	115Uup [288]	116Uuh [292]	117Uus	118Uuo

Figure 1 Detectable elements

An additional instrument is required for the determination of gaseous samples: the LPG/GSS Module (AJ) (2.136.0720) or the GLS Sampler CIC (TEI) (2.0136.0630).




Figure 3 Sequence of the determination



NOTE

Do not mix these up!

- **Absorber solution:** e.g. H_2O_2 solution for oxidizing sulfur oxides.
- **Absorption solution:** Combustion gases dissolved in the absorber solution.

The determination is controlled entirely by MagIC Net. You can find a detailed description of the sequence in the online help for MagIC Net under the title "Combustion sequence".

- ## 920 Absorber Module

Capillary feed-throughs and cable feed-throughs




Figure 5 920 Absorber Module – Front

1 Capillary guides

To arrange capillaries and tubings in the instrument.

2 Capillary guide

For the combustion gas inlet capillary (Combustion Module (AJ)).

4 Capillary feed-throughs

Between instrument and base tray. For feeding capillaries from the front to the rear of the instrument.

6 Cable opening

For guiding Dosino cables out of the instrument.

- 

Capillaries and tubings




Figure 7 Capillaries and tubings in use – Overview

- | | |
|--|---|
| <p>1 PTFE capillary, 0.75 mm ID / 3 m (6.1803.160)
Connection between the 6-port valve and the waste container.</p> | <p>2 PEEK capillary, 0.5 mm ID / 15 cm (6.1831.040)
Connection between the 6-port valve and the 10-port valve.</p> |
| <p>3 PTFE capillary, 0.75 mm ID / 3 m (6.1803.160)
Approx. 80 cm. Connection between the 10-port valve and the bottle with the standard solution.</p> | <p>4 PTFE capillary, 0.75 mm ID / 3 m (6.1803.160)
Approx. 80 cm. Connection between the 10-port valve and the bottle with the check standard.</p> |
| <p>5 PTFE capillary, 0.75 mm ID / 3 m (6.1803.160)
Approx. 60 cm. Connection between the 10-port valve and the absorber tube.</p> | <p>6 PTFE capillary, 0.75 mm ID / 3 m (6.1803.160)
Approx. 5 cm. For aspirating air segments.</p> |
| <p>7 PTFE capillary, 0.75 mm ID / 3 m (6.1803.160)
Connection between the 10-port valve and the waste container.</p> | <p>8 Combustion Module (AJ): PTFE capillary 1/8 in. OD / 1/16 in. ID / 1.5 m Combustion Oven (TEI): PTFE capillary 1/8 in. OD / 1/16 in. ID / 5 m
Connection between the absorber tube and the Combustion Module (AJ) / Combustion Oven (TEI).</p> |

- | | |
|--|--|
| <p>9 Combustion Module (AJ): PTFE capillary 0.5 mm ID / 3 m (6.1803.030) Combustion Oven (TEI): PTFE capillary 0.8 mm ID / 2.5 m (capillary absorber solution)</p> <p>Connection between the Metrosep A Trap 1 - 100/4.0 trap column (6.1014.000) and the Combustion Module (AJ) / Combustion Oven (TEI).</p> | <p>10 PTFE capillary, 0.5 mm ID / 3 m (6.1803.030)</p> <p>Approx. 20 cm. Connection between the Dosino with the absorber solution and the Metrosep A Trap 1 - 100/4.0 trap column (6.1014.000).</p> |
| <p>11 FEP tubing M6 for 920 150 cm (6.1805.030) with FEP aspiration tubing (6.1819.100)</p> <p>Aspiration tubing for the absorber solution. Connection between the Dosino with the absorber solution and the bottle with the absorber solution.</p> | <p>12 FEP tubing M6 for 920, 13 cm (6.1805.010)</p> <p>Connection between the Dosino with the absorber solution and the T connector.</p> |
| <p>13 FEP tubing M6 for 920, 13 cm (6.1805.010)</p> <p>Connection between the Dosino with ultrapure water and the T connector.</p> | <p>14 Combustion Module (AJ): PTFE capillary 0.5 mm ID / 3 m (6.1803.030) Combustion Oven (TEI): PTFE capillary 0.8 mm ID / 2.5 m (6.07311.310)</p> <p>Connection between the Dosino with ultrapure water and the Combustion Module (AJ) / Combustion Oven (TEI).</p> |
| <p>15 FEP tubing M6 for 920 150 cm (6.1805.030) with FEP aspiration tubing (6.1819.100)</p> <p>Aspiration tubing for ultrapure water. Connection between the Dosino with the water bottle.</p> | <p>16 PTFE capillary, 0.5 mm ID / 3 m (6.1803.030)</p> <p>Approx. 30 cm. Connection between the Dosino with ultrapure water and the Metrosep I Trap 1 - 100/4.0 trap column (6.1014.200).</p> |
| <p>17 PTFE capillary, 0.5 mm ID / 3 m (6.1803.030)</p> <p>Approx. 40 cm. Connection between the Metrosep I Trap 1 - 100/4.0 trap column (6.1014.200) and the 6-port valve.</p> | <p>18 PTFE capillary, 0.5 mm ID / 3 m (6.1803.030)</p> <p>Connection between the 6-port valve and the 930 Compact IC Flex.</p> |

Colored sleeves

Different liquids move throughout the 920 Absorber Module. In order to provide you an overview of the different flow paths, we recommend using the provided colored sleeves as follows to code the capillaries that route the distinct liquids:

Table 1 Recommended color codes for the flow paths

Liquid	Colored sleeve
Ultrapure water	blue

Liquid	Colored sleeve
Absorber solution	yellow
Standard solution	black
Check standard	red
Waste	white



NOTE

The colored sleeves are heat-shrink tubing; this means they shrink when they are exposed to heat.

Slip the colored sleeve over the capillary and secure it to the desired spot by heating it, for example by using a hairdryer.

4.3 Dosino for ultrapure water

The Dosino for ultrapure water is installed first.




Figure 8 Installing the Dosino for ultrapure water


- | | |
|--|--|
| <p>13 FEP tubing M6 for 920, 13 cm (6.1805.010)</p> <p>Connection between the Dosino with ultrapure water and the T connector.</p> | <p>14 Combustion Module (AJ): PTFE capillary 0.5 mm ID / 3 m (6.1803.030) Combustion Oven (TEI): PTFE capillary 0.8 mm ID / 2.5 m (6.07311.310)</p> <p>Connection between the Dosino with ultrapure water and the Combustion Module (AJ) / Combustion Oven (TEI).</p> |
| <p>15 FEP tubing M6 for 920 150 cm (6.1805.030) with FEP aspiration tubing (6.1819.100)</p> <p>Aspiration tubing for ultrapure water.</p> <p>Connection between the Dosino with the water bottle.</p> | <p>16 PTFE capillary, 0.5 mm ID / 3 m (6.1803.030)</p> <p>Approx. 30 cm. Connection between the Dosino with ultrapure water and the Metrosep I Trap 1 - 100/4.0 trap column (6.1014.200).</p> |
| <p>17 PTFE capillary, 0.5 mm ID / 3 m (6.1803.030)</p> <p>Approx. 40 cm. Connection between the Metrosep I Trap 1 - 100/4.0 trap column (6.1014.200) and the 6-port valve.</p> | |

The Dosino is equipped with a 5 mL dosing unit and is inserted in the back in the accessory space (4-**10**). It has 3 tasks:

- Dosing the water infeed during combustion.


- Dosing the transfer water for matrix elimination of the sample.
- Liquid handling.

Dosing the water infeed during combustion



Feeding water in at the combustion tube's inlet supports the combustion process and increases the recovery rate, particularly for fluoride. To do so, the Dosino aspirates the ultrapure water from the bottle and doses a very small amount through port 4 to the combustion tube's inlet. The time and quantity of the water infeed are defined in MagIC Net™.


Dosing transfer water for the sample



The transfer water transports the sample to the IC instrument. In order to obtain accurate analysis results, the transfer water has to be free of ionic contamination. To remove all traces of ionic contamination, a Metrosep I Trap 1 – 100/4.0 (6.1014.200) column is used between the Dosino for ultrapure water and the 6-port valve.

The Dosino aspirates ultrapure water from the bottle and starts by dosing it to the Metrosep I Trap 1 – 100/4.0 (6.1014.200) through port 1 and to port 5 of the 6-port valve from there.

Liquid handling



A Dosino can either aspirate or discharge defined quantities of liquids via its 4 ports. In conjunction with the 10-port valve, it is used to transport liquids through the Combustion IC system. This process is called liquid handling.

Port 3 of the Dosino is connected to transfer tubing and is used for liquid handling.

Recommended procedure for installation

1. Assemble the Dosino (*see chapter 4.3.1, page 22*).
2. Guide the capillary for the infeed of water during combustion into the instrument (*see chapter 4.3.2, page 23*).
3. Connect port 2 and port 4 (*see "Connecting port 2 and port 4", page 26*).
4. Connect port 1 and port 3 (*see "Connecting port 1 and port 3", page 31*).
5. Place the Dosino in the instrument (*see chapter 4.3.5, page 33*).
6. Connect the Metrosep I Trap 1 - 100/4.0 trap column and insert it in the instrument (*see chapter 4.3.6, page 34*).
7. Connect the ultrapure water bottle (*see chapter 4.3.7, page 37*).

4.3.1 Assembling the Dosino

Accessories

For this step you need:

- 800 Dosino (2.800.0010)
- 807 Dosing Unit 5 mL (6.3032.150)

Attaching the Dosino to the 807 Dosing Unit

The manual for the 800 Dosino describes how to attach the 800 Dosino correctly to a dosing unit.



Plugging in the connection cable


- Next step

Feed the capillary for the water infeed (8-**14**) into the instrument.

Accessories

For this step you need:

- ## Feeding the capillary for the water infeed into the instrument



The arrows in the figure point in the direction the capillaries have to be pushed.

1 Feeding the capillary for the water infeed into the instrument

- Slide a blue colored sleeve over the capillary far enough that it remains visible outside the instrument.
- Push the capillary into the instrument through one of the 5 capillary openings on the back panel.

Next step

Connect the ports (see chapter 4.3.4, page 25).


4.3.3 Feeding the capillary into the instrument for Combustion Oven (TEI)

Accessories

For this step you need:

- PTFE capillary of the oxygen and ultrapure water connection tubing (6.07311.310)
- Colored sleeves, blue

Feeding the capillary for the water infeed into the instrument



The arrows in the figure point in the direction the capillaries have to be pushed.

1 Feeding the capillary for the water infeed into the instrument

- Slide a blue colored sleeve over the capillary far enough that it remains visible outside the instrument.

- Next step* Connect the ports (see chapter 4.3.4, page 25).

Use of the ports The ports on the Dosino for ultrapure water are used as follows:

Accessories The following accessories have to be connected to the Dosino for ultra-pure water:

For this step you also need:

- 920 Absorber Module




Figure 9 Connect the Dosino for ultrapure water

Connecting port 2 and port 4

Port 2 and port 4 of the Dosino are on the bottom of the instrument. These ports are connected first.

Accessories

The following accessories have to be connected to these 2 ports:

Port	... install this accessory
2	FEP tubing M6 for 920 150 cm (6.1805.030) with FEP aspiration tubing (6.1819.100) (8- 15)
4	- PTFE capillary for Combustion Module (AJ) (6.1803.030) (8-14) - PTFE capillary for Combustion Oven (TEI) (6.07311.310) - Pressure screw (6.2744.350) - Coupling (6.2744.080) - Adapter Dosino Port 4 (6.1808.280)

For this step you also need:

- Thread adapter (6.1618.020)
- Dosino holder (6.2057.230)
- Colored sleeves, blue


Preparing

Before you can connect the capillary and the aspiration tubing, you have to:

- Install the coupling for the connection between the capillary and Dosino.
- Thread the capillary for the water infeed and the aspiration tubing through the thread adapter and the Dosino holder.


1 Mounting the coupling

The ports on the Dosino are intended for tubing with an M6 connector. If capillaries are connected to the ports, a coupling has to be installed at the end.



- Optional: Guide a blue colored sleeve over the end of the capillary.
- During the installation of the Combustion Oven (TEI), the PTFE capillary (6.07311.310) must be shortened to the required length with a capillary cutter.
- Push the pressure screw over the end of the PTFE capillary far enough that a few millimeters of the capillary extend beyond the tip of the pressure screw.
- Attach the coupling on the pressure screw and tighten it by hand.

2 Threading the Dosino holder and thread adapter




- Start by threading the FEP tubing through the thread adapter from below and then through the Dosino holder.
- Start by threading the PTFE capillary through the thread adapter from below and then through the Dosino holder.

The capillary for the water infeed and the water aspiration tubing can now be connected to the Dosino. Proceed as follows:

Connecting port 2

Prerequisites

- FEP tubing M6 to 920 150 cm (6.1805.030) is threaded in through the thread adapter (6.1618.020) and the Dosino holder (6.2057.230).



1

Tighten the FEP tubing on port 2 of the dosing unit.


2

Guide the loose end of the FEP tubing out of the 920 Absorber Module through one of the capillary openings on the back panel of the instrument.

Connecting port 4

Prerequisites:

- The PTFE capillary for the Combustion Module (AJ) (6.1803.030) or the PTFE capillary for the Combustion Oven (TEI) (6.07311.310) is threaded in through the thread adapter (6.1618.020) and the Dosino holder (6.2057.230).
- The coupling is installed on the end of the capillary.



1 Mounting the adapter

Attach the Adapter Dosino port 4 (6.1808.280) onto port 4 of the dosing unit.


2 Mount the capillary

Tighten the PTFE capillary on the Adapter Dosino port 4.

Tightening the Dosino holder

Prerequisites:

The FEP tubing (6.1805.030) and the PTFE capillary for the Combustion Module (AJ) (6.1803.030) or the PTFE capillary for the Combustion Oven (TEI) (6.07311.310) are threaded in through the thread adapter (6.1618.020) and the Dosino holder (6.2057.230) and are connected to the ports 2 and 4.



- Connect the ports 1 and 3.

Port 1 on the left side of the Dosino is used for dosing water, which is used for transferring the sample from the 6-port valve to the analysis instrument. The transfer water is purified using a Metrosep I Trap 1 – 100/4.0 (6.1014.200) trap column so that it does not contaminate the sample with undesired anions and cations.

Accessories

The following accessories have to be connected to these 2 ports:

Port	... install this accessory
1	30 cm of the PTFE capillary (6.1803.030) (8-16) - Pressure screw (6.2744.350) - Coupling (6.2744.080)
3	FEP tubing M6 for 920, 13 cm (6.1805.010) (8- 13)

For this step you also need:


- Capillary cutter (6.2621.080)

Connecting port 1

1 Cutting the PTFE capillary to size


- Cut off a piece of the PTFE capillary that is approx. 30 cm long using the capillary cutter.

2 Mounting the coupling




- Push the pressure screw over the end of the PTFE capillary far enough that a few millimeters of the capillary extend beyond the tip of the pressure screw.
- Attach the coupling on the pressure screw and tighten it by hand.

3 Connecting the PTFE capillary



- ## Connecting port 3




Tighten the FEP tubing on port 3 of the dosing unit.

Place the Dosino for ultrapure water into the instrument.

Space is limited in the 920 Absorber Module.

When you insert the fully connected Dosino into the instrument, ensure that none of the tubing or capillaries are bent or kinked. Turn the Dosino so that the capillaries and tubing can be as straight as possible.

Inserting the Dosino in the instrument



- 1 Place the Dosino holder onto the mounting bolts and fix it in place using the knurled screw.

Next steps

Connect capillary/ tubing at port ...	... to
1	Metrosep I Trap 1 - 100/4.0 (6.1014.200) trap column inlet (see chapter 4.3.6, page 34)
2	Water bottle (see chapter 4.3.7, page 37)
3	T connector (see chapter 4.4.4, page 47)
4	Combustion oven inlet (see chapter 4.11.1, page 76)

4.3.6 Connecting the Metrosep I Trap 1

The Metrosep I Trap 1 – 100/4.0 trap column (6.1014.200) removes disruptive interfering anions and cations from the ultrapure water used for transferring the absorption solution from the absorber tube to the 6-port valve. The use of the trap column ensures that the analysis is not distorted by unwanted contamination.

Accessories


The following accessories have to be connected to the trap column:

Connecting the Metrosep I Trap 1 – 100/4.0



CAUTION

Be absolutely sure to note this flow direction when connecting the capillaries.



1 Threaded stopper

2 Column outlet

3 Flow direction indicator

4 Column inlet

1 Removing the threaded stoppers


Remove the threaded stoppers at the inlet and outlet of the Metrosep I Trap 1 – 100/4.0.

2 Connecting the capillaries

- Screw the loose end of the PTFE capillary connected to port 1 of the Dosino for ultrapure water to the inlet of the trap column using a pressure screw.
- Cut off a piece from the PTFE capillary (6.1803.030) that is approx. 40 cm long using the capillary cutter (6.2621.080) and tighten it onto the outlet of the trap column using a pressure screw.

Inserting the trap column in the instrument

Insert the Metrosep I Trap 1 – 100/4.0 in the instrument as follows:



- 1** Fasten the Metrosep I Trap 1 – 100/4.0 in place in the rear holding bracket (4-4) of the 920 Absorber Module using the column head.

Next step

Tighten the loose end of the PTFE capillary at the outlet of the trap column to port 5 of the 6-port valve (*see chapter 4.6, page 62*).

4.3.7 Connecting the ultrapure water bottle


The ultrapure water bottle has to be connected to port 2 of the Dosino for ultrapure water.

Accessories

For this step you need:

- Clear glass bottle 2 L (6.1608.070) filled with ultrapure water.
- Bottle cap (6.1602.160)
- FEP tubing M6 for 920 150 cm (6.1805.030) with aspiration tubing (6.1819.100) (8-**15**).
- Adsorber tube (6.1609.000)
- Cotton
- Adsorber granulate

Connecting the bottle with ultrapure water



1 Cotton

2 Adsorber granulate

- 1** Screw the bottle cap onto the bottle.
- 2** Push the FEP aspiration tubing (6.1819.100) through the M6 opening and screw it tight in the M6 opening.
- 3** Fill the adsorber tube with a little cotton and adsorber granulate and seal it with the plastic cover.
- 4** Place the adsorber tube into the large opening of the bottle cap and fasten it in place with the SGJ clip.
- 5** Seal the M8 opening with a stopper.

4.4 Dosino for the absorber solution




Figure 10 Installing the Dosino for the absorber solution

- 9 Combustion Module (AJ): PTFE capillary 0.5 mm ID / 3 m (6.1803.030) / Combustion Oven (TEI): PTFE capillary 0.5 mm ID / 3 m (capillary absorber solution)**

Connection between the Metrosep A Trap 1 - 100/4.0 trap column (6.1014.000) and the Combustion Module (AJ) / Combustion Oven (TEI).

- 10 PTFE capillary, 0.5 mm ID / 3 m (6.1803.030)**

Approx. 20 cm. Connection between the Dosino with the absorber solution and the Metrosep A Trap 1 - 100/4.0 trap column (6.1014.000).

- 11 FEP tubing M6 for 920 150 cm (6.1805.030) with FEP aspiration tubing (6.1819.100)**

Aspiration tubing for the absorber solution. Connection between the Dosino with the absorber solution and the bottle with the absorber solution.


- 12 FEP tubing M6 for 920, 13 cm (6.1805.010)**

Connection between the Dosino with the absorber solution and the T connector.

The Dosino for the absorber solution is equipped with a 10 mL dosing unit and has the following tasks:


- Dosing the absorber solution.
- Transferring the absorption solution to the 10-port valve.

Dosing the absorber solution



The absorber solution is dosed to the combustion tube's outlet in the Combustion Module (AJ) / Combustion Oven (TEI). It absorbs the gases from combustion and carries them to the absorber tube.

Transferring the absorption solution to the 10-port valve



Drawing the absorption solution out of the absorber tube.

4.4.1 Assembling the Dosino

Accessories

For this step you need:

- 800 Dosino (2.800.0010)
- Dosing unit 10 mL (6.3032.210)

Attaching the Dosino to the 807 Dosing Unit

The manual for the 800 Dosino describes how to attach the 800 Dosino correctly to a dosing unit.



Plugging in the connection cable

- #### 4.4.2 Connecting the ports – Overview

The ports on the Dosino for the absorber solution are used as follows:

Accessories

The following accessories have to be connected to the Dosino for the absorber solution:

For this step you also need:

- Thread adapter (6.1618.020)
- Dosino holder (6.2057.230)
- Capillary cutter (6.2621.080)




Figure 11 Connect the Dosino for the absorber solution

Connecting port 2

Port 2 is located on the bottom of the Dosino. Port 2 is used for aspirating the absorber solution.

Accessories

The following accessories have to be connected at port 2:

Port	Install this accessory
Port 2	FEP tubing M6 for 920 150 cm (6.1805.030) with FEP aspiration tubing (6.1819.100) (10- 11)

For this step you also need:

- Thread adapter (6.1618.020)
- Dosino holder (6.2057.230)


1 Threading the Dosino holder and thread adapter

- ## 2 Connecting the FEP tubing to port 2

3 Routing the FEP tubing out of the instrument

Tightening the Dosino holder on the Dosino for the absorber solution

FEP tubing (6.1805.030) is threaded in through the thread adapter (6.1618.020) and the Dosino holder (6.2057.230) and is connected to port 2.



- 1** Tighten the Dosino holder (6.2057.230) with the thread adapter (6.1618.020) on the bottom of the dosing unit.

Next step

Connect the ports 1 and 3.


Connecting port 1 and port 3

Accessories

The following accessories have to be connected to these two ports:

Port	Install this accessory
Port 1	FEP tubing M6 for 920, 13 cm (6.1805.010) (10- 12)
Port 3	- PTFE capillary (6.1803.030), approx. 20 cm (10-10) - Pressure screw (6.2744.350) - Coupling (6.2744.080)


1 Mounting the FEP tubing



Tighten the FEP tubing on port 1 of the dosing unit.


Connecting port 3

1 Mounting the coupling



- Push the pressure screw over the end of the PTFE capillary far enough that a few millimeters of the capillary extend beyond the tip of the pressure screw.
- Attach the coupling on the pressure screw and tighten it by hand.

2 Connecting the PTFE capillary



Tighten the coupling to port 3 of the dosing unit.

Next step

Placing the Dosino for the absorber solution into the instrument.

4.4.3 Inserting the Dosino in the instrument




CAUTION

Space is limited in the 920 Absorber Module.

When you insert the fully connected Dosino into the instrument, ensure that none of the tubing or capillaries are bent or kinked. Turn the Dosino so that the capillaries and tubing can be as straight as possible.

Inserting the Dosino in the instrument



- 1 Place the Dosino holder onto the mounting bolts and fix it in place using the knurled screw.

Next steps

Connect capillary/ tubing at port ...	... to
1	T connector (<i>see chapter 4.4.4, page 47</i>)
2	Bottle with the absorber solution (<i>see chapter 4.4.6, page 53</i>)
3	Connect the Metrosep A Trap 1 – 100/4.0 inlet (<i>see chapter 4.4.5, page 49</i>)

4.4.4 Connecting the T connector


The T connector is connected between the Dosino for ultrapure water, the Dosino for the absorber solution and the transfer tubing.

Accessories

At inlet ...	... install this accessory
1	FEP tubing M6 for 920, 13 cm (6.1805.010) from Dosino for absorber solution (7- 12)

At inlet ...	... install this accessory
2	FEP tubing M6 for 920, 13 cm (6.1805.010) from Dosino for ultra-pure water (7- 13)
3	Transfer tubing (4- 2)

Installing the T connector



1 Connecting the FEP tubing

Screw the loose end of the FEP tubing connected to port 3 of the Dosino for ultrapure water to an inlet **(1)** on the T connector.

2 Connecting the transfer tubing

Screw the longer end of the transfer tubing to the outlet **(3)** on the T connector.

3 Connecting the FEP tubing

Screw the loose end of the FEP tubing connected to port 1 of the Dosino for the absorber solution to the last inlet **(2)** on the T connector.

Next step

- Connect the transfer tubing to the 10-port valve (see *chapter 4.5.1, page 56*).

4.4.5 Connecting the Metrosep A Trap 1

The Metrosep A Trap 1 – 100/4.0 (6.1014.000) trap column continues to remove the smallest anionic contaminants in the absorber solution. The use of the trap column ensures that the analysis is not distorted by unwanted contamination.

The Metrosep A Trap 1 – 100/4.0 is inserted between the Dosino for the absorber solution and the T connector CIC (6.7301.210) at the outlet of the Combustion Module (AJ) or the oven outlet (TEI) and the Combustion Oven (TEI).

The ...	... is connected to
Inlet	Dosino for the absorber solution
Outlet	Combustion tube outlet on the Combustion Module (AJ) / Combustion Oven (TEI)

Accessories

Connect the following accessories to the trap column:

At the ...	... install this accessory
Inlet	▪ Approx. 20 cm of the PTFE capillary (6.1803.030) connected to port 3 of the Dosino for the absorber solution (10-10) ▪ Pressure screw (6.2744.010) ▪ Colored sleeve, yellow
Outlet	▪ Combustion Module (AJ): PTFE capillary (6.1803.030) ▪ Combustion Oven (TEI): PTFE capillary (Is labeled with "Absorber" and emerges from the rear of the Combustion Oven (TEI).) ▪ Optional: Colored sleeve, yellow

Feeding the capillary for the absorber solution into the instrument (Combustion Module (AJ))

Feed the capillary for the absorber solution into the instrument as follows:




Figure 12 Feeding the capillary for the absorber solution into the instrument

- 1
 - Slide a yellow colored sleeve over the capillary far enough that it remains visible outside the instrument.
 - Route one end of the PTFE- capillary (6.1803.030) to the front of the instrument from the right side of the instrument through the capillary feed-through between the base tray and the 920 Absorber Module.
 - Thread the end of the capillary through the bottom horizontal capillary guide on the right inner side wall and then through the vertical capillary guide from the bottom to the top.
 - Optional: Fasten another yellow colored sleeve to the end of the capillary.

Feeding the capillary for the absorber solution into the instrument (Combustion Oven (TEI))

Feed the capillary for the absorber solution into the instrument as follows:

- 1
 - Slide a yellow colored sleeve over the capillary far enough that it remains visible outside the instrument.
 - Route one end of the PTFE capillary that emerges from the rear of the Combustion Oven (TEI) and is labeled with "Absorber" to the front of the instrument from the back or from the left side of the instrument through the capillary feed-through between the base tray and the 920 Absorber Module.
 - Optional: Fasten another yellow colored sleeve to the end of the capillary.

Connecting the Metrosep A Trap 1 – 100/4.0



NOTE


The trap column is sealed with threaded stoppers at both ends.
Remove these threaded stoppers before connecting the capillaries.



CAUTION

Flow can only move through the trap column in one direction. This flow direction is marked on the column with an engraved arrow.

Be absolutely sure to note this flow direction when connecting the capillaries.



1 Threaded stopper

2 Column outlet

3 Flow direction indicator

4 Column inlet

1 Removing the threaded stoppers

Remove the threaded stoppers at the inlet and outlet of the Metrosep A Trap 1 – 100/4.0.

2 Connecting the column's inlet

- Screw the loose end of the PTFE capillary connected to port 3 of the Dosino for the absorber solution in place on the inlet of the trap column using a pressure screw (6.2744.010).

3 Connecting the column's outlet to the Combustion Module (AJ)

- Screw the PTFE capillary (6.1803.030) with the yellow colored sleeve (fed in from the right side of the instrument) in place on the trap column's outlet using a pressure screw (6.2744.010).

Connecting the column's outlet to the Combustion Oven (TEI)

- ## Inserting the trap column in the instrument

6.1014.000

- Next step*

Combustion Oven (TEI): The other end is already connected to the oven outlet.


The bottle with the absorber solution is connected to port 2 of Dosino 1.

For this step you need:

- Clear glass bottle 2 L (6.1608.070) filled with the absorber solution
- Bottle cap (6.1602.160)
- FEP tubing M6 for 920, 150 cm (6.1805.030) fastened to port 2 of the Dosino for the absorber solution.
- Adsorber tube (6.1609.000)
- Cotton

- Adsorber granulate

Connecting the bottle with the absorber solution



1 Cotton

2 Adsorber granulate

- 1 Screw the bottle cap onto the bottle.
 - 2 Insert a sealing ring from the 6.1602.160 accessory set into the bottle cap's M6 opening.
- Put the loose end of the FEP tubing through the bottle cap's M6 opening and tighten it.



NOTE

The end of the aspiration tubing has to reach to the bottom of the bottle.

- 3** Fill the adsorber tube with a little cotton and adsorber granulate and seal it with the plastic cover.
- 4** Place the adsorber tube into the large opening of the bottle cap and fasten it in place with the SGJ clip.
- 5** Seal the M8 opening with a stopper.

4.5 10-port-valve

The 10-port valve can be connected after both Dosinos and the two trap columns have been installed.




Figure 13 Connecting the 10-port valve

- 2** PEEK capillary, 0.5 mm ID / 15 cm (6.1831.040)
Connection between the 6-port valve and the 10-port valve.

- 3** PTFE capillary, 0.75 mm ID / 3 m (6.1803.160)
Approx. 80 cm. Connection between the 10-port valve and the bottle with the standard solution.

Connecting the 10-port valve



NOTE

We recommend connecting the central port first.




Figure 14 Connecting the 10-port valve

1 Connecting the central port

- Screw the adapter in place in the central port.
- Screw one end of the transfer tubing in place on the adapter.

2 Connecting port 1

- Screw the PEEK capillary to port 1 using a long pressure screw.

3 Connecting port 2

- Cut off a piece from the PTFE capillary, 0.75 mm ID / 3 m (6.1803.160) that is approx. 80 cm long using the capillary cutter.
- Screw this piece of capillary in place at port 2 using a long pressure screw.

4 Connecting port 3

- Cut off a piece from the PTFE capillary, 0.75 mm ID / 3 m (6.1803.160) that is approx. 80 cm long using the capillary cutter.
- Screw this piece of capillary in place at port 3 using a long pressure screw.

5 Sealing ports 4 through 7

Ports 4 through 7 are not used.

- Seal each of ports 4, 5, 6 and 7 with a threaded stopper.

6 Connecting port 8

- Cut off a piece from the PTFE capillary, 0.75 mm ID / 3 m (6.1803.160) that is approx. 60 cm long using the capillary cutter.
- Screw this piece of capillary in place at port 8 using a long pressure screw.

7 Connecting port 9

- Cut off a piece from the PTFE capillary, 0.75 mm ID / 3 m (6.1803.160) that is approx. 5 cm long using the capillary cutter.
- Screw this piece of capillary in place at port 9 using a long pressure screw.

8 Connecting port 10

- Screw an entire PTFE capillary, 0.75 mm ID / 3 m (6.1803.160) in place at port 10 using a long pressure screw.

Next step

Connect capillary/ tubing at port ...	... to
1	6-port valve (<i>see chapter 4.6, page 62</i>)
2	Bottle with standard solution (<i>see chapter 4.5.2, page 59</i>)
3	Bottle with check standard (<i>see chapter 4.5.3, page 60</i>)
8	Absorber tube (<i>see chapter 4.7, page 68</i>)

Connect ... to
capillary/
tubing at
port ...

10

Waste container

4.5.2 Connecting the bottle with the standard solution


The bottle with the standard solution is connected to the 10-port valve.

Accessories

For this step you need:

- Clear glass bottle 100 mL (6.1608.050) filled with standard solution
- Bottle attachment (6.1602.150), contains 3 pressure screws


Connecting the bottle with the standard solution



1 Connecting the aspirating capillary for the standard solution

- Slide a pressure screw over the loose end of the capillary connected to port 2 of the 10-port valve (13-3).
- Thread the capillary through one of the three openings on the bottle attachment.
- Screw the bottle attachment onto the bottle.
- Push the capillary far enough into the bottle that it reaches the bottom of the bottle. Then fix it in place with a pressure screw.

Inserting the bottle with the standard solution into the instrument



- 1 Insert the bottle with the standard solution into the bottle holder.

4.5.3 Connecting the bottle with the check standard

The bottle with the standard solution is connected to the 10-port valve.

Connecting the bottle with the check standard

- 1 Guide the loose end of the capillary (13-4) (connected to port 3 of the 10-port valve) in through a capillary feed-through between the device and the bottle holder on the side of the instrument.
- 2 Fasten the capillary in the bottle containing the check standard.
Place the bottle on the bottle holder.

4.5.4 Function of the 10-port valve

The 10-port valve is used for sample preparation. For this, the central port (15-2) can be connected to the other ten ports one after the other using the selector (15-1).

The transfer tubing is connected to the central port. Liquids can be aspirated into the transfer tubing by one of the ten ports with the aid of a Dosino and then transported to any other port by reversing the flow.




Figure 15 10-port valve – Ports


1 Selector

Can be rotated for selecting the active port.

2 Central port

Connected to the transfer tubing.

The selector can rotate in a clockwise or counterclockwise direction. Normally, the shortest path is selected when switching over to a new port.



Secured ports

In order to avoid contamination, any given port can be defined as "secured". The secured port is moved to during switching only if it is the target of the switching procedure.

If the secured port lies along the shortest route during a valve switchover, then the longer path will automatically be selected.

Example:

Switching from port 10 to port 7.

The shorter path leads through ports 9 and 8 (see figure 16, left). If however port 8 is defined as the secured port, then the 10-port valve automatically switches via ports 1 – 2 – 3 – 4 – 5 – 6 (see figure 16, left).




Figure 16 10-port valve – Normal valve switchover and with a secured port

4.6 6-port valve

The 6-port valve is connected next.



Figure 17 Connecting the 6-port valve

- | | |
|---|---|
| <p>1 PTFE capillary, 0.75 mm ID / 3 m (6.1803.160)</p> <p>Connection between the 6-port valve and the waste container.</p> | <p>2 PEEK capillary, 0.5 mm ID / 15 cm (6.1831.040)</p> <p>Connection between the 6-port valve and the 10-port valve.</p> |
| <p>17 PTFE capillary, 0.5 mm ID / 3 m (6.1803.030)</p> <p>Approx. 40 cm. Connection between the Metrosep I Trap 1 - 100/4.0 trap column (6.1014.200) and the 6-port valve.</p> | <p>18 PTFE capillary, 0.5 mm ID / 3 m (6.1803.030)</p> <p>Connection between the 6-port valve and the 930 Compact IC Flex.</p> |

4.6.1 Connecting the ports

Connectors

The ports of the 6-port valve are connected as follows:

Port ...	... is used for ...	... is connected to
1	Inlet for sample, standard solution and check standard	10-port valve
2	Outlet for waste	Waste container
3	Sample loop	6-port valve – Port 6
4	Outlet for sample, standard solution and check standard	Injection valve in the IC instrument

Accessories

64 ■■■■■■

Connecting the 6-port valve




Figure 18 Connecting the 6-port valve

1 Connecting port 1

- Tighten the loose end of the PEEK capillary connected to port 1 of the 10-port valve (17-2) to port 1 of the 6-port valve using a pressure screw.

2 Connecting port 2

- Push one end of the PTFE capillary, 0.75 mm ID / 3 m (6.1803.160), (17-1) through a capillary feed-through below the bottle holder from the back towards the front of the instrument.
- Push the pressure screw over the end of the PTFE capillary far enough that approx. 4 mm of the capillary is exposed.
- Screw the pressure screw in place at port 2 of the 6-port valve.
- Shorten the other end of the capillary using a capillary cutter so that it reaches the waste container.
- Fasten the end to the waste container.

3 Connecting port 4

- Push one end of the PTFE capillary, 0.5 mm ID / 3 m (6.1803.030), (17-18) through the cable feed-through at the bottom left on the front of the instrument and then route it up to the 6-port valve through the two capillary guides on the left inside the instrument.
- Push the pressure screw over the end of the capillary far enough that approx. 4 mm of the capillary is exposed.

- Screw the pressure screw in place at port 4 of the 6-port valve.

4 Connecting port 5

- Push the pressure screw over the loose end of the capillary connected to the outlet of the Metrosep I Trap 1 – 100/4.0 (17-**17**) far enough that approx. 4 mm of the capillary is exposed.
- Screw the pressure screw in place at port 5 of the 6-port valve.

5 Connecting the sample loop

- Push the pressure screws far enough over both ends of the sample loop that approx. 4 mm of the loop is exposed.
- Tighten the one pressure screw at port 3.
- Tighten the second pressure screw at port 6.

4.6.2 Mode of operation of the 6-port valve

A sample loop is installed on the 6-port valve; it is used for measuring the sample amount before it is transferred to the analyzer.

Ports

As its name implies, the 6-port valve has six ports.




Figure 19 6-port valve – Ports

1 Port 1

Connected to the 10-port valve. Used as the inlet for the sample and the standard solution and check standard.

2 Port 2

For connecting the capillary to the waste container.

- | | |
|--|---|
| 3 Port 3
For connecting the sample loop. | 4 Port 4
For connecting the capillary to the ion chromatograph. |
| 5 Port 5
For connecting the capillary from the outlet of the Metrosep I trap column. | 6 Port 6
For connecting the sample loop. |

Valve positions

The 6-port valve can have two valve position:




Figure 20 6-port valve – Valve positions

1 Filling

2 Injecting

- **Filling**

In the **Filling** position, ports 2 and 3, 4 and 5, and 6 and 1 are connected to each other.

The sample flows through port 1 into the 6-port valve, through the sample loop and through port 2 out of the 6-port valve into the waste container.

At the same time, the transfer solution flows through port 5 into the 6-port valve and directly back out through port 4.

The sample loop is filled with the sample at this valve position.

- **Injecting**

In the **Injecting** position, ports 1 and 2, 3 and 4, and 5 and 6 are connected to each other.

The transfer solution now flows through port 5 into the 6-port valve, through the sample loop and back out through port 4. At the same time, the absorber solution flows through port 1 into the 6-port valve and directly back out through port 2.

In this valve position, the absorber solution in the sample loop is flushed to the analyzer together with the transfer solution.

4.7 Absorber tube




Figure 21 Connect the absorber tube

- | | |
|---|--|
| <p>5 PTFE capillary, 0.75 mm ID / 3 m (6.1803.160)
Approx. 60 cm. Connection between the 10-port valve and the absorber tube.</p> | <p>8 Combustion Module (AJ): PTFE capillary 1/8 OD / 1/16 ID / 1.5 m Combustion Oven (TEI): PTFE capillary 1/8 OD / 1/16 ID / 5 m
Connection between the absorber tube and the Combustion Module (AJ) / Combustion Oven (TEI).

Combustion Oven (TEI): The capillary is already installed at the oven outlet and is labeled with "Combustion Gas".</p> |
|---|--|

The absorber tube is connected between the Combustion Module (AJ) / Combustion Oven (TEI) and the 10-port valve.

Accessories

The following accessories have to be connected to the absorber tube:

Assembling the absorber tube

- Place the sealing ring in the holder.
- Screw the holder in place on the bottom end of the absorber tube.
- Check whether the stopper is in place at the upper end of the absorber tube.


920 Absorber Module

3

Tighten the loose end of the capillary (21-**5**) (connected to port 8 of the 10-port valve) to the bottom of the absorber tube holder using a pressure screw (6.2744.010).

Inserting the absorber tube into the instrument

Insert the absorber tube into the instrument:



- 1


Next steps

Connect the Capillary/tubing at ...	... to
side	Secure the outlet of the combustion oven (<i>see chapter 4.11.2, page 79</i>)

- Attach the other end of the silicone tubing to the right-side drainage support on the base tray.

Plugging in the leak sensor connection cable

920 Absorber Module



- 1 Pull the leak sensor connection cable out of the base tray as far as needed.
- 2 Plug the plug for the leak sensor connection cable into the leak sensor connection socket (labeled **Leak Sensor**).

4.9 Connecting the instrument to a computer



NOTE

If the instrument is connected to the computer, then it must be switched off.

Accessories

For this step, you need the following accessories:

- USB connecting cable (6.2151.020)

Connecting the USB cable

- 1 Insert the USB cable into the computer connection socket on the rear of the instrument labeled *PC*.
- 2 Insert the other end into a USB port on the computer.

4.10 Connecting the instrument to the power grid



WARNING

Electric shock from electrical potential

Risk of injury by touching live components or through moisture on live parts.

- Never open the housing of the instrument while the power cord is still connected.
- Protect live parts (e.g. power supply unit, power cord, connection sockets) against moisture.
- Unplug the power plug immediately if you suspect that moisture has gotten inside the instrument.
- Only personnel who have been issued Metrohm qualifications may perform service and repair work on electrical and electronic parts.

Connecting the power cord

Accessories

Power cord with the following specifications:

- Length: max. 2 m
- Number of cores: 3, with protective conductor
- Instrument plug: IEC 60320 type C13
- Conductor cross-section 3x min. 0.75 mm² / 18 AWG
- Power plug:
 - according to customer requirement (6.2122.XX0)
 - min. 10 A



NOTE

Do not use a not permitted power cord!

1 Plugging in the power cord

- Plug the power cord into the instrument's power socket.
- Connect the power cord to the power grid.

4.11 Connecting the 920 Absorber Module to the remaining components

Combustion Module (AJ)

Three capillaries are routed out at the right side of the Combustion-IC-System.


- The capillary connected to port 4 of the Dosino for ultrapure water is marked with a blue colored sleeve.
- The capillary connected to the outlet of the Metrosep A Trap 1 – 100/4.0 is marked with a yellow colored sleeve.
- The thicker capillary installed at the inlet of the absorber tube is routed out under the door.

Combustion Oven (TEI)

Three capillaries are routed out at the back and/or left side of the Combustion Oven (TEI).

- The capillary connected to port 4 of the Dosino for ultrapure water is marked with a blue colored sleeve. The capillary is part of 6.07311.310. The other end of the capillary is already connected.
- The capillary connected to the outlet of the Metrosep A Trap 1 – 100/4.0 is marked with a yellow colored sleeve. The other end of the capillary is already connected to the oven outlet.
- The thicker capillary installed at the inlet of the absorber tube is routed out at the rear. The other end of the capillary is already connected to the oven outlet.

4.11.1 Connecting the Combustion Module (AJ) inlet: water infeed



1 Feeding the capillary into the Combustion Module (AJ)

2 Shortening the capillary

3 Connecting the capillary

■■■■■■■■ 77




Figure 22 Combustion tube – Inlet


- | | |
|---|--|
| <p>1 Water infeed inlet
Connection for the capillary with the blue colored sleeve (7-14).</p> | <p>2 Flame sensor connection
Connection for the flame sensor.</p> |
| <p>3 Oxygen infeed
Connection for the oxygen infeed.</p> | <p>4 Sample inlet
Solid and liquid samples: inlet for sample boats.

Gaseous samples: connection of the LPG/GSS module.</p> |

- Plug the capillary into a FAST connector (*see Chapter "Replacing tubing connections" in the Combustion Module operating instructions*).
- Carefully attach the FAST connector to the connector (22-**1**).

You can find information on connecting the flame sensor, the oxygen infeed and the ABD in the Combustion Module (AJ) operating instructions.

There are 2 connections on the CIC T connector in the Combustion Module (AJ):



1 Absorber solution connection
Connection for the capillary with the yellow colored sleeve (7-9).

2 T connector CIC (6.7301.210)
T connector on the combustion tube's outlet. For feeding the absorber solution in and removing the combustion gases.

4 Combustion gases connector
Connector for the PTFE capillary, 1/8 in. / 1/16 in. ID / 1.5 m (6.1803.130) (7-8).

Connecting the capillary for the combustion gases

1 Feeding the capillary into the Combustion Module (AJ)

- Push the loose end of the PTFE capillary 1/8 in. / 1/16 in. ID / 1.5 m (6.1803.130) (7-8) through the opening in the right-side door of the 920 Absorber Module.
- Guide the capillary inside through the round opening in the bottom of the Combustion Module (AJ).

2 Shortening the capillary

Shorten the capillary using the capillary cutter so that it can be connected to the CIC T connector's outlet (23-4) easily.

3 Connecting the capillary

- Insert the capillary into an angled FAST connector (6.7304.000).
- Carefully attach the FAST connector to the CIC T connector's outlet (23-**4**).

Connecting the capillary for the absorber solution

1 Feeding the capillary into the Combustion Module (AJ)

Guide the loose end of the capillary with the yellow colored sleeve inside through the round opening at the top in the left-side wall of the Combustion Module (AJ).

2 Shortening the capillary

Shorten the capillary using the capillary cutter so that it can be connected to the CIC T connector's outlet (23-**1**) easily.

3 Connecting the capillary

- Plug the capillary into a straight FAST connector (6.7304.010) (see Chapter "Replacing tubing connections" in the Combustion Module operating instructions).
- Carefully attach the FAST connector to the CIC T connector's outlet (23-**1**).

4.11.3 Connecting the Combustion Oven (TEI) inlet




Figure 24 Combustion tube – Inlet

- | | |
|---|--|
| <p>1 Argon infeed
Connection for the argon infeed.</p> | <p>2 Oxygen infeed
Connection for the oxygen infeed.</p> |
| <p>3 Oxygen infeed
Connection for the oxygen infeed.</p> | <p>4 Water infeed inlet
Connector for the ultrapure water supply.</p> |

4.11.4 Connecting the Combustion Oven (TEI) outlet




Figure 25 Combustion Oven (TEI) – Outlet


- | | |
|---|---|
| <p>1 Absorber solution capillary
Capillary for the infeed of absorber solution at the tube outlet.</p> | <p>2 Combustion gases capillary
Capillary for transporting the combustion gases from the tube outlet to the absorber tube.</p> |
|---|---|

4.11.5 Connecting the 930 Compact IC Flex

Combustion Module (AJ): A capillary is routed out on the left side of the 920 Absorber Module.

Combustion Oven (TEI): A capillary is routed out on the right side of the 920 Absorber Module.

This capillary is connected to the injection valve of the 930 Compact IC Flex.



Accessories

For this step you need:

- Capillary cutter (6.2621.080)

Connecting the 930 Compact IC Flex

- 1
 - Shorten the end of the capillary using the capillary cutter so that it can be connected to the injection valve of the 930 Compact IC Flex easily.
 - Connect the end of the capillary to port 1 of the injection valve (*see manual for 930 Compact IC Flex*).

5 Start-up

All of the Combustion IC system components have to be completely installed before putting the 920 Absorber Module into operation; this means:

- The Combustion Module including ABD and optional sampler (MMS 5000) and, if necessary, the LPG/GSS Module is installed according to the specifications in the Combustion IC operating instructions from Analytik Jena/TEI.
- The Combustion Oven (TEI) including, if applicable, Liquid Autosampler CIC, Solid Autosampler CIC and/or GLS Autosampler CIC from Trace Elemental Instruments is installed according to the operating instructions.
- The 920 Absorber Module is installed according to the specifications in the manual for the 920 Absorber Module and is connected to the Combustion Module (AJ) / Combustion Oven (TEI).
- The 930 Compact IC Flex (2.930.2560) is installed according to the specifications in the manual for the ion chromatograph and is connected to the 920 Absorber Module.
- The ion chromatograph is connected to a computer where the MagIC Net software and TEIS are installed. TEIS is only necessary for the Combustion Oven (TEI).

Putting the Combustion IC system (AJ) into operation

1 Starting the software

- Switch on the computer and start the software MagIC Net.
- Enter user name and password.

The software starts and shows the **Workplace** program part.

2 Switching on the instruments

- Switch on the Combustion Module (AJ): the Combustion Module (AJ) and the MMS 5000 sampler (AJ) or the LPG/GSS Module (AJ), if present.
- Switch on the Combustion-IC-System.
- Switch on the ion chromatograph.
MagIC Net detects all of the instruments and asks if they are to be saved to the device table.
- Confirm all of the prompts with **[Yes]**.



When an instrument is not automatically recognized, then the "basic initialization" has to be run in MagIC Net. Additional information concerning the basic initialization can be found in the MagIC Net online help or in the MagIC Net manual.

Workplace ► Run

- #### 4 Starting the equilibration

- The instrument is ready to measure as soon as equilibration has finished

Putting the Combustion IC system (TEI) into operation

- Switch on the computer and start the MagIC Net software and TEIS.
- Enter the user name and the password in each software.

TEIS: The software starts and shows the **Task Manager** program part.

- Switch on the Combustion Oven (TEI) and the Liquid Autosampler CIC (TEI), Solid Autosampler CIC (TEI) or the GLS Sampler CIC (TEI), if present.
- Switch on the ion chromatograph and the 920 Absorber Module. MagIC Net detects both instruments and asks if they are to be saved to the device table.
- Confirm all of the prompts with **[Yes]**.

3 Loading the method in MagIC Net

Workplace ► Run

- Select the **Equilibration** tab.
- In the **Method** field, load the **Preparation of Combustion IC** method.

4 Starting the equilibration of the ion chromatograph

- Start the equilibration with **[Start HW]**.

The IC instrument is ready to measure as soon as equilibration has finished.

- ### 5
- Start the equilibration of the Combustion Oven (TEI) and the Liquid Autosampler CIC (TEI), Solid Autosampler CIC (TEI) or the GLS Sampler CIC (TEI), if present, in the TEIS software.

6 Operation and maintenance

6.1 General

The following chapters describe the care and maintenance tasks that can be carried out by the user.

6.2 Door



CAUTION

The door is made of PMMA (poly(methyl methacrylate)). It must never be cleaned with abrasive media or solvents.



CAUTION

Never use the door as a handle.

6.3 Dosino

Information regarding the maintenance of the Dosino can be found in the Dosino manual enclosed in your delivery. Alternatively, you can also download it from the Internet at www.metrohm.com.

6.4 Dosing unit

Information regarding the maintenance of the Dosino can be found in the Dosino manual enclosed in your delivery. Alternatively, you can also download it from the Internet at www.metrohm.com.

6.5

Information regarding the maintenance of the Metrosep I Trap 1 can be found in the leaflet on the trap column enclosed in your delivery. Alternatively, you can also download it from the Internet at www.metrohm.com.

6.6

Information regarding the maintenance of the Metrosep A Trap 1 can be found in the leaflet on the trap column enclosed in your delivery. Alternatively, you can also download it from the Internet at www.metrohm.com.

6.7

Check on a regular basis whether the connections of the capillaries are tight.

If you find leaking connections, loosen and reconnect them.

More information about how to handle pressure screws and capillaries can be found on the Internet at ic-help.metrohm.com.

6.8

Check on a regular basis whether the connections of the capillaries are tight.


If you find leaking connections, loosen and reconnect them.

More information about how to handle pressure screws and capillaries can be found on the Internet at ic-help.metrohm.com/.

6.9

The absorber tube requires maintenance if the check valve and the absorber tube holder clog up due to sooting resulting from deflagration in the combustion tube.

The absorber tube holder (6.2766.000) and the inlet valve (6.2824.170) in the check valve must be exchanged.



Exchanging the inlet valve

1 Disconnecting the inlet capillary

Unscrew the pressure screw (6.2744.170) from the check valve.

2 Unscrewing the check valve holder

Unscrew the check valve holder (6.2766.010) from the absorber tube holder.

3 Replacing the inlet valve

- Disassemble the check valve holder by unscrewing it.
- Remove the inlet valve.
- Insert the new inlet valve (see figure).
- Reassemble by screwing the parts of the check valve together.

Exchanging the absorber tube holder

1 Unscrewing the outlet capillary

On the bottom side, unscrew the pressure screw from the absorber tube holder.

2 Unscrewing the check valve holder

Unscrew the check valve holder (6.2766.010) from the absorber tube holder.

- Unscrew the absorber tube holder from the absorber tube.
- Place a new sealing ring in the new absorber tube holder.
- On the bottom, screw the new absorber tube holder to the absorber tube.

Screw the check valve holder with the new inlet valve to the side of the absorber tube holder.

Screw the inlet capillary (6.1803.130) to the check valve again using the pressure screw (6.2744.170).

Tighten the loose end of the capillary (6.1803.160) (connected to port 8 of the 10-port valve) to the bottom of the absorber tube holder using a pressure screw (6.2744.010).

8 Accessories

Up-to-date information on the scope of delivery and optional accessories for your product can be found on the Internet. You can download this information using the article number as follows:

Downloading the accessories list

- 1 Enter <https://www.metrohm.com/> into your Internet browser.
- 2 Enter the article number (e.g. **2.920.0010**) into the search field.
The search result is displayed.
- 3 Click on the product.
Detailed information regarding the product is shown on various tabs.
- 4 On the **Included parts** tab, click on **Download the PDF**.
The PDF file with the accessories data is created.



NOTE

Once you have received your new product, we recommend downloading the accessories list from the Internet, printing it out and keeping it together with the manual for reference purposes.

IP degree of protection IP 20

Controls

Indicators LED for standby indicator

On/off switch On the rear of the instrument

9.4 Dosino

Please refer to the manual for the Dosino enclosed in the delivery for information regarding the technical specifications of the Dosino. Alternatively, you can also download it from the Internet at www.metrohm.com.

9.5 10-port valve

Connections 1 out of 10

Resistance to solvents No restrictions

Switching time typ. 100 ms

Pressure resistance 2.0 MPa (20 bar)
Valve function prevents damage at overpressure

9.6 6-port valve

Actuator switching time typ. 100 ms

Maximum operating pressure 35 MPa (350 bar)

Material PEEK

Index

6-port valve	62
Technical specifications	95
Valve position: Fill	67
Valve position: Inject	67
10-port valve	55
Technical specifications	95
930 Compact IC Flex	
Connect	82

A

Ambient conditions	94
--------------------------	----

C

CIC	
See also Combustion IC	7
Combustion IC	
Definition	7
Detectable substances	7
Combustion Module	
Connect	76
Computer connection	74
Connect	
Power grid	75
To computer	74

D

Detector	
Interface	96
Determination	
Sequence	8
Dimensions	94
Door	88
Dosino	
For absorber solution	39
For ultrapure water	19
Drainage tubing	
Installation	72

E

Electrostatic charge	5
----------------------------	---

Energy supply	96
---------------------	----

F

Fill	
6-port valve: valve position .	67
Frequency	96

H

Housing	94
Humidity	94

I

Inject	
6-port valve: Valve position .	67
Installation	
6-port valve	62
10-port valve	55
Dosino absorber solution	39
Dosino ultrapure water	19
Drainage tubing	72
Leak sensor	73
Overview	16
Interface	
MSB	96
USB	96
Interfaces	96
Extension Module	96
Leak sensor	96

L

Leak sensor	
Installation	73
Interface	96
Technical specifications	96

M

Material	94
MSB	96

O

Operation	94
Overvoltage category	94

P

Power connection	75
Power consumption	96
Power supply unit	96

R

Reference conditions	94
----------------------------	----

S

Safety instructions	4
Sea level	94
Service	4
Storage	94
Supply voltage	4
System description	8

T

Technical data	
Reference conditions	94
Technical specifications	
6-port valve	95
10-port valve	95
Detector	96
Interfaces	96
Leak sensor	96
Temperature	94

U

USB	96
-----------	----

V

Voltage	96
---------------	----