

# SpecSuite

SPECTROSCOPY SOFTWARE

## SpecSuite Installation and User Manual





Metrohm Spectro, Inc.  
313 Enterprise Dr  
Plainsboro Township, NJ 08536  
(732) 230-1601  
[www.metrohm.com](http://www.metrohm.com)

# SpecSuite

SpecSuite Software Version 1.1.11

Tutorial

Metrohm Spectro, Inc.  
Plainsboro Township, NJ 08536

This documentation is an original document.

This documentation has been prepared with great care. However, errors cannot be eliminated. Please send comments regarding possible errors to the address above.

### **Copyright**

This documentation is protected by copyright. All rights reserved.

### **Disclaimer**

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g., conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

# Contents

1. Overview .....	11
1.1 Introduction.....	11
1.2 Licensing.....	11
Trial license.....	12
2. Installation and Setup .....	12
2.1 System Requirements.....	12
2.2 Installing SpecSuite.....	13
2.3 Setup Local SQL Express Database .....	17
To connect to a different local server .....	18
To create a new local database .....	19
2.4 First Time Login.....	20
2.5 Software Setup .....	20
Import Licenses for SpecSuite and Metrohm Libraries .....	22
Configure Probe (Ratio3).....	22
Import Performance Validation Reference file .....	23
Import Metrohm Libraries (optional).....	23
3. SpecSuite Login Screen .....	25
4. Explore .....	25
4.1 File Ribbon .....	27
Export Spectral File Format Types:.....	28
Import Spectral File Format Types: .....	30
4.2 Acquisition Settings .....	30

4.3	Acquisition Control Panel .....	32
4.4	Acquisition Controls.....	32
4.5	Process Ribbon .....	34
4.5.1	Peak Analysis .....	36
4.6	Camera.....	38
4.7	Explore Chart.....	39
4.8	Data List.....	41
5.	Instrument .....	42
5.1	Device Info .....	43
5.2	Hardware Setup.....	45
	Detector Settings.....	45
	Probe.....	46
	Input & Output .....	48
	Camera.....	49
5.3	Laser Setup.....	50
5.4	Performance .....	53
	PVR File Creation .....	54
	Performing a Performance Validation .....	56
5.5	Calibration.....	57
	Raman Shift Calibration.....	58
	Laser Wavelength Calibration .....	59
	Relative Intensity Correction.....	61
5.6	Instrument Status Messages.....	62
	Performance Status.....	62
	Laser Key and Interlock .....	62

Connected Instruments .....	63
5.7 Virtual Device .....	64
6. OP (Operating Procedure) .....	66
6.1 OP .....	66
6.2 Building an Operating Procedure .....	67
Acquisition .....	67
Processing .....	69
Identification .....	71
Routine Analysis .....	74
Output .....	75
6.3 Method .....	76
6.4 Variables .....	78
Basic Variable .....	79
Model Variable .....	84
7. Identification .....	84
7.1 Online ID .....	85
How to configure the identification parameters .....	85
How to perform an ONLINE ID .....	86
7.2 Offline ID .....	89
7.3 Re Run .....	90
8. Library .....	90
8.1 Library List .....	91
8.2 Creating a New Library .....	92
8.3 Adding Samples .....	94
Import from Library .....	94

Load from File.....	95
Collect New.....	96
Library Samples Info.....	96
<b>8.4</b> Library History .....	97
9. Model.....	98
<b>9.1</b> Select Samples.....	99
Adding Samples.....	99
Input y-data .....	101
Split the Data.....	103
<b>9.2</b> Parameterize the Model .....	106
<b>9.3</b> Build and Validate the Model .....	107
Reading the Loading Plot .....	108
Reading the Score Plot.....	109
Reading the Influence Plot.....	111
<b>9.4</b> Start RMD (Raman Model Development) .....	112
<b>9.5</b> Exporting Model Data .....	112
10. Routine Analysis .....	114
<b>10.1</b> Preparation .....	114
Offline Routine Analysis.....	115
<b>10.2</b> Basic Operations .....	116
Meta Data input .....	116
Start a Run.....	116
<b>10.3</b> Spectra Plot and Trend Plot.....	118
11. Reports.....	120
<b>11.1</b> Acquisition Reports.....	120

11.2	Identification Reports.....	122
11.3	Routine Analysis Reports .....	124
11.4	Performance Validation Reports .....	127
11.5	Calibration records.....	129
11.6	E-Signatures.....	131
11.7	Bulk Export.....	134
12.	Audit Trail.....	135
12.1	Audit Trail Filtering.....	136
12.2	Audit Trail Reporting.....	139
13.	Settings.....	140
13.1	General .....	140
	Workspace .....	140
	Display.....	140
13.2	License.....	142
13.3	Database.....	144
13.4	Backup and Restore .....	145
	Database Size .....	145
	Backup Database.....	146
	Restore Database .....	146
13.5	Advanced .....	147
	Advanced Settings .....	148
	IT Settings.....	148
13.6	User.....	149
	Default Accounts.....	150
	New Accounts .....	151

User Level Permissions.....	152
<b>13.7</b> About .....	158
14. Terminologies.....	159
15. Appendix.....	161
A. Manually creating a server instance (Optional, for IT only).....	161
B. Peak Analysis Definition .....	168
C. SPC data definition .....	174







By default, a server named "ENTERPRISE" is created locally, and a database named 'SpecSuiteDB' is created on the server. IT users can create multiple databases. Each MS SQL Server Express database created on the server has a size limit of 10GB.



For addition information and installation requirements about MS SQL Server 2022 refer to Microsoft website [SQL Server Installation Guide - SQL Server | Microsoft Learn](#)


## 2.2 Installing SpecSuite

1. Ensure the system PC fulfills the minimum system requirements.
2. Ensure Windows administrative rights are granted before installing SpecSuite.
3. Download the compressed file for the release version of SpecSuite. Extract the files, then open the 'SpecSuite\_[version]' folder. Right-click setup.exe and select Run as Administrator. The InstallShield Wizard for SpecSuite will appear.




NOTE: The location of the installation folder is not important, as the installer copies data to a temporary folder on the PC: C:\SQL2022TempFolder. Please consult with your IT department to confirm Write-access at this location.


4. Click Next to continue the installation.




5. When the Metrohm Software License Agreement appears, choose "I accept the terms in the license agreement" then press Next.




NOTE: If SpecSuite is already installed with the same version, the Program Maintenance options will appear. Select Modify or Repair, then Next to continue.



6. The installer will install SpecSuite to the default location: **C:\Metrohm\SpecSuite\**. Use this default location and click Next.




7. Under Setup Type, select Complete, then click Next.



8. Click Install to begin installation. The installation should take approximately 1 minute.








## To connect to a different local server

1. Log in to SpecSuite with an IT account.
2. Go to Settings > Database.
3. In the Server Name field, enter a known server name to connect to.
4. Press Test Connection to see if the Server is available.
5. If successful, press Connection.
  - a. SpecSuite must restart to finish connecting to the new server.
6. Press Yes on the restart prompt.




NOTE: SpecSuite is only compatible with MS SQL Server 2022. Performance on other versions is not guaranteed.



## To create a new local database

In some cases, the user may not want to use the default database connection. SpecSuite allows users to create a new database inside the selected server.

1. Log in with an IT account.
2. Select Settings > Database tab.
3. Select the server instance where the database will be created.
4. Enter a new database name.
5. Press the Create New Database button:



6. Once a new database is created, choose the name from the Choose Database dropdown.
7. Click Connect.





If a password reset is required, the user will receive a prompt to update their password.

- Usernames are unique for all clients connected to the same database.
- A user can only be logged into one instance of SpecSuite.
  - If a user tries to log in to an account while this account is already logged in elsewhere, the new user will receive a notification that the account is already in use.





6. Enter Probe Name and Probe SN.
7. Select the Select button.
8. Locate the Ratio3 TXT file for your instrument and probe configuration.
  - a. Ratio3 files are unique for each instrument and probe configuration. If a probe is changed, a new Ratio3 file must be loaded or created.
9. Select Save Probe.
10. Verify the Probe file is in the list, and the check mark is selected.

### **Import Performance Validation Reference file**

The Performance Validation Reference (PVR) file is used to validate the instrument's performance. Instruments could be delivered with a PVR file designated for the instrument. This file should be imported to run validation.

If there are changes to the instrument or probe, a new PVR file may be required.

#### **Import PVR file**

1. Turn on the instrument and connect it to the computer.
2. Log in to the ADMIN account.
3. Select Instrument.
4. Select Performance.
5. Select Import PVR.
6. Locate the PVR file for your instrument and probe configuration.
  - a. \*.PVR file type
  - b. Instrument SN, Probe SN, and time of creation are part of the file name
  - c. Example: 19540010.26001-230105959-26-01-01.pvr
7. Select Open.
8. Select OK to clear the successful import message.


### **Import Metrohm Libraries (optional)**

A license for the Metrohm Illicit and General Chemicals library is automatically installed with SpecSuite. The Metrohm Illicit and General Chemicals library package must be imported to display and activate it for use in measurements.






### 3. SpecSuite Login Screen



The Login window appears when SpecSuite is launched. After logging in, SpecSuite opens to different modules based on the user's account type. IT and Administrator account logins open to the Settings module. Developer account logins open to the Explore module, and Operator account logins open to the Operating Procedure module.


Refer to Section 2.4 and 2.5 for additional login information.

### 4. Explore



Explore is the most versatile module in SpecSuite, designed to collect and process data freely. Some applications for the explore workspace include:

- Experimenting with different data processing.
- Focusing the laser on a sample if using a probe stand or microscope.
- Viewing and processing existing data.



**1 Explore**

**3 Acquisition controls**

**5 Laser controls**

**7 Data List**


**2 Import and Export**

**4 Display controls**

**6 Spectrum plot**


**8 Acquisition control panel**

## 4.1 File Ribbon



In the top File ribbon, the user may:

- **Import** – Imports a file from the local hard drive of their computer. Imported data are not saved in the database.
  - Allowed file types are: .CSV, .TXTR, .TXT and .SPC.
- **Export** – Exports the selected spectral data file(s).




- **Export file types are:** .TXTR, .CSV (Raw CSV), .CSV (Simple CSV), and .SPC (Simple SPC).






2. Wavelength (nm): When selected, the x-axis data will be interpolated from pixel to wavelength.
3. Wavenumber ( $\text{cm}^{-1}$ ): When selected, the x-axis data will be interpolated from pixel to wavelength, then converted to absolute wavenumber.
4. Raman shift ( $\text{cm}^{-1}$ ): When selected, x-axis data will be interpolated from pixel to wavelength, then converted to absolute wavenumber, then converted to Raman shift based on the instrument's laser wavelength with a default  $4 \text{ cm}^{-1}$  interpolation and default Raman shift range based on the instrument model.




5. Raman Shift Linear ( $\text{cm}^{-1}$ ): When selected, x-axis data will be interpolated from pixel to wavelength, then converted to absolute wavenumber, then converted to Raman shift based on the instrument's laser wavelength, with the start/finish Raman shift and interpolation interval set by the user in the Export Options dialog.

### Export Options



The available Y-axis are:



1. Raw Data: Spectral intensity directly recorded on the detector, without any processing.



2. Dark: Spectral intensity recorded when the laser is turned off (with the same parameters as the laser-on scan) with no additional processing.
3. Dark Subtracted: Raw Data minus the Dark data.
4. Processed: Processed data is the relative intensity corrected spectra data with all the additional processing (smoothing, derivative, normalization, baseline correction) that the user applied on the spectra entry.

## SPC (Simple SPC)

Simple SPC is a simple 2-column standard SPC format that follows the same axis convention as simple CSV.




Simple CSV and Simple SPC files are not available for the i-Raman Duo.

## Import Spectral File Format Types:

SpecSuite supports all file types mentioned in section 4.1.1 and all TXTR, CSV, and SPC file types from BWSpec 4.15. Multiple different file types can be imported simultaneously.

## 4.2 Acquisition Settings


The Explore page allows the user to adjust acquisition parameters to fine-tune data collection.



Acquisition parameters are set in the right panel of the chart. Under Spectrometer Controls, the user can set:

- Integration Time (range depends on connected spectrometer) and Unit (ms, sec, min)
- Laser power (1- 100%)
- Average

When an i-Raman Duo or DIY Extended Range is connected, SpecSuite supports acquisitions with up to 2 lasers and 1 spectrometer to cover a wide Raman shift range. When two lasers are used, settings for both lasers will appear in the acquisition controls. Laser power and integration time may be different for both. When two lasers are enabled, a new column labeled "Cross Point" appears in the data list. The Cross Point is the point at which the spectra from each laser are concatenated together to form one full spectrum.



The units of each axis can be changed in the same panel.


X-axis options:





- Pixel
- Wavelength (nm)
- Wavenumber (cm<sup>-1</sup>)
- Raman Shift (cm<sup>-1</sup>)

Y-axis options:

- Raw Data
- Dark
- Dark Subtracted
- Processed







-  **Scan** – When Scan is pressed, a spectrum is collected using the Integration Time, Laser Power, and Averages set in Spectrometer Controls, or automatically optimized if Auto Time is turned on. The scan appears in the data list below the chart, and a trace of the data appears in the chart area with the selected x- and y-axes units.
-  **Dark Scan** – The Dark Scan button collects a spectrum using the same Integration Time as the current scan record. The laser is automatically turned off during a dark scan. This dark data is subtracted from the Scan data to produce the dark-subtracted spectrum. The dark data will be used for all subsequent acquisitions at the same integration time. When the integration time is changed, a new dark scan is required. In the Acquisition Control Panel, the user can enable Auto Dark, which automatically collects a dark scan after the data scan.
-  **Check Signal** – The Check Signal button helps the user adjust the fiber probe position to maximize the signal. Pressing Check Signal will begin continuously taking data scans, replacing the previous results until the measurement is either paused or stopped. These Check Signal scans are not saved.
-  **Pause** – The pause button appears during a multi-acquisition run. When pressed, it allows the current scan to finish before pausing all scans until the Scan button is pressed





-  **Smoothing** – Sets the type and parameters for smoothing the data. More info about smoothing can be found in section 6.2.



Smoothing Type

Savitzky-Golay (SG)

Savitzky-Golay (SG)


Boxcar

FFT (Fast Fourier Transform)

Polynomial Order (1-5)

3

Apply Revert

-  **Derivative** – Sets the type of differentiation to use on the data and sets the order or other parameters.



Derivative Type

Savitzky-Golay Derivative

Derivative Order (1-2)

1


Width


15

Polynomial Order (1-5)

3

Apply Revert

-  **Normalize** – The intensity of the spectra is scaled so that spectra can be compared more easily. Normalization types are Scaling (0-1), Mean, and SNV.




Normalization Type

Scaling (0-1)

Scaling (0-1)

Mean

SNV


-  **Baseline Correction** – A standard algorithm optimized for Raman spectra can correct the baseline of the data.



Baseline Correction Type

Standard

Apply Revert

-  **Peaks** – Peaks can be marked with their relevant values. The qualifications for peaks can be adjusted to include or exclude based on absolute intensity and other criteria.


Processing actions are configured in the Process section of the right-side panel. Each processing step is recorded in the Process History column of the data table in the order in which they occur.


### 4.5.1 Peak Analysis

Information about a spectrum's peaks can be shown live in the Explore chart. Peak analysis can be performed on online and offline data.

It is highly recommended to apply baseline correction before performing Peak Analysis.

For more information about how the SpecSuite Peak Analysis feature works, see Appendix B.




Selecting  Peaks on the Process ribbon switches the Explore chart to the Peak Analysis view. The table beneath the chart will list the peaks in the selected spectrum. Position, Intensity, and FWHM are the default columns in the peak table. While the Peak Analysis view is open, File actions, Scan, Dark scan, Check Signal, Reference, Scale, Clear all, and Laser control are disabled.

## Perform Peak Analysis

1. Select a spectrum from Data List.
2. Click Peaks to open the Peak Analysis view.
3. Configure Peak analysis settings in the right-side panel.
4. Press Find Peaks button to perform calculations.

Parameters can be adjusted to the right of the chart in the Peak Parameters section. Peak Type is the algorithm used to calculate which peaks are present.

- Peak analysis is applicable to the **ONE** active spectrum only
- Peak Type options include Data Points, Parabolic, and Curve Fit
  - For Data Points: Peak height is required
  - For Parabolic: Peak height and peak width are required
  - For Curve Fit: Peak height and peak width are required. Initial peaks are required. The user must specify the fit type, such as Gaussian, Lorentzian, or Voigt. If Voigt is selected, the user must also specify the Voigt ratio or use the default value 0.5.



Peak Display determines which peak information is displayed for the calculated peaks.


Checking the Display Peak Labels box displays annotations for calculated peaks in the chart.

Right-clicking an unwanted peak in the table offers the option to delete it from the chart and table.


Clicking the Clear button resets the peak analysis view; peak annotations are removed, and the peak table is emptied.

## 4.6 Camera

When a camera is detected and enabled under Instrument -> Hardware setup, a camera icon will appear in the main Explore menu.



Click on the camera icon to open the camera configured in Hardware setup. This will open a window within the chart area that displays a live feed from the camera.






The right-click menu of the camera image allows the user to:


- Save the displayed image as a .PNG, .JPG or .BMP file (Save Image As...),
- Save up to 3 images to the selected spectrum record in the database (Save to Database),  
and
- Use the mouse to position crosshairs on the camera image (Add Crosshair) or remove them.

When a camera is checked under Hardware Setup, another toggle appears in the Acquisition Controls for Illumination LED. This allows the user to turn the microscope light on and off through the software.



## 4.7 Explore Chart

The explore chart provides tools for visualizing spectral data.



- **Reset** – Resets the x and y axes to the scaling set in the Scale window.
- **Clear** – Clears all the spectra from the chart and the table below.
- **Display** – Opens the Display window where the user can change the x and y axis scale to choose the size of the chart view. They can also set chart display options, such as the chart Legend.

Display
✕


Scale

	Start	End
X-Axis	<input type="text" value="97.34"/>	<input type="text" value="2799.98"/>
Y-Axis	<input type="text" value="-2760.45"/>	<input type="text" value="69193.90"/>

Display


Show Legend

Chart Scale will reset if the X-axis unit and Y-axis units are changed. The chart legend will appear in the top-right corner of the screen.



## 4.8 Data List

Name	Date And Time	Integration Time	Laser power	Wavelength	Average	Color
blue goldstone	05/06/2024 11:51:...	14259 ms	100.00%			#FF008E8A
vitamin_c_simplespc	04/15/2024 01:57:...	1915 ms	100.00%			DarkRed
Amoxicillin	10/17/2025 11:25:...	0 us	100.00%	0.00		#FF2872B3
Naproxen Sodium	10/17/2025 11:25:...	0 us	100.00%	0.00		#FF932279
Prednisone	10/17/2025 11:25:...	0 us	100.00%	0.00		#FF371A09
Cellulose	10/17/2025 11:26:...	2487 ms	100.00%	784.98	1	#FF630239
Creatine	10/17/2025 11:26:...	17 sec	100.00%	784.95	1	#FF0814A6
Fibranol	10/17/2025 11:26:...	715485 us	100.00%	784.94	1	#FF4616EE

The Data List is the table of samples beneath the chart. By default, the columns in the data list are Visibility , Name, Date & Time, Integration Time, Laser power, Wavelength, Average, Color, and Process History.




## 5.1 Device Info

Hardware Information	
VID and PID	16A3-2EC8
Manufacture Date	7/1/2025 12:00:00 AM -0400
SysModel	BWS478-1064-HT
SpModel	BTC284N
Spectrometer FN	VIR
Instrument SN	_19540050.25.03001
Pixel Number	512
Min Int. Time	200 $\mu$ s
X-Axis Min	46
X-Axis Max	497
Detector Temp	-20 °C
Chamber Temp	30 °C
Laser Type	I2C
Laser Wavelength	1064.0681 nm

The Device Info tab displays the basic hardware information of the attached system or spectrometer:

- System Model
- SP model (Spectrometer Model)
- Spectrometer FN (Factory Number)
- Instrument SN
- Laser Type and Wavelength



Show Details reveals the entire contents of the EEPROM laid out in an easy-to-read way, including:

- Hardware Information
- Coefficients
- Firmware version – for BTC66xN series

Read Flash reads the EEPROM and refreshes the information displayed in the Device Info table, and resets all system/spectrometer settings to their defaults. The Read Flash button appears only when Detailed device info is shown.

Pressing Hide Details will show only Hardware info and hide the Read Flash option.




During Service, the instrument EEPROM might be adjusted.

It is recommended to click the Read Flash button to refresh the EEPROM information after service.




## 5.2 Hardware Setup



Hardware Setup contains options for Detector Settings, Probe, Input & Output, and Camera.


### Detector Settings



Detector Temperature






3. Select the Ratio3 file created for the instrument and probe combination.
  - a. The contents of the Ratio file are displayed after the file is loaded.
  - b. Instruments will have a Ratio3 file provided by Metrohm.
  - c. If you do not have a Ratio3 file, request one from your local service representative.
  - d. Experienced customers can follow the procedure in section 5.5.3 Relative Intensity Correction to create their own Ratio3 file. This process requires the appropriate Standard Reference Material to generate a Ratio3 relative intensity correction file.
4. Press  Save to add the new probe to the Probe list.
  - a. If any requirements are not met, an error will appear below the Save Probe button with information about why the probe was not saved.

For DIY setups, there is a button for the user to load a default Ratio3 with all 1s to create a probe record for functions that require a Ratio3 file, like Performance. Using the default Ratio3 file allows data to be collected without having intensity correction applied.

For i-Raman systems, this button is not shown, and the user must load a real Ratio3 file.



^ Probe

New Edit Cancel Save Remove

Apply	Probe Name	Probe SN
<input type="checkbox"/>	DIY Probe	123123123
<input checked="" type="checkbox"/>	default probe	111111111

Probe Name

Probe SN

Relative Intensity Correction Ratio File (Ratio3\_XXX.txt):

OR

Pixel	Ratio3
-------	--------

Success : "default probe" was saved to the Database






### Digital Output (controllable)

- Digital Output1 – Illumination LED for system enabled by default (inverse and dark unchecked)
- Digital Outputs 2-5 – User definable for DIY setups


Please consult the spectrometer user manual for your device for information about the specific Input and Output options for your hardware, or reach out to a Metrohm Service representative

## Camera


The camera section allows the user to select a BAC151C microscope camera or an alternate camera for image/video while taking a scan in Explore.



Users can click  to refresh the drop-down list of available camera modules connected to the PC. Users can also click  to open a preview of the camera feed. Once

enabled, a camera icon  appears in the top ribbon of the Explore module.

## 5.3 Laser Setup



The "Laser Setup" tab allows the user to view pre-set laser parameters, make adjustments to the laser configuration, or set up a laser manually. Lasers used in i-Raman systems are automatically detected and configured. **NOTE:** i-Raman systems do not allow another laser to be added or subtracted from a laser setup.

When a standalone laser device is connected, the user can add or remove a laser from this screen.

Add Laser – Adds a new laser to be configured for this instrument.

Remove Laser – Removes the laser from the configuration table for this instrument.

To allow SpecSuite to control a laser, a compatible laser must be used. The user may select the type of laser from the list:

Analog & I2C are for i-Raman systems and are automatically set by the software based on the connected i-Raman model.

- **Analog** (for i-Raman system 532nm)


- I2C (for i-Raman NxG 785 and 1064, and i-Raman Plus 785)

USB and Serial USB/RS232 are reserved for DIY setup

- USB (for DIY – GLM/Cleanlaze laser USB with the BWT driver)
- Serial USB/RS232 (for DIY - IPS serial USB communication through COM port)


For Serial USB/RS232, the user will find the available IPS laser identified in the drop-down if the hardware is connected correctly. If no laser is available, check connections and troubleshoot by confirming operation in the laser user manual.

For an IPS laser, ensure that Laser Enable and PWM Control are selected while the Laser Enable signal is working on the laser device.



Only up to two lasers can be controlled by SpecSuite at one time. Add a second laser and input the laser wavelength to perform Raman Concatenation.

The user can manually set the laser wavelength. The Set button must be pressed for the laser setting to be saved. If the Wavelength has been set correctly, a message will appear saying “Laser Wavelength updated and saved to database”.



The laser can be temporarily disabled but unchecking the box in the top left of the configured laser. This will stop control in the software without losing the configuration setting.


The Manual laser mode is intended to control the laser externally, while still being able to produce calibrated Raman data. Only the Laser Wavelength is required for input.

Other available Laser Settings:

Open Laser Before Scan – This setting is enabled by default. When this box is checked, the laser will automatically turn on before sample data is collected if it is not already on. When the box is unchecked, the laser will not automatically turn on if it is off when the scan starts. In this case, the laser must be manually turned on before the measurement is started.

Close Laser After Scan – This is enabled by default. When this box is checked, the laser will automatically turn off after the sample data is collected. If the box is unchecked, the laser will continue firing even after the data scan is complete and remain so until the laser is manually turned off.

## 5.4 Performance



The Performance tab provides options to:


- Perform a Performance Validation (PV) to ensure the system is functioning correctly
- Create a new Performance Validation Reference (PVR) file if needed

The PVR file is unique for the instrument, probe, and material being sampled. If any of these parts change, a new PVR file should be created.

- If the reference material changes, a new PVR file should be made.
- If the system's hardware is changed or recalibrated, a new PVR file should be made.
- If a new probe is used, a PVR file for the specific instrument and probe combination should be used.

Even if the product is not used in a regulated environment, a performance validation should be performed routinely to ensure the instrument's reliability. i-Raman systems ship with a factory-created PVR file for instrument validation.

## PVR File Creation



### 1 Create PVR

Collect a spectrum and create a PVR file.

### 3 Performance Settings

Settings used for the Performance Test. Default values are loaded based on the Instrument information.

### 5 Reference Material

Test material to measure Raman shift calibration.

### 7 Laser Power & Average

### 2 Reference Spectrum

Displays the spectrum when the PVR file is created.

### 4 Instrument Information


Software automatically detects the connected instrument for this information.

### 6 Probe SN

Select the probe serial number from the list of created probes in Hardware Setup.

## How to create a PVR file

1. In Instrument > Performance, click the 'Create PVR' button, and the PVR creation window will appear.
  - a. To create a PVR file, an instrument must be online.
  - b. The connected system will automatically be identified and selected by the Instrument Type and Instrument Model fields (4).
2. Select the Reference Material used for the Reference file: Polystyrene, Acetaminophen, or Cyclohexane (5).
3. Select Probe by serial number (6).
  - a. A probe must be created in Hardware setup prior to this operation.
4. Average and Laser Power can be set (7).
  - a. Recommended defaults are 100% laser power and 3 averages.
5. All parameters are set to default values, but some can be changed (3).
  - a. A Performance Precheck can be run to assess the instrument performance prior to creating the PVR file. If checked, the precheck will run: Dark1 and Dark2 (dark scans taken with the minimum integration time), and Dark(60s) scans. If not selected, these dark scans will be skipped to shorten the time to complete PVR file creation.
  - b. If the box next to Electronic Bias is checked, the validation test will do a 60-second dark scan to assess the instrument performance. It is recommended to always create a PVR file with all testing options selected to fully validate the instrument's performance and ensure correct results.
6. Place the selected Reference Material in front of the sample probe and open the probe shutter.
7. Press Create PVR (1).
  - a. During data collection, click the Stop button if the acquisition needs to be canceled.
  - b. A progress bar and status message will appear under the Accept/Reject buttons to indicate the status of the PVR creation process.
8. The reference spectrum can be viewed by switching to the Reference Spectrum tab (2).
9. Press 'Accept' to save the PVR file.



## Performing a Performance Validation

The screenshot shows the SpecSuite Performance Validation window. At the top, a green bar indicates 'Performance: Pass' and 'ADMIN'. Below this, a dropdown menu shows 'PVR-19540010.25.05.001-123456-20251201-155903'. Three buttons are visible: 'Export PVR', 'Import PVR', and 'Create PVR'. The main area is divided into two columns. The left column contains settings for 'Reference Material: Polystyrene', 'Spectral Range: 150 - 3400 (cm-1)', 'Integration Time: 233 (ms)', 'Laser Power: 100 (%)', and 'Average: 3'. The right column contains 'Instrument Model: BWS478-532H-HT', 'Instrument SN: 19540010.25.05.001', 'Spectrometer FN: VDG', 'Probe SN: 123456', and 'Created: 12/01/2025 03:59:03 PM'. A 'Probe Name' dropdown is set to '012345678'. Below these settings, a 'Result: Pass' indicator is shown. A table of results is displayed, with columns for 'Item', 'Low Limit', 'High Limit', 'Measured', and 'Result'. The table contains several rows of data, including 'Spectral Range', 'Electronic Bias', 'Readout Noise', 'Dark Current', 'Raman Peak Width @1001.40 (cm-1)', and 'Intensity Variatio...'. To the right of the table, a 'Performance Validation Finished' progress bar shows '100%'. Below the progress bar is a Raman spectrum plot with 'Processed' on the y-axis (0 to 70000) and 'Raman Shift (cm-1)' on the x-axis (0 to 3000). The plot shows a reference spectrum (blue line) and a validation spectrum (red line) overlaid. A status bar at the bottom shows 'Laser Key: ON' and 'Interlock: Closed'.

1 Select PVR

2 Export PVR

3 Import PVR

4 Create PVR

5 PVR settings

Parameters of the selected PVR file

6 Probe Select

Probe used for validation

7 Result area

8 Results table

Detailed results for each parameter

9 Start / Stop

10 Spectra


Shows the reference spectrum with the validation spectrum overlaid

The Pass/Fail result from the Performance Validation will update the Performance Indicator in the top bar of the SpecSuite window.


Pressing the Stop button will abort the Performance Validation (9).

A history of Performance Validations can be found in the Reports module. See section 11.4 Performance Validation Reports.

## How to run a Performance Validation

1. Select a PVR file from the 'Select a PVR File...' dropdown (1).
  - a. By default, the loaded PVR file will be the last used file.
  - b. Only PVR files for the connected instrument will display.
2. Verify the selected PVR file's parameters are correct in the PVR setting section (5).
3. Select the active probe SN (6).
  - a. The same probe used for pvr file creation should be used for performance validation.
4. Set up the instrument to measure the Reference Material.
5. Click the  Start button (9).
  - a. When the PVR file is using the full settings, the Performance Validation will run for approximately 90 seconds.
  - b. A Pass or Fail result will appear in the Result area (7).
  - c. The Results table (8), and the Reference spectrum and Performance spectrum (10) are created after the validation. The Results table contains detailed information about the Performance Validation results

## 5.5 Calibration






from the specification, it is too large for the Raman shift calibration tool to correct, and the instrument may require service.

Once this calibration is successful, the user can choose to accept or reject the new calibration coefficients. Accepted coefficients will be saved to the database with other instrument-related parameters.

### Run Raman Shift Calibration

1. Select Reference Material from the drop-down menu (1).
2. Set up the instrument to measure the Reference Material.
3. Select Calibrate (2).
4. Review the result of the calibration (5).
5. Select Accept or Reject (4).
  - a. Accepted new coefficients will be saved to the database with other instrument-related parameters.



### Laser Wavelength Calibration


Laser wavelength calibration uses cyclohexane as a standard to calibrate the laser wavelength value. The cyclohexane peaks are used to determine the exact wavelength of the laser. Laser wavelength calibration should be done after installing a new laser.

Laser wavelength can also be entered manually in Hardware Setup > Laser Setup tab.


### To run Laser Wavelength Calibration

1. Set up the instrument to sample cyclohexane.

2. Press Calibrate.
3. Select Accept to save the new wavelength value.
  - a. Or Reject to not save the change.



## Relative Intensity Correction



Metrohm's i-Raman systems ship with a factory-created relative intensity correction file.

The reference material required is determined by the laser wavelength of the instrument:

- 785nm: SRM 2241
- 532nm: SRM 2242a
- 1064nm: SRM 2244

The system will calibrate the Instrument and Probe with the standard.

### To run Relative Intensity Correction with a new Probe

1. Set up the instrument to measure the corresponding reference material.
2. Click **New**, then enter a probe name
3. Select Standard Material from drop down menu.
4. Select Calibrate.
5. If the calibration is successful, click **Accept**.
6. Select Export the Ratio3.
7. Save the Ratio3 file.



If the Laser key is off, the laser button will not appear in the top bar, and the laser will not appear in the laser setup. A system cannot be treated or used as a lone spectrometer.

Interlock does not have a delay in turning the laser back on and will not affect the laser's availability.




It is recommended to set the laser key to ON before connecting the i-Raman instrument to the PC.


Turning the laser key on after connecting to the PC may cause the software to become unresponsive. Recycle the instrument power to reconnect and continue using the instrument.

## Connected Instruments

The bottom left corner of the window shows the status of the connected instrument(s).



## 5.7 Virtual Device



NOTE: The Virtual Device tab will only appear when an Administrator or Service-level login is used.

Virtual Device connection enables all SpecSuite features without connecting a physical instrument. To enable a System or Spectrometer, check the Enable Virtual Device box, then select the desired system from the dropdown menu. Press Activate. SpecSuite will display the chosen system as Online (Virtual) in the bottom-left corner of the main window. Changes can be made using the Hardware Setup and Laser Setup tabs.

List of supported Virtual Devices:

- i-Raman NxG BWS478-785H-HT
- i-Raman NxG BWS478-785H-HT
- i-Raman NxG BWS478-785H-HT
- i-Raman Duo BWS478-860\_1064-HT

List of available samples for each virtual device:

- Polystyrene
- Cyclohexane
- Acetone
- Acetaminophen




### To enable a Virtual Device

1. Select the intended system from the dropdown menu.
2. Press Activate.
3. SpecSuite will identify the chosen system as Online with the displayed model ending in VIR in the bottom left of the main window. Changes can be made using the Hardware Setup and Laser Setup tabs.



## 6.2 Building an Operating Procedure


To create a New OP, press the New **+** button. On the right side of the page, the OP information can be viewed and edited.



OPs have two necessary components: Acquisition and Output. Other components can be added: Processing, Identification, and Routine Analysis.

Users can update the OP name at the top of the right section.

### Acquisition



The screenshot shows the SpecSuite software interface. The top bar includes 'SpecSuite', 'Performance : Fail', and 'ADMIN'. The main window is divided into several sections:


- Left Panel:** A vertical toolbar with icons for New, Edit, Save, Archive, and Withdraw. Below it is a table of existing operations.
- Table:** A table with columns 'OP Name' and 'Modified Date'. It lists several operations, including 'New Operation(1)', 'Cutting Agents', 'testsssss', 'Identification\_Dem...', 'Polystyrene', 'Biodiesel', 'New Operation', 'timelinelimit', 'STANDARD', and 'Identification\_Dem...'. The 'New Operation(1)' entry is highlighted.
- Right Panel:** The 'Acquisition' configuration panel for 'New Operation(1)'. It includes:
  - Integration Time:** 1000 ms
  - Laser power:** 100 %
  - Auto Dark:** Enabled (checkbox checked)
  - Auto Time:** Enabled (checkbox checked)
  - Duo:** Disabled (checkbox unchecked)
  - Average:** 1
  - External Trigger:** Disabled (checkbox unchecked)
  - Trigger Delay:** 0
  - Timeline:** Enabled (checkbox checked)
  - Interval:** 5 ms
  - Total Scans:** 10
  - Output:** A section for output configuration.

At the bottom of the interface, there is a status bar showing '(1) BWS478-785H-HT - VIR', 'Laser Key : ON', and 'Interlock : Closed'.

The following scan parameters can be set in the Acquisition section: Integration time, Laser Power, Average, Trigger, and Timeline.



## Processing



The Processing tab lets you apply data-processing algorithms to acquired spectra. These settings only apply to the Explore module visualization.

Processing Algorithms will apply to the selection of the X-axis and Y-axis defined in the operating procedure.

### To add Processing:

1. Click the + button in the OP details panel.
2. Check **Processing** in the pop-up menu.
3. Click **OK**.
4. Configure spectral processing steps by enabling the toggle for each processing step.
5. Adjust the parameters for each enabled processing step.
6. Click "Save" to save the OP.

### X-Axis and Y-Axis Conversion

#### X-Axis Options:

**Pixel:** The pixel index of the spectrometer.





## Smoothing

- Toggle On/Off
- Smoothing Types:
  - Savitzky-Golay (SG): Polynomial smoothing
    - Window Size: 5-99 (odd numbers only)
    - Polynomial Order: 1-5
  - Boxcar: Moving average filter
    - Width: 3-99 (odd numbers only)
  - FFT: Frequency domain filtering
    - Cutoff Frequency: 1-1024

## Derivative

- Enable/Disable: Toggle switch
- Derivative Type: Savitzky-Golay
  - Derivative Order: 1-2
  - Window Size: 5-99 (odd numbers)
  - Polynomial Order: 1-5

## Identification

**To use an OP in the Identification module, there must be an Identification setup included in the OP.** If the OP does not have an identification section, it will not appear in the Identification OP selector. An OP can operate in one of three identification modes: Normal, ST, or Auto.



OPs with more than one ID mode or with multiple large libraries selected for matching take longer to complete acquisitions.

The configurations available in OP – Identification are the same as those in the Identification module's configuration panel.

Here are the three OP identification modes the user can choose to configure.

## Normal Mode

SpecSuite Performance: Fail ADMIN

OP Method Variable

New Edit Save Archive Withdraw

OP Name	Modified Date
New Operation(1)	3/16/2026 11:39 AM
Cutting Agents	3/5/2026 12:56 PM
testsssss	3/3/2026 11:29 AM
Identification_Dem...	2/23/2026 10:45 AM
Polystyrene	3/13/2026 4:36 PM
Biodiesel	3/5/2026 12:58 PM
New Operation	12/17/2025 5:22 PM
timelinelimit	3/4/2026 5:55 PM
STANDARD	12/17/2025 9:39 AM
Identification_Dem...	2/23/2026 8:51 AM

Name New Operation(1)

Acquisition

Processing ×

Identification ×

Routine Analysis ×

Output

Mode: Normal

Spectral Range (Raman Shift): 176 - 2500

Library: Illicit and General Chemicals

Number Of Hits: 3 HQI Threshold: 60

⌨ (1) BWS478-785H-HT - VIR Laser Key: ON Interlock: Closed

In Normal ID mode, the user can change the spectral range over which the data is processed, the spectral library used for matching, the number of hits, and the HQI threshold.

## See-Through Mode

SpecSuite Performance: Fail ADMIN

OP Method Variable

New Edit Save Archive Withdraw

OP Name	Modified Date
New Operation(1)	3/16/2026 11:39 AM
Cutting Agents	3/5/2026 12:56 PM
testsssss	3/3/2026 11:29 AM
Identification_Dem...	2/23/2026 10:45 AM
Polystyrene	3/13/2026 4:36 PM
Biodiesel	3/5/2026 12:58 PM
New Operation	12/17/2025 5:22 PM
timelinelimit	3/4/2026 5:55 PM
STANDARD	12/17/2025 9:39 AM
Identification_Dem...	2/23/2026 8:51 AM

Name New Operation(1)

Acquisition

Processing ×

Identification ×

Routine Analysis ×

Output

Mode: See Through

Spectral Range (Raman Shift): 200 - 2000

ST Library: Library\_regroup\_Illicit

Container Library: STContainer

Number Of Hits: 1 Threshold: 60

⌨ (1) BWS478-785H-HT - VIR Laser Key: ON Interlock: Closed

In See-Through ID mode, the user can configure spectral ranges, ST Library (Matched library), Number of Hits, and Threshold. To use See-Through mode, the default Container Library must be uploaded.


## AUTO See-Through mode

The screenshot displays the SpecSuite software interface. On the left, there is a vertical toolbar with icons for New, Edit, Save, Archive, and Withdraw. Below the toolbar is a table listing operations with columns for 'OP Name' and 'Modified Date'. The main area is divided into sections for Acquisition, Processing, Identification, and Routine Analysis. The 'Identification' section is active, showing settings for 'Normal' and 'ID See Through' modes. The 'ID See Through' settings include Spectral Range (Raman Shift) from 200 to 2000, ST Library set to 'Library\_regroup\_Illicit', Container Library set to 'STContainer', Number Of Hits set to 1, and Threshold set to 60. The bottom status bar shows 'Laser Key : ON' and 'Interlock : Closed'.

OP Name	Modified Date
New Operation(1)	3/16/2026 11:39 AM
Cutting Agents	3/5/2026 12:56 PM
testsssss	3/3/2026 11:29 AM
Identification_Dem...	2/23/2026 10:45 AM
Polystyrene	3/13/2026 4:36 PM
Biodiesel	3/5/2026 12:58 PM
New Operation	12/17/2025 5:22 PM
timelinelimit	3/4/2026 5:55 PM
STANDARD	12/17/2025 9:39 AM
Identification_Dem...	2/23/2026 8:51 AM

In Auto See-Through Mode, user can configure both Normal mode settings and ST mode settings. Users can also configure Auto See-Through threshold, which will trigger the ST matching algorithm automatically by subtracting containers. Additional instructions on how to set up the Auto See-Through mode can be found in section 7.1 How to configure the Identification parameters

## Routine Analysis



The Routine Analysis tab configures quantitative analysis workflows that combine mixture analysis and prediction models.

### Configuration Areas:

#### 1. Prediction Model Selection

- Select a Method to use in Routine Analysis. The Routine analysis will run the model/variable included in the Method.

#### 2. Run Information Fields


Enable or disable these fields to track metadata during each routine analysis scan.

- Product Name
- Product Number
- Batch Number
- Lot Number

- Container Number
- Supplier Name
- User Defined Fields (UDF 1, 2, 3) (where user can define their own metadata)

The Sample Name is automatically generated by combining Product name, Run name, and Run Index. The user can input a maximum of 50 characters in ALL the fields listed above.

## Output



The Output section defines how the data is saved after the OP completes.

Choose which file format the data is exported in. Choose from SPC, simple CSV, or simple SPC.

Choose which X-axis and Y-axis data is saved for simple CSV and simple SPC file types.

Output also allows the user to automatically save data collected with the OP in any module to a specified location.





When the Model Variable is chosen, use the "Select Model" drop-down and the Select Variable drop-down to choose the variable.

Select Model

Bio-diesel

Select Model

Bio-diesel

New Model

Sunflower Oil

Select Model

Bio-diesel

Select Variable

Bio-diesel\_ Ethanol%

Add variable

Click "Add Variable". The selected variable is now added to the Selected Variable list on the right.

Name New Method

Select Model

Bio-diesel

Select Variable

Bio-diesel\_ Ethanol%

Add Variable

Selected Variable

Bio-diesel\_ Ethanol%

X

When Basic Variable is chosen, users can select the available basic variable from the drop-down list.

Add Variable

Basic Variable

Select a basic variable

- Basic variable\_576
- Basic variable\_129
- Polystyrene 1001 FWHM
- New Variable
- Polystyrene 1031 FWHM
- Test Basic
- polystyrene peak

Click **Add** to add the basic variable to the Method.

Name

---

Add Variable

Basic Variable

Select a basic variable


Polystyrene 1001 FWHM

**Add**

Selected Variable

Polystyrene 1001 FWHM

The user can add up to 6 variables to a single method. The user can also click on the remove  icon to remove a variable from the Selected Variable list.

Once all the variables are added to the Selected Variable list, click the "Save" button  to save the method. A message indicating the method has been saved will appear.

### 6.4 Variables

Variables are items to be included in a Method and can be managed at the variable page.



The types of variables are:



- Basic Variable, which can be created, edited, and reviewed in the Explore page or the Variable page.
- Model Variable, which can be created and edited in the Model module by editing the “y1” field.

## Basic Variable

Basic variables are peak analysis variables that can be derived from a specific spectrum. The user can define the processing steps, the spectral ranges used for each spectrum, and set the peak parameters for the Basic variable.

### Creating a Basic Variable

The user can create the Basic variable from the Explore Page:

- Start peak analysis by clicking 
- Adjust the spectral Range 
- Select either Datapoint or Parabolic for initial peak analysis, adjust the min peak height and max peak counts, then click Find Peaks.

Algorithm Type:

Min Peak Height:

Min Peak Width:


Max Count


- (Optional) Switch to the Curve Fit algorithm and click Find Peaks. Note: the Curve Fit algorithm will only run when initial peaks are found through the Datapoint or Parabolic algorithms.



.....

The user can also create a new basic variable from the Variable page:

1. Click on New  to open up the editing panel on the right.
2. Choose the Spectral Range, X-axis units, and Y-axis type



**Choose Spectral Range:**


Start:  End:

Choose X-Axis:

Choose Y-Axis:

3. Select an algorithm from the Choose Peak Analysis Algorithm drop-down menu and provide the peak parameters:

Choosing Datapoint:




**Choose Peak Analysis Algorithm:**

Min Peak Height:

Users will need to provide the minimum peak height

Choosing Parabolic algorithm:




**Choose Peak Analysis Algorithm:**

Min Peak Height:

Min Peak Width:


Users will need to provide the appropriate minimum peak height and minimum peak width. The minimum peak width unit will use the same X-axis units selected in the Choose Spectral Range dialog.





If Data Points or Parabolic is used for peak analysis, the highest-intensity peak will be selected as the output if multiple peaks are found. Users can use the spectral range and Min peak height to ensure only one peak is qualified.

If Curve fit is used for peak analysis, the user needs to identify the specific peak to output by selecting the labeled peak position through the drop-down menu. The available peak positions are derived from the initial peak labeled in the previous step, and the basic variable outputs the calculated peak value based on that label.




5. Click save to save the basic variable settings.



Identification can be performed online via a system connected to SpecSuite or offline using previously acquired spectral data. Data must be complete, meaning it includes both regular data and dark data. Data formats supported for Offline ID include CSV, TXTR, TXT, and SPC file formats.

## 7.1 Online ID





### How to configure the identification parameters

1. **Acquisition:** Use this module to adjust laser power and integration time. Users can also check "Auto Time" to enable auto integration time.
2. **Identification:** Use this module to choose between the types of identification, including Auto, Normal, and See-through.
  - a. **Auto:** When Auto is selected, the user can define an ST Auto threshold. The system will first perform Normal ID. If any container spectrum is identified and its HQI value exceeds the ST Auto Threshold, then See-Through ID is automatically executed.
  - b. **Normal:** Normal ID uses a Hit Quality Index (HQI) to analyze the spectra and provide the highest HQI match result from the selected libraries.





- a. The instrument is connected to SpecSuite.
  - b. Online status is displayed in the lower left corner of the window.
  - c. Performance status in the top right corner of the window is green , in valid standing.
2. Select Identification Module .
  3. Select a dropdown to choose an Operating Procedure.
    - a. The identification Details will populate below the OP selection.
    - b. The ID types will appear according to the chosen OP.
      - i. Normal mode will show Normal ID.
      - ii. ST mode will show See-Through ID.
      - iii. AUTO mode will show Normal and ST.
    - c. The library options will show the available libraries for matching.
  4. Review the ID Details.
    - a. Only Developers and Administrators can adjust ID-related OP settings from the details panel.
    - b. ID parameters may be adjusted.
    - c. Other parts of the OP must be edited in the OP module.
  5. The software is now properly configured to run the ID.
    - a. Verify the sample is correctly positioned for acquisition.
    - b. Ensure the probe shutter is open.
  6. Press the Run button to begin the ID acquisition.
    - a. While scanning, press the Stop button to stop the scan at any point.
    - b. The scan data will appear in the chart. In the Sample List table, the Result column(s) will say "In progress..." while the matching algorithm runs.



## 7. Review your results.

- ID results will appear in the result table in the middle of the screen with applicable matching statistics. Normal and ST results will show the HQI value.
- The spectrum chart will display the acquired spectrum overlaid with the spectrum selected in the result table.
- Switching to the sample info tab will display information for the chemical selected in the results table: name, synonyms, CAS, Classification, Formula, Molecular weight, NFPA Diamond, and GHS information.
- The ID record is saved in the database, and an entry is added to the bottom scan list on screen.
- Selecting any entries from the scan list will update the info shown in the Result table, Sample Info, Safety Info, and the spectrum chart.

All Online ID results and spectra are automatically saved in the connected database. The user can use the Report module to review and e-sign the identification report.

## 7.2 Offline ID

The screenshot displays the SpecSuite software interface. The top bar shows 'Performance: N/A' and 'ADMIN'. The main window is divided into several sections:


- Left Sidebar:** Contains icons for various functions including a graph, a computer monitor, a document, a fingerprint, a magnifying glass, a line graph, a document with a checkmark, a document with a magnifying glass, and a gear.
- Operating Procedure:** A dropdown menu set to 'STANDARD'.
- Acquisition:** Settings for 'Type' (Auto) and 'ST Auto Threshold' (60).
- Identification:** Settings for 'Spectral Range (Raman Shift)' (176 - 2500), 'Library' (Illicit and General Chemicals), 'Number Of Hits' (1), and 'HQI Threshold' (70).
- SeeThrough ID:** Settings for 'Spectral Range (Raman Shift)' (200 - 2000), 'ST Library' (Illicit and General Chemicals), 'Container Library' (STContainer), 'Number Of Hits' (1), and 'Threshold' (70).
- Matched Chemicals Table:**

HQI	Matched Chemical	Library
98.46	Polystyrene	Illicit and General C...
- Sample Info:** Fields for Name (Polystyrene), Synonyms (Styrene), CAS (100-42-5), Formula ((C<sub>8</sub>H<sub>8</sub>)<sub>n</sub>), Molecular Weight, Classification (Polymer), Hazard Comment (Plastic), and Background Color (Green).
- Raman Spectrum Plot:** Shows 'Processed' intensity vs 'Raman Shift (cm<sup>-1</sup>)'. Two traces are visible: 'Compliance\_Testing\_Polystyrene\_124\_15' (red) and 'Polystyrene' (blue).
- Table at the bottom:**

Scan Index	ID Status	ID Result	ID Algorithm	Mixture Result	OP	Date And Time	Re Run on Scan Ind...
1	Match	Polystyrene	Normal	N/A	STANDARD	11/23/2025 03:32:16 PM	N/A

Identification can be run without a system online or on previously collected data. **Offline ID records are NOT saved in the database.**

### How to perform an OFFLINE ID:

1. Select the Identification Module .
2. Select a dropdown to choose an Operating Procedure.
  - a. The identification Details will populate below the OP selection.
  - b. **For Offline ID**, the acquisition settings for the OP will be ignored. Only the library matching settings will be applied to the data.
3. Review the ID Details.
  - a. Only Developers and Administrators can adjust ID related OP settings from the details panel.
  - b. ID parameters may be adjusted.
4. Select the Run Offline button.



The Library module contains all the tools needed to view, create, and share user-created libraries. Metrohm-licensed libraries are also stored here, with limited permissions. The Library module opens to the Library List. The Library module is only included with the ID license.

## 8.1 Library List

Library Name	Version	Created Date	Mode	Status	Origin	Samples Count	Spectral Range
General Chemicals	14.02	11/23/2024 11:33:...	Normal	Unreleased	User	1366	176 - 2500
Illicit and General...	14.03	09/07/2023 05:08:...	Normal	Unreleased	Metrohm	2398	176 - 2500
MCRL Dyes Pigme...	14.02	09/09/2024 03:03:...	Normal	Unreleased	User	235	176 - 2500
MCRL Inorganics a...	14.02	09/09/2024 04:24:...	Normal	Unreleased	User	2191	176 - 2500
MCRL Minerals	14.00	11/23/2024 12:10:...	Normal	Unreleased	User	3850	176 - 2500
MCRL Organics-Ha...	14.02	11/23/2024 08:43:...	Normal	Unreleased	User	4113	176 - 2500
MCRL Organics-Nit...	14.02	11/23/2024 09:46:...	Normal	Unreleased	User	4057	176 - 2500
MCRL Organics-Ox...	14.02	11/23/2024 10:34:...	Normal	Unreleased	User	3331	176 - 2500
MCRL Organics-Sul...	14.02	11/23/2024 07:49:...	Normal	Unreleased	User	1915	176 - 2500
MCRL Personal Car...	14.02	09/09/2024 02:40:...	Normal	Unreleased	User	2893	176 - 2500
MCRL Polymers an...	14.02	09/09/2024 03:39:...	Normal	Unreleased	User	2526	176 - 2500
SpecSuite Library	1.00	03/18/2025 09:25:...	Normal	Unreleased	User	17	176 - 2500
STContainer	1.04	11/09/2023 05:15:...	Container	Unreleased	User	5	176 - 2500
TestB1	1.01	09/19/2024 07:32:...	Normal	Unreleased	User	2	176 - 2500

The library list is a table for viewing the available spectral libraries in the SpecSuite database.

Information about each library is available in the library table. Some examples are Library Name, Version, Created date and time, Mode, Status, Origin, and Sample Count.

In the library ribbon, there are options for the libraries in the library list below:

- **New** – Creates a library
- **Export** – Saves the selected library in a \*.zip file format to make it portable. Library packages are compatible with SpecSuite and Metrohm TacticID-1064 ST.
- **Import** – A user-created library can be imported either from the database or from the computer hard drive in the form of a Library package \*.zip file. Purchased Metrohm Libraries can also be imported using this feature.





- A maximum of 40 characters can be entered as the library name. Once saved, the library name cannot be changed.
- Library names cannot contain the following symbols: / \ : \* ? > < | .
- A space character is not allowed at the start of the library name.
- A library name cannot match one of the reserved Metrohm Library names.
  - Reserved Metrohm Library names include:
    - Illicit and General Chemicals
    - MCRL Inorganics and Organometallics
    - MCRL Organics
    - MCRL Personal Care Products
    - MCRL Polymers and Polymer Additives
    - MCRL Dyes Pigments and Stains
    - MCRL Minerals
    - MCRL Pesticides
    - MCRL Biochemicals
    - MCRL API
    - General Chemicals
    - CWA
    - EXP
  - If a new Metrohm Library is imported and there is already an existing user library with a matching name, the user library will be appended with '-C' to avoid duplicate library names.
    - Example: 'Material' would be updated to 'Material-C'

In the Mode dropdown, users can choose to create a Container library. Container libraries can only be selected as the Container Library when performing SeeThrough ID.

## 8.3 Adding Samples

The screenshot shows the SpecSuite software interface. At the top, there's a toolbar with icons for 'Back to List', 'Collect New', 'Import from Library', 'Load From File', 'Save', and 'Remove'. Below the toolbar, there are tabs for 'Library Info', 'Library Samples', and 'Library History'. The 'Library Samples' tab is active, showing a search bar and a table of samples. The table has columns for 'Sample Name', 'Synonyms', 'CAS', and 'Origin'. Below the table, there's a 'Sample Info' panel for the selected sample, 'Agarose', which includes fields for Name, Synonyms, CAS, Formula, Molecular Weight, Classification, Hazard Comment, and Background Color. To the right of the table, there's a 'Spectrum' plot showing Intensity vs. Raman Shift (cm⁻¹). The plot shows a blue line representing the spectrum, with a peak around 1000 cm⁻¹. The y-axis ranges from 0 to 60000, and the x-axis ranges from 500 to 2500.

Sample Name	Synonyms	CAS	Origin
*Agarose	Sepharose	9012-36-6	Metrohm
*Alanine	H-Ala-OH-L-Alanine	56-41-7	Metrohm
*Aspartic Acid	(S)-(+)-Aminosuccinic acid:L-Asparti...	56-84-8	Metrohm
*Benzene	Benzole	71-43-2	Metrohm
*Benzoic Acid	Dracylic Acid	65-85-0	Metrohm
*Benzyl Alcohol	Phenylcarbinolalpha-Cresol	100-51-6	Metrohm
*Beta-carotene			User
*Bleach			User
*Buprenorphine			User
*Butylparaben	4-Hydroxybenzoic acid butyl ester	94-26-8	Metrohm
*Citric acid (anhydrous)			User
*corn starch_1064			User

On the Library Samples tab, samples can be added from:

- An existing library.
- A local file.
- Collecting a spectrum with an online instrument.

Samples can be removed from the library by:

1. Select the sample(s) to be removed.
2. Press the Remove button.

The process of removing samples is recorded in the Library History.

All samples within a library must have unique names. A warning message will appear if a sample is given a non-unique name. A library cannot be saved when two samples have duplicate names.

To save a library, it must contain at least one sample.

### Import from Library

Samples from a Metrohm licensed library can be used to create a user library. The samples can be edited once regrouped into a user library.



Select All	Sample Name	Origin
<input type="checkbox"/>	n-Tricosane	Metrohm
<input type="checkbox"/>	n-Tridecane	Metrohm
<input checked="" type="checkbox"/>	Naphtha	Metrohm
<input checked="" type="checkbox"/>	Naphthalene-1,4-dicarboxylic acid	Metrohm
<input checked="" type="checkbox"/>	Neopentylbenzene	Metrohm
<input type="checkbox"/>	Ninhydrin	Metrohm
<input type="checkbox"/>	Nonadecanoic acid	Metrohm
<input type="checkbox"/>	Nonadecanophenone	Metrohm
<input type="checkbox"/>	Nonadecylbenzene	Metrohm
<input type="checkbox"/>	Nonanal	Metrohm
<input type="checkbox"/>	Nonanal Diethyl Acetal	Metrohm

Processed

Raman Shift (cm-1)

3 sample(s) selected

Add

1. Press the Import from Library button.
2. From the Library dropdown, select which library to add samples from.
  - a. Select the sample(s) to be added to the library by checking the box next to each sample.
3. Press the Add button.
  - a. Note that samples can only be selected from one library at a time.
  - b. Changing the library selection will clear the current sample selection.

## Load from File


Sample files saved to the computer can be selected to add to a library.

1. Press the Load from File button.
  - a. A file browser window appears.
2. Choose the file(s) to load and press Open.
  - a. The samples will be loaded into the library.

## Collect New

The screenshot shows the SpecSuite software interface. At the top, there is a status bar with 'Performance: Pass | Expires in 20 hour(s)' and a user name 'dale'. Below this is a toolbar with icons for 'Back To List', 'Collect New', 'Import from Library', 'Load From File', 'Save', and 'Remove'. The 'Collect New' button is highlighted with a red circle. The main window is divided into three sections: 'Library Info', 'Library Samples', and 'Library History'. The 'Library Samples' section is active, displaying a table of library entries. The 'Spectrum' section on the right shows a Raman spectrum plot with 'Raman Shift (cm<sup>-1</sup>)' on the x-axis (ranging from 200 to 1600) and 'Intensity' on the y-axis (ranging from 0 to 1400). The 'Sample Info' section on the right contains fields for Name, Synonyms, CAS, Formula, Molecular Weight, Classification, Hazard Comment, and Background Color.

Sample Name	Synonyms	CAS	Origin
*Acquired 3/24/2026 1:19:15 PM			User
*Agarose	Sepharose	9012-36-6	Metrohm
*Alanine	H-Ala-OH;L-Alanine	56-41-7	Metrohm
*Aspartic Acid	(S)-(-)-Aminosuccinic acid(L-Aspart...	56-94-8	Metrohm
*Benzene	Benzole	71-43-2	Metrohm
*Benzoic Acid	Dracylic Acid	65-85-0	Metrohm
*Benzyl Alcohol	Phenylcarbinol/alpha-Cresol	100-51-6	Metrohm
*beta-carotene			User
*bleach			User
*Buprenorphine			User
*Butylparaben	4-Hydroxybenzoic acid butyl ester	94-26-8	Metrohm
*citric acid (anhydrous)			User
*corn starch_1064			User

Press the Collect New  button to acquire new library spectrum from an Online system. The acquisition is collected using the Standard OP, which cannot be changed.

## Library Samples Info

Detailed library sample information can be edited for each user-created library sample entry, such as:

- Sample Name
- Synonyms
- CAS number
- Formula
- Molecular Weight
- Classification
- Hazard Comment
- Background Color

**Sample Info**

Name  
Acetaminophen

Synonyms  
APAP; paracetamol

CAS  
103-90-2

Formula  
C8H9NO2

Molecular Weight  
151.16

Classification  
Pharmaceutical

Hazard Comment  
Over the Counter Drug

Background Color  
Green

## 8.4 Library History

SpecSuite Performance: Pass | Expires in 20 hour(s) dale

Library Samples

Back To List Save

Library Info Library Samples Library History

Date & Time	User Name	Description
3/24/2026 1:20 PM	dale	created Library library; version 1.00
3/24/2026 1:20 PM	dale	imported Agarose sample to Library library (MCRL Biochemicals)
3/24/2026 1:20 PM	dale	imported Alanine sample to Library library (MCRL Biochemicals)
3/24/2026 1:20 PM	dale	imported Aspartic Acid sample to Library library (MCRL Biochemicals)
3/24/2026 1:20 PM	dale	imported Benzene sample to Library library (Comprehensive USP)
3/24/2026 1:20 PM	dale	imported Benzoic Acid sample to Library library (Comprehensive USP)
3/24/2026 1:20 PM	dale	imported Benzyl Alcohol sample to Library library (Comprehensive US)
3/24/2026 1:20 PM	dale	imported Butylparaben sample to Library library (Comprehensive USP)
3/24/2026 1:20 PM	dale	imported beta-carotene sample to Library library (beta-carotene.txt)
3/24/2026 1:20 PM	dale	imported Bleach sample to Library library (Bleach.txt)
3/24/2026 1:20 PM	dale	imported Buprenorphine sample to Library library (Buprenorphine.txt)
3/24/2026 1:20 PM	dale	imported citric acid (anhydrous) sample to Library library (citric acid (a
3/24/2026 1:20 PM	dale	imported corn starch_1064 sample to Library library (corn starch_1064
3/24/2026 1:20 PM	dale	added Acquired 3/24/2026 1:19:15 PM sample to Library library (BWS-

Spectrum

Sample Info

Name  
Acquired 3/24/2026 1:19:15 PM

Synonyms


CAS  
Formula  
Molecular Weight

Classification  
Hazard Comment  
Background Color

(2) BWS465II-785H - VKY Laser Key : ON Interlock : Closed

Each library contains its own history record. The Library History allows the user to see all actions performed on a library and who performed them. The same records are also included in the Audit Trail.

## 9. Model



Quantitative chemometric analysis tools are a part of SpecSuite. The Model module currently supports PLS1 quantification.

The screenshot shows the SpecSuite Model module interface. At the top, there are three buttons: 'New' (plus sign), 'Edit' (pencil), and 'Withdraw' (circle with slash). Below these is a table of models. To the right of the table is a detailed view of the selected model, 'Ethanol Conc. Model'.

Model Name	Type	Modified Date	Status	Version
basic variable	Quantification	3/4/2026 8:20 AM	Being built	1.00
Ethanol Conc. Model	Quantification	12/9/2025 9:08 AM	Enabled	1.00
Ethanol test model	Quantification	3/3/2026 9:36 AM	Being built	1.01
made with Duo	Quantification	12/8/2025 10:51 AM	Enabled	1.02
New Model	Quantification	3/3/2026 4:00 PM	Being built	1.00

Model name: Ethanol Conc. Model  
Status: Enabled  
Version: 1.00  
Created: 12/9/2025 9:08 AM  
Created by: aspen  
Model type: PLS1  
Calibration samples: 32  
Validation samples: 7  
Num LVs: 2  
Spectral Range (cm<sup>-1</sup>): 200 - 2500  
Processing: Derivative: Savitzky-Golay 2nd 15 3  
Normalization: Mean

Click New or Edit to start a model-building or editing session. There are 3 main steps to building a successful model:

1. Select Samples.
2. Parameterize the Model.
3. Build and Validate the Model.

## 9.1 Select Samples



### Adding Samples

The sample table is located in the lower left of the Build screen. Above the sample table, there are 2 icons that are used to add and remove samples to the model:

 adds new spectra to the data set.

- A file browser opens for the user to select files from their local PC.
  - Allowed formats are: CSV, TXTR, TXT, and SPC.


 removes the selected spectra from the data set.

- Selecting a row in the sample table will highlight it on the XY chart, just like in Explore.
- On the spectral chart, selecting a trace on the XY chart will select the sample in the sample table.
  - Data can be selected visually to be classified as an outlier or not to be used.

Database samples can also be used for model building. First, the sample must be moved to Explore. The sample can then be moved to the model-building stage.




## 8. Choose 'Send for Model Building.



a. Selected samples are moved to model.

## Input y-data


In the Sample table, the model variable is listed in the right-most column. The model variable name "y1" in the column header and the reference values for the model listed in this column can be set next.







3. Users can also quickly update the y-values by right-clicking on the table and selecting 'Load Y1 Reference Values'.
  - a. A file browser will open, allowing users to load a CSV file containing the reference values.
  - b. The CSV file should contain a single column of variable values in the same sequence as the samples included in the model.



Where you right-click to import Y1 Reference Values is where the values will begin to populate from. If you want to import values for all y1 variables, make sure to right-click the top row.

For example, when a user clicks in the second row to load the y values, it will put the first value from the file in the y1 value column. .CSV file in the second row, and then it will continue with rows down until there are no more reference values in the file or no more rows to set the values. Check your reference values before moving forward.

### Split the Data




Splits Dataset – This function can automatically split the dataset. In the pop-up panel, the percentage of validation spectra can be defined. There is also a toggle to detect spectral Outliers during the sample selection process.






By default, the Outliers are not included on the histogram, but can be added by right-clicking on the histogram chart and selecting 'Show Outliers'.



After the samples are added and split, the user may continue to the second step, which is to parameterize the Model.

The parameterization step can be skipped by selecting the Raman Model Developer (RMD) option, which automatically generates the model using the recommended settings. **At least 10 spectra must be added to the model before the model can be built.**

## 9.2 Parameterize the Model



Data Processing:

- Range selection
- Smoothing
- Derivative
- Baseline removal
- Normalization: scaling (0-1); mean; SNV

Users can use the toggle button to manually turn processing steps on or off. The exact functions of each processing step are described in [chapter 3.5](#) and [chapter 5.2.2](#).


## 9.3 Build and Validate the Model

After the model parameters are set, the user can click on the 'Validate quantification Model'



Validate quantification model step to generate the model using the custom parameters set in step 2.

The model will be built, and the final Model Validation page appears.



Terms:

SEC: Standard error of calibration

SECV: Standard error of cross-validation

SEP: Standard error of prediction

R<sup>2</sup>: Coefficient of determination

Criteria - Iterates model parameters to find results in a close SEC: SEP ratio and a high R<sup>2</sup> value.

This is the information displayed in the Model validation screen:

Diagrams:

- Figures of merit: SEC, SECV, SEP





- Example from the plot above:  
The  $\sim 950\text{--}1100\text{ cm}^{-1}$  region has large peaks  $\rightarrow$  **this region explains real chemical differences.**
- Big **positive** or **negative** values both mean “important.”  
The sign simply tells the direction of the correlation.

**You care most about these regions for model building.** These are where chemical signals live.

**Flat or very small values = unimportant / noise**

- Regions where the line is flat (near zero) mean:
  - The model sees **no meaningful chemical variation**
  - Likely just **noise**
  - These wavelengths can be **ignored** or de-weighted

**Each LV curve tells you what each factor is doing**


- LV1 usually captures the strongest chemical variation
- LV2 captures the next strongest pattern
- Higher LVs (LV7–LV10) start to look like random wiggles  $\rightarrow$  **mostly noise**

### **Reading the Score Plot**

A PLS **score plot** is a map showing how your samples relate to each other after the model has compressed the spectrum into a few “latent variables” (LVs).



## Reading the Influence Plot



An influence plot combines:


- **Hotelling's  $T^2$**  (x-axis)
  - Hotelling's  $T^2$  represents how far a sample is from the model center (in the score space).
  - High  $T^2$  = sample is an **extreme but still modeled** point.
- **Q-Residuals** (y-axis)
- Q-residuals represent how much of a sample's spectral information is **NOT explained by the model**.  
This is the "leftover noise or unexplained structure."  
High Q-residual = poor model fit.

The Example Influence Plot above shows:

- The points cluster in a horizontal band, with no one point extremely far up or far right.
- Q-residuals are extremely close together → model explains all samples similarly well.
- $T^2$  ranges from ~50 to ~350, but none are extreme.
- No outliers in the model

- Model is stable, and all samples are described by the model without large leftover error
- Consistent measurement/processing conditions with little variation in Q-residuals.


## 9.4 Start RMD (Raman Model Development)



Once the sample variable values are labeled, the user may press Start RMD. This will automatically split the dataset, parameterize the model and generate the model using default optimized parameters.



If custom parameters have been set in Step 2, they will be overwritten when Start RMD is pressed




## 9.5 Exporting Model Data

At all points of the model-building process, all charts and tables (except the Sample table) can be copied, saved, and exported.

- **For charts:** Right-click the chart to choose either to save data, image, or copy to clipboard.
- **For tables:** Use CTRL/SHIFT select function, then right click on the data to save or copy to the clipboard.

### Save Charts:



### Save Tables:

Reference...	Calculated...	Residual
5.00	4.9999	0.0001
5.00	4.9994	0.0006
5.00	5.0705	0.0705
10.00	9.9991	0.0009
10.00	10.0012	0.0012
10.00	9.9983	0.0017

A context menu is overlaid on the table with the following options:

- Save data
- Copy to clipboard

## 10. Routine Analysis



Routine Analysis is a quick way to run qualitative or quantitative analyses of chemicals using an existing model and method. Routine analysis requires the setup of an operating procedure before it can be used. The following are instructions on how to correctly set up and perform a Routine Analysis.

Performance : Fail    aspen

Start    Pause    Stop    Finish    Load

Run Index    Run Total  
1 / 5

Sample Name  
bagel\_breakfast\_5

MeModel  
MeModel\_y1    25.239

Trend plot

Run Index	Sample Name	Date & Time	MeModel MeModel_y1	Product Name	Product Num...	Batch Number	Lot Number	Container Na...	Supplier Name	lox	salt	pepper	Online
1	bagel_breakfast_1	10/22/2025...	25.157	bagel	everything	002	012	212	bagel grill	no thanks	yes	light	Yes
2	bagel_breakfast_2	10/22/2025...	25.112	bagel	everything	002	012	212	bagel grill	no thanks	yes	light	Yes
3	bagel_breakfast_3	10/22/2025...	25.133	bagel	everything	002	012	212	bagel grill	no thanks	yes	light	Yes
4	bagel_breakfast_4	10/22/2025...	25.21	bagel	everything	002	012	212	bagel grill	no thanks	yes	light	Yes
5	bagel_breakfast_5	10/22/2025...	25.239	bagel	everything	002	012	212	bagel grill	no thanks	yes	light	Yes

(1) BWS478-860.1064-HT - VIR    Laser Key : ON    Interlock : Closed

### 10.1 Preparation

Three things need to be set up before using the Routine Analysis module.

Create a **Variable** → Make a **Method** → Build an **OP** → *Routine Analysis*

1. Create a **Variable**



Variables are either created automatically when a model is created or can be created manually as a Basic Variable. Follow the instructions in the Model section for detailed model creation, or the instructions in the Basic Variable section.

### 2. Build a Method

Methods are set up within the OP Module. Refer to the OP Methods section.

### 3. Build an OP

To use Routine Analysis, an OP must be configured with a **Routine Analysis module enabled**. Only one method can be selected in an OP, but each Method can contain multiple variables.

## Offline Routine Analysis

Offline Routine Analysis provides a testing environment for Developer and Administrator-level users to test the Method created. CSV, simple CSV or simple SPC data can be loaded to make predictions. Multiple data files can be loaded at once.

Once loaded, the system will automatically process the data using the selected algorithms and generate the prediction result on screen. Offline routine analysis data is not saved to the database.


The screenshot displays the SpecSuite software interface. On the left is a vertical toolbar with various icons. The main window is divided into several sections:

- Operating Procedure:** Shows 'Offline' as the current mode and 'Biodiesel Ethanol' as the selected method.
- Batch Info:** Includes a 'Batch Quantity' field set to '1'.
- Control Panel:** Features buttons for Start, Pause, Stop, Finish, and Load.
- Run Summary:** Displays 'Run Index 1 / 1' and 'Run Total 1 / 1'.
- Sample Information:** Shows 'Sample Name: Compliance\_Testing\_Polystyrene\_124\_9', 'Bio-diesel: Ethanol%: 24.798', and 'Sunflower Oil: Olei (mg): 0.000E+000'.
- Trend Plot:** A line graph with 'Run index' on the x-axis (ranging from 2 to 14) and 'Variable' on the y-axis (ranging from 0 to 25). It tracks two variables: Ethanol% (blue line, constant at 24.798) and Olei (mg) (green line, constant at 0.000E+000 until run index 8, then drops to 0).
- Table:** A data table with columns: Run Index, Sample Name, Date & Time, Bio-diesel Ethanol%, Sunflower Oil Olei (mg), and Online. It lists 10 runs.


Run Index	Sample Name	Date & Time	Bio-diesel Ethanol%	Sunflower Oil Olei (mg)	Online
1	Compliance_Testing_Polystyrene_124_15	11/23/2025 05:32:24 PM	24.707	10.819	No
2	Compliance_Testing_Polystyrene_124_14	11/23/2025 05:32:24 PM	24.746	10.9	No
3	Compliance_Testing_Polystyrene_124_13	11/23/2025 05:32:24 PM	24.7	10.881	No
4	Compliance_Testing_Polystyrene_124_12	11/23/2025 05:32:24 PM	24.771	10.949	No
5	Compliance_Testing_Polystyrene_124_11	11/23/2025 05:32:24 PM	24.694	10.874	No
6	Compliance_Testing_Polystyrene_124_10	11/23/2025 05:32:24 PM	24.739	11.001	No
7	Compliance_Testing_Polystyrene_124_9	11/23/2025 05:32:24 PM	24.798	10.854	No
8	Compliance_Testing_Polystyrene_124_8	11/23/2025 05:32:24 PM	24.722	10.897	No
9	Compliance_Testing_Polystyrene_124_15	11/23/2025 05:33:09 PM	24.707	0.000E+000	No
10	Compliance_Testing_Polystyrene_124_14	11/23/2025 05:33:00 PM	24.746	0.000E+000	No



If the timeline is enabled in the operating procedure, the next scan will automatically start at the interval defined in the OP. There will also be a timer counting down to the next scan.




If the external trigger is enabled in the OP, the user must press the external trigger on the probe body to start the measurement after clicking Start.



- **Pause** – The pause button will appear during a multi-acquisition run. Pausing will allow the current scan to complete, then all scans will be paused until the Scan button is pressed again. The next scan will immediately begin regardless of any set delay (Interval) between scans.
- **Stop** – The stop button will immediately halt the acquisition being collected. The cancelled scan will not appear on the chart or in the data list. If Stop is pressed during a multi-acquisition, the entire acquisition will be canceled.






The sample list of a specific Routine Analysis run can be found in the table at the bottom of the screen. The routine analysis record cannot be changed once it has been generated. Each specific run will refresh the spectra plot, trend plot, and the sample list.


Run Index	Sample Name	Date & Time	Bio-diesel Ethanol%	Sunflower Oil Olei (mg)	Product Name	Product Number	Batch number	Lot Number	Container Name	Supplier Name	Online
1	Compliance_Testing...	11/23/2025 05:49:0...	34.524	0.000E+000							Yes
2	Compliance_Testing...	11/23/2025 05:49:0...	34.544	0.000E+000							Yes
3	Compliance_Testing...	11/23/2025 05:49:1...	34.545	0.000E+000							Yes
4	Compliance_Testing...	11/23/2025 05:49:1...	34.49	0.000E+000							Yes
5	Compliance_Testing...	11/23/2025 05:49:2...	34.54	0.000E+000							Yes
6	Compliance_Testing...	11/23/2025 05:49:3...	34.649	0.000E+000							Yes
7	Compliance_Testing...	11/23/2025 05:49:3...	34.53	0.000E+000							Yes
8	Compliance_Testing...	11/23/2025 05:49:4...	34.667	0.000E+000							Yes
9	Compliance_Testing...	11/23/2025 05:49:4...	34.499	0.000E+000							Yes



lookup table for users to search all the scanned samples. Users can select a specific entry and view the detailed instrument acquisition info:



The left panel shows the spectrum-acquisition parameters for instruments, including laser power, integration time, and laser wavelength. Users can also use the drop-down field to select different X-axis and Y-axis to view the spectrum on the right.




The Acquisition table is for filtering and displaying only. Users cannot generate or e-sign a PDF report in the Acquisition tab.

## 11.2 Identification Reports


Scan Index	Instrument SN	Instrument Model	Date & Time	OP	Library	ID Result	Tester	Esign Status
103	195500102601002	BWS465II-785H	3/24/2026 1:54 PM	102 OP	Comprehensive US... Illicit and General...	Acetaminophen, 9... Acetaminophen, 9... Phenacetin, 65.22 Phenacetin, 62.96 N-(4-Hydroxyphen... Methyltestosteron...	dale	Not Signed
102	195500102601002	BWS465II-785H	3/24/2026 1:53 PM	102 OP	Comprehensive US... Illicit and General...	Acetaminophen, 9... Acetaminophen, 9... Phenacetin, 65.30 Phenacetin, 63.03	dale	Not Signed
101	195500102601002	BWS465II-785H	3/24/2026 1:53 PM	102 OP	Comprehensive US... Illicit and General...	Acetaminophen, 9... Acetaminophen, 9... Phenacetin, 65.68	dale	Not Signed
100	195500102601002	BWS465II-785H	3/24/2026 1:51 PM	102 OP	MCRL Organics-Ha... MCRL Pesticides, 1... Minerals, 14.12 Comprehensive US... Illicit and General... Library, 1.00 MCRL API, 14.12 MCRL Biochemical... MCRL Organics-Sul... MCRL Personal Car... MCRL Polymers an... normal operation li...	saved as txt: 97.50 Acetaminophen, 9... Acquired 1/14/202...	dale	Not Signed

After an identification run, the record is saved in the database and will appear in the Identification Tab of the Report Module. ID reports can be viewed, signed, printed, and exported.


ID records can be filtered by Scan Index, Instrument SN, Instrument Model, Tester, OP, Library, ID result, and Date and time. Additional Columns can be added, such as Laser power, Integration time, Scan average, laser wavelength, probe SN, ID mode, container library, container spectrum, Tester's E-sign status, Approver's E-sign status, and E-sign Date and time.



Users can select 1 or multiple ID report entries and click View Report to view the details.



Click on "PDF Report Preview" to preview a PDF report, which can be saved or sent to a physical printer through the PDF report widgets.



SpecSuite Performance: Pass | Expires in 19 hour(s)

## Identification Report

### Sample Metadata

Scan Index	103
Identification Algorithm	Identification - N
Library	Comprehensive USP, 14.12 Illicit and General Chemicals, 14.09
Number Of Hits	9
Threshold	40
Spectral Range	(176,2500)
Created	3/24/2026 1:54 PM
GUID	9860d30f-54c9-41cc-71d1-08de89ce5881
Laser Wavelength (nm)	785.00 nm

### Instrument Information

Instrument Model	BWS465II-785H
Instrument SN	195500102601002
Probe SN	123123123
SpecSuite Version	1.1.5
Last PV Date/Time	3/24/2026 9:20:09 AM -04:00
Last PV Status	Pass

(2) BWS465II-785H - VKY Laser Key: ON Interlock: Closed


## 11.3 Routine Analysis Reports

After a Routine Analysis is run in the Routine Tab, the record is saved to the database and appears on the Routine Analysis page of the Report tab. Routine Analysis reports can be viewed, signed, printed, and exported.

Routine Analysis records can be filtered by

- User Name
- Date & Time
- OP (Operating Procedure)
- Method Name
- Run Name
- Sample Name
- Laser Power
- Integration Time
- Number of Averages
- Laser Wavelength
- Probe Serial Number
- Last Validation ID

- Tester E-sign status
- Approver E-sign status
- E-Sign Date and Time
- Method Variable Names and Values



Users can select a routine analysis report using the checkbox, then click View report, and a PDF Routine Analysis report will be displayed. Users can then use the PDF widget to save the report as a PDF file or send it to a printer.

Navigation sidebar with icons for: Back, Home, Reports, Attachments, Comments, and Settings.


# Routine Analysis Report

Model Name	Variable Name	Result
Ethanol in Beer	Ethanol %	-1.44
Basic Variable	Fret Basic	0.00
Basic Variable	polystyrene peak	22008628.46
Basic Variable	Basic variable_576	0.00
Basic Variable	Basic variable_129	23408485.17


## Sample Metadata

Scan Index	66
Created	3/6/2025 4:07 PM
User Name	ADMIN
GUID	41a1db61-60c2-4225-e72d-08de7bc45e7c
Sample Name	111112_1
Run Name	111112
Product Name	N/A
Batch Number	N/A
Lot Number	N/A
Container Name	N/A
Supplier Name	N/A
Product Sample Name	N/A
UDF2	N/A
UDF3	N/A


## 11.4 Performance Validation Reports




Performance Validation Records can be viewed at all user levels. Users can select the specific performance validation records from the table, then click View Report to view the detailed pass/fail results of each item.



Users can then click the Export button to export selected performance validation spectra in the CSV, simple CSV and simple SPC formats.




Click on "PDF Report Preview" to preview PDF reports.



Click "Back" to return to the previous screen.

## 11.5 Calibration records


All previously accepted calibration records will be saved in the "Calibration" tab. Users can select and view the detailed calibration data during Raman Shift Calibration, Laser Wavelength Calibration, and Relative Intensity Calibration.



The screenshot shows the SpecSuite software interface with the "Calibration" tab selected. The table below lists the calibration records:


Scan Index	Instrument SN	Instrument Model	Tester	Date & Time	Probe SN	Calibration Type	Result
83	195500102601002	BWS465II-785H	aspen	3/19/2026 11:41 AM	123123123	Laser Wavelength...	Accepted
80	ZVDR		aspen	3/19/2026 11:11 AM	123123123	Laser Wavelength...	Accepted
73	ZVDR		aspen	3/19/2026 9:38 AM	N/A	Raman Shift Calibr...	Accepted
52	240830309	BWS478-532H-HT	aspen	3/3/2026 8:28 AM	000000001	Raman Shift Calibr...	Accepted
11	19540030.25.05001	BWS478-785H-HT	aspen	12/4/2025 12:33 PM	1234567890	Raman Shift Calibr...	Accepted
3	195400602509002	BWS478-860_1064...	aspen	12/1/2025 2:51 PM	N/A	Laser Wavelength...	Accepted
2	195400602509002	BWS478-860_1064...	aspen	12/1/2025 2:49 PM	N/A	Laser Wavelength...	Accepted
1	195400602509002	BWS478-860_1064...	aspen	12/1/2025 2:48 PM	N/A	Laser Wavelength...	Accepted

Calibration data can be exported in CSV, Simple CSV, or Simple SPC files.





Note: For Calibrations performed on an i-Raman DUO system, Load to Explore will appear with a gap as shown below. This is not an error. DUO system calibrations use a smaller useful range that excludes data in this region.



## 11.6 E-Signatures

User accounts with Developer, Operator, and Lab Managers privileges are allowed to sign Identification, Routine Analysis, and Performance Validation records.

The person who executed the measurement will be shown under the “Tester” column, and the person who approves/rejects the record will be shown under the “Approver” column. Users can right-click the column to view column choosers and enable/disable additional columns.

The screenshot displays a table with columns: Library, ID Result, Tester, Esign Status, and Approver. A context menu is open over the 'Tester' column, and a 'Column Chooser' dialog box is overlaid on the right. The dialog box contains a search field and a list of columns with checkboxes. The 'Tester' column is highlighted in blue in the table.


Library	ID Result	Tester	Esign Status	Approver
Illicit and General C... MCRL Minerals, 14... MCRL Organics-Ha... MCRL Personal Car...	Polystyrene, 98...		Not Signed	
Illicit and General C... MCRL Minerals, 14... MCRL Organics-Ha... MCRL Personal Car...	Polystyrene, 98...		Approved	M...
Illicit and General C... MCRL Minerals, 14... MCRL Organics-Ha... MCRL Personal Car...	Polystyrene, 98.83	Developer	Rejected	O...
Illicit and General C... MCRL Minerals, 14... MCRL Organics-Ha... MCRL Personal Car...	Polystyrene, 98.71	Developer	Not Signed	
Illicit and General C... MCRL Minerals, 14... MCRL Organics-Ha... MCRL Personal Car...	Polystyrene, 98.83	Developer	Tester Signed	
Illicit and General C... MCRL Minerals, 14... MCRL Organics-Ha... MCRL Personal Car...	Polystyrene, 98.81	Developer	Tester Signed	

Column Chooser

Search Columns...

- Instrument SN
- Instrument Model
- Date & Time
- OP
- Library
- ID Result
- Tester
- Esign Status
- Laser power
- Integration Time
- Average
- Laser Wavelength
- Probe SN
- Last Validation ID
- Validation Status
- ID Mode
- Container Library
- Approver
- Approving Date

Laser Key : ON Interlock : Closed

To sign a record, select the record from the list and click on . The tester can sign the record through the pop-up window. The Tester must enter their username and password, then click Sign. The E-sign status will change to “Tester Signed.”

REPORTS

View Report E-Sign Load to Explore Refresh List

Acquisition Identification Performance Calibration Routine Analysis

Search

Scan Index	Instrument SN	Instrument Model	Date	Tester	Esign Status
<input checked="" type="checkbox"/>	25	195400602509002	BWS478-860_1064...	12/2/2025 9:51 AM	98.79 Operator Not Signed
<input type="checkbox"/>	23	195400602509002	BWS478-860_1064...	12/2/2025 9:50 AM	98.43 Operator Approved
<input type="checkbox"/>	19	195400602509002	BWS478-860_1064...	12/2/2025 9:50 AM	98.83 Developer Tester Signed
<input type="checkbox"/>	18	195400602509002	BWS478-860_1064...	12/2/2025 9:50 AM	98.71 Developer Not Signed
<input type="checkbox"/>	16	195400602509002	BWS478-860_1064...	12/2/2025 9:50 AM	98.83 Developer Tester Signed
<input type="checkbox"/>	15	195400602509002	BWS478-860_1064...	12/2/2025 9:50 AM	98.81 Developer Tester Signed

Page 1 of 1

178-860 1064-HT - LYD Laser Key - ON InterLink - Connect

Tested By

Report Type:  
 Scan Index: Scan Indexes: 25

User Name\*

Password\*

Additional Notes

Cancel Sign

Once the Tester signs their own record, another person (typically a lab manager) can Approve or Reject it. During the approval process, the approver can choose either "Approve" or "Reject" and select a reason based on the drop-down list.

Approval

Report Type:  
 Scan Index: Scan Indexes: 19

User Name\*

Password\*

Reason

2-1: Result Verified  
2-2: Approval for Release  
2-3: Approved with deviation noted (add additional notes)  
2-4: No discrepancies found  
2-6: Other (add additional Notes)

Cancel Sign

Approval

Report Type:  
 Scan Index: Scan Indexes: 19

User Name\*

Password\*

Reason

2-4: Result out of spec  
2-5: Data Incomplete  
2-6: Other (add additional Notes)

Cancel Sign

E-sign rules:

1. The same user (Tester) who executed the identification, performance validation or routine analysis scan must sign the record first. Non-testers cannot sign the record first.




2. Users must sign in with the same credentials as their login credentials.
3. Testers cannot approve or reject the record.
4. The same user cannot sign the record twice.
5. Operator-level users cannot approve or reject the report.

Administrators and IT users cannot sign any records.

After the user has signed the report, the E-sign status, Full name, Username and Signing Date and time will be shown in the report PDF.

## Identification Report

**ID Normal**



Intensity (a.u.)

Raman Shift (cm-1)

HQI	Chemical Name	Classification	Library Name
98.43	Polystyrene	Polymer: Plastic	Illicit and General Chemicals; V14.09

Tester Signature	Tester Signed	Approval Status	Approved
Tester Full Name	Max2	Approver Full Name	Max Lab Manager
Tester User Name	Operator	Approver User Name	Max Lab Manager
Date And Time	12/2/2025 12:06 PM	Date And Time	12/2/2025 12:09 PM

Approval Reason Result Verified


Additional Notes:

Sunday, December 7, 2025 3/3

## 11.7 Bulk Export

In all report tabs, users can select 1 or multiple records and then export the spectra file in a different format.

Selecting the record by checking the checkboxes on the left




The screenshot shows the 'Reports' section of the software interface. At the top, there are icons for 'View Report', 'Load to Explore', 'Export', and 'Refresh List'. Below these are tabs for 'Acquisition', 'Identification', 'Performance', 'Calibration', and 'Routine Analysis'. The 'Identification' tab is active, displaying a table of scan results. The table has columns for 'Scan Index', 'Instrument SN', 'Instrument Model', 'Date & Time', 'OP', 'Library', 'ID Result', 'Tester', and 'Esign Status'. Several rows are selected, indicated by checked checkboxes in the first column. The table data is as follows:

Scan Index	Instrument SN	Instrument Model	Date & Time	OP	Library	ID Result	Tester	Esign Status	
<input checked="" type="checkbox"/>	30	195500302602005	BWS465H-785S	3/5/2026 12:27 PM	Cutting Agents	Comprehensive US-Illicit and General C- Library_regroup_III	Polystyrene, 97.80 Polystyrene, 97.80 Norpropoxyphene...	Operator	Not Signed
<input checked="" type="checkbox"/>	29	195500302602005	BWS465H-785S	3/5/2026 11:51 AM	Cutting Agents	Comprehensive US-Illicit and General C- Library_regroup_III	Polystyrene, 97.79 Polystyrene, 97.79 Norpropoxyphene...	SERVICE	Not Signed
<input checked="" type="checkbox"/>	28	195500302602005	BWS465H-785S	3/5/2026 11:51 AM	STANDARD	Illicit and General C-	Polystyrene, 97.79	SERVICE	Not Signed
<input checked="" type="checkbox"/>	27	195500302602005	BWS465H-785S	3/5/2026 11:49 AM	Cutting Agents	Illicit and General C-	Polystyrene, 97.81 Norpropoxyphene... 1-Phenyl-2-phenet...	ADMIN	Not Signed
<input checked="" type="checkbox"/>	26	195500302602005	BWS465H-785S	3/5/2026 11:47 AM	Cutting Agents	Illicit and General C-	Polystyrene, 97.83 Norpropoxyphene... 1-Phenyl-2-phenet...	ADMIN	Not Signed
<input type="checkbox"/>	25	195500302602005	BWS465H-785S	3/5/2026 11:38 AM	STANDARD	Illicit and General C-	Polystyrene, 97.77	ADMIN	Not Signed
<input type="checkbox"/>	21	19540030.25.05001	BWS478-785H-HT	2/23/2026 10:49 AM	Identification_Dem...	Illicit and General C-	Polystyrene, 97.83	ADMIN	Not Signed
<input type="checkbox"/>	20	19540030.25.05001	BWS478-785H-HT	2/23/2026 10:49 AM	Identification_Dem...	Illicit and General C-	Polystyrene, 97.82 Norpropoxyphene... 1-Phenyl-2-phenet...	ADMIN	Not Signed
<input type="checkbox"/>	19	19540030.25.05001	BWS478-785H-HT	2/23/2026 10:48 AM	Identification_Dem...	Illicit and General C-	Polystyrene, 97.87	ADMIN	Not Signed
<input type="checkbox"/>	18	19540030.25.05001	BWS478-785H-HT	2/23/2026 10:44 AM	Identification_Dem...	Illicit and General C- Comprehensive US-	Polystyrene, 97.82	ADMIN	Not Signed
<input type="checkbox"/>	17	19540030.25.05001	BWS478-785H-HT	2/23/2026 10:43 AM	Identification_Dem...	Illicit and General C-	Polystyrene, 97.81	ADMIN	Not Signed
<input type="checkbox"/>	16	19540030.25.05001	BWS478-785H-HT	2/23/2026 10:43 AM	Identification_Dem...	Illicit and General C-	Polystyrene, 97.79 Norpropoxyphene... 1-Phenyl-2-phenet...	ADMIN	Not Signed
<input type="checkbox"/>	15	19540030.25.05001	BWS478-785H-HT	2/23/2026 10:43 AM	Identification_Dem...	Illicit and General C-	Polystyrene, 97.81	ADMIN	Not Signed
<input type="checkbox"/>	14	19540030.25.05001	BWS478-785H-HT	2/23/2026 10:42 AM	Identification_Dem...	Illicit and General C-	Polystyrene, 97.82 Norpropoxyphene... 1-Phenyl-2-phenet... Diphenoxylate, B3.07 Difenpropriololol, 8- 3-Fluorocetanyl H... (S) Desoxy-DPM... Cyclisanyl fentanyl... Fentanyl methylen... 1,4-DBenzyloperid... Bibacortemone, 97.9K	ADMIN	Not Signed



Then click Export

A new dialog will appear, and the user can select the file type of the spectra to export:



In the Identification, Performance, and Routine analysis report tabs, users can also export the PDF report directly by clicking Export. There is no PDF report for export in the Acquisition or Calibration Tab.

## 12. Audit Trail



Date & Time	User Name	Full Name	Category	Audit Code	Description
11/25/2025 8:50 AM	ADMIN		Login and Account	D-2	Failed to log in (Incorrect password).
11/25/2025 8:50 AM	ADMIN		Login and Account	D-2	Failed to log in (Incorrect password).
11/25/2025 8:50 AM	ADMIN		Login and Account	D-1	Logged in.
11/25/2025 8:50 AM	ADMIN		Login and Account	D-4	Changed password.
11/25/2025 8:51 AM	ADMIN		Login and Account	D-3	Logged out.
11/25/2025 8:56 AM	ADMIN		Login and Account	D-1	Logged in.
11/25/2025 8:56 AM	ADMIN		Instrument	10-0	Connected instrument BWS465-532HH ; 19550010RD11001.
11/25/2025 8:58 AM	ADMIN		Instrument	10-4	Probe (BAC102-532E) added to instrument BWS465-532HH ; 19550010RD11001.
11/25/2025 8:58 AM	ADMIN		Instrument	10-10	Probe (BAC102-532E) enabled on instrument BWS465-532HH ; 19550010RD11001.
11/25/2025 9:22 AM	ADMIN		Instrument	10-0	Connected instrument BWS465-532HH ; 19550010RD11001.
11/25/2025 9:24 AM	ADMIN		Instrument	10-10	Probe (BAC102-532E) enabled on instrument BWS465-532H ; 19550010RD11001.
11/25/2025 9:26 AM	ADMIN		Performance	2-0	Created PVR: 19550010RD11001-250303989-20251125_092622 PVR file for BWS465-532H ; 19550010RD11001.
11/25/2025 9:27 AM	ADMIN		Performance	2-7	Failed performance validation for BWS465-532H ; 19550010RD11001.
11/25/2025 11:10 AM	ADMIN		Instrument	10-0	Connected instrument BWS465-532H ; 19550010RD11001.
11/25/2025 11:11 AM	ADMIN		Instrument	10-0	Connected instrument BWS478-860, 1064-HT ; 195400602509002.
11/25/2025 11:11 AM	ADMIN		Instrument	10-14	Updated Cross Point value to 3400 for BWS478-860, 1064-HT ; 195400602509002.
11/25/2025 11:20 AM	ADMIN		Login and Account	D-3	Logged out.
11/25/2025 11:25 AM	ADMIN		Login and Account	D-1	Logged in.
11/25/2025 11:26 AM	ADMIN		Instrument	10-0	Connected instrument BWS478-860, 1064-HT ; 195400602509002.
11/25/2025 11:26 AM	ADMIN		Instrument	10-14	Updated Cross Point value to 3400 for BWS478-860, 1064-HT ; 195400602509002.
11/25/2025 11:31 AM	ADMIN		Library	5-2	Imported Ilicit and General Chemicals library; version 14.09.
11/25/2025 11:31 AM	ADMIN		Library	5-2	Imported MCRL Organics-Halogenated library; version 14.02.
11/25/2025 11:33 AM	ADMIN		Library	5-2	Imported MCRL Personal Care Products library; version 14.02.
11/25/2025 11:33 AM	ADMIN		Library	5-2	Imported MCRL Minerals library; version 14.00.
11/25/2025 11:34 AM	ADMIN		Operating Procedure	1-1	Created New Operation.
11/25/2025 11:36 AM	ADMIN		Measurement	4-0	Took an acquisition on BWS478-860, 1064-HT-195400602509002 (14642 ms; Laser power:100 %; Scan average:1).
11/25/2025 11:41 AM	ADMIN		Operating Procedure	1-8	Updated module Acquisition in New Operation.
11/25/2025 11:41 AM	ADMIN		Operating Procedure	1-6	Added module ID Normal in New Operation.
11/25/2025 11:41 AM	ADMIN		Operating Procedure	1-100	Updated Output settings in New Operation.
11/25/2025 11:41 AM	ADMIN		Operating Procedure	1-30	Added e4ad3087-55e7-4e3a-8b6b-5ca06b4152e7 in New Operation.
11/25/2025 11:41 AM	ADMIN		Operating Procedure	1-8	Updated module Acquisition in New Operation.
11/25/2025 11:41 AM	ADMIN		Operating Procedure	1-33	Updated CI Normal settings in New Operation.
11/25/2025 11:41 AM	ADMIN		Operating Procedure	1-100	Updated Output settings in New Operation.
11/25/2025 11:43 AM	ADMIN		Measurement	4-0	Took an acquisition on BWS478-860, 1064-HT-195400602509002 (14462 ms; Laser power:100 %; Scan average:1).
11/25/2025 11:43 AM	ADMIN		Measurement	4-6	Took Normal identification: Match 99.14 Polystyrene (New Operation).
11/25/2025 11:43 AM	ADMIN		Operating Procedure	1-8	Updated module Acquisition in New Operation.
11/25/2025 11:43 AM	ADMIN		Operating Procedure	1-33	Updated CI Normal settings in New Operation.
11/25/2025 11:43 AM	ADMIN		Operating Procedure	1-100	Updated Output settings in New Operation.
11/25/2025 11:43 AM	ADMIN		Measurement	4-2	Stopped an acquisition on BWS478-860, 1064-HT-195400602509002.

The Audit Trail contains records of all users and system actions. Audit Trail records include date and time, Username, category, description, notes, and device serial number.

➤ **View Report:** Loads offline file(s) to perform a routine on.

Category is a number associated with a function or module:

0: Login and Account


1: Setup

2: Operating Procedure


3: Library




**Date and time filter:** Users can click the funnel icon in the table headers to filter data by specific time ranges. Users can choose a specific date range by first selecting the start date, then drag the date range across the entire month. If dates across multiple months need to be selected, users can navigate to the next month, then drag the date range while holding down “SHIFT” key on the keyboard.




Users can also click on the Month name and the Year name to select dates for a specific month or year.






User Name filter, Full Name filter, Category filter, Audit Code filter can all be configured through the  icon in the table headers. Click the table icon, then a list of available filters can be selected:



Once the filter is selected, the summary of all filter conditions will show up at the bottom of the screen:











Users can change the display language using the drop-down menu. Once the language is changed, the software will remember to use the selected language for display during startup.



## 13.2 License



Some SpecSuite features and settings are controlled with a license that can be purchased from Metrohm. SpecSuite Licenses are bound to an instrument serial number. Users must connect the software to the instrument with the same serial number to activate and load the license.

To load your license, press the Load button and select the license file sent by Metrohm. If the license is valid for the requested computer, the License's settings will appear on this page. A prompt will appear asking the user to confirm that this license is correct, and if so, SpecSuite will be automatically restarted. Pressing cancel will clear the license information, and nothing will be changed.

Once uploaded, the License file can be found in this location on PC:

C:\ProgramData\Metrohm\SpecSuite

In case an unwanted license is loaded to SpecSuite, go to this location and delete the license, then put the correct license at this location named 'license.lic'



There are 4 different license configurations that can be purchased by customers:

- Basic
- ID
- Quant
- 21CFR.


Default configurations are shown below:

<b>Modules</b>	<b>Basic License</b>	<b>Identification License</b>	<b>Quantification License</b>	<b>Regulatory Compliance License</b>
Explore	on	-	-	-
Identification	off	on	-	-
Routine Analysis	off	-	on	-
Compliance (including Audit Trail)	off	-	-	on
Performance	on	-	-	-
Instrument	on	-	-	-
Operating Procedure	on	-	-	-
Library	off	on	-	-
Settings	on	-	-	-
Report	on	-	-	-

A trial license is included with SpecSuite and enables all features. The trial license expires after 30 days. The license will then revert to a Basic license.

All licenses include the current Illicit and General Chemicals library. Additional library packages and licenses must be purchased from Metrohm.

## 13.3 Database




The database tab is reserved for IT users only. Users can:

1. Switch between local and remote databases.
2. Connect to the Server with Server Name or IP address.
3. Create a new database in the connected server instance.
4. Select the database to connect.


Detailed instructions for setting up server/database connections are described in [chapter 1.3](#) and [chapter 1.4](#).

## 13.4 Backup and Restore



### Database Size

SpecSuite will display the Current Size, Used Space and Free Space of the connected database. The database total size is dynamic and will grow automatically as data is added over time. For SQL Server Express, 10GB is the maximum size your database is allowed to grow.











The screenshot shows the SpecSuite application interface. On the left is a vertical navigation menu with icons for various functions. The main content area is divided into two panels. The left panel, titled 'User List', contains a table with the following data:

User Name	Full Name	Account Level	Password Expirati...	Status
ADMIN		Default Administr...	3/5/2026	Enabled
Dev 222	Dev222	Developer	3/1/2026	Enabled
Dev111	New DEV	Developer	3/2/2026	Enabled
IT		IT	3/5/2026	Enabled
Lab222	Lab 222	Lab Manager	3/2/2026	Enabled
LabMan	Lab Man	Lab Manager	3/4/2026	Enabled
LabMe	Lab Me	Lab Manager	3/4/2026	Enabled
OP222	OP 222	Operator	3/2/2026	Enabled

The right panel, titled 'Edit User Account', contains a form with the following fields and controls:

- Full Name:** An empty text input field.
- User Name:** A text input field containing the value 'ADMIN'.
- Account Level:** A dropdown menu with 'Default Administrator' selected.
- Password:** A text input field with masked characters (dots).
- Change Password:** A button.
- Reset Password:** A button.
- Password Expiration:** A text input field containing the date '3/5/2026'.
- Status:** A radio button labeled 'Enabled' which is selected.
- Save:** A button.

User Name, Full Name, and Level are displayed in the Current Users list. When a User is selected from the table, the row is highlighted, and the user's information is populated in the fields on the right side of the screen.

In the Edit User Account panel, enter Full Name, User Name, Account Level, Password Expiration, and status to create a new account. Once a user account is created, only Level, Password, and Status may be changed. Password expiration cannot currently be changed; by default, it is set to 90 days from account creation. After any change, the account must be saved for the changes to take place. The changed password must be different than the last-used password.

Usernames and passwords may only use ASCII characters and the following special characters: !@#\$%^&\*?. The username length is limited to 4 - 20 characters. Username cannot contain space. Username cannot start with a special character. Usernames are unique.

## Default Accounts

Default Account (User level): For these accounts ONLY default passwords are 999999.

- ADMIN (Administrator)
- IT (IT)



## **New Accounts**

Only an Administrator or IT user can create a new user. Only IT users can create other IT users. The 'New user' button is not available to Operator or Developer-level users.

When the New User button is pressed, the Edit User Account panel will become blank and editable.

User Name: A username must be unique. User Name will be used for e-signature signing.

Full Name: Full name must be unique. It will display, in addition to the User Name on Reports, and at the top of the window to indicate the signed-in user.





Module	Sub – module	IT	Administrator	Developer	Operator	Lab Manager
Instrument	Instrument Module – Menu Tab	Not Visible or Editable				Not Visible or Editable
	Device Info Tab	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
	Device Info - Show Details	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
	Device Info - Read Flash	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
	Hardware Setup Tab	Not Visible or Editable	Full Access	Full Access	Visible Only	Not Visible or Editable
	Hardware Setup - Detector Settings	Not Visible or Editable	Full Access	Full Access	Visible Only	Not Visible or Editable
	Hardware Setup - Accessory	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
	Laser Setup Tab	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
	Laser Setup - Laser Settings	Not Visible or Editable	Full Access	Full Access	Visible Only	Not Visible or Editable
	Laser Setup - Laser info	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
	Laser Setup - Add Laser/Remove Laser	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
	Performance Tab	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
	Performance - Select Reference File	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
	Performance - Export	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
	Performance - Import	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
	Performance - Create PVR	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
	Performance - Select Probe	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
	Performance - Run	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
	Performance - Stop	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
	Calibration Tab	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
Calibration - Raman Shift Calibration	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable	






Module	Sub – module	IT	Administrator	Developer	Operator	Lab Manager
<b>Library</b>	Library Module	Not Visible or Editable			Not Visible or Editable	
	Library - New	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
	Library - Edit	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
	Library - Validate	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
	Library - Invalidate	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
	Library - Archive	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
	Library - Export	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
	Library - Import	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
	Library - Withdraw	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
	Library - Sample (Sample Info section)	Not Visible or Editable	Full Access	Full Access	Visible Only On	Not Visible or Editable
	Library - Sample (Edit, Import, Load, Save, etc.)	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
<b>Routine Analysis</b>	Routine Analysis Module	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
	Select Method	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
	Select OP	Not Visible or Editable	Full Access	Full Access	Full Access	Not Visible or Editable
<b>Module</b>	Model Module – Menu Tab	Not Visible or Editable	Full Access	Full Access	Visible Only	Not Visible or Editable
	Model - New	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
	Model - Edit/Save/Withdraw	Not Visible or Editable	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable
<b>Reports and Audit Trail</b>	Reports and Audit Trail Module – Menu Tab	Visible	Visible	Visible	Visible	Visible
	Report - View Report	Full Access	Full Access	Full Access	Full Access	Full Access





Database Tab (Entire Tab)	Full Access	Visible Only	Visible Only	Visible Only	Visible Only
Advanced Tab	Full Access	Full Access	Visible Only	Visible Only	Full Access
Advanced – IT Settings	Full Access	Full Access	Visible Only	Visible Only	Visible Only
License Tab	Full Access	Full Access	Visible Only	Visible Only	Visible Only
License: view	Full Access	Full Access	Full Access	Full Access	Visible Only
License: Import/add	Full Access	Full Access	Not Visible or Editable	Not Visible or Editable	Visible Only
About	Full Access	Full Access	Full Access	Full Access	Full Access

## 13.7 About



About shows the installed version of SpecSuite along with contact information for Metrohm and Metrohm Service.

A list of open-source licenses used to create SpecSuite can be viewed here.



## 14. Terminologies

**Scan** – Any time the detector is recording information. Scan will usually refer to the collection of data while an external light source is firing.

**Dark Scan** – A scan collected while the light source is not firing.

**“Complete” Scan** – Data scan + Dark scan

**Check Signal** – This is known as “Continuous Scan” in BWSpec4. The laser is on and collecting scans consecutively at the set integration time until the Stop command is given.

**Integration Time** – The duration of time that the detector is collecting data.

**Auto-Time** – The optimal integration time is automatically calculated by the software to achieve a maximum signal of 40,000 to 60,000 counts .

**Average** – Used to describe the averaging of 2 or more spectra. The average is done for each y-value.

**Laser Power** – Refers to the percentage of total laser power being used.

**External Trigger** – When the trigger function is enabled, the instrument will wait for a trigger signal to be received at the External Trigger port. As an example, some i-Raman probes have a small trigger button on them which can trigger the scan start.

**Auto Dark** – When enabled, a dark scan of the same integration time as the data scan will be collected immediately after the data scan without the user having to take any additional action.

**cCode** – A 3–6-digit long string of letters used to name a spectrometer. This is a legacy term and will be depreciated for utilizing a Metrohm defined Serial Number instead.


**Ratio3** – Refers to the output file generated by Relative Intensity Correction. It is a means of regulating data between spectrometers. Ratio3 is a .TXT file which includes information like cCode and coefficients. The data is two columns, one of Pixel and one for the multiplication factor of that pixel. The y-value of each pixel is multiplied by the Ratio value









Select 'Perform a new installation of SQL Server 2022' then press next.




Accept the license terms and privacy statement for SQL Server 2022 Express and press Next.








Name the instance something unique, then press Next.




There are no necessary changes to make at this step, press Next.





The installation should be completed successfully.





## B.2 Peak Analysis Algorithms

SpecSuite provides three algorithms for peak analysis, each suited to different use cases:

- Data Points – A local peak value-based method that identifies peaks directly from the measured data points.
- Parabolic – A parabolic fitting method that uses minimum height and width thresholds to identify and characterize peaks.
- CurveFit – A nonlinear curve fitting method supporting Gaussian, Lorentzian, or Voigt peak shapes.




Figure 2: Data Points vs. Parabolic peak detection.

## B.3 Parabolic – FWHM

When using the Parabolic algorithm, the Full Width at Half Maximum (FWHM) is calculated as the difference between the right mid-height point and the left mid-height point of the peak.



## B.4 Curve Fit – Gaussian Distribution

The Gaussian model is often called normal distribution, which is assuming the peak shape follows Gaussian distribution.

$$y = offset + Ae^{-\frac{(x-\mu)^2}{2\sigma^2}}$$

Where:

- $A$ : The maximum amplitude (peak height) of the curve
- $\mu$ : The center position of the peak
- $\sigma$ : The standard deviation, which determines the width or "spread" of the bell curve

In this distribution,

$$\text{Gaussian } FWHM = 2\sqrt{2\ln 2} \cdot \sigma \approx 2.35482 \sigma$$

This model is commonly applicable for atomic emission spectroscopy (e.g., low-pressure Hg) and NIR molecular spectroscopy.




Figure 4: Gaussian distribution – equation parameters (left) and probability density curves (right).

## B.5 Curve Fit – Lorentzian Distribution

The Lorentzian model is also called Cauchy distribution:

$$y = offset + \frac{A}{1 + \left(\frac{x-\mu}{\gamma}\right)^2}$$

Where:



## B.6 Curve Fit – Voigt Distribution

The Voigt model uses a pseudo-Voigt blend (not a full convolution). It combines a Gaussian component and a Lorentzian component using a mixing ratio  $\eta \in (0,1)$ .

$$\eta = \frac{1}{1+e^{-p^4}} \text{ where } \eta \in (0,1)$$

Where:

- Gaussian component:  $g = e^{-\frac{(x-\mu)^2}{2\sigma_G^2}}$ , with  $\sigma_G = \frac{\text{shared FWHM}}{2\sqrt{2\ln 2}}$
- Lorentzian component:  $l = \frac{1}{1+(\frac{x-\mu}{\gamma_L})^2}$ , with  $\gamma_L = \frac{\text{shared FWHM}}{2}$
- Final peak model:  $y = A \cdot (\eta \cdot l + (1 - \eta) \cdot g)$

Both components share a common FWHM. This model is suitable for mixed environments such as high-pressure Hg and Raman mixture analysis.




Figure 6: Voigt distribution – equation parameters (left) and probability density curves (right).

## B.7 Key Points for Good Curve Fitting

Curve fitting uses initial peak parameters to reconstruct the measured curve by modelling it as a combination of fitted peak shapes. The quality of the fit depends on:

- Initial peak parameters: the number of peaks and the starting values for each peak (position, intensity, width). It defines how peaks are expected in the fitting which will eventually contribute to final fitted curve. Incorrect initial parameters include:

1. Peak exists in original data but not included in initial list





9	Mass (m/z)		9	Millivolts	
10	Parts per million (ppm)		10	Log(1/R)	
11	Days		11	Percent	
12	Years		12	Intensity	Processed
13	Raman shift (cm <sup>-1</sup> )	RamanShift	13	Relative intensity	
14	eV		14	Energy	
16	Diode number		16	Decibel	
17	Channel		19	Temperature (°F)	
18	Degrees		20	Temperature (°C)	
19	Temperature (°F)		21	Temperature (K)	
20	Temperature (°C)		22	Index of refraction (N)	
21	Temperature (K)		23	Extinction coefficient (K)	
22	Data points	Pixel	24	Real	
23	Milliseconds		25	Imaginary	
24	Microseconds		26	Complex	
25	Nanoseconds		128	Transmission	Transmission/Reflection
26	Gigahertz (GHz)		129	Reflectance	
27	Centimeters (cm)		130	Arbitrary	
28	Meters (m)		131	Emission	
29	Millimeters (mm)				



30	Hours				
255	Double interferogram				

\*Y axis option Absorbance and Transmission are to be implemented.