



CH-9101 Herisau/Switzerland

Phone ++41 71 53 85 85

Fax ++41 71 53 89 01

736 GP Titrino

Short instructions for use

8.736.1023

95.12 Ti/pr

Table of contents

	Page
1. Getting started	1
1.1 Principle of data input	1
1.2 Development of a method.....	3
1.3 Storing methods	7
1.4 Working with several dosing units.....	10
2. Titration, dosing and measuring modes	12
2.1 Survey of the titration modes	12
2.2 Mode DET and MET, parameters	14
2.3 Mode SET, parameters.....	15
2.4 Mode KFT, parameters.....	16
2.5 Mode STAT, parameters	17
2.6 Mode DOS, parameters.....	19
2.7 Mode DOC, parameters.....	20
2.8 Mode MEAS and CAL, parameters.....	21
2.9 Mode TIP, parameters	22
2.10 Calculations	23
3. Other functions	25
3.1 Method memory.....	25
3.2 Sample data, keys <smpl data> and <silos>	27
3.3 Configuration, key <configuration>	28
Index	29

1 Getting started

This is a step by step introduction into the operation of the Titrino.

You learn how to go round with the Titrino and you get an idea of the possibilities that the Titrino offers you.

1.1 Principle of data input

Entering data, setting the dialog language

<configuration>

>peripheral units

Press <STOP> if the Titrino is busy. It is now in the inactive standby state. Then press the <configuration> key. The display shows:

This is the title of the inquiry group "peripheral units". This group contains various inquiries about peripheral units.

Press the <configuration> key once more. You see the title of a new group of inquiries:

>auxiliaries

This is the title of the "auxiliaries" group. This group contains the inquiry for the dialog language. Note the ">" sign. All titles are prefixed by this sign.

<enter>

dialog: english

Pressing the <enter> key takes you to the different inquiries of the group "auxiliaries". The display shows

<select>

dialog: deutsch

This is the selection of the dialog language. You select the various dialog languages with the <select> key. Note the ":" sign after the word "dialog". It appears if the values can be selected with key <select>. Press <select> repeatedly until "english" appears in the display.

<enter>

Accept the new "value" with <enter> and the next inquiry appears.

date 1997 - 12 - 16

Other inquiries of the "auxiliaries" group are not of interest at the moment. Exit the inquiries with <QUIT>. Press <QUIT> once again to exit the inquiries of the <configuration> key.

2 x <QUIT>

Summary of the most important signs

Signs in the Titrino dialog:

- > Title of inquiry groups, e.g. ">auxiliaries"
- : Values can be chosen by the key <select>, e.g. "dialog:"

Signs that are used in the instructions for use

- <> sign for "key", e.g. <select>

1.2 Development of a method

You learn how to elaborate a method .
 First you select a titration mode. A general survey of the different titration modes is given on page 12.

Selecting the mode

Press <mode> repeatedly until "DET" appears in the display. DET stands for **D**ynamic **E**quivalence point **T**itration. In this titration, the size of the volume increments vary as a function of the slope of the titration curve. After each increment, a pre-set measured value drift (or time) is awaited until the next increment is added (so called "equilibrium titration"). After the titration, the equivalence points are evaluated automatically. Confirm "DET" with <enter>.

<mode>

mode	DET
-------------	------------

 <enter>

DET: U

Now select the measured quantity: Press <select> until "pH" appears in the display.

<select>

DET:	pH
------	----

Confirm the measured quantity "pH" with <enter>.

<enter>

DET pH	D0a	*****
--------	-----	-------

You are now ready to titrate.

For the titration, you need an Exchange Unit with NaOH, $c=0.1$ mol/L. If you have not yet prepared an Exchange Unit, see page 188 ff of the Instructions for Use.

Plug a combined pH glass electrode into measuring input ("Ind I").

Pipette 2 mL HCl, $c =0.1$ mol/L into the titration vessel, dilute it with ca. 20 mL. dist. water, adjust the stirrer and press <START>.

During the titration, the first line of the display shows the current measured value and the volume already dispensed:

pH	2.78	0.351 ml
----	------	----------

As soon as the instrument has found an equivalence point, this is shown on the second line:

	EP1
--	-----

Let the titration continue for a short while, e.g. until ca. pH=11.50. Now stop it with <STOP>.

The second line shows the equivalence point found, e.g.:

EP1	2.083 ml	pH	7.64
-----	----------	----	------

If more than one equivalence point has been found, the second equivalence point can be viewed with <enter>

Calculation of the result, entry of a formula

You may calculate results with the equivalence point.

<def>

Press <def> . The display shows

<enter>

Press <enter> to move on to the formula entry. The display now shows "RS? "

<1>

Press "1", i.e. the first formula.

You can now enter a formula. Note here the top inscription on the keys of the keypad and the numbers.

<enter>

You can use mathematical operations, parentheses as well as the following symbols:

EP# EP's with 1-digit number, e.g. EP1.

RS# Previously calculated results, e.g. RS1 in the second formula.

C## Calculation values, e.g. C01. C00 is reserved for the sample size. The meaning of the different calculation values is given on page 23.

Calculate the content of your hydrochloric acid in g/L:

$$RS1=EP1*C01*C02/C00$$

$$EP*conc(titrant)*molar\ mass/sample\ size$$

Confirm the formula with <enter>.

<enter>

You may enter a text for the result output, see page 7.

<enter>

Enter the desired number of decimal places for the result.

<enter>

Select the unit g/l with <select> or enter a text as unit.

<select>

<enter>

Quit the formula entry by pressing <QUIT> twice.

2 x <QUIT>

Entry of the calculation constants

<C-fmla>

Press <C-fmla>. The constants which have been used in the formula are requested:

0.1 <enter>

C01: Concentration of your titrant = 0.1 mol/L.

C02: Molar mass of HCl = 36.47 g/mol

36.47 <enter>

Entry of the sample size

4 x <smpl data> Press <smpl data> repeatedly until "smpl size" appears in the display.

smpl size	1.0 g
-----------	-------

Enter 2 (= 2 ml).

2 <enter>

smpl unit:	ml
------------	----

Use <select> to select the unit "ml" and confirm the new value with <enter>.

<select>
<enter>

The result is now calculated and can be displayed in place of the equivalence point. If your method already includes a formula at the end of the titration, the calculated result is displayed directly after the titration. As we have entered the formula later, we now have to select the result display:

Display of the result

<select>

>display results

Press <select> repeatedly until ">display results" appears in the display.

<enter>

Press <enter> to move on to the result display.

RS1	3.80 g/l
-----	----------

Selection of the automatic report output

If you have a printer connected, you can have the curve and a result report printed out automatically at the end of the titration.

4 x <def>

>report

Press <def> repeatedly until the display shows ">report".

<enter>

Press <enter> to move to the definition of reports.

report:

With <select> you select the individual report blocks Use a ";" as delimiter between the report blocks. If you wish to print out a curve and a full result report, enter "curve;full".

<select>

report:curve;full

<enter>
<QUIT>

Confirm the entry with <enter> and quit the inquiries with <QUIT>

1.3 Storing methods

You learn how to operate the method memories.

The Titrino offers 2 types of method memories:

- Internal method memory, key <user meth>
- External method memory on the card, key <card>

Storing a method in the internal method memory

```

2 x <user meth>
>store method
<enter>
method name: *******
< < >
method name: ABCDEFGH
< < > + < > >
<enter>
method name: HIJKLMNO
:
method name: HClmnopq
<QUIT>
method name: HCl
<enter>
DET pH    DOa    HCl
    
```

Press key <user meth> repeatedly until ">store method" appears in the display.

Enter an identification for your method, e.g. HCl:

Press key < < >. The display shows a row of letters. The first position of this row is blinking. This is the input position.

You can move this row of letters pressing the < < > and < > >. Press < > > until an "H" appears in the blinking position and confirm it with <enter>.

The next position is now blinking. Enter the next letter here, confirm it with <enter> etc.

If the identifier is complete, quit the text input with <QUIT>.


Confirm the value with <enter>. The method is stored in the memory.

The method runs now under the identifier "HCl".

A list of contents of the method memory can be printed with the key sequence <print><user meth><enter>.

```

'um
736 GP Titrino          02103  736.0010
date 1997-09-13        time  07:49
user methods           bytes
DET pH      Titer_pH DO    190
DET U       Blank DO     178
DET U       Chloride DO   234
MET U       Diazo DO     208
MET Ipo1    Br-Index DO   226
DET U       Ca-Mg DO     278
DET U       EDTA-NTA DO   216
DET U       Metals DO    190
DET U       Perox.No DO   172
DET pH      HCl DO       108
remaining bytes       7330
-----
    
```



Recall a method from the internal method memory

Stored methods can always be recalled into the working memory.

<p><user meth></p> <div style="border: 1px solid black; padding: 2px; display: inline-block;">>recall method</div>	Press the key <user meth>. In the display appears ">recall method".
<p><enter></p> <div style="border: 1px solid black; padding: 2px; display: inline-block;">method name: HCl</div>	You can choose a method with <select> or enter its name directly.
<p><select></p> <p><enter></p>	Recall the method with <enter>.
<div style="border: 1px solid black; padding: 2px; display: inline-block;">DET pH D0a HCl</div>	The method is now ready to work.

Method memory on the card

The card offers additional memory space for methods. It may be used e.g.

- as a backup-medium for the internal method memory
- as an extended method memory, where the methods can be organized in various directories
- if every user wants to store his methods on his own card
- for exchanging methods between various laboratories, apparatus, and/or users

The card has a battery, which should be exchanged periodically, see page 85 of the Instructions for Use.

As an example, perform a backup from the internal method memory to the card. If the card is not formatted, see page 86 of the Instructions for Use.

<p>7 x <card></p> <div style="border: 1px solid black; padding: 2px; display: inline-block;">>backup</div>	Press the key <card> repeatedly until ">backup" appears in the display.
<p><enter></p> <div style="border: 1px solid black; padding: 2px; display: inline-block;">directory:</div>	Confirm this function with the key <enter> and give a name for the backup directory. If you wish to overwrite an existing directory, choose a name with <select>.
<p>name + <enter></p>	The backup will be effected.

Printing a list of contents of the card

The following card reports can be printed:

Current directory
 <print><card><enter>

name of the current directory

methods in this directory

remaining bytes on this card

```
'ad
736 GP Titrino          02103  736.0010
date 1997-09-13      time 07:49
card label: appl.736 dev.label:
KF Titr.                      bytes
KFT Ipol      TarTiter D0      150
KFT Ipol      H2OTiter D0      150
KFT Ipol      Blank_KF D0      132
KFT Ipol      KF-Blank D0      206
KFT Ipol      KF D0            170
                                remaining bytes      121892
                                -----
```

Contents of the whole card

<print><select><enter>

Press <select> repeatedly until "mem card" appears in the display.

bytes, that are being used in the directory "user meth."

the directory "KF Titr." contains 5 methods

```
'cd
736 GP Titrino          02103  736.0010
date 1997-09-13      time 07:50
card label: appl.736 dev.label:

directory: user meth.
Acid      Tit.NaOH      Tit.HCl
W-Liquor
bytes: 620

directory: KF Titr.
TarTiter  H2OTiter      Blank_KF
KF-Blank  KF
bytes: 808
                                remaining bytes on card: 121892
                                -----
```

1.4 Working with several dosing units

Here you will learn how to work with several dosing units.

Several dosing units can be used

- if a working station will be permanently installed for specific analysis, e.g. having a dosing unit specially prepared for Karl Fisher titrations
- if several titrations are being carried out in the same sample. The linking of various methods is possible with TIP (Titration Procedure), see page 81 of the Instructions for Use.
- for addition of auxiliary solutions

Selection of the dosing unit in the titration method

The dosing unit is selected in the method. The method runs always with the selected dosing unit.

Connect the dosing unit D1 and arm it with an Exchange Unit (the Dosino with a Dosing Unit, resp.).

<user meth>

>recall method

<enter>

method name:

<select>, <enter>

DET pH D0 HCl

n x <parameters>

>titration parameters

m x <enter>

dos.element:internal D0

<select>

<enter>

2 x <QUIT>

DET pH D1 HCl

<START>

DET pH D1a HCl

Recall the method "HCl" from the user memory: Press the key <user meth>. In the display appears ">recall method". Press <enter>.

Select the method name "HCl" with the key <select> and recall the method with <enter>.

The dosing unit of this method is D0, i.e. the Titrino itself.

To change the dosing unit, press the key <parameters> repeatedly until ">titration parameters" appears in the display.

Press <enter> until "dos.element:" appears. With <select> you can choose between

internal D0: Dosing unit of the Titrino

external D1: externally connected dosing unit 1, or 2, resp.

Choose e.g. "external D1", confirm the input with <enter> and quit with the key <QUIT>.

The display shows now e.g. "DET pH D1 HCl". D1 means, that in this method the external dosing unit D1 is being used.

Start the method.

After the titration "D1a" appears in the display. "a" means "active", the dosing unit D1 is now active.

Manual dosing

Manual dosing can be used

- to prepare an Exchange Unit (or Dosing Unit, resp.)
- for quick dosing

If you press the key <DOS>, the dosing unit from the method is working. This dosing unit appears in the display.

If you wish to dose with another dosing unit, use the key <prep> to select the dosing unit and press <DOS> to dose.

<prep>

internal D0	prep
-------------	------

Press the key <prep>. In the display appears

"internal D0 prep", i.e. the dosing unit of the Titrino has been chosen.

<prep>

external D1:	prep
--------------	------

Press the key <prep> again. In the display appears

"external D1: prep".

<DOS>

<FILL>

Press the key <DOS> on the Titrino: The external dosing unit D1 doses. If you press the key <STOP/FILL> on the Titrino, it will fill again.

<START>

With <START>, the function "prep" is started. This prepares the dosing units. (The parameters for "prep" are under the key <configuration>, in the group ">prep.dosing elements").

Attention: With the function "prep", liquid will be expelled!

<select>

external D1:	empty
--------------	-------

If the dosing unit D1 is a Dosino, you can empty the dosing unit completely (e.g. for cleaning). Choose the function "empty" with the key <select> and start with <START>.

Attention: With the function "empty", liquid will be expelled!

<START>



<QUIT>

Quit the selection of dosing units with <QUIT>.




2 Titration, dosing and measuring modes

This chapter gives you a survey over the titration, dosing and measuring modes. You will find the parameters of all modes. The standard values are printed in bold face in the following tables.

2.1 Survey of the titration modes

	DET Dynamic Equivalence point Titration	MET Monotonic Equivalence point Titration
Titration	Reagent feeding: Variable volume increments, depending on the slope of the curve.  Drift controlled ("equilibrium titration") and/or after a fixed equilibration time.	Reagent feeding: Constant volume increments, independent of the slope of the curve.  Drift controlled ("equilibrium titration") and/or after a fixed equilibration time.
Evaluation	The evaluation of EP's is based on the zero crossing of the second derivative with a Metrohm correction for the distortion of the curve from superimposed jumps. Can be combined with selectable recognition criteria. Recognition criteria: as for MET	The evaluation of EP's is based on the Fortuin interpolation. Recognition criteria: all EP's only the last EP only the greatest EP EP windows
Applications	Suitable titration mode for most problems. Specially recommended if jumps lie very close together and for very flat jumps. Note: The reagent feeding algorithm is based on measured data. The curve should therefore not deviate markedly from S-shape.	For <ul style="list-style-type: none"> • slow titration reactions (diazotations, coupling reactions) • sluggish electrode response.

If you are looking for more information to the different quantities, you will find the display text as a keyword in the index of the Instructions for Use.

SET, KFT Endpoint Titration	STAT Holding a measured value	DOC Controlled dosing
<p>Titration to preset end-point.</p> 	<p>The reagent addition is controlled such that a preset measured value is adhered to as accurately as possible.</p> 	<p>The dosing function is controlled such that a nominal measured value is adhered to as accurately as possible. The nominal measured value follows a linear ramp.</p> 
<p>Volume that has been dispensed up to the endpoint (EPX in ml).</p>	<p>Evaluation of</p> <ul style="list-style-type: none"> • rates dV/dt in the set time windows and the mean rate over the entire measuring point list. • volumes associated with a preset time (fix V). • times associated with a preset fraction of the total volume (fix time). 	<p>The total volume (C41) and the time needed for dosing (C42) can be used for calculation.</p>
<ul style="list-style-type: none"> • For rapid, quantitative determinations in analytical chemistry. Requirement: EP of the titration reaction is known and does not change during a determination series. • If an excess of titrant must be avoided. • KFT: Determinations of water contents. 	<p>For</p> <ul style="list-style-type: none"> • the controlled dosing to a preset endpoint in the synthesis laboratory. • activity determinations, e.g. of enzymes. • kinetic investigations (first order kinetics). 	<p>Sweeping a measured value gradient for reaction formulations in the syntheses laboratory.</p>

2.2 Mode DET and MET, parameters

Display	Meaning	Input range
>titration parameters meas. pt. density min. incr. [V step] titr. rate signal drift equilibr. time start V: start V factor dos. rate pause dos. element: meas. input: temperature	Measuring point density. Minimum increment for DET. Volume increment for MET. Dosing rate for increments. Drift for measured value acquisition. Waiting time for value acquisition. Type of start volume: absolute, relative. Volume for absolute start volume. Factor for calculation of relative start volume: factor * smpl size. Dosing rate for start volume. Waiting time, e.g. after start volume. Selection of the dosing unit. Measuring input for pH and U. Titration temperature.	0...4...9 0...10.0...999.9 ul 0...0.10...9.999 ml 0.01...150 ml/min, max. pH, U:0.5..50..999 mV/min, OFF 0...26...9999 s, OFF abs., rel., OFF 0...999.99 ml 0...±999 999 0.01...150 ml/min, max. 0...999 999 s internal D0, external D1/D2 1, 2, diff. -170.0...25.0...500.0 °C
>stop conditions stop V: stop V factor stop pH stop EP filling rate	Type of stop volume: absolute or relative. Volume for absolute stop volume. Factor for calculation of relative stop volume: factor * smpl size. Stop at measured value pH, U, I. Stop after a number of EP's have been found. Filling rate after the titration.	abs., rel., OFF 0...99.99...9999.99 ml 0...±999 999 pH: 0.00...±20.00, OFF U: 0...±2000 mV, OFF 1...9, OFF 0.01...150 ml/min, max.
>statistics status: mean n= res. tab: delete n=	Status of statistics calculation. Number n of single values for statistics. Result table for statistics calculation. Delete data from sample number n.	ON, OFF 2...20 original, delete n, delete all 1...20
>evaluation EPC EP recognition: fix EP1 at pH pK/HNP:	Endpoint criterion. Type of EP recognition: all EP's, only greatest EP, only last EP, only EP's in windows, no evaluation. If "window" is selected, lower and upper limits of windows are inquired. Interpolation of volume at given pH, U, or I, resp. Up to 9 fix EP's. Evaluation of pK or half neutralization potential (HNP).	DET 0...5...200 MET pH: 0.10...0.50...9.99 U: 1...30...999 mV all, greatest, last, window, OFF pH: 0.00...±20.00, OFF U: 0...±2000 mV, OFF ON, OFF
>preselections	Preselections for the sequence see DOS, page 19	

2.3 Mode SET, parameters

Display	Meaning	Input range
>SET1 EP at pH dynamics max. rate min. rate stop crit: stop drift t(delay) stop time	Individual parameters for EP1. Preset EP1 at pH, U, or I, resp. Distance from EP where constant dosing stops and controlling begins. Maximum dosing rate. Minimum dosing rate. Type of stop criteria. Titration stops if stop drift is reached. Stop if there is no dosing during t(delay). If t(delay) is "INF": stop after a time.	pH: 0.00...±20.00, OFF (0...±2000 mV, 0.0...±200.0 uA) pH: 0.01...20.00, OFF (1...2000 mV, 0.1...200.0 uA) 0.01...10...150 ml/min, max. 0.01...25.0...999.9 ul/min drift, time 1...20...999 ul/min 0...10...999 s, INF 0...999 999 s, OFF
>SET2	Parameters for EP2. Identical as SET1.	
>titration parameters titr. direction: pause 1 start V: start V factor dos. rate pause 2 extr. time dos. element meas. input: temperature	General titration parameters. +: Titration to higher pH, voltage, or current. Auto: Direction is set automatically. Waiting time before start volume. Type of start volume: absolute or relative. Volume for absolute start volume. Factor for calculation of relative start volume: factor * smpl size. Dosing rate for start volume. Waiting time after start volume. Extraction time. Selection of the dosing unit. Measuring input for pH and U. Titration temperature.	+, -, auto 0...999 999 s abs., rel., OFF 0...999.99 ml 0...±999 999 0.01...150 ml/min, max. 0...999 999 s 0...999 999 s internal D0, external D1/D2 1, 2, diff. -170.0...25.0...500.0 °C
>stop conditions stop V: stop V factor filling rate	Type of stop volume: absolute or relative. Volume for absolute stop volume. Factor for calculation of relative stop volume: factor * smpl size. Filling rate after titration.	abs., rel., OFF 0...99.99...9999.99 ml 0...±999 999 0.01...150 ml/min, max.
>statistics	Statistics calculation see at DET, page 14	
>preselections conditioning: display drift: drift corr: drift value req. ident: req. smpl size: activate pulse:	Automatic conditioning of titration vessel Display of drift during conditioning. Type of drift correction. Value for manual drift correction. Request of identifications after start. Request of sample size after start. Pulse output on I/O line L6.	ON, OFF ON, OFF auto, man., OFF 0.0...99.9 ul/min id1, id1 & 2, all, OFF value, unit, all, OFF first, all, cond., OFF

2.4 Mode KFT, parameters

Display	Meaning	Input range
>control parameters EP at ... dynamics max. rate min. volume incr. stop crit: stop drift t(delay) stop time	Controlling parameters for KFT. Endpoint for I _{pol} . Endpoint for U _{pol} Distance from EP where constant dosing stops and controlling begins. Maximum dosing rate. Minimum volume increment Type of stop criteria. Titration stops if stop drift is reached. Stop if there is no dosing during t(delay). If t(delay) is "INF": stop time.	U: -2000 ... 250...2000 mV I: -200.0... 25.0...200.0 uA U: 1...100...2000 mV I: 0.1...10...200.0 uA 0.01...150 ml/min, max. 0.1...9.9 ul, min drift, time 1...20...999 ul/min 0...10...999 s, INF 0...999 999 s, OFF
>titration parameters titr. direction: pause 1 start V: start V factor dos. rate pause 2 extr. time dos. element I(pol) U(pol) electrode test: temperature	General titration parameters. +: Titration to higher voltage or current. auto: Direction is set automatically. Waiting time before start volume. Type of start volume: absolute or relative. Volume for absolute start volume. Factor for calculation of relative start volume: factor * smpl size. Dosing rate for start volume. Waiting time after start volume. Extraction time. Selection of the dosing unit. Polarization current for ipol or - voltage for upol in 10 mV-intervals. Performing of electrode test Titration temperature.	+, -, auto 0...999 999 s abs., rel., OFF 0...999.99 ml 0...±999 999 0.01...150 ml/min, max. 0...999 999 s 0...999 999 s internal D0, external D1/D2 -127...50...127 uA bzw. -1270...400...1270mV ON, OFF -170.0...25.0...500.0 °C
>stop conditions stop V: stop V factor filling rate	Type of stop volume: absolute or relative. Volume for absolute stop volume. Factor for calculation of relative stop volume: factor * smpl size. Filling rate after the titration.	abs., rel., OFF 0...99.99...9999.99 ml 0...±999 999 0.01...150 ml/min, max.
>statistics	Statistics calculation see at DET, page 14	
>preselections conditioning: display drift: drift corr: drift value req. ident: req. smpl size: activate pulse:	Automatic conditioning of titration vessel. Display of drift during conditioning. Type of drift correction. Value for manual drift correction. Request of identifications after start. Request of sample size after start. Pulse output on I/O line L6.	ON, OFF ON, OFF auto, man., OFF 0.0...99.9 ml/min id1, id1 & 2, all, OFF value, unit, all, OFF first, all, cond., OFF

2.5 Mode STAT, parameters

Display	Meaning	Input range
>control parameters EP at pH dynamics max. rate min. rate	Endpoint, control point: This value is kept constant. Distance from EP where constant dosing stops and controlling begins. Maximum dosing rate Minimum dosing rate	pH: 0.00...±20.00, OFF (0...±2000 mV, 0.0...±200.0 uA) pH: 0.01...1...20.00, OFF (1..60..2000 mV, 0.1.6..200.0uA) 0.01...10...150 ml/min, max. 0.01...25.0...999.9 ul/min
>titration parameters start V: start V factor dos. rate pause start time start pH start rate time interval titr. direction: dos. element meas. input: temperature	General titration parameters. Type of start volume: absolute or relative. Volume for absolute start volume. Factor of relative start volume. Dosing rate for start volume. Waiting time, after start volume. Start time for data acquisition. Start measured value for data acquisition. Start rate for data acquisition. Time interval for data acquisition. +: Titration to higher pH, U or I. auto: Direction is set automatically. Selection of dosing unit. Measuring input for pH and U. Temperature.	abs., rel., OFF 0...999.99 ml 0...±999 999 0.01...150 ml/min, max. 0...999 999 s 0...999 999 s pH: 0.00...±20.00, OFF (0...±2000 mV, 0.0...±200.0 uA) 0.01...150 ml/min, OFF 1...2...999 999 s +, -, auto internal D0, external D1/D2 1, 2, diff. -170.0...25.0...500.0 °C
>stop conditions stop time: stop time factor t(delta) t(delay) stop V: stop V factor stop rate filling rate	Type of stop time: absolute or relative. Time for absolute stop time. Factor for calculation of relative stop time. Time after EP is once reached. Time after last dosing step. Type of stop volume: absolute or relative. Volume for absolute stop volume. Factor of relative stop volume. Stop when the dosing rate is smaller than.. Filling rate.	abs., rel., delta, delay, OFF 0...999 999 s 0...±999 999 0...999 999 s 0...999 999 s 0...999 999 s abs., rel., OFF 0...99.99...9999.99 ml 0...±999 999 0.01...150 ml/min, OFF 0.01...150 ml/min, max.
>statistics	Statistics calculation see at DET, page 14	
>evaluation low lim up lim fix V1 fix time 1	Evaluation of dosing rates within programmed time windows (→ C8X). Max. 9 windows. Interpolation of volumes at fixed times (→ C5X). Up to 9 fix V's. Interpolation of times at fixed volume rations of the end volume (→ C6X). Up to 9 fix times.	0...999 999 s, OFF 0...999 999 s, OFF 0...999 999 s, OFF 0.01...1, OFF

Display	Meaning	Input range
>moni toring meas. val: low lim pH up lim pH action: rate: temperature: assign output: L4 output:	Monitoring of measured values. Limits for measured values. Action, if a limit is exceeded. End: Abort. Wait: Dosing interrupted, until measured value is within limits again, then continue. Monitoring of rate (same actions as above). Monitoring of temperature (same actions as above). Assignment of output L4 for monitoring. Set a signal when values are out of limit.	ON, OFF pH: 0.00...±20.00 (0...±2000 mV, 0.0...±200.0 uA) end, hold, wait, none ON, OFF ON, OFF meas, temp, rate, all, none active, pulse
>preselections req. ident: req. snpl size: display rate: activate pulse:	Request of identifications after start. Request of sample size after start. Display of rate. Pulse output on I/O line L6.	id1, id1 & 2, all, OFF Value, unit, all, OFF ON, OFF ON, OFF

2.6 Mode DOS, parameters

Display	Meaning	Input range
>dosing parameters dispensing type: volume disp. crit: rate disp. time pause time interval dos. element temperature	Volume, time or rate. Inquiry of 2 values. 3rd value is calculated . 1st criterion i.e. volume. Selection of 2nd criterion. 2nd value, i.e. rate or dispensing time. Waiting time before start. Time interval for data acquisition. Selection of the dosing unit Temperature.	volume, time, rate 0...10...99 999.99 ml time, rate 0.001...150 ml/min, max. 1...100...999 999 s 0...999 999 s 1...10...999 999 s internal D0, external D1/D2 -170.0...25.0...500.0 °C
>stop conditions stop V: stop V factor filling rate	Type of stop volume: absolute or relative. Volume for absolute stop volume. Factor of relative stop volume: Filling rate.	abs., rel., OFF 0...99.99...99 999.99 ml 0...±999 999 0.01...150 ml/min, max.
>statistics	Statistics calculation see at DET, page 14	
>monitoring meas. mode: meas. input: low lim pH up lim pH action: temperature: low lim up lim action: assign output: L4 output:	Measured quantity for monitoring. Meas. input for pH and U. Limits for measured values. Action, if a limit is exceeded. End: Abort. Wait: Dosing interrupted, until measured value is within limits again, then continue. Monitoring of temperature. Limits for temperature. Action, if a limit is exceeded . Assignment of output L4 for monitoring. Set a signal when values are out of limit.	pH, U, I(pol), U(pol), OFF 1, 2, diff pH: 0.00...±20.00 (0...±2000 mV, 0.0...±200.0 uA) end, hold, wait, none ON, OFF -170.0...500.0 °C -170.0...500.0 °C end, hold, wait, none meas, temp, all, none active, pulse
>preselections req. ident: req. snpl size: activate pulse:	Request of identifications after start. Request of sample size after start. Pulse output on I/O line L6.	id1, id1 & 2, all, OFF value, unit, all, OFF ON, OFF

2.7 Mode DOC, parameters

Display	Meaning	Input range
>dosing parameters begin at end at sweep time dynamics max. rate min. rate direction: start V: start V factor dos. rate pause time interval dos. element meas. input: temperature	Begin of measured value sweep. End of measured value sweep. Time required for sweep. Distance from EP where constant dosing starts. OFF means no constant dosing. Maximum dosing rate. Minimum dosing rate. Direction of sweep. +: Titration to higher pH, voltage, or current. Auto: Direction is set automatically. Type of start volume: absolute or relative. Volume for absolute start volume. Factor for of relative start volume: Dosing rate for start volume. Waiting time, after start volume. Time interval for data acquisition Selection of dosing unit. Meas. input for pH and U. Temperature.	pH: 0.00...±20.00, init (0...±2000 mV, 0.0...±200.0 uA) pH: 0.00...±20.00, OFF (0...±2000 mV, 0.0...±200.0 uA) 0...300...999 999 s pH: 0.01...0.25...20.00, OFF (1...15...2000 mV, 0.1...6...200 uA) 0.01...150 ml/min, max. 0.01...5...999.9 ul/min +, -, auto abs., rel., OFF 0...999.99 ml 0...±999 999 0.01...150 ml/min, max. 0...999 999 s 1...10...999 999 s internal D0, external D1/D2 1, 2, diff. -170.0...25.0...500.0 °C
>stop conditions stop V: stop V factor filling rate	Type of stop volume: absolute or relative. Volume for absolute stop volume. Factor of relative stop volume. Filling rate.	abs., rel., OFF 0... 99.99 ...9999.99 ml 0...± 999 999 0.01...150 ml/min, max.
>statistics	Statistics calculation see at DET, page 14	
>monitoring meas. val: low lim pH up lim pH action: temperature: low lim up lim action: assign output: L4 output:	Monitoring of measured value. Limits for measured value. Action, if a limit is exceeded. End: Abort. Wait: Dosing interrupted, until measured value is within limits again, then continue. Monitoring of temperature. Limits for temperature. Action, if a limit is exceeded . Assignment of output L4 for monitoring. Set a signal when values are out of limit.	ON, OFF pH: 0.00...±20.00 (0...±2000 mV, 0.0...±200.0 uA) end, hold, wait, none ON, OFF -170.0...500.0 °C -170.0...500.0 °C end, hold, wait, none meas, temp, all, none active, pulse
>preselections req. ident: req. snpl size: display rate: activate pulse:	Request of identifications after start. Request of sample size after start. Display of rate. Pulse output on I/O line L6.	id1, id1 & 2, all, OFF value, unit, all, OFF ON, OFF ON, OFF

2.8 Mode MEAS and CAL, parameters

MEAS

Display	Meaning	Input range
>measuring parameters signal drift equilibr. time meas. input: temperature	Drift for measured value acquisition. Waiting time for meas. value acquisition. Measuring input for pH and U. Measuring temperature.	pH, U: 0.5..999 mV/min, OFF T: 0.5...999 °C/min, OFF 0...9999 s, OFF 1, 2, diff. -170.0...25.0...500.0 °C
>statistics	Statistics calculation see at DET, page 14.	
>preselections req. ident: req. snpl size: activate pulse:	Request of identifications after start. Request of sample size after start. Pulse output on I/O line L6.	id1, id1 & 2, all, OFF value, unit, all, OFF ON, OFF

CAL

Display	Meaning	Input range
>cal. parameters meas. input: cal. temp. buffer 1 pH signal drift equilibr. time electr. id sample changer cal: activate pulse:	Measuring input. Calibration temperature. pH-value for buffer 1. up to 9 buffer. Drift for measured value acquisition. Waiting time for meas. value acquisition. Electrode identification. Calibration with Sample Changer. Pulse output on I/O line L6.	1, 2, diff. -20.0...25.0...120.0 °C 0...±20.00 0.5..2...999 mV/min, OFF 0...110...999 s, OFF up to 8 ASCII characters ON, OFF all, first, OFF
>statistics	Statistics calculation see at DET, page 14.	

2.9 Mode TIP, parameters

With TIP, methods and other functions can be linked to a titration sequence. Press the key <def> to define the titration sequence.

Display	Meaning	Input range
>sequence method: pause L4 output: L6 output: info prep	Method from the internal method memory or from the card. Waiting time. Can be aborted with <QUIT>. Set a signal when on I/O line L4. Set a signal when on I/O line L6. The sequence will be held and an info appears in the display. Preparation of the dosing unit.	method name 0...999 999 s, inf. active, inactive, pulse, OFF active, inactive, pulse, OFF up to 16 characters internal D0, external D1/D2
>statistics	Statistics calculation see at DET, page 14.	
>preselections req. ident: req. snpl size: meas. mode: meas. input: temperature	Request of identifications after start. Request of sample size after start. Selection of the measured quantity for <meas/hold>. Measuring input for pH and U. Temperature for the compensation of the pH values.	id1, id1 & 2, all, OFF value, unit, all, OFF pH, U, I(pol), U(pol), T, OFF 1, 2, diff. -170.0...500.0 °C

2.10 Calculations

The results are calculated with the formulas given in the method. For the calculation, the raw values, determined in the method (EP's and C variables), are available.

The results and the raw values can be assigned for further calculations, e.g. for statistics or as common variables.

Formula input and assignments, key <def>

Display	Meaning	Input range
>formula RS? RS1=EP1*C01/C00 RS1 text RS1 decimal places RS1 unit:	Input of formula. Enter result number. Enter formula by means of 3 rd functions of keyboard. Text for result output. Number of decimal places for result output. Select result unit. Enter values of calculation variables with <C-fmla>.	1...9 RS1 or up to 8 ASCII characters 0...2...5 %, ppm, g/l, mg/ml, mg/pc,mol/l, mmol/l, g, mg, ml, s, ml/min, no unit or up to 6 ASCII characters
>silc calculations C24= C25= mtch id:	Allocations for silc calculation, see page 93 of the Instructions for Use. Enter values to be stored in the silc memory. Indication which id's must coincide for the statistical evaluation in the silc memory.	RSX, EPX, CXX id1, id1 & 2, all, OFF
>common variables C3? =	Allocations of common variables C3X. Enter the number of the common variable and RSX, EPX, CXX, or MNX.	RSX, MNX, EPX, CXX
>report report:	Selection of report blocks for data output. Depends on the selected mode. If you wish several reports, use ";" as separator.	full, short, mplist, curve, derive, comb, meas crv, temp crv, scalc full, scalc srt, calc, param, calib, ff
>mean MN? =	Allocations for statistics. Enter mean number and result RSX, EPX, or CXX.	RSX, EPX, CXX
>temporary variables C7? =	Allocations of temporary variables for calculations under TIP.	RSX, EPX, CXX

Meaning of the calculation variables

C variables are

- determined by the method (C24...C27, C4X, C5X, C6X, C7X, C8X)
- entered as sample specific data (C00, C21...C23)
- given as fixed constants in the methods (C01...C19)

Variable	Meaning
C00	Sample size, key <smpl data>.
C01. . . C19	Method specific calculation values, such as molecular mass, factors, key <C-fmla>.
C21. . . C23	Sample specific calculation values, such as dilution factors, key <smpl data>.
C24, C25	Variables for storing determination results in the silo memory.
C26, C27	Means from silo calculations.
C30. . . C39	Common variables, e.g. for titer.
C40	Initial measured value of the sample.
C41	End volume.
C42	Determination time.
C43	Volume drift for KFT and SET with conditioning.
C44	Temperature.
C45	Dispensed start volume.
C46	Asymmetry pH (pH calibration).
C47	Electrode slope (pH calibration).
C51. . . C59	Fix EP for DET and MET or fix V for STAT.
C61. . . C69	pK/HNP values for DET and MET or fix times for STAT.
C70. . . C79	Temporary variables for calculations in TIP.
C80	Mean rate for STAT: Overall rate, evaluated over all points of the measuring list.
C81. . . C89	Rates, evaluated in preset time windows for STAT.

Sample data, key <smpl data>

- Sample identifications or sample specific calculation values C21...C23
- Sample size C00

For working with the silo memory, see page 27.

Display	Meaning	Input range
id#1 or C21 id#2 or C22 id#3 or C23 snpl size snpl unit:	Inquiries with silo = OFF (LED "silo" is OFF): } Sample identification. Can be used as } sample specific calculation values. Sample size Unit of sample size	up to 8 ASCII characters -999 999...1...999 999 g, mg, ml, ul, pc, no unit or up to 5 ASCII characters


3 Other functions

This chapter gives you a survey over other functions of the Titrimo.

If you are looking for more information to the different quantities, you will find the display text as a keyword in the index of the Instructions for Use.

3.1 Method memory

Methods can be stored in the internal method memory or on the memory card. These method memories are organized as follows:



Internal method memory, key <user meth>

Display	Meaning	Input range
>recall method method name:	Recall method from the internal method memory to the working memory.	up to 8 characters
>store method method name:	Store method from the working memory to the internal method memory.	up to 8 characters
>delete method method name:	Delete method from the internal method memory.	up to 8 characters

Method memory on the card, key <card>

Display	Meaning	Input range
>recall method method name:	Recall method from the card (current directory) to the working memory.	up to 8 characters
>store method method name:	Store method from the working memory to the current directory on the card.	up to 8 characters
>delete method method name:	Delete method from the current directory on the card.	up to 8 characters
>change directory dir. name:	Change current directory on the card.	up to 10 characters
>create directory dir. name:	Create new directory on the card.	up to 10 characters
>delete directory dir. name:	Delete directory on the card together with its methods.	up to 10 characters
>backup dir. name:	Backup of the internal method memory on the card.	up to 10 characters
>reload dir. name:	Reload of the whole directory from the card to the internal method memory.	up to 10 characters
>format card label format:	Card formatting. Card labeling. Confirm formatting.	up to 8 characters ON, OFF
>change battery date	Date for changing the card battery.	YYYY-MM-DD

3.2 Sample data, keys <smpl data> and <silos>

In the silo memory sample data can be stored on reserve. For working without the silo memory see page 24 of the instructions for use.

Press the key <silos> for working with the silo memory, and the key <smpl data> to input data.

Display	Meaning	Input range
>edit silo lines id#1 or C21 id#2 or C22 id#3 or C23 smpl size smpl unit:	Input for the silo memory: } Sample identification. Can be used as } specific calculation values. Sample size Unit of sample size.	up to 8 ASCII characters -999 999...1...999 999 g, mg, ml, ul, pc, no unit or up to 5 ASCII characters
>delete silo lines delete line n	Delete individual silo lines. Line number of the line to be deleted.	1...99, OFF
>delete all silo lines delete all:	Delete all silo lines. Confirmation.	ON, OFF
cycle lines: save lines:	With "ON", worked off silo lines will be copied to the highest line of the silo memory. Determination results will be stored as C24 or C25 in the silo memory according to the allocations in the methods, see page 23.	ON, OFF ON, OFF

3.3 Configuration, key <configuration>

Display	Meaning	Input range
>peripheral units send to: balance: record:	Settings of peripheral units. Selection of printer. Selection of balance. Selection of record at the analog output.	Epson, Seiko, Citizen, HP, IBM Sartorius, Mettler, Mettler AT, AND, Precisa U, dU/dt, V, dV/dt, U(rel), T
>auxiliaries dialog: date time run number auto start start delay dev. label program	General settings. Selection of dialog language. Current run number for result output. Automatic starts of titrations. Waiting time before start of titration. Device label. Program version.	english, deutsch, français, español, portuguese, italiano, svenska YYYY-MM-DD HH:MM 0..9999 1..9999, OFF 0..999 999 s up to 8 characters read only
>RS232- settings baud rate: data bit: stop bit: parity: handshake: RS control:	Settings of RS232 interface. Baud rate. Data bit. Stop bit. Parity. Handshake. Receiving of commands via RS.	300,600,1200, 2400, 4800, 9600 7, 8 1, 2 none, odd, even HWs, HWf, SWchar, SWline, none ON, OFF
>common variables	Common variables C30...C39.	
>prep. dosing elements power ON prep: report: dos. element: warn. interv. dos. drive: volume DX cycle DX dos. rate. DX fill rate. DX output: len. dos. tub. diam dos. tub. len. asp. tub. diam asp. tub. dos. rate DX fill rate DX	Preparation of titration burets. Warning after power ON. Report. Selection of dosing unit. Warning interval. Selection of type of dosing unit. <i>Parameters for the internal dosing unit and for 685 Dosimat:</i> Volume. Number of cycles. Dosing rate. Filling rate. <i>Parameters for the 700 Dosino:</i> Location, where liquid is expelled. Length of the dosing tube. Diameter of the dosing tube. Length of the aspiration tube. Diameter of the aspiration tube. Dosing rate. Filling rate.	ON, OFF ON, OFF internal D0, external D1/D2 5..9999 min, OFF Dosimat, Dosino 0...3.5...99 999.99 ml 1...2...9 0.01...150 ml/min, max. 0.01...150 ml/min, max. tip, flask 1.0...40.0...999.9 cm 0.1...2.0...9.9 mm 1.0...25.0...999.9 cm 0.1... 2.0...9.9 mm 0.01...150 ml/min, max. 0.01...150 ml/min, max.

4 Index

Backup	8, 26	Recall methods	8
Balance	5	Report	5
Burets	10	- Definitions.....	23
C-Variables	24	Result	
CAL.....	21	- Calculation.....	4
Calculations	4, 23	- Printing	5
Calibration.....	21	- View.....	5
Card	8, 25	Sample data.....	24, 27
Configuration.....	28	Sample size	5
Data input.....	1, 7	Selection	
DET.....	14	- Mode	3
Dialog language	1	- Report.....	5
DOC.....	20	SET.....	15
DOS	19	Silo calculation	23
Dosing unit.....	10	Silo memory	27
Formula.....	4, 23	STAT.....	17
Internal method memory	7, 26	Stop criteria.....	6
KFT	16	Store methods.....	7, 25
Linkage	22	Summary.....	2
Manual dosing.....	11	Survey.....	12
MEAS.....	21	Text input	7
MET	14	TIP	22
Method memory	25	Titration modes	12ff
- Printing	9	-Survey	12
Methods		-Selection.....	3
- Backup	8	Titration procedure	22
- Development	3	Titration stop	6
- Recalling.....	8	Variables	24
- Storing.....	7	Weight.....	5
Modes	12ff		
- Selection.....	3		
- Survey	12		
Operands	4, 24		
Parameters	14ff		
pH Stat.....	17		
Printing.....	5		