

-----  
Method parameters  
-----

Method : AB422 Determination Hg.mth  
Title : AB 422\_2 Determination of Hg  
Remark1 : 10 mL sample + 1 mL electrolyte  
Remark2 : 3% formic acid, 0.06 mol/L KCl, 5 mg/L Fe(III)

Calibration : Standard addition  
Technique : Batch  
Addition : Manual

Sample ID : Sample  
Sample amount (mL): 10.000  
Cell volume (mL): 11.000

Voltammetric parameters  
-----

Mode : DP - Differential Pulse

Highest current range : 10 mA  
Lowest current range : 100 nA

Electrode : SSE/RDE  
Stirrer speed (rpm) : 2000

Initial electr. conditioning : No

No. of additions : 2  
No. of replications : 2

Measure blank : No  
Addition purge time (s) : 0

Initial purge time (s) : 0

Conditioning cycles  
Start potential (V) : -0.300  
End potential (V) : 0.850  
No. of cycles : 5

Hydrodynamic (measurement) : No  
Cleaning potential (V) : 0.750  
Cleaning time (s) : 10.000  
Deposition potential (V) : 0.300  
Deposition time (s) : 90.000

Sweep  
Equilibration time (s) : 5.000  
Start potential (V) : 0.300  
End potential (V) : 0.600  
Voltage step (V) : 0.004  
Voltage step time (s) : 0.400  
Sweep rate (V/s) : 0.010  
Pulse amplitude (V) : 0.050  
Pulse time (s) : 0.040

Cell off after measurement : Yes

Peak evaluation

-----

Regression technique : Linear Regression

Peak evaluation : Height

Minimum peak width (V.steps) : 5

Minimum peak height (A) : 1.000e-010

Reverse peaks : No

Smooth factor : 4

Eliminate spikes : Yes

Substances

-----

Hg : 0.440 V +/- 0.050 V

Standard solution : 1 1.000 mg/L

Addition volume (mL) : 0.050

Mercury : Final result (Hg) =

Conc \* (11 / 10) \* (1e+006 / 1) + 0 - 0

Baseline

-----

Substance	Addition	automatic	start (V)	end (V)	type	scope
Hg	Sample	yes	---	---	linear	wholePeak
	Addition 1	yes	---	---	linear	wholePeak
	Addition 2	yes	---	---	linear	wholePeak

-----

Solutions

-----

No.	Content	Predose (mL)
---	-----	-----

Export options

-----

Export final results as ASCII: no

Export final results as CSV: no

Export final results as XML: no

Export determination to AutoDB: no