



## 946 Portable VA Analyzer

### General information

---

Software version: 1.0

Instrument: -

Sensor: scTRACE Gold

User name:

Report: No

Report elements: -

## Method

### General

Method name: AB433 Plating Ag Film.detp

Remarks: 15 mL H<sub>2</sub>O + 1.5 mL plating electrolyte + 1.5 mL c(Ag) = 20 mg/L

Plating electrolyte: c(SSA) = 0.4 mol/L, c(NaOH) = 1 mol/L

### Determination

Sample volume (mL): 18.0

Total cell volume (mL): 18.0

Stirring time (s): 10.0

Stirring rate (1/min): 2000

Measure blank: No

No. of blanks: 0

Blank value correction: No

No. of replications: 0

No. of additions: 0

### Voltammetric

Measuring mode: Linear sweep

Current measuring range: Auto

#### Cyclovoltammetric pretreatment

Start potential (V): -0.5

Vertex potential (V): 0.4

Potential step (V): 0.01

Sweep rate (V/s): 0.1

No. of cycles: 3

#### Potentiostatic pretreatment

Potential 1 (V): -1.0

Waiting time 1 (s): 60.0

Potential 2 (V): 0.0

Waiting time 2 (s): 0.0

Equilibration time (s): 5.0

#### Sweep

Start potential (V): -0.4

End potential (V): -0.1



## 946 Portable VA Analyzer

Potential step (V): 0.006

Sweep rate (V/s): 0.1



# 946 Portable VA Analyzer

## Evaluation

### Data processing

Smoothing: 1

Calibration method: Standard addition

### Peak evaluation

	-
Characteristic potential (V)	0.0
Tolerance (V)	0.05
Min. width (V)	0.05
Max. width (V)	0.5
Min. measured quantity (µA)	0.001
Baseline type	Linear
Base point automatic	Yes
Start base point (V)	0.0
End base point (V)	0.0

### Standard solutions

	-	Volume (mL)
Standard 1	1.0 mg/L	0.1
Standard 2	-	-
Standard 3	-	-
Standard 4	-	-