



License ID 1955606
Client ID 001-nbk-02849
User Metrohm

Program version viva 3.0 - 69
2024-03-27 09:27:18 UTC+1

Method parameters

Method AB445 - Maintenance - Prepare dosing units
Method saving date 2024-03-27 09:19:43 UTC+1
Method version 1
Method group Main group
Method status original
Method saved by (full name) Metrohm International Headquarters
Method saved by (short name) Metrohm

START

Main track

General

Workplace view

Current view on

Track view for live window

Live display 1 Main track

Live display 2 Main track

Electrode check on

Application note

- Specify the number of preparation cycles for each Dosing Unit in the «Run» window in the respective field «Prep Cycles - ...». For one cycle enter 1, for two cycles 2, and so on. For Dosing Units which should not be prepared enter 0. - The draining pump will empty the measuring vessel during preparation of the dosing units.

Sample data variables

Name	Type	Assignment	Fixed value	Comment	Monitoring
ID1	Text	ID1	Preparation	Sample identification 1	off
Prep Cycles -	Number	ID2		Sample identification 2	off
Prep Cycles -	Number	ID3		Sample identification 3	off
Prep Cycles - DU	Number	ID4			off

Name **ID1**
Type Text
Assignment on. ID1
Fixed value on. Preparation Dosing Unit(s)
Check at start off
Comment Sample identification 1

Name **Prep Cycles - Electrolyte neutral**
Type Number
Assignment on. ID2
Fixed value off.
Check at start on
Comment Sample identification 2




License ID 1955606
Client ID 001-nbk-02849
User Metrohm

Program version viva 3.0 - 69
2024-03-27 09:27:18 UTC+1

Variable monitoring off
Lower limit
Upper limit
Message
Display message on
Record message on
Message by e-mail off
 E-mail template
 Subject Message from viva - method 'New method 1' - command 'Main track'
Acoustic signal off
Action off
Cancel determination on
Cancel determination and series off

Name **Prep Cycles - Electrolyte alkaline**
Type Number
Assignment on. ID3
Fixed value off.
Check at start on
Comment Sample identification 3
Variable monitoring off
Lower limit
Upper limit
Message
Display message on
Record message on
Message by e-mail off
 E-mail template
 Subject Message from viva - method 'New method 1' - command 'Main track'
Acoustic signal off
Action off
Cancel determination on
Cancel determination and series off

Name **Prep Cycles - DU WashStation**
Type Number
Assignment on. ID4
Fixed value off.
Check at start on
Comment
Variable monitoring off
Lower limit
Upper limit
Message



Preparation of dosing units «Electrolyte neutral», «Electrolyte alkaline», and «WashStation»

```
Return immediately . . . . . off
Delete old data . . . . . off
```

Call text	Track name	Sample type	Condition
Call waste position	Waste Position	off Sample	off



License ID 1955606
Client ID 001-nbk-02849
User Metrohm

Program version viva 3.0 - 69
2024-03-27 09:27:18 UTC+1

Call text	Track name	Sample type	Condition
Drain measuring vessel	Drain	off Sample	off

LOOP

Cycles DU neutral

Stop criteria

Max. run number on
Max. run number = 'SD.Prepare Cycles - Electrolyte neutral'
Maximum run time off
Signal assessment for DT off
Condition off

CALL

Call Prep DU neutral

Call text	Track name	Sample type	Condition
Prep DU neutral	Prep DU neutral	off Sample	off

TRACK

DU Electrolyte alkaline

Return immediately off
Delete old data off

CALL

Waste alkaline

Call text	Track name	Sample type	Condition
Call waste position	Waste Position	off Sample	off
Drain measuring vessel	Drain	off Sample	off

LOOP

Cycles DU alkaline

Stop criteria

Max. run number on
Max. run number = 'SD.Prepare Cycles - Electrolyte alkaline'
Maximum run time off
Signal assessment for DT off
Condition off

CALL

Call Prep DU alkaline

Call text	Track name	Sample type	Condition
Prep DU alkaline	Prep DU alkaline	off Sample	off

TRACK

DU neutral & alkaline

Return immediately off
Delete old data off

Waste neutral & alkaline



License ID 1955606
Client ID 001-nbk-02849
User Metrohm

Program version viva 3.0 - 69
2024-03-27 09:27:18 UTC+1

CALL

Call text	Track name	Sample type		Condition
Call waste position	Waste Position	off	Sample	off
Drain measuring vessel	Drain	off	Sample	off

LOOP

1_Cycles DU alkaline

Stop criteria

Max. run number on
Max. run number = 'SD.Prepare Cycles - Electrolyte alkaline'
Maximum run time off
Signal assessment for DT off
Condition off

CALL

1_Call Prep DU alkaline

Call text	Track name	Sample type		Condition
Prep DU alkaline	Prep DU alkaline	off	Sample	off

WAIT

WAIT

Wait

Stop track and waiting for [Continue] off
Stop all tracks and waiting for [Continue] off
Waiting time on
Time 0
Unit s

Message

Record message off
Message by e-mail off
Acoustic signal off

LOOP

2_Cycles DU neutral

Stop criteria

Max. run number on
Max. run number = 'SD.Prepare Cycles - Electrolyte neutral'
Maximum run time off
Signal assessment for DT off
Condition off

CALL

2_Call Prep DU neutral

Call text	Track name	Sample type		Condition
Prep DU neutral	Prep DU neutral	off	Sample	off

TRACK

Waste Position

Return immediately off



License ID 1955606
Client ID 001-nbk-02849
User Metrohm

Program version viva 3.0 - 69
2024-03-27 09:27:18 UTC+1

Delete old data off

SWING Swing Waste Position

Device

Device name 858 Sample Processor

Device type 858.0020 Professional Sample Processor

Target

Tower 1

Swing External position

Number 2

Parameters

Swing rate 55 %/s

LIFT Lift Waste Position

Device

Device name 858 Sample Processor

Device type 858.0020 Professional Sample Processor

Target

Tower 1

Lift position Work position mm

Parameters

Lift rate 25 mm/s

TRACK Drain

Return immediately off

Delete old data off

PUMP PUMP Drain

Device

Device name 858 Sample Processor

Device type 858.0020 Professional Sample Processor

Pumps

Tower 1

Pump(s) 1

Action

Switch on on

Switch off off

Duration off

Time 15 s

TRACK DU WashStation

Return immediately off

Delete old data off

LOOP Cycles DU WashStation

Stop criteria

Max. run number on

Max. run number = 'SD.Prep Cycles - DU WashStation'



License ID 1955606
Client ID 001-nbk-02849
User Metrohm

Program version viva 3.0 - 69
2024-03-27 09:27:18 UTC+1

Maximum run time off
Signal assessment for DT off
Condition off

PREP Prep WashStation

Dosing device
Dosing unit WashStation

TRACK Prep DU neutral

Return immediately off
Delete old data off

PREP Prep Electrolyte neutral

Dosing device
Dosing unit Electrolyte neutral

TRACK Prep DU alkaline

Return immediately off
Delete old data off

PREP Prep Electrolyte alkaline

Dosing device
Dosing unit Electrolyte alkaline

TRACK Shut off

Return immediately off
Delete old data off

PUMP PUMPS OFF

Device
Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor
Pumps
Tower 1
Pump(s) 1+2
Action
Switch on off
Switch off on
Duration off
Time 8.0 s

LIFT Lift Home

Device
Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor
Target
Tower 1
Lift position Home position mm
Parameters



License ID 1955606
Client ID 001-nbk-02849
User Metrohm

Program version viva 3.0 - 69

2024-03-27 09:27:18 UTC+1

Lift rate 25 mm/s

EXIT **Exit track**

CALL **Exit shut off**

Call text	Track name	Sample type		Condition
Exit shut off	Shut off	off	Sample	off

ERROR **Error track**

CALL **Error shut off**

Call text	Track name	Sample type		Condition
Error shut off	Shut off	off	Sample	off