Application Note AN-T-173

pHe value of denatured ethanol fuel

Fast and accurate measurement according to ASTM D6423

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SAMPLE AND SAMPLE PREPARATION

The method is demonstrated for denatured ethanol fuel.

Temperature has a strong influence on the obtained

results. Therefore, sample temperatures are first adjusted to be at the same value in order to be able to compare results.

EXPERIMENTAL

This application is performed on a 913 pH Meter equipped with an EtOH Trode, a temperature sensor, and an external stirring plate. The EtOH Trode is conditioned and calibrated prior to use.

A defined amount of sample is poured into a 100 mL beaker and placed on an external stirring plate. The EtOH Trode and the temperature sensor are immersed, and the measurement is started immediately. The value after 30 seconds is considered to be the acid strength of the sample.

Figure 1. 913 pH Meter equipped with a pH electrode. Example setup for the determination of the pHe value.

RESULTS

The analysis demonstrates acceptable and reproducible pHe values. For the tested denatured ethanol fuel, a pHe value of 8.08 (n = 10, SD(rel) =

0.45%) is obtained. An example measurement curve is displayed in **Figure 2**.


Figure 2. Example curve of a pHe measurement in denatured ethanol fuel.

CONCLUSION

With a 913 pH Meter it is possible to assess quick and easily the pHe of an ethanol fuel according to **ASTM D6423**. The used EtOH Trode is specifically mentioned in the standard, and is therefore fully compliant.

By using a 913 pH/DO meter or a 914 pH/Conductometer, it is also possible to assess the dissolved oxygen concentration or the conductivity in parallel, if necessary.

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