Application Note AN-V-231

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Simultaneous determination on Metrohm DropSens Screen-Printed Carbon Electrodes

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SAMPLE

Drinking water, mineral water, sea water

EXPERIMENTAL

Prior to the first determination, the ex-situ mercury film is deposited in a separate step on the screenprinted electrode. The water sample and the supporting electrolyte are pipetted into the measuring vessel. The simultaneous determination of cadmium and lead is carried out with the 884 Professional VA or with the 946 Portable VA Analyzer using the parameters specified in **Table 1**. The concentration of both elements is determined by two additions of a cadmium and lead standard addition solution.

Figure 1. 946 Portable VA Analyzer (SPE)

Figure 2. 884 Professional VA, semiautomated system

Table 1. Parameters

Parameter	Setting
Mode	SQW – Square wave
Deposition potential	-1.3 V
Deposition time	60 s
Start potential	-1.0 V
End potential	-0.4 V
Peak potential Cd	-0.72 V
Peak potential Pb	-0.52 V

ELECTRODES

- Screen-printed carbon electrode (Metrohm DropSens 11L)

RESULTS

With a 30 s deposition time, the limit of detection is around 1 $\mu\text{g/L}$, and the linear

working range is up to 20 $\mu\text{g/L}$ in the measuring solution for both elements.


Figure 3. Determination in mineral water spiked with 2 $\mu\text{g/L}$ cadmium and lead

Table 2. Result

Sample	Cd (μg/L)	Pb (μg/L)
Mineral water spiked with 2 μg/L Cd and Pb	2.04	1.81

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
Screen-Printed Carbon Electrode (Aux.:C;Ref.:Ag/AgCl)

Screen-Printed Carbon Electrode
(Aux.:C;Ref.:Ag/AgCl)