Application Note AN-V-232

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## Simultaneous determination in low $\mu\text{g/L}$ range on the 11L SPEs modified with a Bi film

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
Drinking water, mineral water

## EXPERIMENTAL

Prior to the first determination, an ex-situ bismuth film is deposited from a Bi solution. In the next step, the electrodes are cleaned with ultrapure water and the bismuth solution is removed. The water sample is placed into the measuring vessel. Ammonia / ammonium chloride buffer along with the complexing agent (dimethylglyoxime) are added, and the simultaneous determination of nickel and cobalt is carried out using the parameters specified in **Table 1**. The concentration is determined by two additions of a nickel and cobalt standard addition solution.

**Figure 1.** 946 Portable VA Analyzer (SPE)

**Figure 2.** 884 Professional VA semiautomated

**Table 1.** Parameters

Parameter	Setting
Mode	DP – Differential Pulse
Deposition potential	-0.9 V
Deposition time	30 s
Start potential	-0.9 V
End potential	-1.3 V
Peak potential Ni	-1.05 V
Peak potential Co	-1.175 V


## ELECTRODES

- Screen-printed carbon electrode (Metrohm DropSens 11L)

## RESULTS

With a 30 s deposition time, this method is suitable for the determination of both nickel and

cobalt in water samples in concentrations from  $\beta$  (Ni) = 0.4–5  $\mu$ g/L and  $\beta$  (Co) = 0.2–8  $\mu$ g/L.

**Figure 3.** Determination of nickel and cobalt in tap water (946 Portable VA Analyzer; 30 s deposition time)

**Table 2.** Result

Sample	Ni (µg/L)	Co (µg/L)
Tap water	1.15	<LOD

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
### Screen-Printed Carbon Electrode (Aux.:C;Ref.:Ag/AgCl)

Screen-Printed Carbon Electrode  
(Aux.:C;Ref.:Ag/AgCl)