



Application Note AN-CIC-034

水中可吸附有机素(AOX)的快速分析方法——燃-子色法

根据 DIN 38409-59 准量 AOCl、AOBr、AOI 和 AOF

可吸附有机素(AOX),即可用活性炭吸附的有机素的和,是一个非常的参数,其中很多有机素及其降解物都会境和人健康来危害[1-4]。在水、代染物溯源追踪以及水理工中AOX去除效率价的各程中,我都需要AOX行。一直以来,AOX的定方法一直是使用活性炭水行吸附,然后将活性炭燃后使用微滴定法行(DIN EN ISO 9562 或 EPA 1650) [1,2]。按照相定,AOX包括可吸

附有机((AOCl)、可吸附有机(AOBr)和可吸附有机(AOI),而不包含可吸附有机(AOF)。新布的准 DIN 38409-59 介了通燃-子色法 (CIC)定 AOCl、AOBr、AOI、AOX(Cl)以及 AOF(目前收到广泛注的全和多基物的参数)的方法。本用告介了如何使用瑞士万通燃炉-子色法参照DIN 38409-59 准 AOX 和 AOF 行分析。

EXPERIMENTAL

This application is focused on the experimental approach of AOX and AOF analysis. More detailed information can be found in related Metrohm literature ([WP-078](#), [WP-081](#), [AN-CIC-033](#)). The complete validation dataset of DIN 38409-59 is available on the [Water Chemistry Society](#) webpage.

The overall sample preparation procedure, i.e., preconcentration and adsorption of organically bound halogens, resembles that of DIN EN ISO 9562, as adsorption on activated carbon is a key point for both methods (**Figure 1**). While for AOF

it is crucial that the samples are neutral to avoid adsorption of inorganic fluorine to the activated carbon, sample acidification is mandatory for the other organically bound halogens, similar to DIN EN ISO 9562. For CIC-AOX_(Cl) determination (i.e., AOCI, AOB_r, and AOI), samples need to be acidified with nitric acid to pH <2 prior to preconcentration (**Table 1**). AOF determination is now within the scope of the new DIN 38409-59, however the sample preparation for such samples requires neutralization. This is done by adding sodium nitrate to the samples (**Table 1**).




Figure 1. Schematic of the procedure for AOX and AOF analysis (WP-081). The first step is adsorption performed with the APU sim (Analytik Jena) for semi-automated and standardized adsorption of up to six samples in parallel. After the second step of sample transfer into the combustion boats, the sample is automatically combusted (step 3, combustion module from Analytik Jena consisting of a combustion oven with Auto Boat Drive (ABD) and an autosampler (MMS 5000)). In the fourth step, the volatilized halogens are transported to the absorber solution via gas stream (920 Absorber Module). The last step (5) is the automatic analysis of AOB_r, AOCI, and AOI, or of AOF with the IC (930 Compact IC Flex) including data evaluation. The complete CIC process is fully automated and controlled by MagIC Net software from Metrohm.

The adsorption of the organically bound halogens is handled in a semi-automated manner using the APU sim system from Analytik Jena (**Figure 1**). Two columns filled with activated carbon (at least 50 mg in each column) are connected in series and 100 mL of sample is

passed through. The organically bound halogens adsorb to the activated carbon (using dedicated disposable columns for AOF and AOX determination, **Table 1**), while inorganic halogens are rinsed off (**Figure 1**).

Table 1. Parameters for AOF and AOX sample preparation.

	AOF	AOCI, AOBr, AOI
pH	Neutralized	Acidified to pH <2 with nitric acid
Buffer	0.5 mL 2 mol/L sodium nitrate	0.5 mL 2 mol/L sodium nitrate, acidified with nitric acid
Sample volume		100 mL
		25 mL
Rinsing solution	0.01 mol/L sodium nitrate	0.01 mol/L sodium nitrate, acidified with nitric acid
Absorption columns	Two activated carbon tubes (disposable, from Analytik Jena)	
	402-880.616	402-880.610
Flow rate APU sim		3 mL/min

After the semi-automated sample preparation is finished, the complete content of the two adsorption columns is transferred into one or two separate ceramic boats for CIC analysis. Combustion occurs at temperatures above 950 ° C in the presence of argon and oxygen (**Figure 1**). For pyrohydrolytic combustion, a water stream is essential as it converts the halogens into their hydrogenous forms. Chlorine, bromine, iodine, and fluorine are volatilized in the combustion step, transported into the absorber solution (ultrapure water) with an argon/oxygen gas stream, and transferred into the liquid phase (**Figure 1**). Dosinos guarantee precise automated liquid handling e.g., the

transfer of the aqueous sample into the IC for analysis, or the water stream essential for pyrohydrolytic combustion.

The ion chromatographic separation is achieved on a Metrosep A Supp 5 - 250/4.0 column in combination with the A Supp 5 Guard/4.0. AOF (as F) elutes in under 7 minutes while AOX (i.e., Br, Cl, and I) elutes in less than 25 minutes (**Figure 2**). Automatic system calibration with MiPT (Metrohm intelligent Partial-Loop Injection Technique) is performed using inorganic anion standards for fluoride, chloride, bromide, and iodide (1 g/L standard solutions, TraceCert® from Sigma-Aldrich).




Figure 2. Chromatogram overlay of the blank and a wastewater sample for the determination of AOCI, AOBr, and AOI measured from absorption column #1. To calculate the mass concentration of the individual AOX fractions, blank correction was performed according to Equation 1. No halogens were adsorbed on column #2, revealing the retention efficiency for AOX on column #1.

Performance checks of AOF and AOX determinations and the standard series for LOD determination (**Table 2**) are accomplished using organic reference standard solutions with varying concentrations (4-fluorobenzoic acid, 4-chlorobenzoic acid, 4-bromobenzoic acid, and 4-iodobenzoic acid), treated in the same way as

the samples.

As the procedure for the determination of AOX and AOF is rather complex, dedicated sample boats and charcoal (i.e., fluoride-free materials for AOF, **Table 1**) and blank measurements are essential to guarantee a low background and an appropriate blank correction (**Equation 1**).

$$c(X_{ads}) = \left(c(X^-)_{IC} * \frac{V_{Abs}}{V_{Smpl}} \right) - \left(c(X^-)_{IC} * \frac{V_{AbsBW}}{V_{SmplBW}} \right)$$

Equation 1.

$c(X_{ads})$	Mass concentration of individual adsorbable organically bound halogens (with X = Cl, Br, I, and F) in µg/L
$c(X^-)$	Halogen concentration in the sample's absorption solution in µg/L (with X = Cl, Br, I, and F) in µg/L
V_{Abs}	Final volume of the absorption solution in L

V_{Smpl}	Volume of the sample that was used for adsorption; always 0.1 L
$c(\text{X}^-)_{\text{BW}}$	Halogen concentration in the absorption solution of the blank in $\mu\text{g/L}$
$V_{\text{AbsB}}^{\text{W}}$	Final volume of the absorption solution of the blank in L
$V_{\text{SmplB}}^{\text{W}}$	Volume of the blank solution that was used for adsorption; always 0.1 L

RESULTS

Individual concentrations for AOCl, AOBr, and AOI, as well as for AOF from neutralized samples are calculated according to **Equation 1**. A sum parameter for AOX ($\text{CIC-AOX}_{(\text{Cl})}$) is calculated

using **Equation 2**. However, due to the novelty of this validated approach, $\text{CIC-AOX}_{(\text{Cl})}$ has not yet replaced AOX in water or wastewater regulations.

$$c(\text{CIC-AOX}_{(\text{Cl})}) = c(\text{AOCl}) + c(\text{AOBr}) \cdot 0.4437 + c(\text{AOI}) \cdot 0.2794$$

Equation 2.

$c(\text{CIC-AOX}_{(\text{Cl})})$	Sum concentration of adsorbable organically bound halogens in $\mu\text{g/L}$ as mass concentration based on chloride
-----------------------------------	---

Dedicated materials and the sensitive analysis of the halogens with suppressed conductivity detection results in low blank values. Blank values were only measurable for fluoride and chloride (**Table 2**). The requirements for DIN 38409-59 are fulfilled—in fact, the overall

procedure here is even more sensitive. During the DIN validation process, several water samples were analyzed from different laboratories using similar setups (validation report: wasserchemische-gesellschaft.de).

Table 2. Blank, LOD (limit of detection for full AOF/AOX procedure), and DIN scope for the determination of adsorbable organically bound halogens. LODs are determined according to DIN 32645. For AOB_r and AOI, the LODs are determined using the calibration curve as no blank values were found. For AOF and AOCl, the blank method was applied (DIN 32645).

	Blank (µg/L)	LOD (DIN 32645) (µg/L)	Scope of DIN application (µg/L)
AOF	1.1	0.38	2
AOCl	2.6	1.36	10
AOB _r	0	0.24	1
AOI	0	0.47	1

Using IC, it is now possible to not only determine the sum parameter CIC-AOX_(Cl), but also to measure the fractions contributing to the AOX

contents (**Figure 2**, WP-081) and to assess AOF (AN-CIC-033, WP-078).

CONCLUSION

Overall, the entire validated procedure profits from its easy, straightforward, and standardized handling, the precise determination of the analytes, automatic calculation of results, and a low maintenance, single-manufacturer setup.

A significant advantage of DIN 38409-59 is that it allows the determination of adsorbable organically bound halogens as individual sum parameters (i.e., AOCl, AOB_r, and AOI) and also provides a fast approach to assess total PFASs content using the validated approach for AOF. Automation (e.g., automated eluent production, MiPT, intelligent and logical MagIC Net features)

improves repeatability, accuracy, and reliability of the results, saves valuable laboratory time for the liquid handling, standard, and eluent preparation, and allows 24/7 analysis – from which every laboratory, either research, routine, or governmental lab – can profit.

The world of organically bound halogens is so varied that these sum parameters enable insights about hot spots, transport pathways, but also particularly vulnerable regions in a very simple way, while complex targeted analysis, if at all, can resolve individual organically bound halogens for deeper investigations afterwards.

REFERENCES

1. Xu, R.; Xie, Y.; Tian, J.; et al. Adsorbable Organic Halogens in Contaminated Water Environment: A Review of Sources and Removal Technologies. *J Clean Prod* **2021**, 283.
2. Müller, G. Sense or No-Sense of the Sum Parameter for Water Soluble "Adsorbable Organic Halogens" (AOX) and "Absorbed Organic Halogens" (AOX-S18) for the Assessment of Organohalogens in Sludges and Sediments. *Chemosphere* **2003**, 52 (2), 371–379.
3. Dann, A. B.; Hontela, A. Triclosan: Environmental Exposure, Toxicity and Mechanisms of Action. *J Appl Toxicol* **2011**, 31 (4), 285–311.
4. Xie, Y.; Chen, L.; Liu, R. AOX Contamination Status and Genotoxicity of AOX-Bearing Pharmaceutical Wastewater. *J Environ Sci* **2017**, 52, 170–177.

Internal reference: AW IC CH6-1438-042021

CONTACT

瑞士万通中国
北京市海淀区上地路1号院
1号楼7702
100085 北京

marketing@metrohm.com.cn

配置



Combustion IC Manuell -

瑞士万通 Combustion IC Manuell - 石英款套餐的作用是,通英燃消解(水解)和随后行的例子色定(Combustion IC) 各形式可燃品中的素和硫行分析。其包含了所有必要的件,其中包括 Trace Elemental Instruments 的 Combustion Oven (TEI) (2.0136.0600)、石英燃管 (6.07311.100)、920 吸收器模、930 Compact IC Flex Oven/SeS/PP/Deg 和 MagIC Net 件。瑞士万通 Combustion IC 套餐需要可以充以下 Autosampler 中的一: Solid Autosampler CIC (TEI)、Liquid Autosampler CIC (TEI) 或 GLS Sampler CIC (TEI)。



Metrosep A Supp 5 - 250/4.0

万通的高效分柱以其高基数用于度大的分任。即便是的分也可用 Metrosep A Supp 5 - 250/4.0 松解决,且其解决方法可重。柱的高容量可在不行品前理的情况下定 150 mg/L 化物中的 1 μ g/L 酸。此柱的使用范不止用来定准子。在半体工中或厂的炉水程中控制高准的度要求,Metrosep A Supp 5 - 250/4.0 柱是最佳。



Metrosep A Supp 5 Guard/4.0

Metrosep A Supp 5 Guard/4.0 IC 子柱 Metrosep A Supp 5 和 7 行可靠保,使其不受品或淋洗液的染。它含有与 Metrosep A Supp 5 相同的分材料,也同由 PEEK(聚)制成,并直接以几乎无死点容的方式到相分柱上(«On Column Guard System»)。保柱将延分析用柱的使用寿命,且不会影其色分效率。建使用 A Supp 5 Guard/4.0,其价格惠且操作。



MiPT

用于安装 Dosino 以行局部循的附件。



858 Professional Sample Processor – Pump

858 Professional Sample Processor – Pump 可理体在 500 μ L 至 500 mL 之的品。行品移,既可以使用内置的双向双通道蠕、也可通 800 Dosino 来行。