



Application Note AN-T-181

# 水中的

## Reliable and inexpensive determination by potentiometric titration

Lithium is a soft metal which is used for many applications, such as production of high-temperature lubricants or heat-resistant glass. Furthermore, lithium is used in large quantities for battery production. It is obtained from brines and high-grade lithium ores. Depending on the lithium concentration, extraction may or may not be economically viable. This Application Note demonstrates a method to determine the lithium concentration in brines by

potentiometric titration. Lithium and fluoride precipitate in ethanol as insoluble lithium fluoride. Using ammonium fluoride as titrant and a fluoride ion selective electrode (ISE), determination of lithium via potentiometric titration is possible.

This method is more reliable, faster, and less expensive than the determination of lithium in brine by other more sophisticated techniques such as atomic absorption spectroscopy (AAS).

Two samples with ~20% (m/m) and ~40% (m/m) lithium chloride are analyzed. No further sample


preparation is required.

## EXPERIMENTAL

The analysis is carried out with an automated system consisting of *tiamo*<sup>TM</sup> in combination with a 905 Titrand. A fluoride ion selective electrode (ISE) in combination with a Long Life ISE reference electrode is used for the indication of the titration.

After transferring an appropriate amount of sample into the titration vessel, enough ethanol is added to cover the diaphragm of the electrode. The titration is carried out using ammonium fluoride until after the first equivalence point is reached.

Calcium will interfere with the analysis and has to be analyzed separately.




**Figure 1.** 905 Titrand with tiamo. Example setup for the analysis of lithium in brine.

For both samples this method gave reliable and

reproducible results.

**Table 1.** Results of the lithium determination in brines.

n = 3	Assay of LiCl brine 1 in %	Assay of LiCl in brine 2 in %
Mean	41.3	14.85
SD(abs)	0.7	0.11
SD(rel)	1.6	0.7



**Figure 2.** Example curve of the lithium determination in brines.

## CONCLUSION

Using ethanol as solvent, ammonium fluoride as titrant, and the fluoride ISE for indication, the determination of lithium in brine can be performed reliably and cost-efficiently.

Therefore, this method provides an inexpensive and simple approach to determine if an extraction of lithium from brine is economically feasible or not.

## CONTACT

Metrohm AG  
Ionenstrasse  
9100 Herisau

[info@metrohm.com](mailto:info@metrohm.com)

## CONFIGURATION



### 907 Titrando

Dosino •

- 800 Dosino
- (DET)(MET)(SET), pH-STAT (STAT)•(KFT)
- "ITrode"
- (MEAS CONC)
- ,LQH
- MSB
- USB
- OMNIS-Software *tiamo*- Touch Control
- , GMP/GLP FDA , 21 CFR 11



F-ISE Pt1000

, Pt1000

ISE :

- F- (10<sup>-6</sup> mol/L )

-

-

:c(KCl) = 3 mol/L