

Combined dCa ISE



6.00502.300

Sensor leaflet

8.0109.8012EN / 2020-11-25



Metrohm AG
Ionenstrasse
CH-9100 Herisau
Switzerland
Phone +41 71 353 85 85
Fax +41 71 353 89 01
info@metrohm.com
www.metrohm.com

Combined dCa ISE

6.00502.300

Sensor leaflet

8.0109.8012EN /
2020-11-25

Technical Communication
Metrohm AG
CH-9100 Herisau
techcom@metrohm.com

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

References

1	Overview	1
1.1	Combined dCa ISE – Product description	1
1.2	Combined dCa ISE – Overview	1
2	Functional description	2
2.1	Combined dCa ISE – Functional description	2
3	Transport and storage	3
3.1	Electrode – Checking the delivery	3
3.2	Electrode – Storing the packaging	3
3.3	Unpacking and checking the electrode	3
3.4	Storing the Combined dCa ISE	4
4	Installation	6
4.1	Preparing the Combined dCa ISE	6
4.2	Mounting the electrode	7
5	Operation and control	9
5.1	Combined dCa ISE – Measurement procedures	9
6	Maintenance	11
6.1	Combined dCa ISE – Changing/refilling the electrolyte ...	11
6.2	Cleaning the Combined dCa ISE	11
6.3	Checking the Combined dCa ISE	12
7	Combined dCa ISE – Troubleshooting	13
8	Electrode – Disposal	14
9	Technical specifications	15
9.1	Ambient conditions	15
9.2	Combined dCa ISE – Dimensions	15
9.3	Combined dCa ISE – Housing	15
9.4	Combined dCa ISE – Connectors specifications	16
9.5	Combined dCa ISE – Display specifications	16
9.6	Combined dCa ISE – Measurement specifications	16

10 Additional information

17

1 Overview

1.1 Combined dCa ISE – Product description

The Combined dCa ISE is a calcium-selective, combined polymer membrane electrode with shockproof membrane for titration, direct measurement and standard addition.

1.2 Combined dCa ISE – Overview




Figure 1 Combined dCa ISE

1	Protective cap	2	Electrode head
3	Filler opening	4	SGJ sleeve SGJ 14/15, movable
5	Sensor surface		

2 Functional description

2.1 Combined dCa ISE – Functional description

An ion-selective electrode only responds to a specific ion in the solution and, ideally, it does not show any change in potential with other ions present.

The measuring ions of the sample solution reach the membrane surface of the ion-selective electrode, after a certain amount of time an equilibrium is established. An electrochemical potential is created between the measuring solution and the membrane.

3 Transport and storage

3.1 Electrode – Checking the delivery

Immediately upon arrival of the merchandise, check the shipment to ensure absence of damage.

3.2 Electrode – Storing the packaging

The product is supplied in extremely protective special packaging. Keep this packaging, as only this ensures safe transportation of the product.

3.3 Unpacking and checking the electrode

1 Unpacking the electrode

Remove the electrode with storage vessel from the packaging.

2 Removing the storage vessel




Figure 2 Releasing the electrode from the storage vessel

- Hold the electrode and storage vessel firmly in your hand so that the electrode cannot slip away.
- Position the tool between the storage vessel and SGJ sleeve.
- **Carefully** push the tool to the side to release the electrode.
Do not tip the tool forwards!



3 Checking the electrode for proper function

- ## NOTICE

3.4 Storing the Combined dCa ISE

1 For short periods

- ## NOTICE

2 For longer periods

- 4 ■ ■ ■ ■ ■ ■ ■ ■



We recommend maintaining some residual moisture between the inner pipe and the three protection flaps to keep the electrode ready for use.

- Push the outer ring of the cable connection over the electrode head.
Ensure that the guide lugs in the electrode head are in the grooves of the cable connection.
- Push the cable connection onto the electrode head until it snaps in place.



NOTICE

To remove the cable, first release the outer ring and then carefully pull the cable connection from the electrode head.

When doing so, be sure not to pull on the cable itself but the cable connector instead.

4.2 Mounting the electrode



The electrode must sit securely in the titration head.



NOTICE

For automatic procedures, ensure that the cables have enough room to move.

During the titration or standard addition, it is important that the solution is mixed well. The stirring rate should be high enough to form a small vor-

5 Operation and control

5.1 Combined dCa ISE – Measurement procedures

Titration

Ion-selective electrodes are well-suited for potentiometric titrations. The resulting titration curves are usually S-shaped and can be evaluated well with automatic titrators.

For application advice on working with ion-selective electrodes, please go to www.metrohm.com.

Direct measurement with calibration

The ion activity in the sample is interpolated by means of a calibration curve. The calibration curve is established with standard solutions. The expected ion activity in the sample should lie in the mid-concentration range of the standard solutions.

Since the concentration of an ion is usually of interest (rather than its activity), measurements are performed at a fixed ionic strength. The ionic strength is measured in an ISA (Ionic Strength Adjuster) solution or a TISAB (Total Ionic Strength Adjustment Buffer) solution. ISATISAB solutions have a high ionic strength so that the various contributions of the measuring ion to the ion strength can be ignored.

For calcium, it is preferable to use a 1 mol/L potassium chloride solution.



NOTICE

Measure samples and calibration standards under identical measuring conditions. The temperature of the standard solutions and the sample solutions should, as far as possible, be the same during measurements. Furthermore, the temperature should fluctuate as little as possible during the measurement.

In order to guarantee reliable results, periodically carry out a control measurement with a calibration standard (e.g. daily). A new calibration curve should be established if the deviation is deemed unacceptable.

Standard addition / standard subtraction

In the standard addition method, a defined quantity of the ion to be determined is added to a known volume of the sample (in several increments). Normally, ISATISAB solutions are used in this process. The unknown concentration can be calculated from the resulting potential differences between the sample and the sample with added standard solution. This calculation is performed automatically by modern ion meters.

The volume of the added standard solutions should not exceed 25% of the sample volume, and the concentration of the standard solutions should be as high as possible (in order to be able to rule out errors due to dilution). The potential differences between the increments should be roughly constant and amount to at least 10 mV. Temperature differences between the standard solution and the sample solution should be avoided. In addition, at least three additions should be carried out.

In the standard subtraction method, a solution that eliminates the ion to be determined is added (complexation or precipitation). Apart from that, the same conditions apply as for standard addition. However, this method is rarely used.

6 Maintenance

6.1 Combined dCa ISE – Changing/refilling the electrolyte

- 1 Open the filler opening (1-3) by turning it.
- 2 Use a plastic pipette to empty the electrode.
- 3 Rinse the inside of the electrode with the new electrolyte.
- 4 Fill the electrode with electrolyte up to the filler opening.
- 5 Close the filler opening (1-3).

6.2 Cleaning the Combined dCa ISE

- 1 Rinse the electrode with distilled water after each measurement or titration.



NOTICE

The surface must be kept clean at all times before the measurement.



NOTICE

Do not treat the electrode in ultrasonic baths! The electrode could become damaged.

6.3 Checking the Combined dCa ISE

- 1 Measure $c(\text{Ca}^{2+}) = 10^{-4}$ mol/L standard solution and write down the potential.
- 2 Measure $c(\text{Ca}^{2+}) = 10^{-3}$ mol/L standard solution and write down the potential.
- 3 Calculate the change of potentials with the 2 previously measured potentials:

The value has to be at least 23.7 mV (at 25 °C) (80% of the theoretical slope).



NOTICE

If the value cannot be reached, the electrode needs to be replaced.


7 Combined dCa ISE – Troubleshooting

If additional or other types of faults occur, ensure that the following points have been fulfilled:

- Is the electrode cable screwed on and plugged in correctly?
- Is the electrode cable functional?
- Is the sensor surface clean and intact?
- Is the electrode new?

If the electrode is too old, the membrane may be depleted.

8 Electrode – Disposal



This product is covered by European Directive, WEEE – Waste Electrical and Electronic Equipment.

The correct disposal of your old instrument will help to prevent negative effects on the environment and public health.

Proceed as follows to dispose of the electrode:

1 Draining the electrolyte

Use a plastic pipette to remove the electrolyte from the electrode.

2 Disposing of the electrolyte

Dispose of the electrolyte in accordance with the legal provisions.

3 Disposing of the electrode

Put the electrode in electronic waste recycling.

More details about the disposal of your old product can be obtained from your local authorities, from waste disposal companies or from your local dealer.

10 Additional information

ISA/TISAB solutions

Table 1 ISA/TISAB solutions

Measuring ion	ISA/TISAB	For 100 mL of solution	Remarks
Ca ²⁺	KCl 1 mol/L	7.46 g	

Interfering ions

The concentrations in mol/L of the interfering ions, which generate an analysis error of approximately 10%, are specified in the following table.

Table 2 Interfering ions

Measuring ion	Interferences
Ca ²⁺	c(Na ⁺) < 0.24; c(K ⁺) < 0.4; c(Mg ²⁺) < 18; c(H ⁺) < 0.12; c(OH ⁻) < 0.11; c(Cu ²⁺) < 8·10 ⁻² ; c(Pb ²⁺) < 3.5·10 ⁻² ; c(Zn ²⁺) < 0.22; c(Fe ²⁺) < 0.45