


dPt ring electrode



6.00403.300

Sensor leaflet

8.0109.8009EN / 2020-08-31



Metrohm AG
Ionenstrasse
CH-9100 Herisau
Switzerland
Phone +41 71 353 85 85
Fax +41 71 353 89 01
info@metrohm.com
www.metrohm.com

dPt ring electrode

6.00403.300

Sensor leaflet

Technical Communication
Metrohm AG
CH-9100 Herisau
techcom@metrohm.com

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

RESULTS

1	Overview	1
1.1	dPt ring electrode – Product description	1
1.2	dPt ring electrode – Overview	1
2	Functional description	2
2.1	dPt ring electrode – Functional description	2
3	Transport and storage	3
3.1	Electrode – Checking the delivery	3
3.2	Electrode – Storing the packaging	3
3.3	Unpacking and checking the electrode	3
3.4	Storing the dPt ring electrode	4
4	Installation	5
4.1	Preparing the dPt ring electrode	5
4.2	Mounting the electrode	6
5	Maintenance	8
5.1	dPt ring electrode – Changing/refilling the electrolyte	8
5.2	Checking the dPt ring electrode	8
6	Troubleshooting	10
7	Electrode – Disposal	11
8	Technical specifications	12
8.1	Ambient conditions	12
8.2	dPt ring electrode – Dimensions	12
8.3	dPt ring electrode – Housing	12
8.4	dPt ring electrode – Connectors specifications	12
8.5	dPt ring electrode – Display specifications	13
8.6	dPt ring electrode – Measurement specifications	13

1 Overview

1.1 dPt ring electrode – Product description

The dPt ring electrode is a combined metal electrode for redox titrations with alteration of the pH value.

1.2 dPt ring electrode – Overview




Figure 1 dPt ring electrode

1 Protective cap

3 Filler opening

5 Metal ring

2 Electrode head

4 Diaphragm

2 Functional description

2.1 dPt ring electrode – Functional description

Pt metal electrodes have a bare metal surface that is exposed to the solution.

If redox-active ions are present in the sample solution, a concentration-dependent potential appears on the Pt surface.

This concentration-dependent equilibrium is characterized by a corresponding potential (Galvani potential).

3 Transport and storage

3.1 Electrode – Checking the delivery

Immediately upon arrival of the merchandise, check the shipment to ensure absence of damage.

3.2 Electrode – Storing the packaging

The product is supplied in extremely protective special packaging. Keep this packaging, as only this ensures safe transportation of the product.

3.3 Unpacking and checking the electrode

1 Unpacking the electrode

Remove the electrode with storage vessel from the packaging.

2 Removing the storage vessel




Figure 2 Releasing the electrode from the storage vessel

- Hold the electrode and storage vessel firmly in your hand so that the electrode cannot slip away.
- Position the tool between the storage vessel and SGJ sleeve.
- **Carefully** push the tool to the side to release the electrode.
Do not tip the tool forwards!



3 Checking the electrode for proper function

- ## NOTICE

3.4 Storing the dPt ring electrode

- 1 Screw the protective cap (1-1) onto the electrode head (1-2).
- 2 Store the electrode in the storage vessel. When doing so, ensure that the diaphragm (1-4) is immersed in the corresponding storage solution.



3 Close the filler opening (1-**3**).



4 ■ ■ ■ ■ ■ ■ ■ ■

4 Installation

4.1 Preparing the dPt ring electrode

1 Filling with reference electrolyte

Open the closure of the filler opening (1-3) and, if necessary, fill reference electrolyte up to the filler opening.

2 Cleaning the electrode

- Rinse the electrode with distilled water.
- If there is excess contamination on the metal ring, clean it with a moist towel and toothpaste or with the polishing set (6.2802.000).
- If necessary, use a suitable solvent to degrease the electrode.



NOTICE

The electrode should be rinsed before each measurement.
Frequent abrasive cleaning is not recommended.

3 Connecting the electrode

- Unscrew the protective cap (1-1).
- Position the cable connection on the electrode head such that the slot in the cable connection is on the guide lug of the electrode head.
- Push the socket in the cable connection into the plug inside the electrode head.
- Push the outer ring of the cable connection over the electrode head.
Ensure that the guide lugs in the electrode head are in the grooves of the cable connection.
- Push the cable connection onto the electrode head until it stops and rotate the outer ring until it snaps in place.



NOTICE

To remove the cable, first release the outer ring and then carefully pull the cable connection from the electrode head.

When doing so, be sure not to pull on the cable itself but the cable connector instead.

4.2 Mounting the electrode



The electrode must sit securely in the titration head.



NOTICE

For automatic procedures, ensure that the cables have enough room to move.

During the titration, it is important that the solution is mixed well. The stirring rate should be high enough to form a small vortex. If the stirring rate is too high, then air bubbles will be aspirated. These may result in incorrect measured values. If the stirring rate is too low, then the solution is only mixed slowly and the reaction time or titration time increases accordingly.

In order for the measurement to be taken in a well-mixed solution after the addition of the titrant, the titration tip should be positioned where tur-

bulence is high. Furthermore, the distance between the addition of the titrant and the electrode should be as large as possible. Therefore, take into account the stirring direction (counterclockwise or clockwise) when positioning the electrode and titration tip.





Figure 3 Diagrams showing rod stirrer, electrode and titration tip during a titration. a) clockwise stirring direction, b) counterclockwise stirring direction.

1 Rod stirrer

2 Electrode

3 Titration tip

7 Electrode – Disposal



This product is covered by European Directive, WEEE – Waste Electrical and Electronic Equipment.

The correct disposal of your old instrument will help to prevent negative effects on the environment and public health.

Proceed as follows to dispose of the electrode:

1 Draining the electrolyte

Use a plastic pipette to remove the electrolyte from the electrode.

2 Disposing of the electrolyte

Dispose of the electrolyte in accordance with the legal provisions.

3 Disposing of the electrode

Put the electrode in electronic waste recycling.

More details about the disposal of your old product can be obtained from your local authorities, from waste disposal companies or from your local dealer.



8 Technical specifications

8.1 Ambient conditions

Nominal function range	+5 to +45 °C	at max. 80% relative humidity, non- condensing
Storage	+5 to +45 °C	

8.2 dPt ring electrode – Dimensions

Measurements

<i>Shaft diameter</i>	12 mm
<i>Maximum installation length</i>	125 mm

8.3 dPt ring electrode – Housing

Materials

<i>Shaft material</i>	Glass
-----------------------	-------

8.4 dPt ring electrode – Connectors specifications

Connector	Metrohm plug-in head Q
-----------	---------------------------



Status display	LED	green - red
----------------	-----	-------------

pH range 0 - 14

Temperature range -5 - 80 °C

Minimum immersion depth	20 mm
--------------------------------	-------