


OMNIS Solvent Module



2.1009.0010

Product manual

8.1009.8002EN / 2021-07-23



Metrohm AG
Ionenstrasse
CH-9100 Herisau
Switzerland
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

OMNIS Solvent Module

2.1009.0010

Product manual

8.1009.8002EN /
2021-07-23

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Disclaimer

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g. conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

Table of contents

1	Overview	1
1.1	OMNIS Solvent Module – Product description	1
1.2	OMNIS Solvent Module – Product versions	1
1.3	Symbols and conventions	2
1.4	Further information	2
1.5	Accessories	2
2	Safety	4
2.1	Intended use	4
2.2	Responsibility of the operator	4
2.3	Requirements for operating personnel	5
2.4	Safety instructions	5
2.4.1	Danger from electrical potential	5
2.4.2	Danger from biological and chemical hazardous substances	5
2.4.3	Danger from highly flammable substances	6
2.4.4	Danger from leaking liquids	6
2.4.5	Danger during transport of the product	7
2.5	Design of warning messages	7
2.6	Meaning of warning signs	8
3	Functional description	10
3.1	OMNIS Solvent Module – Overview	10
3.1.1	Adsorption cartridge – Overview	12
3.2	OMNIS Solvent Module – Function	12
3.2.1	Adsorption cartridge – Function	13
3.3	OMNIS Solvent Module – Indicators and controls	13
3.4	System – Signals	14
3.5	OMNIS Solvent Module – Interfaces	15
4	Delivery and packaging	16
4.1	Delivery	16
4.2	Packaging	16
5	Installation	17
5.1	Installation by Metrohm	17
5.2	Setup location	17

5.3	OMNIS Solvent Module – Removing and mounting the adsorption cartridge	18
5.4	OMNIS Solvent Module – Mounting the bottles	19
6	Start-up	25
6.1	Initial start-up by Metrohm	25
7	Operation and control	26
7.1	Operation	26
7.2	Filling and emptying the KF titration cell	26
8	Maintenance	28
8.1	Maintenance	28
8.2	Cleaning the product surface	28
9	Troubleshooting	30
10	Disposal	31
11	Technical specifications	32
11.1	Ambient conditions	32
11.2	OMNIS Solvent Module – Energy supply	32
11.3	OMNIS Solvent Module – Dimensions	32
11.4	OMNIS Solvent Module – Housing	33
11.5	OMNIS Solvent Module – Connectors specifications	33
11.6	Display specifications	33
11.7	OMNIS Solvent Module – Liquid handling specifications	33

1 Overview

1.1 OMNIS Solvent Module – Product description

The OMNIS Solvent Module is a pump module, which is either controlled manually or via the OMNIS Software. The OMNIS Solvent Module is also equipped with features to secure bottles.


1.2 OMNIS Solvent Module – Product versions

The product is available in the following versions:

Table 1 Product versions

Article number	Designation
2.1009.0010	OMNIS Solvent Module

The article number and serial number for identifying the product can be found on the type plate:



1	(01) = Article number in accordance with GS1 standard	2	(21) = Serial number
3	(240) = Metrohm article number	4	Certification
5	Certification	6	Technical specifications

1.3 Symbols and conventions

The following formatting may appear in the documentation:

(5-12)	Cross-reference to figure legend The first number refers to the figure number. The second number refers to the product part in the figure.
1	Instruction step Numbers indicate the order of the instructions steps.
Method	Names of parameters, menu items, tabs and dialog windows
File ► New	Menu path
[Continue]	Button or key

1.4 Further information

The Metrohm Knowledge Base <https://guide.metrohm.com> always provides the current version of this document. Further instructions, leaflets, release notes etc. may be available, depending on the product. You can directly access the required information or the associated PDF document using the full-text search function and filters.

1.5 Accessories

Up-to-date information on the scope of delivery and on optional accessories can be found on the Metrohm website. Download this information as follows:

Downloading the accessories list

- 1 Go to <https://www.metrohm.com>.
- 2 Enter the article number of the product (e.g. **2.1001.0010**) into the search field.

The search result is displayed.
- 3 Click on the product.

Detailed information regarding the product is shown on various tabs.

- 4



Metrohm recommends downloading the accessories list from the Internet and keeping it for reference purposes.

2 Safety

2.1 Intended use

Metrohm products are used for the analysis and handling of chemicals.

Usage therefore requires the user to have basic knowledge and experience in handling chemicals. Knowledge regarding the application of fire prevention measures prescribed for laboratories is also mandatory.

Adherence to this technical documentation and compliance with the maintenance specifications make up an important part of intended use.

Any utilization in excess of or deviating from the intended use is regarded as misuse.

Specifications regarding the operating values and limit values of individual products are contained in the "Technical specifications" section, if relevant.

Exceeding and/or not observing the mentioned limit values during operation puts people and components at risk. The manufacturer assumes no liability for damage due to non-observance of these limit values.

The EU declaration of conformity loses its validity as soon as modifications are carried out on the products and/or the components.

2.2 Responsibility of the operator

The operator must ensure that basic regulations on occupational safety and accident prevention in chemical laboratories are observed. The operator has the following responsibilities:

- Instruct personnel in the safe handling of the product.
- Train personnel in the use of the product according to the user documentation (e.g. install, operate, clean, eliminate faults).
- Train staff on basic occupational safety and accident prevention regulations.
- Provide personal protective equipment (e.g. protective glasses, gloves).
- Provide suitable tools and equipment to carry out the work safely.

The product may be used only when it is in perfect condition. The following measures are required to ensure the safe operation of the product:

- Check the condition of the product before use.
- Remedy defects and malfunctions immediately.
- Maintain and clean the product regularly.

2.3 Requirements for operating personnel

Only qualified personnel may operate the product. Qualified personnel are persons who meet the following requirements:

- Basic regulations on occupational safety and accident prevention for chemical laboratories are known and complied with.
- Knowledge of handling hazardous chemicals is present. Personnel have the ability to recognize and avoid potential dangers.
- Knowledge regarding the application of fire prevention measures for laboratories is available.
- Safety-relevant information is communicated and understood. The personnel can operate the product safely.
- The user documentation has been read and understood. The personnel operate the product according to the instructions in the user documentation.

2.4 Safety instructions

2.4.1 Danger from electrical potential

Contact with electrical potential can cause serious injuries or death. To avoid danger from electrical potential, observe the following:

- Operate the product only if it is in perfect condition. The housing must also be intact.
- Only use the product with the covers fitted. If covers are damaged or missing, disconnect the product from the energy supply and contact the regional Metrohm service representative.
- Protect live components (e.g. power supply unit, power cord, connection sockets) against moisture.
- Always have maintenance work and repairs on electrical components carried out by a regional Metrohm service representative.
- Disconnect the product from the energy supply immediately if at least one of the following cases occurs:
 - The housing is damaged or open.
 - Live parts are damaged.
 - Moisture penetrates.

2.4.2 Danger from biological and chemical hazardous substances

Contact with biological hazardous substances may cause poisoning from toxins or infections from microorganisms. Contact with aggressive chemical substances may cause poisoning or chemical burns. To avoid danger from biological or chemical hazardous substances, observe the following:

- Label the product according to regulations if it is used for substances that have a potential for chemical hazards and are generally subject to the Hazardous Substances Ordinance.
- Wear personal protective equipment (e.g. protective glasses, gloves).
- Use exhaust equipment when working with vaporizing hazardous substances.
- Dispose of hazardous substances in accordance with regulations.
- Clean and disinfect contaminated surfaces.
- Only use detergents that do not cause any unwanted side reactions with the materials to be cleaned.
- Dispose of chemically contaminated materials (e.g. cleaning material) in accordance with regulations.
- Proceed as follows in case of a return shipment to Metrohm AG or a regional Metrohm representative:
 - Decontaminate the product or product component.
 - Remove the labeling for hazardous substances.
 - Create a declaration of decontamination and enclose it with the product.

2.4.3 Danger from highly flammable substances

Using highly flammable substances or gases may cause fires or explosions. To avoid danger from highly flammable substances, observe the following:

- Avoid ignition sources.
- Use protective grounding.
- Use exhaust equipment.

2.4.4 Danger from leaking liquids

Leaking liquids may cause injuries and may damage the product. To avoid danger from leaking liquids, observe the following:

- Check the product and its accessories for leakages and loose connections.
- Replace leaking parts and connecting elements without delay.
- Tighten loose connecting elements.
- Do not loosen tubing connections under pressure.
- Do not remove aspiration tubing under pressure.
- Carefully pull the ends of the tubing out of the containers.
- Carefully let liquids from tubing drain into suitable containers.
- Insert the buret tips completely into the containers.
- Remove and dispose of leaked liquids in accordance with regulations.
- If you suspect that liquid has penetrated the instrument, disconnect the instrument from the energy supply. Then have the instrument checked by a regional Metrohm service representative.

2.4.5 Danger during transport of the product

Chemical or biological substances may be spilled during the transport of the product. Parts of the product may fall down or may be damaged. There is a risk of injury from chemical or biological substances and pieces of broken glass. To ensure safe transport, observe the following:

- Remove loose parts (e.g. sample racks, sample vessels, bottles) before transport.
- Remove liquids.
- Lift and transport the product with both hands on the base plate.
- Lift and transport heavy products only according to instructions.

2.5 Design of warning messages

There are 4 hazard levels for warning messages. The following signal words are used for classifying the hazard levels in warning messages:

- **DANGER** indicates a hazardous situation which, if not avoided, will result in serious injury or death.
- **WARNING** indicates a hazardous situation which, if not avoided, could result in serious injury or death.
- **CAUTION** indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.
- **NOTICE** indicates a hazardous situation which, if not avoided, could result in property damage.

Warning messages differ in design (color and warning sign) depending on the hazard level:



DANGER

Type and source of danger

Consequences when not observing the notice: An irreversible injury that may result in death is very probable.

- Measures to avoid the danger



WARNING

Type or source of danger

Consequences when not observing the notice: A serious injury that may result in death is probable.

- Measures to avoid the danger



CAUTION

Type or source of danger








Consequences when not observing the notice: A minor to moderate injury is probable.





- Measures to avoid the danger

2.6 Meaning of warning signs

This documentation uses the following warning signs:

Table 2 Warning sign according to ISO 7010

Warning sign	Meaning
	General warning sign
	Warning of electrical voltage
	Warning of hand injuries
	Warning of sharp object
	Warning of hot surface
	Warning of biological hazard
	Warning of toxic materials

Warning sign	Meaning
	Warning of flammable materials
	Warning of corrosive substances
	Warning of optical radiation
	Warning of laser beams

Depending on the intended use of the product, the corresponding warning sign stickers must be placed on the product.




Figure 2 OMNIS Solvent Module – Rear

1 MDL connector

MDL = Metrohm Device Link. Connection socket for connecting cables between OMNIS instruments

2 Air vent

Nozzle for inlet air and exhaust air

3.2.1 Adsorption cartridge – Function

When pumping solvent out of the solvent bottle, the following air stream is led through the adsorption cartridge where it is dried.

For this, the adsorption cartridge must be filled with molecular sieve which has to be replaced regularly, see and the adsorption cartridge must be connected to the solvent bottle, see *OMNIS Solvent Module – Mounting the bottles* (see chapter 5.4, page 19).



NOTICE

Make sure

- that the lid of the adsorption cartridge is tightly sealed.
- that the cotton loosely covers the whole base of the adsorption housing.

3.3 OMNIS Solvent Module – Indicators and controls




Figure 4 OMNIS Solvent Module – Indicators and controls

1 Add key

Deliver liquid (Solvent) into the titration cell

2 Aspirate key

Aspirate waste from the titration cell

3 Status display

Multi-colored

3.5 OMNIS Solvent Module – Interfaces



Figure 5 OMNIS Solvent Module – Interfaces and connectors

1 MDL connector

MDL = Metrohm Device Link. Connection socket for connecting cables between OMNIS instruments

2 Air vent

Nozzle for inlet air and exhaust air (e.g. for controlled extraction of methanol vapor with tubing)

4.1 Delivery

- Check the delivery against the delivery note to ensure completeness.
- Check the product for damage.
- If the delivery is incomplete or damaged, contact your regional Mettler representative.

The product and accessories are supplied in protective special packaging. Keep this packaging to ensure safe transportation of the product. If a transport locking device is present, keep this as well for future reuse.

5 Installation

5.1 Installation by Metrohm


As a basic rule, the installation of the system is carried out by the regional Metrohm service representative.

5.2 Setup location

The product is only suitable for operation indoors and may not be used in explosive environments.

The following requirements apply to the setup location:

- The room is well ventilated, protected against direct sunlight and excessive temperature fluctuations.
- The setup space is stable and free of vibrations. The setup space must be suitable for the dimensions and weight of the components (see Technical specifications).
- All cables and connectors are accessible during operation. The cables are safely installed (no tripping hazards).
- The workplace is ergonomically designed and ensures trouble-free operation of the product.



1. Place the cartridge into the seal of the OMNIS Solvent Module and press it down **(1)**. At the same time, tilt it inwards **(2)** until you hear it click into place.
2. Mount the PVC tubing of the solvent bottle onto the tubing olive.

5.4 OMNIS Solvent Module – Mounting the bottles



NOTICE

The reagents used for Karl Fischer titration must be kept as dry as possible, even if they are conveyed from supply bottles that have been open for a long time. The connected adsorption cartridge prevents humidity from entering the reagent bottle.

Setting up the waste bottle (Waste)




Figure 6 Equipping the GL 45 bottle cap for the waste bottle (Waste)

1	M8 PTFE tubing (6.1805.200)	2	Short PTFE cannula (6.1819.050)
3	Ground-joint stopper SGJ 14/M8 (6.1446.090)	4	M6 threaded stopper (6.1446.040)
5	Bottle cap for GL 45 (6.1602.105)	6	Overflow protection (6.1623.000)
7	Tubing olive (6.1808.050)	8	PVC tubing (6.0184.210)

- 1 Place the threaded stopper (6-4) in the M6 connector (smallest opening) of the bottle cap (6-5) and screw it tight.
- 2 Introduce the overflow protection (6-6) from below into the M8 connector (second-smallest opening) of the bottle cap and fasten it.



NOTICE

Make sure that the overflow protection (6-**6**) is connected to the M8 connector from which the PVC tubing (6-**8**) leads to the tubing olive of the OMNIS Solvent Module.

- 3** Place the tubing olive (6-**7**) in the M8 connector of the bottle cap from above and screw it tight.
- 4** Plug a PVC tubing (6-**8**) onto the tubing olive.
- 5** Insert the ground-joint stopper (6-**3**) into the remaining opening of the bottle cap.
- 6** Insert the short PTFE cannula (6-**2**) into the ground-joint stopper (6-**3**) from above and pull it through.
Make sure that the cannula is pulled through as far as it will go.
- 7** Insert an M8 PTFE tubing (6-**1**) into the ground-joint stopper (6-**3**) from above and screw it tight.
- 8** Place the fully equipped bottle cap (6-**5**) on the clear glass bottle (or on a different bottle with GL 45 thread) and screw it tight.

Setting up the reagent bottle (Solvent)



NOTICE

Only use the reagent bottle (Solvent) with Siphon Breaker!
Other bottle caps are not permitted.




Figure 7 Equipping the Siphon Breaker for the reagent bottle (Solvent)

1	PVC tubing (6.01804.210)	2	Tubing olive, integrated
3	Siphon Breaker (6.01600.200)	4	M8 aspiration tubing (6.01805.130)
5	M8 connector, integrated	6	M8 PTFE tubing (6.1805.200)

- 1 Insert the M8 aspiration tubing (7-4) with the screw nipple from below into the M8 connector of the Siphon Breaker (7-3) and screw it tight.



NOTICE

Make sure that the M8 aspiration tubing (7-4) is fastened from below to the Siphon Breaker and the M8 PTFE tubing (7-6) for the KF titration cell is fixed in place from above.

- 2 Insert an M8 PTFE tubing (7-6) into the M8 connector (7-5) of the Siphon Breaker from above and screw it tight.
- 3 Plug a PVC tubing (7-1) onto the tubing olive (7-2).
- 4 Place the fully equipped Siphon Breaker (7-3) on the reagent bottle (Solvent) and screw it tight.

Mounting and connecting the bottles to the OMNIS Solvent Module




Figure 8 Arrangement of the bottles on the OMNIS Solvent Module, the necessary connections and symbols of the tubing connectors – Overview

1 PVC tubings

From the Siphon Breaker to the adsorption cartridge of the OMNIS Solvent Module and from the waste bottle (Waste) to the OMNIS Solvent Module

2 Reagent bottle (Solvent)

3 Bottle holder (2 pieces)

4 Waste bottle (Waste)



NOTICE

Make sure to limit the time the open reagent bottle (Solvent) and the KF titration cell are exposed to humidity.

The KF titration cell must be fully equipped and mounted on the OMNIS Titrator, see , except for the 2 M8 PTFE tubings from the aspiration tip and the dosing tip.

6 Start-up

6.1 Initial start-up by Metrohm

As a basic rule, the initial start-up of the system is carried out by the regional Metrohm service representative.



Addition and aspiration can also be controlled via the OMNIS Software.




Contact with electrical potential can cause serious injuries or death.

- Prerequisite:**

- ### Required accessories:

- 3** Clean the connectors with a dry cloth.

10 Disposal



Properly dispose of chemicals and of the product to reduce negative effects on the environment and public health. Local authorities, waste disposal companies or dealers provide more detailed information on disposal. Observe the WEEE EU directive (WEEE = Waste Electrical and Electronic Equipment) for the proper disposal of waste electronic equipment within the European Union.

11.4 OMNIS Solvent Module – Housing

Materials

<i>Base</i>	1.4301	stainless steel
<i>Enclosure</i>	PBT	poly(butylene terephthalate)
<i>Front foils</i>	PET	poly(ethylene terephthalate), mat

IP degree of protection IP 30

11.5 OMNIS Solvent Module – Connectors specifications

Energy supply via MDL
Socket round plug

MDL Metrohm Device Link

11.6 Display specifications

Status display LED multi-colored

11.7 OMNIS Solvent Module – Liquid handling specifications

Pump

<i>Type</i>		membrane
<i>Number</i>	1	
<i>Flow rate</i>		
Add	600 mL/min	depending on the fill level in the respective bottle
Aspirate	300 - 400 mL/min	depending on the fill level in the respective bottle