

---

Method parameters

---

Method : AB429\_Det of Cu.mth  
Title : Determination of Cu  
Remark1 : 10 mL (diluted) sample + 2 mL electrolyte  
Remark2 : Electrolyte: c(KCl) = 0.3 mol/L, c(HCl) = 0.1 mol/L

Calibration : Standard addition  
Technique : Batch  
Addition : Manual

Sample ID : Sample  
Sample amount (mL): 10.000  
Cell volume (mL): 12.000

Voltammetric parameters

---

Mode : DP - Differential Pulse

Highest current range : 10 mA  
Lowest current range : 100 nA

Electrode : SSE/RDE  
Stirrer speed (rpm) : 2000

Initial electr. conditioning : No

No. of additions : 2  
No. of replications : 2

Measure blank : No  
Addition purge time (s) : 10

Initial purge time (s) : 10

Conditioning cycles  
Start potential (V) : 0.000  
End potential (V) : 0.750  
No. of cycles : 5

Hydrodynamic (measurement) : No  
Cleaning potential (V) : 0.750  
Cleaning time (s) : 10.000  
Deposition potential (V) : -0.300  
Deposition time (s) : 30.000

Sweep  
Equilibration time (s) : 5.000  
Start potential (V) : -0.100  
End potential (V) : 0.600  
Voltage step (V) : 0.006  
Voltage step time (s) : 0.060  
Sweep rate (V/s) : 0.099  
Pulse amplitude (V) : 0.050  
Pulse time (s) : 0.020

Cell off after measurement : Yes

# Peak evaluation

```
-----
Regression technique      : Linear Regression
Peak evaluation           : Height
Minimum peak width (V.steps) : 5
Minimum peak height (A)   : 1.000e-010
Reverse peaks             : No
Smooth factor             : 4
Eliminate spikes          : Yes
```

# Substances

```
-----
Cu                        : 0.250 V   +/- 0.050 V
```

```
Standard solution        : 1      1.000 mg/L
Addition volume (mL)     : 0.050
```

```
Copper                   : Final result (Cu) =
                          Conc * (12 / 10) * (1e+006 / 1) + 0 - 0
```

# Baseline

```
-----
Substance Addition      automatic start (V) end (V) type      scope
-----
Cu      Sample          yes      ---      ---      linear      wholePeak
        Addition 1      yes      ---      ---      linear      wholePeak
        Addition 2      yes      ---      ---      linear      wholePeak
-----
```

# Solutions

```
-----
No. Content                                     Predose (mL)
--- -----
```

# Export options

```
-----
Export final results as ASCII: no
Export final results as CSV:   no
Export final results as XML:   no
Export determination to AutoDB: no
```