## VA Application Note No. V-10

Title:	Zinc, cadmium, lead, copper, iron, nickel and cobalt in NaOH in one run
Summary:	Determination of Zn, Cd, Pb, Cu, Fe, Ni and Co in NaOH 50% in one run
Sample:	NaOH 50%
Sample Preparation:	none

Zinc, cadmium, lead and copper:					
Electrolyte:	Acetate buffer, pH = 4.6 with CH <sub>3</sub> COOH, NH <sub>3</sub> and HCI.				
AE: RE:	Pt Ag/AgCI/KCI 3M				
Parameters:	Zn at the SMDE, $U_{start} = -800 \text{ mV}$ , $U_{end} = -1200 \text{ mV}$ Cd, Pb and Cu at the HMDE, DPASV (+50 mV) $U_{meas} = -700 \text{ mV}$ (90s), $U_{start} = -700 \text{ mV}$ , $U_{end} = +30 \text{ mV}$ Ep (Zn) = -960 mV, Ep (Cd) = -600 mV Ep (Pb) = -380 mV, Ep (Cu) = -120 mV				


Iron determination in the same vessel:				
Electrolyte:	Add catechol crystals and Pipes buffer. pH = 7.0			
AE: RE:	Pt Ag/AgCl/KCl 3M			
Parameters:	DPCSV (-50 mV), HMDE $U_{meas} = -200 \text{ mV}$ (60s), $U_{start} = -200 \text{ mV}$ , $U_{end} = -600 \text{ mV}$ Ep (Fe) = -450 mV			

Nickel and cobalt determination in the same vessel:						
Electrolyte:	Add dimethylglyoxime and $NH_4CI / NH_3$ buffer					
AE: RE:	Pt Ag/AgCI/KCI 3M					
Parameters:	DPCSV (-50 mV), HMDE U <sub>meas</sub> = -700 mV (30s), U <sub>start</sub> = -800 mV, U <sub>end</sub> = -1300 mV Ep (Ni) = -975 mV, Ep (Co) = -1100 mV					

## A Metrohm

Results:	<b>Zn</b>	Cd	<b>Pb</b>	Cu	Fe	Ni	<b>Co</b>
	mg/L	ng/L	µg/L	µg/L	µg/L	μg/L	µg/L
	5.61	204.9	56.3	36.4	608.4	876.7	44

## Determination of cadmium, lead and copper (zinc not shown)



## Determination of iron

Determination of nickel and cobalt

