Application Note AN-T-251

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Determination of carbonate in lithium carbonate by potentiometric titration according to draft standard ISO 10662

Lithium carbonate is one of the main lithium salts in cathode materials used for lithium-ion batteries. Accurate carbonate and overall purity determination in this raw material are essential because even trace impurities can affect key physicochemical properties, ultimately influencing electrochemical behavior in downstream applications. This is especially significant in battery-grade materials, where

strict compositional tolerances are required to ensure consistent manufacturing performance, long-term cell stability, and overall safety. This Application Note describes carbonate content determination in lithium carbonate on an automated OMNIS system with Dis-Cover lids. The ISO 10662 draft method quantifies the main content of lithium carbonate to enable a precise assessment of material purity.

# SAMPLE AND SAMPLE PREPARATION

This method is applicable to a concentration of battery-grade lithium carbonate ( $\text{Li}_2\text{CO}_3$ ) greater than 95.00%. The samples were dried at 250 ° C

for two hours and allowed to cool down before use.

## EXPERIMENTAL

The analysis is carried out on an OMNIS Advanced Titrator with an OMNIS Sample Robot equipped with a dUnitrode and Dis-Cover lids (Figure 1). A blank analysis is performed with a start volume of hydrochloric acid and titrated until after the equivalence point. Samples are titrated with hydrochloric acid until after the second equivalence point.

**Figure 1.** OMNIS Sample Robot and OMNIS Advanced Titrator equipped with a dUnitrode and Dis-Cover lids.

## RESULTS

The samples were tested in duplicate on one day, and once again on the following day.

Results are displayed in Table 1.

**Table 1.** Results of the carbonate determination of lithium carbonate samples.

| Sample (n = 3) | Content mean value (%) | SD(rel) in % |
|----------------|------------------------|--------------|
| A              | 100.37                 | 0.18         |
| B              | 98.62                  | 0.42         |
| C              | 99.87                  | 0.96         |


**Figure 2.** Titration curve of the carbonate content in lithium carbonate.

## CONCLUSION

The OMNIS Titrator and OMNIS Sample Robot with Dis-Cover lids provide a fast and reliable determination of carbonate in lithium carbonate

according to draft standard ISO 10662. The OMNIS system can be expanded to other titration applications as required.

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