

10 years readability guaranteed

Neo's, the new thermal printer generation

With the new Neo's USB printer you can print out and document your analytical results quickly and conveniently right in your lab.

The installation of the Neo's USB printer is extremely easy: connect the printer to the Titrino plus, select the printer type in the configuration menu of the Titrino plus, KF Titrino plus or Dosimat plus and off you go!

The Neo's thermal printer prints out your results, parameters and curves. It works fast, is very quiet and requires a minimum of space. Its footprint is smaller than half a sheet of A4 paper – the ideal size for any laboratory!

Thanks to a new special paper, we can guarantee that the printouts will still be legible after ten years if stored properly! This paper can be ordered from Metrohm.



Small but nice: the new small Neo's USB printer (on the right) is extremely practical and can be connected to a Titrino plus, KF Titrino plus or Dosimat plus without any problems.

The specialist for the automatic rinsing of titration vessels 843 Pump Station

The 843 Pump Station is a versatile pump unit specially designed for laboratory and plant use. It facilitates the automation of many applications in the fields of titration and voltammetry. Each 843 Pump Station contains two pumps that are accommodated in a single housing in a very compact and space-saving manner. Both diaphragm and peristaltic pumps can be selected. Thanks to its two pump heads, the 843 Pump Station can aspirate and rinse simultaneously and in this way, for example, enables the rapid and efficient rinsing of titration vessels. Both pumps can be switched on and off separately via remote signals, which enables them to be used in automated systems. On request, the pumps can also be operated manually as stand-alone pumps, for example for rinsing permanently mounted titration vessels or measuring cells without having to open them.

The Pump Station is available in four versions, with the Rinse / Aspirate versions containing the complete accessories for rinsing:

2.843.0020	843 Pump Station (diaphragm)
2.843.0030	843 Pump Station (diaphragm) – Rinse / Aspirate
2.843.0120	843 Pump Station (peristaltic)
2.843.0130	843 Pump Station (peristaltic) – Rinse / Aspirate



The new 843 Pump Station is equipped with either peristaltic pumps (left) or diaphragm pumps (right).